Leica FS C

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Leica Microsystems Wetzlar GmbH
Ernst-Leitz-Strasse
D-35578 Wetzlar (Germany)

Responsible for contents:
Verantwortlich für den Inhalt:
Responsable du contenu:
Katja Peter, Karin Schwab
Marketing CM, Compound Microscopy, Product Management,
Marketing CM, Compound Microscopy, Productmanagement,
Marketing CM, Compound Microscopy, Product Management ,

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Leica
MICROSYSTEMS

Copyrights

All rights to this documentation are held by Leica Microsystems Wetzlar GmbH. Reproduction of text or illustrations (in whole or in part) by print, photocopy, microfilm or other method (including electronic systems) is not allowed without express written permission from Leica Microsystems Wetzlar GmbH.

The term “Windows” can be used in the following text without further identification. It is a registered trademark of the Microsoft Corporation. Otherwise, no inference with regard to the free usability of product names may be drawn from the use of those names.

The instructions contained in the following documentation reflect state-of-the-art technology and knowledge standards. We have compiled the texts and illustrations as accurately as possible. Nevertheless, no liability of any kind may be assumed for the accuracy of this manual’s contents. Still, we are always grateful for comments and suggestions regarding potential mistakes within this documentation.

The information in this manual is subject to modification at any time and without prior notification.

Contents

1. Important Notes about this Manual	6	8.9 Diaphragm Setting	24
2. Safety Notes	7	8.10 The Smart Move Control Panel	25
2.1. General Safety Notes	7	8.11 Adjusting the Height of the Stand	25
2.2. Electrical Safety	7	8.12 Magnification Calibration	25
2.3 Weights	7	8.13 Oblique Incident Illumination	26
3. Field of Application	9	8.14 Coaxial Incident Illumination	26
4. Installation Site	10	8.15 Transmitted Light Illumination	28
5. Unpacking the Macroscope	11	8.16 UV Illumination	29
6. Overview of the Instrument	13	9. Accessories	30
7. Assembly	14	9.1 Motorized Work Table	30
7.1 Comparison Bridge	14	9.2 Stackable Stages	31
7.2 Illumination	14	9.3 Slide Holder	33
7.3 Eyepieces	18	9.4 Rotating Mechanism for Toolmark Examination	41
7.4 Objectives	18	9.5 Light Panel	41
7.5 Cabling	19	9.6 Filter for Color Contrast	42
8. Operation	21	9.7 Polarization Devices	43
8.1 Switching on the Macroscope	21	9.8 Stage Scales and Calibration Standard	44
8.2 Stages	22	10. Technical Data	45
8.3 Tube	22	11. Care of the Comparison Macroscope	48
8.4 Eyepieces	23	11.1 Dust Cover	48
8.5 Brightness Setting	23	11.2 Cleaning	48
8.6 Focusing	23	11.3 Handling Acids and Alkalines	49
8.7 Objectives	24	12. Wear and Spare Parts	50
8.8 Magnification Changer	24	13. EU Declaration of Conformity	51

1. Important Notes about this Manual

Caution!

This operating manual is an essential component of the Leica FS C Comparison Macroscope, and must be read carefully before the macroscope is put into operation or used.

This operating manual contains important instructions and information for the operational safety and maintenance of the comparison macroscope and accessories. Therefore, it must be kept and taken care of.

Text symbols and their meanings:

(1.2)

Numbers in parentheses such as (1.2), correspond to illustrations (in the example, Figure 1, Item 2).

→ page 20

Numbers with pointer arrows (for example → page 20), point to a certain page of this manual.

Special safety instructions are indicated with the triangle symbol shown here, and have a gray background.

Caution! The macroscope and accessories can be damaged when operated incorrectly.

Danger due to hot surface.

Explanatory note.

Item not contained in all configurations.

2. Safety Notes

2.1 General Safety Notes

This safety class 1 device is constructed and tested in accordance with EN 61010-1/IEC 1010-1, safety regulations for electrical measuring, control, and laboratory devices.

Caution!

In order to maintain this condition and to ensure safe operation, the user must follow the instructions and warnings contained in this operating manual.

Caution!

The devices and accessories described in this operating manual have been tested for safety and potential hazards.

The responsible Leica affiliate or the main plant in Wetzlar must be consulted whenever the device is altered, modified or used in conjunction with non-Leica components that are outside of the scope of this manual.

Unauthorized alterations to the device or noncompliant use shall void all rights to any warranty claims!

2.2. Electrical Safety

General Specifications

For indoor use only.	
Supply voltage:	90-250 V~
Frequency:	50-60 Hz
Power input:	maximum 125 VA
Fuse (power supply):	F 6,3 A L 250 V -UL-
Fuse (power pack):	T 16 A L 250 V
Ambient temperature:	10-36°C
Relative humidity:	maximum 80% to 30°C
Overvoltage category:	II
Pollution degree:	2

2.3 Weights

Comparison bridge:	15 kg
Macro stand with stage:	32 kg
Motorized work table:	50 kg

2. Safety Notes

Caution!

The power plug may only be plugged into an outlet equipped with a grounding contact.

Do not interfere with the grounding function by using an extension cord without a ground wire. Any interruption of the ground wire inside or outside of the device, or release of the ground wire connection, can cause the device to become hazardous. Intentional ground interruption is not permitted!

Caution!

The comparison microscope's electrical accessory components are not protected against water. Water can cause electric shock.

Do not assemble the comparison microscope and its accessory components in the immediate vicinity of a water connection or other areas where the possibility of water exposure exists.

Caution!

By connecting them to ground, ancillary devices providing their own and/or separate power supplies and which are connected to the comparison microscope can be connected to the same ground wire potential. For connections without a ground connector, Leica Service must be consulted.

Caution!

Protect the comparison microscope from excessive temperature fluctuations. Such fluctuations can lead to the accumulation of condensation, which can damage the electrical and optical components.

Caution!

Never use any fuses as replacements other than those of the types and the current ratings listed here. Using patched fuses or bridging the fuse holder is not permitted.

Caution!

Before exchanging the fuses or lamps, be absolutely certain to switch-off the main power switch and remove the power cable.

3. Field of Application

Comparisons of small surfaces or profiles need to be performed with sufficient precision, if the objects of comparison can be observed simultaneously.

If only a macroscope with just one image beam input is available for such comparisons, then, depending on the interruption of observation by a change of objective, high demands are made on the memory of the observer. Therefore, it is impossible to rule out incorrect interpretation in the case of objects with very minor differences in the surface or profile.

The Leica FS C (**F**orensic **S**olution **C**omparison) Macroscope with two image beam inputs is designed for comparing structures of items such as those listed below in split images and composite images in forensic laboratories, document printing facilities, mints and banks, and documenting the results on photographs and videos:

- Evidence on fired ammunition pieces
- Toolmarks
- Manufacturing traces on adhesive tape
- Playback traces on magnetic tape
- Documents
- Writing samples
- Banknotes
- Stamps
- Coins
- etc.

Illumination equipment are available for incident light (including UV light) and transmitted light.

4. Installation Site

Caution:

Ensure that the workplace area is free of oil and chemical fumes. Vibrations, direct sunlight and strong temperature fluctuations interfere with measurements as well as microphotography and image documentation.

Caution:

Electrical components must be assembled at least 10 cm from the wall and away from flammable substances.

The basic prerequisite is a stable instrument table of optimum height (70–80 cm). Combined with an ergonomically correct, adjustable chair, these are the basic requirements for fatigue-free microscopy.

Allowable ambient conditions

Temperature	10–36 °C
Relative humidity	max. 80% up to 30°C

Microscopes in warm and warm-damp climatic zones require special care in order to prevent the build up of fungus.

For additional instructions, see the chapter entitled “Care of the Comparison Microscope” → page 48.

Fig. 1 Leica FS C with adjustable-height work table and ergonomic working chair

5. Unpacking the Macroscope

The stand is delivered in a solid box with packaging inserts.

The delivery includes the following components:

- Stand with base and control panel, as well as the stages, which are already installed. The stand is equipped with two carrying handles.
- Comparison bridge with integrated tube
- One or two cold light sources
- Accessory case with slide holder
- Optical accessories
- Oblique incident illuminator
- Optional accessories

Fig. 2 Stand with carrying handles

- 1 Removable carrying handle
- 2 Transport anchors

When unpacking the equipment, follow these instructions:

- After opening the box, first remove the comparison bridge and then the cardboard boxes, which hold the accessories.
- Lift the stand out of the carton by its carrying handles (2.1).

! **Caution!**

When lifting the macroscope out of the packaging, and when transporting it to its installation site, carry it by its handles only.

Caution:

Note that the stand weighs approx. 32 kg. Two people are required to lift it out of its packaging and carry it to its installation site.

Later, the carrying handles can be unscrewed and removed from the upper part of the macroscope. Then, seal the threaded holes using seal plugs.

5. Unpacking the Macroscope

- Carefully remove all components from the transportation and packaging materials.

Caution!

If at all possible, avoid touching the lens surfaces of the objectives. If fingerprints do appear on the glass surfaces, remove them with a soft leather or linen cloth. Even small traces of finger perspiration can damage the surfaces of optical surfaces in a short time. For additional instructions, see the chapter entitled “Care of the Comparison Macroscope” → page 48.

- Check the delivery carefully against the packing list, delivery note or the invoice. We highly recommend keeping a copy of these documents with the manual in order to have delivery date and delivery package information available when (for example) placing repeat orders or making service calls. Please make sure that no small parts are lost in the packaging materials.

Note:

Keep the packaging materials for storage purposes and transportation of the individual components.

Caution:

Do not yet connect the comparison macroscope and peripherals to the power supply at this point!

6. Overview of the Instrument

Fig. 3 Leica FSC with motorized work table

- 1 Comparison bridge
- 2 Special articulated arms for oblique incident illuminator
- 3 Universal cold light source
- 4 Motorized work table
- 5 Control unit
- 6 Transmitted light illuminator
- 7 Motorized stage
- 8 Objective turret with objectives
- 9 Tube with eyepieces

7. Assembly

For instructions for assembling the stackable stages and various slide holders, please refer to Chapter “9. Accessories” → page 31 ff.

! **Caution!**

The stages and slide holders have to be mounted before screwing in the objectives.

Only a few commonly used screwdrivers and keys are necessary for assembly, which are included in the delivery package.

! **Caution!**

First, remove the transport anchors (4.1).

7.1 Comparison Bridge

- From above, set the comparison bridge onto the dovetail guide (4.2) and push the bridge downwards until the upper end of the dovetail guide is tightly fastened.
- Fasten the comparison bridge using the side clamping screw.

7.2 Illumination

Either one or two cold light sources can be used. We recommend using the same cold light source for both sides based on the illumination method being used. A double-armed, flexible light line is included for this purpose. The second cold light source should then be used for another illumination method.

Fig. 4 Stand with carrying handles

- 1 Transport anchors
- 2 Dovetail guide

Oblique Incident Illumination

- Screw the special articulated arms onto the left and right of the microscope carrier (5.1).
- Screw the light guide with \varnothing 4.5 mm (Fig. 6) onto the cold light illuminator (8.1).
- Place the cold light source next to the stand or onto the shelf that is to be screwed onto the base of the comparison microscope (7.1).
- Insert both light guide ends into the left and right receptacles of the special articulated arms (5.2) and tighten them using an Allen wrench.
- Mount and fasten the focusing optics (6.1).

Fig. 5 Special articulated arms for oblique incident illumination

- 1 Fastening screws for the articulated arm
- 2 Receptacle for light guide

Fig. 6 Light guide \varnothing 4.5 mm

- 1 Focusing optics

Fig. 7 Shelf for the cold light illuminator

- 1 Fastening screws for the shelf


Fig. 8 Cold light illuminator KL2500 LCD

- 1 Connection for light guide

7. Assembly

Coaxial incident light

- Remove the cover plates (9.1) to the rear left and right on the macroscope carriers.
- In their place, screw on the coaxial illuminator (11.1)
- Screw the light guide with \varnothing 9 mm (Fig. 12) onto the cold light illuminator (8.1).
- Place the cold light source next to the stand or onto the shelf that is to be screwed onto the base of the comparison microscope (7.1).
- Insert the polarization filter cube (Fig. 10). To do so, attach the filter cube onto the side and push it into the holder (11.2).

- Insert the two light guide ends into the left and right receptacles on the coaxial illuminator (11.3).

Note:

If, after installation, point-like objects are appear in the image, the light guide must be pulled 1 to 2 mm out of the coaxial illuminator. The point-like objects are caused by the fibers of the light guide.


Fig. 10 Polarizing filter cube
1 From the rear
2 From the front

Fig. 9 Fastening the coaxial illuminator
1 Fastening screws for the cover plate

Fig. 11 Coaxial illuminator
1 Fastening screws for the coaxial illuminator
2 Filter cube inserted
3 Receptacle for light guide

Transmitted Light


The transmitted illuminator is pre-installed in the factory.

- Screw the light guide with \varnothing 9 mm (Fig. 12) onto the cold light illuminator (8.1).
- Place the cold light source next to the stand or onto the shelf that is to be screwed onto the base of the comparison microscope (7.1).

Fig. 12 Light guide \varnothing 9 mm

Fig. 13 Transmitted light illuminator
1 Connection for light guide


- Insert the two light guide ends into the left and right receptacles on the transmitted light illuminator (13.1).
- Connect the cables of the transmitted light illuminator on both sides with the rear side of the stand (17.2, 17.3).

UV Illumination

- Screw the holder for UV lamps (14.1) onto the comparison bridge.
- Insert the UV lamps (14.3) into this holder. Align them with the object and screw them tight.

Fig. 14 Assembly of the UV illuminator

- 1 Holder for UV lamps
- 2 Rod for fastening a UV lamp
- 3 UV lamp

7. Assembly

7.3 Eyepieces

- The eyepieces are inserted into the eyepiece tubes on the tube.

7.4 Objectives

The position of the individual objectives is preset at the factory. The object turret eyes are numbered.

- Screw the objective into the objective turret (Fig. 15).

Eye	Objective
No. 1	0.4x
No. 2	1x
No. 3	2x
No. 4	4x

Fig. 15 Objective turret from below

Fig. 16 Objective turret with objectives

1 Filter receptacle

7.5 Cabling

The connector panel is labeled.

- Connect the left and right stage cables to terminals (17.4) and (17.6).
- Terminal (17.5) on the rear of the stand is connected to terminal (18.1) on the rear side of the comparison bridge. For this purpose, use the 25-pin bridge cable.
- The remote control cables of the cold light sources are connected to terminals (17.9) and (17.10). When using a cold light source, select the corresponding terminal.


- If you have not already done so when installing the transmitted light illuminator, connect the cables for the transmitted light condenser to terminal (17.2) and/or (17.3).
- The optional Smart Move control panel is connected to the serial interface 3 (17.11).
- The RS232C interface (17.7) is intended for connecting a PC.
- Finally, connect the comparison macroscope and the cold light illuminator to the power supply using the power cable.

Fig. 17 Stand connector panel

- 1 Power supply
- 2 Transmitted light condenser lens, right
- 3 Transmitted light condenser lens, left
- 4 Stage, right
- 5 Comparison bridge
- 6 Stage, left
- 7 Serial interface 1 (PC)
- 8 Serial interface 2
- 9 Cold light source, right
- 10 Cold light source, left
- 11 Serial interface 3 (Smart Move)

Fig. 18 Rear side of the comparison bridge

- 1 Connection

7. Assembly

- UV illumination:

The UV lamps are rated for a 230 V/50 Hz power supply.

Caution:

If this voltage supply is not available, a pre-transformer must be used.

For further details, please refer to the special instructions included with the UV lamps.

Caution:

Wear safety goggles when working with UV lamps!

8. Operation

The Leica FS C Comparison Microscope offers the following functions:

- Simultaneous or separate methods for the left and right stage
- Brightness setting
- Left and right-side focussing
- Left and right-side objective turret
- Height-adjustable stand

Additionally, the comparison bridge contains

- Motorized diaphragms for various imaging techniques
- A magnification changer
- A magnification calibration with control LED


8.1 Switching on the Microscope

The switch (19.1) switches on all components integrated into the stand, as well as the cold light illuminator.

Wait until the initialization of the stages is complete.

Fig. 19 Right stand side

1 On/off switch

8. Operation

8.2 Stages

The stages on the left and right side can be moved either independently of each other, or simultaneously in the X and Y directions.

The right stage is positioned using the setting wheels on the front right (20.5) in the X direction and on the right side (20.6) in the Y direction.

For the left stage, the setting wheels are, respectively, at the front left (20.3) and on the left side (20.2).

To move both stages simultaneously, press the button (20.10). The LED underneath the button lights up to indicate that synchronous adjustment is active. In this case, the setting wheels for the X and Y directions, as well as the focusing handwheel, have the same function on the left and right side.

Note:

The travel speed of the stages is adjusted to the magnification setting. When viewed through the observation tube, the speed always appears the same when rotating the setting wheel. In fact, for lower magnification the stage travels faster.

8.3 Tube

Note:


Close any unused tube openings, as otherwise stray light can interfere with observation.

Adjusting the viewing distance

- Adjust the viewing distance of the eyepiece tubes so that a congruent total image is seen.

Fig. 20 Control unit

- 1 Focus (left stage)
- 2 Adjust in Y direction (left stage)
- 3 Adjust in X direction (left stage)
- 4 Adjusting the height of the stand
- 5 Adjust in X direction (right stage)
- 6 Adjust in Y direction (right stage)
- 7 Focus (right stage)
- 8 Right cold light source
- 9 Control LED (right cold light source)
- 10 Synchronous adjustment in X, Y, Z direction
- 11 Control LED (left cold light source)
- 12 Left cold light source

Adjusting the Angle of Vision

- The viewing angle can be adjusted from 5° to 32° by tilting the binocular eyepiece.

8.4 Eyepieces

Note:

The eyepiece's aperture protector must be removed during microscopy while wearing eyeglasses, or at least folded back.

Eyeglasses with multifocal lenses (bifocals and smooth view glasses) must be removed while operating the macroscope.

Correction for Vision Problems

- With your right eye, look through the right eyepiece and bring the specimen into sharp focus.

If both oculars have an adjustable eyelens, the eyelens of the right ocular can be brought to the position for normal vision, which is indicated by an encircling line.

When using a graticule in the right eyepiece, the eyelens must first be focused on the graticule.

- Then, view the same specimen using the left eye and rotate the eyelens of the left eyepiece until the object appears in sharp focus. While doing so, do not use the focus dial.

8.5 Brightness Setting

Cold Light Illumination

The control dials (20.8 and 20.12) adjust the first and (if applicable) the second cold light source. The function is dependent on the cabling on the rear side of the stand → page 19.

We recommend using the same cold light source for both sides based on the illumination method being used. The second cold light source should then be used for another illumination method.

Illumination method → page 26ff.

Note:

The LED (20.9) or (20.11) is green if the corresponding lamp is operated in standard mode. If the LED is red, the bulb life of the lamp is reduced.

UV Illumination

The UV lamps are equipped with a high-pressure, 180 W quartz radiator and UV filter disks for a wavelength of 366 nm. They are rated for a direct 230 V/50 Hz power supply.

For further details → see the special instructions.

8.6 Focusing

You can focus using the handwheels on the left and right side (20.1 and 20.7).

If the button (20.10) is pressed, i.e. if the synchronous adjustment of the stages is activated, focusing is also synchronous.

8. Operation

Note:

The focus, like the movement of the stages, is also adjusted to the magnification setting. When viewed through the observation tube, the speed always appears the same when rotating the setting wheel.

8.7 Objectives

The desired objectives are swiveled into place manually.

All objectives have an iris diaphragm. With the iris diaphragm, the NUM aperture can be adjusted continuously. Five ratchets make it easy to adjust to identical settings on both sides.

The M macro-objectives can be used for coaxial incident light. To avoid reflex, ARM caps must be used in conjunction with the built-in pole filters → page 27ff.

8.8 Magnification Changer

Use the button (21.3) to switch the magnification changer on and off.

When the magnification changer is switched on, the current magnification of both images is increased by a factor of 1.5x. The LED underneath the button (21.3) is illuminated.

The magnification changer controls the visual and photo outputs simultaneously.

8.9 Diaphragm Setting

The motorized diaphragms allow transition from a superimposed image to the separate view of the right and left image, as well as the setting of the dividing line.

Four diaphragm positions are already predefined and can be selected by pressing the corresponding button (21.4):

- Mix** The left and right images are superimposed (composite image method).
- LR** The left and right image appear side-by-side. A line in the middle of the field of vision separates the two images (split-image method).
- L** Only the left image is displayed.
- R** Only the right image is displayed.

The LEDs below the selector buttons (21.4) are illuminated if one of the predefined settings has been selected or set manually. Similarly, the corresponding LED lights up if the dividing line is moved.

The dividing line is moved to the left or to the right by rotating the knob (21.2).

The knob (21.1) is used to change the width of the dividing line.

Fig. 21 Control panel

- 1 Adjusting the width of the dividing line
- 2 Adjusting the dividing line in x direction
- 3 Magnification changer on/off
- 4 Predefined diaphragm settings
- 5 Adjusting the height of the stand

8.10 The Smart Move* Control Panel

Alternatively, you can connect the Smart Move control panel.

The Smart Move control panel allows control of the stage (22.1, 22.2), as well as setting the focus (22.3).

Use the buttons (23.1 and 23.2) to specify whether the settings are to apply to the left or right stage. The synchronous adjustment is activated using the buttons (23.3).

Fig. 22 Smart Move control panel

- 1 Movement in y direction
- 2 Movement in x direction
- 3 Focus setting
- 4 Individual setting of the knob height position


Fig. 23 Smart Move control panel buttons

- 1 Left stage
- 2 Right stage
- 3 Synchronous stage adjustment

8.11 Adjusting the Height of the Stand

The rocker switch (21.5) adjusts the height of the stand individually. This is provided for aligning the height of objects or object holders which vary in height.

8.12 Magnification Calibration

For specialized applications, such as calibrating thermally modified replicas, the magnification of the right beam path can be changed in the range of +/-4% relative to the left beam path.

Use the button (24.1) to make this adjustment. To prevent the setting from being changed accidentally, the button can only be turned while it is pressed down.

In neutral position, the LED (24.2) is green, i.e. magnification calibration is activated. Otherwise, the LED is red.

Caution!

The objective magnifications are accurately preset at the factory. It is not necessary to calibrate the magnification calibration for this purpose.

Fig. 24 Magnification calibration

- 1 Magnification calibration adjusting knob
- 2 Control LED

8. Operation

8.13 Oblique Incident Illumination

The light is guided to the incident light object by the cold light illuminator KL2500 (25.3) and the double-armed, flexible light guide with **4.5 mm diameter**. Insert the light guide ends into the left and right receptacles on the special articulated arm (25.2). If applicable, attach and tighten the focusing optics (25.1).

The light bundle can be rotated and tilted in every direction. The scales on the movable articulated arms can be used to adjust the illumination accurately. The cold light source KL2500 is used for fine adjustment of the light intensity. It can be adjusted electrically and mechanically. For electrical light control, the color temperature is changed. For mechanical light control, the color temperature remains constant.

Fig. 25 Oblique incident illumination

- 1 Focusing optics
- 2 Receptacle for light guide with 4.5 mm diameter
- 3 Cold light illuminator KL 2500

Note:

The LED (20.9) or (20.11) is green if the corresponding lamp is operated in standard mode. If the LED is red, the bulb life of the lamp is reduced.

The KL2500 can be controlled using the rotary knobs (20.8 and 20.12).

Please make sure that the KL2500 is connected with the included cables on the rear side of the stand.

The light and shadow effect of the oblique illumination yields very pronounced surface reliefs of bullet and shell casing traces, toolmarks, as well as other form, material, and textile texture traces. For reflection protection, a polarization filter can be attached to the focusing optics.

8.14 Coaxial Incident Illumination

The coaxial incident illuminator is mounted on the macroscope carrier. It can be retrofitted at any time. The light is guided to the incident light object by the cold light illuminator KL2500 (25.3) and the double-armed, flexible light guide with **9 mm diameter**. Insert the light guide ends into the left and right receptacles of the coaxial illumination (26.3).

To use the coaxial incident illuminator, it must be inserted completely. When switching to another illumination method, the coaxial incident illumination must be pulled out of the beam path (to the first stop).

**Note:**

If, after installation, point-like objects are appear in the image, the light guide must be pulled 1 to 2 mm out of the coaxial illuminator. The point-like objects are caused by the fibers of the light guide.

The specimen is illuminated by the illumination optics, the reflector and the macro-objectives 1x, 2x and 4x. The 0.4x macro-objective is not intended for the coaxial illumination, because the object field cannot be illuminated completely.

The light source is centered and yields a homogenous object field illumination of up to 22 mm field of vision. The cold light source KL2500 is used for fine adjustment of the light intensity.

Fig. 26 Coaxial incident illuminator in assembly position

- 1 Filter cube inserted
- 2 Stop rod
- 3 Receptacles for light guide(s) with 9 mm diameter

It can be adjusted electrically and mechanically. For electrical light control, the color temperature is changed. For mechanical light control, the color temperature remains constant.

**Note:**

The LED (20.9) or (20.11) is green if the corresponding lamp is operated in standard mode. If the LED is red, the bulb life of the lamp is reduced.

The KL2500 can be controlled using the rotary knobs (20.8 and 20.12).

Please make sure that the KL2500 is connected with the included cable on the rear side of the stand.

For reducing remaining internal reflections, attach the anti-reflex cap (27.1) to the objective currently being used. Rotate it until ideal image contrast is achieved. This only works in connec-

Fig. 27

- 1 Anti-reflex cap

8. Operation

tion with the reflector with built-in pol filter. For this purpose, the coaxial illumination unit must be inserted completely into the beam path.

In particular, traces on smooth, highly reflective surfaces appear true-to-detail and contrast-rich using coaxial illumination.

8.15 Transmitted Light Illumination

The transmitted light illuminator is factory pre-installed and must be ordered along with the stand. It cannot be retrofitted later.

The light is guided to the transparent transmitted light object by the cold light illuminator KL2500 and the double-armed, flexible light guide with **9 mm diameter**. Insert the light guide ends into the left and right receptacles of the transmitted light illumination (28.1).

The specimen is illuminated by the collector lens, the diffusion filter, the mirror, and condenser system.

The encoded objective turret recognizes the objective currently being used, and the condenser is automatically adjusted to the current objective magnification.

The light source is centered and yields a homogeneous object field illumination of up to 55 mm field of vision. The cold light source KL2500 is used for fine adjustment of the light intensity.

It can be adjusted electrically and mechanically. For electrical light control, the color temperature is changed. For mechanical light control, the color temperature remains constant.

Note:

The LED (20.9) or (20.11) is green if the corresponding lamp is operated in standard mode. If the LED is red, the bulb life of the lamp is reduced.

The KL2500 can be controlled using the rotary knobs (20.8 and 20.12).


Please make sure that the KL2500 is connected with the included cables on the rear side of the stand.

Likewise, the cable for the condenser lens on the rear side of the stand must be connected.

Examination with transmitted light illumination is a valuable tool for detecting evidence of manipulations and forgeries of documents and certificates, as well as for identifying trace evidence in textile and material textures, foils, adhesive tape, etc.

Fig. 28 Transmitted light illuminator

1 Receptacle for light guide with 9 mm diameter

8.16 UV Illumination

The UV lamps are fastened to the comparison bridge using a special holder (29.1) and are intended for incident illumination of objects with a large surface area with fluorescent parts (e.g. documents, currency, securities). They are equipped with 180W high-pressure quartz radiators and filter cubes for UV light with 366 nm wavelength.

The UV lamps are rated for a 230 V/50 Hz power supply.

**Caution:**

For other voltage supplies, a pre-transformer must be used!


For further details, please refer to the special instructions included with the UV lamps.

**Caution:**

Wear safety goggles when working with UV lamps!

Fig. 29 Assembly of the UV illuminator

- 1 Holder for UV lamps
- 2 Rod for fastening a UV lamp
- 3 UV lamp

9. Accessories

A specially developed work table and correctly sized work chair allow ergonomically optimized work.

In addition to various optical accessories, such as filters and stage scales, a number of different slide holders can be used.

9.1 Motorized Work Table

The work table with motorized height adjustment is designed for optimum ergonomic conditions when working at the comparison microscope.

Its table surface measures 120 cm x 56 cm and it is mounted symmetrically on the height-adjustable column.

Fig. 30 Motorized work table

- 1 Buttons for height adjustment
- 2 Adjustable-height feet

Caution:

When setting up the comparison microscope, make sure that the stand is positioned exactly over the adjustable-height column.

After the power supply is connected, the work table can be continuously adjusted within a range of approx. 62 to 92 cm over floor height. To do so, use the buttons (30.1).

The lifting speed is approx. 12 mm/sec. The lifting capacity is approx. 2000 N.

If the floor on which the work table is installed is not completely level, differences in level can be compensated for by screwing the feet (30.2) in or out.

9.2 Stackable Stages

The macroscopes are equipped with motorized cross-stages.

Stage surface:	160 mm x 220 mm
Glass stage plates:	80 mm x 80 mm
Travel range of the stage:	50 mm x 50 mm

On the stage surface, there are threaded bore holes as receptacles for the adjustable holder (universal holder), rotating stages \varnothing 118 mm, bullet manipulators K2715A, large sample stages and rotating device for toolmark examinations. For transmitted light, removable glass stage plates are also provided (Fig. 31).

Fig. 31 Table with glass stage plate

Rotating stages \varnothing 118 mm

Rotating stages, which can be stacked on the basic equipment cross-stages (Fig. 32), are fastened using two Allen screws. They are suitable for incident light and transmitted light examinations. The rotating stages each have one opening in the center, which is sealed with a removable glass plate (32.3). The threaded bore holes on the stage surfaces are provided for screwing on the adjustable holder. Additionally, there are also holes for inserting the stage clips. The rotating stage surface can be locked into place using a fastening screw (32.2).

In the bottom parts of the rotating stages, there are receptacles for inserting the polarizers (32.1) with λ -plates (red = OK). These can be rotated 90° .

Fig. 32 Rotating stage \varnothing 118 mm

- 1 Filter slot for polarizer
- 2 Fastening screw
- 3 Glass plate

9. Accessories


Large sample stages

The stackable large sample stages can be screwed directly onto the motorized cross-stages. They are primarily used for holding documents. Documents can be securely held in place and stretched using the magnetic clamps provided (33.1).

Of course, the large sample stages can also be used to hold other large-surface carrier material, such as materials bearing toolmarks.

Fig. 33 Large sample stage (assembly)

1 Magnetic clamps

Inclining rotating stages ø 75 mm

The inclining rotating stages are likewise screwed directly onto the motorized cross-stages. After releasing the lever (34.1), you can rotate them horizontally and incline them up to 45° using the built-in ball-and-socket joint (34.2). After positioning the comparison objects, be sure to close the clamping system using the lever.

The surface of the stages (which are 75 mm in diameter) has a ridged surface (34.3) for fastening objects using plasticine. Four bore holes are provided, into which stage clips (34.4) can be inserted.

In the center of the stage plate, there is a conical bore hole (34.5), into which all object receptacles can be inserted (in a similar fashion as the adjustable holder).

Fig. 34 Inclining rotating stage ø 75 mm

1 Clamping lever, 2 Ball-and-socket joint, 3 Ridged surface

4 Bore hole for stage clips, 5 Conical bore hole

9.3 Slide Holder

Adjustable holder with integrated rotating stage

Fasten the adjustable holders (universal holders) (35.1) to the motorized cross-stage using the two screws (35.4).

After the clamping screw (35.3) is released, the universal holder can be moved sideways. This is especially useful for long, horizontally clamped pieces of evidence.

You can adjust the adapter sleeve (35.6) in the segment guide to an angle of up to 90° after releasing the clamping screw (35.2).

The adapter sleeve itself has a rotation device with gradation (35.7) and can be fixed into place using the clamping screw. The receptacles for pieces of evidence listed below can be inserted into the conical bore of the adapter sleeve (35.6).

After the clamping screw (35.5) is unscrewed, the entire unit can be rotated around its own axis.

Fig. 35 Adjustable holder

- 1 Adjustable holder
- 2 Locking device: angle adjustment of the adapter sleeve
- 3 Locking device: lateral adjustment of the holder
- 4 Fastening onto stage
- 5 Locking device: rotation of the holder
- 6 Adapter sleeve
- 7 Rotating device with gradation and clamping screw

9. Accessories

Bullet Holder

To assemble the bullet holder (36.1), slide the adapter sleeve in the segment guide downwards to the "0" marking. Guide the bullet holder vertically to the adapter sleeve and fasten it with the knurled screw (36.2). Adjust the adapter sleeve horizontally again (slide upwards to the "90" marking) and press the rubber-coated pressure plate (37.3) into the conical bore of the adapter sleeve. After unscrewing the clamping screw (37.2), you can move the centering insert (37.2) horizontally in the bullet holder, thus adjusting it to the length of the projectiles currently being compared.

Caution:

Before removing the unit, be sure there is sufficient distance between it and the objective.

Fig. 36 Adjustable holder with bullet holder (assembly)

- 1 Bullet holder
- 2 Fastening screw


Fig. 37 Adjustable holder with bullet holder

- 1 Rubber-coated pressure plate
- 2 Clamping screw
- 3 Centering insert with spring system
- 4 Bullet

Paste-on Cylinder for Bullets

Using adhesive – or using piceine – you can stick (for example) deformed bullets that fit neither into the bullet receptacles nor into the receptacles for deformed pieces of ammunition onto the paste-on cylinder (38.1). There are paste-on cylinders with diameters of 5, 6, 8 and 10 mm.

Receptacles for deformed bullets

Receptacles for deformed bullets (39.1) are available for holding deformed bullets, which cannot be held by the receptacles described above. Their cavities sufficiently hold such objects in place.

Fig. 38 Adjustable holder with paste-on cylinder for bullets

- 1 Paste-on cylinder
- 2 Bullet

Fig. 39 Adjustable holder with receptacle for deformed bullets

- 1 Receptacle for deformed bullets

9. Accessories

Plug Receptacles for Bullet Casings

For fastening undamaged bullet casings of various calibers, plug receptacles (40.1) are provided. Their plastic holding devices are perfectly suited for holding bullet casings (40.2). The various calibers are marked on the rear sides of the plug receptacles. Worn-out holding devices can be replaced. Corresponding replacement parts are included in the delivery.

Plug receptacles for air rifle ammunition

The stainless steel pyramids (41.1) are intended for attaching air rifle ammunition.

Fig. 40 Adjustable holder with plug receptacle for shell casings

- 1 Plug receptacle
- 2 Shell casing

Fig. 41 Adjustable holder with plug receptacles for air rifle ammunition

- 1 Plug receptacle

Shotgun shell casing holder

The shotgun shell casing holder (42.1) has three spring rods rounded at the edges (42.2), to which shotgun shell casings (43.1) of various calibers can be attached. The conical base (43.2) is provided for centering the shell casings.

Fig. 42 Adjustable holder with shotgun shell casing holder

- 1 Shotgun shell casing holder
- 2 Spring rods

Fig. 43 Adjustable holder with shotgun shell casing holder

- 1 Shotgun shell casing
- 2 Conical base

Wire holder

In the wire holder (44.1), round parts between 0.6 mm and 8 mm in diameter can be clamped in order to (for example) examine cutting marks.

Receptacles for locking cylinders

Using an underlay of plasticine, these receptacles (45.1) can be used to securely fasten lock cylinders of safety locks bearing toolmarks.

Fig. 44 Adjustable holder with wire holder

- 1 Wire holder

Fig. 45 Adjustable holder with receptacle for locking cylinders

- 1 Receptacle for locking cylinder

9. Accessories

Stages \varnothing 20 mm

You can fasten pieces of evidence of various dimensions to the ridged surface of the stages with \varnothing 20 mm (46.1). As the adhesive, you can use plasticine or wax.

Stages \varnothing 60 mm

Objects with a circumference of up to 55 mm and a thickness of up to 15 mm can be fastened with clamps (47.2) to the stages \varnothing 60 mm (47.1). To do so, release the screws (47.3), raise the clamps, attach the objects and press the clamps back in before tightening the screws. The additional threading (47.4) on the border of the stages \varnothing 60 mm is for arranging the clamps asymmetrically. If you remove these, the stages \varnothing 60 mm can also be used as large paste-on stages. You can fasten pieces of evidence of various dimensions to the ridged surface of the stages using plasticine.

Fig. 46 Adjustable holder with stage \varnothing 20 mm

- 1 Stage \varnothing 20 mm
- 2 Joint holder

Fig. 47 Adjustable holder with stage \varnothing 60 mm

- 1 Stage \varnothing 60 mm
- 2 Clamp
- 3 Fastening screw
- 4 Thread

Stages ø 25 mm with ball-and-socket joint

The stages ø 25 mm with ball-and-socket joint (48.1) are inserted into the adapter sleeve of the adjustable holder.

After releasing the clamping screw (48.2) on the ball-and-socket joint, the stage surface, which has a ridged surface, can be inclined arbitrarily. After aligning the object, retighten the clamping screw (48.2).

Joint holder

If, for example, you want to compare the casing sleeves of shell casings, the profiles of the strike pin impressions on the percussion caps or the impressions of the extractor one after the other in one clamping, it is necessary to use joint holders (49.1). Its joint is fastened with a knurled screw (49.2).

The joint holders are provided for use as connecting pieces for the above-mentioned receptacles, with the exception of the bullet holder. In addition to the adjustment options provided by the adjustable holders, they offer an additional option for rotating and inclining the objects. The incline of the objects can be locked in place using the knurled screw (49.2).

Fig. 48 Adjustable holder with stage ø 25 mm with ball-and-socket joint

- 1 Stage ø 25 mm with ball-and-socket joint
- 2 Clamping screw

Fig. 49 Adjustable holder with joint holder

- 1 Joint holder
- 2 Clamping screw

9. Accessories

Bullet manipulators K2715A

As an alternative to the previously described object holders, the bullet manipulators K2715A (order No. 581 075) can be used (Fig. 50).

The bullet manipulators are attached directly to the motorized cross-stage. During assembly, observe the markings for "Left" and "Right".

The bullet manipulators are used for the horizontal examination of bullets and shells. However, a vertical bullet casing holder (order No. 581 074) is also available.

Additional accessories:

Mini stage plate (order No. 581 076)
for examining small documents, toolmark specimens, etc.

Rotating mini plate (order No. 581 077)
for examining toolmark specimens such as cut wire, which is embedded in wax and can be inspected from all angles and rotated.

Brush holder set (order No. 581 078)
for examining deformed bullets

Fig. 50 Cross-stage with attached bullet holder K2715A

Fig. 51 Brush set for bullet holder K2715A

9.4 Rotating Mechanism for Toolmark Examination

The rotating mechanism (order No. 581 088) (Fig. 52) consists of a stage plate, which can be mounted directly onto the motorized cross-stage, and an articulated illumination arm, which allows you to rotate the illumination around the specimen.

Either the inclining rotating stage \varnothing 75 mm or the universal object holder can be installed on the rotating plate.

9.5 Light Panel

For homogeneous illumination of larger objects, you can use the light panel (Fig. 53).

It is installed on a separate articulated arm or onto the oblique incident illuminator and connected to the cold light illuminator KL2500 via the 9 mm light guide.

Dimensions

Housing: 139.7 mm x 50.8 mm

Light panel: 81.28 mm x 30.48 mm

Fig. 52 Rotating mechanism for toolmark examination


Fig. 53 Light panel

9.6 Filter for Color Contrast

For the differentiation of color deviations of two objects being compared, the two filters CRA (red) and CGA (green) in sliders (Fig. 54) are provided. They are inserted into the filter receptacles in the macroscope holders (55.1) above the objectives.

For adjusting the superimposed image with color contrast, first set the button on the control panel to **MIX** (for composite imaging).

If the objects of comparison overlap in their essential structures, these only appear in their original color where matching details are present.

Depending on the assignments of the filters to the two beam paths, the deviations appear in either red or green.

This requires that the intensities of both light sources are carefully matched.

Fig. 54 Color filter

Fig. 55 Objective turret with objectives

1 Filter receptacle

9.7 Polarization Devices

Polarizing devices are available for the comparison microscope for both the oblique incident illumination and the transmitted light equipment.

For both polarizing devices, the analyzers in the slider (Fig. 56) are identical. They are inserted into the filter receptacles of the microscope holders (55.1). Their use is not necessary if the device for coaxial illumination is installed with a fixed analyzer between the macro-objectives and the comparison device.

When using the rotating polarizers for the oblique incident illumination, the focusing lenses must be installed on the ends of the light guides first. Then screw the rotating polarizers onto the threads provided in front of the focusing lenses. By rotating the polarizers, you can reduce disruption from points of light generated by shiny object structures.

To use the polarizers with lambda plates for setting up the transmitted light illumination, first assemble the rotating stages \varnothing 118 mm on the cross-stages (\rightarrow page 31). In the bottom of the rotating stages, there is a filter holder (57.1), into which the polarizer with the lambda plate

(Fig. 58) must be inserted. This receptacle allows 90° horizontal rotation of the polarizer with lambda plate for setting the desired polarization contrast.

For activation and deactivation of the lambda plate, the polarizer must be rotated 180°.

Fig. 56 Analyzers

Fig. 57 Rotating stage \varnothing 118 mm

1 Filter slot for polarizer

Fig. 58 Polarization filter for transmitted light

9.8 Stage Scales and Calibration Standard

During assembly of the comparison microscope, the macro-objectives have been chosen so that they match with regard to their optical properties and permissible tolerances. For identical objectives, the image scales generated thus fall within the permissible tolerances.

The stage scale (stage micrometer) 11 519 963 or the calibration standard 11 581 081 can be used to verify this. To do so, adjust the section image for identical object magnifications and check that the image scales on both sides match exactly.

Fig. 59 Stage micrometer

10. Technical Data

Stand

- Stable and torsion-free cast-iron stand with motorized height adjustment for the carrier of the comparison bridge (stroke 255 mm), motorized stages and motorized focusing
- Motorized synchronous displacement of the stages and focusing unit over the entire travel path
- Built-in power supply for all motorized functions

Focusing

- Two motorized focusing drives
- Focusing speed matched to the current magnification
- Stroke: 25 mm

Stages

- Motorized cross-stages
- Stage surface: 160 mm x 220 mm
- Glass stage plates: 80 mm x 80 mm
- Travel range: 50 mm x 50 mm
- Ergonomic operating knobs for cross and synchronized displacement
- Optional: remote control with "Smart Move" control panel
- Travel speed matched to the current magnification
- Bore hole receptacles for the adjustable holder, the rotating stages, the large sample stages or the bullet holder

Stackable rotating stages ø 118 mm

- Openings ø 50 mm
- Removable glass stage plates
- Clamping device for the stage rotation
- Stage carrier with receptacle for rotating polarizer

Inclining stages ø 75 mm

- Can be locked into place
- Inclines in all directions up to 45°
- Surface with ridged surface

Large sample stages

- Stage surface: 210 mm x 300 mm
- Metal plates for thin objects such as documents
- Can be stacked on the cross-stage

Motorized comparison bridge with integrated ergonomic tube

- For superimposed image or split image comparisons with variable dividing line
- Variable width of dividing line
- Combination of split image and superimposed image possible
- Color variation of deviation in traces possible in superimposed image observations
- Magnification calibration of right beam path possible (+/- 4%)
- Receptacles for macro-objectives 0.4x, 1x, 2x, 4x with adjustable iris diaphragms
- Rotating fasteners for the holders for oblique incident illumination
- Slots for inserting filter slides
- Distance between the optical axes: 400 mm

10. Specifications

Tube factor

1x, 1.5x with magnification changer

Field number

FOV 22

Image orientation

Portrait or landscape

Objectives

- Macro-objectives 0.4x, 1x, 2x, 4x
- Micro-objectives

Eyepieces

HC PLAN S 10x/22

		Eyepiece 10x		Camera 2/3" with C-mount 0.63x	
Objective	FWD	FOV	Total magnification	FOV	Total magnification
0.4x	60	55.00	4.00	43.60	10.10
0.6x	60	36.60	6.00	29.00	15.20
1x	60	22.00	10.00	17.44	25.40
1.5x	60	14.60	15.00	11.60	37.90
2x	60	11.00	20.00	8.70	50.60
3x	60	7.30	30.00	5.80	75.90
4x	60	5.50	40.00	4.40	101.30
6x	60	3.60	60.00	2.90	152.00

Magnifications and object fields**Motorized work table**

- Height adjustment range: 619 mm + 300 mm movement
- Lifting capacity: 2000 N
- Lifting speed: approx. 12 mm/s
- Load capacity: 200 kg
- Control system: up and down switches and foot switch
- Stage surface: 1200 mm x 560 mm
- Universal power supply

Dimensions without camera

- Height: 785 mm
(maximum; Z column extended)
- Width: 1035 mm
(oblique incident illumination with maximum extension of articulated arms)
- Depth: 530 mm
(incl. front control knobs)

Weights

- Basic equipment: 45 kg
- Comparison bridge: 15 kg
- Macroscope stand with table: 32 kg
- Motorized work table: 50 kg

General specifications

For indoor use only.

- Supply voltage: 90-250 V~
- Frequency: 50-60 Hz
- Power input: maximum 125 VA
- Fuse (power supply): F 6,3 A L 250 V -UL-
- Fuse (power pack): T 16 A L 250 V
- Ambient temperature: 10-36°C
- Relative humidity: maximum 80% to 30°C
- Overvoltage category: II
- Pollution degree: 2

11. Care of the Comparison Macroscope

Caution:

Unplug the power supply before performing cleaning and maintenance work!
Protect electrical components from moisture!

Macrosopes in warm and warm-damp climatic zones require special care in order to prevent the build up of fungus.

The comparison macroscope should be cleaned after each use, and the macroscope optics should be kept painstakingly clean.

11.1 Dust Cover

Note

To protect against dust, cover the macroscope and accessories with the dust cover after each use.

11.2 Cleaning

Caution:

Residual fiber and dust can create unwanted background fluorescence.

Cleaning Coated Parts

Dust and loose dirt particles can be removed with a soft brush or lint-free cotton cloth.

Clinging dirt can be cleaned as necessary with all commercially available water solutions, benzene or alcohol.

For cleaning coated parts, use a linen or leather cloth that is moistened with one of these substances.

Caution:

Acetone, xylene or distillations containing nitro can harm the macroscope and thus may not be used.

Test cleaning solutions of unknown composition on a less visible area of the unit first. Be sure when doing so that coated or plastic surfaces do not become matted or etched.

Cleaning the Stage

Remove light-colored spots on the stage by rubbing with paraffin oil or acid-free Vaseline.

Cleaning Glass Surfaces

Remove dust on glass surfaces with a fine, dry and grease-free hair brush, by blowing with a bellows or vacuum suction.

Cleaning Objectives

Caution:

The objective may not be unscrewed during cleaning. If damage appears on inner surfaces, the objectives must be sent to your Leica subsidiary for repair. We also advise against cleaning the inside surfaces of the eyepieces.

The front lenses of objectives are cleaned as described under "Cleaning Glass Surfaces". The upper lens is cleaned by being blown off with a pneumatic pump.

Removing Immersion Oil

Caution:

Follow safety instructions for immersion oil!

First, wipe off the immersion oil with a clean cotton cloth, and then re-wipe the surface several times with ethyl alcohol.

11.3 Handling Acids and Alkalines

For examinations using acids or other aggressive chemicals, particular caution must be taken.

Caution:

Be absolutely certain to prevent the optics and mechanical parts from coming into direct contact with these chemicals.

12. Essential Wear and Spare Parts

Order No. Material No.	Name	Used for
<u>Re</u> <u>Fuses:</u> 11 825795	F 6,3 A L 250 V -UL-	Fuse (power supply)

13. EU Declaration of Conformity

We hereby declare that the device described below, both in its basic design and construction and in the version marked by us, conforms to the relevant safety- and health- related requirements of the appropriate EU directives.

This declaration shall cease to be valid if modifications are made to the device without our approval.

Product: FS C
Model: Microscope

Identification No: 11 581 051

EU directives: Low voltage: 73/23/EEC
EMC: 89/336/EEC

Harmonised standards applied: EN 61010-1 : 2001
EN 61000-6-1 : 2001
EN 61000-6-3 : 2001

Wetzlar, 01.08.2003

Dr. L. Remer
Director of R&D and technological design
microscopy

We hereby declare that the device described below, both in its basic design and construction and in the version marked by us, conforms to the relevant safety- and health- related requirements of the appropriate EU directives.

This declaration shall cease to be valid if modifications are made to the device without our approval.

Product: FS 4000 B / FL
Model: Microscope

Identification No: 11 581 052, 11 581 069

EU directives: Low voltage: 73/23/EEC
EMC: 89/336/EEC

Harmonised standards applied: EN 61010-1 : 2001
EN 61000-6-1 : 2001
EN 61000-6-3 : 2001

Wetzlar, 01.08.2003

Dr. L. Remer
Director of R&D and technological design
microscopy

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Copyright © Leica Microsystems Wetzlar GmbH · Ernst-Leitz-Straße · 35578 Wetzlar · Germany 2002 · Tel. (0 64 41) 29-0 · Fax (0 64 41) 29-25 99 LEICA and the Leica logos are registered trademarks of Leica IR GmbH.
Leica Microsystems Wetzlar GmbH Tel. +49 (0) 64 41-290