INSTRUCTION MANUAL

- *Laboratory Model # 44200*
- *Advanced Model # 44202*
- *Professional Model # 44204*
- *Professional Model # 44206*

Table of Contents

Introduction	2
Laboratory Model # 44200	3
Advanced Model # 44202	5
Professional Model # 44204	8
Professional Model # 44206	12
Care and Maintenance	16
Warranty	17

Introduction

Congratulations on your purchase of a Celestron microscope. Your microscope is a precision optical instrument, made of the highest quality materials to ensure durability and long life. It is designed to give you a lifetime of pleasure with a minimal amount of maintenance.

This instruction manual covers four microscope models. Please be sure to locate your specific model of microscope in order to ensure you read the correct information.

Before attempting to use your microscope, please read through the instructions to familiarize yourself with the functions and operations to maximize your enjoyment and usage. See the microscope diagrams to locate the parts discussed in this manual.

The microscopes described in this manual provide powers from 10 up to 67x. These microscope types are ideally suited for larger specimens such as insects or insect parts, animal parts, rocks, fossils, PC boards and other materials, textiles, plant and flower parts, stamps and coins, and other items. At the same time you can view specimen slides.

The final section provides simple care and maintenance tips for you to follow to ensure that your microscope provides you with years of quality performance, usage, and enjoyment.

Laboratory Microscope – Model # 44200


Figure 1

Standard Accessories Included with your Microscope

<ul style="list-style-type: none"> • 10x Huygens Eyepiece – 20 power • 20x Huygens Eyepiece – 40 power 	<ul style="list-style-type: none"> • 2x Objective Lens • Dust Cover & 5 Prepared Slides
--	---

Specifications

Model # 44200	Specifications
Head	45° Inclined Monocular Head; 360° Rotatable
Stage	60mm Dia. black/white stage plate with slide holder clips
Focuser	Rack and Pinion
Objective	2x
Eyepiece	10x Wide Field with 20mm Field of View
	20x Wide Field with 10mm Field of View
Working Distance	Full range 108mm; per each pillar adjustment – 40mm

Assembling Your Microscope

1. Take the Styrofoam container out of the cardboard box.
2. Remove the tape from the Styrofoam container holding the two sections together.
3. Carefully remove the microscope and other parts and set them on a table, desk, or other flat surface.
4. Remove the plastic bag covering the microscope.
5. Install the stage plate (7) by putting it in the circular opening in the base (6).
6. Make sure the head (1) is parallel to the base (6). If it is not, then hold the head (1) with one hand and with the other hand loosen the pillar stand locking knob (4) by turning it counterclockwise, and then rotate the head until it is perpendicular to the base.
7. The objective lens housing (9) should now be centered directly over the stage plate (7). If it is not, adjust it carefully like number 6 above.
8. Remove the plastic plug from the bottom of the objective lens (9).

You are now ready to use your microscope!

Microscope Operation

Viewing a Specimen

To view a slide specimen:

1. Carefully place a specimen slide under the specimen slide holder clips (8) and center the slide directly in the middle of the stage plate (7).
2. Illuminate the specimen slide by using a nearby desk lamp, table lamp, bright overhead light, or use natural light if located near a window.
3. To focus on your specimen slide, turn the focus knob (5) until the specimen comes into focus. If there is not enough travel in the focuser mechanism to reach focus, it may be necessary to loosen the pillar stand locking knob (4) and raise or lower the entire head assembly (1) until the image comes into sharp focus.

Tip: Experiment with different lighting for different specimens. When using downward lighting, the white side of the stage may wash out fine detail whereas the black side will increase contrast but may make the specimen less bright.

To view larger specimens:

1. For larger specimens such as insects or insect parts, animal parts, rocks, fossils, PC boards, other materials, plant and flower parts, and other larger items --- you will have to make adjustments to the working space to obtain focus.
2. Make sure the binocular head (1) is high enough for the specimen to fit on the stage plate. You may have to raise the head (1) by loosening the pillar stand locking knob (4) move the head up the pillar stand (3) to the desired position and then tighten the pillar stand locking knob (4).
3. Illuminate and focus the large specimen similar to items 2 and 3 in the above section.

Changing Eyepieces

Your microscope is supplied with two eyepieces with the 10x eyepiece installed standard. To change to the 20x eyepiece:

1. You will have to remove the small Phillips head screw (1.4mm size) from the eyepiece tube (10).
2. Pull the existing eyepiece (11) upward and out of the eyepiece tube (10).
3. Put the other eyepiece (11) into the eyepiece tube (10) and reinstall the small Phillips head screw.

Advanced Microscope – Model # 44202


Figure 2

Standard Accessories with your Microscope

<ul style="list-style-type: none"> • Two 10x Huygen Eyepieces • 2x Objective Lens (20 power) • 4x Objective Lens (40 power) • Top & Bottom Illumination 	<ul style="list-style-type: none"> • Black/White Stage Plate • Frosted/Opaque Stage Plate • Rubber Eyecups • Dust Cover & 5 Prepared Slides
---	---

Specifications

Model # 44202	Specifications
Head	45° Inclined Binocular Head; 360° Rotatable
Stage	95mm Dia. stage
Focuser	Rack and Pinion
Objective	Achromatic 2x and 4x
Eyepiece	10x Wide Field with 20mm Field of View
Illuminator	Top & Bottom built-in electric lighting Halogen bulbs – 12volt-15watt
Interpupillary Range	55mm to 75mm
Power	115-volt model or 230-volt model
Working Distance	Full range 110mm; per each pillar adjustment-35mm

Assembling Your Microscope

1. Take the Styrofoam container out of the cardboard box.
2. Remove the tape from the Styrofoam container holding the two sections together.
3. Carefully remove the microscope and other parts and set them on a table, desk, or other flat surface.
4. Remove the plastic bag covering the microscope.
5. Install the frosted/opaque stage plate (7) by putting it in the circular opening in the base (8).
6. Make sure the head (3) is parallel to the base (8). If it is not, then hold the head (3) with one hand and with the other hand loosen the pillar stand locking knob (13) by turning it counterclockwise, and then rotate the head until it is perpendicular to the base.
7. The objective lens housing (5) should now be centered directly over the stage plate (7). If it is not, adjust it carefully like number 6 above.
8. Remove the plastic plug from the bottom of the objective lens (5).
9. Install the rubber eyecups over the eyepieces which will help block out extraneously light to make your viewing more enjoyable if you do not wear eyeglasses. Push the rubber eyecups over the eyepieces (1) making sure the rubber lips are facing up. Then rotate the eyepieces so that the rubber lips are facing out to the sides.

You are now ready to use your microscope!

Microscope Operation

Viewing a Specimen

To view a slide specimen without using the electric illuminator:

1. Carefully place a specimen slide under the specimen slide holder clips (6) and center the slide directly in the middle of the stage plate (7). For specimen slides you should use the black/white stage plate.
2. Illuminate the specimen slide by using a nearby desk lamp, table lamp, bright overhead light, or use natural light if located near a window.
3. To focus on your specimen slide, turn the focus knob (12) until the specimen comes into focus. If there is not enough travel in the focuser mechanism to reach focus, it may be necessary to loosen the pillar stand locking knob (13) and raise or lower the entire head assembly (3) until the image comes into sharp focus.

Tip: Experiment with different lighting for different specimens. When using downward lighting, the white side of the stage may wash out fine detail whereas the black side will increase contrast but may make the specimen less bright.

To view larger specimens without using the electric illuminator:

1. For larger specimens such as insects or insect parts, animal parts, rocks, fossils, PC boards, other materials, plant and flower parts, and other larger items --- you will have to make adjustments to the working space to obtain focus.
2. Make sure the binocular head (3) is high enough for the specimen to fit on the stage plate. You may have to raise the head (3) by loosening the pillar stand locking knob (13) move the head up the pillar stand (14) to the desired position and then tighten the pillar stand locking knob (13).
3. Illuminate and focus the large specimen similar to items 2 and 3 in the section above.

To provide higher quality viewing of specimens you will want to use the electric illuminator. Your microscope has a built-in, three setting illuminator. To use the illuminator:

1. Plug the power cord (10) into a wall outlet. Make sure you have the proper model for your power requirement – either the 115-volt model or the 230-volt model.
2. Turn the illuminator switch (9) to the “I” setting to turn on the bottom illuminator light.
3. Turn the illuminator switch (9) to the “II” setting to turn on the top illuminator light.
4. Turn the illuminator switch (9) to the “III” setting to turn on both the top and bottom illuminator lights at the same time. Note that the bottom illuminator light is only beneficial when using the frosted/opaque stage plate.
5. For many specimens, you should experiment with the different illumination options to get one that gives you the best results.

Focusing

To achieve the best focus for your microscope, not only can you make adjustments with the focus knobs, but you can also adjust the left diopter ring (2) of the left eyepiece (1) as well as the interpupillary distance between the eyepieces:

1. Rotate the knurled ring on the objective lens housing (5) to the 2x position (as seen looking straight at the microscope). It is best to begin viewing a specimen at low power, then switch to high power when needed.
2. Turn the focus knob (12) until the specimen comes into focus. If there is not enough travel in the focuser mechanism to reach focus, it may be necessary to loosen the pillar stand locking knob (13) and raise or lower the entire binocular head assembly (3) until the image becomes sharp.
3. While looking through the eyepieces, adjust the interpupillary distance by rotating both eyepieces and the eyepiece holding mechanisms up and down until the whole field of view can be observed through both eyes at the same time. The reason this is necessary is that people have varying distances between their eyes.
4. Close your left eye and focus the image until sharp looking through your right eye. Then, close your right eye and turn the left eyepiece diopter ring (2) until the image is in sharp focus through the eyepiece. Now the interpupillary distance is set and you will only focus with the focus knob (12).
5. If higher power is desired, rotate the objective lens housing (5) to the 4x position (as seen looking straight at the microscope) and refocus until the specimen becomes sharp.

Rotating the Binocular Head

The head of your microscope can be rotated 360°. Therefore, you can view from any position. To change positions of the head, loosen the binocular head tension screw (4) and rotate the head (3) to the position desired and then tighten the binocular head tension screw (4).

Replacing Illuminator Bulbs

Top Illuminator Bulb – to replace the bulb:

1. Remove the binocular head mechanism ((3) from the pillar stand (14) by loosening the pillar stand locking knob (13) but first remove the pillar stand safety screw (15).
2. Lay the binocular head (3) on its side so you can easily change the bulb in the top illuminator (11).
3. Gently remove the illuminator bulb from the bulb holder of the top illuminator – see Figure 2a.
4. Replace the bulb with a new one of the same type – halogen 12volt, 15watt.
5. Put the binocular head (3) back onto the pillar stand (14) and tighten the pillar stand locking knob (13) in the position desired.
6. Replace the pillar stand safety screw (15).

Bottom Illuminator Bulb – to replace the bulb:

1. Remove the stage plate (7) so that you can easily see the bulb holder inside the base (8).
2. Raise the head (3) up by turning the focus knob (12). If this is not high enough, then loosen the pillar stand locking knob (13) and raise the head to a higher position (hold firmly when doing this) and then tighten the pillar stand locking knob (13).
3. Gently remove the illuminator bulb from the bulb holder of the bottom illuminator – see Figure 2b.
4. Replace the bulb with a new one of the same type – halogen 12volt, 15watt.
5. Put the stage plate (7) back in the base (8).

Figure 2a

Figure 2b

Professional Model # 44204


Figure 3

Standard Accessories with your Microscope

<ul style="list-style-type: none"> • Two 10x Eyepieces • Two 15x Eyepieces • 2x Objective Lens • 4x Objective Lens • Top & Bottom Illumination 	<ul style="list-style-type: none"> • Black/White Stage Plate • Frosted/Opaque Stage Plate • Four Rubber Eyecups • Dust Cover & 5 Prepared Slides • Extra Halogen Bulb 12V-15W
---	--

Specifications

Model # 44204	Specifications
Head	45° Inclined Binocular Head; 360° Rotatable
Stage	95mm Dia. stage
Focuser	Rack and Pinion
Objective	Achromatic 2x and 4x
Eyepiece	10x Extra Wide Field with 20mm Field of View 15x Extra Wide Field with 15mm Field of View
Illuminator	Top & Bottom Lighting -- Top is Halogen bulb w/12volt-15watt Bottom is fluorescent 8watt with 115-volt or 230-volt model
Interpupillary Range	55mm to 75mm
Power	115volt model or 230volt model
Working Distance	Full range 180mm; per each pillar adjustment-43mm

Magnification Table

Objective Lens	2x	4x
10x Eyepieces	20 Power	40 Power
15x Eyepieces	30 Power	60 Power

Assembling Your Microscope

1. Take the Styrofoam container out of the cardboard box.
2. Remove the tape from the Styrofoam container holding the two sections together.
3. Please note that one set of eyepieces will be on the top outer (UP) Styrofoam piece with tape over them. So make sure you remove these.
4. Carefully remove the microscope and other parts and set them on a table, desk, or other flat surface (See Figure 3a)
5. Remove the plastic bags from the stage plate (which is installed on the base) and from the binocular head (14).
6. Put the frosted stage plate (10) back into the circular hole in the base (7).
7. Remove the plastic plug cap from the bottom of the objective lens housing (12) by turning counterclockwise.
8. Make sure the pillar locking knob (5) is very tight before proceeding.
9. Back off the head lock screw (13) from the circular opening in the microscope bracket arm (see Figure 3a) so it does not protrude into the circular opening area.
10. Place the objective lens housing (12) down into the circular opening in the bracket arm and make sure it is all the way down where the head (14) rests upon the top of the bracket arm.
11. Then, tighten the head lock screw (13).
12. Make sure the bracket arm and objective lens housing (12) are parallel to the base (7).
13. Rotate the knurled ring at the top of the objective lens housing (12) counterclockwise all the way until it stops so that the 2x objective lens is facing you when looking from the front of the microscope.
14. Remove the plastic plug caps from the eyepiece tubes.
15. Install the two 10x eyepieces into the eyepiece tubes by pushing down on the eyepiece until it stops. The eyepieces are held in place by friction.
16. If you do not wear eyeglasses you install the rubber eyecups (the ones with the largest circular opening) with the lips up and facing outward. The other set of rubber eyecups are for use with the 15x eyepieces. The rubber eyecups will stop extraneous light.
17. Install the top illuminator by lowering it in the square opening in the bracket arm. Start by backing up the illuminator lock screw (4) so that it does not protrude into the opening in the bracket arm.
18. On the body of the illuminator is a threaded insert. Hold this up adjacent to the lock screw (4) and thread the lock screw (4) into the threaded insert of the illuminator until it is tight.
19. Insert the cable of the top illuminator (2) into the socket at the top of the pillar stand (3).

You are now ready to use your microscope!

Figure 3a

Microscope Operation

Viewing a Specimen

To view a slide specimen without using the electric illuminator:

1. Carefully place a specimen slide under the specimen holder clips (11) and center the slide directly in the middle of the stage plate (10). For specimen slides, you should use the black/white stage plate.
2. Illuminate the specimen slide by using a nearby desk lamp, table lamp, bright overhead light, or use natural light if located near a window.
3. To focus on your specimen slide, turn the focus knob (6) until the specimen comes into focus. If there is not enough travel in the focuser mechanism to reach focus, it may be necessary to loosen the pillar stand locking knob (5) and raise or lower the entire head assembly (14) until the image comes into sharp focus.

Tip: Experiment with different lighting for different specimens. When using downward lighting, the white side of the stage may wash out fine detail whereas the black side will increase contrast but may make the specimen less bright.

To view larger specimens without using the electric illuminator:

1. For larger specimens such as insects or insect parts, animal parts, rocks, fossils, PC boards, other materials, plant and flower parts, and other larger items --- you will have to make adjustments to the working space to obtain focus.
2. Make sure the binocular head (14) is high enough for the specimen to fit on the stage plate (10). You may have to raise the head (14) by loosening the pillar stand locking knob (5) and move the head up the pillar stand (3) to the desired position and then tighten the pillar stand locking knob (5).
3. Illuminate and focus the large specimen similar to items 2 and 3 in the section above.

To provide higher quality viewing of specimens you will want to use the electric illuminator. Your microscope has a top illuminator and a bottom illuminator. To use the illuminators:

1. Plug one end of the power cord into the socket on the back of the microscope and the plug prong end of the power cord into a wall outlet – see Figure 3b below. Make sure you have the proper model for your power requirement – either the 115-volt model or the 230-volt model.
2. Turn the switch (located at the back of the microscope) to the on position – see Figure 3b.
3. Turn the top illuminator knob (9) on by rotating it clockwise and increase the brightness by continuing to rotate the knob to the desired intensity.
4. Turn the bottom illuminator knob (to the right of the top illuminator knob) on by rotating it clockwise. The bottom illuminator does not have an adjustment for brightness. The bottom illuminator can only be useful when using the frosted/opaque stage plate.
5. For many specimens, you should experiment with the different illumination options to get one that gives you the best results.

Figure 3b

Rotating the Binocular Head

The head of your microscope can be rotated 360°. Therefore, you can view from any position. To change positions of the head, loosen the binocular head lock screw (13) and rotate the head (14) to the position desired and then tighten the binocular head tension screw (13).

Focusing

To achieve the best focus for your microscope, not only can you make adjustments with the focus knobs, but you can also adjust the left diopter ring (1) of the left eyepiece as well as the interpupillary distance between the eyepieces:

1. Rotate the knurled ring on the objective lens housing (12) to the 2x position (as seen looking straight at the microscope) all the way until it stops. It is best to begin viewing a specimen at low power, then switch to high power when needed.
2. Turn the focus knob (6) until the specimen comes into focus. If there is not enough travel in the focuser mechanism to reach focus, it may be necessary to loosen the pillar stand locking knob (5) and raise or lower the entire binocular head assembly (14) until the image becomes sharp.
3. While looking through the eyepieces, adjust the interpupillary distance by rotating both eyepieces and the eyepiece holding mechanisms up and down until the whole field of view can be observed through both eyes at the same time. The reason this is necessary is that people have varying distances between their eyes.
4. Close your left eye and focus with the right eyepiece until the image is sharp. Then, close your right eye and turn the left eyepiece diopter ring (1) until the image is in sharp focus through the left eyepiece. Now the interpupillary distance is set and you can focus only with the focus knob (6).
5. If higher power is desired, rotate the objective lens housing (12) to the 4x position (as seen looking straight at the microscope) and rotate as far as you can go and then refocus until the specimen becomes sharp.
6. You can also increase power by changing to the 15x eyepieces. The magnification table will show you the power combinations available. To change the eyepieces, see items 15 and 16 in Assembling Your Telescope earlier.
7. You can adjust the tension of the focus knobs (6) to suit your needs. Looking straight at the microscope, hold the left focus knob tight and rotate the right focus knob clockwise to make the tension stiffer. To make the tension looser, you would rotate the right focus knob counterclockwise.

Replacing Illuminator Bulbs

Top Illuminator Bulb – to replace the bulb:

1. Remove the top illuminator from the bracket arm by unscrewing the illuminator lock screw (4).
2. Hold the top piece of the illuminator (the piece with the cable attached to it) with one hand. With the other hand, pull the body away from the top piece by twisting slightly and the bulb can be seen in its socket holder – see Figure 3c.
3. Gently remove the bulb from its socket holder.
4. Replace the bulb with a new one of the same type – halogen 12volt, 15watt.
5. Push the illuminator body up onto the top piece and twist slightly and it will lock in place.
6. Install the illuminator into the bracket arm by screwing the illuminator lock screw (4) into the illuminator threaded insert.

Bottom Illuminator Bulb – to replace the bulb:

1. Remove the stage plate (10) so that you can easily see the fluorescent bulb inside the base (7) – see Figure 3d.
2. Raise the head (14) up so you have more room to work in.
3. Gently remove the illuminator bulb from the bulb holder by pulling outward and moving it left to right gently.
4. Replace the bulb with a new one of the same type – fluorescent 8watt and either 115-volt or 230-volt depending on the model you own. Put the bulb in gently so the pins of the bulb insert properly into the bulb socket holder.


Figure 3c

Figure 3d

Replacing Fuses

Two fuses are located on the back of the microscope (see Figure 3b). They are T1A fuses. If the power goes out or will not come on, it is possible one or both of the fuses needs to be replaced. Make sure the power switch is turned off and then use a flat head screwdriver to remove the fuses. If they are blown out, replace them with new fuses and then install the fuse box back in place.

Professional Model # 44206


Figure 4

Standard Accessories with your Microscope

<ul style="list-style-type: none"> • Two 15x Eyepieces • Zoom Objective Lens • Top & Bottom Illumination • Black/White Stage Plate 	<ul style="list-style-type: none"> • Frosted/Opaque Stage Plate • Rubber Eyecups • Dust Cover & 5 Prepared Slides • Extra Halogen Bulb 12V 15W
--	--

Specifications

Model # 44206	Specifications
Head	45° Inclined Binocular Head; 360° Rotatable
Stage	95mm Dia. stage
Focuser	Rack and Pinion
Objective	Zoom 0.7x to 4.5x
Eyepiece	15x Extra Wide Field with 15mm Field of View
Illuminator	Top & Bottom Lighting -- Top is Halogen bulb w/12volt-15watt
	Bottom is florescent 8watt with 115volt or 230volt model
Interpupillary Range	55mm to 75mm
Power	115-volt model or 230-volt model
Working Distance	Full range 200mm; per each pillar adjustment-43mm

Assembling Your Microscope

1. Take the Styrofoam container out of the cardboard box.
2. Remove the tape from the Styrofoam container holding the two sections together.
3. Carefully remove the microscope and other parts and set them on a table, desk, or other flat surface (See Figure 4a)
4. Remove the plastic bags from the stage plate (which is installed on the base) and from the binocular head (15).
5. Put the frosted stage plate (10) back into the circular hole in the base (7).
6. Remove the plastic plug cap from the bottom of the objective lens housing (12) by turning counterclockwise.
7. Make sure the pillar locking knob (5) is very tight before proceeding.
8. Back off the head lock screw (13) from the circular opening in the microscope bracket arm (see Figure 4a) so it does not protrude into the circular opening area.
9. Place the objective lens housing (12) down into the circular opening in the bracket arm and make sure it is all the way down where the head (15) rests upon the top of the bracket arm.
10. Then, tighten the head lock screw (13).
11. Make sure the bracket arm and objective lens housing (12) are parallel to the base (7).
12. Remove the plastic plug caps from the eyepiece tubes.
13. Install the two eyepieces (16) into the eyepiece tubes by pushing down on the eyepiece until it stops. The eyepieces are held in place by friction.
14. If you do not wear eyeglasses you install the rubber eyecups with the lips up and facing outward. The rubber eyecups will help stop extraneous light.
15. Rotate the zoom power knob (14) until the dial reads 0.7x.
16. Install the top illuminator by lowering it in the square opening in the bracket arm. Start by backing up the illuminator lock screw (4) so that it does not protrude into the opening in the bracket arm.
17. On the body of the illuminator is a threaded insert. Hold this up adjacent to the lock screw (4) and thread the lock screw (4) into the threaded insert of the illuminator until it is tight.
18. Insert the cable of the top illuminator (2) into the socket at the top of the pillar stand (3).

You are now ready to use your microscope!

Bracket Arm

Figure 4a

Microscope Operation

Viewing a Specimen

To view a slide specimen without using the electric illuminator:

1. Carefully place a specimen slide under the specimen holder clips (11) and center the slide directly in the middle of the stage plate (10). For specimen slides, you should use the black/white stage plate.
2. Illuminate the specimen slide by using a nearby desk lamp, table lamp, bright overhead light, or use natural light if located near a window.
3. To focus on your specimen slide, turn the focus knob (6) until the specimen comes into focus. If there is not enough travel in the focuser mechanism to reach focus, it may be necessary to loosen the pillar stand locking knob (5) and raise or lower the entire head assembly (15) until the image comes into sharp focus.

Tip: Experiment with different lighting for different specimens. When using downward lighting, the white side of the stage may wash out fine detail whereas the black side will increase contrast but may make the specimen less bright.

To view larger specimens without using the electric illuminator:

1. For larger specimens such as insects or insect parts, animal parts, rocks, fossils, PC boards, other materials, plant and flower parts, and other larger items --- you will have to make adjustments to the working space to obtain focus.
2. Make sure the binocular head (15) is high enough for the specimen to fit on the stage plate (10). You may have to raise the head (15) by loosening the pillar stand locking knob (5) and move the head up the pillar stand (3) to the desired position and then tighten the pillar stand locking knob (5).
3. Illuminate and focus the large specimen similar to items 2 and 3 in the section above.

To provide higher quality viewing of specimens you will want to use the electric illuminator. Your microscope has a top illuminator and a bottom illuminator. To use the illuminators:

1. Plug one end of the power cord into the socket on the back of the microscope and the plug prong end of the power cord into a wall outlet – see Figure 4b below. Make sure you have the proper model for your power requirement – either the 115-volt model or the 230-volt model.
2. Turn the switch (located at the back of the microscope) to the on position – see Figure 4b.
3. Turn the top illuminator knob (9) on by rotating it clockwise and increase the brightness by continuing to rotate the knob to the desired intensity.
4. Turn the bottom illuminator knob (to the right of the top illuminator knob) on by rotating it clockwise. The bottom illuminator does not have an adjustment for brightness. The bottom illuminator can only be useful when using the frosted/opaque stage plate.
5. For many specimens, you should experiment with the different illumination options to get one that gives you the best results.

Figure 4b

Rotating the Binocular Head

The head of your microscope can be rotated 360°. Therefore, you can view from any position. To change positions of the head, loosen the binocular head lock screw (13) and rotate the head (15) to the position desired and then tighten the binocular head tension screw (13).

Focusing and Changing Power

To achieve the best focus for your microscope, not only can you make adjustments with the focus knobs, but you can also adjust the left diopter ring (1) of the left eyepiece as well as the interpupillary distance between the eyepieces: Start with the zoom power knob (14) at the 0.7x position on the dial which is the lowest power (10x). It is usually best to begin viewing a specimen at low power, then zoom to high powers. The maximum power is 67x (dial at 4.5x).

1. Turn the focus knob (6) until the specimen comes into focus. If there is not enough travel in the focuser mechanism to reach focus, it may be necessary to loosen the pillar stand locking knob (5) and raise or lower the entire binocular head assembly (15) until the image becomes sharp.
2. While looking through the eyepieces, adjust the interpupillary distance by rotating both eyepieces and the eyepiece holding up and down until the whole field of view can be observed through both eyes at the same time. The reason this is necessary is that people have varying distances between their eyes.
3. Close your left eye and focus with the right eyepiece until the image is sharp. Then, close your right eye and turn the left eyepiece diopter ring (1) until the image is in sharp focus through the left eyepiece. Now the interpupillary distance is set and you can focus only with the focus knob (6).
4. If higher power is desired, rotate the zoom power knob (14) to the desired power. As you increase power, you will have to refocus the object.
5. You can adjust the tension of the focus knobs (6) to suit your needs. Looking straight at the microscope, hold the left focus knob tight and rotate the right focus knob clockwise to make the tension stiffer. To make the tension looser, you would rotate the right focus knob counterclockwise.

Replacing Illuminator Bulbs

Top Illuminator Bulb – to replace the bulb:

1. Remove the top illuminator from the bracket arm by unscrewing the illuminator lock screw (4).
2. Hold the top piece of the illuminator (the piece with the cable attached to it) with one hand. With the other hand, pull the body away from the top piece by twisting slightly and the bulb can be seen in its socket holder – see Figure 4c.
3. Gently remove the bulb from its socket holder.
4. Replace the bulb with a new one of the same type – halogen 12volt, 15watt.
5. Push the illuminator body up onto the top piece and twist slightly and it will lock in place.
6. Install the illuminator into the bracket arm by screwing the illuminator lock screw (4) into the illuminator threaded insert.

Bottom Illuminator Bulb – to replace the bulb:

1. Remove the stage plate (10) so that you can easily see the fluorescent bulb inside the base (7) – see Figure 4d.
2. Raise the head (15) up so you have more room to work in.
3. Gently remove the illuminator bulb from the bulb holder by pulling outward and moving it left to right gently.
4. Replace the bulb with a new one of the same type – fluorescent 8watt and either 115-volt or 230-volt depending on the model you have. Put the bulb in gently so the pins of the bulb insert properly into the bulb socket holder.


Figure 4c

Figure 4d

Replacing Fuses

Two fuses are located on the back of the microscope (see Figure 4b). They are T1A fuses. If the power goes out or will not come on, it is possible one or both of the fuses needs to be replaced. Make sure the power switch is turned off and then use a flat head screwdriver to remove the fuses. If they are blown out, replace them with new fuses and then install the fuse box back in place.

Care and Maintenance

Your Celestron microscope is a precision optical instrument and should be treated with care at all times. Follow these care and maintenance suggestions and your microscope will need very little maintenance throughout its lifetime.

- When you are done using your microscope, remove any specimens left on the stage.
- Turn off any electric illuminator knobs or switches.
- Unplug any power cords being used.
- Always place the dust cover over the microscope when not in use or when being stored.
- Store the microscope in a dry and clean place.
- Be very careful if using your microscope in direct sun light to prevent damage to the microscope or your eyes.
- When moving your telescope, carry it by the arm bracket with one hand and your other hand underneath the base and not by the focuser knobs, eyepiece housing, etc.
- Clean the outside surfaces (metal and plastics) with a moist cloth.
- Always unplug any cords before cleaning.
- Never clean optical surfaces with cloth or paper towels as they can scratch optical surfaces easily.
- Blow off dust with a camel's hair brush or an air blower from optical surfaces.
- To clean fingerprints off of optical surfaces, use a lens cleaning agent and lens tissue available at most photo outlets and when cleaning do not rub in circles as this may cause streaks or scratches to occur.
- Never disassemble or clean internal optical surfaces. This should be done by qualified technicians at the factory or other authorized repair facilities.
- When handling glass specimen slides, be careful as the edges can be sharp.

Celestron Two Year Warranty

- A. Celestron warrants this microscope to be free from defects in materials and workmanship for two years. Celestron will repair or replace such product or part thereof which, upon inspection by Celestron, is found to be defective in materials or workmanship. As a condition to the obligation of Celestron to repair or replace such product, the product must be returned to Celestron together with proof-of-purchase satisfactory to Celestron.
- B. The Proper Return Authorization Number must be obtained from Celestron in advance of return. Call Celestron at (310) 328-9560 to receive the number to be displayed on the outside of your shipping container.

All returns must be accompanied by a written statement setting forth the name, address, and daytime telephone number of the owner, together with a brief description of any claimed defects. Parts or product for which replacement is made shall become the property of Celestron.

The customer shall be responsible for all costs of transportation and insurance, both to and from the factory of Celestron, and shall be required to prepay such costs.

Celestron shall use reasonable efforts to repair or replace any microscope covered by this warranty within thirty days of receipt. In the event repair or replacement shall require more than thirty days, Celestron shall notify the customer accordingly. Celestron reserves the right to replace any product which has been discontinued from its product line with a new product of comparable value and function.

This warranty shall be void and of no force of effect in the event a covered product has been modified in design or function, or subjected to abuse, misuse, mishandling or unauthorized repair. Further, product malfunction or deterioration due to normal wear is not covered by this warranty.

CELESTRON DISCLAIMS ANY WARRANTIES, EXPRESS OR IMPLIED, WHETHER OF MERCHANTABILITY OF FITNESS FOR A PARTICULAR USE, EXCEPT AS EXPRESSLY SET FORTH HEREIN.

THE SOLE OBLIGATION OF CELESTRON UNDER THIS LIMITED WARRANTY SHALL BE TO REPAIR OR REPLACE THE COVERED PRODUCT, IN ACCORDANCE WITH THE TERMS SET FORTH HEREIN. CELESTRON EXPRESSLY DISCLAIMS ANY LOST PROFITS, GENERAL, SPECIAL, INDIRECT OR CONSEQUENTIAL DAMAGES WHICH MAY RESULT FROM BREACH OF ANY WARRANTY, OR ARISING OUT OF THE USE OR INABILITY TO USE ANY CELESTRON PRODUCT. ANY WARRANTIES WHICH ARE IMPLIED AND WHICH CANNOT BE DISCLAIMED SHALL BE LIMITED IN DURATION TO A TERM OF TWO YEARS FROM THE DATE OF ORIGINAL RETAIL PURCHASE.

Some states do not allow the exclusion or limitation of incidental or consequential damages or limitation on how long an implied warranty lasts, so the above limitations and exclusions may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.


Celestron reserves the right to modify or discontinue, without prior notice to you, any model or style microscope.

If warranty problems arise, or if you need assistance in using your microscope contact:

Celestron
Technical Support Department
2835 Columbia Street
Torrance, CA 90503 U.S.A.
Tel. (310) 328-9560
Fax. (310) 212-5835
www.celestron.com
Monday-Friday 8AM-4PM PST

This warranty supersedes all other product warranties.

NOTE: This warranty is valid to U.S.A. and Canadian customers who have purchased this product from an Authorized Celestron Dealer in the U.S.A. or Canada. Warranty outside the U.S.A. and Canada is valid only to customers who purchased from a Celestron Distributor or Authorized Celestron Dealer in the specific country and please contact them for any warranty service.

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All rights reserved
Torrance, Ca 90503 U.S.A.