# 7710

Bedienungsanleitung S. 2  
Operating Instructions P.28  
Mode d'emploi P.54

**Inhalt**

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**Content**

CE marking	29	Terminal settings	37
Classification	29	Function mode	38
Safety instructions	29	Keypad	39
Intended use	29	Display	40
Use of additional devices	30	Symbols in the symbol bar	40
Applicator	30	Symbol for standstill	41
Maintenance and service	30	Getting Started	42
Cleaning	30	Power-on	42
Disposal	31	Power-off	43
Note	31	Zero	43
Warranty - Liability	31	x 10 display	43
Explanation of signs and symbols	32	Alphanumeric input	44
Electronic nameplate and calibration counter	32	Weighing and taring	46
Technical data	33	Weighing with manual tare	46
Additional technical information	33	Weighing with manual tare input	47
Installing the scale	33	Dialysis function	47
Operating concept	33	Weighing with a fixed tare value memory	48
Setting mode	34	Error message when taring	49
Keypad in setting mode	34	Possible displays with the Info key	49
Calling up setting mode	35	BMI function	50
Input in setting mode	35	Printing	51
Settings in setting mode	35	Error message (Cause/Remedy/Information)	52
Program settings	36		

Thank you for choosing this Soehnle Professional product. All the features of this product were designed to the state of the art and are optimised for simple and straightforward use. If you have any queries or experience any problems with the device that are not addressed in the operating instructions please contact your Soehnle Professional service partner or visit us on the Internet at [www.soehnle-professional.com](http://www.soehnle-professional.com).

## CE marking

The product carries the CE mark in accordance with the following directives:


89/336/EEC Electromagnetic compatibility

73/23/EEC Low Voltage Directive

90/384/EEC Nonautomatic Weighing Instruments

93/42/EEC Medical Device Directive

(Classification: Class I with measuring function)

## Classification

Class I medical device with measuring function.

Electrical Safety Class II (totally insulated equipment, no earth terminal).

No applicator.

Protection against water to IEC529: IPX2.

## Safety instructions

Before putting the appliance into use, please read with care the information given in the operating instructions. They contain important instructions for installation, proper use and maintenance of the appliance. The manufacturer shall not be liable for damages arising out of failure to heed the following instructions: When using electrical components under increased safety requirements, always comply with the appropriate regulations.

Never perform work on the appliance while the power is switched on.

Improper installation will render the warranty null and void. Ensure the voltage marked on the power supply unit matches your mains power supply. The appliances are designed for use within buildings. Observe the permissible ambient temperatures for use (Technical Information). The device meets the requirements for electromagnetic compatibility. Do not exceed the maximum values specified in the applicable standards. If you have any problems, contact your local Soehnle Professional service partner.

## Intended use

The appliance is designed for weighing people who are lying down. The appliance may only be used for this purpose. Any other use is prohibited. Messages concerning faults that could endanger the patient or the operating staff, as well as faults leading to the invalidation of measuring results are to be communicated in writing to the manufacturer's medical device adviser (see Address). These operating instructions, as well as the operating instructions for the 3020 display terminal, are component parts of the appliance. They must be kept near the appliance. Precise compliance with these instructions is a prerequisite for the correct use and operation of the appliance.

If the appliance fails to function correctly, the appliance may be damaged. It is essential that the appliance is then taken to an authorised Soehnle Professional service station. Only genuine spare parts may be used in the event of repairs being carried out by an authorised service centre.

These original parts are described in the service documentation (Spare Parts List 470.014.032 Scale / 470.630.269 Analyser) along with their order number.

Soehnle Professional will only assume responsibility for the safety of the appliance if these instructions are observed and the appliance is used in conformity with the instructions for use. The appliance is a medicotechnical device and may only be used by people whose training ensures they know how to handle it correctly.

The user must ensure before using the appliance that it is functioning safely and is in proper working order. The operator must be trained in the use of the appliance. The appliance is not designed for use in medical rooms subject to explosion hazards.

Areas can become subject to explosion hazards through the use of flammable anaesthetics, skin cleaning agents and skin disinfectants.

## Use of additional devices

Only additional devices (printer, computer) that comply with EN 60601-1 may be connected to the serial data output or a suitable disconnecting device must be switched between them.

## Applicator

The scale does not have an applicator. No part is intended for direct contact with the patient. The compensating currents were measured on the weighing machine in accordance with the classification for type B applicators.

## Maintenance and service

The measuring inspection test was performed during the initial calibration process in production. Further regular measuring inspection tests (recalibration) must be carried out in accordance with the respective national regulations by the responsible calibration authorities.

Repairs may only be carried out by an authorised Soehnle Professional repair centre using genuine spare parts.

## Cleaning

Before cleaning, disconnect the appliance from the mains by removing the power plug from the wall socket. The appliance may only be cleaned with a damp cloth. On no account may water be allowed to enter into the appliance. Disinfectant may only be applied to the membrane keyboard for the display.

**The following disinfectants are allowed:**

methylated spirit; isopropanol; 2% Kohrsolin; 1% aqueous Sokrena solution; 5% Sagrotan; 5% Gigasept. Spraying the appliance and the power plug is not permissible.

## Disposal

This product must not be treated as normal waste but should be handed in to an electrical/electronic equipment recycling centre. For further information please contact your local authority, municipal waste disposal company, or the company from which you bought the product.

## Note

Interference suppression of the unit corresponds to the valid EU directive 89/336/EEC. The displayed value may be influenced by extreme electrostatic and electromagnetic interference, however, e.g. operation of a radio device or mobile phone in close proximity to the appliance. When the interference has ceased the product can again be used for its intended purpose.

You may need to power up the appliance again. If there is permanent electrostatic interference, we recommend earthing the platform.

The unit is a measuring instrument. Draught, vibrations, rapid changes in temperature and direct sunlight may have an effect on the accuracy of the scale.

The scale is protected to protection class IP 54. Avoid high levels of humidity, vapours, aggressive fluids and heavy contamination.

## Warranty - Liability

If a fault or defect is present on receipt of the supplied unit which is within Soehnle Professional's scope of responsibility, Soehnle Professional shall be entitled either to repair the fault or to supply a replacement unit at its discretion. Replaced parts shall be the property of Soehnle Professional. Should the fault repairs or replacement delivery not be successful, the statutory provisions shall be valid. The warranty period shall be 2 years, beginning on the day of purchase.


Please retain your receipt as proof of purchase. Should your scale require servicing, please contact your dealer or Soehnle Professional Customer Service.

No responsibility shall be accepted in particular for damages caused through any of the following reasons: Unsuitable or improper storage or use, incorrect installation or commissioning by the owner or third parties, natural wear and tear, changes or modifications, incorrect or negligent handling, in particular overuse, chemical, electrochemical or electrical interference or humidity, unless this is attributable to negligence on the part of Soehnle Professional.

If operating, climatic or any other influences lead to a major change in conditions or material quality, the warranty for perfect unit functioning shall be rendered null and void. If Soehnle Professional provides an individual warranty, this means that the unit supplied shall be free of faults for the length of the warranty period.

Always keep the original packaging in case you have to return the scale.

## Explanation of signs and symbols

- ↔1** Weighing ranges
- Max** Maximum load permitted for calibrated weighing
- Min** Minimum load permitted for calibrated weighing
- e=** Calibration value (resolution)
- SNR** Serial number of the scale (scale type, last digit of the year of manufacture, count number)
- CE 06** CE mark of conformity with year of attachment and simultaneous commencement of validity of calibration
- (III)** Precision class
- D06-09-08-PTB** No. of CE approval
- TYPE** Manufacturer's product number
- XXXX** e.g. 0102, official number of the "notified body" (the body that carried out the initial calibration process)
- M** Symbol for EU calibration
- 001** The calibrated counter status indicates how often the scale was calibrated.  
The stored counter status must coincide with the calibrated counter status on the adhesive label on the calibration plug.

## Electronic nameplate and calibration counter

	↔1	↔2	↔3
Max	300,0 kg	500,0 kg	
Min	2,0 kg	4,0 kg	
e=	0,1 kg	0,2 kg	
SNR	7710/06-0001 Eichz: 0001		

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The appliance has an electronic nameplate.  
To view, press the Info key.

The current calibrated counter status can be displayed in setting mode under Scales.  
Setting mode is obtained by first pressing the Print key and then simultaneously the F-key. Then use the arrow keys to select "Scales" and confirm with the Print key. Confirm the connected platform with the Print key. Confirm the calibrated counter status with the Print key.  
Exit the calibrated counter status display by pressing the ->0< key.  
Exit setting mode by pressing the ORG key.

## Technical data

Weighing range:	500 kg	
Resolution:	up to 300 kg e = 100 g / from 300 kg e = 200 g	
Minimum load:	2 kg	
Maximum bed weight:	250 kg	
Approval number:	D06-09-08-PTB (0102)	
Ambient temperature for operation:	-10°C ... +40°C, rapid changes in temperature, e.g. caused by a cold draught, can adversely affect the measurement results.	
Storage temperature:	-40°C ... +70°C	
Humidity:	20% - 85% relative humidity, not condensing	
Mains operation:	100 or 230 V / 50-60 Hz, 400 mA	
Transport conditions:	the appliance should only be transported in the original packaging. The storage temperature must be maintained; avoid shaking and shocks.	
Membrane keyboard:	19 keys, 1 function key	
Organisation data memory:	4 memories, alphanumeric input via the numerical keypad	
Display:	dot matrix 240 x 64 pixels, dimmable backlighting	
Memory capacity:	fixed value memory for 400 fixed values, valid for all programs, no-volt protected, optional alibi memory for 2 million inputs.	
I/O ports:	1 x RS 232, additional ports optional, including buses.	
Material:	programming terminal:	housing in stainless steel
	Weighing machine:	aluminium

## Additional technical information

Diagrams, lists of spare parts, and setting instructions can be requested through the manufacturer's Customer Service.

## Installing the scale

The scale is designed only for use on level floors. For installation sites in Germany, this conditions is satisfied if the floors comply with the tolerances permitted for surface-ready floors in accordance with DIN 18202, Table 3, Line 4.

For other countries the respective applicable national standards should be taken as a basis.

## Operating concept

The appliance is operated by means of the keypad. The keys are assigned to different functions in different modes.


## Setting mode

Here you can enter user-definable settings on the appliance to optimise the scale to your requirements.

Program settings	Terminal	Scales	IT	Service
General* Weighing + taring Totalising* Counting* Checking*	Version Display Keypad Date/time Organisation data User password Reset terminal	Scale 1 Scale 2* Scale 3*	Data transfer* Alibi memory* I/O port 1* I/O port 2* Print template* I/O port settings* Fieldbus* Barcode*	Only for service technicians (protected by password)

\* not relevant for this application

## Keypad in setting mode

- F ↑
**F-key:**  
 Moves the marking arrow upwards in the list or cursor to the left
- T↓
**T key:**  
 Moves the marking arrow downwards in the list or cursor to the right
- ORG
**ORG key:**  
 Exits setting mode.
- ⊙ ↵
**Print/Enter key:**  
 Calls up the marked menu item or accepts an input or change
- 0←
**Zeroing key:**  
 Returns to superordinated level (returns without accepting a change)
- C
**C key:**  
 Clears an item

## Calling up setting mode

Switch the appliance on.

After the power-on routine, the appliance starts the last activated application program automatically.

To call up setting mode first press the print key and then simultaneously the F key.

You are now in the Setting mode selection menu.

## Input in setting mode

In most cases you can select from several different settings available.

With editable values, you can type in values into the numerical and alphanumeric keypad.

## Settings in setting mode

With the parameters in setting mode you can adjust the "Weighing and taring" application program to suit your individual requirements.

To call up setting mode, first press the Print key and then simultaneously function key "F"; "Setting mode" appears in the display. The arrow points to "Programs".

### Key functions in setting mode

The Print key has two functions in setting mode:

- a) It calls up the program settings for the marked application program.
- b) It is used to confirm and save an input.

The Tare key is used to scroll down in the menu, and the F key is used to scroll up.

The Zeroing key is used to exit a menu item without change. It can be used to scroll step by step to the output menu.

The ORG key is used to exit the setting menu.

Data and input confirmed with the Print key is imported.

## Program settings

In the "Program Settings" setting mode you can adjust the application program for your scale to your specific requirements - for operation as a bed scale this is always "Weighing and Taring".

### Calling up Setting mode

First press the Print key and then simultaneously the F key.

### Editing fixed tare values

In setting mode, you will find the fixed tare values via the program settings → General → Fixed tare values.

### Selecting an existing fixed value from the list

Browse with the Tare and F keys, and then call up the required fixed value with the Print key.

### Editing an existing fixed value from the list

Mark the "Designation" parameter for the called up fixed value.

Open the input box using the Print key. Use the alphanumeric keypad to input data.

Press the Print key to accept the input. (Press the 0 key to cancel without accepting the input.)

The tare value parameter is marked.

Open the input box using the Print key. Use the alphanumeric keypad to input data.

Press the Print key to accept. (Press the Zeroing key to cancel without accepting the input.)

### Clearing an existing fixed value

Browse with the Tare and F keys. Clear the displayed fixed value with the C key.

The gap in the consecutive number scheme is closed.

### Inputting and saving a new fixed value

Press the ORG key; a new fixed value memory is created.

The memory location number is assigned by the appliance. The first figure in brackets indicates the number of free memory locations, and the second figure indicates the memory locations assigned with fixed tare values.

The "Designation" parameter is marked.

Open the input box using the Print key. Use the alphanumeric keypad to input data.

Press the Print key to accept the input. (Press the 0 key to cancel without accepting the input.)

The tare value parameter is marked.

Open the input box using the Print key.

Use the alphanumeric keypad to input data.

Press the Print key to accept. (Press the Zeroing key to cancel without accepting the input.)

### Exiting the fixed value list

Exit the list with the Zeroing key.

### Clearing all fixed tare values

Press the Print key and use the Tare key to scroll to "execute". Confirm with the Print key to clear all fixed tare values.

## Terminal settings

In "Terminal" setting mode you will find the following request and setting options.

### Version

AZG software version	Display only
AZG applications software	Display only
AZG serial number	Display only
CPU board serial number	Display only

### Display

Backlight luminosity setting	Input in percent via the numerical keypad 0% is dark, 100% is full luminosity. Factory setting: 80%
Contrast	Contrast setting, input in percent via the numerical keypad 1% is lowest contrast, 100% is full contrast. Factory setting: 50%
Language	Language: German/English/French
Decimal separator	Comma, point or none Factory setting: comma
Thousand separator	Comma, point or none Factory setting: comma

### Keypad

Keypad keylock	Individual, several or all keys on the keypad can be locked. By pressing the key that is to be locked, the status of the key can be changed from open to locked and vice versa. To accept, press the Print key for 3 seconds.
Key tone selection	"Off" or "on". Confirmation tone when a key is pressed. Factory setting: "off"
Acknowledgement tone	selection "off" or "on" Acknowledgement tone when a function has been executed. Factory setting: "off"

### Date/time

Time	Input via the numerical keypad. Press Print key to start clock count down
Date	Input via the numerical keypad
Date format	Selection European or imperial date Factory setting: European

**Organisation data**

Designation ORG A      Input designation of organisation data memory A  
Type digits or plaintext into keypad, e.g. "Article no.", max. 12 characters

Designation ORG B      Input designation of organisation data memory B  
Type digits or plaintext into keypad, max. 12 characters.

Designation ORG C      Input designation of organisation data memory C  
Type digits or plaintext into keypad, max. 12 characters.

Designation ORG D      Input designation of organisation data memory D  
Type digits or plaintext into keypad, max. 12 characters.

**User password**

Create / change      Input password with up to 8 characters via keypad  
After activating, the setting mode is only accessed by entering a password.  
If you have forgotten your password, press CLEAR to release access.

CLEAR                      Clears the password and releases access

**Reset terminal**

Selection                      do not execute / execute  
"Execute" clears all user-defined settings for the terminal and resets to the factory setting.

**Function mode**

The "F" function key is available for calling up functions.  
This assumes different functions depending upon how it is assigned.  
The current function is selected in setting mode.


Setting mode is obtained by first pressing the Print key and then simultaneously the F key.  
To establish the assignment options for the F key, in setting mode, go to Program Settings ' Weighing + Taring → F key.  
Use the F key or Tare key to select the required function and confirm with the Print key.

The following functions are available:

- Provisional tare
- Dialysis function
- Hold mode
- Second unit
- BMI function

Exit setting mode by pressing the ORG key.

## Keypad

The following keys are available for input and control:

Numeric keys to type in numerical values.  
Multiple programming of numeric keys to type in alphabetic characters.

Calls up organisation data.

Clears last number.

Info key, electronic nameplate with calibration data.  
Calls information for application programs.

Tare key / Cursor down / Cursor right

Zeroing key

Print key / Enter key


Function key / Cursor up / Cursor left

On / Off

## Display

The display shows a number of different views depending on the situation. It is divided into display panes:

### Additional display

Shows second units, weight values for counting, differences from setpoint etc. If required, it can be switched on with the "F" key according to the program settings.


### Symbol bar

Shows information about the status of the activated measuring point.

### Display panel

Shows weight value, piece number, unit characters.

## Symbols in the symbol bar

	Alibi memory capacity 100%
	Alibi memory capacity 80%
	Alibi memory capacity empty
	Keylock
Kg lb	Second unit
	Multirange/Multidivision 1
	Multirange/Multidivision 2
	Multirange/Multidivision 3
×10	Resolution X10
	Zero
HLD	Hold function active
Kg g lb	Unit character


## Symbol for standstill

After the weighing item has been positioned, the scale needs time to settle.

Until it comes to a standstill the unit character (kg, lb) remains blank.

Once it has reached a standstill and there is a stable value, this is indicated by the appearance of the unit character.

## Getting Started

Connect the 4 weighing machines (one for each wheel of the hospital bed) by the plugs to the distributor box. Check that they are arranged correctly according to the coloured marking.

Each cable is 6 m long.

When getting started care must be taken to ensure that the cable is not obstructed.


Connect the display unit to the mains.


### Caution:

Before powering on the scale, position the weighing machines under the 4 wheels of the hospital bed. Position the weighing machines under the wheels right up to the limit stop.

The weighing machines are designed to take wheels with a diameter of 150 mm!

## Power-on

The Soehnle Professional logo is displayed during the power-on routine which includes the offset and reference measurements.

The measuring points are briefly listed as they are detected.

The scale is zeroed after the power-on routine

At power-on the last used application program is reactivated.

### Zero range limits at power-on


Calibratable: zero range at power-on is -10% to +60% of the weighing range.

Non calibratable: zero range at power-on is -99% to +99% of the weighing range.

If power-on is outside zero range, the scale shows an error message.

Once the cause of the underload or overload has been eliminated, the scale goes to zero without the need to power-on again and is operative.

## Power-off


If the sensors are strained or you have called up a menu, the scale will not power off until the ON/Off key has been pressed for 3 seconds.

The scale cannot be powered off while waiting for an input for editing in setting mode. To power-off you must exit the editing area.

### OFF display

In the Service/General/OFF display setting mode you can set whether the word "OFF" is to appear in the display when the scale is switched off (medical directive).

## Zero

Zeroing corrects minor deviations from the zero point caused, for example, by soiling of the scale.

The following symbol in the symbol bar indicates zero: >0<.

Zero range calibratable and non calibratable: -1% to +3% of the weighing range.

If zeroing is not possible, the error message


"Above zero range" or "Below zero range" appears for two seconds.

### Zero tracking

Zero tracking automatically corrects minor deviations from the zero display.

## x 10 display

Indicates the weight value with an additional decimal in 10 times higher resolution.

Switch to x10 display by pressing first the Print key and then simultaneously the Zeroing key.

According to the Calibration Law, the x10 function on a calibratable scale must only be activated briefly for checking purposes.

For bed weighing, the x10 function is programmed to switch off again automatically 5 seconds after being activated.

While the ten times higher display is activated, X10 appears in the symbol bar.

When the appliance is switched off, the function is deactivated.

## Alphanumeric input

It is possible to make alphanumeric inputs with the combined numeric and alphabetic character key pad.

### Types of input box

- Display panels, read only
- Display/writing panels with numeric input option
- Display/writing panels with alphanumeric input option

### Activating input boxes

Situation-dependent by means of function keys or direct numerical input.

The Input view appears in the display. Existing values are displayed.

The description of the input box is displayed in the symbol bar.

A unit character is displayed if the input refers to a dimension.

### Input procedure

#### Numerical boxes

You can type in digits 0 to 9 using the numerical keypad.

Input is right-flush.

A cursor flashes under the last digit entered.

The C key clears the digit marked with the cursor.

The Zeroing key (cancel key) ends the input without accepting new input values and existing values are retained.

The Print key (enter key) accepts the input value and reactivates the application program.

The instruction as to whether there is to be a point, comma or no decimal separator is defined in the setting mode

"Display/Decimal Separator".

Only the preconfigured number of characters may be written into the input box; any further input will not be displayed and will not be accepted.

#### Alphanumeric boxes

You can type in digits 0 to 9 using the numerical keypad.

Input is left-flush.

For multiple assignment, upper case alphabetic characters and special characters are also assigned to the numeric keys. The special character keys ", . /" are also available.

If a numeric key is pressed and then pressed again within a second, this activates the characters assigned to the key concerned in rotation. If the key is not pressed again within a second, the character called up is accepted and the cursor moves on one place.

Purely numeric inputs can be made in sequence without waiting provided a particular digit is not entered several times in a row. A cursor flashes under the last item entered.

With the F key (arrow upwards) and the Tare key (arrow downwards) you can move the cursor to the left and to the right.

The C key clears the last character to the left of the cursor or the character beneath which the cursor is flashing.


The Zeroing key (cancel key) ends the input without accepting new input values and existing values are retained.

The Print key (enter key) accepts the input value and reactivates the application program.

Only the preconfigured number of characters may be written into the input box; any further input overwrites the last character.

## Weighing and taring

### Display view

### Function key

The "F" function key is available for controlling functions. It can be assigned various functions. The functions are changed in the setting menu.

Tip: To additionally display the weight in a second unit, in the setting menu you can assign the F key the function of activating a second weight display at the top left of the display.

Specify the second weight unit for each measuring point in setting mode under Setting mode/ Scale/Second unit. In addition, the function "Second unit" needs to be activated for the F key: setting mode: Program Settings → Weighing + Taring → F key → Second unit.

## Weighing with manual tare

Determine and tare an unknown bed weight.

Position the empty bed, including bedding, on the scale.  
Press the Tare key.

The tare value is saved and the weight display goes to zero.  
Now the patient can be laid on the bed.  
The weight is displayed with the message "Net".

Clear tare by overwriting with a new tare or by pressing the Zeroing key.

Display the tare memory by pressing the Info key, then press the Tare key. Clear individual points of the value displayed in the tare memory here using the C key.  
The Zeroing key ends the display of the tare memory without change. If you wish to change the value of the tare memory you can input a new value using the numerical keypad and import it into the tare memory as a manual tare input using the Tare key. Confirm the changes you have made with the Print key.

## Weighing with manual tare input

If the patient is already on the bed and the bed weight is known, then this can be tared using manual input.

Input the known bed weight (tare weight) using the numerical keypad.

The input appears in the display. The input is displayed in the activated weight unit.

Press the Tare key.

The input is imported into the tare memory and the tare weight deducted from the weight display.

The weight is displayed with the message "Net".

Clear tare by overwriting with a new manual tare input or by pressing the Zeroing key.

Display the tare memory by pressing the Info key, then press the Tare key.

Clear individual points of the value displayed in the tare memory here using the C key. The Zeroing key ends the display of the tare memory without change. If you wish to change the value of the tare memory you can input a new value using the numerical keypad and import it into the tare memory as a manual tare input using the Tare key.

Confirm the changes you have made with the Print key.

## Dialysis function

### Setting the mode of operation

- First press the Print key and the F key simultaneously; setting mode appears in the display.
- The arrow points to "Program Settings" . Confirm with Print key.
- Use the Tare key to scroll down to "Weighing + Taring" and confirm with the Print key.
- Confirm the "F key" display by pressing the Print key.
- Use the Tare key (downwards) or F key (upwards) to "Dialysis function" and confirm with the Print key. There will be an acoustic acknowledgement.
- Exit setting mode with the ORG key. The gross value is displayed.

### Application

Press the F key (with the load on the scale). Message: "Position tare".

Position the additional tare. Confirm with the Tare key. The additional tare is imported into the tare memory and the net value remains unchanged. Repeat the procedure for any additional tares. The weight is displayed with the message "Net". If the load is removed from the scale, it shows the tare weight with a negative sign.

Clear tare by pressing the Zeroing key.

Display the tare memory by pressing the Info key, then press the Tare key.

Clear individual points of the value displayed in the tare memory here using the C key.

The Zeroing key ends the display of the tare memory without change. If you wish to change the value of the

tare memory you can input a new value using the numerical keypad and import it into the tare memory as a manual tare input using the Tare key. Confirm the changes you have made with the Print key.

## Weighing with a fixed tare value memory

Known bed weights can be stored as the tare value in a no-volt protected memory and called up as required (400 different bed weights are possible).

Each memory location has a:

- three-digit consecutive number
- description (plaintext, 18 characters)
- tare value

### Inputting tare weights into the fixed tare value memory

- Press the Print key and function key "F" simultaneously to go to setting mode.
- The arrow points to "Program Settings". Confirm with the Print key.
- Confirm the display "General" with the Print key.
- Scroll to "Fixed tare value" using the Tare key and accept with the Print key.
- If no tare memory has yet been assigned, the first free, or otherwise the last selected fixed tare value memory is displayed.
- Press the ORG key and then the Print key: the input box for describing a new tare memory is displayed.
- Input the description with the alphanumeric keypad and confirm with the Print key.
- The arrow points to "Tare value" .
- Press the Print key: the input box for the tare value is displayed.
- Input the tare value and confirm with the Print key.
- Return to tare memory selection with the Zeroing key.
- If further tare memories are to be assigned, select the next free memory location by pressing the ORG key and repeat the procedure described above.
- To return to weighing mode, press the Zeroing key and then the ORG key.

### Changing or overwriting a fixed tare value memory

- Press the Print key and the F key simultaneously; setting mode appears in the display.
- The arrow points to "Program Settings". Confirm with the Print key.
- The arrow points to "General". Confirm with the Print key.
- Scroll to "Fixed tare value" using the Tare key and accept with the Print key.
- The last selected fixed tare value memory is displayed.
- Select the tare memory you wish to change using the F key or the Tare key, or input the memory location number and confirm with the Print key.
- Press the Print key: the input box for describing the tare memory is displayed.
- Change the description with the alphanumeric keypad and confirm with the Print key.
- The arrow points to "Tare value" .
- Press the Print key: the input box for the tare value is displayed.
- Change the tare value using the numeric keypad and confirm with the Print key.
- Return to tare memory selection with the Zeroing key.
- If you wish to change more tare memories, use the F key or Tare key to select the next memory location you wish to change and repeat the procedure described above.
- To return to weighing mode, press the Zeroing key and then the ORG key.

### Calling up the list of fixed tare value memories for taring

The tare value from the fixed tare value memory can be activated when the scale is unladen (before the weighing operation) or when the scale is laden (during the weighing operation).

- A weight value appears in the display panel: e 0.000kg.
- First press the Print key and then simultaneously the Tare key.
- The last used fixed tare value memory appears in the display.
- Scroll to the required fixed tare value memory with the F key or Tare key or input the memory location number and confirm with the Print key.
- Activate the fixed tare value by pressing the Print key.
- The display shows the gross value less the tare value (net).
- Further tare values from the memory are activated by repeating the procedure described above.

### Clearing an existing fixed value

- Press the Print key and the F key simultaneously; setting mode appears in the display.
- The arrow points to "Program Settings". Confirm with the Print key.
- The arrow points to "General". Confirm with the Print key.
- Scroll to "Fixed tare value" using the Tare key and accept with the Print key.
- The last selected fixed tare value memory is displayed.
- Scroll to the required fixed tare value memory with the F key or Tare key or input the memory location number and confirm with the Print key.
- Clear the displayed fixed tare value memory with the C key.
- The gap in the consecutive number scheme is closed. The following fixed tare value memories are given a number reduced by 1.

### Error message when taring

"Error 06 Not possible to tare"

Possible causes:

Scale empty

Scale in overload or underload

Manual tare input with tare value above maximum load

### Possible displays with the Info key

Display the tare memory by pressing the Info key, then press the Tare key. Pressing the Zeroing key ends the display of the tare memory.

## BMI function

The Body Mass Index (BMI) is intended only for weighing people. It is used to assess whether people are within the normal weight range or are overweight.

The BMI is calculated according to the following formula:

$$\text{BMI} = \frac{\text{Weight in kg}}{(\text{Body height in m})^2}$$

The unit for the BMI is therefore kg/m<sup>2</sup>.

BMI classification separated according to gender:

Classification	m	f
Underweight	<20	<19
Normal weight	20-25	19-24
Overweight	25-30	24-30
Adiposity	30-40	30-40
Severe adiposity	>40	>40

### Setting the mode of operation

- Press the Print key and the F key simultaneously; setting mode appears in the display.
- The arrow points to "Program Settings". Confirm with the Print key.
- The arrow points to "General".
- Use the Tare key to scroll down to "Weighing + Taring" and confirm with the Print key.
- The arrow points to "F key". Confirm with the Print key.
- The arrow points to the last activated function.
- Scroll to "BMI function" using the Tare key and accept with the Print key. A beep tone confirms acceptance.
- Exit setting mode with the ORG key. The weight value is displayed.

### Application

Press the F key.

The message "Input body size" appears.

Input your body size in metres.

Confirm with the Print key.

Now, parallel to the person's weight, the BMI value is displayed in the additional display.

The BMI display can be switched off by pressing the F key again.

## Printing

The Program Terminal 3020 can be connected to a printer via a serial port to print out scale chits. Up to a max. of 15 print templates can be stored.

To design print templates with individual contents the Soehnle Professional service centres offer service software for the 3020/3030.

Parameters listed below in the standard print templates can be used to design individual print templates. In addition, it is possible to copy organisation data (identifiers) and individual parameters such as company names.

Please contact your Soehnle Professional service partner about your print template requirements.

Printing is triggered either by means of the Print key or via another print-command condition, such as the F key.

### Configuring the I/O port for the printer

Call up the setting menu using the Print key and F key (press simultaneously).

- The arrow points to "Program Settings" .
- Scroll to "I/O Ports" with the Tare key.
- Press the Print key: the arrow points to "Data transfer fixed values".
- Scroll to "I/O port 1" using the Tare key.
- Press the Print key: the arrow points to "Application".
- Press the Print key: the arrow points to "2795.11".
- Select the printer you wish to connect using the Tare key or F key.
- Confirm with the Print key.

Information: Soehnle Professional printers are stored in the memory as standard. However, it is also possible to store other printers in the memory with the 3020/3030 service software.

Assign further parameters such as baud rate, data bits, parity and XON/XOFF as indicated in the printer manual.

The parameters set at the factory are configured for Soehnle Professional printers.

### Standard print template factory setting

The print-command keys can be assigned using the print templates in the service software. In the factory setting the following print templates are already stored in the memory and assigned. They can therefore be printed.

This standard print template is an example that shows which parameters are available.

Date		01.09.05
Time		09:13:28
Scale no.		02
Gross	W	0,154 kg
Tare	T	0,015 kg

**Print-command conditions**

The print templates are assigned to print-command conditions using the service software. The following conditions can be assigned to print templates:

- Print key
- Automatic print/data record once immediately
- Automatic print/data record once after weight change at standstill
- Automatic print/data record after weight change at standstill
- Automatic print/data record after weight change at standstill above the empty message
- Automatic print/data record after weight change at standstill above the empty message if the empty message was previously active.
- Continuous transmission
- Enter key for input from ORG A
- Enter key for input from ORG B
- Enter key for input from ORG C
- Enter key for input from ORG D
- Plus key (F key for totalising)
- Key combination INFO, F and Print key (setting-dependent)
- Enter key for manual input
- Print alibi memory
- Print reference statistics

Error message	Cause /Remedy / Information
2 Ident chip error	For Customer Service for error diagnosis
4 REF weight too small	Scale unsuitable for part weight, if applicable change setting in setting mode
5 Zeroing not possible	Deviation from zero point too great
6 Taring not possible	Empty scale, scale in underload or overload, manual tare input with tare value above maximum load
7 Not authorised to print	
8 No unit switchover	
9 Print signal buffer full	
10 Load cell signal overflow	
11 No measuring point signal	
12 Measuring point number ? CAN addr	
13 Converter failure	
15 Measuring point write-protected	
16 No measuring point detected	
18 Measuring point not zeroed	
20 Print format empty	
23 Print error	
24 Print template memory corrupt	
25 No print template(s)	

**Error message****Cause /Remedy / Information**

26 Printer update cancelled	
30 Alibi memory full	
33 Alibi memory contents corrupt	
34 Alibi memory empty	
35 No or defective alibi memory	
36 Cannot write alibi	
37 Alibi memory unreadable	
40 General data format error	
41 Input value too small	
42 Input value too big	
44 Power-off not allowed	
50 Value below zero range limit	Zeroing not possible because measured value is outside the set zero range limit
51 Value above zero range limit	Zeroing not possible because measured value is outside the set zero range limit
52 Overload	Measured value above weighing range
53 Underload	Measured value below 0
55 Position weight	Scale empty or no clip off has yet been set according to the set clip off factor.
60 One Wire hardware error	For Customer Service for error diagnosis
61 One Wire device absent	For Customer Service for error diagnosis
62 One Wire protocol error	For Customer Service for error diagnosis
63 One Wire CRC error	For Customer Service for error diagnosis
64 Unknown family code found	For Customer Service for error diagnosis
71 CAN input error	For Customer Service for error diagnosis
72 CAN messages lost error	For Customer Service for error diagnosis
73 CAN timeout error	For Customer Service for error diagnosis
75 Replace backup battery	
80 Abort external command	
81 Timeout external command	
82 Number already exists	The identification assigned to a fixed value already exists
88 System overload	For Customer Service for error diagnosis
90 Fixed value memory assigned / full	A total of 400 memory locations are available and are valid for all programs.
98 Wrong barcode	
99 Global error	For Customer Service for error diagnosis

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**Soehnle Professional GmbH & Co. KG**