

SANYO

INSTRUCTION MANUAL

MDF-136

MDF-236

MDF-436

BIOMEDICAL FREEZER



MDF-236

Note:

- 1. No part of this manual may be reproduced in any form without the expressed written permission of SANYO.**
- 2. The contents of this manual are subject to change without notice.**
- 3. Please contact SANYO if any point in this manual is unclear or if there are any inaccuracies.**

SANYO Electric Biomedical Co., Ltd. All rights reserved. Printed in Japan.

CONTENTS

PRECAUTIONS FOR SAFE OPERATION	P. 2
CAUTIONS FOR USAGE	P. 6
ENVIRONMENTAL CONDITIONS	P. 7
FREEZER COMPONENTS	P. 8
INSTALLATION	P. 11
OPTIONAL COMPONENTS	P. 13
START-UP OF UNIT	P. 15
REMOTE ALARM TERMINAL	P. 15
TEMPERATURE SETTING	P. 16
ALARMS & SAFETY FUNCTIONS	P. 18
SETTING OF ALARM RESUME TIME	P. 19
ROUTINE MAINTENANCE	P. 20
TROUBLESHOOTING	P. 21
DISPOSAL OF UNIT	P. 21
SPECIFICATIONS	P. 22
PERFORMANCE	P. 23
SAFETY CHECK SHEET	P. 24

PRECAUTIONS FOR SAFE OPERATION

It is imperative that the user complies with this manual as it contains important safety advice.

Items and procedures are described so that you can use this unit correctly and safely. If the precautions advised are followed, this will prevent possible injury to the user and any other person.

Precautions are illustrated in the following way:

WARNING

Failure to observe WARNING signs could result in a hazard to personnel possibly resulting in serious injury or death.

CAUTION

Failure to observe CAUTION signs could result in injury to personnel and damage to the unit and associated property.

Symbol shows;

-  this symbol means caution.
-  this symbol means an action is prohibited.
-  this symbol means an instruction must be followed.

Be sure to keep this manual in a place accessible to users of this unit.

< Label on the unit >
 This mark is labeled on the cover in which the electrical components of high voltage are enclosed to prevent the electric shock.
The cover should be removed by a qualified engineer or a service personnel only.

PRECAUTIONS FOR SAFE OPERATION

WARNING

-  **Do not use the unit outdoors.** Current leakage or electric shock may result if the unit is exposed to rain water.
-  **Only qualified engineers or service personnel should install the unit.** The installation by unqualified personnel may cause electric shock or fire.
-  **Install the unit on a sturdy floor.** If the floor is not strong enough or the installation site is not adequate, this may result in injury from the unit falling or tipping over.
-  **Never install the unit in a humid place or a place where it is likely to be splashed by water.** Deterioration of the insulation may result which could cause current leakage or electric shock.
-  **Never install the unit in a flammable or volatile location.** This may cause explosion or fire.
-  **Never install the unit where acid or corrosive gases are present** as current leakage or electric shock may result due to corrosion.
-  **Use a dedicated power source** as indicated on the rating label attached to the unit.
-  **Remove dust from the power supply plug** before inserting in a power source. A dusty plug or improper insertion may pose a hazard.
-  **Use a power supply outlet with ground (earth)** to prevent electric shock. If the power supply outlet is not grounded, it will be necessary to install a ground by qualified engineers.
-  **Never ground the unit through a gas pipe, water main, telephone line or lightning rod.** Such grounding may cause electric shock in the case of an incomplete circuit.
-  **Do not insert metal objects such as a pin or a wire into any vent, gap or any outlet** for inner air circulation. This may cause electric shock or injury by accidental contact with moving parts.
-  **Never store volatile or flammable substances** in this unit. This may cause explosion or fire.
-  **Never store corrosive substances** in this unit. This may lead to damage to the inner components or electric parts.
-  **Use this unit in safe area when treating the poison, harmful or radiate articles.** Improper use may cause bad effect on your health or environment.
-  **Disconnect the power supply to the unit prior to any repair or maintenance** of the unit in order to prevent electric shock or injury.
-  **Ensure you do not inhale or consume medication or aerosols** from around the unit at the time of maintenance. These may be harmful to your health.

PRECAUTIONS FOR SAFE OPERATION

WARNING

-  **Never splash water directly onto the unit** as this may cause electric shock or short circuit.
-  **Never disassemble, repair, or modify the unit yourself.** Any such work carried out by an unauthorized person may result in fire or injury due to a malfunction.
-  **Disconnect the power supply plug if there is something wrong with the unit.** Continued abnormal operation may cause electric shock or fire.
-  If the unit is to be stored unused in an unsupervised area for an extended period, **ensure that children do not have access and that doors cannot be closed completely.**
-  **The disposal of the unit should be accomplished by appropriate personnel. Remove doors** to prevent accidents such as suffocation.
-  **Prepare a safety check sheet** when you request any repair or maintenance for the safety of service personnel.

CAUTION

-  **Select a level and sturdy floor for installation.** This precaution will prevent the unit from tipping. Improper installation may result in water spillage or injury from the unit tipping over.
-  **Connect the unit to a power source as indicated on the rating label attached to the unit.** Use of any other voltage or frequency other than that on the rating label may cause fire or electric shock.
-  **When removing the plug from the power supply outlet, grip the power supply plug,** not the cord. Pulling the cord may result in electric shock or fire by short circuit.
-  **Never damage or break the power supply plug or cord. Do not use the supply plug if its cord is loose.** This may cause fire or electric shock.
-  **Do not touch any electrical parts** such as the power supply plug or any switches with a wet hand. This may cause electric shock.
-  **Do not put a container with water or heavy articles on the unit.** It may cause injury if the articles fall. Current leakage or electric shock may be resulted from the deterioration of insulation by spilled water.

PRECAUTIONS FOR SAFE OPERATION

CAUTION

-  **Do not climb onto the unit or do not put articles on the unit.** This may cause injury by tipping or damage to the unit.
-  **Put on dry gloves when you take out refrigerated articles from the freezer.** Handling frozen contents or the inside walls with naked hands may cause frostbite.
-  **Do not defrost inside walls using a knife or ice pick.** There are pipelines for cooling behind the walls. **Be careful not to damage the lines** as this could cause a breakdown. Also, **do not make a hole in the wall** for installation of attachments.
-  **Hold the handle when closing the door.** This will reduce the likelihood of a trapped finger.
-  **Disconnect the power supply plug** before moving the unit. Take care not to damage the power cord. A damaged cord may cause electric shock or fire.
-  **Be careful not to tip over the unit** during movement to prevent damage or injury.
-  **Disconnect the power plug** when the unit is not used for long periods.
-  **Do not put the packing plastic bag within reach of children** as suffocation may result
-  **Remove dust and clean the power supply plug periodically.** The operation with a dusty plug may cause fire.

CAUTIONS FOR USAGE


1. Sometimes the alarm may not operate at the time of first start-up. This does not mean malfunction. It is due to the complete discharge of incorporated battery. 2-day continuous operation of the freezer is necessary to charge the battery fully.
2. The digital thermometer of medical freezer is designed to display the temperature of the center part of the freezing compartment. Although the thermometer sometimes displays a temperature a little bit higher than the actual temperature of the center part, it gradually approaches the real temperature.
3. An access port to take out the measuring cable in the case is provided on the back wall of the freezer. Be sure to replace the cap and heat insulator after take out cable or, the inside temperature cannot complete down, and frost may accumulate outside the port surroundings.
4. Do not use brushes, acids, benzine, thinners, powdered soap. Polishing powders, or hot water for cleaning, as they can deteriorate the painted surfaces and parts made of plastic and rubber. Be especially careful not to wipe plastic or rubber parts with volatile solvents such as benzine. Wipe up the neutral detergents with a wet cloth when used.
5. Frost will accumulate on the inside walls of the freezer during use. Then remove the frost with the spatula provided or a similar tool. Do not defrost inside walls using a knife or ice pick. There are pipelines for cooling behind the walls. Be careful not to damage the lines as this could cause a breakdown. Also, do not make a hole in the wall for installation of attachments.
6. The freezing temperature is different at each position inside the case of this freezer. The set temperature represents the temperature at the center of the inside case. The refrigerated articles put on the highest position part, are refrigerated higher than the set temperature, while those put on the lowest part are refrigerated lower.
7. Do not put too many worm articles into a freezer compartment before enough operating. Put items in a few at a time after the freezer compartment temperature has cooled to at least -20°C .
8. In the case of high ambient temperature, the cabinet front may heat up after the freezer starts to operate first. However, this does not denote a malfunction. It is due to heater or hot gas piped around the unit frame to prevent frost and ice sticking around the cabinet.

ENVIRONMENTAL CONDITIONS


This equipment is designed to be safe at least under the following conditions (based on the IEC 1010-1):

1. Indoor use;
2. Altitude up to 2000 m;
3. Ambient temperature 5°C to 40°C
4. Maximum relative humidity 80% for temperature up to 31°C decreasing linearly to 50% relative humidity at 40°C;
5. Mains supply voltage fluctuations not to exceed $\pm 10\%$ of the nominal voltage;
6. Other supply voltage fluctuations as stated by the manufacturer;
7. Transient overvoltages according to Installation Categories (Overvoltage Categories) II; For mains supply the minimum and normal category is II;
8. Pollution degree 2 in accordance with IEC 664.

FREEZER COMPONENTS



MDF-236




FREEZER COMPONENTS

- 1. Door:** To open the door, grip the handle.
- 2. Handle:** Always grip this handle to open and close the outer door.
- 3. Door gasket:** This provides a tight door seal and prevents cold air leak. Keep clean.
- 4. Access port:** This is used for leading the measuring cable from the freezing chamber to the outside.
- 5. Lock:** Turn counterclockwise to 180° with a key and the outer door is securely locked.
- 6. Control panel:** Temperature set key, buzzer key and alarm lamp etc. are installed on the panel. Digital temperature display is also provided on it. See page 10 for the details.
- 7. Caster:** 4 casters are provided to facilitate moving of the cabinet. At the time of installation, make sure that the front 2 casters are not contact with the floor, by adjusting the leveling legs.
- 8. Leveling foot:** 2 feet are provided on the front side (right and left). Keep the unit in level by adjusting these legs at the installation.
- 9. Space for temperature recorder:** An automatic temperature recorder (optional part) can be attached here. See page 13 “Temperature recorder”.
- 10. Drain port:** The water accumulated on the bottom of the chamber can be drained through this port.
- 11. Basket:** Used for storing the materials in the chamber.
- 12. Remote alarm terminal (on back side):** This is used to notice an alarm condition of the unit to remote location. Refer to page 15 “Remote alarm terminal”.
- 13. Battery switch:** This is a switch for power failure. Always set the switch in ON position. Be sure to turn off this switch to save the battery if the freezer is not in operating for the long period (more than 1 month).
- 14. Power switch:** This is for turning ON/OFF the power to the unit. ON – “I” OFF – “O”. This has a function as an over-current protection breaker (15A).

FREEZER COMPONENTS

Control panel and keypad



1. Alarm test key (ALARM TEST): Test key for alarm device. By pressing this key, the alarm lamp is flashed, remote alarm is activated and buzzer sound. This means all alarm function operate correctly.

2. Buzzer stop key (BUZZER): To silence the audible alarm, press this key. The remote alarm is also silenced by pressing this key. (The buzzer cannot be stopped during remote alarm testing.)

3. Defrost key (DEF): By pressing this key for 5 seconds, the refrigerating operation is stopped. Pressing this key again after defrosting leads resumption of the operation.

Note: The refrigerating operation never resume automatically after defrosting.

4. Set key (SET): Temperature setting mode is led by pressing this key. Once the key is pressed, the changeable digit is flashed. Pressing this key again after setting desired temperature, the setting is stored into computer memory. If there is no key operation for 90 seconds during the temperature setting mode, the temperature setting mode is invalid automatically. See page 16 for the details.

5. Digit shift key (▶▶): Pressing this key in the setting mode causes the changeable digit to shift. Key lock is available by pressing this key for more than 5 seconds in the temperature display mode. Refer to page 16 for the key lock.

6. Numerical value shift key (▲): Pressing this key in the setting mode causes the numerical value to shift. "ON-OFF" of key lock can be selected by pressing this key in the key lock mode.

7. Alarm lamp (ALARM): This lamp is flashed when the audible alarm is activated.

8. Digital temperature indicator: This indicator shows the present chamber temperature or set temperature.

INSTALLATION

Installation site

To operate this unit properly and to obtain maximum performance, install the unit in a location with the following conditions:

1. A location not subjected to direct sunlight

Installation in a location subjected to direct sunlight may lead to inadequate cooling.

2. A location with adequate ventilation

Leave at least 10 cm around the unit for ventilation. Poor ventilation will result in a reduction of the refrigeration capacity.

3. A location away from heat generating sources

Avoid installing the unit near heat-emitting appliances such as gas ranges or stoves. Heat can cause inefficient refrigeration.

4. A location with a sturdy and level floor



WARNING

Install the unit on a sturdy floor. If the floor is not strong enough or the installation site is not adequate, this may result in injury from the unit falling or tipping over.

Select a level and sturdy floor for installation. This precaution will prevent the unit from tipping. Improper installation may result in water spillage or injury from the unit tipping over.

5. A location without flammable or corrosive gas



WARNING

Never install the unit in a flammable or volatile location. This may cause explosion or fire.

Never install the unit where acid or corrosive gases are present as current leakage or electric shock may result due to corrosion.

6. A location not prone to high humidity



WARNING

Do not use the unit outdoors. Current leakage or electric shock may result if the unit is exposed to rain water.

Never install the unit in a humid place or a place where it is likely to be splashed by water. Deterioration of the insulation may result which could cause current leakage or electric shock.

INSTALLATION

Installation

1. Remove the packaging materials and tapes

Remove all transportation packaging materials and tapes. Open the doors and ventilate the unit. If the outside panels are dirty, clean them with a neutral detergent and wipe it up with a wet cloth.


2. Adjust the leveling feet

Extend the leveling feet by rotating them counterclockwise to contact them to the floor. Ensure the unit is level.

3. Fix the unit

Two fixtures are attached to the rear of the frame. Fix the frame to the wall with these fixtures and rope or chain.

4. Ground (earth)



⚠️ WARNING

Use a power supply outlet with ground (earth) to prevent electric shock. If the power supply outlet is not grounded, it is necessary to install a ground by qualified engineers.

Never ground the unit through a gas pipe, water main, telephone line or lightning rod. Such grounding may cause electric shock in the case of an incomplete circuit.

OPTIONAL COMPONENTS


Temperature recorder

WARNING


Always disconnect the power supply to the unit prior to attachment of a temperature recorder in order to prevent electric shock or injury.

An automatic temperature recorder is available for this freezer as the optional component. The type of the recorder is MTR-G85. For the attachment, the mounting kit (MDF-S740) is necessary. Following shows the attachment procedure.

1. Remove four screws on the cover for the recorder mounting space and take off the cover. (Fig. 1)



2. Fix temperature recorder to the mounting by using 2 screws and fixture enclosed with the mounting kit. (Fig. 2)




3. Connect the power cord of the recorder with the coupler in the unit compartment of the freezer.

4. Remove the caps and insulation covering the access port on the chamber side and unit compartment side. (Fig. 3)

Then pass the recorder sensor through this port to go into the freezer chamber.


5. Insert the recorder sensor into the tube enclosed with the mounting kit.




OPTIONAL COMPONENTS

6. As shown in the Fig. 4, fix the recorder sensor on the chamber wall by using 2 clips with a 10mm spacer located between the wall and the clip. The clips and spacers are enclosed with the mounting kit.

7. Cover the gap of the access port on the chamber side completely by a silicon sealant after pushing the insulation into the gap.



8. After binding the lead wire on the back of the recorder and fix the recorder to the freezer with 4 screws removed in procedure 1. (Fig. 5)



START-UP OF UNIT

Follow the procedures for the initial and consequent operations of the unit.

1. Turn the power switch ON with the chamber empty.
2. Turn ON the battery switch.
3. The alarm sometimes may operate. In this case, silence the alarm by pressing BUZZER key.
4. Set the chamber temperature to a desired value.
5. Allow the unit to achieve the desired chamber temperature.
6. Check that the alarm lamp is flashed and the buzzer is activated by pressing ALARM TEST key. The remote alarm is also operated. E09 is displayed on the control panel if the battery switch is OFF. Make sure to turn on the battery switch.
7. Now you can put articles into the freezer chamber gradually to minimize the temperature rise.

REMOTE ALARM TERMINAL

WARNING

Always disconnect the power supply cord before connecting an alarm device to the remote alarm terminal.


The terminal of the remote alarm is installed at the back of the unit. The alarm is outputted from this terminal. Contact capacity is DC 30V, 2 A.

Contact output: At normal condition "Open"
At abnormal condition "Close"

Note:

The remote alarm is silenced by pressing the BUZZER key as the remote alarm is operated in conjunction with the buzzer except for the power failure alarm and alarm test.

The contact is closed in the case of power failure.



TEMPERATURE SETTING

Chamber temperature

Table 1 shows the basic procedure for setting the chamber temperature. Perform key operations in the sequence indicated in the table. The example in the table is based on the assumption that the desired temperature is -35°C .

Note: The unit is set at the factory that the chamber temperature -30°C .

Table 1. Basic operation sequence (Example: Chamber temperature -35°C)

	Description of operation	Key operated	Indication after operation
1	Switch ON the freezer.	----	The current chamber temperature is displayed. 
2	Press SET key.	SET	The second digit is flashed. 
3	Press ▶▶ key.	▶▶	When pressed, the settable digit moves. 
4	Press ▲ key and scroll the figure to 5.	▲	When pressed, the figure of settable digit increases. 
5	Press SET key.	SET	Set temperature is memorized and the current chamber temperature is displayed. 

Note:

The temperature set mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation.

Although the value of the chamber temperature setting can range between -18°C and -40°C , the guaranteed temperature without load is -35°C at ambient temperature of 35°C .

Key lock function

This unit is provided with the key lock function. When the key lock is ON, change of temperature setting through the key pad is not available. The key lock is set in OFF at the factory.

Display	Mode	Function
	Key lock is OFF	Enable to change of temperature setting
	Key lock is ON	Disable to change of temperature setting

Table 2. Procedure for key lock setting (change from key lock OFF to key lock ON)

	Description of operation	Key operated	Indication after operation
1		----	The current chamber temperature is displayed. 
2	Press ▶▶ key for 5 seconds.	▶▶	The right digit is flashed. 
4	Press ▲ key and scroll the figure to 1.	▲	When pressed, the figure of settable digit increases. 
5	Press SET key.	SET	The key lock is set to ON. The current chamber temperature is displayed. 

TEMPERATURE SETTING

Alarm temperature setting

This unit is provided with both high and low temperature alarms. The temperature at which the alarm is activated may be changed.

The available set range for high temperature alarm is between +5°C and +15°C, and -5°C and -15°C for low temperature alarm against the chamber temperature.

Note: The temperature alarm is set at $\pm 10^\circ\text{C}$ of the set temperature at the factory.

Display	Mode	Function
F01	High temperature alarm set	See Table 5 on page 18
F02	Low temperature alarm set	See Table 5 on page 18

As an example, Table 3 shows the procedure to set the high temperature alarm so that the alarm can activate when the chamber temperature is 5°C higher than the set temperature.

Table 4 shows the procedure to set the low temperature alarm so that the alarm can activate when the chamber temperature is 5°C lower than the set temperature.

Table 3. Procedure for setting high temperature alarm

	Description of operation	Key operated	Indication after operation
1		----	The current chamber temperature is displayed. -30
2	Press ▲ key for 5 seconds.	▲	The first digit is flashed. F00
3	Press ▲ key and scroll the figure to 1.	▲	When pressed, the figure of settable digit increases. F01
4	Press SET key.	SET	The right digit is flashed. 010
5	Set the temperature to 005 with the ▶▶ key and ▲ key.	▶▶	Pressing the key shifts the digit which can be set.
		▲	When pressed, the figure of settable digit increases. 005
6	Press SET key.	SET	Alarm temperature is memorized and the current chamber temperature is displayed. -30

Table 4. Procedure for setting low temperature alarm

	Description of operation	Key operated	Indication after operation
1		----	The current chamber temperature is displayed. -30
2	Press ▲ key for 5 seconds.	▲	The right digit is flashed. F00
3	Press ▲ key and scroll the figure to 2.	▲	When pressed, the figure of settable digit increases. F02
4	Press SET key.	SET	The right digit is flashed. -10
5	Set the temperature to -05 with the ▶▶ key and ▲ key.	▶▶	Pressing the key shifts the digit which can be set.
		▲	When pressed, the figure of settable digit increases. -05
6	Press SET key.	SET	Alarm temperature is memorized and the current chamber temperature is displayed. -30

ALARMS & SAFETY FUNCTIONS

This unit has the alarms and safety functions shown in Table 5, and also self diagnostic functions.

Table 5 Alarms and safety functions

Alarm & Safety	Situation	Indication	Buzzer	Safety operation
High temperature alarm	If the chamber temperature is higher than the temperature at which the high temperature alarm is activated.	ALARM lamp is flashed. Temperature indicator is flashed.	Intermittent tone with 15 minutes delay.	Remote alarm with 15 minutes delay.
Low temperature alarm	If the chamber temperature is lower than the temperature at which the low temperature alarm is activated.	ALARM lamp is flashed. Temperature indicator is flashed.	Intermittent tone with 15 minutes delay.	Remote alarm with 15 minutes delay.
Power failure alarm	In the case of power failure. When power switch is turned OFF. When the power to the unit is disconnected.	ALARM lamp is flashed.	Intermittent tone	Remote alarm.
Auto-return	When there is no key pressing in each setting mode for 90 seconds.	Chamber temperature is displayed.	----	Finishing of each Setting mode.
Key lock	When the key lock is "ON".	----	----	Change of setting is disable.
Thermal sensor abnormality	If the thermal sensor is disconnected.	ALARM lamp is flashed. E01 and chamber temp. are displayed alternately.	Intermittent tone	Remote alarm. Continuous running.
	If the thermal sensor is short-circuited.	ALARM lamp is flashed. E02 and chamber temp. are displayed alternately.	Intermittent tone	Remote alarm. Continuous running.
Thermal sensor abnormality (MDF-436 only)	If the protective sensor for compressor is disconnected.	ALARM lamp is flashed. E03 and chamber temp. are displayed alternately.	Intermittent tone	Remote alarm. Normal operation.
	If the protective sensor for compressor is short-circuited.	ALARM lamp is flashed. E04 and chamber temp. are displayed alternately.	Intermittent tone	Remote alarm. Normal operation.
Battery switch check	When the battery switch is OFF at the time of alarm test.	E09 is flashed.	Intermittent tone	Remote alarm.
Compressor temp. abnormality (MDF-436 only)	In the case of failure of compressor cooling fan motor. In the case of abnormal high temperature due to the dust on the condenser. In the case of abnormal high ambient temperature.	E10 is flashed.	Intermittent tone	Remote alarm. Compressor stops running.

Note:

1. The above power failure alarm is available when the battery switch is ON and the battery is charged. If the battery switch is OFF or the battery is discharged, only the remote alarm is activated.
2. The power failure alarm can be kept about 12 hours with the battery charged completely. 2-day operation of the freezer is needed to charge the battery full.
3. The chamber temperature is displayed for 5 seconds if the BUZZER key is depressed during the power failure alarm. After that, the alarm stops and the alarm lamp blinks.
4. The remote alarm is silenced by pressing BUZZER key as the remote alarm is operated in conjunction with the buzzer, except for the power failure alarm.
5. After power failure, the operation is resumed with the condition before power failure since the temperature setting and alarm temperature setting are memorized in a nonvolatile memory.

SETTING OF ALARM RESUME TIME

The alarm buzzer and remote alarm are silenced by pressing BUZZER key on the control panel during alarm condition. The buzzer and remote alarm will be activated again after certain suspension if the alarm condition is continued. The suspension time can be set by following the procedure shown in the Table 6 below.

The example in the table is based on the assumption that the desired duration is 20 minutes.

Note: The duration is set in 30 minutes at the factory.

Table 6. Setting procedure for alarm resuming time (change from 30 minutes to 20 minutes)

	Description of operation	Key operated	Indication after operation
1		----	The current chamber temperature is displayed. 
2	Press ▲ key for 5 seconds.	▲	The first digit is flashed. 
3	Set the figure to F25 with the ►► key and ▲ key.	►►	Pressing the key shifts the digit which can be set.
		▲	When pressed, the figure of settable digit increases. 
4	Press SET key.	SET	The current reset time is displayed. The middle digit is flashed. 
5	Set the figure to 020 with the ▲ key.	▲	When pressed, the figure of settable digit increases. 
6	Press SET key.	SET	Alarm temperature is memorized and the current chamber temperature is displayed. 

- The settable alarm resume time are 10, 20, 30, 40, 50, or 60 minutes. The buzzer and remote alarm would not reset if the reset time is set in 000.
- The setting of alarm reset time cannot be changed during the defrosting.
- The buzzer and remote alarm during power failure or alarm testing cannot be silenced.
- The set mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation. In this case, any setting before pressing SET key is not memorized.

ROUTINE MAINTENANCE

WARNING

Always disconnect the power supply to the unit prior to any repair or maintenance of the unit in order to prevent electric shock or injury.

Ensure you do not inhale or consume medication or aerosols from around the unit at the time of maintenance. These may be harmful to your health.

CAUTION

Always put on the dry gloves to protect the hands at the time of maintenance. No gloves may cause cut of the finger by the edge or corner.

Cleaning of cabinet

1. Clean the unit once a month. Regular cleaning keeps the unit looking new.
 2. Use a dry cloth to wipe off small amounts of dirt on the outside and inside of the unit and all accessories. If the unit is very dirty, use a neutral detergent.
 3. After cleaning, wipe away the cleaner completely with a cloth washed in clean water.
 4. Never pour water onto or into the unit. Doing so can damage the electric insulation and cause failure.
- The compressor and other mechanical parts are completely sealed. This unit requires absolutely no lubrication.
 - Remove dust from the power supply plug periodically.

Defrosting

This freezer is direct-cooling type and the frost is built on the chamber wall during long term operation. The excessive frost possibly make some gap between the cabinet and door gasket, which may cause poor cooling. Remove the frost inside the chamber once a month. Following shows the procedure for removing the chamber frost.

1. Temporarily move all the contents of freezer chamber to another freezer.
2. Press DEF key for 5 seconds to stop the refrigerating operation. While the refrigerating operation is stopped, the current chamber temperature and dF is displayed on the control panel alternatively.
3. After a several hours, check visually that all defrost was removed completely.
4. Remove the cap of drain port on the bottom of the chamber and drain the accumulated water.
5. Wipe up the water remaining in the freezer chamber and then replace the cap of drain port.
6. Press DEF key so that the refrigerating operation can be started.
7. Once the chamber temperature has dropped to the desired temperature, place the original contents back in the freezer chamber.

Note:

After the defrosting, the refrigerating operation is never resumed automatically. Make sure to press DEF key to start the freezer operation after defrosting.

While the freezer stops refrigerating operation for defrosting, neither high temperature alarm nor low temperature alarm is effective.

TROUBLE SHOOTING

If the unit malfunctions, check out the following before calling for service.

The chamber is not cooled at all

1. The circuit breaker of power source is active.
2. The voltage is too low. (In this case, call an electrician.)
3. The power switch is not ON.
4. The large amount of articles (load) is stored in the chamber at one time.
5. The freezer is in defrost condition.

The cooling is poor

1. The ambient temperature is too high.
2. The door is not closed firmly.
3. The large amount of frost is built on the chamber wall.
4. The air intake vent is blocked.
5. The set temperature is not inputted properly.
6. The freezer is in the direct sunlight.
7. There is any heating source near the freezer.
8. A rubber cap and insulation for the access port are not set correctly.
9. You put too many unfrozen articles into the freezer compartment.

When the unit does not accept changes of set-point temperature

1. The key lock is set in "ON" mode.

Noise

1. The freezer is not installed on the sturdy floor.
2. The freezer is not leveled with the leveling feet.
3. There is anything touching the frame.
4. The freezer is in the status immediately after start up.

The unit sometimes causes a noise when the chamber temperature is high due to the large load. The noise gets less and less accompanying with the cooling of the chamber.

DISPOSAL OF UNIT

WARNING

If the unit is to be stored unused in an unsupervised area for an extended period **ensure that children do not have access and doors cannot be closed completely.**

The disposal of the unit should be accomplished by appropriate personnel. Always remove doors to prevent accidents such as suffocation.

SPECIFICATIONS

Name	Biomedical Freezer		
Model	MDF-136	MDF-236	MDF-436
External dimensions	W640 x D687 x H881 (mm)	W905 x D687 x H881 (mm)	W1265 x D807 x H905 (mm)
Internal dimensions	W525 x D440 x H715 (mm)	W790 x D440 x H715 (mm)	W1140 x D550 x H735 (mm)
Effective capacity	138 L	221 L	426 L
Exterior	Polyester finish baked on zinc galvanized steel		
Interior	Colored aluminum plate		
Door	Polyester finish baked on zinc galvanized steel		
Insulation	Rigid polyurethane foamed-in place (CFC-FREE)		
Baskets	Polyethylene coated steel wire		
Access port	Diameter 17 mm, Right side and bottom left		
Compressor	Hermetic type, 175 W		Hermetic type, 350 W
Condenser	Wire and tube type		Fin & tube + wire & tube type Fan motor 2W
Evaporator	Aluminum tube on sheet type		
Refrigerant	R-404A		
Temperature controller	Electronics controller (between -18 and -40°C)		
Temperature display	Digital display (between -50 and +50°C)		
Temperature sensor	Thermister sensor		
Temperature alarm	Flash of digital indicator and alarm lamp, Buzzer, (Remote alarm)		
Accessories	Key 1 set, Scraper 1, Basket 2 (MDF-136), 3 (MDF-236), 4 (MDF-436)		
Weight	45 kg	54 kg	71 kg
Battery	For power failure alarm, Nickel hydrogen battery, DC 6V, 1100mAh, Automatic charge		
Optional component	Automatic temperature recorder (MTR-G85) Mounting kit for automatic temperature recorder (MDF-S740) Baskets		

Note: Design or specifications will be subject to change without notice.

PERFORMANCE

MDF-136

Cooling performance	-35°C (ambient temperature; 35°C, no load)			
Temperature control range	-20°C to -35°C			
Rated voltage	AC 110 V	AC 115 V	AC 220 V	AC 220 to 240 V
Rated frequency	60 Hz	60 Hz	60 Hz	50 Hz
Rated power consumption	140 W	140 W	150 W	145 W
Noise level	37 dB [A] (background noise; 20 dB)			
Maximum pressure	2570 kPa			

MDF-236

Cooling performance	-35°C (ambient temperature; 35°C, no load)			
Temperature control range	-20°C to -35°C			
Rated voltage	AC 110 V	AC 115 V	AC 220 V	AC 220 to 240 V
Rated frequency	60 Hz	60 Hz	60 Hz	50 Hz
Rated power consumption	140 W	140 W	170 W	160 W
Noise level	36 dB [A] (background noise; 20 dB)			
Maximum pressure	2260 kPa			

MDF-436

Cooling performance	-35°C (ambient temperature; 35°C, no load)			
Temperature control range	-20°C to -35°C			
Rated voltage	AC 110 V	AC 115 V	AC 220 V	AC 220 to 240 V
Rated frequency	60 Hz	60 Hz	60 Hz	50 Hz
Rated power consumption	230 W	230 W	220 W	220 W
Noise level	40 dB [A] (background noise; 20 dB)			
Maximum pressure	2180 kPa			

Note : The unit with CE mark complies with EC directives 89/336/EEC, 93/68/EEC and 73/23/EEC.

⚠ CAUTION

**Please fill in this form before servicing.
Hand over this form to the service engineer to keep for his and your safety.**

Safety check sheet

1. Refrigerator contents : Yes No
Risk of infection: Yes No
Risk of toxicity: Yes No
Risk from radioactive sources: Yes No

(List all potentially hazardous materials that have been stored in this unit.)

Notes :

2. Contamination of the unit
- Unit interior Yes No
No contamination Yes No
Decontaminated Yes No
Contaminated Yes No
Others:

3. Instructions for safe repair/maintenance of the unit
- a) The unit is safe to work on Yes No
b) There is some danger (see below) Yes No
Procedure to be adhered to in order to reduce safety risk indicated in b) below.

Date :

Signature :

Address, Division :

Telephone :

Product name : Biomedical Freezer	Model : MDF-136/236/436	Serial number :	Date of Installation :
--------------------------------------	----------------------------	-----------------	------------------------

Please decontaminate the unit yourself before calling the service engineer.

7FB6P101279004

Recycled paper **CEM**

SANYO

SANYO Electric Biomedical Co., Ltd.

5-5, Keihan-Hondori 2-Chome
Moriguchi City, Osaka 570-8677 Japan

Printed in Japan