





SAFETY PRECAUTIONS

- Babies, young children, or individuals who cannot express their consent are not suitable for blood pressure measurement.
- Individuals with serious arteriosclerosis are not suitable for blood pressure measurement.
- Blood pressure readings may fluctuate during pregnancy. Pregnant women should consult their doctor before taking measurement.
- If there are unusual values, please consult your doctor before taking medical action.
- Handling batteries properly:
 - As soon as old batteries wear out, replace them with new batteries.
 - Do not use old and new batteries together.
 - Align the polarities of batteries correctly.
 - When the unit will not be used for more than 3 months, remove the batteries. Otherwise, batteries may leak and cause damage to the unit.

QUICK-START GUIDE


Avoid smoking, eating or exercising for 30 minutes before taking a measurement.
Do not move or talk during measurement.
Measurement can be taken on either arm; however, if possible, always take measurements on the left arm.

1. Sit upright in a chair with your feet on the floor.
2. Remove tight-fitting or thick clothing from the lower arm.
3. Pull on the end of the cuff until it wraps securely around your wrist. Place your arm on a table so that the cuff will be at the same level as your heart.
4. Press the "O" button. The cuff will automatically inflate and the measurement will start.
5. When the measurement is complete, the cuff will automatically deflate and your systolic and diastolic pressure values and pulse rate will display.



OVERVIEW


Monitor



Features

- One Touch Operation
- Irregular Heartbeat Detector
- Blood Pressure Mode and Clock Mode
- 4 User / 396 Total Memory
- 4 Alarm Clock Settings
- 5 Mode Averaging Function – Hour, Day, Week, Month, and Overall
- WHO Blood Pressure Indicator with 6-Color LED
- Auto-off Power Saving Feature
- Operated by 2 AAA batteries

Advanced Wrist Blood Pressure Monitor Model 1147 Instruction Manual



Read these instructions before using your Blood Pressure Monitor. Save these instructions for future reference.

1147-INS-LAB-RevA09

PREPARATION

2.1 Battery Installation

1. Remove the battery cover in the direction indicated on the cover.



2. Insert 2 AAA batteries with the polarities as indicated.



3. Slide the battery cover back to close it. The cover will click into place when completely closed.




Note: Replace batteries when the battery indicator shows low voltage.

2.2 Setting the Date and Time

- Note:
- Set the correct date and time before taking a measurement.
 - Reset the date and time after replacing the batteries.

1. Press and hold the "AVG" button for a few seconds until the Year display flashes.
2. Press the "MEMO" or "PAGE" button to select the current year.
3. Press the "AVG" button to confirm Year Setting.
4. Follow the same steps to set the Month, Date, Hour and Minute.

5




5. "ON" and "OFF" will then flash. This is to have the clock display when the unit is not in use; you can choose Clock Function "ON" or "OFF". Press "MEMO" or "PAGE" button to select. Press the "AVG" button to confirm.

Note: If you select Clock Function "ON", you will continue at the Alarm Setting section. Otherwise, setting will end.

2.3 Setting the Alarm

- Note:
- You will not be able to set alarm clock if Clock Function is set to "OFF".
 - You can set a maximum of 4 alarms.



6. To disable Alarm, Press the "MEMO" or "PAGE" button to select "OFF" and press the "AVG" button to confirm.
7. To set 1st Alarm, Press the "MEMO" or "PAGE" button to select "ON". Press the "AVG" button to confirm.
8. When Hour digit flashes, press the "MEMO" or "PAGE" button to select Hour. Press the "AVG" button to confirm.

6

TAKING A MEASUREMENT

- Note:
- Avoid smoking, eating or exercising for 30 minutes before taking a measurement.
 - Do not move or talk during measurement.
 - Measurement can be taken on either arm; however, if possible, always take measurements on the left arm.

3.1 Sitting Correctly


1. Relax.
2. Sit upright in a chair with your feet on the floor.
3. Remove thick or tight-fitting clothing from the lower arm. Do not roll up your sleeve if it is too tight.
4. Place your arm on a table so that the cuff will be at the same level as your heart.



7

3.2 Applying the Wrist Cuff

1. Place your elbow on a table with the palm facing upward. Wrap the wrist cuff around the left wrist as shown below.
2. Pull on the end of the cuff until it wraps securely around your upper wrist. Do not over tighten the cuff. Allow 1-2 cm or 1/2 inch between the bottom of the cuff and your wrist joint.
3. Make sure the cuff is at the same level as your heart.



3.3 User Setting

This unit allows a maximum of 4 users to store measurement records in different User Pages.

Press "PAGE" button to switch between User Page 1, 2, 3 and 4.

3.4 Taking A Measurement

1. Press the "O" button. The cuff will automatically inflate and the measurement will start.
 2. When the measurement is complete, the cuff will automatically deflate and your systolic and diastolic pressure values and pulse rate will display.
 3. Remove the wrist cuff and take note of your measurement results.
- Note: You can press the "O" button to stop measurement at any time.

8

3.5 Heartbeat Readings



Besides Systolic and Diastolic values, Pulse will also be displayed after a measurement

Irregular Heartbeat

This unit has a unique feature that alerts user of irregular heartbeat detection during measurement.

Note: An irregular heartbeat is defined as a heartbeat rhythm that has a variation of more than 25% from the average rhythm detected during measurement.



If an irregular rhythm occurs twice or more during the same measurement, the monitor will display the pulse value and "IH" symbol alternately. Example: (90 / IH)



Note: Please consult your doctor if irregular heartbeat occurs often.

Regular but Below Normal Heartbeat

If the monitor detects a pulse rate below 58 beats per minute during measurement, the monitor will display the pulse value and "LO" symbol alternately. Example: (58/LO)



9


Regular but Above Normal Heartbeat

If the monitor detects a pulse rate above 100 beats per minute during measurement, the monitor will display the pulse value and "HI" symbol alternately. Example: (102/HI)



3.6 WHO Light Indicator

This unit features a 6-color LED light indicator to display blood pressure level in accordance with the WHO (World Health Organization) classification.



10

- WHO Classification Color
- Optimal -----Blue
 - Normal -----Sky Blue
 - High-normal -----Green
 - Mild hypertension -----Purple
 - Moderate hypertension -----Yellow
 - Severe hypertension -----Red



WHO Light Indicator

Example: If your measurement reading is 122/75, the LED light will display the color in blue.

MEMORY FUNCTIONS

After measurement, the last set of measured values will be saved automatically after measurement when the unit is turned off by pressing "O" button or by automatic power-off after 3 minutes.

4.1 View Past Readings

1. Press the "MEMO" button to enter memory function.
2. If the User Page is not what you want to view, press the "PAGE" button to select.
3. Press the "MEMO" button to display the next sets of values for the user selected. Past readings are displayed from the most recent to the oldest.
4. This blood pressure monitor can store up to 99 sets of past readings for each user (4 x 99).

11

4.2 View Average Readings

1. Press the "AVG" button to enter Average function.
2. If the User Page is not what you want to view, press the "PAGE" button to select.
3. Press the "AVG" button to view different Average Values. 5 different Average Values are available:

Symbol on Display	Average Values of
1 h	One Hour
1 d	One Day
7 d	One Week
30 d	One Month
ALL	All Readings

4.3 Delete Past Readings

User can delete all the past readings stored in each User Page.

1. Press the "PAGE" button to select the User Page you want to delete.
2. Press and hold the "MEMO" & "PAGE" buttons simultaneously for 5 seconds. Symbol of "ALL EE" will be displayed and all past readings stored in that User Page will be deleted.

12

TROUBLESHOOTING

Symptom	Possible Cause	Solution
No display when you press the Power/Start button	Batteries are worn out Battery polarities are incorrect	Replace with new batteries Re-align batteries in correct polarities
"Err P" displays	Failure to inflate	Ensure the cuff is properly connected
"Err 1" displays	Deflating too rapidly	Return for service
"Err 2" displays	Movement during measurement	Do not move during measurement
	Signal interference	Remove interfering source, e.g. mobile phones, magnets
"Err 3" displays	Incorrect results	Measure again

Note: If your problem cannot be solved by the above, consult a GF Health Products, Inc. representative. DO NOT disassemble the unit.

STORAGE & MAINTENANCE

- Keep away the unit from direct sunlight, extreme temperatures, humidity or moisture.
- Use a dry or slightly water-dampened cloth to clean the unit.
- Do not use alcohol, benzene, thinner or other volatile liquids to clean the unit.

13

- Do not wash the wrist cuff or expose it to liquid.
- Remove the batteries from the unit when it will not be used for more than 3 months.

SPECIFICATIONS

Model No.	1147
Measurement Method	Oscillometric
Measurement Range	Pressure: 20 – 280 mmHg Pulse Rate: 40 – 195 beats/min
Accuracy	Pressure: +/- 3 mmHg Pulse Rate: +/- 5%
Inflation	Pump driven
Pressure Detection	Semi-conductor
Power Supply	2x AAA batteries
Auto Power Off	3 minutes after last operation
Storage Environment	-4°F to 140°F, 10% to 95% RH
Operation Environment	50°F to 104°F, 30% to 85% RH
External Dimensions	3.15" (L) x 3.35" (W) x 2.36" (H) or 80 x 85 x 60 mm
Net Weight	140g (batteries excluded)
Cuff Dimensions	Wrist circumference of 135~220mm
Package Contents	Blood Pressure Monitor, Instruction Manual, Batteries

Note: Specifications are subject to modification without prior notice.

14

LIFETIME LIMITED WARRANTY

GF Health Products, Inc. offers a lifetime limited warranty on defects in material and workmanship for this product, except as noted below:

This product warranty does not cover damage caused by misuse or abuse, the attachment of any unauthorized accessory, alteration to the product, or any other conditions whatsoever that are beyond the control of GF Health Products, Inc. GF Health Products, Inc. shall not be responsible for any type of incidental, consequential, or special damage.

To obtain warranty service on your Lumiscope Blood Pressure Monitor, please contact the Repair Department at 1-800-672-8293. There is a \$15 fee for return shipping and handling. Please make checks payable to GF Health Products, Inc.

Upon receipt, we will repair or replace, as appropriate, the Blood Pressure Monitor and return it to you. This warranty does not include any labor charges incurred in replacement part(s) installation or any associated freight or shipping charges to GF Health Products, Inc.

Warranty is solely through the Repair Department. Service of this product by anyone other than the Repair Department voids warranty.

15

GF Health Products, Inc.
2935 Northeast Parkway
Atlanta, Georgia 30360
telephone: 800-347-5678, 770-447-1609
fax: 800-726-0601, 678-291-3232



GRAHAM-FIELD
www.grahamfield.com

© December 2009 GF Health Products, Inc.
Made in China