**630-010**

<b>SE</b>	<b>Automatisk blodtrycksmätare för handled</b> .....	<b>2</b>
<b>DK</b>	<b>Automatisk blodtryksmåler til håndled</b> .....	<b>4</b>
<b>NO</b>	<b>Automatisk blodtryksmåler til håndleddet</b> .....	<b>6</b>
<b>FI</b>	<b>Automaattinen verenpainemittari</b> .....	<b>8</b>
<b>UK</b>	<b>Automatic Blood Pressure Monitor</b> .....	<b>10</b>
<b>DE</b>	<b>Vollautomatisches Blutdruckmeßgerät</b> .....	<b>12</b>
<b>PL</b>	<b>Automatyczny ciśnieniomierz do mierzenia ciśnienia w nadgarstku</b> .....	<b>14</b>
<b>RU</b>	<b>Аппарат для автоматического измерения кровяного давления на запясть</b> ....	<b>16</b>

## INLEDNING

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## INTRODUCTION

Please read all instructions before operating your new Idé-line blood pressure monitor. They contain important information on the operation, safety and maintenance of the appliance, keep these instructions in a safe place for future reference and show them to other users if necessary.

## KNOWLEDGE ABOUT BLOOD PRESSURE

Blood pressure is the force that blood exerts on the arteries. This force is constantly changing as the heart beats. When the heart contracts, blood pressure reaches its highest value, this is called systolic blood pressure. When the heart relaxes between beats, the value of blood pressure is lower, this is called the diastolic blood pressure.

High blood pressure, is a condition where an individual's blood pressure is high over a long period of time. Untreated, hypertension can cause many serious medical problems including strokes and heart attacks. To control hypertension it is therefore important to check the blood pressure and decrease it by reducing fat and salt intake, maintain proper weight and exercise. Even people who are born with high blood pressure can prevent worsening by doing suitable controls.

## CLASSIFICATION OF BLOOD PRESSURE

- Normal blood pressure should be under 140mmHg (systolic) and under 90mmHg (diastolic).
- The borderline for normal blood pressure lies between 140mmHg and 160mmHg (systolic) and between 90mmHg and 95mmHg (diastolic).

- High blood pressure is over 160mmHg (systolic) and 95mmHg (diastolic).

	<b>Max. Systolic</b>	<b>Max. Diastolic</b>
High blood pressure	Over 160	Or over 95
Borderline	140-159	Or 90-94
Normal blood pressure	101-139	And 61-89
Low blood pressure	Under 90	Or under 50

## QUESTIONS AND ANSWERS ON BLOOD PRESSURE

Why is it not healthy to have high blood pressure ?

The veins will be tighten when you have high blood pressure, and the vein can in some cases break or thicken with time. This can mean that important organs are not receiving the correct amount of blood, which again can cause serious diseases. Important organs which are mostly affected by high blood pressure includes:

### The Brain

High blood pressure can cause sclerosis of the veins in the brain. This can further cause a brain haemorrhage or a stroke.

### The Heart

High blood pressure can in serious cases cause coronary.

### The Kidneys

Sclerosis of the veins can result in a decreased kidney function, which further can result in urine poisoning.

### **How do you take a correct measurement of the blood pressure ?**

The blood pressure changes constantly. Various conditions are affecting the changes in the blood pressure. Therefore it is extremely important that the measurements are taken at the same time, several days and it is also important that the measurements are saved so you can control the development.

### **Are there any symptoms related to high blood pressure ?**

If your blood pressure is too high, there are no visible symptoms, but examples of symptoms are: headache, problems breathing, dizziness, sleeplessness etc.

### **Which conditions are affecting the blood pressure ?**

- Medication
- Exercise
- Sleep
- Diet
- Physical condition
- Working environment

### **Can you avoid / prevent high blood pressure ?**

If you have high blood pressure, you should be aware of the food you eat: eat less fat and less salt. You should also be aware of your working conditions: avoid stress. Furthermore exercise and plenty of rest are also helpful to prevent high blood pressure.


### **IMPORTANT SAFEGUARDS**

- Do not use the blood pressure monitor, before consulting your doctor.
- Keep the unit out of locations with high temperatures or high humidity and keep it out of direct sunlight.
- Do not immerse the appliance in water or any other liquids.
- Store the unit out of childrens reach.
- Do not push the Start/Stop button before the unit is securely placed on your wrist.
- To obtain the best results make sure that the person doing the measurements is relaxed and calm.
- One hour after eating a meal is a perfect time to measure the blood pressure, since you are very relaxed at that time.
- After smoking or drinking alcohol the blood pressure will increase.
- After exercising or after taking a bath it is recommended that you wait 20 minutes before measuring the blood pressure.
- Low temperatures will make the blood pressure rise, a room temperature of approx. 20° is perfect for measuring the blood pressure.
- If you wish to do a measurement after a meal, do the measurement before drinking alcohol or smoking.
- People with high blood pressure, diabetics, people with kidney problems etc. will properly find a large difference in the blood pressure taken on the upper arm (which is the case when you have it done at the doctor) and when the blood pressure is measured on the wrist. If there is any doubt, please contact your doctor.

## CORRECT METHOD OF MEASUREMENT



To obtain the most accurate blood pressure measurement, please follow these important directions.

- Be seated in a chair with back support.
- Rest your arm on a table so the cuff is at the same level as your heart. Use the storage box or similar item to support the forearm.

- Place both feet on the ground.

## LIST OF COMPONENTS


## DESCRIPTION OF DISPLAY MARKS

- ▼ Deflating
- ▲ Inflating
- Err Measurement error
- ♥ Measuring
- 🔋 Low battery

## PREPARATION BEFORE TAKING A MEASUREMENT

Battery installation:

- Pull off the battery cover in the same direction as the arrows are pointing.
- Install 2 pcs. alkaline AAA batteries, noting the proper orientation of positive and negative terminals of the batteries in the battery compartment.
- Replace the battery cover.
- If the low battery symbol appears on the display replace both batteries.
- It is recommended that you use Alkaline batteries since they are functioning for a longer period of time.
- Remove the batteries if the unit is not being used for a long period of time.

## Applying the waste cuff:

- With your palm side up, center monitor on your wrist and apply the cuff so that the cuff's edge should be approximately 1 cm from your palm.
- Do not apply wrist unit on top of a shirt sleeve.

- Wrap the cuff snugly around your wrist and attach using the fastener. The top and bottom edges of the cuff should be tightened evenly around your wrist.
- Any gap between the wrist and cuff will affect the result of the measurement.
- The fabric cuff can be folded down for ease of viewing the display.

## TAKING A MEASUREMENT

- After the cuff has been properly placed on either the right or left wrist, press the Start/Stop button. All symbols on the display will appear.
- After the display reads "0" the "Fuzzy Logic function will automatically inflate to an approximate pressure and will then begin to slowly deflate as it takes the measurement. When the measurement is done the heart symbol will be shown on the display and the systolic and diastolic values be shown in the display after 5 sec. will the pulse be shown in the display. Now will the values be

shown alternately.

- To end a measurement for any reason, just press Start/Stop button to turn off the unit and release the cuff pressure.
- If during a measurement the power drops to an insufficient level, the unit will end the measurement and display the low battery symbol.
- Replace the batteries and repeat the measurement.
- If the monitor is unable to measure a pulse it will end the measurement attempt. Wait a few minutes make sure that the cuff is positioned properly, and try again.
- The display automatically turns off after 3 minutes.

### **SETTING DATE AND TIME**

- To be able to save subsequent measurement with correct date and time, you have to set date and time after proper battery installation or replacement.
- After the batteries have been installed "0:00" and "1 – 1" on the display.
- Press the Mode button, the individual figures will start to flash.
- When the figure you wish to change is flashing, press the Set button to change it.
- To continue to the next figure, press the Mode button, the next figure starts to flash. Be aware that when setting the date, the month is on the left and the date is on the right.
- When the date and time has been installed, press the Start/Stop button to end the setting procedure.

### **MEMORY FUNCTION**

- When a valid measurement is made, the reading data will be stored with data and time into the Memory automatically as soon as the unit is turned off, automatically shuts off, or when another measurement is initiated
- To recall your last readings, press the memory button, the memory number will appear in the upper right hand corner of the display.
- Less memory number represents the newer reading. Up to 30 previous readings can be recalled.  
When "No:A" appears in display will the represents average value of your systolic, diastolic be shown. The blood pressure reading and pulse will then be displayed alternately. These readings are based on the current readings in the memory.
- By pressing the Start/Stop button at any time you would stop your memory recall and the display will return to normal time/date display.

## ERROR INDICATORS

Err  
330

- Indication: Cuff pressure reaches 330mmHg  
Cause: Twisted bladder inside off the cuff  
Correction: Make sure that the bladder lays flat within the cuff.  
Carefully read and repeat the above mentioned steps.

Err

- Indication: Displays "Err", cuff pressure releases and "0" pressure displays.  
Cause: Excessive arm movement or talking during measurement.  
Correction: Press Start/Stop button to measure again while keeping still.

Batteries are used. Please change the two batteries with new ones.

## CARE AND MAINTENANCE

- Do not clean the monitor with volatile liquids. The monitor should be cleaned with a soft, dry cloth.
- Take out the batteries if the unit is not being used for a long period of time.

Model	MS-900VX
Type	Oscillometric
Display	LCD: size: 26 x 26 mm
Measuring range	0 – 300 mmHg / 0-40 kPa (blood pressure) 30 – 190 pulse beats pr. minute
Accuracy	± 3mmHg / ± 0.4 kPa (blood pressure) ± 5% (pulse)
Pressure	Fuzzy logic
Deflation system	Electro-mechanical release valve
Cuff	Velcro with mechanical closing
Power source	2 pcs. AAA batteries – Alkaline (Not incl.)

## ENVIRONMENTAL TIP

An electronic appliance should, when it is no longer capable of functioning be disposed with least possible environmental damage. The appliance should be disposed according to the local regulations in your municipality, in most cases you can discharge the appliance at your local recycling center

## THE WARRENTY DO NOT COVER

- If the above points have not been observed.
- If the appliance has not been properly maintained, if force has been used against it or if it has been damaged in any other way.
- If the appliance has been repaired or modified or changed in any way or by any person not properly authorized.

Owing to our constant development of our products on both functionality and design we reserve the right to change the product without preceding notice.

#### **IMPORTER**

Adexi A/S

Adexi AB

We take reservations for printing errors.

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