OMRON

Instruction Manual 使用說明書

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Table of Contents

Before Using this Monitor

- Notes on Safety3-4
- Name of the Parts5
- Preparation before Use6-8

Correct Use of this Monitor

- How to Measure Blood Pressure9-10
- How to Use the Memory Function11-12

How to Care for this Monitor

- Maintenance13
- Error Indicators14
- Troubleshooting15
- Specifications16

A Few Words about Blood Pressure

- What You should know about Blood Pressure17
- For the Correct Measurement of Blood Pressure18
- Blood Pressure Classification19

Notes on Safety

Please observe the following safety notes

Caution

Self-diagnosis of measured results and treatment is dangerous. Please follow the instructions of your doctor.

- Self-diagnosis may worsen the disease.
- Diabetes, hyperlipemia, or hypertension may lead to arterial sclerosis. If the conditions are not checked, there is an increased danger of cerebral apoplexy or myocardial infarction, or arteriarctia or peripheral circulatory disturbances. In such cases, wrist and brachial blood pressure may differ greatly. (Even for the healthy person, there may be a difference of 20 mmHg if the measuring conditions are not appropriate. Therefore, do not diagnose the measured results by yourself. Consult your doctor and ask for his/her instruction.)
- The difference between wrist and brachial blood pressure values may be due to physiological conditions at the time of measurement. However, blood pressures measured at the wrist and those at the brachium (upper arm) fluctuate in the same manner. Therefore, you can also check fluctuation tendency by measuring blood pressure at your wrist.

Do not use with infant or person who cannot express one's intention.

- This may lead to accident or trouble.

Do not use for purposes other than to measure blood pressure.

- May lead to accident or trouble.

Do not use portable phones near this monitor.

- The monitor may operate erroneously.

Do not disassemble or modify the main unit or the wrist cuff of the blood pressure monitor.

- Correct measurement will be impossible.

Notes on Safety

Please take note

Worn battery may leak and damage the main unit. Please observe the following.

- When you are not going to use this monitor for a long period of time (approximately three months or more), remove the batteries.
- Replace worn batteries with new ones immediately.
- Do not use worn and new batteries together.
- Do not insert batteries with their polarities in the wrong direction.

Do not use force to stretch or bend the wrist cuff.

Do not press the START/STOP button when the cuff is not wrapped on the wrist.

Do not apply shock or drop the monitor.

Do not store the monitor under direct sunshine, where there is high temperature and humidity, dust, contact with corrosive gases or water.


Do not use benzine, thinner, gasoline, or alcohol to clean the monitor.


Names of the Parts

Main unit


Accessories

Two AAA Alkaline batteries (LR03)

Storage case


INTELLISENSE blood pressure monitor is a global brand name of blood pressure monitor equipped with bio-information sensing and high performance fuzzy logic technology, which are the key technologies of OMRON. With this excellent sensing technology, the monitor can achieve accurate measurement.

Preparation before Use

Insert the batteries.


- (1) Pinch the sides of the battery cover.

- (2) Pull out the battery cover.

- (3) Insert the batteries by matching the polarities (+ and -).

- (4) Replace the battery cover by pinching it with finger tips.



Battery life and replacement:

- With the use of high performance alkaline batteries (two AAA batteries), you can measure approximately 400 times if you use this monitor at room temperature of 22°C and inflate it to 170 mmHg, three times a day.
- The battery life varies with the ambient temperature. At low temperature such as during winter, battery life is shortened compared with the normal temperature of 22°C.
- The batteries included are for monitoring only and may not last 400 times.
- If the battery replacement mark  flashes, or if the monitor does not inflate when the START/STOP button is pressed, replace the two batteries with new ones of the same type.
- If you hold down the START/STOP button, batteries will wear out. Keep the monitor in the storage case when carrying or storing to avoid pressing on the START/STOP button.
- Replace the batteries while the power is turned off (when nothing is displayed).

Preparation before Use



Attach the Wrist Cuff.

1. Place the wrist cuff over your left wrist with your left thumb facing upward.

- * Wrap the wrist cuff over bare skin.
- * Be sure not to roll the sleeve of your shirt into the wrist cuff.
- * As shown in the Figure, leave a space of 10 to 15 mm (about the width of an index finger) between the edge of the wrist cuff and the palm.

2. Hold the bottom part of the wrist cuff and wrap it around the wrist while pulling so that it fits snugly.

- * Unless the wrist cuff is wrapped snugly, correct measurement may not be possible.


- * For the slender wrist (wrist circumference of less than 145 mm), the wrist cuff may not be wrapped snugly and give a feeling of loose wrap. In this case, do not force to wrap since the measurement is not affected.

- * Fasten the Velcro tape tightly.

- * The remaining part of the wrist cuff can be conveniently folded back out of the way.

When measuring on the right wrist

Measurement can also be made on the right wrist. Fit the monitor as shown in the Figure.

Preparation before Use

Align the wrist cuff with the level of your heart and gently support your left hand with your right hand.


- 1. Sit on a chair in correct posture with your back straight.**

* Do not bend your body forward.


- 2. Take deep breaths five to six times, then relax.**

- 3. Align the wrist cuff to the level of your heart. (Refer to the Figure to the right.)**

* Fingers touching the shoulder in this position.

• Measuring postures and hand shapes you should avoid

* If the wrist cuff is lowered or you hold your hand too tight, blood pressure values may be measured as higher than usual.


* Do not support the wrist cuff with your right hand. It will cause erroneous measurement.

How to Measure Blood Pressure


When it is ready to measure, take the measurement in the following procedure. To start measurement, simply press the START/STOP button.

1. Press the START/STOP button.

- (1) Press the START/STOP button.
Release the Button when the display turns on.


- (2) Initial display.
All display symbols turn on for about 0.5 seconds.


* The battery replacement mark  shown here does not mean batteries are worn out.

- (3) Deflation display.
A  mark starts to flash.

- (4) The monitor starts to inflate and measure.

* Because measurement is made during inflation, do not move your hand, wrist, or arm after the monitor starts to inflate.


These marks disappear one by one as the measurement proceeds.

How to Measure Blood Pressure

2. Measurement is over and blood pressure values and pulse are displayed.

When the measurement is over, a  mark flashes, and the cuff deflates automatically. Then a  mark turns on and blood pressure and pulse rate are displayed simultaneously.

Display of measurement results

* If a measurement cannot be made or the blood pressure values are displayed abnormally high (low), refer to Page 15.

* If the Velcro tape comes off during the measurement, wrap the wrist cuff snugly as shown in Page 7 and measure again after fastening the Velcro tape.


3. Turn off the power.

Blood pressure measurement is complete now.

Press the START/STOP button to turn off the power.

To make another measurement, repeat the procedure in Step 1 above.

* Even if you forget to turn off the power, the monitor will turn off automatically in approximately 2 minutes.


How to Use the Memory Function

The R5 Blood Pressure Monitor can record up to seven measurements.

1. Press the MEMORY button ("M" button).


The measurement values can be recalled after completion of the measurement or even when the power is off.

The latest measurements are displayed.


Press the "M" button.

Each time the "M" button is pressed, the previous measurement values are displayed.

Press the "M" button.

When all the stored measurement results are displayed, the display of the present measurement result or the first measurement result will resume (if the memory is recalled when the power is off.)

* Immediately after completion of the measurement, the memory may not be called even by pressing the "M" button. In that case, leave the finger from the "M" button and press it again to call the memory.


* You can recall up to seven measurements.

* Memory capacity is maintained in a "first-in-first-out" basis, ie the oldest readings will be deleted first.

How to Use the Memory Function

2. Press the START/STOP button to turn off the power.

After reading the blood pressure values and pulse rate stored in the monitor, press the START/STOP button to turn off the power.

* Even if you forget to turn off the power, the monitor will turn off automatically in approximately 2 minutes.

* The stored measurements will remain in the memory even when the batteries are removed for battery replacement.


• To delete the memory (Everything in the memory will be deleted.)

(1) Press the "M" button while the power is off to call the measured data.

(2) Press the START/STOP button for more than 2 seconds while pressing the "M" button.

* If you press the "M" button first, the power will turn off.


The display appears as shown to the right and the "M" button starts blinking.


(3) When the "M" mark starts blinking, stop pressing the "M" button.

(4) When the "M" mark stops blinking, the contents of the memory are deleted and the display of no memory appears.

(5) Press the START/STOP button to turn off the power.

Memory being deleted

Display of no memory

Maintenance

How to clean the monitor

Wipe the monitor with a cloth, moistened with water or detergent, then wipe it dry with a dry cloth.


* Please be careful so that water do not get into the unit.

Do not use benzine, thinner, or gasoline to clean the monitor.


Do not wash or wet the wrist cuff. Even if the cloth lining inside the wrist cuff becomes fluffy, it does not affect the measuring accuracy.


How to store the monitor

Keep the main unit in the storage case when carrying or storing the monitor.



If you are not going to use the monitor for a long period of time such as for more than three months, remove the batteries from the monitor to prevent leakage.

Error Indicators

If measurement is not made correctly, the following error will be displayed.

Error	Cause	How to correct
  <p>After several seconds</p>	You moved your arm or body, or talked during measurement.	Press the START/STOP button once to turn off the power. Measure your blood pressure again without moving your arm or body. (Refer to Page 9.)
  <p>After several seconds</p>	You moved your arm or body greatly during measurement.	
 <p>Only mark flashes.</p>	Batteries are worn out. *Display of battery replacement mark during the initial display does not mean worn batteries.	Replace the batteries with new ones. (Refer to 6.)

Troubleshooting

If trouble occurs while you are using the monitor, check the following points first.

Trouble	What to check	How to correct
Nothing is displayed when you press the START/STOP button.	Are the batteries worn out?	Replace the batteries with new ones. (Refer to Page 6.)
	Are the polarities of batteries wrong?	Insert the batteries with the polarities in correct direction. (Refer to Page 6.)
Measurement cannot be made. Or, blood pressure values are displayed abnormally high (low).	Did you keep the wrist cuff at the same level as your heart?	Measure your blood pressure in the correct posture. (Refer to Page 8.)
	Is the wrist cuff wrapped snugly?	Wrap the wrist cuff correctly. (Refer to Page 7)
	Are you putting too much stress on your shoulder or arm?	Relax, and measure. (Refer to Page 8.)
	Are you talking or moving your hands during measurement?	Keep still and measure again. (Refer to Page 8.)
The wrist cuff does not inflate to the required pressure and an E is displayed.	Is the wrist cuff leaking air?	Ask for the repair.
The blood pressure values differ from the brachial pressures.	After referring to the section on pages 17 to 18, consult the nearest OMRON dealer.	
Blood pressure value differs each time I measure or the blood pressure value is displayed abnormally high (low).	Blood pressure value varies according to the mental state during measurement or the measuring time. Take a deep breath, relax, then measure again. (Refer to Page 17 to 18.)	

* If you cannot measure your blood pressure correctly even after taking the above mentioned measures, the monitor may have malfunctioned.

* To confirm the problem or request for repair, consult the nearest OMRON dealer.

* In some very rare cases, some error may be caused by the physical natures. In such a case, please consult the doctor.

Specifications

Name	: OMRON Automatic Wrist Blood Pressure Monitor
Model	: R5
Display	: Digital display
Measurement	: Oscillometric method
Measurement Range	: Pressure; 0 to 280 mmHg Pulse rate; 40 to 180 beats/min
Accuracy	: Pressure; Within ± 4 mmHg Pulse rate; Within $\pm 5\%$ of reading
Inflation	: Automatic inflation by pumping
Deflation	: Automatic rapid deflation
Pressure Detection	: Electrostatic capacity semiconductor pressure sensor
Power supply	: Two AAA batteries (LR03)
Battery Life	: Approximately 400 measurements when using alkaline batteries at the room temperature of 22°C and by using three times a day and inflating to 170 mmHg
Operating Temperature and Humidity	: +10°C to +40°C, 30 to 85% RH
Storage temperature and humidity	: -20°C to +60°C, 10 to 95% RH
Measurable circumference of wrist	: 135 to 215 mm
Weight of Main Unit	: Approximately 130 g (not including batteries)
External Dimensions	: 54 (W) x 68 (H) x 32.5 (D) mm (not including the wrist cuff)
Electric Shock Protection	: Internal power supply appliance type B
Accessories	: Storage case, two AAA alkaline batteries, Instruction manual


* Specifications may be changed without prior notice.

What you should know about Blood Pressure

Blood pressure values measured at home tend to be lower than those measured at the doctor's office.

Blood pressure can vary by as much as 30 to 50 mmHg depending on whether you are nervous or relaxed.

Blood pressure fluctuations throughout the day. (Direct measuring method)

Provided by Dr. Osamu Tochikubo of the Second Internal Medicine,
the Faculty of Medicine at Yokohama Municipal University.

Blood pressure values measured at home are sometimes 25 to 30 mmHg lower than those measured at the doctor's office. This is caused by the fact that you may be nervous at the doctor's office, but are relaxed at home. It is important to know your regular normal blood pressure measured at home.

Factors that may cause the blood pressure to fluctuate


- Breathing • Exercise • Mental stress • Concern
- Environmental and/or temperature change • Eating
- Urination and bowel movement • Conversation • Taking bath
- Drinking alcohol • Smoking, etc.

Please understand that blood pressure fluctuates easily as shown above.

For the Correct Measurement of Blood Pressure

Let's Measure Blood Pressure Correctly.

It is important to observe the correct measuring method in order to measure blood pressure accurately.

Learn your blood pressure trend

In order to know the fluctuation trends of your blood pressure, it is important to measure your blood pressure at the same time every day at home and recording it.


Furthermore, do not worry about your blood pressure each time you measure. Ask your doctor to interpret your blood pressure data.

During the measurement, do not stare at the display but look at something at far.

Wrist blood pressure may differ from brachial pressure.

Wrist blood pressure is measured at a different region from the usual brachial blood pressure measurement and may show slightly different value. In many cases, the difference between the wrist and the brachial blood pressures is approximately ± 10 mmHg both for the systolic and the diastolic pressures in the case of healthy people.


People with the following condition should take note. When arterial sclerosis worsens because of diabetics, hyperlipemia, or hypertension, it may cause peripheral circulatory disturbance. In these cases, the wrist and brachial blood pressure values may differ greatly. If you have any question on the blood pressure values measured at wrist, please contact the nearest OMRON dealer.

Blood Pressure Classification

Classification of Blood Pressure by the World Health Organization

The World Health Organization (WHO) and the International Society of Hypertension (ISH) developed the Blood Pressure Classification shown in the Figure below. (This classification is made based on the blood pressure values measured in a sitting position, in the outpatient department of hospital.)

* There is no universally accepted definition of hypotension. However, those having the systolic pressure below 100 mmHg are assumed as hypotension.

According to the blood pressure classification by the WHO/ISH* (revised in 1999)

* ISH: International Society of Hypertension

The OMRON logo is centered on the page. It consists of the word "OMRON" in a bold, sans-serif typeface. The letter "O" is a solid black circle. The letters "M", "R", "O", "N" are formed by thin vertical lines, with the "M" and "R" having horizontal bars at their bases. The "O" is positioned between the "M" and "R".

OMRON

OMRON Corporation

3-4-10, Toranomon, Minato-Ku,
Tokyo 105-0001, Japan.

OMRON HEALTHCARE SINGAPORE PTE LTD.

83 Clemenceau Avenue,
#11-01, UE Square,
Singapore 239920.