TRUEread®

Owner's Booklet

24/7 Customer Consultation
English or Spanish
1-954-677-9201 • 1-800-803-6025
www.niprodiagnostics.com

Manufactured by:
NIPRO
DIAGNOSTICS™

Fort Lauderdale, FL 33309 U.S.A.

© 2010 Nipro Diagnostics, Inc. TRUEread and the Nipro Diagnostics logo are trademarks of Nipro Diagnostics, Inc.
G4NPD03 Rev. 20

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For quick reference only, not intended as a substitute for complete instructions. Please read entire Owner's Booklet and product Instructions for Use before testing.

2 APPLY BLOOD SAMPLE

1 INSERT TEST STRIP

Two simple steps

Welcome to the TRUEread® Blood Glucose Monitoring System

Congratulations on your purchase of the TRUEread Blood Glucose Monitoring System. You now have a very simple and accurate way to test your blood glucose (sugar) level, anytime, anywhere.

Our Commitment to You

Our goal is to provide you with quality healthcare products and dedicated customer service. If you have questions about using TRUEread products, visit our web site at: **www.niprodiagnostics.com**.

CAUTION!

***Please read complete Owner's Booklet and
all product Instructions for Use.***

Expected Results for people without diabetes¹:

Before eating	<110 mg/dL
Two hours after meals	<140 mg/dL
Bedtime	<120 mg/dL

Importance of Blood Glucose Monitoring

The more you know about diabetes, the better you are able to take care of yourself. A Doctor or Diabetes Healthcare Professional will determine your target range for your blood glucose results and how often you should test. Having the majority of your results within your target range helps slow or stop complications of diabetes.

Having most blood glucose results within your target range shows how well a treatment plan is working to control glucose levels. NEVER change your treatment plan without consulting with your Doctor or Diabetes Healthcare Professional.

Use of the TRUEread Blood Glucose Monitoring System in a manner not specified in this Owner's Booklet is not recommended and may affect your ability to determine your blood glucose.

IMPORTANT INFORMATION:

Health Related Information

Inaccurate results may occur in severely hypotensive individuals or patients in shock. Inaccurate results may occur for individuals experiencing a hyperglycemic – hyperosmolar state, with or without ketosis. Critically ill patients should not be tested with a blood glucose meter.²

- Low blood glucose (hypoglycemia) symptoms may be:
 - trembling, sweating, intense hunger, nervousness, weakness, and trouble speaking.
- High blood glucose (hyperglycemia) symptoms may be:
 - intense thirst, a need to urinate often, a dry mouth, vomiting, and headache.

If you have any of these symptoms, check your blood glucose. If your result does not match the way you feel, repeat the test. If your results still do not match the way you feel, call your Doctor or Healthcare Professional.

Important Information

For the most accurate results using TRUEread:

- **Read all instructions** before testing.
- **Check Codes.** Code in Meter Display must match Code on vial of TRUEread Test Strips being used and Code Number printed on Code Chip.
- **Meter displays results as plasma values.**
- **Use only TRUEread Test Strips and TRUEcontrol Glucose Control with TRUEread Meter.**
- **Do not use for diagnosis of diabetes or for testing blood glucose in neonates.**
- **The TRUEread is recommended for testing of human capillary whole blood only.** TRUEread is not recommended for use with venous samples.
- **Perform Check Strip and Control Tests before** performing a blood test for the first time.

***Note:** Three levels of TRUEcontrol® Glucose Control Solution are available for Quality Control Testing. We recommend testing at least 2 levels of Control. Contact your distributor to order different levels of Glucose Control Solution.*

- **The TRUEread is recommended for testing of human whole blood only.**

Health Related Information For Healthcare Professionals

- If comparing results between TRUEread and a lab system, perform a TRUEread blood test within 15 minutes of lab test. Results from TRUEread are considered accurate if within $\pm 20\%$ of lab results. If patient has recently eaten, fingerstick results from TRUEread can be up to 70 mg/dL higher than venous lab results.²

Table of Contents

Fast Test Guide	inside front cover
Introductory Information	1-6
Know Your System	
Meter	7-9
Test Strip	10
Test Strip Vial Label	11
Control Bottle Label	12
Getting Started	
.....	13
Coding	14-15
Quality Control Tests	16
Check Strip Test	16
Glucose Control Test.....	17-20
Testing your Blood	
How to Obtain a Blood Sample	21-22
How to Test Blood Glucose	23-24
System Out of Range Warning Messages	25
Set Up	
Time and Date	26-27
Memory	
Viewing Results in Memory	28
Care and Troubleshooting	
Caring for TRUEread	29-30
Changing Battery	31
Troubleshooting.....	32
Display Messages	33-35
Specifications	36
Chemical Composition.....	37
Meter Warranty	38
References.....	39

Know Your System


Meter

Front of Meter

Top of Meter

- ①  **Button** - Decrease numbers in Set Up; remove Control Symbol; move backward by time/date when viewing results in Memory.
- ②  **Button** - Select settings; view results in Memory.
- ③  **Button** - Increase numbers in Set Up; add Control Symbol; move forward by time/date when viewing results in Memory.
- ④ **Display** - Shows test results, messages, user prompts, and other information.
- ⑤ **Test Port** - Insert TRUEread Test Strip or Check Strip here.


- ⑥ **Data Port** - Contact Customer Service at the number located on back cover for more information.
- ⑦ **Code Chip Port** - Insert the Code Chip here.
- ⑧ **Battery Compartment** - Use one non-rechargeable 3V lithium battery (#CR2032), positive (“+”) side up.

Note: See Changing Battery for details.

- ⑨ **Meter Label** - Contains telephone number for assistance.
- ⑩ **Serial Number Label** - Identifies Meter when calling for assistance.


Display

1. Time, Date, Code.
 2. Result in Memory.
 3. Check Strip Test correct.
 4. Test result.
 5. Battery Symbol (see *Display Messages*).
 6. Temperature Symbol (see *Display Messages*).
 7. Apply Sample Symbol.
 8. Control Symbol (see *Glucose Control Test*).
 9. Units of Measure.
- Note:** (Factory Set. Cannot be changed by user. See Check Strip Test to verify unit of measure.)
10. Time is P.M.


Note: AM does not appear. “AM” or “PM” is not an option if factory set for a 24-hour clock.

Test Strip

- ① **Contact End** - Insert into Meter with Contacts (metallic blocks) facing up *before* adding sample.
- ② **Sample Tip** - Touch sample (blood or Glucose Control) to edge of Tip.

Placement of Sample to Test Strip


Correct

Incorrect

- Allow sample (blood or Glucose Control drop) to be drawn into Sample Tip.
- Do not smear or scrape drop with Strip.
- Do not apply more sample to Strip after removing from drop.

Note: *Do not insert Sample Tip with sample into Meter for testing*

Test Strip Vial Label

①	CODE	C1027	①
②	LOT	EA645A4	
③	EXP	2012/10/31	
	0	31-61 mg/dL	④
	1	87-117 mg/dL	
	2	247-335 mg/dL	
	May 30, 2012		③

- ① **Code** - Match Code in Meter Display with Code on Code Chip ***and*** Code on Test Strip vial label ***before*** testing. All Codes **MUST** match for accurate results.
- ② **Lot Number (LOT)** - Used for identification when calling for assistance
- ③ **Expiration Dates (EXP)** - Write date first opened on vial label. Discard vial and unused Test Strips 120 days after written date or date printed next to EXP, whichever comes first.
- ④ **Control Level Range** - Range of numbers in which Control Test result must fall to assure System is working properly.

Caution! *Use of Test Strips or Glucose Control past the Expiration Dates may give incorrect test results. Discard out-of-date products and test with new products.*

Control Bottle Label

- ③ **Lot Number (LOT)** - Used for identification when calling for assistance.
- ② **Expiration Dates (EXP)** - Write date first opened on bottle label. Discard bottle if either 3 months after opening or date printed next to EXP on bottle label has passed.
- ③ **Control Level (0,1, or 2)** - Three levels of TRUEcontrol Control are available. We recommend testing at least 2 levels of Control. Use the number on the front cover for information on how to obtain different levels of Control.

Getting Started

The Meter turns on automatically when a Check Strip or Test Strip Contact End is inserted into the Test Port or when  Button is pressed (*see Set Up and Memory*).

The Meter turns off automatically after the Check Strip or Test Strip is removed from the Test Port or after 2 minutes of non-use.

Meter comes with pre-set time and date. Before using the Meter for the first time or after a battery change, check the time and date and update as needed (*see Set Up*).

We recommend performing Quality Control Tests (Check Strip Test, Control Test) for practice before using your Meter for the first time to test your blood.

Quality Control Tests should be performed:

- Occasionally as you use the vial of Strips,
- When opening a new vial of Strips,
- When Meter results seem unusually low or high, or
- Meter was dropped.

Coding

Code Chip is packaged with each box of Test Strips.

Note: *After all Test Strips are used from box, remove Code Chip from Meter and discard. Open new box of Test Strips and insert new Code Chip into Meter before testing.*

1. Start with Meter off. Remove and open Code Chip package from box of Test Strips to be used. Match Code printed on Code Chip to Code printed on vial label.

Caution! *If Codes do not match, do not use Test Strips. Call for assistance.*


2. If Codes match, insert Code Chip into Code Chip Port as shown, locking in place.

3. Insert Test Strip into Test Port.
Meter turns on. “Code” and number
appear in Display.

4. Match to Code number
on vial label. If Codes
match, begin testing.

Note: Letter printed before number on the Vial
Label Code will not appear in Meter
Display.


Caution! If “Code ----” shows in Display,
reinsert Code Chip and then reinsert
Test Strip. If “Code ----” still shows in
Display or if Codes do not match, **DO**
NOT use Test Strips. Test results could
be incorrect. Call for assistance.

Quality Control Tests

Check Strip Test

Check Strip is stored in Carrying Case. Check Strip Test is used to check that Meter is working properly.

1. Start with Meter off. Remove Check Strip from Carrying Case. Insert Check Strip Contact End into Test Port with dimple facing up. Meter turns on.
2. Check Display to make sure that all segments are showing (see Display).
3. After a few seconds, the Display shows OK✓, Apply Sample Symbol, and mg/dL.

Caution! *If segments are missing from Display, OK✓, the Apply Sample Symbol, and mg/dL does not appear, do not use Meter for testing. Call for assistance.*

4. Remove Check Strip. Return to Carrying Case. Meter turns off.

Glucose Control Test

Glucose Control Test uses TRUEcontrol Glucose Control to check performance of System.

Three levels of TRUEcontrol Glucose Control Solution are available for Quality Control Testing.

It is important to perform Control Tests with more than one level of Control to assure your System is working properly and your testing technique is good. Contact place of purchase or call for assistance to obtain Controls.

Caution! *Control ranges printed on Test Strip vial label are for Control Test results only and **are not** suggested levels for your blood glucose.*



Only the most current Control Test result is stored in Memory.

How to Test Glucose Control

1. Wash hands with soap and water, dry thoroughly. Allow Control, vial of Test Strips and Meter to adjust to room temperature (68-77°F). Write date first opened on Control Label and/or Test Strip vial label.
2. Check dates on Control label and Test Strip vial label. Do not use Control if 3 months past written opened date *or* past date printed next to EXP. Do not use Test Strips if 120 days past written opened date *or* past date printed next to EXP. Discard out-of-date products and use new products if either date has passed.
3. Swirl or invert bottle gently to mix Control. **DO NOT SHAKE!**
4. Remove Test Strip from vial. Recap vial immediately.

Note: Use Test Strip quickly after removal from vial. Test Strips that have been left out of vial too long will give error message when used. Discard old Test Strip and test with new Test Strip.


5. Insert Test Strip Contact End (Contacts facing up) into Test Port. Meter turns on. Keep Strip in Meter until testing is finished.

6. Match Code in Display with Code on Test Strip vial label being used. If Codes do not match, see *Coding*.

7. Press  . Control Symbol appears in Meter Display.

Note: *If you decide not to perform a Control Test, press  to remove Control Symbol.*


8. Turn Control bottle upside down. Squeeze one drop of Control onto a clean tissue. Wipe off bottle tip.

9. Gently squeeze a drop onto a small piece of unused aluminum foil or clear plastic wrap.

10. With Strip still in Meter, touch edge of Sample Tip to drop of Control and allow drop to be drawn into Test Strip. Remove Test Strip from drop when Meter beeps. Meter begins countdown


Note: *If Meter does not beep or begin countdown after placing Sample Tip to Control drop, discard Test Strip. Repeat test with new Test Strip. If problem persists, see Troubleshooting.*

11. Compare result to Control range printed on Test Strip vial label.

If result is in range,
System can be used for
testing blood. If result
does not fall within range,
repeat test using a new
Test Strip.

0	31-61	mg/dL
1	87-117	mg/dL
2	247-335	mg/dL

Caution! *If Control Test result continues to be outside range, System should not be used for testing blood. Call for assistance.*

12. Remove Test Strip from Meter and discard.
Meter turns off. Recap Control bottle tightly.


How to Obtain a Blood Sample

Refer to Lancing Device Instructions for Use for detailed instructions.

- Never share lancets or lancing device.
- Lancets are for single use only. Do not re-use.

From Fingertip

1. Prepare fingertip by washing hands in warm, soapy water. Rinse well. Dry thoroughly.
2. Place end of Lancing Device against tip of finger. Lance fingertip.
3. Set Lancing Device aside. To help blood drop form, lower hand to waist level, gently massaging finger from palm to fingertip. Allow blood drop to form before attempting to apply to Test Strip.


Always remove and discard used Lancet in appropriate container when testing is complete.

Note: *Used Strips and lancets are considered biohazardous. Dispose used Strips and lancets in approved biohazard container*

From Forearm

1. Select area to be lanced. Wash with soap and warm water, rinse and dry thoroughly.
2. Rub area vigorously or apply a warm compress to increase blood flow.
3. Place end of Lancing Device firmly against forearm. Press trigger button. Apply firm pressure on lancing device for 10 seconds.

Note: *Some lancing devices include a special end cap for alternate site testing. Check lancing device Instructions for Use.*

Important Notes Regarding Forearm Testing³

- Check with your Doctor or Diabetes Healthcare Professional to see if forearm testing is right for you.
- Results from forearm are not always the same as results from fingertip.
- Use fingertip instead of forearm for more accurate results:
 - Within 2 hours of eating, exercise, or taking insulin,
 - If your blood sugar may be rising or falling rapidly or your routine results are often fluctuating,
 - If you are ill or under stress,
 - If your forearm test results do not match how you feel,
 - If your blood sugar may be low or high,
 - If you do not notice symptoms when blood sugar is low or high.

How to Test Blood Glucose


Check supplies. Always check your supplies before using. Check meter for damage (damaged or cracked display, missing or damaged buttons). If any damage is seen, do not use Meter. Call for assistance. Check Strip vials for damage. Discard any vials that appear cracked or broken.

1. Check dates on Test Strip vial being used. Do not use if 120 days past written date or date next to EXP, whichever comes first
2. Wash hands (and forearm for alternate site testing). Rinse well and dry thoroughly.
3. Remove Test Strip from vial. Recap vial immediately.

Note: Use Test Strip quickly after removal from vial. Test Strips that have been left out of vial too long will give error message when used. Discard old Test Strip and use new Test Strip for testing.

4. With Meter off, insert Test Strip Contact End (Contacts facing up) into Test Port.
Meter turns on. Keep Test Strip in Meter until testing is finished.
5. Match Code in Display to Code on Test Strip vial being used. If Codes do not match, see *Coding*. Apply Sample Symbol appears in Display.

Note: Letter printed before the number on Code Chip and vial label code will not appear in Meter Display.

6. Lance fingertip or forearm. Allow drop to form.
7. With Strip still in Meter, touch edge of Sample Tip to blood drop and hold. Allow blood to be drawn into Strip until Meter beeps.

Note: *If Meter did not beep or begin countdown after touching blood drop to Sample Tip, discard Test Strip. Repeat test with new Test Strip and new blood drop. If problem persists, see Troubleshooting.*

8. After countdown, result is displayed with date and time. Remove Test Strip and discard. Result is stored in Memory. Meter turns off.

9. Record result in log book.

Note: *Used Strips and lancets are considered biohazardous. Dispose used Strips and lancets in approved biohazard container.*

System Out of Range Warning Messages

WARNING!

TRUEread measures blood glucose results from 20-600 mg/dL.

If blood glucose result is less than 20 mg/dL, “**LO**” appears in Meter Display.

If blood glucose result is greater than 600 mg/dL, “**HI**” appears in Meter Display.

Always repeat test to confirm Low (“**LO**”) and High (“**HI**”) results. If results still display “**LO**” or “**HI**”, call your Doctor or Diabetes Healthcare Professional ***immediately***.

***Note:** “**LO**” results are included in the Morning Average as 20 mg/dL. “**HI**” results are included as 600 mg/dL.*

Set Up

1. Start with Meter off. Do not insert Test Strip.

Press and hold . Date/Time appears in display. Release after full display appears (*shown at right*) and Meter beeps.


Full Display

Set Time

2. The hour flashes.

To change, press  or  on top of the Meter to select the hour. Like many alarm clocks, to set “PM”, scroll through hour until “PM” appears in Display. Press  to set.

Hour


Note: “PM” does not display if Meter is factory set to 24-hour clock.

AM / PM

3. The minutes flash.


To change, press  or  to select the minutes. Press  to set.

Minutes

Set Date

4. The month (number) flashes.
To change, press **+** or **-** to select the month. Press **S** to set.
5. The day (number) flashes.
To change, press **+** or **-** to select the day. Press **S** to set.
6. The year flashes.
To change, press **+** or **-** to select the year. Press **S** to set.

Exit Set-up Mode

To review settings, press **S**. Make changes as needed, per instructions above.

Press and hold **S** until Display goes blank.

Options are saved. Meter turns off.

Viewing Results in Memory

Memory stores 200 results, which are displayed from most recent to oldest. The oldest result is removed from Memory when Memory is full and a new result is added.

1. Start with Meter off. Press and release . Meter displays most recent result in Memory.
2. Press  and release to advance to the first blood glucose test result in Memory.
3. Press  and release again to view the latest Control Test result. The Control Test result shows the Control Symbol in Display. Only one Control Test result is stored in Memory.
4. Continue to press and release the  to advance through the blood glucose results. Holding  scrolls through the results quickly. Press and release  to reverse through the results.
5. Press and hold the  to turn Meter off. Meter also turns off automatically after 2 minutes of inactivity.

Control
Symbol

Caring for TRUEread

- Store System (Meter, Check Strip, TRUEcontrol Glucose Control, Test Strips, Code Chip) in Carrying Case to protect from liquids, dust and dirt.
- Store System in a dry place at room temperature (36°-86°F). **DO NOT REFRIGERATE OR FREEZE.**

Meter Care

- Wipe Meter with clean, lint-free cloth dampened with one of the following:
 - Mild detergent or mild soap and water
 - 10% household bleach and water
- Do not use alcohol to clean. Cleaning Meter with alcohol will cause damage.
- Never put Meter in liquids or allow any liquids into Test and Data Ports.

Check Strip Care

- Do not bend Check Strip.
- Do not apply blood, Glucose Control or any other liquid to Check Strip.
- To clean, wipe with damp cloth, dry thoroughly.
- Do not write on Check Strip.

TRUEcontrol Glucose Control Care

- Write date opened on Control label. Discard 3 months after opening or after date printed next to EXP, whichever comes first.
- After use, wipe bottle tip clean and recap tightly.
- Swirl bottle to mix, **do not** shake.
- Store at room temperature (36°-86°F). **DO NOT REFRIGERATE OR FREEZE.**

Test Strip Care


- Store Test Strips in original vial only. Do not transfer old Test Strips to new vial.
- Write date opened on Test Strip vial. Discard unused Test Strips from vial 120 days after opening or after date printed next to EXP, whichever comes first. Use of Test Strips past either date may give incorrect results.
- Close vial immediately after removing Test Strip. Use Test Strip quickly after removal from vial. Never store Test Strips outside of vial. Store in a dry place at room temperature below 86°F.
DO NOT REFRIGERATE OR FREEZE.
- Do not reuse Test Strip.
- Do not bend, cut or alter Test Strips in any way.

Code Chip Care

- Code Chip must remain in Meter until all Test Strips from box are used, then discard old Code Chip.
- Insert new Code Chip from new box of Test Strips immediately before using Test Strips.

Changing Battery

Replace battery when Low or Dead Battery Symbols appear in Display or Meter does not turn on.

1. Lift tab on Battery Cover.
2. Turn Meter over, tap gently to loosen and remove battery.
3. Insert new battery, positive (“+”) side up. Close cover.

Note: Use non-rechargeable 3V lithium battery (#CR2032).

4. Discard old battery in appropriate container.
5. Turn Meter on. If Meter will not turn on, check that battery was installed properly. If not, reinsert battery and try again. Call for assistance if problem persists.

Caution! Batteries might explode if mishandled or incorrectly replaced. Do not dispose of battery in fire. Do not disassemble or attempt to recharge battery. Dispose according to local/country specific regulations.

Note: Replacing battery may affect date and time settings. Check date and time by going to Set Up and scrolling through the time/date settings by pressing .

Troubleshooting


1) After inserting Test Strip, Meter does not turn on.

Reason	Action
Test Strip inserted upside down or backwards	Remove Strip. Re-insert correctly.
Test Strip not fully inserted	Remove Strip. Re-insert Strip fully into Meter.
Test Strip Error	Repeat with new Strip.
Dead or no battery	Replace battery.
Battery in backwards	Battery positive (“+”) side must face up.
Meter Error	Call for assistance.

2) After applying sample, test does not start/ Meter does not beep or begin countdown.

Reason	Action
Sample drop too small	Repeat test with new Strip and larger drop.
Sample applied after two minute shut-off	Repeat test with new Strip and apply drop after Meter beeps.
Problem with Test Strip	Repeat with new Strip.
Problem with Meter	Call for assistance.




Messages

<u>Display</u>	<u>Reason</u>	<u>Action</u>
	Temperature Error Temperature change too quick	Wait 10 minutes for Meter and Test Strips to reach room temperature before testing.
	Temperature Too Cold Meter temperature less than 50°F	Move Meter and Test Strips to area between 50°-104° before testing.
	Temperature Too Hot Meter temperature greater than 104°F	Move Meter and Test Strips to area between 50°-104° before testing.
	Sample Not Detected or Using Wrong Test Strip	Retest with new or correct Strip and larger sample.
	Used Test Strip or Test Strip outside of vial too long.	Repeat with new Strip. If error persists, call for assistance.

Messages *(continued)*

<u>Display</u>	<u>Reason</u>	<u>Action</u>
	Meter Error	Call for assistance.
	Test Strip Error	Retest with new Strip. If error persists, call for assistance.
	Test Strip Removed During Test	Retest with new Strip. If error persists, call for assistance.
	Check Strip Error	Repeat test. If error persists DO NOT use Meter for testing. Call for assistance.
	Code Chip Error	Remove Code Chip and reinsert. If error persists, call for assistance.

Messages *(continued)*

<u>Display</u>	<u>Reason</u>	<u>Action</u>
	Low Battery	About 50 tests can be done before battery must be replaced.
	Dead Battery	Replace battery.
	Code Chip Error	Remove Code Chip and re-insert. If error persists, call for assistance.
	Communication Error	Turn Meter off and on again. If error persists, call for assistance.
 	• <u>WARNING!!</u> Out of range - HI Result > 600 mg/dL Out of range - LO Result < 20 mg/dL	• <u>WARNING!!</u> Retest with new Strip. If result is still "HI" or "LO" contact Doctor <i>immediately.</i>

Specifications

Result Range: 20-600 mg/dL

Sample Size: Minimum 1 microliter (1 μ L)

Sample: Human capillary whole blood or
TRUEcontrol Glucose Control

Test Time: 10 seconds

Result Value: Plasma values

Assay Method: Electrochemical

Reference Method: Yellow Springs Instrument (YSI)

Power Supply: One 3V lithium battery
#CR2032 (non-rechargeable)
Total power when active
at full battery = 8.6mW

Battery Life: Approximately 1,100 tests or 1 year

Automatic shut-off: After two minutes of non-use

Weight: 1.66 ounces

Size: 3.52" x 2.15" x 0.67"

Memory Size: 200 blood glucose results,
1 Control Test result

System Operating Range (Meter & Test Strips):

Relative Humidity: 10-90% (Non-condensing)

Temperature: 50°-104°F

Hematocrit: 30-55%

Note: *Use within specified environmental conditions only.*

Chemical Composition

Test Strips: Glucose Oxidase (*Aspergillus sp.*) 2.5 Units, Mediators, Buffers and Stabilizers.

TRUEcontrol Glucose Control: Contents:
Volume: 3 mL. Water - 73%, D-glucose - 0.09 - 0.20%, Viscosity Enhancing Agent - 25%, Inorganic Salts - 1.8%, Amaranth - 0.08%, and Preservatives - 0.03%.

TRUEread Limited Lifetime Warranty

Nipro Diagnostics, Inc. provides the following Warranty to the original retail purchaser of the TRUEread Meter:

- 1) Nipro Diagnostics Inc. warrants this Meter to be free of defects in materials and workmanship at the time of purchase. If the meter is inoperative ever, Nipro Diagnostics, Inc will replace the Meter with an equivalent Meter, at its option, at no cost to the purchaser. Failure of the meter due to abuse or use not in accordance with the instructions for use is not covered by this Warranty.
- 2) This Warranty does not include the battery supplied with the Meter.
- 3) Do not take the Meter apart. This action will void the Warranty and cause the Meter to display false results.
- 4) The duration of any implied Warranty, including any implied Warranty of merchantability or fitness for a particular purpose shall be limited to the lifetime in use with original user in accordance with any state law to the contrary.
- 5) Nipro Diagnostics, Inc. disclaims liability for incidental or consequential damages for breach of any expressed or implied Warranty, including any implied Warranty of merchantability or fitness for a particular use with respect to the Meter. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitations or exclusion may not apply.
- 6) This Warranty gives the user specific legal rights, and the user may also have other rights which vary state to state.

Your Nipro Diagnostics, Inc. Customer Care Representative will be able to provide detailed information regarding procedures for returning your Meter, if necessary.

References

1. Joslin Diabetes Center, *Goals for Blood Glucose Control* [Electronic Version], Retrieved June 2, 2008.
www.joslin.org/Beginners_guide.523.asp.
2. Food & Drug Administration 2007. *Draft Guidance for Industry and Staff. Total Product Life Cycle for Portable Invasive Blood Glucose Monitoring System*. [Electronic Version]. Retrieved February 5, 2007 from www.fda.gov/cdrh/oivd/guidance/1603.pdf
3. Larson-Conn U: *Difference between capillary and venous blood glucose during oral glucose tolerance tests*. Scand J Clin Lab Invest 36:805-808, 1976.
4. Food and Drug Administration 2003. *Diabetes Information - Glucose Meters and Diabetes Management*. [Electronic Version]. Retrieved August 11, 2003 from www.fda.gov/diabetes/glucose.html.

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page. There are approximately 20 lines visible. The paper has a slight shadow on its right side, suggesting it's resting on a surface.

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Nipro Diagnostics, Inc. ofrece la siguiente garantía al comprador minorista original del medidor TRUEread: