

Tour 1008™ User Manual



YOUR IMAGE IS OUR BUSINESS.



a VCA ANTECH company

Important Safety Information

Warning



*There are no user serviceable parts inside.
Refer all servicing to qualified service personnel.*

To reduce the risk of fire or electric shock, do not expose the unit to moisture or water. Do not allow foreign objects to get into the enclosure. If the unit is exposed to moisture, or a foreign object gets into the enclosure, immediately disconnect the power cord from the wall. Take the unit to a qualified service person for inspection and necessary repairs.

Caution



*Read all the instructions before connecting or operating the component.
Keep this manual so you can refer to these safety instructions.*

Heed all warnings and safety information in these instructions and on the product itself. Follow all operating instructions.

Clean the enclosure only with a dry cloth or a vacuum cleaner.

Clean the DR plate with a damp cloth. Use no chemicals or alcohol.

You must allow 10 cm or 4 inches of unobstructed clearance around the unit. Do not place the unit on a bed, sofa, rug, or similar surface that could block the ventilation slots. If the component is placed in a bookcase or cabinet, there must be ventilation of the cabinet to allow proper cooling. Keep the component away from radiators, heat registers, stoves, or any other appliance that produces heat. Keep the component away from flammable materials.

The unit must be connected to a power supply only of the type and voltage specified on the side panel of the unit.

Connect the component to the power outlet only with the supplied power supply cable or an exact equivalent. Do not modify the supplied cable in any way. Do not attempt to defeat grounding and/or polarization provisions. Do not use extension cords.

Do not route the power cord where it will be crushed, pinched, bent at severe angles, exposed to heat, or damaged in any way. Pay particular attention to the power cord at the plug and where it exits the back of the unit.

Important Safety Information *(continued)*

The power cord should be unplugged from the wall outlet if the unit is to be left unused for a long period of time. Immediately stop using the component and have it inspected and/or serviced by a qualified service agency if:

- *The power supply cord or plug has been damaged.*
- *Objects have fallen or liquid has been spilled into the unit.*
- *The unit has been exposed to rain.*
- *The unit shows signs of improper operation*
- *The unit has been dropped or damaged in any way*

Place the unit on a fixed, level surface strong enough to support its weight. Do not place it on a moveable cart that could tip over.

System will become unresponsive when exposed to Electrostatic Discharge (ESD) or Electrical Fast Transient (EFT) conditions.

This is the normal response to these environmental conditions.

Please restart the system to return to normal operation.

Specifications

Power Requirements

USA: 115 Volts, 60 Hz

Europe: 230 Volts, 50 Hz

Power Consumption

200 Watts

Dimensions (W x H x D)

20 x 14x 10½ inches, 508 x 355 x 267 mm

Weight (net)

37 lb., 16.8 kg



CAUTION

Risk of electric shock.
Do not open.

Caution: To reduce the risk of electric shock, do not remove cover. No user-serviceable parts inside. Refer servicing to qualified service personnel.



This symbol is to alert the user to important operating and maintenance (service) instructions in this manual and literature accompanying the product.

Applicable for USA, Canada or where approved for the usage.

Caution: To prevent electric shock, match wide blade of plug to wide slot. Insert fully.



This symbol is to alert the user to the presence of uninsulated dangerous voltages inside the product's enclosure that may constitute a risk of electric shock.

Table of Contents

Tour 1008™ User Manual

Important Safety Information	i
Warning	i
Caution	i
Specifications.....	ii
Chapter 1: Tour 1008™ System Setup.....	1
Opening the Tour 1008™ Case	4
Connecting the Power and Sync Cables	5
Connecting the Generator	9
Powering Up the Tour 1008™ System	11
Chapter 2: Opening the Tour 1008™ Software.....	13
Chapter 3: Entering Patient Information	14
Searching and Adding Patients	14
Editing an Existing Patient's Information	16
Chapter 4: Creating Studies	17
Creating a New Study for an Existing Patient	17
Adding New Images to an Existing Study.....	18
Chapter 5: Acquiring Radiographs	19
Choosing Anatomy and Positioning	19
Acquiring a Radiograph	20
Cropping a Radiograph	22
Chapter 6: Managing Protocols	24
Using Protocols	24
Adding New Protocols.....	27
Editing Existing Protocols	31
Renaming Protocols	32
Deleting Protocols	34
Chapter 7: Reviewing Images	35
Reviewing a Study	35
Chapter 8: Exporting Studies and Images	38
Burning a CD	38
Exporting Images to a Folder.....	41
Transferring Studies to the Storage Station	44
Transferring Individual Studies to the Storage Station.....	46
Confirm Storage of a Study	47
Chapter 9: Adding Users	48
Chapter 10: Adding and Editing a DICOM Location.....	50
Chapter 11: Gain Calibration	53

Chapter 1: Tour 1008™ System Setup

The Tour 1008™ works seamlessly with the X-ray generator to provide the best image quality in a small, lightweight, self-contained system. The Tour 1008™ portable Acquisition Stations are intended to be used by Veterinarians for intermittent mobile diagnostic applications. They will connect the Acquisition laptop in the case to a Varian Amorphous Silicon Digital X-Ray Detector (2520E+).



Chapter 1: Tour 1008™ System *(continued)*

The system comes with several components that integrate the X-ray machine, plate, and the computer. If you need to order replacement parts, please refer to the part numbers in parenthesis.



Tour 1008™ Transit Case (70-309)

Wheel Kit for Tour 1008™ case (30-580):

Carries and protects the computer, plate, and plate cable. Wheel kit allows for easy transport in airports and on smooth surfaces.



2520E+ Plate (91-402):

Captures digital radiographs.



Handswitch (90-412):

Syncs the generator's prep and expose cycle with the computer.



Ethernet Cable (90-403):

Connects the computer to the network to allow for easy image transfers from the Acquisition Station to the Storage Station.

Chapter 1: Tour 1008™ System *(continued)*



ODU Plate Cable (10-835):

Provides power to the plate and syncs the plate and Tour 1008™ Transit System.



Generator Power Sync Cable (10-748):

Provides power to the X-ray machine and syncs the X-ray machine with the computer. It comes standard with a Min-X HF 80 connector. Other connectors (Medison, Poskom, Vet X-Ray Ultra, Ajex 90-15HF, Sternes Icon 90/20, and the Porta 1030 & 8020 generators) are also available.



Power Cable (10-171):

Connects the Tour 1008™ to any standard U.S. 110W power outlet (Nema 5-15).

Chapter 1: Tour 1008™ System *(continued)*

Opening the Tour 1008™ Case

1. Press in on the metal pin on each latch to open the Tour 1008™ Case.
With the pin pressed, lift the latch to unlock the case.



2. Open the plate door by pulling up on the silver tab and lifting the cover back.




Chapter 1: Tour 1008™ System (continued)

Connecting the Power and Sync Cables

1. Remove the Plate and the Plate Cable.



The plate cable should always remain attached to the plate.



2. Close the Plate storage door and place the plate in the Plate Park in the center of the door.



Chapter 1: Tour 1008™ System

Connecting the Power and Sync Cables (*continued*)

3. Twist the knob counter clockwise on the Connector Panel door to reveal the connection ports.



Power Cable

Generator Power Sync Cable

ODU Plate Cable



Chapter 1: Tour 1008™ System

Connecting the Power and Sync Cables (*continued*)

4. Connect the ODU Plate Cable by aligning the red dot and notch on the cable with the red dot and notch on the receptacle.



Chapter 1: Tour 1008™ System

Connecting the Power and Sync Cables (*continued*)

5. Connect plate cable to the plate if it is not already connected.
Twist the lock ring clockwise until it locks in place.



Chapter 1: Tour 1008™ System (continued)

Connecting the Generator



The power sync cable should stay connected to the generator and should be stored in the generator box.

1. Plug the Generator Power Sync Cable into the generator by connecting the power adapter to the power port on the generator.



2. Plug the hand switch end of the Generator Power Sync Cable into the generator's hand switch port.



If using the system Internationally or on a 220 power, any device plugged into the side power must be able to handle 220.

Chapter 1: Tour 1008™ System

Connecting the Generator (continued)

3. Plug the Lemo cable on the Generator Power Sync Cable into the case by lining up the red dot on the Lemo cable to the top of the connector. To release the connector, pull out on the ring.



4. Plug the power cable into the case.



Chapter 1: Tour 1008™ System *(continued)*

Powering Up the Tour 1008™ System

1. Plug the female end of the power cable into the power cable slot on the Tour 1008™ box. Plug the other end into a grounded power outlet (see page ii Power Requirements).



2. With the cables connected, turn on the Tour 1008™ by pressing the white ON switch to "I" (O = off, I = on) located in the upper left corner. The battery light on the laptop will light up, showing power to the unit.



Chapter 1: Tour 1008™ System

Powering Up the Tour 1008™ System (*continued*)

3. Open the computer by sliding the gray clip on the front edge to the right and lifting the lid.



4. Turn on the laptop by pressing the Power button  on the keyboard.



Chapter 2:

Opening the Tour 1008™ Software

Open the Tour 1008™ software to begin the process of taking a radiograph.

1. Double-click the TruDRlx™ icon  on the desktop. This will load the TruDRlx™ software for acquiring digital radiographs.
2. Enter the user name and password when the login screen appears.


The default login is:

User Name: admin

Password: admin



The software will load directly into the Patient Screen if the “Remember Me” box is checked.

A screenshot of the TruDRlx login screen. On the left is the logo for 'tru|DR™ DIGITAL RADIOGRAPHY'. Below the logo, it says 'To begin, enter your user name and password'. On the right, there are two input fields: 'User Name:' with 'tech' entered, and 'Password:' with four dots. To the right of the password field is a 'Login' button. Below the password field is a checkbox labeled 'Remember Me' which is checked.

tru|DR™
DIGITAL RADIOGRAPHY

To begin, enter your user name and password

User Name: tech

Password: ••••

☒ Remember Me

Login

Chapter 3: Entering Patient Information

Searching and Adding Patients


Search for existing patients or add new patients easily by using the filtered search fields.

1. Search for a patient in the Patient Screen by entering the **Patient ID**, the **Patient Name**, or the **Owner Last Name**. The patient list will automatically filter as you type.



Patient list is displayed based on the User Selected Default.

2. If you do not see the patient you are searching for, click the **Add** button  or press the **Enter** key on the keyboard to bring up the Add Patient Screen.



3. Enter the patient's information. **Patient Name** and **Last Name** are required fields.

Chapter 3: Entering Patient Information

Searching and Adding Patients (*continued*)

- Click one of the Save buttons.



Save

Saves the current patient and returns you to the Patient Screen.



Save+ Add

Saves the current patient and returns to the Add Patient Screen.



Save+ Acquire

Saves the current patient then opens the Anatomy/Protocol Screen.



Save+ Protocol

Saves the patient and brings you to the Protocol Screen to select the desired series.

The screenshot shows the Sound-Eklin software interface. At the top, there is a toolbar with five buttons: 'patient', 'save', 'save+add', 'save+acquire', and 'save+protocol'. The 'save' button is currently selected. Below the toolbar, the 'Patient Information' form is displayed. The form contains the following fields:

- Patient Name: Trigger
- Date Of Birth: 01/18/2001
- Breed: Thoroughbred
- Sex: Stud
- Owner Information:
 - Last Name: Smith
 - First Name: Sammie
 - Institution Name: Vandelay Farms
 - Address: 5817 Dryden Place
 - Address2: Suite 101
 - Address3:
 - City: Carlsbad, CA 92008

At the bottom of the form, there is a red asterisk and the text 'Required Information'.

Chapter 3: Entering Patient Information *(continued)*


Editing an Existing Patient's Information

Edit an existing patient's information for future DICOM tagging.



Previous exam information cannot be changed

1. Select a patient to edit from the Patient Screen.
2. Click the **Edit** button. 



Patient ID	Patient	Mod Date	Owner
113009085810	Sammie	11/30/09 08:59 P	Smith

Study Date	Images	Locked	Description	Archive	Sent
11/30/2009 8:59:00 PM	2	☒	Protocol		

3. Edit the information as needed.



If the patient is from an owner previously entered, owner can be selected from owner list on the left hand side.


4. Click the appropriate **Save** button (see page 15).

Chapter 4: Creating Studies

Creating a New Study for an Existing Patient

Create a new study for a patient previously entered in TruDRlx™ .

1. Search for a patient in the Patient Screen by entering the **Patient ID**, the **Patient Name**, or the **Owner Last Name**. The Patient List will automatically filter as you type.
2. Select a patient by clicking on the name in the Patient List.
3. Click the **New** button  to create a new study.
4. Select the anatomy and the view from the Anatomy Screen (see page 19).




Chapter 4: Creating Studies (continued)


Adding New Images to an Existing Study

Add new images to an existing patient's study.

1. Search for a patient in the Patient Screen by entering the **Patient ID**, the **Patient Name**, or the **Owner Last Name**. The Patient List will automatically filter as you type.
2. Click the Acquire button. 



This will take you to the Views Screen using the anatomy you selected from the original study. This is shown in the description field.




Chapter 5: Acquiring Radiographs

Choosing Anatomy and Positioning


Choose the anatomy and positioning to automatically enhance and position the radiograph.

1. Select the Tech and ordering Vet from the drop-down menus.



Tech: Davola, Joe
Vet: Bania, Kenny

2. Move the mouse over the model to display anatomical region of interest and select the appropriate anatomy.



3. Click on the appropriate button to select the view position.
This will bring you to the acquisition screen.



Chapter 5: Acquiring Radiographs (continued)

Acquiring a Radiograph

1. Set the X-ray generator according to the Sound-Eklin® Technique Chart using the appropriate KVP and mAs (depending on the X-ray generator).



Due to variations between X-ray generators, optimal techniques may vary somewhat from the values printed on the Sound-Eklin® Technique Chart.

Users should be fully acquainted with State and local regulations governing radiation protection and the use of diagnostic X-ray equipment as well as the manual for the X-Ray Generator being used.

2. Align the generator 22 inches from the plate using the light beam to fill the entire capture area defined on the plate.



If the generator is aligned perpendicular, the light will form a rectangle with even sides filling the capture area.

The image capture area is highlighted by a white box with the dotted cross-hairs in the middle.



3. From the Acquisition Screen press the hand switch halfway down to prep and all the way down to expose the X-ray machine.



Proper Two-Stage Hand Switch Technique


1. To Prep, push the hand switch halfway down.
2. Wait until the Generator is ready to fire (i.e. prep light on the X-ray machine is illuminated, etc.).
3. To expose, push the hand switch completely down.

Chapter 5: Acquiring Radiographs

Acquiring a Radiograph (*continued*)



4. The image will appear in less than 3 seconds. Use the **Left**, **Right**, **Flip**, and **Reverse** buttons to make orientation adjustments before shooting the next shot.



Crop the images as needed (see Pages 22).

If you need to reshoot your X-ray, press the hand switch and fire again.

5. Click the **Anatomy** button  to change views.
6. Select the next view from the Views Screen. Click the **Anatomy** button from the Views Screen if another anatomy is needed (see page 19).
7. Click the **Patient** button  to return to the Patient Screen or click the **Review** button. 

Chapter 5: Acquiring Radiographs (continued)

Cropping a Radiograph

Remove unwanted portions of the radiograph.


1. Click the Crop button  in the Acquisition Screen.



*The new image will appear lighter after clicking the crop button.
The image will enhance once the cropping is complete.*



2. Frame the desired image by holding the left mouse button in the upper left corner of the image and dragging to the lower right hand corner.



Chapter 5: Acquiring Radiographs

Cropping a Radiograph (continued)

3. The image will enhance once it is cropped.



4. Re-crop the image if necessary before taking the next radiograph by pressing the crop button.




Images must be cropped before shooting the next image.

Chapter 6: Managing Protocols


Using Protocols

Shoot a series of X-rays in a preset order which automatically advances to the next shot after each shot taken.

1. From the Patient Screen, select the **Patient** and click **Protocol**.



Or, from the Patient Screen, click **New** and click **Protocol**.




OR, from the Add Patient or Edit Patient screen, click **Save + Protocol**.
(To add a Patient see page 14).

Chapter 6: Managing Protocols

Using Protocols (*continued*)

2. Select a **Protocol** button to begin the protocol series.



3. The Acquisition Screen now has a Protocol Information bar showing the Protocol shot number. The current Shot Information is located in the upper left portion of the screen.



Chapter 6: Managing Protocols

Using Protocols (*continued*)

4. Shoot the first image. The software will automatically advance to the next shot in the protocol series. The shot number in the Protocol Information Bar keeps track of the shot sequence, while the current Shot Information appears in the upper left.



5. Use the Protocol Information bar to repeat a shot, skip a shot, restart the protocol, move to the last shot in the protocol, or stop the protocol.



Move to the first shot in the protocol series



Repeat the previous shot



Skip to the next shot



Skip to the last shot in the protocol series



Stop the protocol

6. Click the **Review** button  to open the Review screen (see page 35).
7. Click the **Patient** button  to return to the Patient Screen.

Chapter 6: Managing Protocols (continued)


Adding New Protocols

Create new protocols for your individual studies.


1. Click the **Manage** button  in the Patient Screen.



To access the Manage screen, the user needs to have Administrative rights. Changing items in the Manage Screen other than “Users” and “Protocols” may cause your system to function improperly. If you are concerned about making changes in the Manage Screen, please contact Technical Support at 800.268.5354.




2. Click on the **Plus Sign**  next to the Protocols List to show the existing protocols.
3. Select the **Protocols** menu from the left side of the screen to load the Protocol Screen.




Chapter 6: Managing Protocols

Adding New Protocols (*continued*)

4. Enter the name of your new protocol in the **Protocol Name** field.
5. Click the **Save** button.




6. Filter the shots as necessary by selecting the **Filter by Laterality** drop down menu.




Chapter 6: Managing Protocols

Adding New Protocols (*continued*)

7. Select the first shot from the Shots for Species list.







8. Click the Link button **<-Link** to add the shot to the Shots Linked to Protocol box.



9. Repeat step 6 through 8 as necessary.


Chapter 6: Managing Protocols

Adding New Protocols (*continued*)

10. Click the **Unlink** button  to remove any shots from the **Shots Linked to Protocol** box.
11. Select a shot in the **Shots Linked to Protocol** box and change the order using the **Move Up**, **Move Down**, **Move First**, and **Move Last** buttons.
 -  Moves your protocol shot up one shot in the series.
 -  Moves your protocol shot down one shot in the series.
 -  Moves your protocol shot from its location to the first shot in the series.
 -  Moves your protocol shot from its location to the last shot in the series.
12. Click the **Save** button  to save your protocols.
13. Click the **Patient** button  to return to the Patient Screen.



If you have forgotten to save your protocols, a Save dialogue box will appear. Click the Yes button to save your changes and return to the Patient Screen.



Chapter 6: Managing Protocols (continued)


Editing Existing Protocols

Edit current Protocols to customize them to your shot series.


1. Click the **Manage** button in the Patient Screen.



To access the Manage screen, the user needs to have Administrative rights. Changing items in the Manage Screen other than “Users” and “Protocols” may cause your system to function improperly. If you are concerned about making changes in the Manage Screen, please contact Technical Support at 800.268.5354.



2. Click on the Plus Sign  next to the Protocols list to show the existing protocols.
3. Select a Protocol to edit from the left side menu, make changes and click **Save**.



4. Continue with Step 8 on page 29.



Closing the TruDRiX™ software and re-opening is required to save protocol database changes.

Chapter 6: Managing Protocols (continued)


Renaming Protocols

Rename current protocols.


1. Click the **Manage** button in the Patient Screen.



To access the Manage screen, the user needs to have Administrative rights. Changing items in the Manage Screen other than “Users” and “Protocols” may cause your system to function improperly. If you are concerned about making changes in the Manage Screen, please contact Technical Support at 800.268.5354.




2. Click on the **Plus Sign**  next to the Protocols List to show the existing protocols.
3. Select the **Protocol** to rename from the left side menu click the **Rename** button.




Chapter 6: Managing Protocols

Renaming Protocols (*continued*)

4. Rename the protocol in the text box and click the **OK** button when finished.



5. Restart the **TruDRIX™** software by clicking the close button in the upper right hand corner.



Chapter 6: Managing Protocols (continued)


Deleting Protocols

Delete current protocols.


1. Click the Manage button  in the Patient Screen.




To access the Manage screen, the user needs to have Administrative rights. Changing items in the Manage Screen other than “Users” and “Protocols” may cause your system to function improperly. If you are concerned about making changes in the Manage Screen, please contact Technical Support at 800.268.5354.



2. Click on the Plus Sign  next to the Protocols List to show the existing protocols.
3. Select the Protocol to delete from the left side menu click the Delete button.



4. Click the OK button to permanently remove the protocol.




5. Restart the TruDRIx™ software by clicking the close button in the upper right corner.

Chapter 7: Reviewing Images

Reviewing a Study

Change the layout or delete any unwanted images before the study is archived on the Storage Station.

1. Select a patient by clicking on the name in the **Patient List**.
2. Click on the **Study** from the list on the right.
3. Click the **Review** button. 



4. Select a thumbnail in the left-hand pane to load the image into the review pane on the right. The selected image is highlighted by a yellow box.



5. Review the images by using the following techniques:
 - Zoom:** Use the LEFT and RIGHT arrow keys.
Hold the arrow keys down for continuous zoom.
 - Pan:** Click and hold the left touch pad button and move your finger across the touch pad.
 - Window Level:** Click and hold the right touch-pad button then move your finger diagonally across the touch pad.

Chapter 7: Reviewing Images

Reviewing a Study (continued)

6. Use the toolbar for advanced review features.



Patient: Returns to the Patient Screen.



New: Opens the Anatomy/Protocol Screen for creating a new study or protocol.



Acquire: Opens the Views Screen using the anatomy you selected from the original study.



Anatomy: Opens the Anatomy/Protocol Screen to select a new view to add into the same study.



Save: Saves the image with any changes that have altered the original image.



Warning message will appear to choose “Yes” to overwrite (replace) the image or “No” to create a new (additional) image.



Delete: Removes any unwanted images.



Deleting an image permanently removes it from the computer.



Revert: Reverts to the original image.



Left: Rotates the image 90° counterclockwise.



Right: Rotates the image 90° clockwise.



Flip: Flips the image vertically.



Reverse: Flips the image horizontally.



Magnify: Activates the magnifier which will zoom on the selected area.



Actions: Shows a drop-down menu with additional review features.



Win Level: Selects the window level feature (also accessible by using the right mouse button).



Invert: Inverts blacks and whites, replicating a fluoroscopy.



Pan: Pans the image (also accessible by using the left mouse button).



Crop: Trims an image without applying enhancements.



In the Review Screen, cropping does not re-enhance the image based on the captured area. It only trims the image without re-enhancing the image.

Chapter 7: Reviewing Images

Reviewing a Study (*continued*)



Annotations:

Shows a drop-down menu with annotations.



Display:

Displays the annotations.



Hide:

Hides the annotations.



Add:

Adds annotations to the image using the Annotation Toolbar (shown below).



Select:

Selects an annotation button or previous annotation.



Line:

Click and drag to draw a line.



Rectangle:

Click and drag to create a square.



Ellipse:

Click and drag to create a circle.



Text:

Click and drag a box to open for text to be typed.



Ruler:

Click and drag to draw a line measurement.



Poly Ruler:

Create a measurement using multiple lines by clicking at each point.



Protractor:

To measure an angle, click the center of the angle, then click the base line, then click the angle line.



Cross Product:

Click and drag to create two perpendicular lines of equal length.



Layout:

Click a view box, then click the image to load.



1 Up:

1 image displays on the screen.



2 Up Horizontal:

Image windows tile horizontally.



2 Up Vertical:

Image windows tile vertically.



4 Up:

4 images display.



Overlay:

Opens menu for DICOM text display.



Hide:

Hides all DICOM information.



Summary:

Displays patient and shot information.



Detail:

Displays patient, shot, and additional detailed information.


7. Click the **Patient** button  to return to the Patient Screen.

Chapter 8: Exporting Studies and Images

Burning a CD

Burn a copy of a study to a CD.

1. Search for a patient in the Patient Screen by entering the **Patient ID**, the **Patient Name**, or the **Owner Last Name**. The Patient List will automatically filter as you type.
2. Select a patient and a study. The patient and study are highlighted in a blue box.
3. Click the CD/DVD button. 



4. Select either the DICOM or JPEG image format from the Export As: drop-down menu.



5. When selecting the DICOM option from Step 4, Click the Add DICOM viewer check box only if the computer does NOT have a full DICOM viewer already installed.



The options below are unchecked by default.

Annotation Burnin

- Imprints annotations made to saved images

Anonymize Images

- Strips all client and clinic information from the images

Repository Format

- Places images in a repository folder format



Add DICOM Viewer

- Adds a DICOM viewer on the CD to view the images on a computer without a full DICOM viewer installed.

Efilm Lite and Reveal Lite are optional DICOM viewers to add.

Chapter 8: Exporting Studies and Images

Burning a CD (*continued*)

- When selecting the JPEG format from Step 4, select the Overlay Burnin options from the drop-down menu:



The default option is set to None:

- None:** Does not show overlays on the JPEG image
- Detail:** Shows a detailed DICOM overlay on the JPEG image (shows Image Number, Manufacturers Model Number, Study ID, Image Anatomy, View Position, Image Laterality Extended, Columns X Rows, Technique, Institution Name, Physician of Record, Patients Name, Patient ID, Patient Breed Description, Patients Birth Date, Patients Sex, Acquisition Number, Acquisition Date, Acquisition Time.)
- Summary:** Shows a summary DICOM overlay on the JPEG image (shows Image Number, Image Anatomy, View Position, Image Laterality Extended, Institution Name, Physician of Record, Patients Name, Patient Breed Description, Patients Birth Date, Patients Sex, Acquisition Date, Acquisition Time).



- Click the **Select** button on the Export As dialog box.


Chapter 8: Exporting Studies and Images

Burning a CD (*continued*)

8. Insert a blank CD (CD-R or CD-RW) into the computer.
9. Click the **Write Disc** button.



*The CD burn progress is shown in the lower right hand corner.
You will be prompted to remove the disc when the burn is finished.*



10. Click the **OK** button to complete writing the CD and remove the CD from the CD-Rom.



Chapter 8: Exporting Studies and Images (continued)

Exporting Images to a Folder

Export images from a study to a folder to save on a computer, thumb drive, etc.

1. Search for a patient in the Patient Screen by entering the **Patient ID**, the **Patient Name**, or the **Owner Last Name**. The Patient List will automatically filter as you type.
2. Select a patient and a study. The patient and study are highlighted in a blue box.
3. Click the **Folder** button. 



4. Select either the DICOM or JPEG image format from the Export As: drop-down menu.



5. When selecting the DICOM option from Step 4, Click the Add DICOM viewer check box only if the computer does NOT have a full DICOM viewer already installed.



The options below are unchecked by default.

Annotation Burnin - Imprints annotations made to saved images

Anonymize Images - Strips all client and clinic information from the images

Repository Format - Places images in a repository folder format

Create DICOM DIR - Adds a DICOM directory to quickly import all images into a PACS system.

Chapter 8: Exporting Studies and Images

Exporting Images to a Folder (*continued*)

- When selecting the JPEG format from Step 4, select the Overlay Burnin options from the drop-down menu:



The default option is set to None:


- None:** Does not show overlays on the JPEG image
- Detail:** Shows a detailed DICOM overlay on the JPEG image (shows Image Number, Manufacturers Model Number, Study ID, Image Anatomy, View Position, Image Laterality Extended, Columns X Rows, Technique, Institution Name, Physician of Record, Patients Name, Patient ID, Patient Breed Description, Patients Birth Date, Patients Sex, Acquisition Number, Acquisition Date, Acquisition Time.)
- Summary:** Shows a summary DICOM overlay on the JPEG image (shows Image Number, Image Anatomy, View Position, Image Laterality Extended, Institution Name, Physician of Record, Patients Name, Patient Breed Description, Patients Birth Date, Patients Sex, Acquisition Date, Acquisition Time).



Chapter 8: Exporting Studies and Images

Exporting Images to a Folder (*continued*)


7. Click the **Select** button



8. Select a location to save the images.



Images are saved to a selected location in a folder with the animals name and the clients last name (e.g. Trigger-Smith).



Chapter 8: Exporting Studies and Images *(continued)*

Transferring Studies to the Storage Station

Transfer the studies from the Acquisition Station to the Storage/Review Station.




You must have your Tour 1008™ and your Tour 1008™ Review Station networked. Please contact your local IT Person for any networking needs and Sound-Eklin® technical support for configuring DICOM Server for the TruDRlx™ Software.

1. Connect a Cat-5 Internet cable into the Internet port on the computer.



2. Press the F4 key on your keyboard to bring up the Store Files to DICOM Server screen from the Patient Screen.




Chapter 8: Exporting Studies and Images

Transferring Studies to the Storage Station (*continued*)


3. Click the **Start Store** button to begin the DICOM image transfer to the Storage Station.



This will store all the studies that were taken since your last transfer and clear the queue. If any studies are missing, you can transfer them individually.



4. Status box displays transfer results. When the status box states “completed,” your images have successfully transferred to the DICOM server. Click the Red close button in the upper right corner to close the screen.





Chapter 8: Exporting Studies and Images *(continued)*


Transferring Individual Studies to the Storage Station

Transfer the studies from the Acquisition Station to the Storage/Review Station.

1. Search for a patient's study by entering the **Patient ID**, the **Patient Name**, or the **Owner Last Name**. The Patient List will automatically filter as you type.
2. Select a patient by clicking on the patient name in the Patient List field.
3. Click the **Server** button. 



4. Click the **Start Store** button  to begin DICOM image transfer.
5. Status box displays transfer results. When the status box states “completed”, your images have successfully transferred to the DICOM server. Click the Red close button in the upper right corner to close the screen.



Chapter 8: Exporting Studies and Images *(continued)*

Confirm Storage of a Study

Confirm the study successfully sent to the Storage/Review Station.

1. On the Patient Screen, the Archive field will show a server name (e.g. SOUND-SERV1).
2. Hold the mouse over the server name to show a pop-up window with the transfer data.



The pop-up box displays how many images have been transferred and the date and time of the transfer. The information is color coded for easy reference. If only one study appears, it will be highlighted in blue and the color coding will not appear:

Green Box – Image transfer complete

Red Box – Image transfer failed

Description	Locked	Study Date	Images	Archive	Sent
Abdomen	☐	3/9/2009 8:39:00 PM	1	SOUNDSERV1	
Thorax	☐	3/9/2009 8:38:00 PM	1	SOUNDSERV1	

1 of 1 images archived on: 3/9/2009 9:02:00 PM


Chapter 9: Adding Users

Add veterinarians and technicians to TruDRlx™.


1. Click the **Manage** button  in the Acquisition Screen.



To access the Manage screen, the user needs to have Administrative rights. Changing items in the Manage Screen other than “Users” and “Protocols” may cause your system to function improperly. If you are concerned about making changes in the Manage Screen, please contact Technical Support at 800.268.5354.




2. Select the **Users** menu from the left side of the screen to load the Users Screen.



Chapter 9: Adding Users (continued)

- Click the **Add User** button  in the Acquisition Screen.



User ID	First Name	Last Name	Email	Group List
3	AdminFirstName	AdminLastName	Admin@soundtech.com	Admin
1	TechFirstName	TechLastName	Tech@soundtech.com	Tech
4	TechSupportFirstName	TechSupportLastName	TechSupport@soundtec	TechSupport
2	VetFirstName	VetLastName	Vet@soundtech.com	Vet

Selected User Info

User ID: 3

First Name: User Name:

Last Name: Password:

Email:

Groups: ☐ Tech ☐ Vet ☐ Admin ☐ TechSupport

- Enter the User's name in the **First Name** and **Last Name** fields.
- Create a unique user name and password in the **User Name** and **Password** fields.



Suggested user name: First initial of your first name plus last name (e.g. *jsmith* for John Smith). Passwords can be left blank for ease of use.

- Select a Group or Groups for each user by clicking **Tech** or **Vet**.

Each Group allows different permissions.

You may select multiple check boxes per user:


Tech: User is added to the Tech pull down list
has no access to the Manage features.

Vet: User is added to the Vet pull down list,
has no access to the Manage features.

Admin: Has full access to the Manage features.

Tech Support: Should only be used by Sound-Eklin® Technical Support

- Click the **Save** button  to add the User.



User ID	First Name	Last Name	Email	Group List
3	AdminFirstName	AdminLastName	Admin@soundtech.com	Admin
1	TechFirstName	TechLastName	Tech@soundtech.com	Tech
4	TechSupportFirstName	TechSupportLastName	TechSupport@soundtec	TechSupport
2	VetFirstName	VetLastName	Vet@soundtech.com	Vet

- Repeat steps 3 through 7 to add additional Users.

Chapter 10:


Adding and Editing a DICOM Location

Add DICOM server locations for sending images to a Storage Station, Radiologist, etc.


1. Click the **Manage** button  in the Acquisition Screen.



To access the Manage Screen, the user needs to have Admin rights. Changing items in the Manage Screen other than “Users” and “Protocols” may cause your system to function improperly. If you are concerned about making changes in the Manage Screen, please contact Technical Support at 800.268.5354.




2. Select the **Configuration** menu from the left side of the screen to load the Users screen.



Chapter 10: Adding and Editing a DICOM Location *(continued)*

3. Select the **DICOM** button from the left to bring up the DICOM Server Screen.



4. Click the **Add** button to add a new DICOM server. Or if editing select the existing DICOM server from the DICOM Servers list on the Right.
5. Enter the **AE Title**, **IP Address**, and **Port** number for your DICOM server. OR if editing enter the corrected **AE Title**, **IP Address**, and **Port** number for your DICOM server.



Contact Sound-Eklin® Technical Support or the location you are sending the DICOM images for the appropriate DICOM provider information.


Chapter 10: Adding and Editing a DICOM Location *(continued)*


6. Select additional information by choosing Yes from the drop-down menus.



Only select the following options if you are directed to do so by a Technical Support representative:

- Auto Send Server:** Automatically sends images to the selected DICOM location.
- Allow Worklist:** Designates the DICOM location as a worklist modality.
- Telemedicine Server:** Designates the DICOM location as a telemedicine server.
- Archive Server:** Designates the DICOM location as an archive server.

7. Click the **Save** button  to save the current information.
8. Click the **Verify** button  to verify that you can connect to the DICOM Server.
9. A smiley face will appear when you have verified the connection to the DICOM Server.



If a frowning face appears, please contact Sound-Eklin® Technical Support at 800.268.5354.

10. Click the **Patient** button  to return to the Patient Screen.


Chapter 11: Gain Calibration

Calibrate the X-ray Generator and the plate.




Before beginning, make sure the TruDR software is closed and that there is nothing on the digital panel. Remove any X-ray grid that is directly on the panel.

1. Manually set the X-ray machine to the following settings and open the collimator to it's full size:
62 KVP
10 seconds
2. Set the collimator 40 inches (102 cm) from the plate and ensure that the light field covers the entire plate.
3. Open the Viva program from the desktop using **START > Varian > Paxscan > Viva**.
4. Click the **OK** button on the status window when the program opens.



5. Click the **Acquisition** menu from the toolbar menus on top and select **Gain Calibration**.




6. Click the **Continue** button when the pop-up window appears.



The software will automatically begin an offset calibration that will last about one minute.

Chapter 11: Gain Calibration (*continued*)

7. Shoot an X-ray when the screen prompts you to “Initiate X-ray exposure number: 1”.
8. Shoot again when prompted until you have reached a total of 8 shots.
9. Click the Finish button when the status box reads “Finish” OR “Initiate Next Exposure (#9)”.
10. Read the Gain Median in the pop-up window. The number should be between 1000 to 1500.



11. If the Gain Median is not between 1000 and 1500, repeat steps 7-10 using a different technique.




Whether to increase or decrease your technique:

Gain Median is lower than 1000: RAISE the technique

Gain Median is higher than 1500: LOWER the technique



12. Click the Acquisition menu from the toolbar menus on top, and select Acquire Image.



13. When the text box pops up, shoot an X-ray using the same technique you used to determine your gain median.
14. Examine the image to make sure it is a uniform milky-gray color.
15. Close the Viva software and open the TruDRiX™ software to begin shooting.



Dr. Marc Laxineta, DVM



Sound-Eklin™ Technical Support

6359 Paseo Del Lago
Carlsbad . California 92011 USA

toll free: 800.268.5354 option 3
phone: 760.918.9626
fax: 76.918.9620
international: +1.760.918.9626

www.soundeklin.com

© 2012 All Rights Reserved
Sound-Eklin® — a VCA ANTECH company