ADC COMPACT ID SOFTWARE AND ID TABLET

USER MANUAL

TABLE OF CONTENTS

SAFETY INFORMATION	1
Safety information - ID Station console Safety information - ID Station tablet Technical requirements ADC Compact ID Software PC Configuration	1 4 6 6
ADC COMPACT ID SOFTWARE	7
WORKING WITH THE ID STATION	9
Safety precautions	9
General description of the ADC Compact ID Station	10
Functions and features of the ID Station	10
Components of the ID Station	11
Switching the ID Station on and off	12
Starting the ADC Compact ID Software	12
Exiting the ADC Compact ID programme	14

2201A GB 199805

Working with the menu system	15
Selecting items from a menu	15
Filling in information in an entry screen	17
Online help	17
Function keys	18
Standby Screen	20
File Menu	21
Configuration Menu	22
Operator options in the Configuration menu	23
System Manager options in the Configuration menu	27
Change User option in the Configuration menu	30
Help Menu	30
Inserting and removing a cassette	33
Identifying a cassette	34
Clearing information	54
Recalling a patient	55
New Patient	56
Identifying a cassette with data from the history list	57
Link	60
Hold/Unhold	64
Destination Rerouting	65
Reading a cassette	66
Erasing an image plate	71

Safety information

Safety information - ID Station console

General Warnings

- There are no user serviceable parts on the console.
- Pay attention when working near the pedestal: the foot projects slightly and could be tripped on.
- When rolling the pedestal away from its position, take care not to pull or put any strain on the power lead.
- The pedestal shelves should not be used as a footrest!
- If the shelves are re-positioned at a different height, the locking clips should be re-inserted.
- If the user prefers to install a monitor other than the one provided/sold by Agfa, he should ensure that it is locked to the top shelf of the pedestal at all times, according to the method described in the manual. See drawing supplied with the console.
- When moving the pedestal, take care to ensure that it remains stable at all times. In particular, the following considerations should be kept in mind: the nature of the floor and any floor-coverings, carpets, etc.; any obstructions such as cables; and any slopes to be negotiated.
- Only attach appliances with line voltage equal to that of the ID Station.

Installation of the equipment

The following instructions should be observed for safe installation of the equipment.

Install the tie-wraps that are delivered with the pedestal crosswise over the foot of the monitor to ensure that it is securely fixed.

- Tighten the tie-wraps firmly.
- The tie-wraps may be released by the use of a simple screwdriver.

Monitor requirements

- It is not necessary for the user to buy the monitor offered by Agfa. He may use/buy his own monitor. If he does, it should comply with certain criteria to ensure that it can safely be used in combination with the pedestal with regard to stability and securing method. (Agfa monitors are fully compliant.)
- The maximum size of monitor is 17 inches.
- The foot of the monitor should be such that it may be properly secured to the pedestal with the tie-wraps provided. It is the user's responsibility to ensure this.

Compliance

The equipment (both tablet and pedestal) is approved to the following standards:

- UL1950, Third Edition
- CAN/CSA 22.2 No. 950-95, Third Edition (cUL)
- EN60950 (TÜV)
- GS (TÜV)

The equipment complies with all applicable European Directives and bears the CE mark.

Cleaning

Pull the power plug out of the equipment prior to cleaning the equipment.

2201A GB 199805

Disposal

When its service life is over, do not dispose of the ID station without consideration of local waste disposal regulations.

Please consult your local sales organisation.

Safety information - ID Station tablet

Tablet technical data

- Working voltage (input): 100-230 V ac
- Frequency: 50-60 Hz
- Max. rated current: 5 A (with appliance outlets loaded to the maximum)
- Classification: Class I (earthed)
- Operating conditions: temp.: 15 to 30 °C, RH: 15% to 75%

Fuses:	Europe:	5 Amp Slow/Träge/Temporisé/250 V
	North America:	5 Amp Slow/Träge/Temporisé/125 V min

The equipment is suited for continuous operation.

Power outlets

The equipment is provided with 2 power outlets (1 for a PC and 1 for a monitor).

- Output voltage is equal to input voltage
- Maximum current (outlet 1 + outlet 2): 4 A total
- Frequency: 50-60 Hz

When connecting equipment to these outlets, take care not to exceed the maximum total current of 4A.

2201A GB 199805

Replacement of fuses

Before replacing the fuses, remove the power plug from the lead.

The following important reminder is marked on a label on the equipment:

'For continued protection against a risk of fire, replace fuse with same type and rating.'

The fuse type is

- for Europe: 5 Amp Slow/250 V
- for North America: 5 Amp Slow/ 125 V min

Compliance

The equipment is approved to the following standards:

- UL1950, Third Edition
- CAN/CSA 22.2 No. 950-95, Third Edition (cUL)
- EN60950 (TÜV)
- GS (TÜV)

The equipment complies with all applicable European Directives and bears the CE mark.

Cleaning

Pull the power plug out of the equipment prior to cleaning the equipment.

Disposal

When its service life is over, do not dispose of the ID station without consideration of local waste disposal regulations.

Please consult your local sales organisation.

Technical requirements ADC Compact ID Software

PC Configuration

- Pentium 200 Mhz MMX
- 64 MB RAM (total)
- 2 MB video RAM
- 1 CD ROM
- 1.6 GB hard disk
- Multimedia
- 10/100 Base T ethernet
- 17 inch monitor
- Windows ® 95, Windows NT ®, Unix
- ID Tablet

ADC Compact ID Software

The ADC Compact ID Software is used to record patient demographic and examination data onto the memory chip that is embedded in the ADC Compact cassette. Patient demographic and examination data are written onto the ADC Compact cassette memory by contactless radiofrequency transfer when the cassette is inserted into the transmission slot of the ID tablet.

During the cassette readout, the ADC Compact system links the demographic data with the image and uses the examination data to set the appropriate image processing parameters.

The ADC Compact ID Software runs on standard PC platforms. If PCs conforming to the required standards are already available in the department, you only need to install the ADC Compact ID Software and the ADC Compact ID tablet to have a fully operational ADC Compact ID Station. For a full description of the required hardware performance standards we refer to page 6 of this manual. The ADC Compact ID software can also be installed on Unix workstation platforms. As a result, in smaller departments the ID software can run on the ADC Compact Processing server without the need for an additional PC.



Agfa also offers a fully configured ADC Compact ID station including the PC, the ID Software, the ID tablet and the console with storage racks, ready for use.



ADC Compact ID Software and ID Tablet

Working with the ID Station

Safety precautions

The ID Station complies with the UL544 standard for 'Non Patient Care Medical Equipment'. This means that, although it is absolutely safe, patients may not come into direct contact with the equipment.

Therefore the ID Station must be placed outside a radius of 1.5 m around the patient.

There are no user-serviceable parts inside the ID Station. Perform no other operations on the ID Station than the ones described in this manual.

General description of the ADC Compact ID Station

Functions and features of the ID Station

With the ID Station you can record patient demographics and examination data required to obtain an optimally processed image onto the memory chip which is embedded in the ADC Cassette.

You have only to insert the cassette into the transmission slot of the ADC Compact ID tablet. You can then perform the various operations using the interactive menus of the ADC Compact ID Software.

Besides identifying a cassette, you can also perform the following operations with the ID Station:

- Read information from a previously identified cassette.
- Command the erasure of the image plate inside the ADC Cassette in the next Digitizer cycle.
- Create and use a patient list. You can transfer data from the patient list directly onto the cassette and from the identification screen to the patient list, eliminating the need to retype any information.
- Using the ADC Rislink Toolkit, the ADC Compact ID Software can be integrated into your RIS, permitting rapid and error-free communication of patient data from the RIS system.

As there is no physical link between the ID Station and the Digitizer, you can easily employ multiple ID Stations to serve one Digitizer.

2201A GB 199805

Components of the ID Station

The illustration below shows the components of a fully configured ID Station.



- PC with internal disk drive, floppy disk drive and CD-ROM drive
- monitor and keyboard to operate the menu system
- ID tablet to perform operations on ADC Cassettes
- cassette storage rack
- console with monitor support.

Switching the ID Station on and off

The ADC Compact ID Software is an extremely user-friendly Windows-based application. The assumption is that the user is thoroughly familiar with the look and feel of the graphical user interface (GUI) used in the industry-standard Microsoft Windows operating systems and in Windows-based applications, as this graphical user interface is the most intuitive way for people to interact with a computer.

A principal goal of the GUI in Microsoft Windows is to establish visually and functionally consistent communication between the user and the application, i.e. consistent keyboard commands and consistent menus for identical actions. Such consistency helps ensure that a user who can use one Windows-based application can use them all. We therefore assume that the user of the ADC Compact ID Software is fully conversant with windows and window arrangements, sizing, scrolling and splitting controls, menus and dialog boxes, toolbars and rulers, shortcut and access keys, GUI messages and status indicators, and the like.

- Note The default Graphic User Interface (GUI) can be re-configured. An Agfa Service technician or the hospital's System Manager can configure the following GUI elements on request:
 - Whether a specific field is present or not
 - Specifying the location of fields
 - Specifying the font size

Starting the ADC Compact ID Software

To start the ADC Compact ID Software from the Windows main desktop screen, you can proceed in the following ways. You either

- Start the programme by clicking the Windows Start button and then click the ID Software programme in the menu
- Have the ADC Compact ID Software started each time Windows starts
- Start the ADC Compact ID Software by using the Windows Run command
- Click the ADC Compact ID Software icon on the desktop

2201A GB 199805

The ADC Compact ID Software first runs an internal test procedure, then displays the ID Startup screen, an example of which is shown below. After prompting the user to enter the user name and password, the ADC Compact ID programme is ready for use.





The ADC Compact ID Software is supplied virus-free on CD. However, it is always advisable to install a virus protection programme.

Exiting the ADC Compact ID programme

In order to exit the ADC Compact ID application from the Standby screen or from the Identification screen, proceed as follows:

Click the File menu and then click Exit.

ADC2 IDENTIFICATION STATION (IDEN> IDENTIFICATION SCREEN							
File Darigunitari Mode Per	form Help						
Read concette							
Evree cassette	Inknown	Firstname	Unknown				
Show history	00/00/0000	Patient ID	Unknown				
Ext	Nale	Accession number					
Uger Info							
STIDY							
Radiologist	ASPA DEPAID/7	-					
Examination	SYSTEM DIAGNOSIS	 Sub examination Ld 	s/Aofa CB contrast				
INAGE			_				
Department	AGFA						
Patient position	PA 💌						
Cassette orientation	Portrait ·	Esposure Cla	aa 200 💌				
User Info 1							
DESTINATIONS							
Mardcopy unit	None	Number of copies	0 💌				
Processing station	ADC_PS1	Pilm format	BINXLOIN -				
Archive station	None						
Send station	None						
F1 Write F2 C1	sar F3 Recall F4	New patient F? His	tory F10 Cancel				
Ready for identification	1. 190	TO (HOLD)	11/03/1997 08:52:23				

Or click the Close button X at the far right of the title bar.

The system then returns the following screen window asking for confirmation.

₿E	xit					n.			x
Are	you	sure	you	want	to	leave	this	program?	
				[Ye	8	No			

14 2201A GB 199805

Working with the menu system

Selecting items from a menu

The ADC Compact ID application has a menu-driven, or interactive, user interface. It consists of drop-down menu screens, cascading menu screens and entry screens. You can open drop-down menus by selecting and clicking the menu name. The cascading menus or submenus are accessed by selecting the menu commands that have a right-pointing triangle, as is shown in the Configuration menu below.



The entry fields or text boxes in the entry screens are accessed by using the [TAB] key (or Shift [TAB] key for moving back) until the item of your choice appears highlighted on the menu screen. Then [ENTER] is pressed to activate the item. Entry screens sometimes feature drop-down list boxes, showing a list of items the user can select. Clicking the drop-down arrow opens the list box to show the available options. The option is selected by moving the pointer and clicking the relevant item.

For example, if you have to fill in the sex in the patient field of the main identification screen, you simply click the drop-down arrow, which will open the available options, drag the pointer to your relevant choice and release the mouse button. An example of the 'Sex' drop-down list box in the main identification screen is shown below.

BADC2 IDENTIFICATION STATION (IDEN) IDENTIFICATION SCREEN							
File Scrigenton Mode Per	form Help						
PATIENT							
Patient Name		Firstname					
Birth date	00/00/0000	Patient ID					
8ех	Nale	Accession number					
Uwer Info	Nale						
STUDY	Penale						
Radiologist	AGFA DEPAULT						
Exemination	SYSTEM DIAGNOSIS	Sub examination Lds/Agfa CR contrast 💌					
INAGE							
Department	AOFA						
Patient position	PA 🔹						
Cassette orientation	Portrait •	Exposure Class 200 💌					
Uger Info 1							
DESTINATIONS							
Mardcopy unit	⊅13000 ×	Mumber of copies 0 💌					
Processing station	ADC_PS1	Film format 14INX17IN .					
Archive station	None						
Send station	None						
F1 Write F2 C1	sar F3 Recall F4 R	ew patient F? History F10 Cancel					
Ready for identification	1. INVERT	HOLD 17/03/1997 24(42)10					

Alternatively you can also enter the first letter, or digit, of the item of your choice to highlight it, then press [ENTER]. So, in order to select 'Female' in the Sex entry field in the above identification screen, you can simply type the letter 'F'.

In this manual the term 'Select' is used to designate the combined activity of highlighting an item and pressing [ENTER] to activate it.

Note The number of items in a list may depend on the specific configuration of your ADC System.



Filling in information in an entry screen

The ADC Compact ID application provides the following text editing features for entering or editing text within an entry field:

- Press the [TAB] respectively Shift [TAB] keys to move the cursor to within the highlighted field.
- Press [DEL] to remove the character over which the cursor is positioned.
- Typed characters either overwrite those over which the cursor is positioned or are inserted, depending on whether the [INSERT] key has been pressed. Press [INSERT] to toggle between insert and overwrite mode.
- You can press [ENTER] at any time while editing the highlighted fields, to accept what is in the field and move to the next one for editing.
- Empty fields are filled automatically with default values.

Online help

In a number of cases there is a hint message in the bottom of the screen, telling you what to do to fill in the current screen correctly.

Function keys

In the 'Identification' entry screen, you can use and configure function keys to shortcut the identification procedure. Up to 10 Function keys (F1 to F10) can be configured using the configuration tool.

ADC2 IDENTIFICATION STATION <iden> IDENTIFICATION SCREEN</iden>							
File Configuration Mode Perform Help							
PATIENT							
Patient Name		Firstname					
Birth date	00/00/0000	Patient ID					
8ех	Male 🔹	Accession number					
Uwer Info							
STUDY							
Radiologist	ASPA DEPAULT	-					
Exemination	SYSTEM DIAGNOSIS	 Sub examination Lds/Agfa CR contrast . 					
INAGE							
Department	AOFA						
Patient position	PA 💌						
Cassette orientation	Portrait •	Exposure Class 200 -					
User Info 1							
DESTINATIONS							
Mardcopy unit	Bone	Number of copies 0 .					
Processing station	ADC_PS1	Film format BINX10IN					
Archive station	None						
Send station	None						
Fi Write F2 Cl	sar F3 Recall F4	New patient #7 History #10 Cancel					
Ready for identification	1. Jan t	0[H0LD] 11/03/1997 05:35:01					

The table below is a sample list of a number of available functions and the function keys they could have been assigned to, together with their description (the name of the key) and what happens when you press the relevant key. These functions are displayed in the bottommost toolbar of the identification entry screen.

Кеу	Description	Function			
F1	Write	Enters the information from the current fill-in screen into the cassette at the moment it is pressed.			
F2	Clear	Clears the highlighted field.			
F3	Recall	Recalls 15 previous data of the Patient group, User Info fields and Comment field. When recalling an item, the previously identified information of the last identification is called back and is used for the subsequent identification.			
F4	New Patient	Clears all the patient data fields in the Patient group.			
F5	Add Patient	Adds all the patient information to the patient list and saves the list.			
F6	Link	Initiates a link between two cassettes.			
F7	History	Copies patient information from a list with the last identified patients directly to the identification screen after selection.			
F8	Re-route	Re-routs images that have been autorouted by default to certain hardcopy, review, archive or processing destinations, to another destination.			
F9	RIS	Provides access to the RIS system.			
F10	Cancel	Brings you back to the previous menu or entry screen and eventually to the main Standby screen.			

Also other functions, which will be described later in this manual, like 'Internal worklist', 'Hold', 'Full Leg/Full Spine', 'Build Study', and so on can be customized into a function key.

Standby Screen

After having started the ADC Compact ID application from the Windows main desktop screen, the system returns a Startup screen prompting you to enter your user name and password. User names and passwords differ, depending on whether you log in as a User, a System manager or Service. The User and System passwords are configurable.

After entering your user name and password, the system displays the following screen, featuring three active drop-down menus, 'File', 'Configuration' and 'Help'. The 'Mode' and 'Perform' menus cannot be accessed from the Standby screen. They are only active in the Identification screen proper.



File Menu

The functions in the File Menu of the Standby screen, viz. Read cassette, Erase cassette, Patient List and Exit as displayed in the screen below, can be accessed from the identification screen proper and will be discussed in a further section of this manual.





Configuration Menu

When moving the pointer to 'Configuration' in the Standby screen, as a User, the drop-down menu unfolds two active options : 'Operator' and 'Change user', as is shown in the sample screen dump below.



The System Manager and Service options within this menu are inactive and can only be accessed by the System Manager, respectively Agfa service personnel.

Operator options in the Configuration menu

By selecting the 'Operator' option within the Configuration menu, a cascading menu, as is shown below, can be accessed featuring various options:



- Hold status
- Sound
- Automatic text selection
- Reset window size
- Always to front
- Destination rerouting The Destination rerouting option can only be accessed when the ADC Compact Autorouting Software has been installed.

Hold status

The 'Hold' status can be toggled ON or OFF. When it is set ON, the system displays the 'Hold' indication on the status line of the main identification screen, as shown in the sample screen dump below.

ADC2 IDENTIFICATION STATION (IDEN) IDENTIFICATION SCREEN						
File Serigentist: Mode Per	form Help					
PATIENT						
Patient Name	Unknown	Pirstname	Unknown			
Birth date	00/00/0000	Patient ID	Unknown			
8ех	Male 💌	Accession number	r			
Uwer Info						
STUDY						
Radiologist	ASPA DEPAULT	-				
Exemination	SYSTEM DIAGNOSIS	 Sub examination 	ids/Agfa CR contrast 🔳			
INAGE						
Department	AGPA					
Patient position	PA 💌					
Cassette orientation	Fortrait	Esposure C.	Lang 200 💌			
User Info 1						
DESTINATIONS						
Mardcopy unit	None 💌	Number of copie	ur 🛛 💌			
Processing station	ADC_PS1	Film format	BINX10IN -			
Archive station	None 💌					
Send station	None 💌					
Fi Write F2 Cl	ear F3 Recall F4	New patient F? H	istory F10 Cancel			
Ready for identification	i. 19 7	TO (HOLD)	11/03/1997 08:50:21			

When destinations are on 'Hold', the default routing of images to hardcopy units, archive and review stations is halted. In order to override the 'Hold' function and re-initiate the sending procedure, you must press the [PROCEED] button in the Preview screen, provided the ADC Compact Preview Software has been installed on your system. If you fail to press the [PROCEED] button in the Preview screen, the system will automatically, after a set time-out, re-initiate the sending procedure. When pressing [CANCEL] the default routing of images is discontinued. It can only be re-activated from the Processing Station.

The sending of images to their default destinations can be put 'on hold':

- by interactively toggling the Hold/Unhold option in the Perform menu of the Main identification screen or
- by pressing a function key that has previously been configured as a Hold/Unhold button
- the service technician might also assign the 'Hold' attribute to certain examinations by default, for example bed-side thorax, by altering the cpf file. This 'Hold' attribute can be overridden using the Hold/unhold toggle in the Perform menu of the identification screen or by pressing the Hold function button, if any.

Sound

The Sound item in the Operator menu enables respectively disables the beep sound option. In the dialog box as displayed below you can select the sound device, i.e. sound card respectively motherboard.

🖉 Image hiding	×
Active	🖲 On 🔿 Off
Sound devic	e Audio Speaker
Save	e CanAudio

Automatic text selection

When the Automatic text selection option is on, the data of the entry field selected is highlighted. When new characters are keyboarded, the highlighted text is automatically deleted and replaced with the newly entered characters. This feature only affects editable fields.

Reset window size

Selecting the Reset window size option sets the active window size to default.

Always to front



The Always to front option is only relevant in case the ADC Compact Preview Software has been installed on your system. It is a toggle option which enables the user to set the Identification Standby screen always in front of the Preview screen.

With the Always to front option ON, the identification screen remains in front of the Preview screen, or comes to the front whenever an identification procedure is performed. When the Always to front option is OFF, the Identification screen is put behind the Preview screen and remains there until any identification action is performed.

X
CR contrast
2
10 .
1217EN .

System Manager options in the Configuration menu

Having logged in as System Manager, the System Manager option within the Configuration menu, as shown below, can be accessed featuring the following options:

- Set time zone
- Change passwords
- Screensaver
- Make fields RO
- Set department name
- Change language

Cassette detection



Set time zone

Using the Set time zone option the System Manager can set the system's time zone.



Change passwords

The Change passwords option enables the System Manager to set User and System passwords. Whenever a new password is entered, the old password is asked for.

ě					×
Please enter	username,	old	password,	new	password
Username					
Old password			4		
New password					
Verification					
		Canc	el		

2201A GB 199805

ADC Compact ID Software and ID Tablet

The use of 'User Password' can be enabled, respectively disabled by the service technician. If disabled, the user is not asked for a password when accessing the Identification application.

Screensaver

Note

The Screensaver option is used for enabling respectively disabling the use of a screen saver. In case it is enabled, the time-out can be configured to 1, 5, 15, 30 and 60 minutes after any identification activity.

Make fields RO

This option makes fields 'Read-Only'.

Set department name

This option enables the System Manager to set the department name.

🦉 Set department name			>
Department		k	
	Save	Clear contents	Cancel

Change language

The Change language option enables the System Manager to select from a scroll-down list a country language other than the default.

👹 Change language			_ 🗆 ×
Language	basic:	Standard data and time format	
	basic:	Standard data and time format	
	da_DK:	Danish for Denmark	
	de_AT:	German for Austria	
	de_BE:	German for Belgium	
	de_CH:	German for Switzerland	
	de_DE:	German for Germany	
	de_IT:	German for Italy	
	de_LU:	German for Luxemburg	
	en_AU:	English for Australia	
	en_CA:	English for Canada	_
	en_GB:	English for United Kingdom	
	en_US:	English for USA	
	es_ES:	Spanish for Spain	

Cassette detection

This option enables the System Manager to enable respectively disable the cassette detection feature.

Change User option in the Configuration menu

器ADC2 IDENTIFICATI	ON STATION «IDEN> IDENTIFICATION SCREEN	. 🗆 🗆
File Configuration (1993):	Derivers Help	
Operator 🕨		
System Manager 🕨		
Benéze 🕨		
Change user		
		·

This option enables the user to change the User level. The password associated with the newly requested user level is asked for. Entered characters show up as asterisks (*).

8	×
Please ente	er username/password.
Username	system
Password	*****
	Cancel

Help Menu

When clicking the Help Menu, the drop-down list displays the following options.



2201A GB 199805 ·

Info

This file contains a full description of the configuration of the ID Station, listing the following items:

Info on Configuration

ID Station Name	
Locale	
Date format	
Link option	ON/OFF
RIS-Link option	ON/OFF
Auto-Routing option	ON/OFF
Full leg/Spine option	ON/OFF
Screen saver via password	ON/OFF

Info on Hardware

Processor type

Info on System

Operating version system

Info on Network

- PC hostname ip-address subnetmask default router AE-title CPF host URL Port number
- Info on RIS
 - RIS enabled Yes/No RIS .exe path RIS data file path RIS locale Remove risfile ON/OFF

ADC Compact ID Software and ID Tablet

Version

When selecting this menu item the ID Station's name and version number is shown.

ADC2 IDENTI	► FICATION STATION
ID-station name	IDEN
Version	PRID1.1.01.03
Cancel	
Inserting and removing a cassette

When you want to identify a cassette, or when the ADC Compact ID application prompts for insertion of a cassette, insert the cassette in the slot of the ADC Compact ID tablet so that the ID Station can access the RF-tag memory chip embedded in the cassette.

To do so, proceed as follows:

1 Orient the cassette with the labels facing up (tube side down) as shown in the illustration below.



- **2** Gently insert the cassette.
- **3** Leave the cassette in this position to perform the operation of your choice.

Caution If you decide not to perform any operation on the cassette, press [CANCEL], then withdraw the cassette from the ID tablet slot.

- 2201A GB 199805 33

Identifying a cassette

Identifying a cassette consists in keying in the required information like patient demographics, examination data, image features, destinations, etc.

Proceed as follows to identify a cassette:

1 Insert a cassette in the cassette slot of the ADC Compact ID tablet.

The ID Station returns the 'Identification' entry screen:

In case of a first identification after start-up, the upper patient field in the screen is blank and in the 'Patient Name' field the pointer will be blinking, as in the sample screen below:

ADC2 IDENTIFICATION 9	STATION «IDEN> IDENTIFICATION SCREEN
PATTERT	Termin Trade
Patient Name Birth date Sex Vger Info	Firstname 00/02/0000 Patient ID Nale Accession number
STUDY	
Radiologist Examination	AGPA DEFAULT . STSTEM DIAGNOSIS . Sub exemination Lds/Agfa CR contrast .
INAGE	
Department	AOPA
Patient position	PA 💌
Cassette orientation	Portrait Exposure Class 200
Uger Info 1	
DESTINATIONS	
Mardropy unit	None Number of copies 0 .
Processing station	ADC_PS1 Film format SINX10IN
Archive station	None 💌
Send station	None 💌
Fi Write F2 Cl	and F3 Recall F4 New patient F7 History F10 Cancel
Ready for identificatio	ta. (##70) 19/03/1997 16:33:40

Note

The format and layout of the Identification screen in your application might differ from the one displayed above, as they are fully configurable. Also the number of displayed fields can be configured.



2201A GB 199805

The default entries in the lower part of the screens are configured by your application engineer during installation of the system. They are linked to the Radiologist, Examination type and Examination Subtype. You can manually override these settings, e.g. to process a specific examination which was not performed in the usual way. The default settings themselves can only be modified by a trained application engineer.

If it is not the first identification after start-up, the patient name field will be highlighted. All fields in the Identification screen will contain the information stored on the previous cassette's memory chip, as you can see from the sample screen below.

器ADC2 IDENTIFICATION S	TATION (IDEN) IDENTIFI	CATION SCREEN	
File Contiguition: Mode Peri	'orm Help		
PATIENT			
Patient Name	COORER	Firstname	ALAN
Birth date	01/04/1947	Patient ID	GUH58 6
8ех	Male •	Accession number	SP12
User Info			
STUDY			
Radiologist	A/SPA DEPAUL/T	•	
Examination	SYSTEM DIAGNOSIS	 Sub examination Ld 	s/Agfa CR spatial 🔳
INGAGE			
Department	Quology		
Patient position	PA 💌		
Cassette orientation	Portrait •	Exposure Cla	aa 50 💌
Uwer Info 1			
DESTINATIONS			
Mardcopy unit	DRYSTAR	Number of copies	1 💌
Processing station	ADC_PS3	Film format	BINX10IN 💌
Archive station	None		
Send station	None *		
Fi Write F2 Cl	ear F3 Recall F4	New patient F7 His	tory F10 Cancel
Ready for identification	La .	[MITO]	19/03/1997 17:03:03

- 2201A GB 199805 35

This feature saves time as you do not have to re-enter information when several cassettes of the same patient have to be identified. If you want to enter information for a different patient, press the New Patient function key, if any, or select the New Patient option in the Perform menu to remove the identification data and to move the cursor to the first data field.

器ADC2 IDENTIFICATION S	TATION (IDEN) IDENTIF	ICATION SCREEN	
File Servicention: Mode Per	form Help		
PATIENT			
Patient Name	COOPER.	Firstname	ALAN
Birth date	01/04/1947	Patient ID	GUH28 6
8ех	Nale 💌	Accession number	
Uwer Info			
STUDY			
Radiologist	ASPA DEPAULT	•	
Examination	SYSTEM DIAGNOSIS	 Sub examination Ld 	s/Agfa CR spatial 🔳
INAGE			
Department	Ωκοτοδλ		
Patient position	PA 🔹		
Cassette orientation	Portrait •	Exposure Cla	aa 50 💌
User Info 1			
DESTINATIONS			
Mardcopy unit	DRYSTAR 💌	Number of copies	1 💌
Processing station	ADC_PS3	Film format	BINXIOIN 💌
Archive station	None		
Send station	None *		
Fi Write F2 Cl	sar F3 Recall P4	Hew patient 17 His	tory F10 Cancel
Ready for identification	1	INUTO	19/03/1997 17:10:41

器ADC2 IDENTIFICATIO	IN STATION CIDE	EN> IDENTIFIC	ATION SCREEN		- O ×
File Scritganizar: Mode	Perrora Help				
PATIENT	Wile				
Patient Name	Clear		Firstname	ALAI	9
Birth date	Recall 94	17	Patient 1	D GUB	58.6
Sex	New proent		Iconviou	number	
	Hold Minhold		200 C 10 C		
User Into	Cancel	_			
STUDY	Link				
Radiologist	Reputing 11	ul T	•		
Examination	SISTER PL	USNOSIS	 Sub examination 	ation Lds/Agi	a CR spatial 💌
INAGE					
Department	Grology				
Patient position	PA 💌				
Cassette orientati	on Portrait	-	Expos	sure Class 50	•
User Info 1					
DESTINATIONS					
Mardcopy unit	DRYSTAR		Number of	f copies 1 💌	
Processing station	ADC_PS3	*	Film for	ant BINS	10IN -
Archive station	None	*			
Send station	None	*			
Fi Write F2	Claur F3 3	Secall F4	New patient	F? History	F10 Cancel
Ready for identifica	tion.	1	autoj	19.	/03/1997 17:09:23

2 Fill in the patient name (33 characters maximum), then press [ENTER] or [TAB].

Whenever [ENTER] or [TAB] is pressed, the cursor moves to the next field. If you want to move back to the previous field, press Shift [TAB].

- Note You do not have to fill in all the fields in this data entry screen. Empty fields are automatically filled in with default values.
 - **3** Fill in the patient's first name (20 characters maximum), then press [ENTER] or [TAB].

The cursor moves to the next field.

4 Fill in the date of birth of the patient, then press [ENTER] or [TAB].

5 Fill in the patient ID code (15 characters maximum), then press [ENTER] or [TAB].

The cursor moves to the next field.

6 Select the patient's sex in the drop-down list box, then press [ENTER] or [TAB].

The cursor moves to the next field.

7 Enter the patient's accession number, press [ENTER] or [TAB].

The patient's accession number is the RIS identification code that is used for extracting the current examination information from the Radiology Information System (RIS) or in the Hospital Information System (HIS). If you want more information on how to connect your ADC Compact system with your RIS or HIS, refer to your local Agfa representative or consult the ADC Compact Rislink Toolkit manual.

器ADC2 IDENTIFICATION S	TATION (DEN) IDENTIF	ICATION SCREEN	
File Contrantor: Mode Per	form Help		
PATIENT			
Patient Name	COOPER.	Firstness	ALAN
Birth date	01/04/1947	Patient ID	GUB58 6
8ex	Male 💌	Accession number	RIS-6603
User Info			2
STUDY			
Radiologist	ASPA DEPAULT	•	
Examination	SYSTEM DIASNOSIS	 Sub examination Ld 	s/Agfa CR spatial 💽
INAGE			
Department	Ωιοτοδλ		
Patient position	PA ·		
Cassette orientation	Portrait 🔹	Exposure Cla	aa 50 💌
Uger Info 1			
DESTINATIONS			
Mardcopy unit	None	Number of copies	1 💌
Processing station	ADC_PS1 ·	Film format	14INX17IN 💌
Archive station	None		
Send station	None		
Fi Write F2 Cl	enr F3 Recall F4	New patient F7 His	tory F10 Cancel
Ready for identification	1.	[auto]	19/03/1997 17:45:11

After having pressed [ENTER] or [TAB], the cursor moves to the next field in the sample screen as shown below, the User Info field. User Info fields are optional information fields which you may not want to include or which you can configure, assign a different name and place anywhere in the 'Patient', 'Study' or 'Image' sections of the Identification screen.

器ADC2 IDENTIFICATION S	TATION (IDEN) IDENTIF	ICATION SCREEN	_ O ×
File Danigunities: Mode Peri	form Help		
PATIENT			
Patient Name	COOPER.	Firstname	ALAN
Birth date	01/04/1947	Patient ID	GUB58 6
8ех	Male 💌	Accession number	RIS-6603
Uwer Info	Bed N° 776		
STUDY	4		
Radiologist	ASPA DEPAULT	-	
Examination	SYSTEM DIAGNOSIS	 Sub examination Ld 	s/Agfa CR spatial 💌
INAGE			
Department	Οιοτοδλ		
Patient position	PA ·		
Cassette orientation	Portrait 🔹	Exposure Cla	aa 50 💌
Uwer Info 1			
DESTINATIONS			
Mardcopy unit	Bone ·	Number of copies	1 💌
Processing station	ADC_PS1 ·	Film format	14INX17IN 💌
Archive station	None 💌		
Send station	None *		
Fi Write F2 Cl	ear F3 Recall F4	New patient F7 His	story F10 Cancel
Ready for identification	L	190201	19/03/1997 18:21:55

- 2201A GB 199805 **39**

When selecting New Patient, the Patient Name field turns blank, together with the User Info field (in point of case, the bed number in the hospital).

器ADC2 IDENTIFICATION S	STATION «IDEN» IDENTIFICATION SCREEN	
File Scritzentist: Mode Per	form Help	
PATIENT		
Patient Name	Firstness	
Birth date	00/00/0000 Patient ID	1
Зех	Other Accession number	
Uwer Info		
STUDY	h.	
Radiologist	ASPA DEPAULT	
Examination	SYSTEM DIASNOSIS 🔹 Sub examination Lds/Agfa CR spatial	
INAGE		
Department	Urology	
Patient position	PA .	
Cassette orientation	Portrait Exposure Class 50	
Uwer Info 1		
DESTINATIONS		
Mardcopy unit	None Mumber of copies 1 V	
Processing station	ADC_PS1 Film format 14INX17IN .	
Archive station	Bone	
Send station	Bone	
Fi Write F2 Cl	and F3 Recall F4 New patient F7 History F10 Cancel	
Ready for identification	n. (3070) 19/03/1997 18:24:5	52

This concludes the completion of the patient data. After having pressed [ENTER] or [TAB], the cursor moves to the 'Radiologist' field in the 'Study' section of the Identification screen.

Note Filling in the 'Date of birth' field is not mandatory. Empty fields are automatically filled in with 'unknown'.

8 Press [ENTER] or [TAB] in the Radiologist field if you want to select the currently displayed radiologist.

Otherwise, press the first character of the radiologist's name, or press the downward arrow in the list box to display the list of radiologists. A drop-down menu appears, showing a number of radiologists. Select the radiologist in charge of the exposure.

By selecting a specific radiologist you automatically select default parameters adapted to, or preferred by, this radiologist.

After pressing [ENTER] or [TAB], the cursor moves to the Examination field.

9 In the Examination field, press the first character of the examination type, or select the examination type from the drop-down list box that is displayed when clicking the downward arrow. The system returns a screen similar to the sample screen below.

BADC2 IDENTIFICATION 9	STATION «IDEN» IDEN	TIFICATION	SCREEN				- 0
File Conformation Mode Per	form Help						
PATIENT							
Patient Name		3	firstname				
Birth date	00/00/0000	I	Patient I	D			
Sex	Other	• ,	Accession	number			
User Info							
STUDY							
Radiologist	AGFA DEFAULT	*		_			
Examination	SERVICE	K 50	b examina	tion CE	EST PO	RTABL	E
INAGE	Chest	`					
Department	upper extremit.						
Patient position	shoulder	_					
Cassette orientation	spine skull		Expos	nure Cla	aa 200		
User Info 1	pelvis						
DESTINATIONS	aboomen	_					
Herdcopy unit	None	-	Number of	copies	1 •		
Processing station	ADC_P81	-	Film form	nat	8 I M X 1	018	*
Archive station	None	•					
Send station	None	•					
Pl Write F2 Cl	ear F3 Recall	F4 Neu	patient	P7 His	tory	P10	Cancel
Ready for identification	n.	Impl			19/0	13/1997	15:59:3

10 Press the [ENTER] or [TAB] key to move to the Sub-Examination field.

11 Press [ENTER] or [TAB] if you want to select the currently displayed Sub-Examination type. You can also type the first character of the Sub-Examination name, e.g. 's' for 'skull'. You can also make your selection by clicking the downward arrow of the drop-down list box listing all the examination sub-types that are associated with the Examination type you selected in the previous step, as shown in the sample screen below.

器ADC2 IDENTIFICATION S	TATION (IDEN) IDENTIF	FICATION SCREEN	
File Soniganition Mode Per	torm Help		
PATIENT			
Patient Name		Firstname	
Birth date	00/00/0000	Patient ID	
8ех	Other •	Accession number	
Uwer Info			
STUDY			
Badiologist	ASPA DEPAULT	•	
Examination	SERVICE	 Sub examination C 	HEST PORTABLE 📃
INAGE		<u>a</u>	HEST PORTABLE
Department	Ωκοτοδλ	81	KULL GENERAL
Patient position	AP 💌	E E	TOOT
Cassette orientation	Landscape 💌	Exposure CL	ass 200 💌
Uwer Info 1			
DESTINATIONS			
Mardcopy unit	None	Number of copies	z 1 💌
Processing station	ADC_PS1	Film format	BINXIOIN 💌
Archive station	None		
Send station	None		
Fi Write F2 Cl	an F3 Recall F	A New patient F7 Hi	story F10 Cancel
Ready for identification	L.	[3020]	19/03/1997 19:09:28

By selecting the examination sub-type you automatically set default parameters for the following fields in the Image section of the Identification screen, viz. patient position, cassette orientation and exposure class. You will have to modify these parameters only in exceptional circumstances.

After pressing the [ENTER] or [TAB] key for the selection of the examination sub-type, the cursor automatically moves to the 'Department' field in the Image section of the Identification screen.

12 Press [ENTER] or [TAB] to confirm the currently displayed Department or type the Department name you want to select.

The cursor automatically moves to the next field, the Patient Position.

13 Press [ENTER] or [TAB] to select the currently displayed Patient Position or type the first character of the Patient Position acronym, e.g. 'p' for PA, or select the Patient Position from the drop-down list box, using your mouse or the arrow keys.

ADC2 DENTIFICATION S	TATION (IDEN) IDENT	FICATION SCREEN	- O ×
PATIENT			
Patient Name		Firstness	
Birth date	00/00/0000	Patient 3	ID I
8ех	Other ·	Accession	nunber
User Info			
STUDY			
Radiologist	ASPA DEPAULT	-	
Exemination	SERVICE	 Sub examin 	ation SKULL GENERAL
INAGE			
Department	Orology		
Patient position	IL ·		
Cassette orientation	A.P	Ехро	sure Class 200 💌
User Info 1	PA		
DESTINATIONS	RL 3		
Mardcopy unit	None	Number o	f copies 1 💌
Processing station	ADC_PS1 *	Film for	mat BINX10IN .
Archive station	None		
Send station	None		
Fi Write F2 Cl	ear F3 Recall F	4 New patient	F7 History F10 Cancel
Ready for identification	i	190401	19/03/1997 19:37:16

Type [ENTER] or [TAB] to move the cursor to the following Cassette Orientation field.

14 In the Cassette Orientation field, press [ENTER] or [TAB] to accept the default cassette orientation, type the first character L for 'landscape' or P for 'portrait' or select the appropriate cassette orientation from the drop-down menu.

ADC2 DENT FICATION S	TATION «IDEN» IDENTIFIC	ATION SCREEN
PATIENT		
Patient Name		Firstname
Birth date	00/00/0000	Patient ID
8ех	Other •	Accession number
User Info		
STUDY		
Radiologist	ASPA DEPAULT	•
Examination	SERVICE	 Sub examination SKULL GENERAL
DIGAGE		
Department	Ωκοτοδλ	
Patient position	IL ·	
Cassette orientation	Landscape 💌	Exposure Class 200 -
Uwer Info 1	Landscape	
DESTINATIONS	Portrait 👌	
Mardcopy unit	None	Number of copies 1 -
Processing station	ADC_PS1 💌	Film format BINX10IN .
Archive station	None	
Send station	None	
F1 Write F2 C1	sar F3 Recall F4	New patient F? History F10 Cancel
Ready for identification	1. 1	AUTO 19/03/1997 19:52:37

After pressing [ENTER] or [TAB] the cursor automatically moves to the adjacent Exposure Class field.

15 In the Exposure Class field, press [ENTER] or [TAB] to accept the default value, type the desired speed class or select the exposure class by scrolling down the drop-down list box.

ADC2 DENTIFICATION S	TATION «IDEN» IDENTIF	FICATION SCREEN	
PATIENT			
Patient Name		Firstname	
Birth date	00/00/0000	Patient ID	
8ех	Other •	Accession number	
User Info			
STUDY			
Radiologist	ASPA DEPAULT	•	
Examination	SERVICE	 Sub examination SKULL GENERAL 	×
INAGE			
Department	θεοτοβλ		
Patient position	IL ·		
Cassette orientation	Portrait •	Exposure Class 200 💌	
User Info 1		25	
DESTINATIONS		50	
Mardcopy unit	None *	Number of copies 1 100	
Processing station	ADC_PS1	Film format 8:150	
Archive station	None	200	
Send station	None	600 B	
Fi Write F2 Cl	sar F3 Recall F4	Wew patient F7 Histo 800 Can	zel
Ready for identification	1-	(30/TO) 19/03/1997 19:	58132

After confirming your selection, the cursor automatically moves to the User Info field, provided you have configured and placed a User Info field within the Image section of your Identification screen.

16 Fill in any comments (up to 40 alphanumeric characters) that you want to append and press [ENTER] or [TAB]. With each new patient you identify, the User Info field will be cleared.

17 After pressing [ENTER] or [TAB], the cursor moves to the Hardcopy Unit field in the Destinations section of your Identification screen. Type the destination name or select the Hardcopy Unit destination from the scroll-down list, as shown.

器ADC2 IDENTIFICATION S	TATION (DEN) IDEN	TIFICATIO	N SCREEN				
File Scriganitist Mode Per	form Help						
PATIENT							
Patient Name			Firstname				
Birth date	00/00/0000		Patient 1	D			
8ех	Other	-	Accession	numbe	r		
User Info							
STUDY							
Radiologist	ASPA DEPAULT						
Examination	SERVICE	• 81	ab examina	ation	SKAFF 05	NERAL	*
INAGE							
Department	θεοτοδλ						
Patient position	IL ·						
Cassette orientation	Portrait	-	Екрос	sure c	lass 25		
Uwer Info 1							
DESTINATIONS							
Mardcopy unit	D13000	*	Number of	f copi	az 1 💌		
Processing station	ADC_LR1		Film for	awt	14183	17IN	*
Archive station	DRYSTAR						
Send station	P13000 ···						
Fi Write F2 cl	ear F3 Recall	F4 Res	patient	17 1	istory	r10	Cancel
Ready for identification	1.	190701	1		19/0	3/1997	20122139

Press [ENTER] or [TAB]. The cursor moves to the adjacent Number of Copies field.

Note When 'None' is selected as a Hardcopy Unit destination, the next two fields, viz. The Number of Copies field and the Film Format field remain blank.



18 To select the number of hardcopies you want, enter the number of copies you want or select a number from the drop-down list box.

ADC2 DENTIFICATION S	TATION (IDEN) IDENTIFI	CATION SCREEN
PATIENT		
Patient Name		Firstnape
Birth date	00/00/0000	Patient ID
8ех	Other ·	Accession number
User Info		
STUDY		
Radiologist	ASPA DEPAULT	•
Exemination	SERVICE	 Sub examination SKULL GENERAL
INAGE		
Department	Orology	
Patient position	IT .	
Cassette orientation	Portrait 🔹	Exposure Class 25 💌
User Info 1		
DESTINATIONS		
Mardcopy unit	DRYSTAR *	Number of copies 1 -
Processing station	ADC_P51	Film format 0 101N .
Archive station	None	1
Send station	Bone	3
Fi Write F2 Cl	sar F3 Recall F4	New patient F7 His 4 F10 Cancel
Ready for identification	La .	6 03/1997 20131106

Press [ENTER] or [TAB]. The cursor moves to the Processing Station field.

19 In the Processing Station field, enter the Processing Station name or select the destination from the drop-down list box, as shown.

ADC2 IDENTIFICATION S	TATION (IDEN) IDENTIF	ICATION SCREEN	
The Soundganider: Mode Per	form Melp		
PATIENT			
Patient Name	Van Wesenbeeck	Firstname	Eric
Birth date	04/04/1947	Patient ID	040447
8ех	Nale 💌	Accession number	
User Info			
STUDY			
Radiologist	ASPA DEPAULT	•	
Examination	SYSTEM DIAGNOSIS	 Sub examination Lo 	is/Agfa CR contrast 💌
INAGE			
Department	AOFA		
Patient position	PA ·		
Cassette orientation	Portrait 🔹	Exposure Cla	es 200 💌
Uwer Info 1			
DESTINATIONS			
Mardcopy unit	ADC_LR1	Number of copies	
Processing station	ADC_PS1	Film format	BINXIOIN 💌
Archive station	ADC_PS1		
Send station	ADC P52		
Fi Write F2 Cl	sar F3 Recall F4	New patient F7 His	story F10 Cancel
Ready for identification	1.	190201	19/03/1997 21:28:15

Select 'none' if you don't want to send the image to the Processing Station. Press [ENTER] or [TAB] to confirm your choice and to move the cursor to the next field, the Film Format field. **20** Enter the desired film format or select the format from the drop-down list box, as shown in the sample screen below.

ADC2 IDENTIFICATION S	TATION (IDEN) IDENTIF	ICATION SCREEN	×
The Soundgardson Mode Per	form Melp		
PATIENT			
Patient Name		Firstname	
Birth date	00/00/0000	Patient ID	
8ех	Other •	Accession number	
User Info			
STUDY			
Radiologist	ASPA DEPAULT		
Examination	SERVICE	 Sub examination SKULL GENERAL 	•
INAGE			
Department	Ωκοτοδλ		
Patient position	IT .		
Cassette orientation	Portrait -	Exposure Class 25	
Uwer Info 1			
DESTINATIONS			
Mardcopy unit	ADC_LR1	Number of copies 2 -	
Processing station	ADC_P51	Film format SINXIOIN .	
Archive station	None	BINXLOIN	
Send station	None	11INXI4IN	
Fi Write F2 Cl	sar F3 Recall F4	New patient F? Mistory F10 Cancel	J
Ready for identification	1.	(WITO) 19/03/1997 2014515	9

Only the film formats that are supported by the previously selected Hardcopy Unit destination are displayed in the Film Format field. If you select Drystar 3000 as your Hardcopy Unit destination, for example, your options are limited to 11 x 14' and 14 x 17' film formats, and when selecting Drystar 2000, the only available choice is 8 x 10'.

1 ADDRESS INTERACTORS	INTER-USER BIR	FREE FREE GOVERNM	R E
EACTOR IN COLUMN	N. 194		
Fort start, Space		Columna .	
Barris data	41/14/0301	Patient St	
and a second	Oline 1		and an other states of the sta
Terry Tarley			1000 million - 1000 million
TTWIC .	Contract of the second		
And individual of	HIRA DEPUTAT	1	
Rowshill Berger	accesses.	· Ball Sciences	AND DECIDE LOOMERAN
Peace.			
Separates.	thickey?		
PORTAGE APRENTATION	14 E		
(seasons builded and	Portiaut 1	a angener	ee class 15 El
liter Into 1			
MUTTINITION .			
PROKING MUT	8.21 (W)	water of	COLUMN 1 P
PURPHASE STREET.	NC NL	TILL THINK	141801718
Andure eterioe	Rise		111001-000
Deal station	Acce .		111611704
IS NOT N.C.	Line #1 Becald	It Ber pataret	PT Ration Pill Calcel
Posts the search and	de-	100010	

ADC Compact ID Software and ID Tablet

- 2201A GB 199805 49



Press [ENTER] or [TAB]. The cursor moves to the Archive Station destinations field.

21 If an archiving system is connected to your ADC Compact system, select the Archive Station destination of your choice by entering its name or by selecting it from the drop-down list box.



Select 'none' if you do not want to send the image to the Archive Station. Press [ENTER] or [TAB] to move to the Send destination field.

22 If you have any Review Stations connected to your ADC Compact system, select the Review Station you want the image to be sent to, by entering its name or by selecting it from the drop-down list, as shown in the example below.

ADC2 IDENTIFICATION S	STATION (IDEN) IDEN	ATIFICAT	ON SCREEN			- [C] X
File Companyon Mode Per	rtava Help		A DECEMBER OF THE			
PATIENT						24
Patient Name	Van Werenbeerk		Firstness	6	Eric	
Birth date	04/04/1947	11419	Fatiant 1	ED	040447	
Ser	Male	•	Access10	number.		
User Info						100
TOTA						
Rediclogist	AGPA DEPAULT					
Examination	SYSTEM DIAGNOGIS		Sub exertin	ation Ld	is/Agfa CR o	ontrast .
DRAGE						
Department	A98A					
Patient position	7A .					
Cassette orientation	Portrait	-	Expo	sure cla	as 200 -	
User Info 1						
DESTINATIONS						
Hardcopy unit	ADC_181		Burber o	f copies	0 -	
Processing station	ADC_PS2	-	Film for	nat.	8 MX10IN	
Archive station	None					
Send station	None	2				
P1 Write P2 CJ	None	6 30	no patient	¥7 HL	story F10	Cancel
Ready for identificatio	ADC BC8	1 MUT	101		19/03/1997	21:41:14
	ADC_RS1					
	ADC R.S.L2					
	ADC BRS D					
	185 1888					

Select 'none' if you do not want to send the image to the Review Station. Press [ENTER] or [TAB].

23 Press the Write Function key to record all identification data onto the cassette memory chip. If no such Write Function key has been configured, select Write in the Perform drop-down menu, as shown in the sample screen below.

RADC2 IDENTIFICATIO	N STATION (II	DEN> IDENTIEN	ATION SCREEN			
File Contigunition: Mode	Perform Help					
PATIENT	Write					
Patient Name	Gear 🗟	enbeeck.	Firstness	E	ric	
Nicth data	Recal	47	Datiant 1	- D	40447	
Date of the care	New patient		Further of	. –	10111	
sex	History		Accession	number		
User Info	Hold/Unitiolo	_				
STUDY	Link					
Radiologist	Berouting	AULT	•			
Examination -	SISTER P	IAGNOSIS	 Sub examin 	ation Lds/3	lgfa CR co	ntrast 💌
INAGE						
Department	AUFA					
Patient position	PA 💌					
Cassette orientatio	on Portrait	-	Ехро	sure Class	200 -	
Uger Info 1						
DESTINATIONS						
Mardcopy unit	None		Number o	f copies 0	*	
Processing station	ADC_PS1		Film for	nat 14	inxi7in	•
Archive station	Mone					
Send station	None	*				
F1 Write F2	Cienr 73	Recall F4	New patient	17 Histor	ry 110	Cancel
Ready for identificat	:1.0EL-		(AUTO)		19/03/1997	22104111

After you have written the identification data to the memory chip in the cassette the application returns the following message:

Š			_ 🗆 🗵
Please	wait.	I'm	working

24 Remove the cassette from the cassette slot in the ADC Compact ID tablet and place it into the input buffer of the Digitizer.

After you have removed the cassette, the ID Station returns to the main Standby screen.

ADC2 IDENTIFICATION STATION «IDENTIFICATION SCREEN
*
Please enter a new cassette.
and click with your mouse
in this window or
press any key.

- 2201A GB 199805 53

Clearing information

By selecting the Clear option in the Perform menu, as shown in the sample screen below, you can clear any information contained in the field selected in the Identification screen. This command can be shortcut using a Function key. In the sample screen below the configured Clear Function key is [F2 Clear].



Recalling a patient

The Recall function can only be used in the Patient Identification fields, Comment field and User Info fields of the Identification screen. In order to recall some previous recorded information, proceed as follows:

- **1** Move the cursor to one of the patient identification fields, the Comment field or the User Info fields.
- 2 Press the Recall Function key, if you previously configured the Recall function into a Function key button, or select Recall from the Perform menu in the Identification screen, as shown below. Each time Recall is selected, the previously identified data of the last identifications is called back and used for the next identification procedure.

器ADC2 IDENTIFICATION S	STATION «IDEN» IDENTIF	CATION SCREEN	
File Configuration: Mode Per	form Help		
PATIENT	Wile		
Patient Name 0	lear enbeeck	Firstname	Eric
Birth date	atal 947	Patient ID	040447
8ex H	evropeen.	Accession number	
User Info H	kold/Unihold		
STIDY	ancel		
Badiologist	AUD/7	•	
Examination	SISTER PLASMOSIS	• Sub examination Lo	is/Aqfa CR contrast 💌
INAGE			
Department	AOFA		
Patient position	PA .		
Cassette orientation	Portrait •	Exposure Cla	as 200 -
Uwer Info 1			
DESTINATIONS			
Mardcopy unit	None V	Number of copies	0 -
Processing station	ADC_PS1	Film format	BINXLOIN ·
Archive station	None *		
Send station	None		
Fi Write F2 Cl	ear F3 Recall F4	New patient F? Hi	story F10 Cancel
Used date format: 30/MM	/1999	[3070]	20/03/1997 04:08:22

New Patient

The New Patient function clears all patient data fields in the Patient Identification section of the Identification screen. In order to initiate the New Patient function, proceed as follows:

1 Press the New Patient Function key (F4 in the sample screen below), provided you previously configured the New Patient function button, or select the New Patient option in the Perform menu of the Identification screen.

器ADC2 IDENTIFICATION S	STATION (IDEN) IDENTIF	ICATION SCREEN	
File Contrantor: Mode Per	form Help		
PATIENT	htte		
Patient Name C	loar	Firstname	
Birth date	ecal 000	Patient ID	
No. 1	ew patient		
OCA H	istay	Accession hund	
User Into	ancel		
STUDY	rk .	_	
Radiologist R	erouting AUL/I	•	
Examination -	SISTER PLASMOSIS	 Sub examination 	Lds/Agfa CR contrast 💌
INGAGE			
Department	Aofa		
Patient position	PA •		
Cassette orientation	Portrait •	Exposure (1ass 200 -
User Info 1			
DESTINATIONS			
Mardcopy unit	None	Number of copi	ana 🛛 💌
Processing station	ADC_PS1 ·	Film format	SINXLOIN 💌
Archive station	None		
Send station	None ·		
Fi Write F2 Cl	sar F3 Recall F4	New patient 17	History F10 Cancel
Ready for identification	1-	1904001	20/03/1997 04:19:50

2 Enter new patient data by either typing the data in or by using the Recall function or the patient data contained in the History list, as described in the next section of this manual.

Identifying a cassette with data from the history list

The ADC Compact ID application stores a history list with the last identification data. You can identify an ADC Compact cassette directly with patient information from this history list.

In order to use the information stored in the patient history list when identifying a cassette, proceed as follows:

- 1 Insert a cassette into the slot of the ADC Compact ID tablet. Click in the window of the Standby screen or press any key. The system returns the main identification screen.
- 2 Select Show History from the File menu or History from the Perform menu.



ADC Compact ID Software and ID Tablet



You can also press the History function key, provided you have previously configured a function key to that effect.

器ADC2 IDENTIFICATION S	STATION «IDEN» IDENTIFICATION SCREEN	
File Scritganitism Mode Per	rform Help	
PATIENT		
Patient Name	Firstname	
Birth date	00/00/0000 Patient ID	
8ех	Male · Accession number	
User Info		
STUDY		
Radiologist	AGPA DEPAULT	
Exemination	SYSTEM DIAGNOSIS Sub examination Lds/A	igfa CR contrast 💌
INAGE		
Department	AOFA	
Patient position	PA .	
Cassette orientation	Portrait Exposure Class	200 💌
User Info 1		
DESTINATIONS		
Mardcopy unit	None Mumber of copies 0	•
Processing station	ADC_PS1 Film format 81	nxi0in 💌
Archive station	None	
Send station	None	
Fi Write F2 Cl	Lear F3 Recall F4 New patient F7, Mistor	ry F10 Cancel
Ready for identification	G. (MITO)	19/03/1997 23:09:20

Accession number	i MO3A	
Sex	: Main	
Rigth date	: 19470404	
Patient ID	± 040447	
Firstname	: Eris	
Fatient Name	1 Van Wesenbeeck	
Accession number	: 19970317043614	
Sex	: Male	
Birth date	: 00010101	
Patient ID	i Unknown	
Firstname	: Unknown	
Fatient Name	+ Unknown	
Accession number	: 19970319215350	
8ex	i Male	
Birth date	: 19470404	
Fatient ID	+ 040447	
Firstname	: Eric	
Patient Name	: Van Mesenheeok	
Accession number	: 19970319220900	
Sex	+ Male	
Birth date	: 19470404	
Patient ID	a D40447	
Firstname	: Eric	
Patient Name	: Van Wesenbeeck	

The ID Station returns the history list, similar to the sample screen below:

3 Select the patient of your choice, then press the Transfer Item button.

The ADC Compact ID application copies the selected patient data from the list and moves it to the main Identification screen. If you want to remove a patient from the History list, click the Clear Item button. The Clear All button is used to remove all patient files from the History list.

- 4 Press the Close button. The selected patient is displayed in the main Identification screen. Fill in the remaining fields or press the Write Function key straight away if you wish to use the History List values already displayed.
- **5** Remove the cassette.

Link

The link option in the ADC Compact ID Software allows for the printing of two cassettes on one single film. The film layout is determined by the layout of the first cassette in the Link sequence.

The use of the Link feature is confined to 2/0 layouts (two images on one film).

In order to 'link' two cassettes onto one film, proceed as follows:

- **1** Insert a cassette in the transmission slot of the ADC Compact ID tablet.
- **2** Identify the cassette by entering all appropriate names and parameters, as described in this manual on pages 31 and following.
- **3** If you want the Link feature to be effective for a whole sequence of cassettes, select Link Mode in the Mode menu. The Link status is then indicated as [LINK] in the status bar at the bottom of the Identification screen.

ADC2 IDENTIFICATION 8	STATION (IDEN) IDENTIFI	CATION SCREEN	
PATIENT Normal of	ole		
Patient Name Virk nod	senbeeck.	Firstname	Eric
Birth date Romon	1947	Patient ID	040447
8ex 🗌	Nale •	Accession number	
User Info			
STUDY			
Radiologist	AGPA DEPAULT	<u> </u>	
Examination	SYSTEM DIASNOSIS	 Sub examination L 	ds/Agfa CR contrast 🗾
THOREE		_	
Department	AOPA		
Patient position	PA •		
Cassette orientation	Portrait •	Exposure CL	ass 200 <u>-</u>
User Info 1			
DESTINATIONS			
Nardcopy unit	None ·	Number of copie:	
Processing station	ADC_P51 ·	Film format	BINXIOIN
Archive station	None ·		
Send station	None ·		
F1 Write F2 Cl	sar F3 Recall F4	New patient F7 Hi	story F10 Cancel
Ready for identificatio	n	[AUTO]LINK]	19/03/1997 23:44:09

4 In case you want a 'once-only' linking of two cassettes and the selected Mode is Normal Mode, you can Link two cassettes by selecting the Link option in the Perform menu, as shown.

ADC2 IDENTIFICATION 9	STATION IDENS IDENT	FIGATION SCHEEN	10日日
PARTERIE ANOTHER	form field	_	
Patient Base Likenoo	a be pentered	Firstanes	Bric
Birth date	1947	Fatient IS	040447
Sez	Male -	Accession make	
leer Lafe			
STREY.			
adiologist /	ASPA DEFAULT		
Researching	STRTEM BIAGEORIE	· Dido estendariation	lds/Agfa CR contrast 💌
IIIAAE			
Department:	AGFA		
Patient position	PA #		
Cassette orisentation	Portrait .	Exposure C	Lass 200 *
neer rafe t	1		Contraction of the local sectors of the local secto
DESTIMATIONS	8 20		- Hereine
Bardonyy unit	Hone -	Number of orpla	
Processing station	ADC_PS1	Tile fatest	SINCLOIN .
Archive station	Boat #		
Bead station	Bone is		
F1 Meste F2 C3	eas #3 Becall #	M New potiant PT 1	tistory F10 Cancel
Ready for identification	L.	(MITE)	20/03/1807 24:10:04

ADC2 IDENTIFICATION S	TATION ODENS ID	DENTIFICAT	ION SCHEEN	
Pier Cardia beren Marte Ter	ten fiele			
PATIENT	Ada			
Patient Base	and anticent	1	Firstanan	Bric
Birth date	947	Sec. 1	Patient IS	040447
Sex .	tebry		Accession and	es i
lasr Lafe	DRAWING DRAWING	_		and the second sec
STUDY C	NCH .	_		
Saliologist .	ALL PARTY			
Renation -	NEWIAM STANDO	18 *	Diff. examinant to b	lds/Agfa CR contrast .
DIAME				5
be particulat	AGFA			
Patient position	PA #			
Cassette orisentation	Portrait		Reposition 1	class 200 -
neer rafe t		D ATE:		
DESTINATIONS	e · ·	101		
Bardenyy unit	Hone		Mundaer of cop-	188.0 .
Processing station	ADC_PS1		Tile fater	01N2101N .
Archive station	Boag	*		
Bend station	Bone	*		
F1 Meste F2 C3	ear Fl Recal	1 11 1	tw poblant FT	Ristory F1D Cancel
Ready for identification		130	niir mei	18/83/1887 23:45:13



5 After having identified the first cassette in the sequence, press the Function key or select the Write option in the Perform menu. The system then requests you to remove the first cassette, returning the following prompt.



6 Remove the first cassette from the ID tablet and insert the second cassette in the cassette slot of the ADC Compact ID tablet. The system returns the Identification screen for identifying the second cassette.

Enter the data that are relevant for the second cassette and press the Write button or select Write in the Perform menu.

端ADC2 IDENTIFICATION 9	TATION (DEN) DENTI	FICATION SCREEN	. O ×
File Contganitor: Mode Per	form Help		
PATIENT			
Patient Name	Van Wesenbeeck	Firstname	Eric
Birth date	04/04/1947	Patient ID	040447
8ех	Male •	Accession number	
Uwer Info			
STUDY			
Radiologist	ASPA DEPAULT	-	
Examination	SYSTEM DIAGNOS 18	 Sub examination L 	ds/Agfa CR contrast 🔳
INAGE			
Department	AOFA		
Patient position	AP 💌		
Cassette orientation	Landscape 🔹	Exposure Cl	ass 400 -
Uwer Info 1			
DESTINATIONS			
Hardcopy unit	None	Number of copie	
Processing station	ADC_PS1 ·	Film format	BINX10IN 💌
Archive station	None		
Send station	None 💌		
Fi Write F2 Cl	sar F3 Recall F	4 New patient 17 Hi	astory F10 Cancel
Ready for identification	1-	(AUTO (LUNK)	20/03/1997 24:36:41

During the identification of the second cassette, the patient demographic data and all other data except Sub-examination, patient position, cassette orientation, exposure class and user info can for obvious reasons not be altered. Neither can the destinations, as these have been determined during the identification of the first cassette.

The ADC Compact ID application writes all data to the memory chip embedded in the second cassette and returns to the Standby screen.

- 2201A GB 199805 63

Hold/Unhold

For a detailed description of the Hold/Unhold function, we refer to the treatment of the Hold function and Hold status in the description of the Configuration menu of the Standby screen, on page 24 of this manual.

Destination Rerouting

The Destination Rerouting option can only be accessed when the ADC Compact Autorouting Software has been installed on your system. For a full description of the ADC Compact Autorouting Software, we refer to the relevant manual.

Reading a cassette

The 'Read' function enables you to look at the data that was previously entered into the memory chip on the cassette. This is a read only-function. You cannot edit any information from the 'Read' screen. If you wish to edit information, you have to run the 'Identification' procedure again.

You can read an ADC Compact Cassette using the following procedure:

1 Select the Read option from the File menu in the Standby screen as shown in the sample screen below.



or

Select the Read option from the File menu in the Identification screen proper.

ADC2 IDENTIFICATION S	TATION (IDEN) IDENTIFIC form Help	CATION SCREEN	
Read cooseffe End-2 conserve Preferit list Show history Exit Uner: Inifia	00/00/0000 Other	Firstname Patient ID Accession numb	
STUDY			
Radiologist	ASPA DEPAULT	•	
Examination	SYSTEM DIAGNOSIS	 Sub examination 	Lds/Agfa CR contrast 💌
INAGE			
Department	AOFA		
Patient position	PA 💌		
Cassette orientation	Portrait ·	Exposure (1ass 200 -
User Info 1			
DESTINATIONS			
Mardcopy unit	None	Number of cop:	Lease D 💌
Processing station	ADC_PS1	Film format	BINXIOIN -
Archive station	None		
Send station	None		
Fi Write F2 Cl	ear F3 Recall F4	New patient 17	History F10 Cancel
Ready for identification	La I	(2010)	20/03/1997 24:57:19

- 2201A GB 199805 67

🛱 Read cassette data	
File	
Cassette	status
PATIENT	N 100
Patient Name	Firstname
Patient ID	Accession number
Birth date	8ex
STUDY	
Radiologist	Physician
Examination	Sub examination
Study ID	Study Inst. UID
Study Date	Study Time
Total number	Eody Fart
Processing code	Session number
Inage number	Total image number
IRAGE	
Department	
Patient position	Carrette crientation
Consent	User Info 0
4	,

The following screen is returned:

2 Insert the cassette you want to read into the ID tablet and press the Read button. If you want to quit the Reading procedure press the Exit button.
The data that is stored in the memory chip on the cassette is displayed on the screen, as in the sample screen below. The cassette status is shown in the upper cassette status bar. If the Cassette status is 'exposed' then the cassette has been identified but not cycled through the Digitizer.

Ca	ssette status	EXPOSED	
PATIENT		<u></u>	
Patiss/: Name	Van Genenbeeck	Firstnam	BESS
Patient 10	D40447	Accession number	19970320243058
litth date	19470404	llex	Hals
STUDY	1992 - 19		
Rediologist	ASPA DEPAULT	Physician	Unknown
Remainstan	SYNTEM DIAGNOSIS	Sub examination	Lds/Agfa CF contras
Study ID	Second Second	Budy Inst. WID	1.3.51.0.7.26.19970
Study Date	19970320	Study Tues	243058
Total sumber	2	Body Part	A
froessaing code	20	Supples number	26
Dange number	1	Total Lange subbr	1
DKAGE			10 10 10 10 10 10 10 10 10 10 10 10 10 1
Pepartnest	A/22A		12
Patient prostion	PA	Cassette prientation	Portrait
Comment	10 - 10	Heer lafo II	S. 199

The lower part of the screen displays additional initialization information, the number of cycles the plate has made since it was initialized, the type of image plate, the acquisition date, and so on.

Car	antha at store	EVBORED				
Caz	Sette status	EXPOSED				
Experiers Class	200					
DESTINATIONS						
Hardsopy unit	10	Number of moyies				
Processing station	ADC_PS1	Archive station	Concernance of the second seco			
Send station	Sec.	Filsleyout	1000			
West fleg	D					
CRESETTE INFO			as 10			
Connette muller	B97CVV	Initialization date	19970908			
Consette format	35CM243CM	Scan size	35cm243cm			
image plate type	AGFAID	Scan resolution	High resolution			
Acquisition Data	19970320	Acquisition Time	243058			
Usage count	0	Casaette status	2			
ID station mase	TIGH	character set	0			
Series number	1	Info Wersion	100			
Ctrl Version	200	PSDate Version	101			
Exposure dose	5	Sast of study	1			

- **3** Press [EXIT] to return to the identification screen.
- 4 Remove the cassette from the ID tablet slot.

Erasing an image plate

A cassette that has not been used for three days should be erased before being used again.

Caution An accidental or exceptional X-ray overdose excites the phosphor crystals in deeper layers of the image plate. The deeper layers then still retain a latent image after erasure. This image is called a 'ghost image'. This ghost image is caused by rebouncing, cosmic radiation and intrinsic radioactive elements in the phosphor. To remove the ghost image from an image plate you must select the Standard Erasure mode when passing the cassette in the Digitizer. For more details about the Digitizer Erasure modes, we refer to the ADC Compact Digitizer Manual.

The procedure to erase a cassette is twofold. First you must identify the cassette in the ADC Compact ID Station using the Erase function. Secondly, you must insert the cassette into the ADC Compact Digitizer. The Digitizer automatically recognizes this cassette and clears the image plate by exposing it to a powerful source.

Proceed as follows to erase a cassette:

1 Select 'Erase' from the File menu of the Standby screen.



2 You are prompted to insert the cassette into the ID tablet slot. Insert the cassette you want to erase into the ID tablet slot.

- 2201A GB 199805 71

The application returns a warning message, similar to the one displayed below:

≝ Era	ase cas	sette							
Are	you	\mathbf{sure}	that	you	\mathbf{want}	to	erase	\mathbf{the}	cassette?
				0.	к. С	ance	el		

- **3** Click O.K. if you wish to erase the cassette, otherwise select Cancel.
- **4** Remove the cassette from the ID tablet and place it in the input feed buffer of the ADC Compact Digitizer.

The plate of a cassette which is identified as Erase, will not be scanned by the Digitizer, and there will be no data transmission to the printer or to the Processing Station. Instead, the Digitizer will erase the image plate.



