


hiClot
User Manual



version March 19th, 2010

Chapter 1

General information

1.1 Layout of the device

hiClot is a single channel low cost clot meter, very simple to use and capable of performing the most common clot assays.

PT	Prothrombin time
PTT	Partial thromboplastin time
TT	Thrombin time
FIB	Fibrinogen (Clauss)

A piece of software is available¹ for connecting the **hiClot** to an external computer in order to manage an archive of the assay results and to print reports.

It is a small box that has an LCD back lit display, a keyboard and

- 12 positions to preheat the samples
- 2 positions for the reagents
- 1 reading well

¹The most recent version of this package can be downloaded from www.hilab.it



The device can be controlled by means of a three button keyboard and the display on top of it. The power switch and the plug for the external supply on the right side, the connector of the PC interface is on the left side.

1.2 Structure of the instrument


The *hiClot* upper side has the control panel and the thermostatic block.



1.2.1 Thermostatic block

The thermostatic block is made of aluminum and has some wells to preheat the samples and the reagents as well as the reading station. It is covered with a plastic surface to lower its dissipation thus assuring temperature uniformity which is between $36,5^{\circ}C$ and $37,5^{\circ}C$. The central well incorporates the optical reader and is where to insert the reaction cup to

be read.



1.2.2 Control panel


The control panel is composed by an alphanumeric LCD display and a three button keyboard. It is seldom used during the normal work as it serves only to select the assay type and start the assay.

The two upper keys ▲ and ▼ select the items of a list while the lower one ●_{OK} accepts the current selection.

1.3 Working principle

All the clot assay determine chain reactions that end by transforming the Fibrinogen into Fibrin using Thrombin as a catalyst. The Fibrin forms a high turbidity jelly clot that can be monitored by measuring the absorbance augmentation.


Actually, **hiClot** is fundamentally a mono channel blue light colorimeter that can be schematized as follows.



The semitransparent cup is crossed by the light beam generated by the emitter whose intensity variation is measured by the optical sensor on the opposite side.

The assay sequence begins by inserting a cup into the measuring well containing the needed volume of sample. The cup insertion is revealed by the photometer. When the **OK** key is pressed, a three second count down begins at the end of which the operator adds the reagent: from now on, **hiClot** starts monitoring the clot formation.

The typical absorbance graph after the reagent addition is shown in the following plot.



For about ten seconds nothing happens, then after 12 ÷ 13 seconds the sample starts clotting and the absorbance raises as it is clearly seen in the figure². As the Fibrin is being synthesized the absorbance raises until the upper limit that corresponds to the total clot formation.

hiClot follows completely the absorbance curve and builds the straight segments that control the stopwatch. This results in a good repeatability and improves **hiClot** insensibility to external disturbances.

²Please note that the vertical absorbance scale is arbitrary.

Chapter 2

Assays

2.1 Basics

```
hiClot 1.0  
(C) hiLab 2003
```

```
Heating  
>>>>>>>-----
```

```
Insert sample  
[PT ] PTT>
```

When switched on, **hiClot** shows its name and the software version for a while, then automatically starts heating the measurement block up to 37°C . During the heating phase, lasting for about $5 \div 8$ minutes depending on the external temperature, it displays an arrow that moves slowly rightwards as the temperature rises.

When the heating is over, the message `Insert sample` is shown on the upper row and on the second the method list begins with `PT PTT` and continues rightwards as implied by the arrow.

Please note that `PT` is surrounded by square parentheses to show that it is the current method, i.e. the one the instrument is ready to execute.

You can change the current method by means of the \blacktriangle and \blacktriangledown keys.

```
Insert sample  
PT [PTT] TT >
```

Pressing \blacktriangledown the selection goes to the next method: in the middle of the row now the `PTT` method is surrounded by the

parentheses. Pressing \blacktriangledown once more, the selection goes to the subsequent method (now TT is among parentheses).

Please note that the arrow to the right disappears so we are at the end of the list while an arrow is now to the left since the head of the list is no more shown. Pressing \blacktriangle goes back and the previous method is selected. The methods are ordered as follows: PT, PTT, TT e FIB.

In order to start an assay, a cup must be inserted into the measurement well: the photometer reveals its presence and starts waiting for the reagent dispensation. Now the display shows the selected method and Reagent. If the instrument is not calibrated for the chosen method, an asterisk is shown: in such case no other result but the clot formation time will be delivered.

Press OK and the displays shows a three second count-down. Insert the reagent when the countdown ends. the stopwatch starts and it will be stopped by the clot formation. After that the instrument shows the time and, if calibrated, the calculated results depending on the chosen method as indicated by the following table.

```
Insert sample
<PTT [TT ] FIB
```

```
PT Reagente
```

```
PT Attesa 007.1
```

```
PT 15.0 1.23
Inr 1.20 78.4%
```

test	parameter	display	
PT	Prothrombin time	sec	XXX.X
P	Ratio = PT/PT_{pool}		XX.XX
	INR = $Ratio^{ISI}$		XX.XX
	Activity	%	XXX.X
PTT	Partial Prothrombin time	sec	XXX.X
	Ratio = PTT/PTT_{pool}		XX.XX
TT	Thrombin time	sec	XXX.X
	Ratio = TT/TT_{pool}		XX.XX
FIB	Derived Fibrinogen	sec	XXX.X
	Concentration	mg/dl	XXX.

The last column to the right contains the significant figures: for instance XXX.X means that the result will be shown

by means of a number having one decimal figure and up to three figures before the decimal point (example 12.3)

If you take out the cup from the reading well during the clot formation, the measurement is interrupted and **hiClot** sets for a new reaction. The maximum time for the clot formation is 300 *sec* after which the test is anyway interrupted.

For the PT, PTT and TT assays, the ratio computations needs the PT_{pool} , PTT_{pool} and TT_{pool} coefficients that should be determined on a pool of sera that represents well the healthy population of the laboratory.



2.2 PT - Prothrombin time

The Prothrombin or factor II is a glycoprotein synthesized by the liver whose concentration in plasma is $10 \div 16 \text{ mg/dl}$.

The PT test detects the activity of the coagulation factors of the extrinsic pathway specifically the I,II, V, VII ad X some of which are influenced by the oral anticoagulant therapy. For this reason the PT assay is used for monitoring patients under such therapy.

2.2.1 Assay execution

1. Select the PT method¹
2. Preheat the PT reagent in a reagent well.
3. Insert a hollow cup in a preheating well.
4. Put 100 μL of sample in it (plasma)
5. Preheat the sample for at least 2 *min* or the time indicated on the reagent box.
6. Move the cup in the reading well: the system automatically senses its presence and asks the operator to insert

¹The selected method is the one between parentheses on the second row of the display. Select it by means of the  and 

```
Insert sample
[PT ] PTT>
```

```
PT *Reagent
```

the reagent.

- Press **OK** : a three second countdown starts. When it is over, dispense 200 μL of preheated reagent into the sample cup and the timer starts.
- Once the clot is formed, the timer stops and the results are shown.

PT *Wait 006.4

PT * 15.0

2.2.2 Calibration

A four point calibration is needed in order to calculate the PT Activity: the calibrators are obtained from the normal control plasma² diluting it with a salt solution as shown by the following schema:

	Activity %	Dilution	Preparation	
			plasma	solution
A	100	1 : 1	100 μL	-
B	50	1 : 2	50 μL	50 μL
C	25	1 : 4	25 μL	75 μL
D	12.5	1 : 8	12.5 μL	87.5 μL

Preheat the cups with the calibrators for no less than 2 *min* and then:

- Measure calibrator A (100 %) with the procedure described in the last paragraph. When the assay is over, press the **OK** key: CAL is displayed lower right to indicate that the first calibration point has been accepted.
- Measure the B calibrator (50 %) and press **OK** to enter the second point.
- Measure the C calibrator (25 %) and press **OK** to enter the third point.
- Measure the D calibrator (12.5 %) and press **OK** to enter the fourth point.


PT * 12.1 CAL

PT * 19.9 CAL

PT * 35.2 CAL

PT * 64.7 CAL

²which is supposed to have 100% activity



After such procedure the system is calibrated and can calculate the *Ratio* and the *Activity*. The calibration curve is the regression straight line based on the calibration points as shown in the figure with the inverse *Activity* ($1/A$) on the x axis and the clot time in seconds on the y axis.

Verification of the calibration data

```
PT T0 12.1
[PT ] PTT>

PT T1 19.9
[PT ] PTT>

PT T2 35.2
[PT ] PTT>

PT T3 64.7
[PT ] PTT>
```

The calibration data can be shown and changed pressing **OK** when the current test is PT The coagulation time T_0 of the calibrator A appears.. Pressing **OK** con and on, the coagulation time T_1 of the B calibrator is shown, then T_2 of the C calibrator and T_3 of the last calibrator.

2.2.3 Insertion of the ISI

If you press **OK** once more after the last calibration point, the current value of the ISI³ is displayed which, as known, is used in the calculation of the INR⁴. This parameter is nor-

³International Sensitivity Index

⁴International Normalized Ratio

mally found on the reagent box and generally depends on the production lot.

Let us suppose that the current ISI value is 1.00 and it should be 1.08. Proceed as follows

```
ISI 1.00
[PT ] PTT>
```

1. Press ▲ - the previous value disappears and a 1 appears : accept it pressing OK
2. Press ▼ - a decimal point is shown : accept it pressing OK
3. Press ▼ - besides the decimal point the figure 9 appears, pressing ▼ once more it turns into 8. Continue pressing ▼ until 0 is selected : accept it pressing OK
4. Press ▲ five times until besides 0 a 5 appears : accept the figure with OK
5. Now the selected number is correct: enter it pressing OK another time

```
ISI 1
[PT ] PTT>
```

```
ISI 1.
[PT ] PTT>
```

```
ISI 1.0
[PT ] PTT>
```

```
ISI 1.05
[PT ] PTT>
```

Actually the procedure consists in inserting a decimal number selecting each single digit: ▲ selects the next digit, ▼ the previous one. The decimal point lies after 9 and before 0 and the system inhibits to insert more than one.

2.2.4 Delivered results


When the instrument is well calibrated, the assay results are four values organized as follows.

PT	15.0	1.23
Inr	1.20	78.4%

- top left PT: clotting time in seconds (*15.0 in the example*)
- top right Ratio: ratio between the PT of the sample and that of the A calibrator (*1.23 in the shown example*)
- bottom left $INR = Ratio^{ISI}$ (*1.20 in the example*)
- bottom right Activity expressed as a percentage and based on the described calibration procedure (*78.4% in the example*)

2.2.5 Reset of the PT calibration

The aforementioned procedure is aimed at inserting automatically the calibration points from the calibrator assay results.

In order to recalibrate the instrument, it is mandatory to clear the current calibration: switch the instrument off and switch it on again keeping pressed the  button until the software version is disappears.


2.3 PTT - Partial Thromboplastin Time

It is the clotting time of a platelet deficient plasma after the insertion of phospholipids (in place of platelets) and of calcium ions. It verifies the extrinsic coagulation pathway.

2.3.1 Assay execution

1. Select the PTT assay
2. Preheat the $CaCl_2$ in a reagent well
3. Insert a void cup in a sample well

```
Insert sample
PT [PTT] TT >
```

4. Dispense $100\ \mu\text{L}$ of sample (plasma) and incubate for $2\ \text{min}$
5. Dispense $100\ \mu\text{L}$ of PTT reagent and incubate for another $3 \div 5\ \text{min}$
6. Move the sample cup to the reaction well: the system automatically senses its presence and asks the operator to insert the reagent.
7. Press  : a three second countdown starts and when it is over dispense $100\ \mu\text{L}$ of preheated *CaCl₂* into the sample cup. The stopwatch starts counting.
8. When the clot is developed, the stopwatch stops and shows the results



PTT Reagente



PTT Attesa 002.1



PTT 15.6 1.04

The first result is the PTT in seconds, the second is the ratio between this time and the normal PTT value.

2.4 TT - Thrombin time

It is the clotting time of a citrated platelet deficient plasma to which a known volume of Thrombin is dispensed. This assay investigates the final pathway of the coagulation reaction: the transformation of Fibrinogen into Fibrin and its polymerization.

2.4.1 Assay execution

1. Select the TT assay
2. Preheat the TT reagent in a reaction well.
3. Insert a void cup into a sample well
4. Pipet $200\ \mu\text{L}$ of sample (plasma)
5. Preheat the sample for at least $3\ \text{min}$ or the time you can find on the reagent box.



Insert sample
<PTT [TT] FIB

```
TT Reagente
```

```
TT Attesa 001.4
```

```
TT 15.0 1.07
```

6. Move the sample cup into the reading well: the system senses automatically its presence and asks the operator to insert the reagent.
7. Press **OK** : a three seconds countdown starts and, when it ends, pipet 100 μL of preheated reagent into the sample cup. The stopwatch gets started.
8. When the clot is formed, the timer is stopped and the results are shown

The first result is the Thrombin time in seconds, the second is the ratio between this value and the normal TT.

2.5 FIB - Fibrinogen (Clauss method)

The Fibrinogen or factor I is a dimer of three polymeric chains (α , β and γ) with a residual carboxylic terminal whose normal concentration in plasma is $200 \div 400 \text{ mg/dL}$. The Fibrinogen is converted into Fibrin by means of a polymerization catalyzed by Thrombin or other similar enzymes. **hiClot** uses the Clauss method to determine it that is an indirect assay based on a modified Thrombin Time. In presence of high Thrombin concentration, the polymerization time of Fibrin of a diluted citrated plasma is correlated to the concentration of that protein.

This reaction is influenced by additives having antithrombinic activity (FDP, heparin) or by changes in the structure of the Fibrinogen protein and in these cases it is underestimated.

2.5.1 Assay execution

WARNING For this assay the samples must be diluted 1:10 into IBS tampon.

```
Attesa campione  
<TT [FIB]
```

1. Select the FIB assay

2. Preheat the FIB reagent in a reagent well.
3. Insert a void cup into a sample well.
4. Pipet 200 μL of diluted sample (plasma) into it
5. Preheat the sample for at least 5 *min* or the time printed on the reagent box.
6. Move the sample cup into the reading well: the system automatically senses its presence and prompts the operator to insert the reagent.
7. Press **OK** : starts a three seconds countdown at the end of which pipet 100 μL of preheated reagent into the sample cup. The stopwatch gets started.
8. When the clot is formed, the timer stops and shows the results


2.5.2 Calibration

This assay needs a three point calibration with calibrators that are produced from a normal control plasma which is to be diluted with the following schema

	Dilution	Preparation	
		plasma	IBS
A	1 : 5	100 μL	400 μL
B	1 : 10	50 μL	450 μL
C	1 : 20	25 μL	475 μL

Preheat the calibrators for 2 *min* in the sample wells and operate as follows

1. Measure the A calibrator with the procedure described in the previous paragraph. When the assay is finished press **OK** : bottom right appears CAL to indicate that the first point of the calibration curve has been accepted.



2. Measure the B calibrator and press **OK** to accept the second value.
3. Measure the C calibrator and press **OK** to accept the third.

The calibration curve is the regression straight line approximating the calibration points placed on a bi logarithmic plane as indicated in figure where the logarithm of the concentration runs along the x axis (mg/dl) and the clotting time (sec) is on the y axis.

Verification of calibration data


The calibration data can be shown and edited pressing **OK** when the selected assay is FIB. The clotting time T_0 of the A calibrator. Pressing **OK** on and on the T_1 clotting time of the B calibrator appears and then the T_2 clotting time for the C calibrator.

```
<TT  FIB T0  6.6
      [FIB]
```

```
<TT  FIB T1 13.4
      [FIB]
```

```
<TT  FIB T2 27.5
      [FIB]
```

2.5.3 Insertion of the concentration

Pressing  after the last calibration point has been shown, the concentration in mg/dl of the control plasma is displayed. Insert the new values entering each digit as seen for the PT assay.

2.5.4 Delivered results

When the **hiClot** is calibrated for FIB with the procedure described in the previous paragraph, the following results are delivered


top left Clotting time in seconds (*16.3 in the example*)

top right Ratio between the sample time and that of the A calibrator. (*1.21 in the example*)

bottom Fibrinogen concentration in mg/dl (*228 in the example*)

2.5.5 Reset of the FIB calibration

The calibration procedure described permits to insert automatically the calibration values measuring the calibrators.

In order to recalibrate the instrument it is mandatory to reset the current calibration data. This can be done switching the instrument off and switching it on again pressing the  key until the software version disappears.

2.6 General advice

1. Cold reagents stored in the refrigerator should stay at least 5 minutes in the heated reagent cups before usage.
2. The insertion of reagents or starters at the end of the countdown should be fast. If the insertion is too slow, poor mixing of the components turns into wrong assay results.
3. Please check properly the pipettes. Non repeatable results normally are due to poor repeatability of liquid dispensation.
4. Preparation of the calibration samples by means of dilutions is critical. Errors made in this phase cause wrong measurements.
5. This is a manual instrument, so good results depend on proper manipulation of reagents and samples and on good liquid dispensation as well. The operator should be trained in using both the reagents and the pipettes.
6. During the assay, the second row of the display shows the time elapsed to the right, the measured absorbance in an arbitrary scale to the left.
7. Each **hiClot** is identified by a serial code, written on a sticker under the instrument, composed by HC, four digits a dash and four more digits. This code is needed to download the software package to interface the **hiClot** to a PC.



Chapter 3

Security and installation rules

In order to use the **hiClot** at its best, please follow carefully the security rules and installation advices contained in the following paragraph.

3.1 Security rules

1. Read carefully and preserve this manual.
2. Before cleaning, disconnect all the power plugs.
3. For cleaning the instrument, do not use cleaning liquids, aerosols or any other detergent that can penetrate into the it: use only a damp cloth.
4. Do not modify in any way the instrument without an explicit authorization of the builder.
5. Do not store or use the system in a wet or excessively damp environment.
6. Do not place the system on an unstable plane.

7. The back and bottom of the unit has holes aimed at cooling the instrument: keep them clear.
8. Connect the system only to a power supply conforming the rules and able to deliver the needed power.
9. An uninterruptible power supply (UPS) is strongly recommended.
10. A ground connection of the power lines is required for a correct and safe usage of the system.
11. Place the system as near as possible to the power plug.
12. Do not lay anything on the cables.
13. Do not surcharge the power plug using multiple plugs.
14. Do not insert objects into the system.
15. The system has been tested for a typical room temperature of 25°C ; in case the room temperature is much higher than this, air conditioning is strongly recommended.
16. Do not open the system: any repair should be carried out by trained personnel.
17. Disconnect the system from the power net and call the technical assistance in case of danger.
18. Use only original spare parts for reparation.
19. The system must be used only by trained personnel.
20. Do not use flammable matter near the instrument.

3.2 Package content and installation rules

The **hiClot** is forwarded inside a carton box that contains a shockproof soft and elastic plastic frame that itself contains the instrument.

The box contains:

- The **hiClot**
- The power supply

- A flyer containing the basic instructions
- 30 sample cups
- 6 reagent cups



3.3 Declaration of conformity


The builder company hiLab S.r.l. (www.hilab.it), sited in via Piero Aloisi 36, 00158 ROMA, ITALY, hereby declares that the instrument named **hiClot** conforms to the following EC directives:

89/336/CEE 92/31/CEE 93/68/CEE 93/97/CEE
73/23/CEE

It conforms moreover to the product requisites as stated by the following rules

CEI EN 50082-1
CEI EN 50081-1
CEI EN 60335-1

hiLab S.r.l. - l'amministratore unico
ing. Francesco Iovine



Roma, li 5/5/2004

Index

1	General information	3
1.1	Layout of the device	3
1.2	Structure of the instrument	4
1.2.1	Thermostatic block	4
1.2.2	Control panel	5
1.3	Working principle	5
2	Assays	8
2.1	Basics	8
2.2	PT - Prothrombin time	10
2.2.1	Assay execution	10
2.2.2	Calibration	11
2.2.3	Insertion of the ISI	12
2.2.4	Delivered results	13
2.2.5	Reset of the PT calibration	14
2.3	PTT - Partial Thromboplastin Time	14
2.3.1	Assay execution	14
2.4	TT - Thrombin time	15
2.4.1	Assay execution	15
2.5	FIB - Fibrinogen (Clauss method)	16
2.5.1	Assay execution	16
2.5.2	Calibration	17

2.5.3	Insertion of the concentration	19
2.5.4	Delivered results	19
2.5.5	Reset of the FIB calibration	19
2.6	General advice	20
3	Security and installation rules	21
3.1	Security rules	21
3.2	Package content and installation rules	22
3.3	Declaration of conformity	24