

User Manual





Ритес СмвН

Stauffenbergring 21 D-37075 Göttingen Germany

Tel. +49(0)551 205 29 170 Fax +49(0)551 205 29 176 +49(0)2232 152849

Email: info@pivtec.com

Homepage: http://www.pivtec.com

Version 2.4 - 15 April 2006

Copyright © 2001-2006 PIVTEC GmbH. All Rights reserved worldwide. No part of PIVview software or documentation may be reproduced, transmitted, transcribed, disassembled, reverse engineered, stored in a retrieval system, or translated in any form without the express written permission of PIVTEC GmbH.

Warranty

PIVTEC GmbH (PivTec) warrants that the PIVview computer program and documentation will substantially conform to published specifications. PIVTEC also warrants that the media (CD-ROM) used to transfer the software is free from defects in material and workmanship and that the software is free from substantial programming errors for a period of six (6) months from the date of purchase, unless a longer period is required by local law. During this period defective media will be replaced and substantial programming errors in the software will be corrected by PIVTEC with no charge. If PIVTEC is unable to replace defective media or correct substantial programming errors within sixty (60) days after notification of the defect or error, PIVTEC will refund the license fee to the licensee. There are your sole remidies for any breach of warranty.

Disclaimer of Warranties

YOU EXPRESSLY ACKNOWLEDGE AND AGREE THAT USE OF THE PIVTEC SOFTWARE IS AT YOUR SOLE RISK AND THAT THE ENTIRE RISK AS TO SATISFACTORY QUALITY, PERFORMANCE, ACCURACY AND EFFORT IS WITH YOU. EXCEPT FOR THE LIMITED WARRANTY ON MEDIA SET FORTH ABOVE AND TO THE MAXIMUM EXTENT PERMITTED BY APPLICABLE LAW, THE PIVTEC SOFTWARE IS PROVIDED "AS IS", WITH ALL FAULTS AND WITHOUT WARRANTY OF ANY KIND, AND PIVTEC HEREBY DISCLAIMS ALL WARRANTIES AND CONDITIONS WITH RESPECT TO THE PIVTEC SOFTWARE, EITHER EXPRESS, IMPLIED OR STATUTORY, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES AND/OR CONDITIONS OF MERCHANTABILITY, OF SATISFACTORY QUALITY, OF FITNESS FOR A PARTICULAR PURPOSE, OF ACCURACY, AND NON-INFRINGEMENT OF THIRD PARTY RIGHTS. PIVTEC DOES NOT WARRANT, THAT THE FUNCTIONS CONTAINED IN THE PIVTEC SOFTWARE WILL MEET YOUR REQUIREMENTS, THAT THE OPERATION OF THE PIVTEC SOFTWARE WILL BE UNINTERRUPTED OR ERROR-FREE, OR THAT DEFECTS IN THE PIVTEC SOFTWARE WILL BE CORRECTED. NO ORAL OR WRITTEN INFORMATION OR ADVICE GIVEN BY PIVTEC OR A PIVTEC AUTHORIZED REPRESENTATIVE SHALL CREATE A WARRANTY. SHOULD THE PIVTEC SOFTWARE PROVE DEFECTIVE, YOU ASSUME THE ENTIRE COST OF ALL NECESSARY SERVICING, REPAIR OR CORRECTION. SOME JURISDICTIONS DO NOT ALLOW THE EXCLUSION OF IMPLIED WARRANTIES OR LIMITATIONS ON APPLICABLE STATUTORY RIGHTS OF A CONSUMER, SO THE ABOVE EXCLUSION AND LIMITATIONS MAY NOT APPLY TO YOU.

Limitation of Liability

TO THE EXTENT NOT PROHIBITED BY LAW, IN NO EVENT SHALL PIVTEC BE LIABLE FOR PERSONAL INJURY, OR ANY INCIDENTAL, SPECIAL, INDIRECT OR CONSEQUENTIAL DAMAGES WHATSOEVER, INCLUDING, WITHOUT LIMITATION, DAMAGES FOR LOSS OF PROFITS, LOSS OF DATA, BUSINESS INTERRUPTION OR ANY OTHER COMMERCIAL DAMAGES OR LOSSES, ARISING OUT OF OR RELATED TO YOUR USE OR INABILITY TO USE THE PIVTEC SOFTWARE, HOWEVER CAUSED, REGARDLESS OF THE THEORY OF LIABILITY (CONTRACT, TORT OR OTHERWISE) AND EVEN IF PIVTEC HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. SOME JURISDICTIONS DO NOT ALLOW THE LIMITATION OF LIABILITY FOR PERSONAL INJURY, OR OF INCIDENTAL OR YOU YOU FOR ALLOW THE LIMITATION MAY NOT APPLY TO YOU. In no event shall PIVTEC's total liability to you for all damages (other than as may be required by applicable law in cases involving personal injury) exceed the original purchase cost of the PIVTEC software. The foregoing limitations will apply even if the above stated remedy fails of its essential purpose.

PIVTEC RESERVES THE RIGHT TO MAKE CHANGES FROM TIME TO TIME IN THE CONTENTS OF PIVTEC SOFTWARE WITHOUT OBLIGATION OF PIVTEC TO NOTIFY YOU AS A LICENSE HOLDER OF SUCH REVISIONS OR CHANGES.

Acknowledgements

The PIVTEC Software was and continues to be developed in close cooperation with the Institute of Aerodynamics and Flow Technology of the German Aerospace Center (DLR) in Goettingen, Germany.

Copyrights and Trademarks

netCDF — (network Common Data Format) Data Format Software Library and Utilities, © 1993-2005 University Corporation for Atmospheric Research/Unidata. All right reserved. (http://www.unidata.ucar.edu)

Tecplot — Amtec Engineering, Inc.

Microsoft, Visual C++, Windows, Windows Metafile — Microsoft Corporation.

All other product names are trademarks or registered trademarks of their respective owners.

Contents

Background	1
Features	1
Upcoming features	1
What <i>PIVview</i> is not	1
PIVview Versions	2
Credits	2
System Requirements	3
Installation Method A	3
Installation Method B	3
What is installed where ?	3
Removing PIVview	3
Getting acquainted with PIVview	5
Title bar	5
Menu bar	5
Main window	6
Main Window - Reference Vector	6
Main Window - Scroll Bars	6
Status Bar	6
The PIVview Menu Entries	7
The File Menu	7
The Image Menu	7
The PIV menu	8
The <i>Plot</i> menu	9
The <i>Window</i> menu	9
The <i>Help</i> menu	9
The PIVview Toolbars	10
The PIVview Dialogs	11
The Image Pre-Processing Menu	11 11
The Image Pre-Processing Menu Tab General	11 11 11
The PIVview Dialogs The Image Pre-Processing Menu Tab General Tab Filters	11 11 11 11
The PIVview Dialogs The Image Pre-Processing Menu Tab General Tab Filters Tab Background	11 11 11 11 12
The Image Pre-Processing Menu Tab General Tab Filters Tab Background Tab Mask	11 11 11 11 12 12
The Image Pre-Processing Menu Tab General Tab Filters Tab Background Tab Mask	11 11 11 12 12 12
The Image Pre-Processing Menu Tab General Tab Filters Tab Background Tab Mask PIV Evaluation Setup PIV Evaluation Setup – Grid	11 11 11 12 12 13 13
The Image Pre-Processing Menu Tab General Tab Filters Tab Background Tab Mask PIV Evaluation Setup PIV Evaluation Setup – Grid PIV Evaluation Setup – Correlation	11 11 11 12 12 13 13 14
The Image Pre-Processing Menu Tab General Tab Filters Tab Background Tab Mask PIV Evaluation Setup PIV Evaluation Setup – Grid PIV Evaluation Setup – Correlation. Phase-correlation (Phase only correlation)	11 11 11 12 12 13 13 14 14
The Image Pre-Processing Menu Tab General Tab Filters Tab Background Tab Mask PIV Evaluation Setup PIV Evaluation Setup – Grid PIV Evaluation Setup – Correlation PIV Evaluation Setup – Correlation Multiple Correlation	11 11 11 12 12 13 13 14 14 14
The Image Pre-Processing Menu Tab General Tab Filters Tab Background Tab Mask PIV Evaluation Setup PIV Evaluation Setup – Grid PIV Evaluation Setup – Correlation Phase-correlation (Phase only correlation) Multiple Correlation PIV Evaluation Setup – Algorithm	11 11 11 12 12 13 13 13 14 14 14
The Image Pre-Processing Menu	11 11 12 12 12 13 13 13 14 14 14 15 17
The Image Pre-Processing Menu Tab General Tab Filters Tab Background Tab Mask PIV Evaluation Setup PIV Evaluation Setup – Grid PIV Evaluation Setup – Correlation Phase-correlation (Phase only correlation) Multiple Correlation PIV Evaluation Setup – Algorithm Advanced PIV Processing options PIV Evaluation Setup – Peak Detection	11 11 11 12 12 12 13 13 13 14 14 14 15 17
The Image Pre-Processing Menu Tab General Tab Filters Tab Filters Tab Background Tab Mask PIV Evaluation Setup PIV Evaluation Setup – Grid PIV Evaluation Setup – Correlation Phase-correlation (Phase only correlation) Multiple Correlation PIV Evaluation Setup – Algorithm Advanced PIV Processing options PIV Evaluation Setup – Peak Detection PIV Evaluation Setup – Outliers.	11 11 11 12 12 12 13 13 13 13 14 14 14 15 17 18 20
The Proview Dialogs The Image Pre-Processing Menu Tab General Tab Filters Tab Background Tab Mask PIV Evaluation Setup PIV Evaluation Setup – Grid PIV Evaluation Setup – Correlation Plv Evaluation Setup – Correlation Plv Evaluation Setup – Correlation Plv Evaluation Setup – Algorithm Advanced PIV Processing options PIV Evaluation Setup – Peak Detection PIV Evaluation Setup – Conversion	11 11 11 12 12 12 13 13 13 14 14 14 14 14 14 15 17 18 20
The Proview Dialogs The Image Pre-Processing Menu Tab General Tab General Tab Filters Tab Background Tab Mask PIV Evaluation Setup PIV Evaluation Setup – Grid PIV Evaluation Setup – Correlation PIV Evaluation Setup – Correlation PIV Evaluation Setup – Correlation PIV Evaluation Setup – Algorithm Advanced PIV Processing options PIV Evaluation Setup – Peak Detection PIV Evaluation Setup – Outliers PIV Evaluation Setup – Filters etc.	11 11 11 12 12 12 13 13 13 14 14 14 15 17 18 20 22
The Image Pre-Processing Menu Tab General Tab Filters Tab Background Tab Mask PIV Evaluation Setup PIV Evaluation Setup – Grid PIV Evaluation Setup – Correlation. Phase-correlation (Phase only correlation) Multiple Correlation PIV Evaluation Setup – Algorithm Advanced PIV Processing options PIV Evaluation Setup – Peak Detection PIV Evaluation Setup – Outliers PIV Evaluation Setup – Conversion PIV Evaluation Setup – Conversion PIV Evaluation Setup – Filters etc. Region of Interest	11 11 11 12 12 12 13 13 13 13 14 14 15 17 18 20 22 23 24
The Image Pre-Processing Menu Tab General Tab Filters Tab Background Tab Mask PIV Evaluation Setup PIV Evaluation Setup – Grid PIV Evaluation Setup – Correlation. Phase-correlation (Phase only correlation) Multiple Correlation PIV Evaluation Setup – Algorithm Advanced PIV Processing options PIV Evaluation Setup – Peak Detection PIV Evaluation Setup – Outliers. PIV Evaluation Setup - Outliers. PIV Evaluation Setup - Conversion PIV Evaluation Setup – Pilters etc. Region of Interest. PIO Configuration Menu	11 11 11 12 12 12 13 13 13 13 13 13 13 13 13 13 14 14 14 15 20 22 23 24
The Proview Dialogs The Image Pre-Processing Menu. Tab General Tab General Tab Filters Tab Background. Tab Mask PIV Evaluation Setup. PIV Evaluation Setup – Grid PIV Evaluation Setup – Correlation. Phase-correlation (Phase only correlation) Multiple Correlation PIV Evaluation Setup – Algorithm Advanced PIV Processing options PIV Evaluation Setup – Peak Detection PIV Evaluation Setup – Peak Detection PIV Evaluation Setup – Outliers. PIV Evaluation Setup – Conversion. PIV Evaluation Setup – Filters etc. Region of Interest. Plot Configuration Menu. Vector Plotting Parameters.	11 11 11 12 12 12 12 13 13 13 14 14 14 14 14 14 15 20 22 23 23
The Proview Dialogs The Image Pre-Processing Menu Tab General Tab Filters Tab Background Tab Mask PIV Evaluation Setup PIV Evaluation Setup – Grid PIV Evaluation Setup – Correlation. Phase-correlation (Phase only correlation) Multiple Correlation PIV Evaluation Setup – Algorithm Advanced PIV Processing options PIV Evaluation Setup – Peak Detection PIV Evaluation Setup – Outliers PIV Evaluation Setup – Outliers PIV Evaluation Setup – Conversion PIV Evaluation Setup – Filters etc. Region of Interest Plot Configuration Menu Vector Plotting Parameters Outlier display	11 11 11 12 12 12 12 13 13 13 14 14 14 15 17 18 20 23 23 25 25
The Proview Dialogs The Image Pre-Processing Menu Tab General Tab General Tab Filters Tab Background Tab Mask PIV Evaluation Setup PIV Evaluation Setup – Grid PIV Evaluation Setup – Correlation. Phase-correlation (Phase only correlation) Multiple Correlation PIV Evaluation Setup – Algorithm Advanced PIV Processing options PIV Evaluation Setup – Peak Detection PIV Evaluation Setup – Outliers PIV Evaluation Setup – Filters etc. Region of Interest Plot Configuration Menu. Vector Plotting Parameters Outlier display Offset Vector Outlier display	11 11 11 12 12 12 13 13 13 13 14 14 14 14 15 20 22 23 25 25
The Proview Dialogs The Image Pre-Processing Menu Tab General Tab Filters Tab Background Tab Background Tab Mask PIV Evaluation Setup PIV Evaluation Setup – Grid PIV Evaluation Setup – Correlation Phase-correlation (Phase only correlation) Multiple Correlation PIV Evaluation Setup – Algorithm Advanced PIV Processing options PIV Evaluation Setup – Deak Detection PIV Evaluation Setup – Conversion PIV Evaluation Setup – Peak Detection PIV Evaluation Setup – Pilters etc. Region of Interest Plot Configuration Menu. Vector Plotting Parameters. Outlier display Offset Vector Color Coding	11 11 11 12 12 12 13 13 13 13 13 13 13 13 13 13 13 13 13 14 14 15 20 22 25 25 25
The Image Pre-Processing Menu Tab General Tab Filters Tab Background Tab Mask PIV Evaluation Setup PIV Evaluation Setup – Grid PIV Evaluation Setup – Correlation Phase-correlation (Phase only correlation) Multiple Correlation PIV Evaluation Setup – Algorithm Advanced PIV Processing options PIV Evaluation Setup – Peak Detection PIV Evaluation Setup – Outliers PIV Evaluation Setup – Conversion PIV Evaluation Setup – Pitters etc. Region of Interest Plot Configuration Menu Vector Plotting Parameters Outlier display Offset Vector Color Coding Vector Style	11 11 11 12 12 12 12 13 13 13 13 14 14 14 14 14 15 20 22 23 25 25 26 26
The PIVview Dialogs The Image Pre-Processing Menu Tab General Tab Filters Tab Background Tab Mask PIV Evaluation Setup PIV Evaluation Setup – Grid PIV Evaluation Setup – Correlation. Phase-correlation (Phase only correlation) Multiple Correlation PIV Evaluation Setup – Algorithm Advanced PIV Processing options PIV Evaluation Setup – Peak Detection PIV Evaluation Setup – Outliers PIV Evaluation Setup – Conversion PIV Evaluation Setup – Conversion PIV Evaluation Setup – Filters etc. Region of Interest Plot Configuration Menu. Vector Plotting Parameters Outlier display Offset Vector Contour Plotting Parameters	11 11 11 12 12 12 12 13 13 13 14 14 14 14 15 17 18 20 23 25 25 25 26 26
The Image Pre-Processing Menu Tab General Tab Filters Tab Background Tab Mask PIV Evaluation Setup Piterevaluation Setup Piteva	11 11 11 12 12 12 13 13 13 13 14 14 14 14 15 20 23 23 25 25 25 26 27 27

Axis Plotting Parameters	29
Reference Vector	30
Fonts / Colors	30
The Preferences Menu	31
File Extensions	31
Correlation Plane Window	32
Global Histogram Window	33
1-D PDFs and Peak Locking	34
Batch PIV Processing	35
Batch Processing Options	36
Batch Processing Output Options	36
Definition of statistics terms	37
PIVmask - Image Mask Editor	38
File Menu	38
View Menu	38
Draw Menu	39
Stereo PIV Modules	40
Stereo PIV Procedure in <i>PIVview3C</i>	40
Vector reconstruction in stereo PIV	41
Enabling the Stereo Mode	41
Stereo Menu Options & Toolbar	41
Creating a new stereo project with PIVmap3	43
Background on Stereo PIV	43
The program PIVmap3	45
Grid reconstruction in <i>PIVmap</i>	45
Calculation of the mapping function	47
PIVmap3 Input Menu	48
PIVmap3 Calibration Area	50
Helper axis	53
Description of the <i>PIVmap3</i> menu entries	53
Menu View	53
Menu Mapping	54
Menu Tools	54
Creating a new stereo project with StereoWizard	56
StereoGeom - Geometry Definition	58
PIVmapT01 - Projection Module	59
Some notes on obtaining calibration grid images	59
PIVmapT01 Toolbars & Menu Options	60
PIVmapT01 Tab Dialog	61
PIVmap101 Dialog - Detect Lines	62
PIVmap101 Dialog – Mesh & Map	63
PIVmapT01 Dialog – Projection	64
Stereo View Disparity Correction	65
Disparity Correction Tool Buttons:	65
Command-Line Utility Programs	67
	67
	67
	68
	70
	70
	71
	71
	71
SIEKEU	72
Shortcut keys in Pivview	73
MOUSE SHORTGUTS & ACTIONS	73

Supported File Formats	74
Supported Image Formats	74
Image naming and numbering conventions	74
Supported Output File Formats	
References and Further Reading	75
5	

Overview of PIVview

Background

PIVview is a compact program package for the evaluation of particle image velocimetry (PIV) or speckle displacement recordings. Its original and primary intention was to provide a quick-look data retrieval from PIV recordings. Experience has shown that this quick-look capability is essential in the optimization of the acquisition process.

Features

During its development *PIVview* has matured from a simple PIV image probing program to a (nearly) complete PIV evaluation package including:

- advanced processing algorithms:
 - multiple-pass method
 - mulitple-grid (pyramid) method
 - double-correlation method
 - second order, image deformation methods
 - various peak detectors and fitting algorithms)
 - ensemble correlation
- classical and nonlinear peak detection
- image pre-processing (various filters, masking)
- region of interest selection
- interactive calibration
- outlier detection and replacement
- data postprocessing (vorticity, strain)
- vector and contour display / print capability
- data export (TecPlot, ASCII and netCDF)
- PIV data import for further post-processing
- WYSIWYG export using BMP, PNG, and others
- batch processing of image and data sequences
- data set statistics calculation
- stereo PIV processing capabilities
- drag-and-drop capabilities
- mask editing
- command-line programs

Since version 2.0 *PIVview*'s user interface has been making use of the platform-independent QT-library (Trolltech, www.trolltech.com) and can be made available on most Unix systems (e.g., Linux, Solaris).

Upcoming features

Since *PIVview* is a continuously evolving program suite; new features and improvements are constantly added. *PIVTEC* is always open to new suggestions and encourages users to forward these as well as any bug reports.

What PIVview is not

PIVview does not include any form of data acquisition such as PIV camera and laser control, experiment synchronization, etc. This guarantees that the program can be operated on any PC-system running a 32-bit Windows operating system.

Also *PIVview* has no intention of replacing the capabilities of full-featured plotting packages such as *TecPlot* or *Origin*. The plotting capabilities of *PIVview* merely provide the user with visual feedback of the processing results while also providing the capability of quickly obtaining a hardcopy of the calculated data.

PIVview Versions

Two versions of *PIVview* are available: a demo version as well as a licensed version.

The demo version - downloadable from *PIVTEC*'s internet server - provides nearly all features of the licensed version except for data-export and batch-processing. The demo version can be distributed freely and is especially useful for optimizing a PIV experiment or for demonstrating the principles of PIV evaluation.

Full use of all of *PIVview's* features is only possible by obtaining a license through *PIVTEC* GmbH.

Credits

PIVview was conceived and is maintained by C. Willert, *PIVTEC* GmbH, Göttingen (Germany), and developed in close cooperation with the PIV-Group of the German Aerospace Center (DLR) in Göttingen.

System Requirements

PIVview should run on essentially any of today's *Windows* operating systems.

- 32-bit Microsoft Windows (tested on Windows versions 95, 98, NT4, 2000, XP)
- Pentium class processor
- >32 MB of RAM (depending on *Windows* version)
- 256-color SVGA or better graphics display
- hard disk space: program executables: ~5 MB example files: ~30 MB utility programs: ~4 MB

If you merely want to try out *PIVview* you need not install it on your system. The program can be run from the CD-ROM by double clicking the PivView2CDemo.exe file.

Alternatively *PIVview* can be installed from the setup CD-ROM in several ways:

Installation Method A

- 1. Put the *PIVview* CD-ROM into your CD-ROM drive.
- 2. In your file browser double click the file autorun.exe. This starts an installation option dialog.
- 3. From there the setup procedure can be initiated. Alternatively, the documentation or latest release notes can be viewed.
- 4. Follow the instructions in the installation program.

Installation Method B

- 1. Put the *PIVview* CD-ROM into your CD-ROM drive.
- 2. From the Windows Start button click Run.
- Type the following Run command line: D:\PIVVIEW2\SETUP\SETUP.EXE (If necessary substitute the appropriate drive letter for your CD-ROM drive in the above command.)
- 4. Follow the instructions in the installation program.

What is installed where ?

By default the *PIVview* program and its associated files are installed in the system's Programs or Program Files directory. Other destinations can be chosen during setup. All(!) files associated with *PIVview* are copied into this directory tree, that is, no files are copied to the Windows system directories. This also means that the program can be removed efficiently by deleting this directory although the system's uninstall method should be preferred.

Aside from the files, the installation programs installs a few short cuts to the system's start menu.

PIVview makes extensive use of the *Windows* registry to store its current configuration settings. These are located under:

\HKEY_LOCAL_MACHINE\Software\PivTec\PIVview2

Removing *PIVview*

The best way to remove *PIVview* is to use the Uninstaller which is part of the Windows operating system (found under Configuration\Software).

Alternatively the software is completely removed by deleting the *PIVview* installation directory. This however does not remove the registry entries associated with *PIVview*.

Getting acquainted with PIVview

The user interface of *PIVview* is similar to many common *Microsoft Windows* based applications and features:

- a title bar
- a menu bar
- several toolbars
- a main window to view data (workspace)
- a status bar

PIVview is based on a single document software architecure, that is, only one image (pair) can be opened by each *PIVview* instance. If you need to open several PIV images simultaneously you need to start several instances of *PIVview* which can be started through the menu option: **File/New Window**. (Please note that the settings of last *PIVview* instance closed are stored in the registry.)



Main window of PivView with image pair loaded and interactively sampled interrogation vectors. The intensity LUT is displayed on the left. The top-left shows a reference scale vector.

Title bar

The title bar displays the file name(s) of currently loaded image or data file(s) and *PIVview*'s current mode of operation (e.g. cross-correlation, autocorrelation or speckle shift mode).

Menu bar

The menu bar is used to perform any required operation within *PIVview* and contains seven sub-menus:

- File for loading images, saving data, setting preferences or exiting the program;
- *Image* for basic image display related functions (brightness, look-up tables, zoom, etc.) as well as image enhancement;

- *PIV* for performing PIV analysis, defining PIV interrogation parameters or regions of interest;
- *Plot* for configuration of plotting parameters or exporting a plot of the currently displayed data;
- *Windows* to activate additional data display windows such as correlation plane, 2-D and 1-D histograms, or data probing menu;
- Help to display PIVview version information or available program options.

Main window

Once loaded, PIV images are displayed in the main window. Image pairs are displayed by by selecting the brightest pixels from both images.

The image (pair) can be interactively interrogated at any position using the left mouse button. The interrogation parameters for cross-correlation analysis are defined in various setup menus. The results of the correlation analysis are displayed in the status bar at the bottom as well as in the form of a color-coded vector (defined in the Plot/Vector setup menu).

Main window - Lookup table

A lookup table (LUT) is used to convert image intensities to screen intensities or colors. The image above shows the currently selected LUT as a vertical bar on the left — in this case a linear intensity ramp with 0 (black) at the bottom and 255 (white) at the top.

The number at the top of the LUT represents the intensity of brightest display pixel. The image brightness can be changed by using the menu entries **Image/Brigher** or **Image/Darker** or the short-cut keys, b and d, respectively. Pixel intensities exceeding the maximum display intensity are either clipped at the maximum or wrap back to the bottom of the LUT (**Image/Clipping**).

The display of the LUT can be turned off with the menu option Image/Hide LUT.

Main Window - Reference Vector

The reference vector's main purpose is to have a visual reference when a complete data set is plotted or exported. It's magnitude and display is defined via **Plot/PlotConfig./Ref.Vector**.

Main Window - Scroll Bars

The scroll bars on the right side and bottom edge of the main window can be used to pan the image to a different viewport if the entire image does not fit into the display window.

Status Bar

The status bar generally reports the cursor position in both image units [pixels] as well as physical coordinates. Also the image intensity is given. In case of interactive interrogation, the correlation result is shown (e.g. position, displacement, velocity, correlation coefficient).

If the image has been evaluated these fields provide information about the data point closest to the cursor.

Note: Many of the toolbar buttons and menu entries only become active if image or data is available.

The *PIVview* Menu Entries

The File Menu

Evaluation Mode – specifies one of the three *PIVview* operation modes:

- Cross-Correlation Mode the standard mode for PIV analysis based on double-frame, single-exposure images (image pairs).
- Autocorrelation Mode PIV analysis based on single-frame, double-exposure images which are common in photographic PIV.
- Speckle Displacement Mode for the analysis of speckle displacement image pairs. This mode is similar to the cross-correlation mode but instead of converting the data to velocity it is presented as displacements. Also the two required images are loaded separately.

Preferences – opens a dialog which allows to set a variety of *PIVview*-global settings such as file extensions, default data output, etc.

Open image – opens a file dialog for loading a PIV image (pair) or speckle shift image.

Open data - opens a file dialog for loading a PIV or speckle shift data set.

Save data set – opens a menu with output options to save computed PIV data in a file. The **Save** menu option uses the file title (without extension) to create a name for the output file and places it in the same directory with the original image. With the **Save as** option the user can define the name and destination directory of the output file.

Create Tecplot Data – to save the processed data set in a format compatible with *Tecplot*. Style files for *Tecplot* are also generated.

Parameters... (Load, Save, Save as) – stores or retrieves the current *PIVview* parameter configuration in/from a file. This file is formatted in *netCDF* and therefore is platform independent but cannot be read be a standard text editor.

Recent images, Recent files – are listed by their names and are opened by clicking on them. Note that the PIV mode of operation does not automatically change for different image types.

New Window – starts are another instance of *PIVview* by opening additional main window. New windows operate completely independent of one another.

Exit terminates *PIVview* program. The current PIV, image and plotting parameters are stored in the *Windows* registry and are restored when the program is started again.

Tip: Use parameter files if you plan to work with different types of images.

The Image Menu

Image PreProcess – opens a dialog in which a variety of image pre-processing options can be set, such as splitting of the PCO-camera's double image files, image filtering, histogram adjustment, background subtraction, etc.

Edit Mask – starts the mask editor *PIVmask* with the currently loaded image. Completed mask images can be loaded directly from in the image pre-proceesing menu.

Calc. Minimum Image – allows the selection of an image sequence from which a minimum intensity image is computed.

Horizontal Profile, Vertical Profile – shows plots of image intensity profiles for the current



Img Pre mouse cursor position.



F.

<u>CAL</u>

ROI

Zoom out, **Zoom in** – change the image magnification. The scroll bars on the edges can be used to pan through the enlarged image / data set (Shortcut keys are + and -).

Brighter, **Darker** – changes the visual brightness of the currently loaded image. The image data itself is not affected (The shortcut keys are **b** and **d**).

Auto-scale - sets default intensity on image

Change LUT – changes the current image lookup table (LUT) or switches to previous LUT. Clicking on the function/button repeatedly cycles through all possible LUTs. Currently there is no provision for user defined LUTs. (Shortcut key is L).

Clipping – enables or disables image clipping for image intensties exceeding the LUT maximum intensity. With clipping on intensities greater than the maximum are set to the maximum. With clipped off the intensity wraps back to the bottom of the LUT (Shortcut key is c).

Show LUT – controls the display of the LUT bar on the left side of the image. (Shortcut key is CTRL+L).

Show Image – is useful to hide the image to better view the computed PIV data. (Shortcut key is I).

The PIV menu

- **Calibration** puts PIVview into calibration mode, hiding the current PIV data set and displaying a moveable ruler and origin. Within a separate dialog the magnification factor and reference point coordinates can be defined interactivley. After closing this menu PivView reverts back to the interrogation mode. (Shortcut key is **F3**).
- **PIV ROI** puts PIVview into the ROI-mode in which a rectangular portion within the image can be defined for processing. To use this function the *Use ROI* option in the **PIV Setup** dialog has to be enabled (Shortcut key is **F2**).

Process Log - shows a window with processing messages.

- **PIV Configuration** opens a tabbed dialog in which a wide variety of PIV processing parameters can be defined (Shortcut key is **F4** or via context menu by right-clicking the mouse).
- **Evaluate image** performs an evaluation of the loaded image (pair). For large data sets a progress bar is displayed during processing allowing the user to abort the interrogation. (Shortcut key is **F5**).
 - **Batch Processing** opens a dialog through which ensemble PIV processing can be initiated.
 - **Correlation Plane** displays the correlation plane and associated PIV image samples (if chosen) in a separate window. The data in this window is updated each time a new interactive interrogation is performed on the image (by clicking the left button within the image, Shortcut key is **F6**).
 - **2-D Histogram (PDF)** also is only available for evaluated PIV images. It plots the calculated displacement data in a two-dimensional plane and can be used in the vector validation procedure. Additionally one-dimensional histograms (PDFs) show the projection of the 2-D data onto each major axis. (Shortcut key is **F7**).



PIV Data Probe – is a window that is available after a PIV image has been evaluated. It gives detailed information about any calculated vector within the data set. The vector closest

<u>Σu</u>²

-Ū.,

A

to the mouse cursor is selected.

Statistics – provides information about the currently displayed PIV data set, such as mean and RMS values of displacement, velocity, vorticity, etc.

Live Evaluation – if enabled PIVview will display newly computed vectors during processing of the image.

Clear PIV Data – deletes the currently displayed PIV data set. Usually any modification in the **PIV Configuration** menu has the same effect.

The Plot menu

Plot Settings – opens a dialog in which the layout and properties of the plotted PIV data can be controlled. Among these are vector coding and length, scalar image display, axis and axis labeling, reference vectors, and color definitions (Shortcut key is P).

Subtract Mean – displays the vector field with or without the mean displacement of the data set subtracted. This is especially useful for displaying flow structures that might be hidden in a large, nearly uniform, convective flow.

Print Preview, Printer Setup, Print – handles printing within PIVview and includes previewing, layout and printer selection. (Shortcut key for printing is P).

Export as image> – stores the currently displayed data set or image in an image file (=screen capture). Export formats include BMP, JPEG, PNG and others. Current zoom factors are used in the exported image. The larger the zoom factor, the larger the resulting image.

Position objects – puts PIVview into and edit mode where workspace objects (image, luts, reference vector,...) can be freely placed to suit the user's preferences.

Save current LUT / Load LUT - save or load user specified pseudo-color lookup table

Default positions – restores the original positions of the workspace objects.

Save / Load plot configuration – stores/restores plot parameters

The Window menu

Hide Toolbars – hides all toolbars

Individual toolbars can be activated or deactivated one by one.

The Help menu

PIVview options – displays the command line options available for PIVview.

About PivView gives information about the program itself (e.g. release version number) as well as shows the currently registered end user.



Plot>Print – sends the currently displayed data set to the printer.



Plot>Position Objects – allows free position of workspace objects such as lookup-tables, image or reference vector.

The PIVview Dialogs

The Image Pre-Processing Menu

The purpose of this menu is to perform a variety of operations on the PIV image(s) prior to the PIV analysis. This includes basic image re-orientation, image enhancement using filters, background subtraction or mask overlays. Each of these are described in more detail below.

Me Image PreProcessing	? 🗙 Image PreProcessing
General Filters Background Mask	General Filters Background Mask
F Always Split PCO (B16) Images	Convert Images to Floating Point Data
Rotation Mirror	Invert (for scanned negatives)
None None	- 🥅 Median Filter
C 90° C Left <=> Right	Kemel size 3x3 💌
C 180° C Top <=> Bottom	
C 270°	- 🗖 Dynamic Thresholding (Histogram Clipping)
	Minimum 5% 🚔 Maximum 80% 🚆
	U High Pace Filter
	Kernel size 📴 🗹 Use Gaussian Weighting 📔 Liip at Zero
	- 🗖 Binarize (Thresholding) ————————————————————————————————————
	Threshold 5
	,
	- 🔽 Anti-Alias (Smoothing) ————————————————————————————————————
	🔿 Fixed 3x3 Gaussian 💿 Variable Gaussian Kernel [px] 0.5 🚔
	Cancel Apply OK Cancel

Tab General

Always Split PCO Images – PIV image pairs recorded with the PCO acquisition software are generally stored in a common file. The upper half of the image is associated with the first exposure while the lower half corresponds to the second exposure. This operation splits the input file into the two images. Otherwise the software will try to find the matching second image.

Rotation – rotates the image(s) in 90° increments. This operation as well as all others is performed on both images of a pair.

Mirror – can flip the image either left-to-right or top-to-bottom. Both flip operations together are equivalent to a 180° rotation.

Tab Filters

Filters are applied in the order they appear in the menu.

Convert Images to Floating Point Data – prevents truncation of data during filtering procedures (may increase overall precision).

Invert – computes the negative of the image by turning bright pixels into dark ones and vice-versa.

Median Filter – performs median filtering on the data using a 3x3 or 5x5 kernel.

Dynamic Thresholding – performs percentage-based thresholding of the image. This operation determines the minimum or maximum cut-off thresholds by summing up the individual bin-counts of the image histogram. The threshold is value is determined when the

percentage of pixels in the sum exceeds the percentage specified by the user. (Obviously this only works for images which a not binary.)

High Pass Filter – performs a high-pass filtering of the image by subtracting a lowpassfiltered version of the input image. Both the kernel size and filter-shape (Rectangle or Gaussian) of the low-pass filter can be specified. Gaussian filtering somewhat slower but positively influences the final result. Negative values are clipped only if this option is activated.

Binarization – binarizes the the image data by setting pixel values above the specified threshold to white (255) and lower values to black (0).

AntiAlias (Smoothing) – applies a low-pass filter which can be used to damp out spurious noise in the image.

Tip: The combination of the last three filters can significantly help in processing images with low signal to noise ratio.

Image PreProcessing	Image PreProcessing
General Filters Background Mask	General Filters Background Mask
Subtract Background	Use Mask(s)
C:/PIV/PIV_Challenge2005/B/min_img.tif	Choose Mask
🔲 Use seperate bkgd for frame B	E:/PIV/SampleData/Real_BL/WallMask3.bmp
	Use separate mask for image B Choose Mask2
	j invert Mask(s)
	Randomize masked pixels (experimental)

Tab Background

Subtract Background When activated the specified image or image pair is subtracted from all data images upon loading. Background subtraction can improve the data yield in images suffering from strong background artifacts. The image size of this background has to match that of the PIV image. Separate background images for each image of a pair may be specified.

Tab Mask

Use Mask When specified mask images are used to blot out data in regions covered by the mask. Separate mask images for each image of a pair may be specified.

Invert Mask Inverts the pixels of the mask images prior to application.

Randomize masked pixels This is an empirical method to reduce the influence of the masked areas on the correlation signal. Normally masked pixels are set to zero prior to correlation. This may result in edge artifacts in the correlation plane. In this case the masked area is filled with random data that should have roughly the spectral characteristics of the unmasked data and generally reduces the edge artifacts in the correlation plane.

PIV Evaluation Setup

PIV Evaluation Setup – Grid

This dialog can only be opened when an image (pair) is loaded. The image is then overlaid with a grid of crosses indicating the centers of the correlation samples (i.e. interrogation windows).

📅 PIV Evaluation Setup 💦 🗶
Grid Correlation Algorithm Peak Detection Validation Conversion Filters etc.
Horizontal Vertical Grid alignment Bottom edge
Window Size [px] 32 🔽 Grid Anchor Point
💽 Step Size [px] 16 🚔 16 🚔 X [px] 24 🚔 Y [px] 24 🚍
C Overlap 50% ▲ 50% ▲
Image Offset (px) 0 🚔 0
Current Grid Dimensions: X=(16, 1424) -> 89 Nodes Y=(16, 672) -> 42 Nodes => 3738 Nodes total
- 🗂 Use ROI
0 + 0 + ROJ 512 + 480 + 480 +
Apply Cancel OK

Window Size

The menu entry defines the interrogation window size, that is, the rectangular size with which the PIV images are sampled to perform the local cross-correlation analysis. Frequently the interrogation window is also referred to as interrogation spot. The window dimensions can be freely set (e.g. 24x27 is also possible).

Step Size

defines the increment with which the images are sampled. These values define the mesh size of the final data size. Typical values are around 50% of the window size which satisifies the sampling criterion. Greater sample overlap (e.g. smaller step sizes) produces over-sampled data but does not necessarily yield additional information.

Image Offset

Applies a constant integer offset to the interrogation windows by offsetting the images with respect to each other. This offset is applied to the displayed image (with the exception of double-exposure *autocorrelation* images). Interrogation window offset is especially useful in the analysis of images with large convective displacement (e.g. pipe flow, wind tunnel flow, etc.).

Grid alignment

Allow the grid to be aligned either with the bottom or top edge of the image. Alternatively the grid may be locked to any specified point (pixel) the image.

Use ROI

Restricts the evaluation to a defined rectangular region within the image. This region of interest (ROI) can either be set by entering boundary positions or interactively by using the ROI button on the toolbar or menu entry **PIV**>**Select ROI**.



Process Button – closes this menu and starts evaluating the image.

PIV Evaluation Setup – Correlation

PIV Evaluation Setup	?
Grid Correlation Algorithm Peak Detection Validation Conversion F	Filters etc.
Phase-correlation (Phase-only Filtering)	
Nyquist Frequency Filtering	
┌ 🔽 Multiple (repeated) Correlation	
Number of repeated correlations 9	
Image re-sampling	
C correlation re-sampling based on grid spacing X: 16 px Y: 16 px	
user specified correlation re-sampling offset X 4 px ★ Y 0 px ★ Y	
Correlation accumulation mode	
Multiplication of correlation planes Averaging of correlation planes	
🖄 🔀 🔼 🔄 🖄	<u>0</u> K

Phase-correlation (Phase only correlation)

This approach spectrally filters the data such that only the phase information is used in the calculation. The technique, proposed by Wernet (2005), has the potential of obtaining a higher correlation signal.

Nyquist Frequency Filtering

This option removes the highest frequency components in the cross-power spectrum (Fourier transform of the correlation function). This can be of use in the removal of camera noise artifacts (interlacing, multi-channel readout, etc.). The correlation plane window should be used in diagnosing the function's performance.

Multiple Correlation

This algorithm essentially calculates the correlation from at least two slightly spatially separated correlation planes. An approach based on multiplication of correlation planes was first proposed by Hart (2000). By calculating the product or sum of these correlation planes a new correlation plane results which ideally contains only the single displacement peak. Due to the random position of the noise peaks, which are present in all standard correlation planes, the final correlation plane will show no or strongly reduced noise peaks. Hence, this procedure effectively increases the detectability of the correlation peak associated with the particle shift.

The algorithm requires additional processing time depending on the number of additional correlations calculated and does not perform equally well on all images. A careful adjustment of the offset between the sample positions is sometimes necessary. If correlation planes are too far displaced from each other the correlation signal may be lost entirely in high gradient flows. The correlation plane window should be used in diagnosing the function's performance.

When the number of repeated correlations is higher (more than 3) then the averaging mode should be preferred to the multiplication, because this approach is less susceptible to loss of correlation in one of the repeated correlations.

PIV Evaluation Setup – Algorithm

Advanced Correlation Techniques

If none of the advanced algorithms are activated *PVview* uses the classical cross-correlation method as described in various PIV references (e.g. Willert & Gharib, 1991). More advanced algorithms are briefly described in the following.

📅 PIV Evaluation Setup ? 🗙
Grid Correlation Algorithm Peak Detection Validation Conversion Filters etc.
Interrogation method Multi-grid interrogation (Grid refinement) Multiple-pass options Number of iterations
Multi-grid options Initial sampling window 64 x 64 💽 Final window size: 32x32 px
Image interpolation options
Sub-pixel image shifting Enabled, on all passes 💌
Scheme B-Spline interpolation (best) 🔽 Order 3 🚍
Advanced options
No speed optimization
<u>Apply</u> <u>Cancel</u> <u>D</u> K

Interrogation method

Allows the selection of different PIV interrogation schemes:

- **Standard, single pass interrogation** does classical PV interrogation without additional passes to improve the result
- **Multiple-pass interrogation** The interrogation of the image is repeated at least once more. In the following passes the image sample positions are offset by the integer shift determined in the preceding pass. Once the residual shifts are less than one pixel a re-evaluation of the respective point is no longer necessary (e.g. convergence). The maximum number of interrogation passes can be specified. Two to three passes are generally sufficient for convergence.

Multiple pass interrogation not only increases the data yield due to the higher amount of matched particle images, but also reduces the so called *bias*-error (Westerweel et al. 1997). Outlier detection is used to validate intermediate data sets.

Multigrid interrogation uses a pyramid approach by starting off with larger interrogation windows on a coarse grid and refining the windows and grid with each pass. Outlier detection is used to filter the intermediate results.

Multiple-pass options

Number of iterations – specifies the number of passes during multiple-pass PIV interrogation.

Multi-grid options

Initial sampling window – specifies the starting interrogation window for multi-grid PIV interrogation. The final window size is set in the *Grid* menu tab.

Image interpolation options

Sub-pixel image shifting – enables sub-pixel based image shifting (image deformation) on multiple-pass- or multi-grid interrogation. This technique is a second order method and involves a complete deformation of the image data using the displacement data of the previous interrogation passes.

Interploation scheme – provides a choice of different pixel interpolation schemes as well as the order. In general cubic B-splines were found to produce the best results.

Speed optimization options

To increase processing speed intermediate processing step may use faster peak fitting algorithms or simpler image interpolation schemes. This may have a minor influence on the final result (e.g. reduced precision and/or reduced vector detection rate).



Flowchart for coarse-to-fine multi-grid processing with image deformation for second order accurate PIV interrogation.

Advanced PIV Processing options

Store intermediate deformed images – is active during image deformation schemes. If checked the program stores the iteratively deformed images for each intermediate processing step (after image deformation).

Sample-overlap for multi-grid processing – determines the sampling grid for the intermediate steps of multi-grid schemes.

Multi-grid – Final scan – determines the number of repeated interrogations at the final grid resolution. Generally one final scan is sufficient whereas additional scans may increase the precision.

Limit peak search area – is active only for the final passes of multi-grid algorithm at the final grid resolution. In general the iterative schemes should converge on the underlying data within a few pixels (or less) such that the correlation peak search area can be restricted.

Advanced PIV Processing Options	<u>? ×</u>
File prefix deform_	-Multigrid - Final scan
Sample-overlap for multigrid processing	✓ Limit peak search area +/- 6 px ★ (on final pass only)
×	<u>Cancel</u>

📅 PIV Evaluation Setup 🔗 🔀
Grid Correlation Algorithm Peak Detection Validation Conversion Filters etc.
Sub-Pixel Peak Fit Non-linear Gaussian peak fit Peak Fit Size 🔄 🚔
Horiz. Spx 😴 Vert. Spx 🚔 (Distance from Origin)
Maximum distance from (0,0) 4.0
Apply Cancel OK

PIV Evaluation Setup – Peak Detection

Sub-Pixel Peak Fit

Various peak fit algorithms are available which recover the sub-pixel displacement of the correlation peak:

- **3-Point Gauss Fit** This is the most frequently used scheme as proposed by Willert & Gharib (1991). It selects the four closest neighbors of a correlation maximum and fits a 3-point Gaussian curve along each of the two major axis.
- Least Squares Gauss Fit Similar to the 3-point Gauss fit but uses all 8 points surrounding the maximum. Since it uses more data for the location of the peak the noise in the result is generally less than with more simple fits. This peak fit method is recommended for most applications.
- **3-Point Parabolic Fit** Similar to the 3-point Gauss fit except that it is *only* fits a parabola through the three points.
- **Center-of-Mass Fit** This sub-pixel locator uses a specified number of correlation values surrounding the maximum (3x3, 5x5, etc.) and calculates the associated first moment or center of mass. It is generally not well suited for the narrow correlation peaks recovered from typical PIV images. In cases where larger particle images or larger features have to be correlated, the correlation peak will be wider and this fitting algorithm may yield more accurate data.
- Whittaker Reconstruction Reconstructs the correlation function from all correlation values along the horizontal and vertical line crossing the brightest correlation value (peak) and estimates the maximum from this function.
- **Nonlinear Gaussian Peak Fit** Uses up to N² points for a parametric fit to a Gaussianshaped function allowing for peak ellipticity. The peak shape is estimated from the autocorrelation function of the image samples which may reduce overall processing speed.

Limited Peak Search Area

When activated the correlation peak detection is restricted to an area whose size can be selected. For data with a small dynamic range a limited search area can significantly increase the data yield. That is, only the area enclosed by the given rectangle is scanned for the correlation peaks.

Disable Self-Correlation Peak

disables the use of the correlation peak near the origin which may be useful in images with strong (stationary) background artifacts which may hide the actual signal. Correlation peaks falling within the specified radius from the origin are disabled.

PIV Evaluation Setup - Outliers

Outlier Detection Methods

The definition of outliers in experimental data can be very subjective. In general an outlier in PIV data can be considered to be a spurious vector within an otherwise homogenous displacement field (see Sec.6.1, Data Validation, in Raffel et al., 1998 for more information).

When activated several methods for outlier detection are available. By combining the various methods most outliers are detected with very little valid data disabled (known also as optimization of cost).

PIV Evaluation Setup	<u>? ×</u>
Grid Correlation Algorithm Peak Detection	Validation Conversion Filters etc.
✓ Outlier Detection Maximum Displ. [px] Max. Displ. Difference [px] Dynamic Mean Filter Mean 1.0 ★ Var. 1.0 ★ Minimum signal-to-noise ratio 5.0 ★ Ratio Peak1 - Peak 2 Minimum correlation 28% ★ Global Histogram Filter Maximum number of validation passes 5 ★	Outlier Replacement ✓ Try lower order peaks ✓ Interpolation ✓ Re-evaluate with larger sample (multi-grid only)
× <u>*</u>	Apply <u>C</u> ancel <u>O</u> K

- Maximum Displacement Any displacement vector exceeding the specified magnitude is considered to be an outlier.
- Maximum Displacement Difference This method computes the magnitude of the vector difference between the vector in question and its eight surrounding neighbors. When this difference exceeds the specified value for more than 4 neighbors, then the analysed vector is labelled to be an outlier.
- **Normalized Median Filter** applies a median filter to each vector along with its eight neighbors. A vector is an outlier if the difference to the median displacement value is greater than the given threshold.
- **Dynamic Mean Filter** This filter works very effectively even in flows with discontinuities (e.g. shocks). However it requires a careful adjust of the coefficients to work reliably (see Raffel et al., 1998, p.153 for more information).
- **Minimum signal-to-noise ratio** disables all displacement vectors which have a SNR less than the specified value.
- **Ratio Peak 1 to Peak 2** disables all displacement vectors for which the ratio between the strongest and second strongest correlation peak is less than the specified value.
- **Minimum correlation** disables all displacement vectors for which the correlation peak value is less than the specified value.
- **Global Histogram** restricts valid vectors to be within a certain range of displacements. This rectangular or elliptical domain has to be set interactively within the PDF-window using the mouse.

Outlier Replacement Schemes

Once outliers have been detected they can be optionally replaced with new data. Two replacement schemes are available:

Other Correlation Peaks – The program stores the three highest correlation peaks found in each correlation calculation. If the highest correlation peak turns out to be an outlier this options picks the next highest correlation peak and repeats the outlier detection procedure.

Interpolation – Applies a bi-linear interpolation to estimate the replacement vector. When several of the immediate neighbors are also outliers then a Gaussian-weighted interpolation scheme is used.

Re-evaluate with larger sample – If none of the alternative correlation peaks fulfill the validation criteria, the image can be locally re-evaluated with the next largest interrogation sample in the multi-grid processing scheme.

Note: All correlation peaks are first tried out before any data is interpolated.

The **Apply Button** can be used to immediately see the effect of the outlier detection and replacement provided that the image is already evaluated.

PIV Evaluation Setup - Conversion

This dialog is used for setting the magnitude and time conversion factors which are necessary to obtain velocity data in a physical grid. Within *PIVview* the horizontal display axis is always associated with the X-axis, the vertical with the Y-axis. No Z-coordinate is used. If the measured PIV plane is not aligned with the desired physical space coordinates, the user has to perform the required transformation on the *PIVview* output using other software.

When this menu is activated, a ruler and current origin are displayed within the image. The ruler may be used to mark off know lengths in the image. The magnification factor can be interactively set by supplying the associated physical length in the menu.

📅 PIV Evaluation Setup	?×
Grid Correlation Algorithm Peak Detection	Validation Conversion Fil (
Time, Length, Units Pulse Delay 3 µs v Magnification factor in [px/mm] Length Unit mm v Displacement Unit um v Velocity Unit m/s v	Interactive ruler 336.0 px Magn. = 15.8493 mm Reset Ruler
Origin Image Origin X= 378 Y= 322 pixel Origin Offset X= 0 Y= 0 Z= 0 Image Process X Defaults Apple	mm y <u>C</u> ancel <u>D</u> K

Movement of the origin changes the entries in the menu accordingly.

Pulse Delay – specifies the time delay between the two exposures in the image pair. Generally does not apply to speckle displacment data.

Magnification Factor – The scale factor used to conversion the image displacement to a physical displacement. It can also be set interactively using the calibration menu.

Length Unit – The units used for the physical coordinate system. It is also associated with the magnification factor.

Displacement Unit – The unit associated with speckle displacement data. Also the physical particle image displacement prior to its conversion into velocity.

Velocity Unit – The unit associated with the velocity as estimated by PIV.

Origin

The origin is used as the reference point for the physical coordinate system of the data set. It can be translated within the image by holding down the left mouse button.

Origin Offset – Additional offset in physical space units if reference point is not at (0,0,0) in physical space.

Magnification Factor – This field is used for the calculation of the magnification factor by using the pixel length of the ruler within the image (numerator) and the actual distance in physical units (denominator). Only the physical distance can be modified here. The ruler length can only be changed by altering the ruler end points in the image.

The **Apply Button** can be used to immediately see the effect of the conversion factors on the data and grid.

PIV Evaluation Setup – Filters etc.

The menu tab allows the selection of various post-processing options.

PIV Evaluation Setup
Grid Correlation Algorithm Peak Detection Validation Conversion Filters etc.
Median Filtering
Median kemel size 3 x 3 💌
- 🗖 Data smoothing
Smooth kernel width (in nodes) 1 🚊
Differentiation
Calculate cell size for differentiation
C User defined differentiation cell size
Horiz, 3 Nodes 🚔 Vert, 3 Nodes 🚊
Extrapolate to edges
Apply Cancel OK

Median Filtering

The final data set is once median filtered using either a 3x3 or 5x5 kernel.

Data Smoothing

The final data set is once smoothed using a Gaussian weighted kernel of specified width.

Differentiation

Calculate cell size for differentiation The estimation of differential quantities (e.g vorticity, shear strain, etc.) relies on calculating differences between neighboring vectors and can therefore be very noisy especially if the seeding density is low and the data points are closely spaced. Generally it was found that differentials calculated from PIV data sampled at less than 50% overlap are increasingly noisy due to artifacts introduced as a consequence of nearly identical particle image pairs in neighboring sampling positions.

In other words, for correlation planes calculated from two nearly overlapping sampling areas most of the particle image pairs can be found in both samples. Especially in areas of high gradients the data will be biased according to the mean displacement of all pairs. Closely sampled data will then show a similar bias in neighboring points and thereby provide an underestimated differential.

When activated, this option checks the sampling (interrogation) window overlap and from it determines an optimal distance from which to calculate the differentials.

Alternatively, the user can specify the spacing for the calculation of the differentials.

Extrapolate to edges – Generally differentials cannot be estimated at the borders of the data set. However, by extrapolating the vectors beyond the edge using the interior vectors the differential at the edges can be estimated.

Region of Interest

This display mode can only be activated if the ROI-option is checked in the Grid Tab of the PIV-parameter menu.

The region of interest (ROI) is used limit the evaluation grid to a defined rectangular region. This region is shown as a red/white border whose corners and mid points can be translated by holding down the left mouse button.

Yellow +'s inside the ROI indicate the centers of sampling positions for the correlation analysis. Their spacing is defined in the Grid tab of the PIV configuration menu.



Image overlaid with ROI-selection rectangle. Yellow crosses mark the centers of the interrogation windows.

PivView Plotting Parameters	? ×
Vectors Contours Axis Ref.Vector	Fonts / Colors
Show Vector Field	Default
-Vector Size / Spacing	Color Coding
Stretch Factor 3	Validation
Max. Vector [px] 100	Min. 0
Spacing 1 🌻 1 🚔	Max. 9
Outlier display Mark outliers with X 💌	Color Lookup Table
Offset Vector	E Show LUT
No offset vector	Vector color
7 0 [px]	Vector style
219.415 0 m/s	
1097.07 0 um	
	K Cancel

Plot Configuration Menu

Vector Plotting Parameters

Show Vector Field – has to be enabled to show the vectors. If a data set is available, that is, the image is evaluated, then enabling this option will immediately display the vector field.

Vector Size / Spacing

Both the vector head length and width can be specified here. Also there is the choice to fill in the arrow heads.

Vector Stretch Factor – Specifies the amount the displacement is scaled with respect to the image. For example, a vector with length 10 pixels with a scale factor of 2.0 is 20 pixels long when plotted on the image. The stretch factor is always based on the pixel displacement of the data.

Maximum Vector – Can be used to avoid plotting vector whose length exceeds the value specified (e.g. outliers). The magnitude is specified in pixels.

Spacing – This option allows vectors to be skipped and is useful when the data is too dense to be properly displayed.

Outlier display

Determines how outliers should be treated on the screen:

Mark outliers with X – when activated any vector labeled as an outlier will only be indicated by an x.

Hide Outliers – when activated any vector labeled as an outlier will not be shown.

Show Outlier Vectors - when activated any vector labeled as an outlier will be shown.

Offset Vector

Subtract Image Offset – When activated the image offset specified in the evaluation setup menu will be subtracted from each displacement vector.

Subtract Mean Displacement – When activated the mean displacement of all vectors in the field will be subtracted. If disabled a user specified offset can be subtracted (see below).

Specified Offset Vector – The user can specify a given displacement to be subtracted from the vectors. This value is given in pixel displacement but is also displayed in physical units for reference. This field only become active if the Subtract Mean option (above) is disabled.

Color Coding

Magnitude – color-codes the vectors according to their magnitude. Generally a rainbow-type of LUT is used (several are available). The minimum and maximum magnitudes for the color bar can be specified in the fields below.

Signal-to-noise Ratio – color-codes the vectors according to their signal-to-noise ratio in the displacement estimation.

Validation – color-codes the vectors according their status of validation. The color coding can be set in the **Fonts/Colors** tab of the menu.

Uniform Color – alternatively the vector can be plotted in a desired uniform color; that is without color-coding. The color may be chosen through a color selection menu.

Apply Button – any change to the setup can be immediately displayed by pressing this button.

Vector Style

This sub-menu provides further choices on the vector shape, size and positioning.

Vector Style	×
Arrow Heads	_
Head Width 17%	÷
Head Length 34%	÷
Fill Arrow Heads	
- 🗖 Uniform Heads	_
Size 6	
- 🗖 Uniform vectors ————	_
Length 10	
Anchor Vector at MidPoint	
<u>A</u> pply <u>O</u> K	

Contour Plotting Parameters

Within *PIVview* a contour image describes a color-coded representation of a scalar quantity derived from the PIV data. These include velocity magnitude, vorticitiy, strain or single components of velocity.

Show Contour Image – has to be enabled to show the contour image. If a data set is available (i.e. image is evaluated) then enabling this option will immediate display the contour plot. Bi-linear interpolation is used to fill in the area between the individual data points.

PivView Plotting Parameters	? ×
Vectors Contours Axis Ref.Vector Fonts / Colors	
Show Contour Image	Default
Contour Type X-velocity [m/s]	
C Pixel Units C Physical Units (e.g. Velocity)	
Display Range: Min -2 Max 2	Nice Range
No data available	
Color Coding	
Show Color Bar	
Long Rainbow	
☑ Discrete Contours 20	
🗖 Reverse LUT 🔽 Wrap LUT 🦵 Fill Background	
Apply <u>O</u> K <u>C</u> ancel	

Contour Type

Specifies which kind of data is plotted as a contour plot. The contour display range can be adjusted below.

Vorticity

plots the out-of-plane vorticity ω_z as estimated by finite differences:

$$\omega_z = \left(\frac{\partial v}{\partial x} - \frac{\partial u}{\partial y}\right)$$

Shear Strain

plots the in-plane shear strain ε_{xy} as estimated by finite differences:

$$\varepsilon_{xy} = \left(\frac{\partial u}{\partial y} + \frac{\partial v}{\partial x}\right)$$

Normal Strain

plots the sum of the in-plane extensional strains η as estimated by finite differences:

$$\eta = \varepsilon_{xx} + \varepsilon_{yy} = \left(\frac{\partial u}{\partial x} + \frac{\partial v}{\partial y}\right)$$

Displacement Magnitude

plots the pixel displacement vector magnitude.

Velocity Magnitude

plots the velocity vector magnitude.

Physical Displacement Magnitude

plots the actual physical displacement of the particle images or speckle patterns.

Correlation Coefficient

color-codes the image according to the correlation coefficient from the displacement estimation.

Enstrophy

color-codes the squared magnitude of the out-of-plane vorticity ω_z .

Lambda operator

color-coding according to:

$$\lambda = \left(\frac{\partial u}{\partial x} + \frac{\partial v}{\partial y}\right)^2 - 4\left(\frac{\partial u}{\partial x}\frac{\partial v}{\partial y} - \frac{\partial u}{\partial y}\frac{\partial v}{\partial x}\right)$$

Vector Count

color-codes the number of validated data in averaged PIV data sets (only for averaged data).

3-C vector residuals

color-codes the residuals from a 3-C vector reconstruction for stereo PIV (only available for stereo PIV data).

(other options may also be available)

Contour Display Range

Specifies the range of contour values to be shown. These values are also displayed on the lookup-table above the contour plot.

Nice Range

calculates nice (or round) limits for the color bar using the maxima available from the current data set.

Show Color Bar

This option determines whether the reference lookup table is shown. This color bar is placed above the data field.

Discrete Contours

when enabled discrete contour levels are shown. Otherwise a smooth transition from minimum to maximum is used.

Reverse LUT

when enabled reverses LUT

Wrap LUT

starts from the lower end of the look-up table for values exceeding the specified range.

Fill Background

fills the image background with the zero-contour color.

Axis Plotting Parameters

This dialog controls the type of axis and axis labeling on the displayed data. Either image or physical coordinates may be used. The appearance of each axis may be controlled independently.

PivView Plotting Parameters			?×	
Vectors Contours	Axis Ref.Vector	Fonts / Colors		
Show Axes	,		Default	
Axis Labelling				
C Image Coordinates 💿 Physical Coordinates				
-X-Axis		-Y-Axis		
🔽 Show Labels	[pixel]	🔽 Show Labels	[pixel]	
Skip	0 🛓	Skip	0 *	
Tic Spacing	10	Tic Spacing	10	
Tic Anchor	0	Tic Anchor	0	
Minor Tics	0	Minor Tics	0 🛓	
F Show Title		🗖 Show Title		
Title 🗙		Title Y		
	Apply C)K Cance	1	
_				

Show Axes – has to be enabled to place axes on the edge of the image or plot.

Image coordinates – plots the data in image space [pixels].

- **Physical coordinates** plots the data in physical (object) space using the currently active scaling factors.
- Show Labels has to be enabled to place labels on either of the axes.
- Tic Spacing defines the increment of the tics to be labeled (major ticks).
- **Tic Anchor** sets a (virtual) starting point for the major ticks to ensure a nice layout with even label numbers.
- Skip defines which of the major ticks should not be labeled (0 for to label all major ticks, 1 to label every other, etc.).
- Minor Ticks sets the number of minor ticks between two major ticks.
- **Show Title** has to be enabled to place a title on either of the axes. The title is defined in the entry field below.

Reference Vector

/ectors Contours .	Axis Ref.Vector Fo	onts / Colors	
Show Reference Vec	tor		
Type of Ref. Vector			
Displacement	5	pixel	
C Velocity	1	m/s	
O Physical Shift	1	um	

Show Reference Vector – has to show a reference vector in the upper left corner. It is only displayed in the presence of a vector map or during interactive interrogation.

Type of Ref. Vector – three different forms of labeling the reference vector are possible:

Displacement – chooses a image based reference scale in pixels;

Velocity - labels the vector with the currently active velocity unit;

Phys. Shift – uses the physical shift of particle images.

Fonts / Colors

PivView Plotting Parameters	×
Vectors Contours Axis Ref.Vector Fonts / Colors	
Font	
Select Font Arial	
Font Sizes	
Axis Labels 13 pt 🚔 Axis Title 16 pt 🚔	
Color preferences	
Background Axis/Text Mask	L
Verbauele estre	
vector color coaing	
Valid vector Outlier vector Replaced vec. Interpolated	
Edited vector	
Apply <u>Q</u> K <u>C</u> ancel	

Font – allows the selection of the font type and size used for axis and LUT labels.

Color Preferences / **Vector color coding** – defines the color for various plot items and data types (e.g. outliers vectors).
The Preferences Menu

PivView Preferences	PIV View Preferences
General File Extensions	General File Extensions
Recent Image List contains 6 Files 🚔 Recent File List contains 6 Files 🚔	Image Extensions PIV File Extensions PIV Image Pairs ✓ Image Bairs ✓
Show Splash Screen	

- **Recent Image List** determines how many image files should be listed in the File>Recent Images menu option.
- **Recent File List** determines how many PIV data files should be listed in the File>Recent Files menu option.
- Show Splash Screen if enabled, shows the startup splash screen.

File Extensions

This menu associates a variety of image types with certain file extensions. These extensions are then used as file filters in various image selection dialogs. The drop-down selection offer the following image types:

- **PIV Images Pairs** are used for cross-correlation PIV and defines the first image of a pair of single exposed PIV images.
- **Double Exposed PIV Images** are used for the selection of images for autocorrelation PIV.
- **Speckle Reference Images** defines the zero-state (e.g. un-deformed) image used in speckle shift measurements.
- **Speckle Signal Images** associates images where the speckles have been displaced with respect to the zero-state image.
- Mask Images are images that are overlaid and disable any data covered by them.

Correlation Plane Window

The correlation plane data is color-coded using the currently active lookup table - in the example above, a gray-scale LUT. Bright pixels indicate a high correlation value, dark ones a low correlation. The local correlation values are displayed in the status line by moving the mouse across the image.

Red lines – at the center of the correlation plane marks the origin, that is, the position of zero displacement.

Green Rectangle – near the center of the correlation plane indicates the region of interest for the peak detection algorithm. It is only visible if the ROI peak search is enabled in the processing setup.

Image Samples – to the right of the correlation plane the two image samples from which the correlation was calculated are displayed if the toolbar button is enabled.

Note: In case of double-correlation, the two images represent the two correlation planes from which the plane shown on the left was calculated.



Change Image Magnification – these buttons increase or decrease the magnification of the displayed image data.

Change LUT – alters the lookup table (LUT) used to display the correlation / image data. Each time the button is pressed a new LUT is applied.



Saving the Correlation Data – the currently displayed correlation data and image samples can be stored in ASCII files by pressing this button. Once the data is saved the button is deactivated until new correlation data becomes available.



Global Histogram Window

This plot gives a visual representation of all displacement values available for the evaluated image. Each colored point in the two-dimensional plane represents a displacement vector whose coordinates correspond to the actual shift in pixels. The boundary of this plot is defined by the size of the correlation plane with zero shift at the center.

The rectangular and elliptical regions near the center are user-sizeable regions of interest. These ROIs can be used as a form of data validation: any points on their exterior are considered to be outliers. To use this function, the Global Histogram option must be enabled in the validation menu.



Global histogram of PIV data with two ROIs enabled, one rectangular (currently selected) and one elliptical. The one dimensional PDFs are shown for the X-axis below and for the Y-axis on the right.



Detail of the figure to the left. Points outside of the ROIs are labeled as outliers as indicated by their red color.

One-Dimensional Histograms – immediately bordering the 2-D histogram are vertical (right) and horizontal (below) one-dimensional histogram of the two respective displacement components. They are projections of the two-dimensional histogram onto either of the major axes. Their range therefore corresponds to the respective dimension of the correlation plane.

Sub-Pixel Histograms – these histograms, located the furthest from the 2-D histogram, are generated by plotting only the non-integral part of the displacement values. Their range is from -0.5 to 0.5 with a default bin-width of 1/20th a pixel.

Sub-pixel histograms are especially useful to check whether peak-locking effects are present in the data.



Zoom in / Zoom out – these buttons increase or decrease the magnification of the displayed histogram data.

This button enables the global histogram validation. Alternatively the Outlier Configuration menu can be used.



toggles the currently selected ROI between elliptical and rectangular shape.

- Validation with new ROI if global histogram validation is enabled, then a change to the ROI can be tested directly on the data with this button.
- Default ROI this function resets the currently selected ROI such that it nearly covers the entire histogram plane. This is especially useful if the ROI is lost off the histogram plane.

1-D PDFs and Peak Locking

The histograms (PDFs) below illustrate the effect of peak-locking which is a tendency of the peak-fitting algorithm to lock into integer values. The histogram on the top shows the presence of peak-locking in the form of integer-spaced peaks. The corresponding sub-pixel histogram below shows a triangular shape with increased counts near integer values (middle of histogram).

The histograms on the below were obtained from the same image but processed with a different peak fitting algorithm (least-squares fit). No peak-locking can be observed - the sub-pixel histogram has a near uniform distribution.

Peak-locking generally arises through a combination of certain particle image sizes (typically too small) whose narrow peaks in the correlation plane are analysed with the wrong type of peak fitting algorithm. Another cause may be the non-uniform spatial pixel response of the recording device (CCD) in combination with particle images with diameters less than 3 pixels.



One-dimensional PDFs of PIV data containing peak-locking effects.



Same PIV-data as on the left but processed with a different peak-fitting algorithm – no peak locking is noticable.

Batch PIV Processing

Batch processing is a means by which large number of PIV recordings can be evaluated in an efficient manner. A common set of analysis parameters is used for all images. During the processing the program accumulates the mean and RMS displacement and velocity data which are stored in a separate file after completion of the processing sequence.

Batch PIV Processing	?	×		
File Selection Processing Options Ou	Itput Options			
Directory e:/piv/Jet3C/Cam1				
File Type Single-exposed PIV Imag	ges(*.b16;*.ima;*.bmp;*.tif)]		
Filter #1 ×	▼ Filter #2 ×]		
Available Files	Selected Files			
min_img.b16 pivC1_0001.b16 pivC1_0003.b16 pivC1_0004.b16 pivC1_0005.b16 pivC1_0006.b16 pivC1_0008.b16 pivC1_0009.b16 pivC1_0009.b16 pivC1_0010.b16	→ dd All nove All			
Cancel	<u>Start batch</u>			

Available Files – the left column shows the files currently available in the specified directory. Multiple files can be selected by holding down the CTRL or SHIFT keys (just as in the *Windows Explorer*) Using the **Add** button the file selection is transferred to the list of files to be processed.

Selected Files – The right column shows the files currently chosen for batch processing. In analogy to the left column multiple files can be selected by holding down the CTRL or SHIFT keys. Using the Remove button the file selection is removed from the list.

Directory – this field displays the currently selected data (image) directory. It can be changed using the Directory button on the left.

File Type – this field acts as a file filter showing only file which match the chosen extensions (File type association is defined through File>Preferences>File Type.)

Filter #1 - Selects files in data directory with this (prefixed) filter. Allows selection of all odd or all even numbered files.

Filter #2 - A file filter that may be freely modified. Caution: Without the asterisc (*) no files may be displayed.

Batch Processing Options

É	Batch PIV Processing	? ×
	File Selection Processing Options Output Options	
	Perform ensemble correlation	
	Only perform image pre-processing and store resulting images	
	Load processing parameters from file	
	Choose	

Perform ensemble correlation – Rather than calculating the displacement field for each image pair, this method computes the averaged correlation planes from all images before detecting any correlation peaks. This procedure is common in the analysis of μ PIV applications and reduce the effects of underlying Brownian motion.

Load processing parameters from file – When enabled, the images will be processed using the parameters in the specified file. Otherwise the evaluation uses the currently defined parameters.

Batch Processing Output Options

A variety of output options may be defined in this menu. Preferably the data from each processing image should be stored in the platform independent netCDF format. This data can then be exported in any desired form using the same batch-processing menu.

🗟 Batch PIV Processing			
File Selection Processing Options Output Options			
Stub for generated statistics files seq			
Store data for each processed image/file			
NetCDF (Platform independent)			
🗂 Compute Reynolds stresses 🔲 Save deviations from mean			
· 🔽 Other output formats			
Tecplot (TM) ASC-II			
- 🗖 Displacement (raw) data output including			
🔽 Magnitude 🔲 SNR-data 🔲 Correlation coefficients 🗖 All correlation peaks			
· 🔽 Velocity data output including			
🔽 Magnitude 🥅 Vorticity 🥅 Normal strain 🥅 Lambda2			
🗖 Z-Coordinate 🦷 Shear strain 🦷 Enstrophy 🦳 Differentials			
Generate plot for each processed image/file Export format: BMP ✓ Put images/data into separate destination directory ✓ Choose C:/dpiv/tmp			
Cancel Start batch			

Stub for Statistics Files – specifies the beginning letters of the statistics files generated during batch processing. These data sets are computed over all available files and are stored after the batch processing is completed.

For the stub seq the output files will be:

seqAvgPiv.dat	mean pixel displacement
seqAvgVel.dat	mean velocity
seqRmsPiv.dat	RMS pixel displacement
seqRmsVel.dat	RMS velocity

Store data for each processed image/file – When enabled, this option will store the displacement and velocity data for each processed image or file. Otherwise only the mean and RMS data sets computed from all images will be stored.

Generate plot for each processed image/file – when enabled, this option will generate an image file (plot) for each processed image using the currently defined plot configuration. This is useful if animated sequences of the data are to be compiled later on. (Note: These files are quite large and may result in the use of much disk space.)

Put image / data into separate directory – this option can be enabled to store the computed data in a directory different from the image directory. This can be useful when batch-processing images retrieved from a read-only medium (e.g. CD-ROM).

Definition of statistics terms

Statistical quantities such as the mean and standard deviation of velocity are calculated by the following equations:

Quantity	Formula	Notes
Mean value	$\overline{x} = \frac{1}{N} \sum_{i=0}^{N} x_i$	here N is the number of valid measurements x_i (including interpolated values)
Variance	$\sigma_x^2 = \frac{N \sum_{i=0}^{N} x_i^2 - \left(\sum_{i=0}^{N} x_i\right)^2}{N(N-1)}$	
Standard deviation	σ_{x}	defined as square root of the variance
Mean velocity components	$(\overline{u},\overline{v},\overline{w})$	calculated from sequence of PIV data sets using formula for mean value given above
RMS-velocity	$(\sigma_u, \sigma_v, \sigma_w)$	using formula for standard deviation
Fluctuation terms	$u'u' = \frac{1}{N} \sum_{i=1}^{N} (u_i - \overline{u})^2$	Values for v'v', w'w' and v'w' are calculated correspondingly.
	$u'v' = \frac{1}{1} \sum_{i=0}^{N} (u_i - \overline{u})(v_i - \overline{v})$	Only validated and interpolated vectors are used in the calculations.
	$u'w' = \frac{1}{N} \sum_{i=0}^{N} (u_i - \overline{u})(w_i - \overline{w})$ \vdots	A multiplication with density ρ results in the respective Reynolds stresses.

PIVmask - Image Mask Editor

The program PIVmask allows the creation of image masks that can be overlaid on PIV images to prevent processing in these areas. Any PIV image can be opened using the File/Open function. Several simple shapes then are available for the creation of the mask.

Once finished the mask can be stored as a bitmap file. To enable its use by PIVview the appropriate option in the image preprocessing menu must be activated.



PIVmask in edit mode displaying several shapes during the creation of a polygon.

File Menu

- File / Load image loads one of the supported PIV image types.
- File / New resets the image mask by removing all objects.
- File / Save mask image saves the currently displayed mask as a bitmap file (BMP).

View Menu

View / Zoom in, Zoom out - changes the magnification factor of the image on the screen.



View / Darken image, Brighten image - changes the brightness of the displayed image.

View / Change LUT - changes the image display lookup table

View / Toggle pixmap size - displays larger tool buttons.

Draw Menu

5НОЦ МАБН

 \mathbb{R}

K



Draw / Invert mask - calculates the negative of the current mask.

overlaying on the image or switches back to edit mode.

Draw / Move/resize objects - allows selection of objects and as well as their size

Draw / Toggle current mask - either calculates and displays the current mask by

Draw / Move object vertices - allows the shape of polygon objects to be changed

Draw / Draw rectangles - creates a rectangle on the screen by dragging the mouse. Its shape can be later modified using the resize objects entry. Pressing the shift key allows the

Draw / Draw ellipses - creates ellipses. Circles are generated by pressing the shift key.

Draw / Draw centered ellipses - creates ellipse starting from their center point. Circles are generated using the shift key.

Draw / Draw Polygons - allows the creation of irregular polygon shapes. Points are added to the polygon with the left mouse button. The right button closes the current polygon.

Draw / Delete object - deletes the currently selected object.



PIVmask in mask display mode. Objects extending beyond the image dimensions have no influence on the mask.



PIVmask in mask display mode showing an inverted mask.

Stereo PIV Modules

Stereo PIV processing is possible through a separate set of modules for *PIVview*. These modules are tied into *PIVview3C* but can be also be used in stand-alone fashion.

Stereoscopic PIV is considerably more complex than standard 2-C PIV implementations. The use of the stereo modules is not quite as straight forward and requires that the user has a basic understanding of stereo PIV implementations - namely the concept of image registration through back-projection. Aside from reading this brief introduction new users are encouraged to study the stereo PIV example(s) supplied on the program CD.

Stereo PIV Procedure in PIVview3C

The figure at left illustrates the general procedure followed by PIVview for stereo PIV processing. The following files/information are required:

- a set of simultaneous acquired PIV images for both camera views (e.g. 1 PIV image pair per camera)
- images of a calibration grid placed parallel to the light sheet - one image per camera. The grid should also have some clearly identifiable reference marks on it to identify its orientation.
- the position of the cameras in physical space with respect to reference grid (more on this later).



The first step involves determining two sets of coefficients which allow the PIV recordings to be matched to each other. This projection of the images into a common image plane shall be referred to as image back-projection. The grid images are essential for determining the projection coefficients as the spacing between the grid lines is well defined. The modules **PIVmap3** and **PIVmapT01** (formerly named T01) use semi-automated line or dot grid detection algorithm for exactly this purpose.

Once the projection coefficients are known, PIV images for each view can be back-projected. The command-line program *mapimage* can also be used for this purpose. Note that *PIVview* does not allow you to load images as long as the mapping coefficients are unity (e.g. not calculated).

The back-projected PIV images can now be processed by standard (2-C) PIV interrogation. In *PIVview3C* the same interrogation parameters are used for both views to ensure that the resulting data sets exactly match in size and position. Command-line based PIV processing is possible by the programs *dpiv* and *dpiv3C*.

The resulting pair of 2-C PIV data sets (or files) is now merged into a single stereo (3-C) data set by taking into account the stereo viewing geometry. The geometry information can be supplied relatively easily using the **StereoGeom** module. Views of the geometrical

arrangement are available for three different orthogonal directions as a form of visual verification. Data set merging is either performed by *PIVview* internally or by the command-line program stereo.

In order to make file handling more user-friendly all information relevant for stereo PIV processing is contained in a common project file. This file is shared between the modules. Individual modules only update their specific sections. Platform-independence is guaranteed through the use of the *netCDF* data format that is also used for the data sets as well as other parameter files.

Vector reconstruction in stereo PIV

As mentioned above image pairs are first de-warped and processed for the respective view using standard 2-C PIV correlation processing algorithms. The resulting 2-C PIV data points from both views lie on a common grid but have to be combined using the local viewing angles. The figure below illustrates the viewing geometry in the xz-plane.



The system of equations given on the left comprise a once over-determined (3 unknown velocity components vs. 4 known displacements) and are solved in a least-squares fashion. The residuals of this solution are a measure of quality of the vector reconstruction. In practice the residuals should be below 0.5 pixel. The validation tab of the PIV configuration menu allows a vector validation based on these residuals.

Enabling the Stereo Mode

The stereo modules of *PIView* are enabled by selecting the menu entry *File/Stereo*. A check mark next to the option signals that *PIVview* is in stereo PIV mode. The Stereo menu is then enabled.

Once the stereo mode is enabled certain features of standard (2-C) PIVview are disabled:

- The calibration mode is no longer available. All calibration is now supplied through the *PIVmap* image registration module, that is, the physical grid spacing together with the mapping coefficients define the image magnification factor. The origin is also defined through *PIVmap*.
- PIV images can only be loaded if the mapping coefficients are available. This prevents processing of images that have not been registered with respect to each other.

Stereo Menu Options & Toolbar

8-0

New project - allows a step-by-step creation of a new stereo project using the StereoWizard program including the image reconstruction (back-projection) procedure.

Open project - opens an existing stereo project. It loads previously calculated projection coefficients, the viewing geometry, the calibration images used by PIVmap as well as previously loaded PIV images (if available).

Save project - saves the current stereo configuration (coefficents, geometry & parameters) in the current file name.

Save as - allows the stereo configuration to be stored under a different name.

Recent projects - allows the selection of one of the (4) most recently accessed stereo projects.

Geometry - opens the stereo geometry definition window (module *StereoGeom*). It is here where the camera position with respect to the viewing area (light sheet) is specified.

Image registration - opens the *PIVmap* module which is used for determining the backprojection coefficients.

A Show Image A - displays the image as viewed from camera A (if loaded). If the data has been processed, the image will be overlaid with the corresponding (2-C) vector or contour map(s).

B Show Image B - displays the image as viewed from camera B (if loaded).

Show combined data - displays the reconstructed 3-C data if available.

3C

Creating a new stereo project with PIVmap3

The menu entry *Stereo/Image Registration* calls up the program *PIVmap3* which allows a step-by-step calibration of a stereo PIV configuration including the grid registration procedure.

Background on Stereo PIV

Most planar imaging techniques require some sort of calibration to associate image coordinates with physical coordinates. A viewing arrangement in which the object plane (e.g. laser light-sheet) is parallel to the image plane (co-planar arrangement) involves a simple linear transformation consisting of translation, scaling and possibly rotation. This is no longer the case for oblique viewing arrangements or if lens distortion is present. In effect the co-planar (on-axis) viewing arrangement should be considered as a special case of the oblique viewing arrangement. The calibration procedure for the more generalized oblique viewing arrangement typically requires high order mapping functions to transform image pixels to physical space coordinates.

Another calibration step is necessary for multi-component measurements such as DGV or stereoscopic PIV. In this case the positions of the cameras with respect to the image plane are required. The most pragmatic approach is a direct measurement of the camera positions with respect to the object plane (Willert, 1997) which is not necessarily trivial in complex geometries involving windows or change of media (i.e. air-water interface).

Soloff et al. (1997) introduced a mathematical formalism that combines back-projection and reconstruction in a single step. The procedure requires multiple images of the calibration target as it is translated at known displacements along the Z-axis (typical 5-7 recordings). The drawback of this approach is that the camera position is not retrieved directly and therefore unsuited for DGV. In effect the resulting projection equations connect object volume coordinates with planar image coordinates thus simplifying the vector reconstruction in stereo-PIV.

A third approach to camera calibration originates from the field of photogrammetry and image vision. Common here is the use of a pin-hole camera model with added parameters to account for radial distortions as suggested by Tsai (1986). The calibration procedure involves a so-called bundle adjustment of object-space calibrated feature points using a least squares method. In practice this can be achieved by translation of the target as in the previous case. More general approaches can even calibrate on sequences of randomly oriented target recordings. The major drawback of this method is that the camera parameters – so called internal parameters – such as pixel pitch need to be specified. Specification of these intrinsic parameters is especially difficult when imaging fiber bundles are used. Also this model does not seem to account for distortions introduced by tilting of the back-plane (Scheimpflug angle adjustment in stereoscopic PIV). In addition to the estimated camera position the Tsai camera model also yields distortion coefficients and estimated lens focal length which are not necessarily of direct interest in image coordinate calibration.

A slightly different approach to camera calibration is described in the following and has been implemented in the *PIVmap3* program. The basic aim is to provide an accurate image mapping to de-warp images on the hand as well as to retrieve the position for each individual camera without knowledge of the camera characteristics as in the pinhole camera model (i.e. pixel pitch, lens focal length). The image mapping function is calculated from a single recording of a calibration target that is placed in the object plane. The target should contain a known pattern such as a rectangular grid of lines, dots or crosses as shown in figure 4.

Calibration of the camera position is also based on an image series of the Z-translated targets as proposed by Soloff et al. (1997). The approach followed here is to use a functional relationship between the out-of-plane displacement and associated in-plane displacements to estimate the camera position. In practice, this involves solving the projection equations for

the point at which the image coordinates (x_i, y_i) vanish. With the object coordinates (x_o, y_o, z_o) and image coordinates (x_i, y_i) the first order projection equations are:

$$x_{i} = \frac{a_{00}x_{o} + a_{01}y_{o} + a_{02}z_{o} + a_{03}}{a_{30}x_{o} + a_{31}y_{o} + a_{32}z_{o} + a_{33}}$$
$$y_{i} = \frac{a_{10}x_{o} + a_{11}y_{o} + a_{12}z_{o} + a_{13}}{a_{30}x_{o} + a_{31}y_{o} + a_{32}z_{o} + a_{33}}$$
with $a_{33} \equiv 1$

Although these equations do not include higher order distortions, they nonetheless are sufficient for most PIV applications and show deviations on the order of 1/10th of degree with respect to the more precise pinhole camera model by Tsai (1986). These higher order methods are also included in *PIVmap3*.



Figure 1: Coordinate system and viewing arrangements. Left: co-planar viewing arrangement with object plane parallel to image plane. Right: non co-planar (off-axis) arrangement with image plane tilted with respect to object plane.

a) Grid imaged without distortions on co-planar arrangement b) Grid imaged obliquely in a non co-planar arrangement

c) Grid imaged with strong lens distortion

5	_	
•••••		
•••••		
		•••••
•••••		•••••••••••

Figure 2: Examples of imaged grid under different viewing arrangements.



Figure 3: Multiple plane calibration procedure based on translation of calibration target. Camera A and B comprise a single-sided stereo configuration. Camera A and C view the light sheet from two sides. Here the light sheet thickness should be accounted for.



Figure 4: Images of various calibration targets. In all cases the camera is located to the right of the reference marks.

The program PIVmap3

This program is intended to provide accurate image-object calibration of oblique (distorted) camera views as well as estimate the position of the camera if a series of Z-translated images is available. The core of the program is an algorithm to reconstruct a regular grid from images of a calibration target. The following gives a brief overview of the reconstruction algorithm before the programs features are described in detail.

Grid reconstruction in PIVmap

The reconstruction of the grid itself is a three step process:

Step 1 – Image enhancement

A variety of filters can be used to increase contrast of the image if the following two steps fail to produce reliable results.

Step 2 – Feature extraction

This is based on correlation-based pattern matching, that is, a grid of dots is correlated with a template containing a single dot. Similarly a template containing a cross is used to detect a grid of crosses or line grid. As the correlation locks into features brighter than the

background, images may need to be inverted a priori. This is done automatically by specifying the *type of target* in the input menu (e.g. black dots on white background). In any case, the image shown in the calibration menu should contain bright features (dots, crosses or lines) on a dark background. (See sections below for further details.)

Step 3 – Grid reconstruction

This step attempts to assign physical coordinates to the detected feature points. The grid origin is assigned to the point closest to the origin of the manually placed helper coordinate axis (red and yellow arrows). The helper axes are then used to find the 4 neighboring points on the major axes as well as 4 neighboring diagonal grid points. The resulting group of 3x3 points forms the basis for the reconstruction of the remainder of the grid using a predictor-corrector approach. To allow for deviations of the grid nodes from their predicted positions a *reconstruction tolerance* is defined (given in percent).

The reconstructed grid will have grid index coordinates assigned to the valid correspondence points. Points on the major axes are shown in purple while unmatched points are shown in orange.



Reconstructed grid showing labeled grid points. Unmatched points are shown in orange.



Reconstructed grid with X-axis pointing down, while Y-axis is pointing to left.

General strategy on feature extraction

The image sample size and the sampling overlap determine the detection efficiency of the grid. As a rule of thumb, the sample size should be near the node distance in the image. In general one should start with a larger sampling window. If not all grid points were found change to next smaller sample or increase sample overlap. If grid points were detected more than once increase the merging radius to combine neighboring points. Depending on the tolerance the grid reconstruction procedure may ignore these misaligned calibration points.

Effect of sample size on dot detection

The following images illustrate the effect on sample size on the feature extraction. Detected dots are labeled with a cross.



→ some dots detected twice (resolved by increase of merging radius)

Calculation of the mapping function Once the grid is reconstructed the image mapping function may be calculated. The mapping

is a combination of two separate steps as shown in the figure below



The first mapping serves to *de-warp* the image to reconstruct the image space onto a Cartesian grid. At this point the image axes are parallel to the axes in physical space. The second step merely applies a linear transformation to convert from image space to physical space.

Mapping from image recording to reconstructed image

A variety of mapping functions are available for the first mapping step:

polynomials up to degree 4

$$\begin{aligned} x_{m} &= c_{x0} + c_{x1}x_{i} + c_{x2}y_{i} \\ &+ c_{x3}x_{i}^{2} + c_{x4}x_{i}y_{i} + c_{x5}y_{i}^{2} \\ &+ c_{x6}x_{i}^{3} + c_{x7}x_{i}^{2}y_{i} + c_{x8}x_{i}y_{i}^{2} + c_{x9}y_{i}^{2} \\ &+ c_{x10}x_{i}^{4} + c_{x11}x_{i}^{3}y_{i} + c_{x12}x_{i}^{2}y_{i}^{2} + c_{x13}x_{i}y_{i}^{3} + c_{x14}y_{i}^{4} \end{aligned}$$

$$\begin{aligned} y_{m} &= c_{y0} + c_{y1}x_{i} + c_{y2}y_{i} \\ &+ c_{y3}x_{i}^{2} + c_{y4}x_{i}y_{i} + c_{y5}y_{i}^{2} \\ &+ c_{y6}x_{i}^{3} + c_{y7}x_{i}^{2}y_{i} + c_{y8}x_{i}y_{i}^{2} + c_{y9}y_{i}^{2} \\ &+ c_{y10}x_{i}^{4} + c_{y11}x_{i}^{3}y_{i} + c_{y12}x_{i}^{2}y_{i}^{2} + c_{y13}x_{i}y_{i}^{3} + c_{y14}y_{i}^{4} \end{aligned}$$

> projection equations, a ratio of polynomials $(1^{st} \div 1^{st}, 1^{st} \div 2^{nd}, 2^{nd} \div 1^{st}, 2^{nd} \div 2^{nd})$

$$x_{m} = \frac{a_{00}x_{i} + a_{01}y_{i} + a_{03} + a_{04}x_{i}^{2} + a_{05}x_{i}y_{i} + a_{06}y_{i}^{2}}{a_{30}x_{i} + a_{31}y_{i} + a_{33} + a_{34}x_{i}^{2} + a_{35}x_{i}y_{i} + a_{36}y_{i}^{2}}$$
$$y_{m} = \frac{a_{10}x_{i} + a_{11}y_{i} + a_{13} + a_{14}x_{i}^{2} + a_{15}x_{i}y_{i} + a_{16}y_{i}^{2}}{a_{30}x_{i} + a_{31}y_{i} + a_{33} + a_{34}x_{i}^{2} + a_{35}x_{i}y_{i} + a_{36}y_{i}^{2}}$$
with $a_{33} \equiv 1$

The choice of function depends on the degree of distortion in the image. In general, if the grid consists of straight lines a ratio of first order polynomials should be favored. However convergence of the polynomial ratios is not always guaranteed. Table 1 gives an overview of the associated number of coefficients and the required number of calibration points. If the number of calibration points exceeds the minimum count a standard least squares fit is used to calculate the coefficients for the polynomial mapping functions; a Levenberg-Marquardt nonlinear least squares fit is used in the calculation of the polynomial ratios.

Mapping function	Number of coefficients	Minimum number of calibration points
Polynomial, 1 st Order	6	3
Polynomial, 2 nd Order	12	6
Polynomial, 3 rd Order	20	10
Polynomial, 4 th Order	30	15
Ratio, 1 st ÷ 1 st	8	4
Ratio, 1 st ÷ 2 nd	11	6
Ratio, 2 nd ÷ 1 st	14	7
Ratio, 2 nd ÷ 2 nd	17	9

Table 1.: Mapping functions and their coefficients.

Mapping from reconstructed image to object space

Mapping between the reconstructed image and physical space is described by the following linear transformation:

$$x_{Obj} = \frac{x_m}{M} + x_{Offset}$$

where the magnification *M* is given in [pixel/mm]
$$y_{Obj} = \frac{y_m}{M} + y_{Offset}$$

PIVmap3 Input Menu

The main menu – shown Figure 5 – consists of two input areas: The left area is used to specify parameters that are common to all camera views and/or all image planes. The right area is specific to each supplied calibration image which is referred to here as *plane*. In the *Global Settings* area the following parameters are specified:

Configuration data file

stores a summary of all parameter settings within the current instance of the program. The load function can retrieve previous parameter settings. The parameter files are ASCII formatted and may be edited with any text editor.

Type of target

specifies whether the images of the calibration target contain dots, crosses and lines. Further the image characteristics with respect to features (black-on-white or white-on-black) are specified here. This may be specific for each individual view.

File Mapping View Tools Help		
🗅 🖄 🖬 🕅 🏞 🕰 🔽 View A 💌 Plane 0 🚔 📕	F 🎬 🖿 🛲 💈 🛛 🔍 100% 🕑	3 q. † † %, H H
Global Settings Corriguration data file In3C/Steeo/Pt/map3Demo/PpeFlow/3planes.map Target / Views Target / Views Type of target Xtransledd images Stranslend images Grid Spacing on Target X(nm) [2] Y Common grid for all views Output Image Sce Wridth [pol] 500] Height [pol [500]	Calibration Images Dir [Image file I c12:05.TIF 2 c12:05.TIF 3 c1:05.TIF Use mark image Invert mesk.	Zpos (mm) Ref 0.5 0 Ref 0.5 0 S 0 Ref 0.5 0 Ref 0.5 0.5 0.5 0.5 0.5 0.5 0.5 0.5

Figure 5: Main input menu of program *PIVmap*

Number of views

allows up to 4 camera views. Stereoscopic PIV would require only 2 views. A triple camera DGV-setup would require 3 views. A single camera view could be used to calibrate a standard 2-C PIV experiment.

Z-translated images

determines the number of calibration images for each view. At least two are required if the camera position is to be determined.

Grid Spacing on target

determines the physical distance of the calibration points on the target.

Output Image Size

specifies the image size of the reconstructed (de-warped) image. This parameter is common to all camera views.

Calibration input area

Load (Images)

permits the simultaneous selection of all calibration images associated with this camera view. The common directory of these files is specified above the table.

Table Area

The table itself lists the image files, their position in Z and the specified reference plane. Clicking on the image file shows a preview of the corresponding image in the area below the table. The out-of-plane translation (Z-position) of each image can be entered directly in the corresponding cell. The reference plane indicates that the camera position and image mapping is calculated with respect to the given image and is set by clicking on the corresponding cell. Single images can be loaded by clicking the icon on the far right of the table.

Mask image

can be used to mask out areas in the calibration images to improve grid reconstruction performance and also increase overall calculation speed (less correlations required). Mask inversion disables the complementary area. A common mask is used for all images of a single view.

PIVmap3 Calibration Area

Once a sequence of images has been loaded they may be displayed one-by-one in the calibration menu. The display should display the image with bright features (dots, crosses) on a dark background. If not, then the image type was not properly set in the main menu (i.e. black-on-white vs. white-on-black). Also any optional filtering that has been applied to the image will be made visible here. The image will be overlaid with a pair of coordinate axis, herein referred to as *Helper Axis*. These may be adjusted with the mouse and serve to define the principle directions in the grid. Details on this are given in the following section.

The calibration area is accompanied by a tab-menu containing the tabs *Filters*, *Scan*, *Scale* and *Map*. Their purpose is described in the following.



Figure 6: Calibration view of program PIVmap3 with Filters Menu shown.

Filters Tab

To improve the feature detection several filters may applied to the image. They are shown top-down in the order of their application.

Median Filter

applies median filtering on the image. Either a 3x3 or 5x5 kernel may be used.

Histogram Equalization

performs histogram equalization on the image.

High-Pass Filter

applies a Gaussian-weighted high-pass filter to the image. The width of the Gaussian filter function can be specified. Values between 3 and 10 are recommended. The processing speed decreases with increased width of the weight function.

Binarize

converts the image to contain only black (0) or white (255) pixels by setting values above the specified threshold to 255. Lower values are set to 0.

Low-Pass Filter

applies a Gaussian weighted smoothing filter to the image. In combination with the high-pass filter the resulting image will be band-pass filtered. Recommended kernel widths are the range of 0.5 - 1.5 pixels.

Scan Tab

Feature extraction in *PIVmap3* is based on a correlation approach in which a template of either a dot or cross is matched locally to the image. The position of best match indicates the

center of the dot or cross in the locally sampled area. The performance of this matching process depends on the choice of sample size, sampling overlap and a few more parameters that are set in this tab.

Spot/Template radius

specifies the radial dimension of the dot or cross in the template.

Sample size

provides a choice of image sampling areas. Typically this size should roughly correspond to the mean spacing of the features to be extracted.

Sample overlap

specifies the image sampling distance. A higher sampling rate increases the overall performance but decreases processing speed.

Minimum count

As each sample produces an estimate of a feature point location, adjoining samples will detect the same feature point more than once. This value specifies the minimum number of times a feature point needs to be detected to qualify as such. The merging radius given below has a direct influence.

Merging radius

specifies the scan area to combine near-by feature point estimates. Values between 3 and 10 are recommended.

Reconstruction tolerance

Vmap - I	Contro			×
filters 2	Scan	Scale	Map	
Template	radius (px]	4	Ň
		Spot cove	ers 16.7% o	of sample
Scan				
Sample	size (p)	<] 48x48		•
Sample	overlap	80%		•
Minimu	m count		4	A V
Merging	g radius	[px]	7	Ť
Grid Reco	onstruct	ion ——		
Reconstruction Tolerance 10%				
Arrow tips mark neighboring points				
- 🔲 Limit	range c	if nodes -		
×	10	*	Y: 10	*
🗖 Ellip	otic regio	on of intere	est	
E Editio		into monu	-11.	
	iove po	inis manua i	any	
📿 Sca	an			
GRID Buil	ld Grid			

P

is required to reconstruct the grid from the detected feature points. The reconstruction is based on a predictor-corrector approach where valid grid points are used to estimate the position of neighboring grid points. The search area is determined by this tolerance value. Values between 5 and 10 percent are recommended.

Arrow tips mark neighboring points

is used to aid the reconstruction algorithm in finding the first points near the origin.

Limit range of nodes

restricts the range of nodes in the reconstructed grid. This is useful if the image has many invalid feature points on the edges. The use of a mask (specified in the main menu) should however be favored. Circular or elliptic calibration geometries can be treated as well.

Edit/move points manually

In some cases the grid reconstruction may fail, if falsely detected or missing feature points near the grid origin mislead the reconstruction. By manually moving these points the grid may be properly reconstructed nonetheless. However, these movements are lost if the scan operation is repeated.

Scan

performs the feature extraction scan by the localized template matching scheme using the current parameter settings. Any manual change to previously detected feature points are lost.

Build grid

tries to reconstruct the grid using the currently detected feature points.

Scale Tab

Dot Spacing

specifies the actual spacing of the grid points. This value is set in the main input menu of the program and cannot be changed here.

Global physical offset

assigns the coordinates of the calibration origin in physical space. Initially this value is zero and should be set according to the experiment's reference frame.

Camera position

reflects the camera's position with respect to the origin of the grid in the specified reference plane (set in main menu). These values are irrespective of the global physical offset.

Estimate

tries to calculate the camera position using at least two calibrated planes. A valid image mapping (see Map) is required. The result, even for non-convergence is placed in the edit boxes above.

Various estimation schemes are available, some of which require additional input regarding the type of camera (sensor) with which the images were acquired.

PIVmap - Control 🛛
Filters Scan Scale Map
Dot Spacing [mm]
Spacing X 2
Spacing Y 2
optioning i ja
Global physical offset
×[mm] 0
Y [mm] 0
Z [mm] 0
Lamera position
×[mm] 218.68
Y [mm] -0.378627
Z [mm] 220.842
Underlying camera model
Tsai Model, fully optimized 📃
Camera pixel size [µm]
9 9
Kodak ES1 0 / KAL1010
E stimate

Map Tab

Parameters specified in this menu determine image mapping coefficients that perform the back-projection from camera view to an un-distorted image. The choice of mapping function is specified on the top of this menu.

Magnification

determines the final magnification factor for the reconstructed image. This value is constant throughout the image.

Translate image

By default the origin of the grid will be shifted to the lowerleft corner. These parameters allow the origin to be shifted.

Copy to all views

copies the current parameters to the remaining camera views. This is especially of interest when all reconstructed views should coincide after mapping them.

Fit to output image

is useful to approximately match the reconstructed grid to the specified output image size and should be used on a new set of calibration images. An even better match is achieved by manually fine tuning the magnification and shift parameters.

Dewarp this view

calculates the mapping coefficients for each reconstructed grid of this camera view. These mappings are used to back-project the original images which are then displayed. The overall speed of the mapping process depends on the output image size and the number of reconstructed planes.

PI¥map ·	- Control				×
Filters	Scan	Scale	Мар		
Project	ion Map (I	Order 2)	-		
[px/mm]	J3U	J			
Translat	e image –				-
	Shift X [m	m] 21		1	
	Shift Y [m	m] 21		÷	
(1) C	opy to all v	/iews			
Fi	t to output	t image			
D	ewarp this	: view			

Helper axis

Reconstruction of the grid is based on the a priori knowledge of the direction of the major axes as well as an approximate location of the origin. This is achieved by manually specifying a pair of *helper axes* for each calibration image. The figures below show how these axes may be placed in the image using the mouse. The tips of each axis should be placed on the points neighboring the origin to allow a more efficient grid reconstruction.

The function 'Copy helper axis' can be used to copy the coordinates from the previous plane to the currently displayed plane. This speeds up positioning over many calibration images.



Position of axes at beginning

Dragging the X-axis to the origin.

Alignment of Y-axis and arrow tips.

Figure 6: Manual positioning of the *helper axis* in the image.

Description of the PIVmap3 menu entries

🗋 File/New

creates a new PIVmap3 project. All previously loaded images are removed and file names are reset.

File/Open...

loads a previously created PIVmap3 project. All current settings are lost. File/Save

stores the current selection of image and parameters in a user-specified project file. This ASCII file is editable with any text editor (i.e. Notepad). The default extension is 'map'.

File/Save

stores the current selection of image and parameters in a user-specified project file. This ASCII file is editable with any text editor (i.e. Notepad). The default extension is 'map'.

File/Exit

leaves the program after asking the user for confirmation.

Menu View



View/Zoom out/in – changes the image display scale factor



View/Brighter - increases the image brightness on the screen (Short-cut key: B).



View/Darker – decreases the image brightness on the screen (Short-cut key: D).



d2

View/Auto-scale - sets default intensity on image

View/Change color LUT – skips to the next available lookup table (LUT) to convert the grayscale images to pseudo-color.

-

*

Menu Mapping

Mapping/Input menu

changes to the main input area of the program

CAL Mapping/Calibration mode

changes to the calibration area and displays the filtered image of the currently selected view and plane.

Mapping/Scan image

performs the feature extraction scan by the localized template matching scheme using the current parameter settings. Any manual change to previously detected feature points are lost.

Mapping/Create grid

tries to reconstruct the grid using the currently detected feature points.

Mapping/Show mapped image

calculates the mapping coefficients for each reconstructed grid of this camera view. These mappings are used to back-project the original images which are then displayed. The overall speed of the mapping process depends on the output image size and the number reconstructed planes.

Mapping/Process current view

performs grid reconstruction and image back-projection for all images in the currently selected camera view.

Mapping/Overlay mapped images

computes a combined image of all available back-projected views for the currently selected image plane. At least two valid (back-projected) camera views are required. This is especially useful to check the coincidence of the back-projections.

Mapping/Process all images

processes all images for all views including camera position estimation

Mapping/Center helper axis

This utility re-centers the helper axis to the middle of the image.

Mapping/Copy helper axis

This utility copies the helper axis from the previous plane to the current plane.

Mapping/Estimate camera position

tries to calculate the camera position using at least two calibrated planes of the current view. A valid image mapping is required. The result, even for non-convergence is placed in the calibration tab menu.

Mapping/Save calibration points

exports the detected and validated feature points for all available plane of the current view to a user-specified ASCII file.

Mapping/Export for PIVview

stores the current image projection coefficients along with all relevant data (i.e. image size, camera positions, etc.) in a user-specified *PIVview3C* project file. If an existing file is selected the new data is merged into it.

Menu Tools

Tools/Batch image dewarping

opens a menu allowing a series of images to be back-projected using the



current image maps. The image interpolation function and destination directory may be specified.

Tools/Show log

opens a window with past error messages, warnings and comments.

Creating a new stereo project with StereoWizard

The menu entry *Stereo/New project* calls up the program *StereoWizard* which allows a stepby-step definition of a stereo projection including the grid registration procedure.

Step 1 - Calibration image

Here the user is asked to either load a pair of calibration grid images or an existing project he wishes to modify.

Step 2 - Calib.Image - Basic Information

Using the dropdown entries give the actual (physical) spacing of the line in the grid (here 2 mm). Later the actual magnification factor is calculated from this.

Also, some information regarding the image type needs to be given. For proper grid detection it is important for the program to know whether the image contains dark lines on a bright background or bright lines on a dark background.

Step 3 - Calibration Images – Orientation

In case the image needs to be flipped in some way (mirror view), this can be achieved here.

🖀 Create new Stereo-PIV Project 🔹 👔 🗴					
Calibration Images					
A new stereo PIV project requires a pair of calibration images. These should basically consist of grid lines of known spacing similar to the image on the right. Proper processing requires white lines on a dark background. Dark lines on a white background are also fine if specified as such.					
Pick File A C:/DPIV/StereoJet/gridC1_0001.b16					
Pick File B C:/DPIV/StereoJet/gridC2_0001.b16					
Import existing project Show Map(s)					
Help < Back Cancel					
🖀 Create new Stereo-PIV Project 🔹 👔 🗴					
Calibration Images - Basic information					
Grid Spacing Horizontal 2 Vertical 2 Vertical 10 Mits mm					
Image splitting is required for images of the B16 file format which may contain a pair of images. Further, proper calibration image processing requires bright (white) grid lines on a dark background. The inverse - black lines on white - need to be specified as such.					
Image Information B16 image treatment Image contains Image A No image split Image A No image split Image B No image split Image B Black Lines Image A Preview B					
Help < <u>B</u> ack <u>N</u> ext > <u>C</u> ancel					
*** Create new Stereo-PIV Project Calibration					
In order to properly match up the grids the images must have the same orgination					
Below you can flip either of the images about the horizontal and/or vertical axis.					
Flip L<=>R Flip T<=>B Flip L<=>R Flip T<=>B					
Show full scale					
Help < Back Next > Cancel					

Step 4 - Set Reference Point

The buttons *Display A* and *Display B* show the currently loaded grid image. Using the mouse the position of the origin can be selected. The grid detection algorithm later snaps to the nearest grid crossing.

Step 5 - Camera Positions

This dialog allows the user to enter the positions of the camera with respect to the origin. In addition the physical (offset) coordinates of the reference point may be specified. Alternatively the *StereoGeom* program can be used to enter the coordinates.

The furthe	er nrocessin	a will lock	onto the	e neares	t arid cr	nssin	1	
1110 10001				- nearer	, gind of			
	100	T						
	-	-	-	-	-		_ [Display A
	100	1000		830		21		
			1	D				
		L	1	H	26			
		199.9		29				Display E
			-					
Help	1		(Back	1	N	evt s	1	Cano
Help reate ne amera p The stere with respe should pre The coord	ew Stereo ositions oscopic ca ect to a kno eferably be dinate conv	-PIV Pro	< <u>B</u> ack ject equires t nce poin the grid ir such that	hat the p t in the o mage pa	No position calibrati ir. plane c	of the	camer I. This	anc
Help Treate ne amera p The stere with respe should pre The coord sheet. Th below to i -Camera F	w Stereo ositions oscopic ca cot to a kno eferably be dinate conv e Z-axis is n theractively Positions	-PIV Pro own refere visible in vention is normal to v position i	< <u>Back</u> ject equires the nce poin the grid in such that the light s the came	hat the p t in the o mage pa t the XY sheet. U gras.	No cosition calibrati ir. plane c se the s	of the on grid	camer I. This les with	anc
Help The stere with respe should pre The coorr sheet. Th below to i -Camera F	ew Stereo ositions oscopic ca set to a kno eferably be dinate conv e Z-axis is is nteractively ^p ositions — >	-PIV Pro	< Back	hat the p t in the o mage pa t the XY sheet. U gras.	osition alibrati ir. plane c se the s	of the on gric	camer I. This les witi Z	Canc
Help Treate ne omera p The stere with respe should pre The coord sheet. Th below to i -Camera A	ew Stereo ositions oscopic ca sect to a kno eferably be dinate conv e Z-axis is n nteractively nteractively ositions > ositions	-PIV Pro	< Back	hat the p t in the o mage pa t the XY sheet. U gras.	Dosition calibrati ir. plane o se the s	of the on gric tereo	came I. This config Z	anc as is known reference In the light uration tool
Help Treate ne amera p The stere with respe should pre The coord sheet. Th below to i -Camera A Camera B	ew Stereo ositions oscopic ca cet to a kno eferably be dinate conve eferably be dinate conve eferable be dinate conve effabricatione dinate conve effabricatione di conve effabricat	-PIV Pro	< Back	hat the p t in the o mage pa t the XY sheet. U ras. Y	<u>N</u> osition alibrati ir. plane c se the s	ext > of the on gric tereo 400	camel I. This les with config	as is known reference In the light uration tool

Step 6 - Save Project

Finally the project must be given a name and is saved at the chosen location.

To begin with the actual grid detection, start the *PIVmap* module by pressing the button. The next steps show how the grid is recovered and the projection coefficients are determined.

Once the wizard is closed *PIVview* asks the user whether to load the new project file into PIVview3C.

Parameters such as the camera positions may be modified at a later stage from within *PIVview3C*.

Create new Stereo-PI	¥ Project	<u>: ×</u>				
Save Project						
Finally store the previous settings in a configuration file. All information including the projection coefficients will be stored in this file. This file allows the stereo project to be reloaded later in time.						
Project Name jet_flow	v1	Directory				
File Name c:/dpiv/	jet_flow1.cfg					
	Save Project					
The next steps to obt	ain stereo PIV data are :					
1 - Run the PIV-MAP gr	id registration program	Show Map(sì				
2 - load a pair of stereo Pl (a valid projection is requir	V images 📕 ed!!)					
Help	< <u>B</u> ack <u>F</u> inisł	h <u>C</u> ancel				

alul

StereoGeom - Geometry Definition

This stereo module is used to define the stereo viewing arrangement, namely the camera positions with respect to the light sheet. The entered coordinates are specified with respect to the same same reference point used by the *PIVmap* module. Please note that in current implementation of *PIVview* the z=0 plane always defines the light sheet plane, that is, the X and Y-axis lie in the light sheet, while the Z-axis points normal to this plane.

The camera positions can be entered directly in the fields on the right hand side. The positions within the geometry plot are updated immediately and give the user a visual feedback.

 $\downarrow z \neq z \neq \downarrow$ These tool buttons allow the user to view the stereo geometry from three major (orthogonal) directions.

The Apply button is used to transfer the settings to *PIVview* such that the reconstruction of the 3-C velocity map is repeated using the new settings.

StereoGeom may also be used in stand-alone fashion. In this case a file menu allows the user to load or save existing stereo project files.



StereoCfg module showing XZ-projection of the camera / light-sheet geometry. In this example the cameras are placed on either side of the light sheet.

PIVmapT01 - Projection Module

This module is a semi-automated procedure for determining the image back-projection coefficients that are necessary for stereo PIV processing. For input it requires images of a calibration grid (see notes below) for each camera view. The grid reconstruction process can be split up into three major steps - the coefficient calculation tab dialog follows these steps.

- **Step1**: Line detection is performed by interrogating image locally using the Hough transform which transforms the image intensity distribution into a polar coordinate space. Since straight lines show up as distinct bright spots in this transformation line segments are detected for each interrogation area.
- Step 2: The segments line are subsequently assembled into sets of continuous lines spanning the entire image. Separation into two major directions as well as sorting of the lines are then performed to reconstruct a grid (meshing). In order to properly overlay the images, the origin can be shifted in grid intervals using a desired point of reference in the calibration grid.
- **Step 3**: At this point the projection coefficients can be calculated using a desired magnification factor for the back-projected image.



For more information about this module please consult the paper by Ehrenfried (2001).

Some notes on obtaining calibration grid images

The grid must be made up of two sets of continuous, evenly spaced, parallel lines that are orthogonal to each other. The grid should either have bright lines on a dark background or dark (black) lines on a bright (white) background. The latter lends it self especially well to be printed by a laser printer. The grid should be aligned parallel with the light sheet of cover the entire field of view (if possible). The grid image should contain 10-30 lines in each dimension (roughly 30-60 pixel spacing) - this results in 100 to 1000 grid crossings.

Misalignment of the grid with respect to the light sheet may increase the measurement error. Reference marks placed on the grid allow doubtless identification of image orientation and also help in obtaining an image-object space referencing.

Important: Make sure to check the invert option in case your image has dark lines on a bright background. Otherwise *PIVmapT01* will not be able to properly detect the grid lines.

PIVmapT01 Toolbars & Menu Options

- File/Load image opens a file dialog for loading new calibration images.
- File/Load project opens and loads an existing project.
 - File/Save project save current project.
- Tools/Place origin allows rough placement of the origin using the cursor. The grid detection set the origin to the line crossing closest to this point.
- Tools/Complete evaluation performs full reconstruction on loaded image(s) using current parameter settings.
- Tools/Line detection performs line detection on current image using the current Hough transform processing parameters.
 - Tools/Form mesh tries to reconstruct the mesh from the detected set of lines.
 - Tools/Reconstruct grid detects the line crossing in the mesh and determines the major axis.
- Tools/Find optimal map provides an initial guess for the back-projected image by trying to keep the magnification factors near unity. Fine tuning is possible with by changing the parameters in the Projection tab of the parameters.
- Tools/Eval. parameters opens a tab dialog that allows step-by-step processing of the grid images.

Tools/Show Info - gives information about the status of PIVmap.

- Tools/Measure places a ruler into the image allows the user to measure distances between two points.
- **ROI** Tools/Set ROI allows interactive definition of a region of interest in which grid detection is to be performed.
- Tools/Elliptic ROI toggles the ROI shape between rectangular and elliptic.

Tools/Best fit ROI - maximizes the ROI to fit the image size.

- View/Show Grid A displays the calibration image for camera view A.
- **B** View/Show Grid B displays the calibration image for camera view B.
- AB View/Back-projected image(s) dewarps the loaded calibration images using the current sets of projection coefficients (if avaiable).
- The View menu has a number of entries that are identical to those found in PIVview: zoom functionality, image brightness, etc.

쿣



Main window of PIVmapT01 showing calibration image overlaid with reconstructed grid. Red crossings indicate the major axis. The origin itself is marked by the yellow circle. Note that the image is a negative of the original which consists of dark lines on a bright background.

PIVmapT01 Tab Dialog

The tab dialog is accessed with the button an upper portion giving feedback on the current status of *PIVmapT01* such as which image is available and displayed or whether grid reconstruction has been successful. The lower tab dialog (described below) follows the main processing steps of *PIVmapT01*.

The buttons **A** and **B** toggle the display between images A and B. The parameters in the first two tabs are unique to each image a change accordingly. The Projection-tab parameters are identical for both images.

The **Go-Button** starts a complete processing of the loaded images using the current parameter settings. If successful the mapped image pair is displayed in the main window.

PIVmapT01 Dialog - Detect Lines

This tab is devoted to the line detection portion of T01 based on the Hough transform. The following parameters can be set:

- Pixel threshold threshold for pixels to be used in the Hough transform. High value = less lines, but faster. Low value = more lines, slower.
- Line threshold threshold for detection of the maximum in the Hough plane. High value = less lines, low value = more lines.
- Filter size kernel size of filter for smoothing the Hough transform (Can normally be set to 1, e.g. no filtering)
- Scan size theta angular increment used during maximum search. Default of 5 usually does not need to be altered.
- Scan size radius radius increment used during maximum search. Default of 5 usually does not need to be altered.

			•	
TU1 - C:/DP14/	stereotest/	GridA.ti	l Man	<u> </u>
Image L	aded Lines	mesn 4	мар	
A	$\checkmark \checkmark$	\checkmark	-	
В	√ √	\checkmark	\checkmark	
Proj				Go
Detect Lines	Mesh and Map	o Proj	iection	1
-Hough			1 🗖	(Nogativo
Pixel threshold:	80%	<u>*</u>		negauve 1
Line threshold:	50%	<u>+</u>	-	Detect lines
Filter size:	1	-		
Scan size theta:	5	÷		
Scan size radius	5	<u>+</u>		
-Scanning			J	
Window size:	129	-		
	20			
Window step:]:	db	T		
Г	Image desa	ampling		
Done				

- **Window size** size of the sampling window used in the Hough transform. Its size should be larger than the typically distance between two grid lines. This window size also defines the length of the individual line segments. In case of strongly curved grid lines the window size should not be chosen to be too large which might impede the accuracy.
- **Window step** sampling distance between Hough interrogation windows. This distance should be smaller than the window size to ensure overlapping line segments. Small value = slower, but high accuracy High value = faster, less precise.
- **Image desampling** down-samples the image to speed up processing at the price of slightly reduced precision especially useful during optimization of line detection parameters.
- **Invert** turns loaded image into its negative. The Hough line detection algorithm requires bright line data on a dark background.
- **Detect lines button** starts the detection algorithm. This process is the most lengthy step in the grid reconstruction process but can be aborted at any point, available line segments are displayed.

PIVmapT01 Dialog – Mesh & Map

This dialog contains all parameters necessary for reconstruction of the grid from the previously detected line segments.

- **Form mesh** button connects line segments to form continuous grid lines, The result is displayed by overlaying the image with green lines.
- **Mapping** button calculates the back-projection (mapping) coefficents.
- Scan size image search distance between neighboring line segments. Default of 5 usually does not need to be altered.
- **Theta tolerance** maximum possible angle (degrees) allowed between two neighboring line segments. For strongly curved grid lines this value may need to be increased (default = 5°).
- Minimum count minimum number of line segments necessary to form a grid line;

high value = less short lines, less chance of outlier lines,

low value = permits more short lines. In case of outlier lines just increase this value.

Base angle - defines the principle direction of the Xaxis. An angle of 0° corresponds to having the Xaxis run horizontal from left to right. In case of diagonal lines choose a value of either +45° or -45°.

T01 - C:/DPIV/stereotest/GridA.tif	<u>? ×</u>
Image Loaded Lines Mesh Maj	p
	÷
B ✓ √ √ √	,
Proj	Go
Detect Lines Mesh and Map Projection	m
Mesh forming	
Scan size image: 5 🚔	Form mesh
Theta tolerance: 5	Mapping
Minimum count: 5	
Base angle: 0 🚔	
- Mapping	
Iterations: 10 🚔 🗔 Flip L-R	
Ist Order ☐ Flip T-B	
C Mix: 2nd / 0th	
C Mix: 2nd / 1st	
C 2nd Order	
Done	

- **Horiz. Shift** shifts the vertical (Y) axis horizontally across the image. Defines the horizontal coordinate of the image origin.
- **Vert. Shift** shift the horizontal (X) axis vertically across the image starting from the bottom. Defines the vertical coordinate of the origin.
- Flip L-R flips the image left-to-right. Useful in case of using mirrors or viewing the light sheet from two different sides.

Flip T-B - flip the image top-to-bottom.

- **Iterations** maximum number of iterations used for the calculation of the coefficients; low value = faster, higher chance of result not being converged, high value = slower, better chance of convergence.
- 1st Order sets mapping function to be a ratio of polynomials of first order.
- Mix: 2nd/0th mapping function only a second order polynomial (divisor is unity).
- **Mix: 2nd/1th** mapping function = ratio of 2nd order polynomial and first order polynomial (divisor).
- 2nd Order mapping function = ratio of two second order polynomials.

Note: higher order mapping functions may result in more precision but may require more iteration steps.

PIVmapT01 Dialog – Projection

This dialog is used to set the image size of the back-projected images.

- **Calculate** button calculates back-projected images which are displayed in the main window of PIVmap. The images a overlaid with each other.
- **Only A** calculates back-projected image only for view A.
- **Only B** calculates back-projected image only for view B.
- **Origin ox/oy** x-, y-coordinates (pixel) of the image origin within the back-projected image.
- **Grid spacing x, y [pixel]** distance between grid lines in the back-projected image. This value should chosen to be close to the line spacing in the original images.
- **Grid spacing x, y [mm]** distance between grid lines in physical space. The ratio of this value and the pixel grid spacing gives the image magnification factor (also shown in dialog).
- Image Size x,y pixel dimensions of the resulting image which is shown as a red box in the main window.

🚻 T01 - C:/DPI¥,	/stereo	test/(GridA.t	if	<u>? x</u>
Image	Loaded	Lines	Mesh	Мар	
A	\checkmark	\checkmark	\checkmark	\checkmark	÷
В	\checkmark	\checkmark	\checkmark	\checkmark	
Proj	A+B				Go
Detect Lines	Mesh ar	nd Map	Pro	jectior	
–Output Image–					
Origin	ox:	290		÷	Calculate
	oy:	203		- -	Only A
Grid spacing	x (px):	40		ㅋ	Only B
	y [px]:	40		Ð	Best quess
Grid spacing	x [mm]:	5		-	
(physical)	y (mm):	5		-1	
Magn. [px/mm]	Mx=8	My=8			
Image Size	Width:	570		÷	
	Height:	480		÷	
Rotate	deg:	0/100)	÷	
Done					

Stereo View Disparity Correction

In practice it generally is not possible to perfectly align the calibration target with the light sheet. As a result the sampling windows of the PIV interrogation do not coincide between the two views. One way to account for this is to perform a PIV correlation analysis on the particle images acquired at the same instant from the two views, that is, correlating an image from view A with an image from view B. This correlation results in a vector map showing the misalignment between the two camera views. This disparity (misalignment) map can then be used to adjust the projection coefficients.

PIVview3C offers this alignment procedure in the Stereo Menu entry *Disparity Correction*¹, which opens the window shown on the right (provided a stereo project and image set has been loaded).

The alignment window shows an overlayed image of the two camera views. The mouse can be used to interactively correlation-probe the image, similar to the main *PIVview* window.



Main window of the *PIVview3C* stereo camera disparity correction.

Displacement vector map showing the misalignment between the camera views.

Disparity Correction Tool Buttons:

■

The functions of the view tool buttons are equivalent to those of other views so only those unique to this window are described here.

shows the correlation plane window for interactive image probing.

opens the configuration menu shown to the right of image view.

performs a full interrogation of the image pair followed by data validation. The result is displayed as a green vector map.

applies a least squares fit to the validated data and displays a red vector map of the fit.

saves the coefficients of the current least squares fit and prompts the user whether to enable the fit in which case *PIVview* reloads the image set with modified projection coefficients.

Once the alignment map has been activated this window can be used to verify the actual alignment. Further alignment iterations are not possible.

¹ In earlier versions of *PIVview* the misalignment correction between the camera views was referred to as *residual alignment*.

Some notes on performance:

In some instances the particle images viewed by one camera are only weakly distinguishable on the other. This is frequently the case for air flows seeded with small particles utilizing a forward-backward scattering arrangement. For this reason rather larger interrogation windows are necessary. When in doubt, we suggest increasing the interrogation window until consistent displacement fields are produced. The correlation plane window is also useful for this optimization.

The correlation module of alignment procedure actually makes use of all four images of a stereo PIV recording and uses a multiplication of correlation planes to enhance the correlation peak.

The alignment procedure allows a choice on which of the two projection maps are modified; the default is a symmetric adjustment of both views.



Disparity correction window after least-squares fit to offset data.
Command-Line Utility Programs

The *PIVview* program suite includes several command-line based utility programs that can be used for batch processing purposes, debugging, data retrieval, etc. These programs should be run from the command shell (e.g. DOS-box), and their path (.\PIVview\Bin\ or .\PIVview\Utils\) should be added to the system's PATH variable.

These program are useful for batch processing large amounts of PIV-images without using *PIVview* itself (frequently done on Linux/Unix systems).

AVERAGE

This program first builds a list of files from the first and last file provided and then calculates statistics from these. The result is placed in a new file whose name is determined from the common file stub of the input files. It is useful in the calculation of quantities such as mean velocity field or Reynolds stresses.

Calling procedure:

average <dataIN_1> <dataIN_N>

Parameter description

<dataIN_1> (input) first file in sequences to be processed.

<dataIN_N> (input) last file in sequences to be processed.

Command line options

-prefix<stub> overrides the output file stub to the name provided

-rms	also stores the RMS data in a separate file
-reyn	computes Reynolds stresses from velocity data
-fluct	computes fluctuations by removing mean from each data set
-info	provides feedback during processing of data
-tec	store result in Tecplot ASCII file
-plt	store result in <i>Tecplot</i> binary file.

DPIV

The program *dpiv* reads the specified image(s) and parameter file and then performs complete PIV processing including:

- > image pre-processing (filtering, background subtraction, etc.)
- > correlation analysis (particle image shift extraction)
- data validation
- > conversion to velocity or physical shift

The result is stored in a *netCDF* file such that it can be read by other *PIVview* programs. Use ncdf2tec to convert the data to a *Tecplot*-readable format. Unless specified with the -o option the output file name is determined from the image name.

Calling procedure:

dpiv <par_file> <imageA> <imageB>

Parameter description

<par_file> (input) PIVview generated parameter or data file (netCDF-format) with
processing information.

<imagea></imagea>	(input) PIV or speckle image
<imageb></imageb>	(input) Second image of pair (if required)

Command line options

-def	performs second-order accurate image deformation PIV
-mg	performs multigrid (pyramid) PIV
-auto	perform autocorrelation processing; only one image is required.
-noflip	avoids flipping or rotating the image (important for stereo-PIV processing)
-noip	avoids image pre-processing options
-norotate	avoids flipping or rotating the image (important for stereo-PIV processing)
-info	provides feedback during processing of data
-o <fname></fname>	stores output data in file <fname></fname>
-pco	treats the input image as a PCO-double exposure image
-speckle	perform speckle shift processing, where the first image is considered the reference (unshifted) image.
-wx y	specifies the interrogation window size. Overrides the value specified in the parameter file.
-sx y	specifies the interrogation window step size. Overrides the value specified in the parameter file.
-bx y	specifies the global image offset. Overrides the value specified in
-mask <fname></fname>	uses specified mask image. Useful for batch processing with varying masks. Note: mask image must have the same size as PIV images.

DPIV3C

The program *dpiv3c* reads the specified stereo configuration and parameter files to perform complete stereo PIV processing using the given image set. This processing including:

- image pre-processing (filtering, background subtraction, etc.)
- image back-projection using projection maps (including disparity correction)
- correlation analysis for both view (particle image shift extraction)
- data validation
- merging of the two views
- conversion to velocity

The result is stored in a netCDF file such that it can be read by other PIVview programs. Use ncdf2tec to convert the data to a Tecplot-readable format. Unless specified with the -o option the output file name is determined from the first image name (View A).

Calling procedure:

```
dpiv3c <cfg_file> <imageA1> <imageA2> <imageB1> <imageB2>
dpiv3c -pco <cfg_file> <imageA> <imageB>
```

Parameter description

<cfg_file></cfg_file>	(input) <i>PIVview3C</i> generated stereo configuration file (netCDF-format).
<imagea?></imagea?>	(input) PIV image pair for view A
<imageb?></imageb?>	(input) PIV image pair for view B

-align	perform camera disparity correction
-def	performs second-order accurate image deformation PIV
-mg <fact></fact>	performs multi-grid (pyramid) PIV, initial de-sampling factor <fact> (default=3)</fact>
-info	provides feedback during processing of data
-o <fname></fname>	stores output data in file <fname></fname>
-noip	avoids image pre-processing options
-рсо	treats the input images as a PCO-double exposure image (both images must be PCO images)
-wx y	specifies the interrogation window size. Overrides the value specified in the parameter file.
-sx y	specifies the interrogation window step size. Overrides the value specified in the parameter file.
-bx y	specifies the global image offset. Overrides the value specified in
-mask <fname></fname>	uses specified mask image. Useful for batch processing with varying masks. Note: mask image must have the same size as PIV images.

IMAGEPRE

This utility program reads the specified image and performs image pre-processing (filtering, background subtraction, etc.) as specified in the parameter file.

Calling procedure:

```
imagepre <par_file> <imageIN> <imageOUT>
```

Parameter description

<imageout></imageout>	(output) pre-processed output image
<imagein></imagein>	(input) original image
<par_file></par_file>	(input) PIVview generated parameter or data file (netCDF-format) with image pre-processing information.

Command line options

-noflip	avoids flipping or rotating the image (important for stereo-PIV processing)
-norotate	avoids flipping or rotating the image (important for stereo-PIV processing)
-info	provides feedback during processing of data
-pco	treats the input image as a PCO-double exposure image and generates two output images with the extensions .ima, .imb
-pco2	same as the -pco option, but stores the resulting image pair in a single PCO- double expsure image file.

NCDF2TEC

The program reads the specified *PIVview* data file and generates a *Tecplot*-formatted ASCII or binary output file.

Calling procedure:

ncdf2tec <fileIN> <fileOUT>

Parameter description

<filein></filein>	(input) netCDF data file as generated by PIVview

<fileOUT> (output) Tecplot formatted output file

-asc	creates ASC-II data file
-bin	creates Tecplot binary data file
-vel	writes out the velocity data contained in the input file if it is available. This is the default output option.
-disp	writes out the pixel (image) displacement data contained in the input file.
-piv	writes out the pixel (image) displacement data contained in the input file.
-mag	also writes out either velocity magnitude or pixels shift magnitude depending on which option was set.
-raw	skip post-processing of data after loading (i.e. validation)
-diff	stores differentials (implies -vel)
-dvort	stores vorticity data (implies -vel)
-dshear	stores shear strain data (implies -vel)
-dnstrain	stores normal strain data (implies -vel)

-s <num></num>	sets output stride to <num>, that is only every <num>th point is stored.</num></num>
-sx <num></num>	sets horizontal output stride to <num>.</num>
-sy <num></num>	sets vertical output stride to <num>.</num>
- Z	stores Z-coordinate
-z <num></num>	stores user-specified Z-coordinate
-info	provides feedback during processing of data

NCDUMP

The NCDUMP program reads the specified data file and writes the contents onto the command (shell) screen.

Calling procedure:

ncdump <fileIN>

Parameter Description

<fileIN> (input) any netCDF-formatted data file

- -h (option) shows only variable and attribute definitions and their values.
- -help (option) shows all available options for ncdump.

NCGEN

The NCGEN program reads the specified raw CDF data file and generates a netCDF (binary) data file. The makes it possible to convert the ASC-II output obtained from *ncdump* back to *netCDF*. The resulting data file is written to the screen unless a file name is specified with the -o option.

Calling procedure:

ncgen <fileIN>

Parameter Description

- <fileIN> (input) ASC-II formatted data file with variable and attribute definitions
- -o <fileOUT> (option) stores the resultant netCDF data in the specified output file. Otherwise data is written to stdout (screen).
- -help (option) shows all available options for ncdump.

POSTPROC

This program first builds a list of files from the first and last file provided and then performs post-processing (i.e. validation, conversion, filtering, etc.) on each PIV displacement data set. The results are placed in new files whose names are determined from the common file stub of the input files. It is useful when only post-processing of previously calculated PIV displacement data is necessary.

Calling procedure:

postproc <PARfile> <dataIN_1> <dataIN_N>

Parameter description

<PARfile> (input) first file in sequences to be processed.
<dataIN_1> (input) first file in sequences to be processed.
<dataIN_N> (input) last file in sequences to be processed.

-prefix<stub> overrides the output file stub to the name provided

-A	use parameters associated with camera view A, if available
-B	use parameters associated with camera view B, if available
-info	provides feedback during processing of data
-tec	store result in Tecplot ASCII file
-plt	store result in Tecplot binary file.

STEREO

The program reads the specified *PIVview* data files supplied for each camera view and reconstructs the 3-C velocity field using the geometry information specified in the parameter file. Use the *StereoGeom* module to specify the geometry information.

Calling procedure:

stereo <CfgFile> <fileA> <fileB> <fileOut>

then this data is interpolated)

Parameter Description

<CfgFile> (input) Stereo project file as generated by *PIVview3C* and/or stereo modules.

<fileA> (input) PIV data set for view A (netCDF format required)

<fileB> (input) PIV data set for view B (netCDF format required)

<fileOut> (output) PIV data with reconstructed 3-C velocity data (default = *netCDF* format).

-info	(option) gives feedback during processing.
-plt	(option) writes out the data in Tecplot-binary format.
-tec	(option) writes out the data in Tecplot-ASCII format.
-A	use grid from <filea> as output grid (If the grid from view B does not coincide then this data is interpolated)</filea>
-B	use grid from <fileb> as output grid (If the grid from view A does not coincide</fileb>

Shortcut keys in PIVview

- F2 region-of-Interest selection (if enabled)
- F4 PIV parameter menu (also with right-mouse button on image)
- F5 start PIV evaluation
- F6 show correlation plane
- F7 show 2-D scatter plot (histogram)
- F10 redraw current view
- B image brighter
- D image darker
- I show / hide image after evaluation
- L change pseudo-color look-up table
- Shift+L show previous pseudo-color look-up table
- P open plotting parameter menu
- C show / hide scalar image after evaluation
- V show / hide vector field after evaluation
- + increase image magnification (zoom in)
- decrease image magnification (zoom out)
- CTRL+i image pre-processing menu
- CTRL+L hide / show LUT bar
- CTRL+M open mask editor
- CTRL+O open file dialog for loading images
- CTRL+P print currently displayed data or image
- CTRL+Q quit PIVview
- space-bar toggle between images of image pair
- CTRL+Space overlay images of image pair

Mouse ShortCuts & Actions

Left click on image pair

- single-point PIV evaluation using current parameter settings.

- Right click on image pair
 - opens context-sensitive menu
- Double-click on evaluated image
 - opens a data monitoring window showing statistics for each point currently captured by the mouse.

Supported File Formats

Supported Image Formats

Tagged image file format (TIFF)

at 8-bit/pixel and 1-bit/pixel including compression (Not supported are LZWcompressed TIFF images on which Unisys holds a patent. If you plan to use LZW-TIFF images use a image converter such as IrfanView to convert the images to another supported format). Color images are converted to 8bit grayscale.

Windows Bitmap (BMP)

1, 4 and 8 -bit/pixel, color images are converted to 8-bit grayscale.

JPEG Color JPEG's are converted to 8-bit grayscale.

Portable Network Graphics (PNG)

Color images are converted to 8-bit grayscale.

Portable Bitmap Graphics (PBM, PGM, PPM)

This includes various formats from a collection of ASCII based bitmap types: PBM, PGM and PPM. PBM (Portable Bitmap format) is a 1-bit, black and white bitmap type, PGM (Portable Greymap format) is an 8-bit greyscale bitmap type and PPM (Portable Pixelmap format) is a 24-bit full colour bitmap.

PCO-B16 a 16-bit/pixel image format used by PCO's camera acquisition software

B32 an extension of PCO's B16-format to 32-bit/pixel

LaVision Images (IMG)

a 16-bit/pixel integer and 32-bit/pixel floating point image format used by LaVision's DAVIS software. Not supported is LaVision's proprietary, compressed IMX file format.

Note: The command line utilities such as *dpiv* only support the following image formats:

- > TIFF all except LZW-TIFF and 24-bit color
- > BMP 1 and 8 bit/pixel, uncompressed
- ➢ PCO-B16
- LaVision-IMG

Image naming and numbering conventions

Generally PIVview automatically determines the image name of the second image using the name of the first. Several rules are used to determine the second image:

$ima \rightarrow imb$	If the first image has the extension 'ima' then the program will look for a file where the extension changed to 'imb' Example: (1 st image: myImage.ima) → (2 nd image: myImage.imb)
a. → *b.*	If the first image has the letter 'a' before the extension then program will look for a file changed where 'a' is replaced by 'b'

- Example: $(1^{st} \text{ image: myImage}_a.tif) \rightarrow (2^{nd} \text{ image: myImage}_b.tif)$
- (n) \rightarrow (n+1) If the images are sequentially numbered the program will look for the file with the next highest index. Example: (1st image: img099.tif) \rightarrow (2nd image: img100.tif)

Excepted from these rules are images loaded in the speckle or autocorrelation modes, as well as images that are loaded as self-contained image pairs (e.g. PCO-B16 double images).

Supported Output File Formats

- **ASCII** Plain ASCII data can be generated by PIVview which can then be read by other programs without much effort. A short header specifies the number of points in the data set along with a description of each column's content.
- TecplotPIVview data can be exported to be compatible with the TecPlot graphics
package (AMTEC Engineering, Inc., www.amtec.com), Version 7.0 and higher).
Style files for both displacement and velocity data are generated to ease the
data presentation with TecPlot.
- **netCDF** is a self-describing, platform-independent data file format frequently used in the scientific community. NetCDF data generated by PIVview is aimed to be compatible with the layout utilized by PivNet/EuroPIV cooperation projects.

Since NetCDF is an open format, additional variables can be added without interfering with a predefined layout. Further information on NetCDF can be found at Unidata (www.unidata.ucar.edu).

Two separate utilities are redistributed with PivView that allow NetCDF file viewing (ncdump) as well as NetCDF file generation (ncgen). In principle these utilities allow the user to edit a NetCDF file by first converting it to a text file with ncdump and recompiling it after modification using ncgen.

Parameter Files

The current parameter settings are generally stored in the system registry when the program exits. Once the program restarts these parameters are restored. However, the most important parameters, especially those regarding evaluation, conversion and post-processing can also be stored in a separate file which has the same name as the image and the extension .par. These files are formatted in NetCDF (see above) and are compatible with the existing UNIX-based PIV software of the German Aerospace Center (DLR) as well as with the specifications defined by the EuroPIV / PivNet activities.

References and Further Reading

- Raffel, M.; Willert, C. & Kompenhans, J., 1998: Particle Image Velocimetry A Practical Guide, Springer Verlag, Berlin, ISBN 3-540-63683-8.
- Astarita T, Cardone G (2005) Analysis of interpolation schemes for image deformation methods in PIV" Exp Fluids 38:233-243 Hart, D.P., 2000: PIV error correction, Exp Fluids **29**: 13-22.
- Brent RP (1973) Algorithms for minimization without derivatives Prentice Hall, Englewood Cliffs NJ, Ch.5
- Coudert SJM, Schon J-P (2001) Back-projection algorithm with misalignment corrections for 2D3C stereoscopic PIV. Meas Sci Technol **12**:1371–81.
- Meinhart C, Wereley S, Saniago J (2000) A PIV algorithm for estimating time-averaged velocity fields. J.Fluids Eng.-T. ASME **122**(2):285-289
- Pérenne N, Foucaut JM, Savatier (2004) Study of the accuracy of different stereoscopic reconstruction algorithms. In: Stanislas M, Westerweel J, Kompenhans J (eds) Particle image velocimetry: Recent improvements. Proceeding of the EUROPIV 2 workshop, Zaragoza, Spain, March/April 2003. Springer, Berlin Heidelberg NewYork, pp. 375-389.
- Press WH, Flannery BP, Teukolsky SA, Vetterling WT (1988) *Numerical recipes in C the art of scientific computing* Cambridge University Press, Cambridge, Ch.10.2, pp. 299.
- Tsai RY (1986) An efficient and accurate camera calibration technique for 3D machine vision. Proc. IEEE Conference on Computer Vision and Pattern Recognition (CVPR'86), Miami, FL, pp. 364–374.

- Wereley ST, Meinhart CD (2001) Second-order accurate particle image velocimetry. Exp Fluids **31**:258-268
- Wernet MP (2005) Symmetric phase only filtering: a new paradigm for DPIV data processing. Meas Sci Technol **15**:2001-2028
- Westerweel J, Scarano F (2005) Universal outlier detection for PIV data. Exp.Fluids (*Online First* publication).
- Westerweel J, Dabiri D, Gharib M (1997) The effect of a discrete window offset on the accuracy of cross-correlation analysis of digital particle image velocimetry. Exp.Fluids **23**: 20-28.
- Wieneke B (2005) Stereo-PIV using self-calibration on particle images, Exp Fluids **39**, 267-280
- Willert, C., 1997, Stereoscopic digital particle image velocimetry for applications in wind tunnel flows, Meas Sci Technol 8: 1465-1479

Willert, C. & Gharib, M., 1991: Digital particle image velocimetry, Exp Fluids 10: 181-193.