

Caddo 10B
Digital Electrocardiograph

Learning Material
Ver 1.1

An ISO 9001:2008 company



Innovative Technology Ecosystem

94, Electronic Complex, Pardeshipura

Indore - 452 010 India

Tel: 91-731 4211100

Fax: 91-731-2555643

e mail: info@scientechnbz

Websites: www.caddo.bz

www.scientechnbz

www.scientechnbz

Certificate

Standard: **ISO 9001:2008**

Certificate Registr. No. **85 100 001 10182**

TÜV Rheinland India Pvt. Ltd.;

Certificate Holder: **Sciencetech Technologies Pvt. Ltd.**

Unit 1: 94 – 101, Electronics Complex, Pardeshi Pura,
Indore – 452 010, Madhya Pradesh, India

Unit 2: 90 – 91, Electronics Complex, Pardeshi Pura,
Indore – 452 010, Madhya Pradesh, India

Scope: **Design, Manufacture of Electronic Test & Measuring
Instruments, Training Products for Electrical & Electronics
Education and Providing Technology Training**

An audit was performed, Report No. 10182. Proof has been
furnished that the requirements according to ISO 9001:2008
are fulfilled.

The due date for all future audits is 04-10 (dd.mm).

Validity: The certificate is valid from **2010-12-13** until **2013-12-12**.



Bangalore, 2010-12-20

The Certification Body of
TÜV Rheinland India Pvt. Ltd.

The validity of this certificate is subject to timely completion of surveillance audits as agreed

in the Contract. The validity of the Certificate can be verified under WWW.tuv.com with the identification No. 8510000110182

ISO/IEC 17021

www.tuv.com

 **TÜVRheinland®**
Precisely Right.

**Caddo 10B
Digital Electrocardiograph
Table of Contents**

1.	Security Notice	4
2.	Safety Instructions	5
3.	Features	6
4.	Technical Specifications	7
5.	Maintenance Regulation	8
6.	Apparatus Characteristic	9
7.	ECG Caddo 10B Panel Sketch Map	10
8.	Operation Regulation	12
9.	Preparation Before Operation	13
10.	Attention during Operation	14
11.	Recording Paper Loading	15
12.	Electrode Installation	16
13.	Grounding And Power Connection	19
14.	Battery Operation Regulation	16
15.	Keypad And Controls	21
16.	Troubleshooting	24
17.	Maintenance Transportation And Preservation	27
18.	Warranty	28
19.	List of Accessories	28

Security Notice

1. Make sure the instrument grounding properly during installation.
2. If the ground cable is not integrated, please run the device with battery.
3. Please pull out power supply plug before change the fuse.
4. This device must be operated and preserved by professional doctor.
5. The operator must read this user manual carefully before operation and operate the device according to operation regulation strictly.
6. The design of this device with mature consideration of security, but operator should never neglect attention to device state and patient's situation.
7. Please dismantle the battery and pull out power supply plug before cleanout and disinfection of this device.
8. Please don't operate this device in the environment which contains flammable anaesthesia gas.
9. If use this device with cardiac defibrillator or other electric stimulate devices at same time, please use our company's Ag - AgCl chest electrode and ECG lead, if use the electric stimulate device over 55 seconds, please choose one-off chest electrode. We suggest ECG (Caddo 10B) not be used with other electric stimulate device, if it is compulsory, there should be professional technician guided on the scene.
10. When other devices are connected with this ECG instrument, they must be Type I devices which accord with IEC6060 1-1. Because the total amount of leakage current may hurt patients, the monitoring of leakage current is carried out and taken charge by connect devices.

Safety Instructions

Read the following safety instructions carefully before operating the instrument. To avoid any personal injury or damage to the instrument or any product connected to it.

Do not operate the instrument if you suspect any damage within.

The instrument should be serviced by qualified personnel only.

For your safety :

Use proper Mains cord : Use only the mains cord designed for this instrument. Ensure that the mains cord is suitable for your country.

Ground the Instrument : This instrument is grounded through the protective earth conductor of the mains cord. To avoid electric shock the grounding conductor must be connected to the earth ground. Before making connections to the input terminals, ensure that the instrument is properly grounded.

Observe Terminal Ratings : To avoid fire or shock hazards, observe all ratings and marks on the instrument.

Use only the proper Fuse : Use the fuse type and rating specified for this instrument.

Use in proper Atmosphere : Please refer to operating conditions given in the manual.

- 1. Do not operate in wet / damp conditions.**
- 2. Do not operate in an explosive atmosphere.**
- 3. Keep the product dust free, clean and dry.**

Features

- One-touch Operation
- High Resolution Thermal Printing System
- Free from adjustment in operation as well as high frequency response and elimination of artificial errors which are common in stylus recording system
- Capable of printing continuous I channel trace and annotations including lead mark and such parameters as Sensitivity, Paper speed, and Filter operation status
- High Resolution Digital Filter
- Inhibits baseline drift without distorting ECG waveform
- Able to enhance efficiency and relieve operating strain
- Concept of floating input circuit meets IEC safety standard
- Visual status due to light indicated touch keys
- Rechargeable battery supporting more than 50 patient examinations
- Under the automatic mode, pressing key "MODE" continuously can keep recording when abnormal lead is found
- Modular design guarantees upgrade to 3-channel ECG and communication with computer database and network
- The Waveform is displaying synchronously on the LCD with CCFL 240 x 128 dots

Technical Specifications**Normal work environment****Operation**

- (a) **Environment temperature** : +5°C ~ +35°C
- (b) **Relative humidity** : $\leq 80\%$
- (c) **Power supply** : AC 220V, 50Hz (110V, 60Hz)
: DC 7.4V, 3700 mAh rechargeable lithium battery
- d) **Atmosphere pressure** : 86kpa ~ 106kPa

Store and Transportation

- a) **Environment temperature** : 10°C ~ 55°C
- b) **Relative humidity** : $\leq 95\%$
- c) **Atmosphere pressure** : 50kPa ~ 106kPa

Input Circuit : Floating and Defibrillation Protection

Lead : Standard 12 Leads

Patient Current Leakage : $<10\text{mA}$

Input Impedance : $\geq 50\text{m}\Omega$

Frequency Response : 0.05 Hz ~ 150 Hz (-3dB)

Time Constant : >3.2 Seconds

CMRR : $>60\text{ dB}$ $>100\text{dB}$ (add filter)

EMG interface Filter : 35 Hz (-3dB)

Recording way : Thermal Printing system

Specification of recording paper

- 50mm (W) \times 20m (L) high speed thermal paper (Roll)

Paper Speed : 25 mm/s, 50mm/s, error +5%

Sensitivity choice : 5, 10, 20mm/mv, error + 5% Standard sensitivity is 10mm/mv + 0.2 mm/mv

Classification : Class I, CF applied part

Noisy Level : $+15\mu\text{Vp-p}$

Safety Standard : IEC class I, type CF

Skin Voltage Tolerance : $\pm 300\text{mV}$

Fuse : 2 pcs $\phi 5 \times 20\text{m}$ T250mA /250V
(Power supply 220V) AC time lag

Size (mm) : 315 (L) \times 215 (W) \times 77 (H)

Weight : 2.25Kg

Maintenance Regulation

1. Under the condition of normal use according to the user manual and operation notice, if this instrument has any problem, please contact with our customer service department. Our company has the sales record and customer archives for each instrument. The customer has two year's warranty service from the beginning of shipping date according to the below time and condition. To supply all-around and fast maintenance service to our customers, please mail the maintenance card to us in time.
2. Our company may adopt the ways of instruction, mailing to company by courier, visiting customers' company, etc to carry out the maintenance promise.
3. Even in the period of free maintenance, we charge for reparation in the following archives:
4. Faults or damnification caused by misuse because not operate according to user manual and operation notice.
5. Faults or damnification caused by dropping accidentally when users move after purchasing.
6. Faults or damnification caused by preparation, reconstruction, decomposition, etc outside of our company.
7. Faults or damnification caused by natural disasters such as fire, flood, earthquake, etc.
8. Faults or damnification caused by unapt thermal recording paper.
9. The free maintenance period for spare parts and fray parts is half a year. Power cable, recording paper, operation manual and packing material are excluded
10. Our company is not responsible for the faults of other connecting instruments caused by the faults of this device directly or indirectly.
11. The free maintenance service will be cancelled if we find the protection label has been destroyed.
12. For charge maintenance beyond the warranty period, our company advise to continue to use "Maintenance contract regulation". Please consult our customer service department for specific situation.

Apparatus Characteristic

1. Recording system: Thermal-array (8 dots/mm), it needs not be adjusted. Frequency Response: 150Hz (IEC)
2. The device can record exact single ECG waveform and remark. The remark includes: lead sign, sensitivity, paper speed, filter state.
3. Under automatic mode, just press the button once, it starts record procedure, which can enhance your work efficiency.
4. The keyboard is convenient to operate, and the LCD can display the operation state, which is convenient and readable.
5. Classification: Class: I, CF applied part
6. The device can use AC and DC and it includes built-in chargeable lithium battery.
7. This instrument can record 150 pieces of ECG waveform and print 90 minutes continually under the best DC state.
8. The figure of whole device is elegancy and gliding.
9. According to dependence degree of deleterious fluid, this device is belong to common device.
10. The device can't be used in the environment, which contain flammable anaesthesia gas mixed with Air
11. Adopting digital signal which deals with the work filter, the baseline filter and the EMG filter will obtain the higher quality of the ECG.
12. The device can AUTO print the normal ECG, which can lighten the doctor's burden and enhance your work efficiency
13. According the working mode class, this device is not suit to continuous work state.
14. Explanation of some symbols in this device:

~ **AC** AC work mode

OFF Power supply is disconnected

ON Power supply is connected



Equipotential point



Places need to be noticed, please refer to user manual



Device type is CF applied part, which has defibrillation protection function



RS232 connector




USB connector



PATIENT Lead connector

Caddo 10B Panel Sketch Map


Front View



Side view



Rear view



Bottom view



Function button : ON/OFF & Time Display



Function button : plus adjust



Function button : Paper speed adjust



Function button : filter function select



Function button : Pause /on



Function button : switch work mode



Function button : Marker



Function button : print



Function button : system menu



Function button : upwards



Function button : downwards



Function button : leftwards



Function button : rightwards

C. Indicator Definition



The indicator turns green when there is AC power supply, and when the indicator turns green and red same time it is being recharged.



Indicator for instrument when power on.

Operation Regulation

1. You are required to read the operation regulation so as to ensure taking proper operation of the instrument
2. Installation and maintenance of the instrument shall be carried out as the following:
3. There shouldn't have high voltage cable, X radial engine, ultrasound instruments and electrotherapeutics engine around the ECG.
4. Do not install the instrument in the place where it might be affected by bad humidity and ventilation, direct sunlight, as well as air containing dust, salt, and sulphur, etc.
5. The device should be placed in evenness, move gently and should avoid the strong vibration and impact.
6. AC frequency and voltage value should be accorded with the need and the current capacity should be enough
7. Do the instrument grounding properly during installation. Don't put the patients and the lead which connect with patients contact with other conductors, including the ground or the sickbed which ground properly.

Preparation before Operation

1. Check that the instrument properly grounded and that cable connections safe or not.
2. Check the electrode which connected with patient safe or not
3. When power supply is direct current (UPS), please check the voltage of battery before use.
4. The gel should be separated from each other and the chest electrodes shouldn't be contacted with the others, as this operation can avoid short circuit.
5. The AC power supply cable and leads should be separated

www.hik-consulting.pl/edu

Attention during Operation

1. Keep close observation of state of the patient and instrument.
2. Make sure that the patient and device only be connected by leads.
3. The device and patient can't be moved during working.
4. Turn off device after operation.
5. Pull out the power supply plug, then move the leads lightly.
6. Tidy up the devices and accessories for next use.
7. The installation recording paper.



8. This device use high-speed thermal paper whose specification is 50mm (W) x 20m (L).
9. First, open the paper cabinet, take out the paper axis, then put the paper axis in recording paper, then put it in the relevant position of paper cabinet.
10. Cover the paper cabinet with paper cabinet cover, 2cm of the beginning of paper should be left out of the cabinet exit.


Recording Paper Loading

1. If the recording paper is used up during the recording process, the paper record will over, and a notice will be displayed on the LCD screen.
2. There is a line at the verge of paper at the last two meters of the recording paper, this line means the paper is not enough, please change the paper immediately. We suggest you choose our company's print paper, as for its detailed information, please consult with our company or agency.
3. The possible reason which will make the recording paper disable includes: high temperature, humidity and sunshine irradiation. The recording paper which needs long time stock should be deposit in dry, dark and cool environment.
4. The instance which may contaminate the recording paper are Gel, glue and wet diazo compound paper including their organic solvent.
5. The materials that may cause the record wave disappear: the folder contain soft PVC; plastic; the demagnetize ware and tape contain elasticizer; the high-lighter pen, stamp-pad ink, and so on.

Electrode Installation

You'd better install the chest electrode firstly, then Limb electrode.

Chest electrode, as shown in figure



Chest Electrode Positions

The position of installing chest electrodes are as following:

V1: Fourth inter-costal space at right border of sternum.

V2: Fourth inter-costal space at left border of sternum.

V3: Midway between V2 and V 4

V4: Fifth inter-costal space at left mid-clavicular line.

V5: Left anterior auxiliary line at the horizontal lever of V4.

V6: Left mid-auxiliary line at the horizontal lever of V4.

Cleaning the chest skin with alcohol, then put the gel in the diameter about 25mm and the edge of the chest electrodes, press the ball of the chest electrodes, the chest electrodes will be attracted in the position of VI-V6.

Attention: The chest electrodes should be separated from gel coats, this operation can avoid short circuit.

Limb Electrode**Limb electrode Positions**

Clean all the limb electrodes and the positions around to which limb electrodes are to be attached with alcohol before applying ECG cream to them, then firmly attach the electrodes to the positions.

Caddo 10B

Attention: the fix knob should be screwed down tightly after lead connected with main unit.



Check-List for Electrode connection and ECG cable








Electrode Location	Electrode Code	Socket Number
Right Alarm	RA/R	9
Left Alarm	LA/L	10
Left Leg	LL/F	11
Right Leg	RL/N	14
Chest I	VI/C1	12
Chest 2	V2/C2	1
Chest 3	V3/C3	2
Chest 4	V 4/C4	3
Chest 5	VS/CS	4
Chest 6	V6/C6	5


Grounding and Power Connection

Make sure the status of the instrument is power off, and then make the instrument be properly grounded through a 3-prong outlet. When the outlet, a grounding cable may be utilized to connect the grounding terminal of the instrument. Do not use other pipeline. Properly grounding could guarantee the safety and prevent from the interference of AC power and electromagnetic wave.

Battery Operation Regulation

1. This device includes built-in chargeable lithium battery, which doesn't need maintenance. This battery is with perfect automatic charge and discharge monitoring system. When you connect power supply adapter with alternating current, the charge will start automatically. When this device is open, an icon  will be displayed on the top right corner of the LCD screen.  means the battery is charging. The whole charge process needs four hours.
2. When the battery is full, the device can be operated for one hour, when the battery is used as power supply. An icon of battery will be displayed in the LCD screen of the front panel, this icon includes five degrees indicating power of battery. When the battery is power off, the device will turn off automatically, this setting is for avoiding permanent damage on battery caused by excessive discharge.
3. Please charge the battery after power off. When this device is deposited for a long time, the battery should be charged once every six months, this operation will prolong the use-span of battery.
4. The icon of seven different states of power supply as following:

	The alternating current is power supply & the battery is full or no battery
	The battery is only power supply and its power is full
  	The battery is only power supply and its power is not full
	The battery is only power supply and its power is exhausted.
	Charge up

5. If the battery is full, but the power of battery is exhausted within 10 minutes. Please change new battery. If the battery can't be charged, please change new battery.
6. When the icon  displays on screen. Please charge the battery immediately, or the device will turn off.


Caddo 10B


Warning

- Please don't connect the anode and cathode with lead of battery directly, it will cause danger.
- Please don't put the battery on fire. It may cause explosion.
- Please don't disassemble the battery privately.
- The battery should be take gently, please don't strike it with other article.


www.hik-consulting.pl/edu

Keypad and Controls

1. Press power supply key , for several seconds, the device will enter auto-check mode, at this time, the display will be boot-strap menu.
2. After auto-check mode, the display as following:



1. The operation of lead indicate column

Press button  to choose relevant lead, the device will switch to appointed lead check state, it switch among according following order: I, II, III, aVR, aVL, aVF, VI, V2, V3, V4, V5, V6.

2. The operation of system state information column:

Switch by press relevant function key (The function key as following)

Sensitivity: 5mm/mV, 10mm/mV, 20mm/mV, three kinds of sensitivity in all.

Switch Mode: MANU, AUTO.

Under AUTO-MODE, the device will note 12 leads, 3 second ECG signal every lead.

Filter: OFF, 50Hz, 60Hz, 50Hz+, 60Hz+, five filter mode in all.

The mode of 50Hz+ & 60Hz+ mean open 35Hz EMG filter.

Attention: The range of recording R wave will be fallen a little, which caused by attaching the EMG filter

Speed: 25mm/s, 50mm/s two kind of paper speed in all.

3. Leads state indication.

Caddo 10B



When the leads state is "NORM", you can print the ECG.




When the leads state is "OVER", you can't print the ECG, please check whether electrodes are placed well. Stop printing and print date again after collecting the wave.

When the leads state is "SAT", printed ECG is disordered, please check whether electrodes are placed well. Stop printing and print date again after collecting the wave.


When the leads state is "DROP", leads shown on the screen have been off .Please reconnect them.

4. Print operation

Press  under this state, you can start print system setup and ECG wave, press  again the device will be turned off.


Attention: when the paper cabinet is empty, press  or  the device will indicate no paper, please put in the paper then press .

5. Mark operation


Press  you can print a 1mv standard voltage marker, which is helpful to know current sensitivity. Attention: the marking procedure is automatically, after this procedure you need not press any key, the interface will be back automatically.


6. Operation of waveform frozen

Press  you can freeze current waveform in LCD screen, which is helpful for preview.







Press  again, back to previous interface

7. Operation of turning off

Press  for several seconds, the device will be turned off.

System menu

(1) **Operation of menu**

Press  to enter above interface, you can choose relevant item by press  , then you can press   to adjust the content, after setup, press  to be back.

(2) **Introduction of every item**

Backlight: 0-99seconds, back light start time, when choosing 0s, the back light will be turned off, when choosing 99s, the backlight will be turned on for 99s.

Contrast: 00-20, please choose different contrast degree according to different device state.

Demo: ON, OFF, if you need not inspection practice, just choose ON for demo.


About: Software Version

Troubleshooting

Automatic Switch off

1. Please check whether the power of battery is used up. Turn off is for protecting circuit.
2. Please check whether the alternating current voltage is too high, Turn off is for protecting circuit.
3. Please check whether the alternating current disturb too high, whether the fix knob of lead plug too tight. Shut automatically is for protecting circuit when overload.

AC interference




1. Is the ECG device ground cable proper?
2. Is the electrode and leads' ground cable proper?
3. Is the electrode and skin covered with enough Gel?
4. Is the metal bed grounding proper?
5. Does the patient touch the wall or metal sickbed?
6. Does other people touch the patient?
7. Whether there is powerful electric device working beside ECG device? For example: X radial device or B-Ultrasound devices.

EMG interference



1. Whether the patient room is comfortable.
2. Is the patient nervous?
3. Is the sickbed too narrow?

Baseline drift

1. Is the installation of the electrode instability?
2. Is the connection between leads and electrodes credibility?
3. Check the cleaning of electrode and patient skin. Is the electrode and skin covered with enough Gel?
4. Does it cause by the patients' moving or breathing?
5. Is the connection between lead and electrode proper?

Please use filter if still having above-mentioned interference.

Caddo 10B

Troubleshooting List

Phenomenon	Reason	Resolve method
Disturbance too big, the waveform is in disorder	<ol style="list-style-type: none">1. Whether the ground cable proper.2. The connection of leads is not stable.3. Whether there is disturbance from alternating current.4. Patient is nervous	<ol style="list-style-type: none">1. Please check the lead, ground cable and power supply.2. Please dispose the patient in proper state.
Baseline is rough	<ol style="list-style-type: none">1. Disturbance from alternating current too strong2. Patient is nervous and the disturbance EMG too strong	<ol style="list-style-type: none">1. Change a comfortable environment for patient2. If the sickbed is metal, please change it.3. The power line and lead is not parallel or too close.
Wave form is not regular, with too great wave or beeline	<ol style="list-style-type: none">1. The conductivity of electrode is not well.2. Power of battery is used up3. Contact between electrode and skin is not proper4. The plug between lead and main unit is not tight.5. The contact between lead and electrode is not proper.	<ol style="list-style-type: none">1. Use alcohol of high quality.2. Clean the electrode and patient's skin where touch the electrode3. Charge the Battery4. Keep the electrode reed clamping.
Baseline drift	<ol style="list-style-type: none">1. Power of Battery is used up.2. Patient is moving	<ol style="list-style-type: none">1. Charge the battery2. Keep patient hold still
Waveform is not clear	<ol style="list-style-type: none">1. The printer head is dirty2. The paper is not right	<ol style="list-style-type: none">1. Clean the printer head with alcohol when the power is off, use the printer head after the alcohol is volatile.2. Use the appointed thermal print paper.

Maintenance Transportation and Preservation

1. Customer is not permitted to open the instrument, in archive of any electronic shock. Any maintenance or update should execute by the trained and authorised professionals from our company. The maintenance should be done with the original accessories from our company.
2. Please pull out the power supply plug when power off. If the device out of use for long time, please put the device in a shady cool dry place, and the device should be charged once every three months.

www.hik-consulting.pl/edu

Warranty

- 1) We guarantee this product against all manufacturing defects for 12 months from the date of sale by us or through our dealers. Consumables like dry cell etc. are not covered under warranty.
- 2) The guarantee will become void, if
 - a) The product is not operated as per the instruction given in the Learning Material
 - b) The agreed payment terms and other conditions of sale are not followed.
 - c) The customer resells the instrument to another party.
 - d) Any attempt is made to service and modify the instrument.
- 3) The non-working of the product is to be communicated to us immediately giving full details of the complaints and defects noticed specifically mentioning the type, serial number of the product and date of purchase etc.
- 4) The repair work will be carried out, provided the product is dispatched securely packed and insured. The transportation charges shall be borne by the customer.

List of Accessories

1. Learning Material1 No.
2. Suction electrode 6 Nos.
3. Clamp electrode..... 4 Nos.
4. Electrode lead cable1 No.
5. ECG roll paper.....1 No.

B-Ultrasound Diagnostic System

Caddo 12B

Learning Material

Ver 1.1

An ISO 9001:2008 company



Innovative Technology Ecosystem

94, Electronic Complex, Pardeshipura

Indore - 452 010 India

Tel: 91-731 4211100

Fax: 91-731-2555643

e mail: info@scientechnbz

Websites: www.caddo.bz

www.scientechnbz

www.scientechnbz

Certificate

Standard: **ISO 9001:2008**

Certificate Registr. No. **85 100 001 10182**

TÜV Rheinland India Pvt. Ltd.;

Certificate Holder: **Sciencetech Technologies Pvt. Ltd.**

Unit 1: 94 – 101, Electronics Complex, Pardeshi Pura,
Indore – 452 010, Madhya Pradesh, India

Unit 2: 90 – 91, Electronics Complex, Pardeshi Pura,
Indore – 452 010, Madhya Pradesh, India

Scope: **Design, Manufacture of Electronic Test & Measuring
Instruments, Training Products for Electrical & Electronics
Education and Providing Technology Training**

An audit was performed, Report No. 10182. Proof has been
furnished that the requirements according to ISO 9001:2008
are fulfilled.

The due date for all future audits is 04-10 (dd.mm).

Validity: The certificate is valid from **2010-12-13** until **2013-12-12**.



Bangalore, 2010-12-20

The Certification Body of
TÜV Rheinland India Pvt. Ltd.

The validity of this certificate is subject to timely completion of surveillance audits as agreed

in the Contract. The validity of the Certificate can be verified under 00000_Lux2008 with the identification No. 0100027000

000000000000

www.tuv.com

 **TÜVRheinland®**
Precisely Right.

B-Ultrasound Diagnostic System**Caddo 12B****Table of Contents**

Introduction	6
Chapter One Summary	8
1.0 Features1.	
2.0 Technical specification	9
3.0 Panel chart	11
4.0 Basic principle	11
5.0 Equipment constituent	11
6.0 Standard configuration pieces	11
7.0 Optional pieces	12
8.0 Shape	12
Chapter Two Installation	14
2.1 Operating environmental requirements	14
2.2 Unpacking in specification	14
2.3 Installation	14
• Connecting probe and main unit	14
• Installingandremoving battery	15
2.4 Power Supply	16
• Power supply with adapter	16
• Battery Operation	16
2.5 Battery charging	16
• Charging through main unit	16
Chapter Three Key boardand Trackball Operation	18
3.1 Screen display.	18
3.2 Operatingkeyboard	18
3.3 Trackball.	18
Chapter Four Operation Procedures	28
4.1 Start the equipment	28
4.2 Diagnosis	28
4.3 Notation	28

Caddo 12B

4.4	Probe frequency switchover	31
4.5	Body mark	32
4.6	Measurement	34
	• Distance measurement	34
	• Circumference and Area Measurement	35
	• Volume measurement	38
	• Heart rate measurement(only at B/M modes)	39
4.7	Recall Penetration Sampling Line	40
4.8	Cine loop and image storage function	41
	• Saving a image	41
	• View a saved image.	41
	• Review story board mode	42
	• Obstetrics calculation	43
	• Report	46
	• Image processing	47
	• Perrectum examination procedure	49
	• Pervaginam examination procedure	49
	• Image print	50
	• Image up load to computer	50
	• Shutdown the equipment	50
Chapter Five Transportation and Storage.....		51
5.1Environmental Requirements off transportation and Storage.....		51
5.2 Transportation.....		51
5.3 Storage		51
Chapter Six Check and Maintenance. ...		52
6.1 Check		52
6.2 Service life.....		52
6.3 Main unit maintenance.....		52
6.4 Probe maintenance.....		52
6.5 Cleansing.....		41
6.6 Correct usage of probe.....		54

Caddo 12B

6.7 Battery information.....	55
6.8 Instrument test and calibration	57
Chapter Seven Malfunction Examination and Troubleshooting	58
7.1 Examination.....	58
7.2 Troubleshooting.....	58
Appendix A Sound Output Specifications Table.....	59
Appendix B Gestational Table.....	60
Warranty	91
List of Accessories	91

Introduction

Matters need Attention :

To ensure operational safety and long-term stable equipment performance, please read this operation manual closely and understand the device functions, operation and maintenance at all points before operating the device, especially contents of "Warning", "Caution" and "Note".

Mis operation or inobservance of the instructions given by manufacturer or its agents may result in device damage or personal injury.

The following convention works through this manual to lay special emphasis on some information.

"Warning": Stands for neglect of it will cause severe personal injury, death or realized property loss. "Caution": Stands for neglect of it will cause slight personal injury or property damage.

Note": To remind user of installation, operation or maintenance information. These information is very significant but with no risk. Any warning against dangers shall not be contained in note.

Safety labels

Device labels explanation:



Type B device



Note! Refer to accompanying documents



Turn-on (general supply)




Disconnect (general supply)



Signal output



USB port



RS-232







IPX7 Water tight-proof



II type equipment

Caddo 12B

Packing and transportation labels explanation:

	Handle carefully
	Temperature limit
	Upward
	Piling limit
	Keep dry
	Protect against heat

- a. Environment temperature range: + 10°C - + 40°C
- b. Relative humidity range: 30% ~ 75%
- c. Atmosphere pressure range: 70KPa ~ 106KPa

Environment requirements on device storage and transportation:

- a. Environment temperature range: - 40°C ~ + 55°C
- b. Relative humidity range: 10% - 100 %
- c. Atmosphere pressure range: 50KPa ~ 106KPa

Tabu illness:

- This equipment is not suit for contain-gas viscera, such as lung examine.
- Suggesting not to examine the wounded or urgent inflammation part, it could Cause across infect.
- Below patient forbid to use vagina and recta probe:
Vagina colpitis, such as colpitis and sexual illness discoverture people vagina abnormality menses vagina shrink vagina ultrasound exame difficulty person vagina bleeding placenta front sufferer.
- Below sufferer is forbidden to puncture:
High blood pressure, coronary heart disease, cruor obstacle and haemorrhage Separative sufferer.
Implementary criteria (Safety)

Chapter One Summary

1.1 Features

This equipment is high resolution linear/convex ultrasound scanning diagnostic equipment. It adopts microcomputer control and digital scan converter (DSC), digital beam-forming (DBF), real time dynamic aperture (RDA), real time dynamic receiving apodization, real time Dynamic receiving focusing (DRF), Digital frequency Scan (DFS), 8 segments TGC, frame correlation technologies to endue its image with clarity, stability and high resolution.

There are five display modes: B, B+B, B+M, M, 4B And 256 gray scale.

The system can process real time image display, freeze, save, load, zoom, up and down flip, left and right flip, black and white flip, and capacity cine loop Multi-level scanning depth, angle, dynamic range, acoustic power, frame correlation factor regulation and focus number, focal space, focus position, etc. It offers more than 40 body marks.

Date, clock display Name, sex, age, doctor, hospital annotation Distance, circumference, area, volume, heart rate measurement preset two obstetric tables to measure GA, FW ~ EDD. Many probes are optional for clinic diagnosis demands.

PAL-D video output offers connection to external video image printer and big display and other equipments. High speed USB port provides real time image transfer to the PC.

Using the touch pucker keyboard and trackball to operate, it's quick, convenience, agility.

The equipment is jet molding enclosure and potable structure, the usage of non-industrial frequency transformer switching power supply, programmable parts (FPGA) and surface mounting technology (SMT) make the whole unit highly compact.

The equipment is build up by host, probe and adapter. The maid probe is 60R/3.5MHz convex probe, the optional probe is 20R/5.0MHz convex probe, 13R/6.5MHz antrum probe, 7.5MHz High frequency linear probe, 7.5MHz recta probe.

1.2 Range of application

Be the same with the body belly viscera's ultrasound examine.

Caddo 12B


Technical Specifications

Model		Caddo 12B			
Probe		Standard	Option		
		60R/3.5 MHz convex array	20R/5.0 MHz Micro convex probe	7.5MHz HF linear probe 7.5MHz endo-rectal probe	13R/6.5MHz endo-vaginal probe
Display depth (mm)		240 (max), 16 levels adjustable			
Maximal detect depth (mm)		≤ 160	≤80	≤80	≤60
Resolution (mm)	Lateral	≤2 (depth ≤80) ≤3 (80< depth	≤3 (depth ≤ 60)	≤1 (depth ≤60)	≤1 (depth ≤40)
	Axial	≤2 (depth ≤ 80 ≤3 (80 <depth	≤1 (depth ≤ 60	≤1 (depth ≤ 80	≤1 (depth ≤ 40
Blind zone (mm)		≤5	≤8	≤3	≤7
Geometric position precision	Horizontal	≤15	≤15	≤5	≤10
	Vertical	≤10	≤10	≤5	≤5
Monitor size		10.4inch			
Display mode		B, B+B, B+M, M, 4B			
Image gray scale		256 level			
Scan depth		40mm-240mm			
Image storage		32 frame			
Scan angle		50% ~ 100%			
Acoustic power		2 steps			
Dynamic range		100dB ~ 130 dB			
Focus position		Adjustable			
Focal space		5 Level			

Caddo 12B

Cine loop	809 frame (max)
Part magnify	2 times (both real time and freeze)
Depth exaltation	B, B +B mode real time exaltation
Frame related coefficient	B, B +B mode eight level adjustable
Image flip	Up/down , left/right, black/white
Measurement	Distance, circumference, area, volume, heart rate. GA, FW, EDD
Notation	Date, time, name, sex, age, doctor, hospital name, full screen words edit.
Posture output	Two types
Power consumption (MAX)	45VA

Panel chart



1.5 Basic principle

Digital Ultrasonic Diagnostic System works in this following procedure: Different tissues of human body possess different densities and speeds of transmission of ultrasound, i.e. different acoustic impedance (product of media density and sound speed). when piezo-chip (transducer) gets certainly regulated electric impulse, it will produce ultrasound with certain frequency. when this ultrasound (sound energy) is injected into human body, different organ surfaces will produce reflection echo, the different size reflection is received by the transducer which emitted ultrasound and is changed into electric impulse, when this electric impulse is amplified, demodulated, digital scanned, shifted and some other handling, video standard signal is produced and organ cross-sectional images are displayed on the monitor

1.6 Equipment constituent

1.6.1 Standard configuration pieces

- Mainframe (included one HYLB-3 19 battery)
- 60R/3.5MHz convex probe
- SP-2 power adapter
- Adapter power cable
- Video frequency connect cable
- Charge wiring
- High speed USB connect cable
- User manual
- Inspection report
- Packing list

Caddo 12B



1.6.2 Optional pieces

20R/5.0MHz convex probe

- 7.5MHz High frequency linear probe.
- 13R/6.5MHz vagina probe
- Monitor
- Video frequency printer--Sony UP-895MD
- Car charger
- HYLB-319 battery

Warning

Customer should use the prescribe fittings listed above, if Use other fittings, the manufacturer should take no responsibility of the risk caused by this reason.



Mainframe shape size's sketch map




Power Supply port

VGA output

Video out

Probe port

Mainframe back port sketch map



Battery

Charging battery indicator light

Charge/electronics level light and battery sketch map

Radiating hole



Under side sketch map of mainframe

Chapter Two Installation

2.1 Operating environmental requirements

- a) Environment temperature range: +5°C ~+ 40°C
- b) Relative humidity range: ~80%
- c) Atmosphere pressure range: 86KPa-106KPa

When using, avoid strenuous vibration, keep it away from devices with high field, intense magnetic field or high voltage avoid strong sunlight blazing down on the display keep the device well-ventilated, moisture proof and dustproof.

2.2 Unpacking inspection

After unpacking, check the device according to "Packing List" and install it according to requirements and methods described in "Installation" after affirm that there is no shipping damage.

Warning

If there is breakage at unpacking check, it is banned to use the device to ensure security.

2.3 Installation

2.3.1 Connecting probe and main unit

The probe socket is in the rear of the equipment and in shape of letter D. Plug the probe connector and then screw down the two screws on the connector to finish connection of the main unit and probe. There is only one probe socket but it is compatible with all the optional probe types of the equipment.

A reverse process can demount the probe from the main unit.

Warning

Do not by all means unplug or plug the probe connector at state of log on in case the probe and main unit be damaged.

Once the probe is connected with the main unit, do not unplug nor plug it at discretion in case poor contact happen.

Warning

The probe should be protected from felling off or crashing and the manufacturer assumes no responsibility for the kind of hazard.

Warning

Must not touch the contact pin in the probe connector.

Warning

Please handle the equipment carefully.

2.3.2 Installing and removing battery

Install battery: Set the battery into the battery slot in the right direction and then fasten the battery with the fixing

Screw (refer to figure 2-1)



Battery Installation

Remove battery: Loosen the fixing screw, pull the battery out from the slot (refer to figure 2-2)



Battery removal

Caddo 12B

Power Supply

The equipment provides two automatic switch-over modes to supply power: adapter and built-in battery.

2.4.1 Power supply with adapter

- 1 Check the input power cord plug of the adapter to see if it matches the EPS outlet.
- 2 Check the EPS to see if it is in the specified range and the power cord to see if it is connected well.
- 3 Check the adaptor to see if it works well: Plug the power cord into the AC input outlet switch on the power switch of the outlet if the DC output indicating light turns on green, it works well.
- 4 Shut the power switch of the outlet.
- 5 Insert the DC output plug of the adapter into the DC 14V 13A outlet in the rear of the equipment Switch on the outlet
- 6 Turn on the main unit and the equipment is ready for operation now.

2.4.2 Battery Operation

- 1 Install the battery correctly on the main unit.
- 2 Switch on the power switch on the left side of the equipment. When the main unit is power on, the power indicator will turn on.
- 3 The equipment can start operation.

Note:

When the charging/electric quantity indicator flashes, it means the battery is running up and needs charging.

Warning:

It is prohibited to use any other power supply except the standard adapter as the external power supply for the main unit

2.5 Battery charging

There are 3 ways to charge the battery.

2.5.1 Charging through main unit

1. Install the battery correctly into the main unit.
2. Insert the plug of "output DC 14 V /3A" of the adapter into the "DC 14 V /3A" interface on the rear.
3. Connect power cord of "Input AC100 - 240V, 50/60Hz" of the adapter to the AC. EPS.
4. No matter the main unit is power on or shut down, when the "Charging / Electric Quantity" indicator turns on in circulation, the adapter is charging the battery When the four LED turn on at the same time, the battery is fully charged.

Tips:

To prolong the service life of the battery, please do not stop charging until the charging indicator light indicates the battery is fully charged.



Figure 2-3. Charging through main unit

Tips:

1. The input voltage of the auto-charger is DC9~14V/1.5A
2. The output voltage of the auto-charger is DC12.6V/1A
3. The operations and storage environment are the same as those of the main unit.

Chapter Three Keyboard and Trackball Operation

3.1 Screen display




Figure 3-1 Equipment interface illustration

3.2 Operating keyboard



 ~  Alphabet keys

Press these keys at annotation mode to put in characters at the cursor position

Among the alphabet keys, there are some keys within dual functions:

Figure 3-2 Operating keyboard sketch

Caddo 12B

**Alphabet keys**

Press these keys at annotation mode to put in characters at the cursor position.

Among the alphabet keys, there are some keys with dual functions:

SN	Key	Function
1		At freeze, press this key to enter or exit report interface (detail in 5. 7)
2		At real time or freeze, press this key to adjust acoustic, the Equipment has 80 percent and 100 percent, can switch to each other circle
3		At real time or freeze, press this key to clear notation and report information on the image or the characters area
4		Total gain control key. There are four steps: 25, 30, 35, 40. Press this key to set them circularly
5		Under the real time estate (single B mode) press this key T 0 display the puncture lead (detail operation see 4.6.5)
6		In real time estate, press this key to switch the probe frequency (detail operation see 4.4)
7		Press this key to over turn up/down of the image (see chart 3-3)
8		Press this key to over turn left/right of the image(see chart 3-4)
9		Under" B", "BB", "4B" modes, press this key to activate The body mark(detail operation see 4.5)




Figure 3-3. Image up and down flip sketch




Figure 3-4. Image left and right flip sketch

Tips:

Probe scan direction indicator (inside the red circle) is the image left/right overturn mark. The equipment default scan mode is the left chart's scan mode.

SPACE Space key

At annotation, press this key to input space after cursor.



Backspace key .

At annotation, press this key to delete the input characters.

At non-annotation, this key is the fake color control key, the equipment provides five display Colours, when start up, the equipment's default colour is blue (green background), press This key, the five colour cycle display, then choose the colour needed.

The sequence is: blue (green background), yellow (blue background), red (light blue background), black and white (red background), black and white.






0 ~ 9 Number character and direction key

Numbers are used for time, date settings, age notation and function menu and selection.

'when doing the measurement, the direction key is used to move the measure cursor, sampling Line, chart Statistic window, puncture lead, part enlarge sampling window and enlarge window.

Number character includes below double function keys:

Caddo 12B

SN	Key	Function
1		Press this key to move cursor or sampling line left
2		Press this key to move cursor up
3		Press this key to move cursor down
4		Press this key to move cursor or sampling line right
5		Control the cursor moving speed by direction keys

**Clear**

Press this key to clear off the measuring marks, notation data and measurement results on the image (When some menu is displayed, it can not be cleared until the current task is finished).

**Notation**


Press this key to display the notation menu and press the above mentioned numbers or letters for notation.

**OB Table**


Press this key to display the OB table and press the number keys to carry out the corresponding function.

**Image Post-process menu**

Press this key to display image post-process menu and press the number keys to carry out the corresponding function.

**Reset**

Press this key to restart the system when malfunction and mis operation make the system "DEAD"(unable to start the system by pressing any other key).

**Frame-correlation factors regulation /M mode velocity option**

This key has two functions.

At "B/M" mode or "M" mode, enter real time state, press this key to adjust the refresh

Caddo 12B

velocity of "M mode image"

. "B/M" mode equipment has 3.00S, 2.50S, 2.00S, 1.25S, four level reflash speed, "M" mode equipment has 6.00S, 5.00S, 4.00S, 2.50S four level reflash speed, Press this key repeatedly to set the velocity circularly. The current velocity displays on the right of the screen.






Figure 3-5 M speed display sketch map

At "B", "B/B" and "4B" modes enter real time state press this key the "frame correlation factors." will display on the right top of the screen, then press the  or  to adjust the parameters of frame correlation factors. The system provides eight factors: 0.25, 0.35, 0.45, 0.55, 0.65, 0.75, 0.85 and 0.95, Press the continuous key, the eight factors is setting circularly. The current parameter appears on the right top of the screen (see figure 3-1).



Save image

At real time of freeze, press this key to save the current image: refer to "4.5.1 Image Saving" for details




Load Image

Press this key to recall the stored image refer to "4.5.2 image loading" for details






Measurement

For volume, fetus weight measurement and histogram (see 4.4 for details)

**Cine loop**

At freeze state, press this key to enter cine loop, refer to "4.5.3 Cine Loop" for details

**Measurement operation**

Cooperate with    keys to finish the measurement of distance, girth, area, cubage. Detail operation see chapter four.

**Confirmation**

For menu selection when operating by trackball.

Confirm the contents when noting.

At real time or freeze, press this key to display body mark ←

**Cancel key**

Within the annotate operation, cancel the annotate comment

Within the distance measurement, finish the distance's measurement.

**B display model**


At any state, press this key to enter B mode (default mode is single B at starting-up)

**Double B display mode key**

Press this key to enter BB mode. And there will be two B mode images displayed on the screen. One of them is a "Frozen" image and the other is a "real-time" image. Repress this key can make the two images switch between the states of "Frozen" or "real-time". Press "Freeze" key to freeze these two images


**BM/M mode display key**

Press this key to enter B/M mode. B mode and M mode images will be displayed on the screen ("BM" or "B+M" for shortened form). At the same time, B mode real time image is on the left and M mode real time image is on the right.


**4B mode display Key**

Press this key again to enter M mode, then M mode image will be displayed on the screen. (Press this key to convert from BM to M, M to BM.) mode display key


Under the single B or BB mode, press this key to enter 4B mode and display 4B mode images. One is real time, the other three are all frozen images Press this key repeatedly to switch each image between "freeze" and "real time" circularly. Press the freeze key to get four "Frozen" Images




Single B mode




BB Mode



BM mode



M mode



4B mode

Figure 3-6. Five display modes sketch

**Adjusting keys**

Adjusting, keys are used to adjust focus position, dynamic range, scanning angle, zoom, cine loop See" Chapter 4 Operation Procedures")

**Focus number selection**

Press this key to lighten up the" focus number" on the screen right top, press again to change the focus numbers.



There are two focus at most. Press this key to switch between one focus to two.

**Focus space selection**

When set the focus number as two, their distance can be adjusted. There are 5 levels focus space 2, 3, 4, 5, 6 and press this key repeatedly to switch among them circularly.

Caddo 12B

**Scan angle selection key**



Press this key to lighten up the "Angle" on the up right screen, press  or  to change the scanning angles. The angel range is 50%-100%.

Tips: When linear probe or recta/linear probe is used, this function can not be used. Reduce the angle can improve the image's frame rate and increase the cine loop frames

**Dynamic range adjusting key**


Press this key to lighten up "DR" on the screen top right corner. Press dynamic key to adjust the range between 100dB-131dB. 118 dB is recommended

**Depth Adjusting key**

Press  or  keys to adjust scan depth between 40mm-240mm. The current depth displays at bottom screen (see figure 3-1). The depth adjusting range of each probe is given below:

Probe mode	60R/3.5MHz Convex probe	7. 5 MHz Endo-rectal probe 7. 5MHz HF linear probe	20R/5.0MHz Micro convex probe	13R/6. 5MHz Endo-cavity probe
Depth range (mm)	70~240	40~90	80~150	50~120

Table 3-1. Probes and their depth adjusting range



Zoom
Zoom key

At real time or freeze state, press this key to display or close sampling windows and zooming window.

Move "zooming window" and "sampling window" by trackball. Press **Reference** Key to switch the current window





Figure 3-7. Image zooming sketch

Tips:

Change the sizes of "zooming window" and "sampling window" by pressing **↑** or **↓** keys at the same time. The windows change pro rata. The maximum size of zooming window is a quarter of the image display area of the screen.

**Image freeze key**

Overall gain adjusting key Press Δ . key to increase the overall gain and press ∇ to reduce the overall Gain.

**Near field gain adjusting key**

Press Δ . key to increase the Near field gain and press ∇ to reduce the Near field Gain.

**Far field gain adjusting key**

Press Δ . key to increase the Far field gain and press ∇ to reduce far field Gain.

The current near field gain, far field gain and overall gain are displayed respectively after on the upper screen.

Tips:

under real "G" key is the overall gain adjust key, equipment has 25,30,35,40, totally four level parameter, continually press this key, those four level parameter will circle setting.

3.3 Track ball.

Track ball is fast, convenient in operation. In this equipment, the trackball functions show as below:

- Move measuring mark during measurement.
- Select these items in the menu during operation.
- Move the cursor at annotation state.
- Move the sampling line at BM mode.
- Display the puncture guide line's and moves it.
- Control movie replay on single frame show at cine loop mode.
- Move the sampling window and zooming window at zooming state.
- At histogram statistic mode, Move the statistic sampling window.

Tips:

1. Do not press the trackball too heavy when operating the trackball.
2. Keep the trackball surface clean.

Chapter Four Operation Procedures

4.1 Start the equipment

Switch on the equipment, indicator light on the panel turns on and the start interface shows on screen. Press any key to enter the state of scanning. Adjust the display brightness, contrast, overall gain to acquire a satisfying visual effect.

4.2 Diagnosis

Apply an amount of contact agent (medical ultrasound coupling gel) to the diagnostic parts of patient body, and then press the acoustical window of the probe on the part closely. The cross-section acoustic image of tissues will be displayed on the screen, move the probe properly to locate the optimum position for the investigation depth meanwhile, adjust the overall gain, near field gain, far field gain or the focus properly to acquire an optimal cross section acoustic image of the diagnostic part.


Note

1. Do not press the probe too much to avoid breaking the probe or causing discomfort.
2. Do use proper probe and frequency for the target part.

Annotation
ABC


Notation

Press key to activate notation menu. Show as below:



1. NAME
2. AGE
3. SEX
4. COMMENT
5. TIME
6. HASP
7. DOCT
8. LANGUAGE
9. ERASE
A. OBTABLE
B. EXIT

Figure 4-2 Notation menu

Press  key to select" 1. NAME" to enter the patient name, 15 characters maximum

can be put in (letter, number or space), press  to delete the input mistakes.

Caddo 12B

(Same as the following comment) the name input box will appear on screen as follows:

Please Enter Name :

After input, press any key (except letter, number keys) to exit. To cancel input, you can press any key. The name will display after "NAME" on top of the screen. For example: input WANG NING, then press any non-character



Key (e.g.) to finish input.

Press "2" to select "2. AGE" to put in patient age, 3 characters at most. The age input box is given below:

Please Enter Age:

After input, press any key (except letter, number keys) to exit. To cancel input, you can press any key. The age will display after "AGE" on top of the screen. For example: input 27, then press any non-number key.



To finish input.

Press "3" to select "3. SEX" to put in patient sex. The input box is given below:

Please Enter Sex:
1. Male 2. Female

Press 1 to select "1. MALE" and "2" to select "2. FEMALE".

Press "4" to select "4. COMMENT" to enter image notation. The white cursor will display on the screen.

Move the trackball cursor, you can annotate anywhere within the image area. At the same time the bottom of the screen will display the shortcut key, press the corresponding number keys (0-9), the cursor will input the same characters.

When finishing the comment press "Annotation ABC" key on the board to end this

operation or press "Cancel" to cancel the comment.

Tips : At real time or freeze state, Press "Confirm" to approve the location mark ✓,

then press "Reference" to change the mark's direction, use the direction key or trackball to mark the image area or body specimen.

Caddo 12B

Press **5** to select "5. Time to modify the system's time and date. The time and date input box is listed as bellow:

YY-MM-DD
HH-MM-SS

For example, the date and time are: 2006-9-22 9:35:30, then the input should be:

YY-MM-DD
060922
HH-MM-SS
093530

Press **6** to select "6. HOSP" to input hospital name, input letters or numbers maximum 18 characters. The hospital name input box is listed bellow:

Please Enter HOSP :

After input, press any key (except letter and number) to confirm exit. To cancel input, you can press any key. The hospital name will display after "Hospital" at the bottom of the screen. (see figure 3-1 shows)

Press **7** to select "7. DOCT" to put in doctor's name in letter or number maximum 15 characters. The doctor name input box is listed bellow:

Please Enter Doct :

Press **8** to select "8. LANGUAGE", choose the system's interface language.

Pre install function.

Press: **9** to select "9. ERASE", erase image's storage section. The remind box is given bellow:

Erase All Storage
1 Yes 2. No

Press number key **1** to confirm clearance. When processing, a word "ERASING..." displays on the upper left screen to indicate that system is processing clearing and no other operation can be done at this moment. When the reminder disappears, the image storage is cleared.

Press **2** to give up and exit.

Caddo 12B

Press "**A**" to select "A. OBT ABLE" to put in obstetric table. The system is built in two OB tables "TOKYOHADLOCK". Table TOKYO is suitable for Asian, and Table HADLOCK is suitable for European. The remind box is given bellow:

Please Enter Obtable:
1.Tokyo 2.Hadlock

Press "**1**" to select TOKYO, press "**2**" to select HADLOCK.

Tips:

Head circumference (HC). abdomen circumference (A C) only have the Hadlock form, this equipment does Not provide the HC, AC form under the Tokyo formula


Press "**B**" to select "B. EXIT" to exit comment menu.

Tips:

1. When image storage erasing is in process (system reminds "Erasing"), please do not take other operations in case the equipment is damaged.
2. The trackball can be used for operating the menu: when the menu displays, move the trackball up and down to select the items, when the selected item is lightened, press to enter it
3. **N** Key functions:
 - At real time or freeze state, press this key to clear off the measure marks, comments, measurement results, doctor, name, age and sex in image area
 - At report interface, clear off all the information and measurement results except "Hospital" doctor" and "comment" information.

4.4 Probe frequency switchover



At real time estate, press  key, change the working &frequency of the probe. The & frequency of each probe is given bellow:

2.5 MHz, 3.5 MHz, 5.0 MHz-60R/3.5 MHz Convex Array

4.5 MHz, 5.0 MHz, 5.5 MHz-20R/5.0 MHz Micro Convex Probe

5.5 MHz, 6.5 MHz, 7.5 MHz-R13/6.5 MHz Endo-Vaginal Probe


6.5 MHz, 7.5 MHz, 8.5 MHz-7.5 MHz HF linear probe or 7.5 MHz endo-rectal probe

The current probe &frequency displays on the up right of the screen (see figure 3-1).

Tips:

1. The device can automatically identify probes and display the probe's model when it get into the scan working estate (see the top right corner of chart 3-1).
2. Please shut down the system first before replacing probes. Restart the system, it can realize automatic with default setting.

4.5 Body mark


At B, B/B, B/M, 4B modes, press  key to active body mark. There are 40 body marks, press this key repeatedly to display these marks circularly.

1. abdomen

Body front (M)



Body front (F)



Body back



Body left side



Body right side

2. Obstetrics:

Mother (F)



Fetus



Fetus

Caddo 12B


3. Head:



4. Heart:



5. Others:



R



L



4.6 Measurement

Distance, circumference, area and volume can be measured by controlling the direction keys or trackball Unit of distance and circumference is mm, unit of area is mm² and volume is cm³

There are four measuring marks corresponding to measurement display:

+: D1

x: D2

*: D3





*: D4

In circumference, area measurement, there have two measuring marks corresponding to measurement display as below:

+: C1, A1

x: C2, A2

4.6.1 Distance measurement Operation steps:

1. Press “” key, and the first cursor appears on the screen
2. Press  "key, the second cursor appears on the screen Reference
3. Move the cursor to the start point of this measurement by controlling the trackball, press  key again to confirm the start point
4. Move another cursor to the end of this measurement, press" "to fulfill the measurement and exit (Tips: Press"  " key repeatedly to switch between the start and end points of the measurement).
5. Go on measuring just repeat step 1-4 to get at most 4 groups of distance measurement. The measuring results display on the screen right as the following figure shows




Figure 4-3. Distance measurement sketch

The four groups of data are D1, D2, D3, D4, among them

D1/D2 is the ratio of D1 and D2

D3/D4 is the ratio of D3 and D4.

1.6.2 Circumference and Area Measurement

The circumference and area can be measured in two ways.


Press A/C key to display this menu as given below:

Put in :
1. Freehand 2. Ellipse

Option 1 is FREEHAND measuring method and option 2 is ELLIPSE measuring method.

a. Freehand method operation:

1. Press 1 key to select freehand method and call out the measuring cursor, move the cursor to the start point of the measurement.
2. Press Reference key and move the cursor along the edge of the measurement area to the end point by controlling the trackball
3. Press Reference key again to end the measuring of circumference and area.

4. Press  key and repeat the above step 2-3 to process another measurement. Two groups of data can be measured at most. And the measuring results display on the screen right as given below:




Figure 4- 4 Circumference and area measuring sketch (Freehand method)

C1 and A 1 are the circumference and area of the first group data



C2 and A2 are the circumference and area of the second group data

C2 and A2 are the circumference and area of the second group data

C1/C2 is the ratio of the two circumference


A1/A2 is the ratio of the two areas.


b. Ellipse measuring method

1. Press  key to select Ellipse method, an ellipse measuring mark display on the screen, which is named the measured area, can be moved to any place on the image by trackball.
2. Press  key, then move the trackball to change the size of this measured area.

Move trackball left and right, the measured area shrinks or enlarges horizontally
Move trackball up and down to shrink or enlarge the measured area vertically key to finish measuring.

Tips:

Press  key repeatedly to shift between moving the mark and adjust the size by trackball.

3. Press  key, Move trackball to adjust the angle of the measured area.

4. Press  key to finish measuring.


5. Press  key and repeat the above steps 2-4 to process another measurement. Two group's data can be measured at most. And the measuring results display on the screen right as given below:



Figure 4-5 Circumference and area measuring sketch (Ellipse method)

C 1 and A 1 are the circumference and area of the first group data

C2 and A2 are the circumference and area of the second group data

C 1/C2 is the ratio of the two circumferences

A 1 / A2 is the ratio of the two areas.

Attention

Ellipse measure operation can only be done under the freeze estate.

4.6.3 Volume measurement

Volume can be measured in two ways.

1. Volume measurement utilizes 3-axle method to measure 3 groups of distances and then calculate them.

To fulfill volume measurement, first measure three distances and then press

Measure

key to get the volume.

If the measured distance data is less than three groups, pressing Measure does not get volume if four distances information are measured, when pressing

Measure

key, the result is the value of the first three data (D1, D2, D3).

Operation procedure: (take kidney as an example)


- Get the kidney cross section image and profile section image respectively and freeze them.
- Measure the long axis and short axis lengths of kidney cross section with distance measurement method.
- Measure the diameter of kidney profile section with distance measurement method.
- Press Measure key to fulfill volume measurement. The volume cost displays at the lower-right corner of the screen behind "Vml" as the following figure shows



Figure 4-6. volume measurement sketch (3-axil method)

2. Utilizes ellipse measuring method to measure 2 groups of areas and then calculate them.

Operation procedure: (take kidney as an example)

- Get the kidney cross section image and profile section image respectively and freeze them.
- Measure the circumference and area of kidney cross and profile section with circumference and area measuring method.
- Press  key to fulfill volume measurement. The volume value will be calculated automatically and displays on screen right behind "Vml" as the following figure shows:




Figure 4-7. volume measurement sketch (Ellipse method)

Tips:

If only a group of circumference and area are measured, Volume will not be displayed.

4.6.4 Heart rate measurement (only at RIM modes)

1. At B/M modes, move trackball or direction keys to change the sampling line position and freeze it to acquire a satisfying heart beat oscillogram.
2. Measure the distance between the three cardiac cycles wave peaks according to distance measurement method 4 groups data display at screen lower-right corner, they are: heart rate (HR) (unit: B.P.M.), rate of slope (EF) (unit: mm/s), hour (T) (unit: mS), refresh velocity (ST) (unit: s). As following figure shows:




Figure 4-8. Heart rate measurement

Tips:

1. Please identify the functions of **[N]** key and **[Clear]** key: Press ~ key to delete
2. measuring marks and results, press **[N]** key to clear off all annotation notes (including Doctor, name, age and sex) and measuring marks, results as well as report data.
3. You can fulfill measurement by operating **[1][2][4][3]** direction keys as well as trackball.

4.7 Recall Penetration Sampling Line

Under the Real Time Estate (single B mode), press **[P]** key, the screen will display the PT Sampling line, it calls puncture lead line, the trackball or left/right direction keys can move to confirm the position,

if want to exit, please press **[Cancel]** key


Warning


There is risk exist in the ultrasound lead puncture, it has to be Implemented by qualified people, it has to strictly prepare the Before surgent work, including to check the cruor time, blood platedet calculation, cardiogram, blood pressure, puncture bag, puncture probe antisepsis and sign the surgery treatment agreement.


4.8 Cine loop and image storage function

The system provides separate memory area for "storyboard view" and "image save"


4.8.1 Saving a image




When getting a satisfying image, press  to save the current image. At the same time, the current stored image serial number automatically displays at left top corner of the image area. Such as "SAVING 05". When the image's saving is finished, the serial number will disappear. After the saving, the system enters a freeze mode,

Press  key to return to real time state.



The equipment can store 32 images at most, and these images will be automatically numbered by their saving order.

For example: The number 01- 20 images has already been saved, press  key, then the currently saved image will be numbered 21 When the memory is full (I.e.

Already 32 images in it), then pressing  key, the screen will display as below:



Storage is Full.Erase
No.01?
1. Yse 2. No.

Tips replace the stored image number 01 press key  to replace it with current image press key  to exit this operation.


Select "2 No" to abandon saving image. then "Save" again, the system will remind you if you want to cover image number 02, reason by analogy in turn.

Tips:


When the image memory is full, if some saved images are reloaded, then processing the saving operation, the system will remind you if you want to delete this loaded image and save the current there from this serial number.

4.8.2 View a saved image



At real time or freeze, press "  "key to display the following dialog box:


Please put in image number:






Please Put in image number :
|

Caddo 12B

Input the saved image number, e.g. 01, then presses any key (except numbers and



letters. Press , backspace key to clear the mistakes and input again to load the first image). On the screen left bottom, 01/32 means: 01 is the serial number and 32 is the total amount of saved images in the memory. Press   0 to load other images.

Press  key to return to real time or freeze state. Repeat the above steps to load other Images.

Tips :





- **You can not load any image from any empty memory area.**
- **The input number should be from 01 to 32, when the input number is more than 32, the system will not take any operation.**

4.8.3Review storyboard mode

After starting up, system has entered the real time state, first collect cine loop images it will last for 30 seconds.



Switch the equipment into freeze mode. Press  key to review images in circulation.

During review, press  or to enter manual review mode press  to play increases Image by frame, and press  to play decrease by frame. And press  key to return to auto- review mode.

When review at "B", "B/B" mode, switch between "B/B" and "B" to play the images in difference windows.

To exit cine loop, press  key.

Tips:

1. If changing the image scanning mode, probe or display mode, cine loop can not be processed immediately. Normally cine loop can be done 30 seconds after scanning.
2. At cine loop mode, roll trackball can play cine by single frame. When roll trackball up, cine is played from small number to big, roll it down, the reverse. Press "Cine loop" key to go on auto playback.

Caddo 12B

3. Replace probe or change the scan angle, the cine loop frame may change. The standard results depend on the display data on the left bottom of the screen.

4.9 Obstetrics calculation

The equipment can calculate gestation age, expected date of delivery according to 7 parameters such as BPD, FL,

AC, HC, CRL, GS, first day of the last menses. There in BPD, AC, HC can be used to calculate fetus weight.


Press Key **OB Table** to display obstetrics table menu as given below:



TOKYO
1. BPD-FW
2. FL
3. AC-FW
4. HC-FW
5. AC-FW
6. GS
7. LMP
8. Exit

Figure 4-9. Obstetrics table menu sketch

Press **1** to select "1. BPD-FW", then measure the distance according to distance measurement method, the corresponding obstetrical table data display on screen right. As below:



BPD-061mm
G.A = 24W.2D
EDD= 06-03-27
FW = 0490g

BPD- Biparietal diameter

G.A - Gestation sac

EDD-Expected date of delivery

FW - Fetus weight

Press **2** to select "2. FL", then measure the distance according to distance measurement method, the corresponding obstetrical table data display on screen right. As below:

FL = 051 mm
G.A = 28W.0D
EDD= 06-03-27

FL-Femur length.

Press 3 select "3. AC FW", then measure the distances (two distances need measuring) according to distance measurement method, and press Measure.
+. key again, the corresponding obstetrical table data display on

+D1 = 044mm
×D2 = 086mm
AC = 204mm
G.A = 24W.4D
EDD = 06-04-21
FW = 0.68 Kg.

D1--distance 1

D2--distance 2

AC-Abdomeninal circumference

Press 4 to select "4. HC-FW", then measure the distances (two distances need measuring) according to distance measurement method and press Measure.
+. key again, the corresponding obstetrical table data display on screen right. As below:

+D1 = 136mm
×D2 = 121mm
HC = 256mm
G.A = 24W.4D
EDD = 06-04-07
FW = 0.85Kg.

Caddo 12B

HC-Head circumference

Press **5** to select "5. CRL", then measure the distance according to distance measurement method, the corresponding obstetrical table data display on screen right. As below:

CRL= 032mm
G.A = Low.4D
EDD = 06-07-29

CRL-Crown rump length

Press **6** to select "6. GS", then measure the distance according to distance measurement method, the corresponding obstetrical table data display on screen right. As below:

GS=035mm
G.A=08W.0D
EDD=06-08-17

GS-Gestation sac

Press **7** to select "7. First day of last menses", a reminder displays:

Put in the first date of the last menses:
MM-DD
I

This reminder ask to put in the first date of the pregnant woman's last menses in the date format MM-DD. For example: May II, put in 05-11. When the correct date is put, measuring result will displays at "EDD" on screen right in the format of YY-MM-DD, for example 07-02-21 means that the expected date of delivery is on February 2-2007.

Note

In expected date of delivery calculation, the system date must be correct. The default standard pregnancy period in the system is 40 weeks. In last menses method measuring, if the time interval between the input date and the system date exceeds 40 weeks, the inward date will not be accepted and need to be put in again. Under the condition of item 1, if the inward date is bigger than the current system date, it will be regarded as the date of the last year.

After measurement, press "  , key to clear the screen for the next measurement. Otherwise.

Press 8 to select "8. Exit" to exit OB table.

4.10 Report

The information of patient, diagnostic comments, measuring results as well as hospital, date, time, doctor will be saved in the report. System will generate a routine report page automatically according to different measurement items. Therein, distance, circumference, area, volume, etc. will be saved in Routine Report measuring results such as gestation age, expected date of delivery, fetus weight, etc. will be saved in Obstetrics Report At B, BB, BM, 4B mode, when fulfilling patient information comment, diagnoses, and foci distance, circumference, area, volume

measurements, keep the image in freeze state, press R
Report key to switch to Routine Report interface as given bellow:

HOSP:		DOCT:		2007-07-06	
NAME:		AGE:	SEX:	10:46:14	
+D1=		+C1= 223mm			
+D2=		+A1= 3556mm ²			
+D3=		XC2= 144mm			
+D4=		XA2= 1632mm ²			
D1/D2=		C1/C2= 154%			
D3/D4=		A1/A2= 217%			
Um1=		Um1= 92cm ³			
COMMENT:					

Figure 4-10. Routine report sketch

At B, BB 4B mode, when fulfilling patient information comment, diagnoses, gestation age, and expected date of delivery, Fetus weight measurements, keep the

image in freeze state, press R
Report key to display Obstetrics Report inter face as given bellow:

HOSP:				DOCT:			
NAME:				AGE:		SEX:	
1 2 3 AVG				G.A:		EDD:	
BPD:				FL:		FW:	
AC:				HC:		CRL:	
GS:							
FW(BPD.AC):							
RATIO: HC/AC=				FL/AC =			
				FL/BPD=			
COMMENT:							

Figure 4-11. Obstetrics Report sketch

Tips :

At report interface, press Automation
A B C key to add comments but other content can not be amended. At most three times Biparietal diameter (BPD) , femur length (FL), Crown rump length (CRL), Gestation sac (GS) can be measured at most three times, 1,2,3 indicate the measuring times, A VG indicates the average value of these three measurements.

1. Fetus weight FW (BPD, AC) is acquired from calculation on Biparietal diameter (BPD) and circumference of abdomen (A C).

2. Press R
Report key repeatedly to switch between diagnoses interface and report interface.

4.11 Image processing

Press Image
Processing key to display the image processing menu (as given below). It provides image positive-negative flip, GAMA calibration, statistic histogram, edge gain and other image processing functions.

1. Posi-Neg
 2. GAMA calibration
 3. Histogram
 4. Edge enhancement
 5. Exit

Press 1 to select " 1. Posi - Nega" to change the image polarity as given below:







Figure 4-12. polarity sketch

Press **2** to select "2. GAMA calibration" to lighten up "GAMA" on the up right of the screen. There are 4 levels G1, G2, G3 and OFF, press   to select them circularly.

Tips:

This function is set for video output and external devices (e.g. video printer) adjustment and is recommended to set at "OFF" during operation.

Press **3** to select "3. Histogram" to process histogram statistics as the following procedure:

Ai freeze, press "Image Processing" key to display the menu, select "3. Histogram" to call out a rectangle window on the image, press "Measure. 




Figure 4-13. Histogram sketch

Horizontal ordinate indicates gray scale degree, vertical ordinate indicates number.

Caddo 12B

PT indicates total quantity of pixel dots in the rectangular window.

Gm indicates the corresponding gray scale of the curve at the vertical peak.

Pm indicate the quantity of pixel dots at Gm gray scale.

From the above figure, in the rectangular area, the total number of pixel dots is 10000. At gray scale 16, there are

461 dots, the most image pixel dots.

Press **4** to select "4. Edge enhancement" to lighten up the "EDGEEN" on the up right screen. Then edge enhancement is adjustable. There are 9 levels: OFF, EE1-EE8 which can be adjusted by pressing **↑** **↓**.

Tips:

This function is a special item and recommended to set at "OFF" during operation.

Press **5** to select "5. Exit" to exit the menu.

4.12 Per rectum examination procedure

1. Empty the rectum before inspection, after cleansing enema, the patient take left lateral position, put some coupling gel around the endo-rectal probe tip and put on the rubber sleeve (one-time use).
2. During inspection, ask the patient to release the anus or take defecating, then breathe deeply to cooperate doctor's inspection. The first is anus digital examination to know the position, size, texture of the lesion and see if it is painful when pressing it and if there is blood on the finger.
3. Put the probe into anus about 2-3cm slowly, take a close look at the image, then go deep into the recta gradually and check each position. At the same time, press gut gently by putting hand on pubes at lower rump to enhance image definition.

4.13 Per vaginam examination procedure

1. Preparation before inspection: It is different from ultrasonic examination per abdomen, no need to suffuse the bladder. But please tell the patient that it is needed to put the probe into her vagina. Explain the course and advantages of per vaginam ultrasonic examination and release her tension.
2. Patient posture: Take the posture of bladder lithotripsy.
3. Scanning method: Put some coupling gel into the disposable sterilizing plastic jacket or a condom and then cover it on probe, apply some coupling gel on the cover and separate the nymphae with left hand to expose vaginal orifice, hold the probe with right hand and put it into the vagina gently can check. The position of probe in vagina will influence the image definition. In general, put the probe closely against vagina vault or cervix. Sometimes it is depending, so move or rotate the probe for a better known of the cavum pelvis. For the

moment, there is no hard and fast rule for image position. The near field can be the upper as well as the lower. Please identify the image direction. As same as per abdomen ultrasonic examination, first take vertical section scanning, and check the uterine position, size, outline and each stool section. muscle wall and endomembrane echo state as well as cervix vertical section. Wiggle the probe right and left to display the two ovaries' image. Then contra-rotate probe 90° to scan the uterine and two appendices cross section. Since most of per vaginam ultrasonic examinations are taken on average examining table instead of special gynecologic examining table, probe operation might be influenced. Some measures can be taken, for example, for some fatty woman or with anteposition of uterus, ask her to hold fists and put it under the hips to raise the rump, which is better for probe operation and complete observation. Sometimes, the target organs are in a higher position, press the abdominal wall gently to make them approach the probe.

4. Limitation: Endo-vaginal probe adopts high frequency probe mostly. The focal area is within 10cm. Cysts out the range of 10cm can not be observed well thus abdomen examination is necessary to cover the shortage.

4.14 Image print

Connect the VIDEO IN port of the video printer to the VIDEO OUT port on rear panel of the equipment, then operate according to video printer operation manual.

4.15 Image upload to computer

Connect the USB communication port of the equipment to the computer's USB port.

The high speed USB2.0 mouthpiece can upload images to the computer at current time.

In the accompanying disc, there are USB device driver. image preview program Amcap and SDK documents for relevant redevelopment.

Install the USB device driver, open Amcap.EXE, then you can view images according to the instruction files in the disc. Using SDK documents in the disc, you can process redevelopment as required.

4.16 Shut down the equipment

Turn power off.

Tips: It would be better to unplug the A C supply if do not use the equipment for a long time.

Note

Must not plug or unplug the power plug when the equipment is not shut down If it needs to switch on the equipment immediately after just shutdown, please wait for 2 or 3 minutes to avoid breaking the equipment.

Chapter Five Transportation and Storage

5.1 Environmental Requirements of Transportation and Storage

Environmental temperature: - 40°C ~ + 55°C

Relative humidity: 10% ~ 100%

Atmospheric pressure: 50KPa ~ 106KPa

5.2 Transportation

All marks printed on the packing box are in accordance with the requirements of GB/T191-2008 <Packing, Storage: and Transportation Marks>. Simple shockproof installations are made in the packing box, which adapt to voyage, railway, and highway and steamship transportation. Avoid rain, inversion and impact.

5.3 Storage

- When the device is stored for more than 6 months, take out the device from the packing box. Electrify the device for 4 hours, and then put the device into the packing box according to the instructions way on the packing box. Do not crisscross/piling devices or keep the device to floor, wall and ceiling.
- Keep the repository ventilative. Avoid direct sunshine and caustic gas.

Chapter Six Check and Maintenance

6.1 Check

Check the equipment's power line and probe cable regularity, if found any disrepair, Broken phenomena, taboo and replacing it immediately.

6.2 Service life

service life is 6 years.

Note

The Discard the device according to local law.

Do not discard it mixing with other household garbage


Warning

The manufacturer shall not assume the responsibility of risks caused by using the device beyond its service life.

Main unit maintenance

Instrumentation environment should accord with "2.1 operation environmental requirement". If device enclosure needs cleaning, shutdown the device first and then wipe with alcohol sponges. Device should not turn on and off frequently. When the device does not work for a long time, pack the device according to the instructions on the packing. Store it properly in the warehouse. The storage environment should accord with "5.1 Transportation and storage environmental requirements".

6.4 Probe maintenance

Probe is an expensive part and frangible. Never hit it or drop it on floor. When not working, place the probe in the box and press  key to set the probe at freeze state.

Please use medical ultrasound coupling gel during diagnoses. The water prevention level is IPX7, so do not let water exceed the acoustic window (see figure 6-1, 6-2, 6-3, 6-4,6-5). Do daily inspection on the probe enclosure to see if it is cracked and avoid liquid leakage to spoil the inner components.



Figure 6-1 . C3-1/60R/3.5MHz Convex probe water-proof sketch



Figure 6-2. EC1-1 Endocavity probe water-proof sketch




Figure 6-3. L3-1/7.5MHz HF linear probe water proof sketch



Figure 6.4. C1-6/20R/5.0MHz convex probe water- tight sketch



Figure 6-5. EL3-1/7.5MHz rectal probe water-proof sketch

Caddo 12B

Note

Do not press the probe on the patient too long to avoid discomfort.

6.5 Cleansing

- I. When the enclosure needs cleaning, wipe it with soft dry cloths and then wipe gently with sponge dipped with 75% medical alcohol.
- II. When the inner part of the equipment needs cleaning, power off the equipment first and open the enclosure and vacuum it.

Warning

To prevent accidents, please power off the equipment when cleaning it or the probe.

Caution

Please refer to instructions prescribed by the manufacturer closely when using detergents.

Be careful with cleaning of the display, because it is very easy to scratch and spoil. Please wipe it with dry soft cloth.

Please do not clean the inner part of the device.

Please do not place the device in liquid.

Do not leave any detergent on the device surface.

Though there will be no chemical reaction between the device enclosure and most of those detergents, We still suggest no detergent in cleaning lest the device surface is spoiled.

Warning

Must not use extender, ethylene oxide or any other organic solvent which tend to deface the probe's protective foil.

Must not place the probe in liquid or detergent.

Must keep the equipment or probe from any type of liquid's infiltration. Must not clean device or probe by airing or heating.

6.6 Correct usage of probe

In order to extend probe's service life and obtain optimum performance, follow the instructions as below:

1. Periodic inspection on probe cable, socket and acoustic window.
2. Shutdown the device first and then connect or disconnect the probe.
3. Do not drop probe or flint body, and never hit the probe acoustic window, otherwise probe should be damaged.
4. Put the probe in the probe box when it is not in use.

Caddo 12B

5. Never heat the probe.
6. Never bend or pull probe cable, otherwise the internal connection should be broken.
7. Use coupling gel only on probe header and then clean probe.
8. Inspect probe acoustical window, enclosure and cable seriously after probe cleaning. Never use the probe again if any crack or breakage is found.
9. Most of Probe could not be disinfected. If antiseptis needed, please use the asepsis probe head gear.
10. Suggesting to use asepsis probe headgear and asepsis gel during ail antrium probe clinic application. Do not open the probe headgear and gel package.

Warning

Probe headgear may have caoutchouc latex, it may cause some people hype rsusceptibility reaction

Warning

Can not work with high frequency surgery equipment, the applicate part needs to prevent to harmed by the high frequency equipment.

6.7.1 Battery information

1. The equipment is fitted with rechargeable Li-ion battery.
2. For optimum efficiency, the new battery must be charged and discharged (regular service, not enforced discharging) two or three rounds completely.
3. The battery can be charged and discharged for hundred times, but it will be worn-out. When the work time shortens apparently, please replace it with a new one.
4. Be sure to use electricity charger appointed (Le. AC adapter) to charge the battery. Do not connect the battery to the electricity charger (AC adapter) when charging is not needed. Do not connect the battery to the electricity charger (AC adapter) longer than 10 hours otherwise the battery life may shorten. The fully charged battery will discharge by itself if it is long-time out of use.
5. Extreme environment temperature (overcooling or overheating) will influence battery charging effect. Must not charge the battery near the ignition source or under extreme hot condition! Do not use or store battery near source of heat (such as fire or heater)! If find the battery is leaking or smelling, move the battery away from the naked flame immediately.
6. Don't go on using non-serviceable battery and electricity charger (AC adapter).
7. Don't try demounting battery.
8. Don't short circuit the battery.
9. Do not throw the battery into the fire or heat it, otherwise it would trigger an explosion.

Caddo 12B

10. Do not submerge or wet the battery.
11. Do not incorrectly connect the positive and negative polarity.
12. Do not directly connect the battery to wall outlet or car-lit socket.
13. Must not short circuit the positive and negative polarity of the battery with lead or other metal objects. Must not transport nor store the battery with necklace, hair pin or other metal objects.
14. Must not pierce battery shell with nail or other sharp objects, must not hammer nor step on the battery.
15. Must not hit, cast the battery and avoid mechanical shock on it.
16. Must not bend the battery terminals.
17. Must not decompose the battery in any way.
18. Must not place the battery in microwave oven or pressure vessels.
19. Must not combine the battery with primary battery (such as dry battery) or battery with different capacity, models and types.
20. Do not use the battery if it is smelling, heating, straining, discoloured or with other abnormal phenomena and remove it from the current consumer or electricity charger immediately and stop using it any longer.
21. Do deal carefully with the discarded battery according to local related waste handling regulations.

6.8 Instrument test and calibration

1. Check the leakage current of the device annually referring to the following data.

Test Items			Standard Requirements
Continuous leakage current under normal temperature (□A)	Terra creepage current	Normal	≤ 500
		Single malfunction	≤ 1000
	Crust creepage current	Normal	≤ 100
		Single malfunction	≤ 500
	Sufferer creepage current	Normal	≤ 100
		Single malfunction	≤ 500
Dielectric endurance under normal temperature (V)	A-a1		1500V /1 min No flashover No hit pierce
	A-a2		4000V/lmin No flashover No hit pierce
	B-a		4000V/lmin No flashover No hit pierce

2. Test the software of obstetrics area and circumference measurement please refer to Appendix C for detailed data.

Chapter Seven Malfunction Examination and Troubleshooting

7.1 Examination

Examine whether the supply power is in normal state. The power line of the main frame has been properly connected and inserted in the electric socket.

Examine whether the probe has been connected to the device correctly.

7.2 Troubleshooting

Troubleshooting (see the following diagram)

No.	Malfunction phenomena	Troubleshooting
1	Power adapter indicate light does Not work	<ol style="list-style-type: none">1. Check the supply power2. Check the adapter power cable and connector plug.
2	Main unit power light doe not work	Check the main unit supply power Cable and whether the connection Is all right.
3	The screen display disconnected stripe disturb, snowflake shape disturb.	<ol style="list-style-type: none">1. Check the adapter supply power, the Supply power effected by other equipment strike fire distur2. Environment check, effected by Around magnetic field disturb.3. Check the probe connector plugs, and connector jack.
4	Screen image display not clear	<ol style="list-style-type: none">1. Adjust overall gain2. Adjust fake colour control key

Appendix A Sound Output Specifications Table

TABLE: Summary of the acoustic quantities 20R/5.0MHz convex probe B-Mode

No.	Test Item	Requirements and Clauses	Results	Verdict	Remarks
1	Maximum power(m W)	Maximum temporal-average power output. IEC 1157 4.2.2a)	3.9	-	Test frequency : 5.0MHz
2	P-(MPa)	Peak-negative acoustic pressure. IEC 1157 4.2.2b)	0.77	-	
3	Iob(m W /cm2)	Output beam intensity. IEC 1157 4.2.2c)	0.89	-	
4	Ispta(m W /cm2)	Soatial-peak temporal average derived intensity. IEC 1157 4.2.2d)	2.30	-	
5	System settings	Ultrasound instrument console settings. IEC 1157 4.2.2e)	F3 * 1.0	-	
6	Lp(mm)	Distance form the transducer output face to the point of maximum pulse-pressure-squared integral. IEC 1157 4.2.21)	25	-	
7	Wpb6(mm)	-6dB pulse-width. IEC 1157 4.2.2g)	(II): 1.99 3.31 -		
8	pr(KHz)	Pulse repetition rate. IEC 1157 4.2.2h)	33.33	-	
9	srr(Hz)	Scan repetition rate. IEC 1157 4.2.2h)	Not applicable	-	

Caddo 12B

10	Output beam dimensions(cm2)	<p>Dimensions parallel (II)and perpendicular J:L) to the reference direction shall be specific.</p> <p>For scanning modes, these shall refer to the central scan line only. In many cases, especially contact systems, these dimensions may be taken as the geometrical dimensions of the ultrasonic transducer or ultrasonic transducer element group.</p> <p>IEC 1157 4.2.2i)</p>	4.37		
11	fawf(MHz)	<p>Arithmetic-mean acoustic-working frequency.</p> <p>IEC 1157 4.2.2j)</p>	3.87	-	-
12	APF	<p>Acoustic power up fraction .</p> <p>IEC 1157 4.2.2n</p>	B-Mode freeze	-	

Caddo 12B

113	Power-up mode	In system in which the user defines the power- up mode, this shall be stated as either "user defined" or "not applicable" (n/a). IEC 1157 4.2.21)	B-Mode Freeze	-	
1.4	AIF	Acoustic initialization fraction. IEC 1157 4.2.2m)	100%	-	
15	Initialization mode	If appropriate. in systems in which the user defined the initialization mode, this shall be stated as either "user defined " or "not applicable "(n/a) IEC 1157 4.2.2n)	B-Mode	-	
16	Accoustic output freeze	If the system has acoustic output freeze then this shall be stated as "yes", otherwise it shall be stated as "no". IEC 1157 4.2.20)	Yes	-	

TABLE: Summary of the acoustic quantities 20R/5.0MHz convex probe M-Mode

No.	Test Item	Requirements and Clauses	Results	Verdict	Remarks
I	Maximum power(mW)	Maximum temporal-average power output. IEC 1157 4.2.2a)	0.070	-	Test frequency: 5.0MHz
2	P-(MPa)	Peak-negative acoustic pressure. IEC 1157 4.2.2b)	0.71	-	
3	lob(mW/cm2)	Output beam intensity. IEC 1157 4.2.2c)	0.039	-	
4	Ispta(mW/cm2)	Soatial-peak temporal average derived	1.26	-	

Caddo 12B

		intensity. IEC 1157 4.2.2d)			
5	System settings	Ultrasound instrument console settings. IEC 1157 4.2.2e)	F3 * 1.0	-	
6	Lp(mm)	Distance form the transducer output face to the point of maximum pulse-pressure-squared integral. IEC 1157 4.2.2t)	25	-	
7	Wpb6(mm)	-6dB pulse-width. IEC 1157 4.2.2g)	(II):3.30 1 :2.91	-	
8	pr(KHz)	Pulse repetition rate. IEC 1157 4.2.2h)	0.148	-	
9	srr(Hz)	Scan repetition rate. IEC 1157 4.2.2h)	Not applicable	-	
10	Output beam dimensions(cm2)	Dimensions parallel (II) and perpendicular }L) to the reference direction shall be specific . For scanning modes, these shall refer to the central scan line only. In many cases, especially contact systems, these dimensions may be taken as the geometrical dimensions of the ultrasonic transducer or ultrasonic transducer element group.	1.79	-	

Caddo 12B

		IEC 1157 4.2.2i)			
11	fawt<MHz)	Arithmetic-mean acoustic- working frequency. IEC 1157 4.2.2j)	4.22	-	
12	APF	Acoustic power-up fraction. IEC 1157 4.2.2k)	0%	-	
13	Power-up mode	In system in which the user defines the power-up mode, this shall be stated as either "user defined" or "not applicable" (n/a). IEC 1157 4.2.2l)	B-Mode Freeze	-	
14	AIF	Acoustic initialization fraction. IEC 1157 4.2.2m)	108.45%	-	
15	Initialization mode	If appropriate, in systems in which the user defined the initialization mode, this shall be stated as either "user defined" or "not applicable "(n/a) IEC 1157 4.2.2n)	B-Mode	-	
16	Accoustic output freeze	If the system has acoustic output freeze then this shall be stated as "yes", otherwise it shall be stated as "no". IEC 1157 4.2.2o)	Yes	-	

TABLE: Summary of the acoustic quantities 60R/3.5MHz convex probe B-Mode

Index	MI	TIS	TIS	TIS	TIB	TIB	TIC
Mode		- Scanning	Non-scanning	Non-scanning	Scanning	Non-scanning	-
Maximum index value	0.50	0.061			0.061		
Acoustic working frequency (MHz)	2.91	2.91			2.91		
Output power (mw)							
Bounded output Power (mw)		4.42			4.42		
Attenuated output Power (mw)							
Spatial-peak temporal-average intensity(mw/cm2)							
Attenuated spatial-peak temporal-average intensity(mw/cm2)							
Pulse-pressure squared integral(mj/cm2)	0.046						
Attenuated Pulse-pressure squared integral(mj/cm2)	0.031						
Peak-rare fractional acoustic pressure (MPa)	1.03						

Caddo 12B

Attenuated peak-rarefactional acoustic pressure (MPa)	0.85						
-1 2dB output beam area (cm2)							
Equivalent aperture diameter (cm)							
Depth for TIS (cm)							
Depth for TIB (cm)							
Depth at max. attenuated pulse-intensity integral (cm)	5.55						
Supplementary information:							

Caddo 12B

TABLE: Summary of the acoustic quantities 60R/3.5MHz convex probe (B+M)-Mode

Index	MI	TIS	TIS	TIS	TIB	TIB	TIC
Mode		Scanning	Non-scanning Aaprt: Icm2	Non-scanning	Scanning	Non-scanning	
Maximum index value	0.50	0.061			0.061		
Acoustic working frequency(MHz)	2.91	2.91		3.34	2.91	3.34	
Output power (mw)				0.29		0.29	
Bounded output Power(mw)		4.42			4.42		
Attenuated output Power (mw)				0.20		0.20	
Spatial-peak temporal-average intensity(mw/cm2)				2.39		2.39	
Attenuated spatial-peak temporal-average intensity(mw/cm2)				0.45		0.45	
Pulse-pressure squared integral(mj/cm2)	0.046					0.021	
Attenuated Pulse-pressure squared integral(mj/cm2)	0.031					0.0049	
Peak-rarefactional acoustic pressure (MPa)	1.03						

Caddo 12B

Attenuated peak-Rarefactional acoustic Pressure (MPa)	0.85						
-I 2dB output beam area (cm ²)				1.33			
Equivalent aperture diameter (cm)				1.30			
Depth for TIS (cm)				1.96			
Depth for TIB (cm)						1.70	
Depth at max. attenuated pulse-intensity integral (cm)	5.55						
Supplementary information:							

TABLE: Summary of the acoustic quantities 20R/5.0MHz convex probe B-Mode

Index	MI	TIS	TIS	TIS	TIB	TIB	TIC
Mode	-	Scanning	Non-scanning	Non-scanning	Scanning	Non-scanning	-
			Aaprt lcm ²	Aaprt> lcm ²			
Maximum index value	0.35	0.019			0.019		
Acoustic working frequency (MHz)	3.88	3.88			3.88		
Output power (mw)							
Bounded output Power (mw)		1.05			1.05		
Attenuated output Power							

Caddo 12B

(mw)							
Spatial-peak temporal- Average intensity(mw/cm2)							
Attenuated spatial-peak temporal-average intensity(mw/cm2)							
Pulse-pressure squared integral(mj/cm2)	0.0095						
Attenuated Pulse-pressure squared integral(mj/cm2)	0.0077						
Peak-rarefactional acoustic pressure (MPa)	0.77						
Attenuated peak- rarefactional acoustic pressure (MPa)	0.69						
-1 2dB output beam area (cm2)							
Equivalent aperture diameter (cm)							
Depth for TIS (cm)							
Depth for TIB (cm)							
Depth at max. attenuated pulse-intensity integral (cm)	3.05						
Supplementary information:							

TABLE: Summary of the acoustic quantities 20R/5.0MHz convex probe M-Mode

Index	MI	TIS	TIS	TIS	TIB	TIB	TIC
-------	----	-----	-----	-----	-----	-----	-----

Caddo 12B

Mode		Scanning	Non-scanning	Non-scanning	Scanning	Non-scanning	
Maximum index value	0.24		0.0014			0.0030	
Acoustic working frequency (MHz)	4.23		4.23			4.23	
Output power (mw)			0.071			0.071	
Bounded output Power (mw)							
Attenuated output Power (mw)						0.05	
Spatial-peak temporal-average intensity(mw/cm ²)						1.27	
Attenuated spatial-peak temporal-average intensity(mw/cm ²)						0.52	
Pulse-pressure squared integral(mj/cm ²)	0.0085					0.0085	
Attenuated Pulse-pressure squared integral(mj/cm ²)	0.0040					0.0040	
Peak-rarefactional acoustic pressure (MPa)	0.69						
Attenuated peak - rarefactional acoustic Pressure (MPa)	0.47						
-1 2dB output beam area (cm ²)			0.22				

Caddo 12B

Equivalent aperture Diameter (cm)							
Depth for TIS (cm)							
Depth for TIB (cm)						2.30	
Depth at max. attenuated pulse-intensity integral (cm)	2.30						
Supplementary information:							

TABLE: Summary of the acoustic quantities 13R16.5MHz antrum probe B-Mode

Index	MI	TIS	TIS	TIS	TIB	TIB	TIC
Mode	-	Scanning	Non-scanning	Non-scanning	Scanning	Non-scanning	-
			Aaprt Icm2	Aaprt > Icm2			
Maximum index value	0.49	0.021			0.021		
Acoustic working Frequency (MHz)	4.36	4.36			4.36		
Output power (mw)							
Bounded output Power (mw)		1.00			1.00		
Attenuated output Power (mw)							
Spatial-peak temporal-Average intensity(mw/cm2)							
Attenuated spatial-peak temporal-average intensity(mw/cm2)							

Caddo 12B

Pulse-pressure squared integral(mj/cm2)	0.0097						
Attenuated Pulse-pressure squared integral(mj/cm2)	0.0079						
Peak-rarefactional acoustic pressure (MPa)	1.08						
Attenuated peak-rarefactional acoustic pressure (MPa)	0.98						
-I 2dB output beam area (cm2)							
Equivalent aperture diameter (cm)							
Depth for TIS (cm)							
Depth for TIB (cm)							
Depth at max. attenuated pulse-intensity integral (cm)	2.40						
Supplementary information:							

TABLE: Summary of the acoustic quantities 13R/6.5MHz antrum probe M-Mode

Index	MI	TIS	TIS	TIS	TIB	TIB	TIC
Mode	Scanning		Non-scanning	Non-scanning	Scanning	Non-scanning	
Maximum index value	0.24		0.0012			0.0015	

Caddo 12B

Acoustic working frequency (MHz)	4.62		4.62			4.62	
Output power (mw)			0.056			0.056	
Bounded output Power (mw)							
Attenuated output Power (mw)						0.026	
Spatial-peak temporal-average intensity(mw/cm2)						0.48	
Attenuated spatial-peak temporal-average intensity(mw/cm2)						0.22	
Pulse-pressure squared integral(mj/cm2)	0.0029					0.0029	
Attenuated Pulse-pressure squared integral(mj/cm2)	0.0014					0.0014	
Peak-rarefactional acoustic pressure (MPa)	0.70						
Attenuated peak-rarefactional acoustic pressure (MPa)	0.48						
-1 2dB output beam area (cm2)			0.15				

Caddo 12B

Equivalent aperture diameter (cm)							
Depth for TIS (cm)							
Depth for TIB (cm)						2.35	
Depth at max. attenuated pulse-intensity integral (cm)	2.35						
Supplementary information:							

TABLE: Summary of the acoustic quantities 13R/6.5MHz antrum probe (B+M)-Mode

Index	MI	TIS	TIS	TIS	TIB	TIB	TIC
Mode	-	Scanning	Non-scanning Aaprt Icm2	Non-scanning Aaprt> Icm2	Scanning	Non-scanning	-
Maximum index value	0.49	0.022			0.022		
Acoustic working frequency (MHz)	4.36	4.36	4.62		4.36	4.62	
Output power (mw)			0.056			0.056	
Bounded output Power (mw)		1.00			1.00		
Attenuated output Power (mw)						0.026	
Spatial-peak temporal-average intensity(mw/cm2)						0.48	

Caddo 12B

Attenuated spatial-peak temporal-average intensity(mw/cm2)						0.22	
Pulse-pressure squared integral(mj/cm2)	0.0097					0.0029	
Attenuated Pulse-pressure squared integral(mj/cm2)	0.0079					0.0014	
Peak-rarefactional acoustic pressure (MPa)	1.08						
Attenuated peak - Rarefactional acoustic pressure (MPa)	0.98						
-I 2dB output beam area (cm2)			0.15				
Equivalent aperture diameter (cm)							
Depth for TIS (cm)							
Depth for TIB (cm)						2.35	
Depth at max. attenuated pulse-intensity integral (cm)	2.40						
Supplementary information:							

Caddo 12B

TABLE: Summary of the acoustic quantities 7.5MHz High frequency linear probe B-Mode

Index	MI	TIS	TIS	TIS	TIS	TIB	TIC
Mode		Scanning	Non-scanning Aaprt lcm2	Non-scanning Aaprt> lcm2	Scanning	Non-scanning	s
Maximum index value	0.83	0.135			0.135		
Acoustic working frequency (MHz)	5.53	5.53			5.53		
Output power (mw)							
Bounded output Power (mw)		5.12			5.12		
Attenuated output Power (mw)							
Spatial-peak temporal-average intensity(mw/cm2)							
Attenuated spatial-peak temporal-average intensity(mw/cm2)							
Pulse-pressure squared integral(mjlc2)	0.16						
Attenuated Pulse-pressure squared integral(mjlc2)	0.059						
Peak-rarefactional acoustic pressure (MPa)	2.72						

Caddo 12B

Attenuated peak- rarefactional acoustic pressure (MPa)	1.66						
-I 2dB output beam area (em2)							
Equivalent aperture diameter (em)							
Depth for TIS (em)							
Depth for TIS (cm)							
Depth at max. attenuated pulse-intensity integral (cm)	2.60						
Supplementary information:							

Caddo 12B

TABLE: Summary of the acoustic quantities 7.5MHz High frequency linear probe M-Mode

Index	MI	TIS	TIS	TIS	TIB	TIB	TIC
Mode	-	Scanning	Non-scanning Aaprt lcm2	Non-scanning Aaprt> lcm2	Scanning	Non-scanning	-
Maximum index value	0.98		0.0014			0.0052	
Acoustic working Frequency (MHz)	5.37		5.37			5.37	
Output power (mw)			0.060			0.060	
Bounded output Power (mw)							
Attenuated output Power (mw)						0.021	
Spatial-peak temporal-average intensity(mw/cm2)						31.34	
Attenuated spatial-peak temporal-average intensity(mw/cm2)						11.95	
Pulse-pressure squared integral(mjlcms2)	0.19					0.19	
Attenuated Pulse-pressure squared integral(mjlcms2)	0.074					0.074	
Peak-rarefactional acoustic pressure (MPa)	2.83						

Caddo 12B

Attenuated peak - Rarefactional acoustic Pressure (MPa)	1.97						
-I 2dB output beam area (cm2)			0.61				
Equivalent aperture diameter (cm)							
Depth for TIS (cm)							
Depth for TIB (cm)						2.40	
Depth at max. attenuated pulse-intensity integral (cm)	2.55						
Supplementary information:							

TABLE: Summary of the acoustic quantities 7.5MHz High frequency linear probe (B+M)Mode

Index	MI	TIS	TIS	TIS	TIS	TIS	TIC
Mode	-	Scanning	Non-scanning Aaprt lcm2	Non-scanning Aaprt> lcm2	Scanning	Non-scanning	-
Maximum index value	0.98	0.136			0.140		
Acoustic working frequency (MHz)	5.37	5.53	5.37		5.53	5.37	
Output power (mw)			0.060			0.060	
Bounded output Power (mw)		5.12			5.12		

Caddo 12B

Attenuated output Power (mw)						0.021	
Spatial-peak temporal-average intensity(mw/cm2)						3 1.34	
Attenuated spatial-peak temporal-average intensity(mw/cm2)						11.95	
Pulse-pressure squared integral(mj/cm2)	0.19					0.19	
Attenuated Pulse-pressure squared integral(mj/cm2)	0.074					0.074	
Peak-rarefactional acoustic pressure (MPa)	2.83						
Attenuated peak - rarefactional acoustic pressure (MPa)	1.97						
-I 2dB output beam area (cm2)			0.61				
Equivalent aperture diameter (cm)							
Depth for TIS (cm)							
Depth for TIB (cm)						2.40	
Depth at max. attenuated pulse-intensity integral (cm)	2.55						
Supplementary information:							

Appendix B Gestational Table**Head circumference (HC) Table Set1: Hadlock, F.P., et al. Method**

HC mm	WEEKS w. d.	S.D ±days	HC mm	WEEKS w. d.	S.D ± days
56	12.0	8	210	23. 1	11
63	12. 3	8	217	23.5	11
70	12.6	8	224	24.3	15
77	13.2	8	231	25. 1	15
84	13.5	8	238	25.6	15
91	14. 1	8	245	26.4	15
98	14.4	8	252	27. 3	15
105	15.0	8	259	28. 1	15
112	15.3	8	266	29.0	15
119	15.6	8	273	29.6	15
126	16. 3	8	280	30.5	21
133	16.6	8	287	31. 4	21
140	17.2	8	294	32.4	21
147	17.6	8	301	33.3	21
154	18.3	11	308	34.3	21
161	18.6	11	315	35.2	19
168	19.4	11	322	36.3	19
175	20.0	11	329	37.3	19
182	20.4	11	336	38.4	19
189	21. 1	11	343	39.4	19
196	21. 6	11	350	40.5	19
203	22. 3	11			

Femur length (FL)**Table Set1: Tokyo University Method**

FL mm	WEEKS w. d.	S.D ± days	FL mm	WEEKS w. d.	S.D ± days
8	12.3	10	42	24.5	24
10	13.0	10	44	25.2	25
12	13.4	10	46	26.0	25
14	14. 1	10	48	27.0	25
16	14.5	10	50	28.0	25
18	15.2	10	52	29.0	30
20	16.0	10	54	29.5	30
22	16.4	10	56	30.2	30
24	17. 1	10	58	31. 3	32
26	17.6	10	60	33.0	38
28	18.4	14	62	34.0	42
30	19.2	17	64	35.0	46
32	20. 5	17	66	36.5	50
34	21. 5	18	68	38.0	57
36	22.3	19	70	40.0	64
38	23.0	21	72	40.2	64
40	24.0	22			

Femur length (FL)**Table Set2: Hadlock, F. P., et al. Method**

FL nun	WEEKS w. d.	S.D + days	FL nun	WEEKS w. d.	S.D +days
8	12.3	10	46	25.2	15
10	13.0	10	48	26.1	15
12	13.4	10	50	26.6	15
14	14.1	10	52	27.5	15
16	14.5	10	54	28.4	15
18	15.2	10	56	29.4	15
20	16.0	10	58	30.2	21
22	16.4	10	60	31.1	21
24	17.1	10	62	32.1	21
26	17.6	10	64	33.1	21
28	18.4	14	66	34.0	21
30	19.2	14	68	34.6	21
32	20.0	14	70	35.6	21
34	20.5	14	72	36.6	22
36	21.3	14	74	37.6	22
38	22.1	14	76	38.6	22
40	22.6	14	78	39.6	22
42	23.5	14	80	40.6	22
44	24.4	15	82	42.0	22

Crown-rump length (CRL) Table Set1: Tokyo University Method

CRL mm	WEEKS w. d.	S.D ± days	CRL mm	WEEKS w. d.	S.D ± days
6	6.3	7	54	12.4	7
8	6.6	7	56	12.5	7
10	7.1	7	58	13.0	7
12	7.4	7	60	13.1	7
14	7.6	7	62	13.2	7
16	8.1	7	64	13.3	7
18	8.4	7	66	13.4	7
20	9.1	7	68	13.5	7
22	9.2	7	70	13.6	7
24	9.4	7	72	14.0	8
26	9.6	7	74	14.1	8
28	10.2	7	76	14.2	8
30	10.3	7	78	14.3	8
32	10.4	7	80	14.4	8
34	10.6	7	82	14.5	8
36	11.0	7	84	14.6	8
38	11.1	7	86	15.0	8
40	11.3	7	88	15.1	14
42	11.4	7	90	15.2	14
44	11.6	7	92	15.3	14
46	12.0	7	94	15.4	14
48	12.1	7	96	15.5	14
50	12.2	7	98	15.6	14
52	12.3	7	100	16.0	14

Crown-rump length (CRL) Table Set2: Hadlock, F.P.,et al. Method

CRL mm	WEEKS w. d.	S.D +days	CRL mm	WEEKS w. d.	S.D + days
4	6. 1	3	68	13. 1	6
6	6. 3	3	70	13. 1	6
8	6.5	3	72	13. 3	6
10	7. 1	3	74	13.3	6
12	7.3	3	76	13.5	6
14	7.5	3	78	13.6	6
16	8.0	3	80	14.0	6
18	8.2	3	82	14. 1	6
20	8.4	4	84	14.2	6
22	8.6	4	86	14.3	6
24	9. 1	4	88	14.5	6
26	9. 3	4	90	14.6	6
28	9.4	4	92	15. 1	6
30	9.6	4	94	15.2	6
32	10. 1	5	96	15.3	6
34	10.2	5	98	15.4	6
36	10.4	5	100	15.6	6
38	10.5	5	102	16. 1	6
40	10.6	5	104	16.2	6
42	11. 1	5	106	16. 3	7
44	11. 1	5	108	16. 5	7
46	11. 3	5	110	16.6	7
48	11. 4	6	112	17. 1	7
50	11. 5	6	114	17.2	7
52	11. 6	6	116	17.3	7
54	12.0	6	118	17.5	7
56	12. 1	6	120	17.6	7
58	12.2	6			

Caddo 12B

60	12.3	6	
62	12.4	6	
64	12.6	6	
66	12.6	6	

Gestational Sac (GS)**Table Set1: Tokyo University Method**

GS nun	WEEKS w. d.	S.D + days	GS nun	WEEKS w. d.	S.D + days
10	4.0	7	42	9. 1	14
12	4. 1	7	44	9.3	14
14	4.4	7	46	9.4	14
16	5.0	8	48	10.0	15
18	5. 1	8	50	10. 1	15
20	5.4	8	52	10.3	15
22	6.0	11	54	10.4	15
24	6.1	11	56	10.6	15
26	6.6	12	58	11. 1	16
28	7. 1	12	60	11. 3	16
30	7.3	12	62	11. 4	16
32	7.4	12	64	11. 6	16
34	8.0	13	66	11. 8	16
36	8. 1	13	68	12. 1	17
38	8. 3	13			
40	8.6	13			

Gestational Sac (GS)**Table Set2: Hadlock, F.P., et al. Method**

GS mm	WEEKS w. d.	S.D ± days	GS mm	WEEKS w. d.	S.D ± days
6	5.0				
8	6.0				
10	6.0				
12	6.2				
14	6.4				
16	7.0				
18	7.2				
20	7.2				
22	7.4				
24	7.4				
26	8.0				
28	8.2				
30	8.4				
32	8.4				
34	8.6				
36	9.0				
38	9.3				
40	9.3				
42	9.5				
44	10.0				
46	10.0				
48	10.2				
50	10.4				
52	11.0				
54	11.2				
56	11.6				
58	12.0				

Caddo 12B

60	12.4	
62	13.0	
64	13.0	

www.hik-consulting.pl/edu

Table Set1: Tokyo University Method

BPD mm	WEEKS w. d.	S.D ± days	BPD mm	WEEKS w. d.	S.D ± days
16	11.3	7	56	23.0	11
18	11.6	7	58	23.5	11
20	12.0	7	60	24.2	12
22	12.4	7	62	25.0	12
24	13.0	7	64	25.6	12
26	13.6	7	66	26.3	13
28	14.2	7	68	27.3	13
30	14.6	7	70	28.0	13
32	15.2	7	72	29.0	14
34	16.0	8	74	29.5	14
36	16.3	8	76	30.1	15
38	17.0	8	78	31.1	16
40	17.5	8	80	32.1	16
42	18.2	9	82	33.0	18
44	19.0	9	84	34.0	20
46	19.5	10	86	35.5	25
48	20.2	10	88	37.0	25
50	21.0	10	90	39.0	25
52	21.4	10	92	42.0	25
54	22.2	10			

Table Set2: Hadlock, F.P., et al. Method

BPD mm	WEEKS w. d.	S.D ± days	BPD mm	WEEKS w. d.	S.D ±days
14	11. 6	7	60	24.4	16
16	12.2	9	62	25. 1	16
18	12.6	9	64	25.6	16
20	13. 1	9	66	26.4	16
22	13.4	9	68	27.3	16
24	14. 1	9	70	28. 1	16
26	14.4	9	72	28.6	16
28	15.0	9	74	29. 5	16
30	15.4	9	76	30.4	22
32	16.0	9	78	31. 2	22
34	16.4	9	80	32. 1	22
36	17.0	9	82	33.0	22
38	17.4	9	84	33.6	22
40	18. 1	14	86	34.5	22
42	18. 5	14	88	35.4	22
44	19.2	14	90	36.4	24
46	19.6	14	92	37. 3	24
48	20.4	14	94	38.2	24
50	21. 1	14	96	39. 1	24
52	21. 6	14	98	40. 1	24
54	22.3	14	100	41. 1	24
56	23. 1	14			
58	23.6	14			

Abdominal circumference (AC) Table Set1: Hadlock, F.P., et al. Method

AC mm	WEEKS w. d.	95% Conf. Limits Cd)	AC mm	WEEKS w. d.	95% Conf. Limits Cd)
100	15.4	13.7-17.5	235	27.5	25.5-29.9
105	16.1	14.2-18.0	240	28.1	26.0-30.4
110	16.4	14.6-18.4	245	28.5	26.5-30.9
115	16.6	15.0-18.8	250	29.1	27.0-31.4
120	17.2	15.4-19.2	255	29.5	27.5-31.9
125	17.6	15.9-19.7	260	30.1	27.1-33.1
130	18.1	16.2-20.2	265	30.4	27.6-33.6
135	18.4	16.6-20.6	270	31.1	28.1-34.1
140	19.1	17.1-21.1	275	31.4	28.6-34.6
145	19.4	17.5-21.5	280	32.1	29.1-35.1
150	20.0	18.0-22.0	285	32.4	29.6-35.6
155	20.3	18.4-22.4	290	33.1	30.1-36.1
160	20.6	18.8-22.8	295	33.4	30.6-36.6
165	21.2	19.3-23.3	300	34.1	31.1-37.1
170	21.5	19.7-23.7	305	34.4	31.6-37.6
175	22.0	20.2-24.2	310	35.1	32.1-38.1
180	22.4	20.6-24.6	315	35.4	32.6-38.6
185	23.1	21.1-25.1	320	36.1	33.6-38.6
190	23.4	21.6-25.6	325	36.4	34.1-39.1
195	24.0	21.8-26.2	330	37.1	34.6-39.6
200	24.4	22.3-26.7	335	37.4	35.1-40.1
205	24.6	22.7-27.1	340	38.0	35.6-40.6
210	25.3	23.2-27.6	345	38.5	36.2-41.2
215	25.6	23.7-28.1	350	39.1	36.7-41.7
220	26.2	24.1-28.5	355	39.5	37.2-42.2
225	26.6	24.6-29.0			
230	27.2	25.1-29.5			

Warranty

- 1** We guarantee this product against all manufacturing defects for 12 months from the date of sale by us or through our dealers. Consumables like dry cell etc. are not covered under warranty.
- 2** The guarantee will become void, if
 - a.** The product is not operated as per the instruction given in the Learning Material.
 - b.** The agreed payment terms and other conditions of sale are not followed.
 - c.** The customer resells the instrument to another party.
 - d.** Any attempt is made to service and modify the instrument.
- 3.** The non-working of the product is to be communicated to us immediately giving full details of the complaints and defects noticed specifically mentioning the type, serial number of the product and date of purchase etc.
- 4.** The repair work will be carried out, provided the product is dispatched securely packed and insured. The transportation charges shall be borne by the customer.

List of Accessories

- | | | |
|-----------|-----------------------------|-------|
| 1. | 3.5 MHz Convex Probe | 1 No. |
| 2. | Mains Power Cord. | 1 No. |
| 3. | Ultrasound gel..... | 1 No. |
| 4. | Learning Material (CD)..... | 1 No. |