


3Shape Implant Studio™

Implant Planning and Surgical Guide Solution



User Manual

www.3shape.com

Rx only

CE
0086

16/07/2014

IS-2.1.9.0-E-EN

Table of Contents

1	Introduction	1
1.1	Intended Use	2
1.2	Indications for Use	2
1.3	Software Device Performance	2
1.4	The User Interface.....	3
1.5	Workflow Overview.....	6
2	The Workflow - How Do I ?... ..	7
2.1	Create New Order	7
2.2	Load CBCT/CT and Surface Scans.....	9
2.3	Crop the CBCT/CT Scan	11
2.4	Align CBCT/CT to Surface Scan	12
2.5	Setup Panoramic Curve and CBCT Views	14
2.6	Define Mandibular Nerve	14
2.7	Place Virtual Crowns	16
2.8	Place and Optimize Implant Position	17
2.9	Create Surgical Guide	21
2.10	Approve Order and Create Surgical Report and Drill Protocol.....	25
2.11	Send / Export Data Files.....	26
3	Your Tools and Settings	27
3.1	Order Form	27
3.2	Load CBCT/CT and Adjust CBCT/CT Scan Visibility	31
3.3	Load Surface Scan and Optimize Visibility.....	34
3.4	Crop the CBCT/CT Scan	35
3.5	Align the CBCT/CT and the Surface Scans	37
3.6	Visualize Your Settings and Options	43
3.7	Setup Default Views and User Interface	47
3.8	Define Panoramic View.....	50
3.9	Define Nerve Channel	54
3.10	Apply Virtual Crown Tools.....	57

3.11 Apply Implant Planning Tools	62
3.12 Use the Implant Library	70
3.13 Design Surgical Guide	73
3.14 Approve Order and Generate Surgical Report and Drill Protocol	83
3.15 General Settings.....	84
4 Communicate and Export / Import Design Data	97
4.1 Export / Import Your Planning and Guide Design Data	97
4.2 Communicate With My Partner IPS Software.....	98
4.3 3Shape Dental System™	99
4.4 3Shape Communicate™	99
5 System Settings	101
5.1 System Requirements.....	101
5.2 Setup Network Installation	101
6 Manufacturer Information	102
7 Warnings and Precautions.....	103
8 Information on the CBCT/CT Scanning.....	107

1 Introduction

3Shape Implant Studio™ is a pre-operative planning software used to plan the placement of one or more implants based on an imported and aligned Computer Tomography (CT) and 3D surface scan. The implants are provided by a secure implant library. In addition, a surgical guide for a guided surgery can be designed and exported for 3rd party fabrication based on the approved implant position and a specified sleeve file defined by the manufacturer. The use of the software requires having the necessary expertise in implant dentistry.

Implant Studio™ provides a method of importing medical imaging information from radiological imaging systems, such as CT combined with an optional 3D surface scan, to a computer. Visual representations of the imaged anatomical structures (e.g. the jaw with teeth and roots) are derived allowing for a three- and two-dimensional assessment of the patient without patient contact.

Implant Studio™ adds features for pre-operative simulation of dental implant placement and surgical treatment options and allows the user to visualize the placement of dental implants and other planning objects (like mandibular nerve canals and virtual crowns) within the visualized 3D volume data. Thus, dental professionals, like implantologists, are able to precisely plan the positions, orientations, types and sizes of implants, which are to be placed in the patient's mandible/maxilla, together with the related surgical procedures keeping the aesthetic aspect in mind. Build-in safety feature, like collision detection and safety zone, minimize the risk of harm to the patient to the highest possible degree.

The case planning can be used to generate an output (STL-file) for a surgical drill guide, thus transferring the virtual case planning into physical tools enabling the surgical installation according to the virtual case planning.

Surgical drill guides are patient specific templates, which are intended to transfer the pre-operatively determined positioning of the dental implants to the patient intra-operatively, which will assist the surgeon in placing dental implants by guiding and marking drill locations and assist during the implant insertion.

A surgical report allows the user to share the planning case with others and to approve the case, if the planning has been done by a 3rd party. After the final case approval, a drill protocol is provided listing the needed parts to be ordered for the implant surgery and summarizing the planned surgery. The position of the virtual crown can be transferred to 3Shape Dental System Lab software to design the temporary or final abutment(s) and restoration(s) as indicated in the order form.

1.1 Intended Use

3Shape Implant Studio™ is a pre-operative planning software for the placement of dental implant(s) based on imported CT image data, optionally aligned to an optical 3D surface scan. The CT image data are provided from medical scanners such as a CBCT/CT scanner. Virtual Crowns can be used for optimized implant positioning under the prosthetic aspect. Furthermore, a surgical guide for use with guided surgery can be designed based on the planned implant position(s). The surgical guide data can be exported for manufacturing on an external system.

Use of the software requires that the user has the necessary medical training in implantology and surgical dentistry.



Warning! The US Federal law restricts this device to sale by or on the order of a dentist.

1.2 Indications for Use

The 3Shape Implant Studio™ software is indicated for use as a medical front-end software that can be used by medically trained people for the purpose of visualizing gray value images. It is intended for use as a pre-operative software program for simulating/evaluating dental implant placement and surgical treatment options.

Indications of the dental implants do not change with guided surgery compared to conventional surgery.

Contraindications: Pre-operative planning with Implant Studio™ must only be applied if additional radiation dose from 3D radiology (compared to conventional x-ray) can be justified.

1.3 Software Device Performance

Implant Studio™ has been developed with the intention of providing maximum safety and accuracy in regard to planning of the implant positions and the design of surgical guides. The actual performance of the software is limited by the quality, resolution and accuracy of its input data (DICOM CBCT/CT and eventual surface scan) and by the user's decisions on how to define the input data's alignment. The actual implant planning process and surgical guide design can only be performed within the limitations defined by the input data and their alignment. The input data and alignment also defines the achievable degree of accuracy of all measurements.

The software does not apply any modifications, compression, or alterations to the input data (DICOM and eventual surface scan) during their upload, alignment, usage and/or export.


Implant Studio™ DICOM Conformance Statement is available for download here:

<http://www.3shapedental.com/implant-studio-DICOM-conformance-statement>

Overall, the intended use of the software is considered to be safe if used correctly and while taking the necessary precautions. The medical benefits outweigh the overall residual risk.

1.4 The User Interface


LOGIN PAGE



When you launch Implant Studio™, you need to login to start using the system. Click the **Default operator** button to enter. You can setup the default operator or create a new one on the [Configure](#) page.

HOME PAGE


Implant Studio™ Home page displays the **Calendar** with sessions **(1)**, the **Service toolbar** **(2)** with various options and settings and the **Main toolbar** **(3)** which helps you navigate through the application.




- **(1) CALENDAR**


Allows you to create new sessions. Please see chapter [Order Form](#) for details.


- **(2) SERVICE TOOLBAR**




- Contains a list of patients, their sessions and information on every order, including images of scans. Enter patient's name in the Search field to find his sessions.
- Allows you to create a new order for the selected patient - click the **New session** button to open the [Order Form](#).




- Displays a calendar with existing sessions and available time for new sessions.
- Allows you to create a new order - click on the empty time field in **Calendar** or the  button to open the [Order Form](#).




- Opens a list of existing orders and allows you to search orders using various filters.



- Opens an inbox with the incoming system and external messages. When a new message(s) arrives, a red square with a number of unread messages appears on the tab.
- Informs you when your Implant Studio™ system receives an order with the Surface scan - a number of newly arrived orders appear in a green square on the tab. Click on the message to open the order on the Order form page. The order itself is automatically stored under **Cases**.




- Allows you to setup your Implant Studio™ system. See chapter [General Settings](#) for details.



- **Online Support** - launches a screen sharing session for online support.
- **Training Center** - opens a dialog window with training videos.
- **User Manual** - opens Implant Studio™ online help.
- **Safety Guide** - opens Implant Studio™ Safety Guide.
- **Terms of Service** - opens the document with 3SHAPE SOFTWARE GENERAL LICENSE TERMS AND CONDITIONS.
- **Third Party Licenses** - displays information on the external applications used in Implant Studio™.
- **Send feedback** - opens a form for reporting an eventual system malfunction.
- **Close** - closes the *Help* window.


- **Logout** - signs the user out and displays the Login page.
- **Shutdown** - turns off the Implant Studio™ system (a confirmation dialog window appears to make sure you want to shut down Implant Studio™).
- **Cancel** - closes the window and returns to the Main page.



The *Close* and *Help* windows also contain information on the Implant Studio™ version you are currently using.

- **(3) MAIN TOOLBAR**

The Main Toolbar provides easy access to the main software pages. It can appear both at the top or bottom of the page:



[Home page](#)



[Order form](#)




[Workflow page](#)



[Communicate page](#)

WORKFLOW TOOLBAR

The process of implant planning is comprised of a succession of important steps which are reflected in the Workflow toolbar:





To open the Workflow page, click the corresponding button in the Main toolbar. The window will open on the [Cropping](#) page allowing you to crop the CBCT scan. Note that some of the steps will remain disabled until the previous step is completed or in case an error was made during planning (e.g. implant-nerve collision).

1.5 Workflow Overview

The Implant Studio™ typical workflow in combination with other Dental software packages is illustrated below. The steps are described in detail in the corresponding chapters.

- ▶ **Step 1: Load and align open format CBCT scan and surface scan**
- ▶ **Step 2: Place virtual crowns**
- ▶ **Step 3: Plan implant positions**
- ▶ **Step 4: Design surgical guides**
- ▶ **Step 5: Create restorative components in Dental System™**
- ▶ **Step 6: Communicate with a dentist and receive his final approval of the planned case**




2 The Workflow - How Do I ?...

This section helps you to quickly find information about specific tasks:


- 2.1 [Create New Order](#)
- 2.2 [Load CBCT/CT and Surface Scan](#)
- 2.3 [Crop the CBCT/CT Scan](#)
- 2.4 [Align CBCT/CT to Surface Scan](#)
- 2.5 [Setup Panoramic Curve and CBCT/CT Views](#)
- 2.6 [Define Mandibular Nerve](#)
- 2.7 [Place Virtual Crowns](#)
- 2.8 [Place and Optimize Implant Position](#)
- 2.9 [Create Surgical Guide](#)
- 2.10 [Approve Order and Create Surgical Report and Drill Protocol](#)
- 2.11 [Send / Export Data Files](#)

2.1 Create New Order




▶ Step 1: Open a blank order form

Open a [Calendar](#) from the Service toolbar and press an empty time slot to open a new Order Form:



► **Step 2: Choose the required Lab and Template** (if available)



▶ **Step 3: Enter General information for the patient** (and other details if available)

▶ **Step 4: Select Construction details for the restoration**

Choose the relevant teeth in the map and specify the restoration and bridge type, as well as the implant and surgical guide positions.

▶ **Step 5: Use the comment field to add special instructions or comments to the lab**

See chapter [Order Form](#) for details.

2.2 Load CBCT/CT and Surface Scans

Once an order has been created, press the corresponding **Load** buttons in the Order form to load **Surface** (upper/lower) and **CBCT/CT** scans:

Implant Studio Lab Change lab

View Export Reschedule Delete Copy

General information Delivery information Patient images

First Name: John Last Name: Doe Dentist: Default operator

Comments: Delivery: Select a date Patient ID:

Upper jaw scan Load Lower jaw scan Load CT scan Load

Choose desired construction Add pre-preparation scan Open color tool

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16

32 31 30 29 28 27 26 25 24 23 22 21 20 19 18 17

Type Crown zirkon Color -- Not chosen --

See chapters [Load CBCT/CT and Adjust CBCT/CT Scan Visibility](#) and [Load Surface Scan and Optimize Visibility](#) for details.

For the CBCT/CT scans, 3Shape Implant Studio™ supports all open DICOM .dcm file formats, which comply with the 3Shape DICOM Conformance Statement. In case a DICOM .dcm fails to load, it can be due to a corrupted file structure or that the file fails to comply with the quality check during loading (please see hint below for further instructions).



Hint! In case a DICOM .dcm fails to load, please check the [3Shape DICOM Conformance Statement](#) on our website or in chapter [Information on the CBCT/CT Scanning](#) of the manual. Alternatively, please contact the original CT manufacturer for clarification.

3Shape Implant Studio™ allows to import open format STL files and 3shape .dcm surface scans. In case it fails to load, please check the quality of the surface scan and the used triangulation settings.



Warning! Make sure that the quality of the loaded CT scan is sufficient for planning the case and that the relevant areas on the image for the bone analysis are adequate in order to make an accurate and sustainable decision for implant planning and surgery!



Hint! Please use the [visualization settings for the CBCT/CT scan](#) and the ability to adjust the grey values to find the optimum for each individual scan.



Warning! User needs to ensure that the surface scan is of the required quality for planning the case, and that the relevant areas for the implant planning and placing the surgical guide fully exist.



Hint! Please check that the scan data are correctly and completely loaded.



Hint! Please check that the visualization of the CBCT and Surface Scan is working correctly to allow a secure and accurate placement of the implant during the navigated planning and surgical guide design process.




Hint! To assure optimal scan quality, please make sure that the surface scan is taken with a reliably functioning 3D scanner, and the scanner has recently been calibrated.

2.3 Crop the CBCT/CT Scan



Once the order has been created and scans loaded, click the **Workflow page** button in the Main toolbar. The software will automatically bring you to the first step - **Cropping**, where you can reduce the scan size by cropping the CBCT/CT scan and keeping only that part of a scan that will be used for planning.



- (1) Axial plane
- (2) Coronal plane
- (3) Sagittal plane
- (4) Graphics information
- (5) Cropping information

Please see chapter [Crop the CBCT/CT Scan](#) below for details.

2.4 Align CBCT/CT to Surface Scan



Once the cropping of the CBCT/CT was done (optional), you can proceed to the **Scan alignment** step. To align the CBCT/CT and Surface scans, please follow steps below:



Warning! Implant Studio™ must not be used for aligning CT and 3D optical surface scans in edentulous cases! The available alignment tools are not suited for achievement of reliable and accurate alignment results!

▶ Step 1: Place the alignment points

Place 3 points on the **Surface** scan in the left pane and 3 corresponding points on **CBCT/CT** scan in the right pane for correct alignment of the scans. The scans are automatically aligned after placing the points on marking the Confirmation checkbox while the **Surface** and **CBCT/CT** windows get hidden.


▶ Step 2: Check and correct the alignment if necessary

Use the following tools to check and improve the alignment:

- [Region of interest](#) tools
- [Difference map](#)
- [2D mode](#) in the CT pane
- [2D cross section](#)

▶ Step 3: Confirm Alignment

Approve the alignment by marking the **Confirm alignment** checkbox.



See chapter [Align the CBCT/CT and the Surface Scans](#) for more information.



Warning! Before clicking the **Confirm Alignment** check-box, please assure that the CBCT/CT and surface scans are well aligned, particularly in the area(s) where the implant(s) will be placed. Please check these also for the alignment in the 2D cross section views.



Note! You must confirm the alignment with a click on the **Confirm Alignment** in order to be able to proceed to further steps which otherwise appear disabled.



Hint! Place the 3 points in approximately the same place and same order in both scans, and ensure that the area is clearly visible in both scans. The optimal alignment should be reached within the area where the implant(s) are placed later. Please use also the cross-section views to check the alignment!




Hint! Please assure that the visualization of the CBCT/CT and surface scans appear correct and complete.

2.5 Setup Panoramic Curve and CBCT Views



After aligning the two scans, proceed to the next step **Panoramic curve** where the panoramic view for the CBCT scan can be defined:

Using the cursor, click the points along the arch to define the panoramic curve (2):




- (1) Matched surface and CBCT/CT scans.
- (2) Green panoramic curve drawn with the cursor.
- (3) The blue lines indicate the depth of the CT scan shown in the panoramic curve view.
- (4) The orange line indicates the section view along the panoramic curve.
- (5) The panoramic curve view.
- (6) Cross-section view within the blue lines called the **Sagittal** view.
- (7) Section view defined by the orange line called the **Tangential** view.
- (8) Visualization sliders.

See chapter [Define Panoramic View](#) for details.

2.6 Define Mandibular Nerve



To define a nerve channel on the **Nerve definition** step, use the panoramic view, two cross section views and the 3D view to help you assess the CBCT/CT scan from different angles and place the nerve correctly. See chapter [Define Nerve Channel](#) for details.



Warning! Please make sure, that the nerve is correctly placed and that you have chosen the appropriate size of the nerve channel diameter to compensate for potential inaccuracies in the nerve position!



Warning! If the nerve position is not clear due to poor CBCT/CT image quality, the dataset must not be used!




Hint! Please ensure that the nerve channel is clearly visible in the different views after it is placed. Please see also the chapter [Define Nerve Channel](#) to learn more about how to place them correctly.

2.7 Place Virtual Crowns




Once the panoramic curve has been defined, press the **Anatomy design** button in the Workflow toolbar to start placing virtual crowns.


► Step 1: Select desired smile from the library




► Step 2: Shape virtual crowns

Use the necessary **Transformation**, **Morphing** and **Wax Knife** tools to shape the crowns as illustrated below.




 **Hint!** Please assure that the selected crown design correlates with the intended final design.


 **Hint!** The positions of the virtual crowns need to be approximated to the intended final restoration position.

See chapter [Apply Virtual Crown Tools](#) for details.


2.8 Place and Optimize Implant Position

 Once the virtual crowns have been placed, press the **Implant planning** button in the workflow toolbar to start planning implants:


► **Step 1: Select implant manufacturer**



Press the green plus button in the Implant selection menu and select the implant manufacturer from those displayed.







► **Step 2: Select an implant from the manufacturer**

 **Warning!** Be aware that not all displayed implant systems may have regulatory clearance in all countries. Please direct any questions regarding regulatory status in your country towards the implant manufacturer or reseller.

Double-click the implant to select it and close the *Select implant* window.


Once an implant is selected, its model appears in the 3D and 2D views of the main window.

Avinent | Select implant


Platform (P) / Diameter (Ø), mm:

Length, mm	P 3.5		P 4.1				P 5.1	
	Ø 3.3	Ø 3.3	Ø 3.8	Ø 4	Ø 4.2	Ø 4.8	Ø 4.8	
7								
8.5								
10								
11.5								
13								
15								





Warning! Please make sure that the correct implant type incl. platform type/ length and diameter is chosen according to the patient clinical situation. Please consider this in particular when planning an implant case near a main anatomical feature like artery or alveolar nerve!


► **Step 3: Adjust implants**

Adjust implants and optimize their positions using the 2D views at the bottom of the screen, the main 3D modeling window and the Visualization sliders. Use the  **Enable implant aligned plane** button in the **Sagittal** and **Tangential** views to enable a plane centered around an implant.

Please see chapter [Apply Implant Planning Tools](#) for details.



 **Warning!** Please assure the correct minimum size of the safety zone for the implant and mandibular nerve. A collision of the implant with the nerve, another implant or any other main anatomical structure like nerve or artery can cause severe damage, paresthesia and complication to the patient! (see chapter [Apply Implant Planning Tools](#)). For further information on how to change the safety zone size, please see chapter [General Settings - Planner and Guide](#).

 **Warning!** In case of collision between implant and nerve or between implants, the software will disable access to further steps i.e. creation of report, protocols as well as the surgical guide design step.



Warning! Please assure that an implant does not collide with an existing implant or tooth root channel. Also assure that the implant does not penetrate the jawbone border.



Warning! During implant placement, please consider that the bone-density in reality can be different than the visualized bone-density displayed on the screen. Please see section [Load CBCT/CT and Adjust CBCT/CT Scan Visibility](#) for more details.



Hint! Please assure that the chosen implant position corresponds to the intended restoration position and the fixation/mounting method for the prosthesis on the implant.




Warning! The user must be able to recognize red and green colors to detect warnings such as *“Collision between implant and nerve”*.


2.9 Create Surgical Guide





Once implants have been placed and positioned, you can press the **Surgical Guide** button in the Workflow toolbar in case you want to create surgical guides.


► Step 1: Place the splines to mark the surgical guides' outlines

Press the **Start new patch**  button to enable a spline-placing mode and create one or several patches (you can place a few splines if necessary). Make sure you select the correct manufacturer under **Settings**.



 **Warning!** Please assure that the required material thickness and offset parameters are applied when creating the surgical guide – according to the specific material and manufacturing machine requirements you will use.

 **Warning!** Please check for potential collision between the sleeves, the visualized drill cylinders and any existing tooth or dental restoration and consider the consequences for the intended treatment procedure. Note that in case of the sleeve collision, all further steps will be disabled.

 **Warning!** The surgical guide should sit stably and firmly on the teeth, and it must not move during surgery. This is to ensure optimal guidance for the drills and implants during surgery. **Do not use the surgical guide unless you judge the fit to be optimal, safe and completely correct.**




Hint! Please design the surgical guide in such a way that sharp edges are avoided.




Hint! As you design the surgical guide, ensure that there is sufficient support and connection between the main structure and the support for the sleeves.

► Step 2: Create a surgical guide

Press the **Preview** button to create the surgical guide.



► Step 3 (optional): Add bars and windows



See chapter [Design Surgical Guides](#) for details.

3Shape Implant Studio™ will create an open STL-file ready for export after approval. This surgical guide STL file is intended to be used only for production of the guide and not supposed to be modified in any 3rd party solution.



Warning! Before you start using a surgical guide from a new supplier or machine for the first time, we strongly recommend that you validate the physical and functional quality plus the accuracy of the surgical guide! It is the responsibility of the guide manufacturer to ensure that the equipment used is capable of producing surgical guides suitable for guided surgery.



Hint! We recommend that you or your supplier apply a 3D inspection solution and the 3Shape protocol: “Implant Studio_Validation Procedure to Approve Surgical Guides” in order to validate the printed surgical guides. We recommend that you do so both before you or your supplier start manufacturing, and also at frequent intervals during the production run. In this connection, we advise you to use the 3Shape Convince Quality Control System.



Hint! Please make sure that the Surgical guide is adequately disinfected.

2.10 Approve Order and Create Surgical Report and Drill Protocol



- When you proceed to the **Save/Send** step, the software automatically creates a surgical report which is intended to be used for the review and approval procedures.
- A drill protocol, however, will be generated only after you read and agree to the terms of the "[Approval of Planning](#)" form which appears when you enter the step.

► Step 1: Go to the Save/Send step

The "Approval of Planning" form appears on the screen.

► Step 2: Read carefully and agree to the terms of the "Approval of Planning" form

► Step 3: View the created surgical report and a drill protocol by clicking the corresponding buttons

The Surgical report file can be sent out via e-mail or [3Shape Communicate](#).

See a detailed chapter below [Approve Order and Generate Surgical Report and Drill Protocol](#) for more information.



Warning! Please make sure that the presented screenshots and information for the implant positions and implant types as well as, eventually, the surgical guide in the Surgical report are sufficient in order to make a clear and safe decision to approve or modify the implant or surgical guide position and design.



Warning! Please follow the instructions in the drill protocol!



Hint! Please check the Surgical report and Drill Protocol every time before a surgery and make sure that

- the provided information is correct

- you ordered correct implants based on their order number
- the right drills and sleeves are used for the selected implant type and manufacturer
- drills and sleeves are in a good condition and ready to be used for the surgery

2.11 Send / Export Data Files



When you proceed to the **Save/Send** step, the software automatically sends the order to your partner via 3Shape Communicate IF **two conditions** are met:

1. You configured your Implant Studio™ [Connection settings](#) under the **Configure** tab.
2. You read and agreed to the "[Approval of Planning](#)" form which appears when you enter the **Save/Send** step.

After you have agreed to the "Approval of Planning" form, the order is automatically sent. If the software was not set up to send orders via 3Shape Communicate, then the output files can be found in the output folder.

In case a surface scan was used together with a CBCT/CT scan for planning, Implant Studio™ will export/send the approved implant position and surface scan as a 3Shape .dcm file. This file can be imported into any 3Shape dental software application.

See chapter [Communicate and Export / Import of Design Data](#) for details.


3 Your Tools and Settings

Implant Studio™ offers a variety of tools and settings for your implant planning and surgical guide design. See chapters below for details.


3.1 Order Form

OPENING THE ORDER FORM

- To open the Order Form and create a new order, click on the empty time field **(1)** in the **Calendar** or the **+** button **(2)**:




- You can also create a new order for the existing patient through the **Patients** tab by clicking the **New session (3)** button:





ORDER DETAILS

Fill in the provided fields with the patient's information and construction details:

- (1) Choose the required Lab and Template (if available).
- (2) Enter **General information** for the patient (and other details if available).
- (3) Upload the CBCT/CT and Surface scans.
- (4) Choose necessary teeth in the overview map with a left click.
- (5) Specify restoration details: **Type**, **Color**, **Implant Planning** and/or **Surgical Guide**.




 **Note!** The availability of the **Implant Planning** and **Surgical guide** options depends on your license version.

 **Note!** In order to be able to create a surgical guide, you must select the **Surgical guide** option in the order form. Otherwise, the **Surgical Guide** step will not be available.



Note! You must select the **Implant Planning** option for your order to be able to proceed to the next step.


TOOLS AND OPTIONS




Allows you to place markings and annotations on the teeth.



Opens order details on a separate page which can be printed out.




Allows you to export orders to a separate folder.




Opens a calendar and the list showing free time slots. Press the preferred date and time and then **Select** to reschedule the session:



Deletes an entire order.




Creates a Standard or Pre-preparation copy of an order (see the description for each below):



3.2 Load CBCT/CT and Adjust CBCT/CT Scan Visibility


UPLOAD CBCT/CT SCAN

- ▶ **Step 1:** Click the CT scan *Load* button in the Order form
- ▶ **Step 2:** Go to the folder containing the DICOM CBCT/CT scan and select a .dcm scan file
- ▶ **Step 3:** Press *Open*



► Step 4: Preview scans in the CBCT/CT loader

The CBCT/CT loader will assist you with selecting the correct scan by placing the available scans into two categories: **Supported scans** and **Unsupported scans**. Click the scan to read additional information in the window to the right.




► **Step 5: Confirm the selected scan to upload it**


Select the supported scan you want to load and then click the **Confirm** button to upload the CBCT/CT scan.

A **green tick** will appear in the **CT scan** check-box after the successful upload:



ADJUST CBCT/CT SCAN VISIBILITY

 **Hint!** Please use the visualization settings for the CBCT/CT scan and the option to adjust the grey values in order to find the optimum for each individual scan.




► **Step 1: Click the *CT visibility* button below the Visualization sliders**

This will open the CT Element Visibility with the sliders for **Air**, **Skin**, **Bone** and **Enamel** visibility.

► **Step 2: Click and drag the sliders to setup the desired CBCT/CT scan visibility**

- Move sliders to the left/right to change the voxel intensity.
- Move handles up/down to adjust the opacity.



Hint! Please ensure that the points representing bone, skin, enamel and enamel end represent correctly the different body parts and do not cause misinterpretation.


3.3 Load Surface Scan and Optimize Visibility

UPLOAD SURFACE SCAN

► **Step 1: Click the *Load* button for upper/lower jaw in the Order form**

► **Step 2: Select a surface scan file from the selection window that appears (.stl and 3Shape .dcm formats are supported)**


► **Step 3: Click *Open***




A **green tick** will appear in the **Upper/Lower jaw scan** check-boxes after the successful upload:



OPTIMIZE VISIBILITY


Use the **Preparation** slider  to toggle the visibility of the Surface scan. Click on slider's icon to hide the scan completely.

3.4 Crop the CBCT/CT Scan

 Once the order has been created and scans loaded, click the **Workflow page** button in the Main toolbar to start the Implant planning process:



The software will automatically bring you to the first step - **Cropping**, where you can reduce the scan size by cropping the CBCT/CT scan and keeping only that part of a scan that will be used for planning.



- (1) Axial plane
- (2) Coronal plane
- (3) Sagittal plane
- (4) Graphics information
- (5) Cropping information

All three windows provide different views of the - (1) **Axial**, (2) **Coronal** and (3) **Sagittal** and each has a slide bar which you can move sideways to look through the slices of the CBCT scan and better see which part to crop.

Click and drag the blue cropping box inside the windows to select the part of the CBCT/CT scan that you wish to work with (area outside the box will be deleted)

In the (4) **Graphical information** window you can view information on your video card and the size of the original CBCT/CT scan.

The **(5) Cropping information** window shows the size of the new cropped scan and if it is too small, the **Volume region** field displays **Invalid** and you will not be able to continue with Implant planning unless you crop the CBCT/CT scan correctly (in which case the field will display **Valid**).

3.5 Align the CBCT/CT and the Surface Scans




To align the CBCT/CT and Surface scans, please follow steps below:

► Step 1: Place three corresponding points on each scan




Place 3 points on the **CT** scan in the left 3D pane and 3 corresponding points on the **Surface** scan in the right 3D pane - the scans will automatically align after you place the last point. You can move the points in order to correct the alignment, or remove the points all altogether to start over.




Hint! Please use only the teeth area while placing the three points for both CBCT/CT and Surface scans. Please do the same while using the [Region of Interest](#) tool.





When moving the cursor over the **CT/Surface** panes, the following buttons appear:

- The **Clear the points**  button lets you remove incorrectly placed alignment points.
- The **Lock views**  button lets you lock and synchronize the movement of the two scans during rotation.
- The **Switch to 2D mode**  button (CT pane only) becomes available after the scans are aligned:


- o Click the button to use a cross section plane for evaluation of the alignment.
- o Move and rotate the plane to place it in the necessary position.
- o Click **Reset cross section plane** to return the plane to its initial position.
- o Place the measurement points on the CT and Surface scans with a mouse click to measure the distance between scans if necessary.





- o Click the **Enable Threshold highlighting**  button to highlight with green those areas on the CT scan which are used for the base alignment between the CT and surface scans. These areas, which are marked green, are defined by the **Threshold (HU)** slider (see 'Step 2' below).



- ▶ **Step 2: Check and correct the alignment if necessary**
 - **Threshold (HU)**



The **Threshold (HU)** slide bar lets you adjust the Hounsfield scale for the enamel in the CT scan, the effect of adjustment is shown on the images:



- **Region of interest**

Manually improve the alignment by marking with the cursor a specific area on the Surface Scan in the 3D pane to perform a regional alignment optimization for the selected areas.




Add - lets you mark the area of greater importance for the alignment algorithm.



Remove - removes only the part of the marking indicated by the cursor.




Clear selection - clears all your previous marking.



- **Difference Map**


Click the **Show difference map**  button to view the alignment result with the help of a color scale. Please pay special attention to the areas marked red and blue - they might require additional alignment, e.g. with the help of the **Region of interest** tool.




You can also use [2D Cross Section](#) to check the alignment.


► **Step 3: Confirm Alignment**

Approve the alignment by marking the **Confirm alignment** checkbox. The **CT** and **Surface** panes get hidden.



 **Note!** The scans alignment must be performed in order to proceed to further steps which, otherwise, appear disabled.

3.6 Visualize Your Settings and Options




The visualization tools can help you while planning implants and designing surgical guides.

Right-click anywhere in the modeling window to open the visualization tools menu. Three tools are available:

1. **Visualization sliders**
2. **Measurement tools**
3. **2D cross-section**


 **VISUALIZATION SLIDERS**



Use the Visualization sliders to turn on/off the visibility of

- scans
- bone
- enamel
- skin
- opacity
- contrast
- brightness
- sleeve with implant insertion direction
- drill diameter
- etc.

You can pin the slider to the modeling window and set up the default views (see [Setup Default Views and User Interface](#) for details).



Note! Only sliders that are relevant for the current step will be available.



MEASUREMENT TOOLS



Right-click the 3D scan on the **Implant Planning** step and click one of the measurement options to let the software automatically measure the distance from:




Implant to implant – shows the shortest distance to next neighboring implant.




Implant to nerve - shows the shortest distance from implant to nerve.



Implant to surface – shows the shortest distance from implant to surface scan.




2D CROSS-SECTION




Use the cross section to create a 2D cutting plane at a position of your choice. You can also use the tool to define a free definable cutting plane to inspect the implant position in the 2D format.

- Draw a section line across the model by holding down the left mouse button.
- To **rotate** the 2D plane, click and drag a **red** control point or a blue rim (see image below).
- To **move** the 2D plane alongside the scan, click and drag a **blue** control point.



After you have placed the cross section plane, a 2D preview window appears in the bottom left corner.



- Move the mouse cursor over the window and scroll to zoom the sectioned area in/out.
- Press down the mouse wheel to move the cross section in the preview window.
- Move the preview window within the screen by clicking and dragging the green border of a window (mouse cursor changes into a "hand").
- Change the window's size by clicking and dragging its borders (mouse cursor changes into arrows).
- Place two points to measure the distance between objects (see image above).
- Click the green cross button to close the 2D cross section tool.



Hint! To enlarge the 2D cross section window, move the cursor over it and press 'E' on your keyboard.

3.7 Setup Default Views and User Interface


During implant planning you can use the visualization sliders to toggle the visibility of objects in the modeling window, as well setup the default views.


► Step 1: Open the right-click menu

Right-click in the modeling window to open the menu.

See chapter [Visualize Your Settings and Options](#) for details on the right-click menu.


► Step 2: Pin the sliders to the modeling window

For easier access, you can display the visualization sliders in the modeling window by clicking the **Pin**  icon.




► **Step 3: Use sliders to control the visibility of objects**

With the sliders accessible from the modeling window, they can be used for turning the visibility of scans on/off, controlling the brightness, setting contrast etc. Drag the sliders sideways or click their icons to toggle the visibility of items.



► Step 4: Setup the default views



The software can save the default views under the **1 / 2 / 3** buttons **(1)**, so that when you click a number, the corresponding default view is displayed.


To save the default view:

1. Press **Shift** and click one of the numbers.
2. Select the desired check-boxes **(2)** to mark the sliders you want to use for the view.
3. Drag the sliders or click their icons to setup the desired visualization.
4. The software will save the settings under the selected number.

3.8 Define Panoramic View



After aligning the two scans, proceed to the next step where the panoramic view for the CBCT/CT scan can be defined. This will allow you to view the scan from different angles:




- (1) Axial view
- (2) Panoramic view
- (3) Sagittal view (blue)
- (4) Tangential view (orange)

- Every view is shown in its own box.
- The cross section plane is defined by a specific color - **blue** for **Sagittal** and **orange** for **Tangential**.
- You can move along the axis of a view by hovering the cursor over the view and scrolling with the mouse wheel while pressing Ctrl (see image above).
- If you scroll in any view window, the zoom in/out operation will be applied.

(1) Axial view - here you must place a green spline to define the **Panoramic view (2)** of a scan. Prior to placing a spline, use the adjustable plane in the 3D window to easily define the position and angle of the axial plane:


To **move** the plane up/down: click and drag the blue control point in the 3D window or move the cursor to the Axial view, press Ctrl and scroll the mouse wheel.


To **tilt** the plane: click and drag the red control points or the blue rim of the plane:



Then place the spline by left-clicking the control points along the teeth row and double-click on the last point to close the spline. Try to place each point in the center of the tooth:




 **Minimize window** - Press to hide the **Axial view** window.


 **Reset panoramic curve plane position** - Press to return the adjustable plane to its initial position.

 **Clear Panoramic points** - Press to remove the spline.

(2) Panoramic view - projected 2D view of the CBCT/CT scan defined by the green spline in the **Axis view (1)** window:




The **Height** slider is used to scale the height of the **Panoramic, Sagittal** and **Tangential** window views simultaneously.




The **Depth** slider is used to control the thickness of the visualized CBCT/CT scan projected from the buccal to the lingual side of for the panoramic view.


Click the blue plane and drag it sideways to view different parts of the scan in cross section.

The **Sagittal view (3)** and the **Tangential view (4)** show the planes which are sagittal and tangential accordingly to the panoramic spline defined in the **Axis view**:




Hover the cursor over one of the views to visualize the clipping plane options:



Enable clipping plane - cuts the scan sagittal or tangentially depending on the view you are in. Click  to disable the plane.

To move the planes sagittal or tangentially, hover the cursor over one of the corresponding views and scroll the mouse wheel while pressing **Ctrl**:





The Clipping Plane is ON in the Sagittal view (Implant planning step)

3.9 Define Nerve Channel






To define a nerve channel, use the 2D views to help you view the CBCT/CT scan from different angles and place the nerve correctly:




 **Hint!** To enlarge a specific 2D view for a better visualization, move the cursor over the desired view and press 'E' on your keyboard.


- Use the **Sagittal** and **Tangential** views to place the nerve spline: left-click in the views to define the nerve (Pic. 1 and 2 below).
- To change the position of the nerve, move the cursor over it until it turns blue and drag it in the desired direction (Pic. 3 below).
- To change the diameter of the nerve channel, move the cursor over the **Sagittal** or **Tangential** view and scroll the mouse wheel while holding down **Shift**.

 **Hint!** If you are not sure about the correct nerve position, enlarge the diameter of the nerve channel to create a larger safety zone.


- In the **Panoramic** view drag the blue plane sideways to move along the scan and find the best angle to view the channel.
- Use the clipping plane  to cut the scan sagittal or tangentially depending on the view you are in (Pic. 4 below). See more in chapter [Define Panoramic View](#).
- Click  to clear the nerve spline.




Pic 1. Defining the nerve channel - Sagittal view



Pic 2. Defining the nerve channel - Tangential view



Pic 3. Changing the position of the nerve




Pic 4. Tangential view of the nerve channel with clipping plane ON

3.10 Apply Virtual Crown Tools




After you have aligned both scans and defined a nerve channel, you can place your virtual crown:

- ▶ **Step 1: Click the *Choose smile library* button**




- ▶ **Step 2: Choose the required library in the *Select smile library* window**





► **Step 3: Click the relevant place on the scan to position the crown**

The scan needs to be rotated in such a way, that you could place the crown looking from the top:




► **Step 4: Modify your crown using the provided tools**


Select the crown you wish to modify by clicking the corresponding tooth number at the top and choose the desired tool:




The **Global transformations** tool allows you to move, scale and rotate the crown.



- Rotate the model by clicking and dragging the **red** control points. The adjacent arrows illustrate the direction of orientation.
- Drag the **green** points to change the shape of a whole model by moving it sidewise, back/forth or up/down.












The **Morphing** tool helps you shape your crown by dragging its parts.



- The net shows the direction in which the morphing is applied.
- Set up the morphing radius  by dragging the bar. You can also change the radius by holding the **Shift** button and scrolling the mouse wheel.



Use the **Wax knife** tool to add  / remove  material to the crown or smooth  its surface.

- Choose one of the three options and click on the model to apply the operation.
- Select the radius  and amplitude  of the wax knife by using the two slide-bars (or use **Shift** + mouse wheel and **Ctrl** + mouse wheel respectively).
- The buttons numbered **1 - 6** contain the preset values which you can apply while using a Wax knife:
 - **1 - 2** is "Add material" with different sets of the default values for radius and amplitude
 - **3 - 5** is "Remove material" with different sets of the default values for radius and amplitude
 - **6** is "Smooth"




- The **Undo**  button reverses the most recent editing command.
- The **Redo**  button restores the changes back.

Cancelling any operation requires that the operations following are cancelled first.


3.11 Apply Implant Planning Tools



After placing a crown, go to the **Implant planning** step to position your implant:



► Step 1: Select implant manufacturer




Press **+** or use the shortcut keys to select the implant manufacturer and implant type for the selected indications.


See the chapter [Use the Implant Library](#) for details.

► Step 2: Modify your implant's position


2D VIEWS

Use the **2D views** to position your implant correctly:


- Enable/disable the clipping planes  in the **Sagittal** and **Tangential views** - see the chapter [Define Panoramic View](#) for details.
- Use the **Distance measurement** tool to measure distances between objects:









Click the **Insert distance measurement** button to start placing measurement points.




Place two points and view the result. Click the **Delete the active measurement** button to clear your measurement.


- Use the  **Enable implant aligned plane** button in the **Sagittal** and **Tangential** views to enable a plane centered around an implant. When enabled, you can scroll in either of the views and the planes (now of different colors) rotate around an implant:



- Click the **Reset implant position**  button if you wish to manually change the implant's position which initially was automatically pre-positioned based on the Virtual Crown placement.
 - Select the **Reset manually**  button.
 - Click two points in the **Sagittal** or **Tangential** views, making sure that the first point you click is for the occlusal side of the implant.
 - Click the **Reset from crown pose**  button to go back to the automatic placement based on the Virtual Crown.




Manual implant positioning - placing 2 points




Manual implant positioning - result

- Rotate** the implant around apical and occlusal ends using the handlers in the **Sagittal** and **Tangential** views:



- Press and hold the **Ctrl** button on your keyboard and move your implant with a cursor up/down or left/right. This will restrict the movement of an implant along the axis:




Hint! To enlarge a specific 2D view for a better visualization, move the cursor over the desired view and press 'E' on your keyboard.


3D MODELLING WINDOW

Use the 3D modelling window to change the length and the position of your implant:


Moving



Rotating




Move the cursor over the implant until it turns into a "hand". Then left-click the implant and move it. Left-click the red control point and drag it sideways to rotate the implant.




Place the implant in such a way that you view it from the top - this will visualize the red handles for rotating the implant around its axis.

Changing length



Left-click the green control point and drag it up/down to change the length of the implant. The implant diameter will not be changed.

Invalid positioning






If an implant's safety zone border (*) intersects that of another implant or a nerve, the implant turns red and all the following steps will be disabled.


VISUALIZATION SLIDERS

Use the Visualization sliders in the upper right corner of the window to turn the visibility of objects on the screen ON/OFF - see chapter [Visualize Your Settings and Options](#) for details.

For example:

- Use the **Center view clipping** slider  to center view around the selected implant (see image below) and analyze the bone structure around the implant by zooming in and out in the range between 0 and 4mm.

- Use the **Implant insertion directions** slider  to visualize the implant direction and sleeves **(1)**.
- Move the **Drilling cylinders** slider  to visualize the screw hole direction **(2)**:




3.12 Use the Implant Library

On the **Implant planning** step you will be using the Implant Library to choose the desired manufacturer and implant type for the selected indication.

You can open the library **a.** using the shortcut keys or **b.** using the buttons on the Implant Selection box:

a. Shortcut keys


Select the indication with a mouse-click (see image) and use one of the following shortcut keys:





+	Add a new implant for the selected indication.
-	Remove the existing implant for selected indication.
L	Change type of the currently selected implant. This will open the implant library with the previously chosen implant system.
Esc	Close the implant library

b. Implant Selection box buttons




Click  to add an implant to your scan.

Click  to remove the existing implant.

Click  to change the current implant manufacturer or the type of implant:

► Step 1: Select the desired manufacturer from the list


Left click on a manufacturer to select their implant solution. This will open the list of implants by this particular manufacturer.



Note! To add or remove an implant manufacturer to/from this list, please go to the [Configure](#) page.

► Step 2: Select the implant type

Left-click on the desired implant type to select it:



► **Step 3: Select the relevant implant size**

Left-click in the table to choose the necessary **Length** and **Diameter** of the implant. Here you can also select a **Platform** type from the drop-down menu.


Platform (P) / Diameter (Ø), mm:

P	3.5	4.1				5.1	
Ø	3.3	3.3	3.8	4	4.2	4.8	4.8
Length, mm	7						
8.5							
10							
11.5							
13							
15							


► **Step 4: View additional information about the implant**

After selecting an implant, an illustration of the implant together with the additional information appears in the preview window to the right.


←
Avinent | Select implant
✖




Coral HE



Coral HI




Ocean HE



Ocean HI

Platform (P) / Diameter (Ø), mm:



	3.5	4.1				5.1	
	3.3	3.3	3.8	4	4.2	4.8	4.8
Length, mm	7						
	8.5						
	10						
	11.5						
	13						
	15						



Additional information:

Avinent Coral HE
 Implant length = 10 mm
 Implant diameter = 3.3 mm
 Platform: 3.5

4

- Preview the newly selected implant in the 3D modeling window and in the 2D cross section windows to determine if the selected implant is right for the case.
- Navigate between the Manufacturer and Implant list by clicking  in the upper left corner of the form.
- Click  to close the Manufacturer/Implant selection form.


3.13 Design Surgical Guide



After the implants have been planned, proceed to **Surgical Guide** design if applicable:

► **Step 1: Press the *Start new patch*  button**



This will enable the spline placing mode.



► **Step 2: Place the splines to mark the surgical guides' outlines**

Move a cursor over the scan - a green control point appears allowing you to place a spline.

- You can place a spline by **(1)** clicking points one after another or **(2)** drawing the spline with a cursor:




(1)

(2)

- You can also fast edit the spline by drawing or clicking close to the spline to correct its placement:




Fast editing





Result

- To close the spline - place the last control point of the spline on top of the first one or click




- inside the yellow zone around the point
- You can create several patches of splines or create only one spline covering both neighboring teeth and a drill cylinder as illustrated below:



- o Click  to start placing another spline. As a result, you will have several splines, one that is active (green) and another that is inactive (grey). Click on a spline to make it active.


- o Click  to remove an active spline.

- o Click  to clear all the splines.

▶ Step 3: Set Insertion direction

Click **Set From View** to define the insertion direction from the view perpendicular to the screen (the blue arrow on the screen shows the direction).

Click **Auto Detect** to allow the software to automatically set the insertion direction.





► **Step 4: Adjust the sleeve and the drill cylinder**

Sleeve: 30



Change the sleeve type if you are using a Pilot Drill.



If the implant system allows it, you can define the offset between a sleeve and an implant by clicking the upper control point and moving it up and down.

Click and drag the three smaller control points to adjust the supporting material around the sleeve, i.e. the diameter and the support below the sleeve.

► **Step 5: Configure the Surgical guide parameters**

Select a manufacturer with the parameters which have been pre-set in the [Planner and Guide](#) setting under the **Configure** tab.

Manufacturer

Dreve ProDiMed_M12 ▼

Thickness 2.00 mm

Offset from teeth 0.075 mm

Retention amount 0.000 mm

Offset from sleeve 0.030 mm

Thickness - defines the thickness of the surgical guide shell.

Offset from teeth - defines the offset distance from teeth to shell.


Retention amount - defines how tightly the surgical guide will fit on the teeth (see images below).

Offset from sleeve - defines the offset distance from sleeve to shell.




Note! The created surgical guide STL file is intended to be used only for production on the selected machine, material and manufacturer. 3Shape gives no guaranty regarding the performance of a surgical guide produced on machines and materials other than the selected ones.

2D cross-section: A distance measured from a tooth to a surgical guide = 0.99 mm




Retention amount = 0.000




Retention amount = 0.100

► **Step 6: Create surgical guides**




Press the **Preview** button to create a surgical guide.

The software needs a few seconds to create the surgical guide. View the progress bar for status.




► **Step 7 (optional): Add bars and windows**





1. Go to the **Bars and windows** sub-step.
2. Choose to add a support bar or a viewing window from the menu.
3. Place a support bar and/or a viewing window on the guide.
4. Modify the bar and the window by using the handles for rotation and size adjustment. Click and drag bar/window to move it.
5. Click **Preview** to apply the supplements. Press the **Preview** button one more time to go back to the editing mode.

Adding a viewing window - before and after **Preview**:



Adding a support bar - before and after **Preview**:



3.14 Approve Order and Generate Surgical Report and Drill Protocol



When you proceed to the **Save/Send** step, the software automatically creates a surgical report with information about the chosen implant types.

- Click the **Show surgical report** button if you want to send a pre-approval copy of the surgical report to the doctor responsible for the case so he could approve it by signing the report.
- A drill protocol, on the other hand, will be generated only after you read and agree to the terms of the "Approval of Planning" form which appears when you enter the step (see below).

To approve implant planning, click the **Approve & Save** button. If you do not wish to approve it, close the "Approval of Planning" form by clicking on the cross in top right-hand corner of the form (note that in this case the drill protocol will NOT be created and the output file will not be saved).

Show surgical report

Approval of Planning


You have now completed your implant planning using the 3Shape Implant Planning software. By sending the order to your Lab and/or Surgical Guide partner, you also approve the quality of your own planning. Before sending the order, please accept the following items on the checklist:

- I am satisfied with, and approve the medical and clinical aspects of my planning.
- I approve that I performed the planning carefully, chose the most appropriate implant solution and safety zone in relation to any potential existing teeth, dental restorations, the main anatomical structure and the intended final dental restoration solution.
- In case it's applicable I approve that the marking of the mandibular nerve is performed correctly.
- I confirm that the scan data used for the CBCT and surface scan are up-to-date, and that their quality and visualization is sufficient for the planning situation and intended final solution.
- I confirm that the design of the surgical guide was performed correctly and ensures a stable and firm fit during surgery.
- I confirm that I will read the drill protocol before surgery and verify that the provided information corresponds with the intended planned surgery.
- I agree to the fact that neither 3Shape nor its resellers will perform any medical or clinical review of my planning.
- I am fully aware of the technical constraints that govern the manufacturing of the surgical guide.
- I agree that 3Shape and its reseller do not assume any liability for damages or injuries resulting from my planning and treatment.

I have read and agree to all of the above

Approve & Save

After you have agreed to the "Approval of Planning" form, it closes and you enter the following options of the **Save/Send** step become available:




Show surgical report - click to open a Surgical report (the Surgical report after the approval has a status "Approved" written inside it, while the one before the approval has a status "For Approval").

Show drill protocol - click to open a Drill protocol.

Show approval - click to open the "Approval of Planning" form.

Open save folder - click to open the output folder. In case a surface scan was used during planning, Implant Studio™ will create a file of the approved implant position with the surface scan as a 3Shape .dcm file. This file can be imported into any 3Shape dental software application.

3.15 General Settings



Click the **Configure** button in the main toolbar to use various Implant Studio™ settings.

GENERAL

Settings

General

- System
- Sound
- Operator
- Connection
- Planner and Guide
- Implant library
- Lab Connections
- Order
- Order Inbox
- Messages

General

About

Software version: Implant Studio 2014
Software build: 2.1.6.0 (159625)
Dongle ID: 2601772594
Site ID: 59996
Expiry date: 27-Oct-14
Next key renewal date: 26-Dec-13

Check for upgrades

User interface language

English / English

What's new

Renew license online

3Shape support

About - provides information on software version, dongle number etc.

Check for upgrades - checks for software updates allowing you to keep your Implant Studio™ system up-to-date.

User interface language - allows you to select the preferred interface language.

What's new - opens a *What's new in Implant Studio* pdf.

Renew license online - opens a web page for the online renewal of your license.

3Shape support - provides you with the remote support directly from 3Shape. Tell 3Shape the ID and password displayed on the screen for desktop sharing.

SYSTEM

The screenshot shows the 'System' settings page. On the left, a sidebar lists various settings categories: General, System (highlighted in green), Sound, Operator, Connection, Planner and Guide, Implant library, Lab Connections, Order, Order Inbox, and Messages. The main content area is titled 'System' and contains the following items:

- Log**: A button labeled 'Open system log file' is positioned to the right of the 'Log' label.
- Free harddisk space**: Displayed as '4.2 GB'.
- Auto show On-screen Keyboard**: A green checkmark icon indicates this feature is enabled.
- Support application path**: A text box containing the path 'C:\Program Files\3Shape\Implant Stu'.

Log - opens system log file.

Free hard disk space - shows the remaining space on your hard disk.

Sound Volume - controls the system sound volume.

SOUND

Settings

- General
- System
- Sound**
- Operator
- Connection
- Planner and Guide
- Implant library
- Lab Connections
- Order
- Order Inbox
- Messages

Sound

Mute On-screen Keyboard

OPERATOR

The screenshot shows the 'Settings' application with the 'Operator' section selected. The left sidebar lists various settings categories: General, System, Sound, Operator (highlighted), Connection, Planner and Guide, Implant library, Lab Connections, Order, Order Inbox, and Messages. The main 'Operator' panel contains a list of operators with 'Default operator' and 'New operator' visible. To the right of the list are three buttons: 'Add new', 'Modify', and 'Delete'. Below the list, there are input fields for 'Name' (containing 'Default operator') and 'Password'. There are also checkboxes for 'Is administrator' (checked) and 'Automatically login' (unchecked). At the bottom right, there is a 'Photo' field with a placeholder image of a person in a white lab coat.

Add new - creates new Implant Studio operator.

Modify - lets you edit the operator settings.

Delete - removes the selected operator from the list.

Name - operator's name.

Password - operator's password. If you do not want to have a password, leave this field blank.

E-mail address - operator's e-mail address.

Is administrator - allows you to add/delete users when enabled. The first Implant Studio™ operator is always an administrator by default.

Automatically login - automatically logs in the Implant Studio™ operator.


Import photo - allows you to upload your personal photo in .png or .jpg formats. User photos appear next to the operators names in the *Choose user* login screen.

Clear photo - removes the photo.

Save - saves the photo.

Cancel - removes the photo.

CONNECTION




PLANNER AND GUIDE

The screenshot displays the 'Planner and Guide' settings interface. On the left is a sidebar with the following menu items: Settings, General, System, Sound, Operator, Connection, **Planner and Guide** (highlighted), Implant library, Lab Connections, Order, Order Inbox, and Messages. The main content area is titled 'Planner and Guide' and contains two sections:

- Implant planning settings**:
 - Implant safety zone: Radial distance (1.50) and Apical distance (2.00), each with a slider control.
- Surgical guide settings**:
 - Default manufacturer: A dropdown menu showing 'Drevo ProDiMed_M120_FotoDentLED'.
 - Guide thickness: 2
 - Offset from teeth to guide: 0.075
 - Retention amount: 0
 - Offset from sleeve to guide: 0.03

Radial / Apical distance of implant's safety zone (mm) - adjusts the safety zone around the implant. Click and drag the slider bar to change the distance from an implant to the safety zone.



(1) Radial distance

(2) Apical distance

Default manufacturer - select from the list one of the default manufacturers which has predefined minimum and maximum values for each of the settings below. You can modify these settings manually:

Shell thickness - defines the thickness of the surgical guide shell.

Offset from teeth to shell - defines the offset distance from teeth to shell.

Retention amount of shell - defines how tightly the surgical guide will fit on the teeth (see chapter [Design Surgical Guides](#) for images).

Offset from sleeve to shell - defines the offset distance from sleeve to shell.

IMPLANT LIBRARY

You can upload implant libraries you wish to work with into your Implant Studio™.

1. To check for new libraries, click the **Check again** button.
2. Select library in the list.
3. Click **Download**. The library will be added to your list of implant manufacturers.

The screenshot shows the 'Settings' window with the 'Implant library' section selected. The left sidebar lists various settings categories, and the main area displays a list of implant categories and manufacturers, all with checkboxes and 'Up to date' status. At the bottom right, there are 'Check again' and 'Download' buttons.

Category	Manufacturer	Status
Implant systems	Avinent	Up to date
	Bego	Up to date
	Biodenta	Up to date
	Glidewell	Up to date
Sleeves	Avinent	Up to date
	Bego	Up to date
	Biodenta	Up to date
	Glidewell	Up to date
Surgical trays	Full combinations	Up to date

LAB CONNECTIONS

Settings

- General
- System
- Sound
- Operator
- Connection
- Planner and Guide
- Implant library
- Lab Connections**
- Order
- Order Inbox
- Messages












Lab Connections

External Connection Method: Communicate

- Implant Studio (local)

Configure Communicate Refresh

ORDER

Settings	Order
 General	<input type="text" value=""/>
 System	Dentist address
 Sound	Tooth numbering system <input type="text" value="UNN"/>
 Operator	Color system <input type="text" value="Vita color"/>
 Connection	Remember last indication settings <input checked="" type="checkbox"/>
 Planner and Guide	Calendar
 Implant library	Time interval between orders (in minutes) <input type="text" value="60"/>
 Lab Connections	Workday start <input type="text" value="7:00 AM"/>
 Order	Workday end <input type="text" value="7:00 PM"/>
 Order Inbox	Patient names privacy <input type="text" value="Show toggle button"/>
 Messages	

ORDER INBOX

Settings

- General
- System
- Sound
- Operator
- Connection
- Planner and Guide
- Implant library
- Lab Connections
- Order
- Order Inbox**
- Messages

Order Inbox

Direct Connect Folder:

Direct Connect Username:

Direct Connect Password:

MESSAGES

Settings

- General
- System
- Sound
- Operator
- Connection
- Planner and Guide
- Implant library
- Lab Connections
- Order
- Order Inbox
- Messages**

Messages

Delete messages older than Month

Delete only read messages

4 Communicate and Export / Import Design Data

Exchange data with dentists or labs and discuss cases online. See chapters below for information.



Note! To be able to send/receive orders and communicate with your partners, you must first configure your [Connection settings](#) under the **Configure** tab.

4.1 Export / Import Your Planning and Guide Design Data



When you proceed to the **Save/Send** step, the software automatically sends the order to your partner via 3Shape Communicate **IF** the following **conditions** are met:

1. You configured your Implant Studio™ [Connection settings](#) under the **Configure** tab.
2. You read and agreed to the "Approval of Planning" form (see below) which appears when you enter the **Save/Send** step.

To approve implant planning, click the **Approve & Save** button. If you do not wish to approve the order, close the "Approval of Planning" form by clicking on the cross in top right-hand corner of the form (please note that in this case the order will not be sent).

Show surgical report ✕

Approval of Planning

You have now completed your implant planning using the 3Shape Implant Planning software. By sending the order to your Lab and/or Surgical Guide partner, you also approve the quality of your own planning. Before sending the order, please accept the following items on the checklist:


- I am satisfied with, and approve the medical and clinical aspects of my planning.
- I approve that I performed the planning carefully, chose the most appropriate implant solution and safety zone in relation to any potential existing teeth, dental restorations, the main anatomical structure and the intended final dental restoration solution.
- In case it's applicable I approve that the marking of the mandibular nerve is performed correctly.
- I confirm that the scan data used for the CBCT and surface scan are up-to-date, and that their quality and visualization is sufficient for the planning situation and intended final solution.
- I confirm that the design of the surgical guide was performed correctly and ensures a stable and firm sit during surgery.
- I confirm that I will read the drill protocol before surgery and verify that the provided information corresponds with the intended planned surgery.
- I agree to the fact that neither 3Shape nor its resellers will perform any medical or clinical review of my planning.
- I am fully aware of the technical constraints that govern the manufacturing of the surgical guide.
- I agree that 3Shape and its reseller do not assume any liability for damages or injuries resulting from my planning and treatment.

I have read and agree to all of the above Approve & Save

After you have agreed to the "Approval of Planning" form, the order is automatically sent.



If the software was not set up to send orders via 3Shape Communicate, then all output files and data can be found in the output folder by pressing the **Open save folder** button in the next window.



4.2 Communicate With My Partner IPS Software



Communicate with your partner using 3Shape Communicate - access cases and communicate anywhere and at any time, directly from your PC, tablet or smart phone.



4.3 3Shape Dental System™




3Shape Dental System™ provides 3D scanning, CAD modeling, order management, intraoral scanning, dentist-lab communication, and it interfaces to essentially all manufacturing equipment and materials.

Your implant planning and guide design can be sent directly to 3Shape Dental System™ installations in labs via [3Shape Communicate](#) - an advanced communication and visualization tool allowing the labs and dentists discuss patient cases while viewing their 3D models on-line.

4.4 3Shape Communicate™

3Shape Communicate™ is an advanced service that lets labs and dentists share case information, view 3D designs online and exchange messages.



3shape Communicate™ - Login page



Note! To use 3Shape Communicate™, you must register your clinic or your lab on the 3Shape Communicate™ web site - www.3shapecommunicate.com and configure your Implant Studio™ Connections settings. You can use your credentials to access your 3Shape Communicate account from any remote PC, tablet or smart phone.

Cases

Dentist Site

Patient Name	Order No.	Delivery Date		
 Carrew, John	63080_20110907_0952_Techw...	9/8/2011	Preview	Open Case
 Williams, Serena	63080_20110906_1419_Techs...	9/8/2011	Preview	Open Case
 Wosniaki, Caroline	63080_20110906_1407_Techi...	9/8/2011	Close	Open Case

Material: Zirkon
Type : Coping 37

Patient No.: -
Status: WaitingApproval



3shape Communicate™ - Main page

5 System Settings

Please make sure your computer meets the minimum [system requirements](#) prior to using Implant Studio™.

5.1 System Requirements

Item	Minimum Requirements	Recommended
OS:	Windows 7 or 8 64-bit	Windows 7 or 8 64-bit
RAM:	8GB or better	16 GB
Monitor Resolution:	1280x800 or similar	1920x1200 or similar
Video Card Memory:	1GB GeForce	2GB GeForce
Available HDD Space:	250GB	500 GB
CPU:	IntelCore i5 or equivalent	IntelCore i7 or equivalent
Network:	Network Internet connection	Network Internet connection
Mouse:	With the wheel button	With the wheel button



Note! We recommend using the software in combination with a 3DConnexion space ball.



Hint! Please assure you have an appropriate graphic card capable of visualizing the CBCT and surface scan in the required quality and resolution according to the system specifications.

5.2 Setup Network Installation

In order to import/export orders and communicate with your partner, your computer must be connected to the Internet - through Wi-Fi or cable.

6 Manufacturer Information

3Shape is a global company with offices in Europe, Asia and the Americas.

3Shape Implant Studio™ is manufactured by the 3Shape Medical A/S.

3Shape Medical A/S



Europe, Middle East & Africa Sales

Holmens Kanal 7

1060 Copenhagen K

Denmark

Representative offices:

3Shape North America

North American Sales

Somerset Hills Corporate Center

10 Independence Boulevard,
Suite 150

Warren, New Jersey 07059, USA

3Shape (Shanghai) Co., Ltd

Asian Sales

Room 906, Tower A of Eton Place

No. 69, Dongfang Road

200120 Shanghai, China

3Shape Latin America

Latin American and Caribbean
Sales

Carrera 13 # 82-91

Oficina 401

110221 Bogotá, Colombia

3Shape.com - Visit our corporate web-site to learn more about 3Shape and our innovations.

3shapedental.com - Visit our dental web-site to learn what solutions 3Shape has for your clinic or lab.

support.3shape.com - Visit our support site to download a printable version or to request to receive a printed copy of the '[3Shape Implant Studio™ User Manual](#)'.

7 Warnings and Precautions

Load CBCT/CT and Surface Scans



Warning! Make sure that the quality of the loaded CBCT/CT is sufficient for planning the case and that the relevant areas in the CT for the bone analysis are adequate in order to make an accurate and sustainable implant planning and surgery!



Warning! User needs ensure that the surface scan is of the required quality for planning the case, and that the relevant areas for the implant planning and placing the surgical guide fully exist.

Align CBCT/CT to Surface Scan



Warning! Implant Studio™ must not be used for aligning CT and 3D optical surface scans in edentulous cases! The available alignment tools are not suited for achievement of reliable and accurate alignment results!



Warning! Before clicking the **Confirm Alignment** check-box, please assure that the CBCT/CT and surface scans are well aligned, particularly in the area(s) where the implant(s) will be placed. Please check these also for the alignment in the 2D cross section views.

Define Mandibular Nerve



Warning! Please make sure, that the nerve is correctly placed and that you have chosen the appropriate size of the nerve channel diameter to compensate for potential inaccuracies in the nerve position!



Warning! If the nerve position is not clear due to poor CBCT/CT image quality, the dataset must not be used!

Place and Optimize Implant Position



Warning! Be aware that not all displayed implant systems may have regulatory clearance in all countries. Please direct any questions regarding regulatory status in your country towards the implant manufacturer or reseller.



Warning! Please make sure that the correct implant type incl. platform type/ length and diameter is chosen according to the patient clinical situation. Please consider this in particular when planning an implant case near a main anatomical feature like artery or alveolar nerve!



Warning! Please assure the correct minimum size of the safety zone for the implant and mandibular nerve. A collision of the implant with the nerve, another implant or any other main anatomical structure like nerve or artery can cause severe damage, paresthesia and complication to the patient! (see chapter [Apply Implant Planning Tools](#)). For further information on how to change the safety zone size, please see chapter [General Settings - Planner and Guide](#).



Warning! In case of collision between implant and nerve or between implants, the software will disable access to further steps i.e. creation of report, protocols as well as the surgical guide design step.



Warning! Please assure that an implant does not collide with an existing implant or tooth root channel. Also assure that the implant does not penetrate the jawbone border.



Warning! During implant placement, please consider that the bone-density in reality can be different than the visualized bone-density displayed on the screen. Please see section [Load CBCT/CT and Adjust CBCT/CT Scan Visibility](#) for more details.



Warning! The user must be able to recognize red and green colors to detect warnings such as *“Collision between implant and nerve”*.

Create Surgical Guide



Warning! Please assure that the required material thickness and offset parameters are applied when creating the surgical guide – according to the specific material and manufacturing machine requirements you will use.



Warning! Please check for potential collision between the sleeves, the visualized drill cylinders and any existing tooth or dental restoration and consider the consequences for the intended treatment procedure. Note that in case of the sleeve collision, all further steps will be disabled.



Warning! The surgical guide should sit stably and firmly on the teeth, and it must not move during surgery. This is to ensure optimal guidance for the drills and implants during surgery. **Do not use the surgical guide unless you judge the fit to be optimal, safe and completely correct.**



Warning! Before you start using a surgical guide from a new supplier or machine for the first time, we strongly recommend that you validate the physical and functional quality plus the accuracy of the surgical guide! It is the responsibility of the guide manufacturer to ensure that the equipment used is capable of producing surgical guides suitable for guided surgery.



Note! The created surgical guide STL file is intended to be used only for production on the selected machine, material and manufacturer. 3Shape gives no guaranty regarding the performance of a surgical guide produced on machines and materials other than the selected ones.

Approve Order and Create Surgical Report and Drill Protocol



Warning! Please make sure that the presented screenshots and information for the implant positions and implant types as well as, eventually, the surgical guide in the Surgical report are sufficient in order to make a clear and safe decision to approve or modify the implant or surgical guide position and design.



Warning! Please follow the instructions in the drill protocol!

Information on the CBCT/CT Scanning



Warning! Clinicians are responsible for generating a CT scan of optimal quality using as low radiation dose as possible.



Warning! The CT scanner should always be calibrated within the intervals outlined in the manufacturer's calibration guidelines.



Warning! If the alignment process is not conducted on the basis of high quality scans, the foundation for the implant planning procedure is not present.

8 Information on the CBCT/CT Scanning

I. CBCT/CT SCAN INFORMATION

To use 3Shape Implant Studio™ for implant planning, a digital image of the patient is required from an external CT scanner. In cooperation with the 3D optical surface scan, the CT scan forms the basis for clinical diagnostics, planning and evaluation of optimal implant positions. To acquire the CT scan, the following recommendations should be considered.

CBCT/CT Scanner Compatibility Recommendations

3Shape Implant Studio™ is compatible with medical single-slice CT scanners, multi-slice CT scanners and Cone-beam CT scanners, which are exporting DICOM files conforming to 3Shape DICOM Conformance Statement. In addition, these basic recommendations should be followed:

- The field of view for the jaw bone should be large enough to allow a diagnostic evaluation of the clinical situation, as well as ensuring enough information for a reliable alignment to the surface scan.
- The voxel size should not exceed 0.5 mm in any direction.
- The diagnostic image quality should be adequate for clinicians to appropriately read the CT image data.



Warning! Clinicians are responsible for generating a CT scan of optimal quality using as low radiation dose as possible.



Warning! The CT scanner should always be calibrated within the intervals outlined in the manufacturer's calibration guidelines.

II. SCAN PARAMETERS

1. Single-slice and Multi-slice CT

a. Patient scan settings:

The correct settings for tube voltage and effective tube current depend on the CT type, the anatomy and the intended scan of the patient.

- Spiral CT
- No gantry tilt



Hint! For optimal settings of tube voltage (kV) and effective tube current (mAs) consult the manufacturer of the CT scanner to obtain settings optimal for dental CT scanning.

b. Reconstruction settings:

- Reconstruction interval: 0.5 mm or half detector width (typically 0.5 mm or smaller).
- Reconstruction kernel: A sharp bone filter is preferred.

2. Cone-beam CT scan settings

a. Patient scan settings:

- The voxel size should not exceed 0.5 mm in any direction.

b. Reconstruction settings:

- Tilting of axial slices should not be allowed.

III. SCANNING PROCESS

Scanning and positioning of the patient should always be conducted according to the guidelines provided by the manufacturer of the CT scanner.



Hint! Do not scan the patient with a closed bite. Instead, stabilize the mandible relative to the maxilla, e.g. by instructing the patient to bite gently on a bite separator like a bite pad or cotton pads.



Hint! Non-fixed metal parts in both jaws (e.g. dentures) should be removed prior to scanning unless they are needed for diagnostic reasons.

Avoiding potential errors

- Patient positioning: If the patient moves during the scanning process, motion artifacts will occur and discontinuity in the anatomy or double anatomical borders will appear in the image. This degrades the quality of the image data and makes the process of aligning the CT scan to the surface scan difficult or maybe even impossible.
- Metal artifacts: Shapes and anatomy of a tooth restored using metal or ceramic materials, will not be clearly reproduced in a CT scan due to metal artifacts. If such teeth and also adjacent teeth, are clearly impaired by metal artifacts, they should not be used in the alignment process.



Warning! If the alignment process is not conducted on the basis of high quality scans, the foundation for the implant planning procedure is not present.