

X-ray multilayers optic

User's manual

INEL [®] Z.A C.D 405 – 45410 ARTENAY [©] +33.(0)2.38.80.45.45 [©] 0+33.(0)2.38.80.08.14 ⊠ info@inel.fr http://www.inel.fr

Summary

PRESENTATION	1 – 2
EXAMPLES	3 – 4
SETUP	5 -10
$1 - \underline{\text{Description}}$ 2 - Setting the mechanics	6 6 - 10
2a - REMOVING THE MONOCHROMATOR MECHANICS	7
2b – SETTING THE MIRROR/MONOCHROMATOR MECHANICS	7
2c - Differents kinds of setups	8
2d – Adjustements	9
2e – MIRROR CHARACTERISTICS	10

PRESENTATION

Multilayers parabolics mirrors for x-radiation are comprised of nanometric stacks of alternating transparent and absorbent compounds. The mirror captures a wider solid angle of incoming beam than a standard flat germanium or graphite monochromator and its graded multilayers output a concentrated parallel beam. Results have shown that depending upon experiment type, for equivalent resolution and settings, intensities can increase by factors of 3 to 7 times.



Figure 1 : Mirror mounted on the Inel diffractometer, located between the tube shield and monochromator housing.

The parabolic mirror does not offer $K\alpha 1$ and $K\alpha 2$ separation. We therefore have designed our system to be multifunctional. The mirror can be used by itself, or in combination with the flat germanium

monochromator if an intense, pure $K\alpha 1$ beam is desired. (See figure 1.) Parabolics mirrors have graded d-spacing designed to work with specifics wavelengths.

EXAMPLES



Figure 2 : Comparison of silicon standard data (SRM640b) with the direct beam (line 2), the filtered direct beam (line 3), and with the mirror (line 1). The experimental conditions were identical for the 3 measurements.

With the mirror (line 1), the intensity is double that of the unfiltered direct beam (line 2). In line 1 only the K α 1 and K α 2 contributions are present whereas with the unfiltered direct beam, all energies in the copper spectra, K α 1, K α 2, K β give diffraction contributions. Using a nickel filter with the direct beam (line 3) similar resolution is obtained, but the intensity is 5 time less than with the mirror.



Figure 3 : System comparison between setup with a germanium crystal alone vs. germanium crystal plus X-ray mirror in tandem.

Silver behenate diffractograms, 1 mm diameter capillary, with mirror (black curve) and without the X-ray mirror (red curve). Identical experimental conditions (monochromator Ge 111, 40 kV - 30 mA, Cu K α 1, spinning sample, 10 minutes acquisition time). An intensity gain of 3.5 is evidenced in this case.

SETUP

The Inel diffractometer equipped with the X-ray mirror, can beset belongs 3 differents possibilities:

- mirror + monochromator,
- mirror only,
- monochromator only.

The mirror is used as an intensifier of the primary beam, and is able to make parallel the primary beam (with 0.7 mm thickness). However, it is not able to separate $k\alpha 1$ and $k\alpha 2$.

The monochromator is able to separate the k α doublet, to give a monochromatic beam, but decreases the intensity. It is the reason why that for high intensity and monochromatic beam, it is required to use the setup Mirror + Monochromator. When high flux is necessary, then the monochromator can be removed. The third setup is possible to do, but not interesting, because it corresponds to the first one with a low intensity.

<u>1- DESCRIPTION</u>





Figure 5 : Beam path and angles on the Inel diffractometer using the X-ray mirror, for the copper radiation

<u>2- SETTING THE MECHANICS</u>

IMPORTANT : Each time user needs to dismount any mechanics, the power and the water cooling must be off.

2a- Removing the monochromator mechanics

- Remove the cover of the monochromator, holding the crystal and rotation system, by loosing in S1.
- Remove the whole mechanics, by loosing attachment screws on the tube shield in S2. **Be careful**, the system is quite heavy. During the loosing the mechanical safety shutter must be jumped up [3]. **Be careful** to hold it during the procedure.

Only the monochromator housing and the slits will be dismount to set on the mirror assembly.

- Slits are removed by loosing S4
- Monochromator housing is dismounted by loosing the M4 screw (S5) defining the rotation axis, in rear.
- The pipe is dismounted as following scheme.





- S6 allow to loose the "tube shield interface. It is necessary to remove the "rear cover".
- Instead of the "tube shield interface" will be set the "mirror pipe".
- In the next part are described the several possibilities of configuration. In case of it is necessary to remove the crystal itself, loose the S8 screws on the side.

2b- Setting the Mirror/Monochromator mechanics

- Attach "A" the optics holder on the tube shield, with S2 screws. Before tightening the plate, set the "tube shield interface" and the "mechanical safety shutter" [3].
- Attach the "mirror housing" "B", by screwing S9.
- Attach the part "C" (Monochromator housing and slits system), with "S5" in "H1". Take care to introduce the mirror pipe on the "Monochromator housing"
- Fix the slit on the new pipe assembly in S4.





Rear side of "C" part:

- The screw on H2 must fix :
- in "I" if the mirror only is used,
 - or in "II" if the assembly mirror/monochromator is used

2c- DIFFERENT KIND OF SETUP

	 Figure 6 : Mirror + monochromator setup. The angle between the beam coming from the mirror, and the beam exiting to the monochromator is 27.28°. This angle is fixed by using 3 screws (ref.12 on fig.4) : 1 screw for the axis of the monochromator housing 1 screw maintaining the arm 1 handled screw to remove when we wish to switch from one mode to the other. There are only 2 positions allowed : straight or inclined at 27.28°.
--	---





2d- Adjustements

2da- The mirror

The mirror can be aligned by using the screws 1 (S1) and the micro-screw 9 (MS9) (figure 4). We usually need the MS9. S1 are used for adjusting the horizontally. MS9 is used for tilting the mirror, in order to get the reflected beam. To adjust it, it is required to remove the monochromator housing cover (using 11), and replace it with the cover containing the fluo-screen (see tool box). When it is fixed, put X-ray on, and check with a counter that everything is safe. Turn MS9, to see the 2 differents beams : at low incidence (on the top-left of the fluo-screen) the direct beam must be seen. By turning MS9, it must disappeared and the reflected beam must appears at 1.15° of incidence (The print must be displaced to the bottom right). When the print of the reflected beam is selected (figure 9–3), switch off the X-ray, and replace the fluo-screen cover by the Germanium holder. Refers to the monochromator adjustment for the second part.



Figure 9: Prints appearing when MS9 is rotated counter clockwise. (1) is the full direct beam (top), (2) are prints of direct beam (top) and weak reflected beam (bottom), (3) the full reflected beam.

2db- The monochromator

Put the fluo-screen on the fixed sample holder, in order to see the beam coming through the slits. Open in large the vertical slit. Switch on X-rays. Turn the MS8 on the monochromator housing until the beam is found. The different steps of the adjustment are described in figure 10. The best adjustment correspond to the scheme (3).



Figure 10 : Prints appearing when MS8 is rotated counter clockwise. (1) is the beginning of k a 2, (2) is fully k a 2 with a small amount of k a 1, (3) is fully k a 1, and (4) is a ended part of k a 1. Then the beam must be adjusted as (3).

2dc- Centering the beam between the slits

The following procedure is to close gradually the slit (Fv), and check that the beam is not "cut" by the slit. If it is the case, move slightly up and down the pipe in order the center the beam on the slit axis.

2dd- Centering the whole optics on the diffracting axis

Please refers to the chapter : **ALIGNING THE DIFFRACTOMETER** and follows the section "Centering the optics on the diffracting axis "

2e- Mirror characteristics

Mirror characteristics:

Equal distance	100mm
Focal distance	10011111
Dimensions (L x l x e)	40 x 20 x 10mm
Incidence angle	1.15°
Outcoming beam width	0.85mm
Capture angle	0.54°
d-spacing range	32.9 – 40.3 Å
Average reflectivity	65 ± 5 %
FWHM of rocking curve	0.042 à0.075°

