

driving X-ray innovation



PX Scanner™ + SuperNova™

Unique *In-situ* X-ray
and Optical Imager

Highest Quality
& Reliability Protein System



PX Scanner


The PX Scanner is Oxford Diffraction's unique instrument for the X-ray screening of protein crystals, *in-situ* in the multi-well crystallisation plate.

A combined optical imager and X-ray system, the PX Scanner allows the visualisation, identification and X-ray evaluation of protein crystals undisturbed within the growth media.

An unlimited number of crystals can be queued for X-ray analysis, either within one well or over an entire plate. The PX Scanner then automatically processes the queue, leaving the user free to continue with other crystallisation work before reviewing their diffraction results.

By removing the manual handling of potential crystals during screening, the PX Scanner can accelerate the study of crystals at all stages of the crystallisation work flow, including:

- High throughput crystallisation condition screening and
- Crystallisation scale-up for diffraction crystals



Multi-well crystallisation plate and Phi Axis


BENEFITS OF THE PX SCANNER

- Early stage identification of crystalline material
- Discrimination of salt from protein crystal
- Identification of optimal crystal growth conditions
- Comparison of crystal diffraction quality
- Automated crystal screening leading to accelerated results
- Study of crystals too small and fragile to extract

PHI AXIS

Having identified a suitable crystal using *in-situ* diffraction, the crystal can then be extracted from the well and mounted atop the Phi Axis. Accepting protein crystals in either loops or capillaries the Phi Axis attachment is inserted into the front of the PX Scanner and allows full X-ray data collection. The Phi Axis has a rotational range of up to 360 degrees and a resolution limit of up to 2.0Å

APPLICATIONS



Optical and *in-situ* screening

Comparison of crystal diffraction quality

Small, fragile crystals

Structure determination

Screen sensitive, membrane proteins

The complete PX Scanner system including
monitor / keyboard table



PX Scanner



PX Scanner - Loading a crystallisation plate

THE SYSTEM

The PX Scanner is unique as the world's only commercially available system for *in-situ* X-ray screening of protein crystals. Specially designed for this purpose, the PX Scanner is a compact system occupying a footprint of approximately 0.5m² (not including separate monitor table). Contained within a fully interlocked radiation proof cabinet, the PX Scanner consists of an x-y-omega goniometer table and crystallisation plate carriage, a high intensity X-ray source and a 165 mm diameter, large active area CCD detector.

Entirely self contained, the PX Scanner requires only a low power mains electricity and water supply, having a built in control computer, system interface, X-ray generator and compact internal chiller units.

For ease of integration into high throughput research laboratories, the PX Scanner comes complete with a dedicated barcode reader allowing the quick and easy identification of the plate prior to loading.

OPTICAL IMAGING

The x-y-omega goniometer table positions the crystallisation plate beneath a series of digital microscopes and cool LED lights for optical imaging of the droplets within the well. Using Oxford Diffraction proprietary algorithms, a series of optical images are recorded of differing focal depth and lighting. A single optimised composite image results.

IN-SITU X-RAY IMAGING

Housed within the base of the system the high intensity Nova™ X-ray source is orientated to provide a vertical X-ray beam. Having optically imaged the plate and a potential crystal selected for *in-situ* X-ray study, the plate and crystal are accurately positioned directly over the X-ray beam which shoots vertically up through the horizontal plate and crystal. The plate and thereby crystal can be tilted up to 6° using the goniometer, allowing scanned X-ray images to be collected. The diffraction data is collected above the plate by an Oxford Diffraction 165 mm CCD detector.

PLATE COMPATIBILITY

Highly versatile, the PX Scanner accepts virtually all crystallisation plates, micro-fluidic chips and counter-diffusion set-ups which fully conform to SBS format.



View inside the PX Scanner showing the:

- X-ray safety panel
- System interface
- System computer
- X-ray generator
- X-ray and CCD chiller

SuperNova

As the latest protein diffraction system from Oxford Diffraction, the SuperNova is compact and easy to use, providing the highest quality crystallographic data from a system of the highest reliability.

SuperNova is the ideal companion to the PX Scanner. Together the two systems provide all you need for the complete protein crystallography laboratory.

The SuperNova combines Oxford Diffraction's high intensity Nova X-ray micro-source with the 165 mm fast, high performance Titan™ CCD detector. Mounted on a 4-circle goniometer, the copper radiation Nova X-ray source provides up to 3x the intensity of a 5kW rotating anode generator with optics and the data quality of a micro-focus rotating anode.

Occupying a 0.5m² footprint, the SuperNova's stylish design of safety interlocked and radiation proof cabinet is entirely self contained and incorporates the X-ray generator, system interface, system computer and two compact water chiller units. For installation, the SuperNova requires only a low power mains electricity and water supply.



Ease of sample mounting on the SuperNova. Configuration shown includes the optional second wavelength X-ray source

KAPPA PLATFORM

With its 4-circle kappa goniometer, the SuperNova is ideally suited for fast, strategy based data collections as well as high resolution, in-house phasing experiments. The kappa goniometer incorporates a motorised sample to detector distance for the study of large cell proteins and a rotatable beam stop for ease of crystal mounting using cryo-tools.


OPTIONAL CRYO-DEVICES

The SuperNova goniometer has been designed to accept all major open flow cryogenic sample cooling devices, including Oxford Diffraction's liquid nitrogen device (90 - 300K).


APPLICATIONS




Molecular Replacement (MR)




Screening



Ligand binding



SAD Phasing or Heavy Atom Phasing



Small Molecules

Highest Quality
& Reliability



The complete SuperNova system including
monitor / keyboard table and optional cryo-device

SuperNova




View of the SuperNova goniometer

NOVA X-RAY MICRO-SOURCE

The SuperNova has been designed around the low power consumption yet high intensity Nova X-ray micro-source.


Up to 3x more intense than a 5 kW rotating anode and with more than 30 systems installed, the Nova is the best selling, most reliable micro-source from any manufacturer in single crystal X-ray diffraction. The Nova combines a state-of-the-art multi-layer optic with a micro-focus X-ray tube running at 50W. Power to the Nova is provided by a software controlled X-ray generator mounted in the base of the SuperNova system.

- Copper wavelength, 50W X-ray micro-source
- Multi-layer X-ray optic
- Up to 3x more intense than a 5kW rotating anode with optic
- Approximately 300 micron beam size



Nova X-ray source


TITAN CCD DETECTOR



The SuperNova incorporates Oxford Diffraction's high performance, fast CCD detector, Titan. At 165 mm in diameter, the Titan is the largest single chip CCD detector available today and is ideal for high resolution protein data collection using fine slicing.

- Active area Ø165 mm
- Gain, 24 e-/X-ray (Cu)
- Precision, 18 bit
- Speed, 0.28 sec (duty cycle)
- Low noise

ATLAS CCD DETECTOR



For studies requiring higher sensitivity combined with large active area, the SuperNova can be configured with the 135 mm Atlas™ CCD.

- Active area Ø135 mm
- Gain, 65 e-/X-ray (Cu)
- Precision, 18 bit
- Speed, 0.28 sec (duty cycle)
- Low noise

SuperNova CrysAlis^{Pro} Software

The SuperNova comes complete with CrysAlis^{Pro}[™], Oxford Diffraction's intelligent data collection and data processing software for protein crystallography. Designed around an easy to use graphical user interface, CrysAlis^{Pro} can be operated under fully automatic, semi-automatic or completely manual control.

AUTOMATIC CRYSTAL SCREENING


At the heart of CrysAlis^{Pro} are the automatic crystal screening and data collection strategy software modules. For a typical crystal, a pre-experiment lasting less than 5 minutes is recorded to evaluate the crystal quality. From the first frame of data, CrysAlis^{Pro} automatically evaluates the crystal quality and provides the user with information regarding the unit cell, angular resolution and suggested frame exposure time.

FASTEST STRATEGY SOFTWARE

Following the pre-experiment, the sophisticated CrysAlis^{Pro} strategy software automatically calculates the optimal data collection conditions for fast, high quality, complete data. All strategies are calculated based on the specific crystal orientation and unit cell and the user has complete control to optimise the data collection for a wide variety of targets including multiplicity, time and resolution. Strategy calculations are extremely fast and efficient allowing the user to quickly adapt the data collection conditions.

AUTOMATIC AND CONCURRENT DATA REDUCTION

Data reduction and processing initialise automatically with the start of data collection and intelligent routines tune the parameters for the best data quality. Since data reduction is concurrent with data collection, processed data are always available and accompanied by real time on-screen feedback of data quality and completeness.



Main CrysAlis^{Pro} GUI

SPECIALIST CRYSTALLOGRAPHIC TOOLS

In addition to the automatic routines, CrysAlis^{Pro} provides hundreds of specialist tools and functions for dealing with non-standard and problematic crystals, including:

- Single image screening
- Advanced unit cell finding
- Reciprocal space viewer
- Space group determination
- Twinning
- Remote device control (X-ray generators, Cryo-devices)
- Axial photos

COMPATIBILITY AND SOFTWARE LICENCES


CrysAlis^{Pro} automatically outputs data in MTZ format and interfaces seamlessly with the CCP4 suite and third-party data reduction packages including MOSFLM, XDS and HKL2000. CrysAlis^{Pro} is provided with the SuperNova under a multi-site, multi-user licence.

PX Scanner CRYSTALEYES Software

Once the multi-well plate has been loaded, user interaction with the PX Scanner is entirely via the CRYSTALEYES™ software.

EASY TO USE, GRAPHICAL INTERFACE

CRYSTALEYES can be used online to control the PX Scanner or offline for the review of previously recorded data. The easy to use graphical user interface (GUI) uses a schematic of the crystallisation plate for user selection of wells. Dependent upon the type of plate, one or more drops are displayed, together with a list of all data recorded, be it optical or X-ray. The user can then select any part of this data for display as a large image within the GUI.



Main CRYSTALEYES GUI in plate mode

OPTICAL IMAGING

Using CRYSTALEYES the PX Scanner can optically image either the entire plate or individual user selected wells / drops. All tasks are added to an on-screen queue allowing the PX Scanner to be used either in real time manual mode or unattended automatic mode.

MOUSE SELECTION OF OBJECTS FOR X-RAY

Whilst reviewing optical images, objects for *in-situ* X-ray are easily added by simply clicking on the object using the mouse cursor. The software automatically queues the task, centres the crystal on the X-ray beam and collects data.

Once recorded, X-ray images can be analysed within CRYSTALEYES to provide the following information:

- Angular resolution limit of diffraction
- Possible unit cell parameters

SECURE LOGIN AND DATA STORAGE

CRYSTALEYES operates a sophisticated, custom database which uses secure user login to control access to personal and group data.


EASE OF EXPORT OF RESULTS

All results from CRYSTALEYES can be exported as a concise html report with more in depth crystallographic analysis available using Oxford Diffraction's CrysAlis^{Pro} software.

PHI AXIS SOFTWARE

CRYSTALEYES supplies a simple step by step guide to Phi Axis data collections, including:

- Crystal mounting and alignment using the PX Scanner's video microscope
- Crystal screening for evaluation of the unit cell
- Automatic strategy based data collection
- Automatic and concurrent data reduction and finalisation




CRYSTALEYES GUI in Phi Axis mode

System dimensions



Suggested system layout (in mm)



Environmental requirements for each system

Power	240 /110 V (switchable) \pm 10% / 16 A fuse protected
Water	2 l/min flow required / 3 -5 bar gauge / 10 – 20 °C range / drain required
Relative Humidity	<80 % non-condensing
Weight loading	Able to bear system weight of 300 kg
Temperature	18 – 23 °C

Please refer to the system user manual for the most up to date and complete requirements.

For more information on Oxford Diffraction or its products please contact us at: sales@oxford-diffraction.com

Europe UK: +44 (0)1235 443630
 North America: +1 540 443 9272
 Latin America: +1 540 443 9272
 Asia Pacific Hong Kong: +852 2571 9188

For other sales offices and dealers throughout the world - please visit our website at:
www.oxford-diffraction.com or www.varianinc.com



Chromatography • Spectroscopy • Mass Spectrometry • Magnetic Resonance Spectroscopy and Imaging • X-Ray Crystallography • Dissolution • Consumables • Data Systems • Vacuum
 Oxford Diffraction reserves the right to change product specifications without notice.

© Oxford Diffraction Ltd, 2008. All trademarks, copyrights and registrations acknowledged. All rights reserved. Printed by Real Affinity Agency in England.

SI-1617/A/07/08