

ALTRUA™ 60 Pacing System

Specifications

ALTRUA 60: Models S601, S602, S603, and S606

The ALTRUA 60 pacing system from Boston Scientific offers multiple features that can be specifically programmed to tailor therapy for patient needs:

- Multiple AV Delay programming options (fixed, dynamic, and AV Search Hysteresis – now extendable to 400 ms) provide increased flexibility in minimizing unnecessary RV pacing
- Minute Ventilation (MV) Blended Sensor is designed to provide a physiologic response for various levels of work in patients of all ages, and restores chronotropic competence¹
- Stored Onset EGMs provide a snapshot of the patient’s rhythm before, during, and after a triggered event without sacrificing battery longevity
- AutoLifestyle automatically adjusts the blended sensor to ensure appropriate rate response for each individual patient
- Automatic Capture enhances patient safety and projected device longevity over the same model^{2,3}
- Ventricular Rate Regulation (VRR) for atrial arrhythmia management

All models offer atrial arrhythmia management features and leading-edge diagnostics. Ease-of-use tools such as Auto Sense and Quick Check help to streamline patient follow-up visits.



MECHANICAL SPECIFICATIONS

| Model | Type | Size (mm) (H × W × D) | Mass (g) | Volume (cc) | Projected Longevity (yrs) ³ | Connector |
|-------------|------|--------------------------|-------------|----------------|---|-----------------------------|
| S601 | SR | 42 × 42 × 8 | 23.4 | 10.0 | 8.6 | IS-1 |
| S602 | DR | 49 × 43 × 8 | 29.6 | 12.6 | 8.8 | IS-1 compatible / 3.2 mm |
| S603 | DR | 44 × 42 × 8 | 25.4 | 10.8 | 6.5 | IS-1 |
| S606 | DR | 49 × 43 × 8 | 29.6 | 12.1 | 8.8 | IS-1 |

ALL MODELS OF ALTRUA 60

| | |
|-----------------------|---|
| Shape | Modified elliptical |
| Envelope | Hermetically sealed titanium |
| Sensors | Minute-ventilation sensor, integrated circuit accelerometer |
| Power Supply | 2.8-V solid-state lithium-iodine battery |
| Setscrew Style | Preinserted captive setscrews and seal plugs |
| Lead Barrel | Various lead connectors accept IS-1 and 3.2 mm leads (IS-1 refers to the international standard ISO 5841.3:1992.) |

*Select Auto to enable Automatic Capture under V Pulse Amplitude parameter options

BRADY ARRHYTHMIA PACING

| Parameter | Programmable Range (Increments) | Nominal by Device Type | |
|-------------------------------|--|------------------------|-----|
| | | DR | SR |
| Mode | DDD(R), DDI(R), DOO(R), VDD, VVT, VVI(R), VOO(R), AAT, AAI(R), AOO(R), ODO, OOO, SSI(R), SOO(R), SST, OSO (Modes beginning with O are available in temporary mode only) | DDD | SSI |
| Lower Rate Limit (LRL) | 30–50 ppm (5-ppm), 50–90 ppm (1-ppm), 90–150 ppm (5-ppm) (155–180 ppm [5-ppm], 180–300 ppm [10-ppm], 300–380 ppm [20-ppm] in temporary mode only and only in SSI, SOO, VVI, VOO, AAI, AOO modes) | 60 | 60 |
| Maximum Tracking Rate (MTR) | 80–185 ppm (5-ppm) | 130 | 130 |
| Maximum Sensor Rate (MSR) | 80–185 ppm (5-ppm) | 130 | 130 |
| A or V Pulse Width | 0.05 ms, 0.1–1.0 ms (0.1-ms) | 0.4 | 0.4 |
| A Pulse Amplitude | 0.1–3.5 V (0.1-V), 4.0–5.0 V (0.5-V), 6.5 V | 3.5 | 3.5 |
| V Pulse Amplitude | | | |
| Manual | 0.1–3.5 V (0.1-V), 4.0–5.0 V (0.5-V), 6.5 V | 3.5 | 3.5 |
| Automatic Capture | Auto | | |
| AV Delay (Paced) ¹ | 0–400 ms (10-ms) | DYN | ■ |

SENSORS

| Parameter | Programmable Range (Increments) | Nominal by Device Type | |
|--------------------------------------|---|------------------------|-----|
| | | DR | SR |
| Minute Ventilation | | | |
| Single chamber | OFF, ON, 4→ON | ■ | OFF |
| Dual chamber | OFF, ON, 4→ON–A, 4→ON–V | OFF | ■ |
| Response Factor | Passive, 1–16 (1) | 3 | 3 |
| High Rate Response | OFF, 55%, 70%, 85% | 70% | 70% |
| High Rate Break Point | 80–185 ppm (5-ppm) | 110 | 110 |
| Accelerometer | ON, OFF, ATR only | OFF | OFF |
| Activity Threshold | V-low, Low, Med-low, Medium, Med-high, High, V-high | Med | Med |
| Reaction Time | 10–50 sec (10-sec) | 30 | 30 |
| Response Factor | Passive, 1–16 (1) | 8 | 8 |
| Recovery Time | 2–16 minutes (1-minute) | 2 | 2 |
| Time-Dependent Blended Sensor | ON, OFF | OFF | OFF |
| AutoLifestyle | ON, OFF, RESET | ON | ON |
| 4-Minute Fast Walk Within 30 Minutes | YES, NO | YES | YES |

RATE ENHANCEMENTS

| Parameter | Programmable Range (Increments) | Nominal by Device Type | |
|---------------------|------------------------------------|------------------------|-----|
| | | DR | SR |
| Hysteresis Offset | OFF, -5 to -80 ppm (5-ppm) | OFF | OFF |
| Search Hysteresis | OFF, 256–4096 cycles (powers of 2) | OFF | OFF |
| Dynamic AV Delay | ON, OFF | ON | ■ |
| Maximum AV Delay | 20–400 ms (10-ms) | 150 | ■ |
| Minimum AV Delay | 10–290 ms (10-ms) | 80 | ■ |
| Sensed AV Offset | OFF; -100 to -10 ms (10-ms) | -30 | ■ |
| AV Search Interval | OFF; 32–1024 cycles (powers of 2) | OFF | ■ |
| AV Delay Increase | 10–100% (10% increments) | 30 | ■ |
| PVARP (fixed) | 150–500 ms (10-ms) | 250 | ■ |
| Dynamic PVARP | ON, OFF | ON | ■ |
| Maximum PVARP | 160–500 ms (10-ms) | 250 | ■ |
| Minimum PVARP | 150–490 ms (10-ms) | 240 | ■ |
| PVARP after PVC/PAC | OFF, 150–500 ms (50-ms) | 400 | ■ |

ATRIAL ARRHYTHMIA MANAGEMENT

| Parameter | Programmable Range (Increments) | Nominal by Device Type | |
|---|---|------------------------|-----|
| | | DR | SR |
| Ventricular Rate Regulation (VRR) | ON, OFF | ON | OFF |
| (VRR) Maximum Pacing Rate | 60–150 ppm (5-ppm) | 110 | 110 |
| Sudden Bradycardia Response (SBR) | ON, OFF | OFF | ■ |
| SBR Detect Time | 1–15 minutes (1 minute) | 5 | ■ |
| SBR Number of Beats | 1–8 cycles (1) | 4 | ■ |
| SBR Therapy Duration | 1–15 minutes (1 minute) | 10 | ■ |
| SBR Therapy Rate Offset | 5–40 ppm (5-ppm) | 5 | ■ |
| SBR MV Offset | OFF, 10%–50% (10% increments) | OFF | ■ |
| Atrial Tachy Response (ATR) | ON, OFF | ON | ■ |
| Trigger Rate | 100–200 ppm (5-ppm) | 170 | ■ |
| Fallback Mode | VDI(R), DDI(R) | VDI | ■ |
| Duration | 0, 8–2048 cycles (powers of 2) | 8 | ■ |
| Fallback Time | 0–120 sec (5-sec) | 30 | ■ |
| ATR Entry Count | 1–8 cycles (1) | 8 | ■ |
| ATR Exit Count | 1–8 cycles (1) | 8 | ■ |
| ATR Lower Rate Limit | 30–50 ppm (5-ppm), 50–90 (1-ppm), 90–150 (5-ppm). ATR Lower Rate Limit must be equal to or greater than the permanent Lower Rate Limit. | 70 ppm | ■ |
| Atrial Flutter Response | OFF, 130–230 ppm (10-ppm); DDI(R) mode only: 230 ppm | OFF | ■ |
| Rate Smoothing programmable for increase and decrease | OFF; 3%–24% (3% increments). Separately | OFF | OFF |

SENSITIVITY ADJUSTMENT

| Parameter | Programmable Range (Increments) | Nominal by Device Type | |
|-------------------------|--|------------------------|-----|
| | | DR | SR |
| Atrial Sensitivity | Auto, 0.15, 0.25, 0.5, 0.75, 1.0–8.0 mV (0.5-mV), 9.0, 10.0 mV | 0.75 | ■ |
| Ventricular Sensitivity | Auto, 0.25, 0.5, 0.75, 1.0–8.0 mV (0.5-mV), 9.0, 10.0 mV | 2.5 | 2.5 |

LEAD CONFIGURATION

| Parameter | Programmable Range (Increments) | Nominal by Device Type | |
|---------------------------|---------------------------------|------------------------|----|
| | | DR | SR |
| A or V Lead Configuration | Unipolar, Bipolar, Split | BI | BI |

REFRACTORY

| Parameter | Programmable Range (Increments) | Nominal by Device Type | |
|-------------------------|---------------------------------|------------------------|-----|
| | | DR | SR |
| A Refractory Period | 150–500 ms (10-ms) | 300 | ■ |
| V Refractory Period | 200–500 ms (10-ms) | 250 | 250 |
| A Blanking after V Pace | 30–200 ms (10-ms) | 120 | ■ |
| V Blanking after A Pace | 30–200 ms (10-ms) | 40 | ■ |

OTHER FEATURES

| Parameter | Programmable Range (Increments) | Nominal by Device Type | |
|---------------------------|---------------------------------|------------------------|-------|
| | | DR | SR |
| PMT Termination | ON, OFF | ON | ■ |
| Magnet Response | OFF, ASYNC, EGM | ASYNC | ASYNC |
| A or V Lead Safety Switch | ON, OFF, RESET | OFF | OFF |
| Runaway Protection | Not Programmable (ppm) | 210 | 210 |

¹ Chronotropic competence is defined by: Wilkoff BL, Corey J, Blackburn G. A mathematical model of cardiac chronotropic response to exercise. J Electrophysiol. 1989;3(3):176–180. Refer to Physician's System Guide for more information on adaptive-rate therapy. Additional clinical performance was assessed using INSIGNIA Ultra clinical data with the AutoLifestyle feature programmed On. Data on file.

² Back up pulses are 1.5 V above the last threshold with a minimum of 3.5 V and maximum of 4.5 V.

³ Longevity projections as described in device user manual. Settings: 60 ppm, A=2.5 V, V=1.0 V, 500 ohms, 100% paced, **MV Blended Sensor ON, Onset EGMs ON, Automatic Capture ON.**

Pacing Systems from Boston Scientific CRM

Indications

Pacemaker indications include: symptomatic paroxysmal or permanent second- or third-degree AV block; symptomatic bilateral bundle branch block; symptomatic paroxysmal or transient sinus node dysfunction with or without associated AV conduction disorders; bradycardia-tachycardia syndrome, to prevent symptomatic bradycardia or some forms of symptomatic tachyarrhythmias; neurovascular (vaso-vagal) syndromes or hypersensitive carotid sinus syndromes. Adaptive-rate pacing is indicated for patients who may benefit from increased pacing rates concurrent with increases in minute ventilation and/or level of physical activity. Pacemakers' dual-chamber and atrial tracking modes are also indicated for patients who may benefit from maintenance of AV synchrony. Dual-chamber modes are specifically indicated for: conduction disorders that require restoration of AV synchrony, including varying degrees of AV block; VVI intolerance (eg, pacemaker syndrome) in the presence of persistent sinus rhythm.

Contraindications

Pacemakers are contraindicated for the following patients under the circumstances listed: patients with unipolar pacing leads or in MV mode with an implanted ICD because it may cause unwanted delivery or inhibition of ICD therapy; use of the MV sensor in patients with only unipolar leads, because a bipolar lead is required in either the atrium or the ventricle for MV detection; single-chamber atrial pacing in patients with impaired AV nodal conduction; atrial tracking modes for patients with chronic refractory atrial tachyarrhythmias, which might trigger ventricular pacing; dual-chamber and single-chamber atrial pacing in patients with chronic refractory atrial tachyarrhythmias; asynchronous pacing in the presence (or likelihood) of competition between paced and intrinsic rhythms.

Warnings

Read the product labeling thoroughly before implanting the pulse generator to avoid damage to the system. Such damage can result in patient injury or death. Inappropriate sustained high-rate pacing occurred in the PULSAR MAX clinical study in 5 out of 130 patients with MV ON, 4 to 14 days after implant. If sustained high-rate pacing could be of concern, consider programming a reduced Max Sensor Rate or MV to Passive. These programming recommendations are intended to assure that MV calibration is evaluated and, if necessary, recalibrated (4→ON) when the patient and pacing system have stabilized post implant. Continued monitoring of the MV sensor performance should be performed at all follow-up visits until implant stabilization has occurred.

Precautions

For information on precautions, refer to the following sections of the product labeling: clinical considerations; sterilization, storage and handling; lead evaluation and connection; implantation; programming and pacemaker operation; MV initialization; environmental and medical therapy hazards. Advise patients to avoid sources of electric or magnetic interference (EMI). If the pacemaker inhibits or reverts to asynchronous operation at the programmed pacing rate or at the magnet rate while in the presence of the EMI, moving away from the source or turning it off will usually allow the pulse generator to return to its normal mode of operation.

Potential Adverse Events

Potential adverse events from implantation of the pacing system include, but are not limited to, the following: allergic/physical/physiologic reaction, death, erosion/migration, fibrillation or other arrhythmias, lead or accessory breakage (fracture/insulation/lead tip), hematoma/seroma, inappropriate or inability to provide therapy (pacing/sensing), infection, procedure related, and component failure. In rare cases severe complications or device failures can occur.

Refer to the product labeling for specific indications, contraindications, warnings/precautions and adverse events. Rx only.
(Rev. K)

**Boston
Scientific**

Delivering what's next.™

Cardiac Rhythm Management

Boston Scientific Corporation
4100 Hamline Avenue North
St. Paul, MN 55112-5798 USA
Tel: 651.582.4000 Fax: 651.582.4166
Medical Professionals:
1.800.CARDIAC (227.3422)
Patients and Families:
1.866.484.3268
www.bostonscientific.com

Copyright © 2009 by
Boston Scientific Corporation
or its affiliates. All rights reserved.

C2-180-0409