

CARDELL® 9500 HD Veterinary Monitor

User's Guide



For Veterinary Use Only
003-2566-00 Rev. C (4/4/13)

Product Information

Product model: CARDELL® 9500 HD

Product name: Veterinary monitor

Explanation of Symbols

The following symbols appear on the monitor and its packaging.

Table 1 Monitor and Packaging Symbols




Symbol	Description	Symbol	Description
	Keep dry		Caution, consult accompanying documents
	Atmospheric pressure limitation		Temperature limitation
	Keep upright		Humidity limitation
	Recovery		Fragile, handle with care
	Recyclable		Maximum stacking
	Waveform Freeze key		On/Standby key
	Alarm Silence key		Print key
	NIBP history key		NIBP key
	AC Power LED		Charging LED
ECG/RESP	ECG connector	SpO₂	SpO ₂ probe connector
NIBP	NIBP cuff connector	CO₂	CO ₂ connector

Table 1 Monitor and Packaging Symbols (Continued)

Symbol	Description	Symbol	Description
TEMP	TEMP connector		Dangerous Voltage
	CF applied part: including F applied part (float/insulation) and defibrillation-proof function.		BF applied part: including F applied part (float/insulation) and defibrillation-proof function.
	Reserved port		Fuse
	Protective grounding		Serial Number
	Equipotential grounding		Manufacturer address
	Date of manufacture		Configuration number
	Batch code		Consult Instructions for Use
	Catalogue number		Compliance with WEEE standard
IPX1	Protective grade		USB port/ USB Export Menu

Safety Conventions

The manual uses the following conventions for Notes, Cautions, and Warnings.

Note — A **Note** calls attention to an important point in the text.

Caution A **Caution** calls attention to a condition or possible situation that could damage or destroy the product or the user's work.

Warning A **Warning** calls attention to a condition or possible situation that could cause injury to the user and/or patient.

Safety Requirements

Note — The safety indications in this chapter apply to general monitor use. Safety indications in other chapters apply to specific monitor measurements.

Follow the instructions in this user manual when using the monitor. However, conventional medical practices always supersede this document. Significance of safety requirements set forth here in this manual is not in order of reading sequence.

Warning The monitor is not for home use.

The monitor is not for diagnostic or therapeutic use.

To avoid mixture, clear all historical data of the last patient before monitoring a new patient.

Ensure that the veterinary monitor is in good working condition and is placed in proper position before clinical use.

Never rely exclusively on the alarm system for veterinary monitoring. You must periodically check that monitor alarms are working properly. The most reliable method of veterinary monitoring combines close, personal surveillance along with the correct operation of the monitor.

Never use the monitor during MRIs or CT scans.

Never use cables with exposed conductors. Only use undamaged cables and connectors to avoid personal or patient injury.

Never simultaneously touch the patient and the input or output terminal of the monitor.

The monitor can be used with electrosurgical devices. To avoid burns or death, ensure that the electrosurgical circuitry is properly connected.

Ensure that the patient is not connected to any metal, the surface of any conductors, or grounded parts of any devices when a defibrillator is in use. Never touch the animal, table, or instruments during defibrillation.

Never use the monitor in an environment of uncontained anesthetic gases.

Disconnect the AC power supply and remove the battery before disassembling the monitor.

Caution

To avoid personal or patient injury, only use accessories and parts produced or recommended by the manufacturer. Otherwise, damage to the monitor can occur.

Position power cords and other cables to prevent tangling, patient contact, or electrical interference.

Always use a three-wire receptacle in order to properly ground the monitor. Use the monitor on battery power only if you suspect any problem with external grounding or wiring.

The battery will discharge over time when stored inside a monitor with no AC power. Keep the monitor plugged in when storing it and perform a monthly check of battery power.

If multiple devices are connected to a patient, the sum of the leakage currents may cause a dangerous situation. Therefore, qualified service representatives should always perform current leakage tests before setting up the connections. This ensures that leakage is within limits and prevents personal or patient injury and environmental damage. If any question remains, consult the manufacturer for correct usage.

Unknown electromagnetic interference could be caused by radio transmitters or nearby televisions. Move veterinary monitors away from these devices or insert shielding materials around the monitor.

If patient data is lost, closely monitor the patient or replace the monitor immediately until the monitor restores normal function again.

When the monitor is powered off, patient monitoring will be stopped and user configurations will be saved.

To reduce the hazard of burns in high-frequency surgical neutral connections, the leads and connectors must be located away from the surgical site. Keep cables away from other devices.

Periodically check all reusable accessories for damage. Replace and dispose of damaged accessories when necessary according to your local hospital waste disposal regulations.

ECG electrodes are disposable accessories. Always properly dispose of ECG electrodes according to your local hospital waste disposal regulations.

Clean and sterilize the monitor and accessories according to local requirements. Turn off the monitor and disconnect the power cords before cleaning or sterilization.

Keep all monitor packing materials away from children, or dispose of them in accordance with your local environmental regulations.

Always properly dispose of the monitor and all accessories at the end of their service life. Dispose of batteries according to your local regulations. Never incinerate batteries or expose them to high temperatures.

Ensure that no water condenses into or on the monitor. Condensation can occur from changes in temperature or exposure to humidity.

Contents

1 Overview

1.1 Indications for Use - - - - -	1-1
1.2 Intended Use- - - - -	1-1
1.3 CARDELL® 9500 HD Configurations - - - - -	1-2
1.3.1 Main Components - - - - -	1-2

2 Basic Operation

2.1 Before You Begin - - - - -	2-1
2.2 Setting up the Monitor - - - - -	2-1
2.2.1 Power Supply - - - - -	2-2
2.2.2 Inspecting the Monitor - - - - -	2-2
2.2.3 Checking Your Software Release - - - - -	2-2
2.3 Product Description - - - - -	2-2
2.3.1 The Front Panel - - - - -	2-3
2.3.2 The Rear Panel - - - - -	2-4
2.3.3 The Side Panels - - - - -	2-5
2.4 The Screen Display - - - - -	2-6
2.5 Menus - - - - -	2-7
2.6 Configuring the Display - - - - -	2-8
2.6.1 Using the Large Font Screen - - - - -	2-8
2.6.2 Changing the Waveform Display - - - - -	2-8
2.6.3 Changing the Waveform Speed - - - - -	2-8
2.6.4 Freezing the Waveforms - - - - -	2-9
2.7 Changing the Monitor Settings - - - - -	2-9
2.7.1 System Maintenance - - - - -	2-9
2.7.2 Changing the Date and Time - - - - -	2-9
2.8 Changing the Volume and Brightness - - - - -	2-10
2.8.1 Adjusting the HR and PR Volume - - - - -	2-10
2.8.2 Adjusting the Brightness - - - - -	2-10
2.9 Managing Patients - - - - -	2-10
2.10 Exporting ECG Waveform and Trend Data to a USB Disk - - - - -	2-11

3 Alarms

3.1 Alarm Types- - - - -	3-1
3.2 Alarm Levels - - - - -	3-1
3.2.1 High Priority Alarms - - - - -	3-1
3.2.2 Medium Priority Alarms - - - - -	3-1
3.2.3 Low Priority Alarms - - - - -	3-2
3.3 Visual Alarms - - - - -	3-2
3.3.1 LED Flashing Alarms - - - - -	3-2
3.3.2 Flashing Icons - - - - -	3-2
3.3.3 Alarm Messages - - - - -	3-2
3.4 Audible Alarms - - - - -	3-2
3.4.1 Turning on Audible Alarms - - - - -	3-3
3.4.2 Turning off Audible Alarms - - - - -	3-4
3.4.3 Silencing Alarms - - - - -	3-4
3.4.4 Adjusting Alarm Volume- - - - -	3-4
3.5 Setting the Alarm Mode - - - - -	3-5
3.5.1 Automatic - - - - -	3-5
3.5.2 Standard- - - - -	3-5

3.6 Setting Alarm Limits - - - - -	3-5
3.6.1 Setting Individual Alarm Limits - - - - -	3-5
3.6.2 Restoring Default Alarm Limits- - - - -	3-6

4 Monitoring ECG

4.1 Before Monitoring - - - - -	4-1
4.1.1 Skin Preparation for Electrode Placement - - - - -	4-1
4.1.2 Placing the Electrodes - - - - -	4-1
4.1.3 Changing a Lead Set - - - - -	4-3
4.1.4 Lead Recommendations for Surgical Patients - - - - -	4-3
4.1.5 Connecting ECG Cables - - - - -	4-4
4.2 The ECG Display - - - - -	4-4
4.2.1 ECG Waveform and Numeric Pane - - - - -	4-4
4.3 Configuring the ECG Display- - - - -	4-5
4.3.1 Changing the Size of the ECG Waveform - - - - -	4-5
4.3.2 Changing the Speed of the ECG Waveform - - - - -	4-5
4.4 Performing 1 mV Calibration - - - - -	4-5
4.5 Changing the Heart Rate Volume - - - - -	4-6
4.6 Changing the ECG Mode- - - - -	4-6
4.7 Changing HR Alarm Settings - - - - -	4-6
4.8 Configuring HR/PR Priority - - - - -	4-6
4.9 Monitoring Pacemaker Patients - - - - -	4-7
4.10 Safety Information - - - - -	4-8

5 Monitoring Respiration

5.1 Lead Placement for Monitoring Respiration - - - - -	5-1
5.1.1 Selecting a Respiration Lead - - - - -	5-1
5.2 The Respiration Display - - - - -	5-2
5.3 Changing the Respiration Display- - - - -	5-2
5.3.1 Changing the Respiration Waveform Size - - - - -	5-2
5.3.2 Changing the Respiration Waveform Speed - - - - -	5-2
5.3.3 Changing the Respiration Waveform Display Mode - - - - -	5-3
5.4 Changing the Respiration Alarm Settings - - - - -	5-3
5.5 Changing the Apnea Alarm Settings- - - - -	5-3
5.6 Safety Information - - - - -	5-3

6 Monitoring SpO₂

6.1 SpO ₂ Monitoring Theory - - - - -	6-1
6.2 Placing an SpO ₂ - - - - -	6-1
6.3 Connecting SpO ₂ Cables - - - - -	6-2
6.4 Measuring SpO ₂ - - - - -	6-2
6.5 Assessing a Suspicious SpO ₂ Reading- - - - -	6-2
6.6 The SpO ₂ Display - - - - -	6-3
6.7 Changing the SpO ₂ Display- - - - -	6-3
6.7.1 Changing the SpO ₂ Waveform Speed - - - - -	6-3
6.7.2 Changing the SpO ₂ Waveform Display Mode - - - - -	6-3
6.8 Configuring HR/PR Priority - - - - -	6-3
6.9 Changing the SpO ₂ Alarm Settings - - - - -	6-3
6.10 Changing the PR Alarm Settings- - - - -	6-4
6.11 Safety Information - - - - -	6-4

7 Monitoring NIBP

7.1 Selecting and Placing an NIBP Cuff - - - - -	7-1
7.1.1 Cuff Placement for Cat- - - - -	7-1
7.1.2 Cuff Placement for Dog - - - - -	7-2
7.1.3 Cuff Placement for Large Animal- - - - -	7-3
7.1.4 Cuff Size Selection - - - - -	7-3
7.2 Connecting the Cuff - - - - -	7-4
7.3 Measurement Limitations - - - - -	7-4
7.4 Measuring the NIBP - - - - -	7-4
7.5 The NIBP Display - - - - -	7-5
7.6 Changing the NIBP Cuff Size - - - - -	7-5
7.7 Selecting NIBP Measurement Methods - - - - -	7-5
7.8 Selecting the NIBP Interval Time- - - - -	7-6
7.9 Changing the NIBP Unit - - - - -	7-7
7.10 Viewing the NIBP Measurement Data - - - - -	7-7
7.11 Changing the NIBP Alarm Settings - - - - -	7-8
7.12 Calibrating NIBP - - - - -	7-8
7.13 Safety Information - - - - -	7-8

8 Monitoring Temperature

8.1 Selecting and Connecting Temperature Probe - - - - -	8-1
8.2 Influencing Factors - - - - -	8-1
8.3 Monitoring Temperature - - - - -	8-1
8.4 Changing the Temperature Unit - - - - -	8-2
8.5 Changing the Temperature Alarm Settings - - - - -	8-2
8.6 Safety Information - - - - -	8-2

9 Monitoring Carbon Dioxide

9.1 Selecting CO ₂ Accessories - - - - -	9-1
9.2 Zeroing the Sensor- - - - -	9-1
9.3 Measuring Carbon Dioxide- - - - -	9-2
9.3.1 Capnostat® 5 Mainstream CO ₂ Sensor- - - - -	9-2
9.3.2 LoFlo SideStream CO ₂ Sensor- - - - -	9-3
9.4 The CO ₂ Display - - - - -	9-5
9.5 Changing the CO ₂ Display - - - - -	9-5
9.5.1 Changing the CO ₂ Waveform Speed - - - - -	9-5
9.5.2 Changing the CO ₂ Waveform Display Mode - - - - -	9-5
9.6 Changing the CO ₂ Alarm Settings - - - - -	9-5
9.7 Safety Information- - - - -	9-6

10 Managing Data

10.1 Viewing Trends - - - - -	10-1
10.2 Viewing Data List - - - - -	10-2
10.2.1 Viewing All Measurement Data - - - - -	10-2

11 Printing

11.1 About the Printer- - - - -	11-1
11.2 Loading the Printer Paper- - - - -	11-1
11.3 Printing - - - - -	11-2
11.4 Turning on Alarm-triggered Printing - - - - -	11-2
11.5 Setting Timed Printing - - - - -	11-3
11.6 Setting the Print Time- - - - -	11-3
11.7 Setting the Advance Time- - - - -	11-3

11.8 Selecting the Print Mode	- - - - -	11-3
11.9 Changing the Print Speed	- - - - -	11-4
11.10 Sample Print Strip	- - - - -	11-4
11.11 Printer Paper Guidelines	- - - - -	11-4

12 Cleaning and Care

12.1 General Guidelines	- - - - -	12-1
12.2 Cleaning the Monitor	- - - - -	12-1
12.2.1 Cleaning the Display	- - - - -	12-2
12.3 Cleaning the Accessories	- - - - -	12-2
12.3.1 Cleaning the Cables	- - - - -	12-2
12.3.2 Cleaning the SpO ₂ Sensor	- - - - -	12-2
12.3.3 Cleaning the NIBP Cuff	- - - - -	12-2
12.3.4 Cleaning the Temperature Probe	- - - - -	12-3
12.3.5 Cleaning the ECG Clips	- - - - -	12-2
12.3.6 Cleaning the CO ₂ Accessories	- - - - -	12-3
12.4 Cleaning the Print Head	- - - - -	12-4
12.5 Disinfecting the Monitor	- - - - -	12-4
12.6 Sterilizing the Monitor	- - - - -	12-4

13 Maintaining the Battery

13.1 Charging LED	- - - - -	13-1
13.2 Battery Indicator	- - - - -	13-1
13.3 Installing the Battery	- - - - -	13-2
13.4 Recharging the Battery	- - - - -	13-3
13.5 Optimizing the Battery	- - - - -	13-3
13.5.1 Conditioning the Battery	- - - - -	13-3
13.5.2 Checking Battery Performance	- - - - -	13-3
13.6 Disposing of a Battery	- - - - -	13-3
13.7 Safety Information	- - - - -	13-4

14 Maintenance and Troubleshooting

14.1 Inspecting the Equipment and Accessories	- - - - -	14-1
14.2 Maintenance Task and Test Schedule	- - - - -	14-1
14.3 Troubleshooting	- - - - -	14-2
14.4 Disposing of the Monitor	- - - - -	14-5

15 Warranty/Service

15.1 Limited Warranty	- - - - -	15-1
15.2 After-sale Service	- - - - -	15-2

A. Specifications

A.1 ECG Specifications	A-1
A.2 Respiration Specifications	A-2
A.3 SpO ₂ Specifications	A-2
A.4 Temperature Specifications	A-3
A.5 Non-invasive Blood Pressure Specifications-	A-4
A.6 Carbon Dioxide Specifications	A-4
A.7 Power Specifications-	A-5
A.8 Display Specifications	A-5
A.9 Printer Specifications	A-6
A.10 Physical Specifications	A-6
A.11 Environmental Specifications	A-6

B. Alarm Specifications

B.1 Alarm Classifications	B-1
B.2 Physiological Alarms	B-1
B.3 Technical Alarms	B-3
B.4 Default Alarm Limits	B-4

C. Default Settings

C.1 System Default Settings	C-1
C.2 Alarm Default Settings-	C-1
C.3 ECG Default Settings	C-1
C.4 Respiration Default Settings	C-2
C.5 SpO ₂ Default Settings	C-2
C.6 NIBP Default Settings-	C-3
C.7 Temperature Default Settings	C-3
C.8 Carbon Dioxide Default Settings-	C-3
C.9 Printer Default Settings	C-4

D. Glossary

D.1 Terminology Definitions-	D-1
------------------------------	-----

E. Tips

E.1 Blood Pressure	E-1
E.2 ECG-	E-1
E.3 SpO ₂	E-2
E.4 CO ₂	E-2
E.5 Temperature	E-2

F. Accessories

1 Index

1

Overview

The CARDELL® 9500 HD veterinary monitor is used for monitoring, recording, and alarming of multiple physiological parameters of animals. Your monitor stores data in trend database. You can view measurement trend graphs to help you identify changes in the monitored animal's physiological condition.



Figure 1-1 CARDELL® 9500 HD Monitor

This Instructions for Use is an integral part of the monitor and provides detailed introductions on how to use and maintain it. Keep the document next to the monitor for reference. In addition, read and keep any documents that come with any accessories, as these contain important information not in this guide.

This guide describes all of the monitor features and options. Your monitor may not have all of them. What you see on the screen, including the menus, depends on the way it has been configured for your hospital.

1.1 Indications for Use

The CARDELL® 9500 HD veterinary monitor is for use by healthcare professionals whenever there is a need for monitoring the physiological parameters of animals.

1.2 Intended Use

The CARDELL® 9500 HD veterinary monitor is intended to be used in special laboratory and other veterinary hospital or clinic where veterinary monitoring is needed. Additionally, the monitor is also intended for use in transport situations.

1.3 CARDELL® 9500 HD Configurations

The CARDELL® 9500 HD veterinary monitor includes several configurations.

1.3.1 Main Components

The monitor can contain the following components:

- Main unit
- ECG
- RESP
- NIBP
- SpO₂
- TEMP
- CO₂
- Printer

The following table lists the available configurations, the standard features in each monitor, as well as optional features. In the table, a solid circle indicates a standard feature.

Table 1-1 Product Configurations

Parameter	CARDELL® 9500 HD
ECG	•
RESP	•
NIBP	•
SpO ₂	•
TEMP	•
CO ₂	• (with optional sensor)
Printer	•
USB Port	•
Battery	• (One battery)
Display Screen	7"

2

Basic Operation

This chapter describes how to set up and begin using the CARDELL® 9500 HD veterinary monitor.

2.1 Before You Begin

Before veterinary monitoring, inspect the monitor to ensure it is in good condition and used in proper environment.

1. Unpack the equipment and ensure that you have the following:
 - One veterinary monitor
 - Instructions for Use
 - AC power cord (country-specific)
 - Set of accessories
2. Keep the packing case for future transport or storage. If it is damaged, immediately contact the courier company.
3. Examine the monitor for cleanliness and general physical condition. Ensure that:
 - The housing is not cracked or broken.
 - Power plugs and cords are in good condition and plug prongs are not bent.
 - The accessory cables are intact.

Caution If the equipment, accessories or their packaging show signs of damage, do not use the equipment.

2.2 Setting up the Monitor

Place and use the monitor in a location that conforms to the environmental specifications outlined on page A-6. Locate the monitor in an area that:

- is convenient for observation and operation.
- is free from movement, dust, corrosive or explosive gases.
- has a five-centimeter space around to ensure good ventilation and smooth heat dissipation.

You can put the monitor on a flat surface.

Warning Never attempt to open the monitor case. Only qualified personnel should service the monitor.

Never block the monitor's air vent during operation.

2.2.1 Power Supply

You can power the monitor with either an AC power source or its internal battery. You can switch back and forth between supplies without powering off the monitor.

When the monitor is connected to an AC power source, you can turn on and operate the veterinary monitor.

When powered on, the monitor is under normal working condition. It can be used to measure vital signs and charge the battery.

When it is standby, the monitor can only charge the battery.

When AC power is lost, a battery powered monitor continues to run without interruption on battery power. For more information on battery power, see Chapter 13, Maintaining the Battery on page 13-1.

If the monitor is only powered by an AC power source, the monitor will be powered off automatically when the AC power is lost; when powered on again, the monitor is in normal working condition.

2.2.2 Inspecting the Monitor

Before you start using the monitor, ensure that:

1. There is no damage to the monitor or the accessories. Before turning on the power, check the monitor and all the cables, connectors and accessories and ensure that the extension cables are correctly connected to the fittings.
2. The AC power supply is in conformity with standards. Plug one end of the specified power cord into an AC power supply with a protection-grounded wire. Plug the other end into the power connector on the monitor. When the power cord is correctly connected, the AC power LED is yellow. If the monitor is powered by battery, verify that adequate power remains and the battery power gauge is full.
3. The monitor is properly grounded. Connect the small end of the grounding cable to the grounding post on the back of the monitor. The large end of the grounding cable should be connected to any metal surface or copper pipe.
4. The accessories are correctly connected. For further testing of individual measurement alarms, perform the measurement on a simulator (for example, SpO₂). Adjust alarm limits and check that alarm behavior is observed.
5. All settings like cuff size, alarm limits and so on, are correct. Ensure that the monitor is in good working condition and that all functions of the monitor meet the technical requirements mentioned in the Instructions for Use.
6. When you press the **On/Standby** key on the front panel, the main screen appears.

2.2.3 Checking Your Software Release

To check the software release of the monitor:

1. Rotate the navigation knob when no menu displays on the monitoring screen.
2. Rotate the navigation knob to the **System** menu and press the knob. The version number is shown in the top left corner, displaying as Version: Vn.nn.
3. Rotate the navigation knob to **Exit** and press the knob to quit the menu.

2.3 Product Description

The monitor has a color TFT screen that displays waveforms, menus, alarm status and physiological measurements. This section describes all the monitor features and options. Your monitor may not have all of them.

2.3.1 The Front Panel

The following figure shows the front panel of the monitor:




Figure 2-1 CARDELL® 9500 HD Front Panel

The following table describes the controls on the monitor.

Table 2-1 Controls

Item Number	Icon	Meaning	Description
1		On/Standby key	Press and hold for three seconds to turn the monitor on or off.
2		Power LED	Green: the monitor is powered on. Yellow: the monitor is off and connected to an AC power. Off: there is no AC power source.
3		Charging LED	Flashing green: the monitor is connected to AC power, and the battery is charging. Steady green: the battery is fully charged. Off: the monitor is not connected to an AC power or no battery loaded.
4		Waveform Freeze key	Press to freeze/unfreeze the waveform. After 30 seconds, the wave unfreezes (if frozen) and the monitor returns to the previous screen.

Table 2-1 Controls (Continued)

Item Number	Icon	Meaning	Description
5		NIBP key	Press to start or stop an NIBP measurement.
6		Alarm Silence key	Press to turn the alarm audio off for 2 minutes. The monitor resumes alarm audio when the counter reaches 00:00 . Press again when alarms are paused to restart alarm audio.
7		Print key	Press to start or stop printing.
8		NIBP History key	Press to view NIBP history on the monitoring screen.
9		Navigation knob	When a menu is displayed, rotate and press the navigation knob to select and change the settings; When there is no menu on the screen, press to enter the Data Review (history) window.
10		Alarm LED	Flashing red: indicates a high priority alarm or life-threatening situation. Flashing yellow: indicates a medium priority alarm. Steady yellow: indicates a low priority alarm.

2.3.2 The Rear Panel

The following figure shows the rear panel of the monitor. Your monitor may not have all of the following items.




Figure 2-2 CARDELL® 9500 HD Rear Panel

The following table describes each item on the rear panel.

Table 2-2 Rear Panel

Item Number	Item	Name	Description
1		Audio outlet (speaker)	Enables audio output
2		USB port	
3		Reserved port	
4		Signal input/output	Software update port only, DB9 connector
5/6		Fuse socket	T1. 6A, 250V fuse
7		Equipotential grounding lug	Connects to the grounding system
8		External AC power connector	Connects to external AC power supply

2.3.3 The Side Panels

The following figure shows the left and right side panels of the CARDELL® 9500 HD monitor. Your monitor may not have all of the following items.




Figure 2-3 CARDELL® 9500 HD Side Panels

The following table describes each item on the side panels.

Table 2-3 Side Panel

Item Number	Symbol	Description
1	TEMP	TEMP connectors
2	CO₂	CO ₂ connector
3	ECG/RESP	ECG connector
4	SpO₂	SpO ₂ probe connector
5	NIBP	NIBP cuff connector
6		Battery compartment cover
7		Printer

2.4 The Screen Display

The monitor screen displays parameters, waveforms and relevant information about patient measurements. The main screen display might look different, depending on the configuration of your monitor.

The following figure shows a typical screen display.




Figure 2-4 CARDELL® 9500 HD Typical Screen Display

The following table describes the screen display.

Table 2-4 Screen Display Description

Item Number	Description	Item Number	Description
1	Waveform area	2	Numeric panes
3	Menu bar and message area	4	Information bar

2.5 Menus





Main menu and the submenu are two typical menus in the monitor.

When no menu displays on the screen:

1. Rotate the navigation knob to access the **Main** menu.
2. Scroll down the menus to highlight the menu you want and then press the navigation knob to select the menu and make necessary changes.

The following table describes the Main menu indicators:

Table 2-5 Main Menu

Icon	Description	Icon	Description
	Parameter menu		Trend Display menu
	System menu		Volume & Brightness Menu
	Print menu		NIBP History List
	Large Font		USB Export Menu
	Exit		

Submenu is the menu under the **Main** menu.

1. Rotate the navigation knob to access a submenu, then press the navigation knob to a menu item and press the knob.
2. Make your changes then rotate the navigation knob to **Exit** and press the knob to return to the previous menu.

2.6 Configuring the Display

You can configure the display to a variety of layouts.

2.6.1 Using the Large Font Screen



Select the **Large Font** option to enlarge the information on the parameter panes on the screen. This provides a better view of parameters for remote monitoring.

The following figure shows an example of the screen display when **Large Font** is selected:




Figure 2-5 Large Font Screen

To display or exit the Big Number screen:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to highlight **Big Number** and press the knob to confirm.
3. Press the navigation knob again to exit the **Big Number** screen.

Note — The parameter and layout of the **Large Font** screen depends on the parameters monitored.

2.6.2 Changing the Waveform Display

The monitor has three modes of waveform display: single line, bold line and filled. Waveforms display in a single line by default. When **Enhanced** is selected in **Waveform Mode**, waveforms display in a bold line.

To change the waveform display:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob to enter a submenu.
3. Rotate the navigation knob to a parameter (except ECG) and press the knob.
4. Rotate the navigation knob to **Waveform Mode** and press the knob.
5. Rotate the navigation knob to an option (**Enhanced**, **Fill** or **Standard**) and press the knob to confirm the selection.
6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

2.6.3 Changing the Waveform Speed

You can configure the speed at which physiological waveforms are scanned. The unit of waveform speed is **mm/s**. Increasing the waveform speed spreads the waveform scale in the channel for more detailed observation. Reducing

the waveform speed compresses the waveform display so that waveform segments display for longer period of time.

To change the waveform speed:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to a parameter and press the knob.
4. Rotate the navigation knob to **Waveform Speed** and press the knob.
5. Rotate the navigation knob to an option and press the knob to confirm the selection.
6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

For information on the default waveform speed of each parameter, see Chapter C, Default Settings on page C-1.

2.6.4 Freezing the Waveforms



You can freeze the real-time waveforms to check measurements more closely.

When monitoring:

- Press the **Waveform Freeze** key on the front panel to freeze waveforms.

After 30 seconds, the waveforms automatically unfreeze and the monitor returns to the previous screen.

- To unfreeze waveforms, press the **Waveform Freeze** key on the front panel.

2.7 Changing the Monitor Settings

Warning Check for the correct cuff size and alarm settings before changing the monitor settings.

2.7.1 System Maintenance

Use System Maintenance menu to access “Demo” mode for use in training and simulation.
The password is “8727”

2.7.2 Changing the Date and Time

The monitor displays the current date and time in the bottom right corner of the monitoring screen. Once it is set, the monitor saves the information in its internal clock even when powered off. All patient data, including waveforms, trends, and measurements are stored with the preset timestamp.

To change date and time:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **System** menu and press the knob.
3. Rotate the navigation knob to **Date/Time** and press the knob to confirm the selection.
4. Rotate the navigation knob to set a date and time and press the knob to confirm. Date and time is shown as mm-dd-yy.
5. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

Warning Changing the Date and Time deletes all data and affects the storage of trends.

2.8 Changing the Volume and Brightness



When the monitor is being transported or using battery power, you can lower the brightness to save battery power.

2.8.1 Adjusting the HR and PR Volume



You can configure the HR and PR volume through the **Volume & Brightness** menu.

To configure the HR and PR volume:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Volume & Brightness** menu and press the knob to confirm.
3. Rotate the navigation knob to adjust the volume to an option and press the knob to confirm the selection.
4. Press three times the navigation knob to exit the menu.

Note — The sound depends on the **HR/PR Priority** configuration. If it is set to **ECG**, the monitor sounds HR tones. Otherwise, it sounds PR tones.

2.8.2 Adjusting the Brightness



To adjust the monitor brightness:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Volume & Brightness** menu and press the knob to confirm.
3. Press two times the navigation knob to show the **Brightness** icon and rotate the navigation knob to adjust the brightness to an option and press the knob to confirm.

2.9 Managing Patients

For better management and storage of the patient data, the patient information needs to be properly administered and documented. The monitor displays physiological data and stores it in the trends as soon as a patient is connected. This lets you monitor a patient not yet admitted. It is however important to admit patients properly so that you can identify your patient on recordings, reports and so on.

Patient Name: modifies the patient information.

New Patient: adds a new patient but this clears all old patient information.

To enter a patient:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **System** menu and press the knob.
3. Rotate the navigation knob to **Patient Name** and press the knob to enter the text field.
4. To enter data (up to 20 characters), press the navigation knob again and rotate the knob to a letter or number (numbers appear after the last letter) and press the knob to select each letter (or number).

5. Rotate the navigation knob to enter the next letter or number and repeat until you have entered the patient's name.
6. Rotate the navigation knob to the end of the text field to exit text input.
7. Rotate the navigation knob to **Exit** and press the knob to exit the **System** menu.

When the patient information is configured, the related information is shown at the bottom of the monitoring screen.

Note — The patient information is saved when there is a power loss.

2.10 Exporting ECG Waveform and Trend Data to a USB Stick

You can export some patient data to a USB stick by connecting the USB stick to the USB port on the rear panel of the monitor. The data is stored in two files whose names consist of seven numbers and one letter. The file name format is as below.

Take 8101320E for example.

The first number indicates the month. 1 to 9 = January to September; A = October; B = November and C = December. In this case, 8 means August.

The second and third numbers indicate the date. In this case, 10 means the tenth day of August.

The fourth and fifth numbers indicate the hour. In this case, 13 means 1:00 pm.

The sixth and seventh numbers indicate the minute. In this case, 20 means the 20th minute.

The last letter E or T indicates the data is on ECG Waveform or Trend.

To export data to a USB stick:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **USB Export** menu and press the knob.

A confirmation window appears to ensure that a USB stick is connected.

3. Rotate the navigation knob to **OK** and press the knob, the monitor automatically exports data to the USB stick. Select **Cancel** to cancel data export and exit the menu.

3

Alarms

This chapter describes how to configure alarms, including alarm indicators, and alarm limits. The alarm information in this chapter applies to all measurements. Alarms for specific measurements are described in their respective chapters.

Alarms sound when the monitor senses an irregular condition either in monitor operation or alarm limits.

Alarms can include:

- audible alarms
- visual alarms

3.1 Alarm Types

Alarms consist of physiological alarms and technical alarms. For more information on physiological and technical alarms, see Chapter B, Alarm Specifications.

Physiological alarms are triggered by physiological values exceeding the limits or by a patient's abnormalities.

Technical alarms are triggered by equipment malfunction or improper equipment use. When a technical alarm occurs, there might be problems with readability of patient data, however measurement is not interrupted.

3.2 Alarm Levels

Both physiological and technical alarms are classified by priority:

- High
- Medium
- Low

3.2.1 High Priority Alarms

High priority alarms indicate a life-threatening or other serious situation; for example, cardiac arrest. These alarms require immediate action from the clinician.

High priority alarms can be:

- asystole
- low differentiation between SYS and DIA
- apnea
- SpO₂ ultra-low

3.2.2 Medium Priority Alarms

Medium priority alarms indicate a physiological or technical situation that requires timely handling. These alarms can be exceeded values or system irregularity. Most medium priority alarms are triggered by values exceeding alarm limits.

3.2.3 Low Priority Alarms

Low priority alarms are typically equipment malfunctions. These alarms require less immediate handling. For example, a low priority alarm might indicate low battery power.

3.3 Visual Alarms

The monitor uses the following visual alarm indicators: flashing LED, flashing icons and alarm messages.

3.3.1 LED Flashing Alarms

There is an alarm LED on the front panel of the monitor. When the alarm is activated, the alarm LED changes according to the alarm level. The following table describes the alarm LED colors:

Table 3-1 Alarm LED Color

Alarm Priority	Color	Frequency
High	Flashing red	Every 1/2 second
Medium	Flashing yellow	Every 2 second
Low	Steady yellow	-

Note — If more than one alarm condition is active, the monitor signals the most severe.

3.3.2 Flashing Icons

When a parameter value exceeds an alarm limit, the numeric value is underlined and the corresponding alarm icon flashes. When battery power is low, the power indicator flashes every second.

3.3.3 Alarm Messages

When an alarm is triggered, the alarm message appears in the corresponding waveform area and alarm area. NIBP alarm messages display in the NIBP parameter list.

Higher priority alarm messages pre-empt other priority alarm messages for the same parameter.

For more information on alarm messages, see Chapter B, Alarm Specifications.

3.4 Audible Alarms

When audible alarms are turned on, the audible alarm patterns are repeated at certain intervals. If alarms of more than one priority are triggered, the monitor sounds the most severe. The audible alarms of your monitor are configurable. You can turn on or off the alarm sound, pause the alarm sound or adjust the alarm volume.

Note — Changing the alarm sound does not affect visual alarm indicators. You cannot turn off the audible alarms for low battery messages.

The alarm sound depends on the alarm priority. The following table describes the alarm sound in **Auto** alarm mode.

Table 3-2 Alarm Sound

Alarm Priority	Alarm Sound	Frequency
High	Double high volume beeps 5 times.	Every 10 seconds.
Medium	Single medium volume beep 3 times.	Every 25 seconds.
Low	Single low volume beep.	Every 30 seconds.

Warning Do not rely exclusively on the audible alarm system for veterinary monitor. Adjustment of alarm volume to a low level or off during patient monitoring may result in patient danger. Remember that the most reliable method of patient monitoring combines close personal surveillance with correct operation of monitoring equipment.

3.4.1 Turning on Audible Alarms



When all audible alarms are turned on, the alarm-on icon displays in the middle of the screen. To turn on all audible alarms:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **System** menu and press the knob.
3. Rotate the navigation knob to **Alarm ON/OFF** and press the knob to confirm the selection.
4. Rotate the navigation knob to **ON** and press the knob to save the setup.
5. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

To turn on an individual audible alarm:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to a parameter submenu (**ECG**, **SPO2**, **CO2**, **RESP-TEMP** or **NIBP**) and press the knob.
4. Rotate the navigation knob to the corresponding alarm setup item and press the knob to confirm the selection.
5. Rotate the navigation knob to **ON** and press the knob to save the setup.
6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

3.4.2 Turning off Audible Alarms



When all audible alarms are turned off, the alarm-off icon displays in the middle of the screen. When an individual alarm is off, the icon displays beside the parameter value on the numeric pane.

To turn off all audible alarms:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **System** menu and press the knob.
3. Rotate the navigation knob to **Alarm ON/OFF** and press the knob to confirm the selection.
4. Rotate the navigation knob to **OFF** and press the knob to save the setup.
5. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

Note — When all audible alarms are off, pressing once the **Alarm Silence** key on the front panel turns on alarm sound again. The Alarm Off configuration is not saved when the monitor is turned off. When you turn the monitor on again, audible alarms are on by default.

To turn off an individual audible alarm:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to a parameter submenu and press the knob.
4. Rotate the navigation knob to the corresponding alarm setup item and press the knob to confirm the selection.
5. Rotate the navigation knob to **OFF** and press the knob to save the setup.
6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

Warning When all alarms are turned off, there are no audible alarms for any alarm-triggered event. However, visual alarm indicators and messages continue to appear. You cannot turn these off.

3.4.3 Silencing Alarms

To silence all audible alarms, press the **Alarm Silence** key on the front panel.

Audible alarms silence for two minutes and resume when the timeout counter reaches **00:00**. When all audible alarms are silenced, the monitor displays the alarm silenced icon in the middle of the screen.

To activate all audible alarms before the timeout counter reaches **00:00**, press the **Alarm Silence** key again.

Note — Audible alarms reactivate when new alarm-triggered event occur during an alarm pause. These audible alarms sound according to the alarm priority.

3.4.4 Adjusting Alarm Volume



You can adjust the alarm volume ranged from one to seven. The minimum volume is one.

To adjust the alarm volume:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Volume & Brightness** menu and press the knob.

3. Press the navigation knob once to access the alarm volume item and rotate the navigation knob to an option and press the knob to confirm.
4. Press the navigation knob once to exit the menu.

3.5 Setting the Alarm Mode

The monitor has two types of alarm mode: automatic and standard. Changing the alarm mode only affects audible alarms.

To change the alarm mode:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **System** menu and press the knob.
3. Rotate the navigation knob to **Alarm Mode** and press the knob to confirm the selection.
4. Rotate the navigation knob to an option (**Auto** or **Standard**) and press the knob to save the setup.
5. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

3.5.1 Automatic

When the alarm mode is set to **Auto**, and alarms are on, the monitor automatically activates an audible alarm with appropriate messages in the corresponding area for any alarming event.

- If a high priority alarm or limit violation (excluding NIBP limits) occurs, the monitor sounds a high priority audible alarm until the alarm is acknowledged.
- If an alarming event (**Leads Off**, for example) or limit violation for NIBP occurs, the monitor sounds a medium priority audible alarm and the alarm LED flashes until the alarm is acknowledged.

3.5.2 Standard

When the alarm mode is set to **Standard** and alarms are on, the monitor sounds audible alarms with appropriate messages in the corresponding area for any alarming event. The alarm continues until acknowledged.

3.6 Setting Alarm Limits

You can set alarm limits based on an individual patient's vital signs measurements. These limits determine the conditions that trigger alarms based on the high and low limits set. You set alarm limits in the specific parameter alarm setup menu. See the chapters in this guide for setting specific parameter alarm limits.

Warning Always check that the alarm settings are appropriate for your patient before you start monitoring.

3.6.1 Setting Individual Alarm Limits

Warning For certain parameter, for example SpO₂, alarm limits range from 0 to 100. Setting the upper limits to 100 or setting the lower limit to zero is equivalent to no alarm limit, but without the alarm off icon.

To set an individual alarm limit:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob to confirm the selection.
3. Rotate the navigation knob to a parameter and press the knob.
4. Rotate the navigation knob to the upper or lower limit of that parameter and press the knob.
5. Rotate the navigation knob to set an appropriate value for the limit and press the knob.
6. Repeat step 4 and 5 until you complete the alarm limits setup.
7. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

Note — Both alarm limits remain in the system until they are changed, including when you turn the monitor on and off, or after a power loss.

3.6.2 Restoring Default Alarm Limits

You can restore the default alarm limits when you need to restore the factory settings.

To restore default alarm limits:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **System** menu and press the knob.
3. Rotate the navigation knob to **Default Setup** and press the knob. A confirmation window appears.
4. Rotate the navigation knob to **OK** and press the knob to confirm and resume the monitoring screen; rotate the navigation knob to **Cancel**, and press the knob to cancel the operation.

For more information on default alarm limits, see Default Alarm Limits on page B-4.

Monitoring ECG

This chapter describes ECG monitoring on the CARDELL® 9500 HD veterinary monitor.

An electrocardiogram (ECG) monitors the electrical activity of the heart. The veterinary monitor processes these electrical signals and presents an ECG waveform on the screen with its numeric on the numeric pane. The monitor can:

- Display an ECG waveform
- Compute and display heart rate (HR) values
- Detect and filter pacemaker-generated signals

4.1 Prior to Monitoring

4.1.1 Skin Preparation for Electrode Placement

Good electrode-to-skin contact is important for a good ECG signal because the skin is a poor conductor of electricity. To obtain valid signals, the patient skin must be cleaned before applying any electrodes.

To prepare the skin:

1. Select sites with intact skin and no impairment of any kind.
2. Clip or shave hair from sites as necessary.

4.1.2 Placing the Lead Wires

The labels and colors of the ECG lead wires differ according to the standards. The lead wire placement illustrations in this guide use the AAMI and IEC labels and colors. The following table describes the AAMI and IEC labels and lead colors.

Table 4-1 Lead Description

AAMI		IEC	
Lead	Color	Lead	Color
RA	White	R	Red
LA	Black	L	Yellow
LL	Red	F	Green

Table 4-1 Lead Description

AAMI		IEC	
Lead	Color	Lead	Color
RL	Green	N	Black
V	Brown	C	White

When placing lead wires, choose a flat, non-muscular site where the signal will not be interfered with by movement or bones. Correct lead placement is important for accurate diagnosis, especially in the precordial leads, which are close to the heart, QRS morphology can be greatly altered if an electrode is moved away from its correct location.

The following figure shows the 3-lead ECG electrode locations:




Figure 4-1 Three-Lead Locations

The following figure shows the 5-lead ECG electrode locations:




Figure 4-2 Five-Lead Locations

4.1.3 Changing a Lead Set

A pair of lead wires forms a lead (electrical pathway). Each lead provides a different view of the same cardiac activity. You can select a proper lead mode and an ECG lead. The selected lead mode and ECG lead displays on top of the ECG waveform.

To change the ECG lead wire set:

1. Remove the lead wires and then replace them as required.
2. Reselect the **ECG Lead Mode** in the **System** menu and **ECG Lead** in the ECG menu.

Select the appropriate **ECG Lead Mode** based on the number of lead wires you are using. If you remove lead wires and do not select the appropriate **ECG Lead Mode**, the monitor displays a **Leads Off** alarm message. Select the appropriate **ECG Lead Mode** and the alarm message disappears.

To select an ECG lead mode (lead wires):

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **System** menu and press the knob.
3. Rotate the navigation knob to **ECG Lead Mode** and press the knob.
4. Rotate the navigation knob to a lead mode option (**Five Lead** or **Three Lead**) and press the knob.
5. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

To select an ECG lead (electrical pathway):

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to **ECG** and press the knob.
4. Rotate the navigation knob to **ECG Lead** and press the knob.
5. Rotate the navigation knob to a lead set option (see Table 4-2) and press the knob.
6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

Table 4-2 ECG Lead Mode

When Lead Mode is	The options are
Three Lead	I, II, III
Five Lead	I, II, III, V, aVR, aVL, aVF

I, II and III are bipolar extremity leads. aVR, aVL and aVF are voltage-added extremity leads. V is the chest lead.

4.1.4 Lead Recommendations for Surgical Patients

The lead location depends on the type of surgical operation. For example, for chest surgery, the lead wires should be attached to the side or the back of the chest.

When HF electrosurgical equipment is used, patient leads and wires should be away from the surgical operation site and other devices. This reduces the burning risk due to a poor connection of the neutral electrode.

Warning Electrosurgical equipment must be properly grounded to avoid a current inflow hazard. This can cause interference on the ECG signal.

4.1.5 Connecting ECG Cables

To connect ECG cables:

1. Attach the lead wires to the ECG cable. Plug the ECG cable into the ECG connector on the side panel of the monitor (see CARDELL® 9500 HD Side Panels on page 2-5).
2. Attach the alligator clips or snaps to the lead wires before placing them. If you are not using pre-gelled electrodes, apply electrode gel to the alligator clips before placement.
3. Place the alligator clips/lead wires on the patient according to the selected lead placement.

Warning To protect the monitor from damage during defibrillation, to catch accurate ECG information, and to protect against noise and other interference, use only ECG electrodes and cables approved by manufacturer.

4.2 The ECG Display

When ECG is being monitored, the monitor displays the ECG waveform and its numeric HR values.

4.2.1 ECG Waveform and Numeric Pane

The following figure shows a typical ECG waveform and numeric pane:




Figure 4-3 Typical ECG Waveform

- 1-- ECG lead
- 2-- ECG gain
- 3-- ECG waveform speed
- 4-- Pacemaker status
- 5-- ECG mode
- 6-- Heart rate/Pulse rate label(unit)
- 7-- 1 mV calibration bar
- 8-- ECG waveform

4.3 Configuring the ECG Display

By changing the ECG display, you can get a screen effect as needed.

4.3.1 Changing the Size of the ECG Waveform

If the displayed ECG wave is too small or it is too tall, use the **ECG Gain** menu to adjust the size of the waveform. Changing the ECG gain only affects the visual appearance of the ECG wave on the screen. It does not affect the ECG signal analyzed by the monitor.

To change the ECG gain:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to **ECG** and press the knob.
4. Rotate the navigation knob to **ECG Gain** and press the knob.
5. Rotate the navigation knob to an option (**X1/4**, **X1/2**, **X1**, **X2**, **X4** or **Auto**) and press the knob.
6. Rotate the navigation knob to **Exit** and press the knob.

Selecting different gains displays the ECG waveform in different scales. In addition, a 1 mV calibration bar displays on the left of the ECG channel. The following table describes the ECG scale:

Table 4-3 ECG Scale

When the ECG Gain is X	The Scale Height
1/4	is 2.5 mm.
1/2	is 5 mm.
1	is 10 mm.
2	is 20 mm.
4	exceeds the ECG display area.
Auto	varies in step with the ECG gain.

4.3.2 Changing the Speed of the ECG Waveform

The waveform speed setting in the ECG menu determines the speed at which the waveform is drawn across the screen. For more information on changing this setting, see Changing the Waveform Speed on page 2-8.

The ECG waveform speed options are:

- **50 mm/s**
- **25 mm/s**
- **12.5 mm/s**

The selected waveform speed displays on top of the wave.

4.4 Performing 1 mV Calibration



During ECG monitoring, select the **1mV Calibration** button to display a 1 mV square wave on the ECG waveform. Compare the waveform with the square wave to estimate the signal strength.

To perform 1 mV calibration:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to **ECG** and press the knob.
4. Rotate the navigation knob to the **1mV Calibration** button and press the knob.
5. Wait until the calibration completes and the button recovers.
6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

4.5 Changing the Heart Rate Volume

To change the HR volume, see Adjusting the HR and PR Volume on page 2-10 for more information.

4.6 Changing the ECG Mode

The monitor has three ECG modes for different monitoring environments: **Monitor**, **Operation** and **Diagnosis**. The current ECG mode displays on top of the ECG channel.

To change the ECG mode:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to **ECG** and press the knob.
4. Rotate the navigation knob to **ECG Mode** and press the knob.
5. Rotate the navigation knob to an option among **Monitor**, **Operation** and **Diagnosis** and press the knob.

In **Monitor** mode, false alarms are filtered.

In **Operation** mode, the monitor reduces false readings and interference. Use this only in locations (operating room, for example) where there is significant external interference that might cause ECG waveform distortion.

In **Diagnosis** mode, the original (unfiltered) ECG waveform displays.

6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

4.7 Changing HR Alarm Settings

To turn on the HR alarm sounding, see Turning on Audible Alarms on page 3-3.

To turn off the HR alarm sounding, see Turning off Audible Alarms on page 3-4.

For more information on HR alarm settings, see Setting Individual Alarm Limits on page 3-5.

4.8 Configuring HR/PR Priority

Configure the HR/PR priority according to your site's needs. The monitor provides tones, measurement values and color of the heart rate indicator in step with the HR/PR priority configuration.

To configure the HR/PR priority:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to **ECG** (or **SPO2**) and press the knob.
4. Rotate the navigation knob to **HR/PR Priority** and press the knob.
5. Rotate the navigation knob to an option (**ECG** or **SPO2**) and press the knob.
6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

The monitor provides the following tones when the HR/PR priority is set:

- When **ECG** is the first priority and a heart rate is detected, the monitor sounds HR tones; if a heart rate is not detected, the monitor sounds the PR tones of the parameter which is second priority.
- When **SPO2** is the first priority and the SpO₂ pulse rate is detected, the monitor sounds SpO₂ pulse rate tones; if SpO₂ pulse rate is not detected, the monitor sounds HR tones of the parameter which is second priority.

Table 4-4 HR/PR Priority

Priority	HR Source	ECG Numeric Pane Displays	Monitor Provides	HR Indicator and Numeric Color is
ECG	HR detected	HR	HR tone	same as ECG waveform
	PR detected and HR not	PR	PR tone	same as SpO ₂ waveform
	no HR or PR detected	--	None	None
SpO ₂	PR detected	PR	PR tone	same as SpO ₂ waveform
	HR detected and PR not	HR	HR tone	same as ECG waveform
	no PR or HR detected	--	None	None

4.9 Monitoring Patients with Pacemakers



Pacemaker detect function must be switched on for paced patients. The paced status indicator displays on top of the ECG channel

When pacemaker-generated signal is detected, the pace pulse markers (green short lines) display on the ECG wave to suggest the position where the pace pulse is detected.

To select a certain ECG lead for the pacemaker detection:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to **ECG** and press the knob.
4. Rotate the navigation knob to **Pacemaker Det.** and press the knob.
5. Rotate the navigation knob to **OFF** and press the knob to turn the pacemaker detection off.
6. Rotate the navigation knob to an option (**ON I**, **ON II**, **ON III** or **ON V**) to turn the pacemaker detection on.

ON I: use lead I signal for pacemaker detection.

ON II: use lead II signal for pacemaker detection.

ON III: use lead III signal for pacemaker detection.

ON V: use lead V signal for pacemaker detection.

7. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

Warning Some pace pulses can be difficult to detect. When this happens, the pulses are counted as a QRS complex, and could result in an incorrect HR and failure to detect cardiac arrest or some arrhythmias. Keep pacemaker patients under close observation. When changing configurations and at admission/discharge, always check that paced status is correct for the patient.

With a properly placed three-lead connection, use signals from the channel in conformity with the ECG lead as the pacemaker detect signal. With a properly placed 5-lead connection, use leads I, II, III and V with comparatively high amplitude signals as a pacemaker detect signal.

4.10 Safety Information

Warning Use only ECG cables or leads specified by the manufacturer. Other ECG cables and leads can cause improper performance and/or provide inadequate protection during defibrillation.

Using electrodes of dissimilar metal materials can cause over polarization or accelerated polarization.

Ensure that the ECG connectors or electrodes never come into contact with other conductive parts, or with grounding cable. In particular, ensure that all of the ECG electrodes are attached to the patient.

When a defibrillator is used, ensure that the electrodes and patient cables are not in contact with metal or other conductive surface or device grounds.

Do not use non-defibrillation cables to perform defibrillation on a patient.

Never entangle the electrosurgical (ES) cables with ECG cables. When using ES equipment, place the ECG electrodes halfway between the ES grounding plate and the ES knife to avoid burning.

When using electrosurgical equipment, never place ECG electrodes near the grounding plate of the ES devices, as this can cause interference on the ECG signal.

Properly ground electrosurgical equipment to reduce hazards from current inflow which might cause interference to the ECG signals. When using ES equipment, place the patient lead wires and cables far way from the operating room table and other devices.

Note — ECG cables can be damaged when connected to a patient during defibrillation. Check cables that have been connected to a patient during defibrillation for functionality before using them again.

Monitoring with the pace detection feature enabled does not normally affect the monitoring of non-pacemaker patients. However, in some instances, if the patient does not have a pacemaker, it can be desirable to turn the detection function off so that artifacts in the waveform are not mistaken for a pacemaker signal.

When the ECG circuitry is overloaded or any part of the amplifier saturated, the top or the bottom of the ECG waveform appears dashed or becomes a straight, broken-line.

Thorough skin cleaning plays an important role in minimizing electrostatic interference. Before using HF devices, properly clean the skin and attach electrodes correctly to patient site. This reduces the hazard of burns in case of a defective neutral electrode at the HF device.

Line isolation monitor transients may resemble actual cardiac waveforms and thus inhibit heart rate alarms. Such transients may be minimized by proper electrode and cable placement, as specified in this manual and electrode directions for use.

When the monitor is connected to a paced patient or other electric generators, or when the patient has an arrhythmia problem, incorrect HR measurements could occur. In this case, ECG signals might still appear after the patient is dead. Keep the paced patient under close observation in order to make correct diagnoses.

Monitoring Respiration

This chapter describes respiration monitoring on the CARDELL® 9500 HD veterinary monitor. Respiration can be measured using one of the following methods:

- **ECG.** The monitor measures the thoracic impedance between two ECG electrodes on the patient's chest. Changes in the impedance due to thoracic movement produce the Resp waveform on the monitor screen. The monitor counts the waveform cycles to calculate the respiration rate (RR).
- **CO₂.** You can measure a patient's airway respiration rate (awRR). The awRR value is calculated by directly measuring air movement in and out of the patient's airway. This is the preferred method of measurement as it is measuring breaths not chest movement. (see Monitoring CO₂ section 9)

5.1 Lead Placement for Monitoring Respiration with ECG

Correct patient skin preparation techniques for electrode placement are important for respiration measurement. See Before Monitoring on page 4-1.

Cardiac Overlay: Cardiac activity that affects the respiration waveform is called cardiac overlay. Correct electrode placement can help to reduce cardiac overlay: avoid the liver area and the ventricles of the heart in the line between the respiratory electrodes. This is particularly important for small patients.

Lateral Chest Expansion: Some patients, especially small patients, expand their chests laterally. In this case it is best to place the two respiratory electrodes in the right mid axillary and left lateral chest areas at the patient's maximum point of breathing movement to optimize the respiratory wave.

Abdominal Breathing: Some patients with restricted chest movement breathe mainly abdominally. In these cases, you may need to place the left leg electrode on the left abdomen at the point of maximum abdominal expansion to optimize the respiratory wave.

5.1.1 Selecting a Respiration Lead

You can select a respiration lead as necessary. Resp is detected from impedance changes between the selected ECG lead wires.

To select a respiration lead:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to **RESP-TEMP** and press the knob.
4. Rotate the navigation knob to **RESP Lead** and press the knob.
5. Rotate the navigation knob to a lead set option and press the knob.
 - When ECG **Lead Mode** is set to **Three Lead**, the **RESP Lead** options are: **RA-LA**, **RA-LL**.
 - When ECG **Lead Mode** is set to **Five Lead**, the **RESP Lead** options are: **RA-LA**, **RA-LL**, **LA-RL**, **LL-RL**.
6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

5.2 The Respiration Display

The respiration measurement is displayed on the monitor as a continuous wave and a numeric respiration rate. The following figure shows an example of the respiration waveform and numeric pane:




Figure 5-1 Resp Waveform

- 1 -- Respiration lead
- 2 -- Respiration gain
- 3 -- Respiration waveform speed
- 4 -- Respiration waveform
- 5 -- Respiration indicator (unit)
- 6 -- Respiration value

5.3 Changing the Respiration Display

You can configure the monitor respiration display by changing the respiration settings.

5.3.1 Changing the Respiration Waveform Size

You can configure the respiration waveform size by selecting a respiration gain to increase/decrease the waveform height.

To change the respiration gain:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to **RESP-TEMP** and press the knob.
4. Rotate the navigation knob to **RESP Gain** and press the knob.
5. Rotate the navigation knob to an option (**X 1/2**, **X 1**, or **X 2**) and press the knob.
6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

5.3.2 Changing the Respiration Waveform Speed

The speed of waveform scanning on the screen is determined by the configured waveform speed. For more information on changing the speed of waveform, see Changing the Waveform Speed on page 2-8.

The options for waveform speed are: **25 mm/s**, **12.5 mm/s** and **6.25 mm/s**.

5.3.3 Changing the Respiration Waveform Display Mode

To change the respiration waveform mode, see Changing the Waveform Display on page 2-8 for more information.

5.4 Changing the Respiration Alarm Settings

To turn on the respiration alarm, see Turning on Audible Alarms on page 3-3. To turn off the respiration alarm, see Turning off Audible Alarms on page 3-4.

For more information on respiration alarm limits, see Setting Individual Alarm Limits on page 3-5.

5.5 Changing the Apnea Alarm Settings

The apnea alarm is a high priority alarm used to detect apneas. The respiration measurement does not recognize obstructive and mixed apneas. It only indicates an alarm when a user-defined time has elapsed since the last detected breath.

To turn apnea alarm on or off:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to **RESP-TEMP** and press the knob.
4. Rotate the navigation knob to **APN Alarm Time** and press the knob.
5. Rotate the navigation knob to an option and press the knob to save the configuration.

Apnea Alarm Time is the preset time for apnea alarm detection. It ranges from **5** seconds to **120** seconds at 5 seconds intervals. **OFF** indicates that the apnea alarm is turned off.

6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

5.6 Safety Information

Warning Radiated electromagnetic signals from sources external to the patient and might cause inaccurate respiration readings. Keep patients under close surveillance when monitoring respiration.

Excessive patient movement could result in incorrect respiration measurements. Pay close attention to patients who move excessively.

Note — When the ECG lead wires are applied to patient's forelimbs, respiration impedance monitoring is not available.

6

Monitoring SpO₂

This chapter describes how to monitor SpO₂. The monitor measures the functional arterial oxygen saturation. This is the percentage of oxygenated hemoglobin in relation to the sum of oxyhemoglobin and deoxyhemoglobin.

6.1 SpO₂ Monitoring Theory

SpO₂ monitoring is determined by spectrophotometric principles allowing a clinician to measure oxygen saturation without drawing blood. The monitor uses an SpO₂ sensor that clips onto a patient's tongue, for example. The probe (sensor) determines SpO₂ and pulse rate by passing two kinds of light: red light (660 nm) and infrared (940 nm) through the body tissue to a photo detector. The photoelectricity detector converts the red light and infrared into digital signals and the pulse rate of the photoelectricity signal is in step with the PR. Therefore, PR can be calculated together with a signal cycle.

6.2 Placing an SpO₂

Follow the SpO₂ sensor's instructions for use, adhering to all warnings and cautions.

To correctly apply the sensor:

- Ensure that the sensor is applied to smooth skin sites with sufficient arterial blood flow and that no excessive patient movement and/or vibration occurs.
- Ensure that the sensor is not on extremities with an arterial catheter or intra vascular venous infusion line.
- Ensure that the sensor is the appropriate size. The sensor should not be so loose as to fall off, nor should it be too tight.
- Ensure that the light emitter and the photo-detector are directly opposite each other. All light from the emitter must pass through the patient's tissue.

Warning At elevated ambient temperatures be careful with measurement sites that are not well perfused, because this can cause severe burns after prolonged application.

6.3 Connecting SpO₂ Cables

Connect one end of the sensor cable to the SpO₂ probe connector on the side panel of the monitor (see The Side Panels on page 2-5) and apply the other end (or via an extension cable) to patient site.

Warning Position the sensor cable and connector away from power cables to avoid electrical interference.

6.4 Measuring SpO₂

Note —Clean the sensor by gently wiping it with an agent such as 70% Isopropyl Alcohol, before and after use. Never immerse the sensor in liquid.

To start monitoring SpO₂:

1. Connect the SpO₂ cables and ensure that the connection is good.
2. Select appropriate sized clip and apply sensor to patient
3. Confirm proper SpO₂ alarm limits. See Setting Individual Alarm Limits on page 3-5.
4. The monitor displays the SpO₂ waveform and the SpO₂ value when patient SpO₂ is detected.

During measurement, the following factors can influence SpO₂ readings:

- Pigmented skin
- Injected dyes like methylene blue or intra vascular dyshemoglobins (methemoglobin and carboxyhemoglobin).
- Strong ambient lighting and patient motion.
- Misplacing the sensor.
- Poor perfusion.

Warning Inspect the application site every two to three hours to ensure skin quality and correct optical alignment. If the skin quality changes, move the sensor to another site. Change the sensor position at least every hour.

6.5 Assessing a Suspicious SpO₂ Reading

Inaccurate readings may be caused by patient motion or weak signals, especially for small patients. Observe the following guidelines when monitoring small patients:

- There should be a six-second uninterrupted non-artifact SpO₂ waveform.
- A stable SpO₂ value should remain for at least six seconds.

When pulse rate is very low, or serious arrhythmia occurs, the SpO₂/Pleth pulse rate may differ from the heart rate calculated from ECG. However, this does not indicate an inaccurate SpO₂ value.

If you suspect the measured SpO₂ value is inaccurate, use the waveform and perfusion indicator to assess the signal quality.

6.6 The SpO₂ Display

The following figure shows an example of an SpO₂ waveform and numeric pane.




Figure 6-1 Pleth Wave and SpO₂ Numeric Pane

- 1 -- SpO₂ waveform speed
- 2 -- SpO₂ waveform
- 3 -- SpO₂ saturation unit (100%)
- 4 -- SpO₂ perfusion indicator bar
- 5 -- SpO₂ value

6.7 Changing the SpO₂ Display

You can configure the monitor SpO₂ display by changing the SpO₂ settings.

6.7.1 Changing the SpO₂ Waveform Speed

The speed of waveform scanning on the screen is determined by the waveform speed configured. For more information on changing the speed of waveform, see [Changing the Waveform Speed](#) on page 2-8.

The options for waveform speed are: **50 mm/s**, **25 mm/s** and **12.5 mm/s**.

6.7.2 Changing the SpO₂ Waveform Display Mode

To change the SpO₂ waveform mode, see [Changing the Waveform Display](#) on page 2-8 for more information.

6.8 Configuring HR/PR Priority

To configure HR/PR priority, see [Configuring HR/PR Priority](#) on page 4-6 for more information.

6.9 Changing the SpO₂ Alarm Settings

To turn on the SpO₂ alarm, see [Turning on Audible Alarms](#) on page 3-3. To turn off the SpO₂ alarm, see [Turning off Audible Alarms](#) on page 3-4.

For more information on SpO₂ alarm limits, see [Setting Individual Alarm Limits](#) on page 3-5.

6.10 Changing the PR Alarm Settings

To turn on the PR alarm, see Turning on Audible Alarms on page 3-3. To turn off the PR, see Turning off Audible Alarms on page 3-4.

To configure the PR alarm limits:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to **SPO2** and press the knob.
4. Rotate the navigation knob to **PR Upper Limit** and press the knob; rotate the navigation knob to set an alarm limit.
5. Rotate the navigation knob to **PR Lower Limit** and press the knob; rotate the navigation knob to set an alarm limit.
6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

6.11 Safety Information

Warning Continuous monitoring might make skin red, blistered, or cause compression necrosis. This is especially true for small patients or patients with perfusion disorder or unhealthy skin. Periodically inspect the sensor application site and change the application site at least every hour.

Avoid placing the sensor on extremities with an arterial catheter, or intravascular venous infusion line.

If a sensor is too loose, it might compromise the optical alignment or fall off. If it is too tight, for example because the application site is too large or becomes too large due to edema, excessive pressure may be applied. This can result in venous congestion distal from the application site, leading to interstitial edema, hypoxemia and tissue malfunction.

Caution The sensor consists of sensitive electronic components which may be damaged if misused. Keep sharp-edged objects away from the cable.

Improper cleaning of the sensor could compromise the SpO₂ sensor life span.

Note — Use only the specified sensor for your monitor.

Do not use a sensor if you see any signs of damage.

Monitoring NIBP

This chapter describes NIBP monitoring on the CARDELL® 9500 HD veterinary monitor.

The 9500HD uses oscillometric principles to calculate the systolic, diastolic, and mean arterial pressure (MAP) values. The MAP is measured at the lowest cuff pressure that provides the maximum cuff oscillations. Therefore, MAP is the largest signal received and is the most accurate reading using oscillometric methods. Systolic pressure is calculated as the cuff pressure at which an increase in cuff oscillations is perceived. The diastolic pressure is the cuff pressure when oscillations are no longer decreasing as pressure is released from the cuff. It uses the changes of the oscillometric amplitude under different cuff pressures to identify mean pressure and calculate the systolic and diastolic pressure. Special veterinary specific algorithms have been designed to ensure reliable and accurate measurements from kittens to horses.

7.1 Selecting and Placing an NIBP Cuff

Caution Only accessories recommended by Midmark should be used.

Note — Place the patient on a padded surface to provide comfort, and warmth. Shivering will inhibit the monitor from making a determination

7.1.1 Cuff Placement for a Cat


A cat may be left in its owner's lap to keep it calm. Measurements are best done in an area of the hospital away from noise and bright lights. The animal may be held so that the front limbs are free for cuff placement. In conscious patients, the tail may be the most appropriate location for placement of the cuff. Cats may be most comfortable in sternal recumbency making the tail a more preferable site. For the median artery on the foreleg, place the cuff around the forelimb, between the elbow and carpus. It is not necessary to center the cuff over the artery which is on the medial side of the leg because of the fully encircling bladder design. Hair need not be clipped except when heavily matted. In cats less than five (5) pounds when measurements are difficult to obtain, place the cuff around the leg above the elbow to obtain measurements from the brachial artery. Measurements from the coccygeal artery may be used by placing the cuff around the base of the tail but not in anesthetized patients.



7.1.2 Cuff Placement for a Dog

For measurements in dogs, it is preferable to use the right lateral, sternal or dorsal recumbent positions. That is not a problem in anesthetized patients, but it may be difficult to get large dogs to cooperate for proper positioning. If the dog is in a sitting position, place the front paw on the operator's knee and take measurements from the metacarpus. Sites for cuff placement are the metacarpus, metatarsus and anterior tibial. In anesthetized patients, most surgeries are done on the posterior part of the body so the metacarpal area of the forelimb is most convenient. In situations where this is not possible, the cuff should be wrapped around the metatarsus just proximal to the tarsal pad or around the hind leg just distal to the hock. The tail site should not be used for cuff placement during anesthesia. It is not necessary to center the cuff over the artery because of the fully encircling bladder design. If the hair over the artery site is too thick or matted for good contact, it should be clipped.

Figure 7-1 NIBP Cuff Placement



7.1.3 Cuff Placement for Large Animals

A large animal such as a horse should be in a stock, standing still, or lying down. For horses and cows, the cuff can be wrapped around the base of the tail using the coccygeal artery on the ventral surface.

Note — To achieve the most accurate readings, it is important to keep the cuff on a horizontal plane with the heart.

Use the cuff selector to measure limb where cuff will be placed and determine appropriate size cuff for patient.

7.1.4 Cuff Size Selections

The widest cuff that can be placed on the patient, without extending beyond the joint, should be selected. Appropriate sized cuffs may be selected based on published guidelines that cuff width should be 40 – 60% of limb circumference. The cuff should be wrapped for a snug fit. Overlapping the cuff will not affect measurement results. Make sure the hook and loop sections of the cuff are fully engaged when it is wrapped around the limb. If not fully engaged, the cuff will detach during bladder inflation. If that happens, select the next size bigger cuff. Adhesive tape or other material should not be used to secure the cuff. Use the following table as a guide to select the correct size.

Small Animal Cuff Selection

Reorder Number	Bladder Size (Width)	Limb Circumference Range
SV1	2.0 cm	3-6 cm
SV2	2.5 cm	4-8 cm
SV3	3.5 cm	6-11 cm
SV4	4.0 cm	7-13 cm
SV5	5.0 cm	8-15 cm
SV600 (kit)	Includes all of the above	

Large Animal Cuff Selection

Cuff Model Number	Bladder Size (Width)	Limb Circumference Range
SV8	8 cm	13 – 20 cm
SV10	10.2 cm	18 – 26 cm

References:

Pedersen KM, Butler MA, Ersboll AK, Pedersen HD (2002). Evaluation of an oscillometric blood pressure monitor for use in anesthetized cats. JAVMA 221: 646-650.

Sawyer DC, Guikema AH, Siegel EM (2004). Evaluation of a new oscillometric blood pressure monitor in isoflurane anesthetized dogs. Vet Anaesth Analg 31: 27 – 39.

Warning When monitoring over an extended period of time, or at frequent intervals, periodically observe the patient's limb to make sure that the circulation is not impaired for a prolonged period of time.

If the cuff is too small or too tight, the result will be higher than the normal value, and vice versa. Select a proper cuff for different patients and correctly place the cuff to collect reliable results.

For small patients, use only the NIBP cuff recommended by the manufacturer to avoid overloaded inflation pressure.

7.2 Connecting the Cuff

Connect the selected cuff to the hose and the hose to the NIBP cuff connector on the side panel of the monitor, as shown in the CARDELL® 9500 HD Side Panels on page 2-5. The hose clicks into place when seated correctly.

Note — Do not compress the hose or restrict the pressure.

7.3 Measurement Limitations

Measurements are impossible with pulse rate extremes of less than 20 bpm or greater than 265 bpm,

The measurement may be inaccurate or impossible:

- If a regular arterial pressure pulse is hard to detect.
- With cardiac arrhythmias.
- With excessive and continuous patient movement such as shivering or convulsions.
- With rapid blood pressure changes.
- With severe shock or hypothermia that reduces blood flow to the peripheries.
- With obesity, where a thick layer of fat surrounding a limb dampens the oscillations coming from the artery.
- On an edematous extremity.

7.4 Measuring the NIBP



To measure the NIBP:

1. Select an NIBP cuff of proper size and wrap the cuff around the patient's limb at heart level. For detailed information, see *Selecting and Placing an NIBP Cuff* on page 7-1.
2. Connect the selected cuff. See *Connecting the Cuff* on page 7-2.
3. Select a proper NIBP cuff size and a measurement mode and set the appropriate alarm limits.
4. Press the **NIBP** key on the front panel of the monitor to start NIBP measurement. When NIBP is being measured, press the **NIBP** key again to stop NIBP measurement.

7.5 The NIBP Display

The following figure shows an example of the NIBP numeric pane.




Figure 7-2 NIBP Numeric Pane

- 1 -- NIBP label (unit)
- 2 -- Systolic pressure value
- 3 -- Cuff size
- 4 -- Diastolic pressure value
- 5 -- Work mode (Auto mode in this example)
- 6 -- Interval time (7 minutes in this example)
- 7 -- Mean arterial pressure value

7.6 Changing the NIBP Cuff Size

Set an appropriate cuff size before taking an NIBP measurement. The current cuff size displays in the NIBP numeric pane. Use Small for patients using the white vinyl cuffs (SV1-SV5), and Large when using the nylon cuffs (SV8 & SV10).

To change the NIBP cuff size:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to **NIBP** and press the knob.
4. Rotate the navigation knob to **Cuff Size** and press the knob.
5. Rotate the navigation knob to a proper option (**Large** or **Small**) and press the knob.
6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

7.7 Selecting NIBP Measurement Methods

There are three methods of measuring NIBP:

- **Manual**
- **Auto**
- **Stat**

Manual

Selecting **Manual** lets you manually take NIBP readings. When you press the NIBP key on the front panel, the monitor starts/stops an NIBP measurement. When an NIBP measurement is started, the monitor inflates the cuff to a certain pressure and then deflates it to get an NIBP value.

Auto

Selecting **Auto** specifies that the monitor take NIBP measurements at regular intervals. The monitor inflates the cuff at the start of each measurement cycle and collects NIBP values. In Auto mode, if the patient blood pressure is high, the monitor adjusts the maximum inflation pressure by re-inflating the cuff. NIBP measurements automatically stop if NIBP values cannot be correctly collected in Auto mode.

Stat

Selecting **Stat** specifies that the monitor continuously measure patient NIBP for 5 minutes. Use Stat for continuous NIBP monitoring of supervised patients in special cases. The monitor displays 0: 05: 00 in the NIBP parameter area and starts to count down after starting the first NIBP measurement.

When Stat mode times out or Auto mode pauses due to a measurement error or technical malfunction, the monitor switches to Manual mode.

To change the NIBP measurement mode:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to **NIBP** and press the knob.
4. Rotate the navigation knob to **Work Mode** and press the knob.
5. Rotate the navigation knob to an option (**Manual**, **Auto** or **Stat**) and press the knob.
6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

Note — When monitoring NIBP, the initial inflation pressure is adjusted (except in **Manual** mode) according to the previous measurement value. The initial inflation pressure specifies the maximum amount that the cuff will inflate for the first NIBP measurement. The factory default value for initial inflation pressure is 150 mmHg.

In **Auto** mode, the interval time displays in the NIBP numeric pane. The monitor starts to count down after an NIBP measurement completes. When the measurement is paused, it reverts to the preset time display until the next measurement.

During induction and following patient alarms for hypotension and bradycardia it is recommended to cycle NIBP in **Stat** mode until the patient is stabilized

7.8 Selecting the NIBP Interval Time

An interval is measured from the end of one NIBP measurement to the start of the next measurement and displays in the NIBP numeric pane. For example, **0: 01** is displayed as **0: 01: 00** in the numeric pane. You can configure the NIBP interval time only in **Auto** mode. During anesthetic procedures the recommended interval time is 3-5 minutes.

To change the NIBP interval time:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to **NIBP** and press the knob.
4. Rotate the navigation knob to **Interval Time** and press the knob.
5. Rotate the navigation knob to an option and press the knob to confirm the selection. The options are: **0: 01** to **0: 10** at one-minute intervals, **0: 15**, **0: 30**, **1: 00**, **1: 30**, **2: 00**, **3: 00**, **4: 00** and **8: 00**.
6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

Warning In Auto or Stat mode, prolonged measurements could increase friction between the cuff and patient skin which could cause purpura, ischemia, and neuropathy. Inspect the application site regularly to ensure skin quality and inspect the extremity of the cuffed limb for normal color, warmth and sensitivity. If any abnormality occurs, place the cuff on another site or stop NIBP monitoring.

7.9 Changing the NIBP Unit

You can configure the NIBP unit to **mmHg** or **kPa**. The unit displays in the NIBP numeric pane.

To change the NIBP unit:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to **NIBP** and press the knob.
4. Rotate the navigation knob to **Unit** and press the knob.
5. Rotate the navigation knob to an option (**mmHg** or **kPa**) and press the knob to confirm the selection.
6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

Note — When you change the NIBP unit, the change displays wherever pressure parameters are displayed on the monitor.

7.10 Viewing the NIBP Measurement Data

- After a complete NIBP measurement, the NIBP measurement value displays in the NIBP numeric pane.
- The NIBP values measured are listed in the NIBP history. Measurement data is shown in the list in time sequence, with the latest one being on top.

When no menu is shown on the monitoring screen, you can view the NIBP list by pressing the **NIBP History** key on the front panel of the monitor.

The NIBP list appears in the SpO₂ waveform channel and the RESP waveform channel.

The following figure shows a typical NIBP list:

DATE	TIME	SYS	DIA	MAP	PR
06-20	01:11	254	199	215	79
06-20	01:10	253	200	216	80
06-20	01:09	254	201	217	79
06-20	01:08	254	201	217	80
06-20	01:07	255	200	217	80

Figure 7-3 Example of NIBP List

1 -- NIBP list title

2 -- NIBP measurement data

7.11 Changing the NIBP Alarm Settings

To turn on the NIBP alarm, see Turning on Audible Alarms on page 3-3. To turn off the NIBP alarm, see Turning off Audible Alarms on page 3-4.

For more information on NIBP alarm limit settings, see Setting Individual Alarm Limits on page 3-5.

You can configure an NIBP alarm for pressure difference between the systolic value and the diastolic value. When the difference is less than the set value, the monitor sounds an alarm.

To set the alarm:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to **NIBP** and press the knob.
4. Rotate the navigation knob to **SYS-DIA Alarm** and press the knob.
5. Rotate the navigation knob to an option and press the knob to confirm the selection.
6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

7.12 Calibrating NIBP

Calibration of NIBP is not typically necessary during the warranty period. NIBP calibration is not user-configured. Calibration is included as part of the Preventative Maintenance service recommended every two years post warranty.

7.13 Safety Information

Warning Continual NIBP measurements can cause injury to the patient being monitored. Weigh the advantages of frequent measurement and/or use of Stat mode against the risk of injury.

In some cases, rapid, prolonged cycling of an NIBP cuff has been associated with any or all of the following: ischemia, purpura, or neuropathy. Apply the cuff according to the directions and check the cuff site and cuffed extremity regularly when blood pressure is measured at frequent intervals or over extended periods of time.

Check the patient's limb to assure that circulation is not constricted, i.e. no discoloration or ischemia of the extremities. Check the limb at regular intervals based on specific situation.

Never place the cuff on an area where circulation is compromised or has the potential to be compromised.

Never use the NIBP cuff on a limb with an intravenous infusion or arterial catheter in place. This could cause tissue damage around the catheter when the infusion is slowed or blocked during cuff inflation.

Never apply the blood pressure cuff to the same extremity as the one to which an SpO₂ sensor is attached because the cuff inflation disrupts SpO₂ monitoring and leads to nuisance alarms.

Select an appropriate cuff and ensure that the correct cuff size is selected before monitoring. Wrong cuff sizes can overload inflation pressure on small patients.

Use clinical judgment before using NIBP monitoring on patients with serious blood clot disease due to the risk of hematoma in the limb with the cuff.

Monitoring Temperature

This chapter describes how to monitor body temperature. The monitor measures temperature through the sensor temperature and electrical impedance. The monitor can measure one temperature and displays the temperature value in the numeric pane.

8.1 Selecting and Connecting Temperature Probe

Included in your accessory kit is an esophageal/rectal probe.

1. Connect the temperature cable to the TEMP connector (see CARDELL® 9500 HD Side Panels on page 2-5) on the monitor.
2. Place the other end of the cable on the patient site.

The numeric pane of the monitor displays the temperature value collected.

8.2 Influencing Factors

The following factors can affect temperature measurements:

- Improper position of the temperature probe.
- Poor heat balance, heat dissipation.
- Patient temperature fluctuation.
- Site specific influences, including air temp during esophageal use and presence of stool during rectal use.

8.3 Monitoring Temperature

Warning Before taking temperature measurements, ensure that the temperature probe is not in contact with a heat source.

To measure the temperature:

1. Connect the probe to the monitor. See Selecting and Connecting Temperature Probe on page 8-1.
2. Apply the probe to the correct location on the patient.
3. Set the temperature unit and appropriate alarm limits
4. Start temperature monitoring.

The temperature value displays in the numeric pane.

Note — Clean the probe before and after each use.

8.4 Changing the Temperature Unit

You can configure the temperature unit to °C or °F. The unit displays in the temperature numeric pane.

To change the temperature unit:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to **RESP-TEMP** and press the knob.
4. Rotate the navigation knob to **Unit** and press the knob.
5. Rotate the navigation knob to an option (°C or °F) and press the knob to confirm the selection.
6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

8.5 Changing the Temperature Alarm Settings

To turn on the temperature alarm, see Turning on Audible Alarms on page 3-3. To turn off the temperature alarm, see Turning off Audible Alarms on page 3-4.

For more information on temperature alarm limits, see Setting Individual Alarm Limits on page 3-5.

8.6 Safety Information

Warning Use only the specified probes for your monitor.

When the probe is disconnected or loose, the monitor stops measuring the temperature with no audible alarm. To ensure patient safety, check the probe connection and position periodically.

Do not use a probe if you see any signs of damage.

Monitoring Carbon Dioxide

This chapter describes CO₂ monitoring on the CARDELL® 9500 HD veterinary monitor.

The carbon dioxide (CO₂) measurement principle is infrared transmission, where the intensity of infrared light passing the respiratory gas is measured with a photo detector. As some of the infrared light is absorbed by the CO₂ molecules, the amount of light passing the gas probe depends on the concentration of the measured CO₂.

Use the CO₂ measurement to monitor a patient's respiratory and ventilation status. The CO₂ measurement produces:

- a CO₂ waveform.
- an end tidal CO₂ (etCO₂) value: the CO₂ value measured at the end of the expiration phase.
- an inspired minimum CO₂ (inCO₂) value: the smallest value measured during inspiration.
- an airway respiration rate (awRR): the number of breaths per minute, calculated from the CO₂ waveform.

9.1 Selecting CO₂ Accessories

The monitor measures partial pressure of carbon dioxide in a patient's expired gas using the mainstream (Capnostat®) or sidestream (LoFlo®) method.

9.2 Zeroing the Sensor

Zeroing is required when a sensor is used for the first time or whenever you suspect incorrect CO₂ readings. Zeroing the sensor will temporarily interrupt the CO₂ waveform. This does not affect the CO₂ measurement. Generally, after an airway adapter is installed, it is only required to be zeroed once.

To zero the sensor manually:

1. Expose the sensor to room air and keep it away from all sources of CO₂ including the ventilator, the patient's breath and your own.
1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Parameter** menu and press the knob.
3. Rotate the navigation knob to **CO2** and press the knob.
4. Rotate the navigation knob to **Zero Adapter**.
5. Press the navigation knob to start zeroing. When the monitor is not ready, the screen shows **Not Ready Try Later**. If the monitor is ready, when you press the navigation knob the monitor displays **Zero In Progress**. When the zeroing process completes, the monitor shows **Zero Successfully**.
6. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

Note — Other measurement functions are not available during CO₂ zeroing.

9.3 Measuring Carbon Dioxide

To start a CO₂ measurement:


1. Set correct CO₂ alarm limits. See Changing the CO₂ Alarm Settings on page 9-4.
2. Insert the CO₂ sensor connector into the receptacle on the side panel of the monitor (see The Side Panels on page 2-5). Make sure the arrows on the connector are at the top of the connector and line up the two keys of the connector with the receptacle and insert.
3. Choose the appropriate airway adapter/filter line and connect it to the sensor head. The airway adapter/filter line clicks into place when seated correctly.
4. Wait two minutes to allow sensor to warm up.
5. Zero the transducer, if necessary. See Zeroing the Sensor on page 9-1. When the zero calibration is finished you can begin monitoring.
6. Install the airway adapter/filter line at the proximal end of the circuit between the endotracheal tube and the Wye section.

To remove the connector, grasp the body portion of the connector back and remove. Do not remove by pulling cable.

Note — Always connect the airway adapter to the sensor before inserting the airway adapter into the breathing circuit. In reverse, always remove the airway adapter from the breathing circuit before removing the sensor.

The purple airway adapter/filter line is for use with 4.0 mm and smaller ET tubes. The clear airway adapter/filter line is for use with 4.5mm and larger ET tubes.

9.3.1 Capnostat® 5 Mainstream CO₂ Sensor



Warning To prevent stress on the endotracheal tube, support the sensor and airway adapter.

Position sensor cables and tubing carefully to avoid entanglement or potential strangulation. Do not apply excessive tension to any cable.


Replace the airway adapter, if excessive moisture or secretions are observed in the tubing or if the CO₂ waveform changes unexpectedly without a change in patient status.

To avoid infection, use only sterilized, disinfected or disposable airway adapters.


Inspect the airway adapters prior to use. Do not use the airway adapter if you see any damage to the airway adapter.

9.3.2 LoFlo® Sidestream CO₂ Sensor


Connect the airway adapter or sample line to the sensor. It will click into place when seated correctly.



Install the airway adapter at the proximal end of the circuit, between the endotracheal tube and the Wye section.



For intubated patients with an integrated airway adapter in the breathing circuit, connect the mail luer connector on the straight sample line to the female port on the airway adapter.



Warning Position sensor cables and tubing carefully to avoid entanglement or potential strangulation. Do not apply excessive tension to any cable.

Replace the filter line, if excessive moisture or secretions are observed in the tubing or if the CO₂ waveform changes unexpectedly without a change in patient status.


Filter lines should be replaced after 3-4 months of use to keep moisture from entering the sensor.

Inspect the airway adapters prior to use. Do not use the airway adapter if you see any damage to the airway adapter.

Using the Sidestream Sensor Holder:

The sidestream sensor holder can be used to clamp the sensor onto an IV pole or a shelf.

1. Push the sensor into the holder until it clicks into position.
2. Clamp the holder onto an IV pole, a shelf, or another appropriate location.
3. To remove the sensor from the holder, release the clip and pull the sensor out of the holder.



Removing Exhaust Gases from the System:

Use an exhaust tube to remove the sample gas to a scavenging system. Attach it to the sidestream sensor at the outlet connector.

Warning Regarding Anesthetics: When using the sidestream CO₂ measurement on patients who are receiving or have recently received anesthetics, connect the outlet to a scavenging system, to avoid exposing the veterinary staff to anesthetics.

9.4 The CO₂ Display

Note — CO₂ waveform displays in the RESP waveform channel when CO₂ is connected.

The following figure shows an example of the CO₂ waveform and numeric pane.




Figure 9-1 CO₂ Waveform

- 1 -- CO₂ waveform label
- 2 -- CO₂ waveform speed
- 3 -- CO₂ waveform
- 4 -- EtCO₂ indicator
- 5 -- EtCO₂ value
- 6 -- InCO₂ indicator
- 7 -- InCO₂ value

9.5 Changing the CO₂ Display

The CO₂ display can be changed for waveform speed and mode.

9.5.1 Changing the CO₂ Waveform Speed

The speed of waveform scanning on the screen is determined by the waveform speed configured. For more information on changing the speed of waveform, see [Changing the Waveform Speed](#) on page 2-8.

The options for waveform speed are: **25 mm/s**, **12.5 mm/s** and **6.25 mm/s**.

9.5.2 Changing the CO₂ Waveform Display Mode

To change the CO₂ waveform mode, see [Changing the Waveform Display](#) on page 2-8 for more information.

9.6 Changing the CO₂ Alarm Settings

CO₂ alarms include the etCO₂ alarms, the inCO₂ alarms and respiration rate alarms.

To turn on CO₂ audible alarms, see [Turning on Audible Alarms](#) on page 3-3. To turn off CO₂ audible alarms, see [Turning off Audible Alarms](#) on page 3-4.

For more information on CO₂ alarm limits, see [Setting Individual Alarm Limits](#) on page 3-5

9.7 Safety Information

Warning Use only approved accessories to ensure accurate CO₂ measurements.

Never use CO₂ measurements in the presence of flammable uncontained anaesthesia gases mixed with air and/or Oxygen or Nitrous Oxide.

Never measure CO₂ in the presence of pharmaceuticals in aerosols.

EtCO₂ measurement accuracy may decrease temporarily while performing electro-surgery or defibrillation. This does not affect patient or equipment safety.

Leakages in the breathing system or sampling system may cause the displayed etCO₂ values to be too low. Always connect all components securely and check for leaks according to standard clinical procedures.

10

Managing Data

This chapter describes how to manage and view patient data on the CARDELL® 9500 HD veterinary monitor.

The monitor can store measurement data in the trend database. You can view the trend information in graphic form through the trend display.

Measurements (for example, SpO₂ and temperature) are captured and stored in the trend database every 4 minutes or at the time that the measurement is started.

10.1 Viewing Trends

Trends display the patient data in graphic form. The Trend Display window opens the most recent data in the past 48 hours.

To enter the **Trend Display** window:



1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Trend Display** menu and press the knob.
3. Rotate the navigation knob to **Exit** and press the knob.




Figure 10-1 Trend Display

In the window, the parameter values of the trend display on the left of the graphics. The buttons at the bottom of the screen change the displaying time scale and the trend segments.

Parameter values are plotted vertically along the y-axis of the display. A time range (48 hours) displays along the horizontal x-axis. The window displays 24 hours of data per page. The most recent measurements appear on the right side of the graph and the data scrolls to the left so that the oldest data appears on the far left.

In the **Trend Display** window, you can check trend data:

1. Rotate the navigation knob to **Label** and press the knob to confirm the selection. Rotate the navigation knob to check trends and press the knob again.
2. Rotate the navigation knob to **Up/Down** and press the knob to confirm the selection. Rotate the navigation knob to check trends of different parameters and press the knob again.
3. Rotate the navigation knob to **Page** and press the knob to confirm the selection. Rotate the navigation knob to check trends on different pages.

10.2 Viewing Data List


The monitor stores and allows you to check no less than 600 groups of the data in a list.

10.2.1 Viewing All Measurement Data

In the Data List window, you can recall all the measurement data. The monitor shows the stored parameter values in 4-minute intervals. NIBP measurement data is saved in real-time.

The following figure shows an example of the data list window:

Date	Time	SYS	DIA	SPO2	HR/PR	RR	TEMP1	TEMP2	EtCO2	InCO2
12-18	18:45	---	---	---	---	---	---	---	---	---
12-18	18:49	---	---	---	---	---	---	---	---	---
12-18	18:53	---	---	---	---	---	---	---	---	---
12-18	18:57	---	---	---	---	---	---	---	---	---



To check all measurement data:

1. Press the navigation knob when only the **Main** menu icon is on the main screen to enter the data list window.
2. Press the navigation knob again to exit the data list window.

11

Printing

This chapter describes how to print waveforms and data on the CARDELL® 9500 HD veterinary monitor.

The monitor Printer prints real time parameter values and waveforms. The speed, content, time length, timing interval, and mode of printing are all configurable. You can print up to three waveforms at the same time. If the monitor is configured to print when an alarm occurs, the system automatically prints the related alarm data.

11.1 About the Printer

The following figure shows the built-in printer of the monitor.




Figure 11-1 Printer

- 1 -- Outlet slot for printer paper
- 2 -- Printer latch
- 3 -- Printer door

11.2 Loading the Printer Paper

Use only thermal paper with a width of 50 mm.

To load the printer paper:

1. Insert your finger into the printer latch (marked OPEN) and lift it to open the printer door.




Figure 11-2 Opening the Printer Door

2. Insert a new roll on the holder so that the end comes from under the roll, extending beyond the Printer door.



Figure 11-3 Loading Paper

3. Pull at least 25 mm of the paper out of the door, removing any slack.
4. Firmly close the printer door.
5. To test the printer and ensure that the paper is correctly loaded, start printing. If no printout generated, the paper might be loaded backwards. Try reloading the paper.

Caution Use only the printer paper specified by the manufacturer. Using wrong paper can damage the printer.

11.3 Printing

To turn on or off the printer:



- Press the **Print** key on the front panel to turn on recording manually.
- To stop printing, press the **Print** key again when recording is ongoing.

If the alarm-triggered printing or the timed printing is enabled, the Printer starts to print at a specific starting point and stops when the configured print time is reached.

Note — Printing stops if the Printer paper runs out or for technical problems.

11.4 Turning on Alarm-triggered Printing

The alarm-triggered recording only prints waveforms. If the alarm trigger menu is enabled, the Printer automatically prints waveforms when an upper or lower alarm limit is exceeded. Alarm-triggered recordings are based on the **Print Time** setting configured.

To turn on the alarm-triggered recording:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Print** menu and press the knob.
3. Rotate the navigation knob to **Alarm Trigger** menu and press the knob to confirm the selection.

4. Rotate the navigation knob to an option (**ON** or **OFF**) and press the knob.
5. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

11.5 Setting Timed Printing

Note — Timed Printing is only available for waveform printing.

With timed printing enabled, the monitor generates a recording of waveforms when a preset time occurs.

To turn on **Timed Printing**:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Print** menu and press the knob.
3. Rotate the navigation knob to **Timed Printing** menu and press the knob to confirm the selection.
4. Rotate the navigation knob to an option (**5 minutes**, **10 minutes**, **15 minutes**, **30 minutes**, **1 Hour**, **4 Hours**, **8 Hours** or **12 Hours**) to turn on timed printing and set the interval.
5. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

11.6 Setting the Print Time

When you set alarm-triggered or timed printing, the monitor generates printouts according to the preset alarm print time.

To set the alarm print time:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Print** menu and press the knob.
3. Rotate the navigation knob to **Print Time** menu and press the knob to confirm the selection.
4. Rotate the navigation knob to an option (**8s**, **16s** or **32s**) and press the knob to confirm.
5. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

11.7 Setting the Advance Time

When alarm-triggered or timed printing is enabled, you can set an advance time that generates a printout containing the waveforms that occurred from a short period of time (2 seconds or 5 seconds) before the **Print** key was pressed.

11.8 Selecting the Print Mode

The monitor has three print modes: **Data**, **Waveform** and **Timed**.

When it is configured to **Data** mode, the monitor prints data list. When it is configured to **Waveform** or **Timed** mode, the monitor prints waveforms.

Data: the monitor prints the most recent 10 groups of data in the data list window.

Waveform: the monitor prints real-time parameter waveforms for a maximum of 90 seconds, or until you press the **Print** key again.

Timed: the monitor prints according to the length of time specified in the **Print Time** menu. For example, if **Print Time** is set to 8 seconds, the Printer prints eight seconds of waveforms after the **Print** key is pressed.

Note — When **Waveform** print mode is selected, the Printer creates a 90-second recording and stops if no user interaction is taken. This does not include alarm-triggered printing and the timed printing.

11.9 Changing the Print Speed

To change the print speed:

1. Rotate the navigation knob to access the **Main** menu.
2. Rotate the navigation knob to the **Print** menu and press the knob.
3. Rotate the navigation knob to **Print Speed** menu and press the knob.
4. Rotate the navigation knob to an option (**12.5 mm/s**, **25 mm/s** or **50 mm/s**) and press the knob to confirm the selection.
5. Rotate the navigation knob to **Exit** and press the knob to exit the menu.

11.10 Sample Recording Strip

The following figure shows an example of a recording printouts.




Figure 11-4 Sample Recording Printout

11.11 Printer Paper Guidelines

Keep the Printer paper properly for future review. For long term storage, make copies of the printouts.

Observe the following precautions when keeping the printouts:

- Keep the Printer paper in a dry, dark and cool place.
- Avoid exposure to:

- high temperature, humidity, direct sunshine, and fluorescent lighting and radiation.
 - ambient temperatures over 50°C.
 - liquids.
 - starch, adhesive agent, oily paper or carbon paper.
 - ethane, oxide, plastic rubber, adhesive tape, fluorescent ink or red inks.
- Never squeeze, apply friction or scratch the paper.
- Never place two printing surfaces together.

12

Cleaning and Care

This chapter gives general guidelines on the cleaning and care of your monitor and accessories. Use only the approved cleaning methods and agents listed in this chapter. The warranty does not cover damage caused by using unapproved substances. In addition, qualified professionals should perform periodic calibration and maintenance of the monitor.

12.1 General Guidelines

Keep the monitor, cables and accessories free of dust and dirt. After cleaning and disinfecting, check the equipment carefully. Do not use if you see signs of deterioration or damage. If you need to return any equipment to the manufacturer, decontaminate it first.

Observe these general precautions:

- Always dilute cleaning agents according to the instructions in this chapter or use the lowest possible concentration.
- Never allow liquid to enter the case.
- Never immerse any part of the equipment in liquid or allow any liquid to enter electrical contacts.
- Never use abrasive material (such as steel wool or silver polish).
- Never use bleach on electrical contacts or connectors.
- Never use strong solvents such as acetone or acetone-based compounds.
- Never autoclave, steam sterilize, or ultrasonically clean the monitor or cables.

Caution Use the equipment according to instructions accompanied with the sensors.

If you spill liquid on the exterior of the monitor, use a clean cloth to dry the monitor. If you believe the liquid may have gotten inside the monitor, power off the monitor and contact your service representative.

Use the accessories according to instructions accompanied.

12.2 Cleaning the Monitor

To clean the monitor:

1. Turn off the veterinary monitor and unplug the power cord.
2. Clean the monitor at room temperature. Use a soft cloth and any of the following cleaning agents:
 - Isopropyl alcohol (70% solution in water)
 - Mild soap and water
 - Sodium hypochlorite (chlorine bleach), 3% solution in water
3. Dry the cleaned parts with a cloth.

Note — Do not permit any liquid to enter the monitor case and avoid pouring liquids on the monitor while cleaning. Do not allow water or cleaning solution to enter the measurement connectors.

12.2.1 Cleaning the Display

Occasionally, as needed, clean the display window using a soft, lint-free cloth sprayed with an alcohol free glass cleaner. Do not use either isopropyl alcohol or solvent to clean this display. Use of these cleaners can cause damage to the display. The use of paper towels is not recommended as it may scratch the surface.

Caution Use care when cleaning the display. Do not use a paper towel to clean the display as this may cause scratches.

12.3 Cleaning the Accessories

When possible, clean the accessories according to the accessory manufacturer's instructions. Use the general guidelines in this chapter if you do not have specific product cleaning instructions.

12.3.1 Cleaning the Cables

As necessary, the cables should be cleaned with a mild detergent and water to remove excess bioburden. When necessary, the cables may be disinfected using a soft cloth saturated with a 10% (1:10) solution of chlorine bleach in tap water or 70% isopropyl alcohol. Dry the cables with a soft cloth.

12.3.2 Cleaning the SpO₂ Sensor

To clean the SpO₂ sensor:

1. The sensor may be surface-cleaned by wiping it with an agent such as 70% Isopropyl Alcohol.
2. The clip may be cleaned by either wiping it with, or soaking it for ten (10) minutes in, 70% Isopropyl Alcohol. If the clip is soaked, be sure to rinse it with water and air-dry it prior to use on the next patient.
3. Then dry all sensor and cable surfaces with a clean, dry cloth. Set the sensor aside allowing it to dry completely

Note — After each cleaning and prior to each use, inspect the sensor and cable for fraying, cracking, breakage, or other damage. Inspect the clip for cracking or breakage, or loss of spring tension that would allow slippage or movement of the sensor from its proper position.

12.3.3 Cleaning the NIBP Cuff

Prior to each patient use, inspect the blood pressure cuff and its hose for damage.

To clean the NIBP cuff:

1. Clean the cuff with detergent (mild soap water) or dilute bleach solution (1-2%),
2. Wipe down with soft cloth or air-dry the cuff.

Note — Midmark does not recommend submersion of the cuff. Liquid should not be permitted to enter the cuff bladder because instrument damage may occur. The cuff should be allowed to thoroughly dry before use.

12.3.4 Cleaning the Temperature Probe

As necessary, the probe should be cleaned with a mild detergent and water to remove excess bioburden. When necessary, the probe may be disinfected using a soft cloth saturated with a 10% (1:10) solution of chlorine bleach in tap water or 70% isopropyl alcohol. Dry the probe with a soft cloth.

Caution Check the temperature probe and do not use if you see signs of deterioration or damage.

12.3.5 Cleaning the ECG clips

The clips may be cleaned with 70% isopropyl alcohol.

Note —The conductivity of the ECG clips will be reduced if electrolyte gel or other material is allowed to dry and build up on clips

12.3.6 Cleaning the CO₂ Accessories

Cleaning the outside of the Capnostat 5 or LoFlo CO₂ sensor :

1. Use a cloth dampened with isopropyl alcohol 70%, a 10% aqueous solution of sodium hypochlorite (bleach), or mild soap. Take care not to allow moisture to enter the LoFlo sensor.
2. Wipe down with a clean water-dampened cloth to rinse and dry before use.

Caution **Make certain that the sensor windows are clean and dry before reuse. Only use lens cleaner or 70% alcohol with lens paper to clean sensor windows or irreparable damage may occur to the CO₂ windows.**

Cleaning Capnostat airway adapters:

1. Clean by rinsing in a warm soapy solution followed by soaking in a liquid disinfectant such as 70% isopropyl alcohol, a 10% aqueous solution of sodium hypochlorite (bleach), a glutaraldehyde 2.4% solution such as Cidex®, Steris System 1® or ammonia. It should then be rinsed out with sterile water and dried.
2. May be disinfected using the methods listed below:
Immerse and soak in Cidex® or equivalent 2.4 glutaraldehyde solution for a 10 hour soak.
Immerse and soak in Perasafe® or equivalent peracetic acid .26% solution for a 10 minutes.
3. Before reusing the adapter, ensure the windows are dry, residue free and that the adapter has not been damaged during handling or the cleaning/disinfecting process.

Note —Capnostat airway adapters are intended for single patient use in human medicine, but in the veterinary setting, may be reused as long as any moisture is allowed to dry between uses. It is estimated that with fairly regular use, the adapters should be replaced every 60 to 90 days.
Refer to the manufacturer's instructions enclosed with each sensor for more information.

Caution **DO NOT insert any object, such as a brush, into the CAPNOSTAT 5 CO₂ airway adapter. Irreparable damage may occur to the adapter windows.**

Note —The LoFlo filter lines should not be cleaned other than to wipe any debris off of the adapter piece

12.4 Cleaning the Print Head

After long term use, deposits of paper debris may collect on the print head making recordings unevenly fainter in horizontal stripes on the printouts, and shortening the life of the print head and roller.

To clean the print head:

1. Use electrostatic discharge (ESD) precautionary measures.
2. Power off the monitor and open the Printer door to remove the paper roll from the Printer.
3. Soak a cotton swab in alcohol (70%) and gently wipe the surface of the thermal printhead.
4. Clean the roller using the same method above.
5. Reload the paper roller after the alcohol is dried and close the Printer door.

Note — Do not use any hard metal or sharp tool (for example, a screwdriver) to remove debris on the print head.

12.5 Disinfecting the Monitor

To avoid long-term damage to the equipment, disinfect the equipment as determined by hospital's policy.

Warning Never mix disinfecting solutions (such as bleach and ammonia) as hazardous gases may result.

Always shut off the monitor and clean it before performing any disinfection.

Clean the equipment before disinfecting. Wipe the monitor after disinfection.

You can use the following disinfectants:

- Ethanol (up to 70%), 1- and 2- Propanol (up to 70%)
- Glutaraldehyde (up to 3.6%)

Note — Disinfect the accessories according to the accessory manufacturer's instructions. Use the general guidelines in this chapter if you do not have specific product cleaning instructions

12.6 Sterilizing the Monitor

Sterilization is not recommended for this monitor, related products, accessories or supplies unless otherwise indicated in the Instructions for Use that accompany the accessories and supplies.

Maintaining the Battery

This chapter describes how to charge and maintain the battery on the CARDELL® 9500 HD veterinary monitor. You can power the monitor with either an AC power source or its internal battery.

Warning Only use batteries specified by the manufacturer. Other batteries can cause monitor damage or personal injury.

13.1 Charging LED



A battery symbol on the front panel indicates the charging LED.

- Flashing green: the monitor is connected to an AC power source and charging the battery.
- Steady green: the monitor is fully charged.
- Off: the monitor is not connected to an AC power source.

When the monitor is charging the battery, the charging LED flashes. The battery needs 8 hours before it is fully charged. After the charging is complete, the charging LED is steady green.

Note — Always leave the monitor plugged into an AC power source when not in use. This ensures that the battery is charged. Check battery status at least once a month.

13.2 Battery Indicator



When the monitor is powered by the internal battery, a battery power gauge displays on the bottom right corner of the screen and shows the remaining battery capacity. More lines inside the indicator mean a larger capacity remains.

The battery capacity depends on the method and frequency of use. A new, fully charged battery can support a continuously working monitor for at least 1.5 hours. However, NIBP measurements and printing might use more battery power.

When the gauge has only one or no lines, the gauge flashes and the monitor provides an alarm, suggesting low battery. The monitor will shut off in 5 minutes if you do not plug it into an AC power source to charge the battery.

Note — When the battery does not have enough power to support normal monitor use, the monitor shuts off and will not function until you plug it into an AC power source.

13.3 Installing the Battery

Caution Always remove the battery when transporting the monitor.

To remove the battery:

1. Power off the monitor and unplug the power cable or other wires.
2. Open the battery cover.



Figure 13-1 Open the Battery Cover

3. Rotate the battery holder counterclockwise to release the battery from the monitor.



Figure 13-2 Remove the Battery

To install the battery:

1. Check the polarity indicators on the battery cover and insert the new battery according to the indicators.
2. Push the battery into the monitor and rotate the battery holder clockwise to hold it in place.
4. Install the battery cover.

Warning Ensure that the battery is completely inserted into the monitor, the battery holder is in place and that the battery door is firmly closed. An improperly installed battery could result in serious injury to small patients if it falls out of the monitor.

13.4 Recharging the Battery

Before putting the monitor into service, fully charge the battery. You can charge a battery in a monitor used to monitor patients.

To charge the battery:

1. Insert the battery into the monitor. See Installing the Battery on page 13-2.
2. Connect the monitor to an AC power source.

The charging LED flashes, indicating the charging status.

3. Charge the battery at least 8 hours. When the charging LED is steady green and the battery power gauge is filled, continue charging the battery for at least 2 more hours.

13.5 Optimizing the Battery

The battery life depends on the time and frequency of usage. Correctly maintain and use the battery to extend the battery life span.

13.5.1 Conditioning the Battery

Before putting a new battery into use, perform at least two complete reconditioning cycles. This improves the battery life and increases performance specifications. A reconditioning cycle is complete (charging and discharging the battery) when the monitor powers off from lack of battery power. Perform a recondition when battery performance is shortened or when the battery is stored for 2 months or more.

To condition the battery:

1. Disconnect the monitor from the patient and turn it off.
2. Place the battery you want to condition into the monitor battery compartment (see Installing the Battery on page 13-2).
3. Connect the monitor to AC power and continuously charge the battery for more than 8 hours.
4. Disconnect the monitor from AC power and power the monitor with the battery until the monitor turns off.
5. Connect the monitor to AC power and continuously charge the battery for more than 8 hours.

Conditioning is complete.

13.5.2 Checking Battery Performance

Battery performance decreases with increased use. To monitor battery performance:

1. Disconnect the monitor from the patient and turn it off.
2. Connect the monitor to AC power and continuously charge the battery for more than 8 hours.
3. Disconnect the monitor from AC power and power the monitor with the battery until the monitor turns off.
4. Observe how long the battery works. The length of time the battery can power the monitor indicates the current battery performance.

If the battery does not work according to specifications, replace the battery or contact your service representative.

13.6 Disposing of a Battery

When a battery displays signs of damage or when it no longer charges, replace it. Dispose of the used battery in accordance with local laws and regulations.

13.7 Safety Information

Warning Never open a battery or incinerate it. The battery could ignite, explode, leak or heat up, causing personal injury if not disposed of properly.

Never heat the battery above 104°F.

Keep the battery out of the reach of children and in its original package until you are ready to use it.

Under extreme conditions, the battery could leak, causing corrosive liquids to enter eyes and burn skin. If this occurs, flush with water and seek immediate medical attention.

Maintenance and Troubleshooting

This chapter describes how to maintain and troubleshoot the CARDELL® 9500 HD veterinary monitor.

Warning

Failure on implementing a satisfactory maintenance schedule for this equipment could cause undue equipment failure and potential health hazards.

If you discover a problem with any of the equipment, contact your service personnel.

14.1 Inspecting the Equipment and Accessories

Perform a visual inspection before every use.

1. Examine the unit exterior for general physical condition. Make sure that the monitor is clean and the housings and connectors are not cracked or broken.
2. Inspect all accessories (cables, transducers, sensors and so on). If any show signs of damage, do not use.
3. Examine all system cables, the power plug and cord for damage. Make sure that the prongs of the plug do not move in the casing. If damaged, replace it with an appropriate one.
4. Inspect the patient cables, leads and their strain reliefs for general condition. Make sure there are no breaks in the insulation. Make sure that the connectors are properly engaged at each end to prevent rotation or other strain.

14.2 Maintenance Task and Test Schedule

The following tasks are for qualified service professionals only. Ensure that these tasks are carried out as indicated by the monitor’s maintenance schedule, or as specified by local laws. Clean and disinfect equipment to decontaminate it before testing or maintaining it. There is no required maintenance during warranty period.

Table 14-1 Maintenance Task and Test Schedule

Maintenance and Test Schedule	Frequency
Safety checks	Once every two years and after repairs where the unit has been opened (front and back separated) or the monitor has been damaged by impact.
NIBP calibration	Once every two years, or more often if specified by local laws.

Table 14-1 Maintenance Task and Test Schedule (Continued)

Maintenance and Test Schedule	Frequency
CO ₂ calibration and performance test	First calibration at 1,200 hours of continuous use. After first calibration, once every two years or after 4,000 hours-whichever comes first. Perform following any instrument repairs or the replacement of any instrument parts. Replace entire CO ₂ module after 20,000 hours of continuous use.
Performance assurance for any other measurements	Once every two years, or if you suspect the measurement is incorrect.
Battery	See Chapter 13, Maintaining the Battery.

14.3 Troubleshooting

- If you suspect a problem with an individual measurement, read the related section in this guide and double check the measurement settings.
- If you suspect an intermittent, system-wide problem, call Midmark

Use the following tables to troubleshoot problems with the monitor and accessories..

Table 14-2 ECG Problems

Problem	Solution
Inaccurate HR measurements	• Check the patient ECG signal. • Check/adjust lead position. • Check/re-clean skin. • Check/replace ECG electrode. • Check ECG wave amplitude (GAIN).
The lead wires are connected but there is no ECG waveform and the screen displays Leads Off or No Signal .	• Ensure that all alligator clips are in good contact with the patient body. • Check all the extension parts of the ECG lead wires. The circuit between the plugs and pins of the 5 ECG extension cables are open. Replace the lead wire if the circuit is open. • Check the ECG lead wires and replace the lead wires after long-term use. • The ECG might have a problem communicating with the host. Turn off and restart the monitor. If the problem continues, contact your supplier. • The ECG lead mode might be set on 5 instead of 3
The monitor does not show the ECG baseline waveform consistently on the screen and the waveform moves out of the display area at times.	• Ensure that the environment is dry and that the monitor is not wet. If there is dampness in the environment or monitor, run the monitor continuously for 24 hours to dry it out. • Check the quality of the lead wire/clip and whether the patient site where the clip is attached is clean. Replace the lead wire/clip or clean the contact site.

Table 14-3 RESP Problems

Problem	Solution
There is RESP interference or noise	Turn off and restart the monitor. If the problem continues, contact Midmark technical support.
The RESP waveform is too weak for stable observation.	• Check the ECG lead wires/clips for correct placement. • Check the patient site for cleanness. • Select a RESP lead suitable for the patient.

Table 14-4 SpO₂ Problems

Problem	Solution
There is no SpO ₂ waveform or values.	• Check the position of the sensor on the patient. • Check whether the red light of the probe is flashing. If there is no red light flashing, there might be poor contact. Check the extension cable and the connector. • Check whether the patient's limb is under pressure. Never take blood pressure and SpO₂ measurements on the same limb. • Check whether the environmental temperature is too low. Do not expose the patient's limb to cold air since this can affect readings.
The SpO ₂ values turn on and off during SpO ₂ monitoring.	In long term monitoring or operation, patient movement may result in SpO ₂ interruptions. Keep the patient stabilized. Patient limb motion resulting in SpO ₂ interruptions is normal.
No Signal displays in the SpO ₂ waveform channel	The SpO ₂ function has a problem with the host communication. Turn off and restart the monitor. If the problem continues, contact Midmark technical support.
Search Timeout displays in the SpO ₂ waveform channel.	• Check for a loose SpO₂ probe or move the probe to another site. • If the problem remains, reconnect the SpO₂ probe and reboot the monitor. • If the problem remains, contact your Midmark technical support.

Table 14-5 NIBP Problems

Problem	Solution
Air Leak displays in the NIBP channel.	• Check the cuff or hose for leakage.
Measurement Error displays in the NIBP channel.	There is a malfunction during measurement and the monitor cannot process the measurement data. • Check the patient's status and connection. • Replace the cuff and restart the measurement.

Table 14-5 NIBP Problems (Continued)

Problem	Solution
Cuff Loose or Cuff Position Err. displays in the NIBP numeric pane.	• Ensure the proper cuff size is being used. • Check the cuff position. • Check the extension tube for breakage and if the tube is damaged, replace it.
Overpressure displays in the NIBP channel	Check the hose for kinks or jams and restart the measurement.
Pressure range exceeded.	A blood pressure reading has exceeded the pressure range. The monitor cannot take any more NIBP readings until the alarm is acknowledged.
Measurement Error displays in the NIBP numeric pane.	This error can be caused by: system self-testing, excessive patient and/or limb movement, or air leakage. Consult Midmark technical support.

Table 14-6 TEMP Problems

Problem	Solution
There is no temperature displayed.	The TEMP sensor might be off the patient or the monitor. • Check the sensor connection. • If you cannot solve the problem, replace the sensor. • If the problem continues, consult Midmark technical support.

Table 14-7 CO₂ Problems

Problem	Solution
There is no CO ₂ waveform or EtCO ₂ /InCO ₂ displayed.	• Check whether the accessories are connected correctly then reboot the monitor. • If the problem continues, consult Midmark technical support.

Table 14-8 Printer Problems

Problem	Solution
The Printer does not work.	• Make sure the Printer has paper. • If the Printer has paper, reboot the monitor. • If the problem continues, consult your service representative. • Make sure the paper is loaded correctly

Table 14-9 USB Exporting Problems

Problems	Solutions
The USB stick is not detected by the monitor or an error occurs during data exporting.	• Check whether the USB stick is connected correctly. • Remove the USB stick from the monitor and re-connect it. • Format the USB stick to FAT 16 or FAT 32, and then run data exporting. • If the problem continues, try a new USB stick formatted in FAT 16 or FAT 32, and then run data exporting again. • If the problem continues, consult Midmark technical support.

If you suspect an intermittent, system-wide problem, call Midmark technical support.

14.4 Disposing of the Monitor

Warning To avoid contaminating or infecting personnel, the environment or other equipment, make sure you disinfect and decontaminate the monitor appropriately before disposing of it in accordance with your country's laws for equipment containing electrical and electronic parts. For disposal of accessories, where not otherwise specified, follow local regulations regarding disposal of hospital waste.



Do not dispose of waste electrical and electronic equipment as unsorted municipal waste. Collect it separately, so that it can be safely and properly reused, treated, recycled, or recovered.

15

Warranty/Service

15.1 Limited Warranty

REGISTRATION

Register your monitor at midmark.com/vet-register to:

- Log your product warranty with Midmark
- Keep you informed on important warranty information and product updates
- Provide you with faster, more convenient service in the event you experience a problem
- Enhance customer service benefits tailored to your product and account

SCOPE OF WARRANTY

Midmark Corporation (“Midmark”) warrants to the original purchaser its Veterinary products and components (except for components not warranted under “Exclusions”) manufactured by Midmark to be free from defects in material and workmanship under normal interior use and service. Midmark’s obligation under this warranty is limited to the repair or replacement, at Midmark’s option, of the parts or the product defects of which are reported to Midmark within the applicable warranty period and which, upon examination by Midmark, prove to be defective.

APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original user, shall be as follows:

• **The Cardell® 9500HD unit is warranted against defect in material and workmanship for a period of two years from the time of delivery. Monitor Accessories are warranted against defect in material and workmanship for a period indicated below from the time of delivery:**

- C-Stat Mainstream and LoFlo Sidestream CO₂ Modules carry a 2 year warranty
- Battery and temperature probes carry a 1 year warranty
- V-Sat SpO₂ sensors carry a 9 month warranty
- Blood pressure cuffs, CO₂ sidestream sampling lines, CO₂ mainstream adapters and ECG cable/wire sets shall be free of defects at the time of delivery to the customer.

EXCLUSIONS

This warranty does not cover and Midmark shall not be liable for the following; (1) repairs and replacements required because of misuse, abuse, negligence, alteration, accident, freight damage, or tampering; (2) changes in color caused by natural or artificial light; (3) products which are not installed, used, and properly cleaned as required in the Midmark written “User” and or “Quick Reference” Guide applicable for this product; (4) products considered to be of a consumable nature such as light bulbs and surge suppression product; (5) accessories or parts not manufactured by Midmark; (Any warranties on these items are extended directly by the manufacturer of these items to the original purchaser. Information on these manufacturers’ warranties will be enclosed with the applicable products. Also, Midmark will furnish copies of any of the warranties extended by any such manufacturers upon request.) (6) specially manufactured products; (7) damage caused by animals; (8) charges by anyone (including Midmark’s authorized dealers) for adjustments, repairs, replacement parts, installation, or other work performed upon or in connection with such products which is not expressly authorized in writing in advance by Midmark.

EXCLUSIVE REMEDY

Midmark’s only obligation under this warranty is the repair or replacement of defective parts. Midmark shall not be liable for any direct, special, indirect, incidental, exemplary or consequential damages or delay including, but not limited to, damages for loss of profits or loss of use.

NO AUTHORIZATION

No person or firm is authorized to create for Midmark any other obligation or liability in connection with the products.

THIS WARRANTY IS MIDMARK’S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

15.2 After-sale service

To obtain service or product support, please contact Midmark at 800-643-6275 or visit the website at midmark.com. Have the following information available:

- Model and serial number of the equipment
- Date of purchase and distributor name

A

Specifications

This appendix describes the specifications of the CARDELL[®] 9500 HD veterinary monitor.

A.1 ECG Specifications

The following table indicates the ECG specifications.

Table A-1 ECG Specifications

Parameter	Specifications
CMRR (common mode rejection ratio)	≥ 90 db
HR measurement and alarm range	15 bpm - 300 bpm
HR accuracy	±1 bpm or ±1%, whichever is greater
ECG channel	1 channel
ECG leads	3-lead and 5-lead, user selectable
Alarm indication	Visual and audible indicators
QRS indicator	Audible tone and heart symbol flashes
ECG speed	12.5 mm/s, 25mm/s, 50 mm/s
ECG gain	X 1/4, X 1/2, X 1, X 2, X4, Auto
Bandwidth	Monitor mode: 0.5 Hz - 40 Hz Operation mode: 1 Hz - 25 Hz Diagnosis mode: 0.5 Hz - 100 Hz
Heart rate alarm response time	< 12 seconds
Defibrillator protected	Tested with 5 kV
Recovery time following defibrillation	< 5 seconds
Signal delay	≤ 0.5 seconds
Auxiliary current (leads off detection)	Active electrode: < 100 nA Reference electrode: < 900 nA
Response time of heart rate meter to change in heart rate	HR change from 80 bpm to 120 bpm: 8 seconds maximum HR change from 80 bpm to 40 bpm: 10 seconds maximum

Table A-1 ECG Specifications (Continued)

Parameter	Specifications
Time to alarm for tachycardia	Vent tachycardia (1 mVpp, 206 bpm): Gain 0.5, range 5.5 to 8.5 seconds, average 7.2 seconds Gain 1.0, range 6.1 to 7.5 seconds, average 6.5 seconds Gain 2.0, range 5.9 to 6.7 seconds, average 6.3 seconds Vent tachycardia (2 mVpp, 195 bpm): Gain 0.5, range 5.4 to 6.2 seconds, average 5.8 seconds Gain 1.0, range 5.7 to 6.5 seconds, average 6.1 seconds Gain 2.0, range 5.3 to 6.1 seconds, average 5.7 seconds

A.2 Respiration Specifications

The following table indicates the respiration specifications.

Table A-2 Respiration Specifications

Parameter	Specification
Work mode	Thoracic impedance
Respiration rate measurement and alarm range	0 rpm - 150 rpm
Apnea alarm duration	5 seconds - 120 seconds
Respiration accuracy	± 2 rpm
Respiration speed	6.25 mm/s, 12.5 mm/s, 25 mm/s
Respiration gain	X 1/2, X 1, X 2

A.3 SpO₂ Specifications

The update rate for the SpO₂ value and pulse rate is 6 to 7 seconds. Data averaging and other signal processing on the displayed and transmitted data values of SpO₂ and pulse rate is not more than 20 seconds. Depending on the magnitude of difference between the alarm limit and the displayed value, the alarm signal generation delay may be from 1 second to 20 seconds.

Note — It is not necessary to have an SpO₂ calibration when the monitor is in use.

The following table indicates the SpO₂ specifications.

Table A-3 SpO₂ Specifications

Parameter	Specification
SpO ₂ measurement and alarm range	1% - 100%
SpO ₂ average	8 beat average
Accuracy	+/-2% (70-100%), +/-3% (50-69%) (Tested basing on Nellcor V-SAT digital lingual sensor provided with small and large clip)
SpO ₂ pulse rate range	30 bpm -300 bpm
SpO ₂ pulse rate average	8 seconds
SpO ₂ pulse rate accuracy	+/-3% ±1 digit
SpO ₂ pulse rate sound	PR tone
Wavelength range ¹	Red light: 660nm Infrared light: 900nm
Maximum optical output power	15 mW for all specified sensors
Pulse rate display	Numeric

1. Information about wavelength ranges can be useful for clinicians performing photo dynamic therapy.

A.4 Temperature Specifications

The following table indicates the temperature specifications.

Table A-4 Temperature Specifications

Parameter	Specification
Temperature measurement and alarm range	0°C - 50°C (32°F - 122°F)
Unit	°C or °F
Temperature accuracy	± 0.1°C
Resolution	0.1°C
Updating time	Approximately 1 second

A.5 Non-invasive Blood Pressure Specifications

The following table indicates the NIBP specifications.

Table A-5 NIBP Specifications

Parameter	Specification
Measuring method	Oscillometric method
Measuring parameter	SYS, DIA, MAP, PR
Unit	mmHg or kPa
Work mode	Manual, Auto, Stat
Cycle period	1-10, 15, 30, 60, 90, 120, 180, 240, 480 minutes
Measurement range	20 mmHg - 265 mmHg (2.7 kPa - 35.3 kPa)
NIBP accuracy	Max. Std. Deviation: 8 mmHg (1.1 kPa) Max. Mean Error: ± 5 mmHg (± 0.7 kPa)
Cuff pressure range	0 mmHg - 280 mmHg (0 kPa - 37.3 kPa)
Initial inflation pressure	150 mmHg (20 kPa)
Subsequent cuff inflation	30 mmHg above the last measured Systolic value
Measuring time	Manual/Auto: 25 seconds Stat: 20 seconds
PR measuring range	30 bpm - 300 bpm
PR accuracy	$\pm 5\%$ or ± 5 bpm, whichever is greater

A.6 Carbon Dioxide Specifications

The following table indicates the CO₂ specifications.

Table A-6 CO₂ Specifications

Parameter	Specification
CO ₂ measurement and alarm range	0 mmHg - 150 mmHg (0 kPa - 20 kPa)
Unit	mmHg or kPa
CO ₂ accuracy	0 mmHg - 40 mmHg (0 kPa - 5.33 kPa): ± 2 mmHg (± 0.3 kPa) 41 mmHg - 70 mmHg (5.5 kPa - 9.3 kPa): $\pm 5\%$ of reading 71 mmHg - 100 mmHg (9.4 kPa - 13.3 kPa): $\pm 8\%$ of reading 101 mmHg - 150 mmHg (13.4 kPa - 20 kPa): $\pm 10\%$ of reading
CO ₂ resolution	0 mmHg - 69 mmHg (0 kPa - 8 kPa): 0.1 mmHg 70 mmHg - 150 mmHg (7.3 kPa - 20 kPa): 0.25 mmHg
RR range	0 rpm - 150 rpm
RR accuracy	1 rpm

Table A-6 CO₂ Specifications

Parameter	Specification
RR resolution	1 rpm
Waveform speed	6.25 mm/s, 12.5 mm/s, 25 mm/s

A.7 Power Specifications

The following table indicates the power specifications.

Table A-7 Power Specifications

AC Power Parameter	Specification
Mains power input	AC 100 V - 240 V, 50 Hz/60 Hz
Power input	115 VA
Fuse	T1.6 A, 250 V
Battery Power Parameter	Specification
Capacity	Lead acid: 12 V/2.3 Ah per piece
Charging time	≥ 8 hours
Typical operating time	1.5 hours (monitoring ECG, SpO ₂ and NIBP measured every 15 minutes)
Charge mode	Automatic (with charge protection feature) when the monitor is connected to an AC power source.
Discharge protection	When the monitor is powered by battery, the monitor powers off if the battery power is almost depleted.

A.8 Display Specifications

The following table indicates the display specifications.

Table A-8 Display Specifications

Parameter	Specification
Type	Color TFT LCD
Size	7" (15.41 cm x 8.66 cm)
Channel	3 channels

A.9 Printer Specifications

The following table indicates the Printer specifications.

Table A-9 Printer Specifications

Parameter	Specification
Type	Thermal array Printer
Mode	Timed, Data and Waveform
Channel	3 channels
Speed	12.5 mm/s, 25 mm/s, 50 mm/s

A.10 Physical Specifications

The following table indicates the physical specifications.

Table A-10 Physical Specifications

Parameter	Specification
Net weight	4.0 kg (8.8 lbs)
Gross weight	7.0 kg (16.4 lbs)
Monitor dimensions	225 mm (w) X 170 mm (d) X 230 mm (h) (8.86" (w) X 6.69" (d) X 9.05" (h))

A.11 Environmental Specifications

The following table indicates the environmental specifications.

Table A-11 Environmental Specifications

Parameter	Specification
Temperature	Operating: 5 °C - 40 °C (41 °F - 104 °F) Storage: -20 °C - 40 °C (-4 °F - 122 °F)
Relative humidity	Operating: ≤ 90% (with Printer: ≤ 80%)
Barometric pressure	Operating: 70.0 kPa - 106.0 kPa Storage & transportation: 50.0 kPa - 106.0 kPa

B

Alarm Specifications

This appendix describes the physiological and technical alarm messages for the CARDELL[®] 9500 HD veterinary monitor. Each table provides alarm priority and descriptions for each alarm. The appendix also provides the factory default alarm limits.

B.1 Alarm Classifications

All alarms (technical and physiological alarms) are classified into three levels for different priorities.

- High (H)
- Medium (M)
- Low (L)

B.2 Physiological Alarms

Warning When a high priority alarm occurs, an alarm message appears on the screen. Middle and low priority alarms do not display alarm messages on the screen. Instead, the line under the corresponding parameter value and alarm icon flashes.

The following table lists physiological alarms alphabetically of the patient monitor.

Table B-1 Physiological Alarms

Alarm Message	Priority	Cause
HR: too high/too low ¹	M	HR value exceeds the alarm limits. Check the patient.
Asystole ²	H	No QRS is detected. Check the patient, electrodes and leads.
RESP: too high/too low ¹	M	RESP value exceeds the alarm limits. Check the patient.
Apnea ³	H	A user-defined time has elapsed since the last detected breath.
SpO ₂ : too high/too low ¹	M	SpO ₂ value exceeds the alarm limits. Check the patient.
ultra-low ⁴	H	SpO ₂ value drops more than 10 points below the low alarm limit. Check the patient.
SpO ₂ : PR too high/too low ¹	M	SpO ₂ PR value exceeds the alarm limits. Check the patient.
SYS: too high/too low ¹	M	SYS value exceeds the alarm limits. Check the patient.

Table B-1 Physiological Alarms (Continued)

Alarm Message	Priority	Cause
DIA: too high/too low ¹	M	DIA value exceeds the alarm limits. Check the patient.
SYS-DIA: out of range ¹	H	Difference between the Systolic and Diastolic blood pressure values is lower than the preset value. Check the patient.
Press. Rge. Exceed ⁵	M	The blood pressure value exceeds the pressure range limit. Check the tube and the patient.
Temp: too high/too low ¹	M	Temp value exceeds the alarm limits. Check the patient.
EtCO ₂ : too high/too low ¹	M	EtCO ₂ value exceeds the alarm limits. Check the patient.
InCO ₂ : too high/too low ¹	M	InCO ₂ value exceeds the alarm limits. Check the patient.
Apnea ⁶	H	A user-defined time has elapsed since the last detected breath.

1. When the alarm is activated, there is no alarm message, but there is a yellow underline under the corresponding parameter value and the yellow underline flashes.
2. When the alarm is activated, the alarm message will display in the ECG waveform channel.
3. When the alarm is activated, the alarm message will display in the RESP waveform channel.
4. When the alarm is activated, the alarm message will display in the SpO₂ waveform channel.
5. When the alarm is activated, the alarm message will display in the NIBP numeric pane.
6. When the alarm is activated, the alarm message will display in the CO₂ waveform channel.

B.3 Technical Alarms

The following table lists technical alarms alphabetically of the veterinary monitor.

Table B-2 Technical Alarms

Alarm Message	Priority	Cause
No Signal ¹	M	There is an ECG module failure or communication between the monitor and the module has failed.
Leads Off ¹	M	Not all the ECG electrodes are connected.
Leads Off ²	M	Not all the ECG electrodes are connected.
No Signal ³	M	There is a SpO ₂ module failure or communication failure between the monitor and the module.
Probe Off ³	M	The SpO ₂ sensor is not on the patient site. The monitor or the sensor has failed, or an incorrect sensor is being used.
Searching timeout ³	M	The monitor is checking for SpO ₂ data.
Overpressure ⁴	M	The NIBP tube may be blocked.
Cuff Position Err. ⁴	M	The blood pressure cuff is in wrong position.
Air Leak ⁴	M	There is an air leakage.
Measurement Error ⁴	M	There is an error in the NIBP measurement.
Cuff Loose ⁴	M	The blood pressure cuff is not connected or placed correctly.
Battery alarm ⁵	L	The monitor battery is low.

1. When the alarm is activated, the alarm message will display in the ECG waveform channel.
2. When the alarm is activated, the alarm message will display in the RESP waveform channel.
3. When the alarm is activated, the alarm message will display in the SpO₂ waveform channel.
4. When the alarm is activated, the alarm message will display in the NIBP numeric pane.
5. Alarm message not displayed. Only battery indicator flashes.

B.4 Default Alarm Limits

The following table lists default alarm limits of the veterinary monitor.

Table B-3 Default Alarm Limits

Parameter		
HR/PR (bpm)	High	180
	Low	60
SpO ₂ (%)	High	100
	Low	90
SYS (NIBP) (mmHg/kPa)	High	160 (21.3)
	Low	70 (9.3)
DIA (NIBP) (mmHg/kPa)	High	100 (13.3)
	Low	40 (5.3)
SYS-DIA (mmHg/kPa)	Low	20 (2.7)
TEMP (°F/°C)	High	104.0 (40)
	Low	96.8 (36.0)
EtCO ₂ (mmHg/kPa)	High	60 (8.0)
	Low	20 (2.7)
InCO ₂ (mmHg/kPa)	High	10 (1.3)
	Low	0
RR (rpm)	High	55
	Low	5
APN Alarm Time (s)		20

C

Default Settings

This appendix documents the most important default settings of your monitors as it is delivered from the factory. If your monitor has been ordered preconfigured to your requirements, the settings at delivery will be different from those listed here.

C.1 System Default Settings

The following table describes the system default settings.

Table C-1 System Default Settings

Configuration	Default
Patient name	No Patient ID
HR/PR sound	5
Brightness	3

C.2 Alarm Default Settings

The following table describes the alarm default settings.

Table C-2 Alarm Default Settings

Configuration	Default
Alarm volume	5
Alarm mode	Auto
Alarm on/off	On

C.3 ECG Default Settings

The following table describes the ECG default settings.

Table C-3 ECG Default Settings

Configuration	Default
ECG lead	II
ECG gain	X1
Waveform speed	25 mm/s

Table C-3 ECG Default Settings (Continued)

Configuration	Default
Pacemaker detect	Off
Lead mode	Three lead
ECG mode	Operation
HR/PR priority	ECG
ECG alarm	On

C.4 Respiration Default Settings

The following table describes the respiration default settings.

Table C-4 Respiration Default Settings

Configuration	Default
RESP lead	RA-LA
RR alarm	On
RESP gain	X1/2
Waveform speed	12.5 mm/s
Waveform mode	Standard
Apnea alarm time	20 s

C.5 SpO₂ Default Settings

The following table describes the SpO₂ default settings.

Table C-5 SpO₂ Default Settings

Parameter	Default
SpO ₂ alarm	On
Waveform speed	25 mm/s
Waveform mode	Standard
HR/PR priority	ECG

C.6 NIBP Default Settings

The following table describes the NIBP default settings.

Table C-6 NIBP Default Settings

Parameter	Default
Cuff Size	Small
NIBP alarm	On
Work mode	Manual
Unit	mmHg
Calibration	Auto

C.7 Temperature Default Settings

The following table describes the temperature default settings.

Table C-7 TEMP Default Settings

Parameter	Default
Temp alarm	On
Unit	°F

C.8 Carbon Dioxide Default Settings

The following table describes the CO₂ default settings.

Table C-8 CO₂ Default Settings

Parameter	Default
RR alarm	On
Waveform speed	12.5 mm/s
Waveform mode	Standard
Apnea alarm time	20 s

C.9 Printer Default Settings

The following table describes the Printer default settings.

Table C-9 Printer Default Settings

Parameter	Default
Recording	Off
Print speed	12.5 mm/s
Channel 1 waveform	ECG
Channel 2 waveform	SpO ₂
Channel 3 waveform	RESP/CO ₂
Print mode	Timed
Alarm-triggered recording	Off
Alarm print time	8 s
Advance time setup	2 s
Timed printing	Off

D

Glossary

This appendix defines some terminology and abbreviations used in the monitor and this Instructions for Use.

D.1 Terminology Definitions

Table D-1 Terminology Definitions

Term	Definition
AO	Aortic pressure
ART	Arterial pressure
aVL	Left arm augmented lead
aVR	Right arm augmented lead
BSA	Body surface area
CO ₂	Carbon dioxide
DIA	Diastolic
ECG	Electrocardiograph
EtCO ₂	End-tidal carbon dioxide
FA	Femoral artery pressure
HR	Heart rate
InCO ₂	Inspiratory carbon dioxide
LA	Left arm
LL	Left leg
MAP	Mean arterial pressure
NIBP	Non-invasive blood pressure
PR	Pulse rate
RA	Right arm
RESP	Respiration
RL	Right leg
RR	Respiration rate
SpO ₂	Arterial oxygen saturation from pulse oximetry

Table D-1 Terminology Definitions (Continued)

Term	Definition
SYS	Systolic pressure
TEMP	Temperature

E Tips

E.1 Blood Pressure

All white vinyl cuffs are considered small cuffs and are to be used for small animals. The two dark colored nylon cuffs are for larger mammals. A test run can be made by placing a small cuff around the base of your thumb holding it over your heart, but don't be concerned about the numbers. For more accurate measurement, use the largest cuff.

Key points to keep in mind in order to achieve effective blood pressure readings:

- Determine appropriate cuff size by measuring the circumference of the limb. The width of the cuff should be about 40% of the circumference of the limb. The 2.5cm size should work on most adult cats (2.0cm for small cats). For awake cats, the cuff can be placed around the base of the tail or on the forelimb between the carpus and elbow. For measurements during anesthesia, do not use the tail site. Use the Cardell® Cuff Selector™ as a guideline.
- Wrap the cuff snugly. Make sure it doesn't slip during readings. (See section 7 on cuff and site selection.) If during the measurement, the monitor indicates an effort or pulse rate is abnormally low, reposition the cuff and make sure it is securely wrapped.
- If the Velcro fails, check that the cuff is snug, that the monitor is set to the proper size, and if the Velcro is worn, replace the cuff, as using tape can result in inaccurate readings.
- The cuff position should be at the approximate level of the right atrium. Other than in very large dogs, this is not an issue with anesthetized patients. However, when assessing blood pressures in awake animals, the cuff should be on a horizontal plane with the heart. Keep the animal comfortable and minimize movement of the appendage where the cuff is placed.
- Allow the patient to relax for a few minutes with the cuff in place. Then make a couple of measurements so that the animal can become accustomed to the procedure.
- Using a subdued environment- no bright lights, no white coats, and not in a noisy part of the hospital or around other animals, take 4 or 5 readings at 1-2 minute intervals, throw out the "wild one" and average the others.

E.2 ECG

- 3-Lead systems- Use all 3 leads for ECG recordings.
- 5-Lead systems - Remove the exploring lead (brown) for ECG monitoring unless being used for detailed diagnostic measurements.
- Standard right forelimb lead is negative and standard left hind leg is positive in Lead II. If these leads are properly placed and the monitor is set to Lead II, an upright complex should be displayed.
- It is not necessary for hair to be clipped over the lead site. If hair is thickly matted, only lightly clip, but not down to skin.
- Prepare the lead sites with alcohol, and then use the supplied electrode paste/gel for better contact.
- For cats and small dogs, place a thin layer of alcohol-saturated cotton over the site, apply gel and attach leads.
- Avoid allowing the metal clips to come in contact with conductive surfaces such as the table or other parts of the body.
- Other monitors and motors should not be close to the Cardell as can cause interference. Fan motors and circulating water pumps are a common source.
- If interference is encountered, turn off each motor one-by-one to determine the source, and try moving the Cardell to different locations. Depending on the model, you can change the filtering selection to remove some of the baseline "noise" or interference.
- Always be sure to keep the clips clean of hair and debris, and inspect for any breaks in the wire prior to use.

E.3 SPO₂

- Extend the life of the sensor by hanging it (not winding or wrapping it) for storage and be sure to unclip the sensor for removal or before moving the patient.
 - In small patients, values may drop below 90% due to restriction of blood flow to the measuring site by the probe. Reposition the sensor periodically to avoid this problem.
 - Do not use on an extremity with a blood pressure cuff or catheter in place.
 - Cover site with opaque material (like blue surgical towels) to minimize interference from bright lights directed toward the patient.
 - If slippage occurs on the tongue, use one layer of gauze around the tongue and then replace the sensor.
 - Alternative sites include lip, ear, prepuce, vulva, webbing between toes in non-pigmented skin.
 - For long term use, it is best to move sensor to new site every few hours to preserve skin integrity.
 - Clean sensor and clips with 70% isopropyl alcohol before using on another patient.
- Note that some vasoconstrictive drugs (alpha-2s, for instance), may limit the ability to obtain readings on peripheral areas. You may need to reposition the sensor or try another site.

E.4 CO₂

- Wipe clean the sensor/adaptor after each use and during long procedures. However, do not submerge them or clean them with any chlorhexadine or other abrasive cleaners. This may cause the plastic to become clouded in appearance and also may harm the line itself.
- For the LoFlo the CO₂ filter lines should be switched out between patients and hung up to air dry completely. Continuing to use the same line after each patient could result in a build up of condensation and cause inaccurate readings or occlusion alarms. Replace the filter lines approximately every 3 months to ensure safe, accurate and reliable readings.
- Check the sensor by breathing through the opening. Inspired reading should be nearly zero and expired values should be in the 30's. A large expired breath will bring it closer to 40 mmHg.
- When attaching the sensor between the endotracheal tube adaptor and breathing tube of the anesthetic unit, make sure the fitting is tight and secure.
- When evaluating measurements, any rise in inspired carbon dioxide values should be considered abnormal. This may indicate a problem with ET tube placement, rebreathing of CO₂ from the anesthesia circuit, depleted CO₂ absorbent, etc.
- For cats and small dogs, the exotic/neonatal adaptor should be used to minimize dead space imposed by the sensor.

E.5 Temperature

- The temperature probe can be used esophageally or rectally, but be sure to flag the probe with a label indicating which site has been designated to avoid contamination in subsequent uses.
- Use gauze with 70% isopropyl alcohol to clean the probe. Do not saturate or submerge your probe in solutions.
- To store the probe between uses, hang it rather than winding or wrapping it.

F

Accessories

The following items are included in the standard monitor kit and can be reordered from your distributor or directly from Midmark using the associated reorder codes:

Reorder number	Description	Included
SV1	2.0 cm bp cuff, for limb circumference 3-6cm (white vinyl fabric)	1
SV2	2.5 cm bp cuff, for limb circumference 4-8cm (white vinyl fabric)	2
SV3	3.5 cm bp cuff, for limb circumference 6-11cm (white vinyl fabric)	3
SV4	4.0 cm bp cuff, for limb circumference 7-13cm (white vinyl fabric)	3
SV5	5.0 cm bp cuff, for limb circumference 8-15cm (white vinyl fabric)	2
SV8	8.0 cm bp cuff, for limb circumference 13-20cm (green nylon fabric)	1
SV10	10.2 cm bp cuff, for limb circumference 18-26cm (blue nylon fabric)	1
01-03-0162	6' bp hose with quick disconnect	1
V-SAT	6' Nellcor SpO2 lingual sensor and clips (large and small)	1
01-02-0183	10' Nellcor extension cable	1
016-1602-00	ECG Trunk Cable	1
016-1604-00	3-lead ECG wire set	1
ECG-A3	MAX-12/9500 ECG alligator clip	3
590004	MAX-1/MAX-12/9500 Flexible Esophageal/Rectal Temp probe	1
PAPER-4F	Thermal Printer paper (4 rolls/tube)	1
PL-200	Redux gel for ECG use	1
BATTERY-7C	MAX-12/9500 rechargeable sealed lead acid battery	1
01-02-0395	Power cord, domestic 10'	1
015-3091-00	Grounding wire	1

The following are optional accessories for use with the 9500HD series:

Reorder number	Description
C-Stat 5	Capnostat Mainstream CO2 probe (incl. airway adapters 6063 & 6312)
6063-00	CO2 Small animal airway adapter (<5cc dead space)
6312-00	CO2 Exotic animal airway adapter (<1cc dead space)
LoFlo	LoFlo Sidestream CO2 module (incl. sampling lines 3473INF & 3473ADU)
3473INF-00	LoFlo airway adapter w/filter – small/neonate
3473ADU-00	LoFlo airway adapter w/filter – adult/pediatric
3475-00	LoFlo sampling line w/filter and male leur adapter
MaxFast-1	Nellcor MaxFast Reflectance sensor & posey wrap
8008-002	9500HD Rolling Stand w/basket and mounting plate
9A465002	9500 Mount for Matrx VMS Plus anesthesia machine
9A465004	9500 Pole mount
EP-XS	MAX-12/9500 series small esophageal probe (ecg, temp, resp.)
EP-S	MAX-12/9500 series small esophageal probe (ecg, temp, resp.)
EP-L	MAX-12/9500 series small esophageal probe (ecg, temp, resp.)

A

advance time, 11-3
 alarm
 adjust volume, 3-4
 apnea setting, 5-3
 default limit, 3-6
 flashing icons, 3-2
 flashing LED, 3-2
 individual limit, 3-5
 levels, 3-1
 limit range, B-4
 messages, 3-2
 mode, 3-5
 pause, 3-4
 physiological, B-1
 PR setting, 6-4
 technical, B-3
 types, 3-1
 auditory alarm, 3-2
 turn off, 3-4
 turn on, 3-3

B

battery
 charging LED, 13-1
 check performance, 13-3
 condition, 13-3
 dispose, 13-3
 optimize, 13-3
 recharge, 13-3
 big number screen, 2-8
 brightness, 2-10

C

clean
 accessory, 12-2
 cable, 12-2
 guideline, 12-1
 monitor, 12-1
 NIBP cuff, 12-2
 Printer printhead, 12-4
 SpO2 sensor, 12-2
 temperature probe, 12-3
 CO2
 accessory, 9-1
 measure, 9-2
 waveform and numeric, 9-5
 waveform mode, 9-5
 waveform speed, 9-5
 CO2 sensor
 zero, 9-1
 components, main, 1-2
 configurations, monitor, 1-2

D

data list, 10-2
 date and time, 2-9
 default
 alarm setting, C-1
 CO2, C-3

ECG, C-1
 NIBP, C-3
 Printer, C-4
 RESP, C-2
 SpO2, C-2
 system, C-1
 TEMP, C-3
 disinfect
 monitor, 12-4

E

ECG
 1mV calibration, 4-5
 cable, 4-4
 lead mode, 4-3
 mode, 4-6
 waveform and numeric, 4-4
 waveform gain, 4-5
 waveform speed, 4-5
 electrode
 place, 4-1

H

HR/PR priority, 4-6

I

inspect
 monitor and accessory, 14-1

M

main components, 1-2
 maintenance
 task and schedule, 14-1
 manage patient, 2-10
 measure limit, 7-2
 menu, 2-7
 monitor
 dispose, 14-5
 front panel, 2-3
 inspect, 2-2
 power supply, 2-2
 screen display, 2-6
 setting, 2-1
 side panel, 2-5
 version, 2-2
 monitor configurations, 1-2

N

NIBP
 alarm setting, 7-8
 calibrate, 7-8
 cuff, 7-1
 interval time, 7-6
 measure, 7-4
 measurement data, 7-7
 mode, 7-5
 cuff size, 7-3
 unit, 7-7
 waveform and numeric, 7-7

P

- paced patient, 4-7
- print
 - mode, 11-3
- print time, 11-3

R

- Printer paper
 - load, 11-1
- recording
 - alarm triggered, 11-2
 - speed, 11-4
 - strip, 11-4
- RESP
 - abdominal breathing, 5-1
 - cardiac overlay, 5-1
 - lateral chest expansion, 5-1
 - lead, 5-1
 - waveform and numeric, 5-2
 - waveform gain, 5-2
 - waveform mode, 5-3
 - waveform speed, 5-2

S

- skin preparation, 4-1
- specification
 - CO₂, A-4
 - display, A-5
 - ECG, A-1
 - environmental, A-6
 - NIBP, A-4
 - physical, A-6
 - power, A-5
 - Printer, A-6
 - RESP, A-2
 - SpO₂, A-2
 - TEMP, A-3
- SpO₂
 - apply sensor, 6-1
 - cable, 6-2
 - measure, 6-2
 - suspicious reading, 6-2
 - waveform and numeric, 6-3
 - waveform speed, 6-3
- sterilize
 - monitor, 12-4
- system maintenance, 2-9

T

- TEMP
 - influencing factor, 8-1
 - measure, 8-1
 - unit, 8-2
- terminology, D-1
- timed printing
 - turn on, 11-3

- trend, 10-1
- troubleshooting, 14-2

V

- volume
 - HR and PR, 2-10

W

- waveform
 - freeze, 2-9
 - mode, 2-8
 - speed, 2-8



453564356211