solo T4-Vet Test

For vet use only!

EUROLyser
DIAGNOSTICA

In Vitro Diagnostic test kit for determination of T4 with solo analyser.

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Order number: VT 0140
Order number: VT 1400

Indication

solo T4-Vet-Test Kit
solo T4-Vet-Control Kit

Kit size

16 tests/pack
2 x 1 ml (low high range)

Test kit preparation: Allow single test min. 10 minutes to warm to room temperature by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

Thyroxine (T4) is produced by the thyroid. Controlled by a complex regulatory mechanism it is released into the bloodstream. Most of the thyroxine in blood is bound to various carrier proteins such as albumine, prealbumine and other specific binding proteins. Diseases of the thyroid with abnormal hormone concentrations occur more likely in older animals. Elevated or lowered concentrations of the thyroid hormones, e.g. thyroxine (T4), lead towards problems in numerous other systems. Hyperthyreosis are merely found in cats, hypothyreosis are found in all species.

Method

The solo T4-Vet test is a homogenous enzymatic immunoassay. By using an 8-anilino-1-naphthalene acid sulfure (ANS) the bound T4 is separated from its binding proteins. The separated thyroxine released from blood competes afterwards with a thyroxine that is linked to the Glucose-6-phosphat dehydrogenase enzyme to a specified binding area. The GAPDH activity at 340nm is associated with the complete thyroxine concentration in blood. If no T4 is in the blood sample the enzyme activity is inhibited by the binding of a specific antibody.

Sample material

Serum or plasma (EDTA, Heparin, Citrat). Sample volume according to the menu settings on solo and according to the use of a 20µl or 5µl pipette.

Measurement Range

When using a 20µl pipette: 0-8 µg/dl (0-96 nmol/l)
When using a 5µl pipette: 8-40 µg/dl (96-580 nmol/l)

Warnings

This testkit is for in vitro diagnostic use only. Do NOT INGEST. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbin to form explosive compounds. Take the necessary precautions for the use of laboratory reagents.

Stability and Storage

Store at 2 - 8°C. DO NOT FREEZE. Kit can be used until expiry date printed on the labels.

Waste Management

Please refer to local legal requirements.

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English

1. Prepare testsystem

- 1.1 Insert RFID-card into analyser
- 1.2 Place ERS cuvette into solo testkit rack.
- 1.3 Place ERS cap into solo testkit rack.
- 1.4 Press symbol on touch screen, required information is entered by use of the touch display of the analyser

2. For further details please refer to the according section in the solo user manual:

Choose sample type:

- 2.1a 20µl pipette: Possible measuring range:
0-8 µg/dl (0-96 nmol/l) OR...
- 2.1b 5µl pipette: Possible measuring range:
8-40 µg/dl (96-480 nmol/l)
- 2.2 Start of analysis by pressing the start button

WARNING!

Do not wipe pipetting tip after sample collection! Allow test minimum 10 minutes to warm up to room temperature before inserting the test into the analyser!

3. Sample preparation serum or plasma (serum primary tube, lithium-heparin or EDTA plasma)

- 3.1 Pipette 20µl or 5µl serum/plasma from primary tube
- 3.2 Pipette 20µl or 5µl sample into ERS cuvette INTO LIQUID
- 3.3 Place ERS cap firmly on ERS cuvette
- 3.4 Place ERS cartridge into solo analyzer
- 3.5 Start automatic sample processing by pressing the  start button.