

EEG-Z3 (SA7680Z)

User Guide



Thought Technology Ltd.

2180 Belgrave Avenue, Montreal, Quebec H4A 2L8 Canada
Tel: (800) 361-3651 • (514) 489-8251 Fax: (514) 489-8255
E-mail: mail@thoughttechnology.com
Website: www.thoughttechnology.com

CAUTION

U.S. federal law restricts this device to sale by or on the order of a physician or any other practitioner licensed by law of the state in which he or she practices, to use or order the use of this device.

To diminish the risk of spreading communicable diseases, always use good hygiene practices with reusable EEG electrodes, particularly if abrasive substances are used. In all cases, refer to your facility's infection control procedure.

IMPORTANT

Do not connect inputs or outputs to line powered devices except through the fiber optic cable.

INTENDED PURPOSE

Intended for Biofeedback and Muscle Re-education.

Not For Diagnostic Purposes.

TABLE OF CONTENTS

FUNCTIONAL DESCRIPTION	1
SENSOR COMPONENTS.....	1
CONNECTING THE SENSOR	1
ABOUT EEG MODE	2
ABOUT EP / 0.01 HZ MODE (ALSO CALLED "LOW FREQUENCY" MODE).....	2
ABOUT SCP / DC MODE	2
ZEROING THE SENSOR WHEN THE SIGNAL DRIFTS OUT OF RANGE	3
ZEROING THE SENSOR.....	3
ZEROING PROCEDURE:	4
SELECTING ALTERNATIVE BASELINES IN SCP/DC MODE	4
MISCELLANEOUS NOTES.....	4
DEVICE SPECIFICATIONS	5
PLACING ORDERS.....	6
TECHNICAL SUPPORT	6
WARRANTY	6
RETURNING EQUIPMENT.....	7
REPAIR RETURN FORM	8

Manual No. SA7682 rev. 0
© copyright Thought Technology Ltd. -2009

Functional description

The EEG-Z3 can be used on its own for EEG neurofeedback or in combination with the TT-AV Sync sensor (T7670) and a push-button switch (SA7660). When combined, these sensors comprise a system designed for evoked and slow cortical potentials and Go/NoGo motor response protocols.

The EEG-Z3 can be used in any of the following three modes. All three modes support impedance checking.

- Standard mode for regular EEG, similar to the EEG-Z sensor.
- EP/0.01 Hz (0.01 Hz frequency cutoff) for evoked, event-related potentials and slow cortical potentials.
- SCP/DC mode for measuring slow cortical potentials.

A Set/Zero button allows quick switching between modes. LED indicators identify the active mode. When it is powered off, the sensor "remembers" the mode in which it was last used.

In all three modes, if the signal approaches or passes the limit of the range that the sensor can measure, the LED will blink rapidly. This is useful in SCP/DC mode, in which the Set/Zero button can be used to return the signal to baseline.

Sensor components

- EEG-Z3 sensor
- Zeroing cable
- Sensor cable (SA9385). For connection of EEG-Z3 sensor to encoder.
- EEG Electrodes

Note: If you will use the EEG-Z3 sensor in EP/0.01Hz or SCP/DC mode, make sure you are using **sintered Silver/Silver Chloride (Ag/AgCL)** electrodes. For other modes, the use of standard EEG electrodes is adequate.

Connecting the sensor



1. Connect the electrodes to the sensor using color-coded DIN cables. Respect the color coding as indicated in the preceding photograph.

Note: When using standard EEG electrodes, use preferred method to affix to the scalp. When using sintered Silver/Silver Chloride (Ag/AgCl) electrodes, follow the procedure outlined in **About SCP / DC mode**.

2. Connect the EEG-Z3 sensor to the encoder using a sensor cable (T9385).

Note: If the EEG-Z3 sensor is being used with the TT-AV Sync for time-locked EEG averaging (such as evoked and slow cortical potential averaging), both sensors must connect to encoder inputs that have the same sampling rate (either both 256 s/s or both 2048 s/s).

3. To select the sensor mode:
 - a) Press and hold the **Set/Zero** button for about 3 seconds. The LED lights will begin to cycle through the 3 modes.
 - b) Release the **Set/Zero** button when the LED lights up for the mode you want to use.

About EEG mode

This mode requires you to zero the sensor (remove the offset, or non-zero baseline) prior to first time use, using the Calibration Wizard in the BioGraph Infiniti software. The procedure is described below under **Zeroing the sensor**. Note that most applications involving signal processing of EEG do not include the very low frequency content of the signal, and this offset removal may not be necessary. It should be done in all case in which the raw EEG or EEG < 2 Hz is used.

About EP / 0.01 Hz mode (also called "low frequency" mode)

When you first enter this mode, a settling period is imposed during which the signal returns from an out-of-range value to the baseline. This normally requires approximately 100 seconds, as it is based on the time constant of the system in this mode.

Note: This mode is sensitive to both eye and head movement, and care should be taken to minimize both actions during sessions in this mode.

About SCP / DC mode

Used for DC measurements and slow cortical potential assessment and training. As is to be expected with DC measurements, in this mode the EEG-Z3 is subject to the DC drift of the system. Since electrode polarization is a major contributor to drift, the following guidelines should be followed to reduce DC drift:

- Prepare skin well using skin prep. Always wipe the area clean with a dry cloth after preparation.
- Spread a thin coat of electrode paste over a small area of prepared skin (*shown right*). This serves to keep track of the site and to prepare for better contact with the electrode.



- Use sintered Ag/AgCl electrodes. Apply a conservative amount of gel to the electrode. Too much gel can contribute to drift, but too little will make poor contact and cause high impedance. The ideal amount of gel should not protrude from the edge of the electrode (*shown below*), and should form a thin ring around the electrode when it is fixed to the scalp.



- Affix the electrode to the pad of paste on the scalp. Press lightly, just enough to fix the electrode and spread the gel slightly around the perimeter of the disc (*shown right*).



Even with special care to reduce electrode polarization, an initial period of polarization is unavoidable. For this reason, after the initial placement of electrodes on the subject, a "settling time" of 10 to 20 minutes is usually necessary before beginning a session.

Zeroing the sensor when the signal drifts out of range

In SCP/DC mode, when the LED blinks to indicate that the signal is approaching the limits of the sensor's range, you can use the Set/Zero button to bring the signal close to the baseline again. To select a different baseline than the middle of the range, follow the guidelines described in **Selecting alternative baselines in SCP/DC mode** on page 4.

Note: *This mode is sensitive to both eye and head movement, and care should be taken to minimize both actions during sessions in this mode.*

Zeroing the sensor

Note: *In this section, zeroing refers to software zeroing.*


To use the sensor in EEG mode with a zero baseline, the EEG-Z3 sensor **must** be zeroed before its first use on a given channel on a given encoder, and then as required to maintain its accuracy. When zeroing is complete for a specific sensor on a specific channel of the encoder, this information is saved on the system.

Zeroing **must** be repeated if:

- A different EEG-Z3 sensor is connected to the same channel of the encoder.
- The same EEG-Z3 sensor is connected to a different channel of the encoder.
- The same EEG-Z3 sensor is connected to a different encoder.

Zeroing procedure:

1. Connect the EEG-Z3 sensor to the encoder. Connect the zeroing cable (shown right) to the EEG-Z3 sensor.
2. Open a session in **BioGraph Infiniti**.
3. From the *Options* menu of the Recording Screen, select **Sensor Calibration**.
4. From the list, select the channel connected to the EEG-Z3 and click **Zero Now** to start the **Zeroing Wizard**. The wizard guides you through the zeroing procedure.



Selecting alternative baselines in SCP/DC mode

The Zero function in this mode is necessary to operate in the presence of DC drift. Since drift tends to be in one direction, there are three baseline levels that you can choose from.

- If the signal is drifting UP, select **Low**. This will place the baseline in the lower half of the sensor range, toward the bottom.
- If the signal is drifting DOWN, select **High**. This will place the baseline in the upper half of the sensor range, toward the top.
- If signal drift is minimal, OR is up and down, OR you don't know the direction of drift, select **Middle**. This will place the baseline close to 0. This is the default baseline setting.

This selection is made by double clicking the Set/Zero button. With a session running in BioGraph Infiniti, double click the button and see what happens. Ensure that there is enough visible range on the y axis of the graph, or that it is in autoscale mode. Each time you double click, the levels will cycle in the following order: **High** (Up), **Low** (Down), and **Zero**, such that you may select the desired baseline level. Once selected, zeroing the sensor by pressing the Set/Zero button once will return the signal to this baseline.

The sensor resets its baseline to the default level (Middle) when it is powered off.

Miscellaneous notes

1. Since the EEG-Z3 cannot be used in SCP/DC mode with a Legacy protocol, it cannot be used with ProComp2 or ProComp+ encoders, or with ProComp Infiniti or FlexComp Infiniti encoders if they have been set to Legacy mode.
2. The phase response of the EEG-Z3 is not well-matched with that of the EEG-Z at very low frequency (less than 5 Hz, more so around 1 Hz). For this reason, we recommend the following:
 - Do not mix EEG-Z3 and EEG-Z sensors when comparing channels for phase at low frequencies. In this situation, use one type or the other.
 - Consider the effects of phase shift whenever the frequency of interest approaches 1Hz for any phase-sensitive measurement.
 - To control for this issue, use multiple EEG-Z3s in low frequency mode for any measurement where uniform phase response is required near or below 1Hz.

Device Specifications

Specification	Value
Accuracy	$\leq 1\%$, $\pm 0.3\mu V_{RMS}$
Input Impedance	Differential: $100G\Omega$ paralleled with $270pF$ Common-mode: $100G\Omega$ paralleled with $200pF$
Noise	$< 0.5 \mu V_{RMS}$
CMRR (excluding CM signal active cancellation)	$> 100dB$
CM active cancellation effect	$> 40dB$ @10-120Hz
Electrode offset tolerance (Slow AC and DC modes)	$\pm 100mV$
Bandwidth, lower cutoff, 3dB EP/.01Hz mode	0.01 Hz
Bandwidth, lower cutoff, 3dB, EEG mode	1.5 Hz
Bandwidth, upper 3dB (all modes)	1600KHz

Placing Orders

Outside USA

Tel: 1-514-489-8251

Fax: 1-514-489-8255

In USA Toll-Free

Tel:1-800-361-3651

E-Mail: mail@thoughttechnology.com

Or contact your local authorized distributor.

Technical Support

Outside USA

Tel: 1-514-489-8251

Fax: 1-514-489-8255

In USA Toll-Free

Tel:1-800-361-3651

E-Mail: techsupport@thoughttechnology.com

Or contact your local authorized distributor.

Warranty

The EEG-Z3 sensor is guaranteed to be free from defects in material and workmanship for 1 year from the date of purchase.

In the unlikely event that repair is necessary, contact Thought Technology Ltd. to receive a Return Authorization number. Then send the unit back by a traceable method. Thought Technology will not be responsible for items not received. We will repair or replace your unit(s) that are still under warranty free of charge.

This warranty does not apply to damage incurred through accident, alteration, or abuse.

Returning Equipment

**Be sure to contact us for a return authorization number (RA)
before returning any equipment!**

1. Send the unit(s) **postage prepaid** and **insured**, with proof of purchase to one of the addresses below.
2. If you are shipping from outside Canada and the USA to Canada, label the package '**Goods to be repaired - Made in Canada**' to avoid unnecessary customs charges. All customs and duties charges will be billed to you if incurred by sending the unit without this label or to the **wrong** address.
3. Provide a **detailed** description of the problem you are experiencing and your telephone/fax number (see form on the last page of this manual).

In the USA, ship insured to:

Thought Technology Ltd.
Cimetra LLC
20 Gateway Drive
Plattsburgh, New York
12901, USA

In Canada and all other countries, contact your dealer or ship insured to:

Thought Technology Ltd.
2180 Belgrave Avenue
Montreal, Quebec
Canada H4A 2L8

Contact information for service



Canada and International +1 514 489-8251



USA 1-800-361- 3651



service@thoughttechnology.com

Repair Return Form

Be sure to contact us for authorization before returning any equipment!

Remove this sheet and include with returned unit(s).

Include copy of original invoice and return to the address in the Returning Equipment section.

Name _____

Company _____

Address _____

Phone No. _____

Fax No. _____

Date Purchased _____

From Whom _____

Model Name _____

Serial Number _____

Problem _____
