Hoefer SQ33 Sequencer

Includes Fan Sensor Kit

Page finder

Important Information	ii
Waste Electrical and Electronic Equipment (WEEE)	ix
1. SQ33 Function and Description	1
2. SQ33 Specifications	2
3. SQ33 Operating Instructions	4
4. Installation and Operation of the Fan Sensor Kit.....	11
5. SQ33 Care and Maintenance.....	13
6. Customer Service Information.....	14
7. Disposal	15
8. Warranty	16
9. Ordering Information.....	18

Important Information – English

SQ33:

- If this equipment is used in a manner not specified by Hoefer, Inc. the protection provided by the equipment may be impaired.
- This instrument is designed for indoor laboratory use only.
- Only accessories and parts approved or supplied by Hoefer, Inc. may be used for operating, maintaining, and servicing this product.
- Only use a power supply that is CE marked or safety certified by a nationally recognized testing laboratory.
- The safety lid must be in place before connecting the power supply leads to a power supply.
- Turn all power supply controls off and disconnect the power leads before removing the safety lid.
- Circulate only water or 50/50 water/ethylene glycol through the heat exchanger if so equipped. Do not connect the heat exchanger to a water tap or any coolant source where the water pressure is unregulated.
- Never introduce antifreeze or any organic solvent into any part of the instrument. Organic solvents will cause irreparable damage to the unit!
- Do not operate with buffer temperatures above the maximum specified technical specifications. Overheating will cause irreparable damage to the unit!

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- Warning! Because this instrument can develop sufficient voltage and current to produce a lethal shock, care must be exercised in its operation.
- This instrument is designed in accordance with the EN61010-1:2001 electrical safety standard. Nevertheless, it should be used only by properly trained operators. Read this entire manual before using the instrument and use only according to the instructions.
- The instrument must always be used with the earth lead of the power cord correctly grounded to earth at the mains outlet.
- Use only undamaged electrical wire and equipment specific for the voltages you will use. All equipment connected to high voltage should be in accordance with EN61010-1:2001.
- Keep the instrument as dry and clean as possible. Wipe regularly with a soft, damp cloth. Let the instrument dry completely before use.
- Do not operate the instrument in extreme humidity (above 80%). Avoid condensation by letting the unit equilibrate to ambient temperature when taking the instrument from a colder to a warmer environment.
- To permit sufficient cooling, ensure that the vents of the instrument are not covered.

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Waste Electrical and Electronic Equipment (WEEE)

English

This symbol indicates that the waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Please contact an authorized representative of the manufacturer for information concerning the decommissioning of your equipment.

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Warning

These units are capable of delivering potentially lethal voltage when connected to a power supply and are to be operated only by qualified technically trained personnel.

Please read the entire operator's manual thoroughly before operating this unit.

These units comply with the statutory EC safety requirements:

- 98/79/EC In Vitro Diagnostic Directive of October 27, 1998
- 73/23/EEC: LOW VOLTAGE DIRECTIVE: IEC 1010-1:1990 plus AMENDMENT 1:1992
- EN 61010-1:1993/BS EN 61010-1:1993 89/336/EC: THE EMC DIRECTIVE FOR ELECTROMAGNETIC COMPATIBILITY

Hoefer single-plate large format gel electrophoresis units are designed to give long service and reproducible results in your laboratory. A few moments spent reading these instructions will ensure that your expectations are reflected in the successful use of the apparatus.

First check that the apparatus has been received complete and undamaged following shipment. Any faults or losses must be reported to Hoefer immediately. Hoefer cannot accept responsibility for goods returned without prior notification.

Refer to the packing list and check that all components and accessories are present.

Please retain all packaging and materials until the warranty period has expired.

1. SQ33 Function and Description

The Hoefer® SQ33 Sequencer offers a wide gel format for increased sample capacity and Shark's Tooth combs to ensure optimal band resolution. The SQ33 Sequencer is perfect for DNA sequencing, DNase and RNase footprinting, and heteroduplex or oligonucleotide analysis. The removable lower buffer chamber combined with the upper buffer chamber drainage tap make buffer disposal easy and safe.

1.1 Unpacking

Unwrap all packages carefully and compare contents with the packing list, making sure all items arrived. If any part is missing, contact your local sales office. Inspect all components for damage that may have occurred while the unit was in transit. If any part appears damaged, contact the carrier immediately. Be sure to keep all packing material for damage claims or to use should it become necessary to return the unit.

1.2 SQ33 Sequencer Components

Description	Quantity
Running Tank + Lid	1
Plain Glass Plates	2
Notched Glass Plates	2
0.35 mm Melinex Spacers	2
0.35 mm Shark's Tooth Comb, 48 well	1

2. SQ33 Specifications

Unit Dimensions (W × D × H)	63 × 21.5 × 50 cm
Plate Dimensions (W × H × D)	33 × 41 × 0.4 cm
Standard Spacer Dimensions (W × H × D)	10 × 410 × 0.35 mm
Active Gel Dimensions (W × H)	29 × 38.5 cm
Maximum Sample Capacity	1 × 96
Recommended Buffer Volume	Upper Buffer Chamber: 1200 ml Lower Buffer Chamber: 1200 ml
Recommended running conditions for 6% denaturing PAGE gel using passive ventilation at 20 °C ambient temperature:	Power: 45–55 W Gel Temperature: 50–60 °C Time: 2.5–3 h
Recommended running conditions for 6% denaturing PAGE gel using FHS-KIT at 20 °C ambient temperature:	Power: 75 W Gel Temperature: 55 °C Time: 1.75 h

2.1 Fan Sensor Kit Specifications

Temperature Control Unit:

Dimensions (H × W × D)	80 × 96 × 140 mm
Temperature Control Range	0–100 °C
Controller Accuracy	+/- 2%
Weight	1.5 kg

Fan:

Dimensions (H × W × D)	38 × 120 × 120 mm
Operating Speed	2,700 rpm
Flow Rate	> 84 CFM
Noise Rating	36 dBA

2.2 Environmental Conditions

- This apparatus is intended for indoor use only.
- This apparatus can be operated safely up to an altitude of 2,000 m.
- The normal operating temperature range is between 4 °C and 65 °C.
- Maximum relative humidity 80% for temperatures up to 31 °C decreasing linearly to 50% relative humidity at 40 °C.
- The apparatus is rated POLLUTION DEGREE 2 in accordance with IEC 664. POLLUTION DEGREE 2, states that: "Normally only non-conductive pollution occurs. Occasionally, however, a temporary conductivity caused by condensation must be expected".

3. SQ33 Operating Instructions

3.1 Regent Preparation and Gel Volumes

For reproducibility and uniform polyacrylamide cross-linking we recommend deionizing, degassing and filtration of acrylamide gel solutions prior to use. Acrylamide solutions should be stored in a cool, dark environment such as a refrigerator and allowed to reach room temperature prior to pouring. Avoid exposure to heat and sunlight.

Polymerization conditions should be adjusted to effect polymerization within about 5 – 10 minutes. Test a small volume in a vial prior to pouring the gel. As a rough guide 100 ml of degassed 6% acrylamide gel will set in about 5 minutes at room temperature when gently mixed with 450 µl of freshly prepared 10% (w/v) Ammonium persulfate plus 200 µl TEMED. The setting time increases to about 10 minutes if the TEMED volume is reduced to 100 µl and to approximately 15 minutes with 75 µl. The amount of catalysts may need to be reduced under warm conditions. Do not pour under direct sunlight.

The SQ33 Sequencer, using 0.35 mm thick spacers, requires approximately 60 mls of acrylamide solution.

3.2 Gel Pouring

1

Remove any greasy fingerprints or water on the edges of the glass plates by wiping with acetone. Lay the plain glass plate on a flat surface and arrange the spacers flush with the edges of the plate. Carefully lay the notched glass plate on top and clamp the plates together. The gel sandwich is now ready for pouring. Refer to the illustration in Figure 1, page 6.

2

Fill a 60 ml syringe with gel mix being careful not to aerate the solution.

3

Hold the syringe as vertical as possible and using a steady continuous motion expel the gel mix along the notched area moving the syringe nozzle repeatedly from one side of the notch to the other. The gel mix should form a continuous pool along the top of the gel space and by capillary action it will move down between the plates. Do not overfill the notched area—the pool should be approximately half the height of the glass ears. Alternatively do not let the pool run dry or air bubbles will be introduced between the glass plates. The boundary of the gel should migrate as a straight line. The gel plates can be tapped slightly behind the moving gel boundary, using your fingers, to prevent any bubble formation. Do not use any hard or metal objects.

4


When the gel boundary reaches the bottom of the glass plates, remove the surplus gel from the notched area. Failure to do so will result in the gel mix dripping from the bottom end. Insert the required comb.

5


If you are using a comb with square wells insert the teeth without trapping bubbles. If you are using a sharks tooth comb insert the flat face of the comb at a slight angle to prevent bubbles from trapped. Add a few drops of gel mix if necessary. Carefully straighten the comb so that it is parallel to the top of the gel plate and reaches 3–5 mm below the notched area.

Figure 1. Pouring gels horizontally.


- Hold syringe vertically.
- Move across notched area to maintain a pool of gel mix.
- Do not overfill.
- Do not attempt to force gel mix down between plates with the syringe jet.

- When gel boundary has reached the end remove surplus gel mix from notched area. Leave just a small amount to ensure bubbles don't get trapped under the comb.

- The comb area may be clamped with bulldog clips.
- Do not let the clips protrude beyond the line of the comb otherwise get thickness variation and hence lane distortion may occur.
- When gel has set carefully remove comb and insert it teeth first (being careful not to damage points) so it just contacts the gel surface.
- The gel is then ready to be attached to the main unit and loaded with samples.

6

The gel polymerization should occur within 10 minutes but allow the gel to polymerize further for at least 90 minutes. If it is to polymerize overnight, enclose wet tissues under a cling film seal to prevent the ends of the gel from drying out.

If you are using a square toothed comb do not remove it until it is submerged below the running buffer. However, if you use a Shark's tooth comb it should be removed before the addition of buffer. Insert it teeth first into the gel slot so that the teeth just come into contact with the gel and then add the running buffer.

3.3 Unit Assembly

1

Place the lower buffer chamber into the main unit. Push it back until it stops. This ensures that the electrode separation and hence field strength are consistent from one run to another.

2

Remove any clamps and the bottom tape, if used.

3

Clamp the gel plate sandwich into position on the apparatus. *Do not over tighten the gel plate clamping screws* as this will lead to the glass plates breaking and also making the insertion and removal of combs difficult.

4

The gel plate assembly will rest against the aluminum back plate. This has been set during manufacture to be parallel to the main frame. *Do not tighten or adjust the aluminum plate securing screws.* The aluminum plate is cushion mounted to provide good contact with the glass plate.

5

Make sure that the upper buffer tank drainage tube is closed.

6

Attach the electrical connectors to the buffer chambers.

IMPORTANT! Do not fill over the maximum fill lines.

3.4 Run Conditions and Buffer Volumes

1

Into the upper buffer chamber pour 800 ml (approximately).

2

Into the bottom buffer chamber pour 800 ml (approximately).

3

Ensure that there is no buffer leakage.

4

If desired pre-run the gel until the glass plates are warm. Use the settings described below.

5

Prior to loading your samples flush out the wells with running buffer to clear them of urea.

6

Heat the sequencing samples at 95 °C for 3 minutes, place on ice and centrifuge 12,000×g for 3 minutes. Return to ice.

7

Load 1–5 µl of sample using a suitable loading tip. If possible avoid taking the sample from the bottom of the tube — particulate materials may cause streaking or smearing. Try to minimize sample dispersion, load the sample directly onto the bottom of the well and keep it as a thin layer.

8

Replace the safety lid ensuring it is positioned fully down over the electrical connectors.

9

Connect and run the gel at the desired power setting. The leads and electrical connectors are tested to 1,500 Volts and users are advised not to exceed this voltage.

IMPORTANT! Make sure that the bottom tank does not overflow.

10

Typically 45–55 Watts constant power will give good results with passive ventilation. If the fan kit is used the power can be increased up to 75 Watts. Ensure that with the buffer system employed the voltage does not exceed the maximum recommended voltage. The resistance of the gel changes throughout electrophoresis when run at constant power so the starting and finishing voltages will be different.

11

In general the gel temperature should reach 50–60 °C. At temperatures above this range there is the risk of generating smeared bands or damaging the gel plates.

3.5 Ending the Run

1

Isolate the power supply before removing the connectors.

2

Remove the safety lid by gripping the edges and pushing down with your thumbs on the pegs located on the top of the unit.

3

Disconnect the leads to the upper and lower buffer chambers.

4

Open the upper buffer tank drainage tube. The buffer will then flow into the bottom buffer tank.

5

Remove the lower buffer chamber and carefully discard the buffer to waste.

6

Dismantle the gel and separate the plates with a strong, thin, broad blade.

7

Note: Do not pry the glass plates apart at the ears. Concentrated leverage in any one area may damage the plates.

8

After removing the gel, clean the plates thoroughly in neutral detergent and rinse with distilled water. A sheet of clean foam rubber placed in the bottom of the sink serves as a useful support and minimizes the risk of plate damage.

9

Rinse buffer chambers in deionized water.
Never use organic solvents.

10

Dry the electrical connectors thoroughly with a tissue before storing the unit.

4. Installation and Operation of the Fan Sensor Kit

4.1 Installation

The kit consists of:

- 1 × 12 volt fan (+4 securing screws)
- 1 × Self adhesive temperature sensor
- 4 × Adhesive cable ties
- 1 × Temperature control unit

①

To fit the fan, simply unscrew and retain the four screws securing the black ventilation cover on the safety lid. Align the fan and note that the fan draws air up through the sequencing unit. An arrow on the fan indicates the direction of the air flow. Secure the fan into position using the long securing screws provided. Replace the black ventilation cover on top of the fan so that the fan blades are not exposed.

②

Set up the polymerized sequencing gel in the unit. Using the adhesive strip, attach the flat face of the heat sensor to the center of the facing glass plate. It is advisable to clean the area with a tissue soaked in acetone prior to attachment. *Do not allow acetone to come into contact with the plastic components.*

③

Route the heat sensor cable up and over the back corner of the gel unit, securing the cable so that it will not drag on the sensor. The fan controller is now ready for use.

4.2 Using the Fan Controller

1

Attach the fan connector to the back of the temperature control box.

2

Attach the heat sensor cable into the remaining socket on the rear of the temperature control box.

3

Select the temperature you require to run the gel at, noting that while the control unit is factory calibrated there will be a temperature gradient across both the glass and aluminum plate. Therefore, it is necessary to set a temperature on the dial approximately 10 °C *lower* than required at the gel plate, *e.g.* if you require to run the gel at 55 °C set the dial to approximately 45 °C. Attach the mains lead to the control unit and turn on the power. The light on the control panel indicates when the contacts are made for the fan operation. When the light is on the fan will be off. When the fan cuts in the light will go out.

4

The fan will automatically switch on/off to maintain the gel temperature at ± 4 °C of the set temperature. After monitoring the gel temperature increase the dial setting as necessary. A temperature indicator strip on the center of your glass plate will allow you to confirm that your dial setting is correct.

IMPORTANT! Acrylic plastic is *not* resistant to aromatic or halogenated hydrocarbons, ketones, esters, alcohol's (over 25%) and acids (over 25%), they will cause "crazing" especially of the UV transparent plastic and should *not* be used for cleaning. *Do not* use abrasive creams or scourers. Dry components with clean tissues prior to use.

5. SQ33 Care and Maintenance

- To remove the safety lid, push thumbs down on the plastic lugs and lift the lid vertically with your fingers.
- Before use clean and dry the apparatus with *distilled water only*.
- Before use, and then on a monthly basis, check the unit for any leaks at the bonded joints. Place the unit on a sheet of dry tissue and then fill with *distilled water only* to the maximum fill line. Any leakage will be seen on the tissue paper. If any leakage is seen, *do not attempt to repair or use the apparatus*, but notify Hoefer technical support immediately.
- After each use, wash all components with a suitable water-soluble disinfectant, to reduce the risk of infection.
- Ensure that the connectors are clean and dry before usage or storage.
- The cooling plate upper surface cannot be cleaned with acetone. Methanol and ethanol can be used.

5.1 Cleaning and Storage

- After use, thoroughly rinse all components in deionized water and gently dry the gold electrode connectors with a soft tissue.
- Ensure that the connectors are clean and dry before usage or storage.
- *Never use organic solvents.*

5.2 Fan Sensor Kit Care and Maintenance

- The unit can be cleaned as required when the mains supply is isolated.
- Cleaning should be carried out with a cloth moistened with water or tissues impregnated with 70% iso-propyl alcohol. *No other cleaning solutions should be used.*
- There are no internal operator serviceable parts in this unit. If the power should fail, the unit must be returned to the authorized service center.

IMPORTANT! Request a copy of the Hoefer, Inc. "Health and Safety Declaration" form before returning the item. No items can be accepted for servicing or return unless this form is properly completed.

Note: A Return Materials Authorization (RMA) number must be obtained from Hoefer, Inc, before returning any item to Hoefer, Inc.

6. Customer Service Information

6.1 Technical Service and Repair

Hoefer, Inc. offers complete technical support for all of our products. If you have any questions about how to use this product, or would like to arrange to repair it, please call or fax your local Hoefer, Inc. representative.

Check the Hoefer, Inc. website at www.hoeferinc.com for the distributor in your area. Or contact us directly at:

Hoefer, Inc.

84 October Hill Road
Holliston, MA 01746

Toll Free: 1-800-227-4750

Phone: 1-508-893-8999

Fax: 1-508-893-0176

E-mail: support@hoeferinc.com

Web: www.hoeferinc.com

7. Disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that it must not be disposed of with unsorted waste.

Instead it is your responsibility to dispose of your equipment correctly at life cycle end by handing it over to an authorized facility for separate collection and recycling.

It is also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so that the persons involved in the disposal and recycling of the equipment are protected from any potential hazard.

For more information about where you can drop off your waste equipment, please contact your local dealer from whom the equipment was purchased originally.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Thank you.

8. Warranty

Hoefer, Inc. warrants that the unit you have received has been thoroughly tested and meets its published specification.

This unit (excluding all accessories) is warranted against defective material and workmanship for a period of twelve (12) months from the date of shipment ex factory.

Hoefer, Inc. will repair all defective equipment returned during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures outlined in this instruction manual and when not having been subject to accident, alteration, misuse or abuse.

No liability is accepted for loss or damage arising from the incorrect use of this unit. Hoefer, Inc.'s liability is limited to the repair or replacement of the unit or refund of the purchase price, at Hoefer, Inc.'s option. Hoefer, Inc. is not liable for any consequential damages.

Hoefer, Inc. reserves the right to alter the specification of its products without prior notice. This will enable us to implement developments as soon as they arise.

Hoefer, Inc. products are for research use only.

A return authorization must be obtained from Hoefer, Inc. before returning any product for warranty repair on a freight-prepaid basis.

8.1 Warning

- *Do not* attempt to remove the outer casing or make repairs to our electrical range of products, should any unit fail.
- Contact Hoefer, Inc. immediately if the need for repair or servicing should arise.

8.2 Safety Precautions

- Should the sensor develop a fault or become disconnected the fan will remain “ON” so safeguarding against gel overheating.
- The fan utilizes a brushless motor with reverse polarity protection.
- A 500 mAmp 240 Volt antisurge fuse is located on the rear panel in a fuse holder immediately below the mains supply socket.

9. Ordering Information

Product	Number of wells	Thickness (mm)	Well width (mm)	Well depth (mm)	Well volume µl per 1 mm depth
Shark's Tooth Combs					
SQ33-C35-48K	48	0.35	3	10	5
SQ33-C35-60KMC	60*	0.35	2.5	10	4
SQ33-C35-96K	96	0.35	2.5	9	3.5

Standard Combs

SQ33-C35-48	48	0.35	4	5	7
SQ33-C35-60MC	60*	0.35	2	5	3.5

*Microtiter spacing for use with a multichannel pipette.

Additional Accessories and Replacement Parts	Catalog number
Glass Plates, Rectangular, 33 × 41 cm (2 pieces)	SQ33-GPLT
Glass Plates, Notched, 33 × 41 cm (2 pieces)	SQ33-NGPLT
Spacers, 0.35 mm thick (2 pieces)	SQ33-SP
Fan Sensor Kit, 115 VAC	SQ33-HS-115V
Fan Sensor Kit, 230 VAC	SQ33-HS-230V

Hoefer, Inc.

84 October Hill Road
Holliston, MA 01746

Toll Free: 1-800-227-4750

Phone: 1-508-893-8999

Fax: 1-508-893-0176

E-mail: support@hoeferinc.com

Web: www.hoeferinc.com

Hoefer is a registered trademark
of Hoefer, Inc.

© 2013 Hoefer, Inc. —
All rights reserved.

Printed in the USA.