

ABI PRISM[®] 3700 Instrument -Hardware-Raw Data Troubleshooting Module

This module is intended to troubleshoot raw data that are displaying anomalies in the <u>electropherogram display</u> in Data Collection Software V2.0 or the raw data profile in the analysis software.

BEFORE PERFORMING ANY TROUBLESHOOTING WORK ON YOUR 3700 INSTRUMENT, PLEASE READ THE ABI PRISM® 3700 DNA ANALYZER USER'S MANUAL FOR SAFETY AND WARRANTY INFORMATION AND FURTHER DETAILS ON USE OF THE SYSTEM.

Please contact technical support if you have any questions.

Text in blue indicates a link to another portion of the document.

Observations	Possible Causes	Recommended Actions
Two Lines per peak	The preferences to display the confidence bands were turned on	Change the preferences to deselect the display confidence bands.
6 0 0 0 0 0 0 0 0 0 0 0 0 0		
Click to enlarge or view more examples		
Pull Up and/or Pull Down Peaks	Incorrect spectral is selected	Rerun the samples verifying that the correct <u>spectral dye set</u> for your chemistry is being used.
Click to enlarge	Optics have been realigned or replaced Different polymer type is being used. Switched to POP5 or POP6	 Rerun the spectral calibration. A spectral will need to be run for the following changes: Changes to the optics. A different polymer type is being used. New dyes are being used on the system. For example in sequencing moving from BigDye V1.1 to BigDye V3.1 or in fragment analysis changing DS-30 to DS-33.

RawDataTroubleshooting_Final



Observations	Possible Causes	Recommended Actions
	Background fluorescence	Refer to the Array View Colors Troubleshooting Module to identify fluorescence.
	For fragment analysis, the capillary that is borrowed is further than 5 capillaries away.	Select a capillary that is within a 5 capillary distance. If a capillary is unavailable, run another spectral. Refer to page 4-41 of the ABI PRISM® 3700 DNA ANALYZER USER'S MANUAL, "Reviewing and Overriding the Spectral Calibration Profiles."
No Data/Weak data No data	If results are sporadic or random: Sample Transfer problems Bubbles in the sample wells	Spin down samples before placing on the instrument. Inspect the
	Bubbles in the sample transfer lines	Verify that the <u>fittings on the ports</u> above the sample transfer syringes are finger tight.
Operation Operation <thoperation< th=""> Operation <thoperation< th=""> Operation <thoperation< th=""> <thoperation< th=""> <thope< td=""><td></td><td>Verify that the <u>sample transfer syringes</u> are tightened to the pump. Finger tight should be sufficient.</td></thope<></thoperation<></thoperation<></thoperation<></thoperation<>		Verify that the <u>sample transfer syringes</u> are tightened to the pump. Finger tight should be sufficient.
Weak data ²⁰⁰ ¹⁰ ¹⁰⁰		Verify that the fluid line sinkers are at the bottom of the reservoir.
		Verify water levels in the water reservoir. Replenish if necessary.
		Run the <u>Change tips wizard</u> . This will prime the syringes and flush water through the tubing going to the loading tips. During the wizard, monitor both sample transfer syringes for bubbles. If bubbles are present in the syringe but not the tubing from the water reservoir, change the syringe.
	Clogged autoloader tips	Replace the <u>autoloader tips</u> . Refer to the <u>change tips wizard</u> for instructions.
	Sample Issues	Run the Long Read standard to determine if the issues sample related and then contact Technical Support.

 $RawDataTroubleshooting_Final$



Observations	Possible Causes	Recommended Actions
	If results are capillary specific Capillaries are clogged	Run a regenerate array to flush out of the capillary, refer to the <u>Change Array wizard</u> and select to clean the array. This may not always be successful. If the capillary continues to fail, one can block the capillary within the software, refer to page 4-12 of the ABI PRISM® 3700 DNA ANALYZER USER'S MANUAL, "manually overriding the capillary state." Or replace the array
	Capillaries are broken	Replace the <u>array</u> .
	Capillaries are not seated properly	Verify that the <u>array loading-end header</u> is seated correctly, not under tension, and the capillary tips are positioned in their injection wells.
	Inaccurate spatial calibration	Run another spatial calibration.
If results are across the entire an Sample Transfer pumps may be functioning	If results are across the entire array Sample Transfer pumps may not be functioning	Refer to the Verifying Syringe Pumps Module for more details.
	Electrophoresis issues	Refer to the Electrophoresis Module for more details.
If results are have a pattern	Refer to the Missing Lanes Module to for more details.	
Spikes or Bright Events A leak in the system allowing particulates to enter the system.		Inspect polymer lines for any source of bubbles or particles, correct any leaks. Replace the <u>inline filter</u> Replace polymer with a fresh bottle Perform a regenerate array
	Bubbles in the cuvette or the array	Run the <u>Cuvette Flush module</u> . Refer to the <u>Checking</u> for <u>Bubbles</u> Module and the <u>Checking</u> for <u>Leaks</u> Module if the bubbles in the cuvette persist.

 $RawDataTroubleshooting_Final$



Observations	Possible Causes	Recommended Actions
A Martin Martin Martin Carlos	Electrical noise	Use an Uninterruptible Power Supply (UPS).
<section-header></section-header>	Capillary specific	Run a regenerate array to flush out of the capillary, refer to the <u>Change Array wizard</u> and select to clean the array. This may not always be successful depending on the contaminant. If the spike does not go away, one can block the capillary within the software, refer to page 4-12 of the ABI PRISM® 3700 DNA ANALYZER USER'S MANUAL.
	Over tightening the <u>thumbscrews</u> on the detection end of the array. The laser beam may simply be too close to the capillaries resulting in direct illlumination of the tips	Remove the array and re-install such that the thumbscrews are tightened until finger tight. Check the CCD View to confirm the presence of horizontal rows of bright spots near the bottom of the CCD display window, typically referred to as "picket fence." Refer to the Checking CCD Module for more details on checking the CCD view.
Anomalous peaks	Bubbles migrating through the capillary degrade resolution and generate a bright signal upon emerging from the capillary tip within the cuvette.	Run the <u>Cuvette Flush module</u> . Refer to the Checking for Bubbles Module and the Checking for Leaks Module if the bubbles in the cuvette persist.
Click to enlarge Data preceding the spike may have lowered resolution		

Page 4 of 36



Observations	Possible Causes	Recommended Actions
Spikes underneath the baseline	Spatial calibration may have shifted	Run another spatial calibration and run Long Read Standard sample.
Jagged Peaks	Fluid other than polymer is mixing in the cuvette resulting in refractive index changes and deflection of the laser beam	Avoid mixing polymer from different bottles. Run the change polymer wizard for a different lot if the lot numbers are different. Inspect sheath-flow syringe for signs of fluid mixing during filling Install a new inline polymer filter and run the <u>Cuvette Flush module</u> to thoroughly flush cuvette
Symptoms are more evident for longer fragments Dark lines seen across array view	Microbial contamination of water and/or buffer reservoirs	Rinse the water and buffer reservoirs with either 1) hot water or 2) a 1:3 dilution of household bleach in water followed by a thorough rinsing with deionized water
Elevated/Separation of the Baseline	Background Fluorescence	Refer to the Array View Troubleshooting Module for more information.

 $RawDataTroubleshooting_Final$



Observations	Possible Causes	Recommended Actions
200 00 00 100 100 100 100 100 10	Contaminant present but electrophoreses out before data is collected	There may be a contaminant in the system which appears early in the run, before the samples. The software incorporates this contaminant as part of the background level. However, the contaminant then electrophoreses out of the system before the samples appear, therefore the amount of background changes. It may not affect the analyzed data depending on the fluorescent intensity of the contaminant. Refer to the Instrument Maintenance Module.
G or C dye breakdown	Chemical breakdown of BigDye™ Terminator G-dye or C-dye	Minimize exposure to air for samples re-suspended in HiDi Formamide by covering the samples with the recommended foil or a heat seal film and use the robot piercing. Formamide absorbs water in air undergoing hydrolysis to generate Formic Acid and Ammonia. One of these is partially responsible for the degradation of the extension fragments.
Peaks appear as doublets or have bulges or small adjacent extra peaks	Excessive sample has been electrokinetically injected into capillary	The recommended injection solution is HiDi Formamide because it is less prone to evaporation, provides denaturing capabilities and yields the most consistent results. However, water can also be used for resuspension of samples but due to extremely efficient injection it generates the most signal but can also lead to excessive sample being injected. If water is used it may be necessary to increase the resuspension volume. Run the Long Read Standard sample using HiDi Formamide to eliminate instrument issues.
Gradual loss of resolution	Changes in the state of the capillary walls can lead to degradation of performance.	Run a regenerate array to flush out the capillary every 80-100 runs, if maintained properly, the array should provide 300 runs. If the regenerate array procedure does not work, replace the capillary array.



Observations	Possible Causes	Recommended Actions
	Old or expired reagents	Replace the buffer or polymer with fresh supply. Verify that the reagents have not expired and have been on the instrument less than 7 days.
	Poor quality water in buffer solution or incorrectly made buffer	Remake 1X running buffer with fresh distilled, deionized water. Filtered water is adequate if the filtration system produces water purified to 18 Mega Ohms resistance. Use a different water source if available.
	Buffer or polymer levels are depleted.	Prepare fresh 1X running buffer and fill the reservoir. Place a new bottle of polymer on the system.
	Bubbles in the cuvette or capillaries are blocking current flow	Check for bubbles within the CCD view, refer to Checking the CCD window Module. Run the <u>cuvetteflush.mod</u> service module.
	Sample issues	Please contact technical support for sample related issues.
Some samples have data start points that are later than most of the samples of the run. Data in the late-starting lanes	"Overloading" of sample due to the high efficiency of injection from water.	Use HiDi Formamide or increase the resuspension volume. Contact technical support for more information on decreasing sample signal intensity.
suffers early loss of resolution		Run the Long Read Standard to verify instrument is running properly.
Click to enlarge	Old or expired reagents which can affect current	Replace the buffer or polymer with fresh supply. Verify that the reagents have not expired and have been on the instrument less than 7 days.
Low read length on one side Failing samples with signal > 50 RFU	Clumped capillaries	Inspect the <u>array</u> for any capillaries that appear to be clumped together. Carefully clean and separate clumped capillaries in place with damp tissue.

 $RawDataTroubleshooting_Final$

Page 7 of 36



Observations	Possible Causes	Recommended Actions
	Sample issues	Run the long read standard to verify the instrument is running properly, if so please contact technical support for sample related issues.
Failing samples with signal <50 RFU	Broken capillaries	Inspect the capillary <u>array</u> for broken capillaries; replace defective arrays.
	Improperly installed array load header	Confirm that the capillary <u>array load header</u> is properly seated and that all capillaries are in their proper wells.
	Non-optimal cuvette temperature	Choose the cuvette temperature with the most uniform signal strengths across the array for all subsequent run.

If problem persists, please contact technical support.

Contacting Technical SupportBy phone:1-800-831-6844, option 5By email:ABTechnicalSupport@appliedbiosystems.com



Restarting the 3700 DNA Analyzer

IMPORTANT To prevent firmware and software memory problems, we recommend that you restart the instrument and the software once a week.

To shut down and restart the instrument:

Step 1: The instrument should not be running or extracting data.

Step 2: Close the 3700 Data Collection software by selecting Shutdown from the File menu.

Note You cannot use the Close button to exit the software.

Step 3: Close the OrbixWeb Daemon software by right-clicking on its button in the taskbar and selecting Close from the pop-up

menu. If you get a run-time message, click **OK** to close the message.

IMPORTANT Do not shut down the OrbixWeb[™] Daemon until after you have shut

down the Data Collection program.

Step 4:

a. Turn off the instrument using the On/Off button.

The following diagram shows the left front grill.



b. Wait 30 seconds.

c. Turn on the instrument.

Step 5: When the green status light is steady, wait 1 minute.

Step 6: Restart the 3700 Data Collection software.

RawDataTroubleshooting_Final



How is the array view accessed in Data Collection V2.0?

During a run

Note Always exit the Array View when you are finished viewing. Do not leave the window open for extended periods during a run as this may cause unrecoverable screen update problems.

1) Select the **Run Status** tab > select **Array View**.

Below is a normal view of a run.

Array View Sub-Page



RawDataTroubleshooting_Final



After a run

Accesses information within the database

1) Select the **Data Acquisition** menu and select **Display Reduced Data for Run**.

2) This opens the **Select the run to display** dialog box.



3) From the drop-down list box, select the run that you want to display. Click **OK**. This opens the **Retrieving block** box.



4) When the retrieving window disappears go the to **run status** tab and select the **array view**.

RawDataTroubleshooting_Final

Page 11 of 36





RawDataTroubleshooting_Final

Page 12 of 36



What is the electropherogram display?

An electropherogram is a graph of relative dye concentration against time, plotted for each dye. The electropherogram display is plotted using Intensity vs. Time. The plot also appears in the capillary view tab but has been rotated 90 degrees, this shows capillary number vs. run time. This view provides a quick view of capillary specific issues.

To access this view select the Run Status tab > select the Capillary View tab



RawDataTroubleshooting_Final



Normal Data

Below is how normal raw data should look and a general overview of what the plot displays.

Example 1: Data from an empty well



• The Y axis represents the intensity. This scale will vary depending on sample and may be misleading if any offscale data exists, as the sample may appear weak so it is important to note the scale of "good" data. In this example the scale is 0 to 100.

• The X axis represents the time in minutes. Note that the peak morphology may differ depending on the part of the run that is being viewed.





RawDataTroubleshooting_Final

Page 14 of 36



Dye sets

Filter Set	Dye Set	Dyes or kit
E	DS-01	BigDye™ Terminator V1.1
Н		BigDye™ Terminator v3.1
D	DS-30	6-FAM, HEX, NED and ROX
D	DS-31	6-FAM, VIC, NED and ROX
F	DS-32	5-FAM, JOE, NED and ROX
G5	DS-33	6-FAM, VIC, NED, PET and LIZ
E5	DS-02	dR110, dR6G, dTAMRA, dROX and LIZ

 $RawDataTroubleshooting_Final$

Page 15 of 36



Two Lines per peak

Changing the preferences

- 1. In the Data Collection V2.0 Software select the Edit menu > select Preferences > select Run Status tab in the Setting Preferences window.
- 2. In the Data Acquisition Option, deselect the Display Confidence Bands

	Setting Preferences
	Users Instrument Plate Setup Run Setup Run Status Data Analysis Data Extraction
	Status Options
	Refresh Interval (Range 0-30 Seconds):
	2 Seconds Record Interval (Range 0-300 Seconds);
	5 Seconds.
	Data Acquisition Options
$\boldsymbol{\triangleleft}$	Display Confidence Bands
	OK Cancel Apply

The Display Confidence Bands allows for the evaluation of the quality of the spectral. The confidence bands provide a measure of the statistical confidence with which the electropherogram data is known.

A small confidence interval indicates a good spectral calibration for that capillary. A large confidence interval indicates a potential problem with the spectral data. For more details on this please refer to the ABI PRISM® 3700 DNA ANALYZER USER'S MANUAL.

RawDataTroubleshooting_Final

Page 16 of 36



Examples of two lines per peak



Example 3



Example 1 shows a good quality spectral.

Example 2 shows the optimal spectral quality as the peaks appear to be a single peak. Note that this sample is the Long Read Standard run immediately after the spectral calibration and that the spectral file is not borrowed from another capillary. **Example 3** shows a questionable spectral applied to a long read standard. Note the separation of the peaks and also the pull up and pull down of the dyes under the main peak. The spectral should be re-run.

.....ng_Final

Page 17 of 36



Pull Up and/or Pull Down Peaks



This exhibits both pull up and pull down in the baseline. Pull down in green can be seen under the yellow and blue peaks. Pull up in red can be seen under the yellow peaks

♦ If there are consistent patterns in the colors that are pulled up or pulled down for example, green always below blue, this can be attributed to an improper spectral being applied to the sample. The function of the spectral is to subtract any overlap between the dyes.

♦ If the pull up and pull down does not exhibit a consistent pattern then it may be sample related.

RawDataTroubleshooting_Final

Page 18 of 36



Weak Baseline or Flat Baseline

Sequencing standards typically give peaks above 800.

Examples:

No sample injection

The scale is only +/- 100 and the baseline remains near 0.



Weak injection

The scale is 0-200 but the peak heights are at 125 and below.



RawDataTroubleshooting_Final

Page 19 of 36



Spikes or Bright Events

The presence of contaminants or bubbles in the system can cause a brief intense signal when hit by the laser beam. Spikes in the raw data can also be due to voltage spikes from either 1)other instruments on a non-dedicated line or 2) from a non-UPS surge-protected circuit. The spike can be in all four colors or a specific color depending on the cause and contaminant.



The intensity of the spike is near 15,000. The sample is a long read standard which has an intensity of 800-1200 but due to the intensity of the spike the sample appears to be weak. It is important to note the scale of the intensity vs. time graph. Also, note that the spike generates a large negative peak, due to the intensity of the spike the spectral can not adequately compensate for the amount of signal observed. It is common for spikes to also exhibit negative peaks.

RawDataTroubleshooting_Final





The array view sub page. A bright event occurring in one capillary can affect neighboring capillaries. The intensity of the spike will be much less than the capillary of origin.

RawDataTroubleshooting_Final

Page 21 of 36



Anomalous Peaks



Page 22 of 36



Spikes underneath the baseline



 $RawDataTroubleshooting_Final$

Page 23 of 36



Jagged Peaks/Striations





Elevated/Separation of the Baseline



 $RawDataTroubleshooting_Final$

Page 24 of 36



Picket Fence





Module

The picket fence is seen only on the right side and is minimal. No action will need to be performed

The picket fence observed to the left appears across the entire array. If data is affected, it may be necessary to adjust the array at the detection end, tightening the thumbscrews until finger tight. If this does not resolve the issue, replace the array.

RawDataTroubleshooting_Final

Page 25 of 36



G or C Breakdown

Formamide absorbs water in air undergoing hydrolysis to generate Formic Acid and Ammonia. One of these is responsible for the degradation of the dyes. To prevent the occurrence of dye breakdown it is recommended to run the samples immediately after resuspension or to store the samples seal the samples to prevent exposure to air.





Double Peaks



RawDataTroubleshooting_Final

Page 27 of 36



Late Start



RawDataTroubleshooting_Final

Page 28 of 36

Sample Transfer syringe and pump The red arrows identify areas of possible leaks.

Autoloader tips

Tips should be tightened until finger tight.

RawDataTroubleshooting_Final

Module

Page 29 of 36

Array Diagram

Detection end of array

Check the area for any polymer leaks, if present remove the array and clean the area of polymer. Clean with a lint free tissue moistened with DI water. See below to check the o-ring.

Thumbscrews: These should be tightened until finger tight

RawDataTroubleshooting_Final

Page 30 of 36

Load Bar

Verify that the capillaries are sitting in the wells of the load bar and that they are covered with water or buffer.

Inline Filter

Dried polymer will appear in these areas.

Load bar

RawDataTroubleshooting_Final

Module

Page 31 of 36

Reservoirs and Fluid Line Sinkers

The Buffer and Water reservoirs should have enough fluids for the scheduled runs. It is recommended to check the reservoirs daily. Also, verify that the fluid line sinkers are at the bottom of the reservoirs and that the lines are submerged under the fluid.

Running the wizards

Step 1: Select Instrument > Select Wizards > Select the appropriate wizard.

Instrument		
Run	Ctrl+R	
Skip	Ctrl+S	
Stop	Ctrl+T	
Wizards	۰,	Change Array
Utilities	+	Change Syringe
Manual Co	ontrol	Change Polymer
		Change Tips
		Longterm

RawDataTroubleshooting_Final

Page 32 of 36

Running service modules within Data Collection V2.0

Page 33 of 36

Service Mod	ule Location	? ×
Look jn: a) AutoPlatel a) BubbleRei a) CuvetteFlu a) PipetTest3 a) PipetTest3 a) PlenumPu	Service Modules Service Modules Os (C:) AppliedBio Support Files Data Collection Support Programs (E:) (R:)	Anange.mod JseOnly.mod smod hange.mod
File <u>n</u> ame: Files of <u>t</u> ype:	All Files (*.*)	<u>O</u> pen ▼ Cancel

The service module utility should automatically search the Service Modules Folder, if not the folder can be found in the D://AppliedBio/Support Files/Data Collection Support Files/ Service Modules

Service Module Location			? ×	
Look in: 🔄	Service Modules	•		*
 AutoPlateDeckCalibration.mod BubbleRemove.mod CuvetteFlush.mod PipetTest384.mod PipetTest96.mod PlenumPurge.mod 		 PopToWaterExch PostBatchInternal PreBatchInternalU STRSyringePrime. WaterFlush.mod WaterToPopExch 	ange med 1 <mark>Up On</mark> Useoniy.mo IseOnly.mo mod ange.mod	e Level od
File <u>n</u> ame: Files of <u>type</u> : All	Files (*.*)		- [-	<u>O</u> pen Cancel

Double click on the module, in this example CuvetteFlush.mod.

RawDataTroubleshooting_Final

Page 34 of 36

After selecting run method, the steps of the service module will be listed.

Instrum	ent Command Freshonse.
You must	choose a file to run.You selected module file CuvetteFlush.mod
"Begin Cu	vette Flush Module."
"Assume:	s cuvette and lines are filled with polymer."
"Initializing	sensors."
"Flushing	tubing to cuvette with POP."
"Filling los	d trough with water"
"Flushing	tubing to cuvette with POP."
"Flush tw	ice with polymer pump while the sheath pump is flushing."
"Waiting f	or sheath pump."
"Flush po	ymer through cuvette, out valve 11."
"Waiting f	or sheath pump."
"Flush 2.5	ml out valve 11 with sheath pump."
"Flush 4rr	I polymer out valve 11 with polymer pump while sheath pump is running.
"Finished	3ml."
"Waiting f	or sheath pump."
"Allow pr	essure to return to zero."
"Flushing	bubbles from cuvette"
"Pressuri:	ze to 60 psi for 1 minute."
"Opening	cuvette to purge bubbles."
"Flush po	ymer through cuvette with polymer pump, out valve 11, then out valve 12
"Repeat f	usning pupples from cuvette. " to to 50 voi for 4 minuto "
"Pressuri. "Eluch po	ze to oo psi tor n minute." Iuraar through currette with polymer purph, out yolyo 11, then out yolyo 11
"Third tim	ynier through cuvette with polynier pump, out valve 11, then out valve 12 a fluebing bubbles from cuvette. "
"Drossuri	ze to 60 psi for 1 minute "
"Flush no	vero oo psi toi i minute. Wher through cuvette with polymer pump, out value 11, then out value 11
"Flush 1 4	ind notware out valve 12 with sheath numn "
"Flush 2m	I polymer out valve 12 with polymer pump while sheath pump is running '
"Filling los	d trough with water"
a	
	Run Method Select Module

The module selected should be listed.

Run Service Modules			×
Instrument Command / Response.			
You selected module file CuvetteFlush	nod		*
			0
Run Method	Select Module.	-	

 $RawDataTroubleshooting_Final$

🛃 Information					
•	D:\AppliedBio\Support Files\Data Collection Support files\Service Modules\\CuvetteFlush.mod Method Execution Has Successfully Finished.				

When this window appears the module has completed, if there is an error during the module, please contact technical support.

RawDataTroubleshooting_Final

Page 36 of 36