

STUDENT WORKBOOK

QNET HVAC Trainer for NI ELVIS

Developed by Quanser

Curriculum designed by:
Karl Johan Åström, Ph.D., Lund University (Emeritus)
Jacob Apkarian, Ph.D., Quanser
Paul Karam, B.A.SC., Quanser
Michel Lévis, M.A.Sc., Quanser
Jeannie Falcon, Ph.D., National Instruments

Curriculum
complies with*



© 2011 Quanser Inc., All rights reserved.

Quanser Inc.
119 Spy Court
Markham, Ontario
L3R 5H6
Canada
info@quanser.com
Phone: 1-905-940-3575
Fax: 1-905-940-3576

Printed in Markham, Ontario.

For more information on the solutions Quanser Inc. offers, please visit the web site at:
<http://www.quanser.com>

This document and the software described in it are provided subject to a license agreement. Neither the software nor this document may be used or copied except as specified under the terms of that license agreement. All rights are reserved and no part may be reproduced, stored in a retrieval system or transmitted in any form or by any means, electronic, mechanical, photocopying, recording, or otherwise, without the prior written permission of Quanser Inc.

Acknowledgements

Quanser, Inc. would like to thank the following contributors:

Dr. Hakan Gurocak, Washington State University Vancouver, USA, for his help to include embedded outcomes assessment, and

Dr. K. J. Åström, Lund University, Lund, Sweden for his immense contributions to the curriculum content.

Contents

1	Introduction	4
2	On-Off Control	5
2.1	Background	5
2.2	On-Off Control Virtual Instrument	6
2.3	Startup [10 min]	7
2.4	Lab 1: Relay Control [30 min]	7
2.5	Lab 2: Modeling [30 min]	8
3	PI Control	9
3.1	Background	9
3.2	Speed Control Virtual Instrument	10
3.3	Startup [10 min]	10
3.4	Lab 1: Qualitative PI Control [30 min]	11
3.5	Lab 2: Saturation and Windup [30 min]	11
3.6	Lab 3: Set-Point Weight [20 min]	12
3.7	Lab 4: PI Control according to Specifications [30 min]	12
4	System Requirements	14
4.1	Overview of Files	14
4.2	On-Off Control VI	15
4.3	PI Control VI	17
5	Lab Report	19
5.1	Template for Content (On-Off Control)	19
5.2	Template for Content (PI Control)	20
5.3	Tips for Report Format	21

1 INTRODUCTION

The QNET-012: heating and ventilation trainer (HVACT) is shown in Figure 1.1. The system consists of a plexiglass duct, with a heater in one end and a blower in the other end. The heater is a halogen lamp and the blower is a variable-speed fan. There is a thermistor sensor placed inside the duct to measure the temperature of the chamber and another thermistor sensor outside the chamber to measure the room temperature.




Figure 1.1: QNET heating and ventilation trainer (HVACT)

The temperature measured at the thermistor inside the chamber is to be controlled using the heater voltage while the fan is ran at a constant speed. Heat is transferred to the thermistor by radiation from the heater and by convection from the air stream. Radiative heat transfer is highly nonlinear and it is therefore difficult to model the system by first principles. As a result, empirical tuning will be used to control the system. This heat transfer plant is very similar to the systems that are used to control wafer temperature in semiconductor manufacturing.

There are two experiments: on-off control and PI control. The experiments can be performed independently.

Topics Covered

- On-off control
- Modeling
- PI control
- Integrator windup
- Set-point weight

Prerequisites

In order to successfully carry out this laboratory, the user should be familiar with the following:

- Transfer function fundamentals, e.g., obtaining a transfer function from a differential equation.
- Using **LabVIEW®** to run VIs.

2 ON-OFF CONTROL

2.1 Background

2.1.1 Relay Control

On-off control or relay feedback is one of the simplest control strategies. The heater is switched on when the temperature is lower than the desired value, and the heater is switched off when the temperature is higher than the desired value. To avoid rapid switches it is common to introduce a hysteresis in the relay switch. A block diagram of a system with relay feedback is shown in Figure 2.1.




Figure 2.1: Block diagram of the heater system with relay feedback

The error, variable e in Figure 2.1, is the difference between the reference temperature, T_r , and the actual chamber temperature, T_c . The on-off controller is implemented using a relay switch with hysteresis, as shown in Figure 2.2. The heater actuator is represented by a saturation block and the chamber plant is represented by the transfer function $P(s)$.




Figure 2.2: Input and output relation for an on-off controller with hysteresis

The hysteresis width, ΔT_h in Figure 2.2, has to be chosen such that a large measurement noise does not generate

any unintentional switches. As depicted in Figure 2.2, the output control signal voltage of the on-off controller can be adjusted using a mean or offset, $V_{h,off}$, and an amplitude, $V_{h,amp}$.

In the experiment, the behavior of the heater system will be investigated for different values of controller parameters. More specifically, the control signal and the measured temperature will be observed.

2.1.2 Modeling

The on-off control input and the measured temperature output from the experiment have an interesting property that makes it possible to find a simple model for the process. The temperature response is a ramp due a voltage step therefore the temperature is the integral of the voltage. Under certain conditions, the process can be modeled by the simple transfer function

$$P(s) = \frac{K_v}{s} \quad (2.1)$$

where the parameter K_v is the slope of the ramp.

2.2 On-Off Control Virtual Instrument

Tracking a reference temperature using a relay control is first examined in this laboratory. Then, when commanding a fixed reference, the system can be modeled. The LabVIEW virtual instrument for on-off temperature control is shown in Figure 2.3.




Figure 2.3: Virtual instrument for on-off heater control

See Wikipedia for more information on [relay](#), [hysteresis](#), [mathematical model](#), [transfer function](#), and [LTI system theory](#).

2.3 Startup [10 min]

After powering up the ELVIS and before running any of the labs, follow this procedure to calibrate the QNET HVAC system.

1. Open the QNET_HVACT_On_Off_Control.vi is open and configured as described in Section 4.2.
2. Make sure the correct *Device* is chosen.
3. Run the QNET_HVACT_On_Off_Control.vi shown in Figure 2.4.
4. The cooling fan is automatically activated when the *Prototyping Board Power* switch on the ELVIS unit is on. Let the actual temperature, T_c , in the *Temperature (C)* scope settle until it stops decreasing.
5. Adjust the *Temperature (C)* scope scales to see both the reference and actual temperatures (see the QNET HVAC User Manual ([1]) for help).
6. As illustrated in Figure 2.4, calibrate the temperature sensors by clicking on the *Calibrate* button. This will align the chamber temperature, T_c , to the measured ambient temperature, T_a .
7. Activate the control by clicking on the *Heater OFF* button (in the top-right corner of Figure 2.4).




Figure 2.4: Calibrating the temperature in the QNET HVAC On-Off Control VI.

2.4 Lab 1: Relay Control [30 min]

1. Make sure the QNET_HVACT_On_Off_Control.vi is running and has been calibrated as instructed in Section 2.3. When running, the VI should look similar to Figure 2.3.
2. In the *Signal Generator* section set:

- Amplitude = 0 °C
 - Frequency = 0.008 Hz
 - Offset = 0.5 °C
3. Examine the actual temperature (red) and reference temperature (blue) responses in the *Temperature (C)* scope.
 4. Gradually vary the *Offset* in the *Signal Generator* between 0.5 °C and 2 °C. How is the reference temperature, T_r , in the *Temperature (C)* scope is set? Attach a sample temperature response.
 5. Vary the relay amplitude, Vh_amp , in the *Control Parameters* section. Explain how the heater voltage affects the temperature variation and, in particular, observe the frequency and amplitude of the chamber temperature. Attach a representative temperature response.
 6. Explain the effect of changing the relay mean, Vh_off . Attach a temperature response.
 7. Examine the effects of changing the relay width (or hysteresis), DTh , between 0.01 °C and 1.00 °C. Give a short explanation and attach a temperature response with a narrow and wide hysteresis.
 8. Click on the *Stop* button to stop running the VI.

2.5 Lab 2: Modeling [30 min]

1. Make sure the QNET_HVACT_On_Off_Control.vi is running and has been calibrated as instructed in Section 2.3. When running, the VI should look similar to Figure 2.3.
2. In the *Signal Generator* section set:
 - Amplitude = 0 °C
 - Frequency = 0.008 Hz
 - Offset = 0 °C
3. In the *Control Parameters* section set:
 - Vh_amp = 4 V
 - Vh_off = 4 V
 - DTh = 0.50 °C
4. Adjust the *Temperature (C)* scope scales to see both the reference and actual temperatures (see the QNET HVACT User Manual ([1]) for help).
5. Adjust the *Offset* in the *Signal Generator* to obtain a relatively symmetrical oscillation (i.e., the rate of increase and decrease should be similar).
6. Observe the heater voltage and the chamber temperature. As discussed in Section 2.1, this can be modeled by the simple transfer function $P(s) = K_v/s$. Find parameter K_v that would describe the relation between the voltage and the temperature signals. Attach both the temperature and voltage responses used to find K_v and show your calculations.
7. Click on the *Stop* button to stop running the VI.

3 PI CONTROL

3.1 Background

The oscillations that occur with on-off control can be avoided by using a linear proportional and integrating controller. To design such a controller analytically a simple model representing the actual plant is needed. Since the conditions shown in Figure 3.2 are representative for what happens when the temperature is controlled, transfer function 2.1 can be used for the model-based approach to find the controller. The block diagram of the closed loop system is shown in Figure 3.1.




Figure 3.1: Block diagram of heater PI closed-loop system.m

The process transfer function is the transfer function in Equation 2.1 and the input-output relation for a PI controller with set-point weighting is

$$U(s) = k_p (b_{sp} R(s) - Y(s)) + \frac{k_i (R(s) - Y(s))}{s} \quad (3.1)$$

The closed loop transfer function from the relative temperature reference, $\Delta T_r = T_r - T_a$, to the output temperature measured relative to the ambient temperature, $\Delta T_c = T_c - T_a$, is

$$G_{\Delta T_c, \Delta T_r}(s) = K_v \frac{k_p b_{sp} s + k_i}{s^2 + K_v k_p s + K_v k_i} \quad (3.2)$$

The closed-loop system has the characteristic polynomial

$$s^2 + K_v k_p s + K_v k_i \quad (3.3)$$

and the desired closed loop characteristic polynomial is

$$s^2 + 2\zeta\omega_0 s + \omega_0^2, \quad (3.4)$$

where ω_0 is the undamped closed loop frequency and ζ is the damping ratio. The characteristic equation in 3.3 matches equation 3.4 with the proportional control parameter

$$k_p = \frac{2\zeta\omega_0}{K_v} \quad (3.5)$$

and the integral control gain

$$k_i = \frac{\omega_0^2}{K_v} \quad (3.6)$$

Large values of ω_0 give large values of controller gain. This implies noise will create large variations in the control signal. The set-point weight parameter b_{sp} can be used to adjust the overshoot of the response.

The sensor signal is noisy and it is therefore necessary to filter the measured signal. A simple first-order filter has the transfer function

$$T_c = \frac{T_{c,meas}}{T_f s + 1}$$

where $T_{c,meas}$ is the measured temperature from the thermistor and T_f is the transfer function time constant. Increasing T_f decreases the cutoff frequency and minimizes noise in the signal at the expense of changing the shape of the signal.

Temperature control typically admits high controller gains. A consequence of this is that the controller output may saturate and result in integrator windup. The heater is therefore useful to illustrate the usefulness of integrator feedback.

3.2 Speed Control Virtual Instrument

The LabVIEW virtual instrument that implements the heater PI control is shown in Figure 3.2. The control parameters k_p , k_i , b_{sp} , the anti-windup tracking time constant T_r , and the filter time constant T_f can all be adjusted.




Figure 3.2: Virtual instrument PI control for heater

See Wikipedia for more information on [process control](#), [control theory](#), and [PID](#).

3.3 Startup [10 min]

After powering up the ELVIS and before running any of the labs, follow this procedure to calibrate the QNET HVACT system.

1. Make sure the correct *Device* is chosen.
2. Run the QNET_HVACT_PI_Control.vi shown in Figure 3.2.
3. The cooling fan is automatically activated when the *Prototyping Board Power* switch on the ELVIS unit is on. Let the actual temperature, T_c , in the *Temperature (C)* scope settle until it stops decreasing.
4. Adjust the *Temperature (C)* scope scales to see both the reference and actual temperatures (see the QNET HVACT User Manual ([1]) for help).

5. As illustrated in Figure 3.2, calibrate the temperature sensors by clicking on the *Calibrate* button. This will align the chamber temperature, T_c , to the measured ambient temperature, T_a .
6. Activate the control by clicking on the *Heater OFF* button (in the top-right corner of Figure 2.4).
7. Adjust the *Temperature (C)* scope scales to see both the reference and actual temperatures (see the HVACT User Manual ([1]) for help).

3.4 Lab 1: Qualitative PI Control [30 min]

1. Ensure the QNET_HVACT_PI_Control.vi is running and has been configured as described in Section 3.3.
2. In the *Signal Generator* section set:
 - Amplitude = $0.5\text{ }^{\circ}\text{C}$
 - Frequency = 0.02 Hz
 - Offset = $1.5\text{ }^{\circ}\text{C}$
3. In the *Control Parameters* section set:
 - $k_p = 4\text{ V}/^{\circ}\text{C}$
 - $k_i = 0.5\text{ V}/(^{\circ}\text{C}\cdot\text{s})$
 - $b_{sp} = 1$
 - $T_r = 1\text{ s}$
4. Examine the temperature response to the square wave input.
5. Set k_i to $0\text{ V}/(^{\circ}\text{C}\cdot\text{s})$ and change the proportional gain k_p between $2\text{ V}/^{\circ}\text{C}$ and $10\text{ V}/^{\circ}\text{C}$. Explain the effect proportional gain has on the temperature control performance. Attach a temperature response when using a low and high proportional gain.
6. Set k_p to $0.5\text{ V}/(^{\circ}\text{C}\cdot\text{s})$ and change the integral gain k_i between $0.25\text{ V}/(^{\circ}\text{C}\cdot\text{s})$ and $2.0\text{ V}/(^{\circ}\text{C}\cdot\text{s})$ and observe its effect on the temperature control performance. Show the temperature response with a low and high integral gain.
7. Click on the *Stop* button to stop running the VI.

3.5 Lab 2: Saturation and Windup [30 min]

1. Ensure the QNET_HVACT_PI_Control.vi is running and has been configured as described in Section 3.3.
2. In the *Signal Generator* section set:
 - Amplitude = $0.75\text{ }^{\circ}\text{C}$
 - Frequency = 0.02 Hz
 - Offset = $1.5\text{ }^{\circ}\text{C}$
3. In the *Control Parameters* section set:
 - $k_p = 8\text{ V}/^{\circ}\text{C}$
 - $k_i = 4\text{ V}/(^{\circ}\text{C}\cdot\text{s})$
 - $b_{sp} = 1$
 - $T_r = 100\text{ s}$
4. What effect does increasing the anti-windup reset parameter have on the control signal and on the temperature response? Attach a response of the temperature and heater voltage. See the *QNET Control Guide* ([2]) for more information on anti-windup.

5. In the *Control Parameters* section, set $Tr = 1.0$ s.
6. What effect does decreasing Tr have on the control signal and on the temperature response? Capture the temperature response as well as the heater voltage.
7. Click on the *Stop* button to stop running the VI.

3.6 Lab 3: Set-Point Weight [20 min]

1. Ensure the QNET_HVACT_PI_Control.vi is running and has been configured as described in Section 3.3.
2. In the *Signal Generator* section set:
 - Amplitude = $0.5\text{ }^{\circ}\text{C}$
 - Frequency = 0.02 Hz
 - Offset = $1.5\text{ }^{\circ}\text{C}$
3. In the *Control Parameters* section set:
 - $k_p = 8\text{ V}/^{\circ}\text{C}$
 - $k_i = 1\text{ V}/(^{\circ}\text{C}\cdot\text{s})$
 - $bsp = 0$
 - $Tr = 1\text{ s}$
4. Examine the response of the measured temperature in the *Temperature (C)* scope as well as the input heater voltage in the *Voltage (V)* scope. Attach the temperature and heater voltage responses.
5. Try the controller with a set-point weight of 1.
6. Study what effects raising bsp has on the measured temperature signal in the *Temperature (C)* scope and the control signal shown in the *Voltage (V)* scope. Capture the temperature response and its corresponding heater voltage.
7. Click on the *Stop* button to stop running the VI.

3.7 Lab 4: PI Control according to Specifications [30 min]

3.7.1 Pre-Lab Questions

1. Find the proportional and integral gains, i.e., k_p and k_i , needed for the response to satisfy the following specifications:
 - $\zeta = 0.60$
 - $\omega_0 = 0.125\text{ rad/s}$

Use the model gain K_v found previously in Section 2.5 and the design principles outlined in Section 3.1.

3.7.2 In-Lab Exercises

1. Ensure the QNET_HVACT_PI_Control.vi is running and has been configured as described in Section 3.3.
2. In the *Signal Generator* section set:
 - Amplitude = $0.5\text{ }^{\circ}\text{C}$

- Frequency = 0.02 Hz
 - Offset = 1.5 °C
3. Enter the control gains obtained in the pre-lab in the *Control Parameters* section of the VI.
 4. Examine the measured temperature response using your design PI gains. How is the performance of the controller compared to the previous controller? Attach the temperature and the heater voltage responses.
 5. Click on the *Stop* button to stop running the VI.

4 SYSTEM REQUIREMENTS

Required Hardware

- NI ELVIS II
- Quanser QNET Heating and Ventilation Trainer (HVACT). See QNET HVACT User Manual ([1]).

Required Software

- NI LabVIEW® 2011 or later
- NI DAQmx 9.3.5 or later
- NI LabVIEW Control Design and Simulation Module 2011 or later
- *ELVIS II Users*: ELVISmx 4.3 or later (installed from ELVIS II CD).

■ **Caution:** If these are not all installed then the VI will not be able to run! Please make sure all the software and hardware components are installed. If an issue arises, then see the troubleshooting section in the QNET HVACT User Manual ([1]).

4.1 Overview of Files

File Name	Description
QNET HVACT User Manual.pdf	This manual describes the hardware of the QNET Heating and Ventilation Trainer system and how to setup the system on the ELVIS.
QNET HVACT Workbook (Student).pdf	This laboratory guide contains pre-lab questions and lab experiments demonstrating how to design and implement controllers on the QNET HVACT system LabVIEW®.
QNET_HVACT_On_Off_Control.vi	Control temperature using on-off control.
QNET_HVACT_PI_Control.vi	Control temperature using a proportional-integral (PI) regulator.

Table 1: Files supplied with the QNET HVACT Laboratory.

4.2 On-Off Control VI

The HVAC On-Off Control VI implements a relay to control the temperature of the chamber. This VI can also be used to model the dynamics between the heater voltage and the temperature. Table 2 lists and describes the main elements of the QNET-HVACT On-Off Control virtual instrument user interface. Every element is uniquely identified through an ID number and located in Figure 4.1.




Figure 4.1: QNET HVAC On-Off Control VI components.

ID #	Label	Symbol	Description	Unit
1	Chamber Temp	T_c	Temperature inside chamber numeric display.	°C
2	Ambient Temp	$T_{a,m}$	Temperature outside chamber numeric display (i.e. measured room temperature).	°C
3		T_a	Latched ambient temperature that is added to reference temperature from Signal Generator	°C
4	Calibrate		Sets the red latched ambient temperature to the measured ambient temperature.	
5	Signal Type		Type of signal generated for the input voltage signal.	
6	Amplitude		Generated temperature reference signal amplitude input box.	°C
7	Frequency		Generated temperature reference signal frequency input box.	Hz
8	Offset		Generated temperature reference signal offset input box.	°C
9	Vh_amp	$V_{h,amp}$	Heater voltage relay amplitude input box.	V
10	Vh_off	$V_{h,off}$	Heater voltage relay offset input box.	V
11	ΔT_h	ΔT_h	Heater relay hysteresis width.	°C
12	Temperature (C)	T_c, T_r	Scope with reference temperature (in blue) and measured chamber temperature (in red).	°C
13	Voltage (V)	V_h	Scope with applied heater voltage (in red).	V
14	Device		Selects the NI DAQ device.	
15	Sampling Rate		Sets the sampling rate of the VI.	Hz
16	Stop		Stops the LabVIEW VI from running.	
17	Heater OFF		Enables heater when pressed in.	

Table 2: QNET HVAC On-Off Control VI Components

Important: The reference temperature is relative to the latched ambient temperature, ID #3 in Table 2. The reference temperature is equal to the sum of the signal generated from the *Signal Generator* and the latched ambient temperature.

4.3 PI Control VI

In the QNET HVAC PI Control VI, a proportional-integral compensator is used to control the temperature of the chamber. The PI control includes anti-windup and set-point weight strategies. Table 3 lists and describes the main elements of the QNET-HVACT PI Control virtual instrument user interface. Every element is uniquely identified through an ID number and located in Figure 4.2.




Figure 4.2: QNET HVAC PI Control VI components

ID #	Label	Symbol	Description	Unit
1	Chamber Temp	T_c	Temperature inside chamber numeric display.	°C
2	Ambient Temp	$T_{a,m}$	Temperature outside chamber numeric display (i.e. measured room temperature).	°C
3		T_a	Latched ambient temperature that is added to reference temperature from Signal Generator	°C
4	Calibrate		Sets the red latched ambient temperature to the measured ambient temperature.	
5	Signal Type		Type of signal generated for the input voltage signal.	
6	Amplitude		Generated temperature reference signal amplitude input box.	°C
7	Frequency		Generated temperature reference signal frequency input box.	Hz
8	Offset		Generated temperature reference signal offset input box.	°C
9	kp	k_p	Controller proportional gain input box.	V/°C
10	ki	k_i	Controller integral gain input box.	V/(°C s)
11	bsp	b_{sp}	Controller set-point weight input box.	
12	Tr	T_r	Anti-windup tracking time constant.	s
13	Temperature (C)	T_c, T_r	Scope with reference temperature (in blue) and measured chamber temperature (in red).	°C
14	Voltage (V)	V_h	Scope with applied heater voltage (in red).	V
15	Device		Selects the NI DAQ device.	
16	Sampling Rate		Sets the sampling rate of the VI.	Hz
17	Stop		Stops the LabVIEW VI from running.	
18	Heater OFF		Enables heater when pressed in.	

Table 3: QNET HVAC PI Control VI Components

5 LAB REPORT

This laboratory contains three groups of experiments, namely,

1. On-off control, and
2. PI control.

For each experiment, follow the outline corresponding to that experiment to build the *content* of your report. Also, in Section 5.3 you can find some basic tips for the *format* of your report.

5.1 Template for Content (On-Off Control)

I. PROCEDURE

1. Relay Control

- Briefly describe the main goal of the experiment.
- Briefly describe the experiment procedure in Step 4 in Section 2.4.
- Effect of varying relay amplitude on response in Step 5 in Section 2.4.
- Effect of varying relay offset on response in Step 6 in Section 2.4.
- Effect of changing hysteresis width on response in Step 7 in Section 2.4.

2. Modeling

- Briefly describe the main goal of the experiment.
- Briefly describe the modeling procedure in Step 6 in Section 2.5.

II. RESULTS

Do not interpret or analyze the data in this section. Just provide the results.

1. Temperature response from Step 4 in Section 2.4.
2. Temperature response when varying relay amplitude in Step 5 in Section 2.4.
3. Temperature response when varying relay offset in Step 6 in Section 2.4.
4. Temperature response when varying hysteresis width in Step 7 in Section 2.4.
5. Response used for modeling system in Step 6 in Section 2.5.

III. ANALYSIS

Provide details of your calculations (methods used) for analysis for each of the following:

1. Effect of changing signal generator offset in Step 4 in Section 2.4.

IV. CONCLUSIONS

Interpret your results to arrive at logical conclusions for the following:

1. Finding the model gain parameter in Step 6 in Section 2.5.

5.2 Template for Content (PI Control)

I. PROCEDURE

1. *Qualitative PI Control*

- Briefly describe the main goal of the experiment.
- Briefly describe the experimental procedure in Step 5 in Section 3.4.
- Effect of changing proportional gain in Step 5 in Section 3.4.
- Effect of changing integral gain in Step 6 in Section 3.4.

2. *Set-Point Weight*

- Briefly describe the main goal of this experiment.
- Briefly describe the experimental procedure in Step 4 in Section 3.5.
- Effect of changing set-point weight in Step 6 in Section 3.6.

3. *Saturation and Windup*

- Briefly describe the main goal of this experiment.
- Briefly describe the experimental procedure in Step 4 in Section 3.6.

4. *PI Control According to Specifications*

- Briefly describe the main goal of the experiment.
- Briefly describe the experimental procedure in Step 4 in Section 3.7.

II. RESULTS

Do not interpret or analyze the data in this section. Just provide the results.

1. Low and high proportional gain temperature responses in Step 5 in Section 3.4.
2. Low and high integral gain temperature responses in Step 6 in Section 3.4.
3. High reset time temperature and heater voltage response in Step 4 in Section 3.5.
4. Low reset time temperature and heater voltage response in Step 6 in Section 3.5.
5. Temperature and heater voltage response with set-point weight of 0 in Step 4 in Section 3.6.
6. Temperature and heater voltage response with set-point weight of 1 in Step 6 in Section 3.6.
7. Temperature and heater voltage response with designed PI gains in Step 4 in Section 3.7.

III. ANALYSIS

Provide details of your calculations (methods used) for analysis for each of the following:

1. Effect of increasing reset time parameter in Step 4 in Section 3.5.
2. Effect of decreasing reset time parameter in Step 6 in Section 3.5.

IV. CONCLUSIONS

Interpret your results to arrive at logical conclusions for the following:

1. Performance of designed PI control response in Step 4 in Section 3.7.

5.3 Tips for Report Format

PROFESSIONAL APPEARANCE

- Has cover page with all necessary details (title, course, student name(s), etc.)
- Each of the required sections is completed (Procedure, Results, Analysis and Conclusions).
- Typed.
- All grammar/spelling correct.
- Report layout is neat.
- Does not exceed specified maximum page limit, if any.
- Pages are numbered.
- Equations are consecutively numbered.
- Figures are numbered, axes have labels, each figure has a descriptive caption.
- Tables are numbered, they include labels, each table has a descriptive caption.
- Data are presented in a useful format (graphs, numerical, table, charts, diagrams).
- No hand drawn sketches/diagrams.
- References are cited using correct format.

REFERENCES


- [1] Quanser Inc. *QNET Heating-Ventillation Control Trainer User Manual*, 2011.
- [2] Quanser Inc. *QNET Practical Control Guide*, 2011.

Six QNET Trainers to teach introductory controls using NI ELVIS

- **QNET DC Motor Control Trainer**
teaches fundamentals of DC motor control




- **QNET HVAC Trainer**
teaches temperature (process) control




- **QNET Mechatronic Sensors Trainer**
teaches functions of 10 different sensors



- **QNET Rotary Inverted Pendulum Trainer**
teaches classic pendulum control experiment



- **QNET Myoelectric Trainer**
teaches control using principles of electromyography (EMG)



- **QNET VTOL Trainer**
teaches basic flight dynamics and control



Quanser QNET Trainers are plug-in boards for NI ELVIS to teach introductory controls in undergraduate labs. Together they deliver added choice and cost-effective teaching solutions to engineering educators. All six QNET Trainers are offered with comprehensive, ABET*-aligned course materials that have been developed to enhance the student learning experience.

To request a demonstration or quote, please email info@ni.com.

*ABET Inc., is the recognized accreditor for college and university programs in applied science, computing, engineering, and technology. Among the most respected accreditation organizations in the U.S., ABET has provided leadership and quality assurance in higher education for over 75 years.

©2013 Quanser Inc. All rights reserved. LabVIEW™ is a trademark of National Instruments.