

PRINTED COPIES UNCONTROLLED UNLESS MARKED "CONTROLLED" IN RED

STANDARD OPERATING PROCEDURE FOR PFAFF BILIMETER III BILIRUBINOMETER AND ASSOCIATED CENTRIFUGE

Q-PULSE FILE NAME	LI-BIO-Bil3SOPUr
EDITION No	2.0
DATE OF ISSUE	18 th September 2012
REVIEW INTERVAL	2 Years
AUTHORISED BY	N A Smith
SIGNATURE	
AUTHOR	N A Smith
COPY	2 of 3
LOCATION OF COPIES	1. Master file, Quality Officer Shelves 2. Point of Care Testing Shelves 3. Intranet / Point of Care Services Electronic archives

Document review history		
Review date	Reviewed by	Signature
June 2009	N Smith	
June 2010	N Smith	
July 2011	N Smith	
July 2012	N Smith	

PRINTED COPIES UNCONTROLLED UNLESS MARKED "CONTROLLED" IN RED

DOCUMENT AMENDMENT FORM FOR STANDARD OPERATING PROCEDURE FOR UNISTAT BILIRUBINOMETER AND ASSOCIATED CENTRIFUGE				
Record Filename		LI-BIO-Bil3SOPUr		
Number	Date	Page Number	Amendment	Authorised by
1	July 12	All	Updated in line with new LI-BIO-SOP4POCT template; changes throughout	N A Smith
2	July 12	5-6	Updated responsibilities	N A Smith
3	July 12	7	Training and Competency Section 1.1 expanded	N A Smith
4	July 12	13	Included contact details for service engineer	N A Smith
5	July 12	14	Included Sections 3.7, 3.8 & 3.9 Audit, Comparability and Uncertainty	N A Smith
6	July 12	15	Included self assessment competency form	N A Smith
7				
8				
9				
10				
The amendment must be authorised by the POCT Manager to ensure all copies including the electronic version are updated simultaneously.				
Any minor amendment must be handwritten without obscuring the previous text. An asterisk is written in the adjacent margin to highlight the alteration and is signed and dated by the POCT Manager.				
Amendments requiring immediate action must be handwritten without obscuring the previous text. An asterisk is written in the adjacent margin to highlight the alteration and is signed and dated by the POCT Manager.				
Ten or less minor amendments may be recorded before a new edition is issued.				
Major changes must result in the immediate review of the procedure.				
Document amendment does not replace the review process.				

PRINTED COPIES UNCONTROLLED UNLESS MARKED "CONTROLLED" IN RED

HAZARDS AND PRECAUTIONS

For appropriate use of protective clothing, please refer to 'Use of Personal Protective Equipment (PPE) in Pathology and Phlebotomy'.

Substance	Hazards	Precaution	First aid Code*
IQC/EQA Specimens	Danger of infection	Follow standard precautions (see above). Dispose of in orange clinical waste bags.	A1, B, C
Blood	Danger of infection	Follow standard precautions (see above). Dispose of in orange clinical waste bags.	A1, B, C
Used glass capillary tube	Danger of infection	Follow standard precautions (see above). Dispose of in sharps bin	A1, B, C
Blood stained gauze etc	Danger of infection	Follow standard precautions (see above). Dispose of in orange clinical waste bags.	A1, B, C
PDI Sanicloth 70	Irritating to eyes (R36) Highly flammable (R11)	Keep away from sources of ignition. No smoking.	A, B, C

IF IN ANY DOUBT CONSULT A SENIOR MEMBER OF STAFF OR YOUR SAFETY REPRESENTATIVE

Key:

- *A1 Ingestion: wash mouth thoroughly with water and give plenty to drink. In severe cases obtain medical attention.
- B Eye contact: irrigate thoroughly with water. Seek medical help.
- C Skin contact: wash off skin thoroughly with water.
- D Inhalation: remove to fresh air. If severe call a physician.

If First Aid treatment has to be given, contact the nearest first aider but do not delay treatment to the casualty.

If necessary, call the crash team 2222.

Treat all body tissue and waste as potentially infective.

In all cases, fill in an Incident form and report to Occupational Health or Accident Service if out of hours.

**See Also Risk Assessments on Intranet - Departments/Pathology/POCT
Using a Bilirubinometer 98
Blood or Body Fluid Analysis 101**

TABLE OF CONTENTS

0 INTRODUCTION.....	5
0.1 Purpose and Scope	5
0.2 Clinical Indication	5
0.3 Responsibilities	5
0.4 References	6
0.5 Definitions	7
0.6 Related Documents	7
1 PRE-EXAMINATION PROCESS	7
1.1 Training and Competency	7
1.2 Patient Preparation	7
1.3 Primary Blood Specimens	7
1.4 Spillage of Clinical and Non-Clinical Waste	8
1.5 Retention of Clinical Material and Records	8
1.6 Disposal of Clinical and Non-Clinical Waste	8
2 EXAMINATION PROCESS	8
2.1 Analytical Principle	8
2.2 Limitations of Analytical Procedure	8
2.2.1 The Device	8
2.2.2 Cross-reactivity	8
2.2.3 Interferences	8
2.2.4 Specimens Unsuitable for Analysis	8
2.2.5 Safety Notices	9
2.3 Central Database and Connectivity	9
2.4 Equipment	9
2.5 Reagents/Consumables	9
2.6 Calibration	9
2.7 Quality Assurance Programme	9
2.7.1 Internal Quality Control	9
2.7.2 External Quality Assessment	10
2.8 The Test Procedure	10
2.9 Results	11
2.10 Maintenance	11
2.10.1 As required	12
2.10.2 Weekly	12
2.11 Troubleshooting	12
3 POST-EXAMINATION PROCESS	13
3.1 Reference Values	13
3.2 Reporting the Results	13
3.3 Interpretative Comments	13
3.4 Reviewing Previous Results	14
3.5 Performance Criteria	14
3.6 Documentation	14
3.7 Audit of Results and Indication	14
3.8 Comparability	14
3.9 Uncertainty of Results	14
APPENDIX 1 SELF ASSESSMENT COMPETENCY FORM	15

PRINTED COPIES UNCONTROLLED UNLESS MARKED "CONTROLLED" IN RED

0 INTRODUCTION

0.1 Purpose and Scope

The Pfaff Bilimeter III bilirubinometer is used for the determination of bilirubin levels in serum for the investigation of neonatal jaundice. This document describes the procedure to be followed to ensure the provision of reproducible and accurate patient results for the Clinician to act upon.

The bilirubinometers are located on Wards 33 and 32 and are used to process specimens from the respective wards. There is also a Hettich Micro haematocrit Centrifuge for separation of plasma from cells.

Bilirubin may also be analysed on some blood gas analysers within the Trust and on the laboratory analysers.

0.2 Clinical Indication

Neonatal jaundice associated with an elevated serum bilirubin level is a common physiological condition occurring shortly after birth due to immaturity of hepatic systems and usually resolves without treatment. However, a high concentration of bilirubin can cause brain damage (kernicterus) and requires investigation of the underlying cause. A dangerously high level of bilirubin is treated by phototherapy or exchange transfusion. Regular monitoring of serum bilirubin is essential in neonates with jaundice and this is most effectively done using Point of Care Test (POCT) methodology.

When following trends in bilirubin results, users must be made aware that results may differ between POCT analysers or POCT analysers and the laboratory method. Analytical variation is discussed in the user training sessions. The laboratory recommends review of results from one method used in series wherever possible, this reveals trends in results most clearly.

0.3 Responsibilities

POCT Management Group ensures:

- The responsibilities, authority and interrelationships of all personnel involved in POCT are specified and communicated with the organisation;
- Staff performing PCOT receive appropriate training, supervision and competence testing;
- The selection of POCT devices and systems includes their practicability and the comparability of their results with those obtained in the laboratory;
- The reports of the POCT quality assurance programme(s) are reviewed by the group and advice on improvement is provided and implemented.

Senior Pathology Management Staff ensure:

- Departmental policies and procedures are in accordance with the requirements of the MHRA bulletins, Trust Medical Devices Management Policy, and all relevant quality standards to which the organisation subscribes;
- The needs and requirements of the users are reviewed regularly;
- Persistent problems are escalated to the relevant manufacturer and actioned;
- Determining and documenting back up plans if equipment is out of action;
- Responds to serious service loss in liaison with Bleep 555 holder.

Clinical Biochemistry Staff with POCT responsibilities ensure:

- Regular tasks in the management of the POCT service are carried out;
- Persistent problems are notified to senior departmental staff.

Clinical Biochemistry staff out of hours:

- Should not leave the laboratory to trouble-shoot POCT as this compromises the on-call service
- Should instigate back up procedures, where appropriate, and cascade serious service loss as appropriate to the on-call Consultant.

PRINTED COPIES UNCONTROLLED UNLESS MARKED "CONTROLLED" IN RED

Users ensure:

- Their annual training is kept up to date and they are competent to practice;
- The correct specimen is obtained from the correct patient under the correct circumstances and the correct specimen is tested;
- Ensure the specimen is collected with minimal discomfort to the patient;
- The test is carried out according to the manufacturer's instructions;
- The specimen is disposed of according to the Trust Waste Management Policy;
- Clinical biochemistry POCT staff are informed in the event of problems with the device;
- If issued with a log-on or barcode this is kept secure and used only by themselves.
Sharing or photocopying barcodes or log-on details is a disciplinary offence according to Trust Information Technology Security Policy;
- The results are reported as appropriate.

Ward / Department managers ensure:

- Users of the equipment keep their training up to date and maintain their competency;
- Adequate supplies of consumables are maintained;
- Appropriate resources are available;
- Proper records of training and competency are maintained;
- That faulty equipment is taken out of use;
- Notifies the POCT team if equipment becomes obsolete or no longer required.

Infection control team / Nurse Safety representative ensures:

- Performs risk assessment for handling biological material;
- Advises on precautions to be undertaken when handling the biological material including spillages, and safe disposal;
- Advises on safety of products before purchase;
- Updates COSHH details and advises on precautions to be taken when handling all material including spillage, breakage and disposal.

Manufacturers ensure:

- The user is kept up to date with information regarding the device and the associated products that they provide;
- Provide training sessions for the user in order to provide competent use of the instrument;
- Provide regular audits of the equipment;
- Fulfil the criteria agreed in the maintenance and service contract.

0.4 References

- Management of In Vitro Diagnostic Medical Devices, MDA Bulletin DB2002 (02)
PD-GEN-MDAIVDMg
- Management and Use of IVD Point of Care Test Devices, MDA Bulletin DB2010 (02)
PD-GEN-MDAPOCTMg
- Trust Point of Care Testing Policy **MP-GEN-LDHPOCTPo**
- Trust Infection Control Policy **MI-GEN-LDHInfCrl**
- Point of Care Testing Policy in Blood Sciences **LP-BLS-POCTPol**
- ACB National Audit Group and UKNEQAS Clinical Chemistry Specialist Advisory Group for Paediatric Investigations 'Investigation of neonatal jaundice: standards for the measurement of total and conjugated bilirubin and interpretation' **MP-BIO-ACBPdBIGd**

PRINTED COPIES UNCONTROLLED UNLESS MARKED "CONTROLLED" IN RED

0.5 Definitions

User – Any person who handles the device whether it is used directly to produce results or indirectly for maintenance and Quality Assurance procedures. This includes Clinicians, Nursing Staff, Healthcare Scientists and Medical Device Technicians as necessary.

EQAS – External Quality Assurance Scheme

POCT – Point of Care Testing

SOP – Standard Operating Procedure

IQC – Internal Quality Control

MHRA – Medicines and Healthcare Products Regulatory Agency

TPN – Total Parenteral Nutrition

0.6 Related Documents

- Trust Health and Safety Policy **MP-GEN-LDHH&SPol**
- Trust Waste Management Policy **MP-GEN-LDHWsMgP**
- Trust Information Technology Security Policy **MP-GEN-LDHIM&TPO**
- Bilimeter 3 + Prufsonde User Manual **LI-BIO-NPTBil3Ur**
- Hettich Haematocrit 210 Centrifuge Manual **LI-BIO-NPTHctHet**
- Trust Guidelines for the management of neonatal jaundice (CG91)
- Neonatal Jaundice (NICE guideline)
- Patient Results and Quality Control Book
- Laboratory Bilirubin record book
- POCT Annual Self-Assessment Competency Statement for Pfaff III Bilirubinometer (See Appendix 1)

1 PRE-EXAMINATION PROCESS

1.1 Training and Competency

Bilirubin analysis should only be carried out by trained and competent users. Training sessions can be arranged by contacting the POCT team on Ext 7991. Staff working on Wards 32 or 33 can be trained. Training records and self-assessment competency documents are kept by the POCT team. The self-assessment competency forms are available as an Appendix to this document.

1.2 Patient Preparation

The area of skin to be lanced should be warm and relaxed. It should be cleaned with water and dried to avoid interference from contaminants on the skin or dilution from water used for washing.

Care should be taken if babies are on TPN containing lipids, as highly lipaemic specimens (serum appears cloudy) are not suitable for analysis. The specimen should ideally be taken 4 hours after the TPN has been stopped in this case. Measurement of bilirubin by bilirubinometers (direct spectrophotometry) is only possible in neonates

1.3 Primary Blood Specimens

Blood from heel pricks using a Unistik must be collected into a bilimeter capillary. These are glass and have a small green band at one end. Only the glass capillaries supplied by Sartorius are suitable for use with this analyser, any other capillary will give unreliable results. A very small amount of serum is required for bilirubin analysis and the capillary should only be ½ to 2/3 filled with blood.

Bilirubin is photosensitive both to artificial light and direct sunlight and specimens should therefore be analysed immediately. .

Alternatively blood can be taken from indwelling venous or arterial lines.

PRINTED COPIES UNCONTROLLED UNLESS MARKED "CONTROLLED" IN RED

1.4 Spillage of Clinical and Non-Clinical Waste

Spillage of clinical and non-clinical waste must be dealt according to the Trust Health and Safety Policy.

1.5 Retention of Clinical Material and Records

There is no clinical material to keep.

Bilirubinometer glass capillaries should be disposed of into a sharps bin after use.

Used Cristaseal putty strips should be disposed of as clinical waste once completely used.

Patient results must be transcribed into the patient's notes.

Patient results must also be immediately transcribed into the Patient Results and Quality Control Book and the patient's notes. The Patient Results and Quality Control Book must be kept beside the analyser. When full these books must be kept in a secure place on the ward for the lifetime of the instrument (or if required by other guidelines, for the longest time period deemed necessary).

Training and self assessment of competency records are kept by the POCT team.

1.6 Disposal of Clinical and Non-Clinical Waste

Clinical and non-clinical waste must be disposed of according to the Trust Health and Safety Policy.

2 EXAMINATION PROCESS

2.1 Analytical Principle

The Pfaff Bilimeter III bilirubinometer is a dual wavelength direct spectrophotometer. Light beams with wavelengths of 455nm (nanometers) and 575nm (nanometers) pass successively through the serum in the haematocrit capillary tube placed in the specimen holder. Bilirubin has a maximum absorption wavelength between 450 and 460nm (nanometers) and haemoglobin at 575nm (nanometers). The possible effects of haemoglobin present in sera derived from haemolysed specimens are compensated for by deducting the absorbance at 575nm (nanometers) from that at 455nm (nanometers).

2.2 Limitations of Analytical Procedure

2.2.1 The Device

The specimen size required is ~4μL (microlitres).

The manufacturer's quoted linear range is 0-488μmol/L (microlitres) (±3%)

The analyser is operational between 15-30°C at under 70% humidity (with no condensation). The analyser should be kept out of direct sunlight and dust should be avoided.

Opaque plasma is not suitable for analysis due to the presence of lipaemia. Severely haemolysed plasma (red in colour) is also not suitable for analysis.

2.2.2 Cross-reactivity

Any substance, which absorbs appreciably at 460nm (nanometers), is a potential interferent. These include lipochromes and carotenoids and these will affect results. In neonates such substances are usually absent, but are present with increasing age.

2.2.3 Interferences

Drugs such as amines, fatty acids, phenols and sulfonamides may interfere with bilirubinometer measurements. Intravenous fat emulsions have also interfered with direct spectrophotometric measurements.

2.2.4 Specimens Unsuitable for Analysis

The serum of newborn infants does not contain the potential interferents such as carotenoids and other pigments, which may be found in older children and adults. The bilirubinometer is thus only suitable for use in neonates (up to one month of age, or one month from term if premature).

PRINTED COPIES UNCONTROLLED UNLESS MARKED "CONTROLLED" IN RED

2.2.5 Safety Notices

There are no known safety notices

2.3 Central Database and Connectivity

There is currently no connectivity or central database for the Bilimeters.

2.4 Equipment

Pfaff Bilimeter III analysers and Hettich Haematocrit 210 centrifuges are located in the Clinical Room on Ward 32 (Serial number B39052133) and Ward 33 (Serial number B38051939).

The analyser and centrifuge require a power supply and adequate bench space. Access to hand washing facilities and suitable disposal facilities for biological waste are also required. The analyser must be installed on a hard level surface in a well-ventilated area.

The power switch is located on the front of the instrument. Before turning on the analyser pull out the cell holder and confirm that there is no capillary in the cell holder, that there is no dust etc around the inlet for the capillary tube or around the slit. Turn on the power switch with the cell holder inserted (but no capillary). The LCD will display 'Total Bilirubin'. A countdown will start from 30 to 0. When the display shows 0.0µmol/L (micromol/L) the device is ready for use. Quality control checks must be carried out before patient testing. Switch on the centrifuge at the rear of the machine and after a few seconds the 'Open Lid' symbol will light up. The lid can now be opened.

The analysers and centrifuge are covered by a maintenance contract.

2.5 Reagents/Consumables

Bilimeter capillaries, (Part number 105000) and Cristaseal putty (Part Number 106000) are ordered by the ward manager from Sartorius Ltd (01372 737100). The putty must be kept in a sealed bag to prevent it drying out.

Unistix Lancets Neonatal (Owen Mumford) are ordered on a stock requisition order form from supplies. 1 box contains 100 lancets (Part No. AT072)

Control cuvette and light proof box. The control cuvette must be kept in the light proof box at all times apart from the short time when it is tested on the analyser.

PDI Sanicloth 70 (stock requisition).

All the above supplies are kept in the clinical room on Wards 32 and 33 respectively at room temperature. Stock records must be kept by the wards.

Internal Quality Control Liquicheck Paediatric Control Level 2 from Bio-Rad. It is only stable for 1 week and is stored in Clinical Biochemistry Fridge 003. It is the same as that used in the laboratory and is tested by Biochemistry staff weekly.

2.6 Calibration

The analyser is calibrated at manufacture and this is checked yearly as part of the engineer's maintenance visit. The users can not perform calibration processes.

2.7 Quality Assurance Programme

2.7.1 Internal Quality Control

Ward Staff:

Ward staff must carry out quality control using the control cuvette before every patient test. The control cuvette should additionally be processed if the analyser has been moved or switched off.

The control cuvette is contained in its own separate light proof box and **must** be kept in this **closed** box when not in use.

1. Insert the control cuvette and holder into the meter, matching the red dot on the holder with the red dot on the meter.
2. The result will be displayed almost immediately. Do not press any buttons!
3. The result must be within the range:

Ward 32: 264-304µmol/L (micromole/L), target 284 µmol/L

PRINTED COPIES UNCONTROLLED UNLESS MARKED "CONTROLLED" IN RED

Ward 33: 173-199µmol/L (micromole/L), target 186 µmol/L

4. If the result is in range, proceed with patient testing. If the result is outside the range DO NOT USE the instrument. Follow the instructions in section 2.7 and 2.8.
5. Record all control values in the patient's results and quality control book.
6. Return the control cuvette to the light proof box and close the lid.

Laboratory POCT staff:

Internal Quality Control will be processed weekly by Clinical Biochemistry staff and when troubleshooting. The quality control material is Liquicheck Paediatric Control Level 2 from Bio-Rad. It is only stable for 1 week. Results are recorded in the Patient's Results and Quality Control Book and in the laboratory Bilirubin book. Out of range IQC results will be investigated.

2.7.2 External Quality Assessment

UKNEQAS Paediatric Bilirubin scheme, monthly, 3 specimens will be processed by Clinical Biochemistry staff. The results are recorded in the Patients Results and Quality Control Book. Poor performance will be investigated by Clinical Biochemistry Staff and is defined as an A score over 200.

POCT staff record EQA failure on 'Laboratory Non Conformity Form' for later review. (A CAPA template has been set up in Q Pulse for use for management of non conformities).

The Ward Managers will receive a simplified report indicating how their analyser is performing. This should be kept with the equipment records for the analyser.

2.8 The Test Procedure

Each user of the Bilirubinometer and the centrifuges must be trained and be self-assessed as competent (See Appendix 1 – Self assessment competency form).

If the analyser is out of commission, report the problem to Clinical Biochemistry Staff (Tel 7991 Monday to Friday, 9am to 5pm). If the analyser is not available, back up may be arranged between the midwife in charge on wards 32 and 33. If both analysers are out of action, an orange paediatric specimen tube should be sent to the laboratory for bilirubin analysis by the laboratory method.

Gloves must be worn when handling blood specimens, using the bilirubin analyser or centrifuge. Specimens should be centrifuged and analysed immediately because bilirubin is photosensitive both to artificial light and direct sunlight.

Test Procedure:

1. Collect blood into the glass capillary tube, holding the end with the green line. Collect blood to fill only ½ to ⅔ of the capillary. When the required volume of blood has been collected, seal the capillary tip with the Cristaseal putty. Hold the end of the capillary with the green line and insert the capillary into the putty at right angles to the face of the putty. Push in and then twist and pull out the capillary. This should ensure a secure seal. Do not use the Cristaseal if it is dry and hard.
2. To open the lid, turn the blue handle on the front panel to the left. Only do this if the 'Open Lid' symbol is lit up. NEVER attempt to open the centrifuge whilst it is running.
3. Remove the plastic retaining plate with haematocrit scale by pushing down on the grey button in the centre.
4. Ensure that your capillary is thoroughly sealed.
5. The sealed capillary should be placed into the micro-haematocrit centrifuge with the **sealed end facing outwards**. If you place the specimen incorrectly, the specimen will be lost and the centrifuge will require careful cleaning.

PRINTED COPIES UNCONTROLLED UNLESS MARKED "CONTROLLED" IN RED

6. Place specimens in centrifuge in pairs so that the specimens are diametrically opposite each other. If only one specimen is being tested a balance (empty) capillary should be used.
 7. Replace the retaining plate, pushing it down firmly until it clicks into place.
 8. Close the lid of the centrifuge by turning the blue handle on the front to the right and press 'start'. The display is preset for 3 minutes for bilirubin. The centrifuge will spin for 3 minutes.
 9. When running the symbol  will rotate and the time remaining will be displayed with the decimal point flashing.
 10. When the time is up, both displays will flash. It should not be possible to open the lid until the centrifuge has come fully to rest.
 11. Open the lid when indicated by turning the handle to the left, remove the retaining lid by pressing down on the grey button in the centre and remove the specimen carefully from the centrifuge keeping the sealed end down.
 12. Replace the retaining plate and close the centrifuge lid.
 13. Place the capillary in the cell holder; ensure that only serum and no red cells or air are visible through the measurement slit.
 14. **Check that the sides of the capillary are clean and free from fingerprints, plasma, dirt, dust or smudges etc. and that there are no air bubbles or red cells within the capillary in the area through which the measurement will take place i.e. where the slit will be. Clean the outside of the cuvette as necessary with a tissue.**
 15. Place the cell holder into the instrument matching the red dot on the holder with the red dot on the meter. A result in $\mu\text{mol/L}$ (micromol/L) will automatically be displayed within 1 second.
 16. Record the patient's result in the Patient Results and Quality Control Book and into the patient's notes.
 17. Remove the capillary from the holder and discard in a sharps bin and replace the empty holder in the meter. Leave the area clean and tidy.
- N.B if a capillary is left in the instrument light will denature the bilirubin in the serum causing results to be lower than the true initial result.

2.9 Results

Record the results in the Patient Results and Quality Control Book and in the patient's notes. The Patient Results and Quality Control Book must show the following:

- Date
- Time
- Quality Control Result
- Patient's Hospital number or NHS number
- Patient's first and last names in full
- Patient's result
- Comments/Action taken
- User's Signature
- Print name of user

The Patient Results and Quality Control Book must be kept beside the analysers. When full, these books must be kept in a secure place on the Ward for the lifetime of the instrument (or if required by other guidelines, for the longest time period necessary). Any problems, error codes or instrument downtime should also be recorded in the book.

2.10 Maintenance

The operator is responsible for routine maintenance of the analyser.

PRINTED COPIES UNCONTROLLED UNLESS MARKED "CONTROLLED" IN RED

2.10.1 As required

The centrifuge should be inspected when used and cleaned if necessary. Sanicloths should be sufficient for any small blood splashes on the outside or inside.

If the slit needs cleaning because of putty or dust/dirt in the slit, clean this out gently with a cotton bud soaked in water or the edge of a PDI sanicloth 70.

Cleaning the centrifuge:

Wear gloves. If the capillary has broken in the centrifuge, or if the end was not sealed properly, and the centrifuge is full of blood and glass, you **MUST** clear it up.

DO NOT use hypochlorite. Use 10% Virkon, 10% Microsol or other acceptable agents such as Sanicloths and wipe down all surfaces thoroughly and dry them. Remove the curved end caps from the capillary holder if dirty and soak in cleaning agent. Remove any adhering Cristaseal and dried blood with a damp cotton bud. Dry thoroughly and replace in their correct positions. Be careful of glass shards. Dispose of all broken glass in a sharps or glass disposal bin, and mop up any spills on the work surfaces.

The clear plastic platen with haematocrit **must not** be cleaned with any abrasive material or hypochlorite, or the scale will be wiped off. Use detergent only.

If blood has contaminated the inner bowl of the centrifuge please DO NOT attempt to remove the specimen platen. Incorrect handling can cause damage. Contact Clinical Biochemistry (Ext 7991) if there is severe contamination of the platen and inner bowl.

2.10.2 Weekly

The exterior of the analyser and centrifuge should be cleaned with a PDI sanicloth 70. Turn off the analyser before cleaning the instrument. Once the exterior is clean and dry the analyser can be turned on again.

The surrounding area must be decontaminated with a biocidal agent use Sanachlor/Actichlor tablets (in solution) or PDI 70 Sanicloths.

2.11 Troubleshooting

Problem/Error displayed	Probable Cause	Solution
No display when power is switched on	The power cord is not properly connected	Connect the power cord correctly
'Lamp is dead'+	Wrongly inserted cell holder Slit blocked with dust/dirt Specimen severely haemolysed	Insert cell holder correctly Clean the slit Obtain a new patient specimen
'Lamp is dark'	Zero adjustment performed whilst patient specimen inserted Slit blocked with dust/dirt	Remove the capillary tube and re-perform zero adjustment Clean the slit
Display does not read 0.0µmol/L(micromol/L) when no capillary tube is inserted		Perform zero adjustment
'F' is displayed before the measured value	Zero adjustment performed whilst patient specimen inserted If the indication shows 'Lamp is dark', 'F' will be displayed before the measured value. These results must not be used for diagnostic purposes.	Perform zero adjustment and/or clean slit

PRINTED COPIES UNCONTROLLED UNLESS MARKED "CONTROLLED" IN RED

Measured values are too high or too low	Zero adjustment performed whilst patient specimen inserted	Perform zero adjustment
	Slit is not overlying serum	Place the capillary tube in the slit correctly.
	An incorrect capillary tube has been used.	Ensure only the designated capillaries (with green band at one end) are used.
QC out of range		Repeat, if still out of range perform zero adjust and repeat. If still experiencing problems clean the slit and then contact the laboratory (details below).

Performing zero adjust:

1. Place the cell holder into the instrument, but ensure that no capillary is inserted.
2. Press Re-zero.
3. The bilimeter should display 0.0.
4. Re-check the control cuvette.

N.B. If re-zero is pressed whilst a capillary is inserted in the cell holder any patient testing carried out after this will be inaccurate and must not be used for patient testing until a zero adjustment is carried out correctly.

If errors occur repeatedly call Clinical Biochemistry staff on extension 7991, Monday to Friday, 9am-5pm, as there may be a more serious problem. If at any time the bilirubinometer is unavailable for use send specimens to the laboratory who will provide a back up service. The analysers are covered by a service contract with Sartorius who can be contacted on 01372 737102.

3 POST-EXAMINATION PROCESS

3.1 Reference Values

There is no reference range as such for bilirubin in neonates. Treatment levels for phototherapy are determined on the wards in conjunction with the Trust Guidelines for the management of neonatal jaundice (CG91) which includes nomograms with the levels for phototherapy and exchange for babies of varying gestations.

3.2 Reporting the Results

The person performing the test is responsible for the results. If there is an unexpected or very abnormal result, check procedure and repeat IQC. If OK repeat patient test. The requesting doctor should be informed immediately of any unexpected or clinically significant results. If there is concern over the results contact Senior nursing staff or the POCT team for further advice.

3.3 Interpretative Comments

Up to 50% of newborn infants develop jaundice after 48 hours of age; the jaundice being described as physiological, arising from the changeover from foetal to neonatal bilirubin metabolism. Total bilirubin does not usually rise above 200 $\mu\text{mol/L}$ (micromol/L), and in >90% of neonates has fallen to within the adult reference range (6-28 $\mu\text{mol/L}$ {micromol/L}) by the 7-10th day.

PRINTED COPIES UNCONTROLLED UNLESS MARKED "CONTROLLED" IN RED

Jaundice ($>50 \mu\text{mol/L}$ {micromol/L}) persisting beyond the 11-14th day in a term infant is unusual and warrants investigation.

For initial diagnostic specimens with results greater than $300 \mu\text{mol/L}$, a specimen MUST be sent to the Biochemistry laboratory for confirmation before intervention. The Bilimeter III has an increasing negative bias above this level. It is not necessary to send blood to the laboratory on consecutive follow-up specimens to check the response to phototherapy.

Care should be taken when interpreting results as lipaemia can significantly affect results.

In addition to this air bubbles within the specimen or any red blood cells carried over after centrifugation can also affect results.

Contamination of the outside of the bilirubinometer capillary can also affect results e.g. fingerprints, large particles of dust etc.

3.4 Reviewing Previous Results

Previous results can only be reviewed in the patient's notes or from the Patient Results and Quality Control Record Book.

3.5 Performance Criteria

At a target level of $278 \mu\text{mol/L}$ (micromol/L) the acceptable imprecision is $\pm 10\%$.

3.6 Documentation

Patient results are recorded as detailed in Section 2.7

3.7 Audit of Results and Indication

EQA results are used to audit this process. In addition all POCT is subject to periodic horizontal and vertical audits and spot checks.

3.8 Comparability

Comparability with the laboratory method is available by review of EQA results and a limited amount of data is also available from the POCT team.

3.9 Uncertainty of Results

The degree of uncertainty for results is expected to be within acceptable limits provided all calibration and control procedures have been followed correctly and that the sample collection, identification and processing has been followed correctly.

See Section 2.2 Limitation of Analytical Procedure and 3.5 Performance Criteria

PRINTED COPIES UNCONTROLLED UNLESS MARKED "CONTROLLED" IN RED

APPENDIX 1 SELF ASSESSMENT COMPETENCY FORM

POCT Annual Self-Assessment Competency Statement for Pfaff III Bilirubinometer

Do not use this device unless you are competent to do so

Surname:

Forename(s):

Dept/Ward:

Job Title:

Self-verification of competence is undertaken by assessment against the statements below. These statements are designed to help you assess your competency to use this POCT (Point of Care testing) equipment. Responsibility for safe use remains with you, so if you are in any doubt regarding your competence to use the device, you should seek guidance.

You must be able to answer yes to all the questions before considering yourself competent and completing the competency statement in section B overleaf. If you are in any doubt or feel that you are not competent please complete section A over leaf and contact the POCT team (Ext 7991) or your manager to instigate further learning or training.

Ask yourself the following general questions:	Response
1. Have I attended a training session about the device in the last year?	Yes / No
2. Have I performed a patient test under supervision by a competent practitioner?	Yes / No
3. Have I read the Standard Operating Procedure (SOP) for the device? (These can be found on the Intranet under Departments, Pathology, Point of Care Services)	Yes / No
4. Do I know what infection control precautions to take when using the device?	Yes / No
5. Do I know the circumstances when the device should not be used including test interferences and limitations?	Yes / No
6. Do I know what action to take based on the results, including if the result is unexpected or suspicious?	Yes / No
7. Do I know when and where to seek help, or how to report an error or incident?	Yes / No
Ask yourself the following device specific questions:	
8. Do I know how to use the centrifuge safely including balancing and cleaning?	Yes / No
9. Do I know how to perform the quality control check before every test and what to do if the result is not in range?	Yes / No
10. Do I know what to do if an F is seen next to my result on the screen?	Yes / No
11. Do I know how to clean the cell holder if required and the impact if the glass capillary is not clean when a measurement is taken?	Yes / No
12. Do I know what should be recorded in the record book for every test performed?	Yes / No

PRINTED COPIES UNCONTROLLED UNLESS MARKED "CONTROLLED" IN RED

POCT Annual Self-Assessment Competency Statement for Pfaff III Bilirubinometer



Do not use this device unless you are competent to do so

Surname:

Forename(s):

Dept/Ward:

Job Title:

Complete either Section A or B whichever is appropriate to you.

Section A

I require further training before I can use the Pfaff III Bilirubinometer POCT device in a competent manner.

Signed:.....

Date:.....

Section B

I certify that I am aware of my responsibility for continuing professional development and I understand that I am accountable for my actions. With this in mind I make the following statement:

I am competent to use the Pfaff III Bilirubinometer POCT device

Signed:.....

Date:.....

Please return this form to the POCT team, Clinical Biochemistry, Luton & Dunstable Hospital. Telephone Ext 7991. We suggest you also keep a copy for your personal records.