

PRECICE® dCK Screening Assay Kit


For research use only. Not for use in diagnostic procedures

I. Introduction


dCK Assay Kit was specially designed to follow the enzymatic activity of purified deoxycytidine kinase (dCK, EC 2.7.1.74) *in vitro*. This user manual gives the instructions for standard assays in 96-well plate.

The principle of the assay is based on the use of deoxyinosine (dIR) as a substrate of dCK and a coupled reaction involving a highly active IMPDH (Inosine Monophosphate Dehydrogenase, bacterial recombinant) for a direct measurement of the deoxyinosine monophosphate (dIMP) formed by dCK.

- (1) In the presence of dIR and ATP, dCK catalyses the phosphorylation of dIR to form dIMP and ADP:



- (2) The dIMP formed is then oxidized to deoxyxanthosine monophosphate (dXMP) by IMPDH in the presence of NAD, leading to NADH₂ formation.



This coupling reaction is immediate when IMPDH activity is much higher than dCK activity in the assay. The enzymatic activity of dCK, which corresponds to the formation kinetics of dIMP, is then stoichiometrically and directly monitored by the formation kinetics of NADH₂.

The velocity of NADH₂ formation is measured with a spectrophotometer at 340nm (molar extinction coefficient of NADH₂ at 340nm = 6220 M⁻¹.cm⁻¹)

II. Kit Contents

A standard **PRECICE® dCK Screening Assay Kit** (one 96-well plate) contains:

1. one tube " **IMPDH** " (800mU)
2. one tube " **Cofactor 1** " (DTT, 17mg)
3. one tube " **Cofactor 2** " (NAD, 66 mg)
4. one tube " **Cofactor 3** " (ATP, 33 mg)
5. one vial " **Reaction Buffer 5x** ", 4ml
6. one tube " **Deoxyinosine** " (50mM, 1.2mL)
7. one tube " **Human dCK enzyme** ", 100mU
8. Transparent 96-well plate (round-bottom 96-well plate Corning, Costar®, ref. 3797)

III. Equipments required

- 1) Plate agitator
- 2) Plate reader fitted with a filter 340nm.

IV. Storage

PRECICE® dCK Screening Assay Kit must be **stored at -20°C** until used.

IMPORTANT:

The following instructions are given to measure the activity of dCK enzyme *in vitro*, in a range allowing this measurement by spectrophotometry as described here below. NovoCIB does not guarantee the use of its PRECICE® dCK Screening Assay Kit or of one or several of its components, in other conditions than those described in this user manual and/or for other purpose than R&D.

V. Experimental Protocol I - for non-nucleoside dCK inhibitors

V.1. Reconstitute IMPDH enzyme

Add 250µL of deionized water to "IMPDH" tube. Agitate gently until complete dissolution of the powder.

V.2. Preparation of standard reaction buffer (1x)

- i) Add the content of "Reaction Buffer 5x" tube (4ml) to 16ml of deionized water to prepare "Reaction Buffer 1x".
 - ii) Transfer quantitatively the content of 3 tubes with "Cofactor 1", "Cofactor 2", "Cofactor 3" to the tube with "Reaction buffer 1x".
- To do so:
- iii) -pipet 1ml of "Reaction buffer 1x" to each tube with cofactors and mix them by inverting or pipeting up and down until the powder dissolved.
 - iv) - transfer by pipeting the content of three tubes back into a vial "Reaction buffer 1x";
 - v) - repeat to be sure that all reagent and enzymes of the small tubes and vial are recovered. mix by gently inverting until complete dissolution. Avoid bubbles.
 - vi) Transfer quantitatively the content of "IMPDH" tube to "Reaction buffer 1x" with co-factors.
 - vii) Solubilize the content of "Human dCK enzyme" tube by adding 1ml of complete "Reaction buffer 1x" with co-factors and IMPDH, transfer by pipeting the content of the tube back into a vial "Reaction buffer 1x".

Composition of complete reaction buffer: 100mM Tris-HCl, 250mM KCl, 10mM MgCl₂, BSA 0.5mg/ml, 5mM NAD, 2.5mM ATP, IMPDH 40mU/ml; human recombinant dCK 5mU/ml.

VI. Following dCK activity in vitro

VI.1. Pre-incubation phase (15')

- i) Program the plate-reader in a kinetics mode with the measurements done every 1 minutes, absorbance at 340 nm, 37°C, agitation before the kinetics for 1 min, duration time 15min.
- ii) Add 200µL of standard reaction buffer per well.
- iii) Agitate and measure absorbance at 340nm (A₃₄₀). Record this first set of data.

5.2.2. Start the reaction and incubate (40')

- i) Eject the plate from the plate-reader
- ii) Program the plate-reader in a kinetics mode with the measurements done every 1 minutes, absorbance at 340 nm, 37°C, agitation before the kinetics for 1 min, duration time 30-40min.
- iii) Start the reaction by adding 10µL of 50mM of "Deoxyinosine" per well
- iv) Place the plate in the plate-reader and start the measurements. Record second set of data.

VII. Experimental Protocol II; for nucleoside analogues dCK inhibitors

VII.1. Reconstitute IMPDH enzyme

Add 250µL of deionized water to "IMPDH" tube. Agitate gently until complete dissolution of the powder.

VII.2. Preparation of standard reaction buffer (1x)

- i) Add the content of "Reaction Buffer 5x" tube (4ml) to 16ml of deionized water to prepare "Reaction Buffer 1x".
- ii) Transfer quantitatively the content of 2 tubes with "Cofactor 1", "Cofactor 2" to the tube with "Reaction buffer 1x".

Important: Do not add "Cofactor3" (ATP)!

To do so:

- iii) -pipet 1ml of "Reaction buffer 1x" to each tube with cofactors and mix them by inverting or pipeting up and down until the powder dissolved.
- iv) - transfer by pipeting the content of three tubes back into a vial "Reaction buffer 1x";
- v) - repeat to be sure that all reagent and enzymes of the small tubes and vial are recovered. mix by gently inverting until complete dissolution. Avoid bubbles.
- vi) Transfer quantitatively the content of "IMPDH" tube to "Reaction buffer 1x" with co-factors.
- vii) Add 1ml of "Deoxyinosine, 50mM" solution to "Reaction buffer 1x" with cofactors and IMPDH.
- viii) Solubilize the content of "Human dCK enzyme" tube by adding 1ml of complete "Reaction buffer 1x" with co-factors and IMPDH, transfer by pipeting the content of the tube back into a vial "Reaction buffer 1x".

Composition of complete reaction buffer: 100mM Tris-HCl, 250mM KCl, 10mM MgCl₂, BSA 0.5mg/ml, 5mM NAD, 2.5mM inosine, IMPDH 40mU/ml; human recombinant dCK.5mU/ml.

VII.3. Preparation of ATP solution for starting the reaction

Add 1.1ml of water to the tube "Cofactor 3" containing 33mg of ATP powder. Mix until dissolved.

VIII. Following dCK activity in vitro

VI.1. Pre-incubation phase (15')

- iv) Program the plate-reader in a kinetics mode with the measurements done every 1 minutes, absorbance at 340 nm, 37°C, agitation before the kinetics for 1 min, duration time 15min.
- v) Add 200µL of standard reaction buffer per well.
- vi) Agitate and measure absorbance at 340nm (A₃₄₀). Record this first set of data.

5.2.2. Start the reaction and incubate (40')

- v) Eject the plate from the plate-reader
- vi) Program the plate-reader in a kinetics mode with the measurements done every 1 minutes, absorbance at 340 nm, 37°C, agitation before the kinetics for 1 min, duration time 30-40min.
- vii) Start the reaction by adding 10µL of 100mM ATP per well
- viii) Place the plate in the plate-reader and start the measurements. Record second set of data.

Assay validation:




Figure 1.

Nucleosides and nucleoside analogues (2µL of 100mM solution) were added to 96-well plate, followed by the addition of 200µL of standard reaction mix containing DTT, NAD, IMPDH and dCK. 10 min later, the reaction was started by the addition of ATP solution (10µL per well) and the increase in absorbance at 340nm was followed for 30min.

Results:

Of 8 deoxy- and ribonucleosides tested, only deoxycytidine and gemcitabine competitively inhibit deoxyinosine phosphorylation by dCK which totally consistent with previously published data on their substrate properties.