

GelAnalyzer 2010 User's manual

Contents

1. Starting GelAnalyzer	2
2. The main window	2
3. Create a new analysis	2
4. The image window	3
5. Lanes mode	3
5.1 Detect lanes automatically	3
5.2 Delete lane(s)	3
5.3 Add new lane manually	3
5.4 Modify existing lanes	4
6. The profile window	4
7. Bands mode	5
7.1 Automatic band detection	5
7.2 Add new band manually	5
7.3 Delete bands on the selected lane	5
7.4 Modify existing bands on selected lane	5
8. Background subtraction mode	5
8.1 Manual baseline tool	6
8.2 Rolling ball background subtraction	6
8.3 Valley to valley background subtraction	6
9. Rf calibration mode	6
10. Quantity calibration mode	7
11. Molecular Weight (MW) Calibration mode	7
12. The analysis info window	8
13. Saving and opening analysis, exporting results	8
14. Language	8

This tutorial will lead you through the process of analyzing gel images with GelAnalyzer. The elements of the GUI and the functions of the program will be described as they occur in the course of a gel image evaluation.

1. Starting GelAnalyzer

GelAnalyzer does not need to be installed. All you need to do is to extract the downloaded archive to an arbitrary directory. To run the program, you need the Java Runtime Environment (JRE) to be installed on your system. JRE version 1.6 or later is required. In Microsoft Windows systems you may start the program by running „GelAnalyzer.exe“. In other systems you need to run the „GelAnalyzer.jar“ file. For further Information refer to the Java documentation on your system.

2. The main window

When GelAnalyzer has successfully started, you will see the main window. The main window consists of a menu bar, a tool bar, and a large empty space that acts as an inner desktop. The inner desktop will hold the other windows that only appear when you have either created a new analysis or opened an existing one.

The elements of the menu bar are the following:

- File menu: commands to create new analysis, open an existing one, save the current, and copy or export contents of inner windows and results.
- Analysis menu: here you will find the MW standard manager (described in details later), an option to switch the image type (described in the next section), and a command to bring up the analysis settings dialog.
- Windows menu: commands to move the inner windows front in the inner desktop, and to restore the default inner window arrangement.
- Help menu: currently only contains the command to bring up the info about GelAnalyzer.

The elements of the toolbar, except the new analysis and open analysis buttons, are inactive when there is no analysis loaded. The first three buttons serve the purpose of creation, loading and saving analysis. The other six buttons serve for switching between program modes. The program modes will be discussed in details later on. When you switch to a program mode, the second row of the toolbar will show the buttons corresponding to the operations in that mode.



3. Create a new analysis

To create a new analysis, use the command in the "File" menu or the button in the toolbar. After this, a file open dialog will be shown, where you need to select the gel image. Only PNG, JPEG, and GIF image formats are supported. If you wish to analyze an image in a different format, there are many free programs available that can make the conversion for you. When you have selected the gel image, a dialog will be shown where you need to select the image type. Choose the "light on dark" option when the spots on the image are lighter than the background, and the "dark on light" option in the opposite case. If you need to change image type

later on, you may do this by using the option in the "Analysis" menu. However, changing image type will cause the deletion of all bands on every lane, and the analysis results will be lost.

4. The image window

When you've successfully created a new analysis, the inner windows will appear. The upper-left corner that displays the gel image is called the "Image window". Operations in "Lanes mode", "Rf calibration mode", "Quantity calibration mode", and "MW calibration mode" take place in this window. The window has a toolbar with a zoom scale selection option.



Lanes mode

5. Lanes mode

In lanes mode you may detect, add, modify and delete lanes. All these operations take place in the image window. To switch to lanes mode, press the "lanes mode" button in the main window's toolbar. The default submode in lanes mode is lane selection. You may select lanes by clicking on them in the image window. The selected lane's profile will be plotted in the profile window, and the info about the selected lane's bands is shown in the analysis info window. It is important that only vertical lanes are allowed, and the lanes are always scanned from the top to the bottom. If your image's orientation differs, use some image manipulation software to rotate it. Hopefully there will be some basic image manipulation possibilities in the following version.



Define lane detection ROI



Detect lanes

5.1 Detect lanes automatically

When you want to detect lanes automatically you can select the region of interest (ROI) on the image. To do this, press the "define lane detection ROI" button. Left-click on the image where you want the upper-left corner of the area to be placed, and move the mouse to draw the rectangle. Left-click again to finalize the ROI. The rectangle's sides can be grabbed with the mouse, so you can modify it. To delete the ROI, right-click on the image. If no lane detection ROI is defined, the automatic lane detection will affect the whole image. To detect the lanes, press „Detect lanes” button. It is important that the image type is properly set for automatic lane detection.



Delete selected lane



Delete all lanes

5.2 Delete lane(s)

- To delete the selected lane press the "Delete selected lane" button.
- To delete all lanes press the "Delete all lanes button".



Add new lane

5.3 Add new lane manually

To add a new lane, press the "add new lane" button. Now you can define the new lane with 3 clicks on the image. First select the upper point of the lane then select the end point. Now if you move the mouse, you can see that the lanes width is changing. Click with the mouse when you have set the right lane width. The new lane will be added.



Modify
lanes

5.4 Modify existing lanes

To switch to lane modification submenu press "Modify lanes" button. The control points of the lanes will be displayed. Use the mouse to add new control point, delete control points and modify the width of a lane.

- If you move the cursor above the lane's middle line, it will change into crosshair. Click where you want the new control point to be on the middle line. The new point will be shown, and if you don't release the mouse, you can move it. The lane curves are interpolating spline curves.
- To move a control point, drag it with the mouse.
- To delete a control point right-click on it. You can't delete the start and end point of the curve, because this would cause the deletion of the lane. To delete a lane use the delete lane option.
- To modify a lane's width, grab the left or right side with the mouse, and drag it.

When you modify a lane, its profile will be rescanned, and if the lane already had bands defined on its profile, those will be lost.

6. The profile window

The profile window displays the intensity profile of the selected lane. It also shows the bands on the selected lane and the background, if it is defined. The horizontal axis of the graph shows the distance from the lane's start position in pixels, while the vertical axis shows the intensity values. Since all input images are converted into 8 bit grayscale for processing, the maximal intensity value is 255. Also, the original lane image is shown below the graph. The shown profiles values are always scaled to fit the graph area vertically. To zoom on the horizontal axis, you may use the options in the profile window's toolbar's zoom box. These options are:

- Fix zoom values (100%, 200%, and 300%): the length of the profile is scaled with the chosen value. Use the horizontal scroll bar to scroll the profile if it is larger than the graph's width.
- Fit profile mode: the length of the profile is fit to the graph's width. Scrolling is not possible in this case.
- Manual zoom: with this tool, you can select the part of the profile to be displayed. When you select it from the options, the cursor turns into crosshair, and you can select the part of the profile with the mouse by pressing left mouse button and dragging the mouse. Right mouse click before area selection cancels the manual zoom. Maximally 1000% magnification is allowed this way. For a subsequent selection choose the option from the zoom tool again, and repeat the process.



Without
background

The window's toolbar has another button ("Without background") in order to display the profile by subtracting the defined background from the intensity values. Defining backgrounds and background subtraction will be explained later on.



Bands
mode

7. Bands mode

In bands mode you can detect bands automatically on the selected lane or on all lanes at once. You may manipulate the bands on the selected lane, namely, modify their boundaries, delete bands, and manually add new ones. The band modifying operations take place in the profile window. You can select the lanes in bands mode the same way as in lanes mode: in the image window.



Detect bands
on every lane

7.1 Automatic band detection

GelAnalyzer is capable of automatically detecting the bands on lanes. To use this option, press the button “Detect bands on selected lane” or “Detect bands on every lane”. When automatic band detection is applied, previous bands on the affected lane(s) will be deleted.



Detect bands
on selected
lane

7.2 Add new band manually



Add new
band

To add a new band to the selected lane, press the “add new band” button. The cursor in the profile window turns into crosshair, and you can define the boundaries of the new band by two clicks. The first one defines the new band’s start position, and the second one defines the end position. The band’s peak will be automatically

calculated as the position corresponding to the maximal intensity value within the new band’s boundaries.



Delete all
bands on
selected lane

7.3 Delete bands on the selected lane

You can delete an existing band, by right-clicking in it in the profile window. You can delete all the bands on the selected lane at once, by pressing the “Delete all bands on selected lane” button.

7.4 Modify existing bands on selected lane

In the profile window, you can grab the boundaries of the bands, and drag them. Also, the position of the peak sign can be modified the same way. If you alter the position of the band boundaries, the peak sign position will be automatically reset to the maximal intensity position. The peak sign’s position is important, since it marks the band’s Rf value, thus it is the basis of the MW calculation.




Background
subtraction
mode

8. Background subtraction mode

In background subtraction mode you may define the background of the individual profiles manually, or use built in methods to automatically define the background on a single lane or on all lanes at the same time. The defined background will be subtracted from the volumes of the bands. This way the uneven illumination of the image can be compensated, thus a more precise quantity calibration and evaluation can be achieved. The previously mentioned “without background” button in the profile window’s toolbar can be used to switch between the display modes when the background is shown separately, and when the background is subtracted from the profile. The usage of the display modes does not affect the calculations. Currently there are three possibilities in GelAnalyzer to define background. If you would like to define the background of an individual lane’s profile, use the manual baseline tool. For automatic background definition you may use the “valley to valley” and “rolling ball methods”. Every

change in the profiles' background results in the recalculation of the quantity calibration (if there is valid quantity calibration in the analysis).



8.1 Manual baseline tool

By using the manual baseline tool, you can define the background of an individual profile by adding and moving control points of the baseline. To activate the tool, press the “manual baseline” button in the main window’s toolbar. In the profile window the cursor turns to crosshair, and with one left-click you can define a vertical baseline. After this, when the cursor is crosshair, you may add new control points with a single left-click, or grab the existing control points and move them. Existing control points – except the first and the last one – can be deleted by right-clicking on them. The baseline modification only works when the “without background” mode is turned off.




8.2 Rolling ball background subtraction

Rolling ball is an automatic way to define the background of a profile. It can be imagined as when a ball is rolled on the inverted profile. Where the ball touches the profile will be the points of the background. The radius of the ball can be set when using the tool. Also, you can decide whether you want to use on all lanes, or just on the selected one.



8.3 Valley to valley background subtraction

Valley to valley is a rather simple background approximation method. It can only be used when there are bands on the lane. It works such a way that it connects the end and start points of the successive bands.



9. Rf calibration mode



Rf calibration is a very useful tool to repair gel run distortions. Rf calibration is performed by defining fix rf values on the image. This is done by adding and modifying rf curves. The rf curves are vertical spline curves. Each rf curve defines a fix rf value along the image. By default, the rf value of a band is calculated in comparison with the lane’s start and end positions. Using rf calibration, rf values of bands in lanes with different start and end positions become comparable. Enabled rf values are from the interval [0.0, 1.0]. To add a new rf curve press the “add rf curve” button in the main window’s toolbar. The cursor in the image window turn’s into crosshair. Click where you want the new rf curve to be placed. The start (0.0) and end (1.0) curves are automatically added. The new curve’s rf value will be calculated automatically. You may modify the rf value of the curves, except the start and end curves. To modify a curve’s rf value click in the area that displays the value. The rf values of the area’s border will turn into orange, and you may type in the new value. Press ESC to cancel the value modification. The start and end curves can be added separately by pressing the “add default start and end rf curves” button. You may add new control points to the rf curves, by clicking on them where the cursor turns into crosshair. The control points can be grabbed and moved with the mouse. Control points – except the first and the last one – can be deleted by right-clicking on them. Rf curves – except the start and end curves – can be deleted by right-clicking on them on a position where there is no control point present. The start and end curves can only be deleted by deleting the entire rf calibration. Note that

the MW calibration and evaluation is based on the Rf value of the bands. The band's Rf value is calculated from the peak sign's position.



Quantity
calibration
mode

10. Quantity calibration mode

Quantity measurements of the individual bands is based on fitting a curve to known {raw volume, calibrated volume} data pairs. The fitted curve is used to calculate the quantity volume of unknown bands. The bands whose quantity value is known are called calibration bands. To assign quantity value to known bands left-click on the circle that marks the band in the image window. A text field will appear where you can type in the value. Press ENTER to finalize the value, or ESC to cancel the process. When a value has been assigned to a band, its marking circle will turn into orange color, and the value will appear next to it. You may modify the value by left-clicking on the band mark again. To delete a band from the quantity calibration, right-click on the marking circle. GelAnalyzer currently contains four curve types for the fitting process. These are:

- Linear ($y = a * x + b$)
- Quadratic ($y = a * x^2 + b * x + c$)
- Linear log ($\ln(y) = a * x + b$)
- Log ($y = a * e^{b*x} + c$)

You may select the used curve from the analysis settings dialog. To open the analysis settings dialog, open the "Analysis" menu in the main window and click "Settings". To fit a linear or linear log curve you need at least 2 calibration bands. For quadratic and log curves at least 3 calibration bands are necessary. The fitted curve will be displayed in the calibration curves window. The cal. curves window has two tabs, one for the MW calibration curve, and the other one for the quantity calibration. If there are enough calibration bands for the selected curve type, the curve will be fitted, and all bands in the analysis will be evaluated immediately. The evaluation is entirely automatic, i.e. at the moment you perform any modification that involves the changing of the number or value of the calibration bands, the calibration curve will be refitted and all bands will be evaluated. Any modification of the calibration bands' raw volumes, including the modification of the background subtraction will result in the recalculation of the calibration curves, and the evaluation of all bands in the analysis. The cal. curves window's panels display the fitted curve, along with markers for the calibration bands. The fitted curve's equation is also shown here along with the R^2 (coefficient of determination) value of the fit.



MW
calibration
mode

11. Molecular Weight (MW) Calibration mode

MW calibration is entirely analogous to quantity calibration. MW calibration and evaluation is based on the known {Rf value, MW value} pairs of the calibration bands. Assigning MW values to bands, selecting curve fit is the same as in quantity calibration, thus these parts won't be discussed again. Besides assigning MW values to bands, GelAnalyzer provides tools to manage MW standards as unit. You may load previously saved standards to lanes, and save the MW values on a lane's bands as a standard. Also GelAnalyzer has an inbuilt MW standard manager, where you can manually create new standards, add or delete bps, and export it to text files, or import standards from external text files. GelAnalyzer uses its own file

format to store the MW standards. All the standards are stored in the “standards.mws” file. To manage the standards, open the standard manager from the analysis menu in the main window. The usage of the standard manager is rather straightforward, and it won’t be discussed thoroughly here.

12. The analysis info window

The analysis info shows the bands in the analysis along with all their important values. In the window’s “Show” menu you may select a display restrict: show bands on all lanes, show bands on the selected lane, show the mw or calibration bands only. The MW and quantity values of the calibration bands are shown in bold font type.



Save
analysis



Open
analysis

13. Saving and opening analysis, exporting results

You may save the entire analysis in GelAnalyzer’s own format. To save the current analysis, choose the “Save analysis” or “Save analysis as” options from the file menu, or the corresponding button from the main window’s toolbar. The analysis will be saved in a “.gap” file. In contrast with the previous version of GelAnalyzer, the image will be also saved in the file. You also have the possibility to export contents and/or copy the content of the internal windows. You may export and copy the current content of the image window, the profile window, and the cal. curve panels. The content of the analysis info window can be copied to the system clipboard, and pasted to any spreadsheet program for further usage.

14. Language

The language elements are loaded from the “selected.txt” file in the “language” directory. The structure of the file is the following:

- Lines that start with the ‘%’ character are comments, and won’t be used.
- Empty lines are also ignored.
- The structure of other lines is [key] = [value], where [key] identifies the language element, and [value] is the assigned text value.

You may create new language files to adapt GelAnalyzer to your own language. If an error occurs during the loading