




User Manual

QuantiGene[®] FlowRNA Assay

in situ gene expression assay for flow
cytometry

**For research use only.
Not for use in diagnostic procedures.**

Trademarks

Affymetrix® and  are registered trademarks of Affymetrix, Inc
QuantiGene® is a registered trademark exclusively licensed to Affymetrix.

All other trademarks are the property of their respective owners.

Limited License

Subject to the Affymetrix terms and conditions that govern your use of Affymetrix products, Affymetrix grants you a non-exclusive, non-transferable, non-sublicensable license to use this Affymetrix product only in accordance with the manual and written instructions provided by Affymetrix. You understand and agree that, except as expressly set forth in the Affymetrix terms and conditions, no right or license to any patent or other intellectual property owned or licensable by Affymetrix is conveyed or implied by this Affymetrix product. In particular, no right or license is conveyed or implied to use this Affymetrix product in combination with a product not provided, licensed, or specifically recommended by Affymetrix for such use.

Citing QuantiGene® FlowRNA Assay in Publications

When describing a procedure for publication using this product, please refer to it as the QuantiGene® FlowRNA Assay from Affymetrix.

Disclaimer

Affymetrix, Inc. reserves the right to change its products and services at any time to incorporate technological developments. This manual is subject to change without notice.

Although this manual has been prepared with every precaution to ensure accuracy, Affymetrix, Inc. assumes no liability for any errors or omissions, nor for any damages resulting from the application or use of this information.

Copyright

© 2012 Affymetrix Inc. All rights reserved.

Contents

About This User Manual	1
Contacting Technical Support	1
About the QuantiGene FlowRNA Assay Kits	1
How it Works	2
QuantiGene FlowRNA Reagent System Contents and Storage Conditions	3
Required Equipment and Materials Not Supplied	4
Experimental Guidelines	4
Assay Procedure	
I. Fixation and Permeabilization	6
II. Target Probe Hybridization	7
III. Signal Amplification and Detection	9
Appendix	
A1: Troubleshooting	10
A2: Lab setup guide	13
A3: Cytometer Setup	14
A4: Incorporating antibody staining with the QG FlowRNA Assay	17
A5: Expected results	19
A6: Validated cell lines	20
A7: Temperature Validation Procedure for Incubator	21

About This User Manual

This manual is for anyone who has purchased a QuantiGene FlowRNA Assay Kit for the following sample types:

- Cultured mammalian suspension cells
- Peripheral Blood Mononuclear Cells (PBMC) from cryopreserved or freshly collected whole blood and purified by Ficoll®.
- Suspension cells sorted by a flow cytometry sorting instrument, e.g. BD FACSAria®, or purified by antibody-coated magnetic beads

Contacting Technical Support



NOTE: For the most current version of user documentation, go to our website at www.ebioscience.com

For technical support, contact the appropriate resource provided below based on your geographical location. For an updated list of FAQs and product support literature, visit our website at www.ebioscience.com.

Location	Contact Information
North America	1-888-810-6168, tech@ebioscience.com
Europe	tech@ebioscience.com

About the QuantiGene FlowRNA Assay Kits

The QuantiGene FlowRNA RNA Assay System consists of 3 modules (each sold separately).

- *QuantiGene FlowRNA Fixation and Permeabilization Kit*. Contains reagents and buffers to ensure that the cell membranes are properly permeabilized and the RNA is fixed within the cell.
- *QuantiGene FlowRNA Assay Kit*. Contains the reagents for running the assay.
- *QuantiGene ViewRNA Probe Set*. Contains gene specific probe sets for specified RNA targets.

Sample Type	Suspension Cells including PBMC, bone marrow cells and suspension cell lines.
Species	mammalian
Plex level	Up to 3 RNA targets
Assay format	1.5-mL microfuge tube
Detection	fluorescent
Instrumentation for detection	Standard flow cytometer equipped with blue (488 nm) and red (633 nm or similar) lasers, and filter sets for FITC, APC and APC-Cy7


How it Works

The QuantiGene® FlowRNA Assay is an expansion of Affymetrix's QuantiGene ViewRNA Assay. The QuantiGene FlowRNA assay is an *in situ* hybridization assay enabling simultaneous detection of up to three RNA targets with high sensitivity and single-cell resolution for use with a standard flow cytometer. This assay can be combined with immunophenotyping using fluorochrome-conjugated antibodies to allow further discrimination of specific cell subpopulations.

The QuantiGene FlowRNA Assay is based on branched DNA technology to amplify the signal of a RNA target. The first hybridization of the assay allows the oligonucleotide probe set to bind to the target RNA sequence. Subsequent signal amplification is predicated on specific hybridization of adjacent Probe Set pairs. A typical Probe Set contains on average 20 oligonucleotide pairs.

Next signal amplification is achieved via a series of sequential hybridization steps. The PreAmplifier molecules hybridize to their respective pair of bound Probe Set oligonucleotides, then multiple Amplifier molecules hybridize to their respective PreAmplifier. Finally, Label Probe oligonucleotides conjugated to fluorescent dye hybridize to their corresponding Amplifier molecules. A fully assembled signal amplification "tree" has 400 Label Probe binding sites. When all target-specific oligos in the Probe Set bind to the target mRNA transcript, an 8,000 fold amplification is achieved. Affymetrix currently offers three different amplification "tree" structures that allow simultaneous measurement of up to three different RNA targets for multicolor flow cytometric analysis.

Once the cells have been processed by the QuantiGene FlowRNA Assay, the data can be collected and analyzed on a flow cytometer.



QuantiGene FlowRNA Reagent System Contents and Storage Conditions

The components of the QuantiGene FlowRNA Assay Kit and their recommended storage conditions are listed below. Refer to the package insert for quantities of individual components supplied. Kit components have a shelf life of six months from the date of receipt.

Each QuantiGene FlowRNA Assay Kit is supplied in three separate parts, based on storage temperature.

Box 1: Shipped on dry ice and store at -20°C

Component	Description	Storage
100x Label Probe Mix	Fluorescent dye-labeled oligonucleotides in aqueous buffered solution containing Alexa 647-LP1 (uses APC channel), Alexa 488-LP4 (uses FITC channel) and Alexa 750-LP6 (uses APC-Cy7 channel)	-20°C

Box 2: Shipped on blue ice and store at 2-8°C

Component	Description	Storage
Target Probe Diluent	Aqueous solution containing formamide and detergent	2-8°C
PreAmp Mix	DNA (PreAmp1, PreAmp4 and PreAmp6) in aqueous solution containing formamide and detergent	2-8°C
Amp Mix	DNA (Amp1, Amp4 and Amp6) in aqueous solution containing formamide and detergent	2-8°C
Label Probe Diluent	Aqueous solution containing detergent	2-8°C
Wash Buffer	Aqueous buffered solution	2-8°C
Storage Buffer	Aqueous buffered solution	2-8°C
1.5-mL microfuge tube	Clear copolymer tube with a flat cap (RNase free) for low cell binding. DO NOT SUBSTITUTE.	Room Temp

Box 3: Fixation and Permeabilization Kit: Shipped at room temperature and store at 2-8°C

Component	Description	Storage
Solution A	Aqueous buffered solution	2-8°C
Solution B	Aqueous buffered solution	2-8°C
Solution C	Aqueous buffered solution containing detergent	2-8°C
Solution D	Aqueous buffered solution	2-8°C

QuantiGene FlowRNA Probe Sets

In addition to the QuantiGene FlowRNA Assay Kit, Probe Set(s) specific to your RNA target(s) of interest must be purchased separately. Refer to the Product Insert for more information. Our current probe set catalog can be found at <http://www.ebioscience.com/application/flowrna.htm>.

Component	Description	Storage
QuantiGene FlowRNA Probe sets	RNA-specific oligonucleotides designed against a target of interest and are compatible with TYPE 1, 4 or 6 Signal Amplifiers.	-20°C

▲ WARNING: All chemicals should be considered potentially hazardous. We recommend that this product and its components be handled by those trained in laboratory techniques and be used according to the principles of good laboratory practice.

Required Equipment and Materials Not Supplied

Required Equipment/Material	Source	Part Number
16% Paraformaldehyde (PFA) aqueous solution (in 10-mL amber tube)	Electron Microscopy Sciences	15710
Methanol (HPLC-grade, 99.9% pure)	Fisher Scientific	A452-4
QuantiGene ViewRNA Temperature Validation Kit	Affymetrix	QV0523
Flow cytometer with the following components: - Two lasers, blue (488 nm) and red (633 nm or similar) - Detection channels optimized for FITC, APC, and APC-Cy7 fluorophores	MLS ¹ (e.g. BD Biosciences, Beckman Coulter...)	Many (e.g. Calibur from BD Biosciences, FC500 from Beckman Coulter...)
Refrigerated low-speed centrifuge with swinging bucket for 15-mL tubes with inserts for 1.5-mL microfuge tubes	MLS, (e.g. Eppendorf)	Eppendorf 5804 with A-4-44 Rotor (swing bucket)
Incubator (validated to maintain 40 ± 1 °C)	Affymetrix	QS0704 (120V) QS0712 (220V)
Metal tube rack (heat block) for 1.5-mL microfuge tubes (placed inside incubator)	MLS (e.g. VWR, Fisher Scientific...)	VWR 13259-002
Vacuum aspirator, adjustable to 0.5 mL/ sec	MLS, (e.g. Fisher Scientific, VWR...)	Fisher Scientific 07200564
Vortexer	MLS, (e.g. VWR, Fisher Scientific...)	VWR 58816-121
(Optional) 15 mL polypropylene tubes	MLS, (e.g. Fisher Scientific, BD...)	Fisher Scientific 14-959-70C
(Optional) RNasin	Promega	N2611

¹Multiple Lab Supplier

Experimental Guidelines

1. Start with cell culture in good physiological conditions.

Cells should be in active growth phase to preserve RNA integrity and minimize cell lysis in processing. If using freshly isolated primary cells, make sure the cells are healthy after purification.

2. Selection of probe type and fluorescent channel for optimal RNA detection

Type 1 /650 probe set produces the most sensitive detection when compared to Type 4 /488 and Type 6/ 750 probe sets. Type 1/650 channel is about two to five folds more sensitive than the other two channels. We recommend using the following probe type/excitation line for the low, medium and high expression genes, respectively. In addition, background fluorescence in the Type 4/488 channel is higher compared to Type 1/650 and Type 6/750 channels.

Gene Expression Level	Probe Type	Cytometer Laser line	Ex/Em
Low	Type 1	633	650/668
Medium to High	Type 4	488	495/519
Medium to High	Type 6	633	749/775

3. QuantiGene FlowRNA Assay Controls

Below are recommended assay controls that we recommend for first time users to perform.

Control	Gene	Catalog Number for Probe Set
Positive Control	B2M - high expression in PBMC	Human (Type 1) VA1-10611 (Type 4) VA4-13460 (Type 6) VA6-11782 Mouse (Type 1) VB1-11142 (Type 4) VB4-14643 (Type 6) VB6-12836
	GAPDH - medium expression in PBMC	Human (Type 1) VA1-10119 (Type 4) VA4-10641 (Type 6) VA6-10337 Mouse (Type 1) VB1-10150 (Type 4) VB4-10414 (Type 6) VB6-10574
Negative assay control	No Target Probe	Not Applicable
	Bacterial DapB gene	(Type 1) VF1-11712 (Type 4) VF4-10408 (Type 6) VF6-10407

Assay Procedure

This protocol contains the following sections. It can be performed in two days with optional stop points.

Day 1:

I. Fixation and Permeabilization

II. Target Probe Hybridization

Day 2:

III. Signal Amplification and Detection

IMPORTANT: This protocol is optimized for the swinging bucket centrifuge. We do not recommend using a fixed-angle centrifuge as severe cell loss may result.

I. Fixation and Permeabilization

NOTE: The fixation and permeabilization can be performed in bulk if desired. However, polypropylene conical tubes (e.g. Fisher Scientific cat. #14-959-70C) must be used. DO NOT use polystyrene tubes such as FACS tube (e.g. BD Falcon cat. #352052) which would lead to low signal.

Step	Action
Step 1.	(Optional) Antibody staining - see Appendix for details.
Step 2.	<p>Perform the following tasks to prepare enough reagents for 24 samples.</p> <ol style="list-style-type: none"> Pre-warm Target Probe Diluent to 40°C. Equilibrate Solutions A, B, C, and D to room temperature (RT). Prepare a 4% paraformaldehyde (PFA) in Solution A. <p>NOTE: PFA solution must be freshly prepared before use.</p> <ul style="list-style-type: none"> ■ Add 30 mL of Solution A to a 50 mL conical tube and keep at RT in the dark. ■ Open an amber tube containing 10 mL of 16% PFA stock solution and add the solution to the conical tube containing Solution A. Mix well and keep at RT in the dark. <ol style="list-style-type: none"> Ensure incubator and heat block are calibrated to 40°C using a traceable digital thermometer (see QuantiGene ViewRNA Temperature Validation Kit for details). Pre-chill methanol to -20°C. Set the temperature of centrifuge to 4°C.
Step 3.	Hybridization must be performed in the 1.5-mL microfuge tubes provided. For each reaction, start by transferring 1-3 million cells into each tube.
Step 4.	<p>Centrifuge the tubes in a swinging bucket at 400 x g at 4°C for 5 min to pellet the cells. Aspirate the supernatant leaving a residual volume of 100 µL. Vortex tubes to resuspend cells evenly. Add 1 mL Solution A, then invert tubes a few times to mix.</p> <p>NOTE: Always leave a residual volume of 100 µL when aspirating in subsequent steps to avoid disturbing the cell pellet. The lowest mark within the graduations on the microfuge tube is 100 µL.</p>
Step 5.	Centrifuge cells at 400 x g at 4°C for 5 min and aspirate the supernatant leaving a residual volume of 100 µL.
Step 6.	Vortex tubes to resuspend cells evenly. To each tube, add 1 mL of 4% PFA in Solution A prepared in Step 2c. Briefly invert the tubes to mix and incubate at RT for 60 min in the dark, inverting the tubes 1-2 times during incubation.
Step 7.	Centrifuge the tubes at 800 x g at 4°C for 5 min and aspirate the supernatant leaving a residual volume of 100 µL. Briefly vortex tubes to resuspend the cell pellets.

Step	Action
Step 8.	To each tube, add 1 mL of Solution B and invert the tubes to mix. Centrifuge cells at 800 x g at 4°C for 5 min and aspirate the supernatant leaving a residual volume of 100 µL. Briefly vortex tubes to resuspend the cell pellets.
Step 9.	To each tube, add 1 mL of ice-cold methanol solution. Briefly vortex the tubes to mix and incubate on ice for 10 min. NOTE: Optional stop point: For storage, transfer cells to -80°C after adding methanol. Samples should be stable for 3 months. However, we do not recommend sample storage in methanol if antibody staining has been performed as in Step 1.
Step 10.	Perform the following tasks: NOTE: PFA solution must be freshly prepared before use. a. Prepare 4% PFA in Solution C. <ul style="list-style-type: none"> ■ Add 30 mL of Solution C to a 50-mL conical tube. Keep at RT in the dark. ■ Open an amber tube containing 10 mL of 16% PFA stock solution and add the solution to the conical tube containing Solution C. Mix well and keep at RT in the dark. b. Prepare a 1:1 solution of methanol and 4% PFA in Solution C <ul style="list-style-type: none"> ■ Add 13 mL of ice-cold methanol to a 50-mL conical tube. ■ Add 13 mL of 4% PFA in Solution C (from Step 10a). Keep on ice in dark. c. Thaw the target probe sets at RT and then place on ice.
Step 11.	Centrifuge tubes at 800 x g at 4°C for 5 min and aspirate the supernatant leaving a residual volume of 100 µL. Briefly vortex tubes to resuspend the cell pellets.
Step 12.	To each tube, add 1 mL of ice-cold methanol/ 4% PFA in Solution C (1:1) prepared in Step 10b. Invert the tubes to mix and incubate on ice for 5 min.
Step 13.	Centrifuge tubes at 800 x g at 4°C for 5 min and aspirate the supernatant leaving a residual volume of 100 µL. Briefly vortex tubes to resuspend cell pellet.
Step 14.	To each tube, add 1 mL of freshly prepared 4% PFA solution in Solution C prepared in Step 10a. Invert the tubes to mix and incubate at RT for 30 min. Invert tubes once during incubation.
Step 15.	Centrifuge tubes at 800 x g at RT for 5 min and aspirate the supernatant leaving a residual volume of 100 µL. Briefly vortex tubes to resuspend the cell pellets.
Step 16.	To each tube, add 1 mL of Solution D. Invert tubes to mix and centrifuge tubes at 800 x g at RT for 5 min, then aspirate the supernatant leaving a residual volume of 100 µL. Briefly vortex tubes to resuspend the cell pellets.
Step 17.	Repeat wash with Solution D and leave a residual volume of 100 µL.

II. Target Probe Hybridization

Step	Action
Step 1.	Dilute the target probes 1/20 with Target Probe Diluent. Add 100 µL of diluted target probes to each tube from Step 17.
Step 2.	Briefly vortex tubes to mix. Place tubes in the pre-warmed heat block inside the incubator. Incubate at 40°C for 2 hours. Invert tubes to mix after 1 hour of incubation.
Step 3.	To each tube, add 1 mL of Wash Buffer and invert tubes to mix. Centrifuge the tubes at 800 x g at RT for 5 min to pellet the cells, then aspirate leaving a residual volume of 100 µL. Briefly vortex tubes to resuspend the cell pellets.

Step	Action
Step 4.	<p>Repeat wash by adding 1 mL of Wash Buffer (Optional: add 1 μL/mL of RNasin if used as stop point) and invert tubes to mix. Centrifuge the tubes at 800 x <i>g</i> at RT for 5 min to pellet the cells, then aspirate leaving a residual volume of 100 μL. Briefly vortex tubes to resuspend the cell pellets.</p> <hr/> <p>NOTE: Optional stop point: store samples overnight at 4°C in the dark.</p> <hr/>

III. Signal Amplification and Detection

Step																					
Step 1.	Perform the following tasks: a. Pre-warm PreAmp Mix, Amp Mix and Label Probe Diluent to 40°C. b. Equilibrate Wash Buffer, Storage Buffer and the samples to RT.																				
Step 2.	Add 100 µL of pre-warmed PreAmp Mix (with green cap) to each tube and briefly vortex tubes to mix. Place tubes in the pre-warmed heat block inside the incubator. Incubate at 40°C for 1.5 hours.																				
Step 3.	Remove tubes from incubator and add 1 mL of Wash Buffer to each tube. Invert tubes to mix. Centrifuge the tubes at 800 x g at RT for 5 min to pellet the cells, then aspirate leaving a residual volume of 100 µL. Briefly vortex tubes to resuspend the cell pellets.																				
Step 4.	Repeat wash.																				
Step 5.	Add 100 µL of pre-warmed Amp Mix to each tube and briefly vortex tubes to mix. Place tubes in the pre-warmed heat block inside the incubator. Incubate at 40°C for 1.5 hours.																				
Step 6.	During Amp Mix incubation, thaw Label Probe Mix on ice.																				
Step 7.	Dilute Label Probe Mix 1/100 with Label Probe Diluent. Store in the 40°C incubator and protect from light until ready to use.																				
	<table border="1"> <thead> <tr> <th></th> <th>Reagent</th> <th>1 Tube</th> <th>10 Tubes</th> <th>24 Tubes</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Label Probe</td> <td>1.2 µL</td> <td>12 µL</td> <td>28.8 µL</td> </tr> <tr> <td>2</td> <td>Label Probe Diluent</td> <td>118.8 µL</td> <td>1188 µL</td> <td>2851.2 µL</td> </tr> <tr> <td>Total</td> <td></td> <td>120 µL</td> <td>1200 µL</td> <td>2880 µL</td> </tr> </tbody> </table>		Reagent	1 Tube	10 Tubes	24 Tubes	1	Label Probe	1.2 µL	12 µL	28.8 µL	2	Label Probe Diluent	118.8 µL	1188 µL	2851.2 µL	Total		120 µL	1200 µL	2880 µL
	Reagent	1 Tube	10 Tubes	24 Tubes																	
1	Label Probe	1.2 µL	12 µL	28.8 µL																	
2	Label Probe Diluent	118.8 µL	1188 µL	2851.2 µL																	
Total		120 µL	1200 µL	2880 µL																	
Step 8.	Remove tubes from incubator and add 1 mL of Wash Buffer to each tube. Invert tubes to mix. Centrifuge the tubes at 800 x g at RT for 5 min to pellet the cells then aspirate the supernatant leaving a residual volume of 100 µL. Briefly vortex tubes to resuspend the cell pellets.																				
Step 9.	Repeat wash.																				
Step 10.	Add 100 µL of Working Label Probe Solution prepared in Step 7 to each tube and briefly vortex tube to mix. Place tubes in the pre-warmed heat block inside the incubator. Incubate at 40°C for 1 hour.																				
Step 11.	Turn on the flow cytometer to allow the lasers to warm up.																				
Step 12.	Remove tubes from incubator and add 1 mL of Wash Buffer to each tube. Invert tubes to mix. Centrifuge the tubes at 800 x g at RT for 5 min to pellet the cells, then aspirate leaving a residual volume of 100 µL. Briefly vortex tubes to resuspend the cell pellets.																				
Step 13.	Repeat wash.																				
Step 14.	Add 1 mL of Storage Buffer to each tube and invert tubes to mix. Centrifuge the tubes at 800 x g at RT for 5 min to pellet the cells, then aspirate supernatant leaving 100 µL residual volume. Briefly vortex tubes to resuspend the cell pellets. Protect samples from light. NOTE: Optional stop point: store samples at 4°C in the dark.																				
Step 15.	To analyze a sample, pipette to mix. Then transfer 50-100 µL sample to a reading tube containing 250 µL of 1X PBS and analyze on a flow cytometer. For storage, keep samples at 4°C in the dark. Samples are stable for 1 week when stored at 4°C.																				

Appendix

A1: Troubleshooting






Observation	Probable Cause	Recommended Solution
Poor Cell Recovery / Low Cell Count on Cytometer	Cell lysis due to cells in poor physiological condition during processing	<ul style="list-style-type: none"> ■ If using frozen cells, thaw cells carefully following proper cell culture procedures. ■ Check cell viability using trypan blue or other viability dyes before assay.
	Improper centrifuge and tube adaptor	<ul style="list-style-type: none"> ■ Use swinging bucket for centrifugation. Do not use fixed-angled centrifuges. ■ Use correct tube adaptor for proper <i>g</i>-force.
	Improper centrifuge settings	<ul style="list-style-type: none"> ■ Confirm settings are for RCF (<i>xg</i>), not RPM. ■ Set centrifuge to the speed specified in the protocol.
	Improper tubes used	Only use the 1.5-mL microfuge tubes provided with the kit. Tubes from other vendors may result in significant cell loss.
	High aspiration rate by pipetting or vacuum aspirator when removing supernatant	Ensure vacuum setting for aspirator is low. Aspirate no faster than 0.5 mL/sec. Place a 1-10 µL pipette tip on the tip of the aspirator to help reduce aspiration rate. Follow the meniscus while aspirating down to the 100 µL mark on side of tube.
	Rough handling	Avoid excessive vortexing to prevent cell damage. Vortex in pulses. When pipetting to mix, do so gently.
	Clog	Ensure cytometer is not clogging due to clumps of aggregated cells. Try using gentle vortexing and pipette mixing to break up aggregates, or use a 70 µL filter.
	Incorrect cell count by flow cytometer	Verify the flow rate of cytometer using beads of known concentration such as the Flow Cytometry Absolute Count Standard (Bangs Laboratories, cat. # 580).
False Positive Population	Compensation not set up correctly	See Appendix section on how to set up compensation.
Low / No Signal	Cytometer not set up properly	Check that your instrument is properly set up according to the manufacturer's recommendations.
	Old 4% PFA solution was used in fixation.	Prepare fresh 4% PFA solution before use.

Observation	Probable Cause	Recommended Solution								
	Polystyrene tube (e.g. FACS reading tube) was used in fixation and permeabilization.	Use polypropylene tube (e.g. Fisher Scientific 14-959-70C) for fixation and permeabilization.								
	Insufficient number of cells recorded on cytometer	<p>Determine the number of total events to collect to ensure at least 100-1000 positive cells are recorded, for example:</p> <p>Total number of cells to count = number of positives desired / frequency of positives (%)</p> <p>To count 1,000 positive cells,</p> <table border="1"> <thead> <tr> <th>Frequency of Positives</th> <th>Total # of cells to count</th> </tr> </thead> <tbody> <tr> <td>10%</td> <td>10,000</td> </tr> <tr> <td>5%</td> <td>20,000</td> </tr> <tr> <td>1%</td> <td>100,000</td> </tr> </tbody> </table>	Frequency of Positives	Total # of cells to count	10%	10,000	5%	20,000	1%	100,000
Frequency of Positives	Total # of cells to count									
10%	10,000									
5%	20,000									
1%	100,000									
	Incorrect preparation of signal amplification reagents	Make sure PreAmp, Amp and Label Probe are used in the specified order.								
	Incorrect incubator temperature	<ul style="list-style-type: none"> ■ Incubator must be able to hold temperature at 40±1°C. ■ Use 1.5-mL heat block for hybridization. ■ Minimize traffic to incubator. Temperature tolerance is 1-2°C. ■ Ensure incubator has been stable at 40°C for at least 12 hours before starting assay as temperature may drift during initial incubator setup. ■ Measure temperature using QG ViewRNA Temperature Validation Kit from Affymetrix (PN QV0523). 								
	Incorrect Wash Buffer	Do not substitute Wash Buffer with Storage Buffer.								
	Incorrect diluents used	Ensure Target Probe Diluent is used with Target Probe, and Label Probe Diluent used with Label Probe.								
	Incorrect dilution of Target Probe/ Label Probe	<ul style="list-style-type: none"> ■ Ensure Target Probe was added to Target Probe Diluent at a 1/20 dilution. ■ Ensure Label Probe Mix was added to Label Probe Diluent at a 1/100 dilution. 								

Observation	Probable Cause	Recommended Solution
	Gene expressed at very low level	<ul style="list-style-type: none"> ■ Check biological model and confirm gene is expressed in sample by running QuantiGene 2.0 Assay for lysate. ■ Keep in mind that protein expression might not correlate with RNA expression in some cases.
High Background	Incorrect incubator temperature	<ul style="list-style-type: none"> ■ Incubator must be able to hold temperature at 40±1°C. ■ Minimize traffic to incubator. Temperature tolerance is 1-2°C. ■ Ensure incubator has been stable at 40°C for at least 12 hours before starting assay as temperature may drift during initial incubator setup. ■ Measure temperature using QG ViewRNA Temperature Validation Kit from Affymetrix (cat. #QV0523).
	Insufficient washing	<ul style="list-style-type: none"> ■ Ensure the Wash Buffer is used at room temperature. ■ Follow instructions of washing steps in the protocol. ■ Ensure uniform cell resuspension by vortexing.
	Cytometer not set up properly	Reduce voltage settings of your cytometer. Note that this assay may require lower settings than typical antibody staining.
	Over-fixation	Follow the protocol for proper fixation time.
	Excessive Target Probe used	Dilute Target Probe at a 1/20 dilution.

A2. Lab Setup Guide

This guide illustrates the setup of typical equipment and their specifications for the QuantiGene FlowRNA Assay. Consult your equipment manufacturers to make sure the equipment meets the specifications. The major equipment include: incubator, swinging bucket centrifuge, flow cytometer, aspiration system for washing, and temperature validation kit.

Equipment	Specifications
<p data-bbox="326 478 667 506">Incubator with heat block inside</p>  <p data-bbox="467 695 553 716">Metal heat block</p>	<ul style="list-style-type: none"> ■ Validated to maintain 40±1°C. ■ Metal block for 1.5-mL microfuge tube. ■ Example shown: Incubator (Affymetrix QS0704) Metal heat block (VWR 13259-002)
<p data-bbox="326 760 753 787">Refrigerator swinging bucket centrifuge</p>  <p data-bbox="553 932 672 953">Eppendorf 5810R</p> <p data-bbox="440 1079 586 1115">Rotor with adaptor for 15-mL conical tube</p> <p data-bbox="634 1079 753 1115">Adaptor for 1.5-mL microfuge tube</p>	<ul style="list-style-type: none"> ■ Swinging bucket with adaptors for 15-mL conical tubes and 1.5-mL microfuge tubes. ■ With refrigeration to 4°C. ■ Examples shown: Centrifuge (Eppendorf 5810R) Rotor (Eppendorf A-4-44) Adaptor for 15-mL conical tube (Eppendorf 5804 755.006) Adaptor for 1.5-mL microfuge tube (Eppendorf 5804 750.004) NOTE: Fixed-angle centrifuge is NOT recommended.
<p data-bbox="326 1144 493 1171">Flow cytometer</p> 	<ul style="list-style-type: none"> ■ Two lasers - blue (488 nm) and red (633 nm or similar) ■ Detection optics optimized for FITC, APC, and APC-Cy7 fluorophores ■ Example shown: BD LSRFortessa
<p data-bbox="326 1417 646 1444">Aspiration system for washing</p>  <p data-bbox="623 1654 699 1675">Aspiration system</p>	<ul style="list-style-type: none"> ■ Aspiration rate adjusted to 0.5 mL / sec. ■ Can use in-house vacuum line or vacuum pump. ■ Example shown: Vacuum bottle (Argos Technologies EV432) Aspirator (Argos Technologies EV514)
<p data-bbox="326 1711 857 1738">QuantiGene ViewRNA Temperature Validation Kit</p> 	<ul style="list-style-type: none"> ■ NIST traceable thermometer with temperature probe in 1.5-mL microfuge tube (shown on left). ■ Affymetrix, cat.# QV0523


A3. Cytometer Setup

The QuantiGene FlowRNA Assay Kit utilizes up to three fluorescent channels for detection of RNA on a flow cytometer. To ensure optimal detection of RNA, it is important that the cytometer is set up properly. In a multi-color assay, signal from a fluorophore often spills over into the other detection channels due to overlapping emission spectra. The process of subtracting the spillover signal is called compensation. The procedures below describe how to set proper PMT voltage and compensation for multi-color RNA analysis.

Reagents and samples needed
Calibration beads
Unstained cells
Single-color samples: Type 1 probe (Alexa Fluor® 647)-only labeled cells Type 4 probe (Alexa Fluor® 488)-only labeled cells Type 6 probe (Alexa Fluor® 750)-only labeled cells Single-color antibody-labeled cells (needed only if protein staining is performed)

1. Ensure proper alignment following the instrument manufacturer's recommendations. Alternatively, calibration beads, such as Spherotech Rainbow Calibration Particles (8-peak beads, catalog number RCP-30-5A) or Rainbow Fluorescent Particles (1-peak beads, catalog number RFP-30-5A [1 peak]) may be used to check the linearity of the PMT, PMT voltage range, and alignment of the instrument.
2. Setting PMT Voltages
 - A. Create an acquisition template/workspace that includes a forward scatter versus side scatter plot, one single-parameter histogram plot for every fluorescent channel being used, and a two-parameter plot for every pairwise combination of fluorescent channels being used in the experiment.
 - B. Place the unstained cell sample on the cytometer and begin to collect events.
 - 1) Adjust voltages for forward scatter and side scatter so that cells/events of interest are on scale. Set a threshold on forward scatter to eliminate debris from acquisition.
 - 2) Looking at the histogram plots, ensure that the signals are on scale and adjust PMT voltages, if necessary.
 - 3) Stop collecting events and remove sample from the cytometer.
 - C. Place a single-color sample on the cytometer and begin to collect events.
 - 1) Looking at the histogram plot for the fluorochrome being collected, ensure that events are on scale for the fluorescent channel being collected and adjust PMT voltages, if necessary.
 - 2) Stop collecting events and remove sample from the cytometer.
 - 3) Repeat for each of the single-color samples.
 - D. Confirm that the voltage settings are appropriate for the experimental samples.
 - 1) Place an experimental sample on the cytometer and begin collecting events.
 - 2) Looking at each of the histogram plots, confirm that all fluorescent signals are on scale and adjust PMT voltages, if necessary.
 - 3) Stop collecting events and remove sample from the cytometer.
3. Setting compensation
 - A. If an autocompensation feature is available for the cytometer being used, follow the manufacturer's instructions for use.
 - B. If an autocompensation feature is not available or to set compensation manually:

- 1) Place the unstained sample on the cytometer and begin to collect events.
- 2) Looking at the two-parameter plots, create quadrant regions on each plot such that the cells fall within the lower left quadrant (see Figure 1).



- 3) For each two-parameter plot, create statistics windows/tables according to your instrument's acquisition software and include the median for the x- and y-parameters.
 - 4) Stop collecting events and remove sample from the cytometer.
- C. Place a single-color sample on the cytometer and begin to collect events. NOTE: If your single-color sample does not contain any negative events, a small amount of unstained cells may be added to each of the single-color samples just prior to loading the sample onto the cytometer. Alternatively, the MFI of the unstained sample run in Step 3B above must be noted for all parameters.
- 1) Looking at the two-parameter plots for the fluorochrome being collected, adjust the quadrants so that the positive events are fully contained within either the lower right (if the parameter is on the x-axis; Figure 2A, left) or the upper left (if the parameter is on the y-axis; Figure 2B, left) quadrant.
 - 2) From the statistics window, note the y-medians (if the parameter is on the x-axis; Figure 2A, left) or the x-medians (if the parameter is on the y-axis Figure 2B, left) of the two populations.
 - 3) Adjust each of the compensation values for the fluorochrome being collected (subtract its fluorescence out of each of the other channels being used in the experiment) until the medians are the same (Figure 2A and 2B, right).
 - 4) Stop collecting events and remove sample from the cytometer.

5) Repeat for each of the single-color samples.

Figure 2A. CD8 APC-only sample





Figure 2B. CD4 APC-eFluor 780-only sample



D. Confirm that the compensation settings are appropriate for the experimental samples.

- 1) Place an experimental sample on the cytometer and begin collecting events.
- 2) Looking at each of the two-parameter plots, confirm that the compensation is set correctly by looking at the medians.
- 3) Stop collecting events and remove sample from the cytometer.

A4. Incorporating Antibody Staining with the QG FlowRNA Assay

Fluorochrome-conjugated antibodies can be used to stain surface proteins for the purpose of immunophenotyping cells that will be further analyzed for RNA levels using the QuantiGene FlowRNA Assay. Due to the nature of the RNA labeling protocol, additional considerations for staining are warranted.

The following two options are recommended for staining with antibodies:

1. Cells may be stained with fluorochrome-conjugated antibodies before the Fixation and Permeabilization Procedure only if the fluorochromes are resistant to exposure to methanol. Methanol-resistant and methanol-sensitive fluorochromes are listed in the table below. Although some fluorochromes are resistant to methanol, there may still be some loss in signal intensity, therefore, performance of a specific reagent in this assay should be determined empirically.

2. If the antibody is not conjugated to a methanol-resistant fluorochrome, it may be possible to utilize an indirect staining protocol. Cells may be stained with an unconjugated antibody before the Fixation and Permeabilization Procedure. A secondary reagent (fluorochrome-conjugated anti-Ig secondary antibody) can then be used to detect the primary antibody. Staining with the secondary reagent must be done after Signal Amplification and Detection. It is important to note that only one unconjugated antibody from a given host may be used unless the isotypes of the unconjugated antibodies used are different. In this case, isotype-specific secondaries must be used; a pan anti-mouse IgG secondary antibody will bind to all mouse IgG antibodies present on the cells.

Methanol-Resistant and Methanol-Sensitive Fluorochromes

Methanol Resistant Fluorochromes	Methanol Sensitive Fluorochromes
FITC	PE
Alexa Fluor® 488	PE-Cy5
Alexa Fluor® 700	PE-Cy5.5
Alexa Fluor® 647	PE-Cy7
eFluor® 450	PerCP
Pacific Blue™	PerCP-Cy5.5
Pacific Orange™	PerCP-eFluor 710
Brilliant Violet 421™	APC
Brilliant Violet 605™	APC-Cy7
	APC-H7
	APC-Alexa Fluor® 750
	APC-eFluor® 780
	Qdot® reagents
	eFluor® Nanocrystal (NC) reagents

How to validate an antibody for QuantiGene FlowRNA Assay:

Each antibody needs to be validated with the assay protocol to confirm its compatibility.

Experimental conditions:

1. Antibody staining alone

2. Antibody staining followed by RNA detection procedure without target probe

Analyze and compare two samples for:

1. Signal intensity of the protein staining
2. Percentage of positive cells

For a qualified antibody:

1. Sufficient separation between the positive and negative signals should be maintained after the RNA detection procedure to allow identification of desired populations
2. The percentage of desired cell populations should be comparable.


A5. Expected Results

This is an example of typical results of the QuantiGene FlowRNA Assay.

Experimental conditions

Human PBMC were stimulated with 1 $\mu\text{g}/\text{mL}$ of bacterial lipopolysaccharides (LPS) and 2.5 $\mu\text{g}/\text{mL}$ of R-848 for 4 hours. Cells were subjected to analysis by the QuantiGene FlowRNA Assay. Cells were stained with anti-human CD14-FITC to identify monocytes and labeled with probes for the RNA expression of IL-1beta, IL-6, IL-8 and TNF alpha. Cells in the monocyte gate were used for analysis. CD14 staining is shown on the y-axis and RNA labels shown on the x-axis.

Result



A6. Validated Cell Lines

The QuantiGene FlowRNA Assay has been validated with the following cell types:

Primary cells

- Human PBMC (cryopreserved and freshly purified by Ficoll)
- Human T cells (CD3+)*
- Human B cells (CD19+)*
- Human monocytes (CD14+)*
- Mouse bone marrow cells

*Cells were purified by antibody-coated magnetic beads.


Cultured suspension cell lines

- NKT, human natural killer cell line
- Jurkat, human lymphocytic cell line
- U937, human monocytic cell line
- K562, human myelocytic cell line
- M1, mouse myeloid cell line

A7. Temperature Validation Procedure for Incubator

Temperature control is critical for the success of the QuantiGene FlowRNA Assay. The incubator should be validated before use. Improper hybridization temperature will result in high background and/or weak signal.



Required materials

Item	Description
Incubator	Capable of maintaining temperature at $40 \pm 1^\circ\text{C}$, e.g. Affymetrix cat. # QS0704 or QS0712
QuantiGene View Temperature Validation Kit (Affymetrix cat. # QV0523)	NIST calibrated digital thermometer with Type-K beaded probe
1.5-mL microfuge tube	1.5-mL microfuge tube used for assay (included in the QuantiGene FlowRNA Assay Kit)
Metal heat block for 1.5-mL microfuge tube	Metal block to hold the microfuge tube, e.g. Fisher Scientific 11-718-9Q 
Parafilm	Major Laboratory Supplier

Procedure

Step	Action
Step 1 Prepare the incubator	A. Turn on the incubator. B. Set the temperature to 40°C. C. Put the metal heat block into the incubator near the center of the middle incubator. D. Allow the incubator and heat block to equilibrate for overnight.
Step 2 Assemble unit	A. Insert the battery to activate the digital thermometer. B. Use a pointed object, e.g. a ballpoint pen to drill a hole at the lid of the 1.5-mL microfuge tube, along the circular markings. C. Add 0.2 mL deionized water to the tube, and then close the lid. D. Insert the Type-K beaded probe into the digital thermometer, and the other end into the pre-drilled hole of the provided 1.5-mL tube containing 0.2 mL of deionized water. E. Wrap parafilm around the top of the 1.5-mL tube and probe to form a seal. Avoid an excessive amount of parafilm around the sides of the tube; otherwise it may not fit properly into the heat block. F. Turn on the digital thermometer device



Step	Action
<p>Step 3 Measure and adjust incubator temperature</p>	<p>A. Insert the 1.5-mL tube with probe and deionized water into the heat block that has been pre-warmed in the incubator in Step 1.</p>  <p>B. Close the door making sure there is sufficient slack in the wiring. C. Wait 15 min for temperature to equilibrate. D. Record the temperature. E. If necessary, adjust the dry incubator temperature settings so that the digital thermometer reads 40°C. After adjustment, allow the incubator and heat block to equilibrate. Then recheck the temperature. F. Repeat Step E as needed.</p> <p>IMPORTANT: We recommend calibrating the incubator at least once a month to ensure accuracy.</p>
<p>Step 4 Assess incubator temperature uniformity</p>	<p>Repeat step 3 to measure the temperature at multiple positions in the incubator to determine temperature uniformity.</p>  <p>IMPORTANT: The temperature for all positions should be $40 \pm 1^\circ\text{C}$.</p>
<p>Step 5 Assess ramp-up time</p>	<p>A. Open the incubator door for 1 min then close the door. Measure the time needed for the temperature to return to 40°C and monitor the temperature profile during recovery. B. Repeat the Step 5A two more times.</p> <p>IMPORTANT: Do not use the incubator for the assay if it takes more than 5 minutes to return to 40°C or if it overshoots by more than 2°C during recovery.</p>

