


SQUID 2

Rev. 2.50

A User's Manual



Ken Ludwig
Berkeley Geochronology Center
April 12, 2009

Funded by:

Geoscience Australia, U.S. Geological Survey, and Berkeley Geochronology Center

together with

*All-Russian Institute of Geological Research, Australian National University,
Australian Scientific Instruments, Curtin University, Geological Survey of Canada, Na-
tional Institute of Polar Research, Stanford University*

To reference: Ludwig, K. (2009), SQUID 2: A User's Manual, rev. 12 Apr, 2009. Berkeley Geochron. Ctr. Spec. Pub. **5** 110 p.

INTRODUCTION	2
PURPOSE OF SQUID-2	2
REQUIREMENTS	2
<i>Software</i>	2
<i>Hardware</i>	2
<i>Printer Setup</i>	3
<i>Input Data Files</i>	3
INSTALLATION.....	3
DATA REDUCTION FOR U-PB GEOCHRONOLOGY	4
U-Pb/Th-Pb STANDARDS	5
U-Pb/Th-Pb STANDARDS	5
<i>Specifying standards when starting data reduction</i>	5
<i>The Common U-Pb Standards List</i>	5
SPECIFYING THE COMMON-Pb INDEX ISOTOPE (U-Pb GEOCHRON PANEL)	6
SBM NORMALIZATION	6
TRIM MASS CHARTS.....	6
CONDENSED RAW-DATA SHEETS PRODUCED BY SQUID.....	7
REDUCED-DATA SHEETS PRODUCED BY SQUID	8
<i>The StandardData Sheet</i>	8
<i>The SampleData Sheet</i>	9
SAMPLE GROUPING AND AGE EXTRACTION	9
CORRECTING U-Pb GEOCHRONOLOGY DATA FOR SECULAR DRIFT OF AGE-STANDARD SPOTS.....	10
<i>Caveats in using secular drift correction</i>	12
REDUCING GENERAL ISOTOPE DATA.....	13
TASKS.....	14
INTRODUCTION	14
THE <u>TASK NAME</u> PANEL.....	16
THE <u>RUN TABLE</u> PANEL.....	16
<i>Defining the Run Table</i>	16
<i>Specifying Isotope Ratios</i>	19
<i>Specifying Isotope Ratios</i>	19
THE <u>EQUATION</u> PANELS.....	20
<i>Basics</i>	20
<i>Restricting equations to subsets of spots</i>	26
<i>“Special U-Pb” Equations</i>	27
<i>Swapping data columns with U-Pb Geochronology Tasks</i>	28
<i>Referencing data in other workbooks or worksheets</i>	29
<i>Viewing Task Equations in the reduced-data workbook</i>	30
<i>Equation Tips</i>	30
THE AUTOCHARTS PANEL	30
SETTING A “TIME WINDOW” FOR SPOTS	31
SPECIFYING SUBSET-RESTRICTED EQUATIONS AT THE TIME OF DATA REDUCTION.....	32
INVOKING EXCEL’S SOLVER AS PART OF A TASK	32
INVOKING EXCEL’S SOLVER AS PART OF A TASK	33
ISOTOPE-RATIO CHART-INSETS	34

SPOT SCAN-GRAPHICS.....	34
MANUAL REJECTION OF SCAN OUTLIERS VIA SPOT SCAN-GRAPHICS	35
APPENDIX I -- THE SQUID-2 TOOLBAR	36
APPENDIX II: ALGORITHMS USED BY SQUID-2.....	37
WEIGHTED AVERAGES, “MEANS”, INTERNAL AND EXTERNAL ERRORS.....	37
<i>The MSWD parameter</i>	37
<i>Probability of fit.....</i>	37
<i>External errors and Weighted Averages assuming the presence of external errors</i>	38
OUTLIER REJECTION FOR THE INTEGRATIONS COMPRISING A SINGLE SCAN OF A PEAK	39
WITHIN-SPOT ISOTOPE RATIOS.....	39
OUTLIERS OF WITHIN-SPOT INTERPOLATED RATIOS	41
THE TUKEY’S BIWEIGHT ALGORITHM.....	42
WITHIN-SPOT, NUMERICALLY-EVALUATED TASK EQUATIONS	43
OVERCOUNTS ON ²⁰⁴ Pb	43
OUTLIER-RESISTANT LINEAR REGRESSION	43
CONCORDIA AGES	43
CONVERTING TERA-WASSERBURG TO CONVENTIONAL-CONCORDIA SLOPE-INTERCEPT	44
²⁰⁷ Pb/ ²⁰⁶ Pb AGE-ERROR INCLUDING ERRORS IN DECAY CONSTANTS:	44
208-CORRECTED AGES AND RATIOS	45
207-CORRECTED AGES AND RATIOS	46
APPENDIX III: THE SQUID-2 PANELS AND CONTROLS	52
PREFERENCES	52
<i>Show splash-screen at startup.....</i>	52
<i>Retain nonessential rows of PD or XML files (file length permitting).....</i>	52
<i>Place data-reduction parameters in a separate worksheet.....</i>	53
<i>Place auto-charts in a separate worksheet</i>	53
<i>Attach Task sheet to reduced-data workbook.....</i>	53
<i>Freeze header rows and spot name column</i>	53
<i>207-corrected 208Pb/232Th ages.....</i>	53
<i>208-corrected Concordia-plot ratios and 207Pb/206Pb ages</i>	53
<i>Create concordia plots for U-Pb Age Standards (U-Pb only)</i>	53
<i>Always construct within-scan isotope ratio sheet</i>	53
<i>Set Y-axis minimum on isotope-ratio plots to zero.....</i>	53
<i>Maximum permitted mass-difference between measured and true nuclide masses.....</i>	53
<i>Default folder for PD and XML files.....</i>	54
<i>Rebuild Task Catalog file.....</i>	54
<i>Isotope Ratios of Common Pb (U/Pb geochron only).....</i>	54
<i>Calculation of within-spot, scan-by-scan means.....</i>	56
<i>U/Pb calibration constant.....</i>	56
<i>CORRECT FOR SECULAR DRIFT.....</i>	56
<i>U-Pb Standards.....</i>	57
<i>Task.....</i>	58
<i>New PD file</i>	58
<i>Number of characters to show in Standard names.....</i>	58
<i>Pb/U or Pb/Th Age.....</i>	59
<i>Index Isotope for common Pb</i>	59
<i>Normalize ion beams to SBM signal</i>	59
<i>Spot subsets for user parameters</i>	59
<i>Time window</i>	59
<i>Trim mass charts.....</i>	60

<i>Task Editor</i>	60
<i>Preferences</i>	60
<i>Change subsets</i>	60
<i>Change Time-window</i>	60
REDUCE GENERAL ISOTOPE DATA	60
<i>New PD</i>	60
<i>Task</i>	61
<i>Normalize to SBM signal</i>	61
<i>Make trim-mass charts</i>	61
<i>Spot subsets for user parameters</i>	61
<i>Process spots within Time Window?</i>	61
<i>Change Subsets</i>	61
<i>Task Editor</i>	61
<i>Preferences</i>	61
<i>Change Time-window</i>	61
SAMPLE GROUPING	62
<i>Spot Names</i>	62
<i># letters to show in Spot-names dropdown</i>	62
<i>Disregard Case, Spaces, Dashes, Slashes...</i>	62
<i>Extract coherent age group</i>	63
<i>Sort samples by age</i>	63
<i>Group all spots together</i>	63
<i>Type of U-Pb Age</i>	63
<i>Grouping Criteria</i>	63
<i>Construct concordia plot for groups if possible</i>	64
<i>Correction for 204 overcounts on Standard</i>	64
<i>Common Pb</i>	64
TIME WINDOW	65
<i>Limits to specify with Standard spots</i>	65
<i>Starting Name/Time</i>	65
<i>Ending Name/Time</i>	65
SUBSETS	66
COPY OR RENAME A TASK.....	67
<i>Task name</i>	67
<i>Task description</i>	67
<i>Primary mineral</i>	67
TASK EDITOR	68
<i>U-Pb Geochronology</i>	68
<i>Isotope Ratio</i>	68
<i>Existing Tasks</i>	69
<i>Edit/view a panel for the selected Task</i>	69
<i>Refresh Task Catalog</i>	69
<i>OK</i>	69
<i>Cancel</i>	69
<i>Constants</i>	69
<i>Save</i>	69
RUN TABLE	70
<i>Ionic species</i>	70
<i>Clear</i>	70
<i>Specify mass station numerically</i>	70
<i>bkrd</i>	70
<i>ref pk</i>	70
<i>cps col</i>	70
<i>Hide?</i>	71
<i>Periodic table of the elements</i>	71
<i>Next</i>	71

<i>Return</i>	71
<i>Cancel</i>	71
<i>Insert</i>	71
<i>Delete</i>	71
<i>Clear one</i>	71
<i>Clear all</i>	71
THE ISOTOPE SELECTION SUB-PANEL.....	72
<i>How many?</i>	72
<i>Charge</i>	72
<i>Enter</i>	72
<i>Complete</i>	72
ISOTOPE RATIOS	73
<i>Ratios (from mass stations at left)</i>	73
<i>Next</i>	73
<i>Back</i>	73
<i>Cancel</i>	73
U-PB SPECIAL EQUATIONS	74
<i>Equation for 206Pb/238Uconst</i>	75
<i>Equation for 208Pb/232Th const.</i>	75
<i>232Th/238U equation</i>	75
<i>Equation for ppm U</i>	75
<i>Next</i>	75
<i>Back</i>	75
<i>Return</i>	75
<i>Cancel</i>	75
<i>Help</i>	75
<i>Append constant</i>	75
<i>Append LITERAL ratios</i>	75
<i>Append INDEX of ratios</i>	75
GENERAL EQUATIONS	76
<i>Next/previous equation (spinner)</i>	76
<i>ST, SA, SC, LA, FO, NU, HI, AR</i>	76
<i>AR #rows #cols</i>	76
<i>Constants</i>	76
<i>Append column-header reference</i>	76
<i>Append LITERAL equations/ratios</i>	76
<i>Next</i>	76
<i>Back</i>	77
<i>Return</i>	77
<i>Cancel</i>	77
<i>Delete (active equation)</i>	77
<i>Insert (active equation)</i>	77
<i>Clear Active</i>	77
<i>Clear All</i>	77
<i>Help</i>	77
<i>Column headers</i>	77
<i>Axis scaling</i>	77
<i>Calculations</i>	78
<i>Back</i>	78
<i>Return</i>	78
<i>Cancel</i>	78
CONSTANTS	79
<i>Edit</i>	79
<i>Add new</i>	79
<i>Delete</i>	79
<i>Exit</i>	79

Append to Eqn 1	79
EXAMPLE TASKS	80
Zircon	80
²³⁰ Th ages, insignificant ²³² Th	82
²³⁰ Th ages, insignificant ²³² Th	83
Simple Monazite	84
APPENDIX IV: WEIGHTED AVERAGE ALGORITHM FOR U/PB AGE STANDARDS.....	85
APPENDIX V: SOME USEFUL EXCEL AND <i>ISOPLOT</i> FUNCTIONS	86
NATIVE EXCEL FUNCTIONS	86
AVERAGE	86
MEDIAN	86
STDEV	86
LINEST	86
INTERCEPT	86
SLOPE	86
MAX, MIN	87
ISOPLOT AND SQUID FUNCTIONS	88
sqBiWeight	88
ChiSquare	88
Concordia	88
ConcordiaTW	88
DisEq68Age	89
DisEq75Age	89
DisEq75Ratio	89
DisEq76Ratio	89
DisEqPbPbAge	89
Drnd	89
ExternalSigma	89
HalfLife	90
InitU234U238	90
La230, La232, La234, La235, La238	90
Lambda	90
Log10	90
Mad	90
MSWD	90
NukeMass	90
Pb76	90
ProbFit	91
RobReg	91
SingleStagePbT	91
SingleStagePbR	91
sqWtdAv	91
StudentsT	91
Th230Age	92
Th230AgeAndInitial	92
Th230U238AR	92
ThUfromPb86	92
U234age	92
U234ageAndErr	92
U234U238AR	92
Output: The present-day ²³⁴ U/ ²³⁸ U activity ratio	92
YorkInter	92
YorkInterErr95	93

<i>YorkMSWD</i>	93
<i>YorkProb</i>	93
<i>YorkSlope</i>	93
<i>YorkSlopeErr95</i>	93
APPENDIX VI: USEFUL <i>ISO</i> PLOT TOOLS FOR <i>SQUID</i>	94
APPENDIX VII: DEFINITIONS OF TERMS	95
<i>204 corrected</i>	95
<i>207 corrected</i>	95
<i>208 corrected</i>	95
<i>Array function or formula</i>	95
<i>Common Pb</i>	95
<i>Concordance</i>	95
<i>Concordia</i>	95
<i>Concordia age</i>	95
<i>Constants</i>	95
<i>Discordance</i>	95
<i>Dropdown</i>	96
<i>Error correlation</i>	96
<i>External error</i>	96
<i>Grouping</i>	96
<i>Internal error</i>	96
<i>Ionic species</i>	96
<i>Isoplot</i>	96
<i>m/e</i>	96
<i>MSWD</i>	96
<i>Normalizing equation</i>	97
<i>Panel</i>	97
<i>PD file</i>	97
<i>Probability of fit</i>	97
<i>Radiogenic Pb</i>	97
<i>Range</i>	97
<i>Resistant (mean, regression...)</i>	97
<i>Robust (mean, regression...)</i>	97
<i>Sample sheet</i>	97
<i>SBM</i>	97
<i>Scan order</i>	97
<i>Single cell formula</i>	97
<i>Spot</i>	98
<i>Stacey-Kramers</i>	98
<i>Standard sheet</i>	98
<i>Trim mass reference</i>	98
REFERENCES	99
INDEX	100
INDEX	100

Introduction

Purpose of SQUID-2

SQUID-2.5 processes PD and XML files produced by a SHRIMP (Sensitive High Resolution Ion Microprobe) into a variety of parameters, including raw isotope ratios, radiogenic Pb isotope ratios, ages for the U-Th-Pb system. The processed data can be applied to a broad range of run tables and applications, ranging from U-Th-Pb geochronology, $^{230}\text{Th}/\text{U}$ geochronology, light stable-isotope geochemistry, and cosmochemistry. Processed data is placed in an *Excel* workbook that, for most output parameters, makes the formulae used in each data processing step explicitly available to the user.

Requirements

Software

SQUID-2 was developed using *Excel 2003* (service pack 3) and *Windows XP* (service pack 2). Earlier versions of *Excel* and other versions of *Windows* are probably compatible, but have not been tested. *SQUID-2* requires the *Isoplot* add-in (ver. 3.7 or later). *Only English versions of Windows and Excel are compatible with SQUID-2 and Isoplot.*

ISOPLOT is BGC's Visual Basic Add-in for Microsoft's Excel® for data analysis and graphical presentation of geochronology, earth science and other radiogenic isotope data only.

BGC's ISOPLOT is not the Isoplot® for analysis of any measuring system as described at (<http://www.shainin.com/>). BGC is not affiliated with Red X Holdings, LC. Isoplot® is a registered trademark of Red X Holdings, and is licensed to Shainin LLC.

At present, *SQUID* is not compatible with *Excel 2007*.

Hardware

WINDOWS

A mid- or high-range CPU (>1 Ghz non-Celeron, dual-core preferred) and 2 Gb or more RAM is recommended. The screen resolution should be at least 1280x800.

EMULATED WINDOWS UNDER MACINTOSH OS-X

At least two of the several programs enabling running of Windows on an Intel Mac, at least two – *Parallels* and *Fusion* – seem to work very well, with virtually no speed penalty compared to Windows. *Crossover* cannot run *Isoplot* or *SQUID*.

Before installing *SQUID* on your Mac, please ensure that your screen resolution is adequate. For wide-screen desktops, set your screen resolution from Windows (right-click anywhere on the desktop, select Properties/Settings) to about 1920x1200. For laptops or smaller desktop screens, use 1440x830 if possible.

Printer Setup

Upon first installation of *Excel*, a printer must be selected (via Page/Setup) before running *SQUID-2* or *Isoplot*, and the page size must be set to Letter or A4.

Input Data Files

SQUID-2.5 only processes PD or XML files for a single ion-counting collector.

Installation

1. If you're running Windows from an Intel Mac under *Fusion* (*Parallels* is incompatible), you must set Windows screen resolution to somewhat higher than the default setting. Right-click anywhere on the Windows desktop (*not* the *Parallels* or *Fusion* desktop) and select Properties, then Settings . If possible, set the resolution to either 1440 x 830 or 1920 x 1200 pixels
2. Remove or rename **all** existing *SQUID* and *Isoplot* add-in file on your computer. *Excel* tends to copy add-in files to its own folders, so you must to a *complete* search of your hard disks (including Hidden and System files) to make sure that *Excel* will not screw everything up by having multiple, incompatible copies of *Isoplot* and *SQUID* loaded at the same time.
3. Start *Excel* 2003. If you find that *Excel* has automatically loaded either *Isoplot* or *SQUID*, close *Excel* and go back to step 2.
4. Select Tools/Add-ins from the *Excel* toolbar at the top of *Excel*'s window (Figs. 1, 2). Check any boxes labeled *Isoplot* or *SQUID* in the Add-ins list. *Excel* should display a message each time saying that the add-in can't be found, and asking if the *Isoplot* or *SQUID* add-in should be deleted from its list. Click Yes. If such a message does *not* appear, close *Excel* and go back to step 2.
5. Transfer the *SQUID-2.5* folder containing the matched pair of *SQUID-2.5* and *Isoplot* 3.7 add-ins the *SquidUser* folder, and the *SQUID-2.5* manual to your computer. Make sure that the folder in which the *SQUID-2.5* folder resides is one for which you have full read/write permission. This folder should probably not be on a network drive.
6. With the Tools/Add-ins box open, click Browse and navigate to the folder containing *SQUID-2.5* and its compatible *Isoplot*. Select the *Isoplot* add-in, click OK, then OK again from the Add-ins box.




Figure 1




Figure 2

7. Repeat Step 6, but this time select the *SQUID-2.5* add-in. *SQUID-2.5*'s splash screen should appear if installation was successful.
8. Before trying out *SQUID-2*, close, then re-open *Excel*. The compatible pair of *SQUID-2.5* and *Isoplot* should now load automatically whenever *Excel* is opened.


Data Reduction for U-Pb Geochronology

Click  in the *SQUID* toolbar (p. 36) or select Process a Pb-U-Th Run from the *SQUID* dropdown list, then navigate to and select the PD or XML file of interest. The data-reduction setup panel for U-Pb geochronology will appear (Fig. 3).

Two Run Tables will appear in the pink box at the left of the panel, the first being the Run Table for the data-reduction *Task* already in memory, the second the Run Table for the actual PD or XML file. If the two run tables don't match, *SQUID* will refuse to process the file.

Note that when you select a PD or XML raw-datafile, *SQUID* will attempt to find a data-processing *Task* whose Run Table matches that of the PD file. This *Task* may or may not be the one desired. If not, select the proper *Task* – say, *Zircon, canonical ANU +UO2* – from the *Task* dropdown.

Run table box



Reduce data for U-Pb Geochronology

PD file (click to change): OG1_070625_GA5005.pd

Task: Zircon, canonical anu +UO2

STANDARDS

☐ 3 Number of characters to show in Standard names

Pb/U or Pb/Th Age

Name: OG1

Pb/U or Pb/Th Age (ma): 3441

²⁰⁷Pb/²⁰⁶Pb age (same as Pb/U if left blank): 3465.4
for 204-overcrt correction only

U or Th concentration

☐ No U or Th standard

Name: SL1

ppm U: 238

☒ U concentration
☐ Th concentration

Index isotope for common-Pb subtraction: 204Pb

☒ Normalize ion beams to SBM signal

Spot subsets for User Parameters: No Subsets

Time Window: No Time Window -- process all spots

☒ Standard spots only ☐ Make trim-mass charts

version 2.31.09.02.10 K.R. Ludwig, BGC, 10 Feb 2009

Figure 3: Setup panel for U-Pb Geochronology data reduction.

U-Pb/Th-Pb Standards

Specifying standards when starting data reduction

SQUID requires a $^{206}\text{Pb}/^{238}\text{U}$ or $^{208}\text{Pb}/^{232}\text{Th}$ standard to be among the analyzed spots in the PD file. A standard for U or Th concentration is optional, and can be either the same or different from the age standard. To specify a standard from the U-Pb data-reduction panel,

1. Select the number of initial characters to characterize the age-standard spots (e.g. 4 for Standard spots with names such as SL13-1.1, SL13-2.1, SL13-3.1...), then select the name of the age standard from the spot-name fragment dropdown.
2. Enter the $^{206}\text{Pb}/^{238}\text{U}$ (or $^{208}\text{Pb}/^{232}\text{Th}$) age of the Age Standard. *Enter the Age Standard's (apparent) $^{207}\text{Pb}/^{206}\text{Pb}$ age only if different from its Pb/U or Th/U age.*
3. Select the U or Th concentration standard (if any) and enter its U or Th concentration.

The Common U-Pb Standards List

You can specify up to 50 common U-Pb age and/or concentration standards from the Preferences panel (p. 57). When the U-Pb data-reduction panel is being filled out, *SQUID* will look for spot names whose first few characters match those of an age standard. If a match is found, that standard's name, Pb/U age (if defined), Pb/Pb age (if defined), and U or Th concentration (if defined) will entered into the appropriate boxes.

Index #	Name (1st chars)	Pb/U Age (Ma)	Pb/Pb Age (Ma)	U conc (ppm)	Th conc (ppm)
1	FC1	1099			
2	TEM	416.8			
3	R33	419			
4	Q	1842	1851.6		
5	BR266	1041			
6	FC1	1099			
7	CZ3	579		555	
8	AS3	1095	1099		
9	AS57	1095	1099		
10	SL13	572		238	
11	BS-1	508		300	
12	MG-1	490		1050	
13	NY/PK	1010		15000	
14	CAP	275.6		100	
15	8153	513.1	479	2250	
16	OG1	3441	3465.4		
17	44069	424		2500	20000
18	PHAL	2060		100	
19	r3	419	419		
20					

Enter only enough initial characters to distinguish the Standard name from other spot names. Non-alphanumeric characters are ignored.

Leave the "Age" field blank if the standard will be used only for uranium or thorium concentration, and vice versa.

Click an Index # to insert a Standard, Shift-Click to delete.

Figure 4: The *U-Pb Standards* page of the *Preferences* panel.

Specifying the common-Pb index isotope (U-Pb geochron panel)

Select ^{204}Pb , ^{207}Pb , or ^{208}Pb from the Index isotope... dropdown (common-Pb isotope ratios are specified in the Preferences panel; p. 52). The use of common Pb index isotopes is discussed in pages 46, 63, 44, and 95)




Figure 5: Selecting the index isotope for correction of common Pb.

SBM Normalization

Indicate whether or not the secondary ion beam (pp. 34 , 59, 61, 96, 97) is to be normalized to the SBM (secondary beam-monitor).

Trim Mass Charts




Figure 6: Trim-mass charts.

If the Make trim-mass graphics box is checked when starting data reduction, *SQUID* will construct scan-by-scan charts of the trim masses for each of the centered mass stations (for Standard spots only).

Click OK.

Data reduction will typically take a few tens of seconds, and yield an *Excel* workbook containing not only the usual StandardData and SampleData sheets, but also the original PD data, condensed and reformatted for ease of examination (Fig. 8, p. 7). An optional Within-Spot Ratios sheet (Fig. 7), containing the calculated scan-by-scan isotope ratios for each spot-burn can also be requested (from Preferences; p. 52).

Reduced-data sheets produced by SQUID

The StandardData Sheet

The Standard sheet contains processed data for only the Age Standard spots. Among the changes from *SQUID-1* are:

- Additional data-columns containing the results for any *Task* equations (e.g. the *Log UO/U* and *Log Pb/U* columns in the example at right; p. 20).
- Up to eight user-specified, *Task*-linked plot insets, (e.g. *Log UO/U* versus *Log Pb/U*) each linked to one or two of the data columns (p. 30).
- A minimum Pb/U external error can be assigned to the Age Standard in the *Preferences* panel (p. 52).¹
- Automatically-generated Concordia plots for the Age-Standard data can be requested (p. 53).
- A correction can be applied for obvious secular drift of the Standard calibration-constant (p.10).

						4-corr 206Pb /238U calibr. const	% err	Age (Ma)	±1σ
238 /196	% err	Log UO/U	% err	Log Pb/U	% err				
0.103	0.55	0.809	0.22	-0.791	0.42	0.00385	0.67	408	2
0.087	0.52	0.817	0.23	-0.772	0.46	0.00389	0.55	413	2
0.170	0.41	0.818	0.17	-0.767	0.33	0.00393	0.39	417	2
0.279	0.35	0.818	0.13	-0.765	0.22	0.00395	0.30	419	1
0.313	0.35	0.818	0.11	-0.766	0.20	0.00394	0.25	418	1.0
0.160	0.38	0.824	0.16	-0.755	0.28	0.00391	0.36	415	1
0.331	0.31	0.820	0.11	-0.762	0.20	0.00396	0.30	420	1
0.151	0.39	0.820	0.16	-0.766	0.29	0.00389	0.41	413	2
0.110	0.43	0.821	0.19	-0.752	0.35	0.00401	0.42	425	2
0.343	0.38	0.820	0.11	-0.759	0.20	0.00398	0.24	422	1.0
0.165	0.38	0.821	0.16	-0.762	0.28	0.00392	0.35	416	1
0.148	0.39	0.814	0.17	-0.760	0.30	0.00406	0.37	430	2
0.287	0.36	0.812	0.12	-0.767	0.21	0.00404	0.26	428	1
0.232	0.43	0.825	0.17	-0.745	0.25	0.00399	0.39	423	2
0.371	0.38	0.825	0.11	-0.747	0.19	0.00399	0.34	423	1
0.259	0.42	0.825	0.13	-0.753	0.23	0.00394	0.29	418	1

Wtd Mean of Std Pb/U calibr. 0.00395
 1σ error of mean 0.21%
 1σ external spot-to-spot error 0.90%
 MSWD 8.85
 Prob. of fit 0.00
 rejected spot(s) 1, 20
 For visualization & preliminary evaluation only



Figure 9: Portion of StandardData worksheet showing the U/Pb calibration columns.

CORRECTION FOR OVERCOUNTS ON ²⁰⁴Pb

SQUID will attempt to quantitatively assess the presence of non-²⁰⁴Pb counts at the ²⁰⁴Pb mass position by 1) assuming precise concordance of the ²⁰⁶Pb/²³⁸U-²⁰⁷Pb/²³⁵U ages, and also 2) assuming precise concordance of the ²⁰⁶Pb/²³⁸U-²⁰⁸Pb/²³²Th ages. The (resistant) mean value of the calculated 204-overcounts on the age-

		7-corr ²	8-corr ²
		206Pb	206Pb
204Pb	204Pb	/238U	/238U
/206Pb	/206Pb	const.	const.
(fr. 207)	(fr. 208)	delta%	delta%
		204	204
		overcts	overcts
		/sec.	/sec.
		(fr. 207)	(fr. 208)

Figure 10: Column headers for 204 overcounts on the U-Pb age standard.

¹ Assigning a reasonable external error, based on an experienced operator's judgment from past data-sets, is important if the age standard has low enough count-rates that the counting statistics errors mask a (smaller) external variance. In such a case, if the sample spots have a significantly higher count rate, so that the external error is comparable to or greater than their counting-statistics errors, the errors assigned by *SQUID* to the sample spots will be underestimated, so that even well-behaved, age-equivalent spots will not yield a statistically coherent Group.

standard spots can then, if desired, be used to correct ^{204}Pb measurements for all of the Sample spots when the samples are Grouped (p.9). The output-columns related to calculation or ^{204}Pb overcounts are shown in Fig. 10.

The SampleData Sheet

As in *SQUID-1*, the Sample sheet contains the basic raw and processed data for all spots but those of the Age Standard. For U-Pb geochronology, *the Sample Sheet is to be used only as an intermediate step in producing the final, Grouped Sample worksheets.*

Sample Grouping and Age Extraction

Data reduction for U-Pb geochronology is not complete until the sample-spot data has been grouped by spot name. Group names are selected from one or more of six drop-down lists containing all Spot names, trimmed to the number of initial characters shown in # of characters spinner-box (red arrow, fig. 11). The dropdown lists will not show spot-name fragments unless at least 2 spots are thus defined.

To Group samples,

Process Samples as a Group

Samples

Spot Names: B1, D1, SM, ZR, [], []

Clear: []

#characters to show in Spot-names drop down: 2

Extract groups based on spot names, and place in separate worksheets, OR

Put all spots together in a single worksheet

Extract coherent age from group(s): [x]

Sort worksheet rows by age: []

Disregard

Case: [x] Aa

Spaces: [x] " "

Dashes: [x] -

Slashes: [x] /

Periods: [x] .

Commas: [x] ,

Colons: [x] :

Semicolons: [x] ;

Age Groups

Type of age

204-corr. 206Pb/238U

207-corr. 206Pb/238U

208-corr. 206Pb/238U

Total 206Pb/238U

204-corr. 207Pb/206Pb

204-corr. 208Pb/232Th

207-corr. 208Pb/232Th

208-corr. 207Pb/206Pb

Grouping criteria

0.05 Min. probability of fit

Min. fract. of spots: 0.60

Correction for apparent 204-overcounts on Standard

No correction -- use measured ^{204}Pb

Force concordance of 207Pb/235U - 206Pb/238U ages

Force concordance of 208Pb/232Th - 206Pb/238U ages

Common Pb

Use Stacey-Kramers ratios for the age type selected at left

OR

Enter specific ratios

206/204

207/206

208/206

Construct concordia plots for groups if possible: [x]

Compare UO/U (or UO2/U or UO2/UO) of Grouped spots with Std: [x]

Figure 11: The Sample Grouping panel.

- Click  in the *SQUID* toolbar, or press the Group Me button at the top of the SampleData sheet.
- Select the number of initial spot-name characters to define a Group (red arrow, below),
- Indicate whether case is significant in the name fragments, and whether or not to ignore spaces, dashes, slashes ... when identifying name-fragment matches,
- Specify from 1 to 6 spot-name groups,
- Indicate if a statistically-coherent age group is to be extracted from each spot-name group.
- If age grouping is requested, then
 - Indicate whether the data-rows in each grouped-sample worksheet should be sorted by spot age,
 - Specify the type of age to use for age grouping, the criteria that define an age group, and whether to construct a concordia-plot inset for the spot group,
 - Specify the method of correction, if any, for ^{204}Pb overcounts (*No correction* is recommended unless you have a good understanding of the concept),
 - Define the common-Pb ratios for the groups.

Correcting U-Pb Geochronology data for Secular Drift of Age-Standard Spots

Normally, the calibration constant (e.g. $\text{Pb}^+/\text{U}^+ \div \text{UO}^+/\text{O}^+$) measured for the different Pb/U or Pb/Th age-standard is expected to remain, within some statistical scatter, constant during the length of the analytical session, so that extracting the true Pb/U or Pb/Th ratios of samples is a simple matter of normalizing the sample calibration constant to the average value of the Standard spots. If, however, the analyst can plainly see a trend in the Standard calibration constants (whether monotonic or cyclic), it is possible that the Sample spots will have been similarly affected. *SQUID-2* provides a very flexible, tuneable, assumption-independent method for dealing with such cases, enabled by checking the secular drift box (right) in the Preferences/Interpolation panel (p. 52).




Figure 12: Setting parameters for correction of secular drift of the U/Pb calibration constant. (*Preferences* panel).

If the user has selected Specified, the secular drift correction algorithm starts by fitting an outlier-resistant smoothed-spline curve (a LOWESS fit; Cleveland, 1979) to the Standard calibration constants, using a moving window for smoothing comprising the specified number of spots. Instead of the usual weighted-average chart, *SQUID-2* will then construct a chart showing the rela-

tive change of the calibration constant versus spot-time, and the corresponding smoothed-spline curve.




Figure 13: Showing the effect of different smoothing-window widths (yellow boxes, in number of spots) on the secular drift curve of the Age Standard calibration constant relative to its median value (red line). Estimated external (non-counting statistics) error is shown in the blue boxes. Small red circles indicate the Age standard calibration constant, black vertical lines their 2σ internal errors.

If, however, the user has selected Automatic for the smoothing window, *SQUID-2* will automatically find the largest smoothing-window that yields an external scatter no greater than the assigned Minimum acceptable external error for the age Standard.(p. 57). The Min. spots parameter (p. 12) sets the lower bound on the automatic smoothing-window.

The secular drift of Age Standard's calibration constant (relative to the median value) is then used to correct the calibration constants of the Standard spots (i.e. to remove the presumed effects of the Standard's drift).

Note that secular-drift correction will only be applied for Tasks with a single parent-daughter system (Pb/U or Pb/Th) specifically calculated (that is, where calculation of the non-Primary parent-daughter ratio is specified as Indirectly in the U-Pb Special Equations panel (p. 27).

Caveats in using secular drift correction

It can be argued that calibration-constant secular drift should not occur in a properly constructed mount analyzed with a properly functioning SHRIMP, and that when such behavior is observed the best remedy is to solve the instrumental problem rather than trying to cleverly manipulate the data. To the extent that this is the case, it is quite possible that attempting to correct for such drift will degrade rather than improve the accuracy and precision of the resulting U/Pb ages. Of particular concern is the possibility that, though the apparent external scatter of the Standard spots may be reduced, drift-corrected Sample spots of the same true age will have increased scatter (thus making it impossible to find a coherent statistical group without excessive rejection), as well as displaying an more-or-less inverse secular drift from the Standard.

Note also that the algorithm used to infer external scatter from the secular-drift curve may yield under- or over-estimated values, as it has not seen extensive testing with real data.

Reducing General Isotope Data

For *SQUID-2*, all analyses that will not be used for U-Pb geochronology are considered *General Isotope* analyses. For *General Isotope* analyses, no separate *StandardData* worksheet is produced (though normalization to a standard is straightforward), and the only data columns that are automatically produced are those for burn time, cps of the Trim-Mass Reference peak, and the isotope ratios specified for the *Task*. Basic processing such as correction for mass fractionation, correction for isobaric interferences, secular changes in ratios, and normalization to a standard are readily handled by *Task* equations.




Figure 14: Setup panel for General Isotope Ratio data reduction.

The Reduce General Isotope data panel for General Isotope data (Fig. 14) is similar to, though much simpler than that of U/Pb geochronology.

Tasks

Introduction

Tasks are the templates used by *SQUID-2* to process the raw PD files. Each *Task* contains the algorithms by which the data is processed into final isotope ratios, ages, element concentrations, and any other parameter of interest to the user.

The gateway into defining, viewing, and editing *Tasks* is the *Task Editor* panel (Fig. 13 – press  on the *SQUID-2* toolbar). There are two classes of *Tasks*: U-Pb Geochronology and Isotope Ratio. The former is specifically constructed for processing analyses on minerals such as zircon into a complete set of relevant total and radiogenic isotope ratios, together with their associated ages. The latter is for any other type of analysis.

Each *SQUID-2 Task* is defined by an *Excel* workbook (*.xls) file contained by the *SquidUser* folder, located in the same folder as the *SQUID-2* add-in file. Names of *Task* files always begin with *SquidTask_*, followed by the name of the *Task* and name or initials of the *Task*'s creator.

Existing *Tasks* defined in the default *SquidTasks* file are shown in the scrollbox at the lower left of the panel. After you create any additional *Tasks* of your own, all of your interactions with the *Task* editor will be to modify existing *Tasks* via one or more of the *Edit/view* buttons.

Task definition is a four or five step process (General Isotope and U-Pb geochron data, respectively), in which you will:

1. Define mass stations for the PD file's Run Table,
2. Specify isotope ratios to be calculated from the Run Table,
3. Define equations for calculation of U/Pb (or Th/Pb) and $^{232}\text{Th}/^{238}\text{U}$ ratios and for U (or Th) concentration (U-Pb geochron only),
4. Define equations for any additional data-processing, and
5. Specify chart insets for x - y plots, means, or secular trends of any of the final data columns.

The best way of to get an idea of what a *Task* consists of is to go through the process of creating one from scratch, using, for example, the canonical (ANU-RSES) run table and equations.

Task Editor [X]

Purpose

☐ U/Pb Geochronology
☒ Isotope Ratios

Edit Constants
 Refresh Task Catalog

OK Save all changes, then exit
 Cancel Discard all changes, then exit
 Save Save changes to selected Task, don't Exit
 Help Help for this panel

Existing Tasks (from Task Catalog)

New Copy Delete

230Th-U ages
 Common Pb, no 204
 Common Pb, with 204
 Common Pb, with 204 & bkrd
 temp Common Pb, no 204
 Xenotime matrix-corr factors, no 207
 Xenotime matrix-corr factors, no 207
 xxxxxxxxxxxxxxxx

Edit/view a Task panel

Name, descr...
 ↓
 Run Table
 ↓
 Isotope Ratios
 ↓
 General Eqns
 ↓
 Auto Charts

Description

Metal and oxides

opal Mineral
 krl Creator


 version 2.30.08.11.28
 K.R. Ludwig, BGC, 28 Nov 2008

Figure 15: First panel of the Task Editor.

The Task Name Panel

Each *Task* is assigned a unique name (Figure 16) – preferably a reasonably pithy one that fits within the *Task* list-box of the *Task* Editor. The *Task* Description should contain additional information about the purpose of the *Task* and type of data required. The *Primary Mineral* and *Creator name of initials* fields are self explanatory.

Edit Task name, description...

Task name
CanonicalZirconPlus270

Task description
UO2 peak included to test for odd spatial effects, as well as alternate U/Pb calculation

Primary mineral
zircon and baddelyite

Creator name or initials
Alistair Schnuckf III

Next Retain panel changes, go to Run Table

Return Retain panel changes, return to Task Editor

Cancel Discard panel changes, return to Task Editor

Figure 16: The Task name, description, mineral, and creator panel.

The Run Table Panel

Defining the Run Table

From the Existing *Tasks* list of the *Task* Editor (Fig. 15), click the New button, thus invoking a blank Run Table panel. For illustration, figure 17 shows a completed Run Table panel (for U-Pb/zircon) , but of course when defining a new *Task*, all of the entries in the right-hand scrollbox will be empty. The steps for definition of the canonical zircon Run Table are:

Run Table

Change the active nuclide

OK
Cancel
Clear

Numeric mass (m/e)

Specify mass station numerically

Instructions

- Click desired element.
- On the ISOTOPES dropdown, specify number of atoms of the element (e.g. "2" for $^{232}\text{Th}_2$), then the isotope. Click "Complete" if this is the last element/isotope to be added to the active nuclide, or "Enter" if you wish to add additional isotopes.
- OR
- Click the "Specify mass station numerically" button and enter the mass as a number, e.g. "204.15".
- If the nuclide/mass station will be used as a trim-mass reference or background position, check the appropriate box in the "Mass station attributes" area.

cps	col	hide?	ref	pk	bkrd	Ionic species (e.g. $^{90}\text{Zr}_2\text{ }^{16}\text{O}$) or numeric mass	True mass (m/e)	Nominal mass (m/e)
1	☒	☐	☒	☐	☐	$^{90}\text{Zr}_2\text{ }^{16}\text{O}+$	195.804	196
2	☒	☐	☐	☐	☐	$^{204}\text{Pb}+$	203.973	204
3	☒	☐	☐	☐	☐	$^{204.1}+$	204.1	204.1
4	☒	☐	☐	☐	☐	$^{206}\text{Pb}+$	205.974	206
5	☐	☐	☐	☐	☐	$^{207}\text{Pb}+$	206.976	207
6	☐	☐	☐	☐	☐	$^{208}\text{Pb}+$	207.977	208
7	☐	☐	☐	☐	☐	$^{238}\text{U}+$	238.051	238
8	☒	☒	☐	☐	☐	$^{232}\text{Th }^{16}\text{O}+$	248.033	248
9	☒	☒	☐	☐	☐	$^{238}\text{U }^{16}\text{O}+$	254.046	254
10	☐	☐	☐	☐	☐			
11	☐	☐	☐	☐	☐			
12	☐	☐	☐	☐	☐			
13	☐	☐	☐	☐	☐			
14	☐	☐	☐	☐	☐			
15	☐	☐	☐	☐	☐			
16	☐	☐	☐	☐	☐			
17	☐	☐	☐	☐	☐			
18	☐	☐	☐	☐	☐			
19	☐	☐	☐	☐	☐			
20	☐	☐	☐	☐	☐			
21	☐	☐	☐	☐	☐			

Next
Back
Return
Cancel
Help

Nuclides
Delete
Insert
Clear one
Clear all

Figure 17: The Run Table panel.

- Set the active Ionic Species or numeric mass box (the one in the right half of the panel) to Scan Order 1, using either the mouse or the  control. The active Ionic Species... box is tan rather than pink, as for the Ionic species box for scan order 10 in the example panel (figure 17). We now wish to enter $^{90}\text{Zr}_2\text{ }^{16}\text{O}^+$ for the first mass position in the run table.
- Click on the zirconium (Zr) box in the periodic table.
The resulting isotope selection dropdown (figure 18) will show all of the stable and long-lived isotopes of the element, with the most abundant selected by default (relative abundances are shown to the left of the isotopes). In this case, the default of ^{90}Zr is fine (to select any other Zr isotope click on one of the black isotope boxes).
- Select the number of atoms of this isotope in the ionic species (using the How Many spin-box), then click OK. $^{90}\text{Zr}_2+$ will appear in the active Ionic Species box in the right half of the panel.
- We know that the final nuclide will have a charge of +1, so we can leave the Charge box as it is.
- Click the *Enter* button, which will dismiss the *Select Isotope* panel, then click on the yellow *O*(xygen) box in the periodic table. Select ^{16}O , quantity 1, charge

ISOTOPES

%abundance	Isotope
51.46	^{90}Zr
11.23	^{91}Zr
17.11	^{92}Zr
17.4	^{94}Zr
2.8	^{96}Zr

Complete
Enter
Cancel

Add, then proceed to next nuclide
Then select another element to add to the active nuclide

How many 90Zr? 2
Charge? +1

Figure 18: The Select isotope box.

1. Click the *Complete* button. The nuclide (ionic species) box should now contain $^{90}\text{Zr}_2^{16}\text{O} + \text{ref}$ in the ionic species box, with 195.804 as its mass to charge ratio, and the active ionic species box should move down one in the list.
6. Check the small *ref pk* box to the left of the ionic species box (figure 20) so that this nuclide will be used as the reference mass for any Trim Mass graphics
7. Enter the singly-charged ^{204}Pb ion as the next (scan order 2) ionic species/mass station by selecting Pb as the element, and ^{204}Pb as the isotope. Click *Complete* again to move to scan order 3.

In PD files that this *Task* will process, scan number 3 is the background mass of, say, 204.09. Because the background “mass” doesn’t correspond to that of any actual nuclide, you must click the Specify mass station numerically button (just above the periodic table; figure 17), type in 204.09, then click OK. Now check the *Bkrd* box to the left of the ionic species box to indicate that this mass station is to be used for background. Move on to the next mass station as before.

Figure 19: The *Specify mass numerically* box.

And so on.

Note that if the ionic species you are defining is, say, doubly or triply charged, you can specify the charge with the spin-button to the right of the box you used to enter the numeric mass. The mass/charge ratio of the ionic species will automatically adjust.

cps col	hide?	ref pk	bkrd	Ionic species (e.g. 90Zr2 16O) or numeric mass	True mass (m/e)	Nominal mass (m/e)
1	☒	☐	☒	90Zr2 16O +	195.804	196

Figure 20: Completed entry for scan order #1.

If you want to have a column containing the (total) counts/second of a mass station placed on the output worksheets, check the corresponding *cps col* box to the left of the *ref pk* box (figure 20). *Cps col* boxes are automatically checked for reference and background mass-stations, and also for the ^{204}Pb and ^{206}Pb mass stations.

When you’re done, look at the Nominal mass column to the right of the True mass column. The Nominal masses are usually (but not always) an integer, but *SQUID-2*, if necessary, will have added just enough additional decimal places to distinguish the masses of each mass station from the other. Thus true mass of 204.09 that you entered for the background mass-position will be shown as 204.1 to distinguish it from the nominal 204 assigned to the $^{204}\text{Pb}^+$ mass station. When referring to the isotope ratios or mass positions of the Run Table in the Equations panels, always use these nominal masses.

Specifying Isotope Ratios

The next step is to specify which isotope ratios should be calculated. The order of the ratios will be preserved in the output worksheets, but otherwise is unimportant. To enter a ratio,

- Click on a blue ratio-box (figure 21) to activate it (its color will change from blue to orange when activated),
- Click on one of the nominal masses in the green boxes to the left (say 204) to use as the numerator of the ratio,
- Click on another of the green boxes to select the denominator isotope (say 206), giving 204/206 in the ratio box, whose color will have changed back to blue.
- To clear the ratio box and start over, click it a third time.
- The ratio boxes don't have to be all filled in sequence – empty boxes are ignored.

Step 2 -- specify isotope ratios

Click on a blue ratio-box to activate or clear it., then click on the green mass-boxes (numerator then denominator) to create the ratio.

Mass Stations		Ratios (from mass stations at left)	
196	204	A 204/206	B 207/206
204.1	206	C 208/206	D 206/238
207	208	E 208/248	F 254/238
238	248	G 248/254	H 238/196
254		I	J
		K	L
		M	N
		O	P
		Q	R
		S	T
		U	V
		W	X
		Y	Z
		AA	AB

Figure 21: The Isotope Ratios panel.




Figure 22: The General Equations Task panel.

The Equation Panels

Basics

Task equations are the backbone of the built-in flexibility of *SQUID-2*. The variables in a *Task Equation* can be any of the worksheet cells in the processed-data sheet(s), including isotope ratios, and the results of other *Task* equations. There are two Equation panels for U-Pb geochronology, and one for General Isotope data. I will discuss the General Equations panel first, as it exists in *Tasks* for both U-Pb/geochronology and General Isotope Ratios.

The simplest type of equation have no “switches” (p. 23) set, and use only the *Tas*’s isotope ratios and numeric constants as arguments for the algebraic functions. Equation results are placed in the SampleData or StandardData output sheets, in data columns to the right of all of the isotope-ratio output columns. If no switches are set, uncertainties for each of the equation results are calculated numerically by *SQUID-2*, with all of the relevant isotope-ratio errors and error correlations correctly propagated.

For example, consider a U-Pb/zircon *Task* that has two User Equations defined (in addition to the four “special” U-Pb equations), as shown above.

	<u>Equation</u>	<u>Equation Name</u>
(1)	$\text{Ln}(['254/238'])$	Log UO/U
(2)	$\text{Ln}(['^{206}\text{Pb}/^{238}\text{U}'])$	Log Pb/U

The first part of the output worksheet will look something like the example in Figure 23, with the results of the two equations (red arrows) placed in columns to the right of the *Task*'s specified isotope ratios.


REFERENCING AN ISOTOPE RATIO IN THE EQUATIONS PANEL:

To add an isotope-ratio reference as an argument in a *Task* equation, you can:


- Type in the literal ratio, enclosed in square brackets plus a quotation mark, e.g. ["207/206"], or
- Type in the ratio's *index letter* (adjacent to the blue isotope-ratio boxes, as in figure 22, enclosed in square brackets without quotation marks, e.g. [b], or
- Just click on the blue isotope-ratio box for the ratio you want.

Spot Name	Date/Time	Hours	Bkrd		Tot 204		Tot 206		204		206		207		208		209		238		248		254		262		Log UO/U		Log Pb/U	
			cts	cts	cts	cts	err	err	err	err	err	err	err	err	err	err	err	err	err	err	err	err	err	err	err	err	err	err	err	err
r33-1.1	2005-05-02, 10:28	0.0	0.05	0.08	395	8.8E-5	69	0.0575	1.6	0.208	2.1	0.162	0.76	0.00868	1.9	6.45	0.41	0.598	0.32	0.103	0.55	0.809	0.22	-0.791	0.42					
r33-2.1	2005-05-02, 10:44	0.3	0.00	0.11	352	3.2E-4	22	0.0567	1.7	0.192	2.3	0.169	0.81	0.00893	2.1	6.55	0.44	0.553	0.35	0.087	0.52	0.817	0.23	-0.772	0.46					
r33-3.1	2005-05-02, 11:48	1.3	0.10	0.08	670	-2.2E-5	93	0.0551	1.2	0.200	1.6	0.171	0.59	0.00890	1.5	6.58	0.32	0.583	0.25	0.170	0.41	0.818	0.17	-0.767	0.33					
r33-4.1	2005-05-02, 12:53	2.4	0.00	0.00	1108	---	100	0.0553	1.0	0.152	1.4	0.172	0.39	0.00883	1.3	6.58	0.26	0.445	0.26	0.279	0.35	0.818	0.13	-0.765	0.22					
r33-5.1	2005-05-02, 13:59	3.5	0.00	0.05	1256	3.7E-5	85	0.0552	1.0	0.188	1.2	0.172	0.35	0.00864	1.1	6.58	0.22	0.561	0.19	0.313	0.35	0.818	0.11	-0.766	0.20					
r33-6.1	2005-05-02, 15:04	4.6	0.00	0.09	561	1.7E-4	22	0.0558	1.4	0.231	1.6	0.175	0.53	0.00908	1.5	6.63	0.32	0.668	0.27	0.138	0.39	0.822	0.17	-0.757	0.30					
r33-7.1	2005-05-02, 16:12	5.7	0.20	0.13	593	-1.2E-4	12	0.0538	1.4	0.239	1.6	0.173	0.51	0.00864	1.5	6.64	0.31	0.716	0.25	0.148	0.38	0.822	0.16	-0.762	0.29					
r33-8.1	2005-05-02, 17:19	6.9	0.06	0.04	614	-3.8E-5	85	0.0564	1.4	0.211	1.6	0.174	0.50	0.00887	1.5	6.58	0.30	0.629	0.26	0.150	0.47	0.818	0.16	-0.759	0.29					


Figure 23 A reduced zircon SampleData worksheet for a Task having two user-defined equations (red arrows).

If you append a ratio by just clicking its ratio box, the ratio will be appended to the end of the existing expression (if any) in the active equation-box. You can choose which style will be used by toggling between Index (appends the letter-index of the ratio, e.g. [h]) and Literal (appends the ratio itself, e.g. ["238/196"]) with the  button.

ADDING A REFERENCE TO A NAMED NUMERIC CONSTANT:

You can define, view, and insert up to 250 constants constants (such as the decay constant of ^{238}U , natural $^{238}\text{U}/^{235}\text{U}$ ratio, age of a secondary standard...) by clicking  at the bottom the Equation panel. Select a constant from the dropdown (figure 24), then click Append to append the constant to the end of any existing text in the active Equation box. The constant will appear either as its name enclosed in angled brackets, e.g. <Lambda238>.

Constants can also be edited from the main *Task* Editor panel, via the  button..



Constants for Task Equations	
1	BSppmTh 8790
2	Lambda230ka .00919525
3	Lambda232 .000049475
4	Lambda232ka .000000049475
5	Lambda234ka .002834
6	Lambda235 .00098485
7	Lambda238 .000155125
8	Lambda238ka .000000155125
9	MG20638 .079
10	MGage 490
11	MGppmTh 1380
12	MGppmU 900
13	MgREEfix 3
14	MultiTask 127.00

Figure 24: The Constants dropdown.

Named numeric constants are stored both in the User Preferences workbook and with the *Task* itself. If a conflict in the values stored for a named *Task* is detected when the *Task* is loaded, you will be prompted to resolve the conflict, and the selected value will be made consistent in both of the stored locations.

ADDING A REFERENCE TO OTHER TASK EQUATIONS:

Any equation can use the results of any other equation as an argument. To do this, you can

1. Type in the name of the equation (see below), enclosed in square brackets plus a quotation mark, e.g. ["Log UO/U"]. Case, spaces, and slashes (/) are ignored by *SQUID*, so you don't have to type them in when defining an equation.
2. Type in the equation's index number (in the orange equation# box to the left of the equation box itself), enclosed in square brackets without quotation marks, e.g. [13], or
3. Just click the equation-number box to the left of the equation box.

Unless you use the FO switch (see below), you should make only *forward* references to other equations – that is, to higher-numbered equations – so that the results of the equation referred to will have been calculated when needed.

EQUATION NAMES:

Each *Task* equation must be assigned a name in the corresponding Equation Name box (figure 20). The equation name will be used as the column header for the equation's output column the SampleData or StandardData worksheet and, if the equation's SC switch is set, as an *Excel* range name². To avoid overly-wide columns for the equation results, you can use the | character to indicate a linefeed. Thus an equation name entered as

Mean rad. | ppm Pb | (uncorr.)

would appear in the equation's column header as:

For Array Equations (see p. 25) whose output occupies more than a single column, you can specify names for each of the output columns using a double vertical-line (||) as the column-delimiter. Thus the output for a name such as

² *SQUID* will ignore non-alphanumeric characters (underscore excepted) in SC equation names, as only alphanumeric characters are permitted in *Excel* range names. Initial number-characters are also forbidden.

Th230 | raw age | (ka) | | err | | Init | 234/238 | | err

assigned to the array equation

Th230AgeAndInitial([3],[±e],[1],[±f], ,true)

with the ARrc switch marked as 14 (output occupies 1 row, 4 columns) would be something like:

Th230 raw age (ka)	err	Init 234/238	err
104.22	1.44	1.1472	0.0001

EQUATION Switches

To provide more control over how and when the *Task* Equations are evaluated, *SQUID-2* provides up to 7 “switches” that can be turned on or off for each equation. The Switches and their abbreviations are (shown in blue font if restricted to U-Pb/geochron):

Switch	Function
Standard	The equation is placed only in the <i>Standard</i> output sheet (default is both).
Sample	The equation is placed only in the <i>Sample</i> output sheet (default is both).
Single Cell	The equation results occupy a Single Cell rather than a column, and operates on whole columns of data rather than the data in just one row.
LAst	The equation is to be calculated last, after all of the other equations without the LA switch, and after all of the columns produced automatically (e.g. radiogenic isotope ratios, U concentration, apparent ages...) by <i>SQUID-2</i> .

The column-order of calculation for U-Pb geochronology is:

- Isotope ratios as specified in the *Task*.
- *Task* equations for which neither the SC nor LA switches are set.
- Automatically calculated columns (radiogenic isotope ratios, U concentration, apparent ages...).
- *Task* equations whose SC, but not LA switch is set.
- *Task* equations with the LA switch set.

The order of calculation within each category follows the equation index-number.

F ormula	The <i>Formula</i> for the <i>Excel</i> equation is to be placed in the output cell(s), rather than its numeric result.
N umeric	The <i>Numeric value</i> of the <i>Excel</i> formula for the equation is to be placed in the output cell(s), rather than its <i>Excel</i> formula.
H idden	The worksheet column containing the Equation results will be hidden.
A rray formula	The equation is an Excel Array Formula whose output can occupy more than one row and/or column.
A Rrc	The number of rows and columns occupied by the output of the array formula (two numbers separated by a space).

The ST and SA switches (U-Pb/geochron): When examining and calculating data that is relevant to, say, just the age-standard spots, the output of user equations that have been created for this purpose (say to examine the Log(Pb/U) vs. Log(UO/U) slope) would only clutter up the Sample output sheet. With the ST switch *on*, the column headers and output for such equations will be calculated and placed only on the Standard sheet. Conversely, the SA switch restricts the equation to the Sample worksheet. If neither of these switches is *on*, the equation is placed on both sheets.

The SC switch (U-Pb/geochron): This switch specifies an equation that operates on whole columns of data, using *Excel* functions such as *Average*, *Median*, *STDEV*..., and *Isoplot* or *SQUID* functions such as *BiWt*, *RobReg*, *AgePb76*, *Th230age*, *InitU234238*..... An SC-switched equation is placed in only one cell (rather than each cell in a data column). Thus if interested in, say, the median of the measured 204/206 ratios, one would create an equation with the formula

Median(["204/206"]

The LA switch (U-Pb/geochron): If a cell reference in a (numerically evaluated) User Equation hasn't been calculated and placed on the output sheet, the output cell of the equation will either incorrect or an *Excel* error. Unless the Equation's SC switch is *on*, data for isotope ratios and lower-numbered equations will always exist. However, the radiogenic-isotope and age-columns created by *SQUID-2* are placed on the output sheet *after* all equations, *except* for those equations whose LA switch is *on*.

The FO and NU switches: *SQUID* puts the results of its calculations in the output-sheets cells as either the numeric value of the calculation (e.g. 1.234) or as an *Excel* formula (e.g. =SC13/\$F13-SQRT(\$AE13-\$AF13)*["ppmu"]). Advantages of a NUmERIC result are:

- Quickest start-to-finish data- processing,
- Least "fragile" output sheets, in that the equation-result cells are invulnerable to changes or errors in the source cells arising from later changes to the worksheet,
- Automatically-calculated errors are (if the equation references only isotope-ratio columns) rigorously accurate, and placed in the column to the right of the equation column.

Advantages of a FOrmulaic result are:

- The equation's output cells are updated automatically when changes are made to cells referenced by the equation,
- Other equations on whose results are referenced by an FO equation need not have been calculated before the FO equation.
- The user can check the contents of the equation cells to make sure that the equation actually implements the user's intent.

The H switch: If set, this switch causes the column containing the Equation's results to be hidden (useful if the output worksheet is too cluttered with non-essential information). To unhide a hidden column, select the two columns flanking the hidden column, then unhide with **Format/Column/Unhide**. To unhide *all* hidden columns in a worksheet, select all cells by clicking on the blank rectangle at the intersection of the Row and Column indices.

The AR switch: This switch, which requires some knowledge of *Excel's* Array functions, adds another level of flexibility in data processing. You will have to use the *Excel* documentation to learn the general ins and outs of Array Functions, but an important feature is that the output of these functions can occupy not just one cell, but any group of contiguous cells with an arbitrary number of rows and columns. For example, the complete output of the *Isoplot* array function,

230ThAgeAndError()

occupies (up to) 5 contiguous cells in the same row, comprising

Age, age $\pm$, initial $^{234}\text{U}/^{238}\text{U}_{\text{AR}}$, initial $\pm$, age-initial error correlation

Other array functions may occupy, say, a 2-column by 5-row range. You must therefore specify the number of rows and columns in the output sheets that *SQUID-2* will need to reserve. To do this, enter first the number of rows (1 to 9), then the number of columns (1 to 9) in one of the Equation panel's **AR #rows #cols** boxes, without spaces or commas – for example, as 24 to indicate a 2 row by 4 column output range.

If the array function's output range occupies both multiple rows and multiple columns, the output cells can be referred to by other equations by their assigned name plus the characters for the #rows and #columns. For example, the second row, first column of a 3-row by 2-column array-function output named "ThAge" would be named "ThAge21". See p. 22 regarding other naming conventions for Array equations.

If, for an Array equation, you specify fewer rows or columns than the maximum possible for the equation (e.g. 2 rows, 1 column instead of the 3 rows by 2 columns that the Array function in the equation's *can* populate), only the that number of rows/columns will be occupied in *SQUID-2's* reduced-data worksheet.

The above description of how to use Array Functions is admittedly inadequate, but should be enough to allow you to experiment until better documentation is written.

Restricting equations to subsets of spots

The operation of a *Task* equation can be restricted to a specific suite of spots that, like the Age Standard spots, can be grouped by the first few characters of their spot names. Subset-restricted equations can be defined either at the time of data reduction (p. 29) or as part of the *Task* equations. When a *Task* equation is restricted to a subset of spots, cells in the equation's output column will be empty for all spots but those in the subset.

Subset-restricted equations are normally defined along with the data-reduction *Task* (p. 26). You can, however, also restrict any equation to a spot subset at the time of data reduction, which is convenient if the names of the subset spots tends to vary from one PD file to another. Instructions are within the panel itself.

Spot subsets are defined for a *Task* equation by appending the appropriate spot-name fragment to the Equation Name, enclosed in curly braces, e.g. Mean $^{207}\text{Pb}/^{206}\text{Pb}$ rad. {qgng}.

Spot subsets are useful when you want to use certain spots as a normalizing standard for some parameter. For example, suppose that, for U-Pb geochronology data, you wish to correct Pb isotope ratios for mass fractionation, and that you have included in the grain mount a secondary standard, say RRhy, with a known $^{207}\text{Pb}/^{206}\text{Pb}$ (the Age Standard $^{207}\text{Pb}/^{206}\text{Pb}$ being insufficiently precise). To do this, one could:

1. Add an equation consisting simply of the SampleData sheet's radiogenic $^{207}\text{Pb}/^{206}\text{Pb}$ column, e.g. ["4corr207*/206*"] with the SA, FO, and LA switches,
2. Name the equation, e.g. RRrhy rad 7/6 {rrhy} ,
3. Add an equation for the resistant mean of the RRhy $^{207}\text{Pb}/^{206}\text{Pb}$, e.g. BiWt(["rrhyrad7/6"] , named, say, Mean RRrhy rad7/6 , with the SA, SC, FO, and LA switches. Note that neither this equation nor the two that follow are marked as Subset equations, since they operate on all of the spot-data within each column.
4. Add an equation for the mass-fractionation factor for Pb, e.g. ["Mean RRrhy rad7/6"]/0.10736 with the SA, FO, and LA switches, where the number in the denominator is the known $^{207}\text{Pb}/^{206}\text{Pb}$ for the RRhy standard. Name the equation, say, MassFract ,
5. Add an equation (approximate, since the 204/206 is to be left uncorrected) for the fractionation-corrected $^{207}\text{Pb}/^{206}\text{Pb}$ of all of the non-RRrhy spots, e.g. ["4corr207*/206*"]*["MassFract"], named, say, Fract. corr. 207*/206* .

You can also restrict *Task* equations to spot subsets via the data-reduction startup panels (p. 26), to avoid defining separate *Tasks* if only the names of, say, a normalizing standard change from one PD file to another.

REFERENCING ERROR COLUMNS

To reference an error column of an isotope ratio or *Task* equation, include **+-** just after the left-hand square bracket, e.g. $[+ - 3]$, $[+ - \text{LogUO/U}]$, $[+ - b]$, or $[+ - 207/206]$. The reference will then refer to the column immediately *after* that of the reference immediately following.

Example (complex):

$$\text{Sum}([1]/([+ - 1]/100*[1])^2)/\text{Sum}(1/([+ - 1]/100*[1])^2)$$

with switches

$$\text{SC} = \text{TRUE}, \text{AR} = \text{TRUE}, \text{ARc} = 1 \ 1$$

will yield the error-weighted mean of the column-values for *Task* equation 1.

“Special U-Pb” Equations

Figure 25: The U-Pb Special Equations panel.

Tasks for U-Pb geochronology have this extra panel, which pertains only to the particular needs of this type of analysis. In this panel, you will specify:

- Which parent-daughter pair is the one of primary interest (^{206}Pb - ^{238}U or ^{208}Pb - ^{232}Th),

- Whether you want the $^{232}\text{Th}/^{238}\text{U}$ calculated directly, by an equation that you specify, or indirectly, from calculated $^{206}\text{Pb}/^{238}\text{U}$ and $^{208}\text{Pb}/^{232}\text{Th}$ ratios,

As well as:

- The normalizing equation to use for Pb/U (or Pb/Th),
- Either the normalizing equation to use for Pb/Th (direct calculation of Pb/Th) or an equation to calculate Th/U from the isotope ratios defined in the isotope ratio panel, and
- The normalizing equation for [U] (or [Th]) from a concentration standard.

By “normalizing equation” I mean an expression whose numeric value remains constant for the parameter of interest. Thus in the present example,

$$\frac{\text{Pb/U}}{(\text{UO/U})^2}$$

is expected to be constant (during the analytical session) for any given value of Pb/U, and

$$\frac{\text{U/Zr}}{(\text{UO/U})^{0.66}}$$

should be constant for any given U/Zr (the proxy for [U]).

Swapping data columns with U-Pb Geochronology *Tasks*

For U-Pb geochronology data, *SQUID* produces several data-columns having to do with calibration constants, radiogenic-isotope ratios and U-Pb ages that are produced in addition to and independent of any *Task* equations. Most of these data columns reference one or more of the isotope-ratio columns, which cannot normally be directly modified by *Task* equations, thus limiting the ability of the user to compensate for factors such as isobaric interferences and mass fractionation.

However, *Task* equations can be modified using a “column swapping” feature, which permits modification of any of the isotope ratio columns³, and thus most of the otherwise-unmodifiable *SQUID* columns. The general idea is to create a *Task* equation that modifies (corrects) an isotope ratio in some way (say by correcting the ratio for an isobaric interference), then to add

³ The four “Special U-Pb” equations are exceptions, as they are calculated directly from the within-scan peak-height measurements.

“<=>” to the equation directing *SQUID* to exchange the modified value with the original result.

For example, suppose there is an interference on ^{204}Pb that is roughly proportional (say about double) to the counts on mass 203. The correction can be made as follows:

1. Add the *Task* equation $["204/206"] - 2*["203/206"] <=> ["204/206"]$,
2. Name the equation.

When done, the 204/206 column of output-data sheet will contain the raw (measured) 204/206 minus double the measured 203/206, and retain the original 204/206 column header, and the *Task* equation's column will contain the original (measured) 204/206, with the header RAW 204/206.

Referencing data in other workbooks or worksheets

For U-Pb *Tasks*, the same column header can exist in both the SampleData and StandardData worksheets. *SQUID* will first look for the referenced column-header in the worksheet doing the referencing, then (if not found) in the other (StandardData or \SampleData) worksheet.

To reference data in either another workbook or to override the above procedure, you must include the workbook and worksheet name using *Excel*'s syntax (required characters in red), e.g.

To refer to cell A1 in another workbook:

'[WorkbookName]SheetName'!RowColAddress

or

'[WorkbookName]SheetName'! RangeName

for example:

'[MonaziteStdData]CorrectionFactors'!B23

to reference cell B23 in the worksheet “CorrectionFactors”, contained in the workbook “MonaziteStdData”.

To refer to cell A1 in another worksheet in the same workbook as the PD file being processed:

'SampleData'!A1

Viewing Task Equations in the reduced-data workbook

Equations associated with the data-reduction Task used for data reduction are always shown in either the *StandardSpots* or *SampleSpots* (U-Pb geochronology or General Isotope Ratio Tasks, respectively) worksheets via the *Equations Box* (see example at right). If the Attach Task sheet to reduced-data workbook box in *Preferences* is checked, the equations will also be shown in the attached Task sheet.

Task Equations for T2 = 9pk Curtin	
$(238/196)/(254/238)^{0.66}$	U concentration
$(0.03446*(254/238)+0.868)*(248/254)$	232Th/238U
const.	
$(206/238)/(254/238)^2$	206Pb/238U
const.	

Equation Tips

- Use the sqBiWeight (p. 88) or Median functions instead of Excel's Average to avoid worrying about outliers.

If you need to use an error-weighted mean, use the sqWtdAv function (p. 91).

- Use the RobReg function (p. 91) instead of Slope, Linest, Yorkfit, or Yorkslope functions, again to eliminate sensitivity to outliers.
- Use the Min and Max functions to restrict the range of an output, e.g. $\text{Min}(3, \text{Max}(1, \text{Median}([\text{"CalibrSlope"}])))$ to force CalibrSlope within the range of 1 and 3.
- To deal with initial disequilibrium in the ^{238}U or ^{235}U decay chains, use one of the relevant *Isoplot* functions (p. 89),

DisEq68Age	DisEq68Ratio
DisEq75Age	DisEq75Ratio
DisEqPbPbAge	DisEq76Ratio

The Autocharts Panel

The final panel in the definition of a Task is for Auto Charts (figure 26). Auto Charts are simple x - y or weighted average plots of any of *SQUID-2*'s output sheets that are placed as insets in these sheets. Up to 8 *Autocharts* (p. 30) can be requested. You can choose the x and (if not a weighted average plot) y axes from any of the output sheet's column-headers using the x - y dropdown. You may also request that an (outlier resistant) linear regression be calculated and plotted for x - y data, or an (outlier resistant) mean for just the y data. If Y mean is selected, the x axis is always Hours.

For U-Pb geochronology Tasks, Autocharts will be created only for the age-standard spots.

Chart #	Column headers	Axis scaling	Calculations
1	X: Ln[UO/U] Y: Ln[Pb/U]	Autoscale, Starts at 0, Log	None, X-Y regression, Y mean
2	X: Y:	Autoscale, Starts at 0, Log	None, X-Y regression, Y mean
3	X: Y:	Autoscale, Starts at 0, Log	None, X-Y regression, Y mean
4	X: Y:	Autoscale, Starts at 0, Log	None, X-Y regression, Y mean

Buttons on the right: Return (Retain panel changes, return to Task Editor), Back (Discard panel changes, go back to General eqns), Cancel (Discard panel changes, return to Task Editor), Help.

Figure 26: The Auto Charts panel.

In figure 26 (above), for example, the results of the two (user-defined) equations corresponding to $\text{Ln}({}^{206}\text{Pb}/{}^{238}\text{U})$ versus $\text{Ln}(254/238)$ are plotted, as well as the line and numeric results of the x - y regression.

Normally, *Autocharts* will be placed on the reduced-data worksheet of the spots themselves. will be placed on the *Standards* reduced-data sheet

Setting a "Time Window" for spots

Data reduction can be restricted to any sequentially-measured series of spots. For U-Pb geochronology, the series of spots is defined by the bracketing Age-Standards. For General Isotope data, the time window is defined by any pair of bracketing spots.

Standard spots that bracket the window

Limits to specify with Standard spots:

- ☐ Lower (first) only
- ☐ Upper (last) only
- ☒ Both lower and upper
- ☐ No window (reduce all spots)

Starting: Name: r33-1.1, Time: 2005-05-02, 10:28

Ending: Name: r33-24.1, Time: 2005-05-03, 09:04

Figure 27: The Time Window panel.

Specifying Subset-restricted equations at the time of Data Reduction

A “subset-restricted equation” is a user-defined *Task* equation that will operate only on a subset of the spots in the PD file, where the subset is defined by the initial characters of the spot names. This can be very useful, for example if one wishes to know the mean U concentration of just those spots whose names start with “MG”.

Subsets can be assigned either by the data-reduction *Task* (p. 37) or at the time of data reduction by clicking the Subsets button on the data-reduction panel.

Scope of user parameters (spots for which equations will be applied)

To restrict the application of any equation to a subgroup of spots that share the first few letters of their names:

- 1) Specify (using the "# of initial..." spinner) the number of initial spot-name characters that are shared by the subgroup (for example, 2 for all spots whose names start with, say, SW).
- 2) Select (from the "Trimmed spot-names") the first N letters that define the subgroup (say, SW), where N is the number specified in the spinner box above.
- 3) Click on the yellow box corresponding to the name of the equation whose application is to be restricted only those spots whose names start with, say, SW. Click the box again to cancel the action.

SUBSET EQUATIONS DO NOT APPLY TO AGE-STANDARD SPOTS

Name of equation to apply to spot-name subset	Initial spot-name characters
Log UO/U	1
Log Pb/U {q5}	q5
Expo	

Number of initial letters to show in spot-names drop-down: 2

Trimmed spot-names: q5

Figure 28: The Subset panel for specifying subset equations at the time of data reduction.

Invoking Excel's *Solver* as part of a Task

Excel includes a powerful *Solver* add-in that numerically finds the solution to a wide variety of nonlinear equations having more than one unknown. In particular, *Solver* will calculate and place in specified worksheet cells the values for one or more parameters that minimize the value in another cell. Thus one could determine the values of, say, a , b , c , and d that minimize the sums of squared weighted residuals (or, equivalently, the sigma-mean of a column of data), thus providing a least-squares solution to an column of scattered data that you wish to fit to, say, $y=ax^3+bx^2+cx+d$.

To have a *Task* invoke *Solver*, add a *Task* equation into the list of *Task* Equations and enter first the location of the (single) cell whose contents are to be minimized, followed by semicolon, then the location of the parameters whose values you wish to optimize. Assign a name (which will be ignored) to the equation, followed by <<solve>>. Enter two worksheet range-names(or cell addresses) separated by a semicolon. The first of these names specifies the single cell whose contents are to be minimized, the second the range containing the parameters to adjust, for example

SumWtdResids;CubicEquationCoefficients

where SumWtdResids is the range-name of the parameter to be minimized, and CubicEquationCoefficients is the name of the range containing, the trial values for a , b , c , and d . If the parameters to adjust are in separate cells, separate the range names or cell addressed by a comma. For example

SumWtdResids; Coef1, Coef2, Coef3, Coef3, Coef4

If the cell whose value is to be minimized is part of an Array equation, you can use the row-column values to specify which parts of the Array-function output to use. Thus if a *Task* equation named biwthtpbage is defined as sqBiWeight(WeightedResids) with the .FO, .AR switches set and with ARrc = 3 2 i.e. (3 rows, 2 columns), the output cells created by *SQUID-2* might look like:

xeno1	
mean	
8*/2 age	
994.25	Biweight Mean
21.77	Biweight Sigma
14.19	95%-conf. error

If the cell *BiWeight Sigma* contains the parameter to be minimized, the *Task* equation would look like

biwthtpbage(2,1); Coef1, Coef2, Coef3, Coef3, Coef4

indicating the value in the second row, first column of the Array equation's output.

Solver will be invoked after all of the data columns in the output worksheet have been calculated and filled with values.

Isotope-Ratio Chart-Insets

After PD-file data reduction is complete, and if the Construct within-scan... box in the *Preferences* panel (p. 52) is checked, SQUID will construct a *Within-scan isotope ratios/eqns* sheet (figures 7, 29) as part of the data reduction process, so that chart insets of any of the isotope ratios or Equation results can be constructed on demand.

To invoke one or more within-scan isotope-ratio (or within-scan equation) chart, select the cell(s) whose row (s) and column(s) correspond to the spot(s) and isotope ratio/equation(s) of interest, then click the  button in the *SQUID* toolbar.




Figure 29: Some isotope ratio and Equation result charts


The Y-axis of the chart-insets will be auto-scaled to the data plotted, unless the “Set Y-axis minimum on isotope-ratio charts to zero” box in *Preferences* is checked.

Spot Scan-Graphics

You can request spot-scan graphics – by which I mean charts of the total counts at each mass station for each peak jump – for any spot, by selecting one or more cells whose rows correspond to the desired spots (in either the StandardData or SampleData worksheets), then clicking the  button in the *SQUID* toolbar. The resulting charts (figure 30a) – one for each mass station plus one for the SBM – will all be placed on a single new worksheet having the same name as the spot.

Clicking the  button produces spot-scan graphics for *all* of the spots in the raw-data file.

This can take several minutes. For extremely large files or run-tables, it is possible for *Excel* to run out of memory, resulting in some obscure and misleading error message such as “too many fonts”.



Manual Rejection of Scan Outliers Via Spot Scan-Graphics

SQUID's algorithm for within-scan isotope-ratio and equation means is extremely resistant to outliers, so that manual rejection by the user of individual scans (peak-top measurements) should not be needed. If, however, you do wish to reject scans for whatever reason, you can do so from any spot-scan graphics sheet in the following way:

1. Request construction of spot-scan graphics for the desired spot.
2. Click on the chart for the mass station of interest. A bar containing numbered light-blue boxes will appear at the top of the chart (figure 30b).
3. Click on the box corresponding to the scan number to be rejected (figure 28c). The box will turn light red, the rejected data-point will disappear from the chart, the curved line through the data points will be redrawn, and the values in the condensed PD file worksheet will change as shown in figure (to un-reject a scan, click again on the numbered, light-red box).

Now, when you re-process the PD file, *SQUID* will ignore the rejected scan (figure 31).

30	tot cts	mass	30	tot cts	mass
Secs	206.946	±1sig	Secs	206.946	±1sig
67	716.599	33.1	67	716.599	33.1
182	830.449	28.1	182	830.449	28.1
297	444.143	19.7	297	REJ	19.7
411	794.447	35.7	411	794.447	35.7
			67	716.599	33.1
			182	830.449	28.1
			411	794.447	35.7

Source data for plot after rejection.




Figure 30b: Manual rejection of scan 3 of mass 207.

Figure 31: Change in the condensed PD file for the chart in 28b after scan rejection.

Remember that deleted scans will only affect subsequent processing of the PD file, and will have no effect on existing SampleData or StandardData worksheets.

Appendix I -- The *SQUID-2* Toolbar

(to make the toolbar visible, select View/Toolbars/Squid2 from the *Excel* dropdown menus)



Load and process data from a U-Pb geochronology PD file.



Load and process data from a General Isotope PD file.



Group samples from a processed U-Pb geochronology workbook.



Condense and reformat one or more PD files.



Task Editor



Preferences



Construct one or more chart insets of the within-spot isotope ratios for the spot and ratio defined by the selected row and column of the SQUID StandardData or SampleData worksheet. To obtain more than one plot at a time, select more cells, each indicating an isotope ratio.



Construct charts of cps versus time for one or more spots defined by the selected row(s) of the SQUID StandardData or SampleData worksheet, for all of the mass stations.



Construct charts of cps versus time for *all* spots in both the StandardData and SampleData worksheets, for all of the mass stations.



“Freeze” the active worksheet by replacing all cell-formulae with their numeric results.



“Freeze” all of the worksheets in the workbook.



Make a publication-type table from the processed U-Pb geochronology data.



Widen all columns of the worksheet by 1 point.



Narrow all columns of the worksheet by 1 point.



Update all SQUID and Isoplot references to the currently-installed versions.



Quit SQUID 2 (uninstall the SQUID add-in). You can re-install SQUID 2 at any time



Unhide/Rehide toggle for hidden columns in SQUID 2's SampleData and StandardData sheets.



Show this description of the toolbar buttons.

Appendix II: Algorithms Used by *SQUID-2*

Weighted averages, “means”, Internal and External errors

The term “weighted average” is generally used to imply a mean weighted by the inverse variance (variance = σ^2) of the assigned errors to each data point, so for N data points,

$$\bar{x} = \frac{\sum_{i=1}^N \frac{x_i}{\sigma_{x_i}^2}}{\sum_{i=1}^N \frac{1}{\sigma_{x_i}^2}}$$

with the “internal” error on $\bar{x}$ (that is, the error propagated from just the N σ_{x_i}) calculated as

$$\sigma_{\bar{x}}^2 = \frac{1}{\sum_{i=1}^N \frac{1}{\sigma_{x_i}^2}}$$

The *MSWD* parameter

Two extremely useful parameters associated with the error-weighted average are the *MSWD* and *probability of fit*. *MSWD* (mean square of weighted deviates)⁵ is defined as the mean of the N inverse-variance weighted residuals, (residual = measured x minus weighted-mean x), squared. So

$$\frac{\sum_{i=1}^N \left(\frac{x_i - \bar{x}}{\sigma_{x_i}} \right)^2}{N - 1}$$

(note that degrees of freedom is used as the divisor rather than N).

Intuitively, one can see that, if each σ_i is truly the cause of scatter from $\bar{x}$, on average the residuals for the x_i should be σ_i , so that the expected mean of the both the weighted residuals, and their square should be about 1. An *MSWD* much greater than one, therefore, implies that the assigned data-point errors are unlikely to be responsible for the observed residuals (or amount of scatter), and an *MSWD* much less than one implies that the assigned data-point errors are probably too high.

Probability of fit

It turns out that $MSWD \times v$, where v ($=N-1$ for a simple weighted average) follows a chi-squared distribution about v , so that the probability that one will obtain a given *MSWD* or greater – in other words, the *probability of fit*, can be calculated from the *Excel* function as

$$= \text{CHIDIST}(MSWD * \nu, \nu).$$

External errors and Weighted Averages assuming the presence of external errors

If the probability of fit is low – say less than 0.05 – it is likely that some other source of error in addition to the assigned data-point errors is present. If one assumes that this *external error* is the same from data point to data point, so that for each data point the true variance is

$$total \sigma_{x_i}^2 = internal \sigma_{x_i}^2 + external \sigma^2$$

where $\sigma_{internal}$ is the data point's assigned error and that this external error has a Gaussian distribution (both of which seem to be the case for SHRIMP analyses), its (likely) magnitude can be calculated. For U-Pb geochronology SHRIMP analyses, if the age standard analyses indicate that an external error for the U/Pb calibration constant (likely) exists, as indicated by a low probability of fit, the additional variance should be added to the internal variance of the sample (non-standard) spots before taking a weighted average.

To solve for the external error, we first set up the Maximum Likelihood Error (MLE) equation,

$$L(\bar{x}, \sigma_E^2) = \text{constant} - \frac{1}{2} \sum \log(\sigma_{x_i}^2 + \sigma_E^2) - \frac{1}{2} \sum \frac{(x_i - \bar{x})^2}{\sigma_{x_i}^2 + \sigma_E^2}$$

where $\bar{x}$ is the MLE estimate of the true (mean) x and σ_E^2 is the constant, external variance. We then set up the first-derivative constraints which, when solved (numerically) yield the desired MLE $\bar{x}$ and σ_E^2 .

$$\frac{\partial L}{\partial(\sigma_E^2)} = -\frac{1}{2} \sum \frac{1}{\sigma_{x_i}^2 + \sigma_E^2} + \frac{1}{2} \sum \frac{(x_i - \bar{x})^2}{(\sigma_{x_i}^2 + \sigma_E^2)^2} = 0$$

$$\frac{\partial L}{\partial \bar{x}} = + \sum \frac{x_i - \bar{x}}{\sigma_{x_i}^2 + \sigma_E^2} = 0$$

So

$$\bar{x} = \frac{\sum \frac{x_i}{\sigma_{x_i}^2 + \sigma_E^2}}{\sum \frac{1}{\sigma_{x_i}^2 + \sigma_E^2}} \quad \text{and} \quad \sum \frac{(x_i - \bar{x})^2}{(\sigma_{x_i}^2 + \sigma_E^2)^2} = \sum \frac{1}{\sigma_{x_i}^2 + \sigma_E^2}$$

Uncertainties of the two parameters are given by inversion of the negative values of the second derivatives of L , which are

$$\frac{\partial^2 L}{\partial(\sigma_E^2)^2} = +\frac{1}{2} \sum W_i^2 - \sum W_i^3 r_i^2 \quad \frac{\partial^2 L}{\partial \bar{x}^2} = - \sum W_i \quad \frac{\partial^2 L}{\partial \bar{x} \partial(\sigma_E^2)} = - \sum W_i^3 r_i$$

where $W_i = \frac{1}{\sigma_{x_i}^2 + \sigma_E^2}$ and $r_i = x_i - \bar{x}$, and then substituting zero for the r_i (their expected value) to obtain the expected value of the second derivatives,

$$-E \frac{\partial^2 L}{\partial (\sigma_E^2)^2} = -\frac{1}{2} \sum W_i^2 \quad -E \frac{\partial^2 L}{\partial \bar{x} \partial (\sigma_E^2)} = 0 \quad -E \frac{\partial^2 L}{\partial \bar{x}^2} = + \sum W_i$$

The variance-covariance matrix of $\bar{x}$ and σ_E^2 is then obtained by inverting the Fisher Information Matrix,

$$\mathbf{F} = \begin{pmatrix} -\frac{1}{2} \sum W_i^2 & 0 \\ 0 & + \sum W_i \end{pmatrix} = \begin{pmatrix} \sigma_{\bar{x}}^2 & cov(\bar{x}, \sigma_E^2) \\ cov(\bar{x}, \sigma_E^2) & \sigma_{\sigma_E^2}^2 \end{pmatrix}$$

Outlier rejection for the integrations comprising a single scan of a peak

SHRIMP's raw-data PD file reports the secondary ion-beam counts as ten integrations for each mass station, regardless of the time of the integration. If the median of these ten integrations is <100, *SQUID-2* identifies 95%-confidence outliers as being any integration whose value is outside the Poisson distribution for the 10-integration median, and uses the mean of remaining integrations. If the median is ≥ 100 counts, *SQUID-2* uses the Tukey's Biweight algorithm to obtain an outlier-resistant mean.

Within-spot isotope ratios

If either the numerator or denominator peak of a ratio comprises fewer than 32 counts during a burn (that is, the sum of the peak-top integrations for the n scans), isotope ratios are calculated as the ratio of the sum of the peak-top counts for the two isotopes rather than by interpolation.

Otherwise, *SQUID-2* calculates ratios and NU-switched *Task* equations either regressing the double-interpolated ratios (, 1978) versus time and takes the value of the regression at mid-burn time to cancel out any secular trends during the burn. The errors of individual interpolations are calculated by

$$T = \frac{1+f}{2} \quad f_1 = \frac{1-f}{2} \quad f_2 = \frac{1+f}{2}$$

$$A_i = f_1 a_i + f_2 a_{i+1} \quad B_i = f_2 b_i + f_1 b_{i+1}$$

$$\sigma_{r_i}^2 = \frac{1}{B_i^2} \left[f_1^2 a_i^2 \frac{\sigma_{a_i}^2}{a_i^2} + f_2^2 a_{i+1}^2 \frac{\sigma_{a_{i+1}}^2}{a_{i+1}^2} \right] + \frac{A_i^2}{B_i^4} \left[f_2^2 b_i^2 \frac{\sigma_{b_i}^2}{b_i^2} + f_1^2 b_{i+1}^2 \frac{\sigma_{b_{i+1}}^2}{b_{i+1}^2} \right]$$

where

a and b are the two peaks comprising the ratio $r = a/b$,

T is the time of interpolation,

f is the fraction of the time interval between successive integrations of peaks a or b ,

i is the scan-number index,

σ_a and σ_b are the counting-statistic errors of peaks a and b , and

A and B are the interpolated values of a and b .

Because each pair of time-adjacent interpolated ratios shares two of the four peaks used to calculate each ratio, the errors of such ratio pairs are necessarily correlated, which must be taken into account when calculating the error and $MSWD$ of the mean ratios for a spot, using

$$\rho_{r_i, r_{i+1}} = \frac{1}{4} \frac{1-f^2}{B_i B_{i+1} S_{r_i} S_{r_{i+1}}} \left(\frac{\sigma_{a_{i+1}}^2}{r_i r_{i+1}} + \sigma_{b_{i+1}}^2 \right)$$

Good approximations, which assume the absence of major changes in the size of the peaks during the spot burn, are

$$\frac{S_{r_i}^2}{r_i^2} \approx \frac{1+f^2}{2} (S_a^2 + S_b^2)$$

and

$$\rho_{r_i, r_{i+1}} \approx \frac{1-f^2}{2(1+f^2)}$$

For non-negligible error correlations, the general formula for the weighted mean and error must be used,

and

$$\sigma_r^2 = \frac{1}{\sum_{i,j} \Omega^{ij}}$$

where Ω^{-1} is the variance-covariance matrix for the interpolated ratios r_i , and Ω^{ij} indicates the element in the i th row, j th column of Ω . If the r_i are uncorrelated, Ω^{-1} is simply the diagonal matrix for the inverse variances of the r_i . For doubly-interpolated ratios, however, the near-diagonal terms are nonzero, so the variance-covariance matrix is:

$$\Omega^{-1} = \begin{pmatrix} \sigma_{r_1}^2 & \frac{\sum_{i,j} r_i \Omega^{ij}}{\sum_{i,j} \Omega^{ij}} & \dots & 0 \\ \text{cov}(r_1, r_2) & \sigma_{r_2}^2 & \text{cov}(r_3, r_2) & \dots \\ 0 & \text{cov}(r_2, r_3) & \dots & \dots \\ \dots & \dots & \dots & \text{cov}(r_{N-1}, r_N) \\ 0 & \dots & \text{cov}(r_N, r_{N-1}) & \sigma_{r_N}^2 \end{pmatrix}$$

If user Preferences specify that isotope ratios be calculated as the estimated value at mid-burn time via a linear regression versus time, the error-weighted linear regression must again take into account the error correlations between time-adjacent ratios, so that, for

$$r_i = mt_i + b$$

where r_i are the doubly-interpolated isotope ratios and t_i are the corresponding times of interpolation, and m, b is the time-slope and intercept of the linear regression. The best-fit slope and intercept are given by

$$m = \frac{2WP_r - P_t P_r}{4WM_t - P_t^2} \quad b = \frac{P_r - mP_t}{2W}$$

where

$$\begin{aligned} P_r &= \sum \sum (t_i r_j + t_j r_i) \Omega^{ij} & M_t &= \sum \sum t_i t_j \Omega^{ij} \\ P_t &= \sum \sum (t_i + t_j) \Omega^{ij} & P_r &= \sum \sum (r_i + r_j) \Omega^{ij} \\ W &= \sum \sum \Omega^{ij} \end{aligned}$$

and Ω is defined in p. 40.

Error for the slope and intercept are obtained from inversion of

$$\mathbf{F} = \begin{pmatrix} M_t & P_t \\ P_t & W \end{pmatrix}$$

which gives

$$\mathbf{F}^{-1} = \begin{pmatrix} \sigma_m^2 & cov(m, b) \\ cov(m, b) & \sigma_b^2 \end{pmatrix}$$

The error in r at any time t is then obtained from

$$\sigma_r^2 = t^2 \sigma_m^2 + \sigma_b^2 + 2 cov(m, b)$$

Outliers of within-spot interpolated ratios

For calculation of isotope-ratio means by either a simple weighted average or error-weighted linear regression at mid-burn time, the following algorithm is used:

1. the effect on the *MSWD* (of the error-weighted average) of removing each of the $n-1$ ratios is calculated ($n = \#$ scans), and the ratio whose rejection results in the greatest reduction in *MSWD* is determined.
2. If this lowest *MSWD* is lower than the no-rejection *MSWD* by a selected factor⁶, the point is rejected. A second rejection pass is permitted if there are more than 8 or 9 (weighted average and linear regression, respectively) ratios.
3. If the probability-of-fit of the weighted average or linear regression is ≤ 0.05 , the internal sigma-mean of the weighted average is accepted as the 1-sigma error. If not, a $\sqrt{\text{MSWD}}$ multiplier is applied.

The regression is error-weighted (by counting statistics), and removes outliers by sequentially deleting each data point (the isotope ratio or equation result). If any single deletion drops the *MSWD* of the regression by a factor of 3 or more, the point is rejected and the process repeated. The maximum permitted number of rejections is 1 for fewer than 9 ratios or equations, 2 if more. Or,

- Calculates the error-weighted mean of the interpolated ratios. If the probability of fit (calculated from the *MSWD* of the weighted mean) is greater than 0.05, this mean and standard error (propagated from the counting statistics errors of the individual scans) is accepted.
- Otherwise, the mean and standard error of the outlier-resistant Tukey's Biweight algorithm (Hoaglin *et al*, 1983) is used.

The Tukey's Biweight algorithm

The initial estimate of mean ($\bar{x}$) is the median; and of σ is MAD (the median absolute deviation from the median). Schematically, the algorithm is:

$T = \text{tuning const.} \times \sigma$. The tuning constant is either 6 or 9; *SQUID-2* uses 6.

$\Delta = x_i - \bar{x}_j$ where x_i is the i th point of the data-set, and $\bar{x}_j$ is the j th iteration for $\bar{x}$.

Do:

For each data-point ($i=1$ to n),

If $|\Delta| < T$ then:

⁶ 0.15, 0.2, 0.2, 0.25 for 3, 4, 5, 6, and 7 d.f. (=degrees of freedom; $n-1$ for weighted average, $n-2$ for linear regression), respectively, 0.3 for 8 or more d.f. Outlier rejection is forbidden for <3 d.f.

$$\Delta_i = x_i - \bar{x}_j$$

$$u_i = \Delta_i / T$$

End If

$$SA = \sum_i \left[\Delta_i (1 - u_i^2) \right]^2$$

$$SB = \sum_i (1 - u_i^2) (1 - 5u_i^2)$$

$$SC = \sum_i u_i (1 - u_i^2)^2$$

$$\sigma = \sqrt{n \times SA / |SB|}$$

$$\bar{x}_{j+1} = \bar{x}_j + T \times SC / SB$$

Loop until $\sigma_j = \sigma_{j-1}$ and $\bar{x}_j = \bar{x}_{j-1}$

Within-spot, numerically-evaluated Task equations

Same as for isotope ratios, but instead the equation itself is calculated N times ($N = \text{\#scans} - 1$) at the mean time of all of the mass-stations involved in the equation. Errors of each scan's equation are calculated numerically by perturbing all of the peaks involved in the equation.

Overcounts on ^{204}Pb

Apparent overcounts on ^{204}Pb are calculated in two ways: 1) assuming that the true $^{206}\text{Pb}^*/^{238}\text{U}$ and $^{207}\text{Pb}^*/^{235}\text{U}$ ages are concordant, and 2) assuming that the true $^{206}\text{Pb}^*/^{238}\text{U}$ and $^{208}\text{Pb}^*/^{232}\text{Th}$ ages are concordant. The $^{204}\text{Pb}/^{206}\text{Pb}$ ratio required to achieve concordance is calculated, and thus the number of non- ^{204}Pb counts that must be present to explain the apparent discordance.

Outlier-Resistant Linear Regression

This outlier-resistant linear regression algorithm is used by *SQUID-2* for *Task Autocharts*, and when invoked by the user in a *Task* equation. The regression-line slope/intercept is the median of all pairwise slopes/intercepts, with errors based on code in Rock & Duffy (1986), derived from Vugrinovich (1981).

Concordia Ages

See Ludwig (1998) for a detailed explanation of the definition and algorithm for the Concordia Age.

Converting Tera-Wasserburg to Conventional-Concordia Slope-Intercept

With $x = {}^{238}\text{U}/{}^{206}\text{Pb}_{\text{rad.}}$ and $y = ({}^{207}\text{Pb}/{}^{206}\text{Pb})_{\text{rad.}}$, $m = u/b'$ and $b = -m'/b'$, where m, b, u are as before and m', b' are the Tera-Wasserburg slope and intercept, the transformed errors are given by

$$S_m^2 = S_{b'}^2$$

$$S_b^2 = S_{m'}^2 + S_{b'}^2 - 2S_{m'}S_{b'}\rho_{m'b'}$$

$$\rho_{m,b} = \frac{S_{m'}\rho_{m',b'} - S_{b'}}{S_b}$$

where

$$S_m = \frac{\sigma_m}{m}$$

$$\rho_{m,b} = \frac{\text{cov}(m,b)}{\sigma_m\sigma_b}$$

${}^{207}\text{Pb}/{}^{206}\text{Pb}$ age-error including errors in decay constants:

The simple (first-derivative expansion) expression for ${}^{207}\text{Pb}/{}^{206}\text{Pb}$ age-errors is:

$$\sigma_t^2 = \frac{(e^{\lambda_{238}t} - 1)^2 \sigma_r^2 + (ute^{\lambda_{235}t})^2 \sigma_{\lambda_{235}}^2 + (rte^{\lambda_{238}t})^2 \sigma_{\lambda_{238}}^2}{(u\lambda_{235}e^{\lambda_{235}t} - r\lambda_{238}e^{\lambda_{238}t})^2}$$

where $r = {}^{207}\text{Pb}/{}^{206}\text{Pb}$ and $u = {}^{235}\text{U}/{}^{238}\text{U}$

208-Corrected Ages and Ratios

Constrain the $^{206}\text{Pb}^*/^{238}\text{U}$ age to equal the $^{208}\text{Pb}^*/^{232}\text{Th}$ age, so that

$$\frac{1}{\lambda_{238}} \text{Ln} \left(1 + \frac{\alpha_T - \alpha_0}{\alpha_T} r \right) = \frac{1}{\lambda_{235}} \text{Ln} \left(1 + \frac{\beta_T - \beta_0}{\beta_T} R \right)$$

Then define

$$r = \frac{^{206}\text{Pb}_T}{^{238}\text{U}} \quad s = \frac{^{207}\text{Pb}_T}{^{235}\text{U}} \quad f = \frac{^{206}\text{Pb}_0}{^{206}\text{Pb}_T} \quad F = \frac{^{207}\text{Pb}_0}{^{207}\text{Pb}_T}$$

$$\varphi_T = \left(\frac{^{207}\text{Pb}}{^{206}\text{Pb}} \right)_T \quad \varphi_0 = \left(\frac{^{207}\text{Pb}}{^{206}\text{Pb}} \right)_0 \quad u = \frac{^{238}\text{U}}{^{235}\text{U}} \quad p = 1 - f = \frac{r^*}{r}$$

$$r^* = (1 - f)r = pr = \frac{^{206}\text{Pb}^*}{^{238}\text{U}}$$

$$\begin{aligned} s^* &= \left(1 - \frac{\phi_0}{\phi_T} f \right) \phi_T u r \\ &= (\phi_T - f \phi_0) u r \\ &= \frac{^{207}\text{Pb}^*}{^{235}\text{U}} \end{aligned} \quad \begin{aligned} m_1 &= \lambda_{238} [1 + (1 - f)r] \\ &= \lambda_{238} (1 + r^*) = \lambda_{238} e^{\lambda_{238} t} \end{aligned} \quad \begin{aligned} m_2 &= \lambda_{235} [1 + (1 - F)s] = \lambda_{235} [1 + (\phi_T - \phi_0 f) u r] \\ &= \lambda_{235} (1 + s^*) \\ &= \lambda_{235} e^{\lambda_{235} t} \end{aligned}$$

$$(\phi_T - \phi_0 f) u = \frac{s^*}{r} \quad j_1 = \frac{p}{m_1} - \frac{s^*}{r m_2} \quad j_2 = \left(\frac{r}{m_1} - \frac{u r \varphi_0}{m_2} \right)^{-1}$$

where T indicates total, 0 indicates common, and t = age. Find the partial differentials

$$\langle dt \rangle = \frac{p \langle dr \rangle - r \langle df \rangle}{m_1} \quad \langle dr^* \rangle = p \langle dr \rangle - r \langle df \rangle$$

So that

$$\sigma_f^2 = j_2^2 \left[j_1^2 \sigma_r^2 + \left(\frac{u r}{m_2} \right)^2 (\sigma_{\phi_T}^2 + f^2 \sigma_{\phi_0}^2) \right] \quad \text{cov}(r, f) = j_1 j_2 \sigma_r^2 \quad \sigma_{r^*}^2 = m_1^2 \sigma_t^2$$

and

$$\begin{aligned} \sigma_t^2 &= \frac{p^2 \sigma_r^2 + r^2 \sigma_f^2 - 2 p r \text{cov}(r, f)}{m_1^2} \\ &= \frac{(p - 2 r j_1 j_2) p \sigma_r^2 + r^2 \sigma_f^2}{m_1^2} \end{aligned}$$

ignoring any covariance among the measured or common Pb-isotope ratios.

207-Corrected Ages and Ratios

Constrain the $^{206}\text{Pb}^*/^{238}\text{U}$ age to equal the $^{207}\text{Pb}^*/^{235}\text{U}$ age, so that

$$\frac{1}{\lambda_{238}} \text{Ln} \left(1 + \frac{\alpha_T - \alpha_0}{\alpha_T} r \right) = \frac{1}{\lambda_{235}} \text{Ln} \left(1 + \frac{\beta_T - \beta_0}{\beta_T} R \right)$$

Then define

$$\begin{aligned} r &= \frac{^{206}\text{Pb}_T}{^{238}\text{U}} & s &= \frac{^{207}\text{Pb}_T}{^{235}\text{U}} & f &= \frac{^{206}\text{Pb}_0}{^{206}\text{Pb}_T} & F &= \frac{^{207}\text{Pb}_0}{^{207}\text{Pb}_T} \\ \phi_T &= \left(\frac{^{207}\text{Pb}}{^{206}\text{Pb}} \right)_T & \phi_0 &= \left(\frac{^{207}\text{Pb}}{^{206}\text{Pb}} \right)_0 & u &= \frac{^{238}\text{U}}{^{235}\text{U}} & p &= 1 - f = \frac{r^*}{r} \\ r^* &= (1 - f)r = pr = \frac{^{206}\text{Pb}^*}{^{238}\text{U}} & s^* &= \left(1 - \frac{\phi_0}{\phi_T} f \right) \phi_T ur = (\phi_T - f\phi_0)ur = \frac{^{207}\text{Pb}^*}{^{235}\text{U}} \end{aligned}$$

$$\begin{aligned} m_1 &= \lambda_{238} [1 + (1 - f)r] & m_2 &= \lambda_{235} [1 + (1 - F)s] \\ &= \lambda_{238} (1 + r^*) & &= \lambda_{235} [1 + (\phi_T - \phi_0 f)ur] \\ &= \lambda_{238} e^{\lambda_{238} t} & &= \lambda_{235} (1 + s^*) = \lambda_{235} e^{\lambda_{235} t} \end{aligned}$$

where T indicates total, 0 indicates common, and t = age, and find the partial differentials

$$\langle df \rangle = j_2 \left(j_1 \langle dr \rangle - \frac{ur}{m_2} \langle d\phi_T \rangle \right) \quad \langle dr^* \rangle = p \langle dr \rangle - r \langle df \rangle \quad \langle dt \rangle = \frac{p \langle dr \rangle - r \langle df \rangle}{m_1}$$

from which follows

$$j_1 = \frac{p}{m_1} - \frac{s^*}{rm_2} \quad j_2 = \left(\frac{r}{m_1} - \frac{ur\phi_0}{m_2} \right)^{-1} \quad (\phi_T - \phi_0 f)u = \frac{s^*}{r} \quad \sigma_{r^*}^2 = m_1^2 \sigma_t^2$$

and

$$\sigma_t^2 = \frac{(p - j_1 j_2 r)^2 \sigma_r^2 + (uj_2 r^2 / m_2)^2 \sigma_f^2}{m_1^2}$$

ignoring any uncertainty in ϕ_0 and any cov

Equations for Isotope Geochronology

Simple Daughter/Parent ages ($r = e^{\lambda t} - 1$):
$$\sigma_t^2 = \frac{\sigma_r^2}{\lambda^2(1+r)^2} + t^2 \frac{\sigma_\lambda^2}{\lambda^2} .$$

$$\frac{\sigma_t^2}{t^2} = \left[\frac{r}{\lambda t(1+r)} \right] \frac{\sigma_r^2}{r^2} + \frac{\sigma_\lambda^2}{\lambda^2}$$

$^{207}\text{Pb}/^{206}\text{Pb}$ age-error as a function of $^{207}\text{Pb}/^{206}\text{Pb}$ ratio and of errors in decay constants:

$$\sigma_t^2 = \frac{(\mathbf{e}^{\lambda_{238}t} - 1)^2 \sigma_r^2 + (ute^{\lambda_{235}t})^2 \sigma_{\lambda_{235}}^2 + (rte^{\lambda_{238}t})^2 \sigma_{\lambda_{238}}^2}{(u\lambda_{235}\mathbf{e}^{\lambda_{235}t} - r\lambda_{238}\mathbf{e}^{\lambda_{238}t})^2}$$

where $r = ^{207}\text{Pb}/^{206}\text{Pb}$ $u = ^{235}\text{U}/^{238}\text{U}$

Y -error in conventional concordia curve at a specified X :
$$\sigma_Y^2 = \left[\sigma_{\lambda_{238}}^2 + \left(\frac{\lambda_{238}}{\lambda_{235}} \right)^2 \sigma_{\lambda_{235}}^2 \right] (te^{\lambda_{238}t})^2$$

X -error in conventional concordia curve at a specified Y :
$$\sigma_X^2 = \left[\sigma_{\lambda_{235}}^2 + \left(\frac{\lambda_{235}}{\lambda_{238}} \right)^2 \sigma_{\lambda_{238}}^2 \right] (te^{\lambda_{235}t})^2$$

X - and Y -error in conventional concordia curve at a specified t :
$$\sigma_X = te^{\lambda_{235}t} \sigma_{\lambda_{235}}$$

$$\sigma_Y = te^{\lambda_{238}t} \sigma_{\lambda_{238}}$$

Y -error in Tera-Wasserburg concordia curve at a specified X :

$$\sigma_Y^2 = (uXte^{\lambda_{235}t})^2 \left[\sigma_{\lambda_{235}}^2 + \left(\frac{\lambda_{235}}{\lambda_{238}} \right)^2 \sigma_{\lambda_{238}}^2 \right]$$

Equations for errors of common-Pb isotope ratios in zircon, sphene...

a = common ^{206}Pb , b = common ^{207}Pb , $w = ^{204}\text{Pb}$, b = blank, t = total, i = initial

α = common $^{206}\text{Pb}/^{204}\text{Pb}$, β = common $^{206}\text{Pb}/^{204}\text{Pb}$.

$$a_t = \frac{\alpha_b w_b + \alpha_i w_i}{w_t} = \frac{\alpha_b w_b + \alpha_i (w_t - w_b)}{w_t} = \frac{w_b}{w_t} \alpha_b + \alpha_i - \frac{w_b}{w_t} \alpha_i$$

$$d\alpha_t = \frac{w_b}{w_t} d\alpha_b + \frac{\alpha_b}{w_t} dw_b - \frac{\alpha_b w_b}{w_t^2} dw_t + \left(1 - \frac{w_b}{w_t}\right) d\alpha_i - \frac{\alpha_i}{w_t} dw_b + \frac{\alpha_i w_b}{w_t^2} dw_t$$

$$\begin{aligned} \sigma_{\alpha_t}^2 = & \left(\frac{\alpha_b - \alpha_i}{w_t} \right)^2 \sigma_{w_b}^2 + \left[(\alpha_i - \alpha_b) \frac{w_b}{w_t^2} \right]^2 \sigma_{w_t}^2 + \left(\frac{w_b}{w_t} \right)^2 \sigma_{\alpha_b}^2 + \left(1 - \frac{w_b}{w_t} \right)^2 \sigma_{\alpha_i}^2 \\ & + \left(\frac{\alpha_b - \alpha_i}{w_t} \right)^2 \left[\sigma_{w_b}^2 + \left(\frac{w_b}{w_t} \right)^2 \sigma_{w_t}^2 \right] + \left(\frac{w_b}{w_t} \right)^2 \left[\sigma_{\alpha_b}^2 + \left(1 - \frac{w_t}{w_b} \right)^2 \sigma_{\alpha_i}^2 \right] \end{aligned}$$

$$\begin{aligned} \text{cov}(\alpha_t, \beta_t) = & \frac{(\alpha_b - \alpha_i)(\beta_b - \beta_i)}{w_t^2} \sigma_{w_b}^2 + \frac{(\alpha_i - \alpha_b)(\beta_i - \beta_b)}{w_t^2} \left(\frac{w_b}{w_t} \right)^2 \sigma_{w_t}^2 \\ & + 2 \left(\frac{w_b}{w_t} \right)^2 \text{cov}(\alpha_b, \beta_b) + \left(1 - \frac{w_b}{w_t} \right)^2 \text{cov}(\alpha_i, \beta_i) \end{aligned}$$

Equations for Isotope Geochronology

Simple Daughter/Parent ages ($r = e^{\lambda t} - 1$):
$$\sigma_t^2 = \frac{\sigma_r^2}{\lambda^2(1+r)^2} + t^2 \frac{\sigma_\lambda^2}{\lambda^2} .$$

$$\frac{\sigma_t^2}{t^2} = \left[\frac{r}{\lambda t(1+r)} \right] \frac{\sigma_r^2}{r^2} + \frac{\sigma_\lambda^2}{\lambda^2}$$

$^{207}\text{Pb}/^{206}\text{Pb}$ age-error as a function of $^{207}\text{Pb}/^{206}\text{Pb}$ ratio and of errors in decay constants:

$$\sigma_t^2 = \frac{(e^{\lambda_{238}t} - 1)^2 \sigma_r^2 + (ute^{\lambda_{235}t})^2 \sigma_{\lambda_{235}}^2 + (rte^{\lambda_{238}t})^2 \sigma_{\lambda_{238}}^2}{(u\lambda_{235}e^{\lambda_{235}t} - r\lambda_{238}e^{\lambda_{238}t})^2}$$

where $r = ^{207}\text{Pb}/^{206}\text{Pb}$ $u = ^{235}\text{U}/^{238}\text{U}$

Y -error in conventional concordia curve at a specified X :
$$\sigma_Y^2 = \left[\sigma_{\lambda_{238}}^2 + \left(\frac{\lambda_{238}}{\lambda_{235}} \right)^2 \sigma_{\lambda_{235}}^2 \right] (te^{\lambda_{238}t})^2$$

X -error in conventional concordia curve at a specified Y :
$$\sigma_X^2 = \left[\sigma_{\lambda_{235}}^2 + \left(\frac{\lambda_{235}}{\lambda_{238}} \right)^2 \sigma_{\lambda_{238}}^2 \right] (te^{\lambda_{235}t})^2$$

X - and Y -error in conventional concordia curve at a specified t :
$$\sigma_X = te^{\lambda_{235}t} \sigma_{\lambda_{235}}$$

$$\sigma_Y = te^{\lambda_{238}t} \sigma_{\lambda_{238}}$$

Y -error in Tera-Wasserburg concordia curve at a specified X :

$$\sigma_Y^2 = (uXte^{\lambda_{235}t})^2 \left[\sigma_{\lambda_{235}}^2 + \left(\frac{\lambda_{235}}{\lambda_{238}} \right)^2 \sigma_{\lambda_{238}}^2 \right]$$

Miscellaneous Equations for Isotope Geochronology

$^{207}\text{Pb}/^{206}\text{Pb}$ age-error as a function of $^{207}\text{Pb}/^{206}\text{Pb}$ ratio and of errors in decay constants:

$$\sigma_t^2 = \frac{(e^{\lambda_{238}t} - 1)^2 \sigma_r^2 + (ute^{\lambda_{235}t})^2 \sigma_{\lambda_{235}}^2 + (rte^{\lambda_{238}t})^2 \sigma_{\lambda_{238}}^2}{(u\lambda_{235}e^{\lambda_{235}t} - r\lambda_{238}e^{\lambda_{238}t})^2}$$

where $r = ^{207}\text{Pb}/^{206}\text{Pb}$ $u = ^{235}\text{U}/^{238}\text{U}$

Y -error in conventional concordia curve at a specified X :
$$\sigma_Y^2 = \left[\sigma_{\lambda_{238}}^2 + \left(\frac{\lambda_{238}}{\lambda_{235}} \right)^2 \sigma_{\lambda_{235}}^2 \right] (te^{\lambda_{238}t})^2$$

X -error in conventional concordia curve at a specified Y :
$$\sigma_X^2 = \left[\sigma_{\lambda_{235}}^2 + \left(\frac{\lambda_{235}}{\lambda_{238}} \right)^2 \sigma_{\lambda_{238}}^2 \right] (te^{\lambda_{235}t})^2$$

X - and Y -error in conventional concordia curve at a specified t :
$$\sigma_X = te^{\lambda_{235}t} \sigma_{\lambda_{235}}$$

$$\sigma_Y = te^{\lambda_{238}t} \sigma_{\lambda_{238}}$$

Y -error in Tera-Wasserburg concordia curve at a specified X :

$$\sigma_Y^2 = (uXte^{\lambda_{235}t})^2 \left[\sigma_{\lambda_{235}}^2 + \left(\frac{\lambda_{235}}{\lambda_{238}} \right)^2 \sigma_{\lambda_{238}}^2 \right]$$

X -error in Tera-Wasserburg concordia curve at a specified Y :

$$\sigma_X^2 = \left[(uXte^{\lambda_{235}t})^2 \sigma_{\lambda_{235}}^2 + (1 + X\lambda_{238}e^{\lambda_{238}t})^2 \sigma_{\lambda_{238}}^2 \right] (X^2te^{\lambda_{238}t})^2$$

7

Equations for Isotope Geochronology

Simple Daughter/Parent ages ($r = e^{\lambda t} - 1$):
$$\sigma_t^2 = \frac{\sigma_r^2}{\lambda^2(1+r)^2} + t^2 \frac{\sigma_\lambda^2}{\lambda^2} .$$

$$\frac{\sigma_t^2}{t^2} = \left[\frac{r}{\lambda t(1+r)} \right] \frac{\sigma_r^2}{r^2} + \frac{\sigma_\lambda^2}{\lambda^2}$$

$^{207}\text{Pb}/^{206}\text{Pb}$ age-error as a function of $^{207}\text{Pb}/^{206}\text{Pb}$ ratio and of errors in decay constants:

$$\sigma_t^2 = \frac{(e^{\lambda_{238}t} - 1)^2 \sigma_r^2 + (ute^{\lambda_{235}t})^2 \sigma_{\lambda_{235}}^2 + (rte^{\lambda_{238}t})^2 \sigma_{\lambda_{238}}^2}{(u\lambda_{235}e^{\lambda_{235}t} - r\lambda_{238}e^{\lambda_{238}t})^2}$$

where $r = ^{207}\text{Pb}/^{206}\text{Pb}$ $u = ^{235}\text{U}/^{238}\text{U}$

Y -error in conventional concordia curve at a specified X :
$$\sigma_Y^2 = \left[\sigma_{\lambda_{238}}^2 + \left(\frac{\lambda_{238}}{\lambda_{235}} \right)^2 \sigma_{\lambda_{235}}^2 \right] (te^{\lambda_{238}t})^2$$

X -error in conventional concordia curve at a specified Y :
$$\sigma_X^2 = \left[\sigma_{\lambda_{235}}^2 + \left(\frac{\lambda_{235}}{\lambda_{238}} \right)^2 \sigma_{\lambda_{238}}^2 \right] (te^{\lambda_{235}t})^2$$

X - and Y -error in conventional concordia curve at a specified t :
$$\begin{aligned} \sigma_X &= te^{\lambda_{235}t} \sigma_{\lambda_{235}} \\ \sigma_Y &= te^{\lambda_{238}t} \sigma_{\lambda_{238}} \end{aligned}$$

Y -error in Tera-Wasserburg concordia curve at a specified X :
variance among the measured or common Pb-isotope ratios.

Appendix III: The *SQUID-2* Panels and Controls

Preferences (General panel)

The screenshot shows the 'Preferences' dialog box with the 'General' tab selected. The dialog has a title bar with a close button (X). Below the title bar are four tabs: 'General', 'Common Pb', 'Interpolation', and 'U-Pb Standards'. The 'General' tab contains several sections of settings:

- Show splash-screen at at startup:** A checked checkbox.
- Reduced-data workbook:** A group box containing five checkboxes:
 - ☐ Retain nonessential rows of PD or XML files (file length permitting)
 - ☐ Place data-reduction parameters in a separate worksheet
 - ☐ Place Task autocharts in a separate worksheet
 - ☐ Attach Task sheet to reduced-data workbook
 - ☒ Freeze header row and spot name column
- Special output columns (U/Pb only):** A group box containing two checked checkboxes:
 - ☒ 207-corrected 208Pb/232Th ages
 - ☒ 208-corrected Concordia-plot ratios and 207Pb/206Pb ages
- Chart insets:** A group box containing three checked checkboxes:
 - ☒ Create Concordia plots for U-Pb Age Standards (U-Pb geochron only)
 - ☒ Always create a worksheet containing the scan-by-scan isotope ratios
 - ☒ Set Y-axis minimum on isotope-ratio plots to zero
- Maximum permitted mass-difference (ppm) between measured and true nuclide masses to recognize 2 peaks as the same:** A text box containing the value '4000'.
- Default folder for PD and XML files:** A text box containing the path 'C:\Documents and Settings\KRL\Desktop\Squid Desktop Xla'.

On the right side of the dialog, there are four buttons: 'OK', 'Cancel', 'Help' (highlighted in purple), and 'Rebuild Task Catalog file'.

Show splash-screen at startup

Display the *SQUID-2* splash screen when the *SQUID-2* add-in is installed.

Retain nonessential rows of PD or XML files (file length permitting)

If this box is checked, the entire original PD or XML file will be retained in the (hidden) columns **A** of the condensed raw-data worksheet. Because this will slow construction of the condensed sheet somewhat, and slightly increase the possibility of out-of-memory problems under some conditions, do this only if you do not intent to archive the original PD or XML files.

Even if the box checked, however, *SQUID-2* will be unable to retain the original lines of some large XML files, simply because the number of lines of these files exceeds the 65,336-row limit of *Excel* 2003.

Place data-reduction parameters in a separate worksheet

Put data-reduction parameters specified in the PD file and Preferences on a separate worksheet attached to the reduced-data workbook.

Place auto-charts in a separate worksheet

Put auto-charts in a separate worksheet rather than in the reduced-data worksheet.

Attach Task sheet to reduced-data workbook

Attach a copy of the *Task* sheet to the reduced-data workbook.

Freeze header rows and spot name column

Forces the first column (spot names) and the column-header row to always be visible when the output worksheets are scrolled.

207-corrected 208Pb/232Th ages

Check this box if you require an output-sheet column containing $^{208}\text{Pb}/^{232}\text{Th}$ ages corrected for common Pb with the assumption that the $^{206}\text{Pb}/^{235}\text{U}$ - $^{207}\text{Pb}/^{235}\text{U}$ ages are concordant.

208-corrected Concordia-plot ratios and 207Pb/206Pb ages

Check this box if you require output-sheet columns containing the $^{206}\text{Pb}/^{235}\text{U}$ - $^{207}\text{Pb}/^{235}\text{U}$ ages, errors, and error correlation where the common-Pb correction assumes concordance of the $^{206}\text{Pb}/^{235}\text{U}$ - $^{208}\text{Pb}/^{232}\text{Th}$ ages. Checking this box also enables the calculation of 207-corrected $^{207}\text{Pb}/^{206}\text{Pb}$ ages.

Create concordia plots for U-Pb Age Standards (U-Pb only)

If checked, a small conventional concordia plot for the Standard spots will be placed within the Standard worksheet.

Always construct within-scan isotope ratio sheet

A within-scan isotope ratio sheet contains the within-spot (i.e. scan-by-scan) isotope ratio and *Task* equation results for each spot, formatted for clarity of viewing. This worksheet is required by *SQUID* for construction of chart insets of isotope ratios/*Task* equations **Error! Bookmark not defined.** versus time.

Set Y-axis minimum on isotope-ratio plots to zero

If unchecked, the Y-axis (the ratio values) for isotope-ratio chart-insets will be auto-scaled to the plotted data.

Maximum permitted mass-difference between measured and true nuclide masses...

If the apparent (PD file value) and true masses of a mass station are within this difference, *SQUID-2* will accept the measurement. Otherwise, *SQUID-2* will not recognize the mass station as simply a measurement error on the true mass, so that the Setup panel will refuse to reduced that PD file.

Default folder for PD and XML files

The folder *SQUID-2* opens when asking for a raw-data file to process.

Rebuild Task Catalog file

Reconstruct the TaskCatalog file (contains basic information about each valid Task file in the SquidUser folder, used to determine compatibility for data reduction with a raw-data file).

Preferences (Common Pb panel)

Preferences

General Common Pb Interpolation U-Pb Standards

Isotope ratios of default Common Pb

Use Stacey-Kramers ratios for the apparent age of the sample, using the age-type and common-Pb index-isotope specified at the time of data reduction. ☒

OK Cancel Help

Use Stacey-Kramers ratios for the age specified at right ☐ S-K Age (ma)

Use the ratios entered below ☐

206Pb/204Pb

207Pb/204Pb

208Pb/204Pb

207Pb/206Pb

208Pb/206Pb

Isotope Ratios of Common Pb (U/Pb geochron only)

Refers to the Pb subtracted from the total Pb to obtain the isotopic abundances of the spot's radiogenic Pb. These corrections will be applied in the StandardData and SampleData worksheets. The common Pb used for Grouped-sample worksheets is assigned separately.

Use Stacey-Kramers ratios for apparent age of sample

Use the single-stage evolution Stacey-Kramers ratios (Stacey and Kramers, 1975) for the sample's age using the age-type specified for data-reduction..

Use Stacey-Kramers ratios for the age specified at right

Use the Stacey-Kramers ratios corresponding to the specific age entered in the S-K Age= box.

Use the ratios entered below

Enter an arbitrary set of Pb-isotope ratios to use for common-Pb correction.

Preferences (Interpolation panel)

Preferences

General | Common Pb | **Interpolation** | U-Pb Standards

Calculation of within spot, scan-by-scan means

	By weighted average	By linear regression interpolated to mid-time
for Special U/Pb/Th Task equations (recommended)	☐	☒
for other Task Equations	☐	☒
for isotope ratios of the same element	☒	☐
for isotope ratios of different elements	☐	☒

U/Pb calibration constant

Minimum acceptable external error in U/Pb and Th/Pb ratios (1-sigma%)

☒ Inhibit initial auto-rejection of Standard-spot calibr. constants

☒ Correct for secular drift (not applicable to Tasks that calculate calibration constants for both U/Pb and Th/Pb)

Secular drift smoothing window (#spots)

Fixed ☒ #spots in fixed window

Automatic ☐ An automatic smoothing-window finds the smallest window giving an external error larger than the minimum external error

OK Cancel Help

Calculation of within-spot, scan-by-scan means

The mean isotope ratios or Task Equation values assigned to each analyzed spot can be done either by

1. Calculating the outlier-resistant mean of the $N-1$ interpolations for the N peak-scans, or
2. Using an outlier-resistant linear regression of the interpolated values to determine the value of the interpolations at mid-burn time,

Because in many (most?) analyses ratios or equations involving different elements seem to increase or decrease in a pseudo-linear manner during the spot burn, selecting the Linear Interpolation option generally results in better “internal” (within spot) precision. It can be argued that for consistency, other ratios or equations should be treated similarly.

U/Pb calibration constant

MINIMUM ACCEPTABLE EXTERNAL ERROR IN U/Pb AND Th/Pb RATIOS (1-SIGMA%)

SQUID-2 always calculates the most-probable external (that is, in addition to counting-statistics error) spot-to-spot error for the calibration constants of a U/Pb Task. If the Age Standard yields rather low count-rates for Pb, statistical resolution of any external error is poor, so that an apparent external error of zero is possible. If you believe that some external error is always present in SHRIMP U/Pb analyses, you may wish to specify a lower bound for the external error, regardless of the statistical calculation. The external error is included in the error assigned to sample-spot U/Pb ages, so is important to consider carefully.

INHIBIT INITIAL AUTO-REJECTION OF STANDARD-SPOT CALIBRATION CONSTANTS

If this box is unchecked, *SQUID-2* will attempt to reject obvious outliers (subject to later user revision) before calculating the weighted-mean calibration constant. Check the box to inhibit such rejection.

CORRECT FOR SECULAR DRIFT

If the calibration constant displays an *obvious* drift or variation with time (whether monotonic or cyclic), you can correct the sample spots for the observed drift by clicking this box.

SECULAR DRIFT SMOOTHING WINDOW (#SPOTS)

The secular-drift curve is calculated via the LOWESS algorithm (Cleveland, 1979), which fits an outlier-resistant smoothing-spline for a moving spot-window. The degree of smoothing is controlled by the size of the spot-window: the larger the window the greater the smoothing. Large windows are best for, say, a pseudolinear drift, smaller windows for drifts with abrupt changes or cyclic appearance. *SQUID-2* offers two choices for determining the width of the smoothing window: fixed (you specify the window) and automatic (*SQUID-2* determines the best window).

If you select Fixed, the initial smoothing window will be the one you specify here. *You can always change the smoothing window, however, after the data have been reduced.*

If you select Automatic, *SQUID-2* will find the largest smoothing window whose resulting external error is not less than the specified Minimum acceptable external error.

This Secular Drift Correction algorithm is experimental, should be used with care and trusted only after considerable comparison with conventional SQUID-2 data reduction.

Preferences (U-Pb Standards panel)

You can specify up to 50 common U-Pb age and/or concentration standards from this panel. When the U-Pb data-reduction panel is being filled out, *SQUID* will look for spot names whose first few characters match those of an age standard. If a match is found, that standard's name, Pb/U age (if defined), Pb/Pb age (if defined), and U or Th concentration (if defined) will be entered into the appropriate boxes.

Index #	Name (1st chars)	Pb/U Age (Ma)	Pb/Pb Age (Ma)	U conc (ppm)	Th conc (ppm)
1	FC1	1099			
2	TEM	416.8			
3	R33	419			
4	Q	1842	1851.6		
5	BR266	1041			
6	FC1	1099			
7	CZ3	579		555	
8	AS3	1095	1099		
9	AS57	1095	1099		
10	SL13	572		238	
11	BS-1	508		300	
12	MG-1	490		1050	
13	NY/PK	1010		15000	
14	CAP	275.6		100	
15	8153	513.1	479	2250	
16	OG1	3441	3465.4		
17	44069	424		2500	20000
18	PHAL	2060		100	
19	r3	419	419		
20					

Enter only enough initial characters to distinguish the Standard name from other spot names. Non-alphanumeric characters are ignored.

Leave the "Age" field blank if the standard will be used only for uranium or thorium concentration, and vice versa.

Click an Index# to insert a Standard, Shift-Click to delete.

When adding Standards to the *U-Pb Standards* panel, keep in mind that *SQUID* will search the spot names for a Standard in the order that they appear in the panel.

Reduce Data for U-Pb Geochronology

Run table box

PD file (click to change)

__JNA_A225_050502095.45.pd

Task: Zircon, canonical anu

Zircon, canonical anu peaks	PD file peaks
196	196
204	203.97
204.1	204.02
206	206
207	207
208	208
238	238
248	248
254	254

STANDARDS

Number of characters to show in Standard names: 3

Pb/U or Pb/Th Age

Name: r33

Pb/U or Pb/Th Age: 419

207Pb/206Pb age (same as Pb/U if left blank): 419

for 204-overct correction only

U or Th concentration

☐ No U or Th standard

Name: sl1

ppm U: 238

☒ U concentration

☐ Th concentration

Index isotope for common-Pb subtraction: 204Pb

☒ Normalize ion beams to SBM signal

Spot subsets for User Parameters: No Subsets

Time Window: No Time Window -- process all spots

☐ Standard spots only

☐ Make trim-mass charts

version 2.30.08.11.28 K.R. Ludwig, BGC, 28 Nov 2008

Task

Select any of the existing *Tasks* to use for data reduction.

New PD file

Load another PD file without leaving this panel.

Number of characters to show in Standard names

The number of initial spot-name characters that will be shown in the Name dropdowns, and that will be used to categorize a spot as a Standard.

Pb/U or Pb/Th Age

NAME

The initial spot-name characters that define a spot as a Standard.

PB/U OR PB/TH AGE

The *apparent* age of the Standard. Though this is generally the same as the true age, it can be different if the Standard is known, say, to have a consistent amount of Pb losses.

207PB/206PB AGE

By default (i.e. left blank), this is the same as the Pb/U or Pb/Th age. However, if, say, the Standard's Pb/U or Pb/Th age is discordant with its $^{207}\text{Pb}/^{206}\text{Pb}$ age, a different value can be entered here. *The Standard $^{207}\text{Pb}/^{206}\text{Pb}$ age is used only for 204-overcount calculation.*

U OR TH CONCENTRATION

The average concentration of either U or Th in the concentration standard

NAME

Select the initial characters of the U (or Th) concentration standard names. If no such standard was present in the mount, leave blank (the first item in the dropdown list).

PPM U (OR TH)

Enter the average concentration of U (or Th) for the concentration standard. Leave blank if there is no concentration standard

Index Isotope for common Pb

The common-Pb correction for radiogenic U-Pb isotope ratios can be corrected using any of ^{204}Pb (corrects assuming the default common Pb specified in *Preferences*), ^{207}Pb (corrects using the assumption of $^{206}\text{Pb}/^{238}\text{U}$ - $^{207}\text{Pb}/^{206}\text{Pb}$ age-concordance), or ^{208}Pb (corrects using the assumption of $^{206}\text{Pb}/^{238}\text{U}$ - $^{208}\text{Pb}/^{232}\text{Th}$ age-concordance).

Normalize ion beams to SBM signal

Ion beams normalized to the signal of a perfectly reliable SBM (secondary beam monitor) should be unaffected by any instability in the secondary ion beam. If there are doubts as to the functioning of the SBM (i.e. the SBM is not constant over the mass range of the Run Table), leave unchecked.

Spot subsets for user parameters

Shows those spots for which one or more *Task* equations are specifically applied.

Time window

Shows the time window within which spot analyses are selected for data reduction. Other spots are ignored.

Trim mass charts

Construct trim-mass charts for each peak-centered mass station.

Task Editor

Invokes the *Task* Editor panels, returning to this Setup panel when done.

Preferences

Invokes the Preferences panel, returning to this Setup panel when done.

Change subsets

Specify *Task* equation that are to be applied only to a selected group of spots, based in the initial characters of the spot names.

Change Time-window

Define a time window within which spot analyses are selected for data reduction. Other spots are ignored.

Reduce General Isotope Data

230Th-U ages	PD file
peaks	peaks
224	196
234	203.97
238	204.02
246.03	206
246.06	207
248	208
251	238
254	248
	254

New PD

Load another PD file, then return to this panel.

Task

Select the *Task* for data reduction from the dropdown list of all defined *Tasks*.

Normalize to SBM signal

Ion beams normalized to the SBM (secondary beam monitor may be unaffected by any instability in the secondary ion beam. If there are doubts as to the functioning of the SBM (i.e. the SBM is not constant over the mass range of the Run Table), leave unchecked.

Make trim-mass charts

Construct trim-mass charts for each peak-centered mass station.

Spot subsets for user parameters

Shows spots having *Task* equations that are specifically applied.

Process spots within Time Window?

The time window within which spot analyses are selected for data reduction. Other spots are ignored.

Change Subsets

Specify *Task* equation that are to be applied only to a selected group of spots, based in the initial characters of the spot names.

***Task* Editor**

Invokes the *Task* Editor panels, returning to this Setup panel when done.

Preferences

Invokes the Preferences panel, returning to this Setup panel when done.

Change Time-window

Define a time window for spot data reduction.

Sample Grouping

Process Samples as a Group

Samples

Spot Names: B1-, D1-, SM1, ZRD, ,
 Clear: ☐ ☐ ☐ ☐ ☐ ☐
 #characters to show in Spot-names drop down: 3

☒ Extract groups based on spot names, and place in separate worksheets, or:
☐ Put all spots together in a single worksheet
☒ Extract coherent age from group(s)
☐ Sort worksheet rows by age

Disregard

☒ Case Aa
☒ Spaces " "
☒ Dashes -
☒ Slashes /
☒ Periods .
☒ Commas ,
☒ Colons :
☒ Semicolons ;

Age Groups

Type of age

☐ 204-corr. 206Pb/238U
☒ 207-corr. 206Pb/238U
☐ 208-corr. 206Pb/238U
☐ Total 206Pb/238U
☐ 204-corr. 207Pb/206Pb
☐ 204-corr. 208Pb/232Th
☐ 207-corr. 208Pb/232Th
☐ 208-corr. 207Pb/206Pb

Grouping criteria

Min. probability of fit: 0.05
 Min. fract. of spots: 0.60

Correction for apparent 204-overcounts on Standard

☒ No correction -- use measured 204Pb
☐ Force concordance of 207Pb/235U - 206Pb/238U ages
☐ Force concordance of 208Pb/232Th - 206Pb/238U ages

☐ Construct concordia plots for groups if possible
☐ Compare UO/U (or UO2/U or UO2/UO) of Grouped spots with Std

Common Pb

Use Stacey-Kramers ratios for the age type selected at left ☐

OR

Enter specific ratios ☒

18.3 206/204
 0.8536 207/206
 2.093 208/206

OK Cancel Help

Spot Names

From any of the 6 spot-name dropdowns, select from 1 to 6 spot name fragments to use as a selection criterion for each group. Only the initial number of characters shown in the dropdown will be used.

letters to show in Spot-names dropdown

Select how many of the initial spot-name characters will appear in the Spot Names dropdowns.

Disregard Case, Spaces, Dashes, Slashes...

Select characters that SQUID is to ignore when grouping the spot names.

Extract coherent age group

After grouping, extract a statistically coherent age grouping, if possible.

Sort samples by age

After grouping, sort the data alphabetically by spot name.

Group all spots together

Treat all spots in the Sample worksheet as a single group.

Age Groups: Type of U-Pb Age

The type of age to be used for spot age-grouping.

Total $^{206}\text{Pb}/^{238}\text{U}$

The $^{206}\text{Pb}/^{238}\text{U}$ age, with no correction for common Pb.

204-corr. $^{206}\text{Pb}/^{238}\text{U}$

Use the $^{206}\text{Pb}/^{238}\text{U}$ age, using ^{204}Pb as the index isotope for common-Pb correction.

207-corr. $^{206}\text{Pb}/^{238}\text{U}$

Use the $^{206}\text{Pb}/^{238}\text{U}$ age, with ^{207}Pb as the index isotope for common-Pb correction, and assuming $^{206}\text{Pb}/^{238}\text{U}$ - $^{207}\text{Pb}/^{206}\text{Pb}$ age-concordance.

208-corr. $^{206}\text{Pb}/^{238}\text{U}$

Use the $^{206}\text{Pb}/^{238}\text{U}$ age, with ^{208}Pb as the index isotope for common-Pb correction, and assuming $^{206}\text{Pb}/^{238}\text{U}$ - $^{208}\text{Pb}/^{232}\text{Th}$ age-concordance.

204-corr. $^{208}\text{Pb}/^{232}\text{Th}$

Use the $^{208}\text{Pb}/^{232}\text{Th}$ age, with ^{204}Pb as the index isotope for common-Pb correction.

204-corr. $^{207}\text{Pb}/^{206}\text{Pb}$

Use the $^{207}\text{Pb}/^{206}\text{Pb}$ age, with ^{204}Pb as the index isotope for common-Pb correction.

207-corr. $^{208}\text{Pb}/^{232}\text{Th}$

Use the $^{208}\text{Pb}/^{232}\text{Th}$ age, with ^{207}Pb as the index isotope for common-Pb correction.

208-corr. $^{207}\text{Pb}/^{206}\text{Pb}$

Use the $^{207}\text{Pb}/^{206}\text{Pb}$ age, with ^{207}Pb as the index isotope for common-Pb correction.

Age Groups: Grouping Criteria*Min. prob. of fit*

The minimum acceptable probability-of-fit that defines an age group.

Min. fract. of spots

The minimum fraction of spots in the group that defines an age group.

Correction for 204 overcounts on Standard

How to correct for apparent overcounts at the ^{204}Pb mass position, using calculations from the Standard spots as a guide.

No correction – use measured 204Pb

Assume that there are no 204 overcounts

Force concordance of 207Pb/235U - 206Pb/238U ages

Subtract the mean 204 overcts/sec. (from 207) calculated in the Standard worksheet from each of the Sample spot ^{204}Pb measurements.

Force concordance of 208Pb/232Th - 206Pb/238U ages

Subtract the mean 204 overcts/sec. (from 208) calculated in the Standard worksheet from each of the Sample spot ^{204}Pb measurements.

Common Pb*Use Stacey-Kramers ratios for the age type specified at left*

Correct each spot's U-Pb, Pb-Pb, and Th-Pb age assuming common-Pb isotope ratios calculated using the (single stage evolution) Stacey-Kramers method (Stacey and Kramers, 1975), for the age-type specified in the *Type of age* box.

Enter specific ratios

Correct each spot's U-Pb, Pb-Pb, and Th-Pb age with the common-Pb isotope ratios specified in Preferences

Construct concordia plot for groups if possible

If an age group is requested and found, construct a concordia plot inset with the grouped worksheet.

Compare UO/U (or UO2/U or UO2/UO) of Grouped spots with std

Construct a chart-inset in each Grouped Spots sheet showing the $\text{UO}_n / \text{UO}_m$ ratio for the Standard and (in a different color) the Grouped spots. Values of n and m can range from 0 to 2, and depend on what was measured for the spots.

Time Window

Limits to specify with Standard spots

Lower (first) only

Ignore spot-analyses acquired before the Starting spot, process all spots after.

Upper (last) only

Ignore spot-analyses acquired after the Ending spot, process all spots before.

Both lower and upper

Ignore spot-analyses acquired before the Starting spot and after the ending spot.

No window (reduce all spots)

Cancel the Time Window.

Starting Name/Time

The first spot to process. For U-Pb geochronology, restricted to age standards.

Ending Name/Time

The last spot to process. For U-Pb geochronology, restricted to age standards.

Standard spots that bracket the window

Limits to specify with Standard spots:

☐ Lower (first) only
☐ Upper (last) only
☒ Both lower and upper
☐ No window (reduce all spots)

	Name	Time
Starting	r33-1.1	2005-05-02, 10:28
Ending	r33-24.1	2005-05-03, 09:04

Figure 29: The Time Window panel.

Subsets

See instructions in the blue box at the top-left of the form (below).

Scope of user parameters (spots for which equations will be applied) [X]

To restrict the application of any equation to a subgroup of spots that share the first few letters of their names:

- 1) Specify (using the "# of initial..." spinner) the number of initial spot-name characters that are shared by the subgroup (for example, 2 for all spots whose names start with, say, SW).
- 2) Select (from the "Trimmed spot-names") the first N letters that define the subgroup (say, SW), where N is the number specified in the spinner box above.
- 3) Click on the yellow box corresponding to the name of the equation whose application is to be restricted only those spots whose names start with, say, SW. Click the box again to cancel the action.

SUBSET EQUATIONS DO NOT APPLY TO AGE-STANDARD SPOTS

Name of equation to apply to spot-name subset	Initial spot-name characters
Log UO/U	
Log Pb/U {q5}	q5
Expo	

Number of initial letters to show in spot-names drop-down

▲ ▼ 2

Trimmed spot-names

q5 ▼

OK
Clear All
Cancel

Copy or Rename a Task

Rename and Copy a task

Task name
230Th-U ages

Task description
Metal and oxides

Primary mineral
opal

Creator name or initials
krl

Next Retain panel changes, go to Run Table

Return Retain panel changes, return to Task Editor

Cancel Discard panel changes, return to Task Editor

Task name

Enter a name for the *Task* that distinguishes it from other *Tasks*.

Task description

Enter a brief description of the *Task*.

Primary mineral

The mineral that to which the *Task* will most-commonly be assigned.

Task Editor

Task Editor

Purpose

☒ U/Pb Geochronology

☐ Isotope Ratio

Edit Constants

Refresh Task Catalog

OK Save all changes, then exit

Cancel Discard all changes, then exit

Save Save changes to selected Task, don't Exit

Help Help for this panel

Existing Tasks (from Task Catalog)

New **Copy** **Delete**

AC-xenotime, PbUO2-MH
 AC-xenotime-AC9
 Monazite, floating slope, ThO test
 Titanite, floating slope
 Z5 = 9pk GA (207-206 normalised via O
 Zircon, canonical anu +UO2
 Zircon, canonical anu with floating slope
 Zircon, canonical anu

Edit/view a Task panel

Name, descr...
 ↓
 Run Table
 ↓
 Isotope Ratios
 ↓
 Special U-Pb Eqns
 ↓
 General Eqns
 ↓
 Auto Charts

Pb/UO2

Description

xenotime
 Mineral

AJC
 Creator

version 2.21.08.03.16
 K.R. Ludwig, BGC, 16 Mar 2008

U-Pb Geochronology

Edit/view the U/Pb geochronology *Tasks*.

Isotope Ratio

Edit/view the General isotope ratio *Tasks*.

Existing Tasks*NEW*

Define a completely new Task.

RENAME OR COPY

Rename or copy the selected *Task*.

DELETE

Remove the selected *Task* (cannot be undone!).

SINGLE CLICK A TASK NAME

Selects the Task for editing by any of the colored buttons to the right.

DOUBLE CLICK A TASK NAME

Selects the Task and opens its General Equations panel.

Edit/view a panel for the selected Task

Click a colored button to edit or view the specified panel for the selected *Task*.

Refresh Task Catalog

Rebuild the *Task* Catalog by abstracting the name, description, mineral, and creator for each *Task* in the *SquidUser* folder.

OK

Save all changes and exit

Cancel

Discard all changes and exit.

Constants

Edit/view the user-defined constants.

Save

Save all changes but remain in the *Task* Editor.

Ionic species

Clear

Specify mass station numerically

bkrd

ref pk

cps col

Construct a column in the processed-data worksheet for counts/second at this mass station?

Hide?

Hide the *CPS* column in the reduced-data worksheet.

Periodic table of the elements

Commonly-used elements for geochronology are shown in bold font, short-lived elements (not selectable) in gray background, and rare earth elements in green boxes.

Next

Save any changes to the Run Table panel, then go to the next (isotope ratios) *Task* Editor panel.

Return

Save any changes to the Run Table panel, then return to the main *Task* Editor panel.

Cancel

Discard any changes to the Run Table panel, then return to the main *Task* Editor panel.

Insert

Insert a row of blank boxes (nuclide, true mass, nominal mass) before the Active Nuclide.

Delete

Delete the Active Nuclide and associated true- and nominal-mass boxes.

Clear one

Clear the Active Nuclide and associated true- and nominal-mass boxes.

Clear all

Clear all of the nuclide, true mass, and nominal mass boxes.

The Isotope Selection Sub-Panel

How many?

How many atoms of the selected isotope (blue) are in the Active Nuclide? Note that 1) all stable and long-lived isotopes of the selected element are shown, 2) the default isotope will be the most-abundant isotope, and 3) the relative isotopic abundances (summed to 100) of the natural element are shown at the left of the inset.

Charge

Specify the charge of the active nuclide.

Enter

Enter this isotope in the active nuclide box then add another element or isotope.

Complete

Add this isotope to the active nuclide, then move down to the next mass station.

The screenshot shows a sub-panel titled "ISOTOPES". On the left, a list of isotopes for lead (Pb) is displayed with their relative abundances: 1.37% for ²⁰⁴Pb, 25.15% for ²⁰⁶Pb, 21.11% for ²⁰⁷Pb, and 52.38% for ²⁰⁸Pb. Each isotope is highlighted in a blue box. To the right of the list are three buttons: "Enter" (blue), "Complete" (red), and "Cancel!" (grey). Below the list, there are two spinners: one for "How many 206Pb ?" with a value of 1, and another for "Charge?" with a value of +1. To the right of the "Enter" button, there is a text instruction: "Then select another element to add to the active nuclide". To the right of the "Complete" button, there is a text instruction: "Add, then proceed to next nuclide".

%abundance	Isotope
1.37	²⁰⁴ Pb
25.15	²⁰⁶ Pb
21.11	²⁰⁷ Pb
52.38	²⁰⁸ Pb

Buttons: Enter, Complete, Cancel!

Spinners: How many 206Pb ? (1), Charge? (+1)

Instructions: Then select another element to add to the active nuclide; Add, then proceed to next nuclide

Isotope Ratios

Step 2 -- specify isotope ratios

Click on a blue ratio-box to activate or clear it., then click on the green mass-boxes (numerator then denominator) to create the ratio.

Mass Stations		Ratios (from mass stations at left)	
196	204	A 204/206	B 207/206
204.1	206	C 208/206	D 206/238
207	208	E 208/248	F 254/238
238	248	G 248/254	H 238/196
254	270	I	J
		K	L
		M	N
		O	P
		Q	R
		S	T
		U	V
		W	X
		Y	Z
		AA	AB

Next Save changes, go to U-Pb eqns
Back Discard panel changes, go back to Run Table
Return Retain panel changes, return to Task Editor
Cancel Discard panel changes, return to Task Editor
Help

Ratios (from mass stations at left)

To specify the calculation of an isotope ratio, 1) click on one of the blue isotope-ratio boxes to activate it (turns tan), 2) click on one of the green boxes (to left) to specify the numerator isotope of the ratio, 3) click on the denominator isotope. To clear ratio box, click it a second time.

Next

Save any changes to the Isotope Ratios panel, then go to the next panel (equations).

Back

Save any changes to the Isotope Ratios panel, then return to the main *Task Editor* panel.

Cancel

Discard any changes to the Isotope Ratios panel and return to the main *Task Editor* panel.

U-Pb Special Equations

Equations specifically for U-Th-Pb geochronology

Age of primary interest
☐ 206Pb/238U ☒ 208Pb/232Th

Calculation of 208Pb/232Th
☐ Indirectly, from 232Th/238U ☒ Directly, from the age standard

U-Th-Pb equations

Equation for 208Pb/232Th normalization

$$\left[\frac{^{208}\text{Pb}}{^{232}\text{Th}} \right] / \left[\frac{^{208}\text{Pb}}{^{232}\text{Th}} \right]_{\text{age standard}} \cdot \exp(-\lambda_{232}t)$$

Equation for 206Pb/238U normalization

$$\left[\frac{^{206}\text{Pb}}{^{238}\text{U}} \right] / \left[\frac{^{206}\text{Pb}}{^{238}\text{U}} \right]_{\text{age standard}} \cdot \exp(-\lambda_{238}t)$$

Equation for 232Th/238U

Equation for thorium concentration

* cannot refer to equations defined in next panel

Ratios available

A	204/206	B	207/206
C	208/206	D	238/199
E	206/238	F	254/238
G	248/254	H	208/232
I	248/199	J	248/232
K	208/248	L	
M		N	
O		P	
Q		R	
S		T	

Click a blue ratio-box to insert its letter-index into the active equation

AGE OF PRIMARY INTEREST

206Pb/238U

Use $^{206}\text{Pb}/^{238}\text{U}$ ages for normalization to the age standard.

208Pb/232Th

Use $^{208}\text{Pb}/^{232}\text{Th}$ ages for normalization to the age standard.

CALCULATION OF 208PB/232TH

Indirectly, from $^{232}\text{Th}/^{238}\text{U}$

Determine $^{208}\text{Pb}/^{232}\text{Th}$ ratios by multiplying the calculated (age-standard normalized) $^{206}\text{Pb}/^{238}\text{U}$ by the measured $^{208}\text{Pb}/^{206}\text{Pb}$ and the $^{232}\text{Th}/^{238}\text{U}$ calculated from the defined $^{232}\text{Th}/^{238}\text{U}$ equation.

Directly, from the age standard

Determine $^{208}\text{Pb}/^{232}\text{Th}$ ratios in the same way as the $^{206}\text{Pb}/^{238}\text{U}$ ratios (by normalizing the results of the $^{208}\text{Pb}/^{232}\text{Th}$ equation to that of the age standard).

Equation for $^{206}\text{Pb}/^{238}\text{U}$ const.

An expression derived from the measured isotope ratios that is (ideally) constant for a given $^{206}\text{Pb}/^{238}\text{U}$.

Equation for $^{208}\text{Pb}/^{232}\text{Th}$ const.

An expression derived from the measured isotope ratios that is (ideally) constant for a given $^{208}\text{Pb}/^{232}\text{Th}$.

 $^{232}\text{Th}/^{238}\text{U}$ equation

An expression derived from the measured isotope ratios that yields the true $^{232}\text{Th}/^{238}\text{U}$.

Equation for ppm U

An expression derived from the measured isotope ratios that is (ideally) constant for a given uranium concentration.

Next

Save any changes to the U-Pb equations panel, then go to the next panel (equations).

Back

Save any changes to the U-Pb equations panel, then go back to the previous (Isotope Ratios) panel.

Return

Save any changes to the U-Pb equations panel, then return to the main *Task* Editor panel.

Cancel

Discard any changes to the U-Pb equations panel and return to the main *Task* Editor panel.

Help

Display detailed information about *Task* Equations and Switches.

Append constant

Select (and/or edit) a constant from the dropdown list of defined User Constants and append to the active equation.

Append LITERAL ratios

When an isotope-ratio box is clicked, append the literal ratio to the active equation rather than its number.

Append INDEX of ratios

When an isotope-ratio box is clicked, append its index letter to the active equation rather than the literal ratio.

General Equations

User-defined equations

Next/previous equation: [Spinner]

Switches: ST SA SC LA FO NU HI AR AR#

Equation names: [Field]

Index #	Equations	Switches	Equation names
1	$\text{Ln}(["270/254"])$	[Switches]	Ln UO2/UO
2	$\text{Ln}(["206/254"])$	[Switches]	Ln Pb/UO
3	$\text{Ln}(["206/270"])$	[Switches]	Ln Pb/UO2
4	$\text{Ln}(["264/232"])$	[Switches]	Ln ThO2/Th
5	$\text{Ln}(["208/264"])$	[Switches]	Ln Pb/ThO2
6	$\text{RobReg}(["\text{Ln UO2/UO}"], ["\text{Ln Pb/UO2}"])$	[Switches]	Expo_U
7	$\text{robreg}(["\text{lntho2/th}"], ["\text{lnpb/uo2}"])$	[Switches]	Expo_u_2
8	$\text{robreg}(["\text{lnuo2uo}"], ["\text{lnpbtho2}"])$	[Switches]	Expo_Th
9	$\text{RobReg}(["\text{Ln ThO2/Th}"], ["\text{Ln Pb/tho2}"])$	[Switches]	Expo_Th_2
10	$(["208\text{Pb}/206\text{Pb}"] * ["264/254"])/\text{StdRad86fact}$	[Switches]	fThU
11	$\text{average}(["\text{fThU}"])$	[Switches]	meanfThU
12	$((\exp(["\text{age}(\text{Ma})"] * 1.55125e-4) - 1) * ["\text{ppm U}"] * 0.9928 * 0.866$	[Switches]	ppm206/std
13	$(["\text{ppm206std}"] * ["208\text{Pb}/206\text{Pb}"] * 1.0097$	[Switches]	ppm208/std
14	$(["207\text{Pb}/206\text{Pb}"] * ["\text{ppm206std}"] * 1.0049$	[Switches]	ppm207/std
15	$(["\text{ppm206std}"] + ["\text{ppm207std}"] + ["\text{ppm208std}"]$	[Switches]	totPb/std

Isotope ratios:

A	89/203	202/203	B
C	232/203	254/203	D
E	264/203	270/203	F
G	204/206	204.1/206	H
I	207/206	208/206	J
K	264/232	206/254	L
M	264/254	270/254	N
O	208/264	206/270	P
Q	202/270	203/270	R
S	202/89	203/89	T
U			V
W			X
Y			Z
AA			AB
AC			AD

Buttons: Next, Return, Cancel, Help

Equations: Insert, Delete Active, Clear Active, Clear All

To append an Equation reference: Click its index # (left of equation)

Append a Column-Header reference

Append Eqns and Ratios as LITERAL

Append a Constant

To append an Isotope Ratio reference: Click one of the blue boxes above.

To refer to cell A1 in another workbook: {WorkbookName}SheetName!A1

To refer to cell A1 in the Sample worksheet from the Standard worksheet: 'SampleData'!A1

Next/previous equation (spinner)

Click on one of the spinner arrows to move the active equation up or down.

ST, SA, SC, LA, FO, NU, HI, AR

Sets/clears the Standard, Sample, SingleCell, Last, Formula, Numeric, Hidden, or Array formula switches.

AR #rows #cols

If the AR switch is set, enter the number of rows and columns of the formula output, separated by a space.

Constants

Click to invoke the user-defined constants dropdown.

Append column-header reference

Click to invoke the column-headers dropdown.

Append LITERAL equations/ratios

Click to append the equation name (rather than its index#) and literal isotope-ratio (rather than its index letter) when an equation-number or isotope-ratio box is clicked.

Next

Save any changes to the Equations panel, then go to the next panel (Auto charts).

Back

Save any changes to the Equations panel, then go back to the previous (U-Pb equations or Iso-tope Ratios) panel.

Return

Save any changes to the Equations panel, then return to the main *Task* Editor panel.

Cancel

Discard any changes to the Equations panel and return to the main *Task* Editor panel.

Delete (active equation)

Delete the active equation and shift all higher-number equations down.

Insert (active equation)

Insert a blank equation box before the active equation.

Clear Active

Clear the active equation.

Clear All

Clear all equations.

Help

Display detailed information about *Task* Equations and Switches. Auto Charts

Column headers

Select the column headers of the x and y parameters to chart. The x parameter is not needed if a Y-mean calculation is specified.

Axis scaling*Autoscale*

Let *SQUID* determine the bounds and tick intervals of the specified axis.

Starts at 0

Set the lower bounds of the specified axis at zero (*SQUID* determines the upper bound).

Log

Use logarithmic scaling for the specified axis.

Calculations*X-Y REGRESSION*

Determine and plot an outlier-resistant best-fit line to the data.

None

Don't perform any calculations.

Y mean

Calculate and plot the outlier-resistant mean of the y values and make *hours* the x axis.

Back

Save any changes to the Auto Charts panel, then go back to the previous (Equations) panel.


Return

Save any changes to the Auto Charts panel, then return to the main *Task* Editor panel.

Cancel

Discard any changes to the Auto Charts panel and return to the main *Task* Editor panel.

Constants



Edit

Change the name and/or value of the selected constant.

Add new

Add a new constant.

Delete

Delete the selected constant.

Exit

Close the Constants panel.

Append to Eqn 1

Place the constant at the end of the active equation .

Example Tasks

Zircon

CANONICAL ANU, 9 PEAKS

Run Table	
Run order	Name (e.g. 90Zr2 16O)
1	90Zr2 16O + ref
2	204Pb +
3	204.1 + bkrd
4	206Pb +
5	207Pb +
6	208Pb +
7	238U +
8	232Th 16O +
9	238U 16O +

Task Equations

Index#		ST	SA	SC	LA	FO	NU	AR	Equation names
1	$\log([f])$						•		Log UO/U
2	$\log([d])$						•		Log Pb/U

Chart #	Column headers	Axis scaling	Calculations
1	X Log UO/U Y Log Pb/U	☒ Autoscale ☐ Starts at 0 ☐ Log ☒ Autoscale ☐ Starts at 0 ☐ Log	☐ None ☒ X-Y regression ☐ Y mean
Clear			

Auto Charts

ZIRCON U-Pb WITH DIRECT Pb/Th

(other panels as Canonical Zircon U-Pb)

Equations specifically for U-Th-Pb geochronology

Age of primary interest

☒ 206Pb/238U ☐ 208Pb/232Th

Calculation of 208Pb/232Th

☐ Indirectly, from 232Th/238U ☒ Directly, from the the age standard

U-Th-Pb equations

Equation for 206Pb/238U normalization

$[e]/[f]^2$

Equation for 208Pb/232Th normalization

$[h]/[f]$

Equation for 232Th/238U

Equation for uranium concentration

$[d]/[f]^{0.66}$

Ratios available

A	204/206	B	207/206
C	208/206	D	238/196
E	206/238	F	254/238
G	248/254	H	208/248
I		J	
K		L	
M		N	
O		P	
Q		R	

Isotope Ratios and Special U-Pb Equations

ZIRCON WITH “FLOATING SLOPE”

(other panels as Canonical Zircon U-Pb)

Calculation of $^{208}\text{Pb}/^{232}\text{Th}$

☒ Indirectly, from $^{232}\text{Th}/^{238}\text{U}$
☐ Directly, from the the age standard

U-Th-Pb equations

Equation for $^{206}\text{Pb}/^{238}\text{U}$ const.

$[d]/[f]^{["\text{Expo}"]}$

Equation for $^{208}\text{Pb}/^{232}\text{Th}$ const.

$^{232}\text{Th}/^{238}\text{U}$ equation

$[g] \cdot (0.0345 + 0.868 \cdot [f])$

Equation for ppm U

$[h]/[f]^{0.66}$

Ratios available

A	$^{204}\text{Pb}/^{238}\text{U}$	B	$^{207}\text{Pb}/^{235}\text{U}$
C	$^{208}\text{Pb}/^{232}\text{Th}$	D	$^{206}\text{Pb}/^{238}\text{U}$
E	$^{208}\text{Pb}/^{238}\text{U}$	F	$^{254}\text{Pb}/^{238}\text{U}$
G	$^{248}\text{Pb}/^{238}\text{U}$	H	$^{238}\text{Pb}/^{238}\text{U}$
K		L	
M		N	
O		P	
Q		R	

Step 4 -- Define user parameters

Index#	Equations	Switches	Equation names
		ST SA SC LA FO NU AR ARc	
1	$\text{Ln}([f])$	[ST] [SA] [SC] [LA] [FO] [NU] [AR] [ARc]	Log UO/U
2	$\text{Ln}([d])$	[ST] [SA] [SC] [LA] [FO] [NU] [AR] [ARc]	Log Pb/U
3	$\text{min}(3, \text{max}(1, \text{roreg}([1], [2])))$	[ST] [SA] [SC] [LA] [FO] [NU] [AR] [ARc]	Expo

^{230}Th ages, insignificant ^{232}Th

Scan order	Ionic species (e.g. 90Zr2 160) or numeric mass	True mass (m/e)	Nominal mass (m/e)
1	224 + ref	224	224
2	234U +	234.041	234
3	238U +	238.051	238
4	230Th 16O +	246.028	246.03
5	246.08 + bkrd	246.08	246.08
6	232Th 16O +	248.033	248
7	235U 16O +	251.039	251
8	238U 16O +	254.046	254
9			
10			

Next/previous equation

Switches

SC

FO

NU

AR

ARrc

Index#	Equations	Equation names
1	$^{234}\text{Th}/^{238}\text{U} \times \frac{\lambda_{234\text{ka}}}{\lambda_{238\text{ka}}}$	U234/U238 AR
2	$\frac{[^{246}\text{Th}/^{238}\text{U}] \times \lambda_{230\text{ka}}}{\lambda_{232\text{ka}}}$	Th230/Th232 AR
3	$\frac{[^{246}\text{Th}/^{254}\text{U}] \times \lambda_{230\text{ka}}}{\lambda_{238\text{ka}}}$	Th230/U238 AR
4	average($^{234}\text{Th}/^{238}\text{U}$)	Std Mean 234 /238 {SecEquilStd}
5	average($^{246}\text{Th}/^{254}\text{U}$)	Std Mean 230 /238 {SecEquilStd}
6	$\frac{[^{234}\text{Th}/^{238}\text{U}]}{[\text{StdMean}^{234}/^{238}\text{U}]}$	Std-Norm U234/U238 AR
7	$\frac{[^{246}\text{Th}/^{254}\text{U}]}{[\text{StdMean}^{230}/^{238}\text{U}]}$	Std-Norm Th230/U238 AR
8	Th230AgeAndInitial([3],[=e],[1],[=f],,true)	Th230 RawAge = 234/238(0) =
9	Th230AgeAndInitial([7],[=e],[6],[=f],,true)	Th230 NormAge = 234/238(0) =
10	$-100 \times (1/[^{251}\text{Th}/^{254}\text{U}]/137.88-1)/3$	Mass Fract /amu
11		

Isotope ratios

A	254/224	246.03/248	B
C	251/254	254/238	D
E	246.03/254	234/238	F
G			H
I			J
K			L
M			N
O			P
Q			R
S			T
U			V

Simple Monazite

cps col	ref pk	bkrd	Ionic species (e.g. 90Zr2 16O) or numeric mass	True mass (m/e)	Nominal mass (m/e)	
1	☒	☐	☐	140Ce 31P 16O2 +	202.869	203
2	☒	☐	☐	204Pb +	203.973	204
3	☒	☐	☒	204.1 +	204.1	204.1
4	☒	☐	☐	206Pb +	205.974	206
5	☐	☐	☐	207Pb +	206.976	207
6	☐	☐	☐	208Pb +	207.977	208
7	☐	☐	☐	232Th +	232.038	232
8	☐	☐	☐	238U +	238.051	238
9	☐	☐	☐	232Th 16O +	248.033	248
10	☐	☐	☐	238U 16O +	254.046	254

U-Th-Pb equations

Equation for 206Pb/238U normalization

$["206/238"] / ["254/238"]^{["expopbu"]}$

Equation for 208Pb/232Th normalization

$["208/232"] / ["248/232"]^{["expopbth"]}$

Equation for 232Th/238U

Equation for uranium concentration

$["238/203"] / ["254/238"]^{0.66}$

* cannot refer to equations defined in next panel

A	204/206	B	207/206
C	208/206	D	238/203
E	206/238	F	254/238
G	248/254	H	208/248
I	248/232	J	208/232

Equations	Switches								Equation names
	ST	SA	SC	LA	FO	NU	AR	ARc	
$\text{Ln}(["254/238"])$	•				•				Ln UO/U
$\text{Ln}(["206/238"])$	•				•				Ln Pb/U
$\text{Ln}(["248/232"])$	•				•				Ln ThO/Th
$\text{Ln}(["208/232"])$	•				•				Ln Pb/Th
$\text{fixed}(\text{robreg}(["\text{LnUO/U}"], ["\text{LnPb/U}"]), 2)$	•		•		•				Expo PbU
$\text{fixed}(\text{robreg}(["\text{LnThO/Th}"], ["\text{LnPb/Th}"]), 2)$	•		•		•				Expo PbTh

Appendix IV: Weighted Average Algorithm for U/Pb Age Standards

- Reject points whose difference (R_i) between the internal error and median internal error is greater than $5 \times N_{mad}$, where

$$N_{mad} = 1.4826 \left(1 + \frac{5}{N-2} \right) \sqrt{\text{Median}(R)}$$

and N is the number of Age Standard spots.

- Calculate the weighted average of the remaining points in the usual way (weighted by internal errors only), e.g.

$$\bar{x} = \frac{\sum (x_i / {}^{int}\sigma_{x_i})^2}{\sum 1 / {}^{int}\sigma_{x_i}^2}, \quad \sigma_{\bar{x}} = \sqrt{1 / \sum {}^{int}\sigma_{x_i}^2},$$

where ${}^{int}\sigma_{x_i}$ is the "internal" error of the spot calibration constant.

- If the probability-of-fit of the above weighted average is less than 0.05 (as calculated from the χ^2 distribution of the $MSWD$ about $N-2$), calculate another weighted average, but this time assuming that a constant (but unknown) "external" error also exists for each data-point. The Maximum Likelihood Estimator of both the mean calibration constant and the "external" error, $\sigma_{ext.}$, are determined by simultaneously solving

$$\bar{x} = \frac{\sum x_i (\sigma_{x_i}^2 + \sigma_{ext.}^2)^{-1}}{\sum (\sigma_{x_i}^2 + \sigma_{ext.}^2)^{-1}} \quad \text{and} \quad \sum \frac{(x_i - \bar{x})^2}{(\sigma_{x_i}^2 + \sigma_{ext.}^2)^2} = \sum \frac{1}{\sigma_{x_i}^2 + \sigma_{ext.}^2}.$$

Apparent outliers are then rejected if:

The wtd residual of the point is less than $Toler$, where p is the number of the rejection "pass", and $E_{\bar{x}}$ is the 95%-conf. error of $\bar{x}$ as determined via (3), above, and

$$Toler = \left(1 + \frac{p-1}{4} \right) \sqrt{4(\sigma_{int}^2 + \sigma_{ext.}^2) + E_{\bar{x}}^2},$$

- The number of un-rejected points is no fewer than 85% of the original number, and
- The data point is nonzero.

Steps 3-4 are repeated until no rejections occur. Note that with each "pass" through steps 3-4 the bar for rejection ($Toler$) is raised.

Note that the rejection algorithm is *not at all rigorous*. Its intent is merely to offer the user an example of possible true outliers, and that the user should only accept the results of the above rejection algorithm if he or she would select the same rejections independently, preferably on the basis of objective criteria such as $^{204}\text{Pb}/^{206}\text{Pb}$, U or Th concentration, *et cetera*.

Appendix V: Some Useful *Excel* and *Isoplot* Functions

Native Excel Functions

AVERAGE

Argument: Values range

Output: The simple average of the numbers.

MEDIAN

Argument: Values range

Output: The median of the numbers.

STDEV

Argument: Values range

Output: The standard deviation of the population. Divide by $\sqrt{n}$ for the standard error of the mean.

LINEST

Arguments: Y-range, X-range [, include stats]

The x- and y-ranges are usually 2 columns of numbers; add the optional *IncludeStats* if slope- and intercept-errors are desired in the output.

Output: Parameters for a simple y on x linear regression. Typically, from a 1×1 up to 3×2 row/column output-range containing

Slope	Intercept
$\pm$ slope	$\pm$ inter
Rho	

where the errors ($\pm$) are 1 standard error, and rho is the x-y correlation coefficient.

INTERCEPT

Arguments: Y-range, X-range

Output: The intercept of a simple y on x linear regression.

SLOPE

Arguments: Y-range, X-range

Output: The slope of a simple y on x linear regression.

MAX, MIN

Arguments: Values range

Output: The maximum or minimum value of the range.

Isoplot and SQUID Functions

sqBiWeight

Arguments: RangeIn, [, Tuning], [, Mean only = False], [, All comers = False]

RangeIn contains the input values in a single column;

Tuning (optional) must be 6 or 9, controls the degree of outlier resistance, 6 being more outlier-resistant than 9;

Mean only (optional) If False (default) returns the mean and standard error in two (1×2 row/column) cells; If True, returns the mean only, in one cell.

All comers (optional) If False (default), accepts only non-blank cells with numeric contents without the strikethrough font attribute; if True, accepts all non-blank, numeric cells in the input range.

Output: The Tukey's biweight (Hoaglin *et al*, 1973) outlier-resistant mean. More efficient than the median for distributions that are fundamentally Gaussian. For SHRIMP data, BiWt is almost always preferable to the simple average.

ChiSquare

Arguments: MSWD, Degrees of freedom

MSWD is the mean square of weighted deviates of a weighted mean or regression (see p. 37);

Degrees of freedom the degrees of freedom for the calculation (= $N-1$ for a weighted average, $N-2$ for a linear regression).

Output: The probability of fit (probability that the specified *MSWD* or greater will occur for the specified degrees of freedom).

Concordia

Arguments: $^{207}\text{Pb}/^{235}\text{U}$, $^{207}\text{Pb}/^{235}\text{U}$ err, $^{206}\text{Pb}/^{238}\text{U}$, $^{206}\text{Pb}/^{238}\text{U}$ err [, $^{207}\text{Pb}/^{235}\text{U} - ^{206}\text{Pb}/^{238}\text{U}$ err. corr.] [, with decay-const errs = False] [, %errors = False] [, σ -level of errors]

Output: The *Concordia Age* (Ludwig, 1998) for the specified U/Pb ratios and errors.

ConcordiaTW

Arguments: $^{238}\text{U}/^{206}\text{Pb}$, $^{238}\text{U}/^{206}\text{Pb}$ err, $^{207}\text{Pb}/^{206}\text{Pb}$, $^{207}\text{Pb}/^{206}\text{Pb}$ error [, $^{238}\text{U}/^{206}\text{Pb} - ^{207}\text{Pb}/^{206}\text{Pb}$ err. corr.]

Output: The *Concordia Age* (Ludwig, 1998) for the specified U/Pb ratios and errors.

DisEq68Age

Arguments: $^{206}\text{Pb}/^{238}\text{U}$, initial $^{234}\text{U}/^{238}\text{U}$ activity ratio, initial $^{230}\text{Th}/^{238}\text{U}$ activity ratio

Output: The corresponding $^{206}\text{Pb}/^{238}\text{U}$ age.

DisEq68Ratio

Arguments: Age (ma), initial $^{234}\text{U}/^{238}\text{U}$ activity ratio, initial $^{230}\text{Th}/^{238}\text{U}$ activity ratio

Output: The corresponding $^{206}\text{Pb}/^{238}\text{U}$ ratio.

DisEq75Age

Arguments: $^{207}\text{Pb}/^{235}\text{U}$, initial $^{231}\text{Pa}/^{235}\text{U}$ activity ratio

Output: The corresponding $^{207}\text{Pb}/^{235}\text{U}$ age.

DisEq75Ratio

Arguments: Age (ma), $^{207}\text{Pb}/^{235}\text{U}$, initial $^{231}\text{Pa}/^{235}\text{U}$ activity ratio

Output: The corresponding $^{207}\text{Pb}/^{235}\text{U}$ ratio.

DisEq76Ratio

Arguments: Age (ma), initial $^{234}\text{U}/^{238}\text{U}$ activity ratio, initial $^{230}\text{Th}/^{238}\text{U}$ activity ratio, initial $^{231}\text{Pa}/^{235}\text{U}$ activity ratio

Output: The corresponding $^{207}\text{Pb}/^{206}\text{Pb}$ ratio.

DisEqPbPbAge

Arguments: $^{207}\text{Pb}/^{206}\text{Pb}$ ratio, initial $^{234}\text{U}/^{238}\text{U}$ activity ratio, initial $^{230}\text{Th}/^{238}\text{U}$ activity ratio, initial $^{231}\text{Pa}/^{235}\text{U}$ activity ratio.

Output: The corresponding $^{207}\text{Pb}/^{235}\text{U}$ age.

Drnd

Arguments: Number The number to be rounded;
SigFigs The number of significant figures to return.

Output: The number rounded to the specified number of significant figures

ExternalSigma

Arguments: ValuesAndErrors (2-column range), [PercentErrs = False], [SigmaLevel = 1]

Output: The external 1-sigma error required to account for the presence of extra scatter (that is, not accounted for by the assigned errors) of the Values.

HalfLife

Arguments: IntegralNuclideMass [, InMyr]

Output: The half life of the nuclide in years (default) or millions of years.

InitU234U238

Arguments: Age (Ka), meas. $^{234}\text{U}/^{238}\text{U}$ activity ratio

Output: The corresponding initial $^{234}\text{U}/^{238}\text{U}$ activity ratio.

La230, La232, La234, La235, La238

Arguments: None

Output: The decay constant, in years⁻¹, of ^{230}Th , ^{232}Th , ^{234}U , ^{235}U , or ^{238}U , respectively

Lambda

Argument: Integral nuclide mass

Output: The decay constant, in years⁻¹, of the indicated nuclide.

Log10

Argument: Value

Output: The corresponding base-10 logarithm.

Mad

Argument: Numeric range

Output: The median absolute deviation from the median of the range.

MSWD

Arguments: ValuesAndErrors (2-column range), [PercentErrs = False], [SigmaLevel = 1]

Output: The *MSWD* (mean square of weighted deviates) parameter for the weighted average of ValuesAndErrors.

NukeMass

Arguments: Nuclide [, significant figures]

Output: Mass of the nuclide ($^{12}\text{C}=12.000$)

Pb76

Argument: Age (ma)

Output: The radiogenic $^{207}\text{Pb}/^{206}\text{Pb}$ for that age

ProbFit

Arguments: ValuesAndErrors (range), [PercentErrs = False], [SigmaLevel = 1]

Output: The probability that the assigned errors will yield at least the amount of observed scatter of the Values.

RobReg

Arguments: x-range, y-range [, with y-intercept = False] [, include errors = False]

Output: The results of a robust (resistant) linear regression of a suite of x-y points as a 1×2 to 1×4 cell range containing the line's Slope, (optional) Slope error (95%-conf.), (optional) Y-intercept, (optional) Y-intercept error. Output of the optional parameters depends on the values of *with y-intercept* and *include errors*.

SingleStagePbT

Arguments: $^{206}\text{Pb}/^{204}\text{Pb}$, $^{207}\text{Pb}/^{204}\text{Pb}$

Output: The Stacey-Kramers single-stage Pb-growth curve age corresponding to the input values

SingleStagePbR

Arguments: Age (Ma), Which ratio

Output: The Stacey-Kramers single-stage Pb-growth curve isotope ratio corresponding to the input values. The isotope ratio returned is $^{206}\text{Pb}/^{204}\text{Pb}$, $^{207}\text{Pb}/^{204}\text{Pb}$, or $^{208}\text{Pb}/^{204}\text{Pb}$ for *Which ratio* = 0, 1, or 2, respectively.

sqWtdAv

Arguments: Values range, Errors range [, PercentErrsIn=False] [, PercentErrsOut=False] [, CanReject=False]

Output: Up to a 7×2 cell range where the first column contains the weighted mean, 1σ (prob. fit >0.05) or 68.3%-confidence (prob. fit ≤ 0.05) error in the weighted mean, external error for the data points (if prob. fit ≤ 0.05), probability of fit, number of data-points rejected, index-numbers of rejected data-points (if non-zero). The second column contains explanatory labels for the values in the first column.

StudentsT

Arguments: Degrees of freedom [, confidence limit = 95]

Output: The Students's-t value for the specified (or default) confidence limit.

Th230Age

Arguments: $^{230}\text{Th}/^{238}\text{U}$ activity ratio, $^{234}\text{U}/^{238}\text{U}$ activity ratio

Output: The $^{230}\text{Th}/\text{U}$ age (Ka) for the specified activity ratios.

Th230AgeAndInitial

Arguments: $^{230}\text{Th}/^{238}\text{U}$, $^{230}\text{Th}/^{238}\text{U}$ err, $^{234}\text{U}/^{238}\text{U}$, $^{234}\text{U}/^{238}\text{U}$ err [, $^{230}\text{Th}/^{238}\text{U}$ - $^{234}\text{U}/^{238}\text{U}$ err-correl] [, errors as %] [, σ -level] [, with decay-const errs,] [as atomic ratios])

Output: A 1×5 cell range containing the $^{230}\text{Th}/\text{U}$ age (Ka), age error, initial $^{234}\text{U}/^{238}\text{U}$ activity ratio, error in $(^{234}\text{U}/^{238}\text{U})_0$, error correlation between age and initial ratio.

Th230U238AR

Arguments: Age (ka), initial $^{234}\text{U}/^{238}\text{U}$ activity ratio

Output: The $^{230}\text{Th}/^{238}\text{U}$ activity ratio for the specified age and initial $(^{234}\text{U}/^{238}\text{U})$.

ThUfromPb86

Arguments: Age (ma), meas. $^{208}\text{Pb}/^{206}\text{Pb}$, meas. $^{206}\text{Pb}/^{204}\text{Pb}$, initial $^{206}\text{Pb}/^{204}\text{Pb}$, initial $^{208}\text{Pb}/^{204}\text{Pb}$

Output: The $^{232}\text{Th}/^{238}\text{U}$ corresponding to the input ratios, calculated assuming that the system's $^{206}\text{Pb}/^{238}\text{U}$ and $^{208}\text{Pb}/^{232}\text{Th}$ ages are concordant.

U234age

Arguments: Present-day $^{234}\text{U}/^{238}\text{U}$ activity ratio , Initial $^{234}\text{U}/^{238}\text{U}$ activity ratio

Output: The $^{234}\text{U}/^{238}\text{U}$ age, in Ka, for the input ratios

U234ageAndErr

Arguments: $^{234}\text{U}/^{238}\text{U}$ activity ratio, $^{234}\text{U}/^{238}\text{U}$ err, Initial $^{234}\text{U}/^{238}\text{U}$ activity ratio, error [, as %errors]

Output: The $^{234}\text{U}/^{238}\text{U}$ age (Ka) and error, in a 1×2 cell-range.

U234U238AR

Arguments: Age (ka), initial $^{234}\text{U}/^{238}\text{U}$ activity ratio

Output: The present-day $^{234}\text{U}/^{238}\text{U}$ activity ratio

YorkInter

Arguments: xyValues&Errors range [, errors in as %] [, σ -level]

Output: The Y -intercept of a 2-error linear regression. The *xyValues&Errors range* must comprise a 4- or 5-column range, where the columns contain x , x -error, y , y -error, x - y error correlation (the last column is optional, and assumed to be zero if missing).

YorkInterErr95

Arguments: *xyValues&Errors range* [, errors in as %] [, σ -level]

Output: The 95%-conf. Y -intercept of a 2-error linear regression. The *xyValues&Errors range* must comprise a 4- or 5-column range, where the columns contain x , x -error, y , y -error, x - y error correlation (the last column is optional, and assumed to be zero if missing).

YorkMSWD

Arguments: *xyValues&Errors range* [, errors in as %] [, σ -level]

Output: The $MSWD$ of a 2-error linear regression. The *xyValues&Errors range* must comprise a 4- or 5-column range, where the columns contain x , x -error, y , y -error, x - y error correlation (the last column is optional, and assumed to be zero if missing).

YorkProb

Arguments: *xyValues&Errors range* [, errors in as %] [, σ -level]

Output: The probability of fit of a 2-error linear regression. The *xyValues&Errors range* must comprise a 4- or 5-column range, where the columns contain x , x -error, y , y -error, x - y error correlation (the last column is optional, and assumed to be zero if missing).

YorkSlope

Arguments: *xyValues&Errors range* [, errors in as %] [, σ -level]

Output: The slope of a 2-error linear regression. The *xyValues&Errors range* must comprise a 4- or 5-column range, where the columns contain x , x -error, y , y -error, x - y error correlation (the last column is optional, and assumed to be zero if missing).

YorkSlopeErr95

Arguments: *xyValues&Errors range* [, errors in as %] [, σ -level]

Output: The 95%-conf. slope-error of a 2-error linear regression. The *xyValues&Errors range* must comprise a 4- or 5-column range, where the columns contain x , x -error, y , y -error, x - y error correlation (the last column is optional, and assumed to be zero if missing).

Appendix VI: Useful *Isoplot* Tools for *SQUID*

- Concordia plots
- *X-Y* plots
- Cumulative Gaussian plots
- Histograms
- Linearized Gaussian plots
- Multicomponent deconvolution (Sambridge & Compston's *MIX*)
- Weighted averages, Tukey's Biweight means, Medians
- $(^{230}\text{Th}/^{238}\text{U}) - (^{234}\text{U}/^{238}\text{U})$ plots, including evolution curves and isochrons.
- Ages of stacked beds (Bayesian constraints).
- Youngest detrital zircon

Appendix VII: Definitions of Terms

204 corrected

The amount of common Pb to subtract is calculated assuming the isotope ratios assigned in Preferences, using 204 as the index isotope.

207 corrected

The true $^{206}\text{Pb}/^{204}\text{Pb}$ is calculated by assuming that the true radiogenic-Pb isotope ratios yield concordant $^{206}\text{Pb}/^{238}\text{U} - ^{207}\text{Pb}/^{235}\text{U}$ ages, determining the $^{206}\text{Pb}/^{204}\text{Pb}$ that produces this result, then subtracting common-Pb using 204 as the index isotope and the isotope ratios assigned in Preferences.

208 corrected

The true $^{206}\text{Pb}/^{204}\text{Pb}$ is calculated by assuming that the true radiogenic-Pb isotope ratios yield concordant $^{206}\text{Pb}/^{238}\text{U} - ^{208}\text{Pb}/^{232}\text{Th}$ ages, determining the $^{206}\text{Pb}/^{204}\text{Pb}$ that produces this result, then subtracting common-Pb using 204 as the index isotope and the isotope ratios assigned in Preferences.

Array function or formula

See p. 25.

Common Pb

For the purposes of *SQUID*, the Pb present at the time the mineral to be dated formed. Also comprises environmental Pb has contaminated the grain.

Concordance

The state of having equivalent $^{206}\text{Pb}/^{238}\text{U} - ^{207}\text{Pb}/^{235}\text{U} - ^{207}\text{Pb}/^{206}\text{Pb}$ apparent ages.

Concordia

The locus of concordant ages on the $^{206}\text{Pb}_{\text{rad}}/^{238}\text{U} - ^{207}\text{Pb}_{\text{rad}}/^{235}\text{U}$ or $^{238}\text{U}/^{206}\text{Pb} - ^{207}\text{Pb}/^{206}\text{Pb}$ diagram.

Concordia age

The most probable age of a system having a arguably concordant U-Pb ages.

Constants

Constants defined and named by the *SQUID* user and available for use in *Task* equations.

Discordance

The absence of equivalent $^{206}\text{Pb}/^{238}\text{U} - ^{207}\text{Pb}/^{235}\text{U} - ^{207}\text{Pb}/^{206}\text{Pb}$ apparent ages for a U/Pb system.

Dropdown

An *Excel* or *SQUID* menu of items from which the user can choose.

Error correlation

A statistical measure of how much two parameters covary. An error correlation of 1 implies a precisely proportional straight-line relationship between two parameters, zero implies complete independence of one parameter with the other, -1 implies a precisely proportional negative variation of one parameter with respect to another.

The formula for the error correlation is:
$$\frac{\sum (x_i - \bar{x})(y_i - \bar{y})}{\sqrt{\sum (x_i - \bar{x})^2 \sum (y_i - \bar{y})^2}}$$

When referring to errors, the error correlation reflects the likelihood that the error in one parameter's measurement is in the same direction (too high or too low) as that of the other parameter.

External error

A component of measurement error that is constant from sample to sample but not comprised by the assigned analytical error.

Grouping

to be added

Internal error

The analytical error assigned to a measurement.

Ionic species

A particular molecular combination of elements and isotopes with a particular charge.

Isoplot

to be added

m/e

Mass/charge of an ionic species, often loosely substituted by “mass”.

MSWD

Mean Square of Weighted Deviates (see p. 37).

Normalizing equation

An equation that yields a constant value when a parameter of interest (Pb/U , $[U]$, ...) is constant, so that the ratio of normalizing equations for an unknown sample to a standard equals the ratio of the unknown's parameter value to the standard parameter value.

Panel

A *Excel* dialog box that allows the user to make one or more choices.

PD file**Probability of fit**

The probability that the observed weighted-mean or weighted-regression *MSWD* will be exceeded if the assigned data-point errors are entirely responsible for the observed scatter. See p. 37.

Radiogenic Pb

Pb produced *in situ* by the decay of U and Th in a mineral or rock.

Range

A suite of worksheet cells (usually contiguous) in an *Excel* worksheet.

Resistant (mean, regression...)

A mean or regression calculated using an algorithm that is insensitive to outliers in the data.

Robust (mean, regression...)

See above.

Sample sheet

The worksheet containing the reduced SHRIMP data for all spots but the U-Pb age standard, and named "SampleData".

SBM

Stands for Secondary Beam Monitor, a collector that measures the secondary ion beam at a point in the flight tube and ion optics such that, in favorable cases, essentially the entire secondary ion beam is meas..

Scan order

The order of a mass station in the peak-jump (mass-jump) sequence of a SHRIMP analysis.

Single cell formula

A *Task* equation that operates on a columnar range of data, whose output occupies a single cell, rather than a data column.

Spot

A single spot on a mineral grain, used for one SHRIMP analysis and assigned a unique name as part of the analysis.

Stacey-Kramers

Refers to the single-stage Pb isotope growth curve presented by Stacey and Kramers (1975). Predicts the average crustal Pb-isotope ratios at a specified age.

Standard sheet

The worksheet containing the reduced SHRIMP data for the U-Pb age standard, and named "StandardData".

Trim mass reference

The mass station to be used as a reference when constructing charts for the trim masses used for each spot.

References

- Dodson M. H., 1978, A linear method for second-degree interpolation in cyclical data collection: J. Phys. E: Sci. Instrum. **11** 296.
- Cleveland, W.S., 1979, Robust, locally weighted regression and smoothing scatterplots: J. Amer. Statist. Assn. **74**, 829-836.
- Hoaglin, D.C., Mosteller, F., and Tukey, J.W., 1983, *Understanding Robust and Exploratory Data Analysis*: John Wiley and Sons, 345-349.
- Ludwig, K.R., 1998, On the treatment of concordant uranium-lead ages: Geochim. Cosmochim. Acta **62**, 665-676.
- Ludwig, K.R., 2001, SQUID 1.02: *A User's Manual*, Berkeley Geochr. Ctr. Spec. Pub. 2.
- Ludwig, K.R., 2003, Isoplot 3: A Geochronological Toolkit for Microsoft *Excel* (rev. August 27, 2003), Berkeley Geochr. Ctr. Spec. Pub. 4.
- McIntyre, G.A., Brooks, C., Compston, W., and Turek, A., 1966, The statistical assessment of Rb-Sr isochrons: J. Geophys. Res. **71**, 5459-5468
- Stacey, J.S., and Kramers, J.D., 1975, Approximation of terrestrial lead isotope evolution by a two-stage model: Earth Planet. Sci. Lett. **26**, 207-221.
- Titterton, D.M., and Halliday, A.N., 1979, On the fitting of parallel isochrons and the method of maximum likelihood: Chem. Geol. **26** 183-195.
- Vugrinovich, R.G., 1981, J. Math. Geol. **13**, 443-454.

Index

- 204-corr. 206Pb/238U..... 63
- 207-corr. 206Pb/238U..... 63
- ²⁰⁷Pb/²⁰⁶Pb age 5, 59, 63
- 208-corr. 206Pb/238U..... 63
- Age Standard..... 5, 8, 9, 31
- AR switch..... 25, 27, 76
- HI switch 25
- ARrc switch 23, 24, 27
- Auto Charts 30, 77, 78
- background box..... 18
- Calibration constant, U/Pb, correcting for
 - drift of 10
- Celeron..... 2
- charge..... 72
- charge, nuclide, specifying 18
- chart insets 14, 53
- coherent age group..... 9, 63
- column swapping in U-Pb Tasks 28
- common Pb 63
- Common Pb 95
- concordance 59, 63, 64
- Concordance 95
- Concordia..... 8, 88
- counts/second columns, specifying..... 18
- CPS column 18
- CPS columns, hide..... 71
- CPS column 18
- Crossover 2
- decay constants 21
- Drift of U/Pb calibration constant, correcting
 - for..... 10
- Equation panels..... 21, 25
- Equation switches 23, 75, 77
- Equation switches 23
- Equations, Task..... 20
- Excel 2, 6, 14, 22, 24, 25, 29, 30, 34, 36, 37,
 - 86, 96, 97, 99
- external error 8, 11
- FO switch..... 22, 24, 25, 26
- Fusion (Windows emulator for Macintosh) 2
- Fusion program 2
- General Isotope data 13, 20
- Group 9, 63
- grouping 9
- grouping spots by sample name..... 62, 63
- Hidden CPS columns 71
- ionic species box 18
- ionic dpecies 17
- index isotope for common-Pb correction.. 63
- Input Data Files..... 2
- ionic species box 17, 18
- Ionic Species box 70
- Isoplot* 2, 24, 25, 30, 88, 94, 96, 99
- isotope ratios 23
- isotope ratios panel 19
- Isotope ratio/equation charts 34
- Isotope Ratios..... 73, 75, 77
- LA*..... 24
- LA switch 23, 24, 26, 76
- Literal references to ratios or equations.... 21
- Macintosh..... 2
- mass stations 14, 18, 73
- nominal mass 18
- Nonlinear equations 33
- normalizing equation 28
- NU switch 24
- numeric mass 17, 18
- overcounts 9, 64
- page size..... 2
- Parallels (Windows emulator for Macintosh)
 - 2
- Parallels program 2
- PD files..... 18
- PD file 53
- PD files..... 2, 7, 14
- plot insets 8
- Preferences 5, 6, 8, 34, 52, 59, 60, 61, 64, 95
- Printer Setup..... 2
- RAM 2
- Raw data files..... 7
- Reference peak..... 13
- Reference peak, trim mass 13
- regression 30, 31, 78
- rejection, manual..... 35
- Run Editor* 14
- Run Table.... 4, 14, 16, 18, 59, 60, 70, 71, 72
- Sample sheet 6, 9, 97
- SampleData sheet*.... 9, 20, 22, 26, 29, 34, 97
- SBM 6, 60
- SBM (secondary beam monitor)6, 34, 59, 97

SC switch.....	22, 23, 24, 26, 27, 76
scan order.....	18
scan order.....	17
Secular drift, correction for.....	10
smoothing window for secular-drift correction	11
Solver.....	33
sorting by sample age.....	63
Spot names.....	9
spot-scan graphics.....	34, 35
<i>SQUID</i> toolbar	4
<i>SQUID-1</i>	8, 9
ST switch.....	24, 76
Stacey-Kramers.....	64
Standard sheet	8, 24, 98
Standard spots	5, 8, 31, 53, 64, 65
standard, age, U-Pb.....	24
standard, U-Pb.....	5, 13, 21, 26, 28, 59, 74
StandardData sheet.....	20, 22
subsets	61
Subsets	26, 66
Switches	23
Switches, for <i>Task</i> equations.....	75, 77
Task.....	53, 68, 69
204-corr. 206Pb/238U.....	63
207-corr. 206Pb/238U.....	63
²⁰⁷ Pb/ ²⁰⁶ Pb age	5, 59, 63
208-corr. 206Pb/238U.....	63
Age Standard.....	5, 8, 9, 31
AR switch.....	25, 27, 76
HI switch	25
ARrc switch	23, 24, 27
Auto Charts	30, 77, 78
background box.....	18
Calibration constant, U/Pb, correcting for drift of	10
Celeron.....	2
charge.....	72
charge, nuclide, specifying	18
chart insets	14, 53
coherent age group.....	9, 63
column swapping in U-Pb Tasks	28
common Pb	63
Common Pb	95
concordance	59, 63, 64
Task constant	21
Task Editor.....	60, 61
Task Equations.....	60
Task Equation panels	20
Task equations	13, 26, 59
Task Equations.....	20, 21
<i>Task Equations</i>	22
Task, equations	20
time window.....	59, 60, 61
Time Window	31, 65
Toolbar, <i>SQUID</i>	4, 34
trim mass.....	60
Trim mass.....	60, 98
Trim-Mass.....	13
trim-mass charts	60
true mass	18
U-Pb geochronology	4, 9, 13, 20, 23, 26, 27, 28, 30, 31, 38, 65
weighted average plots.....	30
Windows	2
within-scan isotope ratio	53
within-scan ratios/equations.....	34
XML files.....	7
Concordance	95
Concordia.....	8, 88
counts/second columns, specifying.....	18
CPS column	18
CPS columns, hide.....	71
CPS column	18
Crossover	2
decay constants	21
Drift of U/Pb calibration constant, correcting for.....	10
Equation panels	21, 25
Equation switches	23, 75, 77
Equation switches	23
Equations, Task.....	20
Excel	2, 6, 14, 22, 24, 25, 29, 30, 34, 36, 37, 86, 96, 97, 99
external error	8, 11
FO switch.....	22, 24, 25, 26
Fusion (Windows emulator for Macintosh)	2
Fusion program	2
General Isotope data	13, 20

Group	9, 63	SBM	6, 60
grouping	9	SBM (secondary beam monitor).....	6, 34, 59, 97
grouping spots by sample name.....	62, 63	SC switch.....	22, 23, 24, 26, 27, 76
Hidden CPS columns	71	scan order	18
ionic species box	18	scan order	17
ionic dpecies	17	Secular drift, correction for.....	10
index isotope for common-Pb correction..	63	smoothing window for secular-drift	
Input Data Files	2	correction	11
ionic species box	17, 18	Solver	33
Ionic Species box	70	sorting by sample age.....	63
<i>Isoplot</i>	2, 24, 25, 30, 88, 94, 96, 99	Spot names	9
isotope ratios	23	spot-scan graphics	34, 35
isotope ratios panel	19	<i>SQUID</i> toolbar	4
Isotope ratio/equation charts	34	<i>SQUID-1</i>	8, 9
Isotope Ratios.....	73, 75, 77	ST switch.....	24, 76
<i>LA</i>	24	Stacey-Kramers	64
LA switch	23, 24, 26, 76	Standard sheet	8, 24, 98
Literal references to ratios or equations....	21	Standard spots	5, 8, 31, 53, 64, 65
Macintosh.....	2	standard, age, U-Pb	24
mass stations	14, 18, 73	standard, U-Pb.....	5, 13, 21, 26, 28, 59, 74
nominal mass	18	StandardData sheet.....	20, 22
Nonlinear equations	33	subsets	61
normalizing equation	28	Subsets	26, 66
NU switch	24	Switches	23
numeric mass	17, 18	Switches, for <i>Task</i> equations.....	75, 77
overcounts	9, 64	Task.....	53, 68, 69
page size.....	2	Task constant	21
Parallels (Windows emulator for Macintosh)		Task Editor.....	60, 61
.....	2	Task Equations.....	60
Parallels program	2	Task Equation panels	20
PD files.....	18	Task equations	13, 26, 59
PD file	53	Task Equations.....	20, 21
PD files.....	2, 7, 14	<i>Task Equations</i>	22
plot insets	8	Task, equations	20
Preferences 5, 6, 8, 34, 52, 59, 60, 61, 64, 95		time window.....	59, 60, 61
Printer Setup.....	2	Time Window	31, 65
RAM	2	Toolbar, <i>SQUID</i>	4, 34
Raw data files.....	7	trim mass.....	60
Reference peak.....	13	Trim mass.....	60, 98
Reference peak, trim mass	13	Trim-Mass.....	13
regression	30, 31, 78	trim-mass charts	60
rejection, manual.....	35	true mass	18
<i>Run Editor</i>	14	U-Pb geochronology 4, 9, 13, 20, 23, 26, 27,	
Run Table....	4, 14, 16, 18, 59, 60, 70, 71, 72	28, 30, 31, 38, 65	
Sample sheet	6, 9, 97	weighted average plots.....	30
<i>SampleData sheet</i>	9, 20, 22, 26, 29, 34, 97	Windows	2

within-scan isotope ratio	53
within-scan ratios/equations.....	34
204-corr. 206Pb/238U.....	63
207-corr. 206Pb/238U.....	63
²⁰⁷ Pb/ ²⁰⁶ Pb age	5, 59, 63
208-corr. 206Pb/238U.....	63
Age Standard.....	5, 8, 9, 31
AR switch.....	25, 27, 76
HI switch	25
ARrc switch	23, 24, 27
Auto Charts	30, 77, 78
background box.....	18
Calibration constant, U/Pb, correcting for	
drift of	10
Celeron.....	2
charge.....	72
charge, nuclide, specifying	18
chart insets	14, 53
coherent age group.....	9, 63
column swapping in U-Pb Tasks	28
common Pb	63
Common Pb	95
concordance	59, 63, 64
Concordance	95
Concordia.....	8, 88
counts/second columns, specifying.....	18
CPS column	18
CPS columns, hide.....	71
CPS column	18
Crossover	2
decay constants	21
Drift of U/Pb calibration constant, correcting	
for.....	10
Equation panels.....	21, 25
Equation switches	23, 75, 77
Equation switches	23
Equations, Task.....	20
Excel 2, 6, 14, 22, 24, 25, 29, 30, 34, 36, 37,	
86, 96, 97, 99	
external error	8, 11
FO switch.....	22, 24, 25, 26
Fusion (Windows emulator for Macintosh) 2	
Fusion program	2
General Isotope data	13, 20
Group	9, 63

XML files.....	7
grouping	9
grouping spots by sample name.....	62, 63
Hidden CPS columns	71
ionic species box	18
ionic dpecies	17
index isotope for common-Pb correction..	63
Input Data Files.....	2
ionic species box	17, 18
Ionic Species box	70
<i>Isoplot</i>	2, 24, 25, 30, 88, 94, 96, 99
isotope ratios	23
isotope ratios panel	19
Isotope ratio/equation charts	34
Isotope Ratios.....	73, 75, 77
<i>LA</i>	24
LA switch	23, 24, 26, 76
Literal references to ratios or equations....	21
Macintosh.....	2
mass stations	14, 18, 73
nominal mass	18
Nonlinear equations	33
normalizing equation	28
NU switch	24
numeric mass	17, 18
overcounts.....	9, 64
page size.....	2
Parallels (Windows emulator for Macintosh)	
.....	2
Parallels program	2
PD files.....	18
PD file	53
PD files.....	2, 7, 14
plot insets	8
Preferences 5, 6, 8, 34, 52, 59, 60, 61, 64, 95	
Printer Setup.....	2
RAM	2
Raw data files.....	7
Reference peak.....	13
Reference peak, trim mass	13
regression	30, 31, 78
rejection, manual.....	35
<i>Run Editor</i>	14
Run Table....	4, 14, 16, 18, 59, 60, 70, 71, 72

Sample sheet	6, 9, 97	Switches, for <i>Task</i> equations.....	75, 77
<i>SampleData</i> sheet....	9, 20, 22, 26, 29, 34, 97	Task.....	53, 68, 69
SBM	6, 60	Task constant	21
SBM (secondary beam monitor)6, 34, 59, 97		Task Editor.....	60, 61
SC switch.....	22, 23, 24, 26, 27, 76	Task Equations.....	60
scan order	18	Task Equation panels	20
scan order	17	Task equations	13, 26, 59
Secular drift, correction for.....	10	Task Equations.....	20, 21
smoothing window for secular-drift		<i>Task Equations</i>	22
correction	11	Task, equations	20
Solver	33	time window.....	59, 60, 61
sorting by sample age.....	63	Time Window	31, 65
Spot names	9	Toolbar, SQUID.....	4, 34
spot-scan graphics.....	34, 35	trim mass.....	60
<i>SQUID</i> toolbar	4	Trim mass.....	60, 98
<i>SQUID-1</i>	8, 9	Trim-Mass.....	13
ST switch.....	24, 76	trim-mass charts	60
Stacey-Kramers.....	64	true mass	18
Standard sheet	8, 24, 98	U-Pb geochronology 4, 9, 13, 20, 23, 26, 27,	
Standard spots	5, 8, 31, 53, 64, 65	28, 30, 31, 38, 65	
standard, age, U-Pb	24	weighted average plots.....	30
standard, U-Pb.....	5, 13, 21, 26, 28, 59, 74	Windows	2
StandardData sheet.....	20, 22	within-scan isotope ratio	53
subsets	61	within-scan ratios/equations.....	34
Subsets	26, 66	XML files.....	7
Switches	23		
bb			