



Agilent 1200 Series Vacuum Degasser



Service Manual



Agilent Technologies

Notices

© Agilent Technologies, Inc. 2006,
2007-2008

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

G1322-90110

Edition

11/08

Printed in Germany

Agilent Technologies
Hewlett-Packard-Strasse 8
76337 Waldbronn

Research Use Only

Not for use in Diagnostic Procedures.

Warranty

The material contained in this document is provided “as is,” and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

If software is for use in the performance of a U.S. Government prime contract or subcontract, Software is delivered and licensed as “Commercial computer software” as defined in DFAR 252.227-7014 (June 1995), or as a “commercial item” as defined in FAR 2.101(a) or as “Restricted computer software” as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies’ standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will

receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

Contents

1	Introduction	5
	Introduction to the Vacuum Degasser	6
	How the Electronics Work	7
2	Site Requirements and Specifications	11
	Site Requirements	12
	Physical Specifications	15
	Performance Specifications	17
3	Installing the Vacuum Degasser	19
	Unpacking the Vacuum Degasser	20
	Optimizing the Stack Configuration	22
	Installing the Vacuum Degasser	24
	Flow Connections to the Vacuum Degasser	27
	Operational Hints for the Vacuum Degasser	30
	Transporting the Vacuum Degasser	33
4	Using the Vacuum Degasser	35
	When to Use a Vacuum Degasser?	36
	Solvent Information	38
	Prevent Blocking of Solvent Filters	39
	Operation Modes of the Vacuum Degasser	41
5	Troubleshooting and Diagnostics	45
	Overview of the Degasser's Indicators	46
	Status Indicators	47
	Hardware Symptoms	49
6	Maintenance	55
	Introduction to Maintenance and Repair	56

7 Repairs 65

- Exchanging Internal Parts 66
- Overview of the Repair Procedures 69
- Overview of Internal Parts 70

8 Parts and Materials for Maintenance 79

- Cover Parts 80
- Power and Status Light Pipes 81
- Accessory Kit 82

9 Parts for Repair 83

- Vacuum Degasser Parts 84
- Sheet Metal Kit 86

10 Identifying Cables 87

- Cable Overview 88

11 Appendix 97

- General Safety Information 98
- The Waste Electrical and Electronic Equipment (WEEE) Directive (2002/96/EC) 101
- Radio Interference 102
- Sound Emission 103
- Agilent Technologies on Internet 104



1

Introduction

Introduction to the Vacuum Degasser 6

How the Electronics Work 7

Block Diagram 10



Introduction to the Vacuum Degasser

The Agilent 1200 Series vacuum degasser, model G1322A, comprises a 4-channel vacuum container, including 4 tubular plastic membranes, and a vacuum pump. When the vacuum degasser is switched on, the control circuit turns on the vacuum pump which generates a partial vacuum in the vacuum container. The pressure is measured by a pressure sensor. The vacuum degasser maintains the partial vacuum by turning on and off the vacuum pump depending on the signal from the pressure sensor.

The LC pump draws the solvents from their bottles through the special tubular plastic membranes of the vacuum container. As the solvents pass through the vacuum tubes any dissolved gas in the solvents permeates through the membranes into the vacuum container. The solvents will be almost completely degassed when leaving the outlets of the vacuum degasser.




Figure 1 Overview (only one of the four solvent channels is shown)

How the Electronics Work

The vacuum degasser has two different normal operation modes and a continuous mode. In operation mode 1 the vacuum degasser works around a defined set point (115 Torr). Due to environmental conditions it is possible that the vacuum degasser cannot reach the pre-defined set point. Under this condition the operation mode 2 becomes active and the vacuum pump is activated in defined time intervals (vacuum level 115 to 190 Torr). In case of a malfunction of the vacuum degasser (vacuum level above 190 Torr) the instrument is turned into the error mode.




Figure 2 Operation Modes of the Vacuum Degasser

The main function of the vacuum degasser control assembly is to control the vacuum pump and to check the vacuum in the vacuum container.

The power section of the control assembly comprises a switching power supply that generates + 24 V from line voltage. The + 24 V is used to drive the vacuum pump and the solenoid valve. The electronic control circuit uses + 12 V which is generated from the + 24 V.

The pressure sensor is connected to the vacuum chamber and checks for the correct vacuum in the system.

The amplifier and comparator circuit determines the working range of the vacuum that has to be built up. When the vacuum degasser is turned on and the vacuum in the chamber is not within working range (above error limit of 190 Torr), the amplifier and comparator circuit sends a signal to the vacuum pump driver and the timers of the vacuum pump (timer 1) and the solenoid valve (timer 2).

The vacuum pump is turned on immediately while the solenoid valve closing is delayed by 15 seconds. This time delay allows the vacuum pump to start without load before it is connected to the vacuum chamber. The status indicator turns yellow when the pump is activated. The status lamp is turned off as soon as the vacuum is below the error limit.

When the vacuum in the chamber reaches its operation mode 1 (approximately 115 Torr) the amplifier and comparator circuit turns off the solenoid valve and the vacuum pump is turned off by a timer with a delay of 15 seconds.

As soon as the pressure sensor detects that the limit of the operation mode 1 has been exceeded (e.g. when dissolved gas from the solvent diffused into the vacuum chamber) the vacuum pump is started again as described before.

The pressure signal is available at the auxiliary output. It allows to monitor the vacuum system. The upper limit of operation mode 1 is 600 mV. Values below 600 mV on the pressure output indicate sufficient vacuum in the chamber. If the 600 mV are exceeded the vacuum pump will be started to keep the vacuum within its working limit.

The amplifier and comparator circuit also activates the timer 3 when the vacuum in the vacuum chambers is below operation mode 1. The timer 3 is reset when operation mode 1 is reached within a maximum of 8 minutes. If operation mode 1 is not reached and the time (8 minutes) of timer 3 elapsed, the timer mode (operation mode 2) is activated. In this mode the vacuum pump is automatically turned on every 2 minutes for a time frame of 30 seconds. Timers 1 and 2 are activated as described earlier.

The error monitor continuously checks the error limit of the degasser (190 Torr). If the error limit is exceeded (for example, leak in chamber), the error timer is activated and the yellow status indicator lamp is turned on. The vacuum pump is turned on continuously. If the vacuum pump cannot reach either of its two operation modes within 8 minutes (error timer limit) the error driver is activated. The error driver will deactivate the vacuum pump and solenoid valve. The status lamp turns red and the error output on the remote connector will be activated.

The error output at the remote connector provides a *contact closure* (potential free open collector maximum load 35 V DC/50 mA) as long as the error condition is active. The error condition is set (closed) when the status light shows the error condition (red).

The continuous mode overwrites all other operation modes of the degasser. When activated (switch SW1 on the electronic board or via the auxiliary cable) the vacuum pump is forced into continuous mode and is turned on as long as the degasser is switched ON.

Block Diagram




Figure 3 Block Diagram — Vacuum Degasser Control Assembly



2

Site Requirements and Specifications

Site Requirements	12
Power Considerations	12
Power Cords	13
Bench Space	14
Environment	14
Physical Specifications	15
Performance Specifications	17



Site Requirements

A suitable environment is important to ensure optimum performance of the instrument.

Power Considerations

The vacuum degasser power supply has wide-ranging capability (see [Table 1](#) on page 15). It accepts any line voltage in the range mentioned below. Consequently there is no voltage selector in the rear of the vacuum degasser. There are two externally accessible fuses, that protect the power supply.

WARNING

Incorrect line voltage at the instrument

Shock hazard or damage of your instrumentation can result, if the devices are connected to a line voltage higher than specified.

→ Connect your vacuum degasser to the specified line voltage.

CAUTION

Provide access to the power plug

In case of an emergency it must be possible to disconnect the instrument from the power line at any time.

→ Make sure the power connector of the instrument can be easily reached and unplugged.

→ Provide sufficient space behind the power socket of the instrument to unplug the cable.

Power Cords

Different power cords are offered as options with the module. The female end of each of the power cords is identical. It plugs into the power-input socket at the rear of the module. The male end of each of the power cords is different and designed to match the wall socket of a particular country or region.

WARNING

Electric Shock

The absence of ground connection and the use of an unspecified power cord can lead to electric shock or short circuit.

- Never operate your instrumentation from a power outlet that has no ground connection.
 - Never use a power cord other than the Agilent Technologies power cord designed for your region.
-

WARNING

Use of unsupplied cables

The use of cables which haven't been supplied by Agilent Technologies can lead to damage of the electronic components or personal injury.

- Never use cables other than the ones supplied by Agilent Technologies to ensure proper functionality and compliance with safety or EMC regulations.
-

CAUTION

Unaccessible power plug.

In case of emergency it must be possible to disconnect the instrument from the power line at any time.

- Make sure the power connector of the instrument can be easily reached and unplugged.
 - Provide sufficient space behind the power socket of the instrument to unplug the cable.
-

Bench Space

The vacuum degasser dimensions and weight (see [Table 1](#) on page 15) allow to place the vacuum degasser on almost any desk or laboratory bench. It needs an additional 2.5 cm (1.0 inches) of space at either side and approximately 8 cm (3.1 inches) at the rear for the circulation of air and electric connections.

If the bench should carry a complete Agilent 1200 Series system, make sure that the bench is designed to carry the weight of all the modules.

Environment

Your vacuum degasser will work within specifications at ambient temperatures and relative humidity as described in [Table 1](#) on page 15.

CAUTION

Condensation within the vacuum degasser

Condensation will damage the system electronics.

- Do not store, ship or use your vacuum degasser under conditions where temperature fluctuations could cause condensation within the vacuum degasser.
 - If your vacuum degasser was shipped in cold weather, leave it in its box and allow it to warm slowly to room temperature to avoid condensation.
-

Physical Specifications

Table 1 Physical Specifications

Type	Specification	Comments
Weight	7 kg (15.4 lbs)	
Dimensions (width × depth × height)	345 × 435 × 80 mm (13.5 × 17 × 3.1 inches)	
Line Voltage	100 – 120 or 220 – 240 VAC, ± 10 %	Wide-ranging capability
Line Frequency	50 or 60 Hz, ± 5 %	
Power consumption	30 W	Maximum
Ambient Operating Temperature	0 – 55 °C (32 – 131 °F) ¹	
Ambient Non-operating Temperature	-40 – 70 °C (-4 – 158 °F)	
Humidity	< 95 %, at 25 – 40 °C (77 – 104 °F)	Non-condensing
Operating Altitude	Up to 2000 m (6500 ft)	
Non-operating Altitude	Up to 4600 m (14950 ft)	For storing the instrument
Safety Standards: IEC, CSA, UL	Installation Category II, Pollution Degree 2	

¹ This temperature range represents the technical specifications for this instrument. The mentioned temperatures may not be suitable for all applications and all types of solvents.

WARNING

Unspecified Conditions

Operating the instrumentation under conditions other than their intended use might result in a potential safety hazard or might damage the instrumentation.

- Never operate your instrumentation under conditions other than specified by the vendor.

2 Site Requirements and Specifications

Physical Specifications

NOTE

The Agilent 1200 Series vacuum degasser has been tested for evaporation of solvents into the atmosphere by an independent institute with approved methods. The tests were performed with Methanol (BIA Nr. 7810) and Acetonitrile (NIOSH, Nr. 1606). Evaporation of these solvents into the atmosphere when operating the degasser was below the limits of detection.

Performance Specifications

Table 2 Performance Specifications Agilent 1200 Series Vacuum Degasser


Type	Specification
Maximum flow rate	10 ml/min per channel
Number of channels	4
Internal volume per channel	Typically 12 ml per channel
Materials in contact with solvent	PTFE, PEEK
pH range	1 – 14
Analog output (AUX)	For pressure monitoring, range 0 – 3 V

NOTE

The G1322 Vacuum Degasser has been tested for evaporation of solvents into the atmosphere by an independent institute with approved methods. The tests were performed with Methanol (BIA Nr. 7810) and Acetonitrile (NIOSH, Nr. 1606). Evaporation of these solvents into the atmosphere when operating the degasser was below the limits of detection.

2 **Site Requirements and Specifications**

Performance Specifications



3 Installing the Vacuum Degasser

Unpacking the Vacuum Degasser	20
Damaged Packaging	20
Delivery Checklist	20
Accessory Kit Contents	21
Optimizing the Stack Configuration	22
Installing the Vacuum Degasser	24
Flow Connections to the Vacuum Degasser	27
Operational Hints for the Vacuum Degasser	30
Priming the Degasser	30
Priming with a Syringe	31
Priming with the Pump	32
Transporting the Vacuum Degasser	33



Unpacking the Vacuum Degasser

Damaged Packaging

Upon receipt of your vacuum degasser, inspect the shipping containers for any signs of damage. If the containers or cushioning material are damaged, save them until the contents have been checked for completeness and the vacuum degasser has been mechanically and electrically checked. If the shipping container or cushioning material is damaged, notify the carrier and save the shipping material for the carriers inspection.

Delivery Checklist

Ensure all parts and materials have been delivered with the vacuum degasser. The delivery checklist is shown below. To aid in parts identification, please see “[Parts and Materials for Maintenance](#)” on page 79. Please report missing or damaged parts to your local Agilent Technologies Sales and Service Office.

Table 3 Vacuum Degasser Delivery Checklist

Description	Quantity
Vacuum Degasser	1
Power Cable	1
Remote Cable	1
Auxiliary Cable	As ordered
Service Manual	1
Accessory Kit (Table 4 on page 21)	1

Accessory Kit Contents

Table 4 Accessory Kit Contents G1322-68705

Description	Part Number	Quantity
Syringe ¹	5062-8534	1
Syringe Adapter	9301-1337	1
Mounting Tool	0100-1710	1
Waste Tube ²	5062-2463	1
Connecting Tubes labeled A to D	G1322-67300	4

¹ Reorder Number (pack of 10)

² Reorder Number (5 m)

Optimizing the Stack Configuration

If your vacuum degasser is part of a system, you can ensure optimum performance by limiting the configuration of the system stack to the following configuration. This configuration optimizes the system flow path, ensuring minimum delay volume (from point of solvent mixing to head of column) and dead volume (from point of injection to detector outlet).




Figure 4 Recommended Stack Configuration (Front View)




Figure 5 Recommended Stack Configuration (Rear View)

Installing the Vacuum Degasser

Parts required	#	Description
	1	Vacuum degasser
	1	Power cable
	1	Interface cable as ordered, see "Cable Overview" on page 88

Preparations	• Locate bench space • Provide power connections • Unpack the vacuum degasser module
--------------	--

CAUTION

"Defective on arrival" problems

If there are signs of damage, please do not attempt to install the detector. Inspection by Agilent is required to evaluate if the instrument is in good condition or damaged.

- Notify your Agilent sales and service office about the damage.
- An Agilent service representative will inspect the instrument at your site and initiate appropriate actions.

-
- 1 Place the vacuum degasser on the bench.
 - 2 Ensure the power switch on the front of the vacuum degasser is OFF (switch stands out).
 - 3 Connect the power cable to the power connector at the rear of the vacuum degasser.
 - 4 Connect the interface cable to the vacuum degasser. The interface cable (remote cable) is a one way connection to send a not-ready signal from the degasser to the other modules to shut down the whole system after an error condition of the degasser.

NOTE

In an Agilent 1200 Series system, the individual modules are connected through a CAN cable. The Agilent 1200 Series vacuum degasser is an exception. The vacuum degasser can be connected via the APG remote connector to the other modules of the stack. The AUX output allows to monitor the vacuum pressure in the degasser chamber. An Agilent 1200 Series control module can be connected to the CAN bus at any of the modules in the system except for the degasser. The Agilent ChemStation can be connected to the system through one GPIB or LAN cable at any of the modules (except for the degasser), preferably at the detector. For more information about connecting the control module or Agilent ChemStation refer to the respective user manual.




Figure 6 Rear of the Vacuum Degasser




Figure 7 Front of the Vacuum Degasser

3 Installing the Vacuum Degasser

Installing the Vacuum Degasser

- 5 Press in the power switch to turn on the vacuum degasser.

NOTE

The power switch stays pressed in and a green indicator lamp in the power switch is ON when the vacuum degasser is turned ON. When the line power switch stands out and the green light is OFF, the vacuum degasser is turned OFF.

Flow Connections to the Vacuum Degasser

Parts required	#	Description
	1	Solvent cabinet including solvent bottles (filled with solvent) and bottle head assemblies
	1	Vacuum degasser
	1	Solvent outlet tubes
	1	Syringe with adapter

Preparations

- Install the vacuum degasser

WARNING

Toxic and hazardous solvents

The handling of solvents and reagents can hold health risks.

- When opening capillary or tube fittings solvents may leak out.
- Please observe appropriate safety procedures (for example, goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the solvent vendor, especially when toxic or hazardous solvents are used.

- 1 Place solvent cabinet with the bottle(s) on top of the vacuum degasser.
- 2 Remove the front cover by pressing the snap fasteners on both sides.




Figure 8 Removing the Front Cover

3 Installing the Vacuum Degasser

Flow Connections to the Vacuum Degasser

- 3 If the vacuum degasser is not used with an Agilent 1200 Series pump, connect the waste tube from the accessory kit to the waste outlet and place into your waste system.
- 4 Put the bottle head assemblies into solvent bottles containing your mobile phase.
- 5 Connect the solvent tubes from the bottle head assemblies to the inlet connectors A to D (typically the left connection of the channel) of the vacuum degasser. Use the mounting tool shown in [Figure 9](#) on page 28 to fix the tube screw. Fix the tubes in the clips of the vacuum degasser.




Figure 9 Mounting Tool

- 6 Connect the outlet tubes to the output ports (typically right connection of the channel) of the vacuum degasser.
- 7 Prime the degasser before first use (see [“Priming the Degasser”](#) on page 30).

NOTE

Atmospheric gases can diffuse through the tubing and dissolve in the mobile phase solvents. For best chromatographic results, keep the length of tubing between the vacuum degasser and your pump as short as possible.




Figure 10 Flow Connections to the Vacuum Degasser (with Quaternary Pump)

Operational Hints for the Vacuum Degasser

Priming the Degasser

The vacuum degasser can be primed either by drawing solvent through the degasser with a syringe or by pumping with the connected pump.

Priming the vacuum degasser with a syringe is recommended, when:

- vacuum degasser is used for the first time, or vacuum tubes are empty, or
- changing to solvents that are immiscible with the solvent currently in the vacuum tubes.

Priming the vacuum degasser by using the pump at high flow rate (3–5 ml/min) is recommended, when:

- pumping system was turned off for a length of time (for example, overnight) and if volatile solvent mixtures are used, or
- if solvents have been changed.

Priming with a Syringe

Before using a new degasser or new tubings for the first time:

WARNING

Toxic and hazardous solvents

The handling of solvents and reagents can hold health risks.

- When opening capillary or tube fittings solvents may leak out.
- Please observe appropriate safety procedures (for example, goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the solvent vendor, especially when toxic or hazardous solvents are used.

- 1 Prime all tubings with at least 30 ml of iso-propanol no matter whether the channels will be used with organic mobile phase or with water.

If you are changing to a solvent that is immiscible with the solvent currently in the tubing continue as follows:

- 2 Replace the current solvent with iso-propanol, if current solvent is organic or with water, if current solvent is an inorganic buffer or contains salt.
- 3 Disconnect solvent outlet tube of the channel that is supposed to be primed from your pump.
- 4 Connect syringe adapter to solvent outlet tube.
- 5 Push syringe adapter onto syringe.
- 6 Pull syringe plunger to draw at least 30 ml of solvent through degasser and tubing.
- 7 Replace the priming solvent with the new solvent of your choice.
- 8 Pull syringe plunger to draw at least 30 ml of solvent through degasser and tubing.
- 9 Disconnect syringe adapter from solvent tube.
- 10 Connect solvent tube to your pump.
- 11 Repeat step 3 on page 31 to step 10 on page 31 for the other solvent channels.

3 Installing the Vacuum Degasser

Operational Hints for the Vacuum Degasser

NOTE

When priming the vacuum degasser with a syringe the solvent is drawn through the degasser tubes very quickly. The solvent at the degasser outlet will therefore not be fully degassed. Pump for approximately 10 minutes with your selected flow rate before starting any application. This will allow the vacuum degasser to properly degas the solvent in the degasser tubes.

NOTE

The pump should never be used for priming empty tubings (never let the pump run dry). Use the syringe to draw enough solvent for completely filling the tubings to the pump inlet before continueing to prime with the pump.

Priming with the Pump

When the pumping system has been turned off for a certain time (for example, overnight) oxygen will rediffuse into the solvent channels between the vacuum degasser and the pump. Solvents containing volatile ingredients will slightly lose these, if left in the degasser without flow for a prolonged period of time. Therefore priming of the vacuum degasser and the pumping system is required before starting an application.

- 1 Open the purge valve of your pump and set flow rate to 5 ml/min.
- 2 Flush the vacuum degasser and all tubes with at least 30 ml of solvent.
- 3 Set flow to required value of your application and close the purge valve.
- 4 Pump for approximately 10 minutes before starting your application.
- 5 Repeat step 1 on page 32 through step 4 on page 32 with other channels, where needed.

Transporting the Vacuum Degasser

WARNING

Solvents leaking out

Solvents remaining in the solvent channels may leak out during transport. This can possibly cause personal damage.

- Drain any remaining solvents from the degassing channels before transporting the micro vacuum degasser.
-

- 1 Disconnect the solvent tubes from solvent inlets from front panel.
- 2 Disconnect one of the solvent tubes from your pump.
- 3 Connect syringe adapter to solvent tube of this solvent channel.
- 4 Push syringe adapter onto syringe.
- 5 Pull syringe plunger to draw solvent out of vacuum degasser and tubing. Continue to draw solvent into syringe until the solvent channel is completely empty.
- 6 Repeat step 2 on page 33 through step 5 on page 33 for the remaining solvent channels.

3 **Installing the Vacuum Degasser**

Transporting the Vacuum Degasser



4 Using the Vacuum Degasser

When to Use a Vacuum Degasser?	36
Solvent Information	38
Prevent Blocking of Solvent Filters	39
Operation Modes of the Vacuum Degasser	41
Normal Operation Mode 1	42
Normal Operation Mode 2 (Timing Mode)	42
Continuous Mode	43
Error Mode	43



When to Use a Vacuum Degasser?

WARNING

Unspecified Conditions

Operating the instrumentation under conditions other than their intended use might result in a potential safety hazard or might damage the instrumentation.

- Never operate your instrumentation under conditions other than specified by the vendor.
-

Pumps that mix the solvent on the low pressure side like the Agilent 1200 Series quaternary pump do need degassing and must be equipped with a vacuum degasser or alternative degassing systems (for example, helium). Isocratic pumps and high-pressure mixing pumps do not always require degassing. However for the following conditions the vacuum degasser is also recommended for an isocratic or a high pressure mixing pump:

- if your detector is used with maximum sensitivity in the low UV wavelength range,
- if your application requires optimum injection precision, or
- if your application requires highest retention time reproducibility (mandatory at flow rates below 0.5 ml/min),
- if your sample or detection is sensitive to dissolved oxygen in the mobile phase (degradation).

Generally a degasser should be used when negative effects due to dissolved gas in the mobile phase exceed the limits that are acceptable for the user. Negative effects that can be caused by dissolved gas are:


- Unstable flow due to unstable pumping conditions. This may result in a high ripple (unstable pressure at constant flow and with constant mobile phase composition) or high standard deviations of peak retention times and peak areas especially at low flow rates.
- Baseline noise on detectors that are sensitive to changes in the refractive index (e.g. RI detector or UV detector in the low UV range, both at maximum sensitivity).
- Sample degradation.

- Fluorescence Quenching due to dissolved oxygen.
- Baseline drift in electrochemical detectors due to dissolved oxygen especially in reduction mode.

Solvent Information

Always filter solvents through 0.4 µm filters, small particles can permanently block the capillaries. Avoid the use of the following steel-corrosive solvents:

- Solutions of alkali halides and their respective acids (for example, lithium iodide, potassium chloride, and so on).
- High concentrations of inorganic acids like sulfuric acid, especially at higher temperatures (replace, if your chromatography method allows, by phosphoric acid or phosphate buffer which are less corrosive against stainless steel).
- Halogenated solvents or mixtures which form radicals and/or acids, for example:



This reaction, in which stainless steel probably acts as a catalyst, occurs quickly with dried chloroform if the drying process removes the stabilizing alcohol.

- Chromatographic grade ethers, which can contain peroxides (for example, THF, dioxane, di-isopropylether). Such ethers should be filtered through dry aluminium oxide which adsorbs the peroxides.
- Mixtures of carbon tetrachloride with 2-propanol or THF dissolve stainless steel.

Prevent Blocking of Solvent Filters

Contaminated solvents or algae growth in the solvent bottle will reduce the lifetime of the solvent filter and will influence the performance of the pump. This is especially true for aqueous solvents or phosphate buffers (pH 4 to 7). The following suggestions will prolong lifetime of the solvent filter and will maintain the performance of the pump.

- Use sterile, if possible amber solvent bottles to slow down algae growth.
- Filter solvents through filters or membranes that remove algae.
- Exchange solvents every two days or refilter.
- If the application permits add 0.0001–0.001 M sodium azide to the solvent.
- Place a layer of argon on top of your solvent.
- Avoid exposure of the solvent bottles to direct sunlight.

Checking the Solvent Filters

The solvent filters are on the low-pressure side of the pumping system. A blocked filter therefore does not affect the pressure readings of the pump. The pressure readings cannot be used to identify blocked filters. If the solvent cabinet is placed on top of the vacuum degasser the filter condition can be checked in the following way:

Remove the tubing at the inlet port of the vacuum degasser. If the filter is in good condition the solvent will freely drip out of the solvent tube (due to hydrostatic pressure). If the solvent filter is partly blocked no solvent or only very little solvent will drip out of the solvent tube.

Cleaning the Solvent Filters

- Remove the blocked solvent filter from the bottle-head assembly and place it in a beaker with concentrated nitric acid (35%) for one hour.
- Thoroughly flush the filter with bidistilled water (remove all nitric acid).
- Replace the filter.

4 Using the Vacuum Degasser

Prevent Blocking of Solvent Filters

NOTE

Never use the system without solvent filter installed.

Operation Modes of the Vacuum Degasser

The vacuum degasser allows various operation modes. Operation mode 1 and 2 are the normal operation modes of the degasser:

- In operation mode 1 the vacuum degasser works at 115 Torr.
- In operation mode 2 the vacuum degasser works in the range between 115 to 190 Torr.
- The continuous mode can be selected for highest degassing efficiency of the degasser. In this mode the vacuum level is below 115 Torr.
- The error mode is activated in case the degasser cannot achieve a vacuum level of 190 Torr.

(1 Torr = 1.33×10^{-3} bar)

Table 5 Operation Modes of the Vacuum Degasser

	Operation Mode 1	Operation Mode 2	Continuous Mode	Error Mode
Pressure Range [Torr]	115	115 -190	Below 115	Above 190
DC Voltage Readings [mV] (rough values)	Lower Limit: 590 - 600 Upper Limit: 600 - 610	600 - 800	Below 600	Above 800
Operation mode	Pump switched on and off at upper and lower limit during operation	Pump switched on every 2 minutes for 30 seconds	Pump permantly running	8 Minutes of continuous pumping before shutdown
Status Indicator	OFF - below 800 mV YELLOW - above 800 mV	OFF - below 800 mV YELLOW - above 800 mV	OFF - below 800 mV YELLOW - above 800 mV	YELLOW
Failure Actions	Switch to Operation Mode 2	Switch to Error Mode	Switch to Error Mode	Shut down module, switch status indicator to RED

NOTE

The voltage readings given for the various modes are approximate values. All values are set in the factory according to the vacuum level. Depending on the vacuum sensor batch variations this will result in different voltage readings. Regardless of the readings of the various stages no adjustment should be performed on the degasser electronics.

NOTE

See “If the Status Lamp is Yellow and the Vacuum Pump Runs Continuously” on page 52 for setting up test meter connections to the degasser for reading the pressure sensor output.

Normal Operation Mode 1

When the degasser is turned ON, the vacuum pump runs and is connected to ambient through the solenoid valve. The solenoid valve activates about 15 seconds after turning ON the degasser (you can hear it click). The vacuum pump then begins to pull a vacuum on the vacuum chamber. The voltage measurements begin to decrease rapidly. The yellow status LED turns off when the vacuum level reaches 190 Torr (DC voltage readings around 800mV). The normal operation mode vacuum level (115 Torr) is achieved when the DC voltage measures approximately 590 to 600 mV. After achieving the normal operation mode vacuum level, the solenoid valve turns off. The vacuum pump continues to run for a few seconds, then it turns OFF. When the DC voltage measurement rises back to approximately 600 to 610 mV, the turn on process begins again. If the vacuum level of normal operation mode one cannot be achieved within 8 minutes the instrument turns into normal operation mode 2.

Normal Operation Mode 2 (Timing Mode)

Under certain operational conditions (large amount of dissolved gas in mobile phase, high flow rates) the 115 Torr trigger level for operation mode 1 cannot be reached. The vacuum degasser automatically turns into operation mode 2. Normal operation mode 2 is a fixed timing mode. Every 2 minutes the degasser is turned ON for 30 seconds. This assures a pressure level in the range from 115 to 190 Torr. The pressure signal measured with the auxiliary cable is in the range between 600 to 800 mV. If the vacuum level of normal operation mode two cannot be achieved within 8 minutes the instrument turns into error mode.

Continuous Mode

The continuous mode is activated either by switching SW1 on the main board to 1 (removal of cover is required, see “[Removing and Refitting the Top Cover](#)” on page 59, for identifying SW1, see “[Overview of Internal Parts](#)” on page 70), or by connecting pin 1 (white cable) and pin 3 (green cable) of the auxiliary cable to each other. When turning ON the vacuum degasser the vacuum pump will run continuously. This will establish a vacuum level that is below the trigger level (600 mV / 115 Torr) of the normal operation mode 1. If the vacuum level of normal operation mode 2 cannot be achieved within 8 minutes the instrument turns into error mode.

NOTE

When set to continuous mode the life time of the vacuum pump will be significantly reduced.

When to use Continuous Mode


In continuous mode the vacuum pump runs continuously thus generating the highest degassing efficiency of the degasser. This mode is only recommended for very sensitive applications (e.g. RI detection).

Error Mode

The error level for the vacuum degasser is 190 Torr (approximately 800 mV). This level cannot be achieved when there is a failure in the degasser (for example, leaks, etc.). When the error level is exceeded the yellow status indicator lamp is turned on and the vacuum pump runs continuously. If the degasser cannot reach one of the normal operation modes within 8 minutes the status indicator turns red and the vacuum pump is turned OFF.

4 Using the Vacuum Degasser

Operation Modes of the Vacuum Degasser



5 Troubleshooting and Diagnostics

Overview of the Degasser's Indicators 46

Status Indicators 47

Power Supply Lamp 47

Instrument Status Indicator 48

Hardware Symptoms 49

All Lamps are Off and the Vacuum Degasser Appears Dead 49

If the Status Indicator is Red 50

If the Status Indicator is Yellow and the Vacuum Pump is not Running 51

If the Status Lamp is Yellow and the Vacuum Pump Runs Continuously 52



Overview of the Degasser's Indicators

Status Indicators

The vacuum degasser is provided with two status indicators which indicate the operational state (ready, busy, and error states) of the vacuum degasser. The status indicators provide a quick visual check of the operation of the vacuum degasser (see [“Overview of the Degasser's Indicators”](#) on page 46).

Hardware Symptoms

A red status lamp at the vacuum degasser indicates a problem with the vacuum system or with the electronic control. The vacuum degasser generates an error output on the remote lines. The following pages describe hardware symptoms which help you to isolate the cause of a hardware failure (see [“Hardware Symptoms”](#) on page 49).

Status Indicators

Two status lamps are located on the front of the vacuum degasser. The left lamp indicates the power supply status, the right lamp indicates the vacuum degasser status.




Figure 11 Location of Status Lamps

Power Supply Lamp

The power supply lamp is integrated into the main power switch. When the lamp is illuminated (green), the power is ON.

Instrument Status Indicator

The instrument status indicator indicates one of three possible instrument conditions.

- When the lamp is OFF, the vacuum degasser is in *ready* condition (only if the power supply lamp is ON, otherwise, the instrument is switched OFF, or there is a defect in the power supply). A ready condition exists when there is sufficient vacuum in the vacuum degasser.
- A *busy* condition is indicated, when the lamp is yellow. A busy condition exists when the vacuum pump of the vacuum degasser is working to generate or maintain vacuum in the chambers. This is the case when the vacuum degasser is turned ON for the first time or the pressure rises above its limit during normal operation.
- An *error* condition is indicated, when the lamp is red. An error condition exists when the vacuum degasser detects an internal defect that does not allow the vacuum to be built up in a certain time frame (approximately 8 minutes).

CAUTION

Internal Leak or Electronic Failure

If the error LED is on, there is either an internal leak in the vacuum system or an electronic failure. In case of an internal leak it is possible that solvent may enter the vacuum chamber and solvent may leak into the waste drain.

- To prevent any damage of the vacuum degasser, switch off the vacuum degasser and remove the solvent bottles from the solvent cabinet to stop any gravity-caused flow of solvent into the vacuum chamber.

Hardware Symptoms

In case of a problem with the vacuum system or the electronic control the vacuum degasser status lamp will be red. The vacuum degasser will generate an error output on the remote lines. This will shut down other system modules when connected via remote cable, see [“Installing the Vacuum Degasser”](#) on page 24. The vacuum degasser itself will not be able to generate any error messages in the system logbook of the Agilent 1200 Series.

All Lamps are Off and the Vacuum Degasser Appears Dead

If all other modules in the system are ON (power switch lamp is green) and are recognized by the connected user interface (module parameters can be set, module-specific screens appear, and so on), then do the following to determine the problem with the vacuum degasser:

- Ensure the power cable is connected to the degasser, and the power cable is connected to line power.
- Ensure the power switch on the front of the module is ON.
- Ensure the power fuses are OK.

The fuse holders are located on the rear panel of the vacuum degasser and are part of the power socket. Check the fuses, and change if necessary, see [“Checking and Changing the Power Fuses”](#) on page 63.

If changing the fuses did not solve the problem, change the control assembly (see [“Exchanging the Degasser Control Assembly and the Sensor Assembly”](#) on page 76).

If the Status Indicator is Red

Sufficient vacuum is normally built up during the initial start-up and is maintained by turning on the vacuum pump whenever triggered by the vacuum sensor.

If the vacuum cannot be reached, or is lost faster than expected during operation, the vacuum degasser will be forced into an error state. Approximately 8 to 16 minutes after turning on the vacuum pump without reaching its normal operation modes the error LED is turned on and the vacuum pump and the solenoid will be disabled.

The error condition can be reset by turning the vacuum degasser off and on again. If the error condition persists the error will occur again after another 8 to 16 minutes.

The following described failure conditions will appear during the start-up procedure before the red error LED is turned on. Every power cycling of the degasser will give you 8 to 16 minutes of time to diagnose the failure conditions described on the following pages.

If the Status Indicator is Yellow and the Vacuum Pump is not Running

- Remove the top cover (see [“Removing and Refitting the Top Cover”](#) on page 59).
- Check the operation of the vacuum pump. Disconnect the vacuum tube from the sensor assembly to the vacuum chamber at the vacuum chamber (see [“Overview of Internal Parts”](#) on page 70) to adjust the pressure inside the vacuum chamber to ambient conditions. Switch the vacuum degasser OFF and ON again. The vacuum pump should start immediately.

NOTE

As long as there is sufficient vacuum in the chamber the pump will not be turned ON until triggered by the vacuum sensor.

-
- Use a test meter to check for + 24 V DC on the connector CN2 of the control assembly between pin 1 and 2 (see [“Overview of Internal Parts”](#) on page 70). This voltage is needed to operate the pump. If this voltage is low, exchange the control assembly (see [“Exchanging the Degasser Control Assembly and the Sensor Assembly”](#) on page 76).
 - Check the resistance of the motor windings. The windings should have a resistance in the kOhm range (use a test meter to check resistance between red and blue, blue and black, black and red wires from connector CN2 to the pump). If the windings are broken or shortened replace the vacuum pump (see [“Exchanging the Vacuum Pump”](#) on page 73).

If the Status Lamp is Yellow and the Vacuum Pump Runs Continuously

Before starting the following troubleshooting procedures check that all tubes on the vacuum container, the pressure sensor, the vacuum pump and the solenoid valve are tight and connected properly. Then follow the troubleshooting instructions step by step to eliminate the individual parts one after the other as the sources of the problem.

Preparing to Measure the Pressure Sensor Readings

First of all remove the top cover (see [“Removing and Refitting the Top Cover”](#) on page 59). Use the [“Overview of Internal Parts”](#) on page 70 as a guide to indentify the individual parts.

For test purposes and for monitoring the readings of the pressure sensor connect the auxiliary cable to the AUX output at the back of the vacuum degasser. Connect a test meter to wire 1 (white) and 2 (brown) of the auxiliary cable. When the vacuum chamber is at ambient pressure, a DC voltage measurement across these wires will give between 1.5 and 2.7 Volts (for pin connections see, [“Auxiliary Cable”](#) on page 95).

NOTE

The voltage readings given for the various modes are rough values. All values are set in the factory according to the vacuum level. Depending on the vacuum sensor batch variations this will result in different voltage readings. Regardless of the readings of the various stages no adjustment should be performed on the degasser electronics.

Eliminate the vacuum pump as the source of the problem.

- 1 Remove the tube going from the sensor assembly to the vacuum chamber at the vacuum chamber.
- 2 For convenience, temporarily remove the vacuum pump and place it on top of the vacuum chamber. Disconnect the tube which is connected to the vacuum pump inlet fitting.
- 3 Switch degasser into continuous mode by turning SW1 to 1.
- 4 To duplicate the normal operation of the degasser, turn ON the degasser, wait approximately 15 seconds, then connect the sensor assembly directly to the inlet of the vacuum pump.

Results

A pressure reading below 500 mV indicates a good functioning vacuum pump. A pressure reading above 600 mV indicates that one of the pump chambers is contaminated or defective. In this case, disconnect the sensor tube and let the pump run without load for 5 minutes. Reconnect the sensor tube and repeat the measurement. If the correct vacuum is obtained the contamination has been flushed away and the vacuum pump is functioning correctly. If the pressure reading is still above 600 mV, the vacuum pump should be replaced (see [“Exchanging the Vacuum Pump”](#) on page 73).

- 5 Switch the degasser back to normal operation modes by turning SW1 back to 0.
- 6 Properly reconnect the tubings.

Eliminate the solenoid valve as the source of the problem.

- 1 Disconnect the outlet tubing from the solenoid valve (NC).
- 2 Remove the tube going from the sensor assembly to the vacuum chamber at the vacuum chamber.
- 3 For convenience remove the solenoid valve and place it on top of the vacuum chamber.
- 4 Connect the tubing from the sensor assembly directly to the outlet of the solenoid valve (NC).
- 5 Switch degasser into continuous mode by turning SW1 to 1.

Results

If the obtained pressure reading is similar to the result from the previous measurement (see above) when testing the pump (below 500 mV), the solenoid valve is functioning correctly. If the obtained pressure reading is different (above 500 mV) from the results obtained when testing the pump (see above), the solenoid valve is defective and needs replacing (see [“Exchanging the Solenoid Valve”](#) on page 75).

- 6 Switch the degasser back to normal operation modes by turning SW1 back to 0.
- 7 Properly reconnect the tubings and reinstall the solenoid valve.

Eliminate the vacuum chamber as the source of the problem.

- 1 Disconnect the outlet tube from the solenoid valve (NC) at the solenoid valve.
- 2 Disconnect the tube which is connected to the inlet of the vacuum pump.

- 3** Switch degasser into continuous mode by turning SW1 to 1.
- 4** To duplicate the normal operation of the degasser, turn on the degasser, wait approximately 15 seconds, then connect the outlet tube from the solenoid valve directly to the pump inlet.


Results

A pressure reading below 500 mV indicates a good working vacuum chamber. A pressure reading above 600 mV indicates a leak in the vacuum chamber. If necessary exchange the vacuum chamber (see [“Exchanging the Vacuum Chamber”](#) on page 71).

- 5** Switch the degasser back to normal operation modes by turning SW1 back to 0.
- 6** Properly reconnect the tubings and reinstall the pump.

Exchange the electronic control assembly

If the above described troubleshooting steps do not lead to an identification of the problem, replace the electronic control module with the pressure sensor assembly (see [“Exchanging the Degasser Control Assembly and the Sensor Assembly”](#) on page 76).



6 Maintenance

Introduction to Maintenance and Repair	56
Simple Repairs - Maintenance	56
Warnings and Cautions	57
Using the ESD Strap	58
Cleaning the Instrument	58
Removing and Refitting the Top Cover	59
Assembling the Main Cover	62
Checking and Changing the Power Fuses	63
Exchanging the Status Light Pipe	64



Introduction to Maintenance and Repair

Simple Repairs - Maintenance

The vacuum degasser is designed for easy repair. The most frequent repairs such as exchanging power fuses and status light pipes can be performed by the user, but require opening the main cover of the vacuum degasser. These repairs are described in this chapter.

Warnings and Cautions

WARNING

Open main cover

The following procedures require opening the main cover of the vacuum degasser.

- To prevent personal injury, remove the power cable from the vacuum degasser before opening the cover.
 - Do not connect the power cable to the vacuum degasser while the covers are removed.
-

WARNING

Toxic and hazardous solvents

The handling of solvents and reagents can hold health risks.

- When opening capillary or tube fittings solvents may leak out.
 - Please observe appropriate safety procedures (for example, goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the solvent vendor, especially when toxic or hazardous solvents are used.
-

CAUTION

Electronic boards and components are sensitive to electrostatic discharge (ESD). ESD can damage electronic boards and components.

- In order to prevent damage always use ESD protection when handling electronic boards and components.
-

CAUTION

The sheet metal plates of the degasser are very thin.

Although they have been deburred, they are still quite sharp. You may cut your hands or fingers.

- Never slide your fingers along the edges of the enclosure.
-

Using the ESD Strap

Electronic boards are sensitive to electronic discharge (ESD). In order to prevent damage, always use an ESD strap supplied in the standard accessory kit when handling electronic boards and components.

- 1 Unwrap the first two folds of the band and wrap the exposed adhesive side firmly around your wrist.
- 2 Unroll the rest of the band and peel the liner from the copper foil at the opposite end.
- 3 Attach the copper foil to a convenient and exposed electrical ground.




Figure 12 Using the ESD Strap

Cleaning the Instrument

WARNING

Liquid dripping into the electronic compartment of your module.

Liquid in the module electronics can cause shock hazard and damage the module.

- Do not use an excessively damp cloth during cleaning.
- Drain all solvent lines before opening any fittings.

The vacuum degasser case should be kept clean. Cleaning should be done with a soft cloth slightly dampened with water or a solution of water and a mild detergent. Do not use an excessively damp cloth that liquid can drip into the vacuum degasser.

Removing and Refitting the Top Cover

Tools required

- Screwdriver Pozidriv #1

Preparations

- Switch OFF the vacuum degasser at the main power switch
- Disconnect the power cable and remote cable
- Disconnect all solvent tubes from the ports of the vacuum degasser
- Remove solvent cabinet from the vacuum degasser
- Remove vacuum degasser from the stack.


WARNING

Toxic and hazardous solvents


The handling of solvents and reagents can hold health risks.

- When opening capillary or tube fittings solvents may leak out.
- Please observe appropriate safety procedures (for example, goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the solvent vendor, especially when toxic or hazardous solvents are used.

1 Remove the front panel.




2 Unclip the clips on the top cover .




6 Maintenance

Introduction to Maintenance and Repair


3 Lift the cover away.




4 Unscrew the screws on the top plate and remove the plate by lifting its back and then sliding the plate to the front.




5 Place the metal cover onto the housing. Lower the front of the plate first, then slide panel to the back. Assure the four metal tabs of the panel slide underneath the Z-plane. Fix the two holding screws.



6 Replace the top cover. Ensure the clips are seated correctly.



7 Reinstall the front cover.



8 Reinstall the vacuum degasser in your system stack and connect the cables and capillaries and turn on the vacuum degasser.

Assembling the Main Cover

When • If cover is broken.

Parts required	#	Part number	Description
	1	5065-9989	Cover kit (includes base, top, left and right)

CAUTION

Wrong assembly


In case you insert the left or right side in the opposite position, you may not be able to remove the side from the top part.

→ Take care not to mix up left and right side.


NOTE

The cover kit contains all parts, but it is not assembled.

1 Place the top part on the bench and insert the left and right side into the top part.



2 Replace the cover.



Next Steps:

- 3** Replace the vacuum degasser in the stack and reconnect the cables and capillaries.
- 4** Turn ON the vacuum degasser.

Checking and Changing the Power Fuses

When • If the vacuum degasser appears dead

Tools required • Test meter (if available)

Parts required

#	Part number	Description
1	2110-0458	Fuse 500 mA

- 1 Switch OFF the power switch at the front of the instrument.
- 2 Remove the power cable from the power connector at the rear of the instrument.
- 3 Press down the clip of the fuse holder and pull out of the power socket.





Figure 13 Fuse Holder and Clip

- 4 Remove the fuses from the fuse holders.
- 5 Ensure the fuse wires inside the fuses are not broken. If a test meter is available, check the resistance of each fuse. A good fuse shows a low resistance (approximately 0 Ohm).
- 6 If a fuse is defective (wire broken or high resistance), insert a new fuse.
- 7 Reinsert the fuse holders and the power cable.
- 8 Switch ON the power switch.


Exchanging the Status Light Pipe

- When**
- If part is broken
- Tools required**
- Screwdriver Pozidriv#1
- Parts required**
- | # | Part number | Description |
|---|-------------|-------------------|
| 1 | 5041-8384 | Status light pipe |
- Preparations**
- Remove the front cover and top cover, see [“Removing and Refitting the Top Cover”](#) on page 59.

1 The status light pipe is clipped into the top cover.



2 Replace the top cover, see [“Removing and Refitting the Top Cover”](#) on page 59.



7 Repairs

Exchanging Internal Parts	66
Warnings and Cautions	67
Overview of the Repair Procedures	69
Overview of Internal Parts	70
Exchanging the Vacuum Chamber	71
Exchanging the Vacuum Pump	73
Exchanging the Solenoid Valve	75
Exchanging the Degasser Control Assembly and the Sensor Assembly	76



Exchanging Internal Parts

Most of the repairs require exchange of defective internal parts. Exchange of these parts requires removing the vacuum degasser from the stack, removing the covers, and disassembling the vacuum degasser. The security lever at the power input socket prevents that the degasser cover is taken off when line power is still connected.

The procedures in this section describe how to exchange internal parts. You must remove the vacuum degasser from the stack in order to open the main cover.

Table 6 Exchanging Internal Parts

Procedure	Typical Frequency	Notes
“Exchanging the Vacuum Chamber” on page 71	If defective	Yellow status lamp, pressure reading above 600 mV
“Exchanging the Vacuum Pump” on page 73	If electrically defective	Pressure reading above 600 mV or vacuum pump not running
“Exchanging the Solenoid Valve” on page 75	If defective or leaking	Pressure reading above 500 mV
“Exchanging the Degasser Control Assembly and the Sensor Assembly” on page 76	If defective	Low voltage at control assembly

Warnings and Cautions

WARNING

Open main cover

The following procedures require opening the main cover of the instrument.

- Always ensure the instrument is disconnected from the line power when the main cover is removed.
 - The security lever at the power input socket prevents that the instrument cover is taken off when line power is still connected.
-

WARNING

Instrument is partially energized when switched off

The power supply still uses some power, even if the switch on the front panel is turned off.

- To disconnect the vacuum degasser from line, unplug the power cord.
-

WARNING

Toxic and hazardous solvents

The handling of solvents and reagents can hold health risks.

- When opening capillary or tube fittings solvents may leak out.
 - Please observe appropriate safety procedures (for example, goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the solvent vendor, especially when toxic or hazardous solvents are used.
-

CAUTION

Electronic boards and components are sensitive to electrostatic discharge (ESD). ESD can damage electronic boards and components.

- In order to prevent damage always use ESD protection when handling electronic boards and components.
-

NOTE

The electronics of the degasser will not allow operation of the degasser when the top cover and the top foam are removed. A safety light switch on the main board will inhibit the operation of the fan immediately. Voltages for the other electronic components will be turned off after 30 seconds. The status lamp will light up red and an error will be logged into the logbook of the user interface. Always operate the degasser with the top covers in place.

Overview of the Repair Procedures




Figure 14 Overview of Repair Procedures

11	Vacuum Chamber, see “Exchanging the Vacuum Chamber” on page 71
22	Control Assembly and Sensor Assembly, see “Exchanging the Degasser Control Assembly and the Sensor Assembly” on page 76
33	Solenoid Valve, see “Exchanging the Solenoid Valve” on page 75
44	Vacuum Pump, see “Exchanging the Vacuum Pump” on page 73

Overview of Internal Parts

Figure 15 on page 70 shows a top view of the vacuum degasser. It shows the connections of the vacuum tubes as well as the electrical connections at the control assembly.




Figure 15 Top View of Vacuum Degasser

Exchanging the Vacuum Chamber

When • If internal membrane defective or vacuum chamber damaged

Tools required • Screwdriver Pozidriv #1

Parts required	#	Part number	Description
	1	G1322-60001	Vacuum chamber

Preparations • Disconnect the solvent inlet tubes at the degasser inlet ports
• Remove the solvent cabinet from the vacuum degasser

WARNING

Toxic and hazardous solvents

The handling of solvents and reagents can hold health risks.

- When opening capillary or tube fittings solvents may leak out.
- Please observe appropriate safety procedures (for example, goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the solvent vendor, especially when toxic or hazardous solvents are used.

-
- 1 Disconnect a solvent tube of one of the channels from your pump.
 - 2 Connect syringe adapter to the disconnected solvent tube.
 - 3 Push syringe adapter onto syringe.
 - 4 Pull syringe plunger to draw solvent out of degasser and tubing. Continue to draw solvent into syringe until there is no solvent in the tubing.
 - 5 Repeat step 1 on page 71 to step 4 on page 71 for other solvent channels.
 - 6 Disconnect outlet tubes from the vacuum degasser.
 - 7 Remove top cover, see [“Removing and Refitting the Top Cover”](#) on page 59.
 - 8 Loosen the clip of the sensor board and remove it from its slot in the vacuum chamber (for identifying internal parts see [“Overview of Internal Parts”](#) on page 70).

WARNING

Solvents may leak out from the degassing channels.

Risk of personal injury from toxic or hazardous solvents

→ Drain any remaining solvents from the degassing channels before removing the vacuum chamber from the vacuum degasser.

- 9** Remove the two holding screws of the vacuum chamber.
- 10** Slide the vacuum chamber out of the front plate.
- 11** Disconnect the vacuum tubes from the vacuum chamber.
- 12** Reconnect the two vacuum tubes to the new vacuum chamber.
- 13** Slide the new vacuum chamber into the front plate and fix it with the two holding screws.
- 14** Slide the sensor board back into the slots of the vacuum chamber and fix the clip.
- 15** Reconnect all inlet and outlet tubes to the vacuum chamber, see [Figure 16](#) on page 72.
- 16** Refit the cover, see [“Removing and Refitting the Top Cover”](#) on page 59.




Figure 16 Exchanging the Vacuum Chamber

Exchanging the Vacuum Pump

When	• If not generating sufficient vacuum or electrically defective								
Tools required	• Screwdriver Pozidriv #1								
Parts required	<table><tr><th>#</th><th>Part number</th><th>Description</th></tr><tr><td>1</td><td>G1322-60000</td><td>Vacuum pump</td></tr></table>	#	Part number	Description	1	G1322-60000	Vacuum pump	1 Remove Top Cover, see “Removing and Refitting the Top Cover” on page 59. 2 Disconnect vacuum pump cable from connector CN2 at the control assembly (for identifying internal parts see “Overview of Internal Parts” on page 70). 3 Using a Pozidriv screwdriver #1 loosen and remove the two holding screws of the vacuum pump. 4 Remove the fixing plates. 5 Lift the vacuum pump, disconnect the inlet tube (connected to the switching valve) and take the pump out of the unit. 6 Connect the tubing from the switching valve to the exchange vacuum pump inlet. 7 Place the vacuum pump in its position. 8 Place the fixing plates onto the rubber feet of the pump. 9 Insert the holding screws and fix them. 10 Connect the vacuum pump cable to the control assembly (CN2).	
#	Part number	Description							
1	G1322-60000	Vacuum pump							

11 Refit the cover, see “[Removing and Refitting the Top Cover](#)” on page 59.




Figure 17 Exchanging the Vacuum Pump

Exchanging the Solenoid Valve

When • Valve not switching correctly or leaking

Tools required • Screwdriver Pozidriv #1

Parts required

#	Part number	Description
1	G1322-60003	Solenoid valve

- 1 Remove Top Cover, see [“Removing and Refitting the Top Cover”](#) on page 59.
- 2 Disconnect the solenoid connector from the control board (CN3) (for identifying internal parts see [“Overview of Internal Parts”](#) on page 70).
- 3 Loosen and remove the holding screws of the solenoid valve.
- 4 Loosen and remove the two screws of the holding plate at the solenoid.
- 5 Disconnect inlet and outlet tubes from the solenoid valve.
- 6 Fix the holding plate to the new solenoid valve.
- 7 Connect the two vacuum tubes to the solenoid valve. The inlet of the solenoid is labeled IN and is connected to the vacuum pump.
- 8 Fix the solenoid valve with the two screws to its holder.
- 9 Connect the solenoid cable to the control board (CN3).
- 10 Refit the cover, see [“Removing and Refitting the Top Cover”](#) on page 59.




Figure 18 Exchanging the Solenoid Valve

Exchanging the Degasser Control Assembly and the Sensor Assembly

When • No power available, vacuum pump or switching valve not activated

Tools required • Screwdriver Pozidriv #1

Parts required	#	Part number	Description
	1	G1322-66500	Degasser control assembly (includes the sensor assembly)

WARNING

The degasser control assembly does not include any serviceable parts. It is not intended to be opened by an unauthorized persons.

Risk of personal injury!

→ Do not open the degasser control assembly.

NOTE

Sensor assembly and degasser control assembly are a matched pair and must be exchanged together.

- 1 Remove Top Cover, see [“Removing and Refitting the Top Cover”](#) on page 59.
- 2 Disconnect all connectors at the control assembly board (for identifying internal parts see [“Overview of Internal Parts”](#) on page 70).
- 3 Loosen and remove the two holding clips at the control board of the control assembly.
- 4 Unclip the power switch light pipe at the degasser control assembly power switch.
- 5 Remove the four holding screws at the side panel – two at the left side, two on the right side of the vacuum degasser.
- 6 Lift the degasser control assembly out of the unit.
- 7 Remove the coupler from the control assembly and place it onto the power switch in the exchange control assembly.
- 8 Loosen the clip of the sensor board and remove it from its slot in the vacuum chamber.

- 9 Slide the new sensor board into the slots of the vacuum chamber and the fix the clip.
- 10 Place a new degasser control assembly into the unit and fix with the four holding screws to the side panels.
- 11 Clip the power switch light pipe to the degasser control assembly power switch.
- 12 Replace the two holding clips and fix the holding screws.
- 13 Connect all the cables to the control board (sensor assembly cable to CN6, solenoid valve cable to CN3, pump cable to CN2).
- 14 Refit the cover, see “[Removing and Refitting the Top Cover](#)” on page 59.





Figure 19 Exchanging the Degasser Control Assembly

7 Repairs

Overview of Internal Parts



8

Parts and Materials for Maintenance

Cover Parts 80

Power and Status Light Pipes 81

Accessory Kit 82



Cover Parts

Table 7 Cover Parts

Item	Description	Part Number
1	Set of all plastic covers, top-side-base	5065-9989
2	Front cover	5065-9990
3	Logo plate, Agilent 1200	5042-8901
4	Tube clip	5041-8387




Figure 20 Cover Parts

Power and Status Light Pipes

Table 8 Power and Status Light Pipes

Item	Description	Part Number
1	Power switch coupler	5041-8383
2	Light pipe – power switch	5041-8382
3	Power switch button	5041-8381
4	Light pipe – status lamp	5041-8384




Figure 21 Power and Status Light Pipes

Accessory Kit

Table 9 Accessory Kit G1322-68705

Item	Description	Part Number
1	Syringe ¹	5062-8534
2	Syringe adapter	9301-1337
3	Mounting Tool	0100-1710
4	Solvent tubing kit (4 tubes degasser to pump)	G1322-67300
5	Waste tube**	5062-2463

¹ Reorder number (pack of 10)**Reorder number (5m)



9 Parts for Repair

Vacuum Degasser Parts 84

Sheet Metal Kit 86



Vacuum Degasser Parts




Figure 22 Vacuum Degasser Parts

Table 10 Repair Parts: Vacuum Degasser

Item	Description	Part Number
1	Vacuum chamber	G1322-60001
2	Sensor assembly (included in the control assembly)	no PN
3	Vacuum degasser control assembly	G1322-66500
4	Fuse 500 mA	2110-0458
5	Board clip	G1322-43100
6	Solenoid valve	G1322-60003
7	Vacuum tubes	G1322-60002
8	Fixing plate	no PN
9	Vacuum pump	G1322-60000
10	Leak tray	G1322-44100
11	Leak pan, degasser	5042-8567

Sheet Metal Kit

Table 11 Sheet Metal Kit

Item	Description	Part Number
1	Sheet metal kit, includes base and top plate	G1322-68701





Figure 23 Sheet Metal Kit



10 Identifying Cables

Cable Overview	88
Remote Cables	90
Auxiliary Cable	95



Cable Overview


NOTE

Never use cables other than the ones supplied by Agilent Technologies to ensure proper functionality and compliance with safety or EMC regulations.

Type	Description	Part Number
Analog cables	3390/2/3 integrators	01040-60101
	3394/6 integrators	35900-60750
	Agilent 35900A A/D converter	35900-60750
	General purpose (spade lugs)	01046-60105
Remote cables	3390 integrator	01046-60203
	3392/3 integrators	01046-60206
	3394 integrator	01046-60210
	3396A (Series I) integrator	03394-60600
	3396 Series II / 3395A integrator, see details in section “Remote Cables” on page 90	
	3396 Series III / 3395B integrator	03396-61010
	HP 1050 modules / HP 1046A FLD	5061-3378
	HP 1046A FLD	5061-3378
	Agilent 35900A A/D converter	5061-3378
	HP 1040 diode-array detector	01046-60202
	HP 1090 liquid chromatographs	01046-60202
	Signal distribution module	01046-60202
BCD cables	3396 integrator	03396-60560
	General purpose (spade Lugs)	G1351-81600
Auxiliary	Agilent 1100 Series vacuum degasser	G1322-61600


Type	Description	Part Number
<i>CAN cables</i>	Agilent 1100/1200 module to module, 0.5m lg	5181-1516
	Agilent 1100/1200 module to module, 1m lg	5181-1519
<i>External contacts</i>	Agilent 1100/1200 Series interface board to general purpose	G1103-61611
<i>GPIB cable</i>	Agilent 1100/1200 module to ChemStation, 1 m	10833A
	Agilent 1100/1200 module to ChemStation, 2 m	10833B
<i>RS-232 cable</i>	Agilent 1100/1200 module to a computer This kit contains a 9-pin female to 9-pin female Null Modem (printer) cable and one adapter.	34398A
<i>LAN cable</i>	Twisted pair cross over LAN cable, (shielded 3m long) (for point to point connection)	5023-0203
	Twisted pair cross over LAN cable, (shielded 7m long) (for point to point connection)	5023-0202

Remote Cables




One end of these cables provides a Agilent Technologies APG (Analytical Products Group) remote connector to be connected to Agilent 1100/1200 Series modules. The other end depends on the instrument to be connected to.


Agilent 1100/1200 to 3390 Integrators

Connector 01046-60203	Pin 3390	Pin Agilent 1100/1200	Signal Name	Active (TTL)
	2	1 - White	Digital ground	
	NC	2 - Brown	Prepare run	Low
	7	3 - Gray	Start	Low
	NC	4 - Blue	Shut down	Low
	NC	5 - Pink	Not connected	
	NC	6 - Yellow	Power on	High
	NC	7 - Red	Ready	High
	NC	8 - Green	Stop	Low
	NC	9 - Black	Start request	Low

Agilent 1100/1200 to 3392/3 Integrators

Connector01046-60206	Pin 3392/3	Pin Agilent 1100/1200	Signal Name	Active (TTL)
	3	1 - White	Digital ground	
	NC	2 - Brown	Prepare run	Low
	11	3 - Gray	Start	Low
	NC	4 - Blue	Shut down	Low
	NC	5 - Pink	Not connected	
	NC	6 - Yellow	Power on	High
	9	7 - Red	Ready	High
	1	8 - Green	Stop	Low
	NC	9 - Black	Start request	Low


Agilent 1100/1200 to 3394 Integrators

Connector01046-60210	Pin 3394	Pin Agilent 1100/1200	Signal Name	Active (TTL)
	9	1 - White	Digital ground	
	NC	2 - Brown	Prepare run	Low
	3	3 - Gray	Start	Low
	NC	4 - Blue	Shut down	Low
	NC	5 - Pink	Not connected	
	NC	6 - Yellow	Power on	High
	5,14	7 - Red	Ready	High
	6	8 - Green	Stop	Low
	1	9 - Black	Start request	Low
	13, 15		Not connected	

NOTE

START and STOP are connected via diodes to pin 3 of the 3394 connector.


Agilent 1100/1200 to 3396A Integrators

Connector03394-60600	Pin 3394	Pin Agilent 1100/1200	Signal Name	Active (TTL)
	9	1 - White	Digital ground	
	NC	2 - Brown	Prepare run	Low
	3	3 - Gray	Start	Low
	NC	4 - Blue	Shut down	Low
	NC	5 - Pink	Not connected	
	NC	6 - Yellow	Power on	High
	5,14	7 - Red	Ready	High
	1	8 - Green	Stop	Low
	NC	9 - Black	Start request	Low
	13, 15		Not connected	


Agilent 1100/1200 to 3396 Series II / 3395A Integrators

Use the cable **part number: 03394-60600** and cut pin #5 on the integrator side. Otherwise the integrator prints START; not ready.


Agilent 1100/1200 to 3396 Series III / 3395B Integrators

Connector03396-61010	Pin 33XX	Pin Agilent 1100/1200	Signal Name	Active (TTL)
	9	1 - White	Digital ground	
	NC	2 - Brown	Prepare run	Low
	3	3 - Gray	Start	Low
	NC	4 - Blue	Shut down	Low
	NC	5 - Pink	Not connected	
	NC	6 - Yellow	Power on	High
	14	7 - Red	Ready	High
	4	8 - Green	Stop	Low
	NC	9 - Black	Start request	Low
	13, 15		Not connected	


Agilent 1100/1200 to HP 1050, HP 1046A or Agilent 35900 A/D Converters

Connector5061-3378	Pin HP 1050/....	Pin Agilent 1100/1200	Signal Name	Active (TTL)
	1 - White	1 - White	Digital ground	
	2 - Brown	2 - Brown	Prepare run	Low
	3 - Gray	3 - Gray	Start	Low
	4 - Blue	4 - Blue	Shut down	Low
	5 - Pink	5 - Pink	Not connected	
	6 - Yellow	6 - Yellow	Power on	High
	7 - Red	7 - Red	Ready	High
	8 - Green	8 - Green	Stop	Low
	9 - Black	9 - Black	Start request	Low

Agilent 1100/1200 to HP 1090 LC or Signal Distribution Module

Connector01046-60202	Pin HP 1090	Pin Agilent 1100/1200	Signal Name	Active (TTL)
	1	1 - White	Digital ground	
	NC	2 - Brown	Prepare run	Low
	4	3 - Gray	Start	Low
	7	4 - Blue	Shut down	Low
	8	5 - Pink	Not connected	
	NC	6 - Yellow	Power on	High
	3	7 - Red	Ready	High
	6	8 - Green	Stop	Low
	NC	9 - Black	Start request	Low

Agilent 1100/1200 to General Purpose


Connector01046-60201	Pin Universal	Pin Agilent 1100/1200	Signal Name	Active (TTL)
		1 - White	Digital ground	
		2 - Brown	Prepare run	Low
		3 - Gray	Start	Low
		4 - Blue	Shut down	Low
		5 - Pink	Not connected	
		6 - Yellow	Power on	High
		7 - Red	Ready	High
		8 - Green	Stop	Low
		9 - Black	Start request	Low

Auxiliary Cable



One end of this cable provides a modular plug to be connected to the Agilent 1100 Series vacuum degasser. The other end is for general purpose.

Agilent 1100 Series Degasser to general purposes

Connector G1322-81600	Color	Pin Agilent 1100	Signal Name
	White	1	Ground
	Brown	2	Pressure signal
	Green	3	
	Yellow	4	
	Grey	5	DC + 5 V IN
	Pink	6	Vent

10 Identifying Cables
Cable Overview



11 Appendix

General Safety Information [98](#)

The Waste Electrical and Electronic Equipment (WEEE) Directive
(2002/96/EC) [101](#)

Radio Interference [102](#)

Sound Emission [103](#)

Agilent Technologies on Internet [104](#)



General Safety Information

General Safety Information

The following general safety precautions must be observed during all phases of operation, service, and repair of this instrument. Failure to comply with these precautions or with specific warnings elsewhere in this manual violates safety standards of design, manufacture, and intended use of the instrument. Agilent Technologies assumes no liability for the customer's failure to comply with these requirements.

General

This is a Safety Class I instrument (provided with terminal for protective earthing) and has been manufactured and tested according to international safety standards.

Operation

Before applying power, comply with the installation section. Additionally the following must be observed.

Do not remove instrument covers when operating. Before the instrument is switched on, all protective earth terminals, extension cords, auto-transformers, and devices connected to it must be connected to a protective earth via a ground socket. Any interruption of the protective earth grounding will cause a potential shock hazard that could result in serious personal injury. Whenever it is likely that the protection has been impaired, the instrument must be made inoperative and be secured against any intended operation.

Make sure that only fuses with the required rated current and of the specified type (normal blow, time delay, and so on) are used for replacement. The use of repaired fuses and the short-circuiting of fuse holders must be avoided.

CAUTION

Ensure the proper usage of the equipment.

The protection provided by the equipment may be impaired.

→ The operator of this instrument is advised to use the equipment in a manner as specified in this manual.

Some adjustments described in the manual, are made with power supplied to the instrument, and protective covers removed. Energy available at many points may, if contacted, result in personal injury.

Any adjustment, maintenance, and repair of the opened instrument under voltage should be avoided as much as possible. When inevitable, this should be carried out by a skilled person who is aware of the hazard involved. Do not attempt internal service or adjustment unless another person, capable of rendering first aid and resuscitation, is present. Do not replace components with power cable connected.

Do not operate the instrument in the presence of flammable gases or fumes. Operation of any electrical instrument in such an environment constitutes a definite safety hazard.

Do not install substitute parts or make any unauthorized modification to the instrument.

Capacitors inside the instrument may still be charged, even though the instrument has been disconnected from its source of supply. Dangerous voltages, capable of causing serious personal injury, are present in this instrument. Use extreme caution when handling, testing and adjusting.

When working with solvents please observe appropriate safety procedures (e.g. goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet by the solvent vendor, especially when toxic or hazardous solvents are used.

Safety Symbols

Table 12 Safety Symbols

Symbol	Description
	The apparatus is marked with this symbol when the user should refer to the instruction manual in order to protect risk of harm to the operator and to protect the apparatus against damage.
	Indicates dangerous voltages.
	Indicates a protected ground terminal.
	Indicates eye damage may result from directly viewing the light produced by the deuterium lamp used in this product.
	The apparatus is marked with this symbol when hot surfaces are available and the user should not touch it when heated up.

WARNING

A WARNING

alerts you to situations that could cause physical injury or death.

- Do not proceed beyond a warning until you have fully understood and met the indicated conditions.

CAUTION

A CAUTION

alerts you to situations that could cause loss of data, or damage of equipment.

- Do not proceed beyond a caution until you have fully understood and met the indicated conditions.

The Waste Electrical and Electronic Equipment (WEEE) Directive (2002/96/EC)

Abstract

The Waste Electrical and Electronic Equipment (WEEE) Directive (2002/96/EC), adopted by EU Commission on 13 February 2003, is introducing producer responsibility on all Electric and Electronic appliances from 13 August 2005.

NOTE



This product complies with the WEEE Directive (2002/96/EC) marking requirements. The affixed label indicates that you must not discard this electrical/electronic product in domestic household waste.

Product Category: With reference to the equipment types in the WEEE Directive Annex I, this product is classed as a "Monitoring and Control instrumentation" product.

Do not dispose off in domestic household waste

To return unwanted products, contact your local Agilent office, or see www.agilent.com for more information.

Radio Interference

Cables supplied by Agilent Technologies are screened to provide optimized protection against radio interference. All cables are in compliance with safety or EMC regulations.

Test and Measurement

If test and measurement equipment is operated with unscreened cables, or used for measurements on open set-ups, the user has to assure that under operating conditions the radio interference limits are still met within the premises.

Sound Emission

Manufacturer's Declaration

This statement is provided to comply with the requirements of the German Sound Emission Directive of 18 January 1991.

This product has a sound pressure emission (at the operator position) < 70 dB.

- Sound Pressure $L_p < 70$ dB (A)
- At Operator Position
- Normal Operation
- According to ISO 7779:1988/EN 27779/1991 (Type Test)

Agilent Technologies on Internet

For the latest information on products and services visit our worldwide web site on the Internet at:

<http://www.agilent.com>

Select Products/Chemical Analysis

It will provide also the latest firmware of the Agilent 1200 Series modules for download.

Index

A

accessory kit 21, 82
 Agilent on internet 104
 algae growth 39
 ambient non-operating temperature 15
 ambient operating temperature 15
 analog output 17, 17
 analog
 cable 88
 APG remote connector 25
 AUX output 25
 auxiliary cable 52
 auxiliary
 cable 88, 95

B

BCD
 cable 88
 bench space 14
 block diagram 10
 bottle-head assembly 27
 bottle 27
 busy condition 48

C

cable
 analog 88
 auxiliary 20, 52, 88, 95
 BCD 88
 CAN 25
 external contacts 89
 GPIO 25, 89
 interface 24
 LAN 89

power 20, 49, 59
 remote 20, 49, 88, 90
 RS-232 89

cables
 overview 88
 CAN cable 25
 changing solvents 31
 ChemStation 25
 cleaning the instrument 58
 condensation 14
 connecting tube 21
 control assembly 76
 control circuit 6
 cover 80

D

dead volume 22
 delay volume 22
 delivery checklist 20
 dimensions 15, 15

E

electronics
 how does it work 7
 electrostatic discharge (ESD) 57, 67
 environment 12, 14
 error condition 48
 error state 50
 exchanging
 control assembly 76
 sensor assembly 76
 solenoid valve 75
 status light pipe 64
 vacuum pump 73

vacuum chamber 71, 71

external contact
 cable 89

F

flow connections 27, 27
 flow rate 17
 frequency range 15
 front of the vacuum degasser 25
 fuse 49
 fuses 12

G

GPIO cable 25
 GPIO
 cable 89

H

hardware symptoms 49
 highest injection precision 36
 highest retention time reproducibility 36
 humidity 15

I

information
 on sound emission 103
 installation 11, 19, 24, 27
 interface cable 24
 internal leak 48
 internet 104
 introduction to the degasser 6

Index

L

- laboratory bench 14
- lamps
 - off, degasser appears dead 49
 - status indicator is red 50
 - status indicator is yellow, vacuum pump not running 51
 - status indicator is yellow, vacuum pump runs continuously 52

LAN

- cable 89
- line frequency 15
- line voltage 15

M

- main cover, assembling 62
- materials in contact with solvent 17
- maximum flow rate 17
- maximum sensitivity 36
- mounting tool 21
- mounting tool 28, 28

N

- non-operating altitude 15
- number of channels 17

O

- operating altitude 15
- operational hints 30

P

- parts
 - accessory kit 82
 - cover 80
 - damaged 20
 - degasser 84
 - light pipes 81
 - missing 20

- sheet metal kit 86
- performance specifications 17
- pH range 17
- power cable 59
- power consumption 15, 15
- power light pipe 81
- power supply lamp 47
- power supply 12
- power switch 24
- power cable 49
- power-input socket 13
- pressure sensor readings 25
- pressure sensor readings 52
- pressure sensor 6
- priming
 - with a pump 30, 32
 - with a syringe 30
- principle of operation 6

R

- radio interference 102
- ready condition 48
- rear of the vacuum degasser 25
- remote cable 49
- remote
 - cable 88, 90
- RS-232
 - cable 89

S

- safety class I 98
- safety
 - general information 98, 98
- screwdriver pozidriv #1 59, 71
- screwdriver pozidriv#1 64
- sensor assembly 76
- sheet metal kit 86
- shipping container 20

- site requirements 12
- snap fastener 27
- solenoid valve 75
- solvent cabinet 59
- solvent filters
 - checking 39
 - cleaning 39
 - prevent blocking 39
- solvent outlet tubes 27
- solvent bottle 27
- solvent cabinet 27, 27
- solvent information 41
- sound emission 103
- specification
 - performance specifications 17
 - physical 15
- stack configuration 22
- status indicator 47
- status light pipe 81
- status light pipe 64
- status indicator 50, 51, 52
- syringe adapter 21, 31, 33
- syringe 21, 27

T

- tool 28
 - mounting tool 28
- top cover, removing and refitting 59
- transportation 33

V

- vacuum pump 6, 73
- vacuum chamber 71, 71
- voltage range 15
- voltage selector 12
- volume per channel 17

Index

W

wall socket 13
waste tube 21
weight 15, 15
when to use a degasser 43

In This Book

This manual contains service information about the Agilent 1200 Series vacuum degasser. The manual describes the following:

- introduction,
- site requirements and specifications,
- installing the vacuum degasser,
- using the vacuum degasser,
- troubleshooting and diagnostics,
- maintenance,
- repairs,
- parts and materials for maintenance,
- parts for repair,
- identifying cables,
- appendix.

© Agilent Technologies 2006, 2007-2008

Printed in Germany
11/08



G1322-90110



Agilent Technologies