



create the difference

User Manual



phyaction u

© 2004, GymnaUniphy N.V.

All rights reserved. Nothing from this publication may be copied, stored in an automated data file, or made public, in any form or in any way, be it electronically, mechanically, by photocopying, recordings or in any other way, without prior written permission from GymnaUniphy N.V.

User Manual Phyaction U

Device for ultrasound therapy

Manufacturer	GymnaUniphy N.V.
Main office	Pasweg 6A B-3740 BILZEN
Telephone	+(32) (0)89-510.510
Fax	+(32) (0)89-510.511
E-mail	info@gymna-uniphy.com
Website	www.gymna-uniphy.com

Version 1.1
February 2005

Abbreviations

EMC	Electromagnetic Compatibility
ESD	Electrostatic Discharge
HAC	Hospital Antiseptic Concentrate
US	Ultrasound

Symbols on the equipment



Read the manual.

Symbols in the manual



Warning or important information.

TABLE OF CONTENTS

1	SAFETY	7
1.1	PURPOSE	7
1.2	SAFETY INSTRUCTIONS	7
1.3	MEDICAL DEVICES DIRECTIVE	9
1.4	LIABILITY	9
2	INSTALLATION	11
2.1	RECEIPT	11
2.2	PLACING AND CONNECTION	11
2.3	PERFORMING THE FUNCTIONAL TEST	11
2.4	SETTING CONTRAST, LANGUAGE AND STAND-BY TIME	11
2.5	USE IN COMBINATION WITH AN OTHER DEVICE	11
2.6	TRANSPORT AND STORAGE	12
2.7	RESELLING	12
3	DESCRIPTION OF THE EQUIPMENT	13
3.1	PHYACTION U AND STANDARD ACCESSORIES	13
3.2	COMPONENTS OF PHYACTION U	14
3.3	DISPLAY	15
3.4	DISPLAY SYMBOLS	16
3.5	PARAMETER SYMBOLS	16
4	OPERATION	17
4.1	THERAPY SELECTION	17
4.2	SELECTION BY THE THERAPY MENU	17
4.3	SELECTION BY THE GUIDE MENU	18
4.4	PERFORMING THERAPY	20
4.5	ULTRASOUND THERAPY	21
4.6	COMBINATION THERAPY	24
4.7	MEMORY	26
4.8	SYSTEM SETTINGS	28
5	INSPECTIONS AND MAINTENANCE	33
5.1	INSPECTIONS	33
5.2	MAINTENANCE	34
6	MALFUNCTIONS, SERVICE AND GUARANTEE	35
6.1	MALFUNCTIONS	35
6.2	SERVICE	36
6.3	GUARANTEE	36

6.4	TECHNICAL LIFE TIME	37
7	TECHNICAL INFORMATION	39
7.1	GENERAL	39
7.2	ULTRASOUND THERAPY	39
7.3	ENVIRONMENTAL CONDITIONS	40
7.4	TRANSPORT AND STORAGE	41
7.5	STANDARD ACCESSORIES	41
7.6	OPTIONAL ACCESSORIES	42
8	APPENDICES	43
8.1	US HEAD PLACEMENTS	43
8.2	EMC DIRECTIVE	44
8.3	TECHNICAL SAFETY INSPECTION	48
8.4	DISPOSAL	50
9	REFERENCE	51
9.1	FUNCTION OVERVIEW	51
9.2	LITERATURE	53
9.3	TERMINOLOGY	53
10	INDEX	55

1 SAFETY

1.1 Purpose

The Phyaction U is intended solely for medical applications. You can use the Phyaction U for ultrasoundtherapy. The device is suited for continuous use.

1.2 Safety instructions

1.2.1 General



- Only qualified people who are trained in the application of the therapies may use the appliance.
- Only a technician authorised by GymnaUniphy N.V. may open the equipment or the accessories.
- Follow the instructions and directions in these user instructions.
- Place the equipment on a horizontal and stable base.
- Keep the ventilation openings at the bottom and rear of the equipment free.
- Do not place any objects on the equipment.
- Do not place the equipment in the sun or above a heat source.
- Do not use the equipment in a damp area.
- Do not let any liquid flow into the equipment.
- Do not disinfect or sterilise the equipment. Clean the equipment with a dry or moistened cloth. *See §5.*
- Only treat patients with electrical implants (pacemaker) after obtaining medical advice.
- The 'Directive on Medical Devices' from the European Commission (93/42/EEG) requires that safe devices are used. It is recommended to perform a yearly technical safety inspection. *See §5.1.2.*
- For optimum treatment, a patient investigation must first be performed. On the basis of the findings of the investigation, a treatment plan with objectives will be formulated. Follow the treatment plan during the therapy. This will limit possible risks, related to the treatment, to a minimum.
- Always keep these user instructions with the equipment.

1.2.2 **Electrical safety**



- Only use the equipment in an area with facilities that meet the applicable legal regulations.
- Connect the equipment to an outlet with a protective earth terminal. The outlet must meet the locally applicable requirements for medical areas.

1.2.3 **Prevention of explosion**



- Do not use the equipment in an area where combustible gases or vapours are present.
- Switch off the equipment when it is not used.

1.2.4 **Electro Magnetic Compatibility**



- Medical electrical equipment requires special precautions for Electro Magnetic Compatibility (EMC). Follow the instructions for the installation of the equipment. *See §2*
- Do not use mobile telephones or other radio, shortwave, or microwave equipment in the vicinity of the equipment. This kind of equipment can cause disturbances.
- Only use the accompanying accessories that are supplied by GymnaUniphy. *See §7.5 and §7.6.* Other accessories can lead to an increased emission or a reduced immunity.

1.2.5 **Ultrasound therapy**



- Move the US head evenly over the skin during the treatment. This prevents internal burns.
- The US treatment heads are exchangeable. The device detects the characteristics and supplies the right power at the right frequency.
- Handle the US heads carefully. With rough handling, the characteristics can change. Test the US head if it falls on the ground or knocks against something. *See §5.1.1.*
- Check the US head at least once a month. During the check, look for dents, cracks and other damage that could allow liquids to ingress. Check whether the insulation of the cable is still intact. Check whether all pins are present and straight in the connectors. Replace the US head if the head, the cable or the connector is damaged. *See §5.1.*

1.3 Medical Devices Directive

The device complies with the essential requirements of the Medical Device Directive of the European Committee (93/42/EEC) as most recently changed.

1.4 Liability

The manufacturer cannot be held liable for injury to the therapist, the patient or third parties, or for damage to or by the equipment used, if for example:

- an incorrect diagnosis is made;
- the equipment or the accessories are used incorrectly;
- the user instructions are wrongly interpreted or ignored;
- the equipment is badly maintained;
- maintenance or repairs are performed by people or organisations that are not authorised to do so by GymnaUniphy.

Neither the manufacturer nor the local GymnaUniphy dealer can be held liable, in any way whatsoever, for the transfer of infections by accessories.

2 INSTALLATION

2.1 Receipt

1. Check whether the equipment has been damaged during transport.
2. Check whether the accessories are intact and complete. *See §7.5 and §7.6.*
 - Inform your supplier of any damage or defects by no later than within 3 working days after receipt. Report the damage by telephone, fax, e-mail or letter.
 - Do not use the equipment if it is damaged or defective.












2.2 Placing and connection

1. Place the equipment on a horizontal and stable base.
 - Keep the ventilation openings at the bottom and rear of the equipment free.
 - Do not place the equipment in the sun or above a heat source.
 - Do not use the equipment in a wet area.
2. Check whether the mains voltage that is stated on the rear of the equipment corresponds with the voltage of your mains supply. The equipment is suited for a nominal mains voltage from 100 V to 240 VAC / 50-60 Hz.
3. Connect the device to an outlet with protective earth terminal.

2.3 Performing the functional test

1. Switch the equipment on with the switch at the rear of the equipment.
2. When the equipment is switched on, it automatically performs a test.

2.4 Setting contrast, language and stand-by time

1. Press . The **System settings** menu appears. *See §4.8.*
2. Select **Contrast** with the corresponding blue key , 1st key in the row.
3. If necessary, change the contrast with  and .
4. Select **Language** with the corresponding blue key .
5. If necessary, change the language with  and .
6. Select **Stand-by time** with the corresponding blue key .
7. If necessary, change the stand-by time with  and .
8. Press  to return to the **Guide** menu.

2.5 Use in combination with an other device

The Phyaction U can be used in combination with:

- the Phyaction E. *See §4.6.*
- the Phyaction I. *See §4.6.*

2.6 Transport and storage

Take account of the following matters if the equipment has to be transported or stored:


- Transport or store the equipment in the original packaging.
- The maximum period for transport or storage is: 15 weeks.
- Temperature: -20 °C to +60 °C.
- Relative humidity: 10% to 100%.
- Atmospheric pressure: 200 hPa to 1060 hPa.

2.7 Reselling

This medical equipment must be traceable. The equipment, the US head and some other accessories have a unique serial number. Provide the dealer with the name and address of the new owner.

3 DESCRIPTION OF THE EQUIPMENT


3.1 Phyaction U and standard accessories



1. Phyaction U. See §3.2.
2. Power cord


3. Contact gel
4. US head

3.2 Components of Phyaction U



- | | |
|------------------------------------|--|
| 1. Display. See §3.3. | 13. Indicator lamp device on/off |
| 2. Select menu option or parameter | 14. Intensity of channel A |
| 3. Scroll through the list/numbers | 15. Indication: Read manual |
| 4. Increase or set a parameter | 16. Connectors for US head |
| 5. Decrease or set a parameter | 17. Indication: Floating patient circuit |
| 6. Therapy menu | 18. On/off switch |
| 7. Guide menu | 19. Fuse holder |
| 8. Memory menu | 20. Connection to mains supply |
| 9. System settings menu | 21. Type plate |
| 10. Back | 22. Ventilation opening |
| 11. Pause | 23. Combination plug connector |
| 12. Stop | |

3.3 Display



- | | |
|---------------------------------------|---|
| 1. Selected channel | 8. Screen for channel A (here, ultrasound therapy). See §4.5.3. |
| 2. Title of the screen | 9. Remaining treatment time |
| 3. Program number | 10. \hat{I}_{set} |
| 4. Ultrasound therapy | 11. Ppk |
| 5. Type of US head | 12. Contact of the US head |
| 6. Parameters of the selected channel | 13. Numbers, selection with the blue keys ● below display. |
| 7. Explanation or recommendation | |

3.4 Display symbols



Ultrasound therapy



Treatment time

A

Channel A



0:00 Treatment completed

3.5 Parameter symbols



10%

US duty cycle 10%

1:10 ms

US on : period time 10%



20%

US duty cycle 20%

2:10 ms

US on : period time 20%



30%

US duty cycle 30%

3:10 ms

US on : period time 30%



40%

US duty cycle 40%

4:10 ms

US on : period time 40%



50%

US duty cycle 50%

5:10 ms

US on : period time 50%



100%

US duty cycle 100%

10:10 ms

US on : period time 100%

\hat{I}_{set}

Set US intensity



US head, ERA 4 cm²

P_{pk}

Peak US output power



US head, ERA 1 cm²




W/cm^2

Unit of the set US intensity

4 OPERATION


4.1 Therapy selection




You can select a therapy with different keys:

- **Therapy Menu**  : Select a therapy method. See §4.2.
- **Guide Menu**  : Gives access to:
 - **Objectives**: Select a therapy on the basis of an objective. See §4.3.1.
 - **Indication list**: Select a therapy on the basis of a medical indication. See §4.3.2.
 - **Program number**: Select a certain program number. See §4.3.3.
 - **Diagnostic programs**: Perform a diagnosis. See §4.3.4.
 - **Contra indications**: Display an overview with contra indications for the ultrasound therapy. See §4.3.5.
- **Memory Menu**  : Select a saved therapy. See §4.7.

Besides this, you can change the system settings. See §4.8.


4.2 Selection by the Therapy menu


Press . The **Ultrasound** screen appears.

A THERAPY		31	A		
Ultrasound therapy 1 Treatment time 7:30 2 ERA 4.0 cm ² 3 US frequency 1.1 MHz 4 Duty cycle 100 % On time : Period time 10:10 ms 5 Contact detection beep ON			7:30  Iset 0.0W/cm ² Ppk 0.00W <hr/>		
1	2	3	4	5	


4.3 Selection by the Guide menu


4.3.1 Therapy selection via objectives

1. Press  to go to the Guide menu.
2. Select Objectives.
3. Select Ultrasound therapy or Phonophoresis.
4. Follow the on screen options to select the desired treatment.





4.3.2 Therapy selection via indication list

1. Press .
2. Select Indication list.
3. Use ◀ and ▶ to select the following indications. See §9.1.4.
4. Select the desired indication.
 - US: Ultrasound therapy
5. With selection via indication list you can view the placement.
 - Select **US head placement**.
 - If necessary, select the location. You get an advice to place the US head.
 - If available, select a number for the precise anatomic location. See §8.1.




4.3.3 Program number selection


1. Press .
2. Select Program number.
3. Select the desired program with Δ or ▽. See §9.1.
4. Select 1.




4.3.4 Diagnostic program selection


With the diagnostic programs, you can search for stress fractures.

1. Press .
2. Select Diagnostic programs.
3. Select Diagnosis stress fracture.




4.3.5 Contra indication selection


1. Press .
2. Select Contra indications.
3. Scroll through the text with ◀ or ▶.





4.4 Performing therapy


4.4.1 Set and start therapy

1. Press  to go to the Guide menu.
2. Select the desired menu item until the treatment appears.
3. Select the desired parameters. You can only change the pre-numbered parameters.
4. Set the **Treatment time** as follows: Select treatment time once to set the minutes, select treatment time twice to set the seconds.
5. Change the value of the parameter with Δ and ∇ . The setting range of the parameter is shown at the bottom of the screen. You can change the parameter as long as the parameter has a black background.
6. Rotate intensity knob to start the treatment and to set the desired intensity. The set intensity is displayed in the screen.





4.4.2 Opening the intensity screen


1. Set the treatment. See §4.4.1.
2. Rotate intensity knob to start the treatment.
3. Once the treatment is started ($I_{\text{set}} \geq 0.1 \text{ W/cm}^2$) press the therapy key  to go to the intensity screen.
4. Press  to return to the setting menu.



4.4.3 Temporary interruption of treatment

1. Press  during the treatment. The treatment time is stopped. **Pause** appears on the screen. The parameter settings are retained.
2. Press  to restart the treatment. The intensity now increases gradually to the set level and the treatment time continues again.

4.4.4 Immediately stop treatment

1. Press . All active treatments are stopped immediately. **Stop** appears on the screen. The parameter settings are retained.
2. Set the intensity of the channel again to continue the treatment.


4.5 Ultrasound therapy

4.5.1 Performing ultrasound therapy




Move the US head evenly over the skin during the treatment. This prevents internal burns.

1. Connects the US head into one of the two connectors « of the Phyaction U. You can connect two US heads, but only one US head can be in operation at one time. The device detects which US head is connected to the connector «.
2. Select the desired ultrasound therapy. With **Indication list** treatments, the parameter **Head placement** is available.
3. Set the parameter **ERA** to 1 or 4 cm². The corresponding US head is selected, the green indication led on the US head is on.
4. Apply contact gel to the skin to be treated and to the US head.
5. Place the head on the skin.
6. Rotate intensity knob to start the ultrasound therapy.
7. Move the US head evenly over the skin during the treatment. This prevents internal burns.
8. Check the patient's reaction and the effect of the treatment. Repeat this check regularly during the treatment.
9. The equipment stops the treatment and indicates that the treatment is completed.




4.5.2 Phonophoresis

Phonophoresis is used to enhance transdermal transport of several drugs, especially anti-inflammatory NSAID and local anesthetics.

1. Use the drugs (gel ointment) instead of the US contact gel.
2. Press .
3. Select **Objectives**.
4. Select **Phonophoresis**. The frequency is 1 MHz, the duty cycle is 20% and the time is at least 5 minutes.

4.5.3 Read-out values for ultrasound therapy

1. Channel
2. Ultrasound therapy
3. Type of US head
4. Remaining treatment time
5. \hat{I}_{set}
6. P_{pk}
7. Contact of the US head



Contact of the US head

The contact of the US head with the skin:

- [] — — — — : Bad contact, US head switched off (0 W).
- [■] — — — — : Bad contact.
- [■] [■] — — — : Sufficient contact.
- [■] [■] [■] — — : Good contact.
- [■] [■] [■] [■] — : Very good contact

Test the US head if its conduction is bad. See §5.1.1.

\hat{I}_{set} (W/cm²)

The power (W) of the US head per cm².

P_{pk} (W)

The peak power of the US head ($\hat{I}_{set} * ERA$). The peak power delivered therefore depends on the size of the US head and the contact with the skin. This value is 0.0 W if the contact with the skin is bad. In this case, the ultrasound treatment of the equipment is stopped to prevent overheating of the transducer.

4.5.4 Parameters for ultrasound therapy

Treatment time (mm:ss)

The duration of the treatment.

Duty cycle (10, 20, 30, 40, 50%, continuous)

Ratio of the pulse duration to the period duration.

- **Continuous:** Continuous ultrasound (100%).
- 10, 20, 30, 40, 50%: Pulsating ultrasound.

Select a high duty cycle for an intensive treatment. Select a low duty cycle for a mild treatment.

ERA (cm²)

The effective radiating area expressed in cm² of the treatment head connected. This area equals the cross-sectional area of the beam at the treatment surface. The ERA depends on the frequency. This parameter remains empty if no US head is connected.

Head placement

Instructions for placing the US head. This is only available with treatment selection via **Indication list**.

US frequency (MHz)

The frequency of the US head. The absorption at a US frequency of 3 MHz is three times higher and the penetration depth is three times less than at a US frequency of 1 MHz. Use 3 MHz for superficial tissue and 1 MHz for deeper tissue.

4.5.5 Indicator light of the US head

The indicator light of the US head provides the following information.

Indication light

Blinking green

Continuous green

Continuous yellow

Alternating yellow/green

Blinking yellow

Situation

The US head is properly connected.

The US head is selected.

The US-emission is in progress.

Bad contact of the US head with the skin.

End of the treatment

4.6 Combination therapy


4.6.1 Performing combined therapy



With combination therapy, a maximum current density of $2.0 \text{ mA}_{\text{rms}}/\text{cm}^2$ is advised. Exceeding this current density can result in skin irritation and burns. The intensity depends on the surface area of the US head. For US U92 (9 cm^2), the current setting may be a maximum of $18 \text{ mA}_{\text{rms}}$; for US U91 (3 cm^2), a maximum of $6 \text{ mA}_{\text{rms}}$.

1. Select an alternating current shape (TENS or interferential) on the Phyaction E or I.
2. Put the current shape in CV mode.
3. Select a **Ultrasound therapy** on the Phyaction U.
4. Connect the electrode and the US-head. See §4.6.2.
5. Place the electrode on the patient. See the User Manual of the Phyaction E or I.
6. Apply contact gel to the skin to be treated and to the US head.
7. Place the head on the skin.
8. Rotate intensity knob from the Phyaction U to start the ultrasound therapy.
9. Rotate intensity knob A or B from the Phyaction E or I to start the electrotherapy. Set the desired voltage.
10. Check the contact between the US head and the skin. The following indications can indicate a bad contact on the Phyaction U:
 - The treatment stops.
 - The peak power of the ultrasound treatment goes to 0.0 Watt.
11. Check the patient's reaction and the effect of the treatment. Repeat this check regularly during the treatment.
12. The equipment stops the treatment and indicates that the treatment is completed.

4.6.2 Connection for combination therapy





1. Connect the two-ply electrode to the connector **Y**A or **Y**B of the Phyaction E or I.
2. Connect the electrode to the red plug of the two-ply electrode cable.
3. Connect the black plug of the two-ply electrode cable via the 4 mm V/V test connector on the combination therapy cord. *See §7.6.*
4. Connect the combination therapy cord to the combination plug connector of the Phyaction U.
5. Connect the US-head to a US-connector of the Phyaction U.

4.7 Memory


You can save 50 of your own programs for later use: programs 500 up to and including 549. You can modify these programs for much-used settings for a certain patient.

4.7.1 Saving a program

1. Select a therapy. See §4.1.
2. Change the settings for the patient. See §4.4.
3. Press .
4. Select **Save**.
5. Select a free program number or overwrite an existing program number.
If desired, scroll through the list with ◀ or ▶.





6. Enter the name of the program. Use the name or the number of the patient, for example.
 - Select a character with Δ and ∇ .
 - Select **Move cursor to the left/right** to change the cursor position.
7. Select **Ready** and save.




4.7.2 Selecting a saved program


Selecting a program by the list

1. Press .
2. Select **Recall by list**.
3. Select the desired program.
If necessary, scroll through the list with ◀ or ▶.




Selecting a program by the number


1. Press .
2. Select Recall by number.
3. Select the desired program with Δ or ∇ .
4. Select Go to selected number.




4.7.3 Erase a saved program


Erase a program by the list

1. Press .
2. Select Erase by list.
3. Select the desired program.
If necessary, scroll through the list with \blacktriangleleft or \blacktriangleright .
4. Select Erase memory number to delete the program.



Erase a program by the number


1. Press .
2. Select Erase by number.
3. Select the desired program with Δ or ∇ .
4. Select Erase selected number twice to delete the program.



4.8 System settings

With the system settings, you can adapt the Standard settings of the equipment. You cannot change the system settings during a therapy.

4.8.1 Changing the system settings

1. Press . The System settings menu appears.
2. Change the desired system setting.

A SYSTEM SETTINGS		
1 Contrast	10	
2 Language	English	
3 Sound settings		
4 Stand-by time	10:00	
5 Text for start up screen		
6 First screen	GUIDE	
7 System information		
8 Error history		
9 Counter working hours		
10 Reset Menu		
11 Stop time if bad US	ON	
1	2	3 4 5 6 7 ▶

4.8.2 Parameters

Contrast (1 - 20)

The contrast of the display.

Language

The language selection: select the language with which the read-out must work.

Sound settings

Sound settings. See §4.8.3.

Stand-by time (5, 10,15, 20 minutes, off)

If the device is not used during the stand-by time, the device goes to the stand-by mode. Press any key to reactive the device.

Text for start up screen

The text that appears in the top of the start up screen, after the equipment is switched on. See §4.8.6.

First screen (guide menu, therapy menu)

The first screen you see when activating the device.

System information

System information of the equipment

Always have this information available when you contact the technical service department.

Error history

The total number of error reports that the equipment has had and details about the last 10 error reports.

Always have this information available when you contact the technical service department.

Counter working hours (hours, minutes, sec.)

The time that the accessories for electrotherapy or ultrasound therapy have been in use. For this, the output of the channel must have been higher than zero.


Reset menu

- **Reset working hours:** Set the number of working hours of a plate electrode or an US head to zero.
- **Change therapy programs:** Change the settings from the programs in the Therapy menu. See §4.8.5.
- **Erase total memory:** Restores the standard settings of the standard programs and of the edited programs.

Stop time if bad US (on, off)

When there is a bad US contact, the treatment time counter stops. When the contact is restored, the counting continues.

4.8.3 Setting the sound

1. Press .
2. Select **Sound settings**.
3. Change the desired sound setting.

A SYSTEM SETTINGS	
Sound settings	
1 Beep at end of treatment	OFF
2 Beep at pressing of a key	OFF
3 Beep volume	3
4 US bad contact beep	OFF
1	2
3	4

4.8.4 parameters sound settings

End of treatment

On: A sound signal will be heard at the end of the treatment.

Pressing a key

On: A sound signal will be heard every time a key is pressed.

Beep volume (min.1, standard 5, max.10)

The volume of the sound signals.



US bad contact


On: A sound signal will be heard if the US head does not make good contact with the skin.

4.8.5 Change therapy programs

Save new therapy program settings


Change the program to your required settings.

1. Use the **Therapy** menu  to select a program.
2. Make the changes in the program.
3. Press .
4. Select **Reset Menu**.
5. Select **Change therapy programs**.
6. Select **Save new ther. progr. settings** twice to change the program settings.



Restore this therapy program

Change the program back to the manufacture's settings.

1. Use the **Therapy** menu  to select a program.
2. Go to the **Change therapy programs** menu.
3. Select **Restore this therapy program** twice.


Restore all therapy programs


Change all therapy programs back to the manufacture's settings.

1. Go to the **Change therapy programs** menu.
2. Select **Restore ALL therapy program** twice.

4.8.6 Set text for start up screen

You can set your own text for the start up screen. For example, you can put your name or address information here.

1. Press .
2. Select **Text for start up screen**.








3. Enter the name for the start up screen.
 - Select a character with Δ and ∇ .
 - Select **Move cursor to the left/right** to change the cursor position.
 - Select **Move to next line** to enter a line.
4. Select **Ready and save**.

5 INSPECTIONS AND MAINTENANCE

5.1 Inspections

Component	Check	Frequency
US head	Dents, cracks or other damage	At least 1x per month
	Test US head. See §5.1.1.	With bad operation or at least 1x per year
Cable of US head	Damage Pins in connector straight	At least 1x per month
Equipment	Technical safety inspection. See §5.1.2.	At least 1x per year

5.1.1 US head test

Test the US head if its conduction is bad. This is the case when the indication bar for the Ppk value displays      or

    .

1. Select an ultrasound therapy.
2. Place the US head in a bowl with water.
3. Rotate intensity knob to start the treatment.
4. Check in the screen of the channel to see if the Ppk value is increasing.
5. Contact your local GymnaUniphy dealer if the indication bar still displays

     or     .

5.1.2 Technical safety inspection

The 'Directive on Medical Devices' from the European Commission (93/42/EEG) requires that safe devices are used. It is recommended to perform a yearly technical safety inspection. If the legislation in your country or your insurer prescribes a shorter period, you must adhere to this shorter period.



- Only a technician authorised by GymnaUniphy N.V. may open the equipment or the accessories.
- The inspection may only be performed by a suitably qualified person. In some countries this means that the person must be accredited.

Inspection points

The technical safety inspection contains the following tests:

1. Test 1: General: Visual inspection and check on the operating functions
2. Test 2: Ultrasound therapy
3. Test 3: Electrical safety inspection: measurement of the earth leakage current and patient leakage current according to DIN/VDE 0751-1 ed. 2.0.

Inspection result

1. A registration must be maintained of the technical safety inspections. Use the inspection report in the appendix for this purpose. See §8.3.
2. Copy this appendix.
3. Complete the copied appendix.
4. Keep the inspection reports for at least 10 years.

The inspection is successful if all inspection items are passed.

Repair all faults on the equipment before the equipment is put back into operation.

By comparing the registered measurement values with previous measurements, a possible slowly-deteriorating deviation can be ascertained.

5.2 Maintenance

Component	Check	Frequency
US head	Cleaning. See §5.2.1.	After each use



Accessories that come in contact with the body of the patient must be washed with pure water after the disinfection to prevent allergic reactions.

5.2.1 Cleaning the US head

1. Clean the US head with a lightly moistened soft cloth.
2. Disinfect the treatment surface with a cotton bud that is soaked in a 10% HAC solution.
3. Rinse the US head thoroughly with clean water.

6 MALFUNCTIONS, SERVICE AND GUARANTEE

6.1 Malfunctions

Component	Problem	Solution
Phyaction U	Equipment cannot be switched on	See §6.1.1.
	Equipment does not react to commands or a fault report appears	See §6.1.3.
	Foreign language on the screen	Change the language. See §4.8.2.

6.1.1 *Equipment cannot be switched on*

1. Check if the mains voltage has failed.
2. Check if the main switch is switched on ("I").
3. Check if the power cord and the fuses are in order. If necessary, replace the fuse. See §6.1.2.
4. Contact your dealer if the equipment still cannot be switched on.

6.1.2 *Replacing a fuse*

1. Switch the main switch off ("O").
2. Unplug the power cord from the equipment.
3. Pull the fuse holder carefully out of the equipment. If necessary, use a screwdriver.
4. Replace the fuse. If necessary, order new fuses from your dealer.
5. Install the fuse holder and plug in the power cord.
6. Switch the main switch on again ("I").

6.1.3 *Equipment does not react to commands or a fault report appears*

The safety system of the equipment has ascertained a fault. You cannot continue to work. An instruction usually appears on the screen.

1. Disconnect the connection to the patient.
2. Switch the main switch off ("O").
3. Wait 5 seconds and switch the main switch on again ("I").
4. Contact your dealer if the error message reappears.

6.2 Service



- Only a technician authorised by GymnaUniphy N.V. may open the equipment or the accessories to perform repairs. The equipment does not contain any components that may be replaced by the user.
- If possible, open the screen with the system settings before you contact the technical service department. See §4.8.

Service and guarantee are provided by your local GymnaUniphy dealer. The conditions of delivery of your local GymnaUniphy dealer apply.

If you have qualified technical personnel that are authorised by GymnaUniphy to perform repairs, your dealer can provide diagrams, spare parts lists, calibration instructions, spare parts and other information on request, for a fee.

6.3 Guarantee

GymnaUniphy and the local GymnaUniphy dealer declares itself to be solely responsible for the correct operation when:

- all repairs, modifications, extensions or adjustments are performed by authorised people;
- the electrical installation of the relevant area meets the applicable legal regulations;
- the equipment is only used by suitably qualified people, according to these user instructions;
- the equipment is used for the purpose for which it is designed;
- maintenance of the device is regularly performed in the way prescribed. See §5.;
- the technical life time of the equipment and the accessories is not exceeded;
- the legal regulations with regard to the use of the equipment have been observed.

The guarantee period for the equipment is 2 (two) years, beginning on the date of purchase. The date on the purchase invoice acts as proof. This guarantee covers all material and production faults. Consumables, such as sponges, adhesive electrodes and rubber electrodes, do not fall under this guarantee period.

This guarantee does not apply to the repair of defects that are caused:

- by incorrect use of the equipment,
- by an incorrect interpretation or not accurately following these user instructions,
- by carelessness or misuse,
- as a consequence of maintenance or repairs performed by people or organisations that are not authorised to do so by the manufacturer.

6.4 Technical life time

The expected life time of the equipment is 10 years, calculated from the date of manufacture. See the type plate for this information.

In so far as possible, GymnaUniphy will supply service, spare parts and accessories for a period of 10 years from the date of manufacture.

7 TECHNICAL INFORMATION

7.1 General

Dimensions Phyaction U

(w x h x d) 265 x 275 x 122 mm

Weight Phyaction U 3,650 kg

Weight including accessories 4,6 kg

Mains voltage 100 - 240 VAC, 50 - 60 Hz

Maximum power, in operation 85 VA

Safety class Class I (earthed socket required)

Insulation Type BF (floating patient circuit)

Fuses 2 x T2AL250 V

7.2 Ultrasound therapy

7.2.1 General

Insulation classification Type BF

Peak power 0 - 2 W/cm², duty cycle = 100%

0 - 3 W/cm², duty cycle < 100%

Accuracy of intensity ± 10% of maximum at set values above
10% of this maximum

Treatment time 0 - 30 min.

Deviation of time clock < 0,5%

Modulation frequency 100 Hz

Modulation type CW (rectangular on/off)

Repetition period of pulses 10 ms

7.2.2 Modulation and pulse duration

Modulation duty cycle	100	50	40	30	20	10	%
Pulse time	∞	5	4	3	2	1	ms
Ratio of $p_{tm} - p$	1	2	2,50	3,33	5	10	

7.2.3 US heads

US head, model U92			
Acoustic operating frequency	1,1	3,2	MHz
Output power	8,0	9,6	W
Effective intensity of output voltage	2,0	2,0	W/cm ²
Effective Radiating Area (ERA)	4,0	4,8	cm ²
Beam Non-uniform Ratio (BNR)	7,5	7,5	
Maximum intensity of beam	15,0	15,0	W/cm ²
Beam type	Collimated	Collimated	

US head, model U91			
Acoustic operating frequency	1,1	3,2	MHz
Output power	1,2	2,0	W
Effective intensity of output voltage	2,0	2,0	W/cm ²
Effective Radiating Area (ERA)	0,6	1,0	cm ²
Beam Non-uniform Ratio (BNR)	5,0	5,0	
Maximum intensity of beam	10,0	10,0	W/cm ²
Beam type	Divergent	Collimated	

7.3 Environmental conditions




Temperature:	+10 °C to +40 °C
Relative humidity	30% to 75%
Atmospheric pressure	700 hPa to 1060 hPa

7.4 Transport and storage

Transport weight	5,5 kg
Storage temperature	-20 °C to +60 °C
Relative humidity	10% to 100%, including condensation
Atmospheric pressure	200 hPa to 1060 hPa
Transport classification	Single pieces, by post




The transport and storage specifications apply to equipment in the original packaging.

7.5 Standard accessories

	Quantity	Description	Art. no.
	1	US head, 1/3 MHz - ERA 4 cm ² incl. holder	323.584
	1	Contact gel, 500 ml	114.827
	1	Power cord ¹	100.689
	1	User manual Phyaction U	EN: 322.835 NL: 322.879 FR: 322.923 DE: 322.967

¹ This power cord has a CEE 7/7 type plug. For countries with other outlets, a different power cord with the appropriate plug is supplied.

7.6 Optional accessories

	Quantity	Description	Art. no.
	1	US head, multi-frequency, 1/3 MHz - ERA 1 cm ² , incl. holder	323.595
	1	Contact gel, can 5 l	100.019
	1	Pump for can, 5 l	100.020
	1	Combination therapy cord	112.823

Article numbers can change in the course of time. Check the article numbers in the most recent catalogue or ask your dealer.

The drawings are merely indicative, no rights can be derived from them.

8 APPENDICES

8.1 US head placements


Select the therapy via indication list to get information about the placement.
See §4.3.2.

8.1.1 Ultrasound therapy


Select the **US head placement** parameter to show the optimal location for the placement of the US head.

You can select the numbers in the illustration with the blue keys for more information.


- 1 Gives information on the precise anatomic location.
- 2 Numbers with a black background gives specific recommendations.



Relevant bone structures are shown for detailed information on the treated area. The number of points below the US head gives an indication of the dimensions of the treated area. The information in the illustration recommends a treatment technique. This illustration shown an example of the dynamic technique.



If other areas are possible for the US head placement a black area is shown. Select the corresponding number 2 for information on the screen. If the area is on the rear a transparent area is shown.



8.2 EMC directive

Use only US heads that are specified in this manual. See §7. The use of other accessories can have a negative effect on the electromagnetic compatibility of the equipment.


If you use the Phyaction U in the vicinity of other equipment, you must check that the Phyaction U is functioning normally.

The following paragraphs contain information about the EMC properties of the equipment.

8.2.1 Guidance and declarations

Guidance and manufacturer's declaration - electromagnetic emissions		
The Phyaction-series devices are intended for use in the electromagnetic environment specified below. The customer or the user of a Phyaction-series device should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Phyaction-series devices use RF energy only for their internal function. Therefore, their RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
	Class B	
Harmonic emissions IEC 61000-3-3	Class B	The Phyaction-series devices are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic immunity			
The Phyaction-series devices are intended for use in the electromagnetic environment specified below. The customer or the user of a Phyaction-series device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact / ±8 kV air No loss of performance	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity must be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV power / ±1 kV I/O No loss of performance	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV diff. / ±2 kV comm. No loss of performance	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	U_T - 100% (0,5 period) No loss of performance U_T - 60% (5 periods) No loss of performance U_T - 30% (25 periods) No loss of performance U_T - 100% (5 seconds) Device resets to a safe state. (60601-1 § 49.2)	Mains power quality should be that of a typical commercial or hospital environment. If the user of a Phyaction-series device requires continued operation during power mains interruptions, it is recommended that the Phyaction-series device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Not applicable	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level			

Guidance and manufacturer's declaration - electromagnetic immunity			
The Phyaction-series devices are intended for use in the electromagnetic environment specified below. The customer or the user of a Phyaction-series device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V _{rms} AM 1 kHz 80% 150 kHz to 80 MHz	10 V.....0,15-80 Mhz 51 V.....6,78 Mhz 54 V.....13,56 Mhz 50 V.....27,12 Mhz 45 V.....40,68 Mhz	Portable and mobile RF communications equipment should be used no closer to any part of a Phyaction-series device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 0,35\sqrt{P}$ $d = 0,07\sqrt{P}$ $d = 0,06\sqrt{P}$ $d = 0,07\sqrt{P}$ $d = 0,08\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m AM 1 kHz 80% 80 MHz to 2,5 GHz	10 V/m...0,08-1,0 Ghz 26 V/m....1,4-2,0 Ghz 30 V/m...433,92 Mhz 30 V/m.....915 Mhz	$d = 0,35\sqrt{P}$ 80 MHz to 800 MHz $d = 0,70\sqrt{P}$ 800 MHz to 2,5 GHz $d = 0,12\sqrt{P}$ $d = 0,23\sqrt{P}$
Radiated RF ENV 50204	3 V/m CW 200 Hz d.c. 50% 895 MHz to 905 MHz	30 V/m.895-905 Mhz	$d = 0,23\sqrt{P}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 The guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and and mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey can be considered. If the measured field strength in the location in which a Phyaction-series device is used exceeds the applicable RF compliance level above, the Phyaction-series devices should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Phyaction-series device.			
b Over the frequency range 150 kHz to 80 MHz, field strengths must be less than 10 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the Phyaction-series device

The Phyaction-series device is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of a Phyaction-series device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Phyaction-series devices as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 0,35\sqrt{P}$	80 MHz to 800 MHz $d = 0,35\sqrt{P}$	800 MHz to 2,5 GHz $d = 0,70\sqrt{P}$
0,01	0,04	0,04	0,07
0,1	0,11	0,11	0,22
1	0,35	0,35	0,70
10	1,11	1,11	2,21
100	3,50	3,50	7,00

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

8.3 Technical safety inspection

Phyaction U with serial number is / is not ¹ in good working order		
Location:	Inspection performed by: Name:	Owner: Name:
Date:	Initials:	Initials:

¹ Cross out what does not apply.

If a specific test does not apply to this equipment, place a mark in the NA (not applicable) column.

8.3.1 Test 1: General

	Yes	No	NA
1. The results of earlier safety inspections are available.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
2. The logbook is present.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
3. The type plate and the supplier's label are legible.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
4. The housing, adjusting knobs, keys and display are undamaged.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
5. The power connection and power cord are undamaged.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
6. The output connectors are undamaged.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
7. The cables and connectors of the US head(s) are undamaged.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
8. The US head(s) do not display any cracks or other damage that can endanger the insulation.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
9. The automatic self-test at switch-on does not give an error message.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
10. The display does not show any defective points or lines.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

8.3.2 Test 2: Ultrasound

	Yes	No
1. Connect the treatment head and place it in an ultrasound measurement device. Select an ultrasound therapy.		
2. Select 1 MHz, continuous (duty cycle 100%), 2 W/cm ² The measured value is within $\pm 20\%$ of the Ppk value in the channel window.	<input type="checkbox"/>	<input type="checkbox"/>
3. Select 1 MHz, duty cycle 50%, 3 W/cm ² The measured value is within $\pm 20\%$ of half the Ppk value in the channel window.	<input type="checkbox"/>	<input type="checkbox"/>
4. Select 3 MHz, continuous (duty cycle 100%), 2 W/cm ² The measured value is within $\pm 20\%$ of the Ppk value in the channel window.	<input type="checkbox"/>	<input type="checkbox"/>
5. Select 3 MHz, duty cycle 50%, 3 W/cm ² The measured value is $\pm 20\%$ of half the Ppk value in the channel window.	<input type="checkbox"/>	<input type="checkbox"/>
6. Select 3 MHz, duty cycle 50%, 0.5 W/cm ² With a dry treatment surface, the Ppk value becomes 0.	<input type="checkbox"/>	<input type="checkbox"/>
7. Select 1 MHz, duty cycle 50%, 0.5 W/cm ² With a dry treatment surface, the Ppk value becomes 0.	<input type="checkbox"/>	<input type="checkbox"/>

The maximum power transfer takes place at the operating frequencies. If the equipment does not function at the correct frequency, this results in a too low output power. It is therefore not necessary to check the operating frequencies.

8.3.3 Test 3: Electrical safety test (VDE 0751)

	Yes	No
1. The resistance of the safety earth is less than 0.2 Ω	<input type="checkbox"/>	<input type="checkbox"/>
2. The housing leakage current is less than 1000 μA	<input type="checkbox"/>	<input type="checkbox"/>
3. The patient leakage current is less than 5000 μA	<input type="checkbox"/>	<input type="checkbox"/>

Notes:

8.4 Disposal


Take account of the following environmental aspects when disposing of the equipment and the accessories:

- The basic device and the cables fall under small chemical waste (or electronic waste). These components contain lead, tin, copper, iron, various other metals and various plastics, etc. Consult the applicable national regulations.
- Gels contain only organic material and do not require any special processing.
- Packaging materials and manuals can be recycled. Deliver them to the appropriate collection points or include them with the normal household waste. This depends on the local organisation of the waste processing.

9 REFERENCE

9.1 Function overview

9.1.1 Therapy menu

Press .

The numbers refer to the program numbers.

Ultrasound therapy

Ultrasound therapy 31

9.1.2 System settings

Press .

Contrast

Language

Sound settings

Stand-by time

Text start up screen

First screen

System information

Error history

Counter working hours

Reset menu

Stop time if bad US

9.1.3 Objectives

Press  and select **Objectives**.

The numbers refer to the program numbers.

Ultrasound therapy

Improve trophic condition

Tendinitis

Stage 3 or 4 (subacute) 63

Stage 1 or 2 (chronic) 62

Ligament lesions

Subacute 64

Chronic 144

Muscle lesions

Subacute 64

Chronic 144

Osteo-chondral lesions 144

Neurogenic lesions 64

Increase extensibility

Superficial contractures 65

Partial joint contractures 145

Improve cell function

Acute joint lesions 66

Acute muscular lesions 66

Acute neurogenic lesion 66

Fracture healing 67

Phonophoresis

Phonophoresis 66

9.1.4 Indication list

Press  and select Indication list.

US: Ultrasound therapy

The numbers refer to the program numbers.

Arthrosis, US	Posttraum. diseases, US
Subacute.....64	Acute66
Chronic144	Subacute64
Bechterew, US.....62	Scar tissue, US
Bursitis, US62	Acute66
Contractures, US	Subacute65
Superficial65	Sprain, US
Deep62	Acute66
Decubitus, US.....88	Subacute64
Dupuytren, US65	Tendinitis, US
Epicondylitis, US	Subacute63
Subacute.....63	Chronic.....62
Chronic62	Ulcus Cruris, US88
Fractures, US.....67	
Frozen shoulder, US.....145	
Myalgia, US.....144	
Neuropathy, US66	

9.1.5 Diagnostics

Stress fracture search.....112

9.1.6 Contra indication

Ultrasound therapy

General

High fever
Severe cardiovascular problems
Psychological problems
Cancer with tumor metastasis
Generalised tuberculosis

Specific relative for continuous ultrasound

Infections
Acute inflammations
Thrombosis, thrombophlebitis
Varices
Increased risk to haemorrhage
Pacemaker
Epiphyseal disc (children)
Decreased sensibility
Menses
Cement of endoprosthesis
Diabetes mellitus

**Specific relative for pulsing
ultrasound**

Pacemaker

Pregnancy

9.2 Literature

A literature list can be sent on request. Please contact GymnaUniphy.

9.3 Terminology

trophic: The state of nourishment.

10 INDEX

A

Abbreviations 4
Accessories 41
 optional ultrasound therapy 42
 standard 41

C

Cell function
 improve 51
Change therapy program 30
Cleaning 34
Combination therapy 24
Combined therapy
 perform 24
Connection 11
 combination therapy 25
Contra indication 52
 selection 19
Contrast 11, 28
Counter working hours 29

D

Diagnostic program
 selection 19
Diagnostics 52
Display 15
 symbols 16
Disposal 50
Duty cycle 22

E

Electrical safety 8
EMC 8
EMC directive 44
Environmental conditions 40
ERA 23
Erase
 list 27
 number 27
Error history 29

Extensibility
 increase 51

F

First screen 28
Function overview 51
Functional test 11

G

Guarantee 36
Guide menu
 selection 18

I

Indication list 18, 52
Inspections 33
Installation 11
Intensity screen 20
Interruption 20

L

Language 11, 28
Liability 9

M

Maintenance 34
Malfunctions 35
Medical Devices Directive 9
Memory 26

O

Objectives 18, 51
Onset 22

P

Parameter symbols
 ultrasound therapy 16
parameters
 sound settings 29
Phonophoresis 51
Placing 11

Ppk 22
Prevention of explosion 8
Program
 erase 27
 number selection 18
 saving 26
 selecting 26
Purpose 7

R
Recall
 list 26
 number 27
Replacing a fuse 35
Reselling 12
Reset menu 29

S
Safety 7
 instructions 7
 technical inspection 33, 48
Selection
 therapy 17
Service 36
Sound
 Set 29
Sound settings 28
Stand-by time 11, 28
Stop 20
Stop treatment 20
Storage 12
 conditions 41
System information 28
System setting
 parameters 28
System settings 28
 changing 28

T
Technical information 39
Technical life time 37
Terminology 53
Text for start up screen 28
Text start up screen
 set 30
Therapy
 program selection 18
 selection via indication list 18
 selection via objectives 18
 set 20
 Start 20
Therapy menu 51
 selection 17
Transport 12
 conditions 41
Treatment
 interruption 20
 stop 20
Treatment time 22
Trophic 53
 improve condition 51

U
Ultrasound therapy 21
 optional accessories 42
 parameters 22
 Perform 21
 Read-out values 22
 safety 8
 technical information 39
US frequency 23
US head
 cleaning 34
 contact 22
 Indicator light 23
 placement 23
 test 33




Uniphy is a division from GymnaUniphy N.V.

Pasweg 6A
B-3740 Bilzen

Tel.: (+32) (0) 89/510.510
Fax: (+32) (0) 89/510.511

www.gymna-uniphy.com
E-mail: info@gymna-uniphy.com

Your dealer:



CE 0344