

ProTime®
Microcoagulation System



For *in vitro* Diagnostic Use

Operator's Manual

TABLE OF CONTENTS

INTENDED USE.....	2
SUMMARY AND EXPLANATION.....	2
SYSTEM FEATURES.....	3
IMPORTANT LABELS AND SYMBOLS	4
PRINCIPLES OF OPERATION	5
SAFETY FEATURES AND QUALITY CONTROL	5
INSTRUMENT SPECIFICATIONS	7
SERVICE AND MAINTENANCE	7
PREPARING THE INSTRUMENT	9
TEST PROCEDURE	10
TROUBLESHOOTING	14
MAIN MENU OPTIONS	16
SET UP	22
FOR PROFESSIONALS ONLY - PROGRAM MODE.....	27
PERFORMANCE CHARACTERISTICS	31
TECHNICAL ASSISTANCE.....	34
SUGGESTED READING.....	35
SAFETY STANDARDS	36
INDEX	39

This manual is published by International Technidyne Corporation (ITC) for use with the ProTime® Microcoagulation System and Tenderlett® Plus fingerstick blood collection device.

Questions or comments regarding the contents of this manual can be directed to the address at the back of this manual or to your ITC representative.

Please read the instructions before use.

© 2000, 2001, 2002, 2003, 2004, 2005. This document is the copyright of ITC and must not be copied or reproduced in any form without prior consent. ITC reserves the right to make technical improvements to equipment and documentation without prior notice as part of a continuous program of product development.

INTENDED USE

The ProTime Microcoagulation System is a portable, battery-operated instrument with a disposable cuvette for quantitative determination of prothrombin time from fingerstick whole blood or anticoagulant-free venous whole blood. The product is intended for use in the management of patients treated with oral anticoagulants by a healthcare professional and for patient self-testing. Additional information for medical professionals to monitor patients is located at the end of this instruction manual. Materials are available through ITC for professionals to train patient self-testers using the ProTime Microcoagulation System. ProTime instruments intended for patient self-testing are available in the U.S. by prescription only. These instruments include patient-specific product instructions.

For *in vitro* Diagnostic Use.

SUMMARY AND EXPLANATION

What Does The ProTime Microcoagulation System Do?

The ProTime Microcoagulation System is designed for testing **prothrombin time (PT)** and **International Normalized Ratio (INR)**. The reagents for a prothrombin time test are in the cuvette. Whole blood clotting time is converted to INR, then the result is calculated for plasma equivalent PT seconds. This test is done to check the status of patients receiving oral anticoagulation therapy.

What is INR?

The International Normalized Ratio (INR) was developed to help the doctor compare an individual's prothrombin time results from one lab to another. An advantage of reporting an INR is to allow for normalization of comparisons from one lab or instrument to another. The precision of the INR is improved when a reagent with a lower ISI is used.

***Note:** ISI stands for International Sensitivity Index. ProTime uses this number to calculate PT seconds from INR.*


Blood Coagulation Test Methodology



Traditional coagulation tests measure the time required for the formation of a fibrin clot following the addition of a coagulation-activating reagent. Laboratory assays typically use plasma recovered from anticoagulated (citrate) whole blood samples. The clotting time is a measure of the functionality of the patient's hemostatic system. Specific coagulation activating reagents are employed in various clotting tests to measure specific components of the hemostasis system. Clotting times are prolonged in the case of either decreased procoagulant activity or increased anticoagulant activity.

Summary of the Blood Coagulation System

The events leading to the formation of a fibrin clot are simplified in coagulation theory into two coagulation pathways: the intrinsic and the extrinsic pathways. There are twelve clotting factors or proteins involved in this cascade scheme, numbered I through XIII (excluding VI). The prothrombin time test measures the extrinsic coagulation pathway and is sensitive to coagulation Factors VII, X, V, II and Fibrinogen (I). With the exception of Factor V, vitamin K is a required co-factor for biosynthesis of these factors in the liver. The prothrombin time (PT) test uses thromboplastin as the active reagent to initiate the extrinsic pathway. The prothrombin time test will be prolonged in patients with liver disease or vitamin K deficiency. The test is widely used to monitor oral anticoagulant therapy that suppresses the synthesis of vitamin K-dependent clotting factors.
















SYSTEM FEATURES



1. **Tenderlett Plus/Tenderlett Plus LV.** The finger incision and blood collection device used with ProTime and ProTime3 cuvettes.
2. **ProTime/ProTime3 Cuvette.** The cuvette performs the PT test.
3. **Display Panel.** The display panel prompts you through the test procedure, displays results and error messages.
4. **MENU/SCROLL Button.** The  button is used to move through menu screens and to scroll through results.
5. **START/SELECT Button.** The  button is pressed to turn ProTime on and off, to start the test and to select menu items.
6. **Cuvette Entry Port.** The cuvette entry port accepts only ProTime and ProTime3 cuvettes.
7. **AC/DC Power Module.** The AC/DC power module is used to connect the ProTime instrument to the AC power cord to recharge the battery.
8. **Input DC Port.** Plug the DC power cord from the AC/DC power module into this port when connecting the ProTime instrument to the AC/DC power module.
9. **Output Data Port.** Use this port to download records in the ProTime memory.

IMPORTANT LABELS AND SYMBOLS

Before using the ProTime Microcoagulation System, it is essential that the contents of this *Operator's Manual* and any instructions accompanying the ProTime cuvettes and Tenderlett Plus incision devices are read and understood by the operator. These materials make reference to various symbols that are explained below:

	Start/Select
	Menu/Scroll
	Expiration Date of Cuvettes
	Serial Number of Device
	Lot Number of Cuvettes/Tenderlett Plus
	ITC Catalogue Number of Device
	Do Not Reuse – Single Use Only
	Upper and Lower Temperature Limitations (For Storage or Use)
	For <i>in vitro</i> Diagnostic Use
	Attention - Read Accompanying Documentation or Instructions
	Consult Instructions for Use
	Class II Protection Against Electrical Shock
 15.0V $\overline{\text{---}}$ / 2.0A	Input Port for DC Power Cord from AC/DC Power Module - Polarity, VDC and A Input
	Output Port for Data Transfer
	Name and Address of Manufacturer

PRINCIPLES OF OPERATION

The ProTime Microcoagulation System measures the PT using fibrin clot formation and detection. The ProTime cuvette is a self-contained, micro volume reaction cell constructed of precision-molded plastic.

There are two user options within the ProTime Microcoagulation System: the standard ProTime cuvette and the ProTime3 cuvette. These cuvettes differ from each other in the amount of blood that needs to be collected and tested.

The standard ProTime cuvette has five micro-channels, which contain the dried reagents required to perform triplicate testing of the PT assay and two levels of controls. The ProTime3 cuvette has three functional micro-channels. Two micro-channels perform the controls, and one micro-channel performs the PT test. The standard ProTime uses the Tenderlett Plus device for performing the fingerstick, and it is designed to hold 65 μL of blood (approximately 3 drops) needed to fill all five micro-channels. The ProTime3 uses the Tenderlett Plus LV (low volume) device for performing the fingerstick, and it collects 27 μL of blood (approximately 1 large drop) needed to fill the three micro-channels of the ProTime3 cuvette.

The instrument draws the precise volume of blood into the micro-channels of each cuvette, which contain thromboplastin and other reagents. An array of LEDs detects the motion of sample/reagent mixtures as they move through a precision restriction in each channel. The blood is pumped back and forth until a clot forms, obstructing the channel and slowing the flow of blood. The instrument detects the clot when the blood movement decreases below a predetermined rate.

SAFETY FEATURES AND QUALITY CONTROL

Each ProTime and ProTime3 cuvette has two channels for performing the two levels of control which work simultaneous with the channel(s) used for testing the PT assay. The built-in safety features of the instrument and integral reagent controls work together to ensure that the instrument and reagent systems are working properly and that the test procedure was performed correctly. Two levels of quality control are performed with each and every test.

Calibration

The ProTime instrument and cuvettes are pre-calibrated. No additional calibration is required.

Reagents

ProTime cuvettes are pre-loaded with dried thromboplastin, stabilizers and buffers. The thromboplastin has a high sensitivity, measured as an ISI near 1.0. Each cuvette performs the prothrombin time test and, in addition, has one channel for a Level I control and one channel for a Level II control. The controls consist of purified plasma-extracted coagulation factors and anticoagulants.

Operating Precautions

- For *in vitro* Diagnostic Use.
- The ProTime instrument is designed for use only with ProTime cuvettes.
- The ProTime instrument will not perform a test if the cuvette is past its expiration date or has been improperly stored.
- The ProTime instrument is designed to be used for testing in a stationary position. DO NOT perform testing while carrying or holding the instrument.
- In order to charge the ProTime instrument, the AC power cord should be plugged into an electrical service outlet and the AC/DC power module while the DC power cord from the AC/DC power module is plugged into the DC port in the back of the ProTime instrument.
- DO NOT expose the ProTime instrument to extreme temperature (above 35°C, 95°F). Such exposure could affect the performance of any type of electronic instrumentation.
- DO NOT drop the ProTime instrument, and do not use the results if the instrument is dropped during a test.
- DO NOT attempt to open the ProTime instrument other than for battery replacement, as there are no user-serviceable parts.
- DO NOT remove the AC/DC power module from the ProTime instrument by pulling on the cord.

Patient specimens and used cuvettes are potentially infectious. The cuvettes include materials that have been prepared from human plasma or serum that has been tested using US FDA recognized methods and found to be non-reactive for HIV antibody and for hepatitis B surface antigen. Handle with appropriate care and dispose of cuvettes and blood collection materials in accordance with standard methods of biohazard control.

The use of accessory equipment (e.g., printers, etc.) not identified in this manual either in the patient vicinity, or that does not comply with either the equivalent safety requirements of this equipment or UL/IEC 60601-1 or IEC 60601-1-2, may lead to a reduced level of safety with the resulting system.

Limitations

- The ProTime instrument uses only fresh capillary or venous whole blood. Plasma or serum cannot be used. Glass tubes or syringes must not be used to collect venous samples. Use only plastic syringes without anticoagulants to collect venous samples. Poor fingerstick blood collection technique may affect results.
- In clinical trials, no significant difference was observed between fingerstick and venous specimens run on ProTime. During those trials, ProTime measured patients with an INR range of 0.8 to 7.0. ProTime is programmed to calculate and report INR results as follows:

If the calculated INR is:	Then ProTime displays:
< 0.6	INR LOW
0.6 - 0.79	INR < 0.8
0.8 - 7.0	Result
7.1 - 9.9	Result - followed by "*"
10.0 - 12.0	INR > 9.9
> 12.0	INR HIGH

- Results may be affected in patients receiving heparin or who have an abnormal response to heparin.
- Correlation of results reported by the ProTime instrument to laboratory results depends on the precision of the laboratory method and on the ISI value of the laboratory reagent and instrument system.
- Do not disturb instrument while test is in progress.

As with all diagnostic tests, ProTime Microcoagulation System test results should be scrutinized in light of a specific patient's condition and anticoagulant therapy. Any results exhibiting inconsistency with the patient's clinical status should be repeated or supplemented with additional test data.

Follow doctor's instructions if you have difficulty performing the test or receive a result outside of the prescribed therapeutic range.

Reagent Preparation and Storage

- Bring the foil-pouched cuvette to room temperature prior to use.
- ProTime cuvettes are ready-to-use. No additional preparation is required.
- Store the foil-pouched cuvettes in a refrigerator (2°C to 8°C, 36°F to 46°F).
- An unopened cuvette pouch is stable when stored at 2°C to 8°C until the date printed on the pouch. Unopened cuvette pouches may be stored at room temperature for 60 days. Once the pouch has been opened, the cuvette must be used within 16 hours.

Cuvette or Tenderlett Disposal

The ProTime cuvette and Tenderlett Plus are for single use only and are not to be reused. After use they contain human blood and should be disposed of in accordance with local regulations for human blood contaminated waste.

INSTRUMENT SPECIFICATIONS

Size	2.77 (h) x 4.5 (w) x 8.75 (l) inches
Weight	1.6 pounds
Operating Temperature	Room temperature (15°C to 30°C, 59°F to 86°F)
PT Test Temperature	37 ± 1.0°C
Battery Type	Nickel Metal Hydride (NiMH)
Operating Time On Battery	Approximately 2 hours (constant run), or 10 complete test cycles per charge. Tests may also be run while ProTime is plugged into AC/DC power module.
Anticipated Battery Life	500 charges
100 - 240 VAC Power Supply/Charger	Input: 100 – 240 volts AC, 47 - 63 Hz Output: 15 volts DC, 2.0 A

SERVICE AND MAINTENANCE

Routine Maintenance and Cleaning

DO NOT immerse the ProTime instrument or allow fluid to enter the cuvette housing. Inspect and clean the outside of the cuvette slot as required. Remove residual dried blood or other foreign matter on the outside of the instrument using gauze dampened with a 10% dilution of household bleach in water or with gauze dampened with isopropyl alcohol or other disinfectant. DO NOT use other solvents or strong cleaning solutions as they may damage the plastic components of the instrument.

Servicing The ProTime Instrument

Other than replacement of the rechargeable battery as described in this manual, the ProTime instrument is not user serviceable. Should service be required, please contact Technical Support at 1-732-548-5700, 1-800-631-5945, or e-mail us at techsupport@itcmed.com. If you are advised to return the instrument to ITC for service or repair, prior to shipping please clean the instrument, as described above.

Instrument Disposal

If instrument disposal is required, follow local regulations for the disposal of electronic devices. For battery disposal, see the **Battery Replacement** section.

Battery Information

The ProTime Microcoagulation System is designed to run either on AC power supplied by the AC/DC power module or on the rechargeable battery supplied within the unit.

Rechargeable Battery Facts

- Batteries discharge naturally over time (approximately 5% per month).
- Battery capacity (the amount of charge the battery will hold) is lower at low temperatures.
- Batteries that are charged to their maximum capacity will generate heat if they continue to be charged.
- The ProTime instrument uses a rechargeable NiMH (Nickel Metal Hydride) type battery. The maximum capacity of any rechargeable battery will gradually decrease over time. To ensure maximum life of the rechargeable battery, read and follow the information in **Care of the ProTime Battery for Maximum Life** section.

Care of the ProTime Battery for Maximum Life

A new instrument, an instrument that is used infrequently, or an instrument with a new replacement battery, should be plugged in for at least 8 hours before use to ensure that the battery is completely charged. The instrument screen will show CHARGING BATTERY when the AC/DC power module is connected to the AC power cord and the ProTime instrument. The screen will show CHARGE COMPLETE when the battery is fully charged. The AC/DC power module should be disconnected after the CHARGE COMPLETE message is seen. The AC/DC power module that has been supplied by ITC has been selected specifically for use with your ProTime Microcoagulation System. Do not use any other AC/DC power module.

When the battery indicator on the screen shows less than 25% charge remaining, please refer to the **Charging the Battery** section for instructions. To maximize battery life, allow your ProTime instrument to discharge completely before re-charging. Avoid charging the ProTime instrument for frequent, short periods of time (such as charging for a few minutes, removing from the AC/DC power module, and then recharging again). Avoid storing and charging the ProTime instrument at extreme temperatures.

Built-in ProTime Battery Features

The ProTime instrument will shut off after 5 minutes if left unattended in order to preserve the battery charge. If the instrument shuts off automatically, press the **O** button to re-start the instrument. Before the start of each test, the ProTime instrument checks the amount of battery charge. If there is not enough charge in the battery to run a test, the ProTime screen will display PLEASE CHARGE IT. If this occurs, follow the charging instructions located in the **Charging the Battery** section. Whether or not the batteries are charged, you may continue to run tests once the ProTime instrument is plugged into its AC/DC power module.

Battery Replacement

Refer to the instructions provided with the replacement battery. The used battery should be disposed of in accordance with local regulations for NiMH batteries.

SPECIMEN COLLECTION

Fingerstick whole blood is the recommended specimen. The Tenderlett Plus device is to be used with the ProTime cuvette, and the Tenderlett Plus LV (low volume) device is to be used with the ProTime3 cuvette. The Tenderlett Plus device will collect approximately 65 μL of blood (approximately 3 drops), while the Tenderlett Plus LV device will collect approximately 27 μL of blood (approximately 1 large drop). Samples should be analyzed immediately after collection. No additional sample preparation is required.

For venous samples, collect venous whole blood into an anticoagulant-free plastic syringe in place of fingerstick sampling steps 3 and 4 of the **Test Procedure** section. Immediately dispense venous sample into the Tenderlett Plus collection cup, filling the Tenderlett Plus collection cup. Follow instructions from step 5 of the **Test Procedure** section.

Note: Serum, plasma or whole blood collected with any anticoagulants is NOT suitable samples.


PREPARING THE INSTRUMENT

Unpacking and Inspection

1. Remove any protective packaging that may be present around the instrument.
2. Examine the packaging material to be sure that the AC/DC power module, AC power cord (see note below), connecting cables, or other components have been removed. The materials that are provided are listed below.

Note: Inspect each component for damage when unpacking. If damage is observed, contact your ProTime representative.

Materials Provided

Article	Quantity
ProTime Microcoagulation Instrument	1
 ProTime Microcoagulation System Information and Training CD (PROTIMEPRO only)	1
AC/DC power module – ITC Part No. IR5226	1
ProTime Microcoagulation Operator's Manual	1

Note: An AC power cord is provided in the United States and Canada only. For use outside the US and Canada, the customer must obtain a power supply cord that is internationally harmonized and marked "<HAR>", 2-conductor, 0.75 mm² minimum wire, rated 300 V, with a PVC insulated jacket. The cord and plug cap must be suitable for medical use. The cord must have a molded on plug cap rated 250 V, 2.5 A.

Materials Needed But Not Supplied

- ProTime cuvettes
- Tenderlett Plus incision device

Optional Materials

Article
Personal Computer Interface Cable – ITC No. PROCABLE (ITC Part No. IR5313X)
Printer/Label Maker – ITC No. LBLKIT
Replacement Printer/Label Maker Interface Cable – ITC No. LBLCABLE (ITC Part No. IR5314X)
Replacement Battery – ITC No. PROBATTERY

Charging the Battery

The battery in the instrument needs to be charged before the first use.

1. Connect the AC/DC power module to an electrical service outlet, using the AC power cord.
2. Connect the DC power cord from the AC/DC power module into the DC port on the back of the instrument. The instrument screen will show CHARGING BATTERY when the AC/DC power module is connected to the AC power cord and the ProTime instrument. The screen will show CHARGE COMPLETE when the battery is fully charged.
3. Allow the battery to charge for at least eight hours.

TEST PROCEDURE

1. Turn on the ProTime Instrument

Press the **0** button to start. The ProTime instrument performs a self-check procedure that may take up to 60 seconds. ProTime will prompt you through the test. Watch the screen and follow the prompts.




2. Insert a Cuvette

Make sure the ProTime cuvette is brought to room temperature before use. Wait for the prompt. Insert the cuvette into the slot with the printing face up and the bar-code down.



The WARMING screen will appear as follows:



3. **Prepare for Finger Incision**

While the cuvette is warming, prepare the finger. Wait for the prompt before incising the finger and collecting blood.


It is easier to collect blood if the hands are warm. Follow these steps to ensure a good sample:

- Wash the hands in warm water or rub hands together to stimulate blood flow.
- Apply firm pressure to the palm and finger. Massage the hand and push blood into the fingertips.
- Cleanse middle or ring finger and dry. To prevent contamination, do not touch the site after cleansing.



4. **Blood Collection**

When this screen appears, it is time for the finger incision.




CAUTION: Blood collection must be finished within 2:10 minutes to prevent early clotting of the sample. ProTime will keep time. If the 2:10 minute interval has expired and the 0 button has not been pressed during this interval, a TIME OUT error message is displayed. To run another test, repeat procedure from Step 2.

- Place the Tenderlett Plus device firmly against the side of the finger. Place thumb on top of the device as shown, and press the red trigger firmly using the other thumb.




- Wipe away the first trace of blood.
- Gently massage from the base of the finger to force blood to the tip so that a large drop of blood forms.

- Touch the large drop of blood to the collection cup. Keep adding blood until the blood level fills the cup above the line.



- For Tenderlett Plus LV, ensure the cup is filled completely. Ensure the blood extends all the way to the bottom of the cup. **Add another drop if you are not sure you have enough.**



5. Snap Tenderlett Plus to ProTime

- Hold the device at an angle and place the front end of the device into the slot in the instrument.
- Press down to click the Tenderlett Plus in place. You should hear a soft click.



***Note:** Proper engagement of the Tenderlett Plus to the cuvette is critical to prevent a sample error.*

6. Start the Test

- Press the **0** button to start the test. This signals ProTime to draw the sample into the cuvette.





- It takes only a few seconds for ProTime to draw the blood into the cuvette. Watch the screen for the next prompt.



7. Remove Tenderlett Plus

Remove Tenderlett Plus immediately when prompted to do so.

CAUTION: Failure to do so will result in an error message. ProTime allows you six seconds.




- Do not press the **0** button after the Tenderlett Plus is removed from the ProTime instrument while you are testing your INR. This will interrupt the test procedure, and you will have to start over with a new blood sample.
- The instrument then progresses to the test and displays the TESTING screen.



8. Read the Result


In a few minutes, the result is ready.



- Press the **0** button to turn off.
- Press the **◀** button to go to the MAIN MENU if you want to run another test, review the data in memory, print results, transfer results to a computer, or perform set up functions.

Notes:

- The result remains displayed for 5 minutes or until the **0** button or the **◀** button is pressed.
- If the cuvette has not been removed and the instrument is left unattended for 5 minutes, the instrument will display the following messages before shutting down and powering OFF:



What Does the Result Mean?

The result indicates the clotting activity of blood. When ProTime is used as a self-testing instrument, the healthcare provider may program ProTime with the upper and lower limits that are right for the patient. In this case, ProTime will display OTR if results are outside of the therapeutic range. The OTR will not display if results are within the limits. If no limits are set, ProTime will display only the result.

Note: If OTR (Out of Therapeutic Range) appears after the INR result, the result is out of the therapeutic range (TR) limits that have been preset by the professional (see the **For Professionals Only – Program Mode** section).

TROUBLESHOOTING

Most often, an error message indicates a problem with blood collection or a mistake in the test procedure. If an error message appears, read the instructions again and repeat the test with a new cuvette. The **Troubleshooting Guide** section gives the messages for the most common errors with possible causes and solutions.

As with all diagnostic tests, the ProTime Microcoagulation System test results should be scrutinized in light of a specific patient's condition and anticoagulant therapy. Any results exhibiting inconsistency with the patient's clinical status should be repeated or supplemented with additional test data.

If you get a persistent error, write down the error message and call ITC Technical Support at 1-732-548-5700, 1-800-631-5945, or e-mail us at techsupport@itcmed.com.

Please have the ProTime serial number and cuvette lot number ready when you contact ITC Technical Support.

Troubleshooting Guidelines

- Become skilled at obtaining blood with the fingerstick technique. This will help avoid sample errors.
- Follow all of the directions on the screen for optimal operation.
- Do not try to hasten the test process by performing the fingerstick while the WARMING screen is displayed, as blood may clot before it is drawn into the cuvette and an error may occur. FOLLOW THE INSTRUCTIONS ON THE SCREEN.


Troubleshooting Guide

Screen Display	Source of Error	Method of Control
INR HIGH REPEAT TEST IF SAME, CALL MD	Patient's PT result is too high (INR > 12.0)	Repeat the test. If it occurs again, patient should be referred to their physician IMMEDIATELY and re-tested at the lab.
INR LOW REPEAT TEST IF SAME, CALL MD	Patient's PT result is too low (INR < 0.6)	Repeat the test. If it occurs again, patient should be referred to their physician IMMEDIATELY and re-tested at the lab.
ON-BOARD QC OUT OF RANGE TRY AGAIN	Level I or II control is too high or too low	Repeat the test. Check for adequate storage of cuvettes or sampling technique. Confirm results with the lab.
NO CLOT DETECTED REPEAT TEST IF SAME CALL MD	At least one channel did not clot	Repeat the test. If it occurs again, patient should be referred to their physician IMMEDIATELY and re-tested at the lab.
TIME OUT TRY AGAIN	Time ran out before starting test	Test again with new Tenderlett Plus and a fresh fingerstick. You may reuse the cuvette if no blood was in contact with the cuvette.
TEST INTERRUPTED TRY AGAIN	Cuvette removed or operator interrupted test	Review correct procedure. Repeat test.
SAMPLE NOT SEEN TURN OFF TRY AGAIN	Can be caused by small samples, clots, air bubbles, or a seal problem in one or more channels	Check for proper collection technique. Verify cup is completely filled. Use another cuvette.
SAMPLE TOO LARGE TURN OFF TRY AGAIN	Sample oversize	Review directions and repeat the test.
SAMPLE TOO SMALL TURN OFF TRY AGAIN	Sample undersize	Repeat test. Exceed the fill line on the cup to ensure sample size.
SAMPLE ERROR TURN OFF TRY AGAIN	Incorrect fluid movement Air bubble detected in one or more of the channels	Turn off and on again and repeat test with new cuvette and fresh fingerstick.



Screen Display	Source of Error	Method of Control
CUVETTE EXPIRED	Expired cuvette	Check to make sure cuvettes have not expired. The expiration date is located on the packaging and alongside the barcode on the cuvette.
BARCODE ERROR REMOVE CUVETTE TRY AGAIN	Instrument cannot read barcode accurately	Visually inspect barcode. If scratched, discard cuvette. If dirty, wipe clean. If barcode is good, review correct procedure and repeat test.
BATTERY ERROR TURN OFF TRY AGAIN	Power supply battery error	Repeat the test. If it occurs again, replace the battery. If problem persists, call ITC Technical Support.
INSTRUMENT ERROR TURN OFF TRY AGAIN	Instrument set up, data log or communication error	Repeat the test. If problem persists, call ITC Technical Support.
CHARGE ERROR TURN OFF TRY AGAIN	Power supply error	Check the AC/DC power module, repeat the test. If problem persists, call ITC Technical Support.
TEMP ERROR TURN OFF TRY AGAIN	Temperature not in range	Check for proper operating temperature. Repeat the test. If problem persists, call ITC Technical Support.
PHOTO ERROR TURN OFF TRY AGAIN	LED blocked or other photo system error	Repeat the test. If problem persists, call ITC Technical Support.

MAIN MENU OPTIONS

The options presented in the MAIN MENU are:



Each of these options will lead to sub menus within the selection. The following buttons are used to navigate the menu:

- The  button is used to move the highlight bar to select the option.
- The  button is used to select the option that is highlighted.

Run Test

To run the test, select the RUN TEST menu item. The instrument does a SELF CHECK procedure, which may take up to 60 seconds. The following screen is displayed for this period of time:



The test sequence continues as described in the **Test Procedure** section.

When the test completes, the result will appear as one of the following screens, depending on PID/OID selections (see SET UP below).

Note: "X" is used for illustrative purposes only in the following examples.

If PID is ON and OID is ON:




If PID is OFF and OID is ON:



If PID is ON and OID is OFF:



If PID is OFF and OID is OFF:



RUN LQC

When RUN LQC is selected from the MAIN MENU, the following screen is displayed and the user can select NOR/ABN (NORMAL/ABNORMAL).

```
LQC NOR
LQC ABN
MAIN MENU

▽ MOVE
0 OK
```

The instrument performs SELF CHECK procedure which may take up to 60 seconds. The testing proceeds as described above except that no prompt for PID will be shown for a QC test, whether PID is ON or OFF.

When the test completes, the result will be displayed as follows, depending on OID ON/OFF and LQC NOR/ABN selection:

NOTE: The result and ID information will be stored in the database as a QC record.

For LQC Normal with OID on:

```
LIQUID QC
NORMAL
OID = XXXXXX
INR = X.X

▽ MAIN MENU
0 TURN OFF
```

For LQC Normal with OID off:

```
LIQUID QC
NORMAL

INR = X.X

▽ MAIN MENU
0 TURN OFF
```

For LQC Abnormal with OID on:

```
LIQUID QC
ABNORMAL
OID = XXXXXX
INR = X.X

▽ MAIN MENU
0 TURN OFF
```


For LQC Abnormal with OID off:

```
LIQUID QC
ABNORMAL

INR = X.X

▽ MAIN MENU
0 TURN OFF
```

SHOW RESULTS




Press the  button to move the highlight bar to the SHOW RESULTS line on the MAIN MENU.

Press the  button to view the SHOW RESULTS menu and observe the following screen:

```
DATA HISTORY
PRINT DATA
SEND DATA
MAIN MENU

 MOVE
 OK
```

DATA HISTORY

Press the  button to view the most recent result in the DATA HISTORY memory. The instrument memory holds 50 results. Pressing the  button will scroll through individual results. The results are stored in memory from the most recent to the oldest. Press the  button to return to the MAIN MENU.



PATIENT RECORDS

The result and ID information will be stored in the database as one of the following screens depending upon the PID/OID ON/OFF selections:

Note: "MM/DD/YYYY" and "HH:MM" are used for illustrative purposes only in the following examples of date and time screens.



If PID is ON and OID is OFF:

```
MM/DD/YYYY HH:MM
PID = XXXXXXXXXXXX



INR = X.X
PT = XX.X SEC
 NEXT
 MAIN MENU
```

If PID is OFF and OID is OFF:

```
MM/DD/YYYY HH:MM



INR = X.X
PT = XX.X SEC
 NEXT
 MAIN MENU
```

If PID is ON and OID is ON:

```
MM/DD/YYYY HH:MM
PID = XXXXXXXXXXXX
OID = XXXXXX
INR = X.X
PT = XX.X SEC
 NEXT
 MAIN MENU
```

If PID is OFF and OID is ON:

```
MM/DD/YYYY HH:MM

OID = XXXXXX
INR = X.X
PT = XX.X SEC
 NEXT
 MAIN MENU
```

LQC RECORDS

The result and ID information will be stored in the database as a QC record. QC result will be displayed as follows, depending upon the OID ON/OFF selections and the LQC NOR/ABN selections.

For LQC Normal with OID on:

MM/DD/YYYY HH:MM
LIQUID QC
NORMAL
OID = XXXXXX
INR = X.X
▾ NEXT
0 MAIN MENU

For LQC Normal with OID off:

MM/DD/YYYY HH:MM
LIQUID QC
NORMAL

INR = X.X
▾ NEXT
0 MAIN MENU

For LQC Abnormal with OID on:

MM/DD/YYYY HH:MM
LIQUID QC
ABNORMAL
OID = XXXXXX
INR = X.X
▾ NEXT
0 MAIN MENU

For LQC Abnormal with OID off:

MM/DD/YYYY HH:MM
LIQUID QC
ABNORMAL

INR = X.X
▾ NEXT
0 MAIN MENU



PRINT DATA OR SEND DATA

Press the ▾ button to move the highlight bar to the SHOW RESULTS line on the MAIN MENU. Then press the 0 button to view the SHOW RESULTS menu and observe the following screen:

DATA HISTORY
PRINT DATA
SEND DATA
MAIN MENU



▾ MOVE
0 OK




The instrument auto send feature allows transmitting of test results directly to a serial printer or to a computer by using the PROCABLE. Contact ITC Customer Service to order the PROCABLE.



To use the AUTO SEND feature, press the  button on the PRINT DATA or SEND DATA display until AUTO SEND is highlighted. Press the  button and, depending on the current settings in the instrument, one of the following displays will appear.



Example: Printing Data


Press the  button to set the feature to ON or OFF and then press the  button to save the setting.

- Highlight the SHOW RESULTS line in the MAIN MENU.
- Press the  button to view the SHOW RESULTS menu.
- Press the  button to move the highlight bar to the PRINT DATA or SEND DATA line.
- Press the  button to access the PRINT DATA or SEND DATA option and select THIS RESULT or ALL RESULTS option in the PRINT DATA or SEND DATA screen.






Upon selecting THIS RESULT, the last recorded test is sent from the instrument memory to a printer or to a computer. If ALL RESULTS is selected, all results (up to 50) are printed/sent from the instrument memory to a printer or to a computer. If MAIN MENU selection is made, the user is returned to the main menu.

Note: If the AUTO SEND feature is on, results will be automatically sent if the PROCABLE is connected prior to the start of a test. To use a laser printer, the data must first be transferred to an IBM-compatible personal computer and the results can be printed from that computer. The following screen displays until the PRINT DATA or SEND DATA process is complete:




SET UP


Go to the MAIN MENU screen and choose the SET UP option by using the  button to scroll and the  button to select SET UP. The following screen appears:



SET LANGUAGE

Press the  button to select SET LANGUAGE. The following options are displayed when SET LANGUAGE is selected:




Highlight a language selection and press the  button to set it.

Note: English will be highlighted as the default language in a new instrument. After changing the language, the new language becomes the default. The instrument will automatically turn off after changing the language selection.

SET TIME / DAY



Select SET TIME/DAY from the SET UP menu. The following screen is displayed:




12 HOUR - MM/DD will be highlighted as a default selection. Highlight your TIME-DAY preference, then press the  button.

Note: In the case of a 24 HOUR selection, AM/PM acronyms will not appear on the screens.


Changing the Hour

The time is pre-set to Eastern Standard Time (EST). For example, the time is 1:25 PM. Change the hour by pressing the  button until the correct hour appears in the highlight bar. In this example, the correct hour is 10 AM. To change the hour, press the  button until 10 AM appears in the highlight bar. The PM changes to AM at midnight. Once the correct hour appears, press the **0** button to set the hour. The procedure will then advance to the MINUTES screen.




```
TIME = 1:25 PM

HOUR = 1 PM

 CHANGE
0 OK
```


Changing the Minutes


The hour on the top line has changed to the time just set and PM has changed to AM.



```
TIME = 10:25 AM


MINUTE = 25


 CHANGE
0 OK
```

Change the minute by pressing the  button until the correct minute appears.

Once the correct minute appears, press the **0** button to set the minutes. The procedure will then automatically advance to the DAY screen.


Changing the Day

Only the day can be changed. The month and year can be changed only in PROGRAM MODE. Change the day by pressing the  button until the correct day appears. Once the correct day appears, press the **0** button to set the day.




```
FEB 10, 2005


DAY = 10

 CHANGE
0 OK
```

SET PID/OID

The PID/OID selection is indicated only by the location of the highlight (text on the screen is unchanged).



- Use the  button to select PID/OID options:








- Selecting PID/OID ON enables both a patient ID and operator ID to be entered.
- Selecting PID ON enables only a patient ID to be entered.
- Selecting OID ON enables only an operator ID to be entered.
- Selecting OFF disables both a patient ID and operator ID.

Note: The selection is highlighted and remains the default until changed again by the user.


Entering and Changing Numerical Information (PID/OID, PASSWORD)

By using a combination of the  button and the  button and the triangular cursor on the display screen, numeric information can be entered into the instrument for operator identification number (OID), patient identification number (PID) and PASSWORD fields. When entering numeric information into a field, pressing a button has the following actions:




-  **CHANGE** - Pressing the  button once will increment the digit in the current cursor position by one, starting at 0 and going through 9 and back to 0 again.
-  **OK** - Pressing the  button briefly enters the digit displayed in the current cursor position and then moves the cursor right, to the next field position. Pressing the  button for 2 seconds or longer saves the complete numeric field.

EXAMPLE: Entering a PID number:

With the PID feature enabled, the following screen appears after a cuvette is inserted into the instrument:



To enter a PID, the first position in the PID field (starting from the left) is indicated by a triangular cursor.

- Press the  button until the desired digit is displayed.
- Press the  button briefly to enter the digit in that position of the field. The cursor will automatically move to the next field position.
- Repeat these steps until the desired field length is entered.
- To save the entered PID value press and hold the  button for approximately 2 seconds until the second audio beep is heard.


Note: These audio beeps will be heard regardless of BEEP setting (ON/OFF).

After the cuvette is inserted, the user has the option to set an operator ID (OID) and/or a patient ID (PID) if the OID and/or PID have been enabled in the SET UP section of the MAIN MENU. The PID may contain up to twelve digits, the OID can be a value with up to six digits.

Note: In the different scenarios related to SET PID/OID ON/OFF settings, some of the following four screens are not presented to the user. For example, none of these four screens will appear if both PID and OID are set to OFF value. The default numerical PID and OID values are zero.



If OID is set to the ON position, the following screen will appear:





Use the  and  buttons to move the cursor and to enter the numeric ID.

The CONFIRM OID confirmation screen will follow the ENTER OID screen:



If  button is pressed, the user will return to the previous screen, otherwise (if the  button is pressed and PID is set to the ON position) the program will proceed to the following ENTER PID screen:



Use the  and  buttons to move the cursor and to enter the numeric ID.

The CONFIRM PID confirmation screen will follow the ENTER PID screen:



After setting of both requested PID and/or OID, the WARMING screen will appear in the center of the display:



The test sequence then continues as described in TEST PROCEDURE.


SET BEEP

If SET BEEP is selected, the following screen appears:



To turn the BEEP SOUND ON, press the  CHANGE button. The SET BEEP ON screen appears:



Then press the  button to turn the beep sound ON.

***Note:** Regardless of BEEP setting, the beep will sound under the following conditions: power on, power off, and following a key hold.*

FOR PROFESSIONALS ONLY - PROGRAM MODE

Note: The Program Mode feature is only available to professional users. The following information is not contained in the Operator's Manual for patients performing self-testing. For PROTIMEPRO, to access the Professionals Only section of the training CD, enter the following code when prompted by the software: **54321** (the training CD is not provided with PROTIMEINT).

If the PROGRAM MODE option is selected from the SET UP Menu, the following screen appears:


A monochrome screen titled "ENTER PASSWORD". Below the title is a downward-pointing triangle icon. Underneath the triangle is a line of eight "X" characters representing a password input field. At the bottom of the screen are two options: a left-pointing triangle icon next to the word "CHANGE" and a circle icon next to the word "OK".

```
ENTER PASSWORD
▼
XXXXXXXX
△ CHANGE
○ OK
```

Note: The PROGRAM MODE option is password protected. The default password in a newly manufactured instrument is 000000, which the user can then change to a password of their own choosing (see the **Change Password** section). If you forget or lose your password, an emergency password is available from ITC Technical Support by calling 1-732-548-5700, 1-800-631-5945, or e-mail us at techsupport@itcmed.com.


If the entered password matches the password stored in the instrument, the instrument will enter the PROGRAM MODE. If the entered password is incorrect, the user will be given two more chances to enter the correct password in the ENTER PASSWORD screen. After three unsuccessful attempts to enter the correct password, the program will return to the MAIN MENU and the user should contact ITC Technical Support. After performing all programming functions the user will return to the MAIN MENU.


After the correct password has been entered the PROGRAM MODE menu screen appears:

A monochrome screen displaying a list of options. The options are: "SET MM/YYYY", "SET TR LIMITS", "SET ISI", "CHANGE PASSWORD", "CLEAR MEMORY", and "MAIN MENU". At the bottom are two options: a left-pointing triangle icon next to the word "MOVE" and a circle icon next to the word "OK".

```
SET MM/YYYY
SET TR LIMITS
SET ISI
CHANGE PASSWORD
CLEAR MEMORY
MAIN MENU
△ MOVE
○ OK
```

SET MM/YYYY


Selecting SET MM/YYYY (set month and year) will guide you through the following sequence of screens: Change the month by pressing the  button until the correct month appears.

A monochrome screen showing the date "FEB 10, 2005". Below the date is the text "MONTH = FEB". At the bottom are two options: a left-pointing triangle icon next to the word "CHANGE" and a circle icon next to the word "OK".

```
FEB 10, 2005
MONTH = FEB
△ CHANGE
○ OK
```

Once the correct month appears, press the **0** button to set the month. The procedure will then advance to changing the year screen.

Change the year by pressing the **▽** button until the correct year appears.



FEB 10, 2005

YEAR = 2005


▽ CHANGE

0 OK

Once the correct year appears, press the **0** button to set the year. The procedure will then advance to MAIN MENU.

SET TR LIMITS

To set therapeutic range (TR) limits, select the SET TR LIMITS from PROGRAM MODE menu screen. The following screens appear for selection of the upper and lower TR limits.



SET TR


UPPER LIMIT = 9.9

▽ CHANGE

0 OK

Set upper limit by pressing the **▽** button until the desired upper limit is displayed. Select the limit using the **0** button. The upper limit can be set to one of eight values ranging from 2.0 to 5.5 with an incremental step of 0.5 or to the value 9.9. Once the upper limit is set, the lower limit will be available for setting.

Set lower limit by pressing the **▽** button until the desired lower limit is displayed. Select the limit using the **0** button. The lower limit can be set to 0.0, 1.2, or to one of six values ranging from 1.5 to 4.0 with an incremental step of 0.5.



SET TR

LOWER LIMIT = 0.0

▽ CHANGE

0 OK


If limits are not changed, the instrument will use its default limits 9.9 and 0.0. Once the limits are set, the instrument will return to the MAIN MENU.

SET ISI

If SET ISI is selected from the PROGRAM MODE, the following screen is displayed:



To change the ISI used in the calculation, access the PROGRAM MODE from the SET UP on MAIN MENU and enter the ISI of the local laboratory reagent

Displayed plasma equivalent values are indicative of the results obtained using a laboratory reagent with an ISI of 1.0. If the ISI of the reagent used in the local facility is significantly different from 1.0, the operator may wish to select an ISI more closely aligned with their lab reagent. Change the pre-programmed ISI by pressing the  button until the desired ISI is displayed. Press the 0 button to set. The values used for the normal control population change as the theoretical ISI is changed. The normal control population values are determined from clinical data.


REPORTABLE RANGE

In clinical trials, no significant difference was observed between fingerstick and venous specimens run on ProTime. During those trials, ProTime measured patients with an INR range of 0.8 to 7.0. ProTime is programmed to calculate and report INR results as follows:


If the calculated INR is:	Then ProTime displays:
< 0.6	INR LOW
0.6 – 0.79	INR < 0.8
0.8 – 7.0	Result
7.1 – 9.9	Result - followed by “*”
10.0 – 12.0	INR > 9.9
> 12.0	INR HIGH

CHANGE PASSWORD

If CHANGE PASSWORD option is selected from the PROGRAM MODE menu, the user must enter the password (up to 6 digits) as a left justified number without any blanks. The rules for entering a numerical PID or OID are also applicable for entering the 6 digit numerical password (see page 24).

A screenshot of the 'CHANGE PASSWORD' screen. At the top, the text 'CHANGE PASSWORD' is displayed. Below it is a downward-pointing triangle icon. Under the triangle is a row of six 'X' characters, representing a masked password. At the bottom left, there are two options: a left-pointing triangle icon next to the word 'CHANGE' and a circle icon next to the word 'OK'.

After password has been entered, the CONFIRMATION OF PASSWORD screen appears:


A screenshot of the 'CONFIRMATION OF PASSWORD' screen. At the top, the text 'CONFIRMATION OF PASSWORD' is displayed. Below it is a downward-pointing triangle icon. Under the triangle is a row of six 'X' characters, representing a masked password. At the bottom left, there are two options: a left-pointing triangle icon next to the word 'CHANGE' and a circle icon next to the word 'OK'.

If the **0** button is pressed, the newly entered password is accepted by the instrument. If the **◀** button is pressed, the user will repeat the change password procedure by returning to the previous CHANGE PASSWORD screen.

Note: Upon manufacturing, the default password is set to **000000**.

CLEAR MEMORY

Upon this selection from the PROGRAM MODE menu all results will be erased from the instrument memory. To prevent accidental loss of data, the following confirmation screen will precede the data removal from the instrument memory:

A screenshot of the 'CLEAR MEMORY ARE YOU SURE?' screen. At the top, the text 'CLEAR MEMORY' is displayed, followed by 'ARE YOU SURE?' on the next line. At the bottom left, there are two options: a left-pointing triangle icon next to the word 'NO' and a circle icon next to the word 'YES'.

Selecting NO will return the instrument to the MAIN MENU, and no results will be erased from the instrument memory.

PERFORMANCE CHARACTERISTICS

Reporting ProTime Results

The ProTime Microcoagulation System reports a result as International Normalized Ratio (INR) and in PT seconds. The ProTime system calculates the INR directly from whole blood clotting time based on a conversion equation that was established in clinical trials. The result in plasma equivalent seconds is then calculated from the INR result.

Historically, laboratories have reported prothrombin time results in seconds, specifically the time it takes for the plasma sample to clot following the addition of the reagents. Reporting the results in seconds is problematic since the clotting time depends on the sensitivity of the reagent used. Results cannot be compared from one laboratory to another. The INR system was introduced by the World Health Organization (WHO) to standardize PT reporting such that patient results could be uniform across different laboratories. INR is calculated using the formula:

$$\text{INR} = (\text{PT}/\text{PT}_n)^{\text{ISI}}$$

where PT is the time for the patient plasma sample to clot, PT_n is the mean clotting time of a non-anticoagulated population, and ISI is the international sensitivity index that is assigned by the manufacturer of the reagent according to a standard protocol developed by the WHO. The ISI value ranges from 1.0 to approximately 3.0. A low ISI indicates that the reagent has the highest sensitivity to vitamin K-dependent clotting factors. Typically, the variable PT_n is generated by the laboratory with each new reagent lot by running 20 samples and calculating the mean.

The ProTime system enhances clinical utility for the physician and patient by providing results in both formats, so that anticoagulant dosing can be managed using a familiar format. Since results reported in PT seconds depend on the sensitivity (ISI) of the reagent employed, the physician has the option of changing the ISI value in ProTime so that the ProTime results reported in PT seconds closely match the results reported by the hospital laboratory. To change the ISI used in the calculation, the physician simply accesses the programming screens and enters the ISI of the historical laboratory reagent. (Please see the **For Professionals Only – Program Mode** section).

The default parameters used to calculate PT seconds are $\text{ISI} = 1.0$, which is the sensitivity of the thromboplastin reagent used in the ProTime cuvette, and $\text{PT}_n = 13.1$ seconds, which was established in clinical trials. The PT_n value used in the calculation is automatically adjusted when the reporting ISI is changed.

Table 1. Options for programming ISI and associated PT_n

ISI	1.0	1.2	1.4	1.6	1.8	2.0	2.2	2.4	2.6
PT_n	13.1	12.9	12.8	12.6	12.5	12.3	12.2	12.0	11.9

Table 2. Relationship of INR to PT seconds with varied ISI

INR	PT seconds vs. ISI			
	ISI = 1.0	ISI = 1.6	ISI = 2.0	ISI = 2.4
1.0	13.1	12.6	12.3	12.0
1.5	19.6	16.2	15.1	14.2
2.0	26.1	19.4	17.4	16.0
2.5	32.6	22.3	19.5	17.6
3.0	39.2	25.0	21.3	19.0
3.5	45.7	27.6	23.0	20.2
4.0	52.2	30.0	24.6	21.4
4.5	58.7	32.3	26.1	22.5

Table 2 shows examples of how the calculation of PT seconds is affected by changing the ISI. As the ISI increases, the plasma equivalent PT decreases at any given INR. For example, a patient with INR = 3.00 would have PT seconds = 39.2 if the ISI is the default ISI = 1.0 or PT seconds = 21.3 if the ISI is reset to ISI = 2.0.

Expected Values

ProTime measures both normal and therapeutic prothrombin times in fresh whole blood. Results are displayed in plasma equivalent seconds and INR. Expected values for patients taking oral anticoagulants depend on the patient's disease state and the target values established by the physician.

Sample Condition	INR	PT Seconds (ISI = 1.0)
Normal	0.8-1.2	10.4-15.7 sec
Low anticoagulation	1.5-2.0	19.6-26.1 sec
Moderate anticoagulation	2.0-3.0	26.1-39.2 sec
High anticoagulation	2.5-4.0	32.6-52.2 sec

Note: The ProTime instrument converts whole blood clotting time to INR, then calculates PT seconds. The default ISI used in this calculation is 1.0. Please see the **For Professionals Only – Program Mode** section if you wish to select an ISI that is more closely aligned with the local lab reagent.

The ProTime instrument has been tested extensively by doctors, nurses and patients. Eighty-four people, ages 7 to 81, participated in the first home-use trial. The trial was conducted to see how well ProTime results compare to the lab. Comparisons of this type are described by correlation values and a correlation value near 1.0 means the comparison is good. The home tests compared to tests run at the clinic using ProTime with venous samples (421 samples) had a correlation of 0.92. The home test results compared to a reference laboratory results (368 samples) had a correlation of 0.87. These correlation values mean that home testing is comparable to lab testing.

Precision

Precision testing was conducted with two levels of standard control plasma substrate preparations (Level I and Level II).

A. Standard ProTime Cuvette

Control	Precision	n	Mean	SD
Level I	Within day	17	0.9	0.06
	Day-to-day (5 days)	4/day	1.0	0.08
Level II	Within day	19	3.2	0.19
	Day-to-day (5 days)	4/day	3.2	0.12

B. ProTime3 Cuvette

Control	Precision	n	mean	SD
Level I	Within day	18	0.9	0.07
	Day-to-day (5 days)	4/day	0.9	0.12
Level II	Within day	20	4.0	0.19
	Day-to-day (5 days)	4/day	4.2	0.22

Accuracy

INR results generated by the ProTime and ProTime3 cuvettes using venous and fingerstick whole blood samples were compared to INR values obtained using standard laboratory plasma PT methods with samples collected in 3.2% sodium citrate tubes. The following accuracy data was obtained.

A. Standard ProTime Cuvette vs. Lab (Plasma)

Sample Type	Regression Equation	r	n
Fingerstick	$Y = 0.94x + 0.38$	0.95	229
Venous	$Y = 0.91x + 0.44$	0.94	232

B. ProTime3 Cuvette vs. Lab (Plasma)

Sample Type	Regression Equation	r	n
Fingerstick	$y = 1.05x + 0.07$	0.95	229
Venous	$y = 0.97x + 0.19$	0.95	219

Clinical Performance Comparison

A. ProTime3 vs. ProTime

Linear regression containing clinical fingerstick results from three sites yielded a regression equation as below:

$$y = 1.02x - 0.14 \quad r = 0.94 \quad n = 229$$

B. Patient Self-Testing

In a trial of patient self-testing (PST) in the home vs. professional testing in the clinical and reference lab, equivalent ProTime results were obtained.

ProTime Professional vs. Patient Self-Testing

$$y = .94x + 0.13 \quad r = 0.92 \quad n = 421$$

ProTime Patient Self-Testing vs. Reference Lab

$$y = .77x + 0.38 \quad r = 0.87 \quad n = 368$$

Sensitivity

The ProTime instrument is sensitive to deficiencies in vitamin K-dependent coagulation factors known to influence the PT test (i.e., Factors II, VII and X.)

Hematocrit levels between 20% and 60% do not significantly affect test results.

Quality Control

The ProTime instrument has been designed with multiple systems to ensure proper instrument function. The instrument self-check at startup checks temperature and timing functions, battery level, and optical, electrical and mechanical functions. The instrument does not require further calibration. Each ProTime cuvette has two integral reagent controls that ensure assay reliability and performance. Both levels of control produce quantifiable clotting endpoints that are compared to pre-set acceptance limits programmed in the instrument.

Other in-process instrument QC features and the integral reagent controls function together to ensure correct sample size and collection technique, correct test procedure, instrument functionality and reagent integrity. A fault message is displayed instead of PT results when any instrument or reagent quality criterion is not met. When a fault message is displayed, the user should review the product instructions and repeat the test.

Additional external control materials may be used to check instrument function, reagent integrity and user technique. Each institution should establish a quality assurance program for prothrombin time testing that complies with existing local, state and federal regulations as applicable.

As with all diagnostic tests, the ProTime Microcoagulation System test results should be scrutinized in light of a specific patient's condition and anticoagulant therapy. Any results exhibiting inconsistency with the patient's clinical status should be repeated or supplemented with additional test data.

TECHNICAL ASSISTANCE

For further information and technical assistance, contact ITC Technical Support at 1-732-548-5700, 1-800-631-5945, or e-mail us at techsupport@itcmed.com.

SUGGESTED READING

Adcock DM, Kressin DC, Marlar RA. Effect of 3.2% vs 3.8% Sodium Citrate Concentration on Routine Coagulation Testing. *Am J Clin Pathol* 1997;107:105-10.

Brien WF, Baskerville JC, Taberner DA, Crawford L. Calculation vs. Calibration Curve for INR Determination: Results of an Interlaboratory Proficiency Scheme. *Am J Clin Pathol* 1999;111:193-201.

Eckman MH, Levine HJ, Pauker SG. Effect of Laboratory Variation in the Prothrombin-Time Ratio on the Results of Oral Anticoagulant Therapy. *N Engl J Med* 1993;329:696-702.

Fairweather RB, Ansell J, van den Besselaar AMHP, Brandt JT, Bussey HI, Poller L, et al. College of American Pathologists Conference XXXI on Laboratory Monitoring of Anticoagulant Therapy. *Arch Pathol Lab Med* 1998;122:768-81.

Gosselin R, Owings JT, White RH, Hutchinson R, Branch J, Mahackian K, et al. A Comparison of Point-of-Care Instruments Designed for Monitoring Oral Anticoagulation with Standard Laboratory Methods. *Thromb Haemost* 2000;83:698-703.

Hirsh J. Antithrombotic Therapy in Deep Vein Thrombosis and Pulmonary Embolism. *Am Heart J* 1992;123:1115-22.

Hirsh J, Poller L. The International Normalized Ratio. A Guide to Understanding and Correcting its Problems. *Arch Intern Med.* 1994 Feb 14;154(3):282-8.

Hubbard AR, Margetts SM, Weller LJ, Macnab J, Barrowcliffe TW. An International Collaborative Study on the INR Calibration of Freeze-Dried Reference Plasmas. *Br J Haematol.* 1999 Mar;104(3):455-60.

Levine M HJ, Landefeld, Raskob G. Hemorrhagic Complications of Anticoagulant Treatment. *CHEST* 1992;102:352s-63s.

Technical Assistance

For further information and technical assistance, call ITC Technical Support at 1-732-548-5700, 1-800-631-5945, or e-mail us at techsupport@itcmed.com.

Ordering Information

For further information on ordering and supplies, contact your ProTime distributor.

SAFETY STANDARDS

The ProTime instrument complies with the following safety standard requirements and directives:

CSA C22.2. 601.1.	Medical Electrical Equipment –General Requirements for Safety
EN 60601-1/ UL/IEC 60601-1	Medical Electrical Equipment – General Requirements for Safety
EN 60601-1-2 / IEC 60601-1-2	Medical Electrical Equipment – Part 1-2 – General Requirements for Safety – Collateral Standard: Electromagnetic Compatibility – Requirements and Tests

Directives: 89/336/EEC and as amended by 91/263/EEC, 92/31/EEC, 93/68/EEC and 98/13/EC, 98/79/EC

Equipment Classifications As Defined Per UL 60601-1:2003/ IEC60601-1 2nd Edition

- Protection against electrical shock: Class II, Internally Powered Equipment with no applied parts
- Protection against ingress of liquids: Ordinary (no protection as defined by IEC 529)
- Product cleaning and disinfection: Only according to recommendations of the manufacturer's accompanying documentation
- Mode of operation of equipment: Continuous
- Degree of safety of application in the presence of flammable anesthetic mixture with air, oxygen or nitrous oxide: Not Suitable

NOTE: As defined in the above standards, the classification of "Not Suitable" DOES NOT MEAN that the instrument is not suitable for use in an Operating Room (OR) environment. Rather, it is intended to indicate that the instrument is not suitable for use in the direct presence of a flammable anesthetic mixture with air, oxygen or nitrous oxide.

All relevant documentation is kept on file at ITC.



E233358



MEDICAL ELECTRICAL EQUIPMENT
WITH RESPECT TO ELECTRICAL
SHOCK, FIRE AND MECHANICAL
HAZARDS ONLY IN ACCORDANCE
WITH STANDARD NO.
UL60601-1, CAN/CSA C22.2 NO.601.1
58YB

Guidance and Manufacturer's Declaration – Electromagnetic Emissions

The ProTime® Microcoagulation System is intended for use in the electromagnetic environment specified below. The customer or the user of the ProTime® Microcoagulation System should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 1	The ProTime® Microcoagulation System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The ProTime® Microcoagulation System is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Complies	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The ProTime® Microcoagulation System is intended for use in the electromagnetic environment specified below. The customer or the user of the ProTime® Microcoagulation System should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient / Burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical domestic, commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical domestic, commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	$< 5\% U_T$ ($> 95\%$ dip in U_T for 0.5 cycle) $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $< 5\% U_T$ ($> 95\%$ dip in U_T) for 5 sec	$< 5\% U_T$ ($> 95\%$ dip in U_T for 0.5 cycle) $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $< 5\% U_T$ ($> 95\%$ dip in U_T) for 5 sec	Mains power quality should be that of a typical domestic, commercial or hospital environment. If the user of the ProTime® Microcoagulation System requires continued operations during power mains interruptions, it is recommended that the ProTime® Microcoagulation System be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical domestic, commercial or hospital environment.

NOTE: U_T is the AC mains voltage prior to application of the test level.

INDEX

attention label	4	professional guidelines	
battery care	8	reporting results	31
battery information	8	programming mode	27
charging	10	quality control	33
battery replacement	8	reporting	31
blood coagulation systems	2	expected values	32
blood collection	11	reportable range	29
data history	19	safety features	
LQC records	20	calibration	5
patient records	19	cuvette disposal	7
print/send	21	limitations	6
finger incision preparation	11	precautions	6
insert a cuvette	11	reagent preparation	7
intended use	2	reagents	5
introduction	2	storage	7
ISI		service and maintenance	
setting	29	battery facts	7
main menu	16	instrument disposal	7
memory		routine maintenance	7
clear	30	servicing	7
month		set up options	22
setting	27	beep	26
OID		day	23
entering and changing	24	hour	23
setting	25	language	22
password		minutes	23
changing	30	OID	24
entering and changing	24	password	24
performance characteristics		PID	24
accuracy	33	specifications	7
comparisons	33	specimen collection	9
patient self-testing	33	technical assistance	14, 34, 35
precision	32	test methodology	2
sensitivity	33	test procedures	10
PID		therapeutic range	
entering and changing	24	setting	28
setting	25	TR limits	
preparation		setting	28
contents	9	troubleshooting	14
unpacking and inspection	9	year	
		setting	27



International Technidyne Corporation
8 Olsen Avenue • Edison, NJ 08820 USA
tel: 732.548.5700 • fax: 732.248.1928
www.itcmed.com
a subsidiary of Thoratec Corporation

IR5246 1/08