



Operation Instructions

[Panorama and Cephalo Capabilities]



Thank you for purchasing the Veraviewepocs 2D.

For optimum performance and safety, read this manual thoroughly before using the equipment. Pay special attention to the cautionary warnings and safety statements.

Keep this manual in a handy place for ready reference.

CONTENTS

	Page
Prevent Accidents	1
For Safe Operation	
Safety and Accident Prevention for the Operation of Electrical Medical Devices	
•	
Parts Identification and Accessories	
(I) Parts Identification	
(II) Patient Frame and Arm Operation Panels	9
(III) Accessories	13
Operation	14
(I) Preliminary Procedures	
(II) Operation Procedures	
Safety Check	
Panorama Exposures	
(1) Turn main switch on	
Emergency Stop Switch	
(2) Cassette Insertion	19
(3) Panorama Settings	
Auto Exposure (Digital Direct Auto Exposure)	
2) Manual Exposure	
3) No X-ray Setting	
(4) Patient Positioning	
(5) Panorama Exposures	
(6) Patient Egress and Image Transmission	
1) Panorama	39
2) Maxillary Sinus (Mag.: 1.5 ×, throughout)	42
3) I TMJ Quadruple (Mag.: 1.3×, throughout)	44
(8) Removing the Digital Cassette	
Cephalo Exposures (option)	
(1) Turn Main switch On	51
Emergency Stop Switch	52
(2) Cassette Insertion	
(3) LA (Lateral) Exposure	
(4) Patient Positioning	
(5) X-ray emission	
(6) Patient Egress and Image Transmission	
(7) PA (posteroanterior) Exposure	
(8) Removing the Digital Cassette	
(III) After Use	/1
CHERTACHUS SWIICH CHE	/ 1

Page

Maintenance, Parts Replacements, and Storage				
(I)	Regular Maintenance			
(II)	Replacement Parts	72		
(III)	Storage	72		
Regula	ar Inspection	73		
Troub	ole shooting	70		
DDAE	E Verification Procedure	83		
Techn	ical Description	102		
(I)	Technical Description			
(II)	Meaning of the Symbols	113		
(III)	Replacement Parts and List of accessories	115		
Apper	ndix- Electromagnetic declaration	102		

This manual covers a fully equipped model; refer to the sections covering the instruments and functions of your own unit.

Prevent Accidents

ATTENTION CUSTOMERS

Do not fail to receive clear instructions concerning the various ways to use this equipment as described in this accompanying Operator's Manual.

Fill out and sign the warranty and give the dealer from whom you purchased the equipment his copy.

ATTENTION DEALERS

Do not fail to give clear instructions concerning the various ways to use this equipment as described in this accompanying Operator's Manual.

After instructing the customer in the operation of the equipment, have him fill out and sign the warranty. Then fill in your own section of the warranty and give the customer his copy. Do not fail to send the manufacturer's copy to J. Morita Mfg. Corp.

PREVENT ACCIDENTS

Most operation and maintenance problems result from insufficient attention being paid to basic safety precautions and not being able to foresee the possibilities of accidents. Problems and accidents are best avoided by foreseeing the possibility of danger and operating the unit in accordance with the manufacturer's recommendations. First thoroughly read all precautions and instructions pertaining to safety and accident prevention; then, operate the equipment with the utmost caution to prevent either damaging the equipment itself or causing bodily injury.

Note the meaning of the following symbols and expressions:

MARNING

This warns the user of danger of death, serious bodily injury or total equipment damage and failure or fire.

⚠ NOTE

This alerts the user to the risk of light to medium injury or equipment

damage.

Usage Note

This alerts the user of important points concerning operation.

The user (i.e., clinic, hospital etc.) is responsible for the management, maintenance and use of medical devices

Disregarding the information on safety is considered ABNORMAL USE.

Special training is not required for the operation of this unit.

Do not use this equipment for anything other than its specified dental treatment purpose.

 $\label{lem:caution:equal} \textbf{Caution: Federal law restricts this device to sale by or on the order of a dentist (for U.S.A.).}$

INTENDED OPERATOR PROFILE

a) Qualification:

Legally qualified persons such as radiologists, radiologic technicians, physicians or dentists for X-ray device operation (it may differ among countries).

b) Education, Knowledge:

Intended for professional use, with a basic knowledge of X-ray hazards and its protection, anatomy and image diagnosis.

c) Language understanding:

English or German (Intended for professional use as described above).

d) Experience:

Experienced Person with operating X-ray devices.

ATTENTION

The J. Morita Mfg. Corp. will not be responsible for accidents, instrument damage, or bodily injury resulting from


- (1) repairs made by personnel not authorized by the J. Morita Mfg. Corp.
- (2) any changes, modifications, or alterations of its products.
- (3) the use of products or instrument made by other manufacturers, except for those procured by the J. Morita Mfg. Corp.
- (4) maintenance or repairs using parts or components other than those specified by the J. Morita Mfg. Corp. and other than in their original condition.
- (5) operating the instrument in ways other than the operating procedures described in this manual or resulting from the safety precautions and warnings in this manual not being observed.
- (6) workplace conditions and environment or installation conditions which do not conform to those stated in this manual such as improper electrical power supply.
- (7) fires, earthquakes, floods, lightning, natural disasters, or acts of God.

The J. Morita Mfg. Corp. will supply replacement parts and be able to repair the product for a period of 10 years after the manufacture of the product has been discontinued.

For Safe Operation

WARNING

- Do not use this unit with patients who have a pacemaker. This equipment could cause a pacemaker to malfunction.
- Electromagnetic wave interference could cause this device to operate in an abnormal, random and
 possibly dangerous manner. Cellular and PHS telephones, transceivers, remote controls and all other
 devices which transmit electromagnetic waves located inside the building should be turned off.
- The following devices could cause the equipment to operate in a random manner and should not be located within the x-ray protection area, or be turned off while the equipment is in use.
 - 1. Electrical diagnostic, examination or treatment devices.
 - 2. Personal computers
- The Veraviewepocs must be installed in an X-ray shield location. Local regulation for radiation protection must be observed.
- If the Veraviewepocs is not enclosed by an X-ray booth or other protective barrier, everyone except the patient must stay outside the area shown in the illustration during X-ray emission.
 - The X-ray protection area should consist of a wall, floor and ceiling with a minimum of 1.5 mm lead shielding or its equivalent and should have glass windows with a minimum of 1.5 mm lead shielding or its equivalent, through which the operator can observe the patient. A sign should clearly identify the area as an X-ray protection area, and a caution sign should light up during X-ray emission. Observe local regulations.



- The patient must be provided with appropriate X-ray protection gear such as lead-impregnated clothing that conforms to local regulations.
- Proper infection control procedures must be established and maintained for each patient.
- It is strongly suggested that no unauthorized personnel be in the immediate area when the equipment
 is in use.
- This unit is not designed for and must not be used for "fluoroscopic examinations".
- Proper radiation safety precautions must be established in accordance with local, state and governmental
 regulations in regards to operator and patient protection. The ultimate responsibility lies with the
 owner/operator to ensure that the protection requirements of national and local codes are met.
- Judgement and caution should be used in regards to radiographs of pregnant women. The decision should be based on "clinical need of diagnostic information".
- The operator must be able to see the exposure emissions lights and hear the audible signal during operation of the equipment.
- The operator must be able to see and hear the patient during the operation of the equipment.
- MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENTS.
- Portable and mobile RF communications equipment can affect MEDICAL ELECTRICAL EQIPMENT.
- Use of parts other than those accompanied or specified by J. Morita Mfg. Corp. may result in increased EMC emissions or decreased EMC immunity of the EQUIPMENT.
- The EQUIPMENT should not be used adjacent to or stacked with other equipment and that if adjacent
 or stacked use is necessary, the EQUIPMENT should be observed to verify normal operation in the
 configuration in which it will be used.
- Focal Spot and Skin Distance should be as great as possible to keep the absorbed X-ray dosage as low as reasonably achievable.

「**コメント [SK1]:** X700 に無いが、ここには 残す

MARNING

- Do not fail to turn the unit off after use; this will eliminate the risk of electrical leakage and accidents.
- The main switch should be turned off during standard maintenance procedures. It is also recommended that the main breaker or fuse on the main electrical distribution panel be turned off as some electricity reaches the main power terminal inside the unit even when the main switch is off.
- Use special care to explain to the patient proper instructions as to when they can move as there are "multiple movements" with certain exposures.
- Do not press the panel too hard not to damage the panel.


Safety and Accident Prevention for the Operation of Electrical Medical Devices

- 1. Only fully trained and qualified personnel may operate equipment.
- 2. Items to be duly noted when installing equipment.
 - 1) Locate the unit in a place where it will not get wet.
 - Install the unit in a location where it will not be damaged by air pressure, temperature, humidity, direct sunlight, dust, salts, or sulfur compounds.
 - 3) The unit should not be subjected to tilting, excessive vibrations, or shocks (including during shipping and handling).
 - 4) Do not install the unit where chemicals are stored or where gas may be released.
 - 5) Follow all electrical specifications including frequency (Hz), voltage (V), and current capacity (A) (power consumption).
 - 6) The equipment must be properly grounded.
 - (Everything which requires a large power supply such as X-ray units, chairs, air compressors etc. must have a three point power connection.).
- 3. Item to be duly noted before use.
 - 1) Inspect all switch connections, polarity, dial settings, meters etc. to confirm that the equipment will operate properly.
 - 2) Confirm that the ground is connected properly.
 - 3) Confirm that all cords are connected properly.
 - 4) Take into consideration that simultaneous use of more than one instrument or device can create a dangerous situation or lead to a mistake in diagnosis.
 - 5) Reconfirm the safety of external circuits or systems which are connected directly to the patient.
- 4. Item to be duly noted during use.
 - 1) Do not exceed volumes or times necessary for the diagnosis or treatment.
 - 2) Maintain a constant vigilance for abnormal conditions in both the equipment and the patient.
 - 3) Appropriate steps, such as shutting the equipment down, should be devised to protect the safety of the patient in case any abnormalities in the equipment or the patient are observed.
 - 4) Make sure the patient does not handle or manipulate the equipment.
- 5. Item to be duly noted after use.
 - 1) Turn the power off after returning dials, switches etc. back to their original positions in the prescribed order
 - 2) Do not use excessive force or pull the cord itself to disconnect cords.
 - 3) The following items should be considered when storing the equipment:
 - a. The storage area should protect the equipment from getting wet.
 - b. The storage area should protect the equipment from any possible damage due to atmospheric pressure, temperature, humidity, wind, direct sunlight, dust or air containing salts or sulfur.
 - The equipment should be protected from tilting, vibrations, percussive shocks, etc. (including when it
 is being moved).
 - d. The storage area should be free of chemicals and gases.
 - 4) All accessories, cords, guides etc. should be cleaned, properly arranged and carefully put away.
 - 5) Before storage, the equipment should be cleaned so that it is ready to be used again.
- 6. In case of a malfunction or defect, operators should attach a written notice indicating that the equipment is out of order without attempting to repair the equipment themselves; repairs should be referred to a qualified serviceman.
- 7. Equipment should not be modified in any way.
- 8. Maintenance and Inspection.
 - 1) All equipment and components should be inspected regularly.
 - Equipment which has not been used recently should always be inspected to confirm that it functions properly and safely before being put back into use.


Parts Identification and Accessories

(I) Parts Identification


Main Unit



Control box




Patient Frame




7

Cephalo Unit (Option)



(II) Patient Frame and Arm Operation Panels

Patient Frame Operation Panel





1. Ready key

Before the Ready key is pressed, the green LED for it will be blinking on and off. Press the Ready key in order to complete patient positioning. The arm will rotate slightly and the positioning beams will light up. Pressing the Ready key also enables X-ray emission. During image data transmission the LED will go out.

2. Lift Up Key

Hold down to raise the lift.

3. Lift Down Key

Hold down to lower the lift.

4. Image Layer Beam Forward Key

Hold down when unit is in Ready mode to move the beam forward. For panorama, line this beam up with the distal side of the upper, left canine.

5. Image Layer Beam Backward Key

Hold down when unit is in Ready mode to move the beam backward. For panorama, line this beam up with the distal side of the upper, left canine.

6. Beam On/Off Keys

Positioning beams automatically go off after 1 minute. Press either one of these keys to turn them back on. Or press one of them to turn the beams off.

7. Incisal Occlusion Key

This key is used to make a panorama exposure for a patient with standard occlusion biting on a mouthpiece. After adjusting the height of the positioning (AF) sensor, press this key to automatically move the arm backwards or forwards to the optimum position.

8. Natural Occlusion Key

This key is used to make a panorama exposure for a patient with standard occlusion without using a mouthpiece. After adjusting the height of the positioning (AF) sensor, press this key to automatically move the arm backwards or forwards to the optimum position.

9. Protruding Maxillary Occlusion Key

This key is used to make a panorama exposure for a patient with protruding maxillary without using a mouthpiece. After adjusting the height of the positioning (AF) sensor, press this key to automatically move the arm to the optimum position.

10. Protruding Mandibular Occlusion Key

This key is used to make a panorama exposure for a patient with protruding mandible without using a mouthpiece. After adjusting the height of the positioning (AF) sensor, press this key to automatically move the arm to the optimum position.

11. Frame Liquid Crystal Display (LCD)

Display various information including Ready mode, Image layer beam position, exposure mode, exposure conditions, patient positioning guide etc.

12. Frankfort Plane Beam Up/Down Knob

Turn the knob to the right to lower the beam or to the left to raise the beam. For panorama exposures, lines up with the patient's Frankfort plane.


13. Auto Positioning (AF) Sensor Beam Up/Down Knob

Turn the knob to the right to lower the beam or to the left to raise the beam. Adjust the beam so it strikes the center of the mouthpiece in the patient's mouth.

Arm Operation Panel

Usage Note

♦ Do not press down with excessive force on any of the operation panels. Do not press on the panels with any sharp objects like ballpoint pens or fingernails etc.



1. Panorama Exposure Keys

Use these keys to make various settings for panorama exposures. The LED for each key lights up when it is selected.

Region Keys : Dental Arch \(\bigcap \), Maxillary Sinus \(\bigcap \bigcap \), and TMJ \(\bigcap \) : Adult and Child Magnification : $1.3 \times \times 1.3$ and $1.6 \times \times 1.6$

: Standard A, Shadow Reduction A, and Orthoradial

2. Cephalo Exposure Keys

Projections

Use these keys to select the type of cephalo exposure, either LA (lateral) \bigcirc or PA (posteroanterior) .

Press the auto density compensation key $\frac{\text{DENSITY}}{\text{COMP}}$ to turn it on and off.

The LED for each key lights up when it is selected.

3. X-ray Emission Keys

Use the X-ray Emission keys (Auto Level $^{AUTO}_{LEVEL}$, kV , and mA) and Up/Down keys $^{\bullet}$ to set X-ray emission conditions.

4. Auto and Manual Emission Keys

For auto exposures, press LEVEL and confirm that the corresponding LED lights up.

For manual exposures, press $\begin{tabular}{ll} kV \end{tabular}$ or $\begin{tabular}{ll} mA \end{tabular}$ and confirm that the corresponding LED lights up.

5. Up/Down Keys

Press the Up/Down keys to change the value displayed for the kV or mA, whichever is selected. (When the Down key is held down continuously, all displayed values disappear, and the unit is set for the No X-ray mode.)

6. X-ray Emission Display

Shows the Auto Level or the tube voltage (kV) and tube current (mA).

7. Memory Key

Use this key to save the X-ray emission conditions currently set.

8. Exposure Time Display

Shows the expected exposure time before the exposure and the actual time after the exposure.

9. Ready key

Press this key to return the arm its start position.


This also enables X-ray emission.

10. Ready LED

When the power switch is turned on, this LED blinks.

When the Ready key is pressed, it stops blinking and remains lit continuously.

Cephalo Operation Panel (Option)



(III) Accessories

- Panorama (Pan) or Panorama & Cephalo (Pan/Ceph) Cassette (1)
- Mouthpieces (1 Box of 50)
- Chin Rest (1)
- Lip-nose Rest (1)
- Hand X-ray Plate (Option) (1)
- Temple Stabilizer Band (1)



13


Operation

- * If an accident occurs, the equipment must not be used until repairs have been completed by a qualified and trained technician provided by the manufacturer.
- * Have patients remove glasses, necklaces, earrings and other accessories which could interfere with diagnosis.
- * If the unit has not been used for some time, make sure it operates normally and safely before actually using it.

MARNING

In case of lightning, avoid the risk of electrical shock: Stop using the equipment immediately and have the patient move away from it. Do not touch the equipment or the main power cord.

(I) Preliminary Procedures



Turn the computer on and set up the application software to receive the image data read-in.

<For the i-Dixel application>

The i-Dixel application will start up automatically. Open the patient database and select a patient folder where the image will be read-in, or create a new patient folder.



Click Lock Folder for Data Read-in

* For details refer to the user manual for the i-Dixel application.



Check Resolution

Move the mouse cursor to the resolution icon () on the task bar at the lower right to display the presently selected resolution.




Change Resolution

one without a check mark.)

Click the resolution icon () on the task bar at the lower right to select the resolution.

(Either High resolution or Super high resolution will be marked with a check. To change the resolution, click the

- * Fine high-speed mode: High resolution mode
- * Superfine mode: Super high resolution mode







Check Panorama Operation

Turn main switch on.

Insert the PAN or PAN/CEPH cassette into the cassette holder and press the Dental Arch key. Press the Ready key to enable X-ray emission. Hold down the emission button and check that the arm rotates, X-rays are emitted, the X-ray emission LED lights up, and the audible signal sounds. Also check that after the emission time has elapsed, X-ray emission stops and the arm stops rotating. Press the emission button again to return the arm to its patient entry position.

Check Cephalo Operation

Turn main switch on.

Insert the PAN/CEPH cassette into the cephalo cassette holder and press the LA key on the operation panel.

Completely close the temple stabilizers for panorama exposures by turning the adjustment knob.

Press the Ready key to move the X-ray head and the arm into their Cephalo positions. Check that the Ready LED on the Arm Operation Panel or the control box is on. Hold down the emission button and check that X-rays are emitted, the X-ray emission LED lights up, and the audible signal sounds. Also check that after the emission time has elapsed, X-ray emission stops, the emission LED goes out and the audible emission signal stops.

Usage Note


 The unit cannot be turned on if the Emergency Stop switch has been accidentally pressed during cleaning etc. Turn the switch in the direction indicated by the arrow to release it.

(II) Operation Procedures



Safety Check


For safety, keep fingers away from moving parts when they are moved.



\triangle NOTE

Keep fingers away from gaps and openings for moving parts such as the cassette and its holder and the temple stabilizers as well as the holes on the support column for threaded bolts.









Panorama Exposures

(1) Turn main switch on.

Press the top ($\bf l$)of the main switch. The blue main LED will light up to show that the unit is on.

⚠ NOTE

Do not turn the main switch on if the patient is standing near the unit. The arm will move and it could hit the patient.

Emergency Stop Switch


In case of an emergency, press the Emergency Stop Switch. This will stop the arm's rotation, lift movement and X-ray emission. Do not use this switch for anything. If the Emergency Stop Switch is pressed, the lift will stop within 10 mm and the arm rotation within 15 degrees.

After Pressing the Emergency Switch

Guide the patient away from the unit and turn the main switch off. This will return the unit to a safe condition. Turn the switch in the direction indicated by the arrow and restart the computer. Then turn the main unit back on and check that it is set for Panorama Exposure. If the unit cannot be returned to a safe condition or will not operate, contact your local dealer or the J. Morita Corp.


⚠ NOTE

The image will be lost if the Emergency Stop Switch is pressed during its transmission or if the main switch is turned off.



(2) Cassette Insertion


* Insert the PAN or PAN/CEPH cassette.



1. Turn the cassette so that the detector side faces the X-ray head.


Hold the bottom of the cassette with one hand and the grip with the other. Push down the button and carefully slide the cassette into the holder until the rod goes inside.

Then release the button.



2. Slide the cassette all the way in until the rod goes into its hole inside the holder. The rod will make and audible click and the button will pop out. A beep will also sound. After a few seconds the green LED on the cassette will start to blink on and off.

Depending on the LAN set up, it may take longer for the LED to start blinking.





⚠ WARNING

- Do not bump, jiggle, vibrate or use excessive force.
- Use both hands to put the cassette in and take it out. The cassette wieghs about 2 kg, and it could injure your foot if you dropped it. The sensor would also be damaged.
- Do not swing the cassette around. The vibrations could disturb the sensor's adjustment.

<< Digital Cassette>>

Use both hands to slowly and carefully insert or eject the cassette from the main unit while pressing the handle button.



FRAGILE

This cassette is expensive to replace and contains a highly accurate sensor that may be damaged

if handled incorrectly.

Please pay close attention to the following precautions:

- Do not drop or bump it or subject it to excessive vibration or stress.
 Do not store in locations subject to high temperatures, high humidity or direct sunlight.
 Do not store in locations where it may get wet or be exposed to chemical solutions.
 Do not both in locations where it may get wet or be exposed to chemical solutions.
 Do not both the connectors as static electricity may damage the sensor.
 Do not both the connectors as static electricity may damage the sensor.
 Do not both the subject and did may damage the connectors or cause poor electrical contact.
 Use both hands to slowly and carefully insert or eject the cassettle

in the immuniture.

Indie the cassettle gentity as excessive force may damage the connectors, issettle's Weight: 2Kg) and the cassettle using a dry cloth, and make sure moisture does not seep inside it.

<Storage Conditions> Temperature: -10°C ~50°C Humidity: less than 85% (without condensation) Avoid direct sunlight.

↑ NOTE

• Do not store under conditions of high temperature and humidity.

Storage Temperature; -10°C~50°C

Storage humidity: 5~85 % relative humidity with no condensation.

No frequent or continuous exposure to direct sunlight.

- Do not get water, detergents or chemical solutions on the cassette. These could seep inside and damage it.
- Put the cassette in and take it out with great care. Using excessive force could damage the connection components.
- Clean by wiping with a soft, dry cloth. Do not touch the connection components directly with bare fingers etc. Static electricity could damage the sensor.
- Do not allow dust or other debris to collect on the connection components; this could cause poor
- Slide the cassette all the way into the holder and make sure it is locked into place.

Otherwise images could be distorted, transmission could be interrupted, or equipment could be damaged.



LED Color and Condition

Green:

Normal.

Amber On:

For X-ray emission and image transmission.

Red Bllinking: Abnormal. Transmission malfunction. Cannot make exposure.

Off:

Power is off or cassette is not properly inserted.

⚠ NOTE

Slide the cassette all the way into the holder and make sure it is locked into place.


Otherwise images could be distorted, transmission could be interrupted, or equipment could be damaged.

* Never rotate the arm manually.

Rotating the arm by hand too fast can cause the arm motor to generate an electric current, and this could activate the protection circuit for the motor switching power. This, in turn, would cause the motors to stop working.


If the protective circuit has been activated, turn the unit off. Wait about 1 minute and then turn it back on again. This will deactivate the protective circuit.

21



(3) Panorama Settings

When the unit is turned on, the operation panel is set to the factory defaults shown in the photo to the left.



1) Auto Exposure (Digital Direct Auto Exposure)

Auto exposure regulates tube kV and mA for each individual patient.

Press the Auto Level key to use auto exposure.



The current auto level value is shown on the left side of the Auto Level key. This value can be changed. Up and Down keys will also appear to the right; use these to change the settings.

Auto level can be adjusted from +4 to -4 in increments of 1.

* The factory setting for the Auto Level is 0.



2) Manual Exposure

Press the kV or mA key to switch to manual exposure.




Press the kV key to display the tube voltage value. This value can be changed. Use the Up and Down keys to increase or decrease the value.

Press the mA key to display the tube current value. This value can be changed. Use the Up and Down keys to increase or decrease the value.

kV can be set from 60 to 80 kV in increments of 1 kV, and mA from 1 to 10 mA in increments of 1 mA.

* The factory settings are 65 kV and 5 mA.



3) No X-ray Setting

Use this to rotate the arm to make sure it will not hit the patient or for other tests which do not require X-ray emission.

To activate No X-ray mode, hold down the Down key until the Auto Level, kV and mA displays go out.

To restore normal operation, hold down the Up key.

MARNING

To use the No X-ray mode, make sure the Auto Level, kV and mA displays are blank.





(4) Patient Positioning

Press the Ready key. The arm will automatically move into position for patient positioning.

The green Ready LED will stop blinking and stay on. The Frankfort plane, Mid-sagittal, Image layer and AF beams will light up.

Usage Note

◆ Before pressing the Ready key, make sure the temple stabilizers are closed all the way.

⚠ WARNING


The laser beams could damage the eye. Do not look directly into them or let them strike any one directly in the eye.

The beams go off automatically after 1 minute. Press either one of the beam on and off keys to turn them back on again.


* The Image layer beam will not light up if the Ready key has not been pressed.

⚠ NOTE

- Have patients remove glasses, earrings and any other accessories that might interfer with making a good exposure.
- Keep fingers away from gaps and openings for moving parts such as the cassette and its holder and the temple stabilizers as well as the holes on the support column for threaded bolts.
- Do not let patients touch the switches on the Frame Operation Panel.
- Take care that the patient's does not get caught in the cassette or other moving parts.







1. <Initial Patient Positioning>

Put an X-ray protection apron on the patient. Have the patient perform the incisal occlusion while biting on an unused mouthpiece.

Line up the center of the mouthpiece with the center of the upper and lower incisors.

Have the patient stand in front of the chin rest. Look at him from the side, and have him pull in his chin and straighten his back. Then look from behind and make sure his shoulders and back are straight.

⚠ WARNING

- Have the patient wear and X-ray protection apron.
- Always use a new, unused mouthpiece to prevent contagion.
- * Keep mouthpieces in a clean, uncontaminated area.

 Use the Lift Up or Down Key to raise or lower the chin rest and line it up with the patient's chin. Release the key when the chin rest is at the right height.
 The Lift has a Slow Start and Slow Stop. It will stop automatically if there is an excessive load on it.

Usage Note

- ◆ Always use the Up and Down Keys to move the Lift. Never try to force it; this would damage it.
- ♦ The over current protection circuit could be activated after 3 seconds for reasons such as an excessive load or the failure to insert the cassette. The lift will stop and an audible signal will sound. The unit will stop working and one of the following messages will appear in the display on the frame:

LIFT OVERLOAD! PROTECTOR ACTIVE LIFT OVERTIME! PROTECTOR ACTIVE


To restore normal working order, press any one of the occlusal type keys underneath the display on the patient frame.

⚠ WARNING


Take care that the temple stabilizers do not strike the patient in the eye.

⚠ NOTE

- Do not use excessive force to close the temple stabilizers.
 - This could be uncomfortable for the patient or damage the stabilizers.
- Forcing the patient in or out could also damage the stabilizers.



3. Use the knob to open the temple stabilizers. Have the patient move forward without slouching or otherwise changing his posture. Then have him put his chin on the chin rest and lightly grip the patient handles. Make sure his shoulders are lowered.




4. Close the temple stabilizers until they lightly touch the patient's temples.



Make sure the patient's face is straight and then line up the Frankfort plane beam with the patient's frankfort plane using the beam Up or Down key.

⚠ WARNING

The laser beam could cause eye damage. Do not look directly into it or let it strike you or the patient in the eye.



6. Move the patient's head to the left or right until the Mid-sagittal beam lines up with the patient's mid-sagittal plane. Then tighten up the temple stabilizer so that the patient's face will not move.




Head Band

For short people and children, use the Head Band to hold the patient's head in place.

Usage Note

♦ Take care that the patient's hair does not get caught up in the Head Band.



Attachment

To attach the Head Band, hold the stabilizer rod steady with one hand, line up the notch on the end of the Head Band with the rod and firmly press it into place.




Removal

To take the Head Band off, hold the rod steady with one hand and pull it off carefully. Do not simply yank it off forcefully.

Usage Note

• Do not apply excessive force to the rods; they might break.






7.-1 < Auto Positioning>


Use the Auto Positioning (AF) Sensor Beam Up/Down knob to raise or lower the beam so that it strikes the center of the patient's mouthpiece.

Usage Note

- Make sure the mouthpiece is perpendicular and the beam strikes the center of it. If the beam is off center, the image layer will not be properly detected and the image could be spoiled.
- ◆ Condensation on the lens for the AF sensor beam could prevent proper detection of the image layer.

 Use a soft cloth to wipe condensation off the lens.
- * For incisal occlusion using a mouthpiece, the image layer extends to about 20 mm from the end of the mouthpiece or about 8 mm from the edges of the incisors.





Press the Incisal Occlusion Key, which is used for auto positioning with the mouthpiece.

The arm (and the Image layer beam) will move to line the image layer up with the patient, and the value of the image layer will be displayed. Check that the Image layer beam is lined up with the distal side of the upper left canine.

There are three types of auto positioning that do not use the mouthpiece.

For these, the AF sensor beam is lined up with the center of the upper incisors.

- * For natural occlusion, press the Natural Occlusion Key
- * For a protruding maxillary, press the Protruding Maxillary Key $\boxed{\langle}$
- * For a protruding mandible, press the Protruding Mandible Key


In each case, the arm will move to match the image layer with the patient's dentition and the image layer value will be displayed.

Usage Note

◆ If the mouthpiece is not used, the AF sensor beam may not detect the correct image layer position for edentulous patients, patients with twisted upper incisors, or patients with lustrous crowns or full orthodontic bands.

In these cases, use Manual Positioning.

* If a patient is not standing in place or is out of the range for the arm's back and forth movement (+20 mm to -20 mm), a continuous beeping will indicate an error when a patient type key is pressed.



7.-2 Manual Positioning

Use the Backward and Forward keys to line the Image layer beam up with the distal side of the patient's upper, left canine.

⚠ WARNING

The laser beam could cause eye damage. Do not look directly into it or let it strike you or the patient in the eye.

* Depending on the shape of the patient's face, the Image layer beam may not directly strike the distal side of the upper left canine.


In this case, position the beam so that an imaginary extension of it is lined up properly.

The beam's movement range is normally from ± 20 mm to ± 20 mm, but it is less for some types of exposures.


Make sure the Frankfort plane, Mid-sagittal and Image layer beams are lined up properly.



Warn the patient not to move during the X-ray exposure (while the melody is playing). If the patient moves, he could be hit by the arm or the exposure could be a failure.



* However, wheel chairs with a width greater than 480 mm will not fit.






Usage Note

- ♦ If the patient is not properly positioned the image may not be useful for diagnosis. Refer to the examples below to better understand proper patient positioning.
- 7.-3 Examples of Patient Positioning and Image Results

Correct Positioning






Face Lowered








Face Turned Up




Jaw bone partially cut off



Face Turned to the right



Enlarged / Reduced`
Reduced X-ray absorption here
(Excessive saturation)



Face Turned to the Left



Reduced Enlarged
Reduced X-ray absorption here
(Excessive saturation)




Image layer beam too far forward



Enlarged





Image layer beam too far back



Reduced



- * The "Double-Ready Function" reduces the length of time that the emission button needs to be held down.
 - Press the Ready key a second time after completing patient positioning in Ready mode; a two-toned beep will sound and the arm will move to its exposure start position. The Ready LED will light up again and the Image layer beam will go out.
- * In Double Ready mode, X-ray emission starts as soon as the emission button is pressed.
- * In Double Ready mode, other operation keys and switches are disabled except for the Ready key and Lift keys.
- * In Double Ready mode, press the Ready key again to return to the patient positioning Ready mode.

Usage Note

♦ If the patient moves out of position after the unit is set for Double Ready, press the Ready key again to return the unit to its normal Ready setting and re-position the patient.


(5) Panorama Exposures

 Make sure the green Ready LED is on; check the arm, patient frame or control box.









Pick up the handswitch and hold down the emission button. The arm will start to rotate and X-rays will be emitted.


During X-ray emission, the yellow Emission LED on the control box will light up and an audible signal will sound.

* During X-ray emission, the LED on the cassette will be amber.

Usage Note

♦ Error number 18 (Err. 18) appears if the computer is not ready.

If this happens, turn off the main switch. When the computer is ready, turn the main switch back on. If the cassette is not in the panorama cassette holder, insert it and press the Ready key on the arm operation panel again.



3. Keep holding the emission button down. X-ray emission and arm rotation will stop, the Emission LED will go out, and the audible signal will stop. Then the arm will go to its patient egress position. This completes the exposure procedure.

Now release the emission button.

The Ready LED will blink and be amber to show that the image is being transmitted. After transmission is completed, the Ready LEDs on the arm and control box will change to green and blink on and off.

⚠ WARNING

- Always leave the X-ray booth, and press the emission button outside of it.
- In case of an emergency, release the emission button; this will completely stop the unit.

⚠ NOTE

- Warn the patient not to move during emission (while the melody is sounding). If the patient moves, he could be hit by the arm or the image might be ruined.
- Hold the emission button down until the entire procedure is completed; releasing the button will terminate the exposure procedure.
- If an exposure has been terminated before completion, guide the patient away from the unit. Check that the LED on the cassette is green and blinking, Then press the Ready key to move the arm into the patient entry position. Repeat the patient positioning and then make the exposure.



4. Hang the handswitch back in its holder.



(6) Patient Egress and Image Transmission

1. <Patient Egress>

After the exposure, the arm will automatically go to the patient egress position, $90\ degrees$.

Open the temple stabilizers all the way up and guide the patient away from the unit. Then close the stabilizers all the way.

MARNING

- Never use excessive force to move the patient away; this could injure the patient or damage the stabilizers.
- Take care that the stabilizers do not strike the patient in the eye.
- * Take the mouthpiece from the patient and throw it away.



Press the Ready key on the arm operation panel or the patient frame to automatically move the arm to the patient entry position.

The Ready LEDs on the operation panel and the control box will light up green to show that the unit is ready to make another panorama exposure.

⚠ WARNING


- Make sure the patient is clear of the unit before returning the arm to the patient entry position; otherwise it could hit the patient.
- In case of an emergency when the arm is returning to the patient entry position, press the emission button, the Ready key on the operation panel, or the emergency switch to stop it.

⚠ NOTE

If the arm return operation is performed without closing the temple stabilizers or without guiding the patient away from the unit after the exposure, the arm operation panel will display "Frr. 42."

operation panel will display "Err. 42."

If you go ahead and press the Ready key without closing the temple stabilizers or without guiding the patient away from the unit, the arm will rotate to its start position. This can cause the arm to strike the temple stabilizers or the patient.



3. During image transmission the LED on the cassette will light up amber. After transmission is completed, it will change to green and start blinking. ("Transferring image!" will appear in the computer display, and after about 30 seconds the image will appear.)

⚠ NOTE

During image transmission the LED on the cassette will light up amber. After transmission is completed, it will change to green and start blinking.


Do not do any of the following things until this LED has changed to green and starts blinking:

- 1. Do not turn off the main switch
- 2. Do not take the cassette out of its holder


Any of the above actions will cause the image to be lost and the computer to freeze.

Usage Note

- ◆ After image transmission, a two-toned beep will sound and the LED will change to green and start blinking. However, another exposure cannot be made until the image appears in the computer's display.
- ♦ During image transmission, pressing the emission button will result in a series of beeps, but X-rays cannot be emitted.



- * Auto density compensation is applied to create the optimum image. However, if some areas of the image are exceptionally dark, the auto density compensation will tend to make the entire image whiter than usual.
- * For an enlarged digital image, there is a junction line that is not visible when the image is displayed initially. However, if the image is magnified, it appears as a fine, horizontal line through the center of the image. This line is where the two CCD sensors come together.




(7) Exposure Regions

1) Panorama

1. Press the Dental Arch key.



2. Press either the adult or the child key.




Pedodontic Panorama

For a Pedodontic panorama, the arm's angle of rotation and exposure range are reduced; the X-ray dosage is also reduced by from 10% to 15%.

* The pedodontic panorama is for people who have a small jaw bone. If the jaw bone is too large, the TMJ may not appear in the image.

For the entire jaw to appear in the image, the length of a horizontal line from the center of the incisors to the outer ear orifice should be less than 70 mm.




3. Press either 1.3 or 1.6.



Enlarged Panorama

If you select $1.6\times$ magnification, the image will have about 20 % more information for the dental arch than $1.3\times$ magnification. However, the TMJ may be cut off.




4. Select one of the three projections (X-ray beam angles).




Standard Panorama Projection






Shadow Reduction Projection





Orthoradial Projection



(1) Standard Panorama Projection

Very good for making measurements for implants etc.

Press the Standard Projection Key.

(2) Shadow Reduction Projection

Reduces shadows obscuring the mandibular ramus.

Press the Shadow Reduction Key.

(3) Orthoradial Projection


Reduces overlapping of teeth because X-ray beam is at a right angle to the dental arch.


Press the Orthoradial Projection Key.

2) Maxillary Sinus (posterior) (Mag.: 1.5×, throughout)


[To examine the posterior maxillary sinus or facial injuries.]

Press the Maxillary sinus key.





Maxillary Sinus Panorama



Chin Rest Position

Put the chin rest in the lowest position to make a maxillary sinus exposure.

* If, for children or short people, the AF sensor beam does not strike the mouthpiece even when it is at its lowest setting, set the chinrest at its Medium position.



Low (for maxillary sinus)

Line up the Low Groove (see illustration) with the back of the chin rest holder.



Medium (maxillary sinus for children and short people)

Line up the Medium Groove (see illustration) with the back of the chin rest holder.





High (standard panorama)

Line up the High Groove (see illustration) with the back of the chin rest holder.

Usage Note

♦ Be careful when replacing the chin rest with the lip-nose rest; it could break if it is dropped.





Mag.: 1.3x

3) 1 1 TMJ Quadruple (Mag.: 1.3×, throughout)

Press the TMJ key.

Press either the Adult key or the Child key.

Estimated Distance between Joints


Adult : 100mm Child : 90mm Thickness : 10.5mm Length : 54mm

Select the size best for the patient.

Four images will appear in the computer display: one each for the mouth open and closed on both sides.

The X-ray beam angle is optimum for the average distance between the joints and the average length.

The arm rotates twice to make a complete set of images.



Patient Positioning


 Replace the chin rest with the lip-nose rest set at medium height.



Put an X-ray protection apron on the patient and stand in front of the lip-nose rest. Look at him from the side, and have him pull in his chin and straighten his back. Then look from behind and make sure his shoulders and back are straight.

⚠ WARNING

Have the patient wear and X-ray protection apron.



3. Open the temple stabilizers and use the Up or Down Lift key to raise or lower the lip-nose rest to the right height for the patient. Release the key when the rest is at the right height.

№ NOTE

Keep your fingers away from moving parts, gaps between the cassette and its holder and for head stabilization components, and holes in the support column.

Usage Note

- ◆ Always use the Up and Down Keys to move the Lift. Never try to force it; this would damage it.
- ♦ The over current protection circuit could be activated after 3 seconds for reasons such as an excessive load or the failure to insert the cassette. The lift will stop and an audible signal will sound. The unit will stop working and one of the following messages will appear in the display on the frame:

 LIFT OVERLOAD! PROTECTOR ACTIVE

LIFT OVERLOAD: PROTECTOR ACTIVE

To restore normal working order, press any one of the occlusal type keys underneath the display on the patient frame.



4. Use the knob to open the temple stabilizers. Have the patient move forward without slouching or otherwise changing his posture. Then have him put his upper lip on the lip-nose rest and lightly grip the patient handles. Make sure his shoulders are lowered.


⚠ WARNING

Take care that the temple stabilizers do not strike the patient in the eye.



 Keep the patient's Frankfort plane horizontal and align the beam with the patient's external auditory orifices by turning the Frankfort plane beam Up/Down knob. Then hold his head in place by lightly closing the temple stabilizers.

Have the patient move his head left or right until the Mid-sagittal beam lines up with his mid-sagittal plane. Then tight up the temple stabilizer so that the patient's face will not move.



- Use the Backward and Forward keys to line the Image layer beam up with the patient's ear orifice.
 (Arm moves from +20 mm to -16 mm.)
- The image layer is about 12 mm in front of the Image layer beam.
 Have the patient close his mouth.




Head Band

For short people and children, use the Head Band to hold the patient's head in place.

Usage Note

 Take care that the patient's hair does not get caught up in the Head Band.





Attachment

To attach the Head Band, hold the stabilizer rod steady with one hand, line up the notch on the end of the Head Band with the rod and firmly press it into place.

Removal

To take the Head Band off, hold the rod steady with one hand and pull it off carefully. Do not simply yank it off forcefully.

Usage Note

 Do not apply excessive force to the rods; they might break.

⚠ WARNING

The laser beams could cause eye damage. Do not look directly into them or let them strike you or the patient in the eye.

⚠ NOTE

Warn the patient not to move during emission (while the melody is sounding). If the patient moves, he could be hit by the arm or the image might be ruined.

- * The "Double-Ready Function" reduces the length of time that the emission button needs to be held down.
 - Press the Ready key a second time after completing patient positioning in Ready mode; a two-toned beep will sound and the arm will move to its exposure start position. The Ready LED will light up again and the Image layer beam will go out.
- * In Double Ready mode, X-ray emission starts as soon as the emission button is pressed.
- * In Double Ready mode, other operation keys and switches are disabled except for the Ready key and Lift keys.
- * In Double Ready mode, press the Ready key again to return to the patient positioning Ready mode.

Usage Note

◆ If the patient moves out of position after the unit is set for Double Ready, press the Ready key again to return the unit to its normal Ready setting and re-position the patient.

Open and Closed Mouth Exposures

Check that the Ready LED (green) on the patient frame, operation panel, or control box is on.







Emission button Verarieusepocs Emission LED

AUTO LEVEL & W W MEMORY READY



Mouth Closed Exposure (First Exposure)

1. Pick up the handswitch and hold down the emission button. The arm will go to its starting position, start to rotate and exposures will be made of the left and right sides. During X-ray emission, the yellow Emission LED on the control box will light up and an audible signal will sound. Release the emission button when the arm stops rotating.


⚠ WARNING

- Always leave the X-ray booth, and press the emission button outside of it.
- In case of an emergency, release the emission button; this will completely stop the unit.

Have the patient stand still and maintain a good posture.
 Press the Ready key on the patient frame or the arm operation panel, or press the emission button.
 The arm will go back to its starting position.

⚠ NOTE

Warn the patient not to move during emission (while the melody is sounding). If the patient moves, he could be hit by the arm or the image might be ruined.



Open Mouth Exposure (Second Exposure)

- 1. Have the patient open their mouth.
- Pick up the handswitch and hold down the emission button. The arm will start to rotate and exposures will be made of the left and right sides. Release the emission button when the arm stops rotating.
 The Ready LED on the control box will go out to show that the exposure procedure has been completed.

⚠ WARNING


- Always leave the X-ray booth, and press the emission button outside of it.
- In case of an emergency, release the emission button; this will completely stop the unit.





- 3. Open the temple stabilizers all the way up and guide the patient away from the unit.
 - Press the Ready key on the patient frame or the arm operation panel.
 - The arm will go back to the patient entry position.



4. Remove the lip-nose rest and replace it with the chin







(8) Removing the Digital Cassette

- * Make sure the green LED on the cassette is either blinking or out.
- Support the cassette holder with one hand and press the release button. Pull the cassette out a little and then release the button.

Usage Note

- Make sure you press the button down far enough; otherwise, the attachment rod may not come all the way out of its hole, and the cassette will not come out.
- 2. Hold the cassette with both hands and carefully slide it out along the rail on the bottom of the cassette holder.

⚠ WARNING


- Do not bump, jiggle, vibrate or use excessive force.
- Use both hands to put the cassette in and take it out.
 The cassette wieghs about 2 kg, and it could injure your foot if you dropped it. The sensor would also be damaged.
- Do not swing the cassette around. The vibrations could disturb the sensor's adjustment.


^NOTE

- Do not store under conditions of high temperature and humidity.
 - Storage Conditions: Temperature: $-10^{\circ} \sim 50^{\circ} C$. Relative Humidity: $5 \sim 85$ % (without condensation). No frequent or continuous exposure to direct sunlight.
- Do not get water, detergents or chemical solutions on the cassette. These could seep inside and damage it.
- Put the cassette in and take it out with great care.
 Using excessive force could damage the connection components.
- Clean by wiping with a soft, dry cloth. Do not touch the connection components directly with bare fingers etc. Static electricity could damage the sensor.
- Do not allow dust or other debris to collect on the connection components; this could cause poor contact.











Cephalo Exposures (option)

(1) Turn Main switch On

1. Press the top ($\bf I$) of the main switch. The blue main LED will light up to show that the unit is on.

- 2. Press either the LA (Lateral) or PA (posteroanterior) key to set the unit for cephalo exposure.
- 3. Before pressing the Ready key, make sure the temple stabilizers for panorama exposures are completely closed.
- 4. Press the Ready key.

The X-ray head will turn to its cephalo direction, and the arm will rotate.

The arm and head automatically go into and lock in their cephalo positions.

Usage Note

◆ If the temple stabilizers are not closed or the patient has not exited after the exposure, the LCD on the patient frame will display "Err. 42."

\triangle NOTE

Make sure a patient is no longer near the unit, before pressing the Ready key; otherwise he could be hit by the arm.

5. The green Ready LED will light up.

Usage Note

♦ Do not move the arm manually.

The arm may not be set in the proper cephalo position if it is suddenly moved by hand. This is also true if the arm hits the operator's shoulder or something while it is moving. If the arm is accidentally turned by hand or is touched during movement, press the Dental Arch key. Then press either the LA or PA key once again. After this press the Ready key.



Emergency Stop Switch

In case of an emergency, press the Emergency Stop Switch. This will stop the arm's rotation, lift movement and X-ray emission. Do not use this switch for anything. If the Emergency Stop Switch is pressed, the lift will stop within 10 mm and the arm rotation within 15 degrees.

After Pressing the Emergency Switch

Guide the patient away from the unit and turn the main switch off. This will return the unit to a safe condition. Turn the switch in the direction indicated by the arrow and restart the computer. Then turn the main unit back on and check that it is set for Panorama Exposure. If the unit cannot be returned to a safe condition or will not operate, contact your local dealer or the J. Morita Corp.

⚠ NOTE

The image will be lost if the Emergency Stop Switch is pressed during its transmission or if the main switch is turned off.



(2) Cassette Insertion


* Insert the PAN/CEPH cassette, used for both panorama and cephalo exposures. Only the PAN/CEPH cassette can be used to make cephalo exposures.




 Press the Cephalo Start Position Key. Both the cassette holder and the secondary slit plate will move forward.



 Hold the bottom of the cassette with one hand and the grip with the other. Push down the button and carefully slide the cassette into the holder until the rod goes inside. Then release the button.







3. Slide the cassette all the way in until the rod goes into its hole inside the holder. The rod will make and audible click and the button will pop out. A beep will also sound. After a few seconds the green LED on the cassette will start to blink on and off.

⚠ WARNING

- Do not bump, jiggle, vibrate or use excessive force.
- Use both hands to put the cassette in and take it out. The cassette weighs about 2 kg, and it could injure your foot if you dropped it. The sensor would also be damaged.
- Do not swing the cassette around. The vibrations could disturb the sensor's adjustment.

⚠ NOTE

- Do not store under conditions of high temperature and humidity.
 - No frequent or continuous exposure to direct sunlight.
- Do not get water, detergents or chemical solutions on the cassette. These could seep inside and damage it.
- Put the cassette in and take it out with great care.
 Using excessive force could damage the connection components.
- Clean by wiping with a soft, dry cloth. Do not touch the connection components directly with bare fingers etc. Static electricity could damage the
- Do not allow dust or other debris to collect on the connection components; this could cause poor contact.
- Slide the cassette all the way into the holder and make sure it is locked into place.
 Otherwise images could be distorted, transmission could be interrupted, or equipment could be damaged.

Usage Note


- Never slide the cassette in with excessive force. This could break the lock rod or damage the electrical connection components.
- 4. Press the Patient Entry and Egress Key to move the cassette holder and secondary lit plate all the way back.

Usage Note

♦ Do not move the cassette holder or the secondary slit plate manually. Suddenly moving these components by hand could cause damage or a malfunction.







(3) LA (Lateral) Exposure

1. Press the LA key.



2. Turn the DENSITY COMP key on to select both the required soft and hard tissues for making cephalo measurements. (DENSITY COMP: automatic density compensation)

- * 80 kV will appear in the display when the DENSITY COMP key is turned on.
- * Cephalo exposures cannot be made with Auto Exposure.
- * In Cephalo mode, the No X-ray function can be used to check the movement of the cassette holder and the secondary slit plate without emitting X-rays.





(4) Patient Positioning

 Use the Up or Down key to raise or lower the craniostat match the patient's height. Release the key to stop the craniostat moving.


Usage Note

- ◆ The craniostat Up and Down keys will not work if the unit is not set for cephalo and the Ready key has not been pressed to turn the X-ray head in the cephalo direction
- ♦ The over current protection circuit could be activated after 3 seconds for reasons such as an excessive load or the failure to insert the cassette. The craniostat will stop and an audible warning signal will sound. The unit will stop working and one of the following messages will appear in the display on the frame: LIFT OVERLOAD! PROTECTOR ACTIVE LIFT OVERTIME! PROTECTOR ACTIVE To restore normal working order, press any one of the occlusal type keys underneath the display on the patient frame.
- ◆ Always use the Up and Down Keys to move the craniostat. Never try to force it; this would damage it.

⚠ NOTE

Keep your fingers away from moving parts, gaps between the cassette and its holder and craniostat components, and holes in the support column.









2. Grip the base of the ear rod plates and rotate the craniostat to its Lateral position.

\triangle NOTE

The nasion plate must be at the front; a proper exposure cannot be made if it is reversed.


Usage Note

- Rotate the craniostat with the nasion plate is down and inside.
 - Turning the craniostat roughly could cause it to hit the cassette holder and damage it.
- ♦ To move the nasion plate, grip the base of its arm. Otherwise this could be damaged.
- 3. Grip the ear rod plates with both hands and open them up all the way.

4. Make it easy for the patient to take his place by raising the nasion plate and moving it out.




5. Put the ear rods on their studs.




6. Have the patient stand directly underneath the craniostat.

MARNING

Do not fail to have the patient wear an X-ray protection apron.



7. With the Up or Down key, raise or lower the craniostat until the ear rods line up with the patient's outer ear orifice and then release the key.




8. Grip the ear rod plates with both hands and carefully close them until the ear rods go into the patient's ears.

MARNING

Be extremely careful when inserting the ear rods and do not move the craniostat after they have been inserted.

This could seriously injure the patient.





9. Use the key for the Frankfort Plane Beam to line it up with the patient's Frankfort Plane and make sure that it is perpendicular to the patient's mid-saggital plane. Adjust the length of the nasion plate and put it against the patient's nasion.

⚠ WARNING

The laser beam could cause eye damage. Do not look directly into it or let it strike you or the patient in the eye.

Frankfort Plane Beam Key





⚠ WARNING

Never let the patient hold onto the patient guard; his fingers could be pinched between it and the cassette holder causing an injury.



11. The cassette holder and secondary slit plate will move into their starting positions.

\triangle NOTE

- Warn the patient not to move during emission (while the melody is sounding). If the patient moves, he could be hit by the arm or the image might be
- Make the cassette holder or the secondary slit plate will not hit the patient's shoulder.


(5) X-ray emission

1. Check the arm operation panel, patient frame or control box, and make sure the green Ready LED is on.










2. Pick up the handswitch and hold down the emission button.

After a few seconds the secondary slit and cassette holder will start to move and X-rays will be emitted. During X-ray emission, the yellow Emission LED on the control box will light up and an audible signal will sound.

* During X-ray emission the LED on the cassette will light up and be amber.



 Keep holding the emission button down. When the exposure is finished, the cassette holder and secondary slit plate will stop moving and X-ray emission will also stop. The Emission LED will go out, and the audible signal will stop.

Now release the emission button.

The Ready LED will blink and be amber to show that the image is being transmitted. After transmission is completed, the Ready LEDs on the arm and control box will change to green and blink on and off.

⚠ WARNING


- Always leave the X-ray booth, and press the emission button outside of it.
- In case of an emergency, release the emission button; this will completely stop the unit.

⚠ NOTE

- Warn the patient not to move during emission (while the melody is sounding). If the patient moves, he could be hit by the arm or the image might be ruined.
- Hold the emission button down until the entire procedure is completed; releasing the button will terminate the exposure procedure.
- If an exposure has been terminated before completion, guide the patient away from the unit. Check that the LED on the cassette is green and blinking, Then press the Ready key to move the arm into the patient entry position. Repeat the patient positioning and then make the exposure.



- 4. Hang the handswitch back in its holder.
- * Hanging a mirror on the wall where the patient can see it helps to keep the patient from moving during X-ray emission.




(6) Patient Egress and Image Transmission

1. Use both hands to carefully open the ear rod plates and get the ear rods out of the patient's ears.

⚠ WARNING

Open the ear rod plates very carefully and make sure the ear rods are well clear of the patient's ears; otherwise the patient could be injured.



2. Pull the nasion plate out and then raise it up, the guide the patient away from the unit.



3. During image transmission the LED on the cassette will light up amber. After transmission is completed, it will change to green and start blinking. ("Transferring image!" will appear in the computer display, and after about 30 seconds the image will appear.)

⚠ NOTE


During image transmission the LED on the cassette will light up amber. After transmission is completed, it will change to green and start blinking.

Do not do any of the following things until this LED has changed to green and starts blinking:


- 1. Do not turn off the main switch
- 2. Do not take the cassette out of its holder
 Any of the above actions will cause the image to be lost
 and the computer to freeze.

Usage Note

- ◆ After image transmission, a two-toned beep will sound and the LED will change to green and start blinking. However, another exposure cannot be made until the image appears in the computer's display.
- ◆ During image transmission, pressing the emission button will result in a series of beeps, but X-rays cannot be emitted.



- * Auto density compensation is applied to create the optimum image. However, if some areas of the image are exceptionally dark, the auto density compensation will tend to make the entire image whiter than usual.
- * For an enlarged digital image, there are junction lines that are not visible when the image is displayed initially. However, if the image is magnified, they appear as fine, horizontal lines dividing the image in thirds. These lines are where the three CCD sensors come together.



(7) PA (posteroanterior) Exposure

1. Press the PA key.



 Turn the DENSITY COMP key on to select both the required soft and hard tissues for making cephalo measurements. (DENSITY COMP: automatic density compensation)
 (Soft tissues will not be selected if the DENSITY COMP key is turned off.)

 * When the Density Compensation key is turned on, the arm operation panel will display "80 kV."

* Cephalo exposures cannot be made with Auto Exposure.





Patient Positioning

Turn the craniostat into the PA exposure position.

* Otherwise, positioning, X-ray emission, patient egress are all the same as for the LA exposure.

However, the nasion plate should be raised up.


45 Degree Slant Exposure

Set the craniostat at a 45 degree angle to either the right or left.

Set the unit for a PA Exposure.

Hand Exposure


1. Select PA Exposure.














- * 80 kV and 2 mA are rough estimates for a child, but these can be adjusted depending on the patient.
- * Cephalo exposures cannot be made with Auto Exposure.



4. Attach the Hand X-ray Plate to the patient guard.




6. Have the patient place his hand inside the rectangle on the Hand X-ray Plate.





⚠ NOTE

- Make sure nothing other than the patient's hand is inside the rectangle on the Hand X-ray Plate.
- Do not fail to take the Hand X-ray Plate off after completing the exposure. Otherwise, the next cephalo exposure made could be ruined.



- * The X-ray emission procedure is the same as that for cephalo exposures.
- Remove the Hand X-ray Plate after completing the exposure.





(8) Removing the Digital Cassette

- * Make sure the green LED on the cassette is either blinking or out.
- Support the cassette holder with one hand and press the release button. Pull the cassette out a little and then release the button.

Usage Note

- Make sure you press the button down far enough; otherwise, the attachment rod may not come all the way out of its hole, and the cassette will not come out.
- 2. Hold the cassette with both hands and carefully slide it out along the rail on the bottom of the cassette holder.

MARNING

- Do not bump, jiggle, vibrate or use excessive force.
- Use both hands to put the cassette in and take it out.
 The cassette wieghs about 2 kg, and it could injure your foot if you dropped it. The sensor would also be damaged.
- Do not swing the cassette around. The vibrations could disturb the sensor's adjustment.

⚠ NOTE

- Do not store under conditions of high temperature and humidity.
 - Storage Conditions:Temperature: -10°~50°C. Relative Humidity: 5~85 % (without condensation). No frequent or continuous exposure to direct sunlight.
- Do not get water, detergents or chemical solutions on the cassette. These could seep inside and damage it.
- Put the cassette in and take it out with great care.
 Using excessive force could damage the connection components.
- Clean by wiping with a soft, dry cloth. Do not touch the connection components directly with bare fingers etc. Static electricity could damage the sensor.
- Do not allow dust or other debris to collect on the connection components; this could cause poor contact.

Calibrating Digital Cephalo Data for Software Analysis

Digital cephalo data must be calibrated for whatever analysis software you are using.

⚠ NOTE

Measurements will not be correct if the data is not calibrated.

- * Cephalo image data has a resolution of 176 dpi.
- * Refer to the user's manual for your analysis software for instruction on how to calibrate the data.

Calibration Data

* A floppy disk with the calibration data is provided. It is titled "Chart(176dpi).jpg"


How to Use the Calibration Data

- ♦ Import the data into your analysis software.
- Refer to the user's guide for your analysis software to perform the calibration.
- * The distance from the center of one cross to the next in the data image is 27.5 mm. Therefore, the total distance of 4 crosses is 110 mm.
- Based on the above distances, select calibration points and check the distance between them.
 For example, calibrate the software so that the distance from one cross to the fourth cross away from it is 110 mm.
- * Keep the calibration data in a handy, safe place.

 We recommend that you copy the data onto your hard drive in a folder named "C:\Program Files\3dxcom"



Floppy Disk Provided



Chart(176dpi).jpg


(III) After Use

Turn Main switch Off

Press the bottom (\circ) of the main switch to turn it off. The main LED will go out.



Do not fail to turn the unit off after use; this will eliminate the risk of electrical leakage and accidents.





Maintenance, Parts Replacements, and Storage

(I) Regular Maintenance

- Disinfect the Temple Stabilizers, Ear Rods, the Chin Rest, the Nasion Plate, the Hand X-ray Plate, the Lip-nose Rest and the Patient Handles after each patient by wiping them with Ethanol for Disinfection (Ethanol 80 vol%).
- Wipe cassettes carefully with a soft dry cloth.

 (Do not directly touch the connection area or components.)
- Wipe the operation panel with Ethanol for Disinfection (Ethanol 80 vol%).
- Every 6 months, inspect and grease the wire cables for the lift.

⚠ WARNING

Always turn the main switch off before performing maintenance. This will eliminate the risk of shocks, burns, and accidental switch operation which could result in an injury.

Usage Note

- ♦ Use only Ethanol for Disinfection (Ethanol 80 vol%) or a neutral detergent to clean outer surfaces. Never use alkaline or acidic solutions, cresol soap, or other chemical solutions; this could cause discoloration or degrade the materials.
- ◆ If any water, detergent, cleaning solvent, or other chemical solutions get on the outer surfaces, wipe them off immediately with Ethanol for Disinfection (Ethanol 80 vol%).
- ♦ Do not use ozone water to clean the unit. Ozone water could damage to the unit.
- ◆ Do not disinfect the clinic with ozone gas or ultraviolet light. This could damage plastic and rubber components.
- When cleaning the unit with Ethanol for Disinfection (Ethanol 80 vol%), take care that none of
 it seeps inside; this could damage the unit.
- ♦ Do not get the main unit or any of the cassettes wet; this could damage them.

(II) Replacement Parts

- Replace the parts listed in the Regular Inspection List as necessary depending on degree of wear and length of use.
- Order replacement parts from your local dealer or the J. Morita Corp.

(III)Storage

- No frequent or continuous exposure to direct sunlight.
- Store the mouthpieces in an aseptic environment.
- Store the cassettes where they will not get wet and where no chemicals will be spilled on them.
- If the unit has not been used for some time, make sure it works properly before using it again.

Regular Inspection

- Maintenance and inspection are generally consider to be the duty and obligation of the user, but if, for some reason, the user is unable to carry out these duties, he may rely on a qualified medical device serviceman. Contact your local dealer or the J. Morita Corp. for details.
- This unit should be inspected for all the items in the following list once a year.
- At the start and end of each business day, make sure that switching the main switch to ON and OFF turns the equipment on or off without fail.

Regular Inspection List

Power Supply and Physical Stability

1. Power Supply Voltage

Use a digital or analog tester to measure the unit's power supply.

The result must be 220/230/240 VAC±10%.

2. Ground connection

Visually inspection the ground connection to make sure it is securely and properly connected.

3. Floor and base securing bolts

Visually inspect the floor and base securing bolts.

Check that the floor is level and make sure the base bolts have not loosened up.

4. Bolt and screw tightness

Inspect all the bolts and screws on the unit.

Make sure that all bolts are in place and properly secured.

5. Electrical circuitry

Make sure all wiring and connections are intact.

6. LAN cables

Make sure no cables are bent or pinched and that they are all securely connected

7. Outer appearance and labels

Make sure outer covers and panels etc are not damaged, cracked, loose or dirty.

Make sure labels for rating, tube, and lasers are all properly in place and securely attached.

8. Main switch

Turn the main switch on and off and make sure the main LED on the control box lights up.

9. Emergency Stop Switch

Turn on the main switch and then press the emergency stop switch.

Check that the power goes off. (Main LED on the control box goes out.)

10. Patient Auto Positioning (AF)

Put a test object in place. Press the patient type keys and make sure the arm moves forward or backward. Repeat this test 3 times with the test object in a different position each time.

Make sure auto positioning sensor moves up and down smoothly.

Light Beams

Make sure the Mid-sagittal, Frankfort plane and Image layer beams light up and can be turned on and off. Also make sure that the Mid-sagittal, Frankfort plane and Image layer beams move smoothly.

12. Patient Handles

Make sure handles are tight and properly secured.

13. Temple stabilizers and chinrest

Turn the temple stabilizer knob to make sure the stabilizers open and close properly. Make sure the chinrest and lip-nose rest are secure in both their upper and lower positions.

14 Lift Mechanism

Press the up and down keys.

Make sure the lift moves smoothly and stops properly. Repeat this 3 times.

15. Wire Cables

Check the wire cables for broken strands. Make sure the ends are properly secured. Grease the cables with the grease provided by the manufacturer.

16. Operation panel and display

Press all the keys on the operation panel and display and make sure they work.

Press all the keys on the patient frame and make sure they work.

Make sure the LED on the control box lights up when the Ready key is pressed.

17. X-ray Emission and Display

Make sure the Emission LED (yellow) on the control box lights up during X-ray emission and that the audible signal sounds.

18. Digital Cassette

Press the lock button and make sure the lock pin moves smoothly.

Repeat 3 times.

19. Digital Cassette

Visually inspect connector. Make sure it is not cracked and that the pins are not bent.

Panorama Exposures

20. X-ray emission and image read-in

Make an X-ray exposure of a test piece and check the resulting image in the computer monitor.

21. Arm Rotation

Hold down the emission switch to rotate the arm.

Make sure it does not make an abnormal noise or slip and that it stops at the specified point.

22. Arm Emergency Stop

Make sure the arms stops when the emission button is released and when the emergency stop switch is pressed.

23. DDAE Verification

Cephalo Exposure

24. X-ray Head changeover

Close the panorama temple stabilizers. Press either LA or PA cephalo and then the Ready key. Make sure the X-ray head turns around and that it and the arm automatically go into their cephalo positions.

Make the above operation is not performed when the panorama temple stabilizers are open even if the cephalo and Ready keys are pressed.

Repeat each procedure 3 times.

25. X-ray emission and image read-in

Make an X-ray exposure of a test piece and check the resulting image in the computer monitor.

26. Craniostat Movement

Move the craniostat up and down with the up and down keys. Make sure it moves smoothly and stops accurately. Repeat 3 times.

27. Cassette Holder and Secondary Slit

Press the Start Position key and the patient Entrance and Egress key and make sure the cassette holder and secondary slit move smoothly and stop accurately.

Repeat 3 times.

28. Craniostat

Grip the ear rod plates at their base with both hands and make sure they open and close smoothly and go securely into position.

Grip the nasion retainer at the base and make sure it moves to the left and right and up and down smoothly and goes securely into position.

29. Light Beam

Make sure the Frankfort plane cephalo light beam can be turned on and off with its switch.

30. Cephalo Emergency Stop

Release the emission button while the cassette holder and secondary slit are moving and make sure they stop.

Trouble shooting

If the equipment operation does not seem to be normal, check or adjust the following before requesting a repair service.

- * If the equipment does not operate properly after the inspection, adjustment, or parts replacement or if you cannot perform the inspection yourself, contact your local dealer or J. Morita Corp.
- * The inside parts of the equipment are charged with high voltage. Do not attempt to perform maintenance or adjustment that is not described in the troubleshooting table.
- * If an accident occurs, the equipment must not be used until repairs have been completed by a qualified and trained technician provided by the manufacturer.
- * Before conducting the inspection or adjustment, confirm that the Main LED (blue) on the control box is lift
- * Contact your local dealer or the J. Morita Corp. for repairs if the apparatus does not work normally even after performing the steps recommended below.

Main Unit			
Problem Possible Cause		Remedy	
	No power supply	Check the breaker on the distribution panel.	
No power when main switch is turned on.	Emergency Switch has been pressed.	Release emergency switch and turn on main switch.	
Switches don't work Display doesn't work Arm doesn't go to starting point	Some function or process is in progress	Wait a while. If normal operation is not restored, turned the unit off, wait 1 minute, and then turn it on again.	

Exposures and Main Unit			
Problem	Possible Cause	Remedy	
Panorama & Cephalo Image too light Line in image Large white border Partial image Completely black	* Noise Interference	Turn unit off. Make sure of patient and user safety. Turn unit back on and see if it works normally. Use a devoted circuit with a minimum capacity of AC 220/230/240V, 16A and make sure ground is properly connected. Do not use devices which might produce noise during an exposure.	
Unit spontaneously goes back to same condition as when it was turned on at first. Or switches will not work at all	* Short, temporary power cut	Turn unit off. Make sure of patient and user safety. Turn unit back on and see if it works normally. Use a devoted circuit with a minimum capacity of AC 220/230/240V, 16A and make sure ground is properly connected.	

Problem	Possible Cause	Remedy	
Panorama & Cephalo Image has strobe effect	Metal prosthetics can cause excessive feedback in the Auto Exposure system and produce a strobe effect. (AE Strobe)	Check by making exposure without using AE	
Incisor area is blurred		Review patient positioning procedures.	
Center of panorama image is white and left side is contracted.	Poor patient positioning	Patient may have moved after auto positioning was completed. For manual positioning, Image layer beam may not have been lined up with the upper left canine	
Inconsistent density	Application setting	Configure application with tool bar	
Exceptionally dark areas or entire image is too white	Wrong exposure conditions (kV, mA, auto level)	Adjust auto level or kV and mA for manual exposure	
Grainy panorama image Resolution setting		Click the Resolution Icon in the task tray to select resolution with the mouse. Establish communication between the application and the X-ray unit to change the resolution.	
LCD goes out. Weird characters in LCD Error message in LCD Unit spontaneously goes back to same condition as when it was turned on at first. Switches will not work at all	Build up of static electricity	Turn unit off. Make sure of patient and user safety. Turn unit back on and see if it works normally. Make sure ground is properly connected. Maintain room temperature.	
Message in arm display LIFT OVERLOAD! PROTECTOR ACTIVE LIFT OVERTIME! PROTECTOR ACTIVE	Patient is leaning or pulling on frame Cassette is not inserted	Have patient relax and not lean or pull on the frame. Insert cassette.	

Error Numbers

Check the list below if an error number appears in the arm LCD display (Veraviewepocs 3D) or the arm operation panel (Veraviewepocs 2D) and try the suggested remedies.

(The error numbers are the same for both the Veraviewepocs 3D and the Veraviewepocs 2D. For the 2D unit, only the numbers appear in the arm operating panel. Some error numbers are only for the 3D unit.)

Error No.	Possible Cause (appears in LCD)	Remedy
00	Overcurrent protection circuit of CPU1 has been activated.	Turn off, wait 1 minute, turn back on
01	Short pin is not installed on DCN12 of the CPU1 board.	Turn off, wait 1 minute, turn back on
02	X-Ray head overheated! Leave the unit on and wait at least 30 minutes for it to cool off.	Wait at least 30 min. for head to cool off.
03	X Axis did not return to its original position in specified time. Press Ready key again.	Press Ready key and put unit into Ready mode
Y Axis did not return to its original position in specified time. Press Ready key again.		Press Ready key and put unit into Ready mode
05	Arm did not return to its original position in specified time. Press Ready key again.	Press Ready key and put unit into Ready mode
07	Horizontal slit did not return to its original position in specified time. Press Ready key again.	Press Ready key and put unit into Ready mode
08 X-Ray Head did not return to its original angle.		Turn off, wait 1 minute, turn back on
Buzzer does not stop sounding Overcurrent protection circuit of PSU has activated. Turn the unit off and wait at 1 seconds before turning it back on.		Turn off, wait 5 seconds, turn back on
C1 Buzzer does not stop sounding	No communication with the host computer	The C1 error occurs when there is not response from the i-Dixel application. If the application is just starting up, wait for it to finish. If there is not response from the i-Dixel application check the computer and HUB power connections and also check that the i-Dixel application is running. Then turn the main unit off and back on again.

Error No.	Possible Cause (appears in LCD)	Remedy	
10	The digital cassette or the secondary slit on the Cephalo-unit is not working.	Do exposure over again	
11	The patient's head is not in the right position for the intended exposure mode.	Check exposure mode and patient direction, LA or PA	
12	Digital Cassette is not turned on.	Turn off, wait 1 minute, turn back on	
13	Digital Cassette is busy. Remove the cassette and reinsert it	Take cassette out and re-insert	
14	Digital cassette communication failure Remove the cassette and reinsert it	Take cassette out and re-insert	
15	AF function cannot be used in the DOUBLE READY state or during other operations.	Use the AF only when the unit is not engaged in some other operation and not set for Double Ready.	
16	CPU 2 firmware does not support Digital Cephalo.	Turn off, wait 1 minute, turn back on	
17	Install Digital Cassette prior to emission.	Check cassette insertion	
18	No communication with the host computer	Start up computer	
19	CPU communication error.	Turn off, wait 1 minute, turn back on	
20	Failure to send data to CPU2.	Turn off, wait 1 minute, turn back on	
The CPU on the patient frame side was busy and did not receive the acknowledgement signal for the Temple Stabilizer status. To return the arm, have the patient move away, close the temple stabilizers all the way, and then press the Ready key		To return the arm, have the patient move away, close the temple stabilizers all the way, and then press the Ready key	

Error No.	Possible Cause (appears in LCD)	Remedy
22	The analog power for the Digital Cassette is not turned on. Take cassette out and re-insert	Take cassette out and re-insert
24	Exposure with the specified resolution is not available with this cassette.	Turn off, wait 1 minute, turn back on
25	TDI table for this exposure mode cannot be found.	Turn off, wait 1 minute, turn back on
28*	CT scan and 2-Direction scout are not available with this cassette. Replace with CT cassette and press CT key	Replace with CT cassette and press CT key
30*	Cannot establish communication with 3DXD. Check if the application software on the PC is ready. PC may be busy reconstructing CT data or other processes.	Start up i-Dixel if it is not running. Wait for reconstruction to finish if it is under way
32*	Vertical Slit did not return to its original position in specified time. Press Ready key again.	Press Ready key again to put unit into Ready mode
33	Hardware backup timer has been activated.	Turn off, wait 1 minute, turn back on
34	X Axis is not in the ready position for panoramic exposure.	Turn off, wait 1 minute, turn back on
35*	The computer is not receiving the image. Exposure is terminated.	Check LAN cables
36 [*]	No response from the touch panel module.	Turn off, wait 1 minute, turn back on
37*	The touch panel is not connected with its controller.	Turn off, wait 1 minute, turn back on
39	Defect has been detected in High Voltage Circuit. Exposure cannot be performed.	Turn off, wait 1 minute, turn back on
40*	Panoramic exposure cannot be made with this cassette. Insert the proper cassette for Panoramic exposures.	Insert a cassette for making panorama images


Error No.	Possible Cause (appears in LCD)	Remedy
41	Could not establish communication with CPU2. Repeat the operation.	Do procedure again
42	Before pressing the READY key, have the patient exit the unit and close the Temple Stabilizer completely. Failure to do so may result in the arm striking the patient or Temple Stabilizer during rotation.	Have the patient move away, close the temple stabilizers all the way, and then press the Ready key
43	Outside the operational range of the unit. Move the coronal beam backward	Move Image layer beam back
44	Outside the operational range of the unit. Move coronal beam forward	Move Image layer beam forward
45	Arm positioning error	Turn off, wait 1 minute, turn back on
46	X Axis positioning error	Turn off, wait 1 minute, turn back on
47	Y Axis positioning error	Turn off, wait 1 minute, turn back on
48	Vertical Slit positioning error	Turn off, wait 1 minute, turn back on
49	Horizontal Slit positioning error	Turn off, wait 1 minute, turn back on.
51 *	Vertical slit error.	Turn off, wait 1 minute, turn back on
52	Horizontal slit error.	Turn off, wait 1 minute, turn back on
53*	Beyond the maximum vertical range of the horizontal positioning beam.	Review CT exposure procedure

 $^{^{}st}$ These are only for the Veraviewepocs 3D.

Contact your local dealer or the J. Morita Corp. for repairs if the apparatus does not work normally even after performing the steps recommended above.

<Cable Routing Diagram>

Cable connections and routing may differ depending on the model.



DDAE Verification Procedure

1. DDAE Verification

DDAE (Digital Direct Auto Exposure) is verified by this procedure.

1) DDAE Verification Flowchart



2) Warning and Caution

If any errors occur during the verification procedure, turn off Veraviewepocs immediately. Restart the procedure from "Start Verification".

2. Setup

1) Test Piece


DDAE verification uses copper plates attached to Veraviewepocs. Test copper piece consists of three copper plates (1), (2), (3)



2) Set Test Piece

2-2)-1


The chinrest should be removed and the temple stabilizer closed.




Temple Stabilizer

2-2)-2 Set the Test Piece as shown below.

Copper plates (1)+(2)+(3):



Copper plates (1)+(2):



3) Startup Test Program


2-3)-1

Exit the data base application and startup Dxladj.exe


2-3)-2

File→Click the host name setting of connected X-ray device.

Change the IP address as following.
Pan cassette : 192.168.240.17
Pan / Ceph cassette : 192.168.240.25



2-3)-3 Acquire Panorama Imaging Average After the exposure, the "Panoramic Image Average" is display in the status bar.




Verification Procedure

- 1) Exit the database application then run the Dxladj.exe and set the host name.
- 2) Turn on Veraviewepocs.
- 3) Expose X-ray steps 1 through 8 in table 1 and record "Panoramic Image Average" by each step.

Table 1: Exposure Setting of DDAE Verification

	1		
	test piece (Copper plates)	KV mA setting	Radiographic Speed
step 1		Auto 2	Superfine mode (14.9sec) Fine high-Speed mode (7.4sec)
step 2	(1)+(2)+(2)	65KV, 10mA	
step 3	(1)+(2)+(3)	Auto 2	
step 4		65KV, 10mA	
step 5		Auto 2	Superfine mode (14 Occe)
step 6	(1) (2)	65KV, 10mA	Superfine mode (14.9sec)
step 7	(1)+(2)	Auto 2	Fine high Speed made (7.4see)
step 8		65KV, 10mA	Fine high-Speed mode (7.4sec)

4) Use the following flowchart to verify the DDAE.



- *) The step shown in the flowchart means Panoramic Image Average of the step.
- 5) Results

Pass : DDAE operates correctly.
Fail : Please contact your local sales representative.

- 6) Turn off Veraviewepocs.
- 7) Exit Dxladj.exe.

4. Troubleshooting

[Problem] The following error message from Dxladj.exe is displayed:



[Cause] The database application and Dxladj.exe are running at the same time.

[Solution]

Shutdown the database application then restart Dxladj.exe.

Technical Description

(I) Technical Description

X550 Model Type EX-2

Classification

European Directive 93/42/EEC IIb Rule 10

Safety according to IEC-60601-1, IEC 60601-1-1, IEC60601-1-2, IEC60601-1-3, IEC60601-1-4, IEC60601-1-6, IEC60601-2-7, IEC60601-2-28, IEC60601-2-32, ISO14971, ISO10993-1, ISO10993-5 and ISO10993-12 standards

Protection against electric shock Class I Degree of protection Type B Protection against ingress of liquids IPX 0

Disinfection methods:

- Every patient, disinfect the Temple Stabilizers, the Chin Rest, the Nasion Plate, the Hand X-ray Plate, the Lip-nose Rest and the Patient Handles by wiping them with Ethanol for Disinfection (Ethanol 80 vol%).
- Once a day, wipe the operation panel with Ethanol for Disinfection (Ethanol 80 vol%) and wipe the LCD with a dry cloth.
- Single use item; mouth piece.

Mode of operation Intermittent Applied part with no conductive connection to patient.

Product Description

Dental Panoramic X-ray unit with a high frequency switching mode X-ray generator. In addition to panoramic exposure, the unit can also take scanograms. Also cephalometric device is available that uses the panoramic X-ray source.

Intended Use

X550 is used for diagnostics in dentistry by exposing X-ray image receptor and for the use by authorised persons in the practice of dentistry or medicine, and/or associated procedures.

X550 is intended for dental radiographic examination and diagnosis over the whole Dent-maxillofacial area such as teeth, periodontal tissues or chin-bone etc.

Essential Performance (EMC pass/fail criteria)

Performance necessary to achieve freedom from unacceptable RISK.

- No unexpected X-ray irradiation
- No unexpected movement of the equipment

Generator/X-ray Head Assembly

D-051 Tube Focal Spot 0.5 Target Angle Target Material Tungsten

Operating Tube Potential 60 ~ 80kV (±1 kV, 21 steps) Digital Display (accuracy of displayed values ±

Maximum (nominal) 80kV, 10mA

Operating Tube Current 1 mA to 10 mA (accuracy of displayed values ± 10 %)

Maximum Output Power 0.8 kW nominal at 80kV, 10mA

Inherent filtration minimum 2.5 mm Al at 70kV Filtration Beam Quality HVL minimum 3.2 mm AI at 80kV

Primary Protective Shielding Minimum 1.5 mm Pb or equivalent

Outer Shell Temperature 45°C maximum 1:60, 80kV/10mA Duty Cycle

Preheated Filament Rectification Direct Current Oil Cooling Cooling

Maximum Heat Unit of 194.45 kJ (1 HU=1.35 Joule)

X-ray Head Assembly

Maximum Heat Dissipation 2 kJ/min.

of X-ray Head Assembly

Leakage Radiation max. 1.0mGy/h at 1m

Minimum mAs 6.0 mAs Panoramic

Cephalometric 3.5 mAs

Weight of X-ray Head 13.8Kg (with collimators)

Auto Exposure (Not applied for Cephalometric)

Setting parameters +4 - -4 (9 steps) Digital Display

Maximum possible excursion 60kV - 80kV1mA - 10mAMinimum exposure time 4.3 seconds Verification method by test piece

Power requirements

Input Voltage EX-2: AC220/230/240V,

50/60 Hz single phase

Fuse at the distribution panel EX-2:16A, 220/230/240V, slow (220/230/240V)

Power Consumption EX-2: max. 2.0KVA

EX-2: max. 9A/8A/8A (AC220/230/240V, with operation) Rated input amperes

EX-2: max. 1.0A (220/230/240V, stand by)

Power line resistance EX-2: max. 1 Ohm (220/230/240V)

Mechanical parameters

Panoramic SID 520mm Cephalometric SID 1,650mm

Magnification Standard Panoramic: 1.3×,1.6×

Maxillary Sinus Panoramic: 1.5× Pedodontic Panoramic: 1.3×, 1.6× Shadow Reduction Panoramic: 1.3×, 1.6× Orthoradial Panoramic: 1.3×, 1.6×

Scanograms: 1.3×,

Cephalometric: 1.1×

X-ray Field $6 \times 146 \text{ mm (+0/-}30\%)$: Panoramic

 $6 \times 220 \text{ mm (+0/-}30\%)$: Cephalometric

Weight

Panoramic Approximately 184 kg Boxed Cephalometric Approximately 268 kg Boxed

Outer Dimensions

W 1,020 × D 1,330 × H 2,355mm (H 2,185mm Option) Main Unit

(Panoramic, Upgradable Cephalometric)

 \hat{W} 2,000 × D 1,330 × H 2,355mm (H 2,185mm Option) (with

Cephalometric)

Control box W 70 \times D 40 \times H 115mm Vertical Height of Focal Spot 1,055 ~ 1,775mm (Panoramic)

970 ~ 1,605mm (Option)

1,125 ~ 1,775mm (with Cephalometric, Upgradable Cephalometric) 1,040 ~ 1,605mm (Option)

Patient Positioning Auto focus using light sensor for distance measurement and

electrically operated positioning system

89

Patient Positioning Beam Class II Laser. Accordance with 21CFR Part1040.10 and IEC60825-1.

less than 1.7 mmAl Attenuation equivalent of Panoramic Patient head stabilizer Cephalometric less than 1.7 mmAl

Exposure Time

Panoramic max.20.3 seconds Scanograms max.17.2 seconds (8.6x2) Cephalometric max. 5.0 seconds

Accuracy of displayed values ± 10 %

Emission Switch Dead Man Type

Leakage technique factors

80kV, 600mAs/h (80kV, 10mA, duty cycle 1:60, for example 7.4 s exposure per 7 minutes 17 seconds cool-down period)

Measurement bases

The kV is measured by monitoring differentially the current flowing through 450 Mohm, 1% feedback resistors connected between the tube anode and ground.

The mA is measured by monitoring current in the HT return line, which equals the tube current.

Exposure time: Starting point of exposure is determined at the time when the kV value reaches to 75% of average kV value. Termination of exposure is determined at the time when the kV value decreases to 75% of average kV value.

Collimator

Panoramic slit and Cephalo slit.

Digital Cassettes

Digital Panorama cassette and sensor (Option)

Cassette model: XDPC700

External dimensions: $241 \times 280 \times 43 \text{ mm} + /-10\%$

Sensor: Full frame transfer type 2 dimension CCD image sensor

 $\begin{array}{lll} \mbox{Active sensor surface:} & 147.5 \times 6.1 \mbox{ mm} \\ \mbox{Pixel size of the sensor:} & 0.048 \mbox{mm} \\ \mbox{Pixel construction:} & 3,072 \times 128 \mbox{ pixel} \\ \end{array}$

Detail recognition: 0.096 mm pixel size for Superfine mode (Resolution) 0.144 mm pixel size for Fine high-Speed mode

and Digital Cephalometric mode

Resulting image format: Panoramic max. 288 × 147.5 mm

Image size: Panoramic $max.3,000 \times 1,536$ pixel (Superfine mode)

 $max.2,000 \times 1,024$ pixel (Fine high-Speed mode)

Imaging method: Time Delay Integration

Digital Panorama and Cephalo cassette and its sensor (Option)

Cassette model: CDCP700

External dimensions: $241 \times 280 \times 43 \text{ mm +/- } 10\%$

Sensor: Full frame transfer type 2 dimension CCD image sensor

Active sensor surface: $221.2 \times 6.1 \text{ mm}$ Pixel size of the sensor: 0.048 mmPixel construction: $4,608 \times 128 \text{ pixel}$

Detail recognition: 0.096 mm pixel size for Superfine mode (Resolution) 0.144 mm pixel size for Fine high-Speed mode

and Digital Cephalometric mode
Resulting image format: Panoramic max. 288 × 147.5 mm

Cephalometric max. 260×147.5 mm

Image size: Panoramic $\max. 3,000 \times 1,536$ pixel (Superfine mode)

 $max.2,000 \times 1,024$ pixel (Fine high-Speed mode)

Cephalometric $max.1,750 \times 1,537$ pixel

90

Imaging method: Time Delay Integration

SIP/SOP statement

LAN Interface: Unshielded twisted pair cable with RJ-45 plug connections,

Length less than 2 m.

Requirements for computers or other devices connected to the computers

- 1. The Veraviewepocs has been tested and found to comply with the limits for medical devices to the IEC 60601-1-2:2001 for electro magnetic compatibility. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
 - Reorient or relocate the receiving device.
 - Increase the separation between the equipment.
 - Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
 - Consult the nearest J. Morita office, its representative or its dealer for help.
- 2. The following equipment connected to the analog and digital interfaces must be certified according to the respective IEC standards (i.e. IEC 60950 for data processing equipment and IEC 60601-1 for medical equipment). Furthermore all configurations shall comply with the system standard IEC 60601-1-1. Everybody who connects additional equipment to the signal input part or signal output part configures a medical system, and is therefore responsible that the system complies with the requirements of IEC 60601-1-1. If in doubt, consult the nearest J. Morita office, its representative or its dealer for help.
- * Some of the following devices may cause some technical problems with the Veraviewepocs. Ask your nearest J. Morita office for proper selection of equipment and connections.

⚠ NOTE

The following devices may not be located in the X-ray protection area or the patient vicinity except the Hub if the Hub is conformed with IEC60601-1 or IEC60950 and enclosure leakage current is conformed with IEC 60601-1.

* The patient vicinity is the area where intentional or unintentional contact can occur between a patient or a patient's attendant and the above devices, or between a patient or a patient's attendant and other persons touching the above devices. This area extends 1.83 m beyond the perimeter of the bed (examination table, dental chair, treatment booth, and the like) in its intended location, and vertically 2.29 m above the floor.

Other system requirements

Hardware

Windows based Personal Computer (Minimum specifications)

Operating system: Microsoft Windows 2000 with Service Pack 2.

CPU: Intel Pentium IV 1.7GHz or higher, or compatible.

Memory: RAM 512MB

HDD: HDD 20GB or more are recommended.

Video capture board

Video board: resolution of 1024×768 and color depth of 24bit

Network protocol: TCP/IP with static IP address.

Network interface: Universal purpose 10BASE-T Ethernet network interface

board

Port occupied: 69/udp, 2102/tcp, 2102/tcp Others: Network board, CD-ROM drive. Display: 17 inch CRT or TFT LCD

16,000,000 colors

Standard: IEC60950 or IEC60601-1

EMC regulation

Related UL standard (addition to USA) Related C-UL standard (addition to Canada)

Local regulations

Hub

10 Base-T, 100 Base-TX

Standard: IEC60950 if it is used in non patient vicinity

IEC60601-1 or IEC60950 with enclosure leakage current

is conformed with IEC 60601-1.

EMC regulation

Related UL standard (addition to USA) Related C-UL standard (addition to Canada)

Local regulations

Recommended Hub, for example
Manufacturer: Bay Networks
Type: Bay Stack 350T

Storage Device

Patient data can be saved safely. MO or CD-R disk drive is

recommended.

Standard: IEC60950 if it is used in non patient vicinity

EMC regulation

Related UL standard (addition to USA) Related C-UL standard (addition to Canada)

Local regulations

Other equipment connected to PC

Standard: IEC60950 if it is used in non patient vicinity

EMC regulation

Related UL standard (addition to USA)
Related C-UL standard (addition to Canada)

Local regulations

Application Software Application software for image processing or data base is provided by J. Morita.

92

It shall be used with above Windows base computer specifications.

It conforms to 93/42/EEC (in EU), IEC6060-1-4 and 21 CFR (in USA),

Medicaldevice regulations (in Canada).

If another application software is used, it must conform to the above regulations and standards, and must match the Dixel driver from J. Morita Mfg. Corp.

Ask your nearest J. Morita Office for the appropriate interface.

Environmental data

Operating Conditions

Ambient temperature range $+10^{\circ} \sim +35^{\circ}$ C

Relative humidity 20 to 85% no condensation Atmospheric pressure range 700 ~ 1060hPa

Atmospheric pressure range Transport and Storage Conditions

Ambient temperature range Ambient temperature range Relative humidity $5 \sim 85\%$ Atmospheric pressure range $500 \sim 1060\text{hPa}$

Original language

English

Working life

The working-life of this unit is 10 years from the date of shipment provided it is regularly and properly inspected and maintained.

Disposal

The package should be recycled. Metal parts of the equipment are disposed as scrap metal. Synthetic materials, electrical components, and printed circuit boards are disposed as electrical scrap. Material must be disposed according to the relevant national legal regulations. Consult specialized disposal companies for this purpose. Please inquire of the local city/community administrations concerning local disposal companies.



This symbol is affixed to fulfill the requirements of EU Directive 2002/92/ED Article 11.


This equipment cannot be disposed of as unsorted municipal waste within the European Union. Follow local regulations for disposal.

Service


J. Morita products may be repaired and serviced by

- the technicians of J. Morita's subsidiaries all over the world.
- technicians employed by authorized J. Morita dealers and specially trained by J. Morita.
- independent technicians specially trained and authorized by J. Morita.

Tube housing assembly heating curve



Tube housing assembly cooling curve




Tube Rating Chart

Tube rating chart

Maximum Rating Chart (Absolute Maximum Rating Charts)


D-051




Anode Thermal Characteristics

Anode Thermal Characteristics

D-051







Reference axis


X-ray Field Panoramic





X-ray field Cephalometric






(II) Meaning of the Symbols


Patient frame

Rating label * Label contents may differ depending on model.



X-ray head (Arm)

Tube label



Operation Instructions




Manufacturer




Authorised representative in the European Community

Support Column

Main Switch



Control box



On

Mouthpiece Package



Do not reuse

Base



Equipotentiality

Lift Unit



Laser warning label



Emergency Stop Switch



Laser Caution label (EX-2)

Package

Temperature Limitation

(III) Replacement Parts and List of accessories

Replacement Parts

Code No.	Description	Rating	Туре	Qu.
6112473	Main Fuse (EX-2)	F10A250V	Fast-acting, High Breaking Capacity	1

Fuse shall be certified according to IEC127 or 241, or manufactured in proportion to IEC127 or 241.

List of accessories

Code No.	Description	
6270750	Mouth Pieces (100)	
6350210	Chin Rest	
6350207	Lip-nose Rest	
6331000	X550 2D. 3D 3-piece Copper Filter	
6290325	Ear-Rod	
6292325	Hand X-ray Plate	
6211957	Head Band A (with buckle)	
6211958	Head Band B (with buckle)	

- the technicians of J. Morita's subsidiaries all over the world.
- technicians employed by authorized J. Morita dealers and specially trained by J. Morita.
- independent technicians specially trained and authorized by J. Morita.

^{*} For repair or other types of service contact your local dealer or the J. Morita Corp.

Appendix- Electromagnetic declaration

$Guidance\ and\ manufacturer's\ declaration-electromagnetic\ emissions$

The X550 Type EX-1 and EX-2 are intended for use in the electromagnetic environment specified below. The customer or the user of the X550 Type EX-1 and EX-2 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1 Class B	The X550 Type EX-1 and EX-2 use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Group 1 Class B	The X550 Type EX-1 and EX-2 are suitable for use in all establishments, When installed in a specified shielded location, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic
*Harmonic emissions IEC61000-3-2	Class A	purposes.
*Voltage fluctuations/flicker emissions IEC 61000-3-3	Clause 5	

 $[\]ast$ These are only for the X550 Type EX-2.

Guidance and manufacturer's declaration - electromagnetic immunity

The X550 Type EX-1 and EX-2 are intended for use in the electromagnetic environment specified below. The customer or the user of the X550 Type EX-1 and EX-2 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance	
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8kV air	±2, ±4, ±6 kV contact ±2, ±4, ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.	
Electrical fast transients/bursts IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.	
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	$\begin{array}{l} \pm 0.5, \pm 1 \text{ kV} \\ \text{line(s) to line(s)} \\ \pm 0.5, \pm 1, \pm 2 \text{ kV} \\ \text{line(s) to earth} \end{array}$	Mains power quality should be that of a typical commercial or hospital environment.	
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5% U _T (>95% dip in U _T) for 0.5 cycle 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (>95% dip in U _T) for 5 sec	<5% U _T (>95% dip in U _T) for 0.5 cycle 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <0% U _T (>95% dip in U _T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If user of the X550 Type EX-1 and EX-2 requires continued operation during power mains interruptions, it is recommended that the X550 Type EX-1 and EX-2 be powered from an uninterruptible power supply or a battery.	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3.15 A/m	Power frequency magnetic field should be at levels characteristic of a typical location in a typical commercial or hospital environment.	
Note U_T is the a.c. mains voltage prior to application of the test level.				

Guidance and manufacturer's declaration - electromagnetic immunity

The X550 Type EX-1 and EX-2 are intended for use in the electromagnetic environment specified below. The customer or the user of the X550 Type EX-1 and EX-2 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the X550 Type EX-1 and EX-2, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance
			Recommended separation distance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3.15 V 3.6 V/m	$d = 1.2 \sqrt{P}$ $d = 1.0 \sqrt{P}$ $d = 2.0 \sqrt{P}$ $d = 2.0 \sqrt{P}$ $d = 2.0 \sqrt{P}$ $800 \text{MHz to } 2.5 \text{ GHz}$
			Where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected be absorption and reflection from structures, objects and people.

b Over the frequency range 150 kHz to 80MHz, field strengths should be less than 3 V/m.

^a Field strengths from fixed transmitters, such as base stations for ratio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicated theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the X550 Type EX-1 and EX-2 are used exceeds the applicable RF compliance level above, the X550 Type EX-1 and EX-2 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting of relocating the X550 Type EX-1 and EX-2.

Recommended separation distances between portable and mobile RF communications equipment and the X550 Type EX-1 and EX-2.

The X550 Type EX-1 and EX-2 are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the X550 Type EX-1 and EX-2 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the X550 Type EX-1 and EX-2 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output	Separation distance according to frequency of transmitter m			
W	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz	
	$d = 1.2 \sqrt{P}$	$d = 1.2 \sqrt{P}$	$d = 2.0\sqrt{P}$	
0.01	0.12	0.10	0.2	
0.1	0.38	0.32	0.63	
1	1.20	1.00	2.0	
10	3.79	3.16	6.32	
100	12	10	20	

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.



J. MORITA MFG. CORP.

680 Higashihama Minami-cho, Fushimi-ku, Kyoto, 612-8533 Japan www.jmorita-mfg.com

Distributors

J. MORITA CORPORATION

Tokyo Office: 11-15, 2-Chome Ueno, Taito-ku, Tokyo, 110-8513 Japan Osaka Office: 33-18, 3-Chome Tarumi-cho, Suita, Osaka, 564-8650 Japan

J. MORITA USA, Inc.

9 Mason Irvine, CA 92618 U.S.A.

TEL:+1-949-581-9600 FAX: +1-949-465-1095

J. MORITA EUROPE GMBH

Justus-Von-Liebig-Strasse 27A, D-63128 Dietzenbach Germany

TEL: +49-6074-836-0 FAX: +49-6074-836-299

Siamdent Co., Ltd.

71/10 Bangpakong Industrial Park I, Bangna-Trad, KM. 52, Bangpakong

Chachuengsao 24130, Thailand

TEL: +66-38-57-3042 FAX: +66-38-57-3043

J. MORITA CORPORATION Australia and New Zealand

Suite 2.05, Aero 247 Coward Street Mascot NSW 2020, Australia

TEL: +61-2-9667-3555 FAX: +61-2-9667-3577

EU Authorized Representative under the European Directive 93/42/EEC



MEDICAL TECHNOLOGY PROMEDT CONSULTING GMBH

Altenhofstraße 80, 66386 St. Ingbert, Germany

TEL: +49-6894-581020 FAX: +49-6894-581021

The authority granted to the authorized representative, MEDICAL TECHNOLOGY PROMEDT Consulting GmbH, by J. Morita Mfg. Corp. is solely limited to the work of the authorized representative with the requirements of the European Directive 93/42/EEC for product registration and incident report.