

One Year Manufacturer's Limited Warranty

WARRANTY TERMS:

1. This unit is covered for one year from the date of purchase for manufacturer's defects only.
2. The following are not covered:
 - A) damage due to improper handling or use
 - B) normal wear and tear
 - C) accessories: power transformers, cords, etc.
3. If repairs are needed after the warranty period expires or in case of improper use or handling of the item, the customer will be responsible for paying the appropriate parts and labor costs.
4. In case repair is necessary in the future, please record the date of purchase below for your own records:

Date of Purchase: _____



REIZEN, Inc.
P.O BOX 3312
Farmingdale, NY 11735
www.ReizenUSA.com

Reizen Mouse Electronic Magnifier Wireless Version: USB (#807566)



User's Manual

Congratulations on your purchase of the Reizen Mouse Electronic Magnifier. This magnifier was designed to enable those with low vision to read comfortably and easily. By simply rolling the mouse magnifier over reading material, you'll be able to view text magnified on your PC.

With adjustable 2x to 32x magnification, variable viewing modes and an Image Lock (Freeze Frame) feature, you'll have the flexibility and functionality you need to make reading easier than ever. With the USB version you can also save screen images and share them with family and friends via email.

Important Safety Precautions

To keep the unit in good condition and to ensure safe usage, please refer to the following guidelines:

- Read the instructions in this user's guide before starting to use the unit.
- Do not attempt to repair the unit or open its case. This product contains no user serviceable parts.
- Keep the unit away from wet areas, rain, and other sources of liquids or high humidity. Do not submerge the unit.
- Always operate the unit in a place where the temperature is between 10° and 40° C (50° and 104° F).
- Always store where the temperature is between -20° and 65° C (-4° and 149° F).
- Always use the unit on a stable, level, hard surface such as a desk or table.
- Never force connectors into the unit's ports. All appropriate connectors should fit easily into the ports. If the connector does not fit, turn it around, and try again.
- Use only the power cable that came with your unit. Other adapters may damage unit.
- When cleaning the unit, make sure it is unplugged. To clean it, lightly dampen a soft, clean cloth with water or mild detergent.
- Do not bend or twist the cable excessively as this could cause fatigue and lead to unit failure.
- Do not leave the power cable or product cable where someone might trip over it.

Video Sampling Card Functions:

- A. Date & Time display
- B. Camera timer: take picture every 5 seconds (Refer to item D for save direction path)
- C. Photo frame and save: take picture (JPG or BMP format)
- D. Document save direction path: Set up camera timer and document path
- E. About software version information and developer information
- F. Advanced settings: set up video parameters, digital camera focus and stunts
- G. Video image resolution settings
- H. Video record button
- I. Video record settings: select coder and set up video
- J. File save direction path
- K. Video display interface

Note: In order to get the best image quality, please refer to following instructions:

1. Select desired image resolution (refer to G.)
2. Save documents to a direction path that is easy to find (Refer to D)
3. Adjust mouse camera to get the best image display

Technical Specifications:

Item	Description
Display Resolution	Multiple Options
Display Mode	Full Color, B/W, W/B, High Contrast
Image Sensor	CMOS, VGA resolution
Power Output	5V DC; 250ma (Max)
Battery Life	3-4 hours
Battery Type	Lithium
Weight	100 grams adapter excluded
Dimensions: Mouse	5" x 2.75" x 2.50"
Dimensions: Video Card	4" x 1.5" x 0.625"
Magnification Range:	32 times depending on screen size

Mouse Magnifier Button Functions:

1. "M": Mode selection button (Full Color - Black on White - White on Black - High Contrast)
2. "+": Zoom in
3. "-": Zoom out
4. "**": Image Lock (Freeze Frame)
5. Camera On/Off Switch

Installing & Using the Mouse Electronic Magnifier (Wireless USB Version for Computer)

1. First, charge the wireless mouse magnifier using the supplied cable (Disconnect AC cable near yellow Video Out plug and connect to jack on mouse.) Mouse fully charges in 3-4 hours. After mouse is charged, re-connect cable.
2. Unplug any existing camera devices and insert the CD in your computer, then Launch "smartvideo.exe"
3. Click "Next" as necessary and then "Finish"
4. Simply plug your video sampling card into any USB 2.0 port on your computer. The hardware will install automatically and a notification will appear confirming your device is ready to use.
5. To use the magnifier, launch the multi video software (it will be in your programs folder under smartvideo direction path), then move the Mouse Magnifier over the material you wish to read and it will be displayed on the main menu of the software main interface window.
6. Use "Plus" & "Minus" keys to adjust the magnifying power (2x to 32x)
7. Use "M" to choose Color, Black on White, White on Black or High Contrast mode
8. Use the "**" (Image Lock) button to 'freeze' an image for extended viewing

What's in the Box?

- Reizen Mouse Electronic Magnifier
- Video sampling card (Receiver w/USB connector)
- Software installation disc
- Charger & Adapter
- Drawstring Carry/Storage Pouch
- User's Manual

Main Features of Mouse Electronic Magnifier:

- 8 adjustable magnification levels (2x-32x) to accommodate varying user needs
- Multiple viewing modes: full color, black on white, white on black, high contrast
- Built-in LED lighting system allows for operation in different lighting environments
- Ergonomically designed and easy to use
- Stable, reliable wireless receiver and transmitter ensure strong signal quality
- Image Lock/Freeze Frame: Freeze images on screen for extended viewing
- Record and replay video
- Share saved files with friends and family members via email

