Technical center: 5-7, avenue de Paris 94300 Vincennes - France

Tel.: + 33 (0)1 60 31 70 40 - Fax: + 33 (0)1 64 02 31 93 - info@parsys.com



User Manual

TELECARDIA

37-021 V4.0

Date of the first placing on the market of the product: 01/03/2007

CE₀₄₅₉





SUMMARY

1	KIT	CONTENTS	. 3
	1.1 C	COMPONENTS	. 3
		IST OF SYMBOLS USED	
	1.3 V	VARNINGS	. 5
2	EQI	JIPMENT DESCRIPTION	. 6
	2.1)VERVIEW	6
	_	LECTRODES POSITION	_
	2.2.		_
	2.2.	· F · · · · · · · · · · · · · · · · · ·	
		EVICE CONTROL BUTTONS	
		CORPARING	
		CG READING	
3		TALLATION OF TELECARDIALYS PC SOFTWARE	
4	BLU	JETOOTH PAIRING PROCEDURE OF THE TELECARDIA TO A PC	13
	4.1	WINDOWS XP OPERATING SYSTEM	14
	4.2	WINDOWS 7 OPERATING SYSTEM	
	4.2.		
	4.2. 4.3	2 Pairing of Telecardia with a Windows 7 PC	
	4.3.		-
	4.3.		
5	IISI	NG TELECARDIALYS PC SOFTWARE	47
		G CAPTURE	
6			
		-LEAD PERIPHERAL PATIENT CABLE MODEL	
		-LEAD PERIPHERAL PATIENT CABLE MODEL	
7		DATING THE PC SOFTWARE / DEVICE FIRMWARE	
8	TEC	CHNICAL FEATURES	61
	8.1	GENERAL FEATURES	61
		MECHANICAL FEATURES	
		LECTRICAL SPECIFICATIONS	
		UNCTIONAL FEATURES	
_			
9		CTROMAGNETIC EMISSIONS	_
10) MA	NTENANCE	67
	10.1	CLEANING / DISINFECTION	
	10.2	PROCEDURE IN THE EVENT THAT THE INSTRUMENT IS DROPPED OR FALLS	
	10.3	RESET PROCEDURE	-
	10.4 10.5	METROLOGICAL CHECKS	
	10.5	PRODUCT SCRAP TREATMENT	
11		ERS' TRAINING	
12		JLT OR FAILURE	
	12.1	THE TELECARDIA ECG DOES NOT SEND ANY TRACE	
	12.2 12.3	THE ECG IS SUBJECT TO INTERFERENCE	
	12.3	THE TELECARDIA ECG IS NOT CHARGING	
	12.5	THE TELECARDIA ECG HAS LOST IS TIME STAMPING	
13	S DVI	RSYS TELEMEDICINE WARRANTY AND AFTER-SALES SERVICE	72
-	, 1,741	TO I TEENIEDIONIE TIANNAMI I AND ALLENIOAELO GENTIGE	. ~

1 Kit contents

We thank you for having chosen PARSYS Telemedicine equipment. We hope that it will bring you full satisfaction in your daily electrocardiography practice.

1.1 Components



Firstly, please **check the contents** of your **Telecardia** Kit (Ref: 100-004). It must contain the following items:

′1.1 ′1.2
1 /11
,

If any of the items are missing, please contact us at the address specified in this manual.

10/15/2015 Version Page 3 of 72



User Manual Telecardia

Ver 4.0

1.2 List of symbols used



Warning: consult attached documents



Please read the notice



BF type applied part



Do not dispose of in the bin but through the appropriate recycling channels.



Non-ionizing radiation



Product manufactured by PARSYS Telemedecine



Device climatic range of use or storage



Serial number



Medical device European compliance

10/15/2015 Version Page 4 of 72

1.3 Warnings

PARSYS Telemedicine draws the attention of the user on the following points regarding:

- ⚠ The device **must not be used simultaneously** with an electric surgery tool.
- ⚠ The device **is not intended** for cardiac direct use (reserved to CF class devices).
- ⚠ The ECG capture **must not be performed** when the device is charging on its charging base.
- ⚠ When the ECG is capturing, electrodes unplugged on the patient **must not come into contact** with the ground or other exposed metal parts.
- ⚠ When the ECG is capturing, the device's differentiator box **must be in contact** with the patient.
- ⚠ In case of a simultaneous use of the device and another, **it should take care** that the sum of the patient leakage currents to the ground doesn't exceed safety limits.
- ⚠ In case of a simultaneous use of the device with a defibrillator, **do not place** the defibrillator's electrodes into contact with the device's electrodes.

10/15/2015 Version Page 5 of 72

2 **Equipment description**

2.1 Overview

Telecardia represents a major advance in remote cardiac monitoring and portable medical equipment. It enables health professionals to rapidly capture and send a patient's ECG under any circumstances, without the constraints of a conventional device.


The trace captured is a digital, simultaneous 12-lead 12-tracing ECG, in keeping with standard electrocardiographs usually used. It can be sent automatically to a care watch centre or be viewed on a PC next to the patient.

Telecardia has CE marking (CE0459). This is an ECG device intended for temporary use, class IIa.

2.2 Electrodes position

Telecardia is an electrocardiograph:

- simultaneous 12-lead 12-channel,
- with fixed thoracic electrodes,
- with dry electrodes,
- no consumables.



2.2.1 Precordial leads

The shape of the unit enables capture of the **6 precordial leads** collected thanks to the 2 articulated arms. These thoracic electrodes are dry and fixed, from V1 to V6.

The Neutral, usually captured on the right leg, is captured in the thorax via the fixed electrode on the housing (see figure 1.).

10/15/2015 Version Page 6 of 72



2.2.2 Peripheral leads

2.2.2.1 2-lead peripheral Patient cable model

With the 2-lead model, Telecardia contains:

• 1 rolling wire electrode, located under the device left arm for the **AvF** lead.









3. Patient cable connecting

■ 1 cable to plug on the Telecardia box right-side, capturing AvR and AvL with expandable bracelets or limb clips dry electrodes (see figure 4 and 5).



4. Cordon Patient AvR & AvL Patient cable




5. Expandable bracelets or limb clips dry electrodes

10/15/2015 Version Page 7 of 72




2.2.2.2 3-lead peripheral Patient cable model

With the 3-lead model, Telecardia contains a Patient cable to plug into the ECG box right-side. It enables the capture of **AvR**, **AvL** and **AvF** with limb clips dry electrodes (see figure 6.).



6. Peripherals cable AvR, AvL and AvF - Limb clips electrodes



7. Peripheral cable connection
The insertion mark, a red point on the cord, must be facing the device's top

10/15/2015 Version Page 8 of 72




2.3 Device control buttons

To control the machine and start all the functions, Telecardia has 2 buttons:

- 1 button used to release the arms and switch the machine ON / OFF, located on the front of the machine and marked with a green circle crossed by a black line.
- 1 button used to access the machine's advanced functions, located under the top cover and marked "OUT".

In addition, combinations of these buttons are used to select the requested functions.



4. Telecardia buttons

10/15/2015 Version Page 9 of 72



2.4 Charging base

Telecardia is equipped with a separate power supply provided by a rechargeable Lithium/Polymer battery, allowing approximately 50 traces per charge. This battery is charged using the charging base, via 2 metal contacts integrated in the back of the machine unit.

The charging base charges the machine completely in approximately **2.5 hours** using a 220V/50Hz source, via the mains adapter supplied.



ALWAYS USE THE MAINS ADAPTER SUPPLIED TO POWER THE BASE; OTHERWISE, IRREVERSIBLE DAMAGE NOT COVERED BY THE MACHINE WARRANTY IS LIABLE TO OCCUR.

To charge Telecardia:


- Connect your charging base to 220V power supply
- Put Telecardia upon its charging base:





CHECK LIGHT INDICATORS ON THE CHARGING BASE'S SIDE TO CONFIRM THE GOOD CONTACT OF THE DEVICE WITH CHARGING CONTACTS.

Light indicators located on the side of the base indicate the machine charging status, and whether the machine is in contact with the support.



The ON/OFF indicator **flashes BLUE** to indicate that the base **power** is **on**.

The CHARGING INDICATOR does not light up when Telecardia is not correctly inserted in its slot, without contact with the charging contacts, or missing from the slot.

The CHARGING INDICATOR remains **RED** as long as Telecardia is **charging**.

The CHARGING INDICATOR is **GREEN** when Telecardia is **charged**.



Not charging Telecardia for more than 1 year means losing the instrument's internal real-time clock: we recommend to regularly charge Telecardia battery.

10/15/2015 Version Page 10 of 72


2.5 ECG reading

Viewing and/or transmission of the trace on a PC

Following capture, the data is sent to the PC automatically via the Bluetooth radio link which has a free field range of approximately 100 meters. This requires the presence of a Bluetooth USB adapter, supplied in the kit, connected to the PC used to operate Telecardia.

On the operating PC, the Telecardialys software (See the User Manual Telecardialys Software - 37-063 V1.0) is used to:

- dialogue with the equipment via Bluetooth,
- display the captured ECGs,
- print the ECGs in several report modes,
- send the captured ECGs by e-mail,
- archive and view the captured ECGs.



10/15/2015 Version Page 11 of 72

3 Installation of Telecardialys PC Software

See the User Manual Telecardialys Software - 37-063 V1.0

10/15/2015 Version Page 12 of 72

4 Bluetooth pairing procedure of the Telecardia to a PC

Windows XP	Window	s Vista	Windows 7	and above
External adapter	External adapter	Internal module	External adapter	Internal module
Installation CD §xxx	Installation CD §xxx	Nothing to do	Insert adapter + Internet	Nothing to do



The pairing procedure of the Telecardia ECG with a PC running Windows 7 and above is specific: it is essential to follow the appropriate procedure for the operating system of your PC.

10/15/2015 Version Page 13 of 72



4.1 Windows XP Operating System

Pairing of Telecardia with a Windows XP PC



The Bluetooth USB adapter is pre-associated by us with your Telecardia. Use the following procedure in the following cases:

- a. you change your Bluetooth USB adapter,
- b. you do not use the supplied adapter to the benefit of an internal Bluetooth module,
- c. you use another Telecardia with the same adapter (potential After-Sales return).

Insert the USB Bluetooth adapter into an USB port on the PC, then run the installation CD.


If the CD does not run automatically, please run it manually in Windows Explorer.




Do nothing until this window appears. Click on "Install Bluetooth© Device".

10/15/2015 Version Page 14 of 72






Click on "Next >".




Accept the agreement, then click on "Next >".

10/15/2015 Version Page 15 of 72






Click on "Next >".




Click on "Install".

10/15/2015 Version Page 16 of 72





Wait until the installation completes.




Click on "Finish".

10/15/2015 Version Page 17 of 72





Make a "right click" on the Bluetooth icon at the bottom right of the screen.



Click on "Quick Connect", then, on "Bluetooth Serial Port", then, on "Find Devices...".

10/15/2015 Version Page 18 of 72





Open the drop down menu: "Show devices" and select "Show all devices".




Open the Telecardia housing.

Note the ID of the used device: Number starting with a T

10/15/2015 Version Page 19 of 72





Make the connection.

Simultaneously press the **ON** and **OUT** buttons until the device's ignition.



Select the device.

In the list, select your Telecardia and click on "Next >".

If necessary, click on "Search Again".

10/15/2015 Version Page 20 of 72





Pairing.

Pairing starts with Telecardia.




Enter the pairing passkey:

The installer asks for a **passkey**: Type the code: **02101124** and click on "**Next** >".

10/15/2015 Version Page 21 of 72





Click on "Finish".



Create shortcuts.

Click on "Do not display this message again", then on "OK".




To finalize the installation of the Bluetooth device, make a "right click" on the Bluetooth icon in the bottom right of the screen. Then, click on "Bluetooth Configuration".

10/15/2015 Version Page 22 of 72






Go to the "Local Services" tab. Click on "Bluetooth Serial Port". Then, click on "Properties...".




Check "Start up automatically" and uncheck "Secure connection". Click on "Apply", then on "OK".

10/15/2015 Version Page 23 of 72





Go to the "Client Applications" tab. Click on "Bluetooth Serial Port". Then, click on "Properties…".



Uncheck "Secure connection".

Click on "Apply", then on "OK".

10/15/2015 Version Page 24 of 72

4.2 Windows 7 Operating System

4.2.1 Requirements

4.1.1.1 Check if the PC (fixed, portable PC, ultraportable, Tablet PC...) to pair with the ECG includes a preinstalled Bluetooth internal module.

Note: In the case of a fixed PC, there is a good chance that there is no internal module; but, mobile computers often have a preinstalled Bluetooth internal module.

- 4.1.1.2 If no internal module is preinstalled, use the Bluetooth USB adapter supplied with your kit.
 - It is fundamental that your computer has an active Internet connection.
 - Just plug your Bluetooth USB adapter on a USB port, driver installation is automatic via the Internet.



In this case, your computer must be connected to the Internet to install the Bluetooth USB adapter (driver updating).

4.1.1.3 If your computer is equipped with a preinstalled Bluetooth internal module, check the origin of its stack.

<u>Microsoft / WidComm stack</u>: Pairing is automatic between Telecardia and the computer, you have no other

manipulation to perform

Another provider stack: We recommend you to communicate the name of the stack provider to the

PARSYS Telemedicine After-Sales Service before any installation.



In the latter case, we recommend not to try the installation before obtaining confirmation procedure by the PARSYS Telemedicine After-Sales Service.

10/15/2015 Version Page 25 of 72



4.2.2 Pairing of Telecardia with a Windows 7 PC

Insert the USB Bluetooth adapter into an USB port on the PC, then run the installation CD.

If the CD does not run automatically, please run it manually in Windows Explorer.




Do nothing until this window appears. Click on "Install Bluetooth© Device".




10/15/2015 Version Page 26 of 72



Click on "Next >".




Accept the agreement, then click on "Next >".




Click on "Next >".

10/15/2015 Version Page 27 of 72






Click on "Install".



Wait until the installation completes.


10/15/2015 Version Page 28 of 72





Click on "Finish".


Open the "Program Files" folder (32-bit processors) or "Program Files (x86)" (64-bit processors).



Then, open the "Parsys" folder and start the "TCP.EXE" software.

10/15/2015 Version Page 29 of 72





When the software is opened, set the device "Call (Client)" mode.




Open the Telecardia housing.

Note the ID of the used device: Number starting with a T


10/15/2015 Version Page 30 of 72





Make the connection.


Simultaneously press the **ON** and **OUT** buttons until the device's ignition.




The different connectable Bluetooth devices are recognized by the system. You can refine the search by type of equipment (Telecardia, Cardialys, etc.).

10/15/2015 Version Page 31 of 72






If the search is unsuccessful you can click on "**Refresh...**", which will start a new search.




Once the desired device is visible, click its icon (although check the serial number). Then, click on "**Choose**".

10/15/2015 Version Page 32 of 72






The interface window with the medical device is then visible.




Click on "Remote Mac Address".
Then, click on "Local Device".
The PC MAC address will be registered.

10/15/2015 Version Page 33 of 72





Click on "Remote Mac Address".
Then, click on "Write" to save the new MAC address.




In the field "Max Number of ECG" replace the number by "1", if you only want to do an ECG at a time.

This is preferable to avoid any mistake in ECG / Patient assignments.


10/15/2015 Version Page 34 of 72





Click on "Configuration".


Then, click on "Write" to save all the changes.



Click on "Disconnect".

10/15/2015 Version Page 35 of 72





Back to the start window of the "**TCP**" software. Close the window to exit.

10/15/2015 Version Page 36 of 72

Ver 4.0

4.3 Windows 8 / 8.1 Operating System

4.3.1 Requirements

4.3.1.1 Check if the PC (fixed, portable PC, ultraportable, Tablet PC...) to pair with the ECG includes a preinstalled Bluetooth internal module.

Note: In the case of a fixed PC, there is a good chance that there is no internal module; but, mobile computers often have a preinstalled Bluetooth internal module.

- 4.3.1.2 If no internal module is preinstalled, use the Bluetooth USB adapter supplied with your kit.
 - It is fundamental that your computer has an active Internet connection.
 - Just plug your Bluetooth USB adapter on a USB port, driver installation is automatic via the Internet.



In this case, your computer must be connected to the Internet to install the Bluetooth USB adapter (driver updating).

4.3.1.3 If your computer is equipped with a preinstalled Bluetooth internal module, check the origin of its stack.

<u>Microsoft / WidComm stack</u>: Pairing is automatic between Telecardia and the computer, you have no other

manipulation to perform

Another provider stack: We recommend you to communicate the name of the stack provider to the

PARSYS Telemedicine After-Sales Service before any installation.



In the latter case, we recommend not to try the installation before obtaining confirmation procedure by the PARSYS Telemedicine After-Sales Service.

10/15/2015 Version Page 37 of 72




4.3.2 Pairing of Telecardia with a Windows 8 / 8.1 PC

Insert the USB Bluetooth adapter into an USB port on the PC, then run the installation CD.

If the CD does not run automatically, please run it manually in Windows Explorer.



Do nothing until this window appears. Click on "Install Bluetooth@ Device".




Click on "Agree and Install Bluetooth Software".

10/15/2015 Version Page 38 of 72



Wait until the installation completes.




Click on "Finish".


10/15/2015 Version Page 39 of 72



Open the "Program Files" folder (32-bit processors) or "Program Files (x86)" (64-bit processors).



Then, open the "Parsys" folder and start the "TCP.EXE" software.



When the software is opened, set the device "Call (Client)" mode.


10/15/2015 Version Page 40 of 72





Open the Telecardia housing.

Note the ID of the used device: Number starting with a T




Make the connection.


Simultaneously press the ${\bf ON}$ and ${\bf OUT}$ buttons until the device's ignition.

10/15/2015 Version Page 41 of 72






The different connectable Bluetooth devices are recognized by the system. You can refine the search by type of equipment (Telecardia, Cardialys, etc.).




If the search is unsuccessful you can click on "**Refresh...**", which will start a new search.

10/15/2015 Version Page 42 of 72






Once the desired device is visible, click its icon (although check the serial number). Then, click on "**Choose**".




The interface window with the medical device is then visible.

10/15/2015 Version Page 43 of 72






Click on "Remote Mac Address".
Then, click on "Local Device".
The PC MAC address will be registered.



Click on "Remote Mac Address". Then, click on "Write" to save the new MAC address.


10/15/2015 Version Page 44 of 72





In the field "Max Number of ECG" replace the number by "1", if you only want to do an ECG at a time.

This is preferable to avoid any mistake in ECG / Patient assignments.




Click on "Configuration".


Then, click on "Write" to save all the changes.

10/15/2015 Version Page 45 of 72





Click on "Disconnect".



Back to the start window of the "**TCP**" software. Close the window to exit.

10/15/2015 Version Page 46 of 72

User Manual Telecardia Ver 4.0

5 <u>Using Telecardialys PC Software</u>

See the User Manual Telecardialys Software - 37-063 V1.0

10/15/2015 Version Page 47 of 72



6 ECG capture

If you have the Telecardia ECG, please follow the following procedures to capture a simultaneous 12-lead ECG trace on the patient:

- Press the "Start Telecardia (ECG)" button on the software, the message "ECG waiting" appears.
- You can capture the ECG, the trace will be automatically transmitted to the PC.


6.12-lead peripheral Patient cable model


Capture ECG NB: 1 For the sake of clarity, the illustrations show a figure with no shirt. However, the subject does not have to take his/her shirt off. SHIRT OPEN - SLEEVES UP is satisfactory. 2 The ECG can be captured with the subject SITTING or LYING DOWN NEVER STANDING UP.



10/15/2015 Version Page 48 of 72









10/15/2015 Version Page 49 of 72









10/15/2015 Version Page 50 of 72









10/15/2015 Version Page 51 of 72









10/15/2015 Version Page 52 of 72







10/15/2015 Version Page 53 of 72



6.23-lead peripheral Patient cable model

Capture ECG

NB:

- For the sake of clarity, the illustrations show a figure with no shirt.
 However, the subject does not have to take his/her shirt off.
 SHIRT OPEN SLEEVES UP is satisfactory.
- The ECG can be captured with the subject SITTING or LYING DOWN NEVER STANDING UP.

Capture ECG





Connect the ECG cables to the limb clips following the colour codes.

Dampen the limb clips with the water spray.

10/15/2015 Version Page 54 of 72








10/15/2015 Version Page 55 of 72









10/15/2015 Version Page 56 of 72









10/15/2015 Version Page 57 of 72








10/15/2015 Version Page 58 of 72





10/15/2015 Version Page 59 of 72

Ver 4.0

7 Updating the PC software / device firmware

We shall contact you periodically by e-mail concerning download of updated software / firmware versions.

The updates are performed during maintenance.

10/15/2015 Version Page 60 of 72



8 Technical features

8.1 General features

Portable, event-driven, 12-lead temporary electrocardiograph.

Radio data transmission (Bluetooth).

Lithium polymer rechargeable battery power supply.

8.2 Mechanical features

Material: two half-body unit + one hinged hatch made

of flame-retardant ABS

Color: RAL 9002 white with marking

Surface condition: 26 arbours

Electrode on unit:

chest on lower section → VN - V1 - V2 - V3 - V4 - V5 - V6 (7/10μ chromium-plated brass)

Electrodes with cable:

VR, VL and VF peripheral limb clips cable

Connectors:

o "female cylinder connectors with collar" under bushing cover for electrocardiograph connection

- analogue trace output

Protection index:

Folded overall dimensions:
210 mm x 100 mm x 50 mm

■ Weight: < 550 g

8.3 Electrical specifications

Device intended for temporary use.

Device not protected against defibrillation shocks.

Device with internal electricity supply: rechargeable lithium polymer batteries

Battery life: 50 peri-critical ECG captures and Bluetooth

transmission in one battery charge

Battery: Built-in Lithium Polymer battery

3,7V/1300mAh

Battery life: 3 years life cycle

Automatic power off after transmission.

Internal protection using 500mA resettable fuse.

Type B class 2A machine, diagnostic function.

Measurement with internal electrodes, 9 simultaneous leads: VR, VL, VF, V1 - V2 - V3 - V4 - V5 -

V6 with reference to VN.

Reproduction of 12 leads:
 D1, D2, D3, VR, VL, VF, V1, V2, V3, V4, V5,

V6

Input impedance: 10 MΩ.

Bandwidth: 0.16Hz to 70Hz at -3db

(Adjustable to 30Hz by configuration)

Measurement linearity: 0,01%
 Measurement sensitivity: 2264 pt/mV
 Sampling frequency: 271 pt/s

 Backup of 6 (peri-critical) 15-second ECGs (rotating or fixed list backup management with saving maintained even without stack).

Saving of identification and configuration data in encrypted form.

 Real-time clock for ECG dating (autonomy of 1 year from full charging level under normal storage conditions).

Displayed and/or voice task procedure or action messages.

10/15/2015 Version Page 61 of 72



Ver 4.0

task procedure message display

PARSYS

Operating mode setting (computer configuration):

o maximum number of ECGs

2-line, 16-character screen:

- voice prompt
- o rotating stack/locking
- connection type
- Bluetooth type interface for computer link and transmission base link:
 - Bluetooth version 2.0 (EDR)
 - o Class 1 (100 mW)

Average transmission speed: 115.2 Kb/s (depends on medium used)

Operating conditions:

Temperature: +10°C to +35°C
 Relative humidity level: 30% to 75%
 Atmospheric pressure: 700 mb to 1060 mb.

Storage and transport conditions:

Temperature: -10°C to +50°C
 Relative humidity level: 10% to 85%
 Atmospheric pressure: 500 mb to 1060 mb

8.4 Functional features

Front panel:

- start-up button,
- configuration button (under housing cover),
- blue/green LED indicating:
 - o steady green machine power on
 - o steady blue Bluetooth transmission to PC or base
- 1 LCD 2 x 16-character line screen without backlighting.

Curved lower section:

3 chromium-plated brass convex round electrodes

Upper section:

Hatch opening outwards and used to support the unit

Upper section under hatch:

- 10 x 4 mm "female cylindrical" connectors with attachment collar between the two bodies
- 1 signal transfer control button to electrocardiograph (analogue output), or configuration

Rear panel (back):

- Instruction label
- 2 charging contacts


Edges:

- 1 folding articulated arm on either side, supporting electrodes V1, V2 and V5, V6.
- 2-lead model:
 - o Under the left arm, 1 output of the wire VF electrode (iliac), with automatic winding system.
 - o Under the right arm, 1 connecting socket for VR and VL cable.
- 3-lead model:
 - o Under the right arm, 1 connecting socket for VR, VL and VF limb clips cable.

10/15/2015 Version Page 62 of 72



8.5 Functional flow chart



10/15/2015 Version Page 63 of 72



9 Electromagnetic emissions

Table 1: Directives and MANUFACTURER'S statement – ELECTROMAGNETIC EMISSIONS

For all DEVICES and EM SYSTEMS

Directives and manufacturer's statement - electromagnetic emissions

The device is made to be used in the electromagnetic environment specified below. The client or user of the device must ensure that the device is used in such environment.

Emission trials	Conformity	Electromagnetic environment – directives	
RF Emissions CISPR 11	Group 1	The device uses RF energy only through his internal functions. As a consequence, its RF emissions are very weak and are not likely to cause interferences in a nearby electronic device.	
RF Emissions CISPR 11	Class B	The charging newer here (actalogue product) must	
Harmonic Emissions CEI 61000-3-2	Non applicable	The charging power box (catalogue product) must solely be used with the device in sleep mode, disconnected from the patient.	
Pulse fluctuations/ Flicker CEI 61000-3-3	Non applicable		

Table 2 – Directives and MANUFACTURER'S statement – electromagnetic IMMUNITY

For all DEVICES and EM SYSTEMS

Directives and manufacturer's statement – electromagnetic immunity

The device is made to be used in the electromagnetic environment specified below. The client or user of the device must ensure that the device is used in such environment.

	1		
Immunity trials	Trial level CEI 60601	Conformity level	Electromagnetic environment – directives
Electrostatic shocks (ES)	+/- 6 kV in contact +/- 8 kV in the air	Conforms	Floors must be made of wood, concrete or ceramics. If the floors are topped with
CEI 61000-4-2			synthetic materials, it is important that the relative humidity remains at 30% at least.
Magnetic field at the rate of the electric network (50/60 Hz)	3 A/m	Conforms	Magnetic fields at the rate of the electric network must have the distinctive levels of a
CEI 61000-4-8			representative place located in a commercial or hospital environment.

10/15/2015 Version Page 64 of 72



Table 3 - Directives and MANUFACTURER'S statement - electromagnetic IMMUNITY

For all DEVICES and EM SYSTEMS

Directives and manufacturer's statement - electromagnetic emissions

The device is made to be used in the electromagnetic environment specified below. The client or user of the device must ensure that the device is used in such environment.

Emission trials	Trial level CEI 60601	Conformity level	Electromagnetic environment – directives
			The portable and mobile RF communicating devices should be not used nearer of any part of the medical devices, including cables, than the recommended separation distance, calculated from the transmitter frequency applicable equation.
			Recommended separation distance
			d= 1,2√P
RF conducted disturbances CEI 61000-4-6	3 Veff from 150kHz to 80MHz	3 Veff	d= 1,2 √ P from 80MHz to 800MHz d= 2,3 √ P from 800MHz to 2,5GHz
RF radiated disturbances CEI 61000-4-3	3 V/m from 80 to 2,5GHz	3 V/m	where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m) of separation. It should that field strengths from fixed RF transmitters, as determined by an electromagnetic survey on site ^a , should be less than the compliance level in each frequency range ^b . Interferences may occur in the vicinity of equipment marked with the following symbol:

NOTE 1: to 80MHz and to 800MHz, the higher d range frequency applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- ^a Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM / FM radio and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, it should be considered to proceed to an electromagnetic survey on site. If the measured field strength in the location where the medical device is used, exceeds the applicable RF compliance level above, it is necessary to observe the behavior of the medical device, to check that the operation is normal. If abnormal performances are observed, additional measures may be necessary, such as reorienting or repositioning of the medical device on its operation site.
- b In the frequency range of 150kHz to 80MHz, it is appropriate that the field strengths are smaller than 3 V/m.

10/15/2015 Version Page 65 of 72



Recommended separation distances between portable and mobile RF communicating devices and the medical device.

The medical device is intended for use in an electromagnetic environment where radiated RF disturbances are controlled. The customer or the user of the device can help to prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communicating device (transmitter), and the medical device, as recommended below, according to the maximum power emission of the communication device.

Transmitter maximum	Separation distance according to the transmitter frequency (m)			
assigned power output (W)	from 150kHz to 80MHz	from 150kHz to 80MHz	from 150kHz to 80MHz	
, ,	d= 1,2√P	d= 1,2√P	d= 2,3√P	
0,01	0,12	0,12	0,23	
0,1	0,38	0,38	0,73	
1	1,2	1,2	2,3	
10	3,8	3,8	7,3	
100	12	12	23	

For transmitters whose maximum assigned transmission power is not given above, the recommended separation distance **d** in meters (m) can be estimated by using the transmitter frequency applicable equation, where **P** is the characteristic of the transmitter maximum transmission power in watts (W) according to its manufacturer.

NOTE 1: to 80MHz and to 800MHz, the higher d range frequency applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

10/15/2015 Version Page 66 of 72



10 Maintenance

10.1 Cleaning / Disinfection

The Telecardia machine can be cleaned and disinfected by wiping the electrodes in contact with the skin and the plastic unit, except on the self-adhesive label at the rear of the machine, with cotton wool or a wipe soaked in a product such as STERANIOS or equivalent (2% glutaraldehyde solution).

Any other product would not ensure proper cleaning without damaging the Telecardia components.

10.2 Procedure in the event that the instrument is dropped or falls


Mechanical damage to the instrument necessarily implies discontinuation of use of the machine and return to the PARSYS Telemedicine After-Sales Service.

10.3 Reset procedure

If, when using the device, it blocks for more than 1 min:

- Use the reset cable or take a metal object (paper clip) and place it on the charging contacts for 1s,
- Place the equipment back on the charger and wait for the Charger indicator to come on,
- Establish radio connection with the PC in configuration mode in order to update the device's internal real-time clock:

Open the cover as shown on the photo: press and hold the "OUT" button and press on the GREEN ON/OFF button at the same time.



When the machine displays "PARSYS Santé", release both buttons.

The machine then indicates the time and date, before displaying "en cours" (in progress) and "connexion BT" (BT connection).

Restart the device.

10/15/2015 Version Page 67 of 72

User Manual Telecardia

Ver 4.0

10.4 Metrological checks

The device shall undergo metrological checks after 1 year maximum.

Checks should be carried out by the PARSYS Telemedicine technical department, it alone being able to guarantee maintenance of the device's metrological performance.

10.5 Discharged battery

We recommend not to let Telecardia discharged for more than 1 year to keep the device's internal real-time clock.

10.6 Product scrap treatment

According to the directive 2002/96/CE relating to DEEE electronic and electric devices, do not throw away in regular trash. Please bring waste to specialized collection points.



10/15/2015 Version Page 68 of 72

Ver 4.0

11 Users' Training

Users of the Telecardia ECG must:

- have kind of medical knowledge of physician, nurse or first aider,
- hold the associate degree,
 have been trained on this medical device.

10/15/2015 Version Page 69 of 72

Ver 4.0

12 Fault or failure

12.1 The Telecardia ECG does not send any trace

After capturing the ECG, no trace is received by paired PC: the Bluetooth® connection is most certainly lost.



We therefore advise you:

- NEVER ATTEMPT any direct manipulation on the PC;
- Contact the PARSYS Telemedicine After-Sales Service.

12.2 The ECG is subject to interference

Check that the capture protocol described in section 6. has been followed, in particular:

- that the correct color code has been followed for the limb clips: red on the right arm, yellow on the left arm and green on the left leg.
- that the electrodes and the skin have been humidified: using the spray.

The patient and the user must be as relaxed as possible in order to avoid any tense movements when capturing the ECG.

12.3 The Telecardia ECG is not charging

Check that the protocol described in section 2.4 has been followed, in particular:

- that the Telecardia is in the correct position on its charging base,
- that the red LED on the charging base is on during charging,
- that the charging base is correctly plugged into the mains.

Should you note a low battery level after 1 hour's charging, please contact the PARSYS Telemedicine After-Sales Service.

10/15/2015 Version Page 70 of 72



12.4 The Telecardia ECG remains blocked

If, when using the device, it blocks for more than 1 min:

- Use the reset cable or take a metal object (paper clip) and place it on the charging contacts for 1s,
- Place the equipment back on the charger and wait for the Charger indicator to come on,
- Establish radio connection with the PC in configuration mode in order to update the device's internal real-time clock:

Open the cover as shown on the photo: press and hold the "OUT" button and press on the GREEN ON/OFF button at the same time.



When the machine displays "PARSYS Santé", release both buttons.

The machine then indicates the time and date, before displaying "en cours" (in progress) and "connexion BT" (BT connection).

Restart the device.

12.5 The Telecardia ECG has lost is time stamping

If your Telecardia ECG has lost is time stamping:

• Establish radio connection with the PC in configuration mode in order to update the device's internal real-time clock:

Open the cover as shown on the photo: press and hold the "OUT" button and press on the GREEN ON/OFF button at the same time.



When the machine displays "PARSYS Santé", release both buttons.

The machine then indicates the time and date, before displaying "en cours" (in progress) and "connexion BT" (BT connection).

Restart the device.

10/15/2015 Version Page 71 of 72

Ver 4.0

13 PARSYS Telemedicine Warranty and After-Sales Service

The Customer agrees to comply with the Warranty Conditions listed on the Warranty Certificate accompanying the Equipment.

The Equipment is supplied with a one-year (1) warranty during which period the Customer will be able to exchange Equipment found to have a latent defect and subject to the Customer having informed PARSYS Telemedicine thereof in writing and in detail.

In the event of a functional fault or failure of the Equipment, the Customer may contact PARSYS Telemedicine:

- by sending an e-mail to support@parsys.com,
- or by calling the PARSYS Telemedicine After-Sales Service:
 - o Monday to Friday, 10 a.m. 12.30 p.m. and 2 p.m. 6p.m., except for bank holidays,
 - o at +33 (0)1 60 31 51 71
- or by sending a letter by recorded delivery with acknowledgement of receipt to:

PARSYS Telemedicine - After-Sales Service 5-7, avenue de Paris 94300 Vincennes - FRANCE

The PARSYS Telemedicine After-Sales Service identifies the nature of the fault or failure of the Equipment before carrying out any repairs or replacing the defective Equipment.

Beyond the warranty period, the Customer has the possibility of subscribing to a maintenance contract with the PARSYS Telemedicine sales department.

10/15/2015 Version Page 72 of 72