

# Enhanced $\beta$ -Galactosidase™

## Assay Kit (CPRG)



A Division of Gene Therapy Systems, Inc

Cat. No.	Content	Qty	Storage
A10100K	5X Lysis Buffer	55 ml	4°C
	Standard Dilution Buffer	55 ml	4°C
	Substrate Buffer	55 ml	4°C
	Stop Buffer	55 ml	4°C
	10X CPRG Substrate Stock (Chlorophenol red- $\beta$ -D-galactopyranoside)	5 x 1 ml	-20°C
	$\beta$ -gal enzyme standard, 40 units	100 $\mu$ l	-20°C

Shipping Condition	Blue ice.
--------------------	-----------

RELATED PRODUCTS	Catalog #
GenePORTER™ Transfection Reagent, 75 reactions	T201007
GenePORTER™ Transfection Reagent, 150 reactions	T201015
GenePORTER™ Transfection Reagent, 750 reactions	T201075
GenePORTER™ 2 Transfection Reagent, 75 reactions	T202007
GenePORTER™ 2 Transfection Reagent, 150 reactions	T202015
GenePORTER™ 2 Transfection Reagent, 750 reactions	T202075
gWiz™ $\beta$ -galactosidase Expression Vector, 25 $\mu$ g	P010200
$\beta$ -galactosidase Assay Kit (ONPG)	A10200K
X-Gal Staining Kit	A10300K

**INTRODUCTION** *LacZ* is a commonly used reporter gene in transfection experiments because the gene product,  $\beta$ -galactosidase, is very stable, resistant to proteolytic degradation, and easy to assay. Levels of active  $\beta$ -galactosidase expression can be quickly measured by its catalytic hydrolysis of Chlorophenol red- $\beta$ -D-galactopyranoside (CPRG) substrate to a dark red product. The assay kit provides all the required reagents, and offers a rapid, simple and sensitive method to quantify the enzyme expression in transfected cells (e.g., transfected with Genlantis' gWiz  $\beta$ -gal vector). The high sensitivity improves the measurement of  $\beta$ -gal activity when the reporter gene expression is low.

### USAGE

- Transfect cells with a plasmid expressing *LacZ* gene.
- Lyse cells using lysis buffer.
- Transfer the lysate to a fresh tube or a microtiter plate. Dilute the lysate if needed.
- Prepare a  $\beta$ -galactosidase standard curve with standard dilution buffer.
- Add the substrate and monitor the color development at 570-595 nm.
- Calculate the expression levels based on a standard curve.

### EXAMPLE PROTOCOLS

**NOTE:** Dilute 5X Lysis buffer to 1X with distilled deionized water before use. Unused 1X Lysis Buffer may be stored at 4°C for future use.

#### A. Protocol For Adherent Cells

1. Aspirate the growth medium from the culture dish 24-72 hours after transfection (use non-transfected cells for control). Optionally, wash cells 1 time with 1X PBS.
2. Add 1X Lysis Buffer to the culture dish. Use the following recommended volumes depending on your culture dish:

Type of Dish	1X Lysis Buffer Volume( $\mu$ l/well or dish)
96-well plate	50
24-well plate	250
12-well plate	500
6-well plate	1000
60 mm dish	2500
100 mm dish	5000

3. Incubate the dish 10-15 minutes at room temperature by swirling it slowly several times to ensure complete lysis. The culture dishes can be observed under a microscope to confirm that the cells are lysed completely.
4. Proceed to Section C. Below.

**NOTE:** A quick freeze/thaw cycle (freeze 1-2 hours at -20°C or -70°C then thaw at room temperature) of the dish can be performed to improve lysis. Keep dish frozen at -70°C until ready to proceed to the colorimetric assay.

**OPTIONAL:** Before proceeding to the colorimetric assay, the plate or dish can be centrifuged for 2-3 minutes to pellet the insoluble material. Use the supernatant for the assay.

#### B. Protocol For Suspension Cells

1. Aspirate the supernatant 24-72 hours post-transfection after centrifuging cells at 250 x g for 5 minutes. Optionally, wash cells 1 time with 1X PBS.
2. Resuspend the cell pellet in 1X Lysis Buffer. The amount of Lysis Buffer depends on the size of the culture dishes used for transfection (i.e., cell pellet size); we recommend using between 50 to 2000  $\mu$ l.
3. Incubate the cell lysate 10-15 minutes at room temperature by gently swirling the dishes several times to ensure complete lysis. Proceed to the colorimetric assay or freeze the plate at -70°C until ready.

**NOTE:** See NOTE above.

**OPTIONAL:** See OPTIONAL above.

4. Proceed to Section C. Below.

## Enhanced $\beta$ -Galactosidase Assay Kit (CPRG) – Cat#A10100K

### C. Colorimetric CPRG Assay

**NOTE:** Dilute 10X CPRG stock to 1X with Substrate Buffer just before performing the colorimetric assay. Unused 1X CPRG may be stored at  $-20^{\circ}\text{C}$  for future use. We recommend using 1X CPRG solution only 2 times after a freeze/thaw cycle.

**CAUTION:** Wear Gloves for manipulating the CPRG since it will stain exposed skin.

#### 96-well Microtiter Plate Assay\*

1. Thaw the dish, tube or plate of lysed cells at room temperature. If the transfection is performed with a 96-well plate, perform the assay directly on the plate (flat bottom only).
2. Add 50  $\mu\text{l}$  of Standard Dilution Buffer to the wells of a 96-well plate except control wells, which are set aside for a standard curve (see below).
3. Prepare a serial dilution of  $\beta$ -galactosidase (*E. coli*) standards using the Standard Dilution Buffer separately. A 50  $\mu\text{l}$  aliquot of each point on the standard curve is transferred to the control wells of the plate - the highest recommended amount of  $\beta$ -galactosidase is 100 milliunits (100,000–200,000 pg). 2X serial dilution of standard curve consisting of 8 points is recommended. A dilution protocol example is shown in the following table:

$\beta$ -gal Standard (miliunits)	Standard Dilution Buffer Volume	$\beta$ -gal Standard Volume
100	995 $\mu\text{l}$	5 $\mu\text{l}$ of $\beta$ -gal standard stock
50	200 $\mu\text{l}$	200 $\mu\text{l}$ 100 mu $\beta$ -gal standard
25	200 $\mu\text{l}$	200 $\mu\text{l}$ of 50 mu $\beta$ -gal standard
12.5	200 $\mu\text{l}$	200 $\mu\text{l}$ of 25 mu $\beta$ -gal standard
6.25	200 $\mu\text{l}$	200 $\mu\text{l}$ of 12.5 mu $\beta$ -gal standard
3.125	200 $\mu\text{l}$	200 $\mu\text{l}$ of 6.25 mu $\beta$ -gal standard
1.562	200 $\mu\text{l}$	200 $\mu\text{l}$ of 3.125 mu $\beta$ -gal standard
0.78	200 $\mu\text{l}$	200 $\mu\text{l}$ of 1.562 mu $\beta$ -gal standard

**NOTES:** Adjust the standard curve to suit the specific experimental conditions, such as cell type, transfection reagent, or plasmid vector. The dilutions for the standard curve must be prepared fresh each the assay.

4. Add 50  $\mu\text{l}$  of each sample/well.  
**NOTE:** It may be necessary to dilute the cell lysate in 1X Lysis Buffer when transfection efficiency is very high. In contrast, when transfection efficiency is low, reduce the volume of lysis buffer used to harvest the cells (see description above) or use up to 150  $\mu\text{l}$  of cell extract for the colorimetric assay. If the transfection is performed in a 96-well plate, perform the assay directly on the plate.
5. Prepare a blank by adding 50  $\mu\text{l}$  of lysis buffer to a well. Add also 50  $\mu\text{l}$  of cell lysate from non-transfected cells (negative control) to a well to control endogenous  $\beta$ -galactosidase activity.

6. Add 100  $\mu\text{l}$  of 1X CPRG Substrate Solution to each well. Incubate plate at room temp until dark red color develops (~10 minutes to 4 hours depending on cell type).
7. Read the absorbance at 570-595 nm with a microtiter spectrophotometer. Stop solution is not required for this format since all wells are read simultaneously without a time gap. Be sure that there are no bubbles present in the wells while reading because they interfere with the absorbance reading (remove bubbles with a fine gauge needle, tips, or very weak gas flow).
8. Quantify  $\beta$ -galactosidase expression based on a linear standard curve.

#### Macro assay

1. Thaw the cell lysate and transfer 100  $\mu\text{l}$  to a fresh tube, or 50  $\mu\text{l}$  to a 96-well plate. If a 96-well plate is used, follow the protocol described above.  
**NOTE:** It may be necessary to dilute the cell lysate in 1X Lysis Buffer when transfection efficiency is very high. In contrast, when transfection efficiency is low, reduce the volume of lysis buffer used to harvest the cells (see description above) or use up to 150  $\mu\text{l}$  of cell extract for the colorimetric assay.
2. Prepare a blank by adding 100  $\mu\text{l}$  of lysis buffer to a tube. Add also 100  $\mu\text{l}$  of cell lysate from non-transfected cells (negative control) to a tube to control endogenous  $\beta$ -galactosidase activity.
3. Add 50  $\mu\text{l}$  of Standard Dilution Buffer to each tube.
4. Prepare a serial dilution of  $\beta$ -galactosidase (*E. coli*) standards with Standard Dilution Buffer separately. Transfer 50  $\mu\text{l}$  of each standard to a fresh tube containing 100  $\mu\text{l}$  cell lysate from a mock transfection. The highest recommended amount of beta-galactosidase is 200,000 pg. (100 milliunits). Adjust the standard curve to suit the specific experimental conditions, such as cell type, transfection reagent, or plasmid vector. 2X serial dilution of standard curve consisting of 8 points is recommended. A dilution protocol example is shown in the section of 96-well plate assay.
5. Add 300  $\mu\text{l}$  of 1X CPRG Substrate Solution to each tube. Incubate the tubes at room temperature till the red color develops (from approximately 10 minutes to 4 hours depending on the cell type). Add 500  $\mu\text{l}$  of Stop Solution to stop the reaction. Final volume is 950  $\mu\text{l}$ .
5. Read the absorbance at 570-595 nm with a spectrophotometer.
6. Quantify  $\beta$ -galactosidase expression based on a linear standard curve.

\*Felgner, J.H. *et al.* Enhanced gene delivery and mechanism studies with a novel series of cationic lipid formulations. *J. Biol. Chem.* **269**, 2550-2561 (1994).