

Issue Date: Jul 1st, 2012

Bacillus Cereus Real Time PCR Kit User Manual

For In Vitro Diagnostic Use Only



For use with LightCycler 2.0 Instrument



Obelis S.A

Boulevard Général Wahis 53 1030 Brussels, BELGIUM Tel: +(32) 2.732.59.54 Fax: +(32) 2.732.60.03 E-Mail: mail@obelis.net







Shanghai ZJ Bio-Tech Co., Ltd. Tel: +86-21-34680596 www.liferiver.com.cn

trade@liferiver.com.cn Fax: +86-21-34680595 2nd floor,No.15 Building,No.188 Xinjunhuan road, PuJiang Hi-tech Park Shanghai China

1. Intended Use

The Bacillus Cereus real time PCR Kit is a test for the detection of Bacillus Cereus in stool or vomit samples in real time PCR systems.

2. Principle of Real-Time PCR

The principle of the real-time detection is based on the fluorogenic 5'nuclease assay. During the PCR reaction, the DNA polymerase cleaves the probe at the 5° end and separates the reporter dye from the quencher dye only when the probe hybridizes to the target DNA. This cleavage results in the fluorescent signal generated by the cleaved reporter dye, which is monitored real-time by the PCR detection system. The PCR cycle at which an increase in the fluorescence signal is detected initially (Ct) is proportional to the amount of the specific PCR product. Monitoring the fluorescence intensities during Real Time allows the detection of the accumulating product without having to re-open the reaction tube after the amplification.

3. Product Description

Bacillus cereus is an endemic, soil-dwelling, Gram-positive, rod-shaped, beta hemolytic bacteria that causes foodborne illness. It is the cause of "Fried Rice Syndrome". B. cereus bacteria are facultative anaerobes, and like other members of the genus Bacillus can produce protective endospores. B. cereus is responsible for a minority of foodborne illnesses (2-5%). It is known to create heavy nausea, vomiting, and abdominal periods. Generally speaking, Bacillus foodborne illnesses occur due to survival of the bacterial spores when food is improperly cooked. This problem is compounded when food is then improperly refrigerated, allowing the spores to germinate. Bacterial growth results in production of enterotoxin, and ingestion leads to two types of illness, diarrheal and emetic syndrome.

The Bacillus Cereus real time PCR Kit contains a specific ready-to-use system for the detection of the Bacillus Cereus using PCR (polymerase chain reaction) in the real-time PCR system. The master contains reagents and enzymes for the specific amplification of the Bacillus Cereus DNA. Fluorescence is emitted and measured by the real time systems' optical unit during the PCR. The detection of amplified Bacillus Cereus DNA fragment is performed by ROB reaction mix in fluorimeter **channel 530nm** with the fluorescent quencher BHQ1. Toxin types are identified by HBL&NHE reaction mix. In addition, the kit contains a system to identify possible PCR inhibition by measuring the 560nm fluorescence of the internal control (IC).

4. Kit Contents

Ref.	Type of Reagent	Presentation	25rxns
1	DNA Extraction Buffer	2 vials, 1.5ml	
2	ROB Reaction Mix	1 vial, 480µl	
3	HBL& NHE Reaction Mix	1 vial, 480µl	
4	PCR Enzyme Mix	1 vial, 22μl	
5	Molecular Grade Water	1 vial, 400µl	
6	Internal Control (IC)	1 vial, 30µl	
7	Positive Control	1 vial, 60µl	

Analysis sensitivity: 1×104copies/ml

Note: Analysis sensitivity depends on the sample volume, elution volume, nucleic acid extraction methods and other factors .If you use the DNA extraction buffer in the kit, the analysis sensitivity is the same as it declares. However, when the sample volume is dozens or even hundreds of times greater than elution volume by some concentrating method, it can be much higher.

5. Storage

- All reagents should be stored at -20°C. Storage at +4°C is not recommended.
- All reagents can be used until the expiration date indicated on the kit label.
 Repeated thawing and freezing (>3x) should be avoided, as this may reduce the sensitivity of
- Cool all reagents during the working steps.
- Reaction mix should be stored in the dark.

6. Additionally Required Materials and Devices

- · Biological cabinet
- Real time PCR system
- Desktop microcentrifuge for "eppendorf" type tubes (RCF max. 16,000 x g)
- · Vortex mixer
- Real time PCR reaction tubes/plates
- Cryo-container Pipets $(0.5\mu l 1000\mu l)$
- · Sterile filter tips for micro pipets Sterile microtubes
- Disposable gloves, powderlessBiohazard waste container
- Refrigerator and Freezer
- Tube racks

7. AWarnings and Precaution

- Carefully read this instruction before starting the procedure.
- For in vitro diagnostic use only.
- · This assay needs to be carried out by skilled personnel.
- Clinical samples should be regarded as potentially infectious materials and should be prepared in a laminar flow hood.
- This assay needs to be run according to Good Laboratory Practice.
- Do not use the kit after its expiration date.
- Avoid repeated thawing and freezing of the reagents, this may reduce the sensitivity of the
- · Once the reagents have been thawed, vortex and centrifuge briefly the tubes before use.
- Quickly prepare the reaction mix on ice or in the cooling block.
- Set up two separate working areas: 1) Isolation of the RNA/ DNA and 2) Amplification/ detection of amplification products.

 • Pipets, vials and other working materials should not circulate among working units
- Use always sterile pipette tips with filters.
- Wear separate coats and gloves in each area 8. Sample Collection, Storage and transportation

- · Collect samples in sterile tubes;
- Specimens can be extracted immediately or frozen at -20°C to -80°C.
- Transportation of clinical specimens must comply with local regulations for the transport of etiologic agents

9. Procedure

9.1 DNA-Extraction

DNA extraction buffer is supplied in the kit. Please thaw the buffer thoroughly and spin down briefly in the centrifuge before use. It's better to use commercial kits for nucleic acid extraction.

9.1.1 Stool samples

- 1) Take about 50mg samples to a 1.5ml tube; add 1.0ml normal saline then vortex vigorously. Centrifuge the tube at 13000rpm for 2 minutes, carefully remove and discard supernatant from the tube without disturbing the pellet.
- 2) Add 100µl DNA extraction buffer, close the tube then resuspend the pellet with vortex vigorously. Spin down briefly in a table centrifuge.
- 3) Incubate the tube for 10 minutes at 100°C.
 4) Centrifuge the tube at 13000rpm for 5 minutes. The supernatant contains the DNA extracted and can be used for PCR template.

9.1.2 Vomit samples

- 1) Take 1 ml vomit to a tube, Centrifuge the tube at 13000rpm for 2 minutes, carefully remove and discard supernatant from the tube without disturbing the pellet.
- 2) Add 100µl DNA extraction buffer, close the tube then vortex for 10 seconds. Spin down briefly in a table centrifuge.

 3) Incubate the tube for 10 minutes at 100°C.
- 4) Centrifuge the tube at 13000rpm for 5 minutes. The supernatant contains the DNA extracted and can be used for PCR template.

Attention:

A. During the incubation, make sure the tube is not open, as the vapor will

volatilize into the air and may cause contamination if the sample is positive.

B. The extraction sample should be used in 3 hours or store at -20°C for one month.


C. Different DNA extraction kits are available. You may use your own extraction systems or the commercial kit based on the yield. For the DNA extraction, please comply with the manufacturer's instructions.

9.2 Internal Control

It is necessary to add internal control (IC) in the reaction mix. Internal Control (IC) allows the user to determine and control the possibility of PCR inhibition. Add the internal control (IC) $1\mu l/rxn$ and the result will be shown in the 560nm.

9.3 PCR Protocol

The Master Mix volume for each reaction should be pipetted as follows:



1) The volumes of Reaction Mix and Enzyme Mix per reaction multiply with the number of samples, which includes the number of the controls, standards and sample prepared. Molecular Grade Water is used as the negative control. For reasons of unprecise pipetting, always add an extra virtual sample. (n: the number of reaction).Mix completely then spin down briefly in a centrifuge.

Reaction Volume	ROB Master Mix Volume	HBL& NHE Master Mix Volume	
Reaction Mix	17μ l × (n+1)	$18\mu l \times (n+1)$	
Enzyme Mix	$0.4\mu l \times (n+1)$	$0.4\mu l \times (n+1)$	
Internal control (IC)	$1\mu l \times (n+1)$	_	

2) Pipet $18\mu l$ Master Mix with micropipets of sterile filter tips to each *Real time* PCR reaction plate/tube. Then separately add $2\mu l$ DNA sample, positive and negative controls to different reaction plate/tubes. Immediately close the plate/tubes to avoid contamination.

3) Spin down briefly in order to collect the Master Mix in the bottom of the reaction tubes. 4) Perform the following protocol in the instrument:

37°C for 2 94°C for 2 93°C for 5

(Fluoresc

emin 1cycle			Selection of fluorescence channels		
2min 1cycle			Reaction Mix	530nm	560nm
5sec, 60°C for 30sec	40cycles		ROB	Bacillus Cereus	IC
cence measured at 60°C)			HBL&NHE	HBL toxin	NHE toxin

10.Threshold setting: Choose Arithmetic as back ground and none as Noise Band method, then adjust the Noise band just above the maximum level of molecular grade water, and adjust the threshold just under the minimum of the positive control.

11.Quality control: Negative control, positive control and internal control must be performed correctly, otherwise the sample results is invalid.

		Crossing point value	
Control	Reaction Mix	530nm	560nm
Molecular Grade Water	ROB	Blank	25~35 (IC)
	HBL&NHE	Blank	Blank
Positive Control	ROB	≤35	
	HBL&NHE	≤35	≤35
QS (quantitative detection)	In channel 530nm of ROB Reaction Mix, Correlation		
	coefficient of QS curve≤-0.98		

12. Data Analysis and Interpretation

	The following results are possible:							
I		RC)B	HBL&NHE				
		Reaction Mix		Reaction Mix		Result analysis		
		530nm	560nm	530nm	560nm			
ſ	1#	≤35		Blank	Blank	Bacillus Cereus Positive, and it is atoxic.		
ſ	2#	≤35		≤35 Blank		Bacillus Cereus Positive, and it contains HBL		
						toxin.		
	3#	≤35		UNDET	≤35	Bacillus Cereus Positive, and it contains NHE		
						toxin.		
	4#					PCR Inhibition; No diagnosis can be		
						concluded.		

*The crossing point value shows 35~40, please repeat again. If the result still shows 35~40,it can be considered negative.

For further questions or problems, please contact our technical support at trade@liferiver.com.cn