

DENVILLE ***SCIENTIFIC INC.***

a division of Harvard Bioscience, Inc.

The IncuCount



User Manual

Order 800.453.0385

www.densci.com • info@densci.com • Fax 508.429.5732
84 October Hill Road • Holliston, MA 01746

Table of Contents

Parts Checklist, Notes, Service Record	3
Warranty	4
Safeguards and Cautions, Limitations of Use	5
Intended Use	7
Specifications, Contact Information	8
Hardware Installation (Model RS-IC-150)	9
Quick Start	10
Step-By-Step Comprehensive Setup	10
Wi-Fi Setup	12
Automatic Updates	13
Colony Profile Setup	13
Saving a Colony Profile	14
Calibration Image	14
Setup – Lighting, Focusing, Brightness and Contrast Controls	16
Live Video Feed	18
Capturing Images	18
Colony Counting	19
Count Sensitivity	20
Count Region	20
Colony Size	21
Colony Shape (Morphology)	21
Zoom	22
Adding/Removing colonies from count	22
Exclusion Zone	22
Selecting Previous Image	23
Saving Data	23
Export to PDF	23
Antibiotic Reactions – Measurement of Inhibition Zones	24
Help	25
Appendix A-Camera Settings	25

This page intentionally left blank.

Order 800.453.0385

www.densci.com • info@densci.com • Fax 508.429.5732
84 October Hill Road • Holliston, MA 01746

Parts Checklist

- IncuCount
- Carrying case
- Light aperture ring
- Tablet CPU
- Tablet CPU Stand
- USB Belkin Hub
- Power OTG cable
- Detachable power supply w/ barrel connector
- Detachable power supply w/ USB connector
- Manual

Notes & Service Record

For your own future reference and to expedite answers to your questions, please complete the following information:

Date Purchased _____

Serial Number _____

Model Number _____

Purchased From _____

Your Notes:

DENVILLE **SCIENTIFIC INC.**

a division of Harvard Bioscience, Inc.

Warranty Registration

Please fill out this page, or a copy of it, and send it to the address below to validate your warranty. Warranty Registration Card for the Denville Scientific IncuCount.
(Please Print)

Date Purchased: _____ Serial Number: _____

Purchased From: _____

(Please check one) _____ Company _____ School _____ Individual

Name of Department or Institution:

Address:

Intended use:

Where did you hear of our product?

Comments:

Warranty Policy

Your Denville Scientific IncuCount is warranted to be free from defects of material or workmanship under normal use for a period of **two years** from date of purchase. If your IncuCount proves to be defective within this warranty period, APR&D, Inc. will repair or replace your unit at our option. This warranty will be void if malfunction is caused by accident, misuse, or negligence, including tampering, abuse, or damage in transit. This warranty excludes any relief for incidental or consequential damages.

Order 800.453.0385

www.densci.com • info@densci.com • Fax 508.429.5732
84 October Hill Road • Holliston, MA 01746

Safeguards & Cautions

To maximize the use of the IncuCount, please follow the recommendations in this short list before using:

1. Examine carton and contents for possible damage to the product caused during shipping. If damage is present contact carrier at once. Failure to do so immediately may release the carrier from any liability of damage.
2. Fill out the checklist (page of instruction manual). If something is missing please contact us.
3. Carefully read this instruction manual before operating the IncuCount.
4. Fill out warranty registration card and mail it to APR&D, Inc.

Limitations of Use

***Denville Scientific is not responsible for misuse of the IncuCount**

1. ALWAYS pre-clean the IncuCount light panel of residue after use with a damp cloth and appropriate mild detergent solution. Never clean IncuCount with organic solvents like thinner or benzene. It will damage the surfaces.

CAUTION: This is the only recommended method of cleaning.

2. Never attempt to dismantle or modify the IncuCount.
3. Never operate the IncuCount with a damaged USB cord or after the appliance malfunctions or has been damaged in any manner.
4. Accessory attachments not recommended by Revolutionary Science may damage the IncuCount or tablet cpu.
5. Never place the IncuCount or any part of the IncuCount in an environment that exceeds 65 degrees Celsius or below 5 degrees Celsius. Humidity should be at or below 65%.
6. Never use IncuCount for any purpose other than the intended use.


DENVILLE
SCIENTIFIC INC.

a division of Harvard Bioscience, Inc.

7. When used in higher temperatures, the IncuCount may be uncomfortable to the touch. When removing the IncuCount from the incubator, it is recommended to use temperature resistant gloves or turn off the incubator and allow the IncuCount to cool to a comfortable temperature before removing it.
8. Never attempt to modify the plugs or cords. This may cause damage to the IncuCount or tablet cpu.
9. Never immerse the IncuCount, cord, or plug in water or operate on wet surfaces.
10. Never allow the cords to contact hot surfaces beyond 65°C. Any slack in the cord should be wrapped up and held together with a band.

Order 800.453.0385

www.densci.com • info@densci.com • Fax 508.429.5732
84 October Hill Road • Holliston, MA 01746



IncuCount

The Denville Scientific IncuCount can be used for the following purposes:

1. Colony Counting:

- Count colonies on Petri dishes
- Multi color colony recognition

2. Measure the inhibition zones of antibiotic samples

- Manual
- Automatic

3. Monitor inside an incubated environment

- Count colonies as they grow in real time
- Measure inhibition zones as anti-biotic reactions occur in real time
- Email images

Order 800.453.0385

www.densci.com • info@densci.com • Fax 508.429.5732
84 October Hill Road • Holliston, MA 01746

DENVILLE
SCIENTIFIC INC.

a division of Harvard Bioscience, Inc.

Specifications


Power requirements:	100 -240 VAC 60 HZ 1.5 amps
Net weight:	5 lb. (shipping weight 14 lb.)
Dimensions:	7.25 inches tall x 6.50 inches diameter
Lights:	High intensity LED lights, top and bottom
Camera:	High-resolution digital camera, Carl Zeiss Optics
CPU:	Android based operating system Compatible with USB flash drive
Petri dish size capacity:	Up to 100 mm in size, including 3M Petri Film
Certifications:	UL approved LED lights, power supply, & camera Compatible with FDA 21 CFR part 11, export as PDF
Warranty:	2 year, repair or replacement

Order 800.453.0385

www.densci.com • info@densci.com • Fax 508.429.5732
84 October Hill Road • Holliston, MA 01746

Hardware Installation

Install hardware as shown in diagram below. (Model RS-IC-150)



START COUNTING

From out of the box to counting colonies

Quick

- 1) Setup device (as seen on page 9) and turn on.
- 2) Adjust lighting by use of the rotary switches in back.
- 3) Click count.

Comprehensive

- 1) Connect the IncuCount base to the Android tablet as shown on page 9. Make sure all components (USB hub, tablet, IncuCount base, USB OTG cable, power supply) are all connected correctly.
- 2) Turn on the Android tablet by holding in the power button for until the boot screen (Google) appears. The device should boot up straight into the IncuCount Software.
- 3) Check the optics by placing a blank Petri dish into the stage of the IncuCount base. A view of the Petri dish should appear on the screen.
- 4) Change the lighting by adjusting the Lighting knobs (p.16). We recommend using the lower lights on the second setting for most plates.
- 5) Auto focus and exposure are enabled by default. For manual camera adjustments, tap on Camera Controls in the menu bar, then select Manual Controls. Adjust the focus, brightness and contrast by moving the slider bars in the dialog that opened up until the desired contrast and brightness is reached (See Camera Controls, p.17).
- 6) Pause the image by tapping Preview in the menu the selecting pause. This will freeze the current frame in the camera for counting.
- 7) Adjust the counting region by tapping 'Manual' in the menu bar then selecting 'Set Counting Region'. Tap 'Center' in the dialog that opens up and + or— to resize the counting area. The counting area will be shaded in green. Centering works best with a blank or no Petri plate in the viewing area.
- 8) Place your sample Petri plate in the IncuCount. Make sure 'Preview' is selected.
- 9) Move the Sensitivity Slider up or down until the best count is reached. Colonies outlined in green are counted. Those Colonies outlined in red are either too big or too small (See Adjusting Colony Size, p.21)
- 10) Limit the size of detected colonies by Tapping Manual in the menu and selecting 'Set Colony Size'. Adjust colony size by moving the min and max size sliders to desired parameters. Colonies outside of parameters will be outlined in red.

11) Save the desired parameters by tapping 'Save'.

This will save a profile that will be automatically loaded the next time the IncuCount app is used.



12) Tap the count button to count the sample and view the details of the count.

13) If a PDF copy is desired, Tap File then select Save PDF.

14) To save a copy of the image, tap 'File' then select 'Capture Image'. Give the image a name and it will be saved to the IncuCount_images directory on the tablet.

15) To save the count data, tap 'File' then 'Save Data'.

16) To perform a new count, Tap the '>' button to the right of Count. A new image will be taken and counted with the previous saved parameters. If the camera needs adjusting again, Tap 'Pause' then select 'Preview'.




17) To receive automatic updates, turn on WI-FI (see p.12).

Wi-Fi Setup

To enable Wi-Fi (to receive automatic updates to the software), swipe down on the wifi icon (below) in the upper right corner of the tablet screen. Select the “Wi-Fi” option .



Select your network and type in your network key and press “Connect”.
When connected to the internet through wi-fi, the tablet can receive updates as they become available.




Automatic updates

The engineers at Denville Scientific are striving to make the feel and accuracy of the software better and better. To this end we release updates from time to time. As updates become available, they can be installed seamlessly by swiping the notification bar (upper left hand corner of the screen) and clicking “Update Available” option. The update will then download onto the tablet. Install the up-date when download is complete.




Colony Profile Setup

Many of the settings can be programmed in and saved for each agar/colony type being counted. To make a colony profile, simply adjust the camera settings (p.15), sensitivity settings (p.18) and click “Save Profile”.



Saving a Colony Profile

To save the profile, give the profile a name and click the enter button. Any type a change is made to the profile, it needs to be saved.




Making a Control “Calibration” Image

If the lighting isn't even throughout the agar plate, it may be necessary to make a calibration or control image to help filter out an unwanted light gradient. To do this, place a blank agar plate of the same type being counted into the IncuCount. Make sure preview is on and adjust the camera and light settings to give the best image.




Tap “Save in the Profile controls menu and give the image a name. This image is saved in a calibration folder for future use.



Applying the Calibration Image

To apply the calibration image, check the apply checkbox. The image will appear with the light gradient filtered out. Uncheck the box and the calibration image will be removed.



Light Control


The IncuCount comes equipped with an upper and a lower bank of LED lights. To increase or decrease the intensity of these lights, turn knob located in the back of the unit. The left knob controls the lower bank of lights and the right knob operates the upper bank. We recommend using the lower lights on a setting of two for most samples.

Light Aperture

Place the light aperture (black ring) into the unit on the light box.


Note: The purpose of this important mechanism is two fold.

1. Centers the Petri dish in the camera's field of view
2. Blocks out excess light around the Petri dish, providing the best possible view of the sample.




Focusing, Exposure, Brightness and Contrast Controls

Make sure camera mode is in preview (will show as 3rd tab on top menu bar). To adjust the image focus, exposure, brightness or contrast tap on “Camera Controls” and then tap on “Manual controls”. The controls will then appear on the right side of the screen. Simply slide up or down for adjustment.



When “Preview” is shown on top menu bar image displayed is live.


When “Pause” is shown on top menu bar image displayed is fixed. The “Pause” mode is setting which all counting occurs



Capturing Images

To capture an image for counting or measuring an inhibition zone, tap the 'CAPTURE' button under the File dropdown menu.


Note: Before capturing, you may need to adjust the focus and lighting.



Colony Counting

To count an image, tap the 'COUNT' button on the lower right side of screen. The results will appear in the notes box on the right side of the screen. In the image, colonies will be colored green, while those too large or small appear as red (see Sensitivity and Colony Size sections on pages 16 and 17). The IncuCount can differentiate colonies based on color. To turn on this feature, select the color checkboxes that most closely match the colony color. For example, to count purple vs. blue colonies, select the purple and blue checkboxes and click 'Count'. The count of each colony type will appear next to the selected colored checkboxes. If no checkbox is checked, the count will default to the green color.

Note: Before counting, you may need to adjust the sensitivity and/or the counting region.




Adjusting the Count Sensitivity

Locate the blue vertical sensitivity bar (located between the image and notes section) . To increase sensitivity, slide the bar up. To decrease the counting sensitivity slide blue vertical sensitivity bar down. A sensitivity of 0 means the average background color will be counted.

Note: Higher sensitivity is usually needed more with lighter colored colonies and lower sensitivity is usually needed more with darker colored colonies.

Adjusting the Counting Region

To adjust the counting region tap “Set Counting Region” under the manual tab. To recenter, tap ‘Center’ in the dialog that appears. To adjust size, by one unit, tap the ‘+’ or ‘-’ buttons. Tap the '+5x' or '-5x' to increase/decrease the size by 5 units. (A unit is .125mm.)




Adjusting Colony Size Parameters

Tap on the “Manual” tab at the top and then on “Set Colony Size” option. You will then see the sliding Max/Min scale to the right side of the screen. To increase the minimum allowable colony size slide up. To decrease the minimum colony size slide down. To change the maximum allowable size, adjust the “Max” slider.

Adjusting Colony Shape Parameters (morphology)

To include long skinny objects as colonies, adjust the L to W Ratio (Length/Width) slider up. To exclude long skinny objects from the count, adjust the slider down.



Zoom

The IncuCount software uses what is typically known as “pinch and zoom” method. To zoom in on a spot simply touch the screen with two fingers together and then slide them apart. To zoom out on an area touch the screen with two fingers apart and then slide them together.



Zoom in



Zoom out

Manually add colonies

Tap on the “Manual” tab at the top and then on “Add Colonies”. Then tap on the image where you want to add a new counted colony. The new counted colonies will show up as a white dots.

Note: added colonies will show up as a positive figure under “Colonies added” in the notes section on the right side of the screen.

Manually remove colonies

Tap on the “Manual” tab at the top and then on “Remove Colonies”. Then tap on a green counted colony to remove it. This will now show a red mark.

Note: removed colonies will show up as a negative figure under “Colonies added” in the notes section on the right side of the screen.

Exclusion zone

To exclude an area from counting, tap on “Manual” on the top menu and then tap on “Exclusion Zone”. Then tap around the region you want to exclude (four separate taps). The highlighted exclusion zone area will then appear.

To remove an existing exclusion zone double tap inside of the zone.

Selecting a Previous Image/Data

Tap on “File” and then tap on “Load Image”. The Load Image list of previous images will appear. Select desired image and confirm.

Saving Data

Tap on “File” and then tap on “Capture Image”. This brings up a “Saved As” menu. Type in desired name for the file and tap enter.

Export to PDF

Tap on “File” in the top menu and tap on “Save PDF”. This brings up a “Saved As” menu. Type in desired name for the file and tap enter.



Antibiotic Reactions, Measurement of Inhibition Zones

Tap on “Light Colonies” or “Dark Colonies” on the top menu bar. If disks on the screen appear light in color then tap on “Antibiotic (Light)”. If disks on the screen appear dark in color then tap on “Antibiotic (Dark)”. Then tap on the “Count” button in the lower right side of the screen. The software will identify the disks with green inhibition zones. Then tap on “Manual” and then on “Add Well”. This will allow the user to manually adjust each the diameter of each green zone.




Manually Add or Delete Inhibition Zones

To manually add a well tap on the “Manual” tab at the top and then on “Add Well”. Then tap on the image where you want to add a new well. The new well will show up as a green dot.

Note: added wells will show up in the notes section on the right side of the screen.

To manually remove a well tap on the “Manual” tab at the top and then on “Remove Well”. Then tap on the well you want to remove.



Help

The help menu has three items:

About - will give company name and model information


User Guide - will open up this document.

Contact information - for Denville Scientific's Technical support.

Email: info@densci.com

Telephone contact: 800-453-0385

www.densci.com



Appendix A—Camera Settings

Focus

Adjust focus for better clarity

Brightness

Adjusting Brightness lightens or darkens the image.

Contrast

Adjusting contrast heightens or lessens the difference between light and dark areas in the image.