


LabOne Programming Manual



Zurich
Instruments

LabOne Programming Manual

Zurich Instruments AG

Publication date Revision 31421

Copyright © 2008-2015 Zurich Instruments AG

The contents of this document are provided by Zurich Instruments AG (ZI), “as is”. ZI makes no representations or warranties with respect to the accuracy or completeness of the contents of this publication and reserves the right to make changes to specifications and product descriptions at any time without notice.

LabVIEW is a registered trademark of National Instruments Inc. All other trademarks are the property of their respective owners.

Revision History

Revision 31421, 8-Jul-2015:

Update of the LabOne Programming Manual for LabOne Release 15.05.

- The LabOne LabVIEW and C (ziAPI) APIs are now ziCore-based.
- Additions and modifications for MFLI support.

Revision 28870, 18-Mar-2015:

Update of the LabOne Programming Manual for LabOne Release 15.01.

- Added description of API Level 5.

Revision 26206, 01-Oct-2014:

Update of the LabOne Programming Manual for LabOne Release 14.08.

- Added PLL Advisor Module section to ziCore Modules.
- Consistency update of ziCore Module parameters.
- Improvements to plots in Software Trigger ziCore Module section.

Revision 23212, 23-Apr-2014:

- First release of the LabOne Programming Manual.
-

Table of Contents

I. LabOne Programming Concepts	4
1. Introduction	5
1.1. LabOne Software Architecture	6
1.2. Comparison of the LabOne Interfaces	9
1.3. LabOne API Levels	10
1.4. Finding settings: The Node Hierarchy	12
1.5. Obtaining Data from the Instrument	14
1.6. Instrument-Specific Considerations	16
2. ziCore Programming Overview	17
2.1. An Introduction to ziCore-based APIs	18
2.2. Sweeper Module	21
2.3. zoomFFT Module	30
2.4. Software Trigger (Recorder) Module	32
2.5. Device Settings Module	37
2.6. PLL Advisor Module	38
2.7. Tips and Tricks	40
II. LabOne APIs	41
3. Matlab Programming	42
3.1. Installing the LabOne Matlab API	43
3.2. Getting Started with the LabOne Matlab API	46
3.3. LabOne Matlab API Tips and Tricks	50
3.4. Troubleshooting the LabOne Matlab API	52
3.5. LabOne Matlab API (ziDAQ) Command Reference	54
4. Python Programming	70
4.1. Installing the LabOne Python API	71
4.2. Getting Started with the LabOne Python API	75
4.3. LabOne Python API Tips and Tricks	79
4.4. LabOne Python API (ziPython) Command Reference	80
5. LabVIEW Programming	113
5.1. Installing the LabOne LabVIEW API	114
5.2. Getting Started	116
5.3. LabVIEW Programming Tips and Tricks	121
6. C Programming	122
6.1. Getting Started	123
6.2. Module Documentation	125
6.3. Data Structure Documentation	201
6.4. File Documentation	237
Glossary	346
Index	352

Part I. LabOne Programming Concepts

This introduction gives an overview of LabOne programming and deals with generic concepts applicable to any of Zurich Instruments' APIs. It also helps guide the user to an appropriate choice of API.

Refer to:

- [Section 1.1](#) for an overview of the [LabOne Software Architecture](#) .
 - [Section 1.2](#) for a [Comparison of the LabOne Interfaces](#) .
 - [Section 1.3](#) for an explanation of [LabOne API Levels](#) .
 - [Section 1.4](#) for [Finding settings: The Node Hierarchy](#) .
 - [Section 1.5](#) for [Obtaining Data from the Instrument](#) .
 - [Section 1.6](#) for [Instrument-Specific Considerations](#) .
 - [Chapter 2](#) for an overview to working with `ziCore`-based APIs and Modules.
-

Chapter 1. Introduction

This chapter briefly describes the different possibilities to interface with a Zurich Instruments device, other than via the LabOne User Interface or ziControl (HF2 Series only). Zurich Instruments devices are designed with the concept that "the computer is the cockpit"; there are no controls on the front panel of the instrument, instead the user can configure their instrument from and stream data directly to their computer. The aim of this approach is to give the user the freedom to choose where they connect to, and how they control, their instrument.

As an example, the user can either work on a computer directly connected to the instrument via USB or remotely from a different computer on the network, away from their experimental setup. Then, on either computer, the user can configure and retrieve data from their instrument via a number of different software interfaces, i.e. via the web-based LabOne User Interface and/or their own custom programs. In this way the user can decide which connectivity setup and combination of interfaces best suits their experimental setup and data processing needs.

Refer to:

- [Section 1.1](#) for an overview of the [LabOne Software Architecture](#) .
- [Section 1.2](#) for a [Comparison of the LabOne Interfaces](#) .
- [Section 1.3](#) for an explanation of [LabOne API Levels](#) .
- [Section 1.4](#) for [Finding settings: The Node Hierarchy](#) .
- [Section 1.5](#) for [Obtaining Data from the Instrument](#) .
- [Section 1.6](#) for [Instrument-Specific Considerations](#) .

Note

New users could benefit by first familiarizing themselves with the instrument using the LabOne User Interface or ziControl; please refer to the appropriate user manual for your instrument for more details.

Note

The Real-time Option (RTK) for the HF2 Series is not a PC-based interface for controlling an instrument and is documented in the HF2 User Manual.

1.1. LabOne Software Architecture

Zurich Instruments devices uses a server-based connectivity methodology. Server-based means that all communication between the user and the instrument takes place via a computer program called a server, the Data Server. The Data Server recognizes available instruments and manages all communication between the instrument and the host computer on one side, and communication to all the connected clients on the other side. This allows for:

- A multi-client configuration: Multiple interfaces (even from multiple computers on the network) can access the settings and data on an instrument. Settings are synchronized between all interfaces by the single instance of the Data Server.
- A multi-device setup: Any of the Data Server's clients can access multiple devices simultaneously.

This software architecture is organized in layers, see [Figure 1.1](#) for a schematic of the software layers.




Figure 1.1. Software architecture

First, we briefly explain some terminology that is used throughout this manual.

- **Host computer:** The computer where the Data Server is running and that is directly connected to the instrument. Multiple remote computers on a local area network can access the instrument by creating an API connection to the Data Server running on the host computer.
- **Data Server :** A computer program that runs on the host computer and manages settings on, and data transfer to and from instruments by receiving commands from clients. It always has the most up-to-date configuration of the device and ensures that the configuration is synchronized between different clients.
- **ziServer.exe:** The Data Server that handles communication with HF2 Instruments.

- `ziDataServer.exe`: The Data Server that handles communication with UHFLI and MFLI Instruments. Note, in the case of MFLI Instruments the Data Server runs on the instrument itself.
- Remote computer: A computer, available on the same network as the host computer, that can communicate with an instrument via the Data Server program running on the host.
- Client: A computer program that communicates with an instrument via the Data Server. The client can be running either on the host or the remote computer.
- API (Application Programming Interface): a collection of functions and data structures which enable communication between software components. In our case, the various APIs (e.g., LabVIEW, Matlab®) provide functions to configure instruments and receive measured experimental data.
- Interface: Either a client or an API.
- GUI (Graphical User Interface): A computer program that the user can operate via images as opposed to text-based commands.
- LabOne User Interface: The browser-based user interface that connects to the Web Server.
- LabOne Web Server: The program that generates the browser-based LabOne User Interface.
- `ziControl`: The GUI for HF2 Instruments.
- `ziCore`: The core C++ library upon which many APIs are based, see [Part II](#) of this document.
- Modules: `ziCoreSoftware` components that provide a unified interface to APIs to perform high-level common tasks such as sweeping data.

1.1.1. LabOne Port and Hostname Selection

In order for an API client to connect to the correct Data Server, both the correct "hostname" and "port" must be provided when the API session is instantiated.

Firstly, LabOne Programmers must be aware that HF2 Instruments use a different Data Server program than UHFLI and MFLI Instruments. A LabOne API client connects to the correct Data Server for their instrument by specifying the appropriate port; use port 8004 for UHF and MFLI Instruments and 8005 for HF2 Instruments. Both servers can handle connections to multiple devices from multiple software interfaces (APIs or User Interfaces) simultaneously, see [Figure 1.2](#) for an example configuration.

Secondly, MFLI Programmers should be aware that, in contrast to HF2 and UHFLI instruments, the Data Server runs on the instrument itself. Thus, the instrument's hostname must be provided as the input argument when the API session is instantiated. This will be the same hostname that is used to access the LabOne User Interface on the instrument remotely from a web browser, please see the Getting Started chapter of the MFLI User Manual for help on determining the hostname.




Figure 1.2. Server port handling: Use 8004 for UHF Instruments and 8005 for HF2 Instruments.

1.2. Comparison of the LabOne Interfaces

The various software interfaces available in LabOne allow the user to pick a programming environment they are familiar with to achieve fast results. All other things being equal, here is a brief discussion of the merits of each interface.

- The [LabVIEW interface](#) allows for quick and efficient implementation of virtual instruments that run independently. These can easily be integrated in existing experiment control performed in LabVIEW. For most applications that need functions not included by `ziControl`, this is a sensible interface to use. This interface requires a National Instruments LabVIEW license and LabVIEW 2009 (or higher).
- The [Matlab® interface](#) allows the user to directly obtain measurement data within the Matlab programming environment, where they can make use of the many built-in functions available. The user can utilize high-level functionality provided by `ziCore`'s [Modules](#) to perform common measurement tasks. This interface requires a Mathworks Matlab license, but no additional Matlab Toolboxes.
- The [Python interface](#) allows the user to directly obtain measurement data within python. The user can utilize high-level functionality provided by `ziCore`'s [Modules](#) to perform common measurement tasks. Python is available as free and open source software; no license is required to use it.
- The [C API, `ziAPI`](#), is a very versatile interface that will run on most platforms. However, since C is a low-level programming language, the development cycle is slower than with the other programming environments and high-level functionality provided by `ziCore`'s [Modules](#) is not available.
- The text-based interface (HF2 Series only) allows the user to manually connect to the HF2 Data Server in a console via telnet. While this interface is a very useful tool for HF2 programmers to verify instrument configuration set by other interfaces, it is limited in terms of performance and maximum demodulator sample rate. See the HF2 User Manual for more details.

1.3. LabOne API Levels

All of the LabOne APIs depend on the base ziAPI (the C API). Needless to say, we try as hard as possible to make any improvements in ziAPI backwards compatible for the convenience of our users. We take care that existing programs do not need to be changed upon a new software release. Occasionally, however, we do have to make a breaking change in our API by removing some old functionality. This old functionality is, however, phased out over several software releases. First, the functionality is marked as deprecated and the user is informed via a depreciation warning (this can be turned off). This indicator warns that this function may be unsupported in the future. If we have to break some functionality we use a so-called **API level**.

With support of new devices and features we need to break functionality on the ziAPI.h e.g. data returned by poll commands. In order to still support the old functionality we introduced API levels. If a program is only using old functionality the API level 1 (default) can be used. If a user needs new functionality, they need to use a higher API level. This will usually need some changes on the existing code.

Available API levels as of LabOne Software Release 15.01 are:

- API Level 1: HF2 support, basic UHF support.
- API Level 4: UHF support with timestamps and PWA, name clean-up.
- API Level 5: Introduction of scope offset for extended (non-hardware) scope inputs (UHF).

Note that Levels 2 and 3 are used only internally and are not available to the general public.

Note

The HF2 Series only supports API Level 1.

Note

New UHFLI and MFLI API users are recommended to use API Level 5.

1.3.1. API Level 4 Features

The new features in API Level 4 are:

- Timestamps are available for any settings or data node.
- Greatly improved Scope data transfer rates (and new Scope data structure).
- Greatly improved UHF Boxcar and PWA support.

1.3.2. API Level 5 Features

API Level 5 was introduced in LabOne Release 15.01 to accommodate a necessary change in the Scope data structure:

- The Scope data structure was extended with the new field "channeloffset" which contains the offset value that must be added to the scaled wave value in order to obtain the physical value

recorded by the scope. For previous hardware scope "inputselects" there is essentially no change, since their offset is always zero. However, for the extended values of "inputselects", such as PID Out value, (available with the DIG option) the offset is determined by the values of "limitlower" and "limitupper" configured by the user.

1.4. Finding settings: The Node Hierarchy

In order to communicate with an Zurich Instruments device via text-based commands in an API, it is necessary to understand how the settings and measurement data of the instrument are accessed. All the settings and data of the instrument are organized in a file-system-like hierarchical structure. The features of the instrument, such as demodulators, are accessed as branches in this tree and their individual settings are leaves of these branches. It is also possible to browse branches inside the tree as if the user were navigating in a file-system. This hierarchy is used, no matter which interface you use when performing measurements.

An example demonstrating the hierarchy is the representation of the first demodulator on the device, given by the node:

```
/devX/demods/0
```

which, as we've already noted, is very similar to a **path** on a computer's file-system. Note that, the top level of the path is the device that you are connected to. The demodulators are then given as a top-level **node** under your device-node and the node of the first demodulator is indexed by 0. This path represents a branch in the node hierarchy which, in this case, if we explore further, has the following nodes:

```
/devX/demods/0/adselect  
/devX/demods/0/order  
/devX/demods/0/timeconstant  
/devX/demods/0/rate  
/devX/demods/0/trigger  
/devX/demods/0/oscsselect  
/devX/demods/0/harmonic  
/devX/demods/0/phaseshift  
/devX/demods/0/sinc  
/devX/demods/0/sample
```

These nodes are **leaves**, the most bottom-level nodes which represent a setting of an instrument or a field that can be read to retrieve measurement data. For example, `/devX/demods/0/adselect` is the leaf that controls the setting corresponding to the choice of signal input for the first demodulator. To set the index of the signal input the user writes to this node. The leaf `/devX/demods/0/sample` is the leaf where the demodulator's output (timestamp, demodulated x-value, demodulated y-value) are written at the frequency specified by `/devX/demods/0/rate`. In order to obtain the demodulator output you read the values from this node by **polling** this node. Polling a node sends a request from the client to ziServer to obtain the data from the node at that particular point in time.

Note

The numbering on the front panel of the UHFI, MFLI and HF2 Instruments and the block numbering in the graphical user interfaces generally start with 1 (1-based indexing). Note, that when accessing settings and data via a software interface, the numbering starting with 0 (0-based indexing).

Note

A useful method to learn about the nodes of your instrument is to look at the output of the history in the bottom of the graphical user interface. The status line always shows the last applied command and you can view the entire history by clicking the 'Show Log' or 'Show History' button. You will find paths like

```
/devx/sigins/0/ac = 1
```

after you switched on the AC mode for signal input 1, or

```
/devx/demods/1/rate = 7200.000000
```

after changing the readout rate of demodulator 2 to 7.2kHz.

1.5. Obtaining Data from the Instrument

The subscribe and poll commands

The easiest way to obtain data from an instrument is via the `poll` command, available in all of the LabOne API interfaces. The `poll` command is a function for synchronous data recording from specified [nodes](#) of an instrument. Synchronous means that the interface is blocked during execution of the command, see [Section 2.1.4](#) for asynchronous alternatives. `poll` takes two obligatory input arguments **recording time** and **timeout**.

The `subscribe` and `unsubscribe` commands are used to select the [nodes](#) from which data should be recorded. After subscribing to the node, the Data Server's internal data buffer will start filling with data from the subscribed nodes. The `poll` command will return the data that was recorded for the specified recording time (obligatory input argument) and any data that was already in the buffer since the last `poll`. To get rid of the data from earlier measurements it's possible to clear the buffer before polling by using the `flush` command.

In order to avoid losing data (the Data Server has a finite amount of memory available for its data buffers), long recording times (> 20s, depending on sampling rates and available memory) should be avoided. However, since internal data buffering on the Data Server ensures that no data is lost between `poll` commands, it's possible to record for longer periods of time by using the `poll` command inside a loop. In order to check that no data has been lost during a poll, the demodulator sample's `time` flags can be checked, see [Section 1.5.1](#).

If no data was stored in the Data Server's data buffer after issuing a `poll`, the command will wait for the data until the timeout time. If the buffer is empty after timeout time passed, `poll` will throw an error.

Note

If Matlab or Python are available, one of the LabOne `ziCore` Modules could be a more efficient choice for data retrieval than the comparably low-level `poll` command, see the [Section ziCore Modules](#) in [Part II](#).

1.5.1. Demodulator Sample Data Structure

An instrument's demodulator data is returned as a data structure (typically a `struct`) with the following fields (regardless of which API Level is used):

<code>sample.timestamp</code>	The instrument's timestamp of the measured demodulator data <code>uint64</code> . Divide by the instrument's clockbase (<code>/dev123/clockbase</code>) to obtain the time in seconds.
<code>sample.x</code>	The demodulator <code>x</code> value in Volts [<code>double</code>].
<code>sample.y</code>	The demodulator <code>y</code> value in Volts [<code>double</code>].
<code>sample.frequency</code>	The current frequency used by the demodulator in Hertz [<code>double</code>].
<code>sample.phase</code>	The oscillator's phase in Radians (not the demodulator phase) [<code>double</code>].
<code>sample.auxin0</code>	Auxiliary input channel 0 value in Volts [<code>double</code>].
<code>sample.auxin1</code>	Auxiliary input channel 1 value in Volts [<code>double</code>].
<code>sample.time.dataloss</code>	Indicator of sample loss (including block loss) [<code>bool</code>].
<code>sample.time.blockloss</code>	Indication of data block loss over the socket connection. This may be the result of a too long break between subsequent poll commands [<code>bool</code>].

`sample.time.invalidtimestamp` Indication of invalid time stamp data as a result of a sampling rate change during the measurement [bool].

Note

[Chapter 6](#) contains some details of other data structures.

1.6. Instrument-Specific Considerations

This section describes some instrument-specific considerations when programming with the LabOne APIs.

1.6.1. UHF-Specific Considerations

UHF Lock-in Amplifiers perform an automatic calibration 10 minutes after power-up of the Instrument. This internal calibration is necessary to achieve the specifications of the system. However, if necessary, it can be ran manually by setting the device node `/devN/system/calib/calibrate` to 1 and then disabled using the `/devN/system/calib/auto` node.

The calibration routine takes about 200 ms and during that time the transfer of measurement data will be stopped on the Data Server level. If a [ziAPI](#) (LabOne C API) or [LabVIEW](#) client is polling data during this time, the user will experience data loss; ziAPI has no functionality to deal with such a streaming interrupt. Clients polling data from [ziCore-based APIs](#) (i.e. Matlab or Python APIs) will be informed of dataloss, which allows the user to ignore this data.

Please see the UHF User Manual for more information about device calibration.

Chapter 2. ziCore Programming Overview

The LabOne Matlab and Python APIs provide interfaces to configure, acquire data from, and run integral functionality of your Zurich Instruments device in powerful high-level programming environments. These high-level interfaces are, however, just thin application layers based on a shared core API, `ziCore`. This chapter aims to describe the common functionality that's available to any of the interfaces based on `ziCore`.

Refer to:

- [Section 2.1](#) for [An Introduction to ziCore-based APIs](#) .
- [Section 2.2](#) for the [Sweeper Module](#) .
- [Section 2.3](#) for the [zoomFFT Module](#) .
- [Section 2.4](#) for the [Software Trigger \(Recorder\) Module](#) .
- [Section 2.5](#) for the [Device Settings Module](#) .
- [Section 2.6](#) for the [PLL Advisor Module](#) .
- [Section 2.7](#) for some `ziCore` programming [Tips and Tricks](#) .

2.1. An Introduction to ziCore-based APIs

All the ziCore-based APIs share a common structure which provides a uniform interface for programming Zurich Instruments devices. The aim of this section is to familiarize the user with the key ziCore programming concepts.

2.1.1. Software Architecture

Each of the ziCore-based APIs are designed to have a minimal code footprint: They are simply small interface layers that use the functionality derived from ziCore, a central C++ API. The derived API interfaces, Matlab and Python, provide a familiar interface to the user and allow the user to receive and manipulate data from their instrument using the API language's native data types and formats. See [Section 1.1](#) for an overview of the LabOne software architecture.

2.1.2. ziCore Modules

In addition to the usual API commands available for instrument configuration and data retrieval, e.g., `setInt`, `poll`, ziCore-based APIs also provide a number of so-called **Modules**: high-level interfaces that perform common tasks such as sweeping data or performing [FFT](#) s.

The Module's functionality is implemented in ziCore and each derived high-level API simply provides an interface to that module from the API's native environment. This design ensures that the user can expect the same behavior from each module irrespective of which API is being used; if the user is familiar with a module available in one high-level programming API, it is quick and easy to start using the module in a different API. In particular, the LabOne User Interface is also based on ziCore and as such, the user can expect the same behavior using a ziCore-based API that is experienced in the LabOne User Interface, see [Figure 2.1](#).




Figure 2.1. The same results and behavior can be obtained from Modules in any ziCore-based interface; [Sweeper Module](#) results from the LabOne Matlab API (left) and the LabOne User Interface (right) using the same Sweeper and instrument settings.

The modules currently available in ziCore are:

- The [Sweeper Module](#) for obtaining data whilst performing a sweep of one of the instrument's setting, e.g., measuring a frequency response.
- The [zoomFFT Module](#) for calculating the FFT of demodulator output.

- The [Software Trigger \(Recorder\) Module](#) for recording instrument data [asynchronously](#) based upon user-defined triggers.
- The [Device Settings Module](#) for saving and loading instrument settings to and from (XML) files.
- The [PLL Advisor Module](#) for modelling/simulating the PLL (phase-locked loop) incorporated in the instrument (available for UHF Lock-in Amplifiers only).

In addition to providing a unified-interface between APIs, modules also provide a uniform work-flow regardless of the functionality the module performs (e.g., sweeping, recording data), see [Section 2.1.3](#).

An important difference to low-level `ziCore` API commands is that Modules execute their commands **asynchronously**, see [Section 2.1.4](#).

Note

The LabOne User Interface Command Log can be set to store commands in either Matlab or Python formats which can then be used to start writing custom programs, see [Section 2.7](#).

Note

Much of the same functionality is provided in `ziControl`, but `ziControl` UI is not based on `ziCore`.

2.1.3. ziCore Module Work-Flow

Regardless of the Module's function, all `ziCore` Modules follow same work flow in all of the derived interfaces:

- create (instantiate) an instance of the module,
- **set** the module's parameters using `path`, `value` pairs,
- **subscribe** to instrument nodes from which to obtain data,
- **execute** the module,
- wait until the module has **finished** executing; intermediate reading of data is possible,
- **read** the module's data,
- **clear** the module to remove it from memory.

The highlighted words above are commands for all the Modules. For interface-specific concepts when using Modules see the following Sections:

- [Using ziCore Modules in the LabOne Matlab API](#).
- [Using ziCore Modules in the LabOne Python API](#).

2.1.4. Synchronous versus Asynchronous Commands

The low-level API commands such as `setInt` and `poll` are **synchronous** commands, that is the interface will be blocked until that command has finished executing; the user can not run any commands in the meantime. Another feature of `ziCore`'s Modules is that each instantiation of a Module creates a new `Thread` and, as such, the commands executed by a Module are performed **asynchronously**. Asynchronous means that the task is performed in the background and the interface's process is available to perform other tasks in the meantime, i.e., Module commands are non-blocking for the user.

2.2. Sweeper Module

The Sweeper Module allows the user to perform sweeps as in the Sweeper Tab of the LabOne User Interface. In general, the Sweeper can be used to obtain data when measuring a [DUT](#)'s response to varying (or **sweeping**) one instrument setting while other instrument settings are kept constant.

2.2.1. Configuring the Sweeper

In the following we briefly describe how to configure the Sweeper Module. See [Table 2.1](#) for a full list of the Sweeper's input parameters and [Table 2.2](#) for a description of the Sweeper's outputs.

Specifying the Instrument Setting to Sweep

The Sweeper's `sweep/gridnode` parameter, the so-called **sweep parameter**, specifies the instrument's setting to be swept, specified as a path to an instrument's [node](#). This is typically an oscillator frequency in a [Frequency Response Analyzer](#), e.g., `/dev123/oscs/0/freq`, but a wide range of instrument settings can be chosen, such as a signal output amplitude or a PID controller's setpoint.

Specifying the Range of Values for the Sweep Parameter

The Sweeper will change the sweep parameter's value `sweep/samplecount` times within the **range** of values specified by `sweep/start` and `sweep/stop`. The `sweep/xmapping` parameter specifies whether the spacing between two sequential values in the range is linear (`=0`) or logarithmic (`=1`).

Controlling the Scan mode: The Selection of Range Values

The `sweep/scan` parameter defines the **order** that the values in the specified range are written to the sweep parameter. In sequential scan mode (`=0`), the sweep parameter's values change incrementally from smaller to larger values, see [Figure 2.3](#).

In binary scan mode (`=1`) the first sweep parameter's value is taken as the value in the middle of the range, then the range is split into two halves and the next two values for the sweep parameter are the values in the middle of those halves. This process continues until all the values in the range were assigned to the sweep parameter, see [Figure 2.5](#). Binary scan mode ensures that the sweep parameter uses values from the entire range near the beginning of a measurement, which allows the user to get feedback quickly about the measurement's entire range. Since the Sweeper Module is an [asynchronous](#) interface, it's possible to continuously read and plot data whilst the sweep measurement is ongoing and update points in a graph dynamically.

In bidirectional scan mode (`=2`) the sweep parameter's values are first set from smaller to larger values as in sequential mode, but are then set in reverse order from larger to smaller values, see [Figure 2.4](#). This allows for effects in the sweep parameter to be observed that depend on the order of changes in the sweep parameter's values.

Controlling how the Sweeper sets the Demodulator's Time Constant

The `sweep/bandwidthcontrol` parameter specifies which demodulator filter bandwidth (equivalently time constant) the Sweeper should set for the current measurement point. The user can either specify the bandwidth manually (`=0`), in which case the value of the current demodulator filter's bandwidth is simply used for all measurement points; specify a fixed bandwidth (`=1`), specified by `sweep/bandwidth`, for all measurement points; or specify that the Sweeper sets the demodulator's bandwidth automatically (`=2`). Note, to use either Fixed or Manual mode, `sweep/bandwidth` must be set to a value > 0 (even though in manual mode it is ignored).

Specifying the Sweeper's Settling Time

For each change in the sweep parameter that takes effect on the instrument the Sweeper waits before recording measurement data in order to allow the measurement to settle. This behavior is configured with the parameters in the `sweep/settling/` branch.

The `sweep/settling/tc` parameter defines the minimum number of time constants (of the demodulator filter) that will Sweeper module wait after setting the `sweep/gridnode` value before recording measurement data. This corresponds to the settling time required by the lock-in amplifier to measure the response in the system. The `sweep/settling/time` parameter specifies the minimum number of seconds to wait before taking the measurement data, which corresponds to the settling time required by the user's experimental setup before measuring the response in the system. The actual time that the Sweeper waits before recording data is defined in [Equation 2.1](#).

$$t_s = \max(\text{sweep_settling_tc} \times \text{tc}, \text{sweep_settling_time})$$

Equation 2.1. The settling time t_s used by the Sweeper for each measurement point; the amount of time between setting the sweep parameter and recording measurement data is determined by the parameters `sweep/settling/tc` and `sweep/settling/time`.

Note, the Sweeper may change the value of the demodulator filter's bandwidth (equivalently time constant), depending on the value of `sweep/bandwidthcontrol` and this will effect the settling time as specified by `sweep/settling/tc`, see [above, Controlling how the Sweeper sets the Demodulator's Time Constant](#). For a frequency sweep, the `sweep/settling/tc` parameter will tend to influence the settling time at lower frequencies, whereas `sweep/averaging/sample` will influence the settling time at higher frequencies.

The recommended value for the `sweep/settling/tc` setting is 15. For noise measurements a higher value should be used, e.g, 50.

Specifying which Data to Measure

Which measurement data is actually returned by the Sweeper's `read` command is configured by subscribing to [node paths](#) using the Sweeper Module's `subscribe` command.

Specifying how the Measurement Data is Averaged

One Sweeper measurement point is obtained by averaging recorded data which is configured via the parameters in the `sweep/averaging/` branch.

The `sweep/averaging/tc` parameter specifies the minimum time window in factors of demodulator filter time constants during which samples will be recorded in order to average for one returned sweeper measurement point. The `sweep/averaging/sample` parameter specifies the minimum number of data samples that should be recorded and used for the average. The Sweeper takes both these settings into account for the measurement point's average according to [Equation 2.2](#).

$$N = \max(\text{sweep_averaging_tc} \times \text{tc} \times \text{sampling_rate}, \text{sweep_averaging_sample})$$

Equation 2.2. The number of samples N used to average one sweeper measurement point is determined by the parameters `sweep/averaging/tc` and `sweep/averaging/sample`.

Note, the value of the demodulator filter's time constant may be controlled by the Sweeper depending on the value of `sweep/bandwidthcontrol` and `sweep/bandwidth`, see [above, Controlling how the Sweeper sets the Demodulator's Time Constant](#). For a frequency sweep, the `sweep/averaging/tc` parameter will tend to influence the number of samples recorded at lower frequencies, whereas `sweep/averaging/sample` will influence averaging behavior at higher frequencies.

An Explanation of Settling and Averaging Times in a Frequency Sweep

[Figure 2.2](#) shows which demodulator samples are used in order to calculate an averaged measurement point in a frequency sweep. This explanation of the Sweeper's parameters is specific to the following commonly-used Sweeper settings:

- `sweep/gridnode` is set to an oscillator frequency, e.g., `/dev123/oscs/0/freq`.
- `sweep/bandwidthcontrol` is set to 2, corresponding to automatic bandwidth control, i.e., the Sweeper will set the demodulator's filter bandwidth settings optimally for each frequency used.
- `sweep/scan` is set to 0, corresponding to sequential scan mode for the range of frequency values swept, i.e., the frequency is increasing for each measurement point made.

Each one of the three red segments in the demodulator data correspond to the data used to calculate one single Sweeper measurement point. The light blue bars correspond to the time the sweeper should wait as specified by the `sweep/settling/tc` parameter and the purple bars correspond to the time specified by the `sweep/settling/time` parameter. The sweeper will wait for the maximum of these two times according to [Equation 2.1](#). When measuring at lower frequencies the Sweeper sets a smaller demodulator filter bandwidth (due to automatic `sweep/bandwidthcontrol`) corresponding to a larger demodulator filter time constant. Therefore, the `sweep/settling/tc` parameter dominates the settling time used by the Sweeper at low frequencies and at high frequencies the `sweep/settling/time` parameter takes effect. Note, that the light blue bars corresponding to the `sweep/settling/tc` parameter get shorter for each measurement point (larger frequency used → shorter time constant required), whereas the purple bars corresponding to `sweep/settling/time` stay a constant length for each measurement point. Similarly, the `sweep/averaging/tc` parameter (yellow bars) dominates the Sweeper's averaging behavior at low frequencies, whereas `sweep/averaging/samples` (green bars) specifies the behavior at higher frequencies, see also [Equation 2.2](#).




Figure 2.2. Plot demonstrating how the Sweeper records three measurement points from demodulator data when using automatic bandwidth control in a frequency sweep. Please see [An Explanation of Settling and Averaging Times in a Frequency Sweep](#), above, for a detailed explanation.

Table 2.1. Sweeper Input Parameters

Setting/Path	Type	Unit	Description
sweep/device	byte array	-	The device ID to perform the sweep on, e.g., dev123 (compulsory parameter).
sweep/gridnode	byte array	Node	The device parameter (specified by node) to be swept, e.g., 'oscs/0/freq'.
sweep/start	double	Many	The start value of the sweep parameter.
sweep/stop	double	Many	The stop value of the sweep parameter.
sweep/samplecount	uint64	Many	The number of measurement points to set the sweep on.
sweep/endless	bool	-	Enable Endless mode; run the sweeper continuously.
sweep/averaging/sample	uint64	Samples	Sets the number of data samples per sweeper parameter point that is considered in the measurement. The maximum of this value and sweep/averaging/tc is taken as the effective calculation time. See Figure 2.2 .
sweep/averaging/tc	double	Seconds	Sets the effective measurement time per sweeper parameter point that is considered in the measurement. The maximum between of this value and sweep/averaging/sample is taken as the effective calculation time. See Figure 2.2 .
sweep/bandwidthcontrol	uint64	-	Specify how the sweeper should specify the bandwidth of each measurement point, Automatic is recommended, in particular for

Setting/Path	Type	Unit	Description
			logarithmic sweeps and assures the whole spectrum is covered. 0=Manual (the sweeper module leaves the demodulator bandwidth settings entirely untouched); 1=Fixed (use the value from <code>sweep/bandwidth</code>); 2=Automatic. Note, to use either Fixed or Manual mode, <code>sweep/bandwidth</code> must be set to a value > 0 (even though in manual mode it is ignored).
<code>sweep/bandwidth</code>	double	Hz	Defines the measurement bandwidth for Fixed bandwidth selection, and corresponds to the noise equivalent power bandwidth (NEP).
<code>sweep/order</code>	uint64	-	Defines the filter roll off to use in Fixed bandwidth selection. Valid values are between 1 (6 dB/octave) and 8 (48 dB/octave).
<code>sweep/maxbandwidth</code>	double	Hz	Specifies the maximum bandwidth used when in auto bandwidth mode (<code>sweep/bandwidthcontrol=2</code>). The default is 1.25MHz.
<code>sweep/omegasuppression</code>	double	dB	Damping of omega and 2omega components when in auto bandwidth mode (<code>sweep/bandwidthcontrol=2</code>). Default is 40dB in favor of sweep speed. Use a higher value for strong offset values or 3omega measurement methods.
<code>sweep/loopcount</code>	uint64	-	The number of sweeps to perform.
<code>sweep/phaseunwrap</code>	bool	-	Enable unwrapping of slowly changing phase evolutions around the +/-180 degree boundary.
<code>sweep/sincfilter</code>	bool	-	Enables the sinc filter if the sweep frequency is below 50 Hz. This will improve the sweep speed at low frequencies as omega components do not need to be suppressed by the normal low pass filter.
<code>sweep/scan</code>	uint64	-	Selects the scanning type: 0=Sequential (incremental scanning from start to stop value, see Figure 2.3); 1=Binary (Non-sequential sweep continues increase of resolution over entire range, see Figure 2.5), 2=Bidirectional (Sequential sweep from Start to Stop value and back to Start again, Figure 2.4).
<code>sweep/settling/tc</code>	double	TC	Minimum wait time in factors of the time constant (TC) between setting the new sweep parameter value and the start of the measurement. The maximum between this value and <code>sweep/settling/time</code> is taken as effective settling time. See Figure 2.2 .

Setting/Path	Type	Unit	Description
sweep/settling/time	double	Seconds	Minimum wait time in seconds between setting the new sweep parameter value and the start of the measurement. The maximum between this value and sweep/settling/ τ_c is taken as effective settling time. See Figure 2.2 .
sweep/settling/inaccuracy	double	-	Demodulator filter settling inaccuracy defining the wait time between a sweep parameter change and recording of the next sweep point. Typical inaccuracy values: 10m for highest sweep speed for large signals, 100u for precise amplitude measurements, 100n for precise noise measurements. Depending on the order the settling accuracy will define the number of filter time constants the sweeper has to wait. The maximum between this value and the settling time is taken as wait time until the next sweep point is recorded.
sweep/xmapping	uint64	-	Selects the spacing of the grid used by sweep/gridnode (the sweep parameter): 0=linear and 1=logarithmic distribution of sweep parameter values.
sweep/historylength	uint64		Maximum number of entries stored in the measurement history.
sweep/clearhistory	bool	-	Remove all records from the history list.
sweep/filename	byte array	-	This parameter is deprecated. If specified, i.e. not empty, it enables automatic saving of data in single sweep mode (sweep/endless = 0).
sweep/savepath	byte array	-	The directory where files are located when saving sweeper measurements.
sweep/fileformat	byte array	-	The format of the file for saving sweeper measurements. 0=Matlab, 1=CSV.




Figure 2.3. Sweeper scanning modes: Sequential (sweep/scan = 0).




Figure 2.4. Sweeper scanning modes: Bidirectional ($\text{sweep}/\text{scan} = 2$).




Figure 2.5. Sweeper scanning modes: Binary (sweep/scan = 1).

Table 2.2. Sweeper Output Values

Name	Type	Unit	Description
timestamp	uint64	Ticks	A timestamp that gets updated each time a new measurement point has been recorded by the sweeper (divide by the device's clockbase to obtain seconds). It is not part of the sweeper's measurement data and only relevant for intermediate reads of sweeper data (before the current sweep has finished).
frequency	double	Hz	The oscillator frequency for each measurement point (for a frequency sweep this is the same as grid).
x	double	Volts	Demodulator x value.
y	double	Volts	Demodulator y value.
r	double	VoltsRMS	Demodulator R value.

Name	Type	Unit	Description
phase	double	Radians	Demodulator phase value.
auxin0	double	Volts	Auxiliary Input 1 value.
auxin1	double	Volts	Auxiliary Input 2 value.
frequencystddev	double	Hz	Standard deviation of the oscillator frequency.
xstddev	double	Volts	Standard deviation of demodulator x value.
ystddev	double	Volts	Standard deviation of demodulator y value.
rstddev	double	VoltsRMS	Standard deviation of demodulator R value.
phasestddev	double	Radians	Standard deviation of demodulator phase value (phase noise).
auxin0stddev	double	Volts	Standard deviation of Auxiliary Input 1 value.
auxin1stddev	double	Volts	Standard deviation of Auxiliary Input 2 value.
xpwr	double	Volts ²	Average power of demodulator x value.
ypwr	double	Volts ²	Average power of demodulator y value.
rpwr	double	Volts ²	Average power of demodulator R value.
rpwr	double	Volts ²	Average power of demodulator x value.
auxin0pwr	double	Volts ²	Average power of Auxiliary Input 1 value.
auxin1pwr	double	Volts ²	Average power of Auxiliary Input 2 value.
grid	double	Many	Values of sweeping setting (frequency values at which demodulator samples where recorded).
bandwidth	double	Hz	Demodulator filter's bandwidth set for each measurement point (assuming we're performing a frequency sweep).
tc	double	Seconds	Demodulator's filter time constant for each measurement point.
settling	double	Seconds	The waiting time for each measurement point.
setttimestamp	uint64	Ticks	The timestamp at which we verify that the frequency for the current measurement point was set on the device (by reading back demodulator data).
nexttimestamp	uint64	Ticks	The timestamp at which we can obtain the data for that measurement point, i.e., <code>nexttimestamp - setttimestamp</code> corresponds roughly to the demodulator filter settling time.

2.3. zoomFFT Module

The zoomFFT Module corresponds to the Spectrum Tab of the LabOne User Interface. It allows the user to perform Fast Fourier Transforms (FFT) on a specified demodulator's output.

See [Table 2.3](#) for the input parameters to configure the ZoomFFT Module and [Table 2.4](#) for a description of the ZoomFFT's outputs.




Figure 2.6. A plot of an FFT created by one of the LabOne Matlab API examples.

Table 2.3. ZoomFFT Input Parameters

Setting/Path	Type	Unit	Description
zoomFFT/device	byte array	-	The device ID to perform the FFT on, e.g., dev123 (compulsory parameter).
zoomFFT/absolute	bool	-	Shifts the frequencies so that the center frequency becomes the demodulation frequency rather than 0 Hz.
zoomFFT/bit	uint64	-	Number of lines of the FFT spectrum (powers of 2). Increasing the bits increases the frequency resolution of the spectrum.
zoomFFT/endless	bool	-	Enable Endless mode; run the zoomFFT continuously.
zoomFFT/loopcount	uint64	-	The number of zoomFFTs to perform.
zoomFFT/mode	uint64	-	Select the source signal for the FFT. 0=FFT(x+iy), 1=FFT(R), 2=FFT(phase), 3=FFT(Freq)
zoomFFT/overlap	double	-	Overlap of the demodulator data used for the FFT. Use 0 for no overlap and 0.99 for maximal overlap.

Setting/Path	Type	Unit	Description
zoomFFT/settling/tc	double	TC	Minimum wait time in factors of the time constant (TC) before starting the measurement. The maximum between this value and zoomFFT/settling/time is taken as effective settling time.
zoomFFT/settling/time	double	Seconds	Minimum wait time in seconds before starting the measurement. The maximum between this value and zoomFFT/settling/tc is taken as effective settling time.
zoomFFT/window	uint64	-	The type of FFT window to use. 0=Rectangular, 1=Hann, 2=Hamming, 3=Blackman Harris.

Table 2.4. ZoomFFT Output Values

Name	Type	Unit	Description
x	double	Volts	The real part, x, of the complex FFT result.
y	double	Volts	The imaginary part, y, of the complex FFT result.
r	double	VoltsRMS	The absolute value, R, of the complex FFT result.
timestamp	uint64	Ticks	Demodulator timestamp of the measurement (divide by the device's clockbase to obtain seconds)
center	double	Hz	The center frequency (corresponds to the demodulation frequency).
rate	double	-	Sampling rate of the demodulator.
filter	double	-	The filter envelope; the filter compensation value for each gridnode.
bandwidth	double	Hz	The bandwidth of the demodulator
grid	double	Hz	The frequency grid.
nenbw	double	-	The normalized equivalent noise bandwidth.
resolution	double	Hz	FFT resolution: Spectral resolution defined by the reciprocal acquisition time (sample rate, number of samples recorded).
aliasingreject	double	dB	How much damping is present at the border of your spectrum.

2.4. Software Trigger (Recorder) Module

The Recorder Module corresponds to the Software Trigger Tab of the LabOne User Interface. It allows the user to record bursts of instrument data based upon pre-defined trigger criteria similar to that of a laboratory oscilloscope, see [Figure 2.7](#) for an example. The types of trigger available are listed in [Table 2.5](#).

Table 2.5. Overview of the trigger types available in the Software Trigger Module.

Trigger Type	Description	trigger/N/type
Manual	For simple recording.	0
Edge	Edge trigger with level hysteresis and noise rejection, see Figure 2.8 .	1
Digital	Digital trigger with bit masking.	2
Pulse	Pulse width trigger with level hysteresis and noise reduction, see Figure 2.9 and Figure 2.10 .	3
Tracking (edge or pulse)	Level tracking trigger to compensate signal drift, see Figure 2.11 .	4
Hardware Trigger	UHFLI and MFLI only. Trigger on one of the instrument's hardware trigger channels.	6




Figure 2.7. The plot produced by `example_record_edge_trigger.py`, an example distributed with the LabOne Python API. The plot shows 10 bursts of data from a single demodulator; each burst was recorded when the demodulator's R value exceeded a specified threshold using a positive edge trigger. See [Section 4.2.3](#) for help getting started with the Python examples.

See [Table 2.6](#) for the input parameters to configure the Software Trigger's Module. Note that some parameters effect all triggers, e.g., `trigger/endless`, whereas some are configured on a per-trigger basis, e.g., `trigger/N/duration`, where N is the index of the trigger, starting at zero. The data output when using the Software Trigger's read command has the same format as returned by `ziCore's poll` command.




Figure 2.8. Explanation of the Software Trigger Module's parameters for an Edge Trigger.




Figure 2.9. Explanation of the Software Trigger Module's parameters for a positive Pulse Trigger.




Figure 2.10. Explanation of the Software Trigger parameters for a positive or negative Pulse Trigger.




Figure 2.11. Explanation of the Software Trigger Module's parameters for a Tracking Trigger.

Table 2.6. Software Trigger Input Parameters.

Setting/Path	Type	Unit	Description
trigger/device	byte array	-	The device ID to execute the software trigger, e.g., dev123 (compulsory parameter).
trigger/buffersize	double	Seconds	Overwrite the buffersize of the trigger object (set when it was instantiated). Recommended buffer size is $2 * \text{trigger}/N / \text{duration}$.
trigger/endless	uint64	-	Enable endless triggering 1=enable; 0=disable.
trigger/forcetrigger	uint64	-	Force a trigger.

Setting/Path	Type	Unit	Description
trigger/filename	byte array	-	This parameter is deprecated. If specified, i.e. not empty, it enables automatic saving of data in single triggering mode (<code>trigger/endless = 0</code>).
trigger/savepath	byte array	-	The directory where files are saved when saving data.
trigger/fileformat	byte array	-	The format of the file for saving data. 0=Matlab, 1=CSV.
trigger/historylength	uint64	-	Maximum number of entries stored in the measurement history.
trigger/clearhistory	uint64	-	Clear the measurement history
trigger/triggered	uint64	-	Has the software trigger triggered? 1=Yes, 0=No (read only).
trigger/N/bandwidth	double	Hz	Only for Tracking Triggers. The bandwidth used in the calculation of the exponential running average of the source signal.
trigger/N/bitmask	uint64	-	Only for Digital triggers. Specify the bitmask used with <code>trigger/N/bits</code> . The trigger value is bits AND bit mask (bitwise).
trigger/N/bits	uint64	-	Only for Digital triggers. Specify the bits used for the Digital trigger value. The trigger value is bits AND bit mask (bitwise)
trigger/N/count	uint64	-	The number of triggers to save.
trigger/N/delay	uint64	Seconds	The amount of time to record data before the trigger was activated, Delay: Time delay of trigger frame position (left side) relative to the trigger edge. For delays smaller than 0, trigger edge inside trigger frame (pre trigger). For delays greater than 0, trigger edge before trigger frame (post trigger), see Figure 2.8 .
trigger/N/duration	double	Seconds	The length of time to record data for, see Figure 2.8 .
trigger/N/edge	uint64	-	Define on which signal edge to trigger. Triggers when the trigger input signal crosses the trigger level from either low to high (edge=1), high to low (edge=2) or both (edge=3). Used for Trigger Type edge, pulse, tracking edge and tracking pulse. In the case of pulse trigger, the value specifies a positive (edge=1) or negative (edge=2) pulse relative to the trigger level (edge=3 specifies either positive or negative).
trigger/N/findlevel	uint64	-	Automatically find the value of <code>trigger/N/level</code> based on the current signal value.
trigger/N/level	uint64	Many	Specify the main trigger level value.
trigger/N/holdoff/count	uint64	-	The holdoff count, the number of skipped triggers until the next trigger is recorded again.

Setting/Path	Type	Unit	Description
<code>trigger/N/holdoff/time</code>	double	Seconds	The holdoff time, the amount of time until the next trigger is recorded again.
<code>trigger/N/hwtrigsource</code>	uint64	-	UHFLI and MFLI only, For HW Triggers only. Specify which of the four hardware trigger channels to trigger on.
<code>trigger/N/hysteresis</code>	double	Many	Specify the hysteresis value (the trigger is re-armed after the signal exceeds <code>trigger/N/level</code> and then falls below <code>trigger/N/hysteresis</code> , if using positive edge).
<code>trigger/N/path</code>	uint64	-	The path to the demod sample to trigger on, e.g., <code>demods/3/sample</code> , see also <code>trigger/N/source</code>
<code>trigger/N/pulse/max</code>	double	-	Only for Pulse triggers: The maximum pulse width to trigger on. See Figure 2.9 .
<code>trigger/N/pulse/min</code>	double	-	Only for Pulse triggers: The minimum pulse width to trigger on. See Figure 2.9 .
<code>trigger/N/retrigger</code>	uint64	-	1=enable, 0=disable. Enable to allow retriggering within one trigger duration. If enabled continue recording data in one segment if another trigger comes within the previous trigger's duration. If disabled the triggers will be recorded as separate events.
<code>trigger/N/source</code>	uint64	-	The value of <code>trigger/N/path</code> to use (0=X, 1=Y, 2=R, 3=angle, 4=freq, 5=phase, 6=auxin0, 7=auxin1).
<code>trigger/N/type</code>	uint64	-	The trigger type, see Table 2.5

Note

For the pulse trigger type, there is a subtle difference between the way the trigger level and the hysteresis are used for positive/negative pulse triggering (`trigger/N/edge= 1` or `2`) and both (`trigger/N/edge= 3`). The difference can be seen in [Figure 2.9](#) and [Figure 2.10](#).

2.5. Device Settings Module

The Device Settings Module provides functionality for saving and loading device settings to and from file. The file is saved in [XML](#) format.

In general, users are recommended to use the utility functions provided by the APIs instead of using the Device Settings module directly. The Matlab API provides `ziSaveSettings()` and `ziLoadSettings()` and the Python API provides `zhinst.utils.save_settings()` and `zhinst.utils.load_settings`. These are convenient wrappers to the Device Settings module for loading settings asynchronously, i.e., these functions block until loading or saving has completed, the desired behavior in most cases. Advanced users can use the Device Settings module directly if they need to implement loading or saving a synchronously (non-blocking).

See [Table 2.7](#) for the input parameters to configure the Device Settings Module.

Table 2.7. Device Settings Input Parameters

Setting/Path	Type	Description
<code>devicesettings/device</code>	byte array	The device ID to save the settings for, e.g., <code>dev123</code> (compulsory parameter).
<code>devicesettings/command</code>	byte array	The command to issue: 'load' (load settings from file); 'save' (read device settings and save to file) or 'read' (just read the device settings) (compulsory parameter).
<code>devicesettings/filename</code>	byte array	The name of the file to load or save to.
<code>devicesettings/path</code>	byte array	The path containing the file to load from or save to.

Table 2.8. Device Settings Parameters for use only by the LabOne Web Server.

Setting/Path	Type	Description
<code>devicesettings/throwonerror</code>	uint64	Throw an exception if there was error executing the command.
<code>devicesettings/errortext</code>	byte array	The error text used in error messages.
<code>devicesettings/finished</code>	uint64	The status of the command (read-only).

2.6. PLL Advisor Module

The PLL Advisor Module corresponds to the PLL Advisor section of the LabOne User Interface PLL tab. The PLL Advisor is a mathematical model of the PLL incorporated in the instrument and provides a convenient way to tune parameters to obtain an optimal feedback loop performance for the desired application.

Note

Note the PLL Advisor Module is only available for UHF Lock-in Amplifiers.

Table 2.9. PLL Advisor Parameters.

Setting/Path	Type	Unit	Description
<code>pllAdvisor/bode</code>	struct	-	Output parameter. Contains the resulting bode plot of the PLL simulation.
<code>pllAdvisor/calculate</code>	uint64	-	Issues a command for the PLL Advisor to calculate values. Set the value to 1 to start the calculation.
<code>pllAdvisor/center</code>	double	Hz	Center frequency of the PLL oscillator. The PLL frequency shift is relative to this center frequency.
<code>pllAdvisor/d</code>	double	Hz/deg s	The PID differential gain.
<code>pllAdvisor/demodbw</code>	double	Hz	The demodulator bandwidth to use for the PLL loop filter.
<code>pllAdvisor/i</code>	double	Hz/deg/s	The PID integral gain
<code>pllAdvisor/mode</code>	uint64	-	Sets the PLL operating mode. Currently only open-loop mode is supported.
<code>pllAdvisor/order</code>	uint64	-	Demodulator filter order to use for the PLL loop filter.
<code>pllAdvisor/p</code>	double	Hz/deg	The PID proportional gain.
<code>pllAdvisor/pllbw</code>	double	Hz	The demodulator bandwidth to use for the PLL loop filter.
<code>pllAdvisor/pm</code>	double	deg	Output parameter. Simulated phase margin of the PLL with the current settings. The phase margin should be greater than 45 deg and preferably greater than 65 deg for stable conditions.
<code>pllAdvisor/pmfreq</code>	double	-	Output parameter. Simulated phase margin frequency.
<code>pllAdvisor/q</code>	double	-	Quality factor. Currently not used.
<code>pllAdvisor/rate</code>	double	Hz	PLL Advisor sampling rate of the PLL control loop.
<code>pllAdvisor/stable</code>	bool	-	Output parameter. When 1, the PLL Advisor found a stable solution with the given settings. When 0, revise your settings and rerun the PLL Advisor.

Setting/Path	Type	Unit	Description
pllAdvisor/targetbw	double	Hz	Requested PLL bandwidth. Higher frequencies may need manual tuning.
pllAdvisor/targetfail	bool	-	Output parameter. 1 indicates the simulated PLL BW is smaller than the Target BW.

2.7. Tips and Tricks

Use the LabOne User Interface's Command Log to start programming

If you use the LabOne User Interface to perform a measurement, you can obtain the commands sent to your instrument in the "Command Log" by clicking the "Show Log" button in the status bar at the bottom of the User Interface. Be sure to set the "Log Format" of the Command Log in the "User Interface" section of the Config Tab first: The log is available in Matlab and Python formats and can be used as a starting point for your own custom program.

Use the included examples to get started programming

Both the [LabOne Matlab API](#) and the [LabOne Python API](#) come with examples to help you get started programming. In particular, both APIs have at least one example for each of the [ziCore modules](#).

Load LabOne User Interface settings files from the APIs.

The [XML](#) files used for device settings can be loaded and saved from the LabOne User Interface or from any of the ziCore-based APIs. This means that an instrument can be conveniently configured via the LabOne User Interface and then its settings saved to file. This settings file can then be loaded via an API in order to configure an instrument for a script. See the [Section 2.5](#).

Part II. LabOne APIs

This part of the Programming Manual documents language-specific installation and usage for each of the LabOne APIs. For details of common functionality and features that are shared by all the LabOne APIs please refer to [Part I](#).

Refer to:

- [Chapter 3](#) for the [LabOne Matlab API](#) (`ziDAQ`).
- [Chapter 4](#) for the [LabOne Python API](#) (`ziPython`).
- [Chapter 5](#) for the [LabOne LabVIEW API](#).
- [Chapter 6](#) for the [LabOne C API](#) (`ziAPI`).

Chapter 3. Matlab Programming

The Mathwork's numerical computing environment [Matlab®](#) has powerful tools for data analysis and visualization and can be used to create graphical user interfaces or automatically generate reports of experimental results in various formats. LabOne's Matlab API, also known as `ziDAQ`, "Zurich Instruments Data Acquisition", enables the user to stream data from their instrument directly into Matlab allowing them to take full advantage of this powerful environment.

This chapter aims to help you get started using Zurich Instruments LabOne's Matlab API, `ziDAQ`, to control your instrument, please refer to:

- [Section 3.1](#) for help [Installing the LabOne Matlab API](#) .
- [Section 3.2](#) for help [Getting Started with the LabOne Matlab API](#) and [Running the Examples](#) .
- [Section 3.3](#) for some [LabOne Matlab API Tips and Tricks](#) .
- [Section 3.4](#) for help [Troubleshooting the LabOne Matlab API](#) .
- [Section 3.5](#) for [LabOne Matlab API \(ziDAQ\) Command Reference](#) .

Note

This section and the provided examples are no substitute for a Matlab tutorial. See either Mathworks' online [Documentation Center](#) or one of the many online resources, for example, the [Matlab Programming wikibook](#) for help to get started programming with Matlab.

3.1. Installing the LabOne Matlab API

3.1.1. Requirements

To use LabOne's Matlab API, `ziDAQ`, a Matlab 7 R2009b installation (or higher) and license on Windows or Linux is required. No additional Matlab Toolboxes are required for `ziDAQ`. Both 32-bit and 64-bit platforms are supported on both Windows and Linux. On Linux, version 6.0.11 or higher of `libstdc++.so` is required by `ziDAQ`, see [Section 3.1.3](#) for more details.

The LabOne Matlab API `ziDAQ` is included in a standard LabOne installation. Some configuration is required, however, in order to use `ziDAQ` in Matlab. The LabOne installer is available from Zurich Instruments' [download page](#) (login required).

3.1.2. Windows or Linux

Technically, no installation is required to use `ziDAQ` on either Windows or Linux; it's only necessary to add the folder containing LabOne's Matlab API library to Matlab's search path. This is done as following:

1. Start Matlab and set the "Current Folder" (current working directory) to the Matlab API folder in your LabOne installation.

On Windows this is typically:

```
C:\Program Files\Zurich Instruments\LabOne\API\MATLAB2012\
```

On Linux this is the location where you unpacked the LabOne `.tar.gz` file:

```
[PATH]/LabOne64/API/MATLAB2012/
```

2. In the Matlab Command Window, run the Matlab script `ziAddPath` located in the `MATLAB2012` directory:

```
>> ziAddPath;
```

On Windows (similar for Linux) you should see the following output in Matlab's Command Window:

```
Added ziDAQ's Driver, Utilities and Examples directories to Matlab's path  
for this session.
```

To make this configuration persistent across Matlab sessions either:

1. Run the `'pathtool'` command in the Matlab Command Window and add the following paths WITH SUBFOLDERS to the Matlab search path:

```
C:\Program Files\Zurich Instruments\LabOne\API\MATLAB2012\
```

or

2. Add the following line to your Matlab startup.m file:

```
run('C:\Program Files\Zurich Instruments\LabOne\API\MATLAB2012\ziAddPath');
```

This is sufficient configuration if you would only like to use `ziDAQ` in the current Matlab session.

3. To make this configuration persistent between Matlab sessions do either one of the next two steps (as also indicated by the output of `ziAddPath`):
 - a. Run the `pathtool` and click "Add with Subfolders". Browse to the "MATLAB2012" directory that was located above in Step 1 and click "OK".
 - b. Edit your `startup.m` to contain the line indicated in the output from Step 2 above. For more help on Matlab's `startup.m` file, type the following in in Matlab's Command Window:

```
>> docsearch('startup.m')
```
4. On Linux, perform the additional step described in [Section 3.1.3](#).
5. Verify your Matlab configuration as described in [Section 3.1.4](#).

3.1.3. Additional Configuration Required on Linux

To use `ziDAQ` in Matlab on Linux it's necessary to additionally perform the following steps:

1. Tell Matlab the location of the `ziDAQ` MEX-file by appending its location to the `LD_LIBRARY_PATH` environment variable. The MEX-file is located in the `Driver` subdirectory found in `[PATH]/LabOne64/API/MATLAB2012/` (the path used above in [Section 3.1.2](#)).

This can be performed in a terminal before starting Matlab, or set in your shell's configuration file. For the Bash shell, for example, in the `~/.bashrc`:

```
export LD_LIBRARY_PATH=${LD_LIBRARY_PATH}:"[PATH]/LabOne64/API/MATLAB2012/Driver"
```

2. Tell Matlab the location of the shared C++ library (`libstdc++.so`) it should use via the `LD_PRELOAD` environment variable. This is necessary because `ziDAQ` is compiled against a newer version of `libstdc++` than is typically included (and used by default) in a standard Matlab installation. More recent versions of `libstdc++` are readily available in Linux distributions; telling Matlab to use this one instead doesn't have any negative side effects since the standard C++ library is backwards compatible.

We proceed as following:

- a. Locate the version of `libstdc++` available within your Linux distribution, for example, on Ubuntu 12.04 LTS the location and versions are:

- `/usr/lib/i386-linux-gnu/libstdc++.so.6.0.16` (32-bit),
- `/usr/lib/x86_64-linux-gnu/libstdc++.so.6.0.16` (64-bit).

Note: Version 6.0.11 or higher is required by `ziDAQ`.

- b. Append this path to the `LD_PRELOAD` environment variable:

```
export LD_PRELOAD=${LD_PRELOAD}:"/path/to/libstdc++.so.6.0.??"
```

3. Restart Matlab if necessary.

3.1.4. Verifying Successful Matlab Configuration

In order to verify that Matlab is correctly configured to use `ziDAQ` please perform the following steps:

1. Ensure that the correct Data Server is running for your device and that your device is currently connected to the Data Server. The quickest way to check is to start the User Interface for your device, see [Section 1.1](#) for more details.
2. Check that `ziDAQ`'s `connect` command can instantiate an API session to the Data Server by using the following command (with the correct port for your instrument, cf. [Section 1.1.1](#)) in the Matlab command window:

```
■ >> ziDAQ('connect', 'localhost', 8005) % 8005 for HF2 Series
■ >> ziDAQ('connect', 'localhost', 8004) % 8004 for UHFLI
■ >> ziDAQ('connect', mf-hostname, 8004) % 8004 for MFLI (see below)
```

Note, using `'localhost'` above assumes that the Data Server is running on the same computer from which you are using Matlab. See [Section 1.1.1](#) for information about port choice and connecting to the Data Server. For MFLI instruments the hostname/IP address of the MFLI instrument must be provided (the value of `mf-hostname`), see [Section 1.1.1](#) and the Getting Started chapter of the MFLI User Manual for more information.

3. If no error is reported then Matlab is correctly configured to use `ziDAQ` - congratulations! Otherwise, please try the steps listed in [Troubleshooting the LabOne Matlab API](#).

3.2. Getting Started with the LabOne Matlab API

This section introduces the user to the LabOne Matlab API.

3.2.1. Contents of the LabOne Matlab API

Alongside the driver for interfacing with your Zurich Instruments device, the LabOne Matlab API includes many files for documentation, utility functions and examples. See the `Contents.m` file located in a LabOne Matlab API directory (see Step 1 in [Section 3.1.2](#) for its typical location) for a description of the API's sub-folders and files. Run the command:

```
>> doc('Contents')
```

in the Matlab Command Window in the LabOne Matlab API directory to access the following contents interactively in Matlab.

```
% ziDAQ : The LabOne Matlab API for interfacing with Zurich Instruments Devices
%
% FILES
%   ziAddPath - add the LabOne Matlab drivers, utilities and examples to
%               Matlab's Search Path for the current session
%   README.txt - a README briefly describing how to get started with ziDAQ
%
% DIRECTORIES
%   Driver/    - contains Matlab driver for interfacing with Zurich Instruments
%               devices
%   Utils/     - contains some utility functions for common tasks
%   Examples/  - contains examples for performing measurements on Zurich
%               Instruments devices
%
% DRIVER
%   Driver/ziDAQ.m          - ziDAQ command reference documentation.
%   Driver/ziDAQ.mex*      - ziDAQ API driver
%
% UTILS
%   Utils/ziAutoConnect    - create a connection to a Zurich Instruments
%                           server
%   Utils/ziAutoDetect     - return the ID of a connected device (if only one
%                           device is connected)
%   Utils/ziDevices        - return a cell array of connected Zurich Instruments
%                           devices
%   Utils/ziCheckPathInData - check whether a node is present in data and non-empty
%
% EXAMPLES/COMMON - Examples that will run on any Zurich Instruments Device
%   example_connect          - A simple example to demonstrate how to
%                           connect to a Zurich Instruments device
%   example_connect_config  - Connect to and configure a Zurich
%                           Instruments device
%   example_poll             - Record demodulator data using
%                           ziDAQServer's synchronous poll function
%   example_record_async     - Record data asynchronously using ziDAQ's
%                           record module
%   example_record_demod_trigger - Record demodulator data upon a rising
%                           edge trigger via ziDAQ's record module
%   example_record_digital_trigger - Record data using a digital trigger via
%                           ziDAQ's record module
%   example_save_device_settings_simple - Save and load device settings
%                           synchronously using ziDAQ's utility
%                           functions
%   example_save_device_settings_expert - Save and load device settings
%                           asynchronously with ziDAQ's
%                           devicesettings module
%   example_sweeper         - Perform a frequency sweep using ziDAQ's
%                           sweep module
```

```
% example_sweeper_rstddev_fixedbw - Perform a frequency sweep plotting the
%                               stddev in demodulator output R using
%                               ziDAQ's sweep module
% example_sweeper_two_demods - Perform a frequency sweep saving data
%                               from 2 demodulators using ziDAQ's sweep
%                               module
% example_zoomfft - Perform an FFT using ziDAQ's zoomFFT
%                               module
%
% EXAMPLES/UHF - Examples specific to the UHF Series
% uhf_example_boxcar - Record boxcar data using ziDAQServer's
%                     synchronous poll function
% uhf_example_scope - Record scope data using ziDAQServer's
%                     synchronous poll function
% uhf_example_scope_offset - Record digitizer data using ziDAQServer's
%                             synchronous poll function
%
% EXAMPLES/HF2 - Examples specific to the HF2 Series
% hf2_example_autorange - determine and set an appropriate range
%                         for a sign channel
% hf2_example_poll_hardware_trigger - Poll demodulator data in combination
%                                     with a HW trigger
% hf2_example_scope - Record scope data using ziDAQServer's
%                     synchronous poll function
% hf2_example_zsync_poll - Synchronous demodulator sample timestamps
%                           from multiple HF2s via the Zsync feature
```

Note

On Windows the MEX-file is called either `ziDAQ.mexw64` or `ziDAQ.mexw32` for 64-bit and 32-bit platforms respectively and on Linux it is called `ziDAQ.mexa64` or `ziDAQ.mexa32`. If more than one MEX-file is present, Matlab automatically selects the correct MEX-file for the current platform.

3.2.2. Using the Built-in Documentation

To access `ziDAQ`'s documentation within Matlab, type either of the following in the Matlab Command Window:

```
>> help ziDAQ
```

```
>> doc ziDAQ
```

This documentation is located in the file `MATLAB2012/Driver/ziDAQ.m`. See [Section 3.5](#), [LabOne Matlab API \(ziDAQ\) Command Reference](#) for a printer friendly version.

3.2.3. Running the Examples

Prerequisites for running the Matlab examples:

1. Matlab is configured for `ziDAQ` as described above in [Section 3.1](#).
2. The Data Server program is running and the instrument is connected to the Data Server. See [Section 1.1.1](#) and the User Manual specific to the instrument for more help connecting.

- Signal Output 1 of the instrument is connected to Signal Input 1 via a BNC cable; many of the Matlab examples measure on this hardware channel.

See [Section 3.2.1](#) for a list of available examples bundled with the LabOne Matlab API. All the examples follow the same structure and take one input argument: the device ID of the instrument they are to be ran with. For example:

```
>> example_sweeper('dev123');
```

The example should produce some output in the Matlab Command Window, such as:

```
ziDAQ version Jul 7 2015 accessing server localhost 8005.  
Will run the example on `dev123`, an `HF2LI` with options `MFK|PLL|MOD|RTK|PID`.  
Sweep progress 9%  
Sweep progress 19%  
Sweep progress 30%  
Sweep progress 42%  
Sweep progress 52%  
Sweep progress 58%  
Sweep progress 68%  
Sweep progress 79%  
Sweep progress 91%  
Sweep progress 100%  
ziDAQ: AtExit called
```

Most examples will also plot some data in a Matlab figure, see [Figure 3.1](#) for an example. If you encounter an error message please ensure that the [above prerequisites](#) are fulfilled and see [Section 3.4](#) for help troubleshooting the error.




Figure 3.1. The plot produced by the LabOne Matlab API example `example_sweeper.m`; the plots show the instruments demodulator output when performing a frequency sweep over a simple feedback cable.

Note

The examples serve as a starting point for your own measurement needs. However, before editing the m-files, be sure to copy them to your own user space (they could be overwritten upon updating your LabOne installation) and give them a unique name to avoid name conflicts in Matlab.

3.2.4. Using ziCore Modules in the LabOne Matlab API

In the LabOne Matlab API ziCore [Modules](#) are configured and controlled via Matlab "handles". For example, in order to use the [Sweeper Module](#) a handle is created via:

```
>> timeout_milliseconds = 500;  
>> h = ziDAQ('sweep', timeout_milliseconds);
```

and the Module's parameters are configured using the `set` command and specifying the Module's handle with a path, value pair, for example:

```
>> ziDAQ('set', h, 'sweep/start', 1.2e5);
```

The parameters can be read-back using the `get` command, which supports wildcards, for example:

```
>> sweep_params = ziDAQ('get', h, 'sweep/*');
```

The variable `sweep_params` now contains a struct of all the Sweeper's parameters. The other main Module commands are used similarly, e.g., `ziDAQ('execute', h)` to start the sweeper. See [Section 2.1.2](#) for more help with Modules and a description of their parameters.

3.3. LabOne Matlab API Tips and Tricks

In this section some tips and tricks for working with the LabOne Matlab API are provided.

The structure of **ziDAQ** commands.

All LabOne Matlab API commands are based on a call to the Matlab function `ziDAQ()`. The first argument to `ziDAQ()` specifies the API command to be executed and is an obligatory argument. For example, a session is instantiated between the API and the Data Server with the Matlab command `ziDAQ('connect')`. Depending on the type of command specified, optional arguments may be required. For example, to obtain an integer node value, the node path must be specified as a second argument to the `'getInt'` command:

```
s = ziDAQ('getInt','/dev123/sigouts/0/on');
```

where the output argument contains the current value of the specified node.

To set an integer node value, both the node path and the value to be set must be specified as the second and third arguments:

```
ziDAQ('setInt','/dev123/sigouts/0/on', 1);
```

See the [LabOne Matlab API \(ziDAQ\) Command Reference](#) for a list of all available commands.

Data Structures returned by **ziDAQ**.

The output arguments that `ziDAQ` returns are designed to use the native data structures that Matlab users are familiar with and that reflect the data's location in the instruments node hierarchy. For example, when the `poll` command returns data from the instruments fourth demodulator (located in the node hierarchy as `/dev123/demods/3/sample`), the output argument contains a nested `struct` in which the data can be accessed by

```
data = ziDAQ('poll', poll_length, poll_timeout);  
x = data.dev123.demods(4).sample.x;  
y = data.dev123.demods(4).sample.y;
```

The instrument's node tree uses zero-based indexing; Matlab uses one-based indexing.

See the tip [Data Structures returned by ziDAQ](#) : The **fourth** demodulator sample located at `/dev123/demods/3/sample`, is indexed in the data structure returned by `poll` as `data.dev123.demods(4).sample`.

Explicitly convert **uint64** data types to **double**.

Matlab's native data type is double-precision floating point and doesn't support performing calculations with other data types such as 64-bit unsigned integers, for example:

```
>> a = uint64(2); b = uint64(1); a - b  
??? Undefined function or method 'minus' for input arguments of type 'uint64'.
```

Due to this limitation, be sure to convert demodulator timestamps to `double` before performing calculations. For example, in the following, both `clockbase` and `timestamp` (both 64-bit unsigned

integers) need to be converted to double before converting the timestamps from the instrument's native "ticks" to seconds via the instrument's clockbase:

```
data = ziDAQ('poll', 1.0, 500);           % poll data
sample = data.(device).demods(0).sample; % get the sample from the zeroth demod
% convert timestamps from ticks to seconds via the device's clockbase
% (the ADC's sampling rate), specify reference start time via t0.
clockbase = double(ziDAQ('getInt',['/' device '/clockbase']));
t = (double(sample.timestamp) - double(sample.timestamp(1)))/clockbase;
```

Use the utility function **ziCheckPathInData**.

Checking that a sub-structure in the nested data structure returned by `poll` actually exists can be cumbersome and can require multiple nested `if` statements; this can be avoided by using the utility function `ziCheckPathInData`. For example, the code:

```
data = ziDAQ('poll', poll_length, poll_timeout );
if isfield(data,device)
    if isfield(data.(device),'demods')
        if length(data.(device).demods) >= channel
            if ~isempty(data.(device).demods(channel).sample)
                % do something with the demodulator sample...
```

can be replaced by:

```
data = ziDAQ('poll', poll_length, poll_timeout );
if ziCheckPathInData( data, ['/' device '/demods/' demod_c '/sample']);
    % do something with the demodulator sample...
```

3.4. Troubleshooting the LabOne Matlab API

This section intends to solve possible error messages than can occur when using `ziDAQ` in Matlab.

Error message: "Undefined function or method 'ziDAQ' for input arguments of type '*'"

Matlab can not find the LabOne Matlab API library. Check whether the `MATLAB2012/Driver` subfolder of your LabOne installation is in the Matlab Search Path by using the command:

```
>> path
```

and repeating the steps to configure Matlab's search path in [Section 3.1.2](#).

Error message: "Undefined function or method 'example_sweeper'"

Matlab can not find the example. Check whether the `MATLAB2012/Examples/Common` subfolder (respectively `MATLAB2012/Examples/UHF` or `MATLAB2012/Examples/HF2`) of your LabOne installation are in the Matlab Search Path by using the command:

```
>> path
```

and repeating the steps to configure Matlab's search path in [Section 3.1.2](#).

Error message: "Error using: ziDAQ ZIAPIException with status code: 32870. Connection invalid."

The Matlab API can not connect to the Data Server. Please check that the correct port was used; that the correct server is running for your device and that the device is connected to the server, see [Section 1.1.1](#).

Error Message: "Error using: ziAutoConnect at 63 ziAutoConnect(): failed to find a running server or failed to find a connected device..."

The utility function `ziAutoConnect()` located in `MATLAB2012/Utils/` tries to determine which Data Server is running and whether any devices are connected to that Data Server. It is only supported by UHFLI and HF2 Series instruments, MFLI instruments are not supported. Some suggestions to verify the problem:

- Please verify in the User Interface, whether a device is connected to the Data Server running on your computer.
- If the Data Server is running on a different computer, connect manually to the Data Server via `ziDAQ`'s `connect` function:

```
>> ziDAQ('connect', hostname, port);
```

where `hostname` should be replaced by the IP of the computer the Data Server is running on and `port` is specified as in [Section 1.1.1](#).

Error Message: "Error using: ziDAQ ZIAPIException on path /dev123/sigins/0/imp50 with status code: 16387. Value or Node not found"

The API is connected to the Data Server, but the command failed to find the specified node. Please:

- Check whether your instrument is connected to the Data Server in the User Interface; if it is not connected the instruments device node tree, e.g., /dev123/, will not be constructed by the Data Server.
- Check whether the node path is spelt correctly.
- Explore the node tree to verify the node actually exists with the `listNodes` command:

```
>> ziDAQ('listNodes', '/dev123/sigins/0', 3)
```

Error Message: "using: ziDAQ Server not connected. Use 'ziDAQ('connect', ...) first."

A `ziDAQ` command was issued before initializing a connection to the Data Server. First use the `connect` command:

```
>> ziDAQ('connect', hostname, port);
```

where `hostname` should be replaced by the IP address of the computer the Data Server is running on and `port` is specified as in [Section 1.1.1](#). If the Data Server is running on the same computer, use `'localhost'` as the `hostname`.

Error Message: "Attempt to execute SCRIPT ziDAQ as a function: ziDAQ.m"

There could be a problem with your LabOne Matlab API installation. The call to `ziDAQ()` is trying to call the help file `ziDAQ.m` as a function instead of calling the `ziDAQ()` function defined in the MEX-file. In this case you need to ensure that the `ziDAQ` MEX-file is in your search path as described in [Section 3.1](#) and navigate away from from the `Driver` directory. Secondly, ensure that the LabOne Matlab MEX-file is in the `Driver` folder as described in [Section 3.2.1](#).

3.5. LabOne Matlab API (ziDAQ) Command Reference

```
%
% Copyright 2009-2015, Zurich Instruments Ltd, Switzerland
% This software is a preliminary version. Function calls and
% parameters may change without notice.
%
% This version is linked against Matlab 7.9.0.529 (R2009b) libraries.
% In case of incompatibility with the currently used Matlab version. Send
% us the requested version (use 'ver' command).
%
% ziDAQ is a interface for communication with the ZI server.
% Usage: ziDAQ([command], [option1], [option2])
%     [command] = 'clear', 'connect', 'connectDevice',
%                 'disconnectDevice', 'finished', 'flush', 'get',
%                 'getAsEvent', 'getAuxInSample', 'getBytes',
%                 'getDIO', 'getDouble', 'getInt',
%                 'getSample', 'listNodes', 'logOn', 'logOff',
%                 'poll', 'pollEvent', 'programRT', 'progress', 'read',
%                 'record', 'setByte', 'setDouble', 'syncSetDouble',
%                 'setInt', 'syncSetInt', 'subscribe',
%                 'sweep', 'trigger', 'unsubscribe', 'update',
%                 'zoomFFT', 'deviceSettings'
%
% Preconditions: ZI Server must be running (check task manager)
%
%     ziDAQ('connect', [host = '127.0.0.1'], [port = 8005], [apiLevel = 1]);
%     [host] = Server host string (default is the localhost)
%     [port] = Port number (double)
%             Use port 8005 to connect to the HF2 Data Server
%             Use port 8004 to connect to the UHF Data Server
%     [apiLevel] = Compatibility mode of the API interface
%                 Use API level 1 to use code written for HF2.
%                 Higher API levels are currently only supported
%                 for UHF devices. To get full functionality for
%                 an UHF device use API level 5.
%     To disconnect use 'clear ziDAQ'
%
%     result = ziDAQ('getConnectionAPILevel');
%     Returns ziAPI level used for the active connection.
%
%     ziDAQ('connectDevice', [device], [interface]);
%     [device] = Device serial string to connect (e.g. 'DEV2000')
%     [interface] = Interface string e.g. 'USB', '1GbE', '10GbE'
%     Connect with the data server to a specified device over the
%     specified interface. The device must be visible to the server.
%     If the device is already connected the call will be ignored.
%     The function will block until the device is connected and
%     the device is ready to use. This method is useful for UHF
%     devices offering several communication interfaces.
%
%     ziDAQ('disconnectDevice', [device]);
%     [device] = Device serial string of device to disconnect.
%     This function will return immediately. The disconnection of
%     the device may not yet finished.
%
%     result = ziDAQ('listNodes', [path], [flags(int64) = 0]);
%     [path] = Path string
%     [flags] = int64(0) -> ZI_LIST_NONE 0x00
%             The default flag, returning a simple
%             listing if the given node
%             int64(1) -> ZI_LIST_RECURSIVE 0x01
%             Returns the nodes recursively
```

```
%          int64(2) -> ZI_LIST_ABSOLUTE 0x02
%          Returns absolute paths
%          int64(4) -> ZI_LIST_LEAFSONLY 0x04
%          Returns only nodes that are leafs,
%          which means the they are at the
%          outermost level of the tree.
%          int64(8) -> ZI_LIST_SETTINGSONLY 0x08
%          Returns only nodes which are marked
%          as setting
%          Or combinations of flags might be used.
%
% result = ziDAQ('getSample', [path]);
%          [path] = Path string
%          Returns a single demodulator sample (including
%          DIO and AuxIn). For more efficient data recording
%          use subscribe and poll functions.
%
% result = ziDAQ('getAuxInSample', [path]);
%          [path] = Path string
%          Returns a single auxin sample. The auxin data
%          is averaged in contrast to the auxin data embedded
%          in the demodulator sample.
%
% result = ziDAQ('getDIO', [path]);
%          [path] = Path string
%          Returns a single DIO sample.
%
% result = ziDAQ('getDouble', [path]);
%          [path] = Path string
%
% result = ziDAQ('getInt', [path]);
%          [path] = Path string
%
% result = ziDAQ('getBytes', [path]);
%          [path] = Path string
%
%          ziDAQ('setDouble', [path], [value(double)]);
%          [path] = Path string
%          [value] = Setting value
%
%          ziDAQ('syncSetDouble', [path], [value(double)]);
%          [path] = Path string
%          [value] = Setting value
%
%          ziDAQ('setInt', [path], [value(int64)]);
%          [path] = Path string
%          [value] = Setting value
%
%          ziDAQ('syncSetInt', [path], [value(int64)]);
%          [path] = Path string
%          [value] = Setting value
%
%          ziDAQ('setByte', [path], [value(uint8)]);
%          [path] = Path string
%          [value] = Setting value
%
%          ziDAQ('subscribe', [path]);
%          [path] = Path string
%          Subscribe to the specified path to receive streaming data
%          or setting data if changed over the poll or pollEvent
%          commands.
%
%          ziDAQ('unsubscribe', [path]);
%          [path] = Path string
%
%          ziDAQ('getAsEvent', [path]);
%          [path] = Path string
```

```
% Triggers a single event on the path to return the current
% value. The result can be fetched with the poll or pollEvent
% command.
%
% ziDAQ('update');
% Detect HF2 devices connected to the USB. On Windows this
% update is performed automatically.
%
% ziDAQ('get', [path]);
% [path] = Path string
% Gets a structure of the node data from the specified
% branch. High-speed streaming nodes (e.g. /devN/demods/0/sample)
% are not returned. Wildcards (*) may be used, in which case
% read-only nodes are ignored.
%
% ziDAQ('flush');
% Flush all data in the socket connection and API buffers.
% Call this function before a subscribe with subsequent poll
% to get rid of old streaming data that might still be in
% the buffers.
%
% ziDAQ('echoDevice', [device]);
% [device] = device string e.g. 'dev100'
% Sends an echo command to a device and blocks until
% answer is received. This is useful to flush all
% buffers between API and device to enforce that
% further code is only executed after the device executed
% a previous command.
%
% ziDAQ('sync');
% Synchronize all data path. Ensures that get and poll
% commands return data which was recorded after the
% setting changes in front of the sync command. This
% sync command replaces the functionality of all syncSet,
% flush, and echoDevice commands.
%
% ziDAQ('programRT', [device], [filename]);
% [device] = device string e.g. 'dev100'
% [filename] = filename of RT program
% Writes down the RT program. To use this function
% the RT option must be available for the specified
% device.
%
% result = ziDAQ('secondsTimeStamp', [timestamps]);
% [timestamps] = vector of uint64 ticks
% Deprecated. In order to convert timestamps to seconds divide the
% timestamps by the value instrument's clockbase device node,
% e.g., /dev99/clockbase.
% [Converts a timestamp vector of uint64 ticks
% into a double vector of timestamps in seconds (HF2 Series).]
%
% Synchronous Interface
%
% ziDAQ('poll', [duration(double)], [timeout(int64)], [flags(uint32)]);
% [duration] = Recording time in [s]
% [timeout] = Poll timeout in [ms]
% [flags] = Flags that specify data polling properties
% Bit[0] FILL : Fill data loss holes
% Bit[1] ALIGN : Align data of several demodulators
% Bit[2] THROW : Throw if data loss is detected
% Records data for the specified time. This function call
% is blocking. Use the asynchronous interface for long
% recording durations.
%
% result = ziDAQ('pollEvent', [timeout(int64)]);
% [timeout] = Poll timeout in [ms]
% Just execute a single poll command. This is a low-level
```



```

%           function. The poll function is better suited for most
%           cases.
%
% Asynchronous Interface
%
%   Trigger Parameters
%   trigger/buffersize      double Overwrite the buffersize [s] of the trigger
%                               object (set when it was instantiated).
%                               The recommended buffer size is
%                               2*trigger/0/duration.
%   trigger/device          string The device ID to execute the software trigger,
%                               e.g. dev123 (compulsory parameter).
%   trigger/endless         bool  Enable endless triggering 1=enable; 0=disable.
%   trigger/forcetrigger    bool  Force a trigger.
%   trigger/0/path          string The path to the demod sample to trigger on,
%                               e.g. demods/3/sample, see also trigger/0/source
%   trigger/0/source        int   Signal that is used to trigger on.
%                               0 = x [X_SOURCE]
%                               1 = y [Y_SOURCE]
%                               2 = r [R_SOURCE]
%                               3 = angle [ANGLE_SOURCE]
%                               4 = frequency [FREQUENCY_SOURCE]
%                               5 = phase [PHASE_SOURCE]
%                               6 = auxiliary input 0 [AUXIN0_SOURCE]
%                               7 = auxiliary input 1 [AUXIN1_SOURCE]
%   trigger/0/count         int   Number of trigger edges to record.
%   trigger/0/type          int   Trigger type used. Some parameters are
%                               only valid for special trigger types.
%                               0 = trigger off
%                               1 = analog edge trigger on source
%                               2 = digital trigger mode on DIO
%                               3 = analog pulse trigger on source
%                               4 = analog tracking trigger on source
%   trigger/0/edge          int   Trigger edge
%                               1 = rising edge
%                               2 = falling edge
%                               3 = both
%   trigger/0/findlevel     bool  Automatically find the value of trigger/0/level
%                               based on the current signal value.
%   trigger/0/bits          int   Digital trigger condition.
%   trigger/0/bitmask       int   Bit masking for bits used for
%                               triggering. Used for digital trigger.
%   trigger/0/delay         double Trigger frame position [s] (left side)
%                               relative to trigger edge.
%                               delay = 0 -> trigger edge at left border.
%                               delay < 0 -> trigger edge inside trigger
%                               frame (pretrigger).
%                               delay > 0 -> trigger edge before trigger
%                               frame (posttrigger).
%   trigger/0/duration      double Recording frame length [s]
%   trigger/0/level         double Trigger level voltage [V].
%   trigger/0/hysteresis    double Trigger hysteresis [V].
%   trigger/0/retrigger     int   Record more than one trigger in a
%   trigger/triggered       bool  Has the software trigger triggered? 1=Yes, 0=No
%                               (read only).
%   trigger/0/bandwidth     double Filter bandwidth [Hz] for pulse and
%                               tracking triggers.
%   trigger/0/holdoff/count int   Number of skipped triggers until the
%                               next trigger is recorded again.
%   trigger/0/holdoff/time  double Hold off time [s] before the next
%                               trigger is recorded again. A hold off
%                               time smaller than the duration will
%                               produce overlapped trigger frames.
%   trigger/N/hwtrigsources int   Only available for devices that support
%                               hardware triggering. Specify the channel
%                               to trigger on.
%   trigger/0/pulse/min     double Minimal pulse width [s] for the pulse

```

```

%
% trigger/0/pulse/max      double Maximal pulse width [s] for the pulse
%                           trigger.
% trigger/filename        string This parameter is deprecated. If specified,
%                           i.e., not empty, it enables automatic saving of
%                           data in single trigger mode
%                           (trigger/endless = 0).
% trigger/savepath        string The directory where files are saved when saving
%                           data.
% trigger/fileformat      string The format of the file for saving data.
%                           0=Matlab, 1=CSV.
% trigger/historylength   bool   Maximum number of entries stored in the
%                           measurement history.
% triggerclearhistory     bool   Remove all records from the history list.
%
% handle = ziDAQ('record' [duration(double)], [timeout(int64)]);
%                           [duration] = Recording time in [s]
%                           [timeout] = Poll timeout in [ms]
%                           Creates a recorder class. The thread is not yet started.
%                           Before the thread start subscribe and set command have
%                           to be called. To start the real measurement use the
%                           execute function. After that the trigger will start
%                           the recording of a frame.
%
% result = ziDAQ('listNodes', [handle], [path], [flags(int64) = 0]);
%                           [path] = Path string
%                           [flags] = int64(0) -> ZI_LIST_NONE 0x00
%                                   The default flag, returning a simple
%                                   listing if the given node
%                           int64(1) -> ZI_LIST_RECURSIVE 0x01
%                                   Returns the nodes recursively
%                           int64(2) -> ZI_LIST_ABSOLUTE 0x02
%                                   Returns absolute paths
%                           int64(4) -> ZI_LIST_LEAFSONLY 0x04
%                                   Returns only nodes that are leafs,
%                                   which means the they are at the
%                                   outermost level of the tree.
%                           int64(8) -> ZI_LIST_SETTINGSONLY 0x08
%                                   Returns only nodes which are marked
%                                   as setting
%                           Or combinations of flags might be used.
%
% ziDAQ('subscribe', [handle], [path]);
%                           Subscribe to one or several nodes. After subscription
%                           the recording process can be started with the 'execute'
%                           command. During the recording process paths can not be
%                           subscribed or unsubscribed.
%                           [handle] = Reference to the ziDAQRecorder class.
%                           [path] = Path string of the node. Use wild card to
%                           select all. Alternatively also a list of path
%                           strings can be specified.
%
% ziDAQ('unsubscribe', [handle], [path]);
%                           Unsubscribe from one or several nodes. During the
%                           recording process paths can not be subscribed or
%                           unsubscribed.
%                           [handle] = Reference to the ziDAQRecorder class.
%                           [path] = Path string of the node. Use wild card to
%                           select all. Alternatively also a list of path
%                           strings can be specified.
%
% ziDAQ('get', [handle], [path]);
%                           [handle] = Reference to the ziDAQRecorder class.
%                           [path] = Path string of the node.
%
% ziDAQ('set', [handle], [path], [value]);
%                           [handle] = Reference to the ziDAQRecorder class.

```

```

%           [path] = Path string of the node.
%
%           ziDAQ('set', [handle], [path], [value]);
%           [handle] = Reference to the ziDAQRecorder class.
%           [path] = Path string of the node.
%
%           ziDAQ('set', [handle], [path], [value]);
%           [handle] = Reference to the ziDAQRecorder class.
%           [path] = Path string of the node.
%
%           ziDAQ('execute', [handle]);
%           Start the recorder. After that command any trigger will
%           start the measurement. Subscription or unsubscription
%           is no more possible until the recording is finished.
%
%           ziDAQ('trigger', [handle]);
%           [handle] = Handle of the recording session.
%           Triggers the measurement recording.
%
%           result = ziDAQ('finished', [handle]);
%           [handle] = Handle of the recording session.
%           Returns 1 if the recording is finished, otherwise 0.
%
%           result = ziDAQ('read', [handle]);
%           [handle] = Handle of the recording session.
%           Transfer the recorded data to Matlab.
%
%           ziDAQ('finish', [handle]);
%           Stop recording. The recording may be restarted by
%           calling 'execute' again.
%
%           result = ziDAQ('progress', [handle]);
%           Report the progress of the measurement with a number
%           between 0 and 1.
%
%           ziDAQ('clear', [handle]);
%           [handle] = Handle of the recording session.
%           Stop recording.
%
% Sweep Module
%
% Sweep Parameters
%   sweep/device      string  Device that should be used for
%                           the parameter sweep, e.g. 'dev99'.
%   sweep/start       double  Sweep start frequency [Hz]
%   sweep/stop        double  Sweep stop frequency [Hz]
%   sweep/gridnode    string  Path of the node that should be
%                           used for sweeping. For frequency
%                           sweep applications this will be e.g.
%                           'oscs/0/freq'. The device name of
%                           the path can be omitted and is given
%                           by sweep/device.
%   sweep/loopcount   int     Number of sweep loops (default 1)
%   sweep/endless     int     Endless sweeping (default 0)
%                           0 = Use loopcount value
%                           1 = Endless sweeping enabled, ignore
%                           loopcount
%   sweep/samplecount int     Number of samples per sweep
%   sweep/settling/time double Settling time before measurement is
%                           performed, in [s]
%   sweep/settling/tc double  Settling precision
%                           5 ~ low precision
%                           15 ~ medium precision
%                           50 ~ high precision
%   sweep/settling/inaccuracy int Demodulator filter settling inaccuracy defining
%                           the wait time between a sweep parameter change
%                           and recording of the next sweep point. Typical

```

```

%                               inaccuracy values: 10m for highest sweep speed
%                               for large signals, 100u for precise amplitude
%                               measurements, 100n for precise noise
%                               measurements. Depending on the order the
%                               settling accuracy will define the number of
%                               filter time constants the sweeper has to
%                               wait. The maximum between this value and the
%                               settling time is taken as wait time until the
%                               next sweep point is recorded.
%   sweep/xmapping             int    Sweep mode
%                               0 = linear
%                               1 = logarithmic
%   sweep/scan                 int    Scan type
%                               0 = sequential
%                               1 = binary
%                               2 = bidirectional
%                               3 = reverse
%   sweep/bandwidth            double Fixed bandwidth [Hz]
%                               0 = Automatic calculation (obsolete)
%   sweep/bandwidthcontrol    int    Sets the bandwidth control mode (default 2)
%                               0 = Manual (user sets bandwidth and order)
%                               1 = Fixed (uses fixed bandwidth value)
%                               2 = Auto (calculates best bandwidth value)
%                               Equivalent to the obsolete bandwidth = 0
%                               setting
%   sweep/order               int    Defines the filter roll off to use in Fixed
%                               bandwidth selection.
%                               Valid values are between 1 (6 dB/octave)
%                               and 8 (48 dB/octave). An order of 0
%                               triggers a read-out of the order from the
%                               selected demodulator.
%   sweep/maxbandwidth        double Maximal bandwidth used in auto bandwidth
%                               mode in [Hz]. The default is 1.25MHz.
%   sweep/omegasuppression    double Damping in [dB] of omega and 2omega components.
%                               Default is 40dB in favor of sweep speed.
%                               Use higher value for strong offset values or
%                               3omega measurement methods.
%   sweep/averaging/tc        double Min averaging time [tc]
%                               0 = no averaging (see also time!)
%                               5 ~ low precision
%                               15 ~ medium precision
%                               50 ~ high precision
%   sweep/averaging/sample    int    Min samples to average
%                               1 = no averaging (if averaging/tc = 0)
%   sweep/phaseunwrap         bool   Enable unwrapping of slowly changing phase
%                               evolutions around the +/-180 degree boundary.
%   sweep/sincfilter          bool   Enables the sinc filter if the sweep frequency
%                               is below 50 Hz. This will improve the sweep
%                               speed at low frequencies as omega components
%                               do not need to be suppressed by the normal
%                               low pass filter.
%   sweep/filename            string This parameter is deprecated. If specified,
%                               i.e. not empty, it enables automatic saving of
%                               data in single sweep mode (sweep/endless = 0).
%   sweep/savepath            string The directory where files are located when
%                               saving sweeper measurements.
%   sweep/fileformat          string The format of the file for saving sweeper
%                               measurements. 0=Matlab, 1=CSV.
%   sweep/historylength       bool   Maximum number of entries stored in the
%                               measurement history.
%   sweep/clearhistory        bool   Remove all records from the history list.
%
%   Note:
%   Settling time = max(settling.tc * tc, settling.time)
%   Averaging time = max(averaging.tc * tc, averaging.sample / sample-rate)
%
%   handle = ziDAQ('sweep', [timeout(int64)]);

```

```
% [timeout] = Poll timeout in [ms]
% Creates a sweep class. The thread is not yet started.
% Before the thread start subscribe and set command have
% to be called. To start the real measurement use the
% execute function.
%
result = ziDAQ('listNodes', [handle], [path], [flags(int64) = 0]);
% [path] = Path string
% [flags] = int64(0) -> ZI_LIST_NONE 0x00
%           The default flag, returning a simple
%           listing if the given node
% int64(1) -> ZI_LIST_RECURSIVE 0x01
%           Returns the nodes recursively
% int64(2) -> ZI_LIST_ABSOLUTE 0x02
%           Returns absolute paths
% int64(4) -> ZI_LIST_LEAFONLY 0x04
%           Returns only nodes that are leafs,
%           which means the they are at the
%           outermost level of the tree.
% int64(8) -> ZI_LIST_SETTINGSONLY 0x08
%           Returns only nodes which are marked
%           as setting
% Or combinations of flags might be used.
%
ziDAQ('subscribe', [handle], [path]);
% Subscribe to one or several nodes. After subscription
% the recording process can be started with the 'execute'
% command. During the recording process paths can not be
% subscribed or unsubscribed.
% [handle] = Reference to the ziDAQSweeper class.
% [path] = Path string of the node. Use wild card to
% select all. Alternatively also a list of path
% strings can be specified.
%
ziDAQ('unsubscribe', [handle], [path]);
% Unsubscribe from one or several nodes. During the
% recording process paths can not be subscribed or
% unsubscribed.
% [handle] = Reference to the ziDAQSweeper class.
% [path] = Path string of the node. Use wild card to
% select all. Alternatively also a list of path
% strings can be specified.
%
ziDAQ('execute', [handle]);
% Start the sweep. Subscription or unsubscription
% is no more possible until the sweep is finished.
%
result = ziDAQ('finished', [handle]);
% [handle] = Handle of the sweep session.
% Returns 1 if the sweep is finished, otherwise 0.
%
result = ziDAQ('read', [handle]);
% [handle] = Handle of the sweep session.
% Transfer the sweep data to Matlab.
%
result = ziDAQ('progress', [handle]);
% Report the progress of the measurement with a number
% between 0 and 1.
%
ziDAQ('finish', [handle]);
% Stop the sweep. The sweep may be restarted by
% calling 'execute' again.
%
ziDAQ('clear', [handle]);
% [handle] = Handle of the sweep session.
% Stop the current sweep.
%
```

```

%          ziDAQ('save', [handle]);
%          Save the measured data to a file.
%          [handle] = Handle of the sweep session.
%          [filename] = File in which to store the data.
%
% Zoom FFT Module
%
% Zoom FFT Parameters
% zoomFFT/device          string Device that should be used for
%                           the zoom FFT, e.g. 'dev99'.
% zoomFFT/bit             int    Number of FFT points 2^bit
% zoomFFT/mode            int    Zoom FFT mode
%                           0 = Perform FFT on X+iY
%                           1 = Perform FFT on R
%                           2 = Perform FFT on Phase
% zoomFFT/loopcount       int    Number of zoom FFT loops (default 1)
% zoomFFT/endless         int    Perform endless zoom FFT (default 0)
%                           0 = Use loopcount value
%                           1 = Endless zoom FFT enabled, ignore
%                               loopcount
% zoomFFT/overlap         double FFT overlap 0 = none, [0..1]
% zoomFFT/settling/time   double Settling time before measurement is performed
% zoomFFT/settling/tc     double Settling time in time constant units before
%                               the FFT recording is started.
%                               5 ~ low precision
%                               15 ~ medium precision
%                               50 ~ high precision
% zoomFFT/window          int    FFT window (default 1 = Hann)
%                               0 = Rectangular
%                               1 = Hann
%                               2 = Hamming
%                               3 = Blackman Harris 4 term
% zoomFFT/absolute        bool   Shifts the frequencies so that the center
%                               frequency becomes the demodulation frequency
%                               rather than 0 Hz.
%
% handle = ziDAQ('zoomFFT', [timeout(int64)]);
%          [timeout] = Poll timeout in [ms]
%          Creates a zoom FFT class. The thread is not yet started.
%          Before the thread start subscribe and set command have
%          to be called. To start the real measurement use the
%          execute function.
%
% result = ziDAQ('listNodes', [handle], [path], [flags(int64) = 0]);
%          [path] = Path string
%          [flags] = int64(0) -> ZI_LIST_NONE 0x00
%                  The default flag, returning a simple
%                  listing if the given node
%          int64(1) -> ZI_LIST_RECURSIVE 0x01
%                  Returns the nodes recursively
%          int64(2) -> ZI_LIST_ABSOLUTE 0x02
%                  Returns absolute paths
%          int64(4) -> ZI_LIST_LEAFSONLY 0x04
%                  Returns only nodes that are leafs,
%                  which means the they are at the
%                  outermost level of the tree.
%          int64(8) -> ZI_LIST_SETTINGSONLY 0x08
%                  Returns only nodes which are marked
%                  as setting
%          Or combinations of flags might be used.
%
%          ziDAQ('subscribe', [handle], [path]);
%          Subscribe to one or several nodes. After subscription
%          the recording process can be started with the 'execute'
%          command. During the recording process paths can not be
%          subscribed or unsubscribed.
%          [handle] = Reference to the ziDAQZoomFFT class.

```

```

%           [path] = Path string of the node. Use wild card to
%           select all. Alternatively also a list of path
%           strings can be specified.
%
%           ziDAQ('unsubscribe', [handle], [path]);
%           Unsubscribe from one or several nodes. During the
%           recording process paths can not be subscribed or
%           unsubscribed.
%           [handle] = Reference to the ziDAQZoomFFT class.
%           [path] = Path string of the node. Use wild card to
%           select all. Alternatively also a list of path
%           strings can be specified.
%
%           ziDAQ('execute', [handle]);
%           Start the zoom FFT. Subscription or unsubscription
%           is no more possible until the zoomFFT is finished.
%
%           result = ziDAQ('finished', [handle]);
%           [handle] = Handle of the zoom FFT session.
%           Returns 1 if the zoom FFT is finished, otherwise 0.
%
%           result = ziDAQ('read', [handle]);
%           [handle] = Handle of the zoom FFT session.
%           Transfer the zoomFFT data to Matlab.
%
%           result = ziDAQ('progress', [handle]);
%           Report the progress of the measurement with a number
%           between 0 and 1.
%
%           ziDAQ('finish', [handle]);
%           Stop the zoomFFT. The zoom FFT may be restarted by
%           calling 'execute' again.
%
%           ziDAQ('clear', [handle]);
%           [handle] = Handle of the zoom FFT session.
%           Stop the current zoom FFT.
%
% Device Settings Module
%
% Device Settings Parameters
%   devicesettings/device      string  Device whose settings are to be
%                                     saved/loaded, e.g. 'dev99'.
%   devicesettings/path       string  Path where the settings files are to
%                                     be located. If not set, the default
%                                     settings location of the LabOne
%                                     software is used.
%   devicesettings/filename    string  The file to which the settings are to
%                                     be saved/loaded.
%   devicesettings/command     string  The save/load command to execute.
%                                     'save' = Read device settings and save
%                                     to file.
%                                     'load' = Load settings from file and
%                                     write to device.
%                                     'read' = Read device settings only
%                                     (no save).
%
%   handle = ziDAQ('deviceSettings', [timeout(int64)]);
%   [timeout] = Poll timeout in [ms]
%   Creates a device settings class for saving/loading device
%   settings to/from a file. Before the thread start, set the path,
%   filename and command parameters. To run the command, use the
%   execute function.
%
%   result = ziDAQ('listNodes', [handle], [path], [flags(int64) = 0]);
%   [path] = Path string
%   [flags] = int64(0) -> ZI_LIST_NONE 0x00
%               The default flag, returning a simple

```

```

%           listing if the given node
%           int64(1) -> ZI_LIST_RECURSIVE 0x01
%           Returns the nodes recursively
%           int64(2) -> ZI_LIST_ABSOLUTE 0x02
%           Returns absolute paths
%           int64(4) -> ZI_LIST_LEAFSONLY 0x04
%           Returns only nodes that are leafs,
%           which means the they are at the
%           outermost level of the tree.
%           int64(8) -> ZI_LIST_SETTINGSONLY 0x08
%           Returns only nodes which are marked
%           as setting
%           Or combinations of flags might be used.
%
%           ziDAQ('subscribe', [handle], [path]);
%           Not relevant for the device settings module.
%
%           ziDAQ('unsubscribe', [handle], [path]);
%           Not relevant for the device settings module.
%
%           ziDAQ('execute', [handle]);
%           Execute the command.
%
%           result = ziDAQ('finished', [handle]);
%           [handle] = Handle of the device settings session.
%           Returns 1 if the command is finished, otherwise 0.
%
%           result = ziDAQ('read', [handle]);
%           [handle] = Handle of the device settings session.
%           Transfer the device settings to Matlab.
%           Not relevant since device settings are saved to a file.
%
%           result = ziDAQ('progress', [handle]);
%           Report the progress of the command with a number
%           between 0 and 1.
%
%           ziDAQ('finish', [handle]);
%           Stop the device settings module. The module may be restarted by
%           calling 'execute' again.
%
%           ziDAQ('clear', [handle]);
%           [handle] = Handle of the device settings session.
%           End the current device settings thread.
%
% PLL Advisor Module
%
% PLL Advisor Parameters
%   pllAdvisor/bode      struct  Output parameter. Contains the resulting bode
%                               plot of the PLL simulation.
%   pllAdvisor/calculate int     Command to calculate values. Set to 1 to start
%                               the calculation.
%   pllAdvisor/center    double  Center frequency of the PLL oscillator. The PLL
%                               frequency shift is relative to this center
%                               frequency.
%   pllAdvisor/d          int     Differential gain.
%   pllAdvisor/demodbw    int     Demodulator bandwidth used for the PLL loop
%                               filter.
%   pllAdvisor/i          double  Integral gain.
%   pllAdvisor/mode       double  Select PLL Advisor mode. Currently only one mode
%                               (open loop) is supported.
%   pllAdvisor/order      double  Demodulator order used for the PLL loop filter.
%   pllAdvisor/p          int     Proportional gain.
%   pllAdvisor/pllbw      int     Demodulator bandwidth used for the PLL loop
%                               filter.
%   pllAdvisor/pm         int     Output parameter. Simulated phase margin of the
%                               PLL with the current settings. The phase margin
%                               should be greater than 45 deg and preferably

```



```

%
%      pllAdvisor/pmfreq      int      greater than 65 deg for stable conditions.
%                                Output parameter. Simulated phase margin
%                                frequency.
%      pllAdvisor/q          int      Quality factor. Currently not used.
%      pllAdvisor/rate       int      PLL Advisor sampling rate of the PLL control
%                                loop.
%      pllAdvisor/stable     int      Output parameter. When 1, the PLL Advisor found
%                                a stable solution with the given settings. When
%                                0, revise your settings and rerun the PLL
%                                Advisor.
%      pllAdvisor/targetbw   int      Requested PLL bandwidth. Higher frequencies may
%                                need manual tuning.
%      pllAdvisor/targetfail int      Output parameter. 1 indicates the simulated PLL
%                                BW is smaller than the Target BW.
%
%      handle = ziDAQ('pllAdvisor', [timeout(int64)]);
%                                [timeout] = Poll timeout in [ms]
%                                Creates a PLL Advisor class for simulating the PLL in the
%                                device. Before the thread start, set the command parameters,
%                                call execute() and then set the "calculate" parameter to start
%                                the simulation.
%
%      result = ziDAQ('listNodes', [handle], [path], [flags(int64) = 0]);
%                                [path] = Path string
%                                [flags] = int64(0) -> ZI_LIST_NONE 0x00
%                                The default flag, returning a simple
%                                listing if the given node
%                                int64(1) -> ZI_LIST_RECURSIVE 0x01
%                                Returns the nodes recursively
%                                int64(2) -> ZI_LIST_ABSOLUTE 0x02
%                                Returns absolute paths
%                                int64(4) -> ZI_LIST_LEAFSONLY 0x04
%                                Returns only nodes that are leafs,
%                                which means the they are at the
%                                outermost level of the tree.
%                                int64(8) -> ZI_LIST_SETTINGSONLY 0x08
%                                Returns only nodes which are marked
%                                as setting
%                                Or combinations of flags might be used.
%
%      ziDAQ('subscribe', [handle], [path]);
%                                Subscribe to one or several nodes.
%
%      ziDAQ('unsubscribe', [handle], [path]);
%                                Unsubscribe from one or several nodes..
%
%      ziDAQ('execute', [handle]);
%                                Start the PLL Advisor.
%
%      result = ziDAQ('finished', [handle]);
%                                [handle] = Handle of the PLL Advisor session.
%                                Returns 1 if the command is finished, otherwise 0.
%
%      result = ziDAQ('read', [handle]);
%                                [handle] = Handle of the PLL Advisor session.
%                                Read pllAdvisor data. If the simulation is still ongoing only a
%                                subset of the data is returned.
%
%      result = ziDAQ('progress', [handle]);
%                                Report the progress of the command with a number
%                                between 0 and 1.
%
%      ziDAQ('finish', [handle]);
%                                Stop the PLL Advisor module.
%
%      ziDAQ('clear', [handle]);
%                                [handle] = Handle of the PLL Advisor session.

```

```

%                               End the current PLL Advisor thread.
%
% PID Advisor Module
%
% PID Advisor Parameters
%   pidAdvisor/advancedmode    int    Disable automatic calculation of the
%                                     start and stop value.
%   pidAdvisor/auto            int    Automatic response calculation triggered
%                                     by parameter change.
%   pidAdvisor/bode            struct  Output parameter. Contains the resulting
%                                     bode plot of the PID simulation.
%   pidAdvisor/bw              double  Output parameter. Calculated system
%                                     bandwidth.
%   pidAdvisor/calculate       int    In/Out parameter. Command to calculate
%                                     values. Set to 1 to start the
%                                     calculation.
%   pidAdvisor/display/freqstart double Start frequency for Bode plot.
%                                     For disabled advanced mode the start
%                                     value is automatically derived from the
%                                     system properties.
%   pidAdvisor/display/freqstop double Stop frequency for Bode plot.
%   pidAdvisor/display/timestart double Start time for step response.
%   pidAdvisor/display/timestop double Stop time for step response.
%   pidAdvisor/dut/bw          double  Bandwidth of the DUT (device under test).
%   pidAdvisor/dut/damping     double  Damping of the second order
%                                     low pass filter.
%   pidAdvisor/dut/delay       double  IO Delay of the feedback system
%                                     describing the earliest response for
%                                     a step change.
%   pidAdvisor/dut/fcenter     double  Resonant frequency of the of the modelled
%                                     resonator.
%   pidAdvisor/dut/gain        double  Gain of the DUT transfer function.
%   pidAdvisor/dut/q           double  quality factor of the modelled resonator.
%   pidAdvisor/dut/source      int    Type of model used for the external
%                                     device to be controlled by the PID.
%                                     source = 1: Lowpass first order
%                                     source = 2: Lowpass second order
%                                     source = 3: Resonator frequency
%                                     source = 4: Internal PLL
%                                     source = 5: VCO
%                                     source = 6: Resonator amplitude
%   pidAdvisor/impulse         struct  Output parameter. Impulse response
%                                     (not yet supported).
%   pidAdvisor/index           int    PID index for parameter detection.
%   pidAdvisor/pid/autobw      int    Adjusts the demodulator bandwidth to fit
%                                     best to the specified target bandwidth
%                                     of the full system.
%   pidAdvisor/pid/d           double  In/Out parameter. Differential gain.
%   pidAdvisor/pid/dlimittimeconstant double In/Out parameter. Differential filter
%                                     timeconstant.
%   pidAdvisor/pid/i           double  In/Out parameter. Integral gain.
%   pidAdvisor/pid/mode        double  Select PID Advisor mode. Mode value is
%                                     bit coded, bit 0: P, bit 1: I, bit 2: D,
%                                     bit 3: D filter limit.
%   pidAdvisor/pid/p           double  In/Out parameter. Proportional gain.
%   pidAdvisor/pid/rate        double  In/Out parameter. PID Advisor sampling
%                                     rate of the PID control loop.
%   pidAdvisor/pid/targetbw    double  PID system target bandwidth.
%   pidAdvisor/pm              double  Output parameter. Simulated phase margin
%                                     of the PID with the current settings.
%                                     The phase margin should be greater than
%                                     45 deg and preferably greater than 65 deg
%                                     for stable conditions.
%   pidAdvisor/pmfreq          double  Output parameter. Simulated phase margin
%                                     frequency.
%   pidAdvisor/stable          int    Output parameter. When 1, the PID Advisor

```

```

% found a stable solution with the given
% settings. When 0, revise your settings
% and rerun the PID Advisor.
% pidAdvisor/step          struct Output parameter. Contains the resulting
%                           step response plot of the PID simulation.
% pidAdvisor/targetbw     double Requested PID bandwidth. Higher
%                           frequencies may need manual tuning.
% pidAdvisor/targetfail   int    Output parameter. 1 indicates the
%                           simulated PID BW is smaller than the
%                           Target BW.
% pidAdvisor/tf/closedloop int    Switch the response calculation mode
%                           between closed or open loop.
% pidAdvisor/tf/input      int    Start point for the plant response
%                           simulation for open or closed loops.
% pidAdvisor/tf/output     int    End point for the plant response
%                           simulation for open or closed loops.
% pidAdvisor/tune          int    Optimize the PID parameters so that
%                           the noise of the closed-loop
%                           system gets minimized.
%
% handle = ziDAQ('pidAdvisor', [timeout(int64)]);
% [timeout] = Poll timeout in [ms]
% Creates a PID Advisor class for simulating the PID in the
% device. Before the thread start, set the command parameters,
% call execute() and then set the "calculate" parameter to start
% the simulation.
%
% result = ziDAQ('listNodes', [handle], [path], [flags(int64) = 0]);
% [handle] = Handle of the PID Advisor session.
% [path] = Path string
% [flags] = int64(0) -> ZI_LIST_NONE 0x00
%           The default flag, returning a simple
%           listing if the given node
% int64(1) -> ZI_LIST_RECURSIVE 0x01
%           Returns the nodes recursively
% int64(2) -> ZI_LIST_ABSOLUTE 0x02
%           Returns absolute paths
% int64(4) -> ZI_LIST_LEAFSONLY 0x04
%           Returns only nodes that are leafs,
%           which means the they are at the
%           outermost level of the tree.
% int64(8) -> ZI_LIST_SETTINGSONLY 0x08
%           Returns only nodes which are marked
%           as settings
%           Or combinations of flags might be used.
%
% ziDAQ('subscribe', [handle], [path]);
% [handle] = Handle of the PID Advisor session.
% Subscribe to one or several nodes.
%
% ziDAQ('unsubscribe', [handle], [path]);
% [handle] = Handle of the PID Advisor session.
% Unsubscribe from one or several nodes..
%
% ziDAQ('get', [handle], [path]);
% [handle] = Handle of the PID Advisor session.
% [path] = Path string of the node.
%
% ziDAQ('execute', [handle]);
% [handle] = Handle of the PID Advisor session.
% Starts the pidAdvisor if not yet running.
%
% ziDAQ('trigger', [handle]);
% Not applicable to this module.
%
% result = ziDAQ('finished', [handle]);
% [handle] = Handle of the PID Advisor session.

```

```

%           Returns 1 if the command is finished, otherwise 0.
%
% result = ziDAQ('read', [handle]);
%           [handle] = Handle of the PID Advisor session.
%           Read pidAdvisor data. If the simulation is still ongoing only a
%           subset of the data is returned.
%
% result = ziDAQ('progress', [handle]);
%           [handle] = Handle of the PID Advisor session.
%           Report the progress of the command with a number
%           between 0 and 1.
%
%           ziDAQ('finish', [handle]);
%           [handle] = Handle of the PID Advisor session.
%           Stop the PID Advisor module.
%
%           ziDAQ('clear', [handle]);
%           [handle] = Handle of the PID Advisor session.
%           End the current PID Advisor thread.
%
%           ziDAQ('save', [handle]);
%           Save the measured data to a file.
%           [handle] = Handle of the PID Advisor session.
%           [filename] = File name string (without extension)..
%
% Debugging Functions
%
%           ziDAQ('setDebugLevel', [debuglevel]);
%           [debuglevel] = Debug level (trace:0, info:1, debug:2, warning:3,
%           error:4, fatal:5, status:6).
%           Enables debug log and sets the debug level.
%
%           ziDAQ('writeDebugLog', [severity], [message]);
%           [severity] = Severity (trace:0, info:1, debug:2, warning:3,
%           error:4, fatal:5, status:6).
%           [message] = Message to output to the log.
%           Outputs message to the debug log (if enabled).
%
%           ziDAQ('logOn', [flags], [filename], [style]);
%           Flags = LOG_NONE:          0x00000000
%                   LOG_SET_DOUBLE:    0x00000001
%                   LOG_SET_INT:        0x00000002
%                   LOG_SET_BYTE:       0x00000004
%                   LOG_SYNC_SET_DOUBLE: 0x00000010
%                   LOG_SYNC_SET_INT:   0x00000020
%                   LOG_SYNC_SET_BYTE:  0x00000040
%                   LOG_GET_DOUBLE:     0x00000100
%                   LOG_GET_INT:         0x00000200
%                   LOG_GET_BYTE:        0x00000400
%                   LOG_GET_DEMOD:       0x00001000
%                   LOG_GET_DIO:         0x00002000
%                   LOG_GET_AUXIN:       0x00004000
%                   LOG_LISTNODES:       0x00010000
%                   LOG_SUBSCRIBE:       0x00020000
%                   LOG_UNSUBSCRIBE:     0x00040000
%                   LOG_GET_AS_EVENT:    0x00080000
%                   LOG_UPDATE:          0x00100000
%                   LOG_POLL_EVENT:      0x00200000
%                   LOG_POLL:            0x00400000
%                   LOG_ALL :            0xffffffff
%           [filename] = Log file name
%           [style] = LOG_STYLE_TELNET: 0 (default)
%                   LOG_STYLE_MATLAB: 1
%                   LOG_STYLE_PYTHON: 2
%           Log all messages sent to the ziServer. This is useful
%           for debugging.
%

```

```
%      ziDAQ('logOff');  
%      Turn of message logging.  
%
```

Chapter 4. Python Programming

Python is open source software, freely available for download from [Python's official website](#) . Python is a high-level programming language with an extensive standard library renowned for its "batteries included" approach. Combined with the [NumPy](#) package for scientific computing, Python is a powerful computational tool for scientists that does not require expensive software licenses. The Zurich Instruments LabOne Python API, also known as `ziPython` enables the user to configure and stream data from their instrument directly into Python.

This chapter aims to help you get started using the Zurich Instruments LabOne Python API, `ziPython`, to control your instrument, please refer to:

- [Section 4.1](#) for help Installing the LabOne Python API .
- [Section 4.2](#) for help Getting Started with the LabOne Python API and Running the Examples .
- [Section 4.3](#) for LabOne Python API Tips and Tricks .
- [Section 4.4](#) for the LabOne Python API (`ziPython`) Command Reference .

Note

This chapter and the provided examples are not intended to be a Python tutorial. For help getting started with Python itself, see either the [Python Tutorial](#) or one of the many online resources, for example, the [learnpython.org](#) . The [Interactive Python Course](#) is an interesting resource for those already familiar with Python basics.

4.1. Installing the LabOne Python API

4.1.1. Requirements

In order to install and use the LabOne Python API you require:

1. Either a Python 2.6 or a Python 2.7 installation on either Windows or Linux.
2. The [NumPy](#) python package installed for your Python installation.
3. The correct version of [ziPython](#) for your Python version and platform, available from the Zurich Instruments [download page](#) (login required).

Note

Linux users must also ensure they download the version of [ziPython](#) that is Unicode compatible with their Linux distribution's Python installation, see [Section 4.1.4](#) for help determining which version is required.

Note

Important: If you your system already has an existing [ziPython](#) installation older than version 14.08, please be sure to either manually uninstall [ziPython](#) or manually remove the existing [zhinst](#) installation folder. This is due to improvements in the [zhinst](#) package structure in 14.08 (examples for different device classes are now organized in separate module/sub-directories) and the Python installer simply overwrites the existing installation, leading to a duplication of some files. For help locating `[PYTHONROOT]\lib\site-packages\zhinst\` on your system, please see [the section called "Locating the zhinst Installation Folder and Examples"](#).

4.1.2. Recommended Python Packages

The following Python packages can additionally be useful for programming with the LabOne Python API:

1. [Matplotlib](#) - recommended to plot the output from many of [ziPython](#)'s examples.
2. [SciPy](#) - recommended to load data saved from the LabOne UI in binary Matlab format (.mat).

Note

Unofficial pre-compiled 32-bit and 64-bit Windows binaries of NumPy, SciPy and matplotlib are available from Christoph Gohlke's [pythonlibs page](#).

4.1.3. Windows Installation

To install ziPython on Windows execute the .msi installer available from the Zurich Instruments [download page](#) (login required). It will guide you through the installation process as displayed in the following screenshots.




Figure 4.1. Windows ziPython installation: Step 1.

If multiple Python Installations are available on your system, the installer will ask which Python version the ziPython package should be installed in. The ziPython package will be installed in selected versions in the folder `[PYTHONROOT]\lib\site-packages\zhinst\`.




Figure 4.2. Windows ziPython installation: Step 2.




Figure 4.3. Windows ziPython installation: Step 3.

4.1.4. Linux Installation

In addition to the [requirements above](#) and selecting the correct version of ziPython for your Python distribution (2.6 or 2.7) and platform (32-bit or 64-bit), on Linux the correct Unicode version must also be installed. This is because some Python distributions on Linux are compiled to use UCS-2 character encoding, whereas some use UCS-4.

Determining the correct Unicode version of ziPython for your Python distribution

In order to determine which version of Unicode your Python distribution uses, please type the following commands in the interactive shell of your target Python distribution:

```
>>> import sys
>>> print sys.maxunicode
```

If the last command prints:

- 65535, use the UCS-2 version of ziPython,
- 1114111, use the UCS-4 version of ziPython.

Note

The installation needs root access rights. If you do not have these permissions, ask your system administrator for help.

To install ziPython on a Debian-derived distribution such as Ubuntu perform the following steps:

1. If required, install Python, NumPy and matplotlib (with elevated access rights):

```
$ sudo apt-get install python python-numpy python-matplotlib
```

2. Unpack the ziPython software bundle:

```
$ tar xzf ziPython-[version]-[build]-[linux32|linux64].tar.gz
```

3. Change directory into the unpacked folder and run the setup script `setup.py` as following:

```
$ cd ziPython-[version]-[build]-[linux32|linux64]
$ python setup.py build
$ sudo python setup.py install --install-layout=deb # Elevated access rights
```

It's possible to skip the `build` step (`install` will automatically perform this step), but splitting the steps avoids creating a directory in your user space which is owned by `root`.

4.2. Getting Started with the LabOne Python API

This section introduces the user to the LabOne Python API.

4.2.1. Contents of the LabOne Python API

Alongside the driver for interfacing with your Zurich Instruments device, the LabOne Python API includes utility functions and examples. See:

- [Section 4.4.1](#) to see which examples are available in `ziPython`.
- [Section 4.4.2](#) to see which utility functions are available in `ziPython`.

4.2.2. Using the Built-in Documentation

`ziPython`'s built-in documentation can be accessed using the `help` command in a python interactive shell:

- On module level:

```
>>> import zhinst.ziPython as ziPython
>>> help(ziPython)
```

- On class level, for example, for the Sweeper Module:

```
>>> import zhinst.ziPython as ziPython
>>> help(ziPython.ziDAQSweeper)
```

- On function level, for example, for the `ziDAQServer.poll` method:

```
>>> import zhinst.ziPython as ziPython
>>> help(ziPython.ziDAQServer.poll)
```

See [Section 4.4, LabOne Python API \(ziPython\) Command Reference](#) for a printer friendly version of the built-in documentation.

4.2.3. Running the Examples

Prerequisites for running the Python examples:

1. The `zhinst` package is installed as described above in [Section 4.1](#).
2. The Data Server program is running and the instrument is connected to the Data Server. See [Section 1.1.1](#) and the User Manual specific to the instrument for more help connecting.
3. Signal Output 1 of the instrument is connected to Signal Input 1 via a BNC cable; many of the Python examples measure on this hardware channel.

It's also recommended to install the [Matplotlib](#) Python package in order to plot the data obtained in many of the examples, see [Section 4.1.2](#).

The API examples are available in the module `zhinst.examples`, which is organized into sub-modules according to the target Instrument class:

- `zhinst.examples.common`: examples compatible with with any class of instrument,

- `zhinst.examples.uhf`: examples only compatible with the UHF Lock-in Amplifier,
- `zhinst.examples.hf2`: examples only compatible with HF2 Series Instruments.

All the examples follow the same structure and take one input argument: The device ID of the instrument to run the example with. The recommended way to run a `ziPython` example is to import the example's module in an interactive shell and call the `run_example()` function. For example, to run the `zoomFFT` Module example:

```
>>> import zhinst.examples
>>> # Use do_plot=False if matplotlib is unavailable
>>> zhinst.examples.common.example_zoomfft.run_example('dev123', do_plot=True);
```

The example should produce some output in the Python shell, such as:

```
Will perform 1 zoomFFTs
Individual zoomFFT 100.00 complete.
sample contains 1 zoomFFTs
Number of lines in first zoomFFT: 65535
```

Most examples will also plot the retrieved data using `matplotlib`, see [Figure 4.4](#) for an example. If you encounter an error message please ensure that the [above prerequisites](#) are fulfilled.




Figure 4.4. The plot produced by the LabOne Python API example `example_zoomfft.py`; the plots show the results of an FFT performed with `ziCore`'s `zoomFFT` module on demodulator output obtained over a simple feedback cable.

Exploring which Examples are available

Python's help system can be used to see which examples are available for a particular device class; when help is called on the module the available examples are listed under the "Package Contents" section. For example, for the `zhinst.examples.common` package:

```
>>> help('zhinst.examples.common')

Help on package zhinst.examples.common in zhinst.examples:

NAME
    zhinst.examples.common - Zurich Instruments LabOne Python API Examples (for any
    instrument class).

FILE
    /usr/lib/python2.7/dist-packages/zhinst/examples/common/__init__.py

MODULE DOCS
    http://docs.python.org/library/zhinst.examples.common

PACKAGE CONTENTS
    example_connect
    example_connect_config
    example_poll
    example_record_edge_trigger
    example_save_device_settings
    example_save_device_settings_expert
    example_save_device_settings_simple
    example_sweeper
    example_zoomfft

DATA
    __all__ = ['example_connect', 'example_connect', 'example_connect_conf...
```

Locating the **zhinst** Installation Folder and Examples

The examples distributed with the **zhinst** package can serve as a starting point to program your own measurement needs. The example python files, however, are generally not installed in user space. In order to ensure that you have sufficient permission to edit the examples and that your modifications are not overwritten by a later upgrade of the **zhinst** package, please copy them to your own user space before editing them.

The examples are contained in a subfolder of the **zhinst** package installation folder

```
[PYTHONROOT]\lib\site-packages\zhinst\
```

If you are unsure about the location of your **PYTHONROOT**, the `__path__` attribute of the **zhinst** module can be used in order to determine its location, for example,

```
>>> import zhinst
>>> print zhinst.__path__
```

will output something similar to:

```
C:\Python27\lib\site-packages\zhinst
```

4.2.4. Using **ziCore** Modules in the LabOne Python API

In the LabOne Python API **ziCore** [Modules](#) are configured and controlled by instantiating an object of the Module's class. For example, in order to use the [Sweeper Module](#) a sweeper object is created as following:

```
>>> daq = ziPython.ziDAQServer('localhost', 8004, 4) # Create a connection to the
      Data Server
      (running on the same PC as the API client)
>>> timeout_milliseconds = 500
>>> sweeper = daq.sweep(timeout_milliseconds);
```

Note, that since creating a Module object without an API connection to the Data Server does not make sense, the Sweeper object is instantiated via the `sweep` method of the `ziDAQServer` class, not directly from the `ziDAQSweeper` class.

The Module's parameters are configured using the Module's `set` method and specifying a path, value pair, for example:

```
>> sweeper.set('sweep/start', 1.2e5);
```

The parameters can be read-back using the `get` method, which supports wildcards, for example:

```
>> sweep_params = sweeper.get('sweep/*');
```

The variable `sweep_params` now contains a dictionary of all the Sweeper's parameters. The other main Module commands are similarly used, e.g., `sweeper.execute()`, to start the sweeper. See [Section 2.1.2](#) for more help with Modules and a description of their parameters.

4.3. LabOne Python API Tips and Tricks

In this section some tips and tricks for working with the LabOne Python API are provided.

Data Structures returned by **ziPython**.

The output arguments that **ziPython** returns are designed to use the native data structures that Python users are familiar with and that reflect the data's location in the instruments node hierarchy. For example, when the `poll` command returns data from the instruments fourth demodulator (located in the node hierarchy as `/dev123/demods/3/sample`), the output argument contains a tree of nested dictionaries in which the data can be accessed by

```
data = daq.poll( poll_length, poll_timeout);
x = data['dev123']['demods']['4']['sample']['x'];
y = data['dev123']['demods']['4']['sample']['y'];
```

Tell **poll** to return a flat dictionary

By default, the data returned by `poll` is contained in a tree of nested dictionaries that closely mimics the tree structure of the instrument node hierarchy. By setting the optional fifth argument of `poll` to `True`, the data will be a flat dictionary. This can help avoid many nested `if` statements in order to check that the expected data was returned by `poll`. For example:

```
daq.subscribe('/dev123/demods/0/sample')
flat_dictionary_key = False
data = daq.poll(0.1, 200, 1, flat_dictionary_key)
if 'dev123' in data:
    if 'demods' in data['device']:
        if '0' in data['device']['demods']:
            # access the demodulator data:
            x = data['dev123']['demods']['0']['sample']['x']
            y = data['dev123']['demods']['0']['sample']['y']
```

Could be rewritten more concisely as:

```
daq.subscribe('/dev123/demods/0/sample')
flat_dictionary_key = True
data = daq.poll(0.1, 200, 1, flat_dictionary_key)
if '/dev123/demods/0/sample' in data:
    # access the demodulator data:
    x = data['/dev123/demods/0/sample']['x']
    y = data['/dev123/demods/0/sample']['y']
```

Use the Utility Routines to load Data saved from the LabOne UI and ziControl in Python.

The utilities package `zhinst.utils` contains several routines to help loading `.csv` or `.mat` files saved from either the LabOne User Interface or `ziControl` into Python. These functions are generally minimal wrappers around [NumPy](#) (`genfromtxt()`) or [SciPy](#) (`loadmat()`) routines. However, the function `load_labone_demod_csv()` is optimized to load demodulator data saved in `.csv` format by the LabOne UI (since it specifies the `.csv` columns' `dtypes` explicitly) and the function `load_zicontrol_zibin()` can directly load data saved in binary format from `ziControl`. See [Section 4.4.2](#) for reference documentation on these commands.

4.4. LabOne Python API (ziPython) Command Reference

The following reference documentation for ziPython is available in from within a python session using python's help (see [Section 4.2.2](#)) command; It is included here for convenience.

The documentation is grouped by module and class as following:

- [Help for the zhinst Python Package](#)
- [Help for zhinst's Utility Functions](#)
- [Help for ziPython's ziDAQServer class](#)
- [Help for ziPython's ziDeviceSettings class](#)
- [Help for ziPython's ziDAQSweeper class](#)
- [Help for ziPython's ziDAQZoomFFT class](#)
- [Help for ziPython's ziDAQRecorder class](#)
- [Help for ziPython's ziPllAdvisor class](#)
- [Help for ziPython's ziPidAdvisor class](#)

4.4.1. Help for the **zhinst** Python Package

```
>>> help('zhinst')

Help on package zhinst:

NAME
    zhinst - Zurich Instruments LabOne Python API

FILE
    /usr/lib/python2.7/dist-packages/zhinst/__init__.py

MODULE DOCS
    http://docs.python.org/library/zhinst

DESCRIPTION
    Contains the API driver, utility functions and examples for Zurich Instruments
    devices.

PACKAGE CONTENTS
    examples (package)
    utils
    ziPython

DATA
    __all__ = ['ziPython', 'utils']
```

4.4.2. Help for zhinst's Utility Functions

```
>>> help('zhinst.utils')

Help on module zhinst.utils in zhinst:

NAME
    zhinst.utils - Zurich Instruments LabOne Python API Utility Functions.
```


FILE

/usr/lib/python2.7/dist-packages/zhinst/Utils.py

MODULE DOCS

<http://docs.python.org/library/zhinst.Utils>

DESCRIPTION

This module provides basic utility functions for:

- Creating an API session by connecting to an appropriate Data Server.
- Detecting devices.
- Loading and saving device settings.
- Loading data saved by either the Zurich Instruments LabOne User Interface or ziControl into Python as numpy structured arrays.

FUNCTIONS

`autoConnect(default_port=None, api_level=None)`

Try to connect to a Zurich Instruments Data Server with an attached available UHF or HF2 device.

Important: `autoConnect()` does not support MFLI devices.

Args:

`default_port` (int, optional): The default port to use when connecting to the Data Server (specify 8005 for the HF2 Data Server and 8004 for the UHF Data Server).

`api_level` (int, optional): The API level to use, either 1, 4 or 5. HF2 only supports Level 1, Level 5 is recommended for UHF and MFLI devices.

Returns:

`ziDAQServer`: An instance of the `ziPython.ziDAQServer` class that is used for communication to the Data Server.

Raises:

`RuntimeError`: If no running Data Server is found or no device is found that is attached to a Data Server.x

If `default_port` is not specified (`=None`) then first try to connect to a HF2, if no server devices are found then try to connect to an UHF. This behaviour is useful for the API examples. If we cannot connect to a server and/or detect a connected device raise a `RuntimeError`.

If `default_port` is 8004 try to connect to a UHF; if it is 8005 try to connect to an HF2. If no server and device is detected on this port raise a `RuntimeError`.

`autoDetect(daq, exclude=None)`

Return a string containing the first device ID (not in the exclude list) that is attached to the Data Server connected via `daq`, an instance of the `ziPython.ziDAQServer` class.

Args:

`daq` (`ziDAQServer`): An instance of the `ziPython.ziDAQServer` class (representing an API session connected to a Data Server).

`exclude` (list of str, optional): A list of strings specifying devices to exclude. `autoDetect()` will not return the name of a device in this list.

Returns:

A string specifying the first device ID not in exclude.

Raises:

RuntimeError: If no device was found.

RuntimeError: If daq is not an instance of ziPython.ziDAQServer.

Example:

```
zhinst.utils
daq = zhinst.utils.autoConnect()
device = zhinst.utils.autoDetect(daq)
```

check_for_sampleloss(timestamps)

Check whether timestamps are equidistantly spaced, if not, it is an indication that sampleloss has occurred whilst recording the demodulator data.

This function assumes that the timestamps originate from continuously saved demodulator data, during which the demodulator sampling rate was not changed.

Arguments:

timestamp (numpy array): a 1-dimensional array containing demodulator timestamps

Returns:

idx (numpy array): a 1-dimensional array indicating the indices in timestamp where sampleloss has occurred. An empty array is returned if no sampleloss was present.

devices(daq)

Return a list of strings containing the device IDs that are attached to the Data Server connected via daq, an instance of the ziPython.ziDAQServer class. Returns an empty list if no devices are found.

Args:

daq (ziDAQServer): An instance of the ziPython.ziDAQServer class (representing an API session connected to a Data Server).

Returns:

A list of strings of connected device IDs. The list is empty if no devices are detected.

Raises:

RuntimeError: If daq is not an instance of ziPython.ziDAQServer.

Example:

```
import zhinst.utils
daq = zhinst.utils.autoConnect() # autoConnect not supported for MFLI
devices
device = zhinst.utils.autoDetect(daq)
```

get_default_settings_path(daq)

Return the default path used for settings by the ziDeviceSettings module.

Arguments:

`daq` (instance of `ziDAQServer`): A ziPython API session.

Returns:

`settings_path` (str): The default `ziDeviceSettings` path.

`load_labone_csv(fname)`

Load a CSV file containing generic data as saved by the LabOne User Interface into a numpy structured array.

Arguments:

`filename` (str): The filename of the CSV file to load.

Returns:

`sample` (numpy ndarray): A numpy structured array of shape `(num_points,)` whose field names correspond to the column names in the first line of the CSV file. `num_points` is the number of lines in the CSV file - 1.

Example:

```
import zhinst.utils
# Load the CSV file of PID error data (node: /dev2004/pids/0/error)
data = zhinst.utils.load_labone_csv('dev2004_pids_0_error_00000.csv')
import matplotlib.pyplot as plt
# Plot the error
plt.plot(data['timestamp'], data['value'])
```

`load_labone_demod_csv(fname, column_names=('chunk', 'timestamp', 'x', 'y', 'freq', 'phase', 'dio', 'trigger', 'auxin0', 'auxin1'))`

Load a CSV file containing demodulator samples as saved by the LabOne User Interface into a numpy structured array.

Arguments:

`fname` (file or str): The file or filename of the CSV file to load.

`column_names` (list or tuple of str, optional): A list (or tuple) of column names to load from the CSV file. Default is to load all columns.

Returns:

`sample` (numpy ndarray): A numpy structured array of shape `(num_points,)` whose field names correspond to the column names in the first line of the CSV file. `num_points` is the number of lines in the CSV file - 1.

Example:

```
import zhinst.utils
sample =
zhinst.utils.load_labone_demod_csv('dev2004_demods_0_sample_00000.csv',
('timestamp', 'x', 'y'))
import matplotlib.pyplot as plt
import numpy as np
plt.plot(sample['timestamp'], np.abs(sample['x'] + 1j*sample['y']))
```

`load_labone_mat(filename)`

A wrapper function for loading a MAT file as saved by the LabOne User Interface with `scipy.io`'s `loadmat()` function. This function is included mainly to document how to work with the data structure return by `scipy.io.loadmat()`.

Arguments:

`filename` (str): the name of the MAT file to load.

Returns:

data (dict): a nested dictionary containing the instrument data as specified in the LabOne User Interface. The nested structure of ``data`` corresponds to the path of the data's node in the instrument's node hierarchy.

Further comments:

The MAT file saved by the LabOne User Interface (UI) is a Matlab V5.0 data file. The LabOne UI saves the specified data using native Matlab data structures in the same format as are returned by commands in the LabOne Matlab API. More specifically, these data structures are nested Matlab structs, the nested structure of which correspond to the location of the data in the instrument's node hierarchy.

Matlab structs are returned by `scipy.io.loadmat()` as dictionaries, the name of the struct becomes a key in the dictionary. However, as for all objects in MATLAB, structs are in fact arrays of structs, where a single struct is an array of shape (1, 1). This means that each (nested) dictionary that is returned (corresponding to a node in node hierarchy) is loaded by `scipy.io.loadmat` as a 1-by-1 array and must be indexed as such. See the ``Example`` section below.

For more information please refer to the following link:
<http://docs.scipy.org/doc/scipy/reference/tutorial/io.html#matlab-structs>

Example:

```
device = 'dev88'
# See ``Further explanation`` above for a comment on the indexing:
timestamp = data[device][0,0]['demods'][0,0]['sample'][0,0]['timestamp'][0]
x = data[device][0,0]['demods'][0,0]['sample'][0,0]['x'][0]
y = data[device][0,0]['demods'][0,0]['sample'][0,0]['y'][0]
import matplotlib.pyplot as plt
import numpy as np
plt.plot(timestamp, np.abs(x + 1j*y))

# If multiple demodulator's are saved, data from the second demodulator,
# e.g., is accessed as following:
x = data[device][0,0]['demods'][0,1]['sample'][0,0]['x'][0]
```

`load_settings(daq, device, filename)`
Load a LabOne settings file to the specified device. This function is synchronous; it will block until loading the settings has finished.

Arguments:

daq (instance of ziDAQServer): A ziPython API session.

device (str): The device ID specifying where to load the settings, e.g., 'dev123'.

filename (str): The filename of the xml settings file to load. The filename can include a relative or full path.

Raises:

RuntimeError: If loading the settings times out.

Examples:

```
import zhinst.utils as utils
daq = utils.autoConnect()
dev = utils.autoDetect(daq)

# Then, e.g., load settings from a file in the current directory:
```

```
utils.load_settings(daq, dev, 'my_settings.xml')
# Then, e.g., load settings from the default LabOne settings path:
filename = 'default_ui.xml'
path = utils.get_default_settings_path(daq)
utils.load_settings(daq, dev, path + os.sep + filename)

load_zicontrol_csv(filename, column_names=('t', 'x', 'y', 'freq', 'dio',
'auxin0', 'auxin1'))
Load a CSV file containing demodulator samples as saved by the ziControl
User Interface into a numpy structured array.
```

Arguments:

`filename` (str): The file or filename of the CSV file to load.

`column_names` (list or tuple of str, optional): A list (or tuple) of column names (demodulator sample field names) to load from the CSV file. Default is to load all columns.

Returns:

`sample` (numpy ndarray): A numpy structured array of shape `(num_points,)` whose field names correspond to the field names of a ziControl demodulator sample. `num_points` is the number of lines in the CSV file - 1.

Example:

```
import zhinst.utils
sample = zhinst.utils.load_labone_csv('Freq1.csv', ('t', 'x', 'y'))
import matplotlib.pyplot as plt
import numpy as np
plt.plot(sample['t'], np.abs(sample['x'] + 1j*sample['y']))

load_zicontrol_zibin(filename, column_names=('t', 'x', 'y', 'freq', 'dio',
'auxin0', 'auxin1'))
Load a ziBin file containing demodulator samples as saved by the ziControl
User Interface into a numpy structured array. This is for data saved by
ziControl in binary format.
```

Arguments:

`filename` (str): The filename of the .ziBin file to load.

`column_names` (list or tuple of str, optional): A list (or tuple) of column names to load from the CSV file. Default is to load all columns.

Returns:

`sample` (numpy ndarray): A numpy structured array of shape `(num_points,)` whose field names correspond to the field names of a ziControl demodulator sample. `num_points` is the number of sample points saved in the file.

Further comments:

Specifying a fewer names in `column_names` will not result in a speed-up as all data is loaded from the binary file by default.

Example:

```
import zhinst.utils
sample = zhinst.utils.load_zicontrol_zibin('Freq1.ziBin')
import matplotlib.pyplot as plt
import numpy as np
plt.plot(sample['t'], np.abs(sample['x'] + 1j*sample['y']))

save_settings(daq, device, filename)
Save settings from the specified device to a LabOne settings file. This
```

function is synchronous; it will block until saving the settings has finished.

Arguments:

`daq` (instance of `ziDAQServer`): A ziPython API session.

`device` (str): The device ID specifying where to load the settings, e.g., 'dev123'.

`filename` (str): The filename of the LabOne xml settings file. The filename can include a relative or full path.

Raises:

`RuntimeError`: If saving the settings times out.

Examples:

```
import zhinst.utils as utils
daq = utils.autoConnect()
dev = utils.autoDetect(daq)

# Then, e.g., save settings to a file in the current directory:
utils.save_settings(daq, dev, 'my_settings.xml')

# Then, e.g., save settings to the default LabOne settings path:
filename = 'my_settings_example.xml'
path = utils.get_default_settings_path(daq)
utils.save_settings(daq, dev, path + os.sep + filename)
```

DATA

```
LABONE_DEMOD_DTYPE = [('chunk', 'u8'), ('timestamp', 'u8'), ('x', 'f8'...
LABONE_DEMOD_FORMATS = ('u8', 'u8', 'f8', 'f8', 'f8', 'f8', 'u4', 'u4'...
LABONE_DEMOD_NAMES = ('chunk', 'timestamp', 'x', 'y', 'freq', 'phase',...
ZICONTROL_DTYPE = [('t', 'f8'), ('x', 'f8'), ('y', 'f8'), ('freq', 'f8'...
ZICONTROL_FORMATS = ('f8', 'f8', 'f8', 'f8', 'u4', 'f8', 'f8')
ZICONTROL_NAMES = ('t', 'x', 'y', 'freq', 'dio', 'auxin0', 'auxin1')
```

4.4.3. Help for ziPython's `ziDAQServer` class

```
>>> help('zhinst.ziPython.ziDAQServer')
```

Help on class `ziDAQServer` in `zhinst.ziPython`:

```
zhinst.ziPython.ziDAQServer = class ziDAQServer(Boost.Python.instance)
|   Class to connect with the instrument server of Zurich Instruments.
|
|   Method resolution order:
|       ziDAQServer
|       Boost.Python.instance
|       __builtin__.object
|
|   Methods defined here:
|
|   __init__(...)
|       __init__( (object)arg1) -> None
|
|       __init__( (object)arg1, (str)arg2, (int)arg3) -> None :
|           Connect to the server by using host address and port number.
|           arg1: Reference to the ziDAQServer class.
```

```

        arg2: Host string e.g. '127.0.0.1' for localhost.
        arg3: Port number e.g. 8004 for the ziDataServer.

    __init__( (object)arg1, (str)arg2, (int)arg3, (int)arg4) -> None :
        Connect to the server by using host address and port number.
        arg1: Reference to the ziDAQServer class.
        arg2: Host string e.g. '127.0.0.1' for localhost.
        arg3: Port number e.g. 8004 for the ziDataServer.
        arg4: API level number.

    __reduce__ = <unnamed Boost.Python function>(...)

connect(...)
    connect( (ziDAQServer)arg1) -> None

connectDevice(...)
    connectDevice( (ziDAQServer)arg1, (str)arg2, (str)arg3, (str)arg4) -> None :
        Connect with the data server to a specified device over the specified
        interface. The device must be visible to the server. If the device is
        already connected the call will be ignored. The function will block
        until the device is connected and the device is ready to use. This
        method is useful for UHF devices offering several communication
        interfaces.
        arg1: Reference to the ziDAQServer class.
        arg2: Device serial.
        arg3: Device interface.
        arg4: Optional interface parameters string.

    connectDevice( (ziDAQServer)arg1, (str)arg2, (str)arg3) -> None

deviceSettings(...)
    deviceSettings( (ziDAQServer)arg1, (long)arg2) -> ziDeviceSettings :
        Create a deviceSettings class. This will start a thread for running an
        asynchronous deviceSettings.
        arg1: Reference to the ziDAQServer class.
        arg2: Timeout in [ms]. Recommended value is 500ms.

disconnect(...)
    disconnect( (ziDAQServer)arg1) -> None

disconnectDevice(...)
    disconnectDevice( (ziDAQServer)arg1, (str)arg2) -> None :
        Disconnect a device on the data server. This function will return
        immediately. The disconnection of the device may not yet finished.
        arg1: Reference to the ziDAQServer class.
        arg2: Device serial string of device to disconnect.

echoDevice(...)
    echoDevice( (ziDAQServer)arg1, (str)arg2) -> None :
        Sends an echo command to a device and blocks until
        answer is received. This is useful to flush all
        buffers between API and device to enforce that
        further code is only executed after the device executed
        a previous command.
        arg1: Reference to the ziDAQServer class.
        arg2: Device string e.g. 'dev100'.

flush(...)
    flush( (ziDAQServer)arg1) -> None :
        Flush all data in the socket connection and API buffers.
        Call this function before a subscribe with subsequent poll
        to get rid of old streaming data that might still be in
        the buffers.
        arg1: Reference to the ziDAQServer class.

get(...)
    get( (ziDAQServer)arg1, (str)arg2, (bool)arg3) -> object :

```

```

    Return a dict with all nodes from the specified sub-tree.
    High-speed streaming nodes (e.g. /devN/demods/0/sample)
    are not returned. Wildcards (*) may be used, in which case
    read-only nodes are ignored.
    arg1: Reference to the ziDAQServer class.
    arg2: Path string of the node. Use wild card to
    select all.
    arg3[optional]: Specify which type of data structure to return.
    Return data either as a flat dict (True) or as a nested
    dict tree (False). Default = False.

    get( (ziDAQServer)arg1, (str)arg2) -> object

getAsEvent(...)
    getAsEvent( (ziDAQServer)arg1, (str)arg2) -> None :
    Trigger an event on the specified node. The node data is returned by a
    subsequent poll command.
    arg1: Reference to the ziDAQServer class.
    arg2: Path string of the node.

getAuxInSample(...)
    getAuxInSample( (ziDAQServer)arg1, (str)arg2) -> object :
    Returns a single auxin sample. The auxin data is averaged in contrast to
    the auxin data embedded in the demodulator sample.
    arg1: Reference to the ziDAQServer class.
    arg2: Path string

getByte(...)
    getByte( (ziDAQServer)arg1, (str)arg2) -> object :
    Get a byte array (string) value from the specified node.
    arg1: Reference to the ziDAQServer class.
    arg2: Path string of the node.

getConnectionAPILevel(...)
    getConnectionAPILevel( (ziDAQServer)arg1) -> int :
    Returns ziAPI level used for the active connection.

getDIO(...)
    getDIO( (ziDAQServer)arg1, (str)arg2) -> object :
    Returns a single DIO sample.
    arg1: Reference to the ziDAQServer class.
    arg2: Path string

getDouble(...)
    getDouble( (ziDAQServer)arg1, (str)arg2) -> float :
    Get a double value from the specified node.
    arg1: Reference to the ziDAQServer class.
    arg2: Path string of the node.

getInt(...)
    getInt( (ziDAQServer)arg1, (str)arg2) -> int :
    Get a integer value from the specified node.
    arg1: Reference to the ziDAQServer class.
    arg2: Path string of the node.

getList(...)
    getList( (ziDAQServer)arg1, (str)arg2) -> object :
    Return a list with all nodes from the specified sub-tree.
    arg1: Reference to the ziDAQServer class.
    arg2: Path string of the node. Use wild card to
    select all.

getSample(...)
    getSample( (ziDAQServer)arg1, (str)arg2) -> object :
    Returns a single demodulator sample (including DIO and AuxIn). For more
    efficient data recording use subscribe and poll methods.
    arg1: Reference to the ziDAQServer class.

```



```

        arg2: Path string

listNodes(...)
    listNodes( (ziDAQServer)arg1, (str)arg2, (int)arg3) -> list :
        This function returns a list of node names found at the specified path.
        arg1: Reference to the ziDAQRecorder class.
        arg2: Path for which the nodes should be listed. The path may
            contain wildcards so that the returned nodes do not
            necessarily have to have the same parents.
        arg3: Enum that specifies how the selected nodes are listed.
            ziPython.ziListEnum.none -> 0x00
                The default flag, returning a simple
                listing if the given node
            ziPython.ziListEnum.recursive -> 0x01
                Returns the nodes recursively
            ziPython.ziListEnum.absolute -> 0x02
                Returns absolute paths
            ziPython.ziListEnum.leafsonly -> 0x04
                Returns only nodes that are leafs,
                which means the they are at the
                outermost level of the tree.
            ziPython.ziListEnum.settingsonly -> 0x08
                Returns only nodes which are marked
                as setting
            Or combinations of flags can be used.

logOff(...)
    logOff( (ziDAQServer)arg1) -> None :
        Disables logging of commands sent to a server.
        arg1: Reference to the ziDAQServer class.

logOn(...)
    logOn( (ziDAQServer)arg1, (int)arg2, (str)arg3, (int)arg4) -> None :
        Enables logging of commands sent to a server.
        arg1: Reference to the ziDAQServer class.
        arg2: Flags (LOG_NONE:          0x00000000
            LOG_SET_DOUBLE:          0x00000001
            LOG_SET_INT:             0x00000002
            LOG_SET_BYTE:            0x00000004
            LOG_SYNC_SET_DOUBLE:     0x00000010
            LOG_SYNC_SET_INT:        0x00000020
            LOG_SYNC_SET_BYTE:       0x00000040
            LOG_GET_DOUBLE:          0x00000100
            LOG_GET_INT:             0x00000200
            LOG_GET_BYTE:            0x00000400
            LOG_GET_DEMOD:           0x00001000
            LOG_GET_DIO:             0x00002000
            LOG_GET_AUXIN:           0x00004000
            LOG_LISTNODES:           0x00010000
            LOG_SUBSCRIBE:           0x00020000
            LOG_UNSUBSCRIBE:         0x00040000
            LOG_GET_AS_EVENT:        0x00080000
            LOG_UPDATE:              0x00100000
            LOG_POLL_EVENT:          0x00200000
            LOG_POLL:                0x00400000
            LOG_ALL :                0xffffffff)
        arg3: Log file name.
        arg4: Log style (LOG_STYLE_TELNET: 0 (default),
            LOG_STYLE_MATLAB: 1, LOG_STYLE_PYTHON: 2).

    logOn( (ziDAQServer)arg1, (int)arg2, (str)arg3) -> None

pidAdvisor(...)
    pidAdvisor( (ziDAQServer)arg1, (long)arg2) -> ziPidAdvisor :
        Create a pidAdvisor class. This will start a thread for running an
        asynchronous pidAdvisor.
        arg1: Reference to the ziDAQServer class.

```

```

        arg2: Timeout in [ms]. Recommended value is 500ms.

    pllAdvisor(...)
        pllAdvisor( (ziDAQServer)arg1, (long)arg2) -> ziPllAdvisor :
        Create a pllAdvisor class. This will start a thread for running an
        asynchronous pllAdvisor.
        arg1: Reference to the ziDAQServer class.
        arg2: Timeout in [ms]. Recommended value is 500ms.

    poll(...)
    poll( (ziDAQServer)arg1, (float)arg2, (long)arg3, (int)arg4, (bool)arg5) ->
object :
        This function returns subscribed data previously in the API's buffers or
        obtained during the specified time. It returns a dict tree containing
        the recorded data. This function blocks until the recording time is
        elapsed.
        arg1: Reference to the ziDAQServer class.
        arg2: Recording time in [s]. The function will block during that.
        time.
        arg3: Poll timeout in [ms]. Recommended value is 500ms.
        arg4[optional]: Poll flags.
            FILL = 0x0001 : Fill holes.
            ALIGN = 0x0002 : Align data that contains a
                        timestamp.
            THROW = 0x0004 : Throw EOFError exception if sample
                        loss is detected.
        arg5[optional]: Specify which type of data structure to return.
        Return data either as a flat dict (True) or as a nested
        dict tree (False). Default = False.

    poll( (ziDAQServer)arg1, (float)arg2, (long)arg3 [, (int)arg4]) -> object

    pollEvent(...)
    pollEvent( (ziDAQServer)arg1, (long)arg2) -> object :
        Execute a single poll command. Note: only one data packet will be
        fetched. To get all data waiting in the buffers this command should be
        executed continuously until nothing is returned anymore. This is a low
        level command. Use the poll command or asynchronous recording instead.
        arg1: Reference to the ziDAQServer class.
        arg2: Poll timeout in [ms]. Recommended value is 500ms.

    programRT(...)
    programRT( (ziDAQServer)arg1, (str)arg2, (str)arg3) -> None :
        Program RT.
        arg1: Device identifier e.g. 'dev99'.
        arg2: File name of the RT program.

    record(...)
    record( (ziDAQServer)arg1, (float)arg2, (long)arg3, (int)arg4) ->
ziDAQRecorder :
        Create a recording class. This will start a thread for asynchronous
        recording.
        arg1: Reference to the ziDAQServer class.
        arg2: Maximum recording time for single triggers in [s].
        arg3: Timeout in [ms]. Recommended value is 500ms.
        arg4[optional]: Record flags.
            FILL = 0x0001 : Fill holes.
            ALIGN = 0x0002 : Align data that contains a
                        timestamp.
            THROW = 0x0004 : Throw EOFError exception if
                        sample loss is detected.

    record( (ziDAQServer)arg1, (float)arg2, (long)arg3) -> ziDAQRecorder

    revision(...)
    revision( (ziDAQServer)arg1) -> int :
        Get the revision number of the Python interface of Zurich Instruments.

```

```

|         arg1: Reference to the ziDAQServer class.
|
| set(...)
|     set( (ziDAQServer)arg1, (object)arg2) -> None :
|         arg1: Reference to the ziDAQServer class.
|         arg2: A list of path/value pairs.
|
| setByte(...)
|     setByte( (ziDAQServer)arg1, (str)arg2, (object)arg3) -> None :
|         arg1: Reference to the ziDAQServer class.
|         arg2: Path string of the node.
|
| setDebugLevel(...)
|     setDebugLevel( (ziDAQServer)arg1, (int)arg2) -> None :
|         Enables debug log and sets the debug level.
|         arg1: Reference to the ziDAQServer class.
|         arg2: Debug level (trace:0, info:1, debug:2, warning:3, error:4,
|             fatal:5, status:6).
|
| setDouble(...)
|     setDouble( (ziDAQServer)arg1, (str)arg2, (float)arg3) -> None :
|         arg1: Reference to the ziDAQServer class.
|         arg2: Path string of the node.
|
| setInt(...)
|     setInt( (ziDAQServer)arg1, (str)arg2, (int)arg3) -> None :
|         arg1: Reference to the ziDAQServer class.
|         arg2: Path string of the node.
|
| subscribe(...)
|     subscribe( (ziDAQServer)arg1, (object)arg2) -> None :
|         Subscribe to one or several nodes. Fetch data with the poll
|         command. In order to avoid fetching old data that is still in the
|         buffer execute a flush command before subscribing to data streams.
|         arg1: Reference to the ziDAQServer class.
|         arg2: Path string of the node. Use wild card to
|             select all. Alternatively also a list of path
|             strings can be specified.
|
| sweep(...)
|     sweep( (ziDAQServer)arg1, (long)arg2) -> ziDAQSweeper :
|         Create a sweeper class. This will start a thread for asynchronous
|         sweeping.
|         arg1: Reference to the ziDAQServer class.
|         arg2: Timeout in [ms]. Recommended value is 500ms.
|
| sync(...)
|     sync( (ziDAQServer)arg1) -> None :
|         Synchronize all data path. Ensures that get and poll
|         commands return data which was recorded after the
|         setting changes in front of the sync command. This
|         sync command replaces the functionality of all syncSet,
|         flush, and echoDevice commands.
|         arg1: Reference to the ziDAQServer class.
|
| syncSetDouble(...)
|     syncSetDouble( (ziDAQServer)arg1, (str)arg2, (float)arg3) -> float :
|         arg1: Reference to the ziDAQServer class.
|         arg2: Path string of the node.
|
| syncSetInt(...)
|     syncSetInt( (ziDAQServer)arg1, (str)arg2, (int)arg3) -> int :
|         arg1: Reference to the ziDAQServer class.
|         arg2: Path string of the node.
|
| unsubscribe(...)
|     unsubscribe( (ziDAQServer)arg1, (object)arg2) -> None :

```

```
|
|         Unsubscribe data streams. Use this command after recording to avoid
|         buffer overflows that may increase the latency of other command.
|         arg1: Reference to the ziDAQServer class.
|         arg2: Path string of the node. Use wild card to
|               select all. Alternatively also a list of path
|               strings can be specified.
|
| update(...)
|     update( (ziDAQServer)arg1) -> None :
|         Check if additional devices are attached. This function is not needed
|         for servers running under windows as devices will be detected
|         automatically.
|         arg1: Reference to the ziDAQServer class.
|
| version(...)
|     version( (ziDAQServer)arg1) -> str :
|         Get version string of the Python interface of Zurich Instruments.
|         arg1: Reference to the ziDAQServer class.
|
| writeDebugLog(...)
|     writeDebugLog( (ziDAQServer)arg1, (int)arg2, (str)arg3) -> None :
|         Outputs message to the debug log (if enabled).
|         arg1: Reference to the ziDAQServer class.
|         arg2: Severity (trace:0, info:1, debug:2, warning:3, error:4,
|               fatal:5, status:6).
|         arg3: Message to output to the log.
|
| zoomFFT(...)
|     zoomFFT( (ziDAQServer)arg1, (long)arg2) -> ziDAQZoomFFT :
|         Create a zoomFFT class. This will start a thread for running an
|         asynchronous zoomFFT.
|         arg1: Reference to the ziDAQServer class.
|         arg2: Timeout in [ms]. Recommended value is 500ms.
|
| -----
| Data and other attributes defined here:
|
| __instance_size__ = 48
|
| -----
| Data descriptors inherited from Boost.Python.instance:
|
| __dict__
|
| __weakref__
|
| -----
| Data and other attributes inherited from Boost.Python.instance:
|
| __new__ = <built-in method __new__ of Boost.Python.class object>
|         T.__new__(S, ...) -> a new object with type S, a subtype of T
```

4.4.4. Help for ziPython's ziDeviceSettings class

An instance of ziDeviceSettings is initialized using the deviceSettings method from ziDAQServer:

```
>>> help('zhinst.ziPython.ziDAQServer.deviceSettings')
```

Help on method deviceSettings in zhinst.ziPython.ziDAQServer:

```
zhinst.ziPython.ziDAQServer.deviceSettings = deviceSettings(...) unbound
zhinst.ziPython.ziDAQServer method
deviceSettings( (ziDAQServer)arg1, (long)arg2) -> ziDeviceSettings :
    Create a deviceSettings class. This will start a thread for running an
    asynchronous deviceSettings.
    arg1: Reference to the ziDAQServer class.
    arg2: Timeout in [ms]. Recommended value is 500ms.
```

Reference help for the ziDeviceSettings class.

```
>>> help('zhinst.ziPython.ziDeviceSettings')
```

Help on class ziDeviceSettings in zhinst.ziPython:

```
zhinst.ziPython.ziDeviceSettings = class ziDeviceSettings(Boost.Python.instance)
| Method resolution order:
|     ziDeviceSettings
|     Boost.Python.instance
|     __builtin__.object
|
| Methods defined here:
|
|     __reduce__ = <unnamed Boost.Python function>(...)
|
|     clear(...)
|         clear( (ziDeviceSettings)arg1) -> None :
|             End the deviceSettings thread.
|
|     execute(...)
|         execute( (ziDeviceSettings)arg1) -> None :
|             Execute the save/loadcommand.
|
|     finish(...)
|         finish( (ziDeviceSettings)arg1) -> None :
|             Stop the load/save command. The command may be restarted by calling
|             'execute' again.
|
|     finished(...)
|         finished( (ziDeviceSettings)arg1) -> bool :
|             Check if the command execution has finished. Returns True if finished.
|
|     get(...)
|         get( (ziDeviceSettings)arg1, (str)arg2, (bool)arg3) -> object :
|             Return a dict with all nodes from the specified sub-tree.
|             arg1: Reference to the ziDAQDeviceSettings class.
|             arg2: Path string of the node. Use wild card to
|                 select all.
|             arg3[optional]: Specify which type of data structure to return.
|                 Return data either as a flat dict (True) or as a nested
|                 dict tree (False). Default = False.
|
|         get( (ziDeviceSettings)arg1, (str)arg2) -> object
|
|     listNodes(...)
|         listNodes( (ziDeviceSettings)arg1, (str)arg2, (int)arg3) -> list :
|             This function returns a list of node names found at the specified path.
|             arg1: Reference to the pyDeviceSettings class.
|             arg2: Path for which the nodes should be listed. The path may
|                 contain wildcards so that the returned nodes do not
|                 necessarily have to have the same parents.
|             arg3: Enum that specifies how the selected nodes are listed.
|                 ziPython.ziListEnum.none -> 0x00
|                     The default flag, returning a simple
|                     listing if the given node
|                 ziPython.ziListEnum.recursive -> 0x01
|                     Returns the nodes recursively
```

```

        ziPython.ziListEnum.absolute -> 0x02
            Returns absolute paths
        ziPython.ziListEnum.leafsonly -> 0x04
            Returns only nodes that are leafs,
            which means the they are at the
            outermost level of the tree.
        ziPython.ziListEnum.settingsonly -> 0x08
            Returns only nodes which are marked
            as setting
        Or combinations of flags can be used.

progress(...)
    progress( (ziDeviceSettings)arg1) -> object :
        Reports the progress of the command with a number between
        0 and 1.

read(...)
    read( (ziDeviceSettings)arg1, (bool)arg2) -> object :
        Read device settings. Only relevant for the save command.
        arg1[optional]: Specify which type of data structure to return.
        Return data either as a flat dict (True) or as a nested
        dict tree (False). Default = False.

    read( (ziDeviceSettings)arg1) -> object

save(...)
    save( (ziDeviceSettings)arg1, (str)arg2) -> None :
        Not relevant for the deviceSettings module.

set(...)
    set( (ziDeviceSettings)arg1, (str)arg2, (float)arg3) -> None :
        Device Settings Parameters
        Path name      Type      Description
        devicesettings/device    string  Device that should be used for
                                         loading/saving device settings,
                                         e.g. 'dev99'.
        devicesettings/path      string  Directory where settings files should be
                                         located. If not set, the default settings
                                         location of the LabOne software is used.
        devicesettings/filename  string  Name of settings file to use
        devicesettings/command   string  The command to execute
                                         'save' = Read device settings and save to
                                         file.
                                         'load' = Load settings from file and
                                         write to device.
                                         'read' = Read device settings only
                                         (no save).

    set( (ziDeviceSettings)arg1, (str)arg2, (int)arg3) -> None

    set( (ziDeviceSettings)arg1, (str)arg2, (str)arg3) -> None

    set( (ziDeviceSettings)arg1, (object)arg2) -> None :
        arg1: Reference to the pyDeviceSettings class.
        arg2: A list of path/value pairs.

subscribe(...)
    subscribe( (ziDeviceSettings)arg1, (str)arg2) -> None :
        Not relevant for the deviceSettings module.

trigger(...)
    trigger( (ziDeviceSettings)arg1) -> None :
        Not applicable to this module.

unsubscribe(...)
    unsubscribe( (ziDeviceSettings)arg1, (str)arg2) -> None :
        Not relevant for the deviceSettings module.

```

```
| -----  
| Data and other attributes defined here:  
|  
| __init__ = <built-in function __init__>  
|     Raises an exception  
|     This class cannot be instantiated from Python  
|  
| -----  
| Data descriptors inherited from Boost.Python.instance:  
|  
| __dict__  
|  
| __weakref__  
|  
| -----  
| Data and other attributes inherited from Boost.Python.instance:  
|  
| __new__ = <built-in method __new__ of Boost.Python.class object>  
|     T.__new__(S, ...) -> a new object with type S, a subtype of T
```

4.4.5. Help for ziPython's ziDAQSweeper class

An instance of ziDAQSweeper is initialized using the sweep method from ziDAQServer:

```
>>> help('zhinst.ziPython.ziDAQServer.sweep')
```

Help on method sweep in zhinst.ziPython.ziDAQServer:

```
zhinst.ziPython.ziDAQServer.sweep = sweep(...) unbound zhinst.ziPython.ziDAQServer  
method  
    sweep( (ziDAQServer)arg1, (long)arg2) -> ziDAQSweeper :  
        Create a sweeper class. This will start a thread for asynchronous  
        sweeping.  
        arg1: Reference to the ziDAQServer class.  
        arg2: Timeout in [ms]. Recommended value is 500ms.
```

Reference help for the ziDAQSweeper class.

```
>>> help('zhinst.ziPython.ziDAQSweeper')
```

Help on class ziDAQSweeper in zhinst.ziPython:

```
zhinst.ziPython.ziDAQSweeper = class ziDAQSweeper(Boost.Python.instance)  
|   Method resolution order:  
|       ziDAQSweeper  
|       Boost.Python.instance  
|       __builtin__.object  
|  
|   Methods defined here:  
|  
|   __reduce__ = <unnamed Boost.Python function>(...)  
|  
|   clear(...)   
|       clear( (ziDAQSweeper)arg1) -> None :  
|           End the sweeper thread.  
|  
|   execute(...)   
|       execute( (ziDAQSweeper)arg1) -> None :  
|           Start the sweeper. Subscription or unsubscription is no more  
|           possible until the sweep is finished.
```

```

finish(...)
    finish( (ziDAQSweeper)arg1) -> None :
        Stop sweeping. The sweeping may be restarted by calling
        'execute' again.

finished(...)
    finished( (ziDAQSweeper)arg1) -> bool :
        Check if the sweep has finished. Returns True if finished.

get(...)
    get( (ziDAQSweeper)arg1, (str)arg2, (bool)arg3) -> object :
        Return a dict with all nodes from the specified sub-tree.
        arg1: Reference to the ziDAQSweeper class.
        arg2: Path string of the node. Use wild card to
            select all.
        arg3[optional]: Specify which type of data structure to return.
            Return data either as a flat dict (True) or as a nested
            dict tree (False). Default = False.

    get( (ziDAQSweeper)arg1, (str)arg2) -> object

listNodes(...)
    listNodes( (ziDAQSweeper)arg1, (str)arg2, (int)arg3) -> list :
        This function returns a list of node names found at the specified path.
        arg1: Reference to the ziDAQRecorder class.
        arg2: Path for which the nodes should be listed. The path may
            contain wildcards so that the returned nodes do not
            necessarily have to have the same parents.
        arg3: Enum that specifies how the selected nodes are listed.
            ziPython.ziListEnum.none -> 0x00
                The default flag, returning a simple
                listing if the given node
            ziPython.ziListEnum.recursive -> 0x01
                Returns the nodes recursively
            ziPython.ziListEnum.absolute -> 0x02
                Returns absolute paths
            ziPython.ziListEnum.leafsonly -> 0x04
                Returns only nodes that are leafs,
                which means the they are at the
                outermost level of the tree.
            ziPython.ziListEnum.settingsonly -> 0x08
                Returns only nodes which are marked
                as setting
            Or combinations of flags can be used.

progress(...)
    progress( (ziDAQSweeper)arg1) -> object :
        Reports the progress of the measurement with a number between
        0 and 1.

read(...)
    read( (ziDAQSweeper)arg1, (bool)arg2) -> object :
        Read sweep data. If the sweeping is still ongoing only a subset
        of sweep data is returned. If huge data sets
        are recorded call this method to keep memory usage reasonable.
        arg1[optional]: Specify which type of data structure to return.
            Return data either as a flat dict (True) or as a nested
            dict tree (False). Default = False.

    read( (ziDAQSweeper)arg1) -> object

save(...)
    save( (ziDAQSweeper)arg1, (str)arg2) -> None :
        Save sweeper data to file.
        arg1: Reference to the ziDAQSweeper class.
        arg2: File name string (without extension).

```


set(...)			
set((ziDAQSweeper)arg1, (str)arg2, (float)arg3) -> None :			
Sweep Parameters			
Path name	Type	Description	
sweep/device	string	Device that should be used for the parameter sweep, e.g. 'dev99'.	
sweep/start	double	Sweep start frequency [Hz]	
sweep/stop	double	Sweep stop frequency [Hz]	
sweep/gridnode	string	Path of the node that should be used for sweeping. For frequency sweep applications this will be e.g. 'oscs/0/freq'. The device name of the path can be omitted and is given by sweep/device.	
sweep/loopcount	int	Number of sweep loops (default 1)	
sweep/endless	int	Sweep endless (default 0) 0 = endless off, use loopcount, 1 = endless on, ignore loopcount.	
sweep/samplecount	int	Number of samples per sweep.	
sweep/settling/time	double	Settling time before measurement is performed.	
sweep/settling/tc	double	Settling time in time constant units 5 ~ low precision 15 ~ medium precision 50 ~ high precision	
sweep/settling/inaccuracy	int	Demodulator filter settling inaccuracy defining the wait time between a sweep parameter change and recording of the next sweep point. Typical inaccuracy values: 10m for highest sweep speed for large signals, 100u for precise amplitude measurements, 100n for precise noise measurements. Depending on the order the settling accuracy will define the number of filter time constants the sweeper has to wait. The maximum between this value and the settling time is taken as wait time until the next sweep point is recorded.	
sweep/xmapping	int	Sweep mode: 0 = linear, 1 = logarithmic.	
sweep/scan	int	Scan type: 0 = sequential, 1 = binary, 2 = bidirectional, 3 = reverse.	
sweep/bandwidth	double	Fixed bandwidth [Hz], 0 = Automatic calculation.	
sweep/bandwidthcontrol	int	Sets the bandwidth control mode, (default 2): 0 = Manual (user sets bandwidth and order), 1 = Fixed (uses fixed bandwidth value), 2 = Auto (calculates best bandwidth value) Equivalent to the obsolete bandwidth = 0 setting.	
sweep/order	int	Defines the filter roll off to use in Fixed bandwidth selection. Valid values are between 0 (6 dB/octave) and 8 (48 dB/octave). An order of 0 triggers a read-out of the order from the selected demodulator.	
sweep/maxbandwidth	double	Maximal bandwidth used in auto bandwidth mode in [Hz]. The default is 1.25MHz.	
sweep/omegasuppression	double	Damping in [dB] of omega and 2-omega	

		components. Default is 40dB in favor of sweep speed. Use higher value for strong offset values or 3-omega measurement methods.
sweep/averaging/tc	double	Min averaging time [tc] 0 = no averaging (see also time!) 5 ~ low precision 15 ~ medium precision 50 ~ high precision
sweep/averaging/sample	int	Min samples to average 1 = no averaging (if averaging/tc = 0)
sweep/phaseunwrap	bool	Enable unwrapping of slowly changing phase evolutions around the +/-180 degree boundary.
sweep/sincfilter	bool	Enables the sinc filter if the sweep frequency is below 50 Hz. This will improve the sweep speed at low frequencies as omega components do not need to be suppressed by the normal low pass filter.
sweep/filename	string	This parameter is deprecated. If specified, i.e. not empty, it enables automatic saving of data in single sweep mode (sweep/endless = 0).
sweep/savepath	string	The directory where files are located when saving sweeper measurements.
sweep/fileformat	string	The format of the file for saving sweeper measurements: 0 = Matlab, 1 = CSV.
sweep/historylength	bool	Maximum number of entries stored in the measurement history.
sweep/clearhistory	bool	Remove all records from the history list.

set((ziDAQSweeper)arg1, (str)arg2, (int)arg3) -> None

set((ziDAQSweeper)arg1, (str)arg2, (str)arg3) -> None

set((ziDAQSweeper)arg1, (object)arg2) -> None :

arg1: Reference to the ziDAQSweeper class.

arg2: A list of path/value pairs.

subscribe(...)

subscribe((ziDAQSweeper)arg1, (str)arg2) -> None :

Subscribe to one or several nodes. After subscription the sweep
process can be started with the 'execute' command. During the
sweep process paths can not be subscribed or unsubscribed.

arg1: Reference to the ziDAQSweeper class.

arg2: Path string of the node. Use wild card to
select all. Alternatively also a list of path
strings can be specified.

trigger(...)

trigger((ziDAQSweeper)arg1) -> None :

Execute a manual trigger.

unsubscribe(...)

unsubscribe((ziDAQSweeper)arg1, (str)arg2) -> None :

Unsubscribe from one or several nodes. During the
sweep process paths can not be subscribed or unsubscribed.

arg1: Reference to the ziDAQSweeper class.

arg2: Path string of the node. Use wild card to
select all. Alternatively also a list of path
strings can be specified.

Data and other attributes defined here:

```
|  __init__ = <built-in function __init__>
|      Raises an exception
|      This class cannot be instantiated from Python
|
| -----
| Data descriptors inherited from Boost.Python.instance:
|
|  __dict__
|
|  __weakref__
|
| -----
| Data and other attributes inherited from Boost.Python.instance:
|
|  __new__ = <built-in method __new__ of Boost.Python.class object>
|      T.__new__(S, ...) -> a new object with type S, a subtype of T
```

4.4.6. Help for ziPython's **ziDAQZoomFFT** class

An instance of **ziDAQZoomFFT** is initialized using the **zoomFFT** method from **ziDAQServer**:

```
>>> help('zhinst.ziPython.ziDAQServer.zoomFFT')
```

Help on method **zoomFFT** in **zhinst.ziPython.ziDAQServer**:

```
zhinst.ziPython.ziDAQServer.zoomFFT = zoomFFT(...) unbound
zhinst.ziPython.ziDAQServer method
  zoomFFT( (ziDAQServer)arg1, (long)arg2) -> ziDAQZoomFFT :
    Create a zoomFFT class. This will start a thread for running an
    asynchronous zoomFFT.
    arg1: Reference to the ziDAQServer class.
    arg2: Timeout in [ms]. Recommended value is 500ms.
```

Reference help for the **ziDAQZoomFFT** class.

```
>>> help('zhinst.ziPython.ziDAQZoomFFT')
```

Help on class **ziDAQZoomFFT** in **zhinst.ziPython**:

```
zhinst.ziPython.ziDAQZoomFFT = class ziDAQZoomFFT(Boost.Python.instance)
| Method resolution order:
|   ziDAQZoomFFT
|   Boost.Python.instance
|   __builtin__.object
|
| Methods defined here:
|
|   __reduce__ = <unnamed Boost.Python function>(...)
|
|   clear(...)
|       clear( (ziDAQZoomFFT)arg1) -> None :
|           End the zoom FFT thread.
|
|   execute(...)
|       execute( (ziDAQZoomFFT)arg1) -> None :
|           Start the zoom FFT. Subscription or unsubscription is no more
|           possible until the zoom FFT is finished.
|
|   finish(...)
|       finish( (ziDAQZoomFFT)arg1) -> None :
```

```

        Stop the zoom FFT. The zoom FFT may be restarted by calling
        'execute' again.

finished(...)
    finished( (ziDAQZoomFFT)arg1) -> bool :
        Check if the zoom FFT has finished. Returns True if finished.

get(...)
    get( (ziDAQZoomFFT)arg1, (str)arg2, (bool)arg3) -> object :
        Return a dict with all nodes from the specified sub-tree.
        arg1: Reference to the ziDAQZoomFFT class.
        arg2: Path string of the node. Use wild card to
            select all.
        arg3[optional]: Specify which type of data structure to return.
            Return data either as a flat dict (True) or as a nested
            dict tree (False). Default = False.

    get( (ziDAQZoomFFT)arg1, (str)arg2) -> object

listNodes(...)
    listNodes( (ziDAQZoomFFT)arg1, (str)arg2, (int)arg3) -> list :
        This function returns a list of node names found at the specified path.
        arg1: Reference to the pyDAQZoomFFT class.
        arg2: Path for which the nodes should be listed. The path may
            contain wildcards so that the returned nodes do not
            necessarily have to have the same parents.
        arg3: Enum that specifies how the selected nodes are listed.
            ziPython.ziListEnum.none -> 0x00
                The default flag, returning a simple
                listing if the given node
            ziPython.ziListEnum.recursive -> 0x01
                Returns the nodes recursively
            ziPython.ziListEnum.absolute -> 0x02
                Returns absolute paths
            ziPython.ziListEnum.leafsonly -> 0x04
                Returns only nodes that are leafs,
                which means the they are at the
                outermost level of the tree.
            ziPython.ziListEnum.settingsonly -> 0x08
                Returns only nodes which are marked
                as setting
        Or combinations of flags can be used.

progress(...)
    progress( (ziDAQZoomFFT)arg1) -> object :
        Reports the progress of the measurement with a number between
        0 and 1.

read(...)
    read( (ziDAQZoomFFT)arg1, (bool)arg2) -> object :
        Read zoom FFT data. If the zoom FFT is still ongoing only a subset
        of zoom FFT data is returned.
        arg1[optional]: Specify which type of data structure to return.
            Return data either as a flat dict (True) or as a nested
            dict tree (False). Default = False.

    read( (ziDAQZoomFFT)arg1) -> object

save(...)
    save( (ziDAQZoomFFT)arg1, (str)arg2) -> None :
        Save zoom FFT data to file.
        arg1: Reference to the ziDAQZoomFFT class.
        arg2: File name string (without extension).

set(...)
    set( (ziDAQZoomFFT)arg1, (str)arg2, (float)arg3) -> None :
        Zoom FFT Parameters

```

Path name	Type	Description
zoomFFT/device	string	Device that should be used for the zoom FFT, e.g. 'dev99'.
zoomFFT/bit	int	Number of FFT points 2^{bit}
zoomFFT/mode	int	Zoom FFT mode 0 = Perform FFT on X+iY 1 = Perform FFT on R 2 = Perform FFT on Phase
zoomFFT/loopcount	int	Number of zoomFFT loops (default 1)
zoomFFT/endless	int	Perform endless zoomFFT (default 0) 0 = endless off, use loopcount 1 = endless on, ignore loopcount
zoomFFT/overlap	double	FFT overlap 0 = none, [0..1]
zoomFFT/settling/time	double	Settling time before measurement is performed.
zoomFFT/settling/tc	double	Settling time in time constant units before the FFT recording is started. 5 ~ low precision 15 ~ medium precision 50 ~ high precision
zoomFFT/window	int	FFT window (default 1 = Hann) 0 = Rectangular 1 = Hann 2 = Hamming 3 = Blackman Harris 4 term
zoomFFT/absolute	bool	Shifts the frequencies so that the center frequency becomes the demodulation frequency rather than 0 Hz.
set((ziDAQZoomFFT)arg1, (str)arg2, (int)arg3) -> None		
set((ziDAQZoomFFT)arg1, (str)arg2, (str)arg3) -> None		
set((ziDAQZoomFFT)arg1, (object)arg2) -> None :		
arg1: Reference to the ziDAQZoomFFT class.		
arg2: A list of path/value pairs.		
subscribe(...)		
subscribe((ziDAQZoomFFT)arg1, (str)arg2) -> None :		
Subscribe to one or several nodes. After subscription the zoom FFT process can be started with the 'execute' command. During the zoom FFT process paths can not be subscribed or unsubscribed.		
arg1: Reference to the ziDAQZoomFFT class.		
arg2: Path string of the node. Use wild card to select all. Alternatively also a list of path strings can be specified.		
trigger(...)		
trigger((ziDAQZoomFFT)arg1) -> None :		
Execute a manual trigger.		
unsubscribe(...)		
unsubscribe((ziDAQZoomFFT)arg1, (str)arg2) -> None :		
Unsubscribe from one or several nodes. During the zoom FFT process paths can not be subscribed or unsubscribed.		
arg1: Reference to the ziDAQZoomFFT class.		
arg2: Path string of the node. Use wild card to select all. Alternatively also a list of path strings can be specified.		

Data and other attributes defined here:		
__init__ = <built-in function __init__>		
Raises an exception		
This class cannot be instantiated from Python		

```
| -----  
| Data descriptors inherited from Boost.Python.instance:  
|  
|   __dict__  
|   __weakref__  
|  
| -----  
| Data and other attributes inherited from Boost.Python.instance:  
|  
|   __new__ = <built-in method __new__ of Boost.Python.class object>  
|           T.__new__(S, ...) -> a new object with type S, a subtype of T
```

4.4.7. Help for ziPython's **ziDAQRecorder** class

An instance of **ziDAQRecorder** is initialized using the **record** method from **ziDAQServer**:

```
>>> help('zhinst.ziPython.ziDAQServer.record')
```

Help on method **record** in **zhinst.ziPython.ziDAQServer**:

```
zhinst.ziPython.ziDAQServer.record = record(...) unbound zhinst.ziPython.ziDAQServer  
method  
record( (ziDAQServer)arg1, (float)arg2, (long)arg3, (int)arg4) -> ziDAQRecorder :  
    Create a recording class. This will start a thread for asynchronous  
    recording.  
    arg1: Reference to the ziDAQServer class.  
    arg2: Maximum recording time for single triggers in [s].  
    arg3: Timeout in [ms]. Recommended value is 500ms.  
    arg4[optional]: Record flags.  
                   FILL = 0x0001 : Fill holes.  
                   ALIGN = 0x0002 : Align data that contains a  
                                   timestamp.  
                   THROW = 0x0004 : Throw EOFError exception if  
                                   sample loss is detected.  
  
record( (ziDAQServer)arg1, (float)arg2, (long)arg3) -> ziDAQRecorder
```

Reference help for the **ziDAQRecorder** class.

```
>>> help('zhinst.ziPython.ziDAQRecorder')
```

Help on class **ziDAQRecorder** in **zhinst.ziPython**:

```
zhinst.ziPython.ziDAQRecorder = class ziDAQRecorder(Boost.Python.instance)  
| Method resolution order:  
|   ziDAQRecorder  
|   Boost.Python.instance  
|   __builtin__.object  
|  
| Methods defined here:  
|  
|   __reduce__ = <unnamed Boost.Python function>(...)  
|  
|   clear(...)  
|       clear( (ziDAQRecorder)arg1) -> None :  
|           End the recording thread.  
|  
|   execute(...)  
|       execute( (ziDAQRecorder)arg1) -> None :
```

```

        Start the recorder. After that command any trigger will start
        the measurement. Subscription or unsubscription is no more
        possible until the recording is finished.

finish(...)
    finish( (ziDAQRecorder)arg1) -> None :
        Stop recording. The recording may be restarted by calling
        'execute' again.

finished(...)
    finished( (ziDAQRecorder)arg1) -> bool :
        Check if the recording has finished. Returns True if finished.

get(...)
    get( (ziDAQRecorder)arg1, (str)arg2, (bool)arg3) -> object :
        Return a dict with all nodes from the specified sub-tree.
        arg1: Reference to the ziDAQRecorder class.
        arg2: Path string of the node. Use wild card to
            select all.
        arg3[optional]: Specify which type of data structure to return.
            Return data either as a flat dict (True) or as a nested
            dict tree (False). Default = False.

    get( (ziDAQRecorder)arg1, (str)arg2) -> object

listNodes(...)
    listNodes( (ziDAQRecorder)arg1, (str)arg2, (int)arg3) -> list :
        This function returns a list of node names found at the specified path.
        arg1: Reference to the ziDAQRecorder class.
        arg2: Path for which the nodes should be listed. The path may
            contain wildcards so that the returned nodes do not
            necessarily have to have the same parents.
        arg3: Enum that specifies how the selected nodes are listed.
            ziPython.ziListEnum.none -> 0x00
                The default flag, returning a simple
                listing if the given node
            ziPython.ziListEnum.recursive -> 0x01
                Returns the nodes recursively
            ziPython.ziListEnum.absolute -> 0x02
                Returns absolute paths
            ziPython.ziListEnum.leafsonly -> 0x04
                Returns only nodes that are leafs,
                which means the they are at the
                outermost level of the tree.
            ziPython.ziListEnum.settingsonly -> 0x08
                Returns only nodes which are marked
                as setting
        Or combinations of flags can be used.

progress(...)
    progress( (ziDAQRecorder)arg1) -> object :
        Reports the progress of the measurement with a number between
        0 and 1.

read(...)
    read( (ziDAQRecorder)arg1, (bool)arg2) -> object :
        Read recorded data. If the recording is still ongoing only a subset
        of recorded data is returned. If many triggers or huge data sets
        are recorded call this method to keep memory usage reasonable.
        arg1[optional]: Specify which type of data structure to return.
            Return data either as a flat dict (True) or as a nested
            dict tree (False). Default = False.

    read( (ziDAQRecorder)arg1) -> object

save(...)
    save( (ziDAQRecorder)arg1, (str)arg2) -> None :

```

```

        Save trigger data to file.
        arg1: Reference to the ziDAQRecorder class.
        arg2: File name string (without extension).

set(...)
set( (ziDAQRecorder)arg1, (str)arg2, (float)arg3) -> None :
    Trigger Parameters
    Path name
    trigger/buffersize      double  Overwrite the buffersize [s] of the
                                trigger object (set when it was
                                instantiated). Recommended buffer size
                                is 2*trigger/0/duration.
    trigger/device          string  The device ID to execute the software
                                trigger, e.g. dev123 (compulsory
                                parameter).
    trigger/endless         bool    Enable endless triggering:
                                1 = enable,
                                0 = disable.
    trigger/forcetrigger    bool    Force a trigger.
    trigger/0/path          string  The path to the demod sample to trigger
                                on, e.g. demods/3/sample, see also
                                trigger/0/source.
    trigger/0/source        int     Signal that is used to trigger on.
                                0 = x [X_SOURCE]
                                1 = y [Y_SOURCE]
                                2 = r [R_SOURCE]
                                3 = angle [ANGLE_SOURCE]
                                4 = frequency [FREQUENCY_SOURCE]
                                5 = phase [PHASE_SOURCE]
                                6 = auxiliary input 0 [AUXIN0_SOURCE]
                                7 = auxiliary input 1 [AUXIN1_SOURCE]
    trigger/0/count         int     Number of trigger edges to record.
    trigger/0/type          int     Trigger type used. Some parameters are
                                only valid for special trigger types.
                                0 = trigger off
                                1 = analog edge trigger on source
                                2 = digital trigger mode on DIO
                                3 = analog pulse trigger on source
                                4 = analog tracking trigger on source
    trigger/0/edge          int     Trigger edge
                                1 = rising edge
                                2 = falling edge
                                3 = both
    trigger/0/findlevel     bool    Automatically find the value of
                                trigger/0/level based on
                                the current signal value.
    trigger/0/bits          int     Digital trigger condition.
    trigger/0/bitmask       int     Bit masking for bits used for
                                triggering. Used for digital trigger.
    trigger/0/delay         double  Trigger frame position [s] (left side)
                                relative to trigger edge.
                                delay = 0 -> trigger edge at left
                                    border.
                                delay < 0 -> trigger edge inside
                                    trigger frame (pretrigger).
                                delay > 0 -> trigger edge before
                                    trigger frame (posttrigger).
    trigger/0/duration      double  Recording frame length [s].
    trigger/0/level         double  Trigger level voltage [V].
    trigger/0/hysteresis    double  Trigger hysteresis [V].
    trigger/0/retrigger     int     Record more than one trigger in a
                                trigger frame.
    trigger/triggered       bool    Has the software trigger triggered?
                                1=Yes, 0=No (read only).
    trigger/0/bandwidth     double  Filter bandwidth [Hz] for pulse and
                                tracking triggers.
    trigger/0/holdoff/count int     Number of skipped triggers until the

```



```

        trigger/0/holdoff/time    double    Hold off time [s] before the next
                                     trigger is recorded again. A hold off
                                     time smaller than the duration will
                                     produce overlapped trigger frames.
        trigger/N/hwtrigsources  int        Only available for devices that support
                                     hardware triggering. Specify the channel
                                     to trigger on.
        trigger/0/pulse/min      double    Minimal pulse width [s] for the pulse
                                     trigger.
        trigger/0/pulse/max      double    Maximal pulse width [s] for the pulse
                                     trigger.
        trigger/filename         string    This parameter is deprecated. If
                                     specified, i.e. not empty, it enables
                                     automatic saving of data in single
                                     trigger mode (trigger/endless = 0).
        trigger/savepath         string    The directory where files are saved when
                                     saving data.
        trigger/fileformat       string    The format of the file for saving data:
                                     0 = Matlab,
                                     1 = CSV.
        trigger/historylength    bool      Maximum number of entries stored in the
                                     measurement history.
        trigger/clearhistory     bool      Remove all records from the history list.

set( (ziDAQRecorder)arg1, (str)arg2, (int)arg3) -> None

set( (ziDAQRecorder)arg1, (str)arg2, (str)arg3) -> None

set( (ziDAQRecorder)arg1, (object)arg2) -> None :
    arg1: Reference to the ziDAQRecorder class.
    arg2: A list of path/value pairs.

subscribe(...)
    subscribe( (ziDAQRecorder)arg1, (str)arg2) -> None :
        Subscribe to one or several nodes. After subscription the recording
        process can be started with the 'execute' command. During the
        recording process paths can not be subscribed or unsubscribed.
        arg1: Reference to the ziDAQRecorder class.
        arg2: Path string of the node. Use wild card to
        select all. Alternatively also a list of path
        strings can be specified.

trigger(...)
    trigger( (ziDAQRecorder)arg1) -> None :
        Execute a manual trigger.

unsubscribe(...)
    unsubscribe( (ziDAQRecorder)arg1, (str)arg2) -> None :
        Unsubscribe from one or several nodes. During the
        recording process paths can not be subscribed or unsubscribed.
        arg1: Reference to the ziDAQRecorder class.
        arg2: Path string of the node. Use wild card to
        select all. Alternatively also a list of path
        strings can be specified.

-----
Data and other attributes defined here:

__init__ = <built-in function __init__>
    Raises an exception
    This class cannot be instantiated from Python

-----
Data descriptors inherited from Boost.Python.instance:

__dict__

```

```
|  __weakref__  
|  
| -----  
| Data and other attributes inherited from Boost.Python.instance:  
|  
|  __new__ = <built-in method __new__ of Boost.Python.class object>  
|          T.__new__(S, ...) -> a new object with type S, a subtype of T
```

4.4.8. Help for ziPython's **ziPllAdvisor** class

An instance of **ziPllAdvisor** is initialized using the **pllAdvisor** method from **ziDAQServer**:

```
>>> help('zhinst.ziPython.ziDAQServer.pllAdvisor')
```

Help on method **pllAdvisor** in **zhinst.ziPython.ziDAQServer**:

```
zhinst.ziPython.ziDAQServer.pllAdvisor = pllAdvisor(...) unbound  
zhinst.ziPython.ziDAQServer method  
  pllAdvisor( (ziDAQServer)arg1, (long)arg2) -> ziPllAdvisor :  
    Create a pllAdvisor class. This will start a thread for running an  
    asynchronous pllAdvisor.  
    arg1: Reference to the ziDAQServer class.  
    arg2: Timeout in [ms]. Recommended value is 500ms.
```

Reference help for the **ziPllAdvisor** class.

```
>>> help('zhinst.ziPython.ziPllAdvisor')
```

Help on class **ziPllAdvisor** in **zhinst.ziPython**:

```
zhinst.ziPython.ziPllAdvisor = class ziPllAdvisor(Boost.Python.instance)  
| Method resolution order:  
|   ziPllAdvisor  
|   Boost.Python.instance  
|   __builtin__.object  
|  
| Methods defined here:  
|  
|  __reduce__ = <unnamed Boost.Python function>(...)  
|  
|  clear(...)  
|      clear( (ziPllAdvisor)arg1) -> None :  
|          End the pllAdvisor thread.  
|  
|  execute(...)  
|      execute( (ziPllAdvisor)arg1) -> None :  
|          Starts the pllAdvisor if not yet running.  
|  
|  finish(...)  
|      finish( (ziPllAdvisor)arg1) -> None :  
|          Stop the pllAdvisor.  
|  
|  finished(...)  
|      finished( (ziPllAdvisor)arg1) -> bool :  
|          Check if the command execution has finished. Returns True if finished.  
|  
|  get(...)  
|      get( (ziPllAdvisor)arg1, (str)arg2, (bool)arg3) -> object :  
|          Return a dict with all nodes from the specified sub-tree.  
|          arg1: Reference to the ziPllAdvisor class.
```

```

    arg2: Path string of the node. Use wild card to
        select all.
    arg3[optional]: Specify which type of data structure to return.
        Return data either as a flat dict (True) or as a nested
        dict tree (False). Default = False.

    get( (ziPllAdvisor)arg1, (str)arg2) -> object

listNodes(...)
    listNodes( (ziPllAdvisor)arg1, (str)arg2, (int)arg3) -> list :
        This function returns a list of node names found at the specified path.
        arg1: Reference to the ziPllAdvisor class.
        arg2: Path for which the nodes should be listed. The path may
            contain wildcards so that the returned nodes do not
            necessarily have to have the same parents.
        arg3: Enum that specifies how the selected nodes are listed.
            ziPython.ziListEnum.none -> 0x00
                The default flag, returning a simple
                listing if the given node
            ziPython.ziListEnum.recursive -> 0x01
                Returns the nodes recursively
            ziPython.ziListEnum.absolute -> 0x02
                Returns absolute paths
            ziPython.ziListEnum.leafsonly -> 0x04
                Returns only nodes that are leafs,
                which means the they are at the
                outermost level of the tree.
            ziPython.ziListEnum.settingsonly -> 0x08
                Returns only nodes which are marked
                as setting
            Or combinations of flags can be used.

progress(...)
    progress( (ziPllAdvisor)arg1) -> object :
        Reports the progress of the command with a number between
        0 and 1.

read(...)
    read( (ziPllAdvisor)arg1, (bool)arg2) -> object :
        Read pllAdvisor data. If the simulation is still ongoing only a subset
        of the data is returned.
        arg1[optional]: Specify which type of data structure to return.
            Return data either as a flat dict (True) or as a nested
            dict tree (False). Default = False.

    read( (ziPllAdvisor)arg1) -> object

save(...)
    save( (ziPllAdvisor)arg1, (str)arg2) -> None :
        Save PLL advisor data to file.
        arg1: Reference to the ziPllAdvisor class.
        arg2: File name string (without extension).

set(...)
    set( (ziPllAdvisor)arg1, (str)arg2, (float)arg3) -> None :
        PLL Advisor Parameters
        Path name      Type      Description
        pllAdvisor/bode    struct    Output parameter. Contains the resulting
                                   bode plot of the PLL simulation.
        pllAdvisor/calculate int      Command to calculate values. Set to 1 to
                                   start the calculation.
        pllAdvisor/center  double   Center frequency of the PLL oscillator.
                                   The PLL frequency shift is relative to
                                   this center frequency.
        pllAdvisor/d        double   Differential gain.
        pllAdvisor/demodbw  double   Demodulator bandwidth used for the PLL
                                   loop filter

```

pllAdvisor/i	double	Integral gain.
pllAdvisor/mode	double	Select PLL Advisor mode. Currently only one mode (open loop) is supported.
pllAdvisor/order	double	Demodulator order used for the PLL loop filter.
pllAdvisor/p	double	Proportional gain.
pllAdvisor/pllbw	double	Demodulator bandwidth used for the PLL loop filter.
pllAdvisor/pm	double	Output parameter. Simulated phase margin of the PLL with the current settings. The phase margin should be greater than 45 deg and preferably greater than 65 deg for stable conditions.
pllAdvisor/pmfreq	double	Output parameter. Simulated phase margin frequency.
pllAdvisor/q	double	Quality factor. Currently not used.
pllAdvisor/rate	double	PLL Advisor sampling rate of the PLL control loop.
pllAdvisor/stable	int	Output parameter. When 1, the PLL Advisor found a stable solution with the given settings. When 0, revise your settings and rerun the PLL Advisor.
pllAdvisor/targetbw	int	Requested PLL bandwidth. Higher frequencies may need manual tuning.
pllAdvisor/targetfail	int	Output parameter. 1 indicates the simulated PLL BW is smaller than the Target BW.
set((ziPllAdvisor)arg1, (str)arg2, (int)arg3) -> None		
set((ziPllAdvisor)arg1, (str)arg2, (str)arg3) -> None		
set((ziPllAdvisor)arg1, (object)arg2) -> None :		
arg1: Reference to the ziPllAdvisor class.		
arg2: A list of path/value pairs.		
subscribe(...)		
subscribe((ziPllAdvisor)arg1, (str)arg2) -> None :		
Subscribe to one or several nodes.		
trigger(...)		
trigger((ziPllAdvisor)arg1) -> None :		
Not applicable to this module.		
unsubscribe(...)		
unsubscribe((ziPllAdvisor)arg1, (str)arg2) -> None :		
Unsubscribe from one or several nodes.		

Data and other attributes defined here:		
__init__ = <built-in function __init__>		
Raises an exception		
This class cannot be instantiated from Python		

Data descriptors inherited from Boost.Python.instance:		
__dict__		
__weakref__		

Data and other attributes inherited from Boost.Python.instance:		
__new__ = <built-in method __new__ of Boost.Python.class object>		
T.__new__(S, ...) -> a new object with type S, a subtype of T		

4.4.9. Help for ziPython's **ziPidAdvisor** class

An instance of **ziPidAdvisor** is initialized using the **pidAdvisor** method from **ziDAQServer**:

```
>>> help('zhinst.ziPython.ziDAQServer.pidAdvisor')
```

Help on method **pidAdvisor** in **zhinst.ziPython.ziDAQServer**:

```
zhinst.ziPython.ziDAQServer.pidAdvisor = pidAdvisor(...) unbound
zhinst.ziPython.ziDAQServer method
  pidAdvisor( (ziDAQServer)arg1, (long)arg2) -> ziPidAdvisor :
    Create a pidAdvisor class. This will start a thread for running an
    asynchronous pidAdvisor.
    arg1: Reference to the ziDAQServer class.
    arg2: Timeout in [ms]. Recommended value is 500ms.
```

Reference help for the **ziPidAdvisor** class.

```
>>> help('zhinst.ziPython.ziPidAdvisor')
```

Help on class **ziPidAdvisor** in **zhinst.ziPython**:

```
zhinst.ziPython.ziPidAdvisor = class ziPidAdvisor(Boost.Python.instance)
| Method resolution order:
|   ziPidAdvisor
|   Boost.Python.instance
|   __builtin__.object
|
| Methods defined here:
|
|   __reduce__ = <unnamed Boost.Python function>(...)
|
|   clear(...)
|       clear( (ziPidAdvisor)arg1) -> None :
|           End the pidAdvisor thread.
|
|   execute(...)
|       execute( (ziPidAdvisor)arg1) -> None :
|           Starts the pidAdvisor if not yet running.
|
|   finish(...)
|       finish( (ziPidAdvisor)arg1) -> None :
|           Stop the pidAdvisor.
|
|   finished(...)
|       finished( (ziPidAdvisor)arg1) -> bool :
|           Check if the command execution has finished. Returns True if finished.
|
|   get(...)
|       get( (ziPidAdvisor)arg1, (str)arg2, (bool)arg3) -> object :
|           Return a dict with all nodes from the specified sub-tree.
|           arg1: Reference to the ziPidAdvisor class.
|           arg2: Path string of the node. Use wild card to
|               select all.
|           arg3[optional]: Specify which type of data structure to return.
|               Return data either as a flat dict (True) or as a nested
|               dict tree (False). Default = False.
|
|       get( (ziPidAdvisor)arg1, (str)arg2) -> object
```

```

listNodes(...)
    listNodes( (ziPidAdvisor)arg1, (str)arg2, (int)arg3) -> list :
        This function returns a list of node names found at the specified path.
        arg1: Reference to the ziPidAdvisor class.
        arg2: Path for which the nodes should be listed. The path may
            contain wildcards so that the returned nodes do not
            necessarily have to have the same parents.
        arg3: Enum that specifies how the selected nodes are listed.
            ziPython.ziListEnum.none -> 0x00
                The default flag, returning a simple
                listing of the given node
            ziPython.ziListEnum.recursive -> 0x01
                Returns the nodes recursively
            ziPython.ziListEnum.absolute -> 0x02
                Returns absolute paths
            ziPython.ziListEnum.leafonly -> 0x04
                Returns only leaf nodes,
                which means the they are at the
                outermost level of the tree.
            ziPython.ziListEnum.settingsonly -> 0x08
                Returns only nodes which are marked
                as settings
            Or combinations of flags can be used.

progress(...)
    progress( (ziPidAdvisor)arg1) -> object :
        Reports the progress of the command with a number between
        0 and 1.

read(...)
    read( (ziPidAdvisor)arg1, (bool)arg2) -> object :
        Read pidAdvisor data. If the simulation is still ongoing, only a subset
        of the data is returned.
        arg1[optional]: Specify which type of data structure to return.
            Return data either as a flat dict (True) or as a nested
            dict tree (False). Default = False.

    read( (ziPidAdvisor)arg1) -> object

save(...)
    save( (ziPidAdvisor)arg1, (str)arg2) -> None :
        Save PID advisor data to file.
        arg1: Reference to the ziPidAdvisor class.
        arg2: File name string (without extension).

set(...)
    set( (ziPidAdvisor)arg1, (str)arg2, (float)arg3) -> None :
        PID Advisor Parameters
        Path name          Type    Description
        pidAdvisor/advancedmode  int    Disable automatic calculation of
                                         the start and stop value.
        pidAdvisor/auto          int    Automatic response calculation
                                         triggered by parameter change.
        pidAdvisor/bode          struct  Output parameter. Contains the
                                         resulting bode plot of the PID
                                         simulation.
        pidAdvisor/bw           double  Output parameter. Calculated system
                                         bandwidth.
        pidAdvisor/calculate     int    In/Out parameter. Command to
                                         calculate values. Set to 1 to start
                                         the calculation.
        pidAdvisor/display/freqstart double  Start frequency for Bode plot.
                                         For disabled advanced mode the
                                         start value is automatically
                                         derived from the system properties.
        pidAdvisor/display/freqstop double  Stop frequency for Bode plot.
        pidAdvisor/display/timestart double  Start time for step response.

```

pidAdvisor/display/timestop	double	Stop time for step response.
pidAdvisor/dut/bw	double	Bandwith of the DUT (device under test).
pidAdvisor/dut/damping	double	Damping of the second order low pass filter.
pidAdvisor/dut/delay	double	IO Delay of the feedback system describing the earliest response for a step change.
pidAdvisor/dut/fcenter	double	Resonant frequency of the of the modelled resonator.
pidAdvisor/dut/gain	double	Gain of the DUT transfer function.
pidAdvisor/dut/q	double	quality factor of the modelled resonator.
pidAdvisor/dut/source	int	Type of model used for the external device to be controlled by the PID. source = 1: Lowpass first order source = 2: Lowpass second order source = 3: Resonator frequency source = 4: Internal PLL source = 5: VCO source = 6: Resonator amplitude
pidAdvisor/impulse	struct	Output parameter. Impulse response (not yet supported).
pidAdvisor/index	int	PID index for parameter detection.
pidAdvisor/pid/autobw	int	Adjusts the demodulator bandwidth to fit best to the specified target bandwidth of the full system.
pidAdvisor/pid/d	double	In/Out parameter. Differential gain.
pidAdvisor/pid/dlimittimeconstant	double	In/Out parameter. Differential filter timeconstant.
pidAdvisor/pid/i	double	In/Out parameter. Integral gain.
pidAdvisor/pid/mode	double	Select PID Advisor mode. Mode value is bit coded, bit 0: P, bit 1: I, bit 2: D, bit 3: D filter limit.
pidAdvisor/pid/p	double	In/Out parameter. Proportional gain.
pidAdvisor/pid/rate	double	In/Out parameter. PID Advisor sampling rate of the PID control loop.
pidAdvisor/pid/targetbw	double	PID system target bandwidth.
pidAdvisor/pm	double	Output parameter. Simulated phase margin of the PID with the current settings. The phase margin should be greater than 45 deg and preferably greater than 65 deg for stable conditions.
pidAdvisor/pmfreq	double	Output parameter. Simulated phase margin frequency.
pidAdvisor/stable	int	Output parameter. When 1, the PID Advisor found a stable solution with the given settings. When 0, revise your settings and rerun the PID Advisor.
pidAdvisor/step	struct	Output parameter. Contains the resulting step response plot of the PID simulation.
pidAdvisor/targetbw	double	Requested PID bandwidth. Higher frequencies may need manual tuning.
pidAdvisor/targetfail	int	Output parameter. 1 indicates the simulated PID BW is smaller than the Target BW.
pidAdvisor/tf/closedloop	int	Switch the response calculation mode between closed or open loop.
pidAdvisor/tf/input	int	Start point for the plant response simulation for open or closed

```

        pidAdvisor/tf/output      int      loops.
                                         End point for the plant response
                                         simulation for open or closed
                                         loops.
        pidAdvisor/tune           int      Optimize the PID parameters so that
                                         the noise of the closed-loop system
                                         gets minimized.

set( (ziPidAdvisor)arg1, (str)arg2, (int)arg3) -> None

set( (ziPidAdvisor)arg1, (str)arg2, (str)arg3) -> None

set( (ziPidAdvisor)arg1, (object)arg2) -> None :
    arg1: Reference to the ziPidAdvisor class.
    arg2: A list of path/value pairs.

subscribe(...)
    subscribe( (ziPidAdvisor)arg1, (str)arg2) -> None :
        Subscribe to one or several nodes.

trigger(...)
    trigger( (ziPidAdvisor)arg1) -> None :
        Not applicable to this module.

unsubscribe(...)
    unsubscribe( (ziPidAdvisor)arg1, (str)arg2) -> None :
        Unsubscribe from one or several nodes.

-----
Data and other attributes defined here:

__init__ = <built-in function __init__>
    Raises an exception
    This class cannot be instantiated from Python

-----
Data descriptors inherited from Boost.Python.instance:

__dict__
__weakref__

-----
Data and other attributes inherited from Boost.Python.instance:

__new__ = <built-in method __new__ of Boost.Python.class object>
    T.__new__(S, ...) -> a new object with type S, a subtype of T

```

Chapter 5. LabVIEW Programming

Interfacing with your Zurich Instruments device via National Instruments' [LabVIEW®](#) is an efficient choice in terms of development time and run-time performance. LabVIEW is a graphical programming language designed to interface with laboratory equipment via so-called VIs ("virtual instruments"), whose key strength is the ease of displaying dynamic signals obtained from your instrument.

This chapter aims to help you get started using the Zurich Instruments LabOne LabVIEW API to control your instrument, please refer to:

- [Section 5.1](#) for help [Installing the LabOne LabVIEW API](#) .
- [Section 5.2.2](#) for an introduction to [LabOne LabVIEW Programming Concepts](#) .
- [Section 5.2.3](#) for [Finding help for the LabOne VIs from within LabVIEW](#) .
- [Section 5.2.4](#) for help [Finding the LabOne LabVIEW API Examples](#) .
- [Section 5.2.5](#) for help [Running the LabOne Example VIs](#) .
- [Section 5.3](#) for some [LabVIEW Programming Tips and Tricks](#) .

Note

This section and the provided examples are no substitute for a general LabVIEW tutorial. See, for example, the National Instruments [website](#) for help to get started programming with LabVIEW.

5.1. Installing the LabOne LabVIEW API

5.1.1. Requirements

A LabVIEW 2009 (or higher) installation is required on either Windows or Linux in order to use the LabOne LabVIEW API.

The LabOne LabVIEW API is included in a standard LabOne installation. However, in order to make the LabOne LabVIEW API available for use within LabVIEW, a directory needs to be copied to a specific directory of your LabVIEW installation. The LabOne installer is available from Zurich Instruments' [download page](#) (login required).

5.1.2. Windows

1. Locate the `instr.lib` directory in your LabVIEW installation and delete any previous Zurich Instruments API directories. The `instr.lib` directory is typically located at:

```
C:\Program Files\National Instruments\LabVIEW 201x\instr.lib\
```

Previous Zurich Instruments installations will be directories located in the `instr.lib` directory that are named either:

- Zurich Instruments HF2, or
- Zurich Instruments LabOne.

These folders may simply be deleted (administrator rights required).

2. On Windows, navigate to the `API\LabVIEW` subdirectory of your LabOne installation, typically found at:

```
C:\Program Files\Zurich Instruments\LabOne\API\LabVIEW\
```

and copy the subdirectory located there

```
Zurich Instruments LabOne
```

to the `instr.lib` directory in your LabVIEW installation as located in Step 1. Note, you will need administrator rights to copy to this directory.

3. Restart LabVIEW and verify your installation as described in [Section 5.1.4](#).

5.1.3. Linux

1. Locate the `instr.lib` directory in your LabVIEW installation and remove any previous Zurich Instruments API installations. The `instr.lib` directory is typically located at:

```
/usr/local/natinst/LabVIEW-201x/instr.lib/
```

Previous Zurich Instruments installations will be folders located in the `instr.lib` directory that are named either:

- Zurich Instruments HF2, or
- Zurich Instruments LabOne.

These folders may simply be deleted (administrator rights required).

2. Navigate to the path where you unpacked LabOne and locate the subdirectory

`[PATH]/LabOneLinux64/API/LabVIEW/`

3. Copy the directory

`[PATH]/LabOneLinux64/API/LabVIEW/Zurich Instruments LabOne/`

to the `instr.lib` directory in your LabVIEW installation as located in [Step 1](#). Note, you will need administrator rights to copy to this directory.

4. Restart LabVIEW and verify your installation as described in [Section 5.1.4](#).

5.1.4. Verifying your Installation

If the LabOne LabVIEW API palette can be accessed from within LabVIEW, the LabOne LabVIEW API is correctly installed. See [Section 5.2.1](#) for help finding the palette.

5.2. Getting Started

5.2.1. Locating the LabOne LabVIEW VI Palette

In order to locate the LabOne LabVIEW VIs start LabVIEW and create a new VI. In the VI's "Block Diagram" (Ctrl-e) you can to access the LabOne LabVIEW API palette with a mouse right-click and browsing the tree under "Instrument I/O" → "Instr. Drivers", see [Figure 5.1](#).




Figure 5.1. Locating the LabOne LabVIEW Palette

5.2.2. LabOne LabVIEW Programming Concepts

As described in [Section 1.1](#) a LabVIEW program communicates to a Zurich Instrument device via a software program running on the PC called the data server. In general, the outline of the instruction flow for a LabVIEW virtual instrument is as following:

1. Initialization: Open a connection from the API to the data server program.
2. Configuration: Perform the instrument's settings. For example, using the virtual instrument `ziSetValueDouble.vi`.
3. Data: Read data from the instrument.
4. Utility: Perform data analysis on the read data, potentially repeating Step 2 and/or Step 3.
5. Close: Terminate the API's connection to the data server program.

The `VI Tree.vi` included the LabOne LabVIEW API demonstrates this flow and lists common VIs used for working with a Zurich Instruments device, see [Figure 5.2](#). The `VI Tree.vi` can be found either via the LabOne VI palette, see [Section 5.2.1](#), or by opening the file in the `Public` folder of your LabOne LabVIEW installation, typically located at:

C:\Program Files\National Instruments\LabVIEW 2012\instr.lib\Zurich Instruments LabOne\Public\VI Tree.vi

Modules (e.g. Sweeper,...) enable additional functionality of your Zurich instrument device in Labview. The outline of the instruction flow for a LabVIEW Module is as following:

1. Initialization: Create a ziModHandle from a ziHandle ziModStart.vi.
2. Configuration: Perform the module's settings. For example, using the virtual instrument ziModSetValue.vi.
3. Subscribe: Define the recorded data node ziModSubscribe.vi.
4. Execute: Start the operation of the module ziModExecute.vi.
5. Data: Read data from the module. For example, using the ziModGetNextNode.vi and ziModGetData.vi.
6. Utility: Perform data analysis on the read data, potentially repeating Step 2, Step 3 and/or Step 4.
7. Clear: Terminate the API's connection to the module ziModClear.vi.




Figure 5.2. An overview of the LabOne LabVIEW VIs is given in VI Tree.vi. Press Ctrl-h after selecting one of the VIs to obtain help.

5.2.3. Finding help for the LabOne VIs from within LabVIEW

As is customary for LabVIEW, built-in help for LabOne's VIs can be obtained by selecting the VI with the mouse in a block diagram and pressing Ctrl-h to view the VI's context help. See [Figure 5.2](#) for an example.

5.2.4. Finding the LabOne LabVIEW API Examples

Many examples come bundled with the LabOne LabVIEW API which demonstrate the most important concepts of working with Zurich Instrument devices. The easiest way to browse the list of available examples is via the NI Example Finder: In LabVIEW select "Find Examples..." from the "Help" menu-bar and search for "LabOne", see [Figure 5.3](#).

The examples are located in the directory `instr.lib/Zurich Instruments LabOne/Examples` found in LabVIEW installation directory. In order to modify an example for your needs, please copy it to your local workspace.




Figure 5.3. Search for "LabOne" in NI's Example Finder to find examples to run with your instrument.

5.2.5. Running the LabOne Example VIs

This section describes how to run a LabOne LabVIEW example on your instrument.

Note

Please ensure that the example you would like to run is supported by your instrument class and its options set. For example, examples for HF2 Instruments can be found in the Example Finder

(see [Section 5.2.4](#)) by searching for "HF2", examples for the UHFLI by searching for "UHFLI" and examples for the MFLI by searching for "MFLI".

Device Connection

After opening one of the LabOne LabVIEW examples, please ensure that the example is configured to run on the desired instrument type. `ziOpen.vi` establishes a connection to a Data Server. The address is of the format `{<host>}{:<port>}::<Device ID>}`. Usually it is sufficient to provide the Device ID only highlighted in [Figure 5.4](#). The host and port are then determined by network discovery. Should the discovery not work, prepend `<host>:<port>::` to the Device ID. Examples are "myhf2.company.com:8004::dev466" or "myhf2.company.com:8004". In the latter case the first found instrument on the data server listening on "myhf2.company.com:8004" will be selected.




Figure 5.4. LabOne LabVIEW Example Poll Data: Device selection.

Running the VI and Block Diagram

The example can be ran as any LabVIEW program; by clicking the "Run" icon in the icon bar. Be sure to check the example's code and explanation by pressing Ctrl-e to view the example's block diagram, see [Figure 5.5](#).




Figure 5.5. LabOne LabVIEW Example Poll Data: Block Diagram.

5.3. LabVIEW Programming Tips and Tricks

Use the User Interface's command log or Server's text interface while programming with LabVIEW

As with all other interfaces, LabVIEW uses the "path" and "nodes" concept to address settings on an instrument, see [Section 1.1](#). In order to learn about or verify the nodes available it can be very helpful to view the command log in the User Interface (see the bar in the bottom of the screen) to see which node has been configured during a previous setting change. The text interface (HF2 Series) provides a convenient way to explore the node hierarchy.

Always close ziHandles and ziModHandles or LabVIEW runs out of memory

If you use the "Abort Execution" button of LabVIEW, your LabVIEW program will not close any existing connections to the ziServer. Any open connection inside of LabVIEW will persist and continue to consume about 12 MB of RAM so that with time you will run out of memory. Completely exit LabVIEW in order to release the memory again.

Use shift registers

The structure of efficient LabVIEW code is distinguished by signals being "piped through" by use of shift registers in loops and by the absence of object replication. Using shift registers in LabVIEW avoids copying of data and, more important, running the garbage collector frequently.




Figure 5.6. Examples of simple LabVIEW programs for the Zurich Instruments HF2 Series. Left: A well implemented loop, Right: An example for-loop gone wrong.

Chapter 6. C Programming

The LabOne C API, also known as ziAPI, provides a simple and robust way to communicate with the Data Server. It enables you to get or set parameters and receive streaming data.

6.1. Getting Started

After installing the LabOne software package and relevant drivers for your instrument you are ready start programming with ziAPI. All you need is a C compiler, linker and editor.

The structure of a program using ziAPI can be split into three parts: initialization/connection, data manipulation and disconnection/cleanup. The basic object that is always used is the ziConnection data structure. First, ziConnection is has to be initialized by calling [ziAPIInit](#) . After initialization ziConnection is ready to connect to a ziServer by calling [ziAPIConnect](#) . Then ziConnection is ready to be used for getting and setting parameters and streaming data. When ziConnection is not needed anymore the established connection to the ziServer has to be hung up using [ziAPIDisconnect](#) before cleaning it up by calling [ziAPIDestroy](#) .

6.1.1. Example

Below you find a simple program, which sets the demodulator rate of all demods for all devices.

```
#include <stdlib.h>
#include <stdio.h>

#include "                ziAPI.h

int main()
{

    ZIResult_enum RetVal;
    char* ErrBuffer;

    ZIConnection Conn;

    //initialize the ZIConnection
    if( ( RetVal =
        ziAPIInit
        ( &Conn ) ) !=
        ZI_INFO_SUCCESS
        )
    {

        ziAPIGetError
        ( RetVal, &ErrBuffer, NULL );
        fprintf( stderr, "Can't init Connection: %s\n", ErrBuffer );
        return 1;
    }

    //connect to the ziServer running on localhost
    //using the port 8005 (default)
    if( ( RetVal =
        ziAPIConnect
        ( Conn,
          "localhost",
          8005 ) ) !=
        ZI_INFO_SUCCESS
        )
    {

        ziAPIGetError
        ( RetVal, &ErrBuffer, NULL );
        fprintf( stderr, "Can't connect: %s\n", ErrBuffer );
    }
    else
    {

```

```
//set all demodulator rates of all devices to 150Hz
if( ( RetVal =
    ziAPISetValueD
    ( Conn,
        "*/demods/*/rate",
        150 ) ) !=
    ZI_INFO_SUCCESS
    )
{
    ziAPIGetError
    ( RetVal, &ErrBuffer, NULL );
    fprintf( stderr, "Can't set parameter: %s\n", ErrBuffer );
}

//disconnect from the ziServer
//since ZIAPIDisconnect always returns ZI_INFO_SUCCESS
//no error handling is required

    ziAPIDisconnect
    ( Conn );

}

//destroy the ZIConnection
//since ZIAPIDestroy always returns ZI_INFO_SUCCESS
//no error handling is required

    ziAPIDestroy
    ( Conn );

return 0;
}
```

6.2. Module Documentation

6.2.1. Connecting to Data Server

This section describes how to initialize the `ZIConnection` and establish a connection to Data Server as well as how to disconnect after all data handling is done and cleanup the `ZIConnection`.

Typedefs

- `typedef ZIConnection`
The `ZIConnection` is a connection reference; it holds information and helper variables about a connection to the Data Server. There is nothing in this reference which the user may use, so it is hidden and instead a dummy pointer is used. See [ziAPIInit](#) for how to create a `ZIConnection`.

Functions

- `ZIResult_enum ziAPIInit (ZIConnection * conn)`
Initializes a `ZIConnection` structure.
- `ZIResult_enum ziAPIDestroy (ZIConnection conn)`
Destroys a `ZIConnection` structure.
- `ZIResult_enum ziAPIConnect (ZIConnection conn, const char* hostname, uint16_t port)`
Connects the `ZIConnection` to Data Server.
- `ZIResult_enum ziAPIDisconnect (ZIConnection conn)`
Disconnects an established connection.
- `ZIResult_enum ziAPIListImplementations (char* implementations, uint32_t bufferSize)`
Returns the list of supported implementations.
- `ZIResult_enum ziAPIConnectEx (ZIConnection conn, const char* hostname, uint16_t port, ZIAPIVersion_enum apiLevel, const char* implementation)`
Connects to Data Server and enables extended `ziAPI`.
- `ZIResult_enum ziAPIGetConnectionAPILevel (ZIConnection conn, ZIAPIVersion_enum* apiLevel)`
Returns `ziAPI` level used for the connection `conn`.
- `ZIResult_enum ziAPIGetRevision (unsigned int* revision)`
Retrieves the revision of `ziAPI`.

Detailed Description

```
#include <stdio.h>

#include "
                                ziAPI.h
                                "

int main( )
{

    ZIResult_enum RetVal;
    char* ErrBuffer;

    ZIConnection Conn;

    if( ( RetVal =
                                ziAPIInit
                                ( &Conn ) ) !=
                                ZI_INFO_SUCCESS
                                )
    {

        ziAPIGetError
        ( RetVal, &ErrBuffer, NULL );
        fprintf( stderr, "Can't init Connection: %s\n", ErrBuffer );
        return 1;
    }

    //connect to the ziServer set as default using the default port
    if( ( RetVal =
                                ziAPIConnect
                                ( Conn, NULL, 0 ) ) !=
                                ZI_INFO_SUCCESS
                                )
    {

        ziAPIGetError
        ( RetVal, &ErrBuffer, NULL );
        fprintf( stderr, "Can't connect to ziServer: %s\n", ErrBuffer );
    }
    else
    {
        /*
         * do something using the ZIConnection here
         */

        //since ZIAPIDisconnect always returns ZI_INFO_SUCCESS
        //no error handling is required

        ziAPIDisconnect
        ( Conn );

    }

    //since ZIAPIDestroy always returns ZI_INFO_SUCCESS
    //no error handling is required

    ziAPIDestroy
    ( Conn );

    return 0;
}
```


Function Documentation

ziAPIInit

ZIResult_enum ziAPIInit (**ZIConnection** * conn)

Initializes a **ZIConnection** structure.

This function initializes the structure so that it is ready to connect to Data Server. It allocates memory and sets up the infrastructure needed.

Parameters:

[out] conn
Pointer to **ZIConnection** that is to be initialized

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_MALLOC on memory allocation failure

See Also:

[ziAPIDestroy](#) , [ziAPIConnect](#) , [ziAPIDisconnect](#)

See [Connection](#) for an example

ziAPIDestroy

ZIResult_enum ziAPIDestroy (ZIConnection conn)

Destroys a [ZIConnection](#) structure.

This function frees all memory that has been allocated by [ziAPIInit](#) . If it is called with an uninitialized [ZIConnection](#) struct it may result in segmentation faults as well when it is called with a struct for which ZIAPIDestroy already has been called.

Parameters:

[in] conn

Pointer to [ZIConnection](#) struct that has to be destroyed

Returns:

- ZI_INFO_SUCCESS

See Also:

[ziAPIInit](#) , [ziAPIConnect](#) , [ziAPIDisconnect](#)

See [Connection](#) for an example

ziAPIConnect

ZIResult_enum ziAPIConnect (**ZIConnection** conn, const char* hostname, uint16_t port)

Connects the ZIConnection to Data Server.

Connects to Data Server using a **ZIConnection** and prepares for data exchange. For most cases it is enough to just give a reference to the connection and give NULL for hostname and 0 for the port, so it connects to localhost on the default port.

Parameters:

[in] conn

Pointer to **ZIConnection** with which the connection should be established

[in] hostname

Name of the Host to which it should be connected, if NULL "localhost" will be used as default

[in] port

The Number of the port to connect to. If 0, default port of the local Data Server will be used (8005)

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_HOSTNAME if the given host name could not be found
- ZI_ERROR_SOCKET_CONNECT if no connection could be established
- ZI_ERROR_OVERFLOW when a FIFO overflow occurred
- ZI_ERROR_SOCKET_INIT if initialization of the socket failed
- ZI_ERROR_CONNECTION when the Data Server didn't return the correct answer
- ZI_ERROR_TIMEOUT when initial communication timed out

See Also:

[ziAPIDisconnect](#) , [ziAPIInit](#) , [ziAPIDestroy](#)

See [Connection](#) for an example

ziAPIDisconnect

ZIResult_enum `ziAPIDisconnect (ZIConnection conn)`

Disconnects an established connection.

Disconnects from Data Server. If the connection has not been established and the function is called it returns without doing anything.

Parameters:

[in] conn
Pointer to ZIConnection to be disconnected

Returns:

■ ZI_INFO_SUCCESS

See Also:

[ziAPIConnect](#), [ziAPIInit](#), [ziAPIDestroy](#)

See [Connection](#) for an example

ziAPIListImplementations

ZIResult_enum **ziAPIListImplementations** (char* implementations, uint32_t bufferSize)

Returns the list of supported implementations.

Returned names are defined by implementations in the linked library and may change depending on software version.

Parameters:

[out] implementations

Pointer to a buffer receiving a newline-delimited list of the names of all the supported ziAPI implementations. The string is zero-terminated.

[in] bufferSize

The size of the buffer assigned to the implementations parameter

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_LENGTH if the length of the char-buffer given by MaxLen is too small for all elements

See Also:

[ziAPIConnectEx](#)

ziAPIConnectEx

ZIResult_enum ziAPIConnectEx (**ZIConnection** conn, const char* hostname, uint16_t port, **ZIAPIVersion_enum** apiLevel, const char* implementation)

Connects to Data Server and enables extended ziAPI.

With apiLevel=ZI_API_VERSION_1 and implementation=NULL, this call is equivalent to plain [ziAPIConnect](#) . With other version and implementation values enables corresponding ziAPI extension and connection using different implementation.

Parameters:

[in] conn

Pointer to the ZIConnection with which the connection should be established

[in] hostname

Name of the host to which it should be connected, if NULL "localhost" will be used as default

[in] port

The number of the port to connect to. If 0 the port of the local Data Server will be used

[in] apiLevel

Specifies the ziAPI compatibility level to use for this connection (1 or 4).

[in] implementation

Specifies implementation to use for a connection, must be one of the returned by [ziAPIListImplementations](#) or NULL to select default implementation

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_HOSTNAME if the given host name could not be found
- ZI_ERROR_SOCKET_CONNECT if no connection could be established
- ZI_ERROR_OVERFLOW when a FIFO overflow occurred
- ZI_ERROR_SOCKET_INIT if initialization of the socket failed
- ZI_ERROR_CONNECTION when the Data Server didn't return the correct answer or requested implementation is not found or doesn't support requested ziAPI level
- ZI_ERROR_TIMEOUT when initial communication timed out

See Also:

[ziAPIListImplementations](#) , [ziAPIConnect](#) , [ziAPIDisconnect](#) , [ziAPIInit](#) , [ziAPIDestroy](#) , [ziAPIGetConnectionVersion](#)

See [Connection](#) for an example

ziAPIGetConnectionAPILevel

ZIResult_enum ziAPIGetConnectionAPILevel (**ZIConnection** conn,
ZIAPIVersion_enum* apiLevel)

Returns ziAPI level used for the connection conn.

Parameters:

[in] conn

Pointer to ZIConnection

[out] apiLevel

Pointer to preallocated ZIAPIVersion_enum, receiving the ziAPI level

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_CONNECTION if level can not be determined due to conn is not connected

See Also:

[ziAPIConnectEx](#), [ziAPIGetVersion](#)

ziAPIGetRevision

ZIResult_enum ziAPIGetRevision (unsigned int* revision)

Retrieves the revision of ziAPI.

Sets an unsigned int with the revision (build number) of the ziAPI you are using.

Parameters:

[in] revision

Pointer to an unsigned int to fill up with the revision.

Returns:

- ZI_INFO_SUCCESS

6.2.2. Tree

All parameters and streams are organized in a tree. You can list the whole tree, parts of it or single items using [ziAPIListNodes](#) or you may update the tree with nodes of newly connected devices by using [ziAPIUpdateDevices](#).

Enumerations

- `enum ZIListNodes_enum { ZI_LIST_NODES_NONE, ZI_LIST_NODES_RECURSIVE, ZI_LIST_NODES_ABSOLUTE, ZI_LIST_NODES_LEAFONLY, ZI_LIST_NODES_SETTINGSONLY, ZI_LIST_NONE, ZI_LIST_RECURSIVE, ZI_LIST_ABSOLUTE, ZI_LIST_LEAFONLY, ZI_LIST_SETTINGSONLY }`

Defines the values of the flags used in [ziAPIListNodes](#).

Functions

- `ZIResult_enum ziAPIListNodes (ZIConnection conn, const char* path, char* nodes, int bufferSize, int flags)`
Returns all child nodes found at the specified path.
- `ZIResult_enum ziAPIUpdateDevices (ZIConnection conn)`
Search for the newly connected devices and update the tree.
- `ZIResult_enum ziAPIConnectDevice (ZIConnection conn, const char* deviceSerial, const char* deviceInterface, const char* interfaceParams)`
Connect a device to the server.
- `ZIResult_enum ziAPIDisconnectDevice (ZIConnection conn, const char* deviceSerial)`
Disconnect a device from the server.

Detailed Description

```
#include <stdio.h>

#include "                ziAPI.h
                        "

void PrintChildren( ZIConnection Conn,
                  char* Path )
{
    ZIResult_enum RetVal;
    char* ErrBuffer;
```

```
char NodesBuffer[ 8192 ];

if( ( RetVal =
    ziAPIListNodes
    ( Conn,
      Path,
      NodesBuffer,
      8192,

      ZI_LIST_NODES_NONE
    ) ) !=
    ZI_INFO_SUCCESS
  )
{
    ziAPIGetError
    ( RetVal, &ErrBuffer, NULL );
    fprintf( stderr, "Can't List Nodes: %s\n", ErrBuffer );
}
else
{
    char* Ptr = NodesBuffer;
    char* LastPtr = Ptr;

    //print out each node on a separate line with dash as prefix
    for( ; *Ptr != 0; Ptr++ )
    {
        if( *Ptr == '\n' )
        {
            *Ptr = 0;
            printf("- %s\n", LastPtr);
            LastPtr = Ptr + 1;
        }

    }

    //print out the last node
    if( Ptr != LastPtr )
        printf("- %s\n", LastPtr);
}
}
```

Enumeration Type Documentation

enum ZIListNodes_enum

Defines the values of the flags used in [ziAPIListNodes](#) .

Enumerator:

- **ZI_LIST_NODES_NONE**
Default, return a simple listing of the given node immediate descendants.
- **ZI_LIST_NODES_RECURSIVE**
List the nodes recursively.
- **ZI_LIST_NODES_ABSOLUTE**
Return absolute paths.
- **ZI_LIST_NODES_LEAFONLY**
Return only leaf nodes, which means the nodes at the outermost level of the tree.
- **ZI_LIST_NODES_SETTINGSONLY**
Return only nodes which are marked as setting.
- **ZI_LIST_NONE**
Default, return a simple listing of the given node immediate descendants.
- **ZI_LIST_RECURSIVE**
List the nodes recursively.
- **ZI_LIST_ABSOLUTE**
Return absolute paths.
- **ZI_LIST_LEAFONLY**
Return only leaf nodes, which means the nodes at the outermost level of the tree.
- **ZI_LIST_SETTINGSONLY**
Return only nodes which are marked as setting.

Function Documentation

ziAPIListNodes

ZIResult_enum ziAPIListNodes (**ZIConnection** conn, const char* path, char* nodes, int bufferSize, int flags)

Returns all child nodes found at the specified path.

This function returns a list of node names found at the specified path. The path may contain wildcards so that the returned nodes do not necessarily have to have the same parents. The list is returned in a null-terminated char-buffer, each element delimited by a newline. If the maximum length of the buffer (bufferSize) is not sufficient for all elements, nothing will be returned and the return value will be [ZI_LENGTH](#) .

Parameters:

[in] conn

Pointer to the ZIConnection for which the node names should be retrieved.

[in] path

Path for which all children will be returned. The path may contain wildcard characters.

[out] nodes

Upon call filled with newline-delimited list of the names of all the children found. The string is zero-terminated.

[in] bufferSize

The length of the buffer used for the nodes output parameter.

[in] flags

A combination of flags (applied bitwise)as defined in [ZIListNodes_enum](#) .

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_ERROR_LENGTH if the path's length exceeds [MAX_PATH_LEN](#) or the length of the char-buffer for the nodes given by bufferSize is too small for all elements
- ZI_ERROR_OVERFLOW when a FIFO overflow occurred
- ZI_ERROR_COMMAND on an incorrect answer of the server
- ZI_ERROR_SERVER_INTERNAL if an internal error occurred in Data Server
- ZI_ERROR_NOTFOUND if the given path could not be resolved
- ZI_ERROR_TIMEOUT when communication timed out

See [Tree Listing](#) for an example

See Also:

ziAPIUpdate

ziAPIUpdateDevices

ZIResult_enum ziAPIUpdateDevices (**ZIConnection** conn)

Search for the newly connected devices and update the tree.

This function forces the Data Server to search for newly connected devices and to connect to run them

Parameters:

[in] conn
Pointer to ZIConnection

Returns:

■ ZI_INFO_SUCCESS

See Also:

[ziAPIListNodes](#)

ziAPIConnectDevice

ZIResult_enum ziAPIConnectDevice (**ZIConnection** conn, const char* deviceSerial, const char* deviceInterface, const char* interfaceParams)

Connect a device to the server.

This function connects a device with deviceSerial via the specified deviceInterface for use with the server.

Parameters:

[in] conn

Pointer to the ZIConnection with which the connection should be established

[in] deviceSerial

The serial of the device to connect to, e.g., dev2100

[in] deviceInterface

The interface to use for the connection, e.g., USB|1GbE

[in] interfaceParams

Parameters for interface configuration

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_TIMEOUT when communication timed out

See Also:

[ziAPIDisconnectDevice](#) , [ziAPIConnect](#) , [ziAPIDisconnect](#) , [ziAPIInit](#)

ziAPIDisconnectDevice

ZIResult_enum **ziAPIDisconnectDevice** (**ZIConnection** conn, const char* deviceSerial)

Disconnect a device from the server.

This function disconnects a device specified by deviceSerial from the server.

Parameters:

[in] conn

Pointer to the ZIConnection with which the connection should be established

[in] deviceSerial

The serial of the device to connect to, e.g., dev2100

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_TIMEOUT when communication timed out

See Also:

[ziAPIConnectDevice](#) , [ziAPIConnect](#) , [ziAPIDisconnect](#) , [ziAPIInit](#)

6.2.3. Set and Get Parameters

This section describes several functions for getting and setting parameters of different datatypes.

Functions

- `ZIResult_enum ziAPIGetValueD (ZIConnection conn, const char* path, ZIDoubleData* value)`
gets the double-type value of the specified node
- `ZIResult_enum ziAPIGetValueI (ZIConnection conn, const char* path, ZIIntegerData* value)`
gets the integer-type value of the specified node
- `ZIResult_enum ziAPIGetDemodSample (ZIConnection conn, const char* path, ZIDemodSample * value)`
Gets the demodulator sample value of the specified node.
- `ZIResult_enum ziAPIGetDIOSample (ZIConnection conn, const char* path, ZIDIOSample * value)`
Gets the Digital I/O sample of the specified node.
- `ZIResult_enum ziAPIGetAuxInSample (ZIConnection conn, const char* path, ZIAuxInSample * value)`
gets the AuxIn sample of the specified node
- `ZIResult_enum ziAPIGetValueB (ZIConnection conn, const char* path, unsigned char* buffer, unsigned int* length, unsigned int bufferSize)`
gets the Bytearray value of the specified node
- `ZIResult_enum ziAPISetValueD (ZIConnection conn, const char* path, ZIDoubleData value)`
asynchronously sets a double-type value to one or more nodes specified in the path
- `ZIResult_enum ziAPISetValueI (ZIConnection conn, const char* path, ZIIntegerData value)`
asynchronously sets an integer-type value to one or more nodes specified in a path
- `ZIResult_enum ziAPISetValueB (ZIConnection conn, const char* path, unsigned char* buffer, unsigned int length)`
asynchronously sets the binary-type value of one ore more nodes specified in the path
- `ZIResult_enum ziAPISyncSetValueD (ZIConnection conn, const char* path, ZIDoubleData* value)`
synchronously sets a double-type value to one or more nodes specified in the path
- `ZIResult_enum ziAPISyncSetValueI (ZIConnection conn, const char* path, ZIIntegerData* value)`
synchronously sets an integer-type value to one or more nodes specified in a path

- `ZIResult_enum ziAPISyncSetValueB (ZIConnection conn, const char* path, uint8_t* buffer, uint32_t* length, uint32_t bufferSize)`
Synchronously sets the binary-type value of one ore more nodes specified in the path.
- `ZIResult_enum ziAPISync (ZIConnection conn)`
Synchronizes the session by dropping all pending data.
- `ZIResult_enum ziAPIEchoDevice (ZIConnection conn, const char* deviceSerial)`
Sends an echo command to a device and blocks until answer is received.
- `ZIResult_enum ziAPIAsyncSetDoubleData (ZIConnection conn, const char* path, ZIDoubleData value)`
- `ZIResult_enum ziAPIAsyncSetIntegerData (ZIConnection conn, const char* path, ZIIntegerData value)`
- `ZIResult_enum ziAPIAsyncSetByteArray (ZIConnection conn, const char* path, uint8_t* buffer, uint32_t length)`
- `ZIResult_enum ziAPIGetValueS (ZIConnection conn, char* path, DemodSample * value)`
- `ZIResult_enum ziAPIGetValueDIO (ZIConnection conn, char* path, DIOSample * value)`
- `ZIResult_enum ziAPIGetValueAuxIn (ZIConnection conn, char* path, AuxInSample * value)`

Function Documentation

ziAPIGetValueD

ZIResult_enum ziAPIGetValueD (**ZIConnection** conn, const char* path, ZIDoubleData* value)

gets the double-type value of the specified node

This function retrieves the numerical value of the specified node as an double-type value. The value first found is returned if more than one value is available (a wildcard is used in the path).

Parameters:

[in] conn

Pointer to ZIConnection with which the value should be retrieved

[in] path

Path to the node holding the value

[out] value

Pointer to a double in which the value should be written

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_ERROR_LENGTH if the path's length exceeds MAX_PATH_LEN
- ZI_ERROR_OVERFLOW when a FIFO overflow occurred
- ZI_ERROR_COMMAND on an incorrect answer of the server
- ZI_ERROR_SERVER_INTERNAL if an internal error occurred in Data Server
- ZI_ERROR_NOTFOUND if the given path could not be resolved or no value is attached to the node
- ZI_ERROR_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>
```

```
#include "
                                ziAPI.h
"
```

```
void UpdateValue( ZIConnection Conn )
{
```

```
    ZIResult_enum RetVal;
    char* ErrBuffer;
    ZIDoubleData ValueD;
```

```
    if( ( RetVal =
```

```
        ziAPISetValueI
        ( Conn,
          "DEV1046/demods/*/rate",
          100 ) ) !=
        ZI_INFO_SUCCESS
    )
{

    ziAPIGetError
    ( RetVal, &ErrBuffer, NULL );
    fprintf( stderr, "Can't set Parameter: %s\n", ErrBuffer );
}

if( ( RetVal =

    ziAPIGetValueD
    ( Conn,
      "DEV1046/demods/0/rate",
      &ValueD ) ) !=
    ZI_INFO_SUCCESS
    )
{

    ziAPIGetError
    ( RetVal, &ErrBuffer, NULL );
    fprintf( stderr, "Can't get Parameter: %s\n", ErrBuffer );
}
else
{
    printf( "Value = %f\n", ValueD );
}
}
```

See Also:

[ziAPISetValueD](#), [ziAPIGetValueAsPollData](#)

ziAPIGetValueI

ZIResult_enum ziAPIGetValueI (**ZIConnection** conn, const char* path, ZIIntegerData* value)

gets the integer-type value of the specified node

This function retrieves the numerical value of the specified node as an integer-type value. The value first found is returned if more than one value is available (a wildcard is used in the path).

Parameters:

[in] conn

Pointer to ZIConnection with which the value should be retrieved

[in] path

Path to the node holding the value

[out] value

Pointer to an 64bit integer in which the value should be written

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_ERROR_LENGTH if the path's length exceeds MAX_PATH_LEN
- ZI_ERROR_OVERFLOW when a FIFO overflow occurred
- ZI_ERROR_COMMAND on an incorrect answer of the server
- ZI_ERROR_SERVER_INTERNAL if an internal error occurred in Data Server
- ZI_ERROR_NOTFOUND if the given path could not be resolved or no value is attached to the node
- ZI_ERROR_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>

#include "
                                ziAPI.h
                                "

void UpdateValue( ZIConnection Conn )
{
    ZIResult_enum RetVal;
    char* ErrBuffer;
    ZIIntegerData ValueI;

    if( ( RetVal =
                                ziAPISetValueD
                                ( Conn,
```

```
        "DEV1046/demods/*/rate",
        5.53 ) ) !=
        ZI_INFO_SUCCESS
    )
{
    ziAPIGetError
    ( RetVal, &ErrBuffer, NULL );
    fprintf( stderr, "Can't set Parameter: %s\n", ErrBuffer );
}

if( ( RetVal =
    ziAPIGetValueI
    ( Conn,
        "DEV1046/demods/0/rate",
        &ValueI ) ) !=
    ZI_INFO_SUCCESS
    )
{
    ziAPIGetError
    ( RetVal, &ErrBuffer, NULL );
    fprintf( stderr, "Can't get Parameter: %s\n", ErrBuffer );
}
else
{
    printf( "Value = %f\n", (float)ValueI );
}
}
```

See Also:

[ziAPISetValueI](#), [ziAPIGetValueAsPollData](#)

ziAPIGetDemodSample

ZIResult_enum ziAPIGetDemodSample (**ZIConnection** conn, const char* path, **ZIDemodSample** * value)

Gets the demodulator sample value of the specified node.

This function retrieves the value of the specified node as an **DemodSample** struct. The value first found is returned if more than one value is available (a wildcard is used in the path). This function is only applicable to paths matching DEMODS/[0-9]+/SAMPLE.

Parameters:

[in] conn

Pointer to ZIConnection with which the value should be retrieved

[in] path

Path to the node holding the value

[out] value

Pointer to a **ZIDemodSample** struct in which the value should be written

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_ERROR_LENGTH if the path's length exceeds MAX_PATH_LEN
- ZI_ERROR_OVERFLOW when a FIFO overflow occurred
- ZI_ERROR_COMMAND on an incorrect answer of the server
- ZI_ERROR_SERVER_INTERNAL if an internal error occurred in Data Server
- ZI_ERROR_NOTFOUND if the given path could not be resolved or no value is attached to the node
- ZI_ERROR_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>

#include "
                                ziAPI.h
                                "

void GetSample( ZIConnection Conn )
{
    ZIResult_enum RetVal;
    char* ErrBuffer;

    ZIDemodSample
```

```
DemodSample;

if( ( RetVal =
    ziAPIGetDemodSample
    ( Conn,
      "DEV1046/demods/0/sample",
      &DemodSample ) ) !=
    ZI_INFO_SUCCESS
    )
{
    ziAPIGetError
    ( RetVal, &ErrBuffer, NULL );
    fprintf( stderr, "Can't get Parameter: %s\n", ErrBuffer );
}
else
{
    printf( "TS = %f, X=%f, Y=%f\n",
            (float)DemodSample.
                timeStamp
            ,
            DemodSample.
                x
            ,
            DemodSample.
                y
            );
}
}
```

See Also:

[ziAPIGetValueAsPollData](#)

ziAPIGetDIOSample

ZIResult_enum ziAPIGetDIOSample (**ZIConnection** conn, const char* path, **ZIDIOSample** * value)

Gets the Digital I/O sample of the specified node.

This function retrieves the newest available DIO sample from the specified node. The value first found is returned if more than one value is available (a wildcard is used in the path). This function is only applicable to nodes ending in "/DIOS/[0-9]+/INPUT".

Parameters:

[in] conn

Pointer to the ZIConnection with which the value should be retrieved

[in] path

Path to the node holding the value

[out] value

Pointer to a **ZIDIOSample** struct in which the value should be written

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_ERROR_LENGTH if the Path's Length exceeds MAX_PATH_LEN or the length of the char-buffer for the nodes given by MaxLen is too small for all elements
- ZI_ERROR_OVERFLOW when a FIFO overflow occurred
- ZI_ERROR_COMMAND on an incorrect answer of the server
- ZI_ERROR_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_ERROR_NOTFOUND if the given path could not be resolved or no value is attached to the node
- ZI_ERROR_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>

#include "
                                ziAPI.h
                                "

void GetSample( ZIConnection Conn )
{
    ZIResult_enum RetVal;
    char* ErrBuffer;
```

```
        zIDIOSample
        DIOSample;

    if( ( RetVal =

        ziAPIGetDIOSample
        ( Conn,
            "DEV1046/dios/0/output",
            &DIOSample ) ) !=
            ZI\_INFO\_SUCCESS
        )
    {

        ziAPIGetError
        ( RetVal, &ErrBuffer, NULL );
        fprintf( stderr, "Can't get Parameter: %s\n", ErrBuffer );
    }
    else
    {
        printf( "TS = %f, bits=%08x\n",
            (float) (DIOSample.
                timeStamp
            ),
            DIOSample.
                bits
            );
    }
}
```

See Also:[ziAPIGetValueAsPollData](#)

ziAPIGetAuxInSample

ZIResult_enum ziAPIGetAuxInSample (**ZIConnection** conn, const char* path, **ZIAuxInSample** * value)

gets the AuxIn sample of the specified node

This function retrieves the newest available AuxIn sample from the specified node. The value first found is returned if more than one value is available (a wildcard is used in the path). This function is only applicable to nodes ending in "/AUXINS/[0-9]+/SAMPLE".

Parameters:

[in] conn

Pointer to the ziConnection with which the Value should be retrieved

[in] path

Path to the Node holding the value

[out] value

Pointer to an **ZIAuxInSample** struct in which the value should be written

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN or the length of the char-buffer for the nodes given by MaxLen is too small for all elements
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no value is attached to the node
- ZI_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>

#include "
                                ziAPI.h
                                "

void GetSample( ZIConnection Conn )
{
    ZIResult_enum RetVal;
    char* ErrBuffer;

    ZIAuxInSample
```

```
        AuxInSample;

    if( ( RetVal =
        ziAPIGetAuxInSample
        ( Conn,
            "DEV1046/auxins/0/sample",
            &AuxInSample ) ) !=
        ZI_INFO_SUCCESS
        )
    {
        ziAPIGetError
        ( RetVal, &ErrBuffer, NULL );
        fprintf( stderr, "Can't get Parameter: %s\n", ErrBuffer );
    }
    else
    {
        printf( "TS = %f, ch0=%f, ch1=%f\n",
            (float)AuxInSample.
                timeStamp
            ,
            AuxInSample.
                ch0
            ,
            AuxInSample.
                ch1
            );
    }
}
```

See Also:

[ziAPIGetValueAsPollData](#)

ziAPIGetValueB

ZIResult_enum ziAPIGetValueB (**ZIConnection** conn, const char* path, unsigned char* buffer, unsigned int* length, unsigned int bufferSize)

gets the Bytearray value of the specified node

This function retrieves the newest available DIO sample from the specified node. The value first found is returned if more than one value is available (a wildcard is used in the path).

Parameters:

[in] conn

Pointer to the ziConnection with which the value should be retrieved

[in] path

Path to the Node holding the value

[out] buffer

Pointer to a buffer to store the retrieved data in

[out] length

Pointer to an unsigned int to store the length of data in. if an error occurred or the length of the passed buffer doesn't reach a zero will be returned

[in] bufferSize

The length of the passed buffer

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN or the length of the char-buffer for the nodes given by MaxLen is too small for all elements
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no value is attached to the node
- ZI_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>
```

```
#include "
                                ziAPI.h
"
```

```
void PrintVersion( ZIConnection Conn )
{

    ZIResult_enum RetVal;
    char* ErrBuffer;

    const char* Path = "ZI/ABOUT/VERSION";
    unsigned char Buffer[ 0xff ];
    unsigned int Length;

    if( ( RetVal =
        ziAPIGetValueB
        ( Conn,
          Path,
          Buffer,
          &Length,
          sizeof(Buffer) - 1 ) ) !=
        ZI_INFO_SUCCESS
      )
    {

        ziAPIGetError
        ( RetVal, &ErrBuffer, NULL );
        fprintf( stderr, "Can't get value: %s\n", ErrBuffer );
    }
    else
    {
        Buffer[ Length ] = 0;
        printf( "%s=\"%s\"\n", Path, Buffer );
    }
}
```

See Also:

[ziAPISetValueB](#), [ziAPIGetValueAsPollData](#)

ziAPISetValueD

ZIResult_enum ziAPISetValueD (**ZIConnection** conn, const char* path, ZIDoubleData value)

asynchronously sets a double-type value to one or more nodes specified in the path

This function sets the values of the nodes specified in path to Value. More than one value can be set if a wildcard is used. The function sets the value asynchronously which means that after the function returns you have no security to which value it is finally set nor at what point in time it is set.

Parameters:

[in] conn

Pointer to the ziConnection for which the value(s) will be set

[in] path

Path to the Node(s) for which the value(s) will be set to Value

[in] value

the double-type value that will be written to the node(s)

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_READONLY on attempt to set a read-only node
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no node given by path is able to hold values
- ZI_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>
```

```
#include "
                ziAPI.h
"
```

```
void UpdateValue( ZIConnection Conn )
{
```

```
    ZIResult_enum RetVal;
    char* ErrBuffer;
    ZIIntegerData ValueI;
```

```
if( ( RetVal =  
    ziAPISetValueD  
    ( Conn,  
      "DEV1046/demods/*/rate",  
      5.53 ) ) !=  
    ZI_INFO_SUCCESS  
    )  
{  
    ziAPIGetError  
    ( RetVal, &ErrBuffer, NULL );  
    fprintf( stderr, "Can't set Parameter: %s\n", ErrBuffer );  
}  
  
if( ( RetVal =  
    ziAPIGetValueI  
    ( Conn,  
      "DEV1046/demods/0/rate",  
      &ValueI ) ) !=  
    ZI_INFO_SUCCESS  
    )  
{  
    ziAPIGetError  
    ( RetVal, &ErrBuffer, NULL );  
    fprintf( stderr, "Can't get Parameter: %s\n", ErrBuffer );  
}  
else  
{  
    printf( "Value = %f\n", (float)ValueI );  
}  
}
```

See Also:

[ziAPIGetValueD](#). [ziAPISyncSetValueD](#)

ziAPISetValueI

ZIResult_enum ziAPISetValueI (**ZIConnection** conn, const char* path, ZIIntegerData value)

asynchronously sets an integer-type value to one or more nodes specified in a path

This function sets the values of the nodes specified in path to Value. More than one value can be set if a wildcard is used. The function sets the value asynchronously which means that after the function returns you have no security to which value it is finally set nor at what point in time it is set.

Parameters:

[in] conn

Pointer to the ziConnection for which the value(s) will be set

[in] path

Path to the Node(s) for which the value(s) will be set

[in] value

the int-type value that will be written to the node(s)

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_READONLY on attempt to set a read-only node
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no node given by path is able to hold values
- ZI_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>

#include "
                                ziAPI.h
                                "

void UpdateValue( ZIConnection Conn )
{
    ZIResult_enum RetVal;
    char* ErrBuffer;
    ZIDoubleData ValueD;
```

```
if( ( RetVal =  
    ziAPISetValueI  
    ( Conn,  
      "DEV1046/demods/*/rate",  
      100 ) ) !=  
    ZI_INFO_SUCCESS  
    )  
{  
    ziAPIGetError  
    ( RetVal, &ErrBuffer, NULL );  
    fprintf( stderr, "Can't set Parameter: %s\n", ErrBuffer );  
}  
  
if( ( RetVal =  
    ziAPIGetValueD  
    ( Conn,  
      "DEV1046/demods/0/rate",  
      &ValueD ) ) !=  
    ZI_INFO_SUCCESS  
    )  
{  
    ziAPIGetError  
    ( RetVal, &ErrBuffer, NULL );  
    fprintf( stderr, "Can't get Parameter: %s\n", ErrBuffer );  
}  
else  
{  
    printf( "Value = %f\n", ValueD );  
}  
}
```

See Also:

[ziAPIGetValueI](#) . [ziAPISyncSetValueI](#)

ziAPISetValueB

ZIResult_enum ziAPISetValueB (**ZIConnection** conn, const char* path, unsigned char* buffer, unsigned int length)

asynchronously sets the binary-type value of one or more nodes specified in the path

This function sets the values at the nodes specified in a path. More than one value can be set if a wildcard is used. The function sets the value asynchronously which means that after the function returns you have no security to which value it is finally set nor at what point in time it is set.

Parameters:

[in] conn

Pointer to the ziConnection for which the value(s) will be set

[in] path

Path to the Node(s) for which the value(s) will be set

[in] buffer

Pointer to the byte array with the data

[in] length

Length of the data in the buffer

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_READONLY on attempt to set a read-only node
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no node given by path is able to hold values
- ZI_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>
```

```
#include "
                ziAPI.h
"
```

```
void ProgramCPU( ZIConnection Conn,
                unsigned char* Buffer,
                int Len )
{
```

```
ZIResult_enum RetVal;
char* ErrBuffer;

if( ( RetVal =
    ziAPISetValueB
    ( Conn,
      "DEV1046/cpus/0/program",
      Buffer,
      Len ) ) !=
    ZI_INFO_SUCCESS
    )
{
    ziAPIGetError
    ( RetVal, &ErrBuffer, NULL );
    fprintf( stderr, "Can't set Parameter: %s\n", ErrBuffer );
}
}
```

See Also:

[ziAPIGetValueB](#). [ziAPISyncSetValueB](#)

ziAPISyncSetValueD

ZIResult_enum ziAPISyncSetValueD (**ZIConnection** conn, const char* path, **ZIDoubleData*** value)

synchronously sets a double-type value to one or more nodes specified in the path

This function sets the values of the nodes specified in path to Value. More than one value can be set if a wildcard is used. The function sets the value synchronously. After returning you know that it is set and to which value it is set.

Parameters:

[in] conn

Pointer to the ziConnection for which the value(s) will be set

[in] path

Path to the Node(s) for which the value(s) will be set to value

[in] value

Pointer to a double-type containing the value to be written. When the function returns value holds the effectively written value.

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_READONLY on attempt to set a read-only node
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no node given by path is able to hold values
- ZI_TIMEOUT when communication timed out

See Also:

[ziAPIGetValueD](#) , [ziAPISetValueD](#)

ziAPISyncSetValueI

ZIResult_enum ziAPISyncSetValueI (**ZIConnection** conn, const char* path, **ZIIntegerData*** value)

synchronously sets an integer-type value to one or more nodes specified in a path

This function sets the values of the nodes specified in path to value. More than one value can be set if a wildcard is used. The function sets the value synchronously. After returning you know that it is set and to which value it is set.

Parameters:

[in] conn

Pointer to the ziConnection for which the value(s) will be set

[in] path

Path to the node(s) for which the value(s) will be set

[in] value

Pointer to a int-type containing then value to be written. when the function returns value holds the effectively written value.

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_READONLY on attempt to set a read-only node
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no node given by path is able to hold values
- ZI_TIMEOUT when communication timed out

See Also:

[ziAPIGetValueI](#) , [ziAPISetValueI](#)

ziAPISyncSetValueB

ZIResult_enum ziAPISyncSetValueB (**ZIConnection** conn, const char* path, uint8_t* buffer, uint32_t* length, uint32_t bufferSize)

Synchronously sets the binary-type value of one ore more nodes specified in the path.

This function sets the values at the nodes specified in a path. More than one value can be set if a wildcard is used. This function sets the value synchronously. After returning you know that it is set and to which value it is set.

Parameters:

[in] conn

Pointer to the ziConnection for which the value(s) will be set

[in] path

Path to the Node(s) for which the value(s) will be set

[in] buffer

Pointer to the byte array with the data

[in] length

Length of the data in the buffer

[in] bufferSize

Length of the data in the buffer

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_READONLY on attempt to set a read-only node
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no node given by path is able to hold values
- ZI_TIMEOUT when communication timed out

See Also:

[ziAPIGetValueB](#) , [ziAPISetValueB](#)

ziAPISync

ZIResult_enum ziAPISync (ZIConnection conn)

Synchronizes the session by dropping all pending data.

This function drops any data that is pending for transfer. Any data (including poll data) retrieved afterwards is guaranteed to be produced not earlier than the call to ziAPISync. This ensures in particular that any settings made prior to the call to ziAPISync have been propagated to the device, and the data retrieved afterwards is produced with the new settings already set to the hardware. Note, however, that this does not include any required settling time.

Parameters:

[in] conn

Pointer to the ZIConnection that is to be synchronized

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_TIMEOUT when communication timed out

ziAPIEchoDevice

ZIResult_enum ziAPIEchoDevice (**ZIConnection** conn, const char* deviceSerial)

Sends an echo command to a device and blocks until answer is received.

This is useful to flush all buffers between API and device to enforce that further code is only executed after the device executed a previous command. Per device echo is only implemented for HF2. For other device types it is a synonym to ziAPISync, and deviceSerial parameter is ignored.

Parameters:

[in] conn

Pointer to the ZIConnection that is to be synchronized

[in] deviceSerial

The serial of the device to get the echo from, e.g., dev2100

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_TIMEOUT when communication timed out

ziAPIAsyncSetDoubleData

ZIResult_enum ziAPIAsyncSetDoubleData (**ZIConnection** conn, const char* path, ZIDoubleData value)

ziAPIAsyncSetIntegerData

[ZIResult_enum](#) ziAPIAsyncSetIntegerData ([ZIConnection](#) conn, const char* path, ZIntegerData value)

ziAPIAsyncSetByteArray

[ZlResult_enum](#) ziAPIAsyncSetByteArray ([ZlConnection](#) conn, const char* path, uint8_t* buffer, uint32_t length)

ziAPIGetValueS

ZIResult_enum ziAPIGetValueS (**ZIConnection** conn, char* path, **DemodSample *** value)

ziAPIGetValueDIO

ZIResult_enum ziAPIGetValueDIO (**ZIConnection** conn, char* path, **DIOSample *** value)

ziAPIGetValueAuxIn

ZIResult_enum ziAPIGetValueAuxIn (**ZIConnection** conn, char* path, **AuxInSample *** value)

6.2.4. Data Streaming

This section describes how to perform data streaming. It allows for recording at high data rates without sample loss.

Data Structures

- struct [ZEvent](#)
This struct holds event data forwarded by the Data Server.
- struct [ziEvent](#)
This struct holds event data forwarded by the Data Server.
Deprecated: See [ZEvent](#) .

Functions

- [ZEvent *](#) [ziAPIAllocateEventEx\(\)](#)
Allocates [ZEvent](#) structure and returns the pointer to it.
Attention!!! It is the client code responsibility to deallocate the structure by calling [ziAPIDeallocateEventEx\(\)](#)
- void [ziAPIDeallocateEventEx\(\)](#) ([ZEvent *](#) ev)
Deallocates [ZEvent](#) structure created with [ziAPIAllocateEventEx\(\)](#) .
- [ZResult_enum](#) [ziAPISubscribe\(\)](#) ([ZConnection](#) conn, const char* path)
subscribes the nodes given by path for [ziAPIPollDataEx](#)
- [ZResult_enum](#) [ziAPIUnSubscribe\(\)](#) ([ZConnection](#) conn, const char* path)
unsubscribes to the nodes given by path
- [ZResult_enum](#) [ziAPIPollDataEx\(\)](#) ([ZConnection](#) conn, [ZEvent *](#) ev, uint32_t timeOutMilliseconds)
checks if an event is available to read
- [ZResult_enum](#) [ziAPIGetValueAsPollData\(\)](#) ([ZConnection](#) conn, const char* path)
triggers a value request, which will be given back on the poll event queue
- [ZResult_enum](#) [ziAPIPollData\(\)](#) ([ZConnection](#) conn, [ZEvent *](#) ev, int timeOut)
Checks if an event is available to read. Deprecated: See [ziAPIPollDataEx\(\)](#) .

Detailed Description

```
#include <stdio.h>
```



```
#include <stdlib.h>

#include "
                                ziAPI.h
                                "

void EventLoop( ZIConnection Conn )
{
    ZIResult_enum RetVal;
    char* ErrBuffer;

                                ZIEvent
                                * Event;

    unsigned int Cnt = 0;

    /*
        allocate ZIEvent in heap memory instead of
        getting it from stack will secure
        against stack overflows especially in windows
    */
    if( ( Event =
                                ziAPIAllocateEventEx
                                ( ) ) == NULL )
    {
        fprintf( stderr, "Can't allocate memory\n" );
        return;
    }

    //subscribe to all nodes
    if( ( RetVal =
                                ziAPISubscribe
                                ( Conn, "*" ) ) !=
                                ZI_INFO_SUCCESS
                                )
    {
        ziAPIGetError
        ( RetVal, &ErrBuffer, NULL );
        fprintf( stderr, "Can't subscribe: %s\n", ErrBuffer );

        ziAPIDeallocateEventEx
        ( Event );

        return;
    }

    //loop 1000 times
    while( Cnt < 1000 )
    {
        //get all demod rates from all devices every 10th cycle
        if( ++Cnt % 10 == 0 )
        {
            if( ( RetVal =
                                ziAPIGetValueAsPollData
                                (
                                    Conn, "*/demods/*/rate" ) ) !=
                                    ZI_INFO_SUCCESS
                                    )
            {

```

```
        ziAPIGetError
        ( RetVal, &ErrBuffer, NULL );
fprintf( stderr, "Can't get value as poll data: %s\n",
        ErrBuffer );

        break;
    }
}

//poll data until no more data is available
while( 1 )
{
    if( ( RetVal =
        ziAPIPollDataEx
        (
            Conn, Event, 0 ) ) !=
        ZI_INFO_SUCCESS
        )
    {
        ziAPIGetError
        ( RetVal, &ErrBuffer, NULL );
fprintf( stderr, "Can't poll data: %s\n", ErrBuffer );

        break;
    }
    else
    {
        //The field Count of the Event struct is zero when
        //no data has been polled
        if( Event->
            valueType
            !=
            ZI_VALUE_TYPE_NONE
            && Event->
            count
            > 0 )
        {
            /*
            process the received event here
            */

        }
        else
        {
            //no more data is available so go on
            break;
        }
    }
}

if(
    ziAPIUnSubscribe
    ( Conn, "*" ) !=
    ZI_INFO_SUCCESS
    )
{
```

```
        ziAPIGetError
        ( RetVal, &ErrBuffer, NULL );
    fprintf( stderr, "Can't unsubscribe: %s\n", ErrBuffer );
}

        ziAPIDeallocateEventEx
        ( Event );

}
```

Data Structure Documentation

struct ZIEvent

This struct holds event data forwarded by the Data Server.

```
#include "ziAPI.h"

typedef struct ZIEvent {

    ZIValueType_enum
        valueType
    ;
    uint32_t
        count
    ;
    uint8_t
        path
    [256];
    void*
        untyped
    ;
    ZIDoubleData*
        doubleData
    ;
    ZIDoubleDataTS
        *doubleDataTS
    ;
    ZIIntegerData*
        integerData
    ;
    ZIIntegerDataTS
        *integerDataTS
    ;
    ZIByteArray
        *byteArray
    ;
    ZIByteArrayTS
        *byteArrayTS
    ;
    ZITreeChangeData
        *treeChangeData
    ;
    TreeChange
        *treeChangeDataOld
    ;
    ZIDemodSample
        *demodSample
    ;
};
```

```
        ZIAuxInSample
    *
    auxInSample
    ;

    ZIDIOSample
    *
    dioSample
    ;

    ZIScopeWave
    *
    scopeWave
    ;

    ScopeWave
    *
    scopeWaveOld
    ;

    ZIPWAWave
    *
    pwaWave
    ;

    union ZIEvent::@1
    {
        value
    }

    uint8_t
    data
    [0x400000];
} ZIEvent;
```

Data Fields

- [ZIValueType_enum](#) valueType
Specifies the type of the data held by the [ZIEvent](#).
- uint32_t count
Number of values available in this event.
- uint8_t path
The path to the node from which the event originates.
- void* untyped
For convenience. The void field doesn't have a corresponding data type.
- ZIDoubleData* doubleData
when valueType == ZI_VALUE_TYPE_DOUBLE_DATA
- [ZIDoubleDataTS](#)* doubleDataTS
when valueType == ZI_VALUE_TYPE_DOUBLE_DATA_TS
- ZIIntegerData* integerData
when valueType == ZI_VALUE_TYPE_INTEGER_DATA
- [ZIIntegerDataTS](#)* integerDataTS
when valueType == ZI_VALUE_TYPE_INTEGER_DATA_TS
- [ZIByteArray](#)* byteArray
when valueType == ZI_VALUE_TYPE_BYTE_ARRAY

- [ZIByteArrayTS](#)* byteArrayTS
when valueType == ZI_VALUE_TYPE_BYTE_ARRAY_TS
- [ZITreeChangeData](#)* treeChangeData
when valueType == ZI_VALUE_TYPE_TREE_CHANGE_DATA
- [TreeChange](#)* treeChangeDataOld
when valueType ==
ZI_VALUE_TYPE_TREE_CHANGE_DATA_OLD
- [ZIDemodSample](#)* demodSample
when valueType == ZI_VALUE_TYPE_DEMOD_SAMPLE
- [ZIAuxInSample](#)* auxInSample
when valueType == ZI_VALUE_TYPE_AUXIN_SAMPLE
- [ZIDIOSample](#)* dioSample
when valueType == ZI_VALUE_TYPE_DIO_SAMPLE
- [ZIScopeWave](#)* scopeWave
when valueType == ZI_VALUE_TYPE_SCOPE_WAVE
- [ScopeWave](#)* scopeWaveOld
when valueType == ZI_VALUE_TYPE_SCOPE_WAVE_OLD
- [ZIPWAWave](#)* pwaWave
when valueType == ZI_VALUE_TYPE_PWA_WAVE
- union ZIEvent::@1 value
Convenience pointer to allow for access to the first entry in
Data using the correct type according to [ZIEvent.valueType](#)
field.
- uint8_t data
The raw value data.

Detailed Description

[ZIEvent](#) is used to give out events like value changes or errors to the user. Event handling functionality is provided by [ziAPISubscribe](#) and [ziAPIUnSubscribe](#) as well as [ziAPIPollDataEx](#).

```
#include <stdio.h>

#include "                                ziAPI.h
                                "

void ProcessEvent(
                                ZIEvent
                                * Event )
{
    unsigned int j;

    switch( Event->
                                valueType
```

```
        )
    {
    case
        ZI_VALUE_TYPE_DOUBLE_DATA
        :

        printf( "%u elements of double data %s\n",
            Event->
                count
            ,
            Event->
                path
            );

        for( j = 0; j < Event->
            count
            ; j++ )
            printf( "%f\n", Event->
                value
                .
                doubleData
                [j]);

        break;

    case
        ZI_VALUE_TYPE_INTEGER_DATA
        :

        printf( "%u elements of integer data %s\n",
            Event->
                count
            ,
            Event->
                path
            );

        for( j = 0; j < Event->
            count
            ; j++ )
            printf( "%f\n", (float)Event->
                value
                .
                integerData
                [j] );

        break;

    case
        ZI_VALUE_TYPE_DEMOD_SAMPLE
        :

        printf( "%u elements of sample data %s\n",
            Event->
                count
            ,
            Event->
                path
            );

        for( j = 0; j < Event->
            count
            ; j++ )
            printf( "TS=%f, X=%f, Y=%f\n",
                (float)Event->
                    value
```

```
        .
        demodSample
    [j].
        timeStamp
    Event->
        '
        value
        .
        demodSample
    [j].
        x
    Event->
        '
        value
        .
        demodSample
    [j].
        y
    );

    break;

case
    ZI_VALUE_TYPE_TREE_CHANGE_DATA
    :

    printf( "%u elements of tree-changed data %s\n",
        Event->
            count
        Event->
            '
            path
        );

    for( j = 0; j < Event->
        count
        ; j++ ){

        switch( Event->
            value
            .
            treeChangeDataOld
        [j].
            Action
        )
        {
        case
            ZI_TREE_ACTION_REMOVE
            :
            printf( "Tree removed: %s\n",
                Event->
                    value
                    .
                    treeChangeDataOld
                [j].
                    Name
            );

            break;

        case
            ZI_TREE_ACTION_ADD
            :
            printf( "treeChangeDataOld added: %s\n",
                Event->
                    value
                    .
                    treeChangeDataOld
```



```
                [j].
                  Name
                );
        break;

    case
        ZI_TREE_ACTION_CHANGE
        :
        printf( "treeChangeDataOld changed: %s\n",
            Event->
                value
            .
                treeChangeDataOld
            [j].
                Name
            );
        break;
    }

}

break;

default:

    printf( "Unexpected event value type: %d\n", Event->
        valueType
    );

    break;

}

}
```

See Also:

[ziAPISubscribe](#) , [ziAPIUnSubscribe](#) , [ziAPIPollDataEx](#)

struct `ziEvent`

This struct holds event data forwarded by the Data Server. Deprecated: See [ZIEvent](#).

```
#include "ziAPI.h"

typedef struct ziEvent {

    ziAPIDataType

    Type

    unsigned int    ;

    Count

    unsigned char   ;

    Path
    [256];
    union ziEvent::Val
    Val

    unsigned char   ;

    Data
    [0x400000];
} ziEvent;
```

Data Structures

- union `ziEvent::Val`

Data Fields

- `ziAPIDataType` Type
- unsigned int Count
- unsigned char Path
- union `ziEvent::Val` Val
- unsigned char Data

Detailed Description

`ziEvent` is used to give out events like value changes or errors to the user. Event handling functionality is provided by `ziAPISubscribe` and `ziAPIUnSubscribe` as well as `ziAPIPollDataEx`.

See Also:[ziAPISubscribe](#), [ziAPIUnSubscribe](#), [ziAPIPollDataEx](#)

```
#include <stdio.h>

#include "                ziAPI.h
                        "

void ProcessEvent(
                        ZIEvent
                        * Event )
{
    unsigned int j;

    switch( Event->
                        valueType
                        )
    {
        case
                        ZI_VALUE_TYPE_DOUBLE_DATA
                        :

            printf( "%u elements of double data %s\n",
                    Event->
                        count
                        ,
                    Event->
                        path
                        );

            for( j = 0; j < Event->
                    count
                    ; j++ )
                printf( "%f\n", Event->
                    value
                    .
                    doubleData
                    [j]);

            break;

        case
                        ZI_VALUE_TYPE_INTEGER_DATA
                        :

            printf( "%u elements of integer data %s\n",
                    Event->
                        count
                        ,
                    Event->
                        path
                        );

            for( j = 0; j < Event->
                    count
                    ; j++ )
                printf( "%f\n", (float)Event->
                    value
                    .
                    integerData
                    [j] );

            break;
```

```
case
    ZI_VALUE_TYPE_DEMOD_SAMPLE
    :
    printf( "%u elements of sample data %s\n",
        Event->
            count
        ,
        Event->
            path
        );
    for( j = 0; j < Event->
        count
        ; j++ )
        printf( "TS=%f, X=%f, Y=%f\n",
            (float)Event->
                value
            ,
            demodSample
            [j].
            timeStamp
        ,
        Event->
            value
        ,
            demodSample
            [j].
            x
        ,
        Event->
            value
        ,
            demodSample
            [j].
            y
        );
    break;
case
    ZI_VALUE_TYPE_TREE_CHANGE_DATA
    :
    printf( "%u elements of tree-changed data %s\n",
        Event->
            count
        ,
        Event->
            path
        );
    for( j = 0; j < Event->
        count
        ; j++ ){
        switch( Event->
            value
            ,
            treeChangeDataOld
            [j].
            Action
        )
        {
            case
                ZI_TREE_ACTION_REMOVE
            :

```

```
        printf( "Tree removed: %s\n",
                Event->
                    value
                    .
                    treeChangeDataOld
                    [j].
                    Name
                );
        break;

    case
        ZI_TREE_ACTION_ADD
        :
        printf( "treeChangeDataOld added: %s\n",
                Event->
                    value
                    .
                    treeChangeDataOld
                    [j].
                    Name
                );
        break;

    case
        ZI_TREE_ACTION_CHANGE
        :
        printf( "treeChangeDataOld changed: %s\n",
                Event->
                    value
                    .
                    treeChangeDataOld
                    [j].
                    Name
                );
        break;
    }

    break;

default:
    printf( "Unexpected event value type: %d\n", Event->
            valueType
        );
    break;
}

}
```

Data Structure Documentation

union ziEvent::Val

```
typedef union ziEvent::Val {  
    void*  
  
    Void  
    ;  
  
    DemodSample  
    *  
    SampleDemod  
    ;  
  
    AuxInSample  
    *  
    SampleAuxIn  
    ;  
  
    DIOSample  
    *  
    SampleDIO  
    ;  
  
    ziDoubleType*  
    Double  
    ;  
  
    ziIntegerType*  
    Integer  
    ;  
  
    TreeChange  
    *  
    Tree  
    ;  
  
    ByteArrayData  
    *  
    ByteArray  
    ;  
  
    ScopeWave  
    *  
    Wave  
    ;  
} ziEvent::Val;
```

Data Fields

- void* Void
- DemodSample* SampleDemod
- AuxInSample* SampleAuxIn
- DIOSample* SampleDIO
- ziDoubleType* Double
- ziIntegerType* Integer

- [TreeChange](#)* Tree
- [ByteArrayData](#)* ByteArray
- [ScopeWave](#)* Wave

Function Documentation

ziAPIAllocateEventEx

ZIEvent * ziAPIAllocateEventEx ()

Allocates [ZIEvent](#) structure and returns the pointer to it. Attention!!! It is the client code responsibility to deallocate the structure by calling [ziAPIDeallocateEventEx](#)!

This function allocates a [ZIEvent](#) structure and returns the pointer to it. Free the memory using [ziAPIDeallocateEventEx](#).

See Also:

[ziAPIDeallocateEventEx](#)

ziAPIDeallocateEventEx

void ziAPIDeallocateEventEx ([ZIEvent](#) * ev)

Deallocates [ZIEvent](#) structure created with [ziAPIAllocateEventEx\(\)](#) .

Parameters:

[in] ev
Pointer to [ZIEvent](#) structure to be deallocated..

See Also:

[ziAPIAllocateEventEx](#)

This function is the compliment to [ziAPIAllocateEventEx\(\)](#)

ziAPISubscribe

ZIResult_enum **ziAPISubscribe** (**ZIConnection** conn, const char* path)

subscribes the nodes given by path for [ziAPIPollDataEx](#)

This function subscribes to nodes so that whenever the value of the node changes the new value can be polled using [ziAPIPollDataEx](#) . By using wildcards or by using a path that is not a leaf node but contains sub nodes, more than one leaf can be subscribed to with one function call.

Parameters:

[in] conn

Pointer to the ziConnection for which to subscribe for

[in] path

Path to the nodes to subscribe

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no node given by path is able to hold values
- ZI_TIMEOUT when communication timed out

See [Data Handling](#) for an example

See Also:

[ziAPIUnSubscribe](#) , [ziAPIPollDataEx](#) , [ziAPIGetValueAsPollData](#)

ziAPIUnSubscribe

ZIResult_enum **ziAPIUnSubscribe** (**ZIConnection** conn, const char* path)

unsubscribes to the nodes given by path

This function is the complement to [ziAPISubscribe](#) . By using wildcards or by using a path that is not a leaf node but contains sub nodes, more than one node can be unsubscribed with one function call.

Parameters:

[in] conn

Pointer to the ziConnection for which to unsubscribe for

[in] path

Path to the Nodes to unsubscribe

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no node given by path is able to hold values
- ZI_TIMEOUT when communication timed out

See [Data Handling](#) for an example

See Also:

[ziAPISubscribe](#) , [ziAPIPollDataEx](#) , [ziAPIGetValueAsPollData](#)

ziAPIPollDataEx

ZIResult_enum ziAPIPollDataEx (**ZIConnection** conn, **ZIEvent** * ev, uint32_t
timeOutMilliseconds)

checks if an event is available to read

This function returns immediately if an event is pending. Otherwise it waits for an event for up to timeOutMilliseconds. All value changes that occur in nodes that have been subscribed to or in children of nodes that have been subscribed to are sent from the Data Server to the ziAPI session. For a description of how the data are available in the struct, refer to the documentation of struct [ziEvent](#). When no event was available within timeOutMilliseconds, the ziEvent::Type field will be ZI_DATA_NONE and the ziEvent::Count field will be zero. Otherwise these fields hold the values corresponding to the event that occurred.

Parameters:

[in] conn

Pointer to the [ZIConnection](#) for which events should be received

[out] ev

Pointer to a [ZIEvent](#) struct in which the received event will be written

[in] timeOutMilliseconds

Time to wait for an event in milliseconds. If -1 it will wait forever, if 0 the function returns immediately.

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_OVERFLOW when a FIFO overflow occurred

See [Data Handling](#) for an example

See Also:

[ziAPISubscribe](#) , [ziAPIUnSubscribe](#) , [ziAPIGetValueAsPollData](#) , [ziEvent](#)

ziAPIGetValueAsPollData

ZIResult_enum **ziAPIGetValueAsPollData** (**ZIConnection** conn, const char* path)

triggers a value request, which will be given back on the poll event queue

Use this function to receive the value of one or more nodes as one or more events using [ziAPIPollDataEx](#), even when the node is not subscribed or no value change has occurred.

Parameters:

[in] conn

Pointer to the [ZIConnection](#) with which the value should be retrieved

[in] path

Path to the Node holding the value

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN or the length of the char-buffer for the nodes given by MaxLen is too small for all elements
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no value is attached to the node
- ZI_TIMEOUT when communication timed out

See [Data Handling](#) for an example

See Also:

[ziAPISubscribe](#) , [ziAPIUnSubscribe](#) , [ziAPIPollDataEx](#)

ziAPIPollData

ZIResult_enum ziAPIPollData (**ZIConnection** conn, **ziEvent** * ev, int timeOut)

Checks if an event is available to read. Deprecated: See [ziAPIPollDataEx\(\)](#) .

Parameters:

[in] conn

Pointer to the [ZIConnection](#) for which events should be received

[out] ev

Pointer to a [ziEvent](#) struct in which the received event will be written

[in] timeOut

Time to wait for an event in milliseconds. If -1 it will wait forever, if 0 the function returns immediately.

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_OVERFLOW when a FIFO overflow occurred

See [Data Handling](#) for an example

See Also:

[ziAPISubscribe](#) , [ziAPIUnSubscribe](#) , [ziAPIGetValueAsPollData](#) , [ziEvent](#)

6.2.5. Error Handling

all functions return [ZIResult_enum](#) . To retrieve a full-text description of such a status one function is provided: [ziAPIGetError](#) .

Functions

- [ZIResult_enum](#) [ziAPIGetError](#) ([ZIResult_enum](#) result, char** buffer, int* base)
returns a description and the severity for a [ZIResult_enum](#)

Function Documentation

ziAPIGetError

ZIResult_enum ziAPIGetError (**ZIResult_enum** result, char** buffer, int* base)

returns a description and the severity for a **ZIResult_enum**

This function returns a static char pointer to a description string for the given **ZIResult_enum** error code. It also provides a parameter returning the severity (info, warning, error). If the given error code does not exist a description for an unknown error and the base for an error will be returned. If a description or the base is not needed NULL may be passed.

Parameters:

[in] result

A **ZIResult_enum** for which the description or base will be returned

[out] buffer

A pointer to a char array to return the description. May be NULL if no description is needed.

[out] base

The severity for the provided Status parameter:

- ZI_INFO_BASE for infos
- ZI_WARNING_BASE for warnings
- ZI_ERROR_BASE for errors

Returns:

- ZI_SUCCESS

6.3. Data Structure Documentation

6.3.1. struct AuxInSample

The [AuxInSample](#) struct holds data for the ZI_DATA_AUXINSAMPLE data type. Deprecated: See [ZIAuxInSample](#).

```
#include "ziAPI.h"

typedef struct AuxInSample {
    ziTimeStampType      TimeStamp
    double                ;
    double                Ch0
    double                ;
    double                Ch1
    ;
} AuxInSample;
```

Data Fields

- ziTimeStampType TimeStamp
- double Ch0
- double Ch1

6.3.2. struct ByteArrayData

The `ByteArrayData` struct holds data for the `ZI_DATA_BYTEARRAY` data type. Deprecated: See `ZIByteArray`.

```
#include "ziAPI.h"

typedef struct ByteArrayData {
    unsigned int
        Len
    ;
    unsigned char
        Bytes
    [0];
} ByteArrayData;
```

Data Fields

- unsigned int Len
- unsigned char Bytes

6.3.3. struct DemodSample

The [DemodSample](#) struct holds data for the ZI_DATA_DEMODSAMPLE data type. Deprecated:
See [ZIDemodSample](#).

```
#include "ziAPI.h"

typedef struct DemodSample {
    ziTimeStampType      TimeStamp
    double                X
    double                Y
    double                Frequency
    double                Phase
    unsigned int          DIOBits
    unsigned int          Reserved
    double                AuxIn0
    double                AuxIn1
} DemodSample;
```

Data Fields

- ziTimeStampType TimeStamp
- double X
- double Y
- double Frequency
- double Phase
- unsigned int DIOBits
- unsigned int Reserved
- double AuxIn0

- double AuxIn1

6.3.4. struct DIOSample

The `DIOSample` struct holds data for the `ZI_DATA_DIOSAMPLE` data type. Deprecated: See `ZIDIOSample`.

```
#include "ziAPI.h"

typedef struct DIOSample {
    ziTimeStampType      TimeStamp
    unsigned int         ;
                        Bits
    unsigned int         ;
                        Reserved
} DIOSample;
```

Data Fields

- `ziTimeStampType TimeStamp`
- `unsigned int Bits`
- `unsigned int Reserved`

6.3.5. struct ScopeWave

The structure used to hold a single scope shot (API Level 1). If the client is connected to the Data Server using API Level 4 (recommended if supported by your device class) please see [ZIScopeWave](#) instead.

```
#include "ziAPI.h"

typedef struct ScopeWave {
    double
        dt
    ;
    unsigned int
        ScopeChannel
    ;
    unsigned int
        TriggerChannel
    ;
    unsigned int
        BWLimit
    ;
    unsigned int
        Count
    ;
    short
        Data
        [0];
} ScopeWave;
```

Data Fields

- double dt
Time difference between samples.
- unsigned int ScopeChannel
Scope channel of the represented data.
- unsigned int TriggerChannel
Trigger channel of the represented data.
- unsigned int BWLimit
Bandwidth-limit flag.
- unsigned int Count
Count of samples.
- short Data
First wave data.

6.3.6. struct TreeChange

The structure used to hold info about added or removed nodes. This is the version without timestamp used in API v1 compatibility mode.

```
#include "ziAPI.h"

typedef struct TreeChange {
    uint32_t          Action

    char              ;
                    Name
                    [32];
} TreeChange;
```

Data Fields

- `uint32_t Action`
field indicating which action occurred on the tree. A value of the `ZITreeAction_enum` (`TREE_ACTION`) enum.
- `char Name`
Name of the Path that has been added, removed or changed.

6.3.7. union ziEvent::Val

```
typedef union ziEvent::Val {  
    void*  
        Void  
        ;  
        DemodSample  
        *  
        SampleDemod  
        ;  
        AuxInSample  
        *  
        SampleAuxIn  
        ;  
        DIOSample  
        *  
        SampleDIO  
        ;  
    ziDoubleType*  
        Double  
        ;  
    ziIntegerType*  
        Integer  
        ;  
        TreeChange  
        *  
        Tree  
        ;  
        ByteArrayData  
        *  
        ByteArray  
        ;  
        ScopeWave  
        *  
        Wave  
        ;  
} ziEvent::Val;
```

Data Fields

- void* Void
- DemodSample* SampleDemod
- AuxInSample* SampleAuxIn
- DIOSample* SampleDIO
- ziDoubleType* Double
- ziIntegerType* Integer

- [TreeChange](#)* Tree
- [ByteArrayData](#)* ByteArray
- [ScopeWave](#)* Wave

6.3.8. struct ZIAuxInSample

The structure used to hold data for a single auxiliary inputs sample.

```
#include "ziAPI.h"

typedef struct ZIAuxInSample {
    ZITimeStamp          timeStamp
    double               ;
    double               ch0
    double               ;
    double               ch1
    ;
} ZIAuxInSample;
```

Data Fields

- ZITimeStamp timeStamp
The timestamp at which the values have been measured.
- double ch0
Channel 0 voltage.
- double ch1
Channel 1 voltage.

6.3.9. struct ZByteArray

The structure used to hold an arbitrary array of bytes. This is the version without timestamp used in API Level 1 compatibility mode.

```
#include "ziAPI.h"

typedef struct ZByteArray {
    uint32_t          length

    ;

    uint8_t          bytes
                    [0];
} ZByteArray;
```

Data Fields

- `uint32_t length`
Length of the data readable from the Bytes field.
- `uint8_t bytes`
The data itself. The array has the size given in length.

6.3.10. struct ZIByteArrayTS

The structure used to hold an arbitrary array of bytes. This is the same as [ZIByteArray](#) , but with timestamp.

```
#include "ziAPI.h"

typedef struct ZIByteArrayTS {
    ZITimeStamp
        timeStamp
    ;
    uint32_t
        length
    ;
    uint8_t
        bytes
    [0];
} ZIByteArrayTS;
```

Data Fields

- ZITimeStamp timeStamp
Time stamp at which the data was updated.
- uint32_t length
length of the data readable from the bytes field
- uint8_t bytes
the data itself. The array has the size given in length

6.3.11. struct ZIDemodSample

The structure used to hold data for a single demodulator sample.

```
#include "ziAPI.h"

typedef struct ZIDemodSample {
    ZITimeStamp      timeStamp

    double           ;
                    x
    double           ;
                    y
    double           ;
                    frequency
    double           ;
                    phase
    uint32_t         ;
                    dioBits
    uint32_t         ;
                    trigger
    double           ;
                    auxIn0
    double           ;
                    auxIn1
} ZIDemodSample;
```

Data Fields

- ZITimeStamp timeStamp
The timestamp at which the sample has been measured.
- double x
X part of the sample.
- double y
Y part of the sample.
- double frequency
Frequency at that sample.
- double phase
Phase at that sample.
- uint32_t dioBits
the current bits of the DIO.
- uint32_t trigger
trigger bits
- double auxIn0

- value of Aux input 0.
- double auxIn1
value of Aux input 1.

6.3.12. struct ZIDIOSample

The structure used to hold data for a single digital I/O sample.

```
#include "ziAPI.h"

typedef struct ZIDIOSample {
    ZITimeStamp          timeStamp
    uint32_t             bits
    uint32_t             reserved
} ZIDIOSample;
```

Data Fields

- ZITimeStamp timeStamp
The timestamp at which the values have been measured.
- uint32_t bits
The digital I/O values.
- uint32_t reserved
Filler to keep 8 bytes alignment in the array of [ZIDIOSample](#) structures.

6.3.13. struct ZIDoubleDataTS

The structure used to hold a single IEEE double value. Same as ZIDoubleData, but with timestamp.

```
#include "ziAPI.h"

typedef struct ZIDoubleDataTS {
    ZITimeStamp
        timeStamp
    ;
    ZIDoubleData
        value
    ;
} ZIDoubleDataTS;
```

Data Fields

- ZITimeStamp timeStamp
Time stamp at which the value has changed.
- ZIDoubleData value

6.3.14. struct ziEvent

This struct holds event data forwarded by the Data Server. Deprecated: See [ZIEvent](#).

```
#include "ziAPI.h"

typedef struct ziEvent {

    ziAPIDataType
        Type
    ;
    unsigned int
        Count
    ;
    unsigned char
        Path
        [256];
    union ziEvent::Val
        Val
    ;
    unsigned char
        Data
        [0x400000];
} ziEvent;
```

Data Structures

- union [ziEvent::Val](#)

Data Fields

- [ziAPIDataType](#) Type
- unsigned int Count
- unsigned char Path
- union [ziEvent::Val](#) Val
- unsigned char Data

Detailed Description

[ziEvent](#) is used to give out events like value changes or errors to the user. Event handling functionality is provided by [ziAPISubscribe](#) and [ziAPIUnSubscribe](#) as well as [ziAPIPollDataEx](#).

See Also:

[ziAPISubscribe](#), [ziAPIUnSubscribe](#), [ziAPIPollDataEx](#)

```
#include <stdio.h>

#include "
                                ziAPI.h
                                "

void ProcessEvent(
                                ziEvent
                                * Event )
{
    unsigned int j;
    switch( Event->
                                valueType
                                )
    {
        case
                                ZI_VALUE_TYPE_DOUBLE_DATA
                                :
            printf( "%u elements of double data %s\n",
                    Event->
                                count
                                ,
                    Event->
                                path
                                );
            for( j = 0; j < Event->
                                count
                                ; j++ )
                printf( "%f\n", Event->
                                value
                                .
                                doubleData
                                [j]);
            break;
        case
                                ZI_VALUE_TYPE_INTEGER_DATA
                                :
            printf( "%u elements of integer data %s\n",
                    Event->
                                count
                                ,
                    Event->
                                path
                                );
            for( j = 0; j < Event->
                                count
                                ; j++ )
                printf( "%f\n", (float)Event->
                                value
                                .
                                integerData
                                );
    }
```

```
        [j] );

    break;

case
    ZI_VALUE_TYPE_DEMOD_SAMPLE
    :

    printf( "%u elements of sample data %s\n",
        Event->
            count
        ,
        Event->
            path
        );

    for( j = 0; j < Event->
        count
        ; j++ )
        printf( "TS=%f, X=%f, Y=%f\n",
            (float)Event->
                value
            ,
            demodSample
            [j].
            timeStamp
        ,
        Event->
            value
        ,
            demodSample
            [j].
            x
        ,
        Event->
            value
        ,
            demodSample
            [j].
            y
        );

    break;

case
    ZI_VALUE_TYPE_TREE_CHANGE_DATA
    :

    printf( "%u elements of tree-changed data %s\n",
        Event->
            count
        ,
        Event->
            path
        );

    for( j = 0; j < Event->
        count
        ; j++ ){

        switch( Event->
            value
            ,
            treeChangeDataOld
            [j].
            Action
        )
```

```
{
    case
        ZI_TREE_ACTION_REMOVE
        :
        printf( "Tree removed: %s\n",
            Event->
                value
                .
                treeChangeDataOld
            [j].
                Name
            );
        break;

    case
        ZI_TREE_ACTION_ADD
        :
        printf( "treeChangeDataOld added: %s\n",
            Event->
                value
                .
                treeChangeDataOld
            [j].
                Name
            );
        break;

    case
        ZI_TREE_ACTION_CHANGE
        :
        printf( "treeChangeDataOld changed: %s\n",
            Event->
                value
                .
                treeChangeDataOld
            [j].
                Name
            );
        break;
    }

    break;

default:

    printf( "Unexpected event value type: %d\n", Event->
        valueType
    );

    break;
}

}
```

Data Structure Documentation

union ziEvent::Val

```
typedef union ziEvent::Val {  
    void*  
        ; Void  
        ;  
        DemodSample  
        * SampleDemod  
        ;  
        AuxInSample  
        * SampleAuxIn  
        ;  
        DIOSample  
        * SampleDIO  
        ;  
    ziDoubleType* Double  
        ;  
    ziIntegerType* Integer  
        ;  
        TreeChange  
        * Tree  
        ;  
        ByteArrayData  
        * ByteArray  
        ;  
        ScopeWave  
        * Wave  
        ;  
} ziEvent::Val;
```

Data Fields

- void* Void
- DemodSample* SampleDemod
- AuxInSample* SampleAuxIn
- DIOSample* SampleDIO
- ziDoubleType* Double

- `ziIntegerType*` Integer
- `TreeChange*` Tree
- `ByteArrayData*` ByteArray
- `ScopeWave*` Wave

6.3.15. struct ZIEvent

This struct holds event data forwarded by the Data Server.

```
#include "ziAPI.h"

typedef struct ZIEvent {

    ZIValueType_enum
        valueType
    ;
    uint32_t
        count
    ;
    uint8_t
        path
        [256];
    void*
        untyped
    ;
    ZIDoubleData*
        doubleData
    ;
    ZIDoubleDataTS
        *
        doubleDataTS
    ;
    ZIIntegerData*
        integerData
    ;
    ZIIntegerDataTS
        *
        integerDataTS
    ;
    ZIByteArray
        *
        byteArray
    ;
    ZIByteArrayTS
        *
        byteArrayTS
    ;
    ZITreeChangeData
        *
        treeChangeData
    ;
    TreeChange
        *
        treeChangeDataOld
    ;
    ZIDemodSample
        *
        demodSample
    ;
    ZIAuxInSample
        *
```

```
        auxInSample
    ;

    ZIDIOSample
    *
    dioSample
    ;

    ZIScopeWave
    *
    scopeWave
    ;

    ScopeWave
    *
    scopeWaveOld
    ;

    ZIPWAWave
    *
    pwaWave
    ;
union ZIEvent::@1
    value
    ;
    uint8_t
        data
        [0x400000];
} ZIEvent;
```

Data Fields

- [ZIValueType_enum](#) valueType
Specifies the type of the data held by the [ZIEvent](#).
- `uint32_t` count
Number of values available in this event.
- `uint8_t` path
The path to the node from which the event originates.
- `void*` untyped
For convenience. The void field doesn't have a corresponding data type.
- `ZIDoubleData*` doubleData
when valueType == `ZI_VALUE_TYPE_DOUBLE_DATA`
- [ZIDoubleDataTS*](#) doubleDataTS
when valueType == `ZI_VALUE_TYPE_DOUBLE_DATA_TS`
- `ZIIntegerData*` integerData
when valueType == `ZI_VALUE_TYPE_INTEGER_DATA`
- [ZIIntegerDataTS*](#) integerDataTS
when valueType == `ZI_VALUE_TYPE_INTEGER_DATA_TS`
- [ZIByteArray*](#) byteArray
when valueType == `ZI_VALUE_TYPE_BYTE_ARRAY`
- [ZIByteArrayTS*](#) byteArrayTS

- when valueType == ZI_VALUE_TYPE_BYTE_ARRAY_TS
- [ZITreeChangeData](#)* treeChangeData
when valueType == ZI_VALUE_TYPE_TREE_CHANGE_DATA
- [TreeChange](#)* treeChangeDataOld
when valueType ==
ZI_VALUE_TYPE_TREE_CHANGE_DATA_OLD
- [ZIDemodSample](#)* demodSample
when valueType == ZI_VALUE_TYPE_DEMOD_SAMPLE
- [ZIAuxInSample](#)* auxInSample
when valueType == ZI_VALUE_TYPE_AUXIN_SAMPLE
- [ZIDIOSample](#)* dioSample
when valueType == ZI_VALUE_TYPE_DIO_SAMPLE
- [ZIScopeWave](#)* scopeWave
when valueType == ZI_VALUE_TYPE_SCOPE_WAVE
- [ScopeWave](#)* scopeWaveOld
when valueType == ZI_VALUE_TYPE_SCOPE_WAVE_OLD
- [ZIPWAWave](#)* pwaWave
when valueType == ZI_VALUE_TYPE_PWA_WAVE
- union ZIEvent::@1 value
Convenience pointer to allow for access to the first entry in
Data using the correct type according to [ZIEvent.valueType](#)
field.
- uint8_t data
The raw value data.

Detailed Description

[ZIEvent](#) is used to give out events like value changes or errors to the user. Event handling functionality is provided by [ziAPISubscribe](#) and [ziAPIUnSubscribe](#) as well as [ziAPIPollDataEx](#).

```
#include <stdio.h>

#include "                ziAPI.h
                        "

void ProcessEvent(
                        ZIEvent
                        * Event )
{
    unsigned int j;

    switch( Event->
                        valueType
                    )
```

```
{  
  
case  
    ZI_VALUE_TYPE_DOUBLE_DATA  
    :  
  
    printf( "%u elements of double data %s\n",  
        Event->  
            count  
        ,  
        Event->  
            path  
        );  
  
    for( j = 0; j < Event->  
        count  
        ; j++ )  
        printf( "%f\n", Event->  
            value  
            .  
            doubleData  
            [j]);  
  
    break;  
  
case  
    ZI_VALUE_TYPE_INTEGER_DATA  
    :  
  
    printf( "%u elements of integer data %s\n",  
        Event->  
            count  
        ,  
        Event->  
            path  
        );  
  
    for( j = 0; j < Event->  
        count  
        ; j++ )  
        printf( "%f\n", (float)Event->  
            value  
            .  
            integerData  
            [j] );  
  
    break;  
  
case  
    ZI_VALUE_TYPE_DEMOD_SAMPLE  
    :  
  
    printf( "%u elements of sample data %s\n",  
        Event->  
            count  
        ,  
        Event->  
            path  
        );  
  
    for( j = 0; j < Event->  
        count  
        ; j++ )  
        printf( "TS=%f, X=%f, Y=%f\n",  
            (float)Event->  
                value  
                .  
                demodSampleData  
                [j] );  
  
    break;  
}
```

```
                demodSample
            [j].
            timeStamp
        Event->
            value
            .
            demodSample
            [j].
            x
        Event->
            value
            .
            demodSample
            [j].
            y
        );

    break;

case
    ZI_VALUE_TYPE_TREE_CHANGE_DATA
    :

    printf( "%u elements of tree-changed data %s\n",
        Event->
            count
        ,
        Event->
            path
        );

    for( j = 0; j < Event->
        count
        ; j++ ){

        switch( Event->
            value
            .
            treeChangeDataOld
            [j].
            Action
            )
        {
        case
            ZI_TREE_ACTION_REMOVE
            :
            printf( "Tree removed: %s\n",
                Event->
                    value
                    .
                    treeChangeDataOld
                    [j].
                    Name
                );
            break;

        case
            ZI_TREE_ACTION_ADD
            :
            printf( "treeChangeDataOld added: %s\n",
                Event->
                    value
                    .
                    treeChangeDataOld
                    [j].
```

```
                Name
            );
        break;

    case
        ZI_TREE_ACTION_CHANGE
        :
        printf( "treeChangeDataOld changed: %s\n",
            Event->
                value
            .
                treeChangeDataOld
            [j].
                Name
            );
        break;
    }

}

break;

default:

    printf( "Unexpected event value type: %d\n", Event->
        valueType
    );

    break;

}

}
```

See Also:

[ziAPISubscribe](#) , [ziAPIUnSubscribe](#) , [ziAPIPollDataEx](#)

6.3.16. struct ZIntegerDataTS

The structure used to hold a single 64bit signed integer value. Same as ZIntegerData, but with timestamp.

```
#include "ziAPI.h"

typedef struct ZIntegerDataTS {
    ZIntegerData
    ZITimeStamp    timeStamp
};
    value
} ZIntegerDataTS;
```

Data Fields

- ZITimeStamp timeStamp
Time stamp at which the value has changed.
- ZIntegerData value

6.3.17. struct ZIPWASample

Single PWA sample value.

```
#include "ziAPI.h"

typedef struct ZIPWASample {
    double
        binPhase
    ;
    double
        x
    ;
    double
        y
    ;
    uint32_t
        countBin
    ;
    uint32_t
        reserved
    ;
} ZIPWASample;
```

Data Fields

- double binPhase
Phase position of each bin.
- double x
Real PWA result or X component of a demod PWA.
- double y
Y component of the demod PWA.
- uint32_t countBin
Number of events per bin.
- uint32_t reserved
Reserved.

6.3.18. struct ZIPWAWave

PWA Wave.

```
#include "ziAPI.h"

typedef struct ZIPWAWave {
    ZITimeStamp
        timeStamp
        ;
    uint64_t
        sampleCount
        ;
    uint32_t
        inputSelect
        ;
    uint32_t
        oscSelect
        ;
    uint32_t
        harmonic
        ;
    uint32_t
        binCount
        ;
    double
        frequency
        ;
    uint8_t
        pwaType
        ;
    uint8_t
        mode
        ;
    uint8_t
        overflow
        ;
    uint8_t
        commensurable
        ;
    uint32_t
        reservedUInt
        ;
    ZIPWASample
        data
        [0];
} ZIPWAWave;
```

Data Fields

- ZITimeStamp timeStamp
Time stamp at which the data was updated.
- uint64_t sampleCount
Total sample count considered for PWA.
- uint32_t inputSelect
Input selection used for the PWA.
- uint32_t oscSelect

Oscillator used for the PWA.

- `uint32_t harmonic`
Harmonic setting.
- `uint32_t binCount`
Bin count of the PWA.
- `double frequency`
Frequency during PWA accumulation.
- `uint8_t pwaType`
Type of the PWA.
- `uint8_t mode`
PWA Mode [0: zoom PWA, 1: harmonic PWA].
- `uint8_t overflow`
Overflow indicators. `overflow[0]`: Data accumulator overflow, `overflow[1]`: Counter at limit, `overflow[6..2]`: Reserved, `overflow[7]`: Invalid (missing frames).
- `uint8_t commensurable`
Commensurability of the data.
- `uint32_t reservedUInt`
Reserved unsigned int.
- [ZIPWASample](#) data
PWA data vector.

6.3.19. struct ZIScopeWave

The structure used to hold scope data. The data may be formatted differently, depending on settings. See the description of the structure members for details.

```
#include "ziAPI.h"

typedef struct ZIScopeWave {
    ZITimeStamp
        timeStamp
    ;
    ZITimeStamp
        triggerTimeStamp
    ;
    double
        dt
    ;
    uint8_t
        channelEnable
    [4];
    uint8_t
        channelInput
    [4];
    uint8_t
        triggerEnable
    ;
    uint8_t
        triggerInput
    ;
    uint8_t
        reserved0
    [2];
    uint8_t
        channelBWLimit
    [4];
    uint8_t
        channelMath
    [4];
    float
        channelScaling
    [4];
    uint32_t
        sequenceNumber
    ;
    uint32_t
        segmentNumber
    ;
    uint32_t
        blockNumber
    ;
    uint64_t
        totalSamples
    ;
    uint8_t
        dataTransferMode
    ;
    uint8_t
        blockMarker
    ;
    uint8_t
        flags
    ;
    uint8_t
        sampleFormat
    ;
    uint32_t
```

```
        sampleCount
    ;
    int16_t
        dataInt16
    [0];
    int32_t
        dataInt32
    [0];
    float
        dataFloat
    [0];
    union ZIScopeWave: {0
        data
    ;
} ZIScopeWave;
```

Data Fields

- ZITimeStamp timeStamp
Time stamp of the first sample in this data block.
- ZITimeStamp triggerTimeStamp
Time stamp of the trigger (may also fall between samples)
- double dt
Time difference between samples in seconds.
- uint8_t channelEnable
Up to four channels: if channel is enabled, corresponding element is non-zero.
- uint8_t channelInput
Specifies the input source for each of the scope four channels: 0 = Signal Input 1, 1 = Signal Input 2, 2 = Trigger Input 1, 3 = Trigger Input 2, 4 = Aux Output 1, 5 = Aux Output 2, 6 = Aux Output 3, 7 = Aux Output 4, 8 = Aux Input 1, 9 = Aux Input 2.
- uint8_t triggerEnable
Non-zero if trigger is enabled: Bit(0): rising edge trigger off = 0, on = 1. Bit(1): falling edge trigger off = 0, on = 1.
- uint8_t triggerInput
Trigger source (same values as for channel input)
- uint8_t reserved0
- uint8_t channelBWLimit
Bandwidth-limit flag, per channel. Bit(0): off = 0, on = 1
Bit(7...1): Reserved.
- uint8_t channelMath
Math Operation (e.g averaging) Bit (7..0): Reserved.
- float channelScaling
Data scaling factors for up to 4 channels.

- `uint32_t` `sequenceNumber`
Current scope shot sequence number. Identifies a scope shot.
- `uint32_t` `segmentNumber`
Current segment number.
- `uint32_t` `blockNumber`
Current block number from the beginning of a scope shot.
Large scope shots are split into blocks, which need to be concatenated to obtain the complete scope shot.
- `uint64_t` `totalSamples`
Total number of samples in one channel in the current scope shot, same for all channels.
- `uint8_t` `dataTransferMode`
Data transfer mode `SingleTransfer` = 0, `BlockTransfer` = 1, `ContinuousTransfer` = 3, `FFTSingleTransfer` = 4. Other values are reserved.
- `uint8_t` `blockMarker`
Block marker: Bit (0): 1 = End marker for continuous or multi-block transfer Bit (7..0): Reserved.
- `uint8_t` `flags`
Indicator Flags. Bit (0): 1 = Data loss detected (samples are 0), Bit (1): 1 = Missed trigger, Bit (2): 1 = Transfer failure (corrupted data).
- `uint8_t` `sampleFormat`
Data format of samples: `Int16` = 0, `Int32` = 1, `Float` = 2, `Int16Interleaved` = 4, `Int32Interleaved` = 5, `FloatInterleaved` = 6.
- `uint32_t` `sampleCount`
Number of samples in one channel in the current block, same for all channels.
- `int16_t` `dataInt16`
- `int32_t` `dataInt32`
- `float` `dataFloat`
- `union ZIScopeWave::@0` `data`
First wave data.

6.3.20. struct ZITreeChangeData

The struct is holding info about added or removed nodes.

```
#include "ziAPI.h"

typedef struct ZITreeChangeData {
    ZITimeStamp          timeStamp
                        ;
    uint32_t             action
                        ;
    char                 name
                        [32];
} ZITreeChangeData;
```

Data Fields

- ZITimeStamp timeStamp
Time stamp at which the data was updated.
- uint32_t action
field indicating which action occurred on the tree. A value of the ZITreeAction_enum.
- char name
Name of the Path that has been added, removed or changed.

6.4. File Documentation

6.4.1. File ziAPI.h

Header File for the Zurich Instruments AG C/C++ API v4 doing Communication with Data Server.

Data Structures

- struct [ZIDoubleDataTS](#)
The structure used to hold a single IEEE double value. Same as ZIDoubleData, but with timestamp.
- struct [ZIIntegerDataTS](#)
The structure used to hold a single 64bit signed integer value. Same as ZIIntegerData, but with timestamp.
- struct [ZITreeChangeData](#)
The struct is holding info about added or removed nodes.
- struct [TreeChange](#)
The structure used to hold info about added or removed nodes. This is the version without timestamp used in API v1 compatibility mode.
- struct [ZIDemodSample](#)
The structure used to hold data for a single demodulator sample.
- struct [ZIAuxInSample](#)
The structure used to hold data for a single auxiliary inputs sample.
- struct [ZIDIOSample](#)
The structure used to hold data for a single digital I/O sample.
- struct [ZIByteArray](#)
The structure used to hold an arbitrary array of bytes. This is the version without timestamp used in API Level 1 compatibility mode.
- struct [ZIByteArrayTS](#)
The structure used to hold an arbitrary array of bytes. This is the same as [ZIByteArray](#) , but with timestamp.
- struct [ScopeWave](#)
The structure used to hold a single scope shot (API Level 1). If the client is connected to the Data Server using API Level 4 (recommended if supported by your device class) please see [ZIScopeWave](#) instead.
- struct [ZIScopeWave](#)
The structure used to hold scope data. The data may be formatted differently, depending on settings. See the description of the structure members for details.

- struct [ZIPWASample](#)
Single PWA sample value.
- struct [ZIPWAWave](#)
PWA Wave.
- struct [ZIEvent](#)
This struct holds event data forwarded by the Data Server.
- struct [DemodSample](#)
The [DemodSample](#) struct holds data for the ZI_DATA_DEMODSAMPLE data type. Deprecated: See [ZIDemodSample](#) .
- struct [AuxInSample](#)
The [AuxInSample](#) struct holds data for the ZI_DATA_AUXINSAMPLE data type. Deprecated: See [ZIAuxInSample](#) .
- struct [DIOSample](#)
The [DIOSample](#) struct holds data for the ZI_DATA_DIOSAMPLE data type. Deprecated: See [ZIDIOSample](#) .
- struct [ByteArrayData](#)
The [ByteArrayData](#) struct holds data for the ZI_DATA_BYTEARRAY data type. Deprecated: See [ZIByteArray](#) .
- struct [ziEvent](#)
This struct holds event data forwarded by the Data Server. Deprecated: See [ZIEvent](#) .
- union [ziEvent::Val](#)

Defines

- #define MAX_PATH_LEN 256
The maximum length that has to be used for passing paths to functions (including terminating zero)
- #define MAX_EVENT_SIZE 0x400000
The maximum size of an event's data block.

Typedefs

- typedef ZIConnection

The `ZIConnection` is a connection reference; it holds information and helper variables about a connection to the Data Server. There is nothing in this reference which the user may use, so it is hidden and instead a dummy pointer is used. See [ziAPIInit](#) for how to create a `ZIConnection`.

Enumerations

- enum `ZIResult_enum` { `ZI_INFO_BASE`, `ZI_INFO_SUCCESS`, `ZI_INFO_MAX`, `ZI_WARNING_BASE`, `ZI_WARNING_GENERAL`, `ZI_WARNING_UNDERRUN`, `ZI_WARNING_OVERFLOW`, `ZI_WARNING_NOTFOUND`, `ZI_WARNING_MAX`, `ZI_ERROR_BASE`, `ZI_ERROR_GENERAL`, `ZI_ERROR_USB`, `ZI_ERROR_MALLOC`, `ZI_ERROR_MUTEX_INIT`, `ZI_ERROR_MUTEX_DESTROY`, `ZI_ERROR_MUTEX_LOCK`, `ZI_ERROR_MUTEX_UNLOCK`, `ZI_ERROR_THREAD_START`, `ZI_ERROR_THREAD_JOIN`, `ZI_ERROR_SOCKET_INIT`, `ZI_ERROR_SOCKET_CONNECT`, `ZI_ERROR_HOSTNAME`, `ZI_ERROR_CONNECTION`, `ZI_ERROR_TIMEOUT`, `ZI_ERROR_COMMAND`, `ZI_ERROR_SERVER_INTERNAL`, `ZI_ERROR_LENGTH`, `ZI_ERROR_FILE`, `ZI_ERROR_DUPLICATE`, `ZI_ERROR_READONLY`, `ZI_ERROR_DEVICE_NOT_VISIBLE`, `ZI_ERROR_DEVICE_IN_USE`, `ZI_ERROR_DEVICE_INTERFACE`, `ZI_ERROR_DEVICE_CONNECTION_TIMEOUT`, `ZI_ERROR_DEVICE_DIFFERENT_INTERFACE`, `ZI_ERROR_DEVICE_NEEDS_FW_UPGRADE`, `ZI_ERROR_ZIEVENT_DATATYPE_MISMATCH`, `ZI_ERROR_MAX`, `ZI_SUCCESS`, `ZI_MAX_INFO`, `ZI_WARNING`, `ZI_UNDERRUN`, `ZI_OVERFLOW`, `ZI_NOTFOUND`, `ZI_MAX_WARNING`, `ZI_ERROR`, `ZI_USB`, `ZI_MALLOC`, `ZI_MUTEX_INIT`, `ZI_MUTEX_DESTROY`, `ZI_MUTEX_LOCK`, `ZI_MUTEX_UNLOCK`, `ZI_THREAD_START`, `ZI_THREAD_JOIN`, `ZI_SOCKET_INIT`, `ZI_SOCKET_CONNECT`, `ZI_HOSTNAME`, `ZI_CONNECTION`, `ZI_TIMEOUT`, `ZI_COMMAND`, `ZI_SERVER_INTERNAL`, `ZI_LENGTH`, `ZI_FILE`, `ZI_DUPLICATE`, `ZI_READONLY`, `ZI_MAX_ERROR` }

Defines return value for all `ziAPI` functions. Divided into 3 regions: info, warning and error.

- enum `ZIValueType_enum` { `ZI_VALUE_TYPE_NONE`, `ZI_VALUE_TYPE_DOUBLE_DATA`, `ZI_VALUE_TYPE_DOUBLE_DATA_TS`, `ZI_VALUE_TYPE_INTEGER_DATA`, `ZI_VALUE_TYPE_INTEGER_DATA_TS`, `ZI_VALUE_TYPE_DEMOD_SAMPLE`, `ZI_VALUE_TYPE_AUXIN_SAMPLE`, `ZI_VALUE_TYPE_DIO_SAMPLE`, `ZI_VALUE_TYPE_BYTE_ARRAY`, `ZI_VALUE_TYPE_BYTE_ARRAY_TS`, `ZI_VALUE_TYPE_TREE_CHANGE_DATA`, `ZI_VALUE_TYPE_TREE_CHANGE_DATA_OLD`, `ZI_VALUE_TYPE_SCOPE_WAVE`, `ZI_VALUE_TYPE_SCOPE_WAVE_OLD`, }

ZI_VALUE_TYPE_PWA_WAVE, ZI_DATA_NONE,
ZI_DATA_DOUBLE, ZI_DATA_INTEGER,
ZI_DATA_DEMODSAMPLE, ZI_DATA_SCOPEWAVE,
ZI_DATA_AUXINSAMPLE, ZI_DATA_DIOSAMPLE,
ZI_DATA_BYTEARRAY, ZI_DATA_TREE_CHANGED }

Enumerates all types that data in a [ZIEvent](#) may have.

- enum [ZITreeAction_enum](#) { ZI_TREE_ACTION_REMOVE,
ZI_TREE_ACTION_ADD, ZI_TREE_ACTION_CHANGE }

Defines the actions that are performed on a tree,
as returned in the [ZITreeChangeData::action](#) or
[ZITreeChangeDataOld::action](#).

- enum [ZIAPIVersion_enum](#) { ZI_API_VERSION_1,
ZI_API_VERSION_4 }

- enum [ZIListNodes_enum](#) { ZI_LIST_NODES_NONE,
ZI_LIST_NODES_RECURSIVE, ZI_LIST_NODES_ABSOLUTE,
ZI_LIST_NODES_LEAFONLY,
ZI_LIST_NODES_SETTINGSONLY, ZI_LIST_NONE,
ZI_LIST_RECURSIVE, ZI_LIST_ABSOLUTE,
ZI_LIST_LEAFONLY, ZI_LIST_SETTINGSONLY }

Defines the values of the flags used in [ziAPIListNodes](#) .

- enum [TREE_ACTION](#) { TREE_ACTION_REMOVE,
TREE_ACTION_ADD, TREE_ACTION_CHANGE }

TREE_ACTION defines the values for the [TreeChange::Action](#)
Variable.

Functions

- enum [ZIResult_enum](#) [ziAPIInit](#) ([ZIConnection](#) * conn)
Initializes a [ZIConnection](#) structure.
- enum [ZIResult_enum](#) [ziAPIDestroy](#) ([ZIConnection](#) conn)
Destroys a [ZIConnection](#) structure.
- enum [ZIResult_enum](#) [ziAPIConnect](#) ([ZIConnection](#) conn, const
char* hostname, uint16_t port)
Connects the [ZIConnection](#) to Data Server.
- enum [ZIResult_enum](#) [ziAPIDisconnect](#) ([ZIConnection](#) conn)
Disconnects an established connection.
- enum [ZIResult_enum](#) [ziAPIListImplementations](#) (char*
implementations, uint32_t bufferSize)
Returns the list of supported implementations.
- enum [ZIResult_enum](#) [ziAPIConnectEx](#) ([ZIConnection](#) conn, const
char* hostname, uint16_t port, [ZIAPIVersion_enum](#) apiLevel,
const char* implementation)
Connects to Data Server and enables extended [ziAPI](#).

- `ZIResult_enum ziAPIGetConnectionAPILevel (ZIConnection conn, ZIAPIVersion_enum* apiLevel)`
Returns ziAPI level used for the connection conn.
- `ZIResult_enum ziAPIGetRevision (unsigned int* revision)`
Retrieves the revision of ziAPI.
- `ZIResult_enum ziAPIListNodes (ZIConnection conn, const char* path, char* nodes, int bufferSize, int flags)`
Returns all child nodes found at the specified path.
- `ZIResult_enum ziAPIUpdateDevices (ZIConnection conn)`
Search for the newly connected devices and update the tree.
- `ZIResult_enum ziAPIConnectDevice (ZIConnection conn, const char* deviceSerial, const char* deviceInterface, const char* interfaceParams)`
Connect a device to the server.
- `ZIResult_enum ziAPIDisconnectDevice (ZIConnection conn, const char* deviceSerial)`
Disconnect a device from the server.
- `ZIResult_enum ziAPIGetValueD (ZIConnection conn, const char* path, ZIDoubleData* value)`
gets the double-type value of the specified node
- `ZIResult_enum ziAPIGetValueI (ZIConnection conn, const char* path, ZIIntegerData* value)`
gets the integer-type value of the specified node
- `ZIResult_enum ziAPIGetDemodSample (ZIConnection conn, const char* path, ZIDemodSample * value)`
Gets the demodulator sample value of the specified node.
- `ZIResult_enum ziAPIGetDIOSample (ZIConnection conn, const char* path, ZIDIOSample * value)`
Gets the Digital I/O sample of the specified node.
- `ZIResult_enum ziAPIGetAuxInSample (ZIConnection conn, const char* path, ZIAuxInSample * value)`
gets the AuxIn sample of the specified node
- `ZIResult_enum ziAPIGetValueB (ZIConnection conn, const char* path, unsigned char* buffer, unsigned int* length, unsigned int bufferSize)`
gets the Bytearray value of the specified node
- `ZIResult_enum ziAPISetValueD (ZIConnection conn, const char* path, ZIDoubleData value)`
asynchronously sets a double-type value to one or more nodes specified in the path
- `ZIResult_enum ziAPISetValueI (ZIConnection conn, const char* path, ZIIntegerData value)`

asynchronously sets an integer-type value to one or more nodes specified in a path

- `ZIResult_enum ziAPISetValueB (ZIConnection conn, const char* path, unsigned char* buffer, unsigned int length)`
asynchronously sets the binary-type value of one ore more nodes specified in the path
- `ZIResult_enum ziAPISyncSetValueD (ZIConnection conn, const char* path, ZIDoubleData* value)`
synchronously sets a double-type value to one or more nodes specified in the path
- `ZIResult_enum ziAPISyncSetValueI (ZIConnection conn, const char* path, ZIIntegerData* value)`
synchronously sets an integer-type value to one or more nodes specified in a path
- `ZIResult_enum ziAPISyncSetValueB (ZIConnection conn, const char* path, uint8_t* buffer, uint32_t* length, uint32_t bufferSize)`
Synchronously sets the binary-type value of one ore more nodes specified in the path.
- `ZIResult_enum ziAPISync (ZIConnection conn)`
Synchronizes the session by dropping all pending data.
- `ZIResult_enum ziAPIEchoDevice (ZIConnection conn, const char* deviceSerial)`
Sends an echo command to a device and blocks until answer is received.
- `ZIResult_enum ziAPIAsyncSetDoubleData (ZIConnection conn, const char* path, ZIDoubleData value)`
- `ZIResult_enum ziAPIAsyncSetIntegerData (ZIConnection conn, const char* path, ZIIntegerData value)`
- `ZIResult_enum ziAPIAsyncSetByteArray (ZIConnection conn, const char* path, uint8_t* buffer, uint32_t length)`
- `ZIEvent * ziAPIAllocateEventEx ()`
Allocates `ZIEvent` structure and returns the pointer to it. Attention!!! It is the client code responsibility to deallocate the structure by calling `ziAPIDeallocateEventEx`!
- `void ziAPIDeallocateEventEx (ZIEvent * ev)`
Deallocates `ZIEvent` structure created with `ziAPIAllocateEventEx()` .
- `ZIResult_enum ziAPISubscribe (ZIConnection conn, const char* path)`
subscribes the nodes given by path for `ziAPIPollDataEx`

- `ZIResult_enum` `ziAPIUnSubscribe (ZIConnection conn, const char* path)`
unsubscribes to the nodes given by path
- `ZIResult_enum` `ziAPIPollDataEx (ZIConnection conn, ZIEvent * ev, uint32_t timeOutMilliseconds)`
checks if an event is available to read
- `ZIResult_enum` `ziAPIGetValueAsPollData (ZIConnection conn, const char* path)`
triggers a value request, which will be given back on the poll event queue
- `ZIResult_enum` `ziAPIGetError (ZIResult_enum result, char** buffer, int* base)`
returns a description and the severity for a `ZIResult_enum`
- `ZIResult_enum` `ReadMEMFile (const char* filename, char* buffer, int32_t bufferSize, int32_t* bytesUsed)`

- `ziEvent *` `ziAPIAllocateEvent ()`
Deprecated: See `ziAPIAllocateEventEx()` .
- `void` `ziAPIDeallocateEvent (ziEvent * ev)`
Deprecated: See `ziAPIDeallocateEventEx()` .
- `ZIResult_enum` `ziAPIPollData (ZIConnection conn, ziEvent * ev, int timeOut)`
Checks if an event is available to read. Deprecated: See `ziAPIPollDataEx()` .
- `ZIResult_enum` `ziAPIGetValueS (ZIConnection conn, char* path, DemodSample * value)`

- `ZIResult_enum` `ziAPIGetValueDIO (ZIConnection conn, char* path, DIOSample * value)`

- `ZIResult_enum` `ziAPIGetValueAuxIn (ZIConnection conn, char* path, AuxInSample * value)`

- `double` `ziAPISecondsTimeStamp (ziTimeStampType TS)`

Detailed Description

ziAPI provides all functionality to establish a connection with the Data Server and to communicate with it. It has functions for setting and getting parameters in a single call as well as an event-

framework with which the user may subscribe the parameter tree and receive the events which occur when values change.

- All functions do not check passed pointers if they're NULL pointers. In that case a segmentation fault will occur.
- The ZIConnection is not thread-safe. One connection can only be used in one thread. If you want to use the ziAPI in a multi-threaded program you will have to use one ZIConnection for each thread that is communicating or implement a mutual exclusion.
- The Data Server is able to handle connections from threads simultaneously. The Data Server takes over the synchronization.

Data Structure Documentation

struct ZIDoubleDataTS

The structure used to hold a single IEEE double value. Same as ZIDoubleData, but with timestamp.

```
#include "ziAPI.h"

typedef struct ZIDoubleDataTS {
    ZITimeStamp
                                timeStamp
    ZIDoubleData                ;
                                value
    ;
} ZIDoubleDataTS;
```

Data Fields

- ZITimeStamp timeStamp
Time stamp at which the value has changed.
- ZIDoubleData value

struct ZIntegerDataTS

The structure used to hold a single 64bit signed integer value. Same as ZIntegerData, but with timestamp.

```
#include "ziAPI.h"

typedef struct ZIntegerDataTS {
    ZITimeStamp          timeStamp
    ZIntegerData          value
} ZIntegerDataTS;
```

Data Fields

- ZITimeStamp timeStamp
Time stamp at which the value has changed.
- ZIntegerData value

struct ZITreeChangeData

The struct is holding info about added or removed nodes.

```
#include "ziAPI.h"

typedef struct ZITreeChangeData {
    ZITimeStamp          timeStamp
    uint32_t              ;
                        action
    char                  ;
                        name
                        [32];
} ZITreeChangeData;
```

Data Fields

- ZITimeStamp timeStamp
Time stamp at which the data was updated.
- uint32_t action
field indicating which action occurred on the tree. A value of the ZITreeAction_enum.
- char name
Name of the Path that has been added, removed or changed.

struct TreeChange

The structure used to hold info about added or removed nodes. This is the version without timestamp used in API v1 compatibility mode.

```
#include "ziAPI.h"

typedef struct TreeChange {
    uint32_t
        Action
    char
        Name
    [32];
} TreeChange;
```

Data Fields

- `uint32_t Action`
field indicating which action occurred on the tree. A value of the `ZITreeAction_enum` (`TREE_ACTION`) enum.
- `char Name`
Name of the Path that has been added, removed or changed.

struct ZIDemodSample

The structure used to hold data for a single demodulator sample.

```
#include "ziAPI.h"

typedef struct ZIDemodSample {
    ZITimeStamp
        timeStamp
        ;
    double
        x
        ;
    double
        y
        ;
    double
        frequency
        ;
    double
        phase
        ;
    uint32_t
        dioBits
        ;
    uint32_t
        trigger
        ;
    double
        auxIn0
        ;
    double
        auxIn1
        ;
} ZIDemodSample;
```

Data Fields

- ZITimeStamp timeStamp
The timestamp at which the sample has been measured.
- double x
X part of the sample.
- double y
Y part of the sample.
- double frequency
Frequency at that sample.
- double phase
Phase at that sample.
- uint32_t dioBits
the current bits of the DIO.
- uint32_t trigger
trigger bits
- double auxIn0
value of Aux input 0.

- `double auxIn1`
value of Aux input 1.

struct ZIAuxInSample

The structure used to hold data for a single auxiliary inputs sample.

```
#include "ziAPI.h"

typedef struct ZIAuxInSample {
    ZITimeStamp          timeStamp
    double               ;
                        ch0
    double               ;
                        ch1
    } ZIAuxInSample;
```

Data Fields

- ZITimeStamp timeStamp
The timestamp at which the values have been measured.
- double ch0
Channel 0 voltage.
- double ch1
Channel 1 voltage.

struct ZIDIOSample

The structure used to hold data for a single digital I/O sample.

```
#include "ziAPI.h"

typedef struct ZIDIOSample {
    ZITimeStamp          timeStamp
    uint32_t             ;
                        bits
    uint32_t             ;
                        reserved
} ZIDIOSample;
```

Data Fields

- ZITimeStamp timeStamp
The timestamp at which the values have been measured.
- uint32_t bits
The digital I/O values.
- uint32_t reserved
Filler to keep 8 bytes alignment in the array of [ZIDIOSample](#) structures.

struct ZIByteArray

The structure used to hold an arbitrary array of bytes. This is the version without timestamp used in API Level 1 compatibility mode.

```
#include "ziAPI.h"

typedef struct ZIByteArray {
    uint32_t          length
    uint8_t           bytes
} ZIByteArray;
```

Data Fields

- `uint32_t length`
Length of the data readable from the Bytes field.
- `uint8_t bytes`
The data itself. The array has the size given in length.

struct ZIByteArrayTS

The structure used to hold an arbitrary array of bytes. This is the same as [ZIByteArray](#) , but with timestamp.

```
#include "ziAPI.h"

typedef struct ZIByteArrayTS {
    ZITimeStamp
        ;
    uint32_t
        length
        ;
    uint8_t
        bytes
        [0];
} ZIByteArrayTS;
```

Data Fields

- ZITimeStamp timeStamp
Time stamp at which the data was updated.
- uint32_t length
length of the data readable from the bytes field
- uint8_t bytes
the data itself. The array has the size given in length

struct ScopeWave

The structure used to hold a single scope shot (API Level 1). If the client is connected to the Data Server using API Level 4 (recommended if supported by your device class) please see [ZIScopeWave](#) instead.

```
#include "ziAPI.h"

typedef struct ScopeWave {
    double
        dt
    unsigned int
        ScopeChannel
    unsigned int
        TriggerChannel
    unsigned int
        BWLimit
    unsigned int
        Count
    short
        Data
    [0];
} ScopeWave;
```

Data Fields

- double dt
Time difference between samples.
- unsigned int ScopeChannel
Scope channel of the represented data.
- unsigned int TriggerChannel
Trigger channel of the represented data.
- unsigned int BWLimit
Bandwidth-limit flag.
- unsigned int Count
Count of samples.
- short Data
First wave data.

struct ZIScopeWave

The structure used to hold scope data. The data may be formatted differently, depending on settings. See the description of the structure members for details.

```
#include "ziAPI.h"

typedef struct ZIScopeWave {
    ZITimeStamp
        timeStamp
    ;
    ZITimeStamp
        triggerTimeStamp
    ;
    double
        dt
    ;
    uint8_t
        channelEnable
    [4];
    uint8_t
        channelInput
    [4];
    uint8_t
        triggerEnable
    ;
    uint8_t
        triggerInput
    ;
    uint8_t
        reserved0
    [2];
    uint8_t
        channelBWLimit
    [4];
    uint8_t
        channelMath
    [4];
    float
        channelScaling
    [4];
    uint32_t
        sequenceNumber
    ;
    uint32_t
        segmentNumber
    ;
    uint32_t
        blockNumber
    ;
    uint64_t
        totalSamples
    ;
    uint8_t
        dataTransferMode
    ;
    uint8_t
        blockMarker
    ;
    uint8_t
        flags
    ;
    uint8_t
        sampleFormat
    ;
    uint32_t
```



```
                                sampleCount
                                ;
int16_t                        ;
                                dataInt16
                                [0];
int32_t                        ;
                                dataInt32
                                [0];
float                          ;
                                dataFloat
                                [0];
union ZIScopeWave::@0
                                data
                                ;
} ZIScopeWave;
```

Data Fields

- ZITimeStamp timeStamp
Time stamp of the first sample in this data block.
- ZITimeStamp triggerTimeStamp
Time stamp of the trigger (may also fall between samples)
- double dt
Time difference between samples in seconds.
- uint8_t channelEnable
Up to four channels: if channel is enabled, corresponding element is non-zero.
- uint8_t channelInput
Specifies the input source for each of the scope four channels: 0 = Signal Input 1, 1 = Signal Input 2, 2 = Trigger Input 1, 3 = Trigger Input 2, 4 = Aux Output 1, 5 = Aux Output 2, 6 = Aux Output 3, 7 = Aux Output 4, 8 = Aux Input 1, 9 = Aux Input 2.
- uint8_t triggerEnable
Non-zero if trigger is enabled: Bit(0): rising edge trigger off = 0, on = 1. Bit(1): falling edge trigger off = 0, on = 1.
- uint8_t triggerInput
Trigger source (same values as for channel input)
- uint8_t reserved0
- uint8_t channelBWLimit
Bandwidth-limit flag, per channel. Bit(0): off = 0, on = 1
Bit(7...1): Reserved.
- uint8_t channelMath
Math Operation (e.g averaging) Bit (7..0): Reserved.
- float channelScaling
Data scaling factors for up to 4 channels.
- uint32_t sequenceNumber

Current scope shot sequence number. Identifies a scope shot.

- `uint32_t segmentNumber`
Current segment number.
- `uint32_t blockNumber`
Current block number from the beginning of a scope shot. Large scope shots are split into blocks, which need to be concatenated to obtain the complete scope shot.
- `uint64_t totalSamples`
Total number of samples in one channel in the current scope shot, same for all channels.
- `uint8_t dataTransferMode`
Data transfer mode `SingleTransfer = 0`, `BlockTransfer = 1`, `ContinuousTransfer = 3`, `FFTSingleTransfer = 4`. Other values are reserved.
- `uint8_t blockMarker`
Block marker: Bit (0): 1 = End marker for continuous or multi-block transfer Bit (7..0): Reserved.
- `uint8_t flags`
Indicator Flags. Bit (0): 1 = Data loss detected (samples are 0), Bit (1): 1 = Missed trigger, Bit (2): 1 = Transfer failure (corrupted data).
- `uint8_t sampleFormat`
Data format of samples: `Int16 = 0`, `Int32 = 1`, `Float = 2`, `Int16Interleaved = 4`, `Int32Interleaved = 5`, `FloatInterleaved = 6`.
- `uint32_t sampleCount`
Number of samples in one channel in the current block, same for all channels.
- `int16_t dataInt16`
- `int32_t dataInt32`
- `float dataFloat`
- `union ZIScopeWave::@0 data`
First wave data.

struct ZIPWASample

Single PWA sample value.

```
#include "ziAPI.h"

typedef struct ZIPWASample {
    double
        binPhase
        ;
    double
        x
        ;
    double
        y
        ;
    uint32_t
        countBin
        ;
    uint32_t
        reserved
        ;
} ZIPWASample;
```

Data Fields

- double binPhase
Phase position of each bin.
- double x
Real PWA result or X component of a demod PWA.
- double y
Y component of the demod PWA.
- uint32_t countBin
Number of events per bin.
- uint32_t reserved
Reserved.

struct ZIPWAWave

PWA Wave.

```
#include "ziAPI.h"

typedef struct ZIPWAWave {
    ZITimeStamp
        timeStamp
        ;
    uint64_t
        sampleCount
        ;
    uint32_t
        inputSelect
        ;
    uint32_t
        oscSelect
        ;
    uint32_t
        harmonic
        ;
    uint32_t
        binCount
        ;
    double
        frequency
        ;
    uint8_t
        pwaType
        ;
    uint8_t
        mode
        ;
    uint8_t
        overflow
        ;
    uint8_t
        commensurable
        ;
    uint32_t
        reservedUInt
        ;
    ZIPWASample
        data
        [0];
} ZIPWAWave;
```

Data Fields

- ZITimeStamp timeStamp
Time stamp at which the data was updated.
- uint64_t sampleCount
Total sample count considered for PWA.
- uint32_t inputSelect
Input selection used for the PWA.
- uint32_t oscSelect
Oscillator used for the PWA.

- `uint32_t harmonic`
Harmonic setting.
- `uint32_t binCount`
Bin count of the PWA.
- `double frequency`
Frequency during PWA accumulation.
- `uint8_t pwaType`
Type of the PWA.
- `uint8_t mode`
PWA Mode [0: zoom PWA, 1: harmonic PWA].
- `uint8_t overflow`
Overflow indicators. `overflow[0]`: Data accumulator overflow, `overflow[1]`: Counter at limit, `overflow[6..2]`: Reserved, `overflow[7]`: Invalid (missing frames).
- `uint8_t commensurable`
Commensurability of the data.
- `uint32_t reservedUInt`
Reserved unsigned int.
- [ZIPWASample](#) data
PWA data vector.

struct ZIEvent

This struct holds event data forwarded by the Data Server.

```
#include "ziAPI.h"

typedef struct ZIEvent {

    ZIValueType_enum
    valueType
    ;
    uint32_t
    count
    ;
    uint8_t
    path
    [256];
    void*
    untyped
    ;
    ZIDoubleData*
    doubleData
    ;
    ZIDoubleDataTS
    *
    doubleDataTS
    ;
    ZIIntegerData*
    integerData
    ;
    ZIIntegerDataTS
    *
    integerDataTS
    ;
    ZIByteArray
    *
    byteArray
    ;
    ZIByteArrayTS
    *
    byteArrayTS
    ;
    ZITreeChangeData
    *
    treeChangeData
    ;
    TreeChange
    *
    treeChangeDataOld
    ;
    ZIDemodSample
    *
    demodSample
    ;
    ZIAuxInSample
    *
    auxInSample
    ;
};
```

```
        ;  
        ZIDIOSample  
        *  
        dioSample  
        ;  
        ZIScopeWave  
        *  
        scopeWave  
        ;  
        ScopeWave  
        *  
        scopeWaveOld  
        ;  
        ZIPWAWave  
        *  
        pwaWave  
        ;  
union ZIEvent::@1  
        value  
        ;  
        uint8_t  
        data  
        [0x400000];  
} ZIEvent;
```

Data Fields

- [ZIValueType_enum](#) valueType
Specifies the type of the data held by the [ZIEvent](#).
- uint32_t count
Number of values available in this event.
- uint8_t path
The path to the node from which the event originates.
- void* untyped
For convenience. The void field doesn't have a corresponding data type.
- ZIDoubleData* doubleData
when valueType == ZI_VALUE_TYPE_DOUBLE_DATA
- [ZIDoubleDataTS](#)* doubleDataTS
when valueType == ZI_VALUE_TYPE_DOUBLE_DATA_TS
- ZIntegerData* integerData
when valueType == ZI_VALUE_TYPE_INTEGER_DATA
- [ZIntegerDataTS](#)* integerDataTS
when valueType == ZI_VALUE_TYPE_INTEGER_DATA_TS
- [ZIByteArray](#)* byteArray
when valueType == ZI_VALUE_TYPE_BYTE_ARRAY
- [ZIByteArrayTS](#)* byteArrayTS
when valueType == ZI_VALUE_TYPE_BYTE_ARRAY_TS

- [ZITreeChangeData](#)* treeChangeData
when valueType == ZI_VALUE_TYPE_TREE_CHANGE_DATA
- [TreeChange](#)* treeChangeDataOld
when valueType ==
ZI_VALUE_TYPE_TREE_CHANGE_DATA_OLD
- [ZIDemodSample](#)* demodSample
when valueType == ZI_VALUE_TYPE_DEMOD_SAMPLE
- [ZIAuxInSample](#)* auxInSample
when valueType == ZI_VALUE_TYPE_AUXIN_SAMPLE
- [ZIDIOSample](#)* dioSample
when valueType == ZI_VALUE_TYPE_DIO_SAMPLE
- [ZIScopeWave](#)* scopeWave
when valueType == ZI_VALUE_TYPE_SCOPE_WAVE
- [ScopeWave](#)* scopeWaveOld
when valueType == ZI_VALUE_TYPE_SCOPE_WAVE_OLD
- [ZIPWAWave](#)* pwaWave
when valueType == ZI_VALUE_TYPE_PWA_WAVE
- union ZIEvent::@1 value
Convenience pointer to allow for access to the first entry in
Data using the correct type according to [ZIEvent.valueType](#)
field.
- uint8_t data
The raw value data.

Detailed Description

[ZIEvent](#) is used to give out events like value changes or errors to the user. Event handling functionality is provided by [ziAPISubscribe](#) and [ziAPIUnSubscribe](#) as well as [ziAPIPollDataEx](#).

```
#include <stdio.h>

#include "                                ziAPI.h
                                "

void ProcessEvent(
                                ZIEvent
                                * Event )
{
    unsigned int j;

    switch( Event->
                                valueType
                                )
    {
        case
```



```
                ZI_VALUE_TYPE_DOUBLE_DATA
                :

printf( "%u elements of double data %s\n",
        Event->
            count
        ,
        Event->
            path
        );

for( j = 0; j < Event->
        count
        ; j++ )
    printf( "%f\n", Event->
        value
        .
        doubleData
        [j]);

break;

case
                ZI_VALUE_TYPE_INTEGER_DATA
                :

printf( "%u elements of integer data %s\n",
        Event->
            count
        ,
        Event->
            path
        );

for( j = 0; j < Event->
        count
        ; j++ )
    printf( "%f\n", (float)Event->
        value
        .
        integerData
        [j] );

break;

case
                ZI_VALUE_TYPE_DEMOD_SAMPLE
                :

printf( "%u elements of sample data %s\n",
        Event->
            count
        ,
        Event->
            path
        );

for( j = 0; j < Event->
        count
        ; j++ )
    printf( "TS=%f, X=%f, Y=%f\n",
        (float)Event->
            value
        .
        demodSample
        [j].
        timeStamp
```

```
        Event->
            '
            value
            .
            demodSample
[j].
            x
        Event->
            '
            value
            .
            demodSample
[j].
            y
            );

    break;

case
    ZI_VALUE_TYPE_TREE_CHANGE_DATA
    :

    printf( "%u elements of tree-changed data %s\n",
        Event->
            count
        Event->
            '
            path
            );

    for( j = 0; j < Event->
        count
        ; j++ ){

        switch( Event->
            value
            .
            treeChangeDataOld
[j].
            Action
            )
        {
            case
                ZI_TREE_ACTION_REMOVE
                :
                printf( "Tree removed: %s\n",
                    Event->
                        value
                        .
                        treeChangeDataOld
[j].
                        Name
                        );

                break;

            case
                ZI_TREE_ACTION_ADD
                :
                printf( "treeChangeDataOld added: %s\n",
                    Event->
                        value
                        .
                        treeChangeDataOld
[j].
                        Name
                        );

                break;
        }
    }
}
```

```
        case
            :
                printf( "treeChangeDataOld changed: %s\n",
                    Event->
                        value
                    .
                        treeChangeDataOld
                    [j].
                        Name
                    );
            break;
        }

    }

    break;

default:

    printf( "Unexpected event value type: %d\n", Event->
        valueType
    );

    break;

}

}
```

See Also:

[ziAPISubscribe](#) , [ziAPIUnSubscribe](#) , [ziAPIPollDataEx](#)

struct DemodSample

The [DemodSample](#) struct holds data for the ZI_DATA_DEMODSAMPLE data type. Deprecated: See [ZIDemodSample](#).

```
#include "ziAPI.h"

typedef struct DemodSample {
    ziTimeStampType
        TimeStamp
    double
        X
    double
        Y
    double
        Frequency
    double
        Phase
    unsigned int
        DIOBits
    unsigned int
        Reserved
    double
        AuxIn0
    double
        AuxIn1
} DemodSample;
```

Data Fields

- ziTimeStampType TimeStamp
- double X
- double Y
- double Frequency
- double Phase
- unsigned int DIOBits
- unsigned int Reserved
- double AuxIn0

- double AuxIn1

struct AuxInSample

The [AuxInSample](#) struct holds data for the ZI_DATA_AUXINSAMPLE data type. Deprecated: See [ZIAuxInSample](#).

```
#include "ziAPI.h"

typedef struct AuxInSample {
    ziTimeStampType      TimeStamp
    double               ;
    double               Ch0
    double               ;
    double               Ch1
    ;
} AuxInSample;
```

Data Fields

- ziTimeStampType TimeStamp
- double Ch0
- double Ch1

struct DIOSample

The [DIOSample](#) struct holds data for the ZI_DATA_DIOSAMPLE data type. Deprecated: See [ZIDIOSample](#).

```
#include "ziAPI.h"

typedef struct DIOSample {
    ziTimeStampType          TimeStamp
                                ;
    unsigned int             Bits
                                ;
    unsigned int             Reserved
                                ;
} DIOSample;
```

Data Fields

- `ziTimeStampType TimeStamp`
- `unsigned int Bits`
- `unsigned int Reserved`

struct ByteArrayData

The [ByteArrayData](#) struct holds data for the ZI_DATA_BYTEARRAY data type. Deprecated: See [ZIByteArray](#).

```
#include "ziAPI.h"

typedef struct ByteArrayData {
    unsigned int
                                Len
                                ;
    unsigned char
                                Bytes
                                [0];
} ByteArrayData;
```

Data Fields

- unsigned int [Len](#)
- unsigned char [Bytes](#)

struct `ziEvent`

This struct holds event data forwarded by the Data Server. Deprecated: See [ZIEvent](#).

```
#include "ziAPI.h"

typedef struct ziEvent {

    ziAPIDataType

    Type

    unsigned int    ;

    Count

    unsigned char   ;

    Path
    [256];
    union ziEvent::Val
    Val
    unsigned char   ;

    Data
    [0x400000];
} ziEvent;
```

Data Structures

- union `ziEvent::Val`

Data Fields

- `ziAPIDataType` Type
- unsigned int Count
- unsigned char Path
- union `ziEvent::Val` Val
- unsigned char Data

Detailed Description

`ziEvent` is used to give out events like value changes or errors to the user. Event handling functionality is provided by `ziAPISubscribe` and `ziAPIUnSubscribe` as well as `ziAPIPollDataEx`.

See Also:[ziAPISubscribe](#) , [ziAPIUnSubscribe](#) , [ziAPIPollDataEx](#)

```
#include <stdio.h>

#include "                ziAPI.h
                        "

void ProcessEvent(
                        ZIEvent
                        * Event )
{
    unsigned int j;

    switch( Event->
                        valueType
                        )
    {
        case
                        ZI_VALUE_TYPE_DOUBLE_DATA
                        :

            printf( "%u elements of double data %s\n",
                    Event->
                        count
                    ,
                    Event->
                        path
                    );

            for( j = 0; j < Event->
                    count
                ; j++ )
                printf( "%f\n", Event->
                    value
                    .
                    doubleData
                    [j]);

            break;

        case
                        ZI_VALUE_TYPE_INTEGER_DATA
                        :

            printf( "%u elements of integer data %s\n",
                    Event->
                        count
                    ,
                    Event->
                        path
                    );

            for( j = 0; j < Event->
                    count
                ; j++ )
                printf( "%f\n", (float)Event->
                    value
                    .
                    integerData
                    [j] );

            break;
```

```
case
    ZI_VALUE_TYPE_DEMOD_SAMPLE
    :
    printf( "%u elements of sample data %s\n",
        Event->
            count
        ,
        Event->
            path
        );
    for( j = 0; j < Event->
        count
        ; j++ )
        printf( "TS=%f, X=%f, Y=%f\n",
            (float)Event->
                value
            ,
            demodSample
            [j].
            timeStamp
        ,
        Event->
            value
        ,
            demodSample
            [j].
            x
        ,
        Event->
            value
        ,
            demodSample
            [j].
            y
        );
    break;
case
    ZI_VALUE_TYPE_TREE_CHANGE_DATA
    :
    printf( "%u elements of tree-changed data %s\n",
        Event->
            count
        ,
        Event->
            path
        );
    for( j = 0; j < Event->
        count
        ; j++ ){
        switch( Event->
            value
            ,
            treeChangeDataOld
            [j].
            Action
        )
        {
            case
                ZI_TREE_ACTION_REMOVE
            :

```

```
        printf( "Tree removed: %s\n",
                Event->
                    value
                    .
                    treeChangeDataOld
                    [j].
                    Name
                );
        break;

    case
        ZI_TREE_ACTION_ADD
        :
        printf( "treeChangeDataOld added: %s\n",
                Event->
                    value
                    .
                    treeChangeDataOld
                    [j].
                    Name
                );
        break;

    case
        ZI_TREE_ACTION_CHANGE
        :
        printf( "treeChangeDataOld changed: %s\n",
                Event->
                    value
                    .
                    treeChangeDataOld
                    [j].
                    Name
                );
        break;
    }

}

break;

default:

    printf( "Unexpected event value type: %d\n", Event->
            valueType
        );

    break;

}

}
```

Data Structure Documentation

union ziEvent::Val

```
typedef union ziEvent::Val {  
    void*  
  
    Void  
    ;  
  
    DemodSample  
    *  
    SampleDemod  
    ;  
  
    AuxInSample  
    *  
    SampleAuxIn  
    ;  
  
    DIOSample  
    *  
    SampleDIO  
    ;  
  
    ziDoubleType*  
    Double  
    ;  
  
    ziIntegerType*  
    Integer  
    ;  
  
    TreeChange  
    *  
    Tree  
    ;  
  
    ByteArrayData  
    *  
    ByteArray  
    ;  
  
    ScopeWave  
    *  
    Wave  
    ;  
} ziEvent::Val;
```

Data Fields

- void* Void
- DemodSample* SampleDemod
- AuxInSample* SampleAuxIn
- DIOSample* SampleDIO
- ziDoubleType* Double
- ziIntegerType* Integer

- [TreeChange](#)* Tree
- [ByteArrayData](#)* ByteArray
- [ScopeWave](#)* Wave

union ziEvent::Val

```
typedef union ziEvent::Val {  
    void*  
          
        Void  
        ;  
          
        DemodSample  
        *  
        SampleDemod  
        ;  
          
        AuxInSample  
        *  
        SampleAuxIn  
        ;  
          
        DIOSample  
        *  
        SampleDIO  
        ;  
    ziDoubleType*  
        Double  
        ;  
    ziIntegerType*  
        Integer  
        ;  
          
        TreeChange  
        *  
        Tree  
        ;  
          
        ByteArrayData  
        *  
        ByteArray  
        ;  
          
        ScopeWave  
        *  
        Wave  
        ;  
} ziEvent::Val;
```

Data Fields

- void* Void
- [DemodSample](#)* SampleDemod
- [AuxInSample](#)* SampleAuxIn
- [DIOSample](#)* SampleDIO
- ziDoubleType* Double
- ziIntegerType* Integer

- `TreeChange*` Tree
- `ByteArrayData*` ByteArray
- `ScopeWave*` Wave

Enumeration Type Documentation

enum ZIResult_enum

Defines return value for all ziAPI functions. Divided into 3 regions: info, warning and error.

Enumerator:

- ZI_INFO_BASE
- ZI_INFO_SUCCESS
Success (no error)
- ZI_INFO_MAX
- ZI_WARNING_BASE
- ZI_WARNING_GENERAL
Warning (general);.
- ZI_WARNING_UNDERRUN
FIFO Underrun.
- ZI_WARNING_OVERFLOW
FIFO Overflow.
- ZI_WARNING_NOTFOUND
Value or Node not found.
- ZI_WARNING_MAX
- ZI_ERROR_BASE
- ZI_ERROR_GENERAL
Error (general)
- ZI_ERROR_USB
USB Communication failed.
- ZI_ERROR_MALLOC
Memory allocation failed.
- ZI_ERROR_MUTEX_INIT
Unable to initialise mutex.
- ZI_ERROR_MUTEX_DESTROY
Unable to destroy mutex.
- ZI_ERROR_MUTEX_LOCK
Unable to lock mutex.
- ZI_ERROR_MUTEX_UNLOCK
Unable to unlock mutex.
- ZI_ERROR_THREAD_START
Unable to start thread.

- `ZI_ERROR_THREAD_JOIN`
Unable to join thread.
- `ZI_ERROR_SOCKET_INIT`
Can't initialize socket.
- `ZI_ERROR_SOCKET_CONNECT`
Unable to connect socket.
- `ZI_ERROR_HOSTNAME`
Hostname not found.
- `ZI_ERROR_CONNECTION`
Connection invalid.
- `ZI_ERROR_TIMEOUT`
Command timed out.
- `ZI_ERROR_COMMAND`
Command internally failed.
- `ZI_ERROR_SERVER_INTERNAL`
Command failed in server.
- `ZI_ERROR_LENGTH`
Provided Buffer length is too small.
- `ZI_ERROR_FILE`
Can't open file or read from it.
- `ZI_ERROR_DUPLICATE`
There is already a similar entry.
- `ZI_ERROR_READONLY`
Attempt to set a read-only node.
- `ZI_ERROR_DEVICE_NOT_VISIBLE`
Device is not visible to the server.
- `ZI_ERROR_DEVICE_IN_USE`
Device is already connected by a different server.
- `ZI_ERROR_DEVICE_INTERFACE`
Device does currently not support the specified interface.
- `ZI_ERROR_DEVICE_CONNECTION_TIMEOUT`
Device connection timeout.
- `ZI_ERROR_DEVICE_DIFFERENT_INTERFACE`
Device already connected over a different Interface.
- `ZI_ERROR_DEVICE_NEEDS_FW_UPGRADE`
Device needs FW upgrade.
- `ZI_ERROR_ZIEVENT_DATATYPE_MISMATCH`

Trying to get data from a poll event with wrong target data type.

- ZI_ERROR_MAX
- ZI_SUCCESS
Success (no error)
- ZI_MAX_INFO
- ZI_WARNING
Warning (general);.
- ZI_UNDERRUN
FIFO Underrun.
- ZI_OVERFLOW
FIFO Overflow.
- ZI_NOTFOUND
Value or Node not found.
- ZI_MAX_WARNING
- ZI_ERROR
Error (general)
- ZI_USB
USB Communication failed.
- ZI_MALLOC
Memory allocation failed.
- ZI_MUTEX_INIT
Unable to initialize mutex.
- ZI_MUTEX_DESTROY
Unable to destroy mutex.
- ZI_MUTEX_LOCK
Unable to lock mutex.
- ZI_MUTEX_UNLOCK
Unable to unlock mutex.
- ZI_THREAD_START
Unable to start thread.
- ZI_THREAD_JOIN
Unable to join thread.
- ZI_SOCKET_INIT
Can't initialize socket.
- ZI_SOCKET_CONNECT
Unable to connect socket.

- ZI_HOSTNAME
Hostname not found.
- ZI_CONNECTION
Connection invalid.
- ZI_TIMEOUT
Command timed out.
- ZI_COMMAND
Command internally failed.
- ZI_SERVER_INTERNAL
Command failed in server.
- ZI_LENGTH
Provided Buffer length doesn't reach.
- ZI_FILE
Can't open file or read from it.
- ZI_DUPLICATE
There is already a similar entry.
- ZI_READONLY
Attempt to set a read-only node.
- ZI_MAX_ERROR

enum ZIValueType_enum

Enumerates all types that data in a [ZIEvent](#) may have.

Enumerator:

- `ZI_VALUE_TYPE_NONE`
No data type, event is invalid.
- `ZI_VALUE_TYPE_DOUBLE_DATA`
[ZIDoubleData](#) type. Use the `ZIEvent.value.doubleData` pointer to read the data of the event.
- `ZI_VALUE_TYPE_DOUBLE_DATA_TS`
[ZIDoubleDataTS](#) type. Use the `ZIEvent.value.doubleDataTS` pointer to read the data of the event.
- `ZI_VALUE_TYPE_INTEGER_DATA`
[ZIIntegerData](#) type. Use the `ZIEvent.value.integerData` pointer to read the data of the event.
- `ZI_VALUE_TYPE_INTEGER_DATA_TS`
[ZIIntegerDataTS](#) type. Use the `ZIEvent.value.integerDataTS` pointer to read the data of the event.
- `ZI_VALUE_TYPE_DEMOD_SAMPLE`
[ZIDemodSample](#) type. Use the `ZIEvent.value.demodSample` pointer to read the data of the event.
- `ZI_VALUE_TYPE_AUXIN_SAMPLE`
[ZIAuxInSample](#) type. Use the `ZIEvent.value.auxInSample` pointer to read the data of the event.
- `ZI_VALUE_TYPE_DIO_SAMPLE`
[ZIDIOSample](#) type. Use the `ZIEvent.value.dioSample` pointer to read the data of the event.
- `ZI_VALUE_TYPE_BYTE_ARRAY`
[ZIByteArray](#) type. Use the `ZIEvent.value.byteArray` pointer to read the data of the event.
- `ZI_VALUE_TYPE_BYTE_ARRAY_TS`
[ZIByteArrayTS](#) type. Use the `ZIEvent.value.byteArrayTS` pointer to read the data of the event.
- `ZI_VALUE_TYPE_TREE_CHANGE_DATA`
[ZITreeChangeData](#) type - a list of added or removed nodes. Use the `ZIEvent.value.treeChangeData` pointer to read the data of the event.
- `ZI_VALUE_TYPE_TREE_CHANGE_DATA_OLD`
[TreeChange](#) type - a list of added or removed nodes, used in v1 compatibility mode. Use the `ZIEvent.value.treeChangeDataOld` pointer to read the data of the event.

- `ZI_VALUE_TYPE_SCOPE_WAVE`
`ZIScopeWave` type. Use the `ZIEvent.value.scopeWave` pointer to read the data of the event.
- `ZI_VALUE_TYPE_SCOPE_WAVE_OLD`
`ScopeWave` type, used in v1 compatibility mode. use the `ZIEvent.value.scopeWaveOld` pointer to read the data of the event.
- `ZI_VALUE_TYPE_PWA_WAVE`
`ZIPWAWave` type. Use the `ZIEvent.value.pwaWave` pointer to read the data of the event.
- `ZI_DATA_NONE`
no data type. the `ziEvent` is invalid.
- `ZI_DATA_DOUBLE`
double data type. use the `ziEvent::Val.Double` Pointer to read the data of the event.
- `ZI_DATA_INTEGER`
integer data type. use the `ziEvent::Val.Integer` Pointer to read the data of the event.
- `ZI_DATA_DEMODSAMPLE`
`DemodSample` data type. use the `ziEvent::Val.Sample` Pointer to read the data of the event.
- `ZI_DATA_SCOPEWAVE`
`ScopeWave` data type. use the `ziEvent::Val.Wave` Pointer to read the data of the event.
- `ZI_DATA_AUXINSAMPLE`
`MiscADValue` data type. use the `ziEvent::Val.ADValue` Pointer to read the data of the event.
- `ZI_DATA_DIOSAMPLE`
`DIOValue` data type. use the `ziEvent::Val.DIOValue` Pointer to read the data of the event.
- `ZI_DATA_BYTEARRAY`
`ByteArray` data type. use the `ziEvent::Val.ByteArray` Pointer to read the data of the event.
- `ZI_DATA_TREE_CHANGED`
a list of added or removed trees. use the `ziEvent::Val.Tree` Pointer to read the data of the event.

enum ZITreeAction_enum

Defines the actions that are performed on a tree, as returned in the [ZITreeChangeData::action](#) or [ZITreeChangeDataOld::action](#).

Enumerator:

- `ZI_TREE_ACTION_REMOVE`
A node has been removed.
- `ZI_TREE_ACTION_ADD`
A node has been added.
- `ZI_TREE_ACTION_CHANGE`
A node has been changed.

enum ZIAPIVersion_enum

Enumerator:

- ZI_API_VERSION_1
- ZI_API_VERSION_4

enum ZIListNodes_enum

Defines the values of the flags used in [ziAPIListNodes](#) .

Enumerator:

- `ZI_LIST_NODES_NONE`
Default, return a simple listing of the given node immediate descendants.
- `ZI_LIST_NODES_RECURSIVE`
List the nodes recursively.
- `ZI_LIST_NODES_ABSOLUTE`
Return absolute paths.
- `ZI_LIST_NODES_LEAFONLY`
Return only leaf nodes, which means the nodes at the outermost level of the tree.
- `ZI_LIST_NODES_SETTINGSONLY`
Return only nodes which are marked as setting.
- `ZI_LIST_NONE`
Default, return a simple listing of the given node immediate descendants.
- `ZI_LIST_RECURSIVE`
List the nodes recursively.
- `ZI_LIST_ABSOLUTE`
Return absolute paths.
- `ZI_LIST_LEAFONLY`
Return only leaf nodes, which means the nodes at the outermost level of the tree.
- `ZI_LIST_SETTINGSONLY`
Return only nodes which are marked as setting.

enum TREE_ACTION

TREE_ACTION defines the values for the [TreeChange::Action](#) Variable.

Enumerator:

- TREE_ACTION_REMOVE
a tree has been removed
- TREE_ACTION_ADD
a tree has been added
- TREE_ACTION_CHANGE
a tree has changed

Function Documentation

ziAPIInit

ZIResult_enum **ziAPIInit** (**ZIConnection** * conn)

Initializes a **ZIConnection** structure.

This function initializes the structure so that it is ready to connect to Data Server. It allocates memory and sets up the infrastructure needed.

Parameters:

[out] conn
Pointer to **ZIConnection** that is to be initialized

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_MALLOC on memory allocation failure

See Also:

[ziAPIDestroy](#) , [ziAPIConnect](#) , [ziAPIDisconnect](#)

See [Connection](#) for an example

ziAPIDestroy

ZIResult_enum ziAPIDestroy (ZIConnection conn)

Destroys a [ZIConnection](#) structure.

This function frees all memory that has been allocated by [ziAPIInit](#) . If it is called with an uninitialized [ZIConnection](#) struct it may result in segmentation faults as well when it is called with a struct for which ZIAPIDestroy already has been called.

Parameters:

[in] conn

Pointer to [ZIConnection](#) struct that has to be destroyed

Returns:

- ZI_INFO_SUCCESS

See Also:

[ziAPIInit](#) , [ziAPIConnect](#) , [ziAPIDisconnect](#)

See [Connection](#) for an example

ziAPIConnect

ZIResult_enum ziAPIConnect (**ZIConnection** conn, const char* hostname, uint16_t port)

Connects the ZIConnection to Data Server.

Connects to Data Server using a **ZIConnection** and prepares for data exchange. For most cases it is enough to just give a reference to the connection and give NULL for hostname and 0 for the port, so it connects to localhost on the default port.

Parameters:

[in] conn

Pointer to **ZIConnection** with which the connection should be established

[in] hostname

Name of the Host to which it should be connected, if NULL "localhost" will be used as default

[in] port

The Number of the port to connect to. If 0, default port of the local Data Server will be used (8005)

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_HOSTNAME if the given host name could not be found
- ZI_ERROR_SOCKET_CONNECT if no connection could be established
- ZI_ERROR_OVERFLOW when a FIFO overflow occurred
- ZI_ERROR_SOCKET_INIT if initialization of the socket failed
- ZI_ERROR_CONNECTION when the Data Server didn't return the correct answer
- ZI_ERROR_TIMEOUT when initial communication timed out

See Also:

[ziAPIDisconnect](#) , [ziAPIInit](#) , [ziAPIDestroy](#)

See [Connection](#) for an example

ziAPIDisconnect

ZIResult_enum **ziAPIDisconnect** (**ZIConnection** conn)

Disconnects an established connection.

Disconnects from Data Server. If the connection has not been established and the function is called it returns without doing anything.

Parameters:

[in] conn
Pointer to ZIConnection to be disconnected

Returns:

■ ZI_INFO_SUCCESS

See Also:

[ziAPIConnect](#), [ziAPIInit](#), [ziAPIDestroy](#)

See [Connection](#) for an example

ziAPIListImplementations

ZIResult_enum **ziAPIListImplementations** (char* implementations, uint32_t bufferSize)

Returns the list of supported implementations.

Returned names are defined by implementations in the linked library and may change depending on software version.

Parameters:

[out] implementations

Pointer to a buffer receiving a newline-delimited list of the names of all the supported ziAPI implementations. The string is zero-terminated.

[in] bufferSize

The size of the buffer assigned to the implementations parameter

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_LENGTH if the length of the char-buffer given by MaxLen is too small for all elements

See Also:

[ziAPIConnectEx](#)

ziAPIConnectEx

ZIResult_enum ziAPIConnectEx (**ZIConnection** conn, const char* hostname, uint16_t port, ZIAPIVersion_enum apiLevel, const char* implementation)

Connects to Data Server and enables extended ziAPI.

With apiLevel=ZI_API_VERSION_1 and implementation=NULL, this call is equivalent to plain [ziAPIConnect](#) . With other version and implementation values enables corresponding ziAPI extension and connection using different implementation.

Parameters:

[in] conn

Pointer to the ZIConnection with which the connection should be established

[in] hostname

Name of the host to which it should be connected, if NULL "localhost" will be used as default

[in] port

The number of the port to connect to. If 0 the port of the local Data Server will be used

[in] apiLevel

Specifies the ziAPI compatibility level to use for this connection (1 or 4).

[in] implementation

Specifies implementation to use for a connection, must be one of the returned by ziAPIListImplementations or NULL to select default implementation

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_HOSTNAME if the given host name could not be found
- ZI_ERROR_SOCKET_CONNECT if no connection could be established
- ZI_ERROR_OVERFLOW when a FIFO overflow occurred
- ZI_ERROR_SOCKET_INIT if initialization of the socket failed
- ZI_ERROR_CONNECTION when the Data Server didn't return the correct answer or requested implementation is not found or doesn't support requested ziAPI level
- ZI_ERROR_TIMEOUT when initial communication timed out

See Also:

[ziAPIListImplementations](#) , [ziAPIConnect](#) , [ziAPIDisconnect](#) , [ziAPIInit](#) , [ziAPIDestroy](#) , [ziAPIGetConnectionVersion](#)

See [Connection](#) for an example

ziAPIGetConnectionAPILevel

ZIResult_enum ziAPIGetConnectionAPILevel (**ZIConnection** conn,
ZIAPIVersion_enum* apiLevel)

Returns ziAPI level used for the connection conn.

Parameters:

[in] conn

Pointer to ZIConnection

[out] apiLevel

Pointer to preallocated ZIAPIVersion_enum, receiving the ziAPI level

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_CONNECTION if level can not be determined due to conn is not connected

See Also:

[ziAPIConnectEx](#), [ziAPIGetVersion](#)

ziAPIGetRevision

ZIResult_enum ziAPIGetRevision (unsigned int* revision)

Retrieves the revision of ziAPI.

Sets an unsigned int with the revision (build number) of the ziAPI you are using.

Parameters:

[in] revision

Pointer to an unsigned int to fill up with the revision.

Returns:

- ZI_INFO_SUCCESS

ziAPIListNodes

ZIResult_enum ziAPIListNodes (**ZIConnection** conn, const char* path, char* nodes, int bufferSize, int flags)

Returns all child nodes found at the specified path.

This function returns a list of node names found at the specified path. The path may contain wildcards so that the returned nodes do not necessarily have to have the same parents. The list is returned in a null-terminated char-buffer, each element delimited by a newline. If the maximum length of the buffer (bufferSize) is not sufficient for all elements, nothing will be returned and the return value will be [ZI_LENGTH](#).

Parameters:

[in] conn

Pointer to the ZIConnection for which the node names should be retrieved.

[in] path

Path for which all children will be returned. The path may contain wildcard characters.

[out] nodes

Upon call filled with newline-delimited list of the names of all the children found. The string is zero-terminated.

[in] bufferSize

The length of the buffer used for the nodes output parameter.

[in] flags

A combination of flags (applied bitwise) as defined in [ZIListNodes_enum](#).

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_ERROR_LENGTH if the path's length exceeds [MAX_PATH_LEN](#) or the length of the char-buffer for the nodes given by bufferSize is too small for all elements
- ZI_ERROR_OVERFLOW when a FIFO overflow occurred
- ZI_ERROR_COMMAND on an incorrect answer of the server
- ZI_ERROR_SERVER_INTERNAL if an internal error occurred in Data Server
- ZI_ERROR_NOTFOUND if the given path could not be resolved
- ZI_ERROR_TIMEOUT when communication timed out

See [Tree Listing](#) for an example

See Also:

ziAPIUpdate

ziAPIUpdateDevices

ZIResult_enum `ziAPIUpdateDevices (ZIConnection conn)`

Search for the newly connected devices and update the tree.

This function forces the Data Server to search for newly connected devices and to connect to run them

Parameters:

[in] conn
Pointer to ZIConnection

Returns:

■ ZI_INFO_SUCCESS

See Also:

[ziAPIListNodes](#)

ziAPIConnectDevice

ZIResult_enum ziAPIConnectDevice (**ZIConnection** conn, const char* deviceSerial, const char* deviceInterface, const char* interfaceParams)

Connect a device to the server.

This function connects a device with deviceSerial via the specified deviceInterface for use with the server.

Parameters:

[in] conn

Pointer to the ZIConnection with which the connection should be established

[in] deviceSerial

The serial of the device to connect to, e.g., dev2100

[in] deviceInterface

The interface to use for the connection, e.g., USB|1GbE

[in] interfaceParams

Parameters for interface configuration

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_TIMEOUT when communication timed out

See Also:

[ziAPIDisconnectDevice](#) , [ziAPIConnect](#) , [ziAPIDisconnect](#) , [ziAPIInit](#)

ziAPIDisconnectDevice

ZIResult_enum **ziAPIDisconnectDevice** (**ZIConnection** conn, const char* deviceSerial)

Disconnect a device from the server.

This function disconnects a device specified by deviceSerial from the server.

Parameters:

[in] conn

Pointer to the ZIConnection with which the connection should be established

[in] deviceSerial

The serial of the device to connect to, e.g., dev2100

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_TIMEOUT when communication timed out

See Also:

[ziAPIConnectDevice](#) , [ziAPIConnect](#) , [ziAPIDisconnect](#) , [ziAPIInit](#)

ziAPIGetValueD

ZIResult_enum ziAPIGetValueD (**ZIConnection** conn, const char* path, **ZIDoubleData*** value)

gets the double-type value of the specified node

This function retrieves the numerical value of the specified node as an double-type value. The value first found is returned if more than one value is available (a wildcard is used in the path).

Parameters:

[in] conn

Pointer to ZIConnection with which the value should be retrieved

[in] path

Path to the node holding the value

[out] value

Pointer to a double in which the value should be written

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_ERROR_LENGTH if the path's length exceeds MAX_PATH_LEN
- ZI_ERROR_OVERFLOW when a FIFO overflow occurred
- ZI_ERROR_COMMAND on an incorrect answer of the server
- ZI_ERROR_SERVER_INTERNAL if an internal error occurred in Data Server
- ZI_ERROR_NOTFOUND if the given path could not be resolved or no value is attached to the node
- ZI_ERROR_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>

#include "
                                ziAPI.h
                                "

void UpdateValue( ZIConnection Conn )
{
    ZIResult_enum RetVal;
    char* ErrBuffer;
    ZIDoubleData ValueD;

    if( ( RetVal =
                                ziAPISetValueI
                                ( Conn,
```

```
        "DEV1046/demods/*/rate",
        100 ) ) !=
        ZI_INFO_SUCCESS
    )
{
    ziAPIGetError
    ( RetVal, &ErrBuffer, NULL );
    fprintf( stderr, "Can't set Parameter: %s\n", ErrBuffer );
}

if( ( RetVal =
    ziAPIGetValueD
    ( Conn,
        "DEV1046/demods/0/rate",
        &ValueD ) ) !=
    ZI_INFO_SUCCESS
    )
{
    ziAPIGetError
    ( RetVal, &ErrBuffer, NULL );
    fprintf( stderr, "Can't get Parameter: %s\n", ErrBuffer );
}
else
{
    printf( "Value = %f\n", ValueD );
}
}
```

See Also:

[ziAPISetValueD](#), [ziAPIGetValueAsPollData](#)

ziAPIGetValueI

ZIResult_enum ziAPIGetValueI (**ZIConnection** conn, const char* path, ZIIntegerData* value)

gets the integer-type value of the specified node

This function retrieves the numerical value of the specified node as an integer-type value. The value first found is returned if more than one value is available (a wildcard is used in the path).

Parameters:

[in] conn

Pointer to ZIConnection with which the value should be retrieved

[in] path

Path to the node holding the value

[out] value

Pointer to an 64bit integer in which the value should be written

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_ERROR_LENGTH if the path's length exceeds MAX_PATH_LEN
- ZI_ERROR_OVERFLOW when a FIFO overflow occurred
- ZI_ERROR_COMMAND on an incorrect answer of the server
- ZI_ERROR_SERVER_INTERNAL if an internal error occurred in Data Server
- ZI_ERROR_NOTFOUND if the given path could not be resolved or no value is attached to the node
- ZI_ERROR_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>

#include "
                                ziAPI.h
                                "

void UpdateValue( ZIConnection Conn )
{
    ZIResult_enum RetVal;
    char* ErrBuffer;
    ZIIntegerData ValueI;

    if( ( RetVal =
                                ziAPISetValueD
                                ( Conn,
```

```
        "DEV1046/demods/*/rate",
        5.53 ) ) !=
        ZI_INFO_SUCCESS
    )
{
    ziAPIGetError
    ( RetVal, &ErrBuffer, NULL );
    fprintf( stderr, "Can't set Parameter: %s\n", ErrBuffer );
}

if( ( RetVal =
    ziAPIGetValueI
    ( Conn,
        "DEV1046/demods/0/rate",
        &ValueI ) ) !=
        ZI_INFO_SUCCESS
    )
{
    ziAPIGetError
    ( RetVal, &ErrBuffer, NULL );
    fprintf( stderr, "Can't get Parameter: %s\n", ErrBuffer );
}
else
{
    printf( "Value = %f\n", (float)ValueI );
}
}
```

See Also:

[ziAPISetValueI](#) , [ziAPIGetValueAsPollData](#)

ziAPIGetDemodSample

ZIResult_enum ziAPIGetDemodSample (**ZIConnection** conn, const char* path, **ZIDemodSample** * value)

Gets the demodulator sample value of the specified node.

This function retrieves the value of the specified node as an **DemodSample** struct. The value first found is returned if more than one value is available (a wildcard is used in the path). This function is only applicable to paths matching DEMODS/[0-9]+/SAMPLE.

Parameters:

[in] conn

Pointer to ZIConnection with which the value should be retrieved

[in] path

Path to the node holding the value

[out] value

Pointer to a **ZIDemodSample** struct in which the value should be written

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_ERROR_LENGTH if the path's length exceeds MAX_PATH_LEN
- ZI_ERROR_OVERFLOW when a FIFO overflow occurred
- ZI_ERROR_COMMAND on an incorrect answer of the server
- ZI_ERROR_SERVER_INTERNAL if an internal error occurred in Data Server
- ZI_ERROR_NOTFOUND if the given path could not be resolved or no value is attached to the node
- ZI_ERROR_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>
```

```
#include "
                                ziAPI.h
"
```

```
void GetSample( ZIConnection Conn )
{
```

```
    ZIResult_enum RetVal;
    char* ErrBuffer;
```

```
    ZIDemodSample
```

```
        DemodSample;

    if( ( RetVal =
        ziAPIGetDemodSample
        ( Conn,
          "DEV1046/demods/0/sample",
          &DemodSample ) ) !=
        ZI_INFO_SUCCESS
        )
    {
        ziAPIGetError
        ( RetVal, &ErrBuffer, NULL );
        fprintf( stderr, "Can't get Parameter: %s\n", ErrBuffer );
    }
    else
    {
        printf( "TS = %f, X=%f, Y=%f\n",
            (float)DemodSample.
                timeStamp
            ,
            DemodSample.
                x
            ,
            DemodSample.
                y
            );
    }
}
```

See Also:

[ziAPIGetValueAsPollData](#)

ziAPIGetDIOSample

ZIResult_enum ziAPIGetDIOSample (**ZIConnection** conn, const char* path, **ZIDIOSample** * value)

Gets the Digital I/O sample of the specified node.

This function retrieves the newest available DIO sample from the specified node. The value first found is returned if more than one value is available (a wildcard is used in the path). This function is only applicable to nodes ending in "/DIOS/[0-9]+/INPUT".

Parameters:

[in] conn

Pointer to the ZIConnection with which the value should be retrieved

[in] path

Path to the node holding the value

[out] value

Pointer to a **ZIDIOSample** struct in which the value should be written

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_ERROR_LENGTH if the Path's Length exceeds MAX_PATH_LEN or the length of the char-buffer for the nodes given by MaxLen is too small for all elements
- ZI_ERROR_OVERFLOW when a FIFO overflow occurred
- ZI_ERROR_COMMAND on an incorrect answer of the server
- ZI_ERROR_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_ERROR_NOTFOUND if the given path could not be resolved or no value is attached to the node
- ZI_ERROR_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>

#include "
                                ziAPI.h
                                "

void GetSample( ZIConnection Conn )
{
    ZIResult_enum RetVal;
    char* ErrBuffer;
```

```
        zIDIOSample
        DIOSample;

    if( ( RetVal =

        ziAPIGetDIOSample
        ( Conn,
            "DEV1046/dios/0/output",
            &DIOSample ) ) !=
        ZI\_INFO\_SUCCESS
        )
    {

        ziAPIGetError
        ( RetVal, &ErrBuffer, NULL );
        fprintf( stderr, "Can't get Parameter: %s\n", ErrBuffer );
    }
    else
    {
        printf( "TS = %f, bits=%08x\n",
            (float) (DIOSample.
                timestamp
            ),
            DIOSample.
                bits
            );
    }
}
```

See Also:[ziAPIGetValueAsPollData](#)

ziAPIGetAuxInSample

ZIResult_enum ziAPIGetAuxInSample (**ZIConnection** conn, const char* path, **ZIAuxInSample** * value)

gets the AuxIn sample of the specified node

This function retrieves the newest available AuxIn sample from the specified node. The value first found is returned if more than one value is available (a wildcard is used in the path). This function is only applicable to nodes ending in "/AUXINS/[0-9]+/SAMPLE".

Parameters:

[in] conn

Pointer to the ziConnection with which the Value should be retrieved

[in] path

Path to the Node holding the value

[out] value

Pointer to an **ZIAuxInSample** struct in which the value should be written

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN or the length of the char-buffer for the nodes given by MaxLen is too small for all elements
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no value is attached to the node
- ZI_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>

#include "
                                ziAPI.h
                                "

void GetSample( ZIConnection Conn )
{
    ZIResult_enum RetVal;
    char* ErrBuffer;

    ZIAuxInSample
```

```
        AuxInSample;

    if( ( RetVal =
        ziAPIGetAuxInSample
        ( Conn,
            "DEV1046/auxins/0/sample",
            &AuxInSample ) ) !=
        ZI_INFO_SUCCESS
        )
    {
        ziAPIGetError
        ( RetVal, &ErrBuffer, NULL );
        fprintf( stderr, "Can't get Parameter: %s\n", ErrBuffer );
    }
    else
    {
        printf( "TS = %f, ch0=%f, ch1=%f\n",
            (float)AuxInSample.
                timeStamp
            ,
            AuxInSample.
                ch0
            ,
            AuxInSample.
                ch1
            );
    }
}
```

See Also:

[ziAPIGetValueAsPollData](#)

ziAPIGetValueB

ZIResult_enum ziAPIGetValueB (**ZIConnection** conn, const char* path, unsigned char* buffer, unsigned int* length, unsigned int bufferSize)

gets the Bytearray value of the specified node

This function retrieves the newest available DIO sample from the specified node. The value first found is returned if more than one value is available (a wildcard is used in the path).

Parameters:

[in] conn

Pointer to the ziConnection with which the value should be retrieved

[in] path

Path to the Node holding the value

[out] buffer

Pointer to a buffer to store the retrieved data in

[out] length

Pointer to an unsigned int to store the length of data in. if an error occurred or the length of the passed buffer doesn't reach a zero will be returned

[in] bufferSize

The length of the passed buffer

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN or the length of the char-buffer for the nodes given by MaxLen is too small for all elements
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no value is attached to the node
- ZI_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>
```

```
#include "
```

```
    ziAPI.h
```

```
"
```

```
void PrintVersion( ZIConnection Conn )
{

    ZIResult_enum RetVal;
    char* ErrBuffer;

    const char* Path = "ZI/ABOUT/VERSION";
    unsigned char Buffer[ 0xff ];
    unsigned int Length;

    if( ( RetVal =
        ziAPIGetValueB
        ( Conn,
          Path,
          Buffer,
          &Length,
          sizeof(Buffer) - 1 ) ) !=
        ZI_INFO_SUCCESS
      )
    {

        ziAPIGetError
        ( RetVal, &ErrBuffer, NULL );
        fprintf( stderr, "Can't get value: %s\n", ErrBuffer );
    }
    else
    {
        Buffer[ Length ] = 0;
        printf( "%s=\"%s\"\n", Path, Buffer );
    }
}
```

See Also:

[ziAPISetValueB](#), [ziAPIGetValueAsPollData](#)

ziAPISetValueD

ZIResult_enum ziAPISetValueD (**ZIConnection** conn, const char* path, ZIDoubleData value)

asynchronously sets a double-type value to one or more nodes specified in the path

This function sets the values of the nodes specified in path to Value. More than one value can be set if a wildcard is used. The function sets the value asynchronously which means that after the function returns you have no security to which value it is finally set nor at what point in time it is set.

Parameters:

[in] conn

Pointer to the ziConnection for which the value(s) will be set

[in] path

Path to the Node(s) for which the value(s) will be set to Value

[in] value

the double-type value that will be written to the node(s)

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_READONLY on attempt to set a read-only node
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no node given by path is able to hold values
- ZI_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>
```

```
#include "
                ziAPI.h
```

```
void UpdateValue( ZIConnection Conn )
{
```

```
    ZIResult_enum RetVal;
    char* ErrBuffer;
    ZIIntegerData ValueI;
```

```
if( ( RetVal =  
    ziAPISetValueD  
    ( Conn,  
      "DEV1046/demods/*/rate",  
      5.53 ) ) !=  
    ZI_INFO_SUCCESS  
    )  
{  
    ziAPIGetError  
    ( RetVal, &ErrBuffer, NULL );  
    fprintf( stderr, "Can't set Parameter: %s\n", ErrBuffer );  
}  
  
if( ( RetVal =  
    ziAPIGetValueI  
    ( Conn,  
      "DEV1046/demods/0/rate",  
      &ValueI ) ) !=  
    ZI_INFO_SUCCESS  
    )  
{  
    ziAPIGetError  
    ( RetVal, &ErrBuffer, NULL );  
    fprintf( stderr, "Can't get Parameter: %s\n", ErrBuffer );  
}  
else  
{  
    printf( "Value = %f\n", (float)ValueI );  
}  
}
```

See Also:

[ziAPIGetValueD](#). [ziAPISyncSetValueD](#)

ziAPISetValueI

ZIResult_enum ziAPISetValueI (**ZIConnection** conn, const char* path, ZIIntegerData value)

asynchronously sets an integer-type value to one or more nodes specified in a path

This function sets the values of the nodes specified in path to Value. More than one value can be set if a wildcard is used. The function sets the value asynchronously which means that after the function returns you have no security to which value it is finally set nor at what point in time it is set.

Parameters:

[in] conn

Pointer to the ziConnection for which the value(s) will be set

[in] path

Path to the Node(s) for which the value(s) will be set

[in] value

the int-type value that will be written to the node(s)

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_READONLY on attempt to set a read-only node
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no node given by path is able to hold values
- ZI_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>
```

```
#include "
                ziAPI.h
```

```
void UpdateValue( ZIConnection Conn )
{
```

```
    ZIResult_enum RetVal;
    char* ErrBuffer;
    ZIDoubleData ValueD;
```

```
if( ( RetVal =  
    ziAPISetValueI  
    ( Conn,  
      "DEV1046/demods/*/rate",  
      100 ) ) !=  
    ZI_INFO_SUCCESS  
    )  
{  
    ziAPIGetError  
    ( RetVal, &ErrBuffer, NULL );  
    fprintf( stderr, "Can't set Parameter: %s\n", ErrBuffer );  
}  
  
if( ( RetVal =  
    ziAPIGetValueD  
    ( Conn,  
      "DEV1046/demods/0/rate",  
      &ValueD ) ) !=  
    ZI_INFO_SUCCESS  
    )  
{  
    ziAPIGetError  
    ( RetVal, &ErrBuffer, NULL );  
    fprintf( stderr, "Can't get Parameter: %s\n", ErrBuffer );  
}  
else  
{  
    printf( "Value = %f\n", ValueD );  
}  
}
```

See Also:

[ziAPIGetValueI](#) . [ziAPISyncSetValueI](#)

ziAPISetValueB

ZIResult_enum ziAPISetValueB (**ZIConnection** conn, const char* path, unsigned char* buffer, unsigned int length)

asynchronously sets the binary-type value of one or more nodes specified in the path

This function sets the values at the nodes specified in a path. More than one value can be set if a wildcard is used. The function sets the value asynchronously which means that after the function returns you have no security to which value it is finally set nor at what point in time it is set.

Parameters:

[in] conn

Pointer to the ziConnection for which the value(s) will be set

[in] path

Path to the Node(s) for which the value(s) will be set

[in] buffer

Pointer to the byte array with the data

[in] length

Length of the data in the buffer

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_READONLY on attempt to set a read-only node
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no node given by path is able to hold values
- ZI_TIMEOUT when communication timed out

```
#include <stdlib.h>
#include <stdio.h>
```

```
#include "
                ziAPI.h
"
```

```
void ProgramCPU( ZIConnection Conn,
                unsigned char* Buffer,
                int Len )
{
```

```
ZIResult_enum RetVal;
char* ErrBuffer;

if( ( RetVal =
    ziAPISetValueB
    ( Conn,
      "DEV1046/cpus/0/program",
      Buffer,
      Len ) ) !=
    ZI_INFO_SUCCESS
    )
{
    ziAPIGetError
    ( RetVal, &ErrBuffer, NULL );
    fprintf( stderr, "Can't set Parameter: %s\n", ErrBuffer );
}
}
```

See Also:

[ziAPIGetValueB](#). [ziAPISyncSetValueB](#)

ziAPISyncSetValueD

ZIResult_enum ziAPISyncSetValueD (**ZIConnection** conn, const char* path, **ZIDoubleData*** value)

synchronously sets a double-type value to one or more nodes specified in the path

This function sets the values of the nodes specified in path to Value. More than one value can be set if a wildcard is used. The function sets the value synchronously. After returning you know that it is set and to which value it is set.

Parameters:

[in] conn

Pointer to the ziConnection for which the value(s) will be set

[in] path

Path to the Node(s) for which the value(s) will be set to value

[in] value

Pointer to a double-type containing the value to be written. When the function returns value holds the effectively written value.

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_READONLY on attempt to set a read-only node
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no node given by path is able to hold values
- ZI_TIMEOUT when communication timed out

See Also:

[ziAPIGetValueD](#) , [ziAPISetValueD](#)

ziAPISyncSetValueI

ZIResult_enum ziAPISyncSetValueI (**ZIConnection** conn, const char* path, **ZIIntegerData*** value)

synchronously sets an integer-type value to one or more nodes specified in a path

This function sets the values of the nodes specified in path to value. More than one value can be set if a wildcard is used. The function sets the value synchronously. After returning you know that it is set and to which value it is set.

Parameters:

[in] conn

Pointer to the ziConnection for which the value(s) will be set

[in] path

Path to the node(s) for which the value(s) will be set

[in] value

Pointer to a int-type containing then value to be written. when the function returns value holds the effectively written value.

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_READONLY on attempt to set a read-only node
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no node given by path is able to hold values
- ZI_TIMEOUT when communication timed out

See Also:

[ziAPIGetValueI](#) , [ziAPISetValueI](#)

ziAPISyncSetValueB

ZIResult_enum ziAPISyncSetValueB (**ZIConnection** conn, const char* path, uint8_t* buffer, uint32_t* length, uint32_t bufferSize)

Synchronously sets the binary-type value of one ore more nodes specified in the path.

This function sets the values at the nodes specified in a path. More than one value can be set if a wildcard is used. This function sets the value synchronously. After returning you know that it is set and to which value it is set.

Parameters:

[in] conn

Pointer to the ziConnection for which the value(s) will be set

[in] path

Path to the Node(s) for which the value(s) will be set

[in] buffer

Pointer to the byte array with the data

[in] length

Length of the data in the buffer

[in] bufferSize

Length of the data in the buffer

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_READONLY on attempt to set a read-only node
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no node given by path is able to hold values
- ZI_TIMEOUT when communication timed out

See Also:

[ziAPIGetValueB](#) , [ziAPISetValueB](#)

ziAPISync

ZIResult_enum ziAPISync (ZIConnection conn)

Synchronizes the session by dropping all pending data.

This function drops any data that is pending for transfer. Any data (including poll data) retrieved afterwards is guaranteed to be produced not earlier than the call to ziAPISync. This ensures in particular that any settings made prior to the call to ziAPISync have been propagated to the device, and the data retrieved afterwards is produced with the new settings already set to the hardware. Note, however, that this does not include any required settling time.

Parameters:

[in] conn

Pointer to the ZIConnection that is to be synchronized

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_TIMEOUT when communication timed out

ziAPIEchoDevice

ZIResult_enum ziAPIEchoDevice (**ZIConnection** conn, const char* deviceSerial)

Sends an echo command to a device and blocks until answer is received.

This is useful to flush all buffers between API and device to enforce that further code is only executed after the device executed a previous command. Per device echo is only implemented for HF2. For other device types it is a synonym to ziAPISync, and deviceSerial parameter is ignored.

Parameters:

[in] conn

Pointer to the ZIConnection that is to be synchronized

[in] deviceSerial

The serial of the device to get the echo from, e.g., dev2100

Returns:

- ZI_INFO_SUCCESS on success
- ZI_ERROR_TIMEOUT when communication timed out

ziAPIAsyncSetDoubleData

ZIResult_enum ziAPIAsyncSetDoubleData (**ZIConnection** conn, const char* path, ZIDoubleData value)

ziAPIAsyncSetIntegerData

ZIResult_enum ziAPIAsyncSetIntegerData (**ZIConnection** conn, const char* path, ZIIntegerData value)

ziAPIAsyncSetByteArray

[ZlResult_enum](#) ziAPIAsyncSetByteArray ([ZlConnection](#) conn, const char* path, uint8_t* buffer, uint32_t length)

ziAPIAllocateEventEx

ZIEvent * ziAPIAllocateEventEx ()

Allocates [ZIEvent](#) structure and returns the pointer to it. Attention!!! It is the client code responsibility to deallocate the structure by calling [ziAPIDeallocateEventEx](#)!

This function allocates a [ZIEvent](#) structure and returns the pointer to it. Free the memory using [ziAPIDeallocateEventEx](#).

See Also:

[ziAPIDeallocateEventEx](#)

ziAPIDeallocateEventEx

void ziAPIDeallocateEventEx ([ZIEvent](#) * ev)

Deallocates [ZIEvent](#) structure created with [ziAPIAllocateEventEx\(\)](#) .

Parameters:

[in] ev
Pointer to [ZIEvent](#) structure to be deallocated..

See Also:

[ziAPIAllocateEventEx](#)

This function is the compliment to [ziAPIAllocateEventEx\(\)](#)

ziAPISubscribe

ZIResult_enum **ziAPISubscribe** (**ZIConnection** conn, const char* path)

subscribes the nodes given by path for [ziAPIPollDataEx](#)

This function subscribes to nodes so that whenever the value of the node changes the new value can be polled using [ziAPIPollDataEx](#) . By using wildcards or by using a path that is not a leaf node but contains sub nodes, more than one leaf can be subscribed to with one function call.

Parameters:

[in] conn

Pointer to the ziConnection for which to subscribe for

[in] path

Path to the nodes to subscribe

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no node given by path is able to hold values
- ZI_TIMEOUT when communication timed out

See [Data Handling](#) for an example

See Also:

[ziAPIUnSubscribe](#) , [ziAPIPollDataEx](#) , [ziAPIGetValueAsPollData](#)

ziAPIUnSubscribe

ZIResult_enum **ziAPIUnSubscribe** (**ZIConnection** conn, const char* path)

unsubscribes to the nodes given by path

This function is the complement to [ziAPISubscribe](#) . By using wildcards or by using a path that is not a leaf node but contains sub nodes, more than one node can be unsubscribed with one function call.

Parameters:

[in] conn

Pointer to the ziConnection for which to unsubscribe for

[in] path

Path to the Nodes to unsubscribe

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no node given by path is able to hold values
- ZI_TIMEOUT when communication timed out

See [Data Handling](#) for an example

See Also:

[ziAPISubscribe](#) , [ziAPIPollDataEx](#) , [ziAPIGetValueAsPollData](#)

ziAPIPollDataEx

ZIResult_enum ziAPIPollDataEx (**ZIConnection** conn, **ZIEvent** * ev, uint32_t
timeOutMilliseconds)

checks if an event is available to read

This function returns immediately if an event is pending. Otherwise it waits for an event for up to timeOutMilliseconds. All value changes that occur in nodes that have been subscribed to or in children of nodes that have been subscribed to are sent from the Data Server to the ziAPI session. For a description of how the data are available in the struct, refer to the documentation of struct [ziEvent](#). When no event was available within timeOutMilliseconds, the ziEvent::Type field will be ZI_DATA_NONE and the ziEvent::Count field will be zero. Otherwise these fields hold the values corresponding to the event that occurred.

Parameters:

[in] conn

Pointer to the [ZIConnection](#) for which events should be received

[out] ev

Pointer to a [ZIEvent](#) struct in which the received event will be written

[in] timeOutMilliseconds

Time to wait for an event in milliseconds. If -1 it will wait forever, if 0 the function returns immediately.

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_OVERFLOW when a FIFO overflow occurred

See [Data Handling](#) for an example

See Also:

[ziAPISubscribe](#) , [ziAPIUnSubscribe](#) , [ziAPIGetValueAsPollData](#) , [ziEvent](#)

ziAPIGetValueAsPollData

ZIResult_enum ziAPIGetValueAsPollData (**ZIConnection** conn, const char* path)

triggers a value request, which will be given back on the poll event queue

Use this function to receive the value of one or more nodes as one or more events using [ziAPIPollDataEx](#), even when the node is not subscribed or no value change has occurred.

Parameters:

[in] conn

Pointer to the [ZIConnection](#) with which the value should be retrieved

[in] path

Path to the Node holding the value

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_LENGTH if the Path's Length exceeds MAX_PATH_LEN or the length of the char-buffer for the nodes given by MaxLen is too small for all elements
- ZI_OVERFLOW when a FIFO overflow occurred
- ZI_COMMAND on an incorrect answer of the server
- ZI_SERVER_INTERNAL if an internal error occurred in the Data Server
- ZI_NOTFOUND if the given path could not be resolved or no value is attached to the node
- ZI_TIMEOUT when communication timed out

See [Data Handling](#) for an example

See Also:

[ziAPISubscribe](#) , [ziAPIUnSubscribe](#) , [ziAPIPollDataEx](#)

ziAPIGetError

ZIResult_enum ziAPIGetError (**ZIResult_enum** result, char** buffer, int* base)

returns a description and the severity for a **ZIResult_enum**

This function returns a static char pointer to a description string for the given **ZIResult_enum** error code. It also provides a parameter returning the severity (info, warning, error). If the given error code does not exist a description for an unknown error and the base for an error will be returned. If a description or the base is not needed NULL may be passed.

Parameters:

[in] result

A **ZIResult_enum** for which the description or base will be returned

[out] buffer

A pointer to a char array to return the description. May be NULL if no description is needed.

[out] base

The severity for the provided Status parameter:

- ZI_INFO_BASE for infos
- ZI_WARNING_BASE for warnings
- ZI_ERROR_BASE for errors

Returns:

- ZI_SUCCESS

ReadMEMFile

[ZlResult_enum](#) ReadMEMFile (const char* filename, char* buffer, int32_t bufferSize, int32_t* bytesUsed)

ziAPIAllocateEvent

ziEvent * **ziAPIAllocateEvent** ()

Deprecated: See [ziAPIAllocateEventEx\(\)](#) .

ziAPIDeallocateEvent

void ziAPIDeallocateEvent ([ziEvent](#) * ev)

Deprecated: See [ziAPIDeallocateEventEx\(\)](#) .

ziAPIPollData

ZIResult_enum ziAPIPollData (**ZIConnection** conn, **ziEvent** * ev, int timeOut)

Checks if an event is available to read. Deprecated: See [ziAPIPollDataEx\(\)](#) .

Parameters:

[in] conn

Pointer to the [ZIConnection](#) for which events should be received

[out] ev

Pointer to a [ziEvent](#) struct in which the received event will be written

[in] timeOut

Time to wait for an event in milliseconds. If -1 it will wait forever, if 0 the function returns immediately.

Returns:

- ZI_SUCCESS on success
- ZI_CONNECTION when the connection is invalid (not connected) or when a communication error occurred
- ZI_OVERFLOW when a FIFO overflow occurred

See [Data Handling](#) for an example

See Also:

[ziAPISubscribe](#) , [ziAPIUnSubscribe](#) , [ziAPIGetValueAsPollData](#) , [ziEvent](#)

ziAPIGetValueS

ZIResult_enum ziAPIGetValueS (**ZIConnection** conn, char* path, **DemodSample *** value)

ziAPIGetValueDIO

ZIResult_enum ziAPIGetValueDIO (**ZIConnection** conn, char* path, **DIOSample *** value)

ziAPIGetValueAuxIn

ZIResult_enum ziAPIGetValueAuxIn (**ZIConnection** conn, char* path, **AuxInSample *** value)

ziAPISecondsTimeStamp

double ziAPISecondsTimeStamp (ziTimeStampType TS)

Deprecated: timestamps should instead be converted to seconds by dividing by the instrument's "clockbase". This is available as an leaf under the instruments's root "device" branch in the node hierarchy, e.g., /dev2001/clockbase.

Parameters:

[in] TS
the timestamp to convert to seconds

Returns:

The timestamp in seconds as a double

Glossary

This glossary provides easy to understand descriptions for many terms related to measurement instrumentation including the abbreviations used inside this user manual.

A

A/D	Analog to Digital See Also ADC .
AC	Alternate Current
ADC	Analog to Digital Converter
AM	Amplitude Modulation
Amplitude Modulated AFM (AM-AFM)	AFM mode where the amplitude change between drive and measured signal encodes the topography or the measured AFM variable. See Also Atomic Force Microscope .
API	Application Programming Interface
ASCII	American Standard Code for Information Interchange
Atomic Force Microscope (AFM)	Microscope that scans surfaces by means an oscillating mechanical structure (e.g. cantilever, tuning fork) whose oscillating tip gets so close to the surface to enter in interaction because of electrostatic, chemical, magnetic or other forces. With an AFM it is possible to produce images with atomic resolution. See Also Amplitude Modulated AFM , Frequency Modulated AFM , Phase modulation AFM .
AVAR	Allen Variance

B

Bandwidth (BW)	<p>The signal bandwidth represents the highest frequency components of interest in a signal. For filters the signal bandwidth is the cut-off point, where the transfer function of a system shows 3 dB attenuation versus DC. In this context the bandwidth is a synonym of cut-off frequency $f_{\text{cut-off}}$ or 3dB frequency $f_{-3\text{dB}}$. The concept of bandwidth is used when the dynamic behavior of a signal is important or separation of different signals is required.</p> <p>In the context of a open-loop or closed-loop system, the bandwidth can be used to indicate the fastest speed of the system, or the highest signal update change rate that is possible with the system.</p> <p>Sometimes the term bandwidth is erroneously used as synonym of frequency range. See Also Noise Equivalent Power Bandwidth.</p>
BNC	Bayonet Neill-Concelman Connector

C

CF	Clock Fail (internal processor clock missing)
Common Mode Rejection Ratio (CMRR)	Specification of a differential amplifier (or other device) indicating the ability of an amplifier to obtain the difference between two inputs while rejecting the components that do not differ from the signal (common mode). A high CMRR is important in applications where the signal of interest is represented by a small voltage fluctuation superimposed on a (possibly large) voltage offset, or when relevant information is contained in the voltage difference between two signals. The simplest mathematical definition of common-mode rejection ratio is: $CMRR = 20 * \log(\text{differential gain} / \text{common mode gain})$.
CSV	Comma Separated Values
D	
D/A	Digital to Analog
DAC	Digital to Analog Converter
DC	Direct Current
DDS	Direct Digital Synthesis
DHCP	Dynamic Host Configuration Protocol
DIO	Digital Input/Output
DNS	Domain Name Server
DSP	Digital Signal Processor
DUT	Device Under Test
Dynamic Reserve (DR)	The measure of a lock-in amplifier's capability to withstand the disturbing signals and noise at non-reference frequencies, while maintaining the specified measurement accuracy within the signal bandwidth.
E	
XML	Extensible Markup Language. See Also XML .
F	
FFT	Fast Fourier Transform
FIFO	First In First Out
FM	Frequency Modulation
Frequency Accuracy (FA)	Measure of an instrument's ability to faithfully indicate the correct frequency versus a traceable standard.
Frequency Modulated AFM (FM-AFM)	AFM mode where the frequency change between drive and measured signal encodes the topography or the measured AFM variable. See Also Atomic Force Microscope .

Frequency Response Analyzer (FRA)	Instrument capable to stimulate a device under test and plot the frequency response over a selectable frequency range with a fine granularity.
-----------------------------------	--

Frequency Sweeper	See Also Frequency Response Analyzer .
-------------------	--

G

Gain Phase Meter	See Also Vector Network Analyzer .
------------------	--

GPIO	General Purpose Interface Bus
------	-------------------------------

GUI	Graphical User Interface
-----	--------------------------

I

I/O	Input / Output
-----	----------------

Impedance Spectroscopy (IS)	Instrument suited to stimulate a device under test and to measure the impedance (by means of a current measurement) at a selectable frequency and its amplitude and phase change over time. The output is both amplitude and phase information referred to the stimulus signal.
-----------------------------	---

Input Amplitude Accuracy (IAA)	Measure of instrument's capability to faithfully indicate the signal amplitude at the input channel versus a traceable standard.
--------------------------------	--

Input voltage noise (IVN)	Total noise generated by the instrument and referred to the signal input, thus expressed as additional source of noise for the measured signal.
---------------------------	---

IP	Internet Protocol
----	-------------------

L

LAN	Local Area Network
-----	--------------------

LED	Light Emitting Diode
-----	----------------------

Lock-in Amplifier (LI, LIA)	Instrument suited for the acquisition of small signals in noisy environments, or quickly changing signal with good signal to noise ratio - lock-in amplifiers recover the signal of interest knowing the frequency of the signal by demodulation with the suited reference frequency - the result of the demodulation are amplitude and phase of the signal compared to the reference: these are value pairs in the complex plane (X, Y), (R, Θ).
-----------------------------	---

M

Media Access Control address (MAC address)	Refers to the unique identifier assigned to network adapters for physical network communication.
--	--

Multi-frequency (MF)	Refers to the simultaneous measurement of signals modulated at arbitrary frequencies. The objective of multi-frequency is to increase the information that can be derived from a measurement which is particularly important for one-time, non-repeating events, and to increase the speed of a measurement since different frequencies do not have to be applied one after the other. See Also Multi-harmonic .
----------------------	---

Multi-harmonic (MH)	Refers to the simultaneous measurement of modulated signals at various harmonic frequencies. The objective of multi-frequency is to increase the information that can be derived from a measurement which is particularly important for one-time, non-repeating events, and to increase the speed of a measurement since different frequencies do not have to be applied one after the other. See Also Multi-frequency .
---------------------	---

N

Noise Equivalent Power Bandwidth (NEPBW)	Effective bandwidth considering the area below the transfer function of a low-pass filter in the frequency spectrum. NEPBW is used when the amount of power within a certain bandwidth is important, such as noise measurements. This unit corresponds to a perfect filter with infinite steepness at the equivalent frequency. See Also Bandwidth .
Nyquist Frequency (NF)	For sampled analog signals, the Nyquist frequency corresponds to two times the highest frequency component that is being correctly represented after the signal conversion.

O

Output Amplitude Accuracy (OAA)	Measure of an instrument's ability to faithfully output a set voltage at a given frequency versus a traceable standard.
OV	Over Volt (signal input saturation and clipping of signal)

P

PC	Personal Computer
PD	Phase Detector
Phase-locked Loop (PLL)	Electronic circuit that serves to track and control a defined frequency. For this purpose a copy of the external signal is generated such that it is in phase with the original signal, but with usually better spectral characteristics. It can act as frequency stabilization, frequency multiplication, or as frequency recovery. In both analog and digital implementations it consists of a phase detector, a loop filter, a controller, and an oscillator.
Phase modulation AFM (PM-AFM)	AFM mode where the phase between drive and measured signal encodes the topography or the measured AFM variable. See Also Atomic Force Microscope .
PID	Proportional-Integral-Derivative
PL	Packet Loss (loss of packets of data between the instruments and the host computer)

R

RISC	Reduced Instruction Set Computer
Root Mean Square (RMS)	Statistical measure of the magnitude of a varying quantity. It is especially useful when variates are positive and negative, e.g., sinusoids, sawtooth,

square waves. For a sine wave the following relation holds between the amplitude and the RMS value: $U_{\text{RMS}} = U_{\text{PK}} / \sqrt{2} = U_{\text{PK}} / 1.41$. The RMS is also called quadratic mean.

RT Real-time

S

Scalar Network Analyzer (SNA) Instrument that measures the voltage of an analog input signal providing just the amplitude (gain) information.
See Also [Spectrum Analyzer](#), [Vector Network Analyzer](#).

SL Sample Loss (loss of samples between the instrument and the host computer)

Spectrum Analyzer (SA) Instrument that measures the voltage of an analog input signal providing just the amplitude (gain) information over a defined spectrum.
See Also [Scalar Network Analyzer](#).

SSH Secure Shell

T

TC Time Constant

TCP/IP Transmission Control Protocol / Internet Protocol

Thread An independent sequence of instructions to be executed by a processor.

Total Harmonic Distortion (THD) Measure of the non-linearity of signal channels (input and output)

TTL Transistor to Transistor Logic level

U

UHF Ultra-High Frequency

UHS Ultra-High Stability

USB Universal Serial Bus

V

VCO Voltage Controlled Oscillator

Vector Network Analyzer (VNA) Instrument that measures the network parameters of electrical networks, commonly expressed as s-parameters. For this purpose it measures the voltage of an input signal providing both amplitude (gain) and phase information. For this characteristic an older name was gain phase meter.
See Also [Gain Phase Meter](#), [Scalar Network Analyzer](#).

X

XML Extensible Markup Language: Markup language that defines a set of rules for encoding documents in a format that is both human-readable and machine-readable.

Z

ZCtrl	Zurich Instruments Control bus
ZoomFFT	This technique performs FFT processing on demodulated samples, for instance after a lock-in amplifier. Since the resolution of an FFT depends on the number of point acquired and the spanned time (not the sample rate), it is possible to obtain very highly resolution spectral analysis.
ZSync	Zurich Instruments Synchronization bus

Index

A

API

- Compatibility, 10
- Levels, 10
- Versions, 10

Asynchronous commands, 20

C

C API (see ziAPI)

C Programming Language (see ziAPI)

Comparison of LabOne APIs, 9

D

Data Server, 6

Device Settings Module, 37- 37

L

LabOne

- API overview, 9
- Comparison of APIs, 9

LabVIEW, 113- 121

- Comparison to other interfaces, 9
- Concepts, 116
- Examples, finding, 118
- Examples, running, 118
- Finding examples, 118
- Finding help, 117
- Getting started, 116
- Installing the API, 114
 - Linux, 114
 - Windows, 114
- LabOne VI Palette, 116
- Palette, LabOne, 116
- Running examples, 118
- Tips and tricks, 121
- VI Palette, 116

M

Matlab, 42- 69

- Built-in help, 47
- Command reference, 54
- Comparison to other interfaces, 9
- Contents of the API package, 46
- Examples, running, 47
- Getting started, 46
- Help, accessing , 47
- Installing the API, 43
- List of Examples, 46
- List of Utility functions, 46
- Modules, 49
- Modules, configuring, 49
- Reference, 54
- Requirements, 43

Running examples, 47

Tips and tricks, 50

Troubleshooting, 52

Verifying correct configuration, 45

N

Node

- Concept, 12
- Leaf, 12
- Node hierarchy, 12

P

PLL Advisor Module, 38

Polling Data

- Concept, 12

Python, 70- 112

- Built-in help, 75
- Command reference, 80
 - zhinst package, 80
 - zhinst's utility functions, 80
 - ziDAQRecorder class, 102
 - ziDAQServer class, 86
 - ziDAQSweeper class, 95
 - ziDAQZoomFFT class, 99
 - ziDeviceSettings class, 92
 - ziPidAdvisor class, 109
 - ziPLLAdvisor class, 106
- Comparison to other interfaces, 9
- Contents of the API package, 75
- Examples, running, 75
- Exploring the available examples, 76
- Getting started, 75
- Help, accessing , 75
- Installing the API, 71
 - Linux, 73
 - Windows, 72
- Loading data in Python, 79
- Locating the zhinst installation, 77
- Modules, 77
- Modules, configuring, 77
- Recommended python packages for ziPython, 71
- Reference, 80
- Requirements for using ziPython, 71
- Running examples, 75
- Tips and tricks, 79

S

Software Trigger Module, 32- 36

Sweeper Module, 21- 29

- Bandwidth control, 22
- Measurement data, 22
- Measurement data, averaging, 23
- Scanning mode, 21
- Settling time, 22
- Settling time, definition, 22
- Sweep parameter, 21

Sweep range, 21
Synchronous commands, 20

U

UHF
Leaf, 16

Z

ziAPI

Comparison to other interfaces, 9
ziAPI, C API functions and data types
Action (see TreeChange)
action (see ZITreeChangeData)
AuxIn0 (see DemodSample)
auxIn0 (see ZIDemodSample)
AuxIn1 (see DemodSample)
auxIn1 (see ZIDemodSample)
auxInSample (see ZIEvent)
binCount (see ZIPWAWave)
binPhase (see ZIPWASample)
Bits (see DIOsample)
bits (see ZIDIOSample)
blockMarker (see ZIScopeWave)
blockNumber (see ZIScopeWave)
BWLlimit (see ScopeWave)
byteArray (see ZIEvent)
ByteArray (see ziEvent::Val)
byteArrayTS (see ZIEvent)
Bytes (see ByteArrayData)
bytes (see ZIByteArray) (see ZIByteArrayTS)
Ch0 (see AuxInSample)
ch0 (see ZIAuxInSample)
Ch1 (see AuxInSample)
ch1 (see ZIAuxInSample)
channelBWLlimit (see ZIScopeWave)
channelEnable (see ZIScopeWave)
channelInput (see ZIScopeWave)
channelMath (see ZIScopeWave)
channelScaling (see ZIScopeWave)
commensurable (see ZIPWAWave)
count (see ZIEvent)
Count (see ScopeWave) (see ziEvent)
countBin (see ZIPWASample)
data (see ZIEvent) (see ZIPWAWave) (see ZIScopeWave)
Data (see ScopeWave) (see ziEvent)
dataFloat (see ZIScopeWave)
dataInt16 (see ZIScopeWave)
dataInt32 (see ZIScopeWave)
dataTransferMode (see ZIScopeWave)
demodSample (see ZIEvent)
DIOBits (see DemodSample)
dioBits (see ZIDemodSample)
dioSample (see ZIEvent)
Double (see ziEvent::Val)
doubleData (see ZIEvent)

doubleDataTS (see ZIEvent)
dt (see ScopeWave) (see ZIScopeWave)
flags (see ZIScopeWave)
Frequency (see DemodSample)
frequency (see ZIDemodSample) (see ZIPWAWave)
harmonic (see ZIPWAWave)
inputSelect (see ZIPWAWave)
Integer (see ziEvent::Val)
integerData (see ZIEvent)
integerDataTS (see ZIEvent)
Len (see ByteArrayData)
length (see ZIByteArray) (see ZIByteArrayTS)
MAX_EVENT_SIZE, 238
MAX_PATH_LEN, 238
mode (see ZIPWAWave)
Name (see TreeChange)
name (see ZITreeChangeData)
oscSelect (see ZIPWAWave)
overflow (see ZIPWAWave)
path (see ZIEvent)
Path (see ziEvent)
Phase (see DemodSample)
phase (see ZIDemodSample)
pwaType (see ZIPWAWave)
pwaWave (see ZIEvent)
ReadMEMFile, 338
Reserved (see DemodSample) (see DIOsample)
reserved (see ZIDIOSample) (see ZIPWASample)
reserved0 (see ZIScopeWave)
reservedUInt (see ZIPWAWave)
SampleAuxIn (see ziEvent::Val)
sampleCount (see ZIPWAWave) (see ZIScopeWave)
SampleDemod (see ziEvent::Val)
SampleDIO (see ziEvent::Val)
sampleFormat (see ZIScopeWave)
ScopeChannel (see ScopeWave)
scopeWave (see ZIEvent)
scopeWaveOld (see ZIEvent)
segmentNumber (see ZIScopeWave)
sequenceNumber (see ZIScopeWave)
TimeStamp (see AuxInSample) (see DemodSample) (see DIOsample)
timeStamp (see ZIAuxInSample) (see ZIByteArrayTS) (see ZIDemodSample) (see ZIDIOSample) (see ZIDoubleDataTS) (see ZIIntegerDataTS) (see ZIPWAWave) (see ZIScopeWave) (see ZITreeChangeData)
totalSamples (see ZIScopeWave)
Tree (see ziEvent::Val)
treeChangeData (see ZIEvent)
treeChangeDataOld (see ZIEvent)
TREE_ACTION, 290
TREE_ACTION_ADD (see TREE_ACTION)
TREE_ACTION_CHANGE (see TREE_ACTION)
TREE_ACTION_REMOVE (see TREE_ACTION)

trigger (see ZIDemodSample)
 TriggerChannel (see ScopeWave)
 triggerEnable (see ZIScopeWave)
 triggerInput (see ZIScopeWave)
 triggerTimeStamp (see ZIScopeWave)
 Type (see ziEvent)
 untyped (see ZIEvent)
 Val (see ziEvent)
 value (see ZIDoubleDataTS) (see ZIEvent) (see ZIIntegerDataTS)
 valueType (see ZIEvent)
 Void (see ziEvent::Val)
 Wave (see ziEvent::Val)
 X (see DemodSample)
 x (see ZIDemodSample) (see ZIPWASample)
 Y (see DemodSample)
 y (see ZIDemodSample) (see ZIPWASample)
 ziAPIAllocateEvent, 339
 ziAPIAllocateEventEx, 192, 331
 ziAPIAsyncSetByteArray, 172, 330
 ziAPIAsyncSetDoubleData, 170, 328
 ziAPIAsyncSetIntegerData, 171, 329
 ziAPIConnect, 130, 293
 ziAPIConnectDevice, 143, 303
 ziAPIConnectEx, 133, 296
 ziAPIDeallocateEvent, 340
 ziAPIDeallocateEventEx, 193, 332
 ziAPIDestroy, 129, 292
 ziAPIDisconnect, 131, 294
 ziAPIDisconnectDevice, 144, 304
 ziAPIEchoDevice, 169, 327
 ziAPIGetAuxInSample, 155, 313
 ziAPIGetConnectionAPILevel, 135, 298
 ziAPIGetDemodSample, 151, 309
 ziAPIGetDIOSample, 153, 311
 ziAPIGetError, 200, 337
 ziAPIGetRevision, 136, 299
 ziAPIGetValueAsPollData, 197, 336
 ziAPIGetValueAuxIn, 175, 344
 ziAPIGetValueB, 157, 315
 ziAPIGetValueD, 147, 305
 ziAPIGetValueDIO, 174, 343
 ziAPIGetValueI, 149, 307
 ziAPIGetValueS, 173, 342
 ziAPIInit, 128, 291
 ziAPIListImplementations, 132, 295
 ziAPIListNodes, 140, 300
 ziAPIPollData, 198, 341
 ziAPIPollDataEx, 196, 335
 ziAPISecondsTimeStamp, 345
 ziAPISetValueB, 163, 321
 ziAPISetValueD, 159, 317
 ziAPISetValueI, 161, 319
 ziAPISubscribe, 194, 333
 ziAPISync, 168, 326
 ziAPISyncSetValueB, 167, 325
 ziAPISyncSetValueD, 165, 323
 ziAPISyncSetValueI, 166, 324
 ziAPIUnSubscribe, 195, 334
 ziAPIUpdateDevices, 142, 302
 ZIAPIVersion_enum, 288
 ZIConnection, 125, 238
 ZIListNodes_enum, 139, 289
 ZIResult_enum, 281
 ZITreeAction_enum, 287
 ZIValueType_enum, 285
 ZI_API_VERSION_1 (see ZIAPIVersion_enum)
 ZI_API_VERSION_4 (see ZIAPIVersion_enum)
 ZI_COMMAND (see ZIResult_enum)
 ZI_CONNECTION (see ZIResult_enum)
 ZI_DATA_AUXINSAMPLE (see ZIValueType_enum)
 ZI_DATA_BYTEARRAY (see ZIValueType_enum)
 ZI_DATA_DEMODSAMPLE (see ZIValueType_enum)
 ZI_DATA_DIOSAMPLE (see ZIValueType_enum)
 ZI_DATA_DOUBLE (see ZIValueType_enum)
 ZI_DATA_INTEGER (see ZIValueType_enum)
 ZI_DATA_NONE (see ZIValueType_enum)
 ZI_DATA_SCOPEWAVE (see ZIValueType_enum)
 ZI_DATA_TREE_CHANGED (see ZIValueType_enum)
 ZI_DUPLICATE (see ZIResult_enum)
 ZI_ERROR (see ZIResult_enum)
 ZI_ERROR_BASE (see ZIResult_enum)
 ZI_ERROR_COMMAND (see ZIResult_enum)
 ZI_ERROR_CONNECTION (see ZIResult_enum)
 ZI_ERROR_DEVICE_CONNECTION_TIMEOUT (see ZIResult_enum)
 ZI_ERROR_DEVICE_DIFFERENT_INTERFACE (see ZIResult_enum)
 ZI_ERROR_DEVICE_INTERFACE (see ZIResult_enum)
 ZI_ERROR_DEVICE_IN_USE (see ZIResult_enum)
 ZI_ERROR_DEVICE_NEEDS_FW_UPGRADE (see ZIResult_enum)
 ZI_ERROR_DEVICE_NOT_VISIBLE (see ZIResult_enum)
 ZI_ERROR_DUPLICATE (see ZIResult_enum)
 ZI_ERROR_FILE (see ZIResult_enum)
 ZI_ERROR_GENERAL (see ZIResult_enum)
 ZI_ERROR_HOSTNAME (see ZIResult_enum)
 ZI_ERROR_LENGTH (see ZIResult_enum)
 ZI_ERROR_MALLOC (see ZIResult_enum)
 ZI_ERROR_MAX (see ZIResult_enum)
 ZI_ERROR_MUTEX_DESTROY (see ZIResult_enum)
 ZI_ERROR_MUTEX_INIT (see ZIResult_enum)
 ZI_ERROR_MUTEX_LOCK (see ZIResult_enum)
 ZI_ERROR_MUTEX_UNLOCK (see ZIResult_enum)
 ZI_ERROR_READONLY (see ZIResult_enum)
 ZI_ERROR_SERVER_INTERNAL (see ZIResult_enum)

ZI_ERROR_SOCKET_CONNECT (see ZIResult_enum)	(see	ZI_VALUE_TYPE_BYTE_ARRAY (see ZIValueType_enum)	(see
ZI_ERROR_SOCKET_INIT (see ZIResult_enum)		ZI_VALUE_TYPE_BYTE_ARRAY_TS (see ZIValueType_enum)	(see
ZI_ERROR_THREAD_JOIN (see ZIResult_enum)		ZI_VALUE_TYPE_DEMOD_SAMPLE (see ZIValueType_enum)	(see
ZI_ERROR_THREAD_START (see ZIResult_enum)		ZI_VALUE_TYPE_DIO_SAMPLE (see ZIValueType_enum)	(see
ZI_ERROR_TIMEOUT (see ZIResult_enum)		ZI_VALUE_TYPE_DOUBLE_DATA (see ZIValueType_enum)	(see
ZI_ERROR_USB (see ZIResult_enum)		ZI_VALUE_TYPE_DOUBLE_DATA_TS (see ZIValueType_enum)	(see
ZI_ERROR_ZIEVENT_DATATYPE_MISMATCH (see ZIResult_enum)		ZI_VALUE_TYPE_INTEGER_DATA (see ZIValueType_enum)	(see
ZI_FILE (see ZIResult_enum)		ZI_VALUE_TYPE_INTEGER_DATA_TS (see ZIValueType_enum)	(see
ZI_HOSTNAME (see ZIResult_enum)		ZI_VALUE_TYPE_NONE (see ZIValueType_enum)	
ZI_INFO_BASE (see ZIResult_enum)		ZI_VALUE_TYPE_PWA_WAVE (see ZIValueType_enum)	(see
ZI_INFO_MAX (see ZIResult_enum)		ZI_VALUE_TYPE_SCOPE_WAVE (see ZIValueType_enum)	(see
ZI_INFO_SUCCESS (see ZIResult_enum)		ZI_VALUE_TYPE_SCOPE_WAVE_OLD (see ZIValueType_enum)	(see
ZI_LENGTH (see ZIResult_enum)		ZI_VALUE_TYPE_TREE_CHANGE_DATA (see ZIValueType_enum)	(see
ZI_LIST_ABSOLUTE (see ZIListNodes_enum)		ZI_VALUE_TYPE_TREE_CHANGE_DATA_OLD (see ZIValueType_enum)	(see
ZI_LIST_LEAFSONLY (see ZIListNodes_enum)		ZI_WARNING (see ZIResult_enum)	
ZI_LIST_NODES_ABSOLUTE (see ZIListNodes_enum)	(see	ZI_WARNING_BASE (see ZIResult_enum)	
ZI_LIST_NODES_LEAFSONLY (see ZIListNodes_enum)	(see	ZI_WARNING_GENERAL (see ZIResult_enum)	
ZI_LIST_NODES_NONE (see ZIListNodes_enum)		ZI_WARNING_MAX (see ZIResult_enum)	
ZI_LIST_NODES_RECURSIVE (see ZIListNodes_enum)	(see	ZI_WARNING_NOTFOUND (see ZIResult_enum)	
ZI_LIST_NODES_SETTINGSONLY (see ZIListNodes_enum)	(see	ZI_WARNING_OVERFLOW (see ZIResult_enum)	
ZI_LIST_NONE (see ZIListNodes_enum)		ZI_WARNING_UNDERRUN (see ZIResult_enum)	
ZI_LIST_RECURSIVE (see ZIListNodes_enum)		ziCore	
ZI_LIST_SETTINGSONLY (see ZIListNodes_enum)		Tips and Tricks, 40- 40	
ZI_MALLOC (see ZIResult_enum)		zoomFFT Module, 30- 31	
ZI_MAX_ERROR (see ZIResult_enum)			
ZI_MAX_INFO (see ZIResult_enum)			
ZI_MAX_WARNING (see ZIResult_enum)			
ZI_MUTEX_DESTROY (see ZIResult_enum)			
ZI_MUTEX_INIT (see ZIResult_enum)			
ZI_MUTEX_LOCK (see ZIResult_enum)			
ZI_MUTEX_UNLOCK (see ZIResult_enum)			
ZI_NOTFOUND (see ZIResult_enum)			
ZI_OVERFLOW (see ZIResult_enum)			
ZI_READONLY (see ZIResult_enum)			
ZI_SERVER_INTERNAL (see ZIResult_enum)			
ZI_SOCKET_CONNECT (see ZIResult_enum)			
ZI_SOCKET_INIT (see ZIResult_enum)			
ZI_SUCCESS (see ZIResult_enum)			
ZI_THREAD_JOIN (see ZIResult_enum)			
ZI_THREAD_START (see ZIResult_enum)			
ZI_TIMEOUT (see ZIResult_enum)			
ZI_TREE_ACTION_ADD (see ZITreeAction_enum)			
ZI_TREE_ACTION_CHANGE (see ZITreeAction_enum)	(see		
ZI_TREE_ACTION_REMOVE (see ZITreeAction_enum)	(see		
ZI_UNDERRUN (see ZIResult_enum)			
ZI_USB (see ZIResult_enum)			
ZI_VALUE_TYPE_AUXIN_SAMPLE (see ZIValueType_enum)	(see		