

**PINILAB® WASHERS**  
**KV40 / PW41**

PW40 : 85499  
PW41 : 85508

Manuel d'Utilisation  
User Manual

Code du manuel - Manual code : **87050**  
N° lot - Lot number: 9807 - 0003

Revision du manuel: Juillet 1998 - Logiciel versions 2.xx.2Y pour PW40 et PW41 -  
Français

Manual Revision : July 1998 - Software versions 2.XX.2Y for PW40 and PW41 -  
English

Impression : 03/00  
Printed : 03/00

**BIO-RAD**

## TABLE OF CONTENTS

About this manual	4	3.7.3	Programming Synopsis	3/14
SAFETY INFORMATION	5	3.7.4	ADD	3/15
SAFETY SYMBOLS	6	3.7.4.1	ADD : KIT	3/15
		3.7.4.2	ERRORS on KIT programming	3/16
		3.7.4.3	ADD : METHOD	3/17
		3.7.4.4	ERROR on	
<b>1• FEATURES AND SPECIFICATIONS</b>	<b>1/1</b>	3.7.4.5	METHOD programming	3/18
PW40 Basic features	1/1	3.7.5	ADD : PLATE	3/19
PW40 Outside dimensions [mmV]	1/2	3.7.5	INSERT	3/23
PW40 Hydraulic circuit	1/3	3.7.5.1	INSERT : KIT	3/23
PW41 Basic features	1/4	3.7.5.2	TEMPORARY KIT	3/23
PW41 Outside dimensions [mm]	1/5	3.7.5.3	INSERT : METHOD	3/23
PW41 Hydraulic circuit	1/6	3.7.5.4	INSERT : PLATE	3/24
Technical features	1/7	3.7.6	DELETE	3/25
Electrical Data	1/7	3.7.6.1	DELETE : KIT	3/25
Physical data	1/7	3.7.6.2	DELETE : METHOD	3/25
Hardware specifications	1/7	3.7.6.3	DELETE: PLATE	3/26
Software specifications	1/8	3.7.7	EDIT	3/27
Software parameters for Wash protocols	1/9	3.7.7.1	EDIT : KIT	3/27
WASH protocol definition :	1/9	3.7.7.2	EDIT : PLATE	3/29
Methods	1/10	3.7.8	COPY	3/30
Parameter's definitions	1/12	3.7.8.1	COPY : KIT	3/30
		3.7.8.2	ERROR on KIT copy	3/30
		3.7.8.3	COPY : PLATE	3/31
		3.7.8.4	ERROR on PLATE copy	3/31
		3.7.9	Plate Test	3/32
		3.7.10	RAM initialization	3/33
<b>2» INSTALLATION</b>	<b>2/1</b>			
Unpacking	2/1			
<b>3» OPERATING INSTRUCTIONS</b>	<b>3/1</b>	<b>4» MAINTENANCE &amp; TROUBLESHOOTING</b>	<b>4/1</b>	
<b>Run Mode</b>	<b>3/1</b>	4.1	Routine maintenance	4/1
How to use the Keyboard?	3/1		Procedures	4/1
Microplate loading	3/2	4.1.1	Daily Maintenance	4/1
Select one KIT	3/3	4.1.2	Washer storage conditions	4/1
WASH Interruption or Power Failure	3/5	4.2	Cleaning and disinfection	4/2
WASH Interruption	3/5	4.2.1	Cleaning the manifold	4/2
POWER Failure	3/6	4.2.2	Disinfection	4/3
<b>RINSE</b>	3/6	4.3	Troubleshooting	4/4
RINSE (with PW40)	3/6	4.3.1	Hardware and software version	4/4
Please RINSE (with PW40)	3/6			
RINSE (with PW41)	3/7			
Automatic Rinsing (with PW41)	3/7			
<b>Configuration</b>	3/8			
Main menu	3/8			
Washer's configuration	3/9			
<b>PROGRAMMING</b>	<b>3/12</b>			
Kit Structure	3/12			
Access to PROGRAMMING mode	3/13			
		<b>Appendix</b>	<b>A/1</b>	
		Accessories and spare part list	A/1	
		PW40 Wash protocols	A/2	
		PW41 Wash protocols	A/6	
		Recommendations for optimized washing conditions	A/10	
		Revision Upgrade Description	A/11	
		Index	A/12	

**IMPORTANT**  
Please read carefully before using the PLATE WASHERS PW40 or PW41

This document is the copyright of SANOFI DIAGNOSTICS PASTEUR and must not be copied or reproduced in any form without prior consent.

SANOFI DIAGNOSTICS PASTEUR reserves the right to make technical improvements to this equipment and documentation without prior notice as part of a continuous program or product development. This manual supersedes all previous editions.

## About this manual

This manual has been written for the laboratory technician and describes how to use the PW40 and PW41 plate washers.

This manual is designed to give you the information you need to :

- install the PW40 or the PW41
- operate with the PW40 or the PW41 to suit your particular needs
- perform the PW40 or the PW41 basic maintenance procedures, including disinfection.

This manual also describes all the features and specifications of the PW40 and PW41 hardware and configuration software.

You will find explanations of all the error messages and a problem solving guide in the PW troubleshooting guide.

---

## SAFETY INFORMATION

---

**WARNING :**  
**ELECTRIC**  
**SHOCK**  
**HAZARD**

**To operate, this instrument must be mandatory connected to a wall socket equipped with earth.**

**Furthermore the electrical line must be protected by an adapted fuse and a differential circuit breaker.**

**The access to the mains supply input plug, located on instrument rear panel, must be kept free, in order to use this plug, if necessary, as cut-off device.**

Although this instrument is fully insulated and grounded, it is important for all users to be aware of the potential hazards of using liquids in close proximity to an electrical supply. In the event of any spillage, the instrument must be immediately isolated from the electrical supply and drained. Spilt liquid should be removed as soon as possible.

Serious damage to the instrument may occur if fluid is allowed to penetrate the case. The electrical supply **MUST NOT** be reconnected until the instrument has been fully inspected by a qualified technician or service engineer.

**SERVICING MUST ONLY BE CARRIED OUT BY A QUALIFIED TECHNICIAN OR SERVICE ENGINEER.**

**WARNING :**  
**SPILLAGES**

In some operating modes (e.g. dispensing fluids or purging the system, particularly under automatic control) it is possible to overflow the wells of a microplate. This can result in spillage of liquids and/or contamination with hazardous substances. Care should be taken to ensure that overflowing does not occur.

If environmentally hazardous fluid is spilled, the instrument must first be isolated from the electrical supply and physically isolated from the working environment by placing it in a fume cupboard or similar location. Contact the nearest service centre for assistance, advising them of the health hazards associated with the spilled materials.

**WARNING :**  
**INCORRECT**  
**OPERATION**

**OPERATING THIS EQUIPMENT IN WAYS OTHER THAN DETAILED IN THIS MANUAL MAY IMPAIR THE PROTECTION PROVIDED BY THE INSTRUMENT.**

**DO NOT** operate this equipment with the covers removed - potentially dangerous voltages are contained within.

**DO NOT** operate this equipment with the safety earth (ground) disconnected.

**CAUTION**

**DO NOT** poke fingers, pencils or other implements into plate carrier entry slot.

## SAFETY INFORMATION

BE PREPARED for unexpected movement of the plate carrier and noise from the Washer pump when the Washer is controlled by external computer.

ENSURE that the main power cable is correctly wired. Colour codes are as follows:


Europe	United States
BROWN	LIVE
BLUE	NEUTRAL
GREEN/YELLOW	EARTH (Ground)
	BLACK LIVE
	WHITE NEUTRAL
	GREEN GROUND

## SAFETY SYMBOLS

The following alarm and information symbols can be found in several places in the washer. Only a qualified personal aware of the security procedures is allowed to repair this instrument.  
Please read this manual carefully before using this instrument.

 Direct current

 Alternating current


 Direct and alternating current

 Earth (ground) terminal

 Protective earth (ground) terminal

 On

 Off

 Caution (see enclosed documents)

ENSURE that all connections are tight and all blocking elements (foams,...) have been removed before switching on the Washer.

If there is any doubt or concern about the safety of the instrument contact the nearest approved service centre.

# 1 • FEATURES AND SPECIFICATIONS

## PW40 1.1 PW40 Basic features





Fig. 1



## 1.2 PW40 Outside dimensions [mm]

### PW40 1.2 PW40 Outside dimensions [mm]



Washer's height: 243 mm (without bottles)

Fig. 2

1«

FEATURES AND SPECIFICATIONS

PW40 1.3 PW40 Hydraulic circuit




Fig. 3

PW40 Hydraulic Circuit

1»

FEATURES AND SPECIFICATIONS

PW41 1.4 PW41 Basic features




Fig. 4

1•

FEATURES AND SPECIFICATIONS

PW41 15 PW41 Outside dimensions [mm]





Fig. 5

FEATURES AND SPECIFICATIONS

PW41

1.6 PW41 Hydraulic circuit




Fig. 6

1 N,σ sul = Q it

## » FEATURES AND SPECIFICATIONS

### 1.7 Technical features

#### 1.7.1 Electrical Data

Voltage	240/100 VAC, 50/60 Hz
Consumption	100 VA max
Fuses	2,5 AT (2 fuses, 5x20 mm)
Power Cord	on CEI socket

#### 1.7.2 Physical data

Dimensions, including bottles	
Width	551 mm, PW40 822 mm, PW41
Length	530 mm, Plate carrier out, PW40 & PW41.
Height	300 mm, PW40 395 mm, PW41.
Weight (without accessories)	11,0 Kg, PW40 11,6 Kg, PW41
Operating conditions	15-30 °C /15 - 85 %R.H.

#### 1.7.3 Hardware specifications

Manifold	PW40: 8 or 12 channel (Default: 8 channel) PW41: 8, 2*8, 3*8, 12 or 2*12 channel (Default: 2*8)
Vacuum power	PW40 : 1 integrated vacuum pump (maxi 8 L/min) PW41 : 2 integrated vacuum pumps (maxi 16 L/min).

#### Volume of WASTE bottle

PW40 :	2 litres
PW41 :	8 litres

#### User interface:

Flat keyboard with 5 diaphragm keys  
2x20 characters LCD screen , backlighted  
onPW41.

1•

FEATURES AND SPECIFICATIONS

**1.7 Technical features (continued)**

**1.7.4 Software specifications  
Kits**

Up to 75 washing kits programmable.  
Each kit is saved under its own name (up to 15 characters) and is a combination of methods repeated or not, with or without soaking times, in strip mode or in plate mode.

**Method**

6 single-cycle washing methods :  
Wash, Aspiration, Dispensing, Bottom Washing, Bottom Aspiration, Agitation.  
4 two-cycle washing methods :  
Wash + Aspiration, Wash + Bottom Aspiration, Bottom Wash + Aspiration, Bottom Wash + Bottom Aspiration.

**Soak time**

0 to 9,9 sec Strip mode  
0 to 59 min 59 sec in Plate mode.

**Repetition**

Every method can be repeated from 1 up to 9 times.

**Wash mode**

Strip mode or Plate mode.

**Omit strip**

Strips can be omitted at kit programming.

**Plates**

Up to 10 plate parameter sets programmable.  
Accepts flat or curved bottom microplates.  
Vertical and horizontal speeds, vertical and horizontal positions of the aspirating needle in relation with the wells are programmable.

**Other features**

- Priming sequence of the hydraulic system at change of Wash solution.
- Disinfection program of the hydraulic circuit.
- User programming can be inhibited.

1•

FEATURES AND SPECIFICATIONS

1.8 Software parameters for Wash protocols


This section describes the Wash protocols, their methods and their software parameters. It is useful to be read prior starting to programme new wash protocols or new kits.

1.8.1 WASH protocol definition :

A WASH protocol, also called a KIT, is made of a succession of Wash Methods.

There are 10 different wash methods available and every method is programmed with its own batch of specific parameters.

Also, the Wash protocol is associated to a batch of plate parameters -Positions- that are set in respect to plate geometry (flat, U or V shaped well,...) as well as Speeds needed for the wash process



The available methods are :

Method	Method Name
Single-cycle Method	
nW	Wash
nA	Aspiration
nD	Dispensing
nW	Bottom Washing
na	Bottom aspiration
nAg	Plate Agitation

Two-cycle Method

nW + A Wash + Aspiration  
 nW + a Wash + Bottom Aspiration  
 nW + A Bottom Wash + Aspiration  
 nW + a Bottom Washing + Bottom Aspiration

Fig. 7



1»

## FEATURES AND SPECIFICATIONS

### 1.8 Software parameters for Wash protocols (continued)

#### 1.8.2 Methods

A method is a completely defined module made of one or two elementary cycles. Every kit is made of a succession of methods.

In a method, elementary cycles can be applied n times (n =1 to 9).

For method using two elementary cycles, only the first cycle is applied "n" time.

When repeated several times, the waiting time between elementary cycles is called Soaking Time.

The waiting time between two methods is called Method Interval.

Whatever the selected mode, STRIP or PLATE (see definition below), a method is always completely finished on the whole plate prior starting the next method.

#### Available methods :

Method	Abbreviated Name	Method Name
	on LCD Display	

#### Single-cycle Method

nW	WASH	Wash
nA	ASPIRATION	Aspiration
nD	DISPENSING	Dispensing
nW	BOTTOM WASH.	Bottom Washing
na	BOTTOM ASP.	Bottom aspiration
nAg	AGITATION	Plate Agitation



#### Two-cycle Method

nW + A	WASH + ASP	Wash + Aspiration
nW + a	WASH+BOT.ASP	Wash + Bottom Aspiration
nW + A	BOT.WASH+ASP	Bottom Wash + Aspiration
nW + a	B.WASH+B.ASP	Bottom Washing + Bottom Aspiration

#### Kit parameters associated to method :

##### •MODE

STRIP or PLATE

 <p>STRIP mode: the total method is applied on the strip prior processing the next strip.</p>	 <p>PLATE mode: the elementary cycle is applied successively on all strips prior applying next elementary cycle.</p>
<p>Example:</p> <p>STRIP Mode</p> <p>3W + A on 1st strip 3W + A on 2nd strip etc... 3W + A on last strip</p>	<p>PLATE Mode</p> <p>W on the whole plate W on the whole plate W on the whole plate A on the whole plate</p>

---

>

## FEATURES AND SPECIFICATIONS

### 1.8 Software parameters for Wash protocols (continued)

#### 1.8.2 Methods (continued)

**N° OF CYCLES**                      Number of Cycles = n  
1 to 9, step 1  
number of elementary cycles that will be applied within the method.

• **SOAKING**                      Soak time  
Omn Osec to 59mn Osec; step 1 sec                      in PLATE mode.  
Omn Osec to 9,9sec; step 0,1 sec                      in STRIP mode.  
Waiting time between two consecutive elementary cycle within one method.

• **MET. INTER**                      Method Interval  
Omn Osec to 59mn Osec; step 1 sec  
Waiting time between two consecutive methods.

## FEATURES AND SPECIFICATIONS

### 1.8 Software parameters for Wash protocols (continued) 1.8.3 Parameter's definitions

Parameters directly depending on KIT or Wash protocol (such as dispensing, aspirating times, liquids, etc... are called "KIT Parameters", Parameters depending on the microplate dimensions (mainly POSITION and SPEED parameters) are called "PLATE parameters".  
In the washer, the elementary cycles are as follows :

**A**

#### ASPIRATION



The aspirating needle goes down and aspirates until reaching the well bottom. The needle is shifted close to wall (*Aspiration Horizontal Position*) for flat bottom well;  
it remains centered for curved bottom well.



Programmed parameters :

**Kit parameters :**

- **CROSW. ASP.** Crosswise aspiration (CWA)  
(for flat bottom only)

YES or NO

the aspiration needle is shifted close to the wall (*Aspiration Horizontal Position*) while going down to well bottom and aspirates; once the aspiration time has elapsed, the aspiration needle is going up to approx middle of well height and is shifted to opposite end of wall (no scratch on the bottom); the needle goes down to bottom and a second aspiration sequence takes place.



• **ASP. TIME**

Aspiration time

0,1 to 9,9 sec, step 0,1 sec

Time during which the aspirating needle is kept at the well bottom and is aspirating.

**Plate parameters :**

- **ASP. HOR. POS.** Aspiration Horizontal Position  
(for flat bottom only)

0,0 (centered) to 2,0 mm, step 0,1 mm  
(see also fig. 8 at the end of §1.8)

• **ASP. VERT. POS.**

Aspiration Vertical Position

0,1 to 15,0 mm, step 0,1 mm (see also fig. 8 at the end of §1.8)

0,1 is the highest needle position

15,0 is the lowest needle position

- **ASP. DOWNW. SPEED** Aspiration Downward Speed  
0 to 9, step 1. (speed mini = 0; speed max=9)


1•

FEATURES AND SPECIFICATIONS

- 1.8 Software parameters for Wash protocols (continued)  
1.8.3 Parameter's definitions (continued)

D

**DISPENSING**



The dispensing needle is placed over the well and dispenses a wash solution in the well.  
When the dispensed volume is higher than well capacity, the aspiration needle aspirates the overflow to prevent liquid spilling on next wells.

Programmed parameters ;

**Kit parameters :**

- **VOLUME** (for one well)  
50 to 3000 l<sub>il</sub>, step 50 l<sub>il</sub>  
(the well volume is approx. 370  $\mu$ L)

• **OVERFLOW**

1,0 to 12,9 mm, step 0,1 mm (see also fig. 8 at the end of §1.8).  
This parameter sets the height of the aspirating needle and thus the height of the liquid in the well.



1,0 is the highest needle position  
12,9 is the lowest needle position (for the overflow).  
It is recommended to dispense a volume slightly higher than wished and have an overflow phase in order to obtain an equal distribution of volume in each well; the slight flow rate variations among dispensing needles are compensated by the overflow phase.

• **LIQUID**

Wash W1 to Wash W9  
When launching a kit in the RUN mode, the display will remind you to connect correct Wash bottle prior starting.

• **FLOW**

flow rate compensation

-5 to +5, step 1.

The flow rate of the dispensing pump is adapted to the operating manifold. However the flow rate can be slightly modified with this parameter (lower flow rate with negative value, higher flow rate with positive value).

**Plate parameters :**

- **DISP. UPWARD SPEED** Dispensing Upward Speed  
0 to 9, step 1. (speed mini = 0; speed max=9)

1»


## FEATURES AND SPECIFICATIONS

- 1.8 Software parameters for Wash protocols (continued)  
1.8.3 Parameter's definitions (continued)

### W

#### WASH

Aspiration sequence followed by Dispensing sequence.




Programmed parameters:

**Kit parameters:**

- **CROSW. ASP.**  
YES or NO

Crosswise aspiration  
(for flat bottom only)

See details page 1/10.



- **ASP. TIME**  
0,1 to 9,9 sec, step 0,1 sec.

Aspiration time  
See details page 1/10.

- **VOLUME**  
50 to 3000 nL, step 50 iL  
(the well volume is approx. 370 iL)

(for one well)

- **OVERFLOW**  
1,0 to 12,9 mm, step 0,1 mm  
1,0 is the highest needle position  
12,9 is the lowest needle position (for the overflow).

See details page 1/11.

- **LIQUID**  
Wash W1 to Wash W9

See details page 1/11.

- **FLOW**  
-5 to +5, step 1.

flow rate compensation  
See details page 1/11.

#### Plate parameters

- **ASP. HOR. POS.**

Aspiration Horizontal Position  
(for flat bottom only)

- 0,0 to 2,0 mm, step 0,1 mm  
(see also fig. 8, at the end of §1.8)

- **ASP.VERT.POS.**

Aspiration Vertical Position

- 0,1 to 15,0 mm, step 0,1 mm  
(see also fig. 8, at the end of §1.8)  
0,1 is the highest needle position  
15,0 is the lowest needle position

- **ASP. DOWNW. SPEED**  
0 to 9, step 1. (speed mini = 0; speed max=9)

Aspiration Downward Speed

- **DISP. UPW. SPEED**  
0 to 9, step 1. (speed mini = 0; speed max=9)

Dispensing Upward Speed

## FEATURES AND SPECIFICATIONS

- 1.8 Software parameters for Wash protocols (continued)  
1.8.3 Parameter's definitions (continued)

### **W BOTTOM WASH.**



#### **Bottom Washing**

the bottom of the well is thoroughly washed, followed by a normal Wash (W) sequence.

Following sequences are applied :

- Downward aspirating movement of the needle into the well down to *Aspiration Vertical Position*.
- Aspiration at *Aspiration Vertical Position* during *Aspiration Time*.
- One or two upward dispensing movements up to *Bottom Wash Vertical Position*, dispensing at this position during a *Bottom Wash Time*
- ...
- ... followed by downward aspirating movement down to *Aspiration Vertical Position*.
- Aspiration at *Aspiration Vertical Position* during *Aspiration Time*.
- Dispensing up to *Overflow Position*.

In case of crosswise aspiration (for flat bottom only):  
every aspiration sequence (close to the wall) is immediately followed by an other aspiration sequence to the opposite end of wall as described on page 1/10.

Programmed parameters :

**Kit parameters :**

- **CROSSW. ASP.** Crosswise aspiration  
YES or NO for flat bottom only  
See details page 1/10.

- **ASP. TIME** Aspiration time  
0,1 to 9,9 sec, step 0,1 sec. See details page 1/10.

- **VOLUME** (for one well)  
50 to 3000  $\mu$ L, step 50  $\mu$ L  
(the well volume is approx. 370  $\mu$ l)  
Volume is applied at the last Dispensing sequence.

- **OVERFLOW**  
1,0 to 12,9 mm, step 0,1 mm  
1,0 is the highest needle position  
12,9 is the lowest needle position (for the overflow).  
Overflow is applied at the last Dispensing sequence.  
See details page 1/11.

- **LIQUID** Wash W1 to Wash W9  
See details page 1/11.

- **FLOW** flow rate compensation  
-5 to +5, step 1. See details page 1/11.


## K

### FEATURES AND SPECIFICATIONS

1.8 Software parameters for Wash protocols (continued)  
1.8.3 Parameter's definitions (continued)

w **Bottom Washing** (continued)

**BOTTOM WASH.** **BOT. WASH NUMBER** Bottom Wash Number  
1 or 2  
(continued) Number of Bottom Wash Occurrence.



**BOTTOM TIME** Bottom Wash Time  
0,1 to 9,9 sec, step 0,1 sec.  
Time when dispensing at the Bottom Wash Vertical Position.

**Plate parameters :**

**ASP. HOR. POS.** Aspiration Horizontal Position  
(for flat bottom only)

0,0 to 2,0 mm, step 0,1 mm (see also fig. 8 at the end of §1.8)  
Aspiration Horizontal Position is applied for all Aspiration, Bottom Wash and Bottom Aspiration sequences.

**ASP. VERT. POS.** Aspiration Vertical Position  
0,1 to 15,0 mm, step 0,1 mm (see also fig. 8 at the end of §1.8)  
0,1 is the highest needle position  
15,0 is the lowest needle position  
Aspiration Vertical Position is applied for all Aspiration sequences.

**B.W. VERT. POS.** Bottom Wash Vertical Position  
0,1 to 15,0mm, step 0,1 mm (see also fig. 8 at the end of §1.8)  
0,1 is the highest needle position  
15,0 is the lowest needle position  
Bottom Wash Vertical Position is the height of the aspirating needle during the dispensing sequence of the bottom wash.

**ASP. DOWNW. SPEED** Aspiration Downward Speed  
0 to 9, step 1. (speed mini = 0; speed max=9)  
This Speed is applied only at the first Aspiration sequence.

**DISP. UPW. SPEED** Dispensing Upward Speed  
0 to 9, step 1. (speed mini = 0; speed max=9)  
Dispensing Upward Speed is applied at the last Dispensing sequence.

**BOT. DOWNW. SPEED** Bottom Downward Speed  
0 to 9, step 1. (speed mini = 0; speed max=9)  
Speed of all downward movements for all bottom sequences (wash and aspiration).

**BOT. UPWARD SPEED** Bottom Upward Speed  
0 to 9, step 1. (speed mini = 0; speed max=9)  
Speed of all upward movements for all bottom sequences (wash and aspiration).



## FEATURES AND SPECIFICATIONS

### 1.8 Software parameters for Wash protocols (continued)

#### 1.8.3 Parameter's definitions (continued)

##### **2 Bottom aspiration**

**BOTTOM ASP.** following sequences are applied :

- Downward aspirating movement of the needle into the well down to *Aspiration Vertical Position*.
- Aspiration at *Aspiration Vertical Position* during *Aspiration Time*.
- One or two upward aspirating movements up to *Bottom (aspiration) Position* followed by downward aspirating movement down to *Aspiration Vertical Position*.
- Aspiration at *Aspiration Vertical Position* during *Aspiration Time*.



In case of crosswise aspiration (for flat bottom only) :  
every aspiration sequence (close to the wall) is immediately followed by an other aspiration sequence to the opposite end of wall as described on page 1/10.

Programmed parameters :

**Kit parameters :**

- **CROSW. ASP.**  
YES or NO

Crosswise aspiration  
for flat bottom only  
See details page 1/10

- **ASP. TIME**  
0,1 to 9,9 sec, step 0,1 sec.

Aspiration time  
See details page 1/10.

- **BOT. ASP. NUMBER**  
1 or 2  
Number of Bottom Aspiration Occurrence.

Bottom Aspiration Number

**Plate parameters :**

- **ASP. HOR. POS.**  
Aspiration Horizontal Position  
(for flat bottom only)

0,0 to 2,0 mm, step 0,1 mm (see also fig. 8 at the end of §1.8)  
Aspiration Horizontal Position is applied for all Aspiration, Bottom Wash and Bottom Aspiration sequences.

- **ASP.VERT.POS.**  
Aspiration Vertical Position  
0,1 to 15,0 mm, step 0,1 mm (see also fig. 8 at the end of §1.8)  
0,1 is the highest needle position  
15,0 is the lowest needle position  
Aspiration Vertical Position is applied for all Aspiration sequences.



## FEATURES AND SPECIFICATIONS

- 1.8 Software parameters for Wash protocols (continued)  
1.8.3 Parameter's definitions (continued)

**Bottom aspiration** (continued)

**BOTTOM ASP.**  
(continued)

**Plate parameters** (continued):

• **BOT.VERT.POS.**

Bottom (Aspiration) Vertical Position

0,1 to 15,0mm, step 0,1mm (see also fig. 8 at the end of \$1.8)

0,1 is the highest needle position

15,0 is the lowest needle position.

Bottom (Aspiration) Vertical Position is the height of the aspirating needle during the aspiration sequence of the bottom aspiration.



• **ASP. DOWNW. SPEED**

Aspiration Downward Speed

0 to 9, step 1. (speed mini = 0; speed max=9)

Aspiration Downward Speed is applied at the first Aspiration sequence.

• **BOT. DOWNW. SPEED**

Bottom Downward Speed

0 to 9, step 1. (speed mini = 0; speed max=9)

Speed of all downward movements for all bottom sequences (wash and aspiration).

• **BOT. UPWARD SPEED**

Bottom Upward Speed

0 to 9, step 1. (speed mini = 0; speed max=9)

Speed of all upward movements for all bottom sequences (wash and aspiration).

---


## FEATURES AND SPECIFICATIONS

### 1.8 Software parameters for Wash protocols (continued)

#### 1.8.3 Parameter's definitions (continued)

##### Ag **Plate Agitation**

The microplate is shaken horizontally during the *Agitation T/m* with an *Agitation Amplitude* and an *Agitation Speed*.



Programmed parameters :

**Kit parameters :**

- **SHAKE TIME**  
0,1 to 59,9 sec, step 0,1 sec.

Agitation time

**Plate parameters :**

- **SHAKING AMPLITUDE**  
0 to 9, step 1
- **SHAKING SPEED**  
0 to 9, step 1

Agitation Amplitude

Agitation speed

1•

FEATURES AND SPECIFICATIONS

1.8 Software parameters for Wash protocols (continued)

1.8.3 Parameter's definitions (continued)

Identification and range of parameters related to the location of the needle in relation to the microplate well.

Graduations in this sketch are only for better understanding. Set these parameters in real conditions with real plate; well dimension may vary from one plate type to another one.




Fig. 8

## 2.

## INSTALLATION

### 2.1 Unpacking

- Carefully unpack accessories and washer.

**IMPORTANT:**

Do not lift up the washer  
by pulling on the foam frames!

- First remove accessories at the top and beside the washer.
- Then, hold the washer (not the foams) and lift up.
- Remove the foams, remove plastic bag and place washer on a stable, horizontal surface.
- Install bottles beside the washer as shown in user manual on §1.1 for PW40 or § 1.4 for PW41.




Fig.9




Fig.10

## 2.1 Unpacking (continued)




Fig.11

- **In washer working area (Fig. 11, 12 & 13):**
- Remove adhesive tape blocking aerosol cover.
- Remove adhesive tape blocking plate carrier to bottom plate (fig. 11).

- **Pull and push plate carrier to check that it moves freely.**

### • Manifold

- Remove the tie blocking the manifold on the left side (fig. 12) : push the tie to the left. Remove the tie from the vertical shaft.
- Turn blocking device as shown on fig. 12 and ensure to not hinder the manifold Up/Down motion.

### • On rear panel :

Connect tubes on rear panel (match colours).

blue	=	Wash 1
gray	=	Wash 2
green	=	Wash 3
transparent	=	Rinse
red	=	to Waste bottle (liquid)
Yellow	=	Vacuum to Waste bottle (air)




Fig. 12

### • Vacuum pump protection :

Cut yellow tube at halfway between WASTE bottle and washer and place hydrophobic filter (supplied in maintenance kit).

### • Connect power cord.

Fuse installed : 2,5 A Slow Blow.

Plug Power Cord to mains socket; washer is equipped with universal input switch mode power supply (85 VAC to 264 VAC; 47 Hz to 440 Hz).




Fig.13

### • Turn the washer ON (ON/OFF switch on rear panel)

Whereas plate carrier and manifold are moving to their home position, following message appears on display :

VERSION : 2.XXX.X.Y  
PW41

After initialization :  
SELECT : RUN  
IN YES OUT





- **Fill one RINSE bottle with demineralized water and prime hydraulic circuit with at least 5 RINSE operations** (see RINSE mode, §3.5).

- The washer is ready to use.

## 3 OPERATING INSTRUCTIONS

### 3.1 Run Mode

#### 3.1.1 How to use the Keyboard?




Press "T" or "I" Keysto  
scroll main menu  
See para 3.6.1 Main  
Menu, page 3/8.

Press "OUT" Key to set  
or remove microplate :  
plate carrier moves out  
and is locked for 10 se-  
conds to allow micropla-  
te to be placed in carrier.

Press "YES" Key to  
enter selected mode

Press "IN" key to move  
Plate Carrier In



The function of the key is displayed  
on the screen.  
Depending on displayed menu,  
• some keys are activated,  
• some others have no function and  
are therefore deactivated

### 3.2 Microplate loading

Place the 96-well microplate on the plate carrier with "A 1" corner resting against metal spring on the right side (fig. 13).





Fig. 13

3•

OPERATING INSTRUCTIONS

3.3 Select one KIT



(B) : See §3.4.1

\* : may appear as  
"CONNECT THE WASH R9".




3•

OPERATING INSTRUCTIONS

3.3 Select one KIT (continued)

ERR:21 with PW41 :



\* Error 21 indicates that Wash buffer W9 (or R9) requested in "KIT 75" is not connected.  
Connect Wash buffer W9 (or R9) at input 1, 2 or 3.  
Then enter washer configuration to modify input assignment accordingly (see §3.6.2).

\* : may appear as  
"ERR:21 WASH R9 NOT  
CONNECTED".


3»

## OPERATING INSTRUCTIONS

### 3.4 WASH interruption or Power Failure


#### 3.4.1 WASH interruption

When the "ESC" key is pressed whereas a kit is running, the following message is displayed :



#### 3.4.2 POWER Failure

When a power failure occurs whereas a kit is running, the following message will be displayed when power is back :




3\*

## OPERATING INSTRUCTIONS

### 3.5 RINSE

#### 3.5.1 RINSE (with PW40)



#### 3.5.2 Please RINSE (with PW40)

The message

PLEASE RINSE  
YES

is displayed to incite user to perform a RINSE operation when the washer has not been used during 10 minutes after a WASH operation.  
Press "YES" Key to go back to main menu and launch a RINSE.

#### **IMPORTANT**


The washer should never be switched off with WASH solution inside the tubing and manifold.  
Solid crystals from dried WASH solution could clog the dispense needles and clamp electrovalve tube.

3\*

OPERATING INSTRUCTIONS

3.5 RINSE (continued)

3.5.3 RINSE (with PW41)



3.5.4 Automatic Rinsing (with PW41)

One automatic rinsing sequence will take place as soon as the time set in "RINSE PARAM." (see §3.6.2. Washer's configuration) has elapsed after the last wash.

**Note :**

When the washer won't be used for several days, it is recommended to connect all inputs to rinse bottle and rinse every input as described above in order to clean the pinch valve tube of every input. See also §4.1.2 Storage condition.

**IMPORTANT**

The washer should never be switched off with WASH solution inside tubing and manifold.


Solid crystals from dried WASH solution could clog the dispense needles and clamp electrovalve tube.

3»

## OPERATING INSTRUCTIONS

### 3.6 Configuration

#### 3.6.1 Main menu

 SELECT : RUN YES 0 UT  
T / IN YES 0 UT  
RUN  
See §3.1 Run Mode

SELECT : PRIME/RINSE  
T / IN YES 0 UT  
PRIME/RINSE  
See §3.5 RINSE

SELECT : DISINFECTION  
T / IN YES 0 UT  
DISINFECTION  
See §4.2.2 DISINFECTION

SELECT : CONFIGURATION  
T 4 IN YES 0 UT  
CONFIGURATION  
See §3.6.2 Washer's configuration

SELECT : SERVICE/VER  
T / IN YES 0 UT  
SERVICE/VER  
See §4.3.1 Hardware and Software  
version




3<

## OPERATING INSTRUCTIONS

### 3.6 Configuration (continued)

#### 3.6.2 Washer's configuration




\* : may appear as

"INPUT 1 SOL.: WASHR1".

### 3.6 Configuration (continued)

#### 3.6.2 Washer's configuration (continued)

(continuation from page 3/9)



(3otopage3/10


3•

OPERATING INSTRUCTIONS

3.6 Configuration (continued)

3.6.2 Washer's configuration (continued)

(continuation from page 3/10)



(A) : Back to main menu  
(See §3.3).



3•

OPERATING INSTRUCTIONS

3.7 PROGRAMMING

3.7.1 Kit Structure

**KIT PROGRAMME**

Parameters depending  
on Kit (reagents)

**PLATE PROGRAMME**

Parameters depending  
on Plate dimensions




Fig. 15

As seen on §1.8.3 Parameter's definitions (KIT and PLATE), parameters directly depending on KIT (such as dispensing, aspirating times, liquids, etc... are called "KIT Parameters".  
Parameters depending on microplate dimensions (mainly POSITION and SPEED parameters) are called "PLATE parameters".

At programming, each kit is mandatory linked to one plate with the "Plate number" parameter (Fig. 15).

This implies that the PLATE must be already existing or must be programmed **prior** the associated KIT is programmed.

Up to 10 PLATES can be programmed.

Up to 75 KITS can be programmed.

3•

OPERATING INSTRUCTIONS

3.7 PROGRAMMING (continued)

3.7.2 Access to PROGRAMMING mode

SELECT :	RUN	YES	OUT
t	I N		



Enter programming mode : press simultaneously  
"IN" and "OUT" keys.

PRG :	ADD	YES
T	i	

### 3.7 PROGRAMMING (continued)

#### 3.7.3 Programming Synopsis




Fig.16

3•

OPERATING INSTRUCTIONS

3.7 PROGRAMMING (continued)

3.7.4ADD 3.7.4.1 ADD: KIT

The kit is added at the end of the existing kit list

- Use "T" and "/" keys to scroll the available menu.
- Press "YES" key to initiate the displayed selection.
- Press "ESC" key to escape and go back to former prompt.

NAME : KIT Nr2  
t 4 -> YES ESC

Enter kit name : select character with "T" and "4" keys, then press "-" key to move to next character, (up to 15 characters).

PLATE : PLATE01  
t j YES ESC

Select plate to be associated to kit.

MANIFOLD : 8  
T j YES ESC

Select manifold type:  
PW40 : 8- or 12-channels  
PW41 : 8-, 2\*8-, 3\*8- or 12-, 2\*12-channels.

STRIP : 1 X 3 4 5 6  
T 4- -> YES ESC

Omit strip : place cursor on strip with "-" key, then press "T" or "/" key to have strip omitted ("X"). In RUN mode, the omitted strips will never be processed.

STRIP : 7 8 9 10 X 12  
t j -> YES ESC

METHOD : WASH + ASP.  
T 4 YES ESC

Select first method of kit.  
See method parameter selection on §3.7.4.3.

t j Y T A \* 1

Answer "NO" in order to add another method.

END OF KIT : NO  
t j YES ESC

Set "Method interval" with "T" or "j" and "-" keys.  
(from 0s up to 59mn 0s).

MET . INTER : 0MN45 S  
T YES ESC

METHOD : DISPENSING  
T j YES ESC

t j YES ESC

Answer "YES" when no more methods have to be added.

END OF KIT : YES  
T / YES ESC


Go to page 3/16)

### 3.7 PROGRAMMING (continued)


#### 3.7.4 ADD (continued)

##### 3.7.4.1 ADD : KIT (continued)

(continuation from page 3/15)



#### 3.7.4.2 ERRORS on KIT programming



### 3.7 PROGRAMMING (continued)

#### 3.7.4ADD (continued)

3.7.4.3 ADD : METHOD The method is added after the last method of the selected kit

```
ADD      : METHOD
↑
T        : YES ESC
```

•Use "T" and "/" keys to scroll the available menu.  
•Press "YES" key to initiate the displayed selection.  
•Press "ESC" key to escape and go back to former prompt.

```
KIT : KIT N123
t   4, YES ESC
```

Select kit to which you want to add a method with "t" and "/" keys

```
MET. INTER: OMN OS
T          : YES ESC
```

Set "Method interval" between last existing method and added method with "T" or "4" and "-" keys, (from Omn O&up to 59mn Os).

```
METHOD : WASH + ASP.
T        : YES ESC
```

Select method to be added

```
MODE : PLATE
t     : YES ESC
```

Set kit parameters associated to selected method.  
See table next page and kit parameter definition §1.8.3.

```
CROSW. ASP.: NO
t         : YES ESC
```

```
ASP. TIME : 0.5S
t         : YES ESC
```

```
VOLUME : 750ul
t      4 : YES ESC
```

```
OVERFLOW : 3.0mm
t        : YES ESC
```

```
LIQUID : WASH W9
t      : YES ESC
```

```
FLOW : -3
t    : YES ESC
```

```
Nr OF CYCLES : 2
t           : YES ESC
```

```
SOAKING : 1MN OS
t       : YES ESC
```

Only the first cycle is played "n" times (here twice) for a "two-cycle" method such as Wash + Asp(nW+A).

```
SELECT : RUN
t      : IN YES OUT
```

### 3.7 PROGRAMMING (continued)

#### 3.7.4ADD (continued)

#### 3.7.4.3 ADD : METHOD (continued)

Summary of kit parameters associated to method.

Kit parameter	Value	See def.	nW+A	nW+a	nW	nA	nB	nW	na	nAg	
METHOD :			WASH + ASP.	WASH+BOT.ASP.	BOT.WASH	B.WASH	ASP.	DI=SPENSING	BOTTOM WASH.	BOTTOM ASP.	AGITATION
Method Code :											
Kit parameter	Value	See def.									
MODE :	STRIP or PLATE	p. 1/10	.	.	.	.	.	.	.	.	
CROS. ASP. :	YES or NO	p. 1/12	.	.	.	.	.	.	.	.	
ASP. TIME :	0,1 to 9,9 sec	p. 1/12	.	.	.	.	.	.	.	.	
VOLUME :	50 to 3000 lL	p. 1/13	.	.	.	.	.	.	.	.	
OVERFLOW :	1,0 to 12,9ram	p. 1/13	.	.	.	.	.	.	.	.	
LIQUID :	WASH W1 to WASH W9	p. 1/13	.	.	.	.	.	.	.	.	
FLOW :	-5 to +5	p. 1/13	.	.	.	.	.	.	.	.	
BOT. WASH NUMBER :	1 or 2	p. 1/16	.	.	.	.	.	.	.	.	
BOTTOM TIME :	0,1 to 9,9. sec	p. 1/16	.	.	.	.	.	.	.	.	
BOT. ASP. NUMBER :	1 or 2	p. 1/17	.	.	.	.	.	.	.	.	
SHAKE TIME :	0,1 to 59,9 sec	p. 1/19	.	.	.	.	.	.	.	.	
NT OF CYCLES :	1 to 9	p. 1/11	.	.	.	.	.	.	.	.	
SOAKING :	0 to 59 mn in PLATE mode 0 to 9,9 sec i STRIP mode	p. 1/11	.	.	.	.	.	.	.	.	

#### 3.7.4.4 ERROR on METHOD programming

```

ERR:25 NOT ENOUGH
MEMORY YES
SELECT:RUN IN YES OUT
t i IN YES OUT
  
```


Memory is full because one or several kits are containing a large number of methods.  
You must first delete an existing kit or existing method prior programming a new one.

### 3.7 PROGRAMMING (continued)

#### 3.7.4 ADD (continued)

#### 3.7.4.5 ADD : PLATE

The plate is added at the end of the existing plate list





### 3.7 PROGRAMMING (continued)

#### 3.7.4ADD (continued)

#### 3.7.4.5 ADD : PLATE (continued)


(continuation from page 3/19)

BOT.VERT.POS. : 11.8mm  
t I YESSC



Bottom (aspiration) vertical position.  
See definition on para 1.8.3.

B.W.VERT.POS. : 11.8mm  
t I YESSC



Bottom Wash Vertical Position.  
See definition on para 1.8.3.


HORIZONTAL SPEED : 6  
t I YESSC

Horizontal speed :  
In / Out plate carrier speed during non-work movements (movement to home position, return movement from last to first strip, etc.).  
0 to 9, step 1.  
(speed mini = 0; speed max=9)

VERTICAL SPEED : 6  
t I YESSC

Vertical speed :  
Up / Down manifold speed during non-work movements (movement to home position, etc.).  
0 to 9, step 1.  
(speed mini = 0; speed max=9)

ASP. DOWNW. SPEED : 6  
t I YESSC



Aspiration downward speed.  
See definition on para 1.8.3.


(J30 to page 3/21P)

### 3.7 PROGRAMMING (continued)

#### 3.7.4 ADD (continued)

#### 3.7.4.5 ADD : PLATE (continued)

(continuation from page 3/20)



3•

OPERATING INSTRUCTIONS

3.7 PROGRAMMING (continued)

3.7.4ADD (continued)

3.7.4.6 ERRORS on PLATE programming

```
ERR:26  ALREADY  
10 PLATES!!  YES
```

10 plates are already programmed! You must first delete an existing plate prior programming a new one.

```
SELECT:RUN  
↑ 4 IN YES OUT
```

```
ERR:33  PLATE NAME  
IMPOSSIBLE  YES
```

You have entered no name or the plate name you have entered is already existing. You are prompted to enter a new name to continue plate programming.

```
NAME :  
↑ / → YES ESC
```

3•


OPERATING INSTRUCTIONS

3.7 PROGRAMMING (continued)

3.7.5 INSERT

3.7.5.1 INSERT: KIT

The new kit is inserted **before** another kit



3.7.5.2 TEMPORARY KIT


```

RUN : TEMPORARY KIT
+-----+
|         |
| 4-     |
| YES    |
| ESC    |
+-----+
    
```

Power failure or erroneous OFF switching during KIT insertion will cause the kit to be lost unless the first method has been completely programmed. In this case the completely programmed methods can be found in a kit named "TEMPORARY KIT". This kit is always inserted at the first place of the kit list. It is then possible to COPY this Temporary Kit and write its correct name, then continue with programming of next methods and kit repetition parameters. Once saved under its correct name, it is highly recommended to DELETE the TEMPORARY KIT.

3.7.5.3 INSERT: METHOD

The new method is inserted **before** another method



**3.**


**OPERATING INSTRUCTIONS**

**3.7 PROGRAMMING (continued)**

**3.7.5 INSERT (continued)**

**3.7.5.4 INSERT : PLATE**

The new plate is inserted **before** another plate




3\*

OPERATING INSTRUCTIONS

3.7 PROGRAMMING (continued)

3.7.6 DELETE


3.7.6.1 DELETE: KIT



Select kit to be deleted.

Press "YES" key to confirm kit is deleted or press "ESC" key to escape.

3.7.6.2 DELETE: METHOD



Select kit to which you want to delete a method.

Select method to be deleted.

Press "YES" key to confirm method is deleted or press "ESC" key to escape.

**Note:** deleting the method of a "one-method" kit is equivalent to delete this kit.


3•

OPERATING INSTRUCTIONS


3.7 PROGRAMMING (continued)

3.7.6 DELETE (continued)

3.7.6.3 DELETE : PLATE




3.7.6.4 ERROR on PLATE deletion



### 3.7 PROGRAMMING (continued)

#### 3.7.7 EDIT

##### 3.7.7.1 EDIT: KIT






### 3.7 PROGRAMMING (continued)

#### 3.7.7 EDIT (continued)

##### 3.7.7.1 EDIT : KIT (continued)




3\*

OPERATING INSTRUCTIONS

3.7 PROGRAMMING (continued)

3.7.7 EDIT (continued)


3.7.7.2 EDIT PLATE




### 3.7 PROGRAMMING (continued)

#### 3.7.8 COPY

##### 3.7.8.1 COPY : KIT



##### 3.7.8.2 ERROR on KIT copy




3•

OPERATING INSTRUCTIONS


3.7 PROGRAMMING (continued)

3.7.8 COPY (continued)

3.7.8.3 COPY: PLATE



3.7.8.4 ERROR on PLATE copy




### 3.7 PROGRAMMING (continued)

#### 3.7.9 Plate Test

##### Step-by-Step RUN

The programmed kit can be used in a step-by-step mode i.e. manifold and plate carrier are stopped at each programmed positions (Aspiration horizontal position, Aspiration vertical position, Overflow position, etc..) in order to check if these positions are correct. If not correct, user must edit (see §3.7.7) the tested kit or its assigned plate in order to modify position parameter(s).




### 3.7 PROGRAMMING (continued)

#### 3.7.10 RAM initialization


Every programmed kit and programmed plate are saved on a RAM (Random Access Memory).  
The contents of this MEMORY can be erased by a "RAM initialization" thus deleting the whole set of kits and plates that have been programmed by the user.  
Instrument software and DEMO kits and plates are not deleted as they are saved in an EPPROM not affected by RAM initialization.

#### How to perform the RAM initialization?

- Turn the unit ON while keeping those keys pressed :



Following message is displayed :





## 4\* MAINTENANCE & TROUBLESHOOTING

### 4.1 Routine maintenance Procedures

#### 4.1.1 Daily Maintenance

##### **Prior starting to wash :**

It is recommended to pour an appropriate disinfection solution in the bottom of the WASTE bottle prior connecting it to the washer. Thus, the infected liquid coming from the washer and going into the WASTE bottle will be deactivated.

##### **IMPORTANT: Avoid bleach**

Fill the rinse bottle with deionized water. Connect it to the washer. Launch a Rinsing sequence and check that pinch valve tube and

manifold are not clogged.

Otherwise, change pinch valve tube, clean manifold (see § 4.2.1 hereafter).

##### **When washing is terminated :**

Rinse the washer . Use only deionized water for rinsing. Turn the washer off. Tubes and manifold will be kept wet with water that will prevent them from clogging.

Empty WASH bottle, rinse it thoroughly with deionized water.

#### 4.1.2 Washer storage conditions

The washer should never be stored with Wash solution inside dispensing tubes and manifold otherwise the liquid pump, the pinch valve tube and the manifold might be clamped if completely dry, thus preventing any

further dispensing of liquid.

In this case priming must be made with the help of a syringe filled with deionized water connected at the Rinse input of the washer.



4-

MAINTENANCE & TROUBLESHOOTING

## 4.2 Cleaning and disinfection

### 4.2.1 Cleaning the manifold

#### **The manifold is infectious.**

#### **Disinfect it prior any further cleaning.**

See disinfection procedure on next § 4.2.2.

#### **Aspirate needles and chamber**

- Remove the manifold.
- Using a duster or a non scratching utensil, thoroughly clean the outside of the dispense and aspirate needles.
- Remove the silicon top seal of the manifold.
- Using the big pin supplied in maintenance kit, clean the big aspirating needles.
- Rinse out the aspirate chamber with deionized or distilled water. Thoroughly clean the walls of the aspirating cavity.
- Look through the aspirate needles to ensure that cleaning is adequate.
- Put the silicon top seal back in its place.

#### **Dispense needle and chamber**

- Push the silicon cap from the chamber of the manifold using a stem or a screwdriver entered from the opposite end of the cylindrical chamber.
- Using one of the small pin supplied in maintenance kit, clean the small dispense needle.
- Using the cylindrical brush supplied in maintenance kit, thoroughly clean the small dispense chamber. Thoroughly rinse dispense chamber with deionized water. Ensure that all impurities are removed from the manifold.
- Plug dispense chamber of the manifold with the silicon plug. Place it as close as possible to the first needle however without blocking it. Extra silicon caps are supplied in the maintenance kit.
- Remount the manifold, replace the tubing and reprime.


4•

MAINTENANCE & TROUBLESHOOTING

4.2 Cleaning and disinfection (continued)

4.2.2 Disinfection

The disinfection program applies to internal surfaces of tubing, connectors, manifold and manifold needles.



**IMPORTANT:**


- Avoid bleach for disinfecting the washer as a thorough rinsing of the tubes will not prevent traces of bleach to remain inside tubing and manifold. This can be enough to corrupt reagents of next washing.
- Disinfect the washer when it has not been used for more than one day. When used daily, disinfect it at least once a week.
- Disinfect the washer prior to any servicing or handling action, shipment or change of site.

4•

MAINTENANCE & TROUBLESHOOTING

### 4.3 Troubleshooting

#### 4.3.1 Hardware and software version



## Appendix

### Accessories and spare part list

Part Number	Designation
82386	PW4X/LP35 WASTE BOTTLE 2L W/O CAP
82387	PW4X/LP35 COMPLETE WASTE BOTTLE 2L
82388	PW4X/LP35 WASTE 2L RUBBER CAP+TUBE
87042	PW4x COMPLETE WASH BOTTLE
87049	PW4x RINSE 2L SCREW CAP+TUBE
87053	PW4x WASH BOTTLE WITHOUT CAP
87054	PW40 WASH 2L SCREW CAP+TUBE
87058	PW4x RINSE BOTTLE WITHOUT CAP
87059	PW4x COMPLETE RINSE BOTTLE
87063	PW41 COMPLETE WASTE BOTTLE+DETECTOR
87086	PW41 SCREW CAP WASH BOTTLE+LEVEL DETECTION
87087	PW41 COMPLETE WASH BOTTLE+TUBE+LEVEL DETECTION
87088	PW41 COMPLETE WASTE BOTTLE 8L+LEVEL DETECTION
87067	PW41 WASTE 2L RUBBER CAP+LONG TUBE
87040	PW4x 12 CHANNEL MANIFOLD
87041	PW4x 8 CHANNEL MANIFOLD
87045	PW4x MANIFOLD CONNECTOR
87046	PW4x MANIFOLD CONNECTOR RINGS
87051	PW4x RUBBER SEAL FOR MANIFOLD
87055	PW41 2x12 CHANNEL MANIFOLD
87056	PW41 2x8 CHANNEL MANIFOLD
87057	PW41 3x8 CHANNEL MANIFOLD
87085	PW41 MANIFOLD CONNECTOR VACUUM SENSING
85509	PW40/LP35 COMPLETE TRAP VIAL
87082	PW40/LP35 TRAP VIAL CAP+TUBE
87083	PW40/LP35 TRAP VIAL WITHOUT CAP
87782	PW41 MOTOR SELECTOR TUBING SET
87094	PW4X/LP35 WASTE TUBING RED & YELLOW
87043	PW4x DISPENSING TUBE BLUE
87044	PW4x EXTERNAL TUBING SET
87047	PW4x PINCH VALVE TUBE
87060	PW4x DISPENSING TUBE RINSE TRANSPARENT
87061	PW41 DISPENSING TUBE WASH2 WHITE
87062	PW41 DISPENSING TUBE WASH3 GREEN
87064	PW41 4 WAY VALVE TUBE
87090	PW41 LINK TUBE 4WAY VALVEUNIT
87091	PW41 SET TUBE MULTI WASH SELECTOR
87048	PW4x REMOVABLE uPLATE CARRIER
87081	PW4x HOLDER BOTTLE ON TOP WASHER
87089	PW41 HOLDER 2L BOTTLES BESIDE WASHER
87004	PW41 FUSE 2-5AT 5X20 250V
85274	PW4X/LP35 HYDROPHOBIC FILTER 11mm
87092	PW41 RS232 CORD
87052	PW40 STANDARD MAINTENANCE KIT
87093	PW MAINTENANCE KIT
87050	PW USER MANUAL F/E

**PW40 Kit sheet**  
**eioidas**

xippeddy

EKJdAS : NOISIA3H  
96/60ZL

DATE :

Kit name NAME	PLATE	MANIFOLD	strip omitted STRIP	MODE	METHOD	CROSSWISE ASPIRATION CROSW. ASP	Aspiration time ASP. TIME	Dispensed volume VOLUME	Overflow position OVERFLOW	Selected liquid LIQUID	Flow compensation FLOW	Bottom Wash Number BOT. WASH NUMBER	Bottom Wash Time BOT. WASH TIME	Bottom Aspiration Nbr BOT. ASP. NUMBER	SHAKE TIME	N° OF CYCLES	SOAKING	Method Interval MET. INTER.	N° OF KITS	Kit Interval KIT INTER	Date		
																					Year	Month	
T01 N1 500uL	Flat 01	8	1	a	Plate	No	0,3s	500µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N1 800uL	Flat 01	8	1	w	Plate	No	0,3s	800µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N1 800uL C	Curved02	8	1	W+a	Plate	No	0,3s	800µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N2 500uL	Flat 01	8	1	w	Plate	No	0,3s	500µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N2 800uL	Flat 01	8	1	a	Plate	No	0,3s	800µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N2 S2 500uL	Flat 01	8	1	w	Plate	No	0,3s	500µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N2 S2 800uL	Flat 01	8	1	w	Plate	No	0,3s	800µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N2 S3 800uL	Flat 01	8	1	w	Plate	No	0,3s	800µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N3 800uL C	Curved02	8	1	W+a	Plate	No	0,3s	800µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N3 800uL	Flat 01	8	1	w	Plate	No	0,3s	800µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N3 1000uL	Flat 01	8	1	w	Plate	No	0,3s	1000µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N3 800uL C	Curved02	8	1	W+a	Plate	No	0,3s	800µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N4 800uL	Flat 01	8	1	w	Plate	No	0,3s	800µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N4 1000uL	Flat 01	8	1	w	Plate	No	0,3s	1000µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N4 800uL C	Curved02	8	1	W+a	Plate	No	0,3s	800µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N4 800uL	Flat 01	8	1	w	Plate	No	0,3s	800µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N5 800uL	Flat 01	8	1	w	Plate	No	0,3s	800µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N5 S10s 800	Flat 01	8	1	w	Plate	No	0,3s	800µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N5 800uL C	Curved02	8	1	W+a	Plate	No	0,3s	800µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N6 800uL	Flat 01	8	1	w	Plate	No	0,3s	800µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1
T01 N6 500uL	Flat 01	8	1	w	Plate	No	0,3s	500µL	2,5	W1	0	1	1	1	1	1	1	0s	1	1	1	1	1

Kit name NAME	Flat 01	8	MANIFOLD	strip oriented STRIP	1	MODE		Crosswise Aspiration CROSW, ASP	Aspiration time ASP. TIME	Dispensed volume VOLUME	Overflow position OVERFLOW	Selected liquid LIQUID	Flow compensation FLOW	Bottom Wash Number BOT. WASH NUMBER	Bottom Wash Time BOTTOM TIME	Bottom Aspiration Nbr BOT. ASP. NUMBER	SHAKE TIME	SOAKING	Method Interval MET. INTER.	Nr. OF KITS	Kit Interval KITINTER	
						W	Plate															
• T01 N6 800µL	Flat 01	8	MANIFOLD	strip oriented STRIP	1	W	Plate	No	0,3s	800µL	2,5	W1	0	1	1	1	1	1	6	1	1	1
						a	Plate															
• T02 N3B1,0s 800	Flat 01	8	MANIFOLD	strip oriented STRIP	1	W	Plate	No	0,5s	800µL	2,5	W1	0	2	1,0s	1	1	1	3	1	1	1
						a	Plate															
• T02 NSB0,5s 800	Flat 01	8	MANIFOLD	strip oriented STRIP	1	W	Plate	Yes	0,5s	800µL	2,5	W1	0	2	0,4s	1	1	1	5	1	1	1
						a	Plate															
• T03 N3 800µL	Flat 01	8	MANIFOLD	strip oriented STRIP	1	W+A	Strip	Yes	0,5s	800µL	2,5	W1	0	1	0,4s	1	1	1	3	1	1	1
						a	Plate															
• T04 N3 800µL	Flat 01	8	MANIFOLD	strip oriented STRIP	1	W+A	Strip	Yes	0,5s	800µL	2,5	W1	0	2	0,4s	1	1	1	3	1	1	1
						a	Plate															
• T05 N3 800µL	Flat 01	8	MANIFOLD	strip oriented STRIP	1	W+A	Strip	Yes	0,5s	800µL	2,5	W1	0	2	0,4s	1	1	1	3	1	1	1
						a	Plate															
• T06 N3 800µL	Flat 01	8	MANIFOLD	strip oriented STRIP	1	W	Strip	Yes	0,5s	800µL	2,5	W1	0	1	1	1	1	1	3	1	1	1
						a	Plate															
• T07 N5 800µL	Flat 01	8	MANIFOLD	strip oriented STRIP	1	W+A	Strip	Yes	0,5s	800µL	2,5	W1	0	1	1	1	1	1	5	1	1	1
						a	Plate															
• T07 N5 S1s 800	Flat 01	8	MANIFOLD	strip oriented STRIP	1	W+A	Strip	Yes	0,5s	800µL	2,5	W1	0	1	1	1	1	1	5	1s	1	1
						a	Plate															
• T08 N3 800µL	Flat 01	8	MANIFOLD	strip oriented STRIP	1	W	Plate	No	0,5s	800µL	2,5	W1	0	1	1	1	1	1	3	1	1	1
						A	Plate															
• T08 N4 800µL	Flat 01	8	MANIFOLD	strip oriented STRIP	1	W	Plate	No	0,5s	800µL	2,5	W1	0	1	1	1	1	1	4	1	1	1
						A	Plate															
• T09 N3 S4s 500	Flat 04	8	MANIFOLD	strip oriented STRIP	1	W+A	Plate	Yes	0,3s	500µL	2,7	W1	0	1	0,8s	1	1	1	3	40s	1	1
						a	Plate															
• T09 N6 S4s 500	Flat 04	8	MANIFOLD	strip oriented STRIP	1	W+A	Plate	Yes	0,3s	500µL	2,7	W1	0	1	0,8s	1	1	1	5	40s	1	1
						a	Plate															
• T10 N2 800µL	Flat 01	8	MANIFOLD	strip oriented STRIP	1	W	Plate	No	0,5s	800µL	2,5	W1	0	1	0,4s	1	1	1	2	1	1	1
						a	Plate															
• T10 N3 800µL	Flat 01	8	MANIFOLD	strip oriented STRIP	1	W	Plate	No	0,5s	800µL	2,5	W1	0	1	0,4s	1	1	1	3	1	1	1
						a	Plate															
• T10 N5 S1s 800	Flat 01	8	MANIFOLD	strip oriented STRIP	1	W	Plate	No	0,5s	800µL	2,5	W1	0	1	0,4s	1	1	1	5	1	1	1
						a	Plate															
• T10 N6 800µL	Flat 01	8	MANIFOLD	strip oriented STRIP	1	W	Plate	No	0,5s	800µL	2,5	W1	0	1	0,4s	1	1	1	6	1	1	1
						a	Plate															

**PW40 Kit sheet  
SDP013**

SIDODS :NOISIA3U  
96/60/A '31VO

Appendix

Kit name NAME	Flat 01	8	MANIFOLD	strip omitted STRIP	MODE	METHOD		CROSSWISE ASPIRATION CROSW. ASP		Aspiration time ASP. TIME	Dispensed volume VOLUME	Overflow position OVERFLOW	Selected liquid LIQUID	Flow compensation FLOW	Bottom Wash Number BOT. WASH NUMBER	Bottom Wash Time BOTTOM TIME	Bottom Aspiration Nbr BOT. ASP. NUMBER	SHAKE TIME	SOAKING	Method Interval MET. INTER.	Nr. OF CYCLES	Np. OF KITS	Kit Interval KIT INTER	
						w	A	No	Yes															
T11 N3 800µL	Flat 01	8				w	A	No	Yes	0.3s	800µL	2.5	W1		0						3	1	1	1
T12 N3 800µL	Flat 01	8				w	a	No	Yes	0.5s	800µL	2.5	W1		0	2					3	0s	1	1
T13 N3 800µL	Flat 01	8				w	Plate	No	Yes	0.5s	800µL	2.5	W1		0						3	0s	1	1
T14 N3 800µL	Flat 01	8				w	Plate	No	Yes	0.5s	800µL	2.5	W1		0						3	0s	1	1
T15 N1 Stimm 800	Flat 01	8				w	A	No	Yes	0.5s	800µL	2.5	W1		0						1	1mn	1	1
T16 N1 Stimm 800	Flat 01	8				w	Plate	No	Yes	0.3s	800µL	2.5	W1		0	1					1	1mn	1	1
T17 N1 Stimm 800	Flat 01	8				w	Plate	No	Yes	0.3s	800µL	2.5	W1		0	1					1	1mn	1	1
T18 N1 Stimm 800	Flat 01	8				w	Plate	No	Yes	0.3s	800µL	2.5	W1		0	2					1	1mn	1	1
T19 N1 Stimm 800	Flat 01	8				w	Plate	No	Yes	0.3s	800µL	2.5	W1		0	2					1	1mn	1	1
T20 N3 800	Flat 01	8				w	Plate	No	Yes	0.3s	800µL	2.5	W1		0						3	0s	1	1
T21 N1 Stimm 800	Flat 01	8				w	Plate	No	Yes	0.3s	800µL	2.5	W1		0						1	1mn	1	1
COCKIT	Flat 01	8				w	Ag	Yes	Yes	0.5s	50µL	10	W1		0						1	0s	1	1

xipusdov

310D05 :NOISIA3H  
96/60/ZISLODas  
kit sheet  
Kit Interval  
KIT INTER

		10	10	10	10	10	10	10	10
	Flat04	Flat04	Flat04	Flat04	Flat04	Flat04	Flat04	Flat04	Flat04
	Curved	Curved	Curved	Curved	Curved	Curved	Curved	Curved	Curved
	Flat01	Flat01	Flat01	Flat01	Flat01	Flat01	Flat01	Flat01	Flat01
	Plate name NAME								
	Bottom shape BOT. SHAPE	Flat	Curved	Flat	Flat	Flat	Flat	Flat	Flat
	Aspiration horizontal position ASP. HOR. POS.	1.4	1	1.7	0.5	0.5	14.5	9.5	9.5
	Plate centering CENTERING	0.3	0.3	0.3	120	120	14.5	9.5	9.5
	Aspiration vertical position ASP. VERT. POS.	18	18	18	9.5	9.5	14.5	9.5	9.5
	Bottom aspiration vertical position BOT. VERT. POS.	9.5	9.5	9.5	9.5	9.5	9.5	9.5	9.5
	Bottom Wash vertical position B.W. VERT. POS.	10	10	10	9.5	9.5	9.5	9.5	9.5
	HORIZONTAL SPEED	0	0	0	20	20	20	20	20
	VERTICAL SPEED	0	0	0	0	0	0	0	0
	Aspiration downward speed	0	0	0	4	4	4	4	4
	ASP. DOWNW. SPEED	0	0	0	10	10	10	10	10
	Dispensing upward speed	10	10	10	10	10	10	10	10
	DISP. UPW. SPEED	10	10	10	10	10	10	10	10
	Bottom Wash downward speed	0	0	0	4	4	4	4	4
	BOT. DOWNW. SPEED	0	0	0	10	10	10	10	10
	Bottom Wash upward speed	10	10	10	10	10	10	10	10
	BOT. UPWARD SPEED	10	10	10	10	10	10	10	10
	SHAKING AMPLITUDE	1	1	1	1	1	1	1	1
	SHAKING SPEED	10	10	10	10	10	10	10	10

elodcis  
plate sheet  
id OfMd

St-Opds : NOISIA3U  
DATE : 9/6/07U

xipueddv







appendix

PW41 Kit sheet  
 SDP014

96/60/81. :31VQ  
 Hddp : N0SISIAJU

Kit name NAME	PLATE	MANIFOLD	strip omitted STRIP	METHOD		MODE		Crosswise Aspiration CROSW. ASP		Aspiration time ASP. TIME		Dispensed volume VOLUME		Overflow position OVERFLOW		Selected liquid LIQUID		Flow compensation FLOW		Bottom Wash Number BOT. WASH NUMBER		Bottom Wash Time BOTTOM TIME		Bottom Aspiration Nbr BOT. ASP. NUMBER		SHAKE TIME		Nr OF CYCLES		SOAKING		Method Interval MET. INTER.		Nr OF KITS		Kit Interval KIT INTER		
				W	A	W	A	No	Yes	0.5s	800µL	2.5	WT	5	2	0.4s	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1		
T11 N3 800µL	Flat 01	16	1	W	Plate	W	Plate	No	Yes	0.5s	800µL	2.5	WT	5	2	0.4s	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1		
T12 N3 800µL	Flat 01	16	1	W	Plate	W	Plate	No	Yes	0.5s	800µL	2.5	WT	5	2	0.4s	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s
T13 N3 800µL	Flat 01	16	1	W	Plate	W	Plate	No	Yes	0.5s	800µL	2.5	WT	5	2	0.4s	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s
T14 N3 800µL	Flat 01	16	1	W	Plate	W	Plate	3	Yes	0.5s	800µL	2.5	WT	5	2	0.4s	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s
T15 N1 S1m 800	Flat 01	16	1	*	Plate	*	Plate	No	Yes	0.5s	800µL	2.5	WT	5	2	0.4s	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s
T16 N1 S1m 800	Flat 01	16	1	S	Plate	S	Plate	No	Yes	0.3s	800µL	2.5	WT	5	2	0.4s	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s
T17 N1 S1m 800	Flat 01	16	1	S	Plate	S	Plate	1.0	Yes	0.3s	800µL	2.5	WT	5	2	0.4s	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s
T18 N1 S1m 800	Flat 01	16	1	W	Plate	W	Plate	No	Yes	0.3s	800µL	2.5	WT	5	2	0.4s	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s
T19 N1 S1m 800	Flat 01	16	1	S	Plate	S	Plate	No	Yes	0.3s	800µL	2.5	WT	5	2	0.4s	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s
T20 N3 800	Flat 01	16	1	W	Plate	W	Plate	No	Yes	0.5s	800µL	2.5	WT	5	2	0.4s	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s
T21 N1 S1m 800	Flat 01	16	1	S	Plate	S	Plate	No	Yes	0.3s	800µL	2.5	WT	5	2	0.4s	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s
T22 N1 S1m 800	Flat 01	16	1	S	Plate	S	Plate	1.0	Yes	0.5s	800µL	2.5	WT	5	2	0.4s	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s
QCKIT	Flat 01	18	1	S	Plate	S	Plate	1.0	Yes	0.5s	800µL	2.5	WT	5	2	0.4s	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s	1	1	3	0s





## Appendix

### Revision

### Upgrade Description

**October 1996**

First official release of the user manual compatible for both PW40 and PW41 washers.

**July 1997**

Implementation of the new motorized multiwash selector

Modified pages

- page 1/6 § 1.6 LP41 Hydraulic circuit
- page 3/9 § 3.6.2 Washer's Configuration

Implementation of the selection between the new multiwash selector (MOTOR) and the old one (PINCH valve).  
Page 3/10 modified accordingly and page 3/11 created.

**September 1997**

Implementation of the new, internal design with integrated blocking device of the manifold.

Modified pages

- page 2/2 § 2. Installation

**April 1998**

Modified warning message for electrical shock hazard.

Modified pages

- page 5

Bottles are removed from top and placed beside the washer.

Modified pages

- page 1/1, 1/2, 1/4, 1/5, 1/7 § 1 « Features and Specifications.
- page 2/1 § 2.1 Unpacking.

The blocking the manifold on the left side is replaced by adhesive tape (reason : the vertical guide is shorter.)

Implemented in serial product in May 98.

Modified pages

- page 2/2 § 2. Installation.

**July 1998**

"WASH Wx" may appear as "WASH Rx" on display in several firmware versions. In the manual, the affected displays are marked with an "\*" and the alternative message is displayed on a foot note. The text associated to the displays mention both "Wx" or "Rx" possibilities.

Modified pages

- page 3/3, 3/4 & 3/9

- page 3/11 (transfer of foot note "A: back to main menu" from page 3/9) § 3 \* Operating Instructions.

Upgrade according to the configuration of the vacuum sensing (adjustable vacuum threshold).

Modified page

- page 3/10 § 3.6.2. Washer's configuration.

The blocking the manifold on the left side is still used.

Modified page

- page 2/2 § 2. Installation.

Upgrade according to the new PW40 trap vial (Yellow tube is now transparent and orange tube is now yellow).

Modified page

- page 1/3 § 1 \* Features and Specifications.

Update of the accessories and spare part list.

Modified page

- page A/1 § Appendix.

## Appendix

### Index

#### Symbols

2\*8-Channel 1/4  
4 liquid inputs 1/4  
8-Channel 1/1

#### A

Access to PROGRAMMING 3/13  
ADD : KIT 3/15  
ADD : METHOD 3/17  
ADD : PLATE 3/19  
Aerosol Protection Cover 1/1, 1/4  
AGITATION 1/19  
Agitation Amplitude 1/19  
Agitation speed 1/19  
Agitation time 1/19  
ASP. DOWNW. SPEED 1/12  
ASP. HOR. POS. 1/12  
ASP. TIME 1/12  
ASP. VERT. POS. 1/12  
ASPIRATION 1/12  
Aspiration Downward Speed 1/12  
Aspiration Horizontal Position 1/12  
Aspiration time 1/12  
Aspiration Vertical Position 1/12  
Autoclavable 1/1, 1/4

#### B

B.W.VERT.POS. 1/16  
B.WASH+B.ASP 1/9, 1/10  
Beep 3/10  
BOT.ASP. NUMBER 1/17  
BOT.DOWNW. SPEED 1/16  
BOT.SHAPE 3/19  
BOT.UPWARD SPEED 1/16  
BOT.VERT.POS. 1/18  
BOT.WASH NUMBER 1/16  
BOT.WASH+ASP 1/10  
Bottom (Aspiration) Vertical Position 1/18  
BOTTOM.ASP. 1/17, 1/18  
Bottom Downward Speed 1/16  
Bottom shape 3/19  
BOTTOM TIME 1/16  
Bottom Upward Speed 1/16  
Bottom Wash + Aspiration 1/10  
Bottom Wash Vertical Position 1/16  
BOTTOM.WASH. 1/15  
Bottom Washing + Bottom Aspiration 1/9, 1/10

#### C

Centering 3/19  
Clean manifold 4/2  
Compensation 1/13  
Contents 2  
COPY : KIT 3/30  
Crosswise aspiration 1/12  
CROSSW.ASP. 1/12  
CWA 1/12

#### D

Daily Maintenance 4/1  
DELETE 3/25  
DELETE : KIT 3/25  
DELETE: METHOD 3/25  
DELETE : PLATE 3/26  
Disinfection 4/3  
DISP. UPW. SPEED 1/13  
DISPENSING 1/13  
Dispensing Upward Speed 1/13

#### E

EDIT : KIT 3/27  
EDIT : METHOD 3/28  
EDIT : PLATE 3/29  
Electrical Data 1/7  
Electric Shock Hazard 5  
End of Edit 3/27  
ERR:24 3/16  
ERR:25 3/16, 3/18  
ERR:26 3/22  
ERR:27 3/30  
ERR:28 3/31  
ERR:29 3/26  
ERR:32 3/76  
ERR:33 3/22  
ERROR on METHOD programming 3/18  
ERROR on PLATE deletion 3/26  
ERRORS on KIT programming 3/16  
ERRORS on PLATE programming 3/22

#### F

FLOW 1/13  
Flow rate compensation 1/13

#### H

Hardware specifications 1/7  
Hardware version 4/4  
HORIZONTAL SPEED 3/20  
Hydrophobic filter 1/3

#### I

IN 3/1  
INSERT 3/23  
INSERT: KIT 3/23  
INSERT: METHOD 3/23  
INSERT : PLATE 3/24

#### K

Keyboard 3/1  
KIT Parameter 1/12  
Kit Structure 3/12

#### L

LIQUID 1/13  
Liquid level detection in bottle 1/6

## Index (continued)

### M

MAIN PARAMETERS 3/27  
Maintenance 4/1  
Manifold 1/1  
MET. INTER 1/11  
Method 1/10, 1/11  
Method Interval 1/11  
Microplate Loading 3/2  
Mini Priming 3/11

### O

Ornit strip 3/15  
OUT 3/1  
OVERFLOW 1/13  
Overflow Position 1/20

### P

Physical data 1/7  
Plate Agitation 1/19  
Plate Carrier 1/1, 1/4  
PLATE mode 1/10  
PLATE parameter 1/12  
PLATE TEST 3/32  
Please RINSE 3/6  
Prime 3/3  
Programming Synopsis 3/14  
PUMP PRIMING 3/3  
PW40 Basic features 1/1  
PW40 Hydraulic circuit 1/3  
PW40 Outside dimensions 1/2  
PW41 Basic features 1/4  
PW41 Hydraulic circuit 1/6  
PW41 Outside dimensions 1/5

### R

RAM initialization 3/33  
Routine maintenance 4/1  
RS232 serial port for connection to PC 1/4

### S

SAFETY INFORMATION 5  
Secured Stop 3/11  
SHAKE TIME 1/19  
SHAKING AMPLITUDE 1/19  
SHAKING SPEED 1/19  
Soak time 1/11  
SOAKING 1/11  
Software parameters for Wash protocols 1/9  
Software specifications 1/8  
Software version 4/4  
STRIP mode 1/10

### T

Technical features 1/7  
Temporary Kit 3/23  
Trap Vial 1/3  
Troubleshooting 4/4

### U

Unpacking 2/1

### V

Vacuum sensing 1/6  
Vertical speed 3/20  
VOLUME 1/13

### W

WARNING 5  
WASH 1/14  
WASH + ASP 1/9, 1/10  
Wash + Aspiration 1/10  
Wash + Bottom Aspiration 1/10  
WASH+BOT.ASP 1/9, 1/10

## Appendix