## User's Manual

→ Electrophoresis Power Supplies.....P. 2

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**PS 202**

**PS 3003**

Please retain all packaging materials until  
the warranty period has expired.

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**Please retain all packaging materials until  
the warranty period has expired.**

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## WARRANTY

APELEX warrants apparatus of its manufacture against defects in materials and workmanship, under normal service; for two years from the date of receipt by the purchaser. This warranty excluded damages resulting from shipping, misuse, carelessness, or neglect.

APELEX liability under the warranty is limited to the repair of such defects or the replacement of the product, at its option, and is subject to receipt of reasonable proof by the customer that the defect is embraced within the terms of the warranty. All claims made under this warranty must be presented to APELEX within 3 years following the date of delivery of the product to the customer.

This warranty is in lieu of any other warranties or guarantees, expressed or implied, arising by law or otherwise. APELEX makes no other warranty, expressed or implied, including warranties of merchantability of fitness for a particular purpose. Under no circumstances shall.

APELEX be liable for damages either consequential, compensatory, incidental or special, sounding in negligence, strict liability, breach of warranty or any other theory, arising out of the use of the product listed herein.

APELEX reserves the right to make improvements in design, construction, and appearance without notice.

### → Declaration of conformity and CE mark

The information outlined in this section applies only to customers located in the European Union (EU).

This laboratory apparatus is identified with the CE mark. This mark indicates that the product complies to the following EU Directives and Standards: see declaration of conformity on page 64.

## SAFETY PRECAUTIONS

**The power supplies are capable of delivering potentially dangerous voltage and is to be operated only by qualified technically trained personnel. Please read the entire user's manual thoroughly before operating this unit.**

**Take care as the mode of operation of the unit is continuous.**

If the power supply is used in a manner not specified by APELEX, then the protection systems of the equipment may be impaired.

For additional information, please call the APELEX or your distributor Technical Resources Department.

Never attempt to remove the outer casing or make any repairs to the unit.

Contact APELEX distributor immediately if the need for repair or servicing should arise.

**The unit must be earthed. Use only the line cord supplied with the unit for safe operation. The use of a line cord other than this or one supplied by APELEX may result in user hazard.**

For UK users, check the mains plug of the line cord to make sure it is equipped with a protection fuse not exceeding 3 A. Connect the line cord directly into a properly rated, 210/250 VAC 50/60 Hz or 105/125VAC 50/60Hz three pins wall outlet.

**For connection between the power supply and the electrophoresis equipment, use only safety output power cords equipped with Black and Red recessed plugs.**

Check the power cords and the black and red recessed safety jacks periodically to make sure that they are in good condition. Do not use cords which are cracked, nicked or in otherwise poor condition.

Always make all connections between the power supply and the electrophoresis equipment prior to start-up of the output.

Never place any objects other than high voltage connectors rated to 1000V into the output jacks.

Site the unit such that the rear panel has at least 10 cm of clearance to provide for adequate unit ventilation.

**The power supply must only be connected to electrophoresis equipment manufactured with built-in safety protected male plugs.**

**It is also recommended to use electrophoresis equipment that can only be connected when the protection lids are closed.**

## CLEANING

The power supply may be cleaned as required when the main supply is isolated.

Cleaning should be carried out with a cloth moistened with water or with tissues impregnated with 70% Iso-propyl alcohol.

No other cleaning solutions should be used.

## MAINTENANCE

There are no internal operator serviceable parts in this Power supply.

If the power should fail, the unit must be returned to the authorised Service centre.

See troubleshooting guide (page 14).

## SHIPMENT

When shipment or transport of the power supply is required, use the packaging supplied with the unit and add the decontamination form clearly complied.

Store the unit in the packaging and in a dry area.

## DECONTAMINATION

Decontamination is required when the unit has to be returned for servicing. One declaration of contamination status form is attached to this manual and must be filled up prior to inspection servicing repair or return.

**If necessary, ask for this form by fax to +33 1 69 30 71 94.**

### → Declaration of conformity and CE mark

The information outlined in this section applies only to customers located in the European Union (EU).

This laboratory apparatus is identified with the CE mark. This mark indicates that the product complies to the following EU Directives and Standards: see declaration of conformity on page 64.

## USER'S MANUAL

### Electrophoresis Power Supplies

#### PS 608


600 Volts - 800mA

#### PS1006

1000 Volts - 600mA

#### PS202

200 Volts - 2Amps

5

## WARNING

- THE POWER SUPPLIES ARE CAPABLE OF DELIVERING POTENTIALLY LETHAL VOLTAGE AND IS TO BE OPERATED ONLY BY QUALIFIED TECHNICALLY TRAINED PERSONNEL.
- PLEASE READ THE ENTIRE OPERATOR'S MANUAL THOROUGHLY BEFORE OPERATING THIS UNIT.
- TAKE CARE AS THE MODE OF OPERATION OF THIS UNIT IS CONTINUOUS

## DESCRIPTION

### Microprocessor controlled switching → Power Supplies

✓ These power supplies are CE marked and complies with the following requirements : EMC directive 89/336/CEE. It is also certified and manufactured according to the IEC 1010-1 standard and test regulation.

✓ They are equipped with an automatic restart system in case of mains failure.  
When the power returns, an audible alarm sounds for 10 seconds and the power automatically restarts with the previous set values.

If during the power failure the electrophoresis unit is disconnected from the power supply, the alarm still sounds when the mains supply returns. Because no load is connected the power supply will immediately shut-down and set the output to zero.

✓ Direct reading of programmed set values and actual values before and during the cycle.

✓ Volt, mA and Timer adjustable during a cycle.

✓ Timer for automatic shutdown of output.

✓ Stabilisation and automatic crossover between the parameters according to the set limitation values and when output limits are reached. 2 Red LED indicate the constant mode.

✓ 2 operating modes : constant voltage - constant current.

✓ Settings by tactile switches on moisture resistant membrane panel.

✓ Battery-backed memory feature “save” last output set values in the event of a power failure or when the run is terminated and the unit is turned off.

✓ 3 Red LED Displays for set and output Volt, Current and time values.

✓ 3 recessed safety output jacks allow simultaneous operation of 3 electrophoresis units.

✓ Fault detection and 500µA current leakage detection or no load automatically shuts down output and indicates the fault by yellow LED blinking.

## SPECIFICATIONS

<b>PS 608</b>	Voltage range mA range Power Timer range	1 - 600 Volts ; 1Volt step 1 - 800 mA; 1mA step 0,6 - 300Watts 1 - 9999 minutes; 1minute step
<b>PS 1006</b>	Voltage range mA range Power Timer range	1 - 998 Volts ; 1Volt step 1 - 600 mA; 1mA step 0,6 - 300Watts 1 - 9999 minutes; 1minute step
<b>PS 202</b>	Voltage range mA range Power Timer range	1 - 200 Volts ; 1Volt step 1 - 1998 mA; 2mA step 0,6 - 250Watts 1 - 9999 minutes; 1minute step
	Mode of operation	Continuous
	Switching frequency	23 kHz
	Output regulation stability	$\pm 0,2\%$ FS $\pm 1$ digit
	Minimum output	5 Volts; 15 $\mu$ A; 0,6 Watts
	Value display accuracy	$\pm 0,2$ FS $\pm 1$ digit
	Mains failure during the run	Buzzer, then automatic restart with previous set values
	Fault detection	Output supply stop, buzzer Yellow LED blinking
	Fault status	Output to earth leakage Output open circuit Output short circuit No regulation (Overheating, power circuit fault)
	Earth leakage detection level	500 $\mu$ A
	Output to earth impedance	10M $\Omega$ min bypassed by 1nF max
	Size	34 cm x 27 cm x 11 cm ( D x W x H )
	Weight	4 kg
	Environmental conditions	Indoor use, Altitude up to 2000m Temperature 10°C - 40°C
	Mains supply, 110V model	90 - 130V; 50 - 60Hz; T4A fuses
	Mains supply, 220V model	180 - 260V; 50 - 60Hz; T2A fuses
	Fuse value in the mains plug (UK)	3 A
	Rated input power/current	400 VA
	Mains supply VAC fluctuations not to exceed $\pm 10\%$ of the normal voltage	

**Maximum relative RH 80% for temperature up to 31°C decreasing linearly  
to 50% RH at 40°C.**

## OPERATION

The power supply automatically RESTART when the power is operating again after a power failure or repetitive micro failures during a cycle.

In this situation, the elapsed time until the failure is saved and memorised.

When the mains is operating again after the failure, an audible alarm rings and the STOP LED blinks for 10 seconds.

Then the power supply RESTART automatically with respect to the set values.

### → Unit power up

1. Connect the AC line cord to a grounded, 3-prong wall outlet.
2. Connect the power supply to an electrophoresis device using the power cords supplied.
3. Press the main POWER switch to turn the power supply on. The STOP LED (red) will illuminate, and the output Leds will display zeros.

### → Adjusting Output Set Values

When the main power switch is turned on, default output values are displayed. To change the output settings, use the Volt, mA or minutes dedicated tactile switches to increase or decrease the voltage, current or time settings.

As soon as one of these switches is activated, the output Leds will display output set values rather than actual output values. When they are released, the output LED will display output set values for 3 seconds, and then will display actual output (All zeros since output has not yet be started).

#### **Maximum output values are as follow:**

##### **PS 608**

for 600 Volts - maximum current 500 mA (300 Watts)  
for 800 mA - maximum voltage 380 Volts (300 Watts)


##### **PS 1006**

for 1000 Volts - maximum current 300 mA (300 Watts)  
for 600 mA - maximum voltage 500 Volts (300 Watts)

##### **PS 202**

for 200 Volts - maximum current 1500 mA (300 Watts)  
for 1998 mA - maximum voltage 150 Volts (300 Watts)

## FRONT PANEL DIAGRAM

### N° DESCRIPTION

### FUNCTION

1	Main Power switch	Turns on unit
2	Voltage output Led	Display actual or set output in Volts
3	Current output Led	Display actual or set output in mA
4	Time in minutes	Display elapsed or set time in minutes
5	STOP switch	Terminates output
6	START switch	Begins output
7	STOP led	Indicates STOP mode
8	START led	Indicates START mode
9	FAULT led	Indicates ground leakage, no load or automatic restart situation
10	Volt mode Led	Indicates constant voltage mode
11	Milliamps mode Led	Indicates constant current mode
12	Timer mode Led	Indicates timer mode ON/OFF
13	Voltage increase/decrease	Increases/decreases voltage set point
14	Current increase/decrease	Increases/decreases current set point
15	Time increase/decrease	Increases/decreases time set mode
16	ON/OFF timer switch	Put timer ON/OFF mode
17	CLEAR switch	Cancel the elapsed time of a previous running
18	Output Jacks	Red = positive; Black = negative
19	Reading switch	Display output set values

## INITIATING OUTPUT WITHOUT TIMER MODE

1. If the green LED (N° 12) located on the left of the TIMER is alighted, depress the ON/OFF switch of the TIMER on OFF position: the green LED located on the left of the display must be off.
2. Depress the START switch to initiate output. The START LED (green) will illuminate, and the output Leds will display the actual values simultaneously. Depending on the output set values one of the two modes Leds (VOLTS, MILLIAMPS) will illuminate, indicating the parameter controlling output. Slight increases or decreases in output readings for those parameters not "limiting" will occur as the experiment progresses. These changes accurately represent changes in the resistive load (electrophoresis unit) due to changes in temperature, buffer capacity, etc. This is to be expected.
3. To view the output set values during the run, depress the READING switch. The Leds will display the output set values as long as READING switch is depressed. Once the READING switch is released, the Leds will display the output set values for three seconds, and then switch back to displaying the actual output values.
4. It is possible to change the values during the run without depressing the STOP switch.  
To change the output set values during the run, depress the increase or decrease tactile switch until the appropriate output set value is reached.  
In START status, when the reading switch is depressed, the Leds will display the output set values during three seconds and then switch back to displaying the actual output values.
5. To establish the limiting (constant) mode for the particular experiment, set the controlling parameter to the output desired, and increase the other output set value until the appropriate MODE LED (VOLT, MILLIAMPS) illuminates. If the non-controlling output set values is reached during the course of the run, the power supply will automatically crossover to the new mode and control output relative to that mode. The appropriate MODE led will illuminate. If automatic crossover is desired during the run, adjust the output set value of the second controlling parameter to the maximum setting desired. When actual output relative to the second controlling parameter equals its output set value, the output will cross over from the first controlling parameter to the second.
6. When the run has been completed, depress the STOP switch to cease power output. Wait one minute before disconnecting the power cords from the gel unit. Turn the main power switch off when the unit is not in use.

## INITIATING OUTPUT WITH TIMER MODE

1. If the green LED (N° 12) located on the left of the TIMER is not alighted, depress the ON/OFF switch to activate the TIMER: then the green Led located on the left of the display illuminates.
2. Depress the CLEAR switch to cancel the previous elapsed values.
3. Set the new value in minutes by depressing the increasing and decreasing switches.  
To change the output settings, use the V - mA dedicated tactile switches to increase or decrease the voltage or current settings.
4. Depress the START switch to initiate output. The START LED (green) will illuminate, and the output Leds will display the actual values simultaneously. Depending on the output set values one of the two modes Leds (VOLTS, MILLIAMPS) will illuminate, indicating the parameter controlling output.  
Slight increases or decreases in output readings for those parameters not "limiting" will occur as the experiment progresses. These changes accurately represent changes in the resistive load ( electrophoresis unit ) due to changes in temperature, buffer capacity, etc. This is to be expected.
5. To view the output set values during the run, depress the READING switch. The leds will display the output values as long as READING switch is depressed. Once the READING switch is released, the Leds will display the output set values for three seconds, and then switch back to displaying the actual output values.
6. It is possible to change the set values during the run without depressing the STOP switch.  
To change the output set values during the run, depress the increase or decrease tactile switch until the appropriate output set value is reached.  
In START status, when the set switch is moved, the leds will display the output set values during three seconds and then switch back to displaying the actual output values.
7. To establish the limiting (constant) mode for the particular experiment, set the controlling parameter to the output desired, and increase the other output set value until the appropriate MODE LED (VOLT, MILLIAMPS) illuminates. If the non-controlling output set values is reached during the course of the run, the power supply will automatically crossover to the new mode and control output relative to that mode. The appropriate MODE LED will illuminate. If automatic crossover is desired during the run, adjust the output set value of the second controlling parameter to the maximum setting desired. When actual output relative to the second controlling parameter equals its output set value, the output will cross over from the first controlling parameter to the second.

## USER'S PROTECTION AND SAFETY

8. When the set time is elapsed, the output is automatically shutdown.

Volt and mA display zero and timer shows the set time.

Wait one minute before disconnecting the power cords from the gel unit.

Turn the main power switch off when the unit is not in use.

A Yellow LED indicates a fault situation and cut off the power supply:

1. Overload
2. Output to earth leakage
3. Output open circuit
4. Chamber lead(s) disconnected or defective
5. Output short circuit

Press STOP to resume and look for the wrong situation. Press START to run again.

## TROUBLESHOOTING GUIDE

CONDITION	PROBABLE CAUSE	REMEDY
Display fails to illuminate when the POWER switch is put on.	Fuses have blown.	<b>See Warning below.</b>
The desired MODE is not flashing.	One of the other parameters is limiting output.	Increase the output set value of the parameter controlling output until the desired output mode is controlling.
Two different modes are blinking alternatively.	Settings for both parameters are too close to the actual output.	Increase the set value for the mode you do not wish to be limiting.
The time does not elapse.	TIMER mode is OFF.	Put the TIMER ON by depressing the ON/OFF switch located under the timer display.
The TIMER mode is ON but after a mains failure the set values are not saved and other set values are displayed.	The new set values were not memorised.	Set again the new values and depress the START switch to memorise the new ues.
TIMER mode is ON, all other parameters are correctly set but the unit returns on FAULT/STOP mode when pressing the START switch.	<ol style="list-style-type: none"> <li>1. Settings are wrong</li> <li>2. Fault situation</li> <li>3. Timer is not cleared.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check the settings V/ mA.</li> <li>2. Check the connections on the electrophoresis chamber.</li> <li>3. Depress the CLEAR switch.</li> </ol>

### WARNING

**Never attempt to remove the outer casing or make any unit repairs.**

**Contact APELEX if the need for repair or servicing should arise.**

**Should the power supply fail, DO NOT remove the outer case of the unit and attempt any repairs.**

**Contact APELEX if the need for servicing the unit should arise.**

## USER'S MANUAL

### Electrophoresis Power Supplies

#### PS 9009 TX

3000 Volts - 300mA - 300W

#### PS 1006 P

1000 Volts - 600mA - 300W

13

## WARNING

- THE POWER SUPPLIES ARE CAPABLE OF DELIVERING POTENTIALLY LETHAL VOLTAGE AND IS TO BE OPERATED ONLY BY QUALIFIED TECHNICALLY TRAINED PERSONNEL.
- PLEASE READ THE ENTIRE OPERATOR'S MANUAL THOROUGHLY BEFORE OPERATING THIS UNIT.
- TAKE CARE AS THE MODE OF OPERATION OF THIS UNIT IS CONTINUOUS.

## PS 9009 TX

### Features.

Maximum Voltage	3000 V
Range of setting	10 to 3000 V by steps of 10 V
Maximum current	300 mA
Range of setting	1 to 300 mA by steps of 1 mA
Maximum power	300 W
Range of setting	1 to 300 W by steps of 1 Watt
Timer	1 to 9999 minutes by steps of 1 minute
Integrator	1 to 99999 Volts/hour by steps of 1 Volts/hour
Temperature	0 to 99°C by steps of 1°C (Temperature probe in option)
Programmability	8 programs each with 10 steps 16 single step programs 1 low current program with 10 steps
Regulation	< 1%
Accuracy	± 1,5 % full scale for each display
Number of output Terminals	2 recessed sets of 4 mm output connectors
Safety interlock	Load sensing shut-down on disconnect Key activation necessary to begin voltage generation. In the event of shut-down due to power interruption, automatic restart is provided
Ground leakage	Leakage of 500 µA will interrupt the generation of high voltage
Dimensions	27 x 34 x 11 cm (W x d x h)
Net Weight	5 kg
Power requirement	110/220 VAC, 50/60 Hz, 750 W
Main Fuse rating	T4.0A/250 V for 115 VAC Mains power T2.0A/250 V for 230 VAC Mains power

## PS 1006 P

### Features.

Maximum Voltage	1000 V
Range of setting	10 to 1000 V by steps of 10 V
Maximum current	600 mA
Range of setting	1 to 600 mA by steps of 1 mA
Maximum power	300 W
Range of setting	1 to 300 W by steps of 1 Watt
Timer	1 to 9999 minutes by steps of 1 minute
Integrator	1 to 99999 Volts/hour by steps of 1 Volts/hour
Temperature	0 to 99°C by steps of 1°C (Temperature probe in option)
Programmability	8 programs each with 10 steps 16 single step programs 1 low current program with 10 steps
Regulation	< 1%
Accuracy	± 1,5 % full scale for each display
Number of output Terminals	2 recessed sets of 4 mm output connectors
Safety interlock	Load sensing shut-down on disconnect. Key activation necessary to begin voltage generation. In the event of shut-down due to power interruption, automatic restart is provided
Ground leakage	Leakage of 500 µA will interrupt the generation of high voltage
Dimensions	27 x 34 x 11 cm (Wxdxh)
Net Weight	5 kg
Power requirement	110/220 VAC, 50/60 Hz, 750 W
Main Fuse rating	T4.0A/250 V for 115 VAC Mains power T2.0A/250 V for 230 VAC Mains power

# TABLE OF CONTENTS

<b>1. Introduction</b> .....	16
1.1 Notice to Customer .....	16
1.2 Warnings.....	16
<b>2. Overview</b> .....	16
2.1 Description .....	16
2.2 Set -Up .....	16
<b>3. Operating Instructions</b> .....	18
3.1 Programming and Operating a Single Step Program .....	18
3.2 Programming and Operating a multiple Step Program .....	18
3.3 Programming and Operating the Low Current Program.....	18
3.4 Temperature Control Program .....	19
3.5 Setting Safe Operating Limits .....	19
3.6 Relationships Between Volts, milliamps, Watts and Chamber Resistance .....	19
3.7 Running Multiple Chambers .....	20
3.8 Using Automatic Crossover.....	20
3.9 Integration .....	20
<b>4. Troubleshooting Guide</b> .....	21
4.1 Display Messages .....	21
4.1.1 Start Up Message .....	21
4.1.2 End of Run Display .....	21
4.2 Fault Detection Messages.....	21
4.2.1 No Load Detection Message.....	21
4.2.2 Ground Leakage Message.....	21
4.2.3 Power Fault Message .....	21
4.2.4 Memory Lost Message .....	21
4.2.5 Power Restart Message.....	21
4.3 General Troubleshooting .....	22
4.4 Interface RS 232 .....	22
<b>5. Additional Information</b>	
Specifications PS 9009 TX .....	14
Specifications PS 1006 P .....	14
Warranty.....	3
<b>Figure</b>	
Model PS 9009 TX.....	17

## 1. Introduction

### 1.1 Notice to Customer

This product is authorized for laboratory research use only. The product has not been qualified or found safe and effective for any human or animal diagnostic or therapeutic application. Uses for other than the labeled intended use may be a violation of applicable law.

### 1.2 Warnings

1. **DANGER ! HIGH VOLTAGE !** This power supply has been designed for use as a source of DC power for electrophoresis. It is capable of generating lethal currents. This unit should always be operated with extreme caution. Careless handling could result in electrical shock.
2. Never operate damaged equipment. Do not use the unit without the cover in place or with any possible short circuit. If the power supply emits smoke or blows the main fuses, turn off the power supply and disconnect the AC line and power cord.
3. Do not operate with connecting cables that have exposed live wires.
4. Always turn off the power supply before connecting or removing the power cords or moving an electrophoresis apparatus.
5. Handle one power cord at a time, with one hand only, when connecting or removing DC power cords at the power output terminals.
6. Do not immerse this unit in water.
7. Do not operate the unit in a damp, humid atmosphere or in a fashion where condensed moisture may short out electrical components.
8. Do not connect the output to earth ground.
9. The AC power cord has a three-pronged plug which must be connected to a grounded line voltage receptacle. Do not use a two-wire receptacle with an adapter. This could create a serious electrical hazard for persons using the unit.
10. Use the correct voltage AC power outlet for the power supply. Be sure to always plug a 110 VAC power supply into a 110-V line, and use a 240-V line with a 240 VAC unit.

## 2. Overview

### 2.1 Description

The model PS 9009 TX Programmable Electrophoresis Power Supply is intended to be used with electrophoretic devices that operate below 3000 V and up to 300 W, and/or 300 mA DC.

The unit includes two sets of output sockets which operate in parallel to provide output in constant voltage, constant milliamps, or constant power (see figure 1). This power supply has the ability to restart itself in the event of a power failure.

### —> Declaration of conformity and CE mark

The information outlined in this section applies only to customers located in the European Union (EU).

This laboratory apparatus is identified with the CE mark. This mark indicates that the product complies to the following EU Directives and Standards: see declaration of conformity on page 64.

### 2.2 Set-Up


**Unpacking the Unit.** Unpack and inspect the power supply unit carefully for any damage.

Do not use the unit if it is damaged. If damage is found, save the packing material and report the problem to APELEX Customer Service.

**Location.** Make sure the unit is set up in a location where it is protected from physical damage moisture, corrosive agents, and extreme temperatures. Make sure the fan at the rear is not obstructed ; leave about 20 cm of space behind the unit. The unit should be readily accessible for safe operation.

**Connection with the AC Mains.** Connect the unit to the AC mains carrying the appropriate specified voltage (V) in accordance with the rating label affixed to the back of the unit. Ensure that both the main receptacle and the power supply plug have the proper three-wire (grounded or earthed) connections.

## FRONT PANEL DIAGRAM

N°	DESCRIPTION	FUNCTION
1	AC Mains Power ON/OFF Switch	Turns on the main power
2	Voltage Output Display	Displays actual or set Volts
3	Current Output Display	Displays actual or set mA
4	Power Output Display	Displays actual or set Watts
5	Volt/Hour Display	Displays actual or set VHR
6	Timer Display	Displays actual or set MIN
7	Program Number	Displays actual or set program
8	STEP Number	Displays actual or set step
9	Resume Switch	Resumes timer and VHR values
10	Start/Stop Switch	Starts or stops High Voltage Output
11	Start LED	Green LED indicates High Voltage on
12	Stop LED	Red LED indicates High Voltage off
13	Volt Selector	Permits entry of volt set point
14	mA Selector	Permits entry of current set point
15	Watt or °C Selector	Permits entry of Watt or temperatures set point
16	VHR Selector	Permits the entry of Volt/Hour end point
17	MIN Selector	Permits the entry of Time end point
18	Program Selector	Permits the entry of program number
19	Previous Step Selector	Previous selection of steps
20	Next Step Selector	Next selection of steps
21	High Voltage Output Jacks	RED=POSITIVE ; BLACK=NEGATIVE
22	Numeric Keypad	Permits entry of all parameters
23	Temperature Probe Output	Connection of the external surface probe

## → 3. Operating instructions

The touch key pads on the front of the Model PS 9009 TX Programmable Electrophoresis Power supply (figure 1) are used to program and run the unit. The power supply can be operated in any of three different modes: Constant Voltage, Constant Current, or Constant Power with automatic crossover. The power supply can also be run with temperature regulation. The power supply stores up to 24 different programs which are arranged into three different types. Programs coded P01 to P08 and P11 to P18 are one time single step programs. Programs coded P21 to P27 are multiple step programs and can have up to 10 steps. The program coded P28 is designed for low current operation. This is a multiple step program for operation < 1 mA.

### 3.1 Programming and Operating a Single STEP Program

There are 16 available single step programs, coded P01 to P08 and P11 to P18, that can be set. These programs include automatic shutoff and are identified by S # on the display panel. This symbol indicates that for each program there is only a single programmed step. Each program is maintained in the memory of the power supply. Electrophoresing a vertical polyacrylamide sequencing gel at 60 W (constant power) for 3 h.

1. Switch the AC Mains Power ON/OFF Switch to the ON position. An alarm sounds to signify the start up of power.
2. Press the Start/Stop keypad to silence the alarm. The last program used is displayed.

```
0000 V  000 mA  000 W  P01
00000 VHR  0000 min  S#
```

3. Press the Program keypad and enter the program number using the numerical keypad.
4. To enter the setting parameter, press the keypad of the parameter to be changed, and then using the numerical keypad, enter the desired value. Values which are out of limits are not accepted.
5. Enter 3000 V, 300 mA, 60 W (entered 060), 99999 VHR and 180 min (entered 0180). The V, mA and VHR functions are set at maximum so they are not limiting factors. After each value has been entered, the display remains for 5 s and then changes to display all zeros. Press any keypad to return to the program that is being entered.

```
3000 V  300 mA  060 W  P01
99999 VHR  0180 min  S#
```

6. Press the Start/Stop keypad to start electrophoresis. The model PS 9009 TX Power Supply has automatic crossover with the limiting parameter being indicated by the blinking display. Electrophoresis stops at the end of the programmed time or can be manually stopped by pressing the Start/Stop keypad. Press the Start/Stop keypad to resume electrophoresis at the point where it was stopped.

### 3.2 Programming and operating a multiple Step Program

There are 7 available multiple step programs, coded P21 to P27, that can be set. Each program can have up to 10 individual steps which are identified by S0 to S9 on the display panel.

Running a program to electrophorese an SDS-PAGE gel first at a constant current of 15 mA for 30 min and then at a constant current of 30 mA for 90 min.

1. Switch the AC Mains Power ON/OFF Switch to the ON position. An alarm sounds to signify the start up of power.
2. Press the Start/Stop keypad to silence the alarm. The last program used is displayed.
3. Press the Program keypad and enter the program number using the numerical keypad.
4. Using the next and previous keypads select the step to be entered (S0 to S9).

```
0000 V  000mA  000 W  P21
00000 VHR  0000 min  SO
```

5. To enter the setting parameters, press the keypad of the parameter to be changed, and then using the numerical keypad, enter the desired value. Values which are out of limits are not accepted.
6. Enter 3000 V, 15 ma, 300 W, 99999 VHR and 15 min on step S0. The V, W, and VHR functions are set at maximum so they are not limiting factors.

```
3000 V  015 mA  300 W  P21
99999 VHR  0015 min  SO
```

7. Press the Next keypad to switch to the next step to be entered.
8. Enter 3000 V, 30 mA, 300 W, 99999 VHR, and 90 min on step S1. The V, W, and VHR functions are set at maximum so they are not limiting factors.

```
3000 V  030 mA  300 W  P21
99999 VHR  0090 min  S1
```

9. Press the Start/Stop keypad to start electrophoresis.
10. At the end of the 15 min first step in the program, the power supply switches automatically to the next step. At the end of the last programmed step, the power supply automatically shuts off the voltage and current to the electro-phoretic chamber.

### 3.3 Programming and Operating the Low Current Program

Program coded P28 is a dedicated low current (<1 mA) program that is operated and programmed in the same fashion as the other multiple step programs. It is not used for sequencing or other applications requiring the use of high current.

### 3.4 Temperature Control Program

The Model PS 9009 TX Power Supply is designed for automatic temperature regulation as well as that for Voltage, Current and Power. The optional temperature probe maintains constant temperature during electrophoresis by reducing or stopping the voltage and current output when the set temperature is surpassed.

Once the temperature drops below the set point, both voltage and current output is resumed at the set levels. The use of the temperature regulation mode prevents the activation of the power (Watt) regulation mode.

1. Switch the AC Mains Power ON/OFF switch to the ON position. An alarm sounds to signify the start up of power.
2. Press the Start/Stop keypad to silence the alarm. The last program used is displayed.
3. Press the Program keypad and enter the program number using the numerical keypad.
4. Program the desired parameters for electrophoresis. Instead of programming Watts, press the W/°C keypad a second time. The display shows - - - °C instead of --W. Enter the desired upper temperature limit for the gel. The range goes from 00° to 99°C. For example 40°C would be entered as 040

0000 V 000 mA 040°C P01
00000 MinR 0000 MinT S#

The display Min R is a timer that is dedicated for increments of minutes while the model PS 9009 TX does not supply current or voltage. This timer works only when the voltage and current are turned off in the temperature regulation mode. The main timer (MinT) remains operational at all times during the operation of the power supply and provides the automatic termination of the voltage and current output. After the electrophoresis has automatically stopped, it may be necessary to continue electrophoresis for the time that the voltage was discontinued.

### 3.5 Setting Safe Operating Limits

Electrophoretic chambers are generally designed for a relatively specific purpose. For example, submarine chambers use agarose gels to separate DNA or RNA fragments, while a DNA sequencing chamber is almost always used to separate DNA in a denaturing polyacrylamide gel. In each case, the voltage, milliamp, and wattage requirements are well defined within a reasonable range of values. In cases such as these, the user can safely assume that the manufacturer has designed the chamber to withstand the voltage and heat energy necessary to perform the electrophoretic separation when standard protocols are followed. Some types of electrophoretic chambers are specifically designed to be multipurpose devices. For example, a vertical slab gel chamber could be used for anything from DNA sequencing to isoelectric focusing, depending on the gel type and buffer system used. Choosing safe operating limits for a chamber of this type requires a higher degree of caution. Chamber manufacturers normally rate their product for maximum voltage and/or maximum wattage. Whenever

possible, contact the manufacturer and request this information. Whether you use special purpose or multipurpose chambers, the maximum operating temperature the chambers will withstand is a critical aspect of safe operation. Most electrophoretic chambers (with a few notable exceptions) are made of acrylic and plastic and should operate below a temperature of 35°C. If there is a lack of information about the capabilities of the chamber being used, regular monitoring of the operating temperature is recommended. This should be accomplished without coming into physical contact with the chamber when voltage is applied.

The following procedure illustrates how normal operating conditions can be determined and how this information can be used to choose safe operating limits for agarose gel electrophoresis in a submarine chamber at a constant voltage of 125V.

1. Adjust the operating limits to 125V, 300mA, 300W
2. Start electrophoresis and note the mA and W values when the voltage has reached 125 V and begins to blink.
3. Adjust the mA and W settings so that they exceed the actual values by 10 % or 10 mA and 5W, whichever is greater. Thus, the SET and ACTUAL values may look like this:

SET	=	0125 V	055 mA	011 W
ACTUAL	=	0125 V	045 mA	006 W

In the procedure listed above, we first determined the actual output wattage and milliamps for operation at 125 V and then chose operating limits that were at a slightly higher level than those indicated for mA and W. This approach ensures that the maximum output from the power supply will never exceed the normal operating conditions by more than 10 %.

Some users prefer to set up the chamber, switch the power supply's output ON, and adjust the operating limits before the samples are actually loaded. Once the user has confirmed that the chamber and the power supply are functioning properly and determined the safe operating limits, the power supply's output is switched OFF, and the samples are loaded. Since the power supply will remember the settings of the operating limits, the output can be switched back ON without further adjustments.

### 3.6 relationships Between Volts, Milliamps, Watts and Chamber Resistance

Three fundamental concepts form the basis for understanding the relationship between volts, milliamps and chamber resistance. When combined with the power formula they also define watts.

1. A movement of free electrons from atom to atom forms and electric current which is measured in milliamps (mA) or amps (A).
2. Electrostatic lines of force between two different charges produce a pressure that can move electrons (measured in volts).
3. All substances oppose the movement of electrons to some extent and are said to have resistance (measured in ohms). These three factors are always present in any operating

electric circuit. It is possible to incorporate them into one inclusive statement :

Ohm's Law

the value of the current that will flow in any circuit will be directly proportional to the value of the voltage applied and inversely proportional to the value of the resistance.

Or

amps = volts / resistance  
combined with The power formula :  
volts x amps = watts  
(where 1 amp = 1,000 mA)

Together, these two formulas define the relationship between volts, milliamps, watts and chamber resistance.

### 3.7 Running Multiple Chambers

This power supply is equipped with two sets of 4-mm output connectors which are connected in parallel because:

1. The voltage is applied equally to all branch paths in a parallel circuit.
2. The current flow in the branch paths of a parallel circuit is determined by the resistance of the individual paths.
3. The sum of the currents entering the branch paths of a parallel circuit is equal to the sum of the currents leaving the branch paths of a parallel circuit. A practical example of this is described as follows :

The power supply is connected to two identical horizontal submarine electrophoresis chambers (cells A and B). The power supply output is adjusted to 100 V, at constant voltage, and the current display indicates 60 mA. By applying the three rules for parallel circuits, we can determine the following information:

1. The voltage applied to both Cell « A » and Cell « B » is 100 volts (Rule 1).
2. The sum of the currents flowing through Cell « A » and Cell « B » is equal to 60 mA (Rule 3).

Switch OFF the power supply and momentarily disconnect Cell « B ». Switch the power supply back ON and note how the output current reading drops to 35 mA. From this, the following information can be derived.

1. The current flow through Cell « B » is equal to 60 mA - 35 mA [a net value of 25mA (Rule 3)].
2. The reason Cell « A » and Cell « B » have different current readings is due to the difference in resistance between Cells « A » and « B » (Rule 2).

### 3.8 Using Automatic Crossover

Certain electrophoretic techniques require the careful adjustment of operating limits and the utilization of the automatic crossover feature of this power supply.

Automatic crossover involves a transition from one mode of operation (e.g., constant current) to another mode of operation such as constant voltage.

Semi-dry electroblotting exemplifies the utility of this feature. Semi-dry transfer chambers contain two closely spaced parallel electrode plates. A "sandwich" consisting of buffer-saturated filter paper sheets on the outside and a gel and a charged transfer membrane on the inside is assembled and placed between the electrode plates.

Typical protocols suggest that the transfer should be carried out at a constant current. As the transfer process progresses, the buffer in the two filter paper layers heats up and begins to break down. This breakdown leads to an increase in the overall resistance between the two plates. In the constant current mode, the increase in resistance leads to a voltage increase. Left unchecked, the increasing voltage can eventually lead to arcing which would damage the electrode plates, the gel and the transfer membrane.

To eliminate the arcing problem, the voltage should be set at an operating limit which is below the threshold needed for the arc. As the voltage increases, it will eventually reach the predetermined operating limit. At this point, the power supply will automatically cross over from the constant current mode to the constant voltage mode of operation. As the transfer is completed, the current will gradually diminish.

### 3.9 Integration

The conductivity of the buffer in the electrophoresis chamber can change during the course of a run and result in changes in the voltage, current or wattage being applied. Integration of the output over a given time period takes such changes into account by monitoring volt-hours or milliamp-hours or watt-hours.

The vh integration method is used most frequently since the rate of migration of molecules through gels is proportional to the voltage that is being applied. The way this feature can be used is described as follows. If we run an experiment at constant voltage, at an output of 110 V for 3 h, we will accrue 330 vh. We know this is true since the voltage was held constant for 3h. During the course of this same experiment, the output current varied from a maximum of 50 mA to a minimum of 25 mA.

Now suppose we wish to repeat the experiment at constant current, but we want the same effect as 3 h at a constant voltage of 110 V.

At constant current, we know the output voltage must vary since the current varied when we performed the experiment at constant voltage. If we select a vh limit of 330 vh, the power supply automatically stops when the total for the vh equals 330, regardless of the fluctuations in the voltage.

## → 4. Troubleshooting Guide

### 4.1 Display Messages

#### 4.1.1 Start Up Message

press STOP to resume

This is the start up message when the power supply is first turned on. An alarm will also sound. Press the start/stop keypad to display the last used program.

#### 4.1.2 End of Electrophoresis Display

Regulation END  
press STOP to resume

Power output has stopped at the end of the last programmed step. Change the step or program to continue.

### 4.2 Fault Detection Messages

#### **WARNING!**

**This notice and the messages below alert you to potentially dangerous situations.**

#### **CAUTION!**

**This notice means serious damage may occur to your power supply or chamber.**

#### 4.2.1 No Load Detection Message

LOAD FAULT  
press STOP to resume

The LOAD Fault Detection message and the continuous tone alarm indicate that the circuit between the power supply and the electrophoresis apparatus is not closed. Check that all of the leads are connected properly and press the start/stop keypad to resume.

An open connection could be caused by any of the following conditions:

1. a broken electrode within the electrophoretic chamber,
2. a broken wire within the connecting cord,
3. insufficient buffer levels in the electrophoretic chamber, or
4. electrochemical energy stored within the electrophoretic chamber.

#### 4.2.2 Ground Leakage Message

LEAKAGE fault  
press STOP to resume

The LEAKAGE Fault Detection message indicates that some of the current output is flowing to earth ground. This situation can be hazardous. Switch the AC Mains Power ON/OFF Switch to the OFF position. Check the leads and connections and make sure that no buffer is leaking from the electrophoresis chamber onto the bench before proceeding. Contact the TECH-LINE if the problem persists.

#### 4.2.3 Power Fault Message

POWER fault  
press STOP to resume

The POWER Fault Detection message indicates that there is a POWER CIRCUIT FAILURE. This may be caused by internal excess current in the power supply.

Reset the power supply by switching off the AC Mains Power ON/OFF Switch to the OFF position and then ON again. Contact the Tech-Line if the problem persists.

#### 4.2.4 Memory Lost Message

MEMORY lost  
press STOP to resume

The Memory Lost message indicates that the permanent memory has been lost due to a circuit failure or a problem in the software. Reset the power supply by switching off the AC Mains Power ON/OFF Switch to the OFF position and then ON again.

**Contact the TECH-LINE if the problem persists.**

#### 4.2.5 Power Restart Message

POWER IS COMING ON  
press STOP to resume

The Power Restart message indicates that there has been a power failure of the AC Mains. The automatic restart of the unit will occur in 10 s at which time the audible alarm will stop. The power supply will restart with the set values.

### 4.3 General Troubleshooting

Many procedural and operational problems can be solved by carefully following the instructions in this manual. Some suggestions for equipment and procedural troubleshooting are given below. Should these suggestions not resolve the problem, or if you have questions regarding procedures, call the TECH-LINE (numbers listed on the back cover of this manual). If you need to return the unit for repair, contact the Customer Service Department or your local distributor for shipping instructions. Include a full description of the problem.

**PROBLEM****COMMENTS**

Display fails to illuminate when the power supply is turned on	Make sure that the AC Mains power cord is plugged into a functioning outlet. Check the AC Mains fuse on the back of the power supply. Replace if blown.
Desired MODE not flashing	Increase the output setting for the MODE that is limiting until the desired MODE is flashing.
Two different MODES are blinking alternatively	Settings for both parameters are too close to actual output. Increase the set value for the MODE that is not limiting.
Under constant power MODE (W), output displays are not stable.	Load and gel resistance are changing. Disregard slight fluctuations. Display voltage and current will stabilize as gel equilibrates. Increase or decrease W set value by 1 or 2 Watts to stabilize.
When using the temperature control	When the temperature limit is reached, the power supply shuts off both voltage and current output. Check the MinR display and then electrophorese the gel for the total time displayed on this timer. Make sure that the correct voltage limiting value was set.
Temperature Display shows symbols instead of temperature	Make sure that the temperature probe is properly connected to power supply.
Voltage is shut off during electrophoresis	Increase the temperature limit value. Decrease the electrophoresis parameters.
Display actual temperature is not stable	Make sure that the probe is in complete contact with the glass plates.

## POWER SUPPLY PS 9009 TX + INTERFACE RS 232

Connect the PS 9009TX on 230VAC.

Connect the RS232 lead between computer serial output and power supply output located on the back panel.

Switch on the power supply first then the computer. Follow this procedure for a proper connection between the two devices.

The unit is automatically on STOP mode at power up.

Install the diskette software on the hard disc of the computer or directly open the diskette.

The computer screen is similar to the power supply display and below, the command parameters are listed.

### IMPORTANT

**To make possible the keyboard control, you must be on SHIFT mode.  
Before starting the run just press C (CLEAR) to cancel the previous elapsed time and VHR.**

e.g To change voltage, press V and enter new setting  
To change minutes, press M and enter new setting  
Etc.

When settings are OK, press S to start the run.

The PS and computer displays show the previous preset values for a few seconds and then come back on actual values.

The blinking parameter is the constant parameter.

To change any parameter during the run, you must stop the run:

**Press S, change the values and press S to start again the run**

If you work on timed mode, the unit will stop automatically and will display the following message :

**Regulation end  
Press stop to resume**

Press S to cease the audible alarm and press C to cancel the elapsed time and VHR

### IMPORTANT


**If any displaying problem appears on the computer, press X to come back on A:\).  
Then enter APELEX.**

## USER'S MANUAL

### Electrophoresis Power Supplies

## PS 3003

3000 Volts - 300mA - 300W

23

## WARNING

- THE POWER SUPPLY IS CAPABLE OF DELIVERING POTENTIALLY LETHAL VOLTAGE AND IS TO BE OPERATED ONLY BY QUALIFIED TECHNICALLY TRAINED PERSONNEL.
- PLEASE READ THE ENTIRE OPERATOR'S MANUAL THOROUGHLY BEFORE OPERATING THIS UNIT.
- TAKE CARE AS THE MODE OF OPERATION OF THIS UNIT IS CONTINUOUS.

## DESCRIPTION

### Microprocessor controlled switching Power Supplies

- ✓ This power supplies is CE marked and complies with the following requirements: EMC directive 89/336/CEE. It is also certified and manufactured according to the IEC 1010-1 standard and test regulation.
- ✓ It is equipped with an automatic restart system in case of mains failure.  
When the power returns, an audible alarm sounds for 10 seconds and the power automatically restarts with the previous set values.  
If during the power failure the electrophoresis unit is disconnected from the power supply, the alarm still sounds when the mains supply returns. Because no load is connected the power supply will immediately shut-down and set the output to zero.
- ✓ Direct reading of programmed set values and actual values before and during the cycle.
- ✓ Volt, mA and Watts adjustable during a cycle.
- ✓ Stabilisation and automatic crossover between the parameters according to the set limitation values and when output limits are reached. 3 Red LED indicate the constant mode.
- 24 ✓ 3 operating modes : constant voltage - constant current – constant power.
- ✓ Settings by optical encoders and moisture resistant membrane panel.
- ✓ Battery-backed memory feature “save” last output set values in the event of a power failure or when the run is terminated and the unit is turned off.
- ✓ 3 Red LED Displays for set and output Volt, Current and power values.
- ✓ 2 recessed safety output jacks allow simultaneous operation of 3 electrophoresis units.
- ✓ Fault detection and 500µA current leakage detection or no load automatically shuts down output and indicates the fault by yellow LED blinking.

### → Declaration of conformity and CE mark

The information outlined in this section applies only to customers located in the European Union (EU).

This laboratory apparatus is identified with the CE mark. This mark indicates that the product complies to the following EU Directives and Standards: see declaration of conformity on page 64.

## SAFETY PRECAUTIONS

**THE POWER SUPPLIES ARE CAPABLE OF DELIVERING POTENTIALLY LETHAL VOLTAGE AND IS TO BE OPERATED ONLY BY QUALIFIED TECHNICALLY TRAINED PERSONNEL. PLEASE READ THE ENTIRE OPERATOR'S MANUAL THOROUGHLY BEFORE OPERATING THIS UNIT.**

### **TAKE CARE AS THE MODE OF OPERATION OF THE UNIT IS CONTINUOUS**

If the power supply is used in a manner not specified by APELEX, then the protection systems of the equipment may be impaired.

For additional information, please call the APELEX or your distributor Technical Resources Department.

Never attempt to remove the outer casing or make any repairs to the unit. Contact us immediately if the need for repair or servicing should arise.

### **The unit must be earthed. Use only the line cord supplied with the unit for safe operation**

Check the mains plug of the line cord to make sure it is equipped with a protection fuse not exceeding 3 A. The use of a line cord other than this or one supplied by APELEX may result in user hazard.

Connect the line cord directly into a properly rated, 210/250 VAC 50/60 Hz or 105/125VAC 50/60Hz three pins wall outlet.

### **For connection between the power supply and the electrophoresis equipment, use only the safety output power cords equipped with Black and Red recessed plugs as supplied with the unit.**

Check the power cords and the black and red recessed safety jacks periodically to make sure that they are in good condition. Do not use cords which are cracked, nicked or in otherwise poor condition.

Always make all connections between the power supply and the electrophoresis equipment prior to start-up of the output. Never place any objects other than high voltage connectors rated to 3000V into the output jacks.

Site the unit such that the rear panel has at least 20cm of clearance to provide for adequate unit ventilation.


**The power supply must only be connected to electrophoresis equipment manufactured with built-in safety protected male plugs. It is also recommended to use electrophoresis equipment that can only be connected when the protection lids are closed.**

## SPECIFICATIONS

V - mA – W setting range	10 - 3000 Volts ; 10Volt step 1 - 300 mA; 1mA step 1 – 300 Watts; 1 watt step
V - mA – Watts display range	0 - 3000 Volts ; 10 Volt step 0 - 300 mA; 1mA step 0 - 300Watts
Mode of operation	Continuous
Switching frequency	23 kHz
Output regulation stability	$\pm 0,2\%$ FS $\pm 1/2$ digit
Minimum no regulated output	5 Volts; 15 $\mu$ A; 0,6 Watts
Minimum regulated output	10 Volts; 1mA; 1 Watts (load >10Kohms)
Value display accuracy	$\pm 0,2$ FS $\pm 1/2$ digit
Mains failure during the run :	Buzzer, then automatic restart with previous set values
Fault detection	Output supply stop, buzzer Yellow LED blinking
Fault status	Output to earth leakage Output open circuit Output short circuit No regulation (Overheating, power circuit fault)
Earth leakage detection level	500 $\mu$ A
Output to earth impedance	10M $\Omega$ min bypassed by 1nF max
Size	34 cm x 27 cm x 11 cm ( D x W x H )
Weight	5 kg
Environmental conditions	Indoor use, Altitude up to 2000m Temperature 10°C - 40°C
Mains supply, 110V model	90 - 130V; 50 - 60Hz; T4A fuses
Mains supply, 220V model	180 - 260V; 50 - 60Hz; T2A fuses
Fuse value in the mains plug (UK)	3 A
Rated input power/current	400 VA
Mains supply voltage fluctuations not to exceed $\pm 10\%$ of the normal voltage	

**Maximum relative RH 80% for temperature up to 31°C decreasing linearly to 50% RH at 40°C.**

## FRONT PANEL DIAGRAM

N°	DESCRIPTION	FUNCTION
1	Main power switch	Turns on unit (located on the right side)
2	Voltage output LED	Displays actual output in Volts
3	Current output LED	Displays actual output in mA
4	Power output LED	Displays actual output in Watts
5	STOP switch	Terminates output
6	START switch	Begins output
7	STOP LED	Indicated STOP output
8	START LED	Indicated START output
9	Fault LED	Indicated ground leakage, no load or short circuit
10	Volt LED	Indicates constant voltage mode
11	mA LED	Indicated constant current mode
12	Watt LED	Indicated constant power mode
13	Voltage Switches	Increase/decrease voltage set point
14	mA Switches	Increase/decrease current set point
15	Watt Switches	Increase/decrease power set point
16	High Voltage output jacks	Red = positive; Black = negative
17	PRESET switch	Displays output preset values

## OPERATION

The power supply automatically RESTART when the power is operating again after a power failure or repetitive micro failures during a cycle.

When the mains is operating again after the failure, an audible alarm rings and the STOP LED blinks for 10 seconds.

Then the power supply RESTART automatically with respect to the set values.

### → Unit power up

1. Connect the AC line cord to a grounded, 3-prong wall outlet.
2. Connect the power supply to an electrophoresis device using the power cords supplied.
3. Press the main POWER switch to turn the power supply on. The STOP LED (red) will illuminate, and the output leds will display zeros.

### → Adjusting Output Set Values

When the main power switch is turned on, default output values are displayed.

To change the output settings, use the Volt, mA or Watts dedicated adjustment switches to increase or decrease the voltage, current or power settings.

As soon as one of these adjustment switches are activated, the output Leds will display output set values rather than actual output values. When they are released, the output LED will display output set values for 3 seconds, and then will display actual output (All zeros since output has not yet be started).

**Maximum output values are as follow:**

**for 3000 Volts - maximum current 100 mA (300 Watts)**

**for 300 mA - maximum voltage 1000 Volts (300 Watts)**

### → Initiating Output

1. Depress the START switch to initiate output. The START LED (green) will illuminate, and the output Leds will display the actual values simultaneously. Depending on the output set values one of the three modes led (VOLTS, MILLIAMPS, POWER) will illuminate, indicating the parameter controlling output.
2. Slight increases or decreases in output readings for those parameters not "limiting" will occur as the experiment progresses. These changes accurately represent changes in the resistive load ( electrophoresis unit ) due to changes in temperature, buffer capacity, etc. This is to be expected.
3. To view the output set values during the run, depress the PRESET switch. The Leds will display the output set values as long as this switch is depressed. Once the PRESET switch is released, the led will display the output set values for three seconds, and then switch back to displaying the actual output values.

4. It is possible to change the values during the run without depressing the STOP switch. During the run, as soon as one of these adjustment switches is activated, the output Leds will display output set values.

To change the output set values during the run, depress the STOP switch. Increase or decrease the adjustment switches until the appropriate output set value is reached.

In START status, when the PRESET switch is depressed, the Leds will display the output set values during three seconds and then switch back to displaying the actual output values.

5. To establish the limiting (constant) mode for the particular experiment, set the controlling parameter to the output desired, and increase the other output set value on their maximum values.

If the non-controlling output set values is reached during the course of the run, the power supply will automatically crossover to the new mode and control output relative to that mode. The appropriate MODE led will illuminate.

If automatic crossover is desired during the run, adjust the output set value of the second controlling parameter to the desired setting. When actual output relative to the second controlling parameter equals its output set value, the output will cross over from the first controlling parameter to the second.

6. When the run has been completed, depress the STOP switch to cease power output. Wait one minute before disconnecting the power cords from the gel unit. Turn the main power switch off when the unit is not in use

### → User's protection and safety

A Yellow LED indicates a fault situation and cut off the power supply :

1. Overload.
2. Output to earth leakage
3. Output open circuit
4. Chamber lead (s) disconnected or defective.
5. Output short circuit.

**Press STOP to resume and look for the wrong situation.**

**Press START to run again.**

### → Power failure during the run

The power supply automatically RESTART when the power is operating again after a power failure or repetitive micro failures during a cycle.

When the mains is operating again after the failure, an audible alarm rings and the STOP LED blinks for 10 seconds.

Then the power supply RESTART automatically with respect to the set values.

## TROUBLESHOOTING GUIDE

CONDITION	PROBABLE CAUSE	REMEDY
Display fails to illuminate when the POWER switch is put on.	Fuses have blown.	See Warning below.
The desired MODE is not flashing.	One of the other parameters is limiting output.	Increase the output set value of the parameter controlling output until the desired output mode is controlling.
Two different modes are blinking alternatively.	Settings for both parameters are too close to the actual output.	Increase the set value for the mode you do not wish to be limiting.
Under constant power mode Voltage and current output. Displays are not stable.	Load, gel resistance are changing. This usually happens when pre-running the gel.	Disregard slight fluctuations. Displays will stabilize as gel equilibrates. Increase or decrease W set value by one or two watts.

### WARNING

**Never attempt to remove the outer casing or make any unit repairs.**

**Contact APELEX if the need for repair or servicing should arise.**

**Should the power supply fail, DO NOT remove the outer case of the unit and attempt any repairs.**


**Contact APELEX if the need for servicing the unit should arise.**

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