How to setup the SmartNode?

What do the icons mean?

How to paste the ThermoSensor?

Download the latest User Manual from: www.cadisense.com Support Hotline: +65 6276 2676





(€ 0434



Medical Device Safety Service GmbH Schiffgraben 41 30175 Hannover Germany

Complies with IDA Standards DB 102490



How to setup the SmartNode?

Before you start...

Please write down the info shown on your SmartNode plug here:

\$\$ID (e.g., CADI_XXXX): _____

PWD (8 chars):

SmartNode ID:

Step 1

Plug in SmartNode in the room where the patient will be resting. Then switch it on.

(Note: SmartNode takes about 30 secs to start up.)

Step 2

Turn on WiFi on your mobile device. Then connect it to SmartNode using the above SSID and password (PWD, 8 chars).

Step 3

Go to web browser and enter:

http://10.10.10.10

You'll see this >

Quick setup is done!



If you'd like to access Cadi. Sense outside your home, please continue:

Step 4

Click the icon , then click "Connect to WiFi Router".

You'll see this →

Click on your router's SSID, key in your router's password, then click "Apply".

(Please be patient as it will take 2 mins.)



Step 5

Go to **www.cadisense.com**, then click "Access Cadi. Sense". Key in your SmartNode ID (see before Step 1), password is "**cadi**", default port is "**50000**", then click "Loain".

or

If you are using a smartphone, download **Cadi.Sense app** from iPhone App Store or Android Market. Click "+" to create a new profile for your child. Then select your child, click another "+" and select "Cadi.Sense Monitoring". Key in your ThermoSensor ID (shown on the sensor), select "Monitoring" and you're done!

What do the icons mean?





Main Page Sensor Selection

Home Page Temperature Monitoring



Temperature Chart

- 1. Go to Home Page (Temperature Monitoring)
- 2. Refresh
- 3. Advanced Settings
- 4. Set Timezone, Date & Time
- 5. Add / Remove Sensor
- 6. Click to go Home Page
- 7. Go to Main Page (Sensor Selection)
- 8. Sensor Settings (enter name, choose ⁹C / F, set temperature threshold, alert message)
- Temperature Chart / History (note: chart will appear only after 15 mins; "Lo" means below 32.1 °C; "Hi" means above 42.0 °C)
- 10. Toggle Sound Alert (only for iPhone / iPad)
- 11. Help
- 12. Sensor Battery Status
- 13. Export Temperature History Data (only for laptop)
- 14. Shift Forward / Backward 12 Hours
- 15. Jump To Date

How to paste the ThermoSensor?

Step 1

Apply a small amount of body lotion onto the lower abdomen area

Step 2

Take a Stick On and peel off the round paper liner to expose the adhesive.



Step 3

Paste the ThermoSensor on the adhesive, with the metal part facing outward and away from the adhesive.



Step 4

Remove the rectangular backing liner.



Step 5

Paste the sensor by sticking the adhesive part to the inside of the diaper/undergarment, so that the sensor's metal part will touch the lower abdomen (left or right) when the diaper/ undergarment is properly worn.



(Note: ThermoSensor takes 10 mins to warm up, temperature is wirelessly sent every 30 secs.)

Keeping the sensor

Please remember to keep your sensor back into the storage case properly to switch off the sensor. This will save battery power and make your sensor last longer!