

Incucenter ICN52/118

User Manual



swiss made +

salvisLAB

Copyright © SalvisLab / Edition 180311

CE DECLARATON OF CONFORMITY	4
IMPORTANT INFORMATION	5
INFORMATION FOR SERVICE	5
TECHNICAL DATA	6
INTRODUCTION	7
OVERVIEW	7
APPLICATIONS	7
CONSTRUCTION	7
CONTROLLER.....	7
SAFETY / POWER CUT-OFF	7
GETTING STARTED.....	8
PARTS DELIVERED	8
INSTALLATION REQUIREMENTS.....	8
INSTALLING	8
CLEANING	8
SYSTEM COMPONENTS	9
CONTROLLER	10
KEYPAD & DISPLAY	10
OPERATING.....	11
HOW TO INTERPRET DISPLAYS MESSAGES DESCRIBED IN THIS MANUAL	11
MAIN MENU - OVERVIEW.....	12
1 MAIN MENU - TEMP & OPTIONS.....	13
2 MAIN MENU - PROGRAM.....	14
2.1 MENU PROGRAM - START	15
2.2 MENU PROGRAM - NEW.....	16
2.3 MENU PROGRAM- EDIT	17
2.4 MENU PROGRAM - DELETE.....	19
2.5 MENU PROGRAM - PRINT.....	19
3 MAIN MENU - CONFIGURATION	20
4 MENU SERVICE MODE	22
STATUS DISPLAYS.....	23
GENERAL	23
5 STATUS DISPLAY: MANUAL MODE – DEFINED START DATE	23
6 STATUS DISPLAY: MANUAL MODE – PROCESS RUNNING.....	23
7 STATUS DISPLAY: PROGRAM MODE - DEFINED PROGRAM START DATE	24
8 STATUS DISPLAY: PROGRAM MODE - PROGRAM RUNNING.....	24
9 CANCEL A RUNNING PROCESS BY ESC BUTTON	25
10 MESSAGES AND ERRORS	25
PRINTER OPERATION	26
CONNECTING A PRINTER	26
EXAMPLES FOR DIFFERENT PRINTOUTS.....	26
APPENDIX A	27
WIRING DIAGRAM ICN.....	27
APPENDIX B	28
DRAWING ICN52	28
APPENDIX C	29
DRAWING ICN118	29

APPENDIX D **30**
DRAWING SPARE PARTS 30

APPENDIX E **31**
SPARE PART NUMBERS 31

CE DECLARATON OF CONFORMITY



Declaration of Conformity

Wir

We

Nous

Renggli AG / SalvisLab

(Name des Anbieters) (supplier's name) (nom du fournisseur)

Birkenstrasse 31, CH-6343 Rotkreuz

(Anschrift) (address) (adresse)

erklären in alleiniger Verantwortung, dass das Produkt
declare under our sole responsibility that the product
déclarons sous notre seule responsabilité que le produit

Incubator

ICN 52 / 118

Year of Construction

2009

(Bezeichnung Typ oder Modell, Los-, Chargen- oder Seriennummer, möglichst Herkunft und Stückzahl)
(name, type or model, lot, batch or serial number, possibly sources and numbers of items)
(nom, type ou modèle, no de lot, d'échantillon ou de série, éventuellement sources et nombre d'exemplaires)

auf das sich diese Erklärung bezieht, mit der/den folgenden Norm(en) oder normativen Dokument(en)
übereinstimmt
to which this declaration relates is in conformity with the following standard(s) or other normative
document(s)
auquel se réfère cette déclaration est conforme à la (aux) norme(s) ou autre(s) document(s) normatif(s)

(Titel und/oder Nummer sowie Ausgabedatum der Norm(en) oder der anderen normativen (Dokumente)
(title and/or number and date of issue of the standard(s) or other normative document(s)
(titre et/ou no et date de publication de la (des) norme(s) ou autre(s) document(s) normatif(s))

Gemäss den Bestimmungen der Richtlinien; following the provisions of Directives; conformément aux
dispositions des Directives; (falls zutreffend) (if applicable) (le cas échéant)

Safety EN 61010-1:2001
Low voltage directive 2006/95/EC

Safety EN 61010-2:2003
EMC directive 2004/108/EC

(Ort und Datum der Ausstellung)
(Place and date of issue)
(Lieu et date)

(Name/Unterschrift oder Kennzeichnung des Befugten)
(name and signature or equivalent marking of authorised person)
(nom et signature du signataire autorisé)

Rotkreuz, 12.10.2009

Marcel Käppeli
Technical Manager

Important Information

Information for Service

Please fill out all necessary information about your Vacucenter

SERIAL NUMBER	
TYPE	
PHONE NUMBER SERVICE	
PURCHASE DATE	
SOFTWARE VERSION (See Display during Power ON)	

Technical Data

		ICN 52	ICN 118
Outer Dimension			
Width	mm	460	570
Height (incl. feet)	mm	525	650
Depth	mm	590	720
Clearance distance from back wall	mm	100	100
Clearance distance from side wall	mm	100	100
Inner Dimension			
Width	mm	340	450
Height	mm	370	500
Depth	mm	410	530
Internal Volume	L	52	118
Wire shelf	Standard/max	2/8	2/8
Wire shelf dimension (WxD)	mm	330x365	437x495
Distance between shelf	mm	40	40
Max. load per shelf	kg	10	10
Max. load per oven	kg	40	40
Weight from oven (empty)	kg	45	55
Temperature data			
Temperature Range approx. > 5°C oRt to	°C	80	80
Temperature Variation ¹⁾	at 37°C ± °C	0.5	0.5
	at 50°C ± °C	0.7	0.8
Temperature fluctuation ²⁾	at 37°C ± °C	0.1	0.1
	at 50°C ± °C	0.1	0.2
Heating up time ³⁾	to 37°C Min	25	41
	to 50°C Min	50	70
Recovering time after 30 sec door opening ³⁾	at 37°C Min	6	7
	at 50°C Min	8	9
Electrical data			
Power supply (± 10%) 50/60 Hz	Volt	230	230
Nominal Wattage	Watt	450	550
Energy consumption	at 37°C Watt	125	130
Housing protection (EN 60529)		IP20	IP20
Recommended ambient values			
	°C	15 – 25	15 – 25
	RhF%	20 - 60	20 - 60
Equipment			
Microprocessor -Temperature Controller LCD Display		Yes	Yes
Timer	Hours / Min.	999 / 59	999 / 59
RS232 Interface		Yes	Yes
Automatic Restart after Power Cut-Off		Yes	Yes
Adjustable Print Interval		Yes	Yes
Programming	Program / Step	50 / 15	50 / 15
Ramp function		Yes	Yes

- 1) Measured with 3 temperature probes on horizontal level / divided in 1/3 of the chamber size
- 2) maximum temperature deviation in time for one temperature probe
- 3) to 98% of set temperature

All technical specification are specified for units with standard equipment at an ambient temperature of 25 °C (77 °F) and a voltage fluctuation of ± 10 %. The temperature data are evaluated in accordance to following DIN 12880, part 2 respecting the recommended wall clearances of 10 % of the height, width and depth of the inner chamber. All indications are average values, typical for units produced in series. Differing ambient temperatures or variances in the design of individual equipment may produce different performance. We reserve the right to alter technical specifications at all times without prior notice.

Introduction

Overview

The Incucenter ICN is a Incubator with **natural convection** and a Microprocessor Controller with enhanced Fuzzy Logic. It allows precise ramping of temperature as well as an excellent reproduction of temperature distribution in the chamber.

Special Insulation - Less heat loss. Saves energy and costs.

Robust Swiss quality design – Made for scientific applications

Work Chamber made of stainless steel and provided with wire shelf

Exterior is of textured powder coated mild steel

Applications

The Incucenter is designed for all purposes of incubation in a variety of laboratory fields.

All Systems have a controller with alphanumeric display and programming capabilities.

Temperature range up to 80°C.

Note: The Incucenter is not built to use as ovens for drying substances which are explosive or let free explosive gases during the drying process.

Construction

Extremely compact construction. Saves valuable space in the lab.

Inner chamber made of stainless steel.

For max. 8 wire shelves.

Door seal can easily be removed for cleaning or replacing.

The door handle is integrated in the door. No hot surfaces. There is no risk of injury or burn.

Controller

Fuzzy Logic Microprocessor controller with digital alphanumeric LCD-Display and real time clock.

Wide range of temperature ramping functions. More user application. In combination with Fuzzy Logic gives you an excellent stability of temperature distribution and accuracy of programmed values.

Brilliant LCD Display for user dialog and easy to operate keypad for fast programming and operating.

User dialog with controller is displaying your local language. Up to five languages can be selected.

Easy to operate and programming with EasyMenu

It allows the storage of 50 programs with 15 program steps (step = 1 ramp, 1 temperature and 1 Hold Time). The programs remain stored in the memory even without external power

Holding Time 0 - 999h 59m

The real time clock allows a process to be started at any time – i.e.: on December 6, 2009 / 10:30 PM

RS232 interface. All data can be achieved with a printer or computer.

Door switch interrupts heating during open door.

Safety / Power Cut-Off

DIN 12880 class 3.1

In case of over temperature, a built in safety controller as a back-up circuit takes over the control of the heating and will shutdown the oven.

There is also an additional mechanical over temperature device which shuts down the oven.

After a Power Cut-Off the Controller restarts automatically and the last stored values will be maintained again (only for Menu Temp & Options)

High quality accurate PT 100 temperature sensor.

Superior "Swiss Made" manufacturing quality.

Getting Started

Parts Delivered

Your System will be delivered with following Parts:

- 1 System Unit
- 2 Wire Shelf
- 1 Power Cord
- 1 User Manual

Installation Requirements

Electric power connection as per Serial Number plate on the inner door

Installing

Place shelf in appropriate position

Plug cord

Close door

Switch power on

Display shows current Software version during power on sequence

To start oven see chapter Operating Menu

Cleaning

To clean the System use mild detergents. No Acid or similar detergents should be used.

System Components

- 1 Air Exhaust
- 2 Flap for Exhaust
- 3 Air Inlet
- 4 Door Handle
- 5 Controller
- 6 RS232 Interface
- 7 Shelf Glider
- 8 Wire Shelf
- 9 Spring loaded door lock
- 10 Door Switch to indicate a open door (alarm switch)



Controller



To control the system the controller has few “easy to use” buttons. Simply turn / press desired button. All information are displayed on a Liquid Crystal Display (LCD).

LED indicators for status of Power, Heating and Program are used to indicate the main process status.

Buzzer is indicating audible Status or Alarm situation


Keypad & Display

Buttons




Down Scroll through the menu structure



Plus Increase a value in flashing input fields , select multiple choices; Input value must flash to make  button operable



Minus Decrease a flashing value, select multiple choices; Input value must flash to make  button operable



Enter Confirm a menu selection, a value input or an answer choice



ESC

Escape Cancel actual operation, goes back in Menu Structure, quit actual state



ON/OFF Switch ON or OFF your System. When ON LED is lit

LED Indicators



When **lit**: Heat Power is on



When **lit**: You are running a program
When **flashing**: Delayed start for a program is entered

How to interpret Displays Messages described in this Manual

Power ON Sequence



Salvis Lab
Revision X.XX

By pressing the power ON button, the display will show the Software version. All data's and stored information's are loaded during this process. After a while display will show Main Menu

Input Field



Set Temperature
37.0°C

A input value which is highlighted with yellow (grey Background) means this value is flashing on the real display.

Multiple Input Fields



Start Date
DDMMYY 01.01.10

If you reach a multiple input display, the first part (... of 3 in this example) of the input field is flashing to indicate input here. Flashing input fields are changed with ↶ or ↷

Multiple Choice Fields



Start
Now Laeter

In a multiple choice field the last actual saved value is flashing.
Change with ↶ or ↷ and confirm with ↵

Definitions of terms

What is a Set Temperature?

A Set Temperature is the target temperature you want operate the oven with.

What is a Gradient?

A Gradient is the slope of the heating up process to the specific Set Temperature. It is indicated as °C / Minute. Negative Gradients are not allowed. The maximal value of a gradient is system depending and has a range and is pre-defined by the factory setup.

What is a Hold Time?

A timer is used to specify how long a set temperature has to be hold. The timer starts counting back when 98% of the set temperature is reached. The maximal time you can set is: 999 hours and 59 minutes. This is equal to 41 days, 15 hours and 59 minutes.

What is a Start Date/Time

If you are using this feature you will be able to set a date/time to start a process or a program in the future.

Main Menu - Overview

General operation buttons

In general you can scroll through the menu points with the \curvearrowright button.

Select the desired menu point with \leftarrow

1 Temperature & Options



Manual Operation: Define a set temperature and options like Gradient and Hold Time

\curvearrowright Change Menu

\leftarrow Accept Menu \rightarrow 1.1

ESC returns to \rightarrow 1

2 Program



Programmed Operation: The menu Program is divided in menus for creating, editing, deleting, starting and printing programs

\curvearrowright Change Menu

\leftarrow Accept Menu \rightarrow 2.1

ESC returns to \rightarrow 1

3 Configuration



This menu point allows you to configure the system

\curvearrowright Change Menu

\leftarrow Accept Menu \rightarrow 3.1

ESC returns to \rightarrow 1

4 Service Mode



This menu point is protected by a access code and has limited access for the user

\curvearrowright Change Menu

\leftarrow Accept Menu \rightarrow 4.1

ESC returns to \rightarrow 1

1 Main Menu - Temp & Options



Manual Operation: Define a set temperature and options like Gradient, Hold Time (Dwell Time) and Start date/time
← Accept menu → 1.1
ESC returns to → 1

1.1 Set Temperature



↶ ↷ Change value
← Accept value → 1.2
ESC restores the old value or returns to → 1

1.2 Select Quick Start or Start with Options



↶ ↷ Change answer
← Accept answer
If **Now** selected: System will start immediately → 6
If **Option** selected → 1.3
ESC cancels and returns to → 1

1.3 Set Gradient



↶ ↷ Change value
← Accept value → 1.4
ESC restores the old value or returns to → 1
Note: Value of 0.0 means maximum heating power up to set temperature!

1.4 Set Holding Time (dwell time)



↶ ↷ Change value
← Accept value and skips to the next input field (HH→MM) or stores the time and go to → 1.5
ESC restores the old value and skips back one input field (MM→HH) or goes back to → 1
Note: Value of 00:00 means endless Holding Time

1.5 Set Start-Date



↶ ↷ Change value
← Accept value and skips to the next input field (DD→MM→YY) or stores the date and go to → 1.6
ESC restores the old value and skips back one input field (YY→MM→DD) or goes back to → 1

1.6 Set Start Time



↶ ↷ Change value
← Accept value and skips to the next input field (HH→MM) If Start Date/Time is in the **past** you will see the operating display → 6
If Start Date/Time is in the **future** you will see the status display → 5
ESC restores the old value and skips back one input field (MM→HH)

2 Main Menu - Program



Main Menu
Programs

Program Operation: The menu Program is divided in menus for creating, editing, deleting, starting and printing programs.
⤴ Scroll through the menu
⤵ Accept menu → 2.1
ESC returns to → 1

2.1 Menu Program - Start



Menu Program
Start

Start a stored program
⤴ Scroll through the menu
⤵ Accept menu → 2.1.1
ESC return to → 2

2.2 Menu Program - New



Menu Program
New

Create a new program
⤴ Scroll through the menu
⤵ Accept menu → 2.2.1
ESC return to → 2

2.3 Menu Program - Edit



Menu Program
Edit

Edit a stored existing program
⤴ Scroll through the menu
⤵ Accept menu → 2.3.1
ESC return to → 2

2.4 Menu Program - Delete



Menu Program
Delete

Delete a stored program
⤴ Scroll through the menu
⤵ Accept menu → 2.4.1
ESC return to → 2

2.5 Menu Program - Print



Menu Program
Print

Print a stored program
⤴ Scroll through the menu
⤵ Accept menu → 2.5.1
ESC return to → 2

2.1 Menu Program - Start



Start an existing Program
← Accept menu → 2.1.1
ESC returns to → 2

2.1.1 Select Program



↶ ↷ Change value
← Accept value → 2.1.2
ESC returns to → 2.1

Note: Only stored Program Numbers will appear in the display. If no Program exist a message will be displayed.

2.1.2 Choose type of program start



↶ ↷ Change answer
← Accept answer

If **Now** selected: Program starts immediately → 6

If **Later** selected : → 2.1.3

ESC returns to → 2.1

2.1.3 Set Start Date



↶ ↷ Change value
← Accept value and skips to the next input field

(DD → MM → YY) or stores the date and go to → 2.1.4

ESC restores the old value and skips back one input field (YY → MM → DD) or goes back to → 2.1

2.1.4 Set Start Time



↶ ↷ Change value
← Accept value and skips to the next input field

(HH → MM) If the Start Date/Time is in the **past**, the display will return back to → 2.1.2

If your Start Date/Time is in the **future** you will see the operating display → 7

ESC restores the old value and skips back one input field (MM → HH) or goes back to → 2.1

2.2 Menu Program - New

The storage capacity is 50 Programs with 15 Steps/Program. Each Step contains Set Temperature, Gradient, and Hold Time.

The following sample assumes creating a Program #4 with 2 Steps.



Create a new Program
↵ Accept menu → 2.2.1
ESC returns to → 2

2.2.1 Create a new program



↶ ↷ Change value
↵ Accept value → 2.2.2
ESC returns to → 2.2

Note: Only free Program Numbers will appear in the display

2.2.2 Set Temperature – Step 1



↶ ↷ Change value
↵ Accept value → 2.2.3

ESC restores the old value or returns to → 2.2 and the step 1 as well as selected program number is not stored!

2.2.3 Set Gradient – Step 1



↶ ↷ Change value
↵ Accept value → 2.2.4

ESC restores the old value or returns to → 2.2 and the step as well as selected program number is not stored

Note: A value of 0.0 means maximal heating power

2.2.4 Set holding time (dwell time) – Step 1



↶ ↷ Change value
↵ Accept value and skips to the next input field
(HH→MM) → 2.2.5

ESC restores the old value and skips back one input field (MM→HH) or goes back to → 2.2 and the step 1 as well as selected program number is not stored!

Note: A value of 00:00 means endless holding time and program will not go to next step

2.2.5 Choose if a additional step is required



↶ ↷ Change answer
↵ Accept answer

If **Yes** selected: step number will increment with 1 → 2.2.6

If **No** selected: → 2.2.1

2.2.6 Set Temperature – Step 2



↶ ↷ Change value
↵ Accept value → 2.2.7

ESC restores the old value or **if in step 2 and higher** returns to → 2.2.5 but the actual step will not be saved!

2.2.7 Set Gradient – Step 2



↺ ↻ Change value
← Accept value → 2.2.8
ESC restores the old value or if in step 2 and higher returns to → 2.2.6 but the actual step will not be saved!

2.2.8 Set Holding Time – Step 2



↺ ↻ Change value
← Accept value and skips to the next input field (HH→MM) or stores the time and go to → 2.2.9
ESC restores the old value and skips back one input field (MM→HH) or if in step 2 and higher returns to → 2.2.7 but the actual step will not be saved!
Note: A value of 00:00 means endless holding time and program will not go to next step

2.2.9 Choose if a additional step is required



↺ ↻ Change answer
← Accept answer
If **Yes** selected: step number will increment with 1 → 2.2.2
If **No** selected: → 2.2.10

2.2.10 End of programming sequence



↺ ↻ Change answer
← Accept answer
If **Yes** selected: → 2.2.11
If **No** selected: → 2.2.9

2.2.11 Confirming & Saving the new program



Confirmation that the new program has been stored
After a few seconds Display will return to → 2.2

2.3 Menu Program- Edit

The example assumes to edit the program #4 with 2 steps



Edit an existing program
← Accept menu → 2.3.1
ESC returns to → 2

2.3.1 Choose program / step to edit



↺ ↻ Change value
← Accept value and skips to the next field (P→S) or → 2.3.2 (to the selected Step Number respective)
ESC returns to → 2.3
Note: Only stored programs/steps will appear in the display

2.3.2 Set Temperature – Step 1



- ↶ ↷ Change value
- ← Accept value → 2.3.3
- ESC restores the old value or returns to → 2.3 and the step 1 as well as selected program number is not stored!

2.3.3 Set Gradient – Step 1



- ↶ ↷ Change value
- ← Accept value → 2.3.4
- ESC restores the old value or returns to → 2.2 and the step as well as selected program number is not stored
- Note:** A value of 0.0 means maximal gradient

2.3.4 Set Hold Time (Dwell Time) – Step 1



- ↶ ↷ Change value
- ← Accept value and skips to the next input field (HH→MM) or stores the time and go to → 2.3.5
- ESC restores the old value and skips back one input field or goes back to → 2.3 and step 1 is not stored!
- Note:** Value of 00:00 means endless holding time and program will not go to next step

2.3.5 Set Temperature – Step 2



- ↶ ↷ Change value
- ← Accept value → 2.3.6
- ESC restores the old value or if in step 2 and higher returns to → 2.3 but changes to the actual step will not be saved!

2.3.6 Set Gradient – Step 2



- ↶ ↷ Change value
- ← Accept value → 2.3.7
- ESC restores the old value or if in step 2 and higher returns to → 2.3 but changes to the actual step will not be saved!

2.3.7 Set Holding Time – Step 2



- ↶ ↷ Change value
- ← Accept value and skips to the next input field (HH→MM) or stores the time and go to → 2.3.8 (MM→HH) or if in step 2 and higher returns to → 2.3 but changes to the actual step will not be saved !
- Note:** Value of 00:00 means endless holding time and program will proceed to next step

2.3.8 Choose if a additional step is required



- ↶ ↷ Change answer
- ← Accept answer
- If **Yes** selected: step number will increment with 1 → 2.3.2
- If **No** selected: → 2.3.11

2.3.9 End of programming sequence



↻ ↻ Change answer
← Accept answer
If **Yes** selected: → 2.3.11
If **No** selected: → 2.3.9

2.3.10 Confirming & Saving the new program



Confirmation that the edited program has been stored
After a few seconds Display will return to → 2.3

2.4 Menu Program - Delete



Delete an existing program
← Accept menu → 2.4.1
ESC returns to → 2.4

2.4.1 Choose program # to be deleted



↻ ↻ Change value
← Accept value → 2.4.2
ESC returns to → 2.4

Note: Only stored program numbers will appear in the display

2.4.2 Deletion selection will be displayed



↻ ↻ Change answer
← Accept answer
If **Yes** selected: → 2.4.3
If **No** selected: → 2.4

2.4.3 Deletion confirmation will be displayed



Confirmation that the selected program has been deleted
After a few seconds Display will return to → 2.4

2.5 Menu Program - Print



Print a program
← Accept menu → 2.5.1
ESC return to → 2

2.5.1 Choose program # to be printed



↻ ↻ Change value
← Accept value → 2.5.2
ESC return to → 2

Note: Only used program numbers will appear in the display

2.5.2 Displaying print in progress



Confirmation that the program has been printed
After a few seconds Display will return to → 2.5

3 Main Menu - Configuration

Configuration of the system by the user



In this menu point you can define and set system options
← Accept menu → 3.1

3.1 Select language



↶ ↷ Change language
← Accept language → 3.2

Attention:

After confirmation subsequent dialogs are displayed in the selected language

3.2 Set Actual Date for the Real Time Clock



↶ ↷ Change value
← Accept value and skips to the next input field (DD→MM→YY) or saves the date and goes to → 3.3
ESC restores the old value and/or skips back one input field (YY→MM→DD)

3.3 Set Actual Time for the Real Time Clock



↶ ↷ Change value
← Accept value and skips to the next input field (HH→MM) or saves the time and goes to → 3.4
ESC restores the old value and/or skips back one input field (MM→HH)

3.4 Set allowed max. Temperature



Set the maximal possible temperature value for manual operation (Main Menu / Temp & Option)

↶ ↷ Change value
← Accept value → 3.5
ESC restores value

3.5 Set Print interval for Printer Log via RS232 Interface



Set the Print Interval time. A value of 00:00 will disable printout of operating values

↶ ↷ Change value
← Accept value and skips to the next input field (HH→MM→SS) and goes to → 3.6
ESC restores the old value and/or skips back one input field (SS→MM→HH)

3.6 Set automatic interval to scroll status displays



Select if operation display will switch automatically instead of manual switching by ↶

↶ ↷ Change value
← Accept value → 3.7

3.7 Select Baud Rate for RS232 Interface



Available Baud Rate are 4800/9600/57600/115200
↻ ↻ Change value
← Accept value → **3.8**

3.8 Set Program End Buzzer



Buzzer will give an audio signal when program has finished
↻ ↻ Change answer
← Accept answer → **3.9**

3.9 Set Safety Alarm Buzzer



In a case of an over temperature alarm situation, the Buzzer will give an audio signal
↻ ↻ Change answer
← Accept answer → **3.10**

3.10 Set Offset between internal PT100 Sensor and actual Display



Offsets the internal PT100 sensor with the actual displayed temperature. Calibrate with an external temperature sensor.
↻ ↻ Change value
← Accept value → **3.11**
ESC restores value

3.11 Confirmation of Storing Entered Values



Confirmation that the Configuration has been Stored.
After a few seconds Display will return to → **3**

4 Menu Service Mode

<pre>Main Menu Service Mode</pre>	<p>↓</p>	<p>End user has only access to the History Data's Accept menu</p>
<pre>Enter Code 888</pre>	<p>↺ ↻ ↓ ↓</p>	<p>Change value to 8 Accept value Repeat until Display shows 888 Accept last value</p>
<pre>History 11</pre>	<p>↓</p>	<p>The number (11) indicates how many entries are stored Accept menu</p>
<pre>26.03.10 12:00 F14 ↓</pre>	<p>↺</p>	<p>Display shows in the first line Date and Time when the entry was stored. Second line shows the code for the Stored value/message Change to next entry</p>
<pre>22.03.10 08:12 F15 ↓</pre>	<p>↺ ↻ ↓</p>	<p>Change to next/previous entry Back to Main Display History</p>
<pre>History 11</pre>	<p>↺</p>	<p>Change to next input field</p>
<pre>Print Histories Yes No</pre>	<p>↺ ↻ ↓</p>	<p>Change answer Accept answer After a few seconds Display will return to → 1</p>

Status Displays

General

Change the display with \cup \cup or set Display Interval to yes in Menu Configuration to let change the display automatically.

The display will differ according to the installed options

5 Status Display: Manual Mode – Defined Start Date

If the programmed start date/time has achieved the systems starts. The display switches to \rightarrow 6

Start Date 01.01.10 14:59	Start Date and Time \cup to next set of operating display
Set Temp 37.0°C Gradient 0.0°C/M	Set Temperature and Gradient \cup to next set of operating display
Hold Time 02:00	Hold Time (Dwell Time) \cup returns to first set of operating display

6 Status Display: Manual Mode – Process Running

Set Temp 37.0°C Act Temp 37.0°C	Set Temperature and Actual Temperature \cup to next set of operating display
Gradient 0.0°C/M Hold Time 02:00	Gradient and Hold Time \cup to next set of operating display
01.01.10 14:59	Actual Date and Time \cup returns to first set of operating display

7 Status Display: Program mode - Defined Program Start Date

If the programmed Start Date/Time has achieved the system is running the selected program and the display switches to → 8




The LED „Prog“ is flashing as long as the defined start date/time has not been reached.

```
Start Date   P14
01.01.10    12:00
```

Start Date and Program number
ESC Cancels the Start Date → 9

8 Status Display: Program mode - Program Running



The LED „Prog“ is on during a programmed operation

```
Set Temp  37.0° C
Act Temp  37.0° C
```

Set Temperature and Actual Temperature
↻ to next set of operating display

```
Gradient  0.0° C/M
```

Gradient
↻ to next set of operating display

```
Hold Time  00:59
01.01.10   12:00
```

Hold Time and actual Date and Time
↻ to next set of operating display

```
Program    04
Step       01
```

Program Number and Step Number
↻ returns to first set of operating display

If the program has finished a buzzer will sound and the following message will be shown:

```
Message:
Program finished
```

← Accept message

9 Cancel a Running Process by ESC Button

Pressing ESC stops a running heating cycle / program

9.1 Safety question when stopping a running system



- ↶ ↷ Change answer
- ← Accept answer
- If **Yes** selected: → 9.2
- If **No** selected: → Back to running status

9.2 Conformation of cancelling a process



Confirmation that the running process has been stopped
After a few seconds display will return to Main Menu → 1

10 Messages and Errors

The messages and Errors are announced by a Buzzer and can be confirmed by pressing ←

Messages



Indicates that no program is in memory
← Accept message



Indicates that the program memory is full
← Accept message



Indicates that the door is open (manual process or starting a program)
← Accept message

Errors



PT100 Sensor or cable defect
← Accept message
Call Service!



Temperature exceeded security range level
← Accept message
Cool down oven
If error persist call Service!



← Accept message
Call Service!


Printer Operation

Connecting a Printer

Printer Requirements

In order to connect a printer with the serial RS232 Interface of the System, following requirements must be fulfilled : Serial RS232 Interface, Min. 4800 Baud Transfer Rate.

Pin Layout RS232 DB9 Connector System



Used Pins: 2:TxD , 3:RxD and 5:Signal Ground
Data format: 8 Data Bits, 1 Stop Bit, No Parity
No specific protocol

Examples for Different Printouts

Printout of a stored program

To print the content of stored programs go to Menu "Print Program" → 2.5

Program Nr::	1
Step::	1
Set Temperature:	37.0°C
Gradient:	0.0°C/Min
Hold time:	1:00HH:MM

Printout during a manual / program cycle

In order to get printouts for logging you must set the print interval (see "Print Interval" → 3.5) to a value greater than 00:00, otherwise it is disabled


Manual Operation	
Date:	01.01.10
Time:	12:00
Set Temperature:	37.0°C
Act Temperature:	37.0°C
Gradient:	0.0°C/Min

Programmed Operation	
Date:	01.01.10
Time:	12:00
Program:	4
Step:	2
Set Temperature:	37.0°C
Act Temperature:	37.0°C
Gradient:	0.0°C/Min
Hold time:	1:00HH:MM


Wiring Diagram ICN

- A1 = Regelplatine/PCB board
- A2 = Display Platine/Display board
- E1-5 = Heizung/Heating element
- F1 = Hauptsicherung/Main fuse
- F2 = Uebertemp Sich./Over temp fuse
- F3 = Uebertemp Sich./Over temp fuse
- F4 = Sicherung Regelplatine / PCB fuse
- S2 = Türschalter/Door switch
- X1 = Netzanschluss/Main power clamp
- X2 = Strom Regelplatine/Main power PCB
- X3 = RS232/RS232
- X4 = RS232 Ext./RS232 Ext.
- X5 = P1100/P1100
- X6 = Leer/Empty
- X7 = Türechscher/Door switch
- X8 = Heizung/Heating element
- X9 = Leer/Empty
- X11 = Heizung/Heating element
- X13 = Display/Display


Elektroschema ICN 52/118
Wiring diagram ICN 52/118




Drawing ICN52



Drawing ICN118



Drawing Spare Parts



Spare Part Numbers

Position	Part. No.	Description
1	31W04141005	Panel New Version
2	31W04144155	Display / Touch Panel 230V
3	31W04144154	Main PCB 230V / 10A
4	31W04144153	Main PCB 115V / 10A
5	31W04142002	Air Inlet Protection Cover
6	31W04144014	Door Handle
7	31W04144013	Cable Relieve
8	31W04144018	Exhaust Knob
9	31W04140003	Door Outer Case ICN52
10	31W04140001	Door Outer Case ICN118
11	31W04144017	Insulation for Exhaust
12	31W04144016	Closing Bolt
13	31W04960700	Door Switch
14	31W04140002	Door Inner Case ICN52
15	31W04140000	Door Inner Case ICN118
16	31W04144123	Temperature Probe PT100
17	31W04144002	Door Hinge
18	31W04962507	Over Temperature Fuse
19	31W04870018	Heating Element 230V/110W
20	31W04870019	Heating Element 230V/133W
21	31W04962923	Heating Element 115V/110W
22	31W04962924	Heating Element 115V/133W
23	31W04812003	Air Distribution Plate ICN52
24	31W04812020	Air Distribution Plate ICN118
25	31W04943203	Door Seal ICN52 1,5 Meter
26	31W04943203	Door Seal ICN118 2,0 Meter