



TNF (Human) ELISA Kit

Catalog Number KA0648

96 assays

Version: 02

Intended for research use only

www.abnova.com

Table of Contents

Introduction	3
Intended Use	3
Background	3
Principle of the Assay	3
General Information	4
Materials Supplied	4
Storage Instruction	4
Materials Required but Not Supplied	5
Precautions for Use	5
Assay Protocol	7
Reagent Preparation	7
Sample Preparation	8
Assay Procedure	9
Data Analysis	12
Calculation of Results	12
Performance Characteristics	12
Resources	13
Troubleshooting	13
References	14
Plate Layout	15

Introduction

Intended Use

For quantitative detection of human TNF- α in cell culture supernatants, human plasma (EDTA, heparin and citrate), serum, cerebrospinal fluid, urine, synovial fluid or other body fluids.

The assay will recognize both natural and recombinant Hu TNF- α .

Background

Tumor Necrosis Factor- α (TNF- α) is a non-glycosylated 17.5 kDa, 157 amino acid protein. TNF- α is a potent lymphoid factor and exerts cytotoxic effects on a wide range of tumor cells and other target cells. It is secreted by macrophages, monocytes, neutrophils, T-cells, and NK-cells following their stimulation by bacterial lipopolysaccharides. TNF- α has been suggested to play a pro-inflammatory role and has been detected in synovial fluid of patients with rheumatoid arthritis. Various pathological conditions are associated with the production of high levels of TNF- α . These include septic shock, cachexia (e.g. HIV, tuberculosis, cancer), autoimmune diseases, hepatitis, leukemia, myocardial ischaemia, organ transplantation rejection, multiple sclerosis, rheumatoid arthritis, and meningococcal septicemia. Annually, many people die from septic shock syndrome, triggered by TNF- α following complications from an infectious disease. In many cases elevated TNF- α serum levels predict a higher mortality.

Principle of the Assay

The TNF- α ELISA kit is an in vitro enzyme-linked immunosorbent assay for the quantitative measurement of human TNF- α in cell culture supernatants, serum, plasma, cerebrospinal fluid, urine, synovial fluid and other body fluids. This assay employs an antibody specific for human TNF- α coated onto a 96-well plate. Standards, samples and biotinylated anti-human TNF- α are pipetted into the wells. TNF- α present in a sample is captured by the antibody immobilized to the wells and by the biotinylated TNF- α specific detection antibody. After washing away unbound biotinylated antibody, HRP-conjugated streptavidin is pipetted to the wells. The wells are again washed. Following the second wash step, TMB substrate solution is added to the wells, resulting in color development proportional to the amount of TNF- α bound. The Stop Solution changes the color from blue to yellow, and the intensity of the color is measured at 450 nm.

General Information

Materials Supplied

List of component

Component	State	Amount
96 Well Plate with 12 Strips Break apart microtiter test strips each with TNF- α antibody coated single wells.	ready to use	1 frame
TNF- α Standard: Lyophilized & stabilized human TNF- α , reconstitute with Sample Diluent volume shown on the label.	Lyophilized	4 vials
Biotinylated TNF- α antibody	ready to use	10 ml
HRP-Conjugated Avidin	ready to use	12 ml
20x Wash solution concentrate (sufficient for 1000 ml): Dilute 1:20.	concentrated	50 ml
Dilution buffer	ready to use	100 ml
Stop solution: 0,9 N H ₂ SO ₄	ready to use	8 ml
TMB-Substrate	ready to use	8 ml

Storage Instruction

Reagent	Storage	Stability
TNF- α antibody coated 96 well plates with 12 strips. Break apart microtiter test strips each with 8 antibody coated single wells	Store at 2-8°C in closed aluminum bag with desiccant Strips which are not used must be stored in the re-sealable aluminum bag in humidity free and airtight conditions	3 months after opening
TNF- α Standard Lyophilized	Store at 2-8°C	Until date of kit expiry in lyophilized format. Unstable. Use immediately after dissolving. Keep on ice if not used within 1 hr after dissolving
Biotinylated antibody. Ready for use.	Store at 2-8°C Avoid contamination (Use clean sterile tips)	3 months after opening
HRP-Conjugated Avidin. Ready for use.	Store at 2-8°C Avoid contamination (Use clean sterile tips)	3 months after opening
Sample Diluent	Store at 2-8°C Avoid contamination (Use clean sterile tips or pipettes)	3 months after opening
20x Concentrated Wash Buffer	Store at 2-8°C To avoid crystal formation, wash buffer concentrate, may also be stored at Room Temperature.	Until expiry date

Diluted Wash Buffer	1x working dilution Bottles used for the working dilution should be cleaned regularly, discard cloudy solutions	1 week at room temperature or one month at 2-8 °C
TMB-Substrate Solution	Store at 2-8 °C, protected from light! Avoid contamination (Use clean sterile tips)	Until expiry date
Stop Solution	Store at 2-8 °C. May also be stored at Room Temperature	Until expiry date at room temperature

Materials Required but Not Supplied

- ✓ Micro plate reader capable of measuring absorbance at 450 nm.
- ✓ Precision pipettes to deliver 2 ul to 1 ml volumes.
- ✓ Multi-channel pipet (25 ul to 350 ul).
- ✓ Adjustable 1-25 ml pipettes for reagent preparation.
- ✓ 100 ml and 1 liter graduated cylinders.
- ✓ Absorbent paper.
- ✓ Distilled or de-ionized water.
- ✓ Log-log graph paper or computer and software for ELISA data analysis.
- ✓ Tubes to prepare standard or sample dilutions.
- ✓ Timer

Precautions for Use

- ✓ This kit has been configured for research use only and is not for diagnostic and clinical use.
- ✓ Caution: TMB substrate (Tetramethylbenzidine) and the Stop solution (H_2SO_4) are toxic or corrosive and should be handled with care. Use gloves during handling.
- Procedure Note
 - ✓ When not in use, kit components should be refrigerated. All reagents should be warmed to room temperature before use.
 - ✓ Microtiter plates should be allowed to come to room temperature before opening the foil pouches.
 - ✓ Once the desired number of strips has been removed, immediately reseal the pouch and store at 2 - 8 °C to maintain plate integrity. Protect from humidity.
 - ✓ Samples should be collected in pyrogen/endotoxin-free tubes.
 - ✓ Samples should be frozen if not analyzed shortly after collection. Avoid multiple freeze-thaw cycles of frozen samples. Thaw completely and mix well prior to analysis.

- ✓ When possible, avoid use of badly hemolyzed or lipemic sera. If large amounts of particulate matter are present, centrifuge or filter prior to analysis.
- ✓ It is recommended that all standards, controls and samples be run in duplicate.
- ✓ Samples that are >400 pg/ml should be diluted with Sample Diluent.
- ✓ When pipetting reagents, maintain a consistent order of addition from well-to-well. This ensures equal incubation times for all wells.
- ✓ Cover or cap all reagents when not in use.
- ✓ Do not use reagents after the kit expiration date.
- ✓ Read absorbances within 20 minutes of assay completion.
- ✓ In-house controls should be run with every assay. If control values fall outside pre-established ranges, the accuracy of the assay is suspect.
- ✓ All residual wash liquid must be drained from the wells by efficient aspiration or by decantation followed by tapping the plate forcefully on absorbent paper. Never insert absorbent paper directly into the wells.
- ✓ Because TMB Chromogen is light sensitive, avoid prolonged exposure to light. Also avoid contact between Stabilized Chromogen and metal, or color may develop.

Assay Protocol

Bring all reagents and samples to room temperature (18-25 °C) before use.


Reagent Preparation

- Antibody coated plate: Before opening the plastic pouch, determine the number of strips required to test the desired number of samples, plus 16 wells needed for running standards and blanks in duplicate. Remove non-used strips from the plate-frame and return them to the foil pouch containing the desiccant for up to 3 month at 2-8 °C.
- Dilution of test standard: Dissolve the lyophilised TNF- α standard with Sample Diluent volume shown on the label. TNF- α standard is unstable after dissolving. Use immediately or keep on ice if not used within 1 hr after dissolving.

To obtain a standard curve dilute it as follows:

- ✓ Add 300 μ l of TNF- α standard from kit standard tube containing 1000 pg/ml of TNF- α (Standard tube 1.)
- ✓ Add 150 μ l of Sample Diluent to all other 6 dilution tubes. Take 150 μ l from the first tube (1000 pg/ml) and start 2-fold serial dilutions in dilution tubes as described in the figure by mixing several times with the pipet in each tube (Total of 7 dilution tubes).
- ✓ 150 μ l of sample Diluent in tube 8 serves as zero standard (0 pg/ml).

300 μ l of TNF- α Standard from dissolved stock (1000 pg/ml)



- Amounts of the reagents needed to perform the test

	Reagents				
No of strips used (8 well each)	Biotinylated antibody 50 µl/well	Avidin-HRP 50 µl/well	TMB substrate 50 µl/well	Stop Solution 25 µl/well	Wash Buffer 300 µl/well
1 (8 wells)	500 µl	500 µl	500 µl	300 µl	30 ml
2 (16 wells)	1 ml	1 ml	1 ml	600 µl	55 ml
4 (32 wells)	2 ml	2 ml	2 ml	1.2 ml	110 ml
6 (48 wells)	3 ml	3 ml	3 ml	1.8 ml	165 ml
8 (64 wells)	4 ml	4 ml	4 ml	2.4 ml	220 ml
12 (96 wells)	6 ml	6 ml	6 ml	4 ml	350 ml

Sample Preparation

- Sample preparation and dilution: Dilution of samples is not required for initial screening. Samples that exceed the measuring range should be diluted in sample diluent serially 1:2, 1:4, or further if necessary, and measured again. The dilution factor must be taken in account when calculating the results.
Dilute and store all samples in tubes or plates made of material with low binding surface, such as polypropylene.
- Sample collection and storage: Serum, EDTA, heparin or citrate anti-coagulated plasmas, cerebrospinal fluid, urine, synovial fluid, other body fluids and cell culture supernatants are suitable for use in the assay (caution: separate plasma/serum and blood cells within 4 hours after collection, non-separated samples must be kept at 2-8°C). Do not use grossly haemolysed or lipemic specimens. If samples are to be run within 24 hours, they may be stored at 2-8°C; otherwise samples should be stored frozen (at least between -18 to -32°C, but preferably < -70°C). Up to 3 freeze-thaw cycles have no effect on the TNF-α levels of samples. Nonetheless, excessive freeze-thaw cycles should be avoided. Prior to the assay, frozen samples should be thawed as quickly as possible in tap water (18-25°C), do not use 37°C or 56°C water bath for this purpose.
- Preparation of reagents:
 - ✓ Wash Buffer: If the 20x concentrated Wash Buffer contains visible crystals, warm it at 37°C and mix gently until dissolved. Dilute 1:20 with de-ionized or distilled water (e.g. 25 ml of Wash Buffer Concentrate and 475 ml distilled water to yield 500 ml of 1x Wash Buffer). Check the pH of the diluted wash buffer and adjust to 7.4 if necessary.
 - ✓ Vortex mix Biotinylated antibody solution gently before use.
 - ✓ Vortex mix peroxidase (HRP) labeled avidin gently before use.

Caution: TMB substrate (Tetramethylbenzidine) and the Stop solution (H₂SO₄) are toxic or corrosive and should be handled with care. Use gloves during handling.

Assay Procedure

1. Bring all reagents and samples to room temperature (18 - 25°C) before use. It is recommended that all standards and samples be run at least in duplicate. Leave some wells as a reagent blank (2 to 4 wells).

FIRST STEP: STANDARD, SAMPLES AND BLANK + BIOTINYLATED ANTIBODY

2. Pipette 50 µl of sample and 50 µl of each diluted standard starting from 1000 pg/ml into appropriate wells. Pipette 50 µl of sample diluent to the wells which will be used as a blank. Incubate 1 hr at room temperature without shaking.

SECOND STEP: BIOTINYLATED ANTIBODY

3. Wash 5x with 1x Wash Solution (300 µl each)

To wash manually: Empty plate contents. Use a multi-channel pipette to fill each well with 300 µl of diluted wash buffer, then empty plate contents again. Repeat procedure 4 additional times for a total of FIVE washes. Gently blot plate onto paper towels or other absorbent material. Never let reaction wells dry. Continue to the next step without delay or interruption.

For automated washing: Aspirate all wells and wash 5 times with 300 µl diluted wash buffer. Blot plate onto paper towels or other absorbent material. Never let reaction wells dry. Continue to the next step without delay or interruption.

4. Promptly add 50 µl of green colored Biotinylated TNF-alpha detection antibody to all wells. Tap the plate gently by hand to homogenize your mixture. Avoid touching to the reaction wells with the pipette tip. Incubate at room temperature for 30 minutes without shaking.

THIRD STEP: HRP-CONJUGATED AVIDIN

5. Wash 5 times 5x as described in Step 3.

Add 50 µl of prepared HRP-conjugated avidin solution (ready to use) to each well. Incubate for 30 minutes at room temperature.

FOURTH STEP: TMB SUBSTRATE


6. Wash 5 times as described in Step 3.
7. Using a multichannel pipette, promptly add 50 µl of TMB ready to use substrate reagent to each well. Incubate for 20 minutes at room temperature in the dark.
8. Add 25 µl of Stop Solution to each well. Read at 450 nm within 15 minutes.
Correcting for optical imperfections in the microplates by subtracting A_{630} nm is recommended, but not an essential procedure.

FIFTH STEP: READING AND CALCULATION

9. Calculate the mean of reagent blank absorbance values and subtract it from all test well values (standard and test samples). Mean reagent blank absorbance value at 450 nm should be less than 0.200.
10. Calculate your results against standard curve.


- Summary
- 1. Prepare all reagents, samples and standards. Dilution of samples not required at initial screening.
- 2. Add 50 µl standard (starting from 1000 pg/ml), test samples and sample diluent as a blank into the appropriate wells of the strips. Incubate 1 hour at room temperature. Wash 5x.
- 3. Add 50 µl ready for use biotin antibody promptly to each well. Incubate 30 min. at room temperature. Wash 5x.
- 4. Add 50 µl ready for use HRP-Streptavidin solution. Incubate 30 minutes at room temperature. Wash 5x.
- 5. Add 50 µl TMB Substrate Solution to each well. Incubate 20 minutes at room temperature.
- 6. Add 25 µl Stop Solution to each well. Read at 450 nm against *630 nm immediately.
Subtract blank values from values for standards and test samples.

*Correcting for optical imperfections in the microplates by subtracting $A_{630\text{ nm}}$ is recommended, but not an essential procedure.




STEP 1

Add 50 μ L of standards and test samples to test well and incubate 1 hr at Room Temperature




STEP 2

Add 50 μ L of Biotinylated TNF- α antibody to test well and incubate 30 min. at Room Temperature





STEP 3

Add 50 μ L of HRP-Streptavidin to test well and incubate 30 minutes at Room Temperature



STEP 4

Add 50 μ L of TMB substrate to test well incubate 20 minutes at Room Temperature



Data Analysis

Calculation of Results

The standard curve must be determined individually for each experiment. Correct the absorbance values of all standards by subtracting from them the mean O.D. value of the reagent blank (B1 = only sample diluent). Calculate the mean absorbance value for each standard from the duplicates.

The standard curve is used to determine the amount of TNF alpha in an unknown sample. The standard curve is generated by plotting the average O.D. (450 nm) obtained for each of the standard concentrations on the vertical (Y) axis versus the corresponding TNF alpha concentration (pg/mL) on the horizontal (X) axis.

Construct the standard curve using graph paper or statistical software.

If samples generate values higher than the highest standard, dilute the samples with sample diluent and repeat the assay. Note that the concentration read from the standard curve must be multiplied by the dilution factor.

Performance Characteristics

	TNF- α
Assay range	0-1000 pg/ml
Standard curve points	1000, 500, 250, 125, 62.5, 31.25, 15.6 and 0 pg/ml.
Intra-Assay-Precision	$\leq 6\%$
Inter-Assay-Precision	$\leq 4\%$
Inter-Lot-Precision	$\leq 8\%$
Cross-Reactivity	No cross reactivity was observed with the following recombinant human proteins: IL-1 β , IL-1 α , IL-2, IL-3, IL-4, IL-5, IL-6, IL-7, IL-8, IL-9, IL-10, IL-12, IL-13, TARC
Interferences	No interferences to bilirubin up to 0.3 mg/mL, haemoglobin up to 8.0 mg/mL and triglycerides up to 5.0 mg/mL
Specificity	Recognizes both natural and recombinant human TNF- α .
Sensitivity	<15 pg/ml

TNF- α levels may vary greatly between different study groups and sample types (such as serum samples, cell culture supernatant, cell extracts or other biological samples). Each research study should include a proper control group (age, sex, locality or geographical region matched) to establish more precise TNF- α values. Disease status or the use of drugs or TNF- α stimulating agents may interfere with the TNF- α levels and should be taken into careful consideration in all studies.

Resources

Troubleshooting

Problem	Cause	Solution
Poor standard Curve	1. Inaccurate pipetting or pipetting error 2. Improper standard dilution	Check pipettes and calibrate regularly. Vortex the stock before use and dilute carefully in an eppendorf tube.
Low signal	1. Shorter incubation than recommended 2. Inadequate reagent volumes or improper dilution or pipetting error	Ensure sufficient incubation time; Check pipettes and ensure correct performance.
Large CV	Inaccurate pipetting and drying of wells during test procedure.	Check pipettes Fill the wells promptly with wash buffer and reagents.
High background	1. Plate is insufficiently washed 2. Contaminated wash Buffer 3. Wash buffer volume is less than advised	Review the manual for proper wash. If using a plate washer, check that all ports are unobstructed and clean. Make a fresh wash buffer Use 300µl per well
Low sensitivity	1. Improper storage of the ELISA kit 2. Stop solution 3. Contamination of reagents	Store test kit components as advised in this user manual. Keep substrate solution protected from light. Stop solution should be added to each well before measure. Use clean sterile tips. Discard contaminated reagents.

References

1. Seriola B, Paolino S, Sulli A, Fasciolo D, Cutolo M.(2006). Effects of anti-TNF- α treatment on lipid profile in patients with active rheumatoid arthritis. *Ann N Y Acad Sci.* 1069:414-9.
2. Intiso D, Zarrelli MM, Lagioia G, Di Rienzo F, Checchia De Ambrosio C, Simone P, Tonali P, Cioffi Dagger RP. (2004). Tumor necrosis factor α serum levels and inflammatory response in acute ischemic stroke patients. *Neurol Sci.* 24:390-396.
3. Beutler, B. et al. (1987) Cachectin: more than a tumor necrosis factor. *N. Engl. J. Med.* 316:379-385.
4. Tracey, K.J. et al. (1987) Anti cachectin/TNF monoclonal antibodies prevent septic shock during lethal bacteraemia. *Nature* 330:662-664.
5. Piguet, P.F. et al. (1987) Tumor necrosis factor/cachectin is an effector of skin and gut lesions of the acute phase of graft-versus-host
6. disease. *J. Exp. Med.* 166:1280-1289.
7. Aukrust, P. et al. (1994) Serum levels of tumor necrosis factor- α (TNF- α) and soluble TNF receptors in human immunodeficiency virus type 1 infection-correlations to clinical immunologic, and virologic parameters. *J. Inf. Dis.* 169:420-424.
8. Waage, A. et al. (1987) Association between tumor necrosis factor in serum and fatal outcome in patients with meningococcal disease. *Lancet* 1:355-357.
9. Noraz, N. et al. (1997) Human cytomegalovirus-associated immunosuppression is mediated through IFN- α . *Blood* 89(7):2443-2452.
10. Lekutis, C. et al. (1998) HIV-1 env DNA vaccine administered to rhesus monkey elicits MHC class II-restricted CD4+ T helper cells that secrete IFN- γ and TNF- α . *J. Immunol.* 158:4471-4477.
11. Neuman, M.G. et al. (1998) Role of cytokines in ethanol-induced cytotoxicity in HepG2 cells. *Gastroenterology* 114(7):157-169.
12. Dong, Z. et al. (1998) Activation of cytokine production, tumoricidal properties, and tyrosine phosphorylation of MAPKs in human monocytes by a new synthetic lipopeptide, JBT3002. *J. Leukocyte Biol.* 63:766-774.
13. Murato, P.A. et al. (1997) Immunodominance of a low-affinity major histocompatibility complex-binding myelin basic protein epitope (Residues 111-129) in HLA-DR4 (B1*0401) Subjects is associated with a restricted T cell receptor repertoire. *J. Clin. Invest.* 100(2):339-349.
14. Ludviksson, B.R. et al. (1998) Active Wegener's granulomatosis is associated with HLA-DR+ CD4+ T cells exhibiting an unbalanced Th1-type T cell cytokine pattern: Reversal with IL-10. *J. Immunol.* 160:3602-3609.

Plate Layout

[illegible]