


Home 	»
Frontline Workhorse	
Frontline Finesse	
Extra Support	
Specialty	
Guide Wire Technology	
Support/Tip Stiffness	
Ordering Information	
<b>Important Safety Information</b>	
referenced within	

**Indications:**

This HI-TORQUE Guide Wire is intended to facilitate the placement of balloon dilatation catheters during Percutaneous Transluminal Coronary Angioplasty (PTCA) and Percutaneous Transluminal Angioplasty (PTA). This guide wire may also be used with compatible stent devices during therapeutic procedures.




# Leading the Way in Guide Wire Technology



Comprehensive  
portfolio of guide wires

Advanced  
Technology

## Abbott Vascular Guide Wires



**HI-TORQUE**  
Confidence in your hands.





Home 

Frontline Workhorse

Frontline Finesse

Extra Support

Specialty

Guide Wire Technology

Support/Tip Stiffness

Ordering Information

**Important Safety Information**  
referenced within

**R<sub>ONLY</sub> HI-TORQUE Guide Wires for PTCA, PTA, and Stents**

**INDICATIONS FOR USE**

This HI-TORQUE Guide Wire is intended to facilitate the placement of balloon dilatation catheters during Percutaneous Transluminal Coronary Angioplasty (PTCA) and Percutaneous Transluminal Angioplasty (PTA). This guide wire may also be used with compatible stent devices during therapeutic procedures.

**CONTRAINDICATIONS**

Not intended for use in the cerebral vasculature or with atherectomy devices.

**WARNINGS**

**This device is designed and intended for ONE-TIME USE ONLY. Do not resterilize and / or reuse.**

**Carefully observe the instructions under “Do Not” and “Do” below. Failure to do so may result in vessel trauma, guide wire damage, guide wire tip separation, or stent damage. If resistance is observed at any time, determine the cause under fluoroscopy and take remedial action as needed.**

**DO NOT:**

- Push, auger, withdraw, or torque a guide wire that meets resistance.
- Torque a guide wire if the tip becomes entrapped within the vasculature.
- Allow the guide wire tip to remain in a prolapsed condition.
- Deploy a stent such that it will entrap the wire between the vessel wall and the stent

**DO:**

- Advance or withdraw the guide wire slowly.
- Use the radiopaque marker of the interventional device to confirm position.
- Examine the tip movement under fluoroscopy before manipulating, moving, or torquing the guide wire.
- Observe the wire under fluoroscopy for tip buckling, which is a sign of resistance.
- Maintain continuous flush while removing and reinserting the guide wire to prevent air from entering the catheter system. Perform exchanges slowly to prevent air entry and / or trauma.
- When reintroducing the guide wire, confirm that the interventional device tip is free within the vessel lumen and that the tip is parallel with the vessel wall.
- Use extreme caution when moving a guide wire through a non-endothelialized stent, or through stent struts into a bifurcated vessel. Use of this technique carries additional patient risks, including the risk that the wire may become caught on the stent strut.

**For PROGRESS family only: The PROGRESS family of guide wires have distal ends of varying stiffness. Operate these guide wires carefully so as to not injure the blood vessel, observing the information in these instructions. The higher torque performance, stiffer distal ends, and / or higher advancement force may present a higher risk of perforation or injury than a guide wire with a more**

**pliable distal end. Therefore, use the guide wire with the least stiff distal end that will treat the lesion, and use extreme care to minimize the risk of perforation or other damage to blood vessels. Use the most suitable guide wire for the lesion being treated.**

**PRECAUTIONS**

Guide wires are delicate instruments and should be handled carefully. Prior to use and when possible during the procedure, inspect the guide wire carefully for bends, kinks, or other damage. Do not use damaged guide wires. Using a damaged guide wire may result in vessel damage and / or inaccurate torque response.

Confirm the compatibility of the guide wire diameter with the interventional device before actual use.

Free movement of the guide wire within the interventional device is an important feature of a steerable guide wire system because it gives the user valuable tactile information. Test the system for any resistance prior to use. Adjust or replace the hemostatic valve with an adjustable valve if it is found to inhibit guide wire movement.

Never attach the torque device to the modified portion of the proximal end of the extendable guide wire; otherwise, guide wire damage may occur, preventing the ability to attach the DOC Guide Wire Extension.

HI-TORQUE Guide Wires with Hydrophilic Coating: Avoid abrasion of the hydrophilic coating. Do not withdraw or manipulate the hydrophilic-coated wire in a metal cannula or sharp-edged object.

**R<sub>ONLY</sub> HI-TORQUE® Guide Wires**

**INDICATIONS**

Refer to the device label for any additional product-specific indications that may apply.

**CONTRAINDICATIONS**

HI-TORQUE Guide Wires are not intended for use in the cerebral vasculature.

Refer to the device label for any additional product-specific contraindications that may apply.

**WARNINGS**

This device is designed and intended for ONE-TIME USE ONLY. Do not resterilize and/or reuse.

**Carefully observe the instructions under “Do Not” and “Do” below. Failure to do so may result in vessel trauma, guide wire damage, guide wire tip separation, or stent damage. If resistance is observed at any time, determine the cause under fluoroscopy and take remedial action as needed.**

**Do Not:**

- Push, auger, withdraw or torque a guide wire that meets resistance.
- Torque a guide wire if the tip becomes entrapped within the vasculature.
- Allow the guide wire tip to remain in a prolapsed condition.

**Do:**

- Advance or withdraw the guide wire slowly.
- Use the radiopaque marker of the interventional device to confirm position.

- Examine the tip movement under fluoroscopy before manipulating, moving or torquing the guide wire.
- Observe the wire under fluoroscopy for tip buckling, which is a sign of resistance.
- Maintain continuous flush while removing and reinserting the guide wire to prevent air from entering the catheter system. Perform exchanges slowly to prevent air entry and / or trauma.
- When reintroducing the guide wire, confirm that the interventional device tip is free within the vessel lumen and that the tip is parallel with the vessel wall.
- Use extreme caution when moving a guide wire through a non-endothelialized stent, or through stent struts into a bifurcated vessel. Use of this technique carries additional patient risks, including the risk that the wire may become caught on the stent strut.
- Consider that if a secondary wire is placed in a bifurcation branch, this wire may need to be retracted prior to stent deployment because there is additional risk that the secondary wire may become entrapped between the vessel wall and the stent.

**PRECAUTIONS**

Guide wires are delicate instruments and should be handled carefully. Prior to use and when possible during the procedure, inspect the guide wire carefully for bends, kinks, or other damage. Do not use damaged guide wires. Using a damaged guide wire may result in vessel damage and / or inaccurate torque response.

Confirm the compatibility of the guide wire diameter with the interventional device before actual use.

Free movement of the guide wire within the interventional device is an important feature of a steerable guide wire system because it gives the user valuable tactile information. Test the system for any resistance prior to use. Adjust or replace the hemostatic valve with an adjustable valve if it is found to inhibit guide wire movement.

Never attach the torque device to the modified portion of the proximal end of the extendable guide wire; otherwise, guide wire damage may occur, preventing the ability to attach the DOC Guide Wire Extension.

HI-TORQUE Guide Wires with Hydrophilic Coating: Avoid abrasion of the hydrophilic coating. Do not withdraw or manipulate the hydrophilic-coated wire in a metal cannula or sharp-edged object

**R<sub>ONLY</sub> ASAHI® PTCA Guide Wires**

**INDICATIONS FOR USE**

ASAHI PTCA guide wires are intended to facilitate the placement of balloon dilatation catheters during percutaneous transluminal coronary angioplasty (PTCA) and percutaneous transluminal angioplasty (PTA).

The ASAHI PTCA guide wire is not to be used in the cerebral blood vessels.

**WARNINGS**

- This guide wire is for single use only. Do not reuse or resterilize. If reused or resterilized, the performance or quality of the guide wire may be compromised and there is a risk of complications, including infection.

- Do not use the guide wire after the expiration date indicated on the label. Discard any guide wire that exceeds the expiration date.
- This guide wire must be used only by a physician who is fully trained in PTCA/PTA treatment.
- The coil section is especially fragile, so do not bend or pull it more than necessary. Otherwise, the guide wire may be damaged.
- Do not use a damaged guide wire. Using a damaged guide wire may result in blood vessel damage and/or inaccurate torque response. Injury to the patient may result.
- Always advance and withdraw the guide wire slowly.
- Observe guide wire movement in the vessels. Before a guide wire is moved or torqued, the tip movement should be examined and monitored under fluoroscopy. Do not move or torque a guide wire without observing corresponding movement of the tip; otherwise, the guide wire may be damaged and/or vessel trauma may occur. In addition, ensure that the distal guide wire tip and its location in the vessel are visible during wire manipulations.
- Never push, auger, withdraw, or torque a guide wire that meets resistance. Torquing or pushing a guide wire against resistance may cause guide wire damage and/or guide wire tip separation or direct damage to a vessel. Resistance may be felt and/or observed under fluoroscopy by noting any buckling of the guide wire tip. If guide wire tip prolapse is observed, do not allow the tip to remain in a prolapsed position; otherwise damage to the guide wire may occur. Determine the cause of resistance under fluoroscopy and take any necessary remedial action.
- If any resistance is felt due to spasm or the guide wire being bent or trapped while operating the guide wire in the blood vessel or removing it, do not move or torque the guide wire. Stop the procedure. Determine the cause of resistance under fluoroscopy and take appropriate remedial action. If the guide wire is moved excessively, it may break or become damaged, which may cause blood vessel injury or result in fragments being left inside the vessel.
- When torquing this guide wire inside the blood vessel, do not torque continuously in the same direction. This may cause the guide wire to become damaged or break apart, causing injury to the blood vessel or leaving fragments inside the vessel. When torquing the guide wire, rotate it clockwise and counterclockwise alternately. Do not exceed two rotations (720°) in the same direction.





Home 

Frontline Workhorse

Frontline Finesse

Extra Support

Specialty

Guide Wire Technology

Support/Tip Stiffness

Ordering Information

**Important Safety Information**  
 referenced within

- Do not push the guide wire more than is necessary to advance the tip through the narrowed part of the vessel. (For example, do not push the guide wire when the distal tip of the guide wire is bent by the force of manipulation.) After crossing the targeted area, do not roughly twist, push or pull the guide wire. If the guide wire is moved excessively, it may be damaged or break apart, which may injure the blood vessel or leave fragments inside the vessel.
- Use proper technique to ensure and verify that no air enters the interventional device when pulling the guide wire from the interventional device or re-inserting it. Otherwise air embolism could occur.
- Flush the guide wire with heparinized saline or other suitable solution while removing and reinserting it to prevent air from entering the interventional device. Perform guide wire exchanges carefully to prevent air entry and/or trauma. When reintroducing the guide wire, confirm that the interventional device tip is free within the vessel lumen and is not against the vessel wall. Failure to do so may result in vessel trauma when the guide wire is removed. Use the radiopaque marker of the interventional device to confirm position.
- Free movement of the guide wire within the interventional device is an important feature of a steerable guide wire system because it gives the user valuable tactile information. Test the system for any resistance prior to use. Adjust or replace the hemostatic valve with an adjustable valve if it is found to inhibit guide wire movement.
- Do not use in areas of vessel that are not or cannot be visualized.
- **For ASAHI CONFianza and ASAHI MIRACLEbros Series Only: ASAHI CONFianza and ASAHI MIRACLEbros Series have stiff distal ends. Operate these guide wires carefully so as not to injure the blood vessel, observing information in these instructions. The higher torque performance, stiffer ends, and/or higher advancement force may present a higher risk of perforation or injury than if using a more flexible guide wire. Therefore, use the most flexible guide wire that will treat the lesion (i.e., the guide wire with the smallest flexibility number that will treat the lesion), and take due care to minimize the risk of perforation or other damage to blood vessels.**

#### FOR ALL GUIDE WIRES

Use the most suitable guide wire that will treat the lesion.

There are patient risks when using any guide wire including those that may result from damage to, or breakage of, the guide wire. If guide wire damage or breakage occurs, it may cause damage to the vessel and injury to the patient, or death. Accordingly, care should be taken that all persons who operate guide wires are properly trained in their use, that they observe proper technique, and that guide wires are used carefully in accordance with the Instructions for Use.

#### PRECAUTIONS

- If the package is opened or damaged, do not use the product. Do not to open the package until just prior to use. Use aseptic technique in handling and using the guide wire.
- Contraindications, warnings, precautions, and intended uses of interventional devices that are compatible with ASAHI PTCA guide wires are described in the user manuals supplied with the respective interventional devices. Before using an ASAHI PTCA guide wire with other interventional devices, read the user manual of the other devices to ensure the other devices are compatible with ASAHI PTCA guide wire. Ensure you choose the correct ASAHI PTCA guide wire and that its use is consistent with the contraindications, warnings, precautions, and Instructions for Use of both the other devices and ASAHI PTCA guide wire.
- Guide wires are delicate instruments and should be handled carefully. When taking the guide wire out of the holder tube, do not handle the guide wire roughly or pull it out abruptly.
- Inspect the guide wire carefully for bends, kinks, or other damage prior to use and whenever possible during the procedure.
- Take due care when using the guide wire to prevent bending or kinking, and stay within standard practice when using the guide wire.
- When shaping the distal end, use the minimum force needed so that the coil is not damaged. Inspect the coil and guide wire for damage after shaping and before using.
- Verify which is the distal end before insertion and be sure to insert the distal end (coiled end).
- Do not wipe this guide wire using an organic solution such as alcohol.

#### ADVERSE EVENTS OF GUIDE WIRE USE INCLUDE, BUT ARE NOT LIMITED TO:

- Failure to cross a lesion • Separation or breakage of the guide wire • Damage to a vessel, including possible vessel perforation

- Coronary artery dissection • Cardiac tamponade due to coronary artery perforation
- Air embolism • Infection • Coronary artery spasm • Coronary artery thrombus
- Hemotoma at puncture site • Cardiac perforation

#### HI-TORQUE ADVANCE™ ONLY Guide Wires

CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. OBSERVE ALL WARNINGS AND PRECAUTIONS NOTED THROUGHOUT THESE INSTRUCTIONS. FAILURE TO DO SO MAY RESULT IN COMPLICATIONS.

#### INDICATIONS FOR USE

The HI-TORQUE ADVANCE™ Guide Wires are intended to facilitate the placement of interventional percutaneous transluminal coronary angioplasty (PTCA) and percutaneous transluminal angioplasty (PTA) catheters, and other interventional devices including: intravascular stents, intravascular ultrasound devices and intravascular drug eluting stents.

#### CONTRAINDICATIONS

The HI-TORQUE ADVANCE™ Guide Wires are not intended for use in the cerebral vasculature.

#### WARNINGS

- This device is designed and intended for ONE TIME USE ONLY. DO NOT RSTERILIZE AND/OR REUSE.
- **Carefully observe the instructions under “Do Not” and “Do” below. Failure to do so may result in vessel trauma, guide wire damage, guide wire tip separation, or stent damage. If resistance is observed at any time, determine the cause under fluoroscopy and take remedial action as needed.**  
**Do Not:**
  - Push, auger, withdraw or torque a guide wire that meets resistance.
  - Torque a guide wire if the tip becomes entrapped within the vasculature.
  - Allow the guide wire tip to remain in a prolapsed condition.**Do:**
  - Advance or withdraw the guide wire slowly.
  - Use the radiopaque marker of the interventional device to confirm position.
  - Examine the tip movement under fluoroscopy before manipulating, moving or torquing the guide wire.

- Observe the wire under fluoroscopy for tip buckling, which is a sign of resistance.
- Maintain continuous flush while removing and reinserting the guide wire to prevent air from entering the catheter system. Perform exchanges slowly to prevent air entry and/or trauma.
  - When reintroducing the guide wire, confirm that the interventional device tip is free within the vessel lumen and that the tip is parallel with the vessel wall.
- Use extreme caution when moving a guide wire through a non-endothelialized stent, or through stent struts into a bifurcated vessel. Use of this technique carries additional patient risks, including the risk that the wire may become caught on the stent strut.
- Consider that if a secondary wire is placed in a bifurcation branch, this wire may need to be retracted prior to stent deployment because there is additional risk that the secondary wire may become entrapped between the vessel wall and the stent.

#### PRECAUTIONS

- Guide wires are delicate instruments and should be handled carefully. Prior to use and when possible during the procedure, inspect the guide wire carefully for bends, kinks, or other damage. Do not use damaged guide wires. Using a damaged guide wire may result in vessel damage and/or inaccurate torque response.
- Confirm the compatibility of the guide wire diameter with the interventional device before actual use.
- Free movement of the guide wire within the interventional device is an important feature of a steerable guide wire system because it gives the user valuable tactile information. Test the system for any resistance prior to use. Adjust or replace the hemostatic valve with an adjustable valve if it is found to inhibit guide wire movement.
- Never attach the torque device to the modified portion of the proximal end of the extendable guide wire; otherwise, guide wire damage may occur, preventing the ability to attach the DOC® Guide Wire Extension.
- Avoid abrasion of the hydrophilic coating.
- Do not withdraw or manipulate the hydrophilic-coated wire in a metal cannula or sharp-edged object.

Abbott Vascular  
 3200 Lakeside Dr., Santa Clara, CA 95054 USA, Tel: 1.800.227.9902

Caution: These products are intended for use by or under the direction of a physician. Prior to use, it is important to read the package insert thoroughly for Instructions for Use, Warnings and Potential Complications associated with use of these devices.

All tests performed by and data on file at Abbott Vascular. All tip stiffness values are averages, based on a sample of 5 guide wires, tested with an automated machine in a controlled environment. All Device Support Level values are test results on a single guide wire under the same ambient conditions. All illustrations included are artist's renditions. Not drawn to scale.

www.AbbottVascular.com  
 ©2012 Abbott. All rights reserved. AP2937584-US Rev. A 10/12


Technical Data are on file with ASAHI. Manufactured by:  
 ASAHI Intecc Co., Ltd. 3-100 Akatsu-cho, Seto, Aichi 489-0071, Japan.  
 ASAHI is a registered trademark of ASAHI Intecc Co., LTD.

 ASAHI INTECC

 **Abbott**  
Vascular



Home
Frontline Workhorse
Frontline Finesse
Extra Support
Specialty
Guide Wire Technology
Support/Tip Stiffness
Ordering Information
Important Safety Information referenced within



**Frontline Workhorse** guide wires are designed to be flexible, have a soft tip, and provide enough support to deliver most interventional devices

NITINOL (ELASTINITE®)

HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL II®

TIP STIFFNESS	CORE	RADIOPACITY	COATING
0.7 g	ELASTINITE®	3 cm	HYDROPHILIC

- Exposed tip coils
- Soft shaping ribbon tip
  - Intermediate polymer sleeve for excellent device interaction



HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL

TIP STIFFNESS	CORE	RADIOPACITY	COATING
0.6 g	ELASTINITE®	3 cm	HYDROPHILIC

- Exposed tip coils
- Soft shaping ribbon tip
  - Intermediate polymer sleeve for excellent device interaction



HI-TORQUE BALANCE MIDDLEWEIGHT


TIP STIFFNESS	CORE	RADIOPACITY	COATING
0.7 g	ELASTINITE®	3 cm	HYDROPHILIC OR HYDROPHOBIC

- Soft shaping ribbon tip





Home
Frontline Workhorse
Frontline Finesse
Extra Support
Specialty
Guide Wire Technology
Support/Tip Stiffness
Ordering Information
Important Safety Information referenced within



**Frontline Workhorse** guide wires are designed to be flexible, have a soft tip, and provide enough support to deliver most interventional devices

STAINLESS STEEL (DURASTEEL or TruTorq®)

HI-TORQUE ADVANCE®

TIP STIFFNESS	CORE	RADIOPACITY	COATING
1.0 g	DURASTEEL	3 cm	HYDROPHILIC
• Core-to-tip for precise steering		• RESPONSEASE™ transitionless grind	



ASAHI PROWATER®

TIP STIFFNESS	CORE	RADIOPACITY	COATING
1.1 g	TruTorq® Steel	3 cm	HYBRID
• Core-to-tip for precise steering			



17 cm  
Hydrophilic  
Coating

3 cm  
Hydrophobic  
Coating

ASAHI PROWATERFLEX®

TIP STIFFNESS	CORE	RADIOPACITY	COATING
0.8 g	TruTorq® Steel	3 cm	HYBRID
• Core-to-tip for precise steering			



17 cm  
Hydrophilic  
Coating

3 cm  
Hydrophobic  
Coating


ASAHI SOFT

TIP STIFFNESS	CORE	RADIOPACITY	COATING
1.0 g	TruTorq® Steel	3 cm	HYDROPHOBIC
• Core-to-tip for precise steering			




HI-TORQUE FLOPPY™ II


TIP STIFFNESS	CORE	RADIOPACITY	COATING
0.4 g	Stainless Steel	2 cm	HYDROPHILIC
• Soft shaping ribbon tip			





[Home](#) 

[Frontline Workhorse](#)

[Frontline Finesse](#) 

[Extra Support](#)

[Specialty](#)

[Guide Wire Technology](#)

[Support/Tip Stiffness](#)

[Ordering Information](#)

[Important Safety Information](#)  
referenced within



**Frontline Finesse** guide wires are designed to be flexible, with a polymer cover and soft tip, designed for outstanding deliverability and distal access in tortuous anatomy

**HI-TORQUE WHISPER® MS**

TIP STIFFNESS	CORE	RADIOPACITY	COATING
1.0 g	DURASTEEL	3 cm	HYDROPHILIC

- Full polymer cover
- RESPONSEASE™ transitionless core
  - Core-to-tip for steerability



**HI-TORQUE WHISPER® EXTRA SUPPORT**

TIP STIFFNESS	CORE	RADIOPACITY	COATING
1.2 g	DURASTEEL	3 cm	HYDROPHILIC

- Full polymer cover
  - Core-to-tip for steerability
- RESPONSEASE™ transitionless core
  - Enhanced core profile for increased support as compared with HT Whisper® MS



**HI-TORQUE PILOT® 50**

TIP STIFFNESS	CORE	RADIOPACITY	COATING
1.5 g	DURASTEEL	3 cm	HYDROPHILIC

- Full polymer cover
- RESPONSEASE™ transitionless core
  - Core-to-tip for steerability



**ASAHI FIELDER®**

TIP STIFFNESS	CORE	RADIOPACITY	COATING
3.7 g	TruTorq® Steel	3 cm	HYDROPHILIC

- Full polymer cover
- Core-to-tip for steerability





Home	Home
Frontline Workhorse	
Frontline Finesse	
Extra Support	>>
Specialty	
Guide Wire Technology	
Support/Tip Stiffness	
Ordering Information	
Important Safety Information	referenced within



**Extra Support** guide wires are designed to provide additional support for delivery of bulky devices

HI-TORQUE BALANCE HEAVYWEIGHT

TIP STIFFNESS	CORE	RADIOPACITY	COATING
0.7 g	ELASTINITE®	4.5 cm	HYDROPHILIC

- Soft shaping ribbon tip



HI-TORQUE ALL STAR

TIP STIFFNESS	CORE	RADIOPACITY	COATING
0.7 g	Stainless Steel	3 cm	HYDROPHOBIC

- Core-to-tip design
- Jointless spring coil
- Intermediate polymer sleeve for excellent device interaction



ASAHI GRAND SLAM®

TIP STIFFNESS	CORE	RADIOPACITY	COATING
0.8 g	TruTorq® Steel	4 cm	HYDROPHOBIC

- Core-to-tip design
- Extremely supportive for bulky device delivery




HI-TORQUE IRON MAN

TIP STIFFNESS	CORE	RADIOPACITY	COATING
1.0 g	Stainless Steel	3 cm	HYDROPHOBIC

- Core-to-tip design






[Home](#) 

[Frontline Workhorse](#)

[Frontline Finesse](#)

[Extra Support](#)

[Specialty](#) 

[Guide Wire Technology](#)

[Support/Tip Stiffness](#)

[Ordering Information](#)

[Important Safety Information](#)  
referenced within



**Specialty** guide wires are designed to provide increasing tip stiffness and excellent torque response for superb crossing performance

FULL POLYMER COVER

HI-TORQUE PILOT® 200

TIP STIFFNESS	CORE	RADIOPACITY	COATING
4.1 g	DURASTEEL	3 cm	HYDROPHILIC


- Full polymer cover
- RESPONSEASE™ transitionless core



ASAHI FIELDER® XT

TIP STIFFNESS	CORE	RADIOPACITY	COATING
1.2 g	TruTorq® Steel	16 cm	HYDROPHILIC

- Full polymer cover
- Tapered .009" tip diameter




TAPERED TIP

ASAHI CONFIANZA PRO® 9

TIP STIFFNESS	CORE	RADIOPACITY	COATING
9.3 g	TruTorq® Steel	20 cm	HYBRID


- Tapered .009" tip diameter



ASAHI CONFIANZA PRO® 12

TIP STIFFNESS	CORE	RADIOPACITY	COATING
12.4 g	TruTorq® Steel	20 cm	HYBRID

- Tapered .009" tip diameter





Home

Frontline Workhorse

Frontline Finesse

Extra Support


Specialty

Guide Wire Technology

Support/Tip Stiffness

Ordering Information

Important Safety Information  
referenced within



# Specialty

guide wires are designed to provide increasing tip stiffness and excellent torque response for superb crossing performance

STRAIGHT TIP

ASAHI MIRACLEBROS® 3

TIP STIFFNESS	CORE	RADIOPACITY	COATING
3.9 g	TruTorq® Steel	11 cm	HYDROPHOBIC



ASAHI MIRACLEBROS® 4.5

TIP STIFFNESS	CORE	RADIOPACITY	COATING
4.4 g	TruTorq® Steel	11 cm	HYDROPHOBIC



ASAHI MIRACLEBROS® 6

TIP STIFFNESS	CORE	RADIOPACITY	COATING
8.8 g	TruTorq® Steel	11 cm	HYDROPHOBIC



ASAHI MIRACLEBROS® 12

TIP STIFFNESS	CORE	RADIOPACITY	COATING
13.0 g	TruTorq® Steel	11 cm	HYDROPHOBIC





Home



Frontline Workhorse

Frontline Finesse

Extra Support

Specialty

Guide Wire Technology



Support/Tip Stiffness

Ordering Information


**Important Safety Information**  
 referenced within

# Celebrating 30 Years of Pioneering Guide Wire Technology for Interventional Success

## Core Materials

3 different core material options let you choose the level of strength, flexibility and durability you need

- ASAHI TruTorq® Steel for excellent durability and torque response
- DURASTEEL stronger than conventional stainless steel for outstanding durability
- ELASTINITE® Nitinol known for its flexible and durable properties giving you exceptional push and torque control



TruTorq® Steel



DURASTEEL



ELASTINITE®

## Core Tapers

Our core taper designs provide excellent guide wire tracking

- Longer core tapers yield superb wire tracking and a lower propensity to prolapse
- Shorter core tapers yield longer segments of consistent support but provide increased propensity to prolapse
- Proprietary RESPONSEASE™ transitionless core taper is designed to maximize tracking while providing progressive support

## Tip Styles

Offering 2 different designs, each providing different advantages

- Shaping ribbon for tip shape retention and softness
- Core-to-tip design provides tactile feedback and tip control, enabling exceptional 1:1 torque



Shaping ribbon



Core-to-tip design

## Spring Coil Design

For outstanding torque transmission and control

- Abbott's spring coil design contributes to the wire's tip shapeability, tip shape retention and tactile feedback
- Platinum alloy tip coils are radiopaque and stainless steel coils are radiolucent
- Asahi uses a proprietary method of fusing stainless steel to platinum spring coils that yields excellent torque response, delivery and device interaction



## Polymer Covers

Designed for deliverability in tortuous anatomy and excellent distal access

- Full polymer covers over tip coils for outstanding deliverability
- Intermediate polymer sleeves maintain exposed gentle tip coils while providing a lubricious surface for device delivery




Full Polymer Cover



Intermediate Polymer Sleeve







# Celebrating 30 Years of Pioneering Guide Wire Technology for Interventional Success

## Coatings

Precision-engineered coatings designed to reduce surface friction, improve device interaction and guide wire tracking

- Hydrophilic coatings (Hydrocoat®) attract water to create a slippery “gel-like” surface for improved lubricity and smooth device delivery
- Hydrophobic coatings (Microglide® and SmoothGlide™) repel water to create a “wax-like” surface for outstanding tactile feedback and smooth device interaction
- Hybrid coatings combine hydrophobic tip coils for tactile feedback and tip control with hydrophilic intermediate coils for smooth device delivery

## Support

- Support is a measure of a guide wire’s resistance to a bending force
- A more supportive wire can aid in bulky device delivery
- A less supportive wire is designed to be flexible and can aid in access through tortuous anatomy

## Tip Stiffness

- Sometimes referred to as “tip load,” tip stiffness is a measure of the buckling force exerted by the tip of the wire on a surface
- A high tip load is designed to help cross a challenging lesion
- A low tip load makes the tip very soft and gentle
- Abbott Vascular reports tip loads measured on the distal 10 mm of the tip





Home



Frontline Workhorse

Frontline Finesse

Extra Support

Specialty

Guide Wire Technology


Support/Tip Stiffness





Ordering Information

**Important Safety Information**  
referenced within



## The Abbott Vascular Family of Guide Wires: HI-TORQUE® and ASAHI®




### Frontline Workhorse Guide Wires




### Frontline Finesse Guide Wires



### Extra Support Guide Wires




### Specialty Guide Wires





Home



Frontline Workhorse

Frontline Finesse

Extra Support

Specialty

Guide Wire Technology

Support/Tip Stiffness

Ordering Information

»

Important Safety Information  
referenced within

Product Ordering Information for HI-TORQUE® and ASAHI® Guide Wires

		STOCK NUMBER	TIP RADIOCAPACITY (cm)	TIP SHAPE	WIRE SIZE (in)	WIRE LENGTH (cm)	
Frontline Workhorse	Nitinol / Elastinite®	HT BALANCE	1009664	3	Straight	.014	190
		MIDDLEWEIGHT	1009664J	3	J	.014	190
		UNIVERSAL II	1009665	3	Straight	.014	300
			1009665J	3	J	.014	300
		HT BALANCE	1009660	3	Straight	.014	190
		MIDDLEWEIGHT	1009660J	3	J	.014	190
		UNIVERSAL	1009661	3	Straight	.014	300
			1009661J	3	J	.014	300
		HT BALANCE	1001780-HC	3	Straight	.014	190
		MIDDLEWEIGHT	1001780J-HC	3	J	.014	190
	Stainless Steel		1001782-HC	3	Straight	.014	300
			1001782J-HC	3	J	.014	300
			1001780	3	Straight	.014	190
			1001780J	3	J	.014	190
			1001782	3	Straight	.014	300
			1001782J	3	J	.014	300
		HT ADVANCE	1044588	3	Straight	.014	190
			1044588J	3	J	.014	190
			1044589	3	Straight	.014	300
			1044589J	3	J	.014	300
ASAHI PROWATER®		12776-01	3	Straight	.014	180	
		12776-02	3	J	.014	180	
		14935-01	3	Straight	.014	300	
		14935-02	3	J	.014	300	
	ASAHI PROWATERFLEX®	82358-01	3	Straight	.014	180	
		82358-11	3	J	.014	180	
		82358-02	3	Straight	.014	300	
		82358-12	3	J	.014	300	
	ASAHI SOFT	12780-01	3	Straight	.014	180	
		12780-02	3	J	.014	180	
	14939-01	3	Straight	.014	300		
	14939-02	3	J	.014	300		
HT FLOPPY II		22339M	2	Straight	.014	190	
		22339MJ	2	J	.014	190	
		22359M	2	Straight	.014	300	
		22359MJ	2	J	.014	300	
		22359MJ-903	30	Straight	.014	300	
	HT WHISPER MS		1005357H	3	Straight	.014	190
			1005357HJ	3	J	.014	190
		1005359H	3	Straight	.014	300	
		1005359HJ	3	J	.014	300	
HT WHISPER EXTRA SUPPORT (ES)			1011834H	3	Straight	.014	190
		1011834HJ	3	J	.014	190	
		1011835H	3	Straight	.014	300	
		1011835HJ	3	J	.014	300	
	HT PILOT 50		1010480-H	3	Straight	.014	190
		1010480-HJ	3	J	.014	190	
		1010483-H	3	Straight	.014	300	
		1010483-HJ	3	J	.014	300	
ASAHI FIELDER®			82359-01	3	Straight	.014	180
		82359-02	3	Straight	.014	300	
		82359-11	3	J	.014	180	
		82359-12	3	J	.014	300	
	HT BALANCE HEAVYWEIGHT		1000462H	4.5	Straight	.014	190
		1000462HJ	4.5	J	.014	190	
		1000463H	4.5	Straight	.014	300	
		1000463HJ	4.5	J	.014	300	
HT ALL STAR		1001740	3	Straight	.014	190	
		1001740J	3	J	.014	190	
		1001741	3	Straight	.014	300	
		1001741J	3	J	.014	300	
	ASAHI GRAND SLAM®		12781-01	4	Straight	.014	180
		12781-02	4	J	.014	180	
		14940-01	4	Straight	.014	300	
		14940-02	4	J	.014	300	
HT IRON MAN			1001309	3	Straight	.014	190
		1001309J	3	J	.014	190	
Specialty	HT PILOT 200		1010482-H	3	Straight	.014	190
			1010482-HJ	3	J	.014	190
			1010485-H	3	Straight	.014	300
			1010485-HJ	3	J	.014	300
	ASAHI FIELDER® XT		AGP140002	16	Straight	.014	190
			AGP140302	16	Straight	.014	300
	ASAHI MIRACLEBROS® 3		12778-01	11	Straight	.014	180
			14937-01	11	Straight	.014	300
	ASAHI MIRACLEBROS® 4.5		12777-01	11	Straight	.014	180
			14936-01	11	Straight	.014	300
	ASAHI MIRACLEBROS® 6		12779-01	11	Straight	.014	180
			14938-01	11	Straight	.014	300
	ASAHI MIRACLEBROS® 12		82903-01	11	Straight	.014	180
			82903-02	11	Straight	.014	300
	ASAHI CONFIANZA PRO® 9		20629-01	20	Straight	.014	180
			20629-02	20	Straight	.014	300
	ASAHI CONFIANZA PRO® 12		82902-01	20	Straight	.014	180
			82902-02	20	Straight	.014	300