

MEKANT

User Manual

Ventilator System



MEKICS Co., Ltd.

5F, A-BLDG, WooLimLion, #144-3, Sangdaewon-dong,
Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea 462-725

TEL : +82-70-7119-2500

FAX : +82-31-735-2761

URL : <http://www.mek-ics.com>



Definitions

This operator's manual uses three special indicators to convey information of a specific nature. They include :



Warning

Indicates a condition that can endanger the patient or the ventilator operator.



Caution

Indicates a condition that can damage the equipment.



NOTE :

Indicates points of particular emphasis that make operation of the ventilator more efficient or convenient.

Operator's Guide

1 Introduction	9
1.1 Warnings and Cautions	9
1.1.1 Warning	9
1.1.2 Caution	11
1.2 General description	15
1.2.1 Front picture and the explanation	16
1.2.2 Back side picture and explanation	18
1.2.3 Left side picture and explanation	20
1.2.4 Right side picture and explanation	21
1.3 Specifications	22
1.3.1 Features	22
1.3.2 Electrical Specifications	22
1.3.3 Functional Specifications	23
1.3.4 Pneumatic Specifications	29
1.3.5 Physical data	29
1.3.6 Environmental data	30
1.4 Standards and safety requirements	30
1.4.1 Standards and safety requirements	30
1.5 Symbols	31
1.5.1 Description of Symbols	31
1.6 Maintenance	35
2. Installation and set up	37
2.1. Installation	37
2.2. External electrical supply	37
2.2.1 AC description	37
2.2.2 External battery	38

2.3 Internal battery	40
2.3.1 Internal battery Indicator	41
2.3.2 How to change internal battery	42
2.4 Oxygen and Air supply	43
2.5 Filters	44
2.6 Transportation	45
2.7 Fuse replacement instruction	45
2. Operating instructions	47
3.1. Front panel	47
3.1.1 Rotational Control Knob	47
3.1.2 Alarm lamp	47
3.1.3 Keypad	48
3.1.4 LED Indicator	51
3.2 Operating.....	52
3.2.1 Operating procedure and how to use.....	52
3.3 Display	54
3.3.1 Status information	55
3.3.2 Measurement information	56
3.3.3 Waveform and alarm status	57
3.3.4 Setting information	58
3.4 Menu	59
3.4.1. VENT Mode	60
3.4.1.1 VACV	62
3.4.1.2 VSIMV	65
3.4.1.3 SPONT	68
3.4.1.4 PACV	71
3.4.1.5 PSIMV	74
3.4.1.6 TBILEVEL	77
3.4.1.7 NIV	81
3.4.1.8 PRVC	83
3.4.1.9 AUTO	85

3.4.2 Alarm set	86
3.4.3 Trend Mode	88
3.4.4 SYSTEM SETUP	90
3.4.4.1 SOUND.....	91
3.4.4.2 NEBULIZER	92
3.4.4.3 Graph SETUP	93
3.4.4.4 User SETUP.....	95
3.4.4.5 SYSTEM SETUP.....	104
4. Alarms	111
4.1 General Information	111
4.2 Alarm Massage.....	112
Appendix	113
A. Operating theory	114
B. Communications	117
C. Parts accessories	120
Service Calls	122

Introduction

SECTION 1

This operator's manual explains the way of installation and usage for MEKANT. Notice for use explains general information to be careful while using the product and other notices are written with each function. Users are recommended to be aware of notice for use and all menus before using the product.

This operator's manual is allowed to use by doctors or authorized persons who are permitted to use, or designed to use by qualified and trained persons for using this products.

Service, test, correction and main continuous working will be written in service manual of MEKANT Ventilator. However, service manual is allowed to be possessed by authorized persons by MEKICS.

1.1 Warnings and Cautions

1.1.1 Warning

Please make sure to use MEKANT by authorized persons and read the manual before use. Check the setting figure of various sensors and devices before use. Users shall check out whether the environment is OK before use.

Please be aware of the manual provided by the manufacturer when it is used or kept.

Do not use it near inflammable anesthetic to avoid any possible explosion.

Do not disassemble or open the main body to avoid electrical impulse. Disassembly of the equipment shall be allowed by service staffs approved by MEKICS Co., Ltd.

Users shall be careful with circuit joint to be connected with patients.

Always shall be careful with leakage of electrical current and consult with hospital technical engineer if any doubt.

In case of abnormal working or any damage in the equipment, it shall not be used for patients.

Do not touch patient, table or the equipment if a defibrillator is used.

Disposable stuffs and disposable circuit for patients shall not be reused.

Reusable stuffs used for patients shall be sterilized or washed for reuse.

Check out the level before setting alarm. Alarm limitation shall be done for the current selected group. Alarm setting can be changed by level change.


Do not correct calibration while used for patient. Please correct it after disconnection from patient if any doubt.

Insulate power source in case that the equipment is wet and call a hospital technical engineer.

Please connect earth cable to earth connection port to protect patients and the product.

Do not disassemble the product in any case.

Do not install the various kinds of sensor to patients too strongly and not to be entangled.

Do not reuse it if you find this symbol  on the accessory, which means disposable stuff.

MEKANT shall be tested regularly if it is working in normal condition.

Accessories are aseptic products, so the package shall not be torn before use and products shall be free from any defects or cracks.

1.1.2 Caution

MEKANT as respiration assistant equipment for adult helps for normal breathing by setting breathing way, breathing volume, breathing times, oxygen density, breathing rate etc in order to solve hypoxia to be appeared from emphysema of the lungs patients or respiratory disease patients. It is a respiration assistant equipment with various alarm function and safety device for safe respiration management, but above all, should be used together with other clinical information for patients' safety and diagnosis.

Avoid the place where there is moisture.

Do not handle the equipment with wet hands.

Avoid the place under direct rays.

Use the equipment within the limits of temperature 10 ~ 40 centigrade and humidity 30 ~ 75%

Do not place the equipment near electronic heater.

Do not place the equipment under high humidity and bad ventilation.

Do not place the equipment where it is possible to take vibration or sudden impulse.

Do not place the equipment where to be damaged by explosive gas or chemicals.

Protect the equipment from dust or metal scrap.

Pull out power code with plug grasped, instead of cable.

Check out the main equipment and sensor at the same time after cleaning the equipment. Do not use if there is any damage in the equipment.

Make the cable and sensor clean after use. Clean them by towel with alcohol soaked and do not use benzene or gasoline.

Do not use the equipment in dangerous place. Always be careful with the possibility of bad influence by environment.

It is possible to be flew more current than allowed when the assistant device is used. Consult with hospital engineer or authorized local dealer.

Check out the basic auxiliary device regularly. Do not use a damaged or electrical unstable device. Consult with hospital engineer if any doubt.

Do not use the goods beyond standard certified one. This product shall be used under standard certification. Manufacturer does not take responsibility for the problem occurred while using the product beyond standard certification.

The product can be disassembled by educated and authorized staffs. If a user violates this item, after sales service is not provided.

Do not install the various kinds of sensor to patients too strongly and not to be entangled.

Do not place the product where it is possible to fall to a patient.

Do not move the equipment when power code or high pressure hose cage is connected. It will be harmed for patient or the equipment has an error by itself

If a measured figure is turned out not to be correct, consult with hospital engineer or authorized local dealer.

Connect the equipment isolated from other devices electronically for patients' electrical safety.

Accessories of MEKANT shall be used with the ones provided by MEKICS or MEKICS's designated dealer or listed in the manual. MEKICS or MEKICS's designated dealer provide standard certified products only.

MEKANT shall be tested regularly if it is working well.

Turn off main switch located backside of the product after use.

Place the product in stable position.

Use the sensor and cable provided and designated by manufacturer, and consult with hospital engineer or authorized local dealer if you want to use other products.

Replace the disposable sensor after use.

Sep up the equipment when it is needed.

Minimize patient's move when it is used.

Make high pressure hose and cable clean after use. Clean them by towel with alcohol soaked and do not use benzene or gasoline.

The Oxygen Air pressure which is impressed pneumatic block should be between 12PSI and 14PSI. Otherwise, it would affect the equipment.

Move the equipment with fastened to avoid any possible damage to other products or people by the monitor or cable.

Basically, MEKANT is designed not to have functional disorder by surrounded electronic wave. But, it is possible to be occurred any disorder due to high frequency of surrounding. It may be possible to calculate wrong figure, lose present figure or show noisy wave. Please ask MEKICS technical service center or local distributor if you need help.

Please turn off the power to recharge the battery when it is not used.


MEKANT ventilator system is a product conforming to Class II b.

1.2 General Description

MEKANT as respiration assistant equipment for adult helps for normal breathing by setting breathing way, breathing volume, breathing times, oxygen density, breathing rate etc in order to solve hypoxia to be appeared from emphysema of the lungs patients or respiratory disease patients.

It is a respiration assistant equipment with various alarm function and safety device for safe respiration management, but above all, should be used together with other clinical information for patients' safety and diagnosis

1.2.1 Front picture and the explanation




<i>No.</i>	<i>Part</i>	<i>Description</i>
Ⓐ	Circuit Arm	<i>Install Circuit Arm for patient.</i>
Ⓑ	LCD Monitor	<i>Show the condition of various ventilators and patients.</i>
Ⓒ	Key Pad	<i>Key button and encoder.</i>
Ⓓ	Inspiration Connector	<i>Joint hole to connect to patient.</i>
Ⓔ	Expiration Connector	<i>Joint hole to connect from patient.</i>
Ⓕ	Expiration Lever	<i>Lock or unlock lever for Expiration Assembly.</i>
Ⓖ	Carriage Grip	
Ⓗ	Nebulizer Port	<i>Connection to spray medicines</i>



Warning

Be careful because heated exhalation assembly due to moisture during exhalation.

1.2.2 Back side picture and the explanation



 \
Instruction

No.	Part	Description
Ⓐ	Speaker	<i>For alarm and sound.</i>
Ⓑ	VGA COMM To Pneumatic	<i>VGA cable connection hole to show patient condition on the other Monitor.</i> <i>Cable connection hole to connect PC or central connection.</i> <i>Connection hole to connect to the main body.</i>
Ⓒ	AC Inlet	<i>Connection hole for AC power input.</i>
Ⓓ	External Battery	<i>Connection hole for external battery.</i>
Ⓔ	Power Switch	<i>Power on/off switch for common use.</i>
Ⓕ	O₂ Inlet	<i>Connection hole to connect to main body using high pressure hose to approve O₂ gas.</i>
Ⓖ	Air Inlet	<i>Connection hole to connect to the main body, using high pressure hose to approve air gas. Device to filter the moisture in the Air.</i>
Ⓗ	Tank Relief	<i>Outlet to exhaust the over pressed air when the pressure of air tank is over than 40 PSI.</i>
Ⓘ	Battery Replacement Door	<i>Door to replace the internal battery.</i>
Ⓙ	To Monitor	<i>Connection for monitor.</i>
Ⓚ	COMM	<i>Connection hole for pneumatic communication.</i>

**Warning**

You should be careful with battery's polarity when it is exchanged. Products can be damaged when polarity is connected wrongly.

It should be careful of battery's terminal to be short on the metal part


1.2.3 Left side picture and the explanation



① **LOCKER**

*Locker to disassemble main body with cart.
Pull to unlock and push to lock.*

1.2.4 Right side picture and the explanation



① **O₂ Exchange Door**

Door to be used for exchanging O₂ Sensor.

② **Inspiration Check Valve**

Safety valve to exhaust by software, hardware and mechanically when the equipment is out of order.



Warning

Do not stop up safety valve to exhaust by software, hardware and mechanically when the equipment is out of order.

1.3 Specifications

1.3.1 Features

User friendly POP-up Menu with rotary and push/pull Encoder.

100%O₂, Inspiration Pause, Expiration Pause, Key-Lock ,

Manual and/or Auto Nebulizer

Dual Alarm (Speaker / Buzzer) system

External battery (for Mobile application)

1.3.2 Electrical Specifications

<i>AC input voltage</i>	<i>100 - 120 VAC, 1.0A 50/60 Hz 200 - 240 VAC, 4A 50/60 Hz</i>
<i>AC input FUSE Rating</i>	<i>250V/6.3A</i>
<i>Power consumption</i>	<i>14W max</i>
<i>External DC input voltage</i>	<i>12V @ 1A Pb battery</i>
<i>Internal battery</i>	<i>Battery ; 12V / 1A Pb battery Operating time: 3 hours minimum . Recharge time; 4 Hours max. Please change it after 2 years use. It can be used for 30 minutes normally. On the other hand, it will be more than 30 minutes if power consumption is small.</i>
<i>Communication</i>	<i>Central interface ; RS-232 serial interface Baud rate 115,200 bps, Upgrade (Main board, Pneumatic board)</i>
<i>Display</i>	<i>10.4" LCD monitor, 640*480,</i>

1.3.3 Functional Specifications

MODE

<i>V-ACV</i>	<i>Volume Assist Controlled Ventilation</i> <i>Controlling patient's mechanical artificial respiration based on air volume and adjusting inspiration and exhalation time by force.</i>
<i>V-SIMV</i>	<i>Volume based Synchronized Intermittent Mandatory Ventilation</i> <i>Controlling artificial respiration based on air volume in accordance with patient's voluntary respiration.</i>
<i>SPONT</i>	<i>Spontaneous Ventilation</i> <i>Giving pressure support centering around voluntary respiration.</i>
<i>P-ACV</i>	<i>Pressure based Assistance Controlled Ventilation</i> <i>Controlling patient's mechanical artificial respiration based on pressure.</i>
<i>P-SIMV</i>	<i>Pressure based Synchronized Intermittent Mandatory Ventilation</i> <i>Controlling artificial respiration based on pressure in accordance with patient's voluntary respiration.</i>
<i>T-BILEVEL</i>	<i>Same as SPONT and exist Upper PEEP (Refer to SPONT mode)</i>
<i>NIV</i>	<i>Non Invasive Ventilation</i> <i>During artificial respiration by mask type, not tubing type, PACV working in allowance of air leakage. (Refer to PACV mode)</i>
<i>PRVC</i>	<i>Pressure Regulated Volume Controlled Ventilation</i> <i>Breaths are delivered mandatorily with a constant flow to assure present volumes.</i>
<i>AUTO</i>	<i>Quick start according to body weight. (Working in PACV mode)</i>

SETTING PARAMETER

	<u>Setting range</u>	<u>Condition</u>	<u>Default</u>
BW	3 kg ~ 150 kg		10 ml/kg
Pressure	5 ~ 18 cmH ₂ O	±(1.1 + 4% of the actual reading) cmH ₂ O	5 cmH ₂ O
Plateau (pause) time	0 ~ 2 sec		0.1 sec
PEEP	0 ~ 45 cmH ₂ O	±(1.1 + 4% of the actual reading) cmH ₂ O	5 cmH ₂ O
High PEEP	5 ~ 45 cmH ₂ O	±(1.1 + 4% of the actual reading) cmH ₂ O	+15 cmH ₂ O
Low PEEP	0 ~ 45 cmH ₂ O	±(1.1 + 4% of the actual reading) cmH ₂ O	+5 cmH ₂ O
High PEEP time	0.2 ~ 24 sec	±10%	1 sec
Low PEEP time	Infinite		
	<u>Setting range</u>	<u>Condition</u>	<u>Default</u>
Inspiration time	0.2 ~ 9.9 sec	±10%	1 sec
Fend	25 ~ 100 %		100 %
Rise	SLOW,MID,FAST		MID
E-sense %	10 ~ 60 %		30 %
Pressure Support	0 ~ 45 cmH ₂ O	±(1.1 + 4% of the actual reading) cmH ₂ O	0 cmH ₂ O
Trigger sensitivity Pressure	0.5 ~ 20 cmH ₂ O	±(1.1 + 4% of the actual reading) cmH ₂ O	5 cmH ₂ O
Trigger sensitivity Flow	0.5 ~ 20 L/min	±15% (V _{MIN} > 11pm) ±0.4pm (V _{MIN} ≤ 11pm)	5 lpm
Volume	5 ml ~ 2500 ml	±15% (V _{TIDAL} > 100ml) ±20ml (V _{TIDAL} ≤ 100ml)	10 ml/kg
Rate	2 ~ 120 bpm	±2 bpm	20 bpm
O ₂ %	21 ~ 100 %	±2%	30 %

	<u>Setting range</u>	<u>Default</u>
Nebulizer	10 ~ 180 min Limited 200 ml Flow trigger base flow can not use	off
Inspiration pause	Resistance measurement	off
Expiration pause	Compliance measurement	off
Manual Inspiration	Setting tidal volume	off
100% O ₂	2 minutes	off
Key lock	Function key lock	off

QUICK RUN (Emergency Start Ventilation)

	<u>Setting range</u>
<i>I : E</i>	1 : 2
<i>O₂</i>	60 %
<i>Mode</i>	PACV
<i>Inspiration Pressure</i>	22 cmH ₂ O
<i>Inspiration time</i>	1 sec
<i>Expiration time</i>	2 sec
<i>Respiration rate</i>	20 bpm
<i>PEEP</i>	5 cmH ₂ O
<i>Trigger type</i>	Pressure type
<i>Trigger sensitivity</i>	5 cmH ₂ O
<i>Rise %</i>	40 %
<i>E-sense</i>	30 %

AlarmsSetting range

<i>High pressure</i>	$4\text{ cmH}_2\text{O} \sim 12\text{ cmH}_2\text{O}$
<i>Low pressure</i>	$0\text{ cmH}_2\text{O} \sim 4\text{ cmH}_2\text{O}$
<i>High tidal volume</i>	$20\text{ ml} \sim 2500\text{ ml}, \text{OFF}$
<i>Low tidal volume</i>	$0\text{ ml} \sim 2500\text{ ml}$
<i>High minute volume</i>	$1 \sim 99\text{L}$
<i>Low minute volume</i>	$0 \sim 98\text{L}$
<i>High respiration rate</i>	$2 \sim 120\text{ bpm}$
<i>Leak volume</i>	$10 \sim 500\text{ ml}$
<i>Apnea</i>	$\text{OFF}, 10 \sim 60\text{ sec}$
<i>O₂ % high</i>	$21 \sim 100\%, \text{OFF}$
<i>Power loss</i>	AC/DC power loss
<i>Low Air/O₂ pressure</i>	
<i>Obstructed tube</i>	
<i>Airway leak</i>	
<i>Vent in-operation</i>	
<i>Open Circuit</i>	

DISPLAY PARAMETER● **V_{ACV} MODE, V_{SIMV} MODE**

No.	Item	Range	Unit	Condition
Ⓐ	VTIDAL	0~4000	ml	±10% (VTIDAL>100ml), ±20ml (VTIDAL≤100ml)
Ⓑ	VITIDAL	0~2500	ml	±10% (VITIDAL>100ml), ±20ml (VITIDAL≤100ml)
Ⓒ	VMIN	0~110	LPM	±10% (VMIN>1LPM), ±0.2LPM (VMIN≤1LPM)
Ⓓ	RRESP	0~150	BPM	±1.0%
Ⓔ	PPEAK	0~15	cmH ₂ O	±(1.7+ 4% of the actual reading) cmH ₂ O
Ⓕ	PMEAN	0~15	cmH ₂ O	±(1.7+ 4% of the actual reading) cmH ₂ O
Ⓖ	PPAUSE	0~15	cmH ₂ O	±(1.7+ 4% of the actual reading) cmH ₂ O
Ⓗ	PEEP	0~5	cmH ₂ O	±(1.7+ 4% of the actual reading) cmH ₂ O
Ⓘ	O ₂	21~100	%	±2.0%
Ⓝ	I:E RATIO	1:1 ~ 1:99		±2.0%
Ⓚ	sVmin	0~110	LPM	±10% (VMIN>1LPM), ±0.2LPM (VMIN≤1LPM)
Ⓛ	sRR	0~150	BPM	±1.0%

● **SPONT MODE**

No.	Item	Range	Unit	Condition
Ⓐ	VTIDAL	0~4000	ml	±10% (VTIDAL>100ml), ±20ml (VTIDAL≤100ml)
Ⓑ	VMIN	0~110	LPM	±10% (VMIN>1LPM), ±0.2LPM (VMIN≤1LPM)
Ⓒ	RRESP	0~150	BPM	±1.0%
Ⓓ	PPEAK	0~15	cmH ₂ O	±(1.7+ 4% of the actual reading) cmH ₂ O
Ⓔ	PMEAN	0~15	cmH ₂ O	±(1.7+ 4% of the actual reading) cmH ₂ O
Ⓕ	FPEAK	0~110	LPM	±10% (VMIN>1LPM), ±0.2LPM (VMIN≤1LPM)
Ⓖ	PEEP	0~5	cmH ₂ O	±(1.7+ 4% of the actual reading) cmH ₂ O
Ⓗ	O ₂	21~100	%	±2.0%
Ⓘ	sVmin	0~110	LPM	±10% (VMIN>1LPM), ±0.2LPM (VMIN≤1LPM)
Ⓛ	sRR	0~150	BPM	±1.0%

● **P_{ACV}, P_{SIMV}, NIV, PRVC, Auto Mode**

No.	Item	Range	Unit	Condition
Ⓐ	VTIDAL	~ 4000	ml	±10% (VTIDAL > 100ml), ±20ml (VTIDAL ≤ 100ml)
Ⓑ	VMIN	~ 10	LPM	±10% (VMIN > 1LPM), ±0.4LPM (VMIN ≤ 1LPM)
Ⓒ	RRESP	~ 100	BPM	±1.0%
Ⓓ	PPEAK	~ 100	cmH ₂ O	±(1.0 + 5% of the actual reading) cmH ₂ O
Ⓔ	PMEAN	~ 100	cmH ₂ O	±(1.0 + 5% of the actual reading) cmH ₂ O
Ⓕ	FPEAK	~ 10	LPM	±10% (VMIN > 1LPM), ±0.4LPM (VMIN ≤ 1LPM)
Ⓖ	PEEP	~ 50	cmH ₂ O	±(1.0 + 5% of the actual reading) cmH ₂ O
Ⓗ	O ₂	21 ~ 100	%	±2.0%
Ⓘ	I:E RATIO	1:1 ~ 1:99		±2.0%
Ⓙ	sVmin	~ 10	LPM	±10% (VMIN > 1LPM), ±0.4LPM (VMIN ≤ 1LPM)
Ⓚ	sRR	~ 100	BPM	±1.0%

● **T_{BILEVEL} MODE**

No.	Item	Range	Unit	Condition
Ⓐ	VTIDAL	~ 4000	ml	±10% (VTIDAL > 100ml), ±20ml (VTIDAL ≤ 100ml)
Ⓑ	VMIN	~ 10	LPM	±10% (VMIN > 1LPM), ±0.4LPM (VMIN ≤ 1LPM)
Ⓒ	RRESP	~ 100	BPM	±1.0%
Ⓓ	PPEAK	~ 100	cmH ₂ O	±(1.0 + 5% of the actual reading) cmH ₂ O
Ⓔ	PMEAN	~ 100	cmH ₂ O	±(1.0 + 5% of the actual reading) cmH ₂ O
Ⓕ	FPEAK	~ 10	LPM	±10% (VMIN > 1LPM), ±0.4LPM (VMIN ≤ 1LPM)
Ⓖ	PEEPH	~ 50	cmH ₂ O	±(1.0 + 5% of the actual reading) cmH ₂ O
Ⓗ	PEEPL	~ 50	cmH ₂ O	±(1.0 + 5% of the actual reading) cmH ₂ O
Ⓘ	O ₂	21 ~ 100	%	±2.0%
Ⓙ	sVmin	~ 10	LPM	±10% (VMIN > 1LPM), ±0.4LPM (VMIN ≤ 1LPM)
Ⓚ	sRR	~ 100	BPM	±1.0%

1.3.4 Pneumatic Specifications

<i>O₂ and air supply Pressure range</i>	1.26 ~ 6.5 Kg/cm ² or 1.2 ~ 6.2 bar or 1.2 ~ 90 psi
<i>Flow</i>	12 L/min (STPD, dry required)
<i>100% O₂ Suction button – auto calibration</i>	
<i>Oxygen sensor life</i>	10,000 hours or 1 Years of use, normal
<i>Safety pressure Pneumatic</i>	1.2 psi +/- 10%
<i>Internal pressure Limit</i>	
<i>Inspiration Pressure Limit</i>	12 cmH ₂ O +/- 10%; mechanical release valve
<i>Double safety mechanism for over pressure</i>	

1.3.5 Physical Data

<i>Dimensions of device</i>	38 W x 43 D x 69 H cm (Width x depth x height of device)
<i>Weight of device</i>	1.2 Kg (with Battery)
<i>Dimensions of Cart</i>	49 W x 56 D x 77 H cm (Width x depth x height of device)
<i>Weight of Cart</i>	28 kg
<i>Connectors</i>	
<i>Inspiratory limb connector</i>	ISO 15mm conical male
<i>Expiratory limb connector</i>	ISO 15mm conical male
<i>Air and Oxygen inlets</i>	DISS male / female

1.3.6 Environmental Data

<i>Operating environment</i>	<i>Temperature : 10 to 40 °C (50 ~ 104 °F)</i>
<i>Atmosphere</i>	<i>700 to 1060 mbar (10.2 ~ 15.4 psi)</i>
<i>Altitude</i>	<i>- 443 ~ 3280 m (- 1450 ~ 10,763 ft)</i>
<i>Storage environment</i>	
<i>Temperature</i>	<i>- 10 to 40 °C</i>
<i>Relative humidity</i>	<i>10 to 85 %</i>
<i>Atmosphere</i>	<i>500 to 1060 mbar (7.3 ~ 15.4 psi)</i>
<i>Altitude</i>	<i>up to 4560 m (15,000 ft)</i>
<i>Oxygen/Air inlet supply pressure</i>	<i>2.56 ~ 7.3 Kg/cm² / 2.4 ~ 6.9 bar / 35 ~ 100 psi</i>

1.4 Standard and Safety requirements

MEKANT is in conformity with the following tables. Obtainment of license and required testing has been completed.







1.4.1 Standards and Safety requirements












<i>EN60601-1:1990</i>	<i>General requirements for Safety</i>
<i>EN60601-1-1:2001</i>	<i>Electromechanical Safety</i>
<i>EN60601-1-2:2001</i>	<i>Electromagnetic Compatibility Requirement and tests</i>
<i>IEC 60601-2-12:2001</i>	<i>Particular requirements for the safety of lung ventilators for medical use</i>

1.5 Symbols and Labels

Descriptions for marking on MEKANT ventilator and its operational manuals.










1.5.1 Description of Symbol

No	Name	Symbol	Explanation
1	POWER ON/OFF		Provision of power by switching on and off the system. IEC 601-1
2	Equipotentiality		It is indicated terminal which do not need ground electric potential to make equipotentiality by connecting with each part of system or equipment. IEC 601-1, IEC 417-5.21
3	Attention, Consult accompanying documents		Indicates the necessity to refer to the operational manual prior to equipment operation. IEC 601-1
4	stand-by		Indicates the stand-by mode of the system. IEC 417-5.09
5	Start (of action)		Indicate the activation of the equipment. IEC 417-51.4
6	Direct current		Indicates appropriateness for the equipment against direct current and designate the connecting part. IEC 601-1, IEC 417-5.31

No	Name	Symbol	Explanation		
17	Battery check		Battery capacity indicator		
			 FULL Charging	 2/3 Charging	 1/2 Charging
			 Low Charging	 Discharging	
18	AC Power plug		Indicates the AC plug without ground earth. IEC 60320-1-1		
19	Battery		It should be indicated rating sticker that equipment is suitable for only direct current (D.C). Also, indicate it is related connecting part. IEC 601-1, IEC 60730-1		
20	EXT. BATT	EXT. BATT	Indicates external battery for vehicle use.		
21	KEY LOCK		Locks KEY BUTTON		
22	Manual control		Indicates manual operation of the apparatus. ISO 9001-1996		
23	CLEAR		Clears the menu tree from the display.		

1

Instruction

No	Name	Symbol	Explanation
14	Dangerous voltage		Warns the high voltage. IEC601-1, IEC617-5.036, ISO 3864
15	Respiration Mode		Indicate a respiration mode in VENT Mode. (VACV. VSIMV. SPONT. PACV. PSIMV. TBILEVEL. NIPPV. AUTO)
16	100% O ₂ MODE		Indicate 100% O ₂ MODE. This mode can be ON / OFF by O ₂ 100% KEY
17	FREEZE MODE		Indicate FREEZE MODE. This mode can be ON / OFF by FREEZE KEY
19	NEBULIZER MODE		Indicate NEBULIZER MODE. This mode can be ON / OFF by NEBULIZER KEY, and the operating time is set from 10 to 180 minutes.
20	CALIBRATION MODE		Indicate CALIBRATION MODE. This mode is displayed when activating CALIBRATION in SET UP menu.
21	MANUAL INSPIRATION KEY		Supply the respiration to the patient when pressing the MANUAL INSPIRATION KEY.
22	INSPIRATION TIME		Increase INSPIRATION TIME to 1 second.
23	EXPIRATION TIME		Increase EXPIRATION TIME to 1 second.



No	Name	Symbol	Explanation
r4	ALARM SOUND OFF Or MUTE		Indicate ALARM SOUND OFF or MUTE. This function can be ON / OFF by pressing √ MIN KEY.
r5	TREND MODE		Indicate TREND MODE.

Table 1-1. Description of Symbols

1.6 Maintenance

The following table shows the periodic check point and maintenance method for MEKANT ventilator system and the circuit system.

Period	Subject	Check point
Daily	Bacteria Filters - Inspiratory Filter - Expiratory Filter	Check if all filters are in good position.
	Patient Circuit System	Check the connection, leakage, Blockage and moisture.
	Collector Vial, Gas Supply Water Traps	Check if there is moisture and remove it if necessary.
Minimum Every half month or when changed the patient	Patient Circuit System, all Bacteria Filters and Humidification Chamber	Exchange everything, and cleaning and sterilizing used one.
Every Half year	All system	Conduct EST
Every year or After 100 times autoclaving	All Bacteria Filters	Exchange to new filter and scrap used one
Between 1 year and 2 year	Oxygen Sensor	Exchange to new one
When necessary	Oxygen Sensor	Conduct Oxygen sensor calibration by pressing 100% O ₂ /CAL.
Every 10,000 hours	Various components	Exchange appropriate components using 10,000 hours PM kit.

Table 1-2. Check and Maintenance Period

**Caution**

User should check the installation status and problem of the filters and patient circuit system with naked eyes every day. In case of any problem, replace or exchange the parts in problem right away.

When patient is changed or one patient uses it for a long period, exchange the patient circuit system including bacteria filter in every 15 days, and clean and sterilize the used one.

All the components of patient circuit system excluding bacteria filter can be sterilized by any way among ETO(Ethylene Oxide) gas sterilization, Steam Autoclave and Chemical sterilization. But, only Steam Autoclave is available for bacteria filter.

Sterilization shall be done for about 20 minutes at 132 degree centigrade. Exchange to new bacteria filter after one year use or 100 times of Autoclave use.

Oxygen sensor is valid for 2 years, but it can be shortened in case of using at the condition of high FiO₂ or high temperature.

Installation and Set up

SECTION 2

2.1 Installation

When installing and moving the equipment, it is required to stay shock free.

2.2 External electrical supply

MEKANT is capable to receive both AC and DC power supply.

2.2.1 AC description

External AC connector is located on the rear side of the apparatus. Refer to Figure 1-1.




Figure 1-1. External AC and DC Power Supply Sockets



Warning

Prior to applying AC or DC to the ventilator, the line frequency has to be adjusted to the local standard otherwise it might cause damages to the apparatus.

Power Cord has to be plugged in to the concentric plug on the ground.

Replace covered wire in case of wear and tear.

Please apply the ground cable in order to protect the patient and the apparatus.

Maintain the same level of input & output voltage.




Caution

Cable Clips are used in order to affix the AC code and also to prevent disconnection.

When connected to DC or AC, LED display will be turned on to indicate the connection.

2.2.2 External battery



Pin	Name	Direction	Explanation
①	Not connected		Not used
②	GND	—	Pins for the ground.
③	+VCC	←	Pins to connect external battery
④	Not connected		Not used

Figure 2-2. Pin Setting of MV2000 Ventilator External Battery



Warning

Watch for polarity when connecting the battery.

Note that when external battery is used, internal battery will not be charged.

It should be used for rating voltage, 12VDC 7A and Electric current should be more than 1A.

**Caution**

Cover the Port Pin by putting the cap when the external battery is not in use.

Do not make pin to be shorted because internal Battery is damaged when EXT. BATT is opened, or each pin is shorted.

Please put a cap next to terminal when EXT. BATT is not used.

**NOTE :**

Fuse is built in the machine.

ATO[®] Fuse (Fast-Acting Type), Catalog Number 257-015, Ampere Rating (A) 15A, Body color code : blue

Vehicle outlet power supply is used.

2**Installation and Set up****2.3 Internal battery**

Internal battery for MV2000 is used for the provision of power when moving the patients. Recharge of the battery is done while AC power connected as the apparatus simultaneously utilizes the power from the same source of AC. In case of power failure, the operation will be done by internal battery.

**Warning**

When replacing the battery, watch for polarity. Red for +, Black for -. Inappropriate connection might cause damages to the apparatus.

Short-circuit might cause explosion of the battery or harm human body due to high voltage.

It should be used for rating voltage, 110V ~ 240V, 7A and Electric current should be more than 1A.

We recommend that you should contact the distributor or CS team if you find any problems in the internal battery.

**Caution**

When difficulties of power supply arise or when transferring the patients, the battery needs to be charged. Replace the battery if it is in deficiency.

In case of the lowest battery voltage, there may be an error in the measured value.

In case of low battery, use the machine after connection to AC common electric power.

**NOTE :**

Fuse is built in the machine.

ATO[®] Fuse (Fast-Acting Type), Catalog Number 257 015, Ampere Rating (A) 15A, Body color code : blue

Internal Battery may be affected by the number of times for discharge, surrounding temperature, charging voltage. Generally it can be used more than 2 years.

2.3.1 Internal battery indicator

It displays battery capability. This is displayed on the upper right hand corner of the LCD display. The battery capacity will be displayed in 5 main statuses.






Battery Indicator	Level Status	Explanation
	Full Charging	
	2/3 Charging	
	1/2 Charging	
	Low Charging	
	Discharging	In case of discharging status, the system will be automatically turn off to protect the battery 2 minutes later.

Table 2-1. Internal battery Indicator

2.3.2 How to change internal battery

Loose 4 screws on the lower rear cover that is for Battery.



You can see the internal battery. After changing battery, assemble the battery cover.



Warning

When replacing the battery, watch for polarity.

Red for +. Black for -.

Inappropriate connection might cause damages to the apparatus.

Short-circuit might cause explosion of the battery or harm human body due to high voltage.



NOTE :

All of bolt is locked when turned clockwise. All of bolt is loosed when turned counter clockwise.

2.4 Oxygen and Air supply

MEKANT is the system to provide oxygen and air to patients. Oxygen and Air are provided through the adaptor end at the rear by pressure (2.56 ~ 7.2Kgf/cm² / 2.4 ~ 6.7bar / 35 ~ 90 PSI) and get mixed to have appropriate concentration .



Warning

Check the oxygen and air connector to be free from oil to prevent from explosion.

When you supply or stop oxygen and air, be careful when it's assembled and disassembled because high pressure hose can harm to human body due to inside pressure.



Caution

Check and see if the Oxygen and Air pressure is set at exact pressure range prior to ventilator application.

When high pressure hose is used or moved, pay attention to wiring not to be tripped over persons.

Connection way of high pressure hose

O₂ consists of green adaptor and green high pressure hose.

Air input consists of a yellow color.



2.5 Filters

Use √ filters between MEKANT Product and patient for the safety of patients. Before supplying for patient INS. Filter is used, and exhalation filter is used at the place for exhaled air from patient's body in order not to supply polluted air to patient.



Warning

Contaminated gas might be in exposure to patients if the air was provided without applying filter.

Do not reuse the accessories when ⊗ symbol is found since it represents its category as disposable.

Replace the disposable sensor with the new one after single use.

Reusable materials which are used for patients shall be sterilized and washed for reuse.



Caution

Filters must be MEKICS genuine or designated products. If you wish to use other makes, please consult biomedical engineers at the hospitals, or local authorized dealer.

Please use quality certified product, otherwise manufacturer will not be responsible for any problems caused by this.

All accessories are aseptic, therefore the wrapping has to be intact before opening and use.

2.6 Transportation



Warning

When moving the apparatus, please make sure that the LCD Monitor is fixed properly since this might cause damage to the Device.
Please move the product when components and power cords are in order.



Caution


Before moving the system, unlock the locker on the caster.

Move the product by handle.



2.7 Fuse replacement instruction

Inlet Fuse



Pull the upper part of inlet toward you.

Fuses are located where two upper arrows are and those are and those lower ↴ arrows indicate the direction when inserted.

Remove the existing fuse and insert the new fuse as shown in the picture after checking the rating.

Close the inlet by doing the steps backward.


Operating instructions

SECTION 3


3.1 Front panel

3.1.1 Rotational control knob

Rotating and Pressing the Encoder allows a user to navigate through and make changes to the display element.



3.1.2 Alarm lamp



STATUS	COLOR	DESCRIPTION
NORMAL	GREEN	No alarm and Ventilator is working within Setting limits.
WARNING	RED	Ventilator does not work or alarm is occurred.

Figure 3-1. Alarm lamp



Warning

Red lamp is on in case of emergency such as Ventilator does not work. Though emergency is clear, it will be still on before pressing RESET.

3

Operating Instructions

3.1.3 Keypad












Figure ƒ-2. Keypad

No	Picture	Description
a		START SWITCH – If you press START switch for 0.5 sec, power would be on. In case of STOP OFF you want , switch should be pressed for longer than 10 seconds in Vent operating or switch should be pressed for longer than ƒ seconds in Vent Ready status.
b		MANUAL INSPIRATION : This function which is for a respiration according to Ventilator's default is operated once by pressing the key.
c		100% O₂ : If you press this button, 100% O ₂ is provided for ƒ minutes.

3

Operating Instructions

No	Picture	Description
d		NEBULIZER : It is acted once as setting value. This allows to nebulize during only inspiration and is stopped by pressing once again.
e		INS. PAUSE : Operating of Ventilator is stopped for a regular second after Inspiration Time.
f		EXP. PAUSE : Operating of Ventilator is stopped for a regular second after Expiration Time
g		VENT SETTING : This changes Respiration mode & each parameter. But if RUN key or button is not pressed, modified value would be canceled.
h		ALARM SETTING : It is set alarm limit of each parameter.
i		FREEZE : This is for stop of graph on the display. Stop graph is measured after its FREEZE
j		ALARM 2MIN : This makes ALARM remove for 2 minutes.
k		ALARM RESET : This actives when alarm is reset. And previous alarm lamp is removed and the status of NORMAL will be.






No	Picture	Description
<p>①</p> <p>②</p> <p>③</p>		<p>ACCEPT RUN : It is used when Ventilator is working as its final value. If ACCEPT key is not pressed, setting value would not be accepted on the change of MODE and VENT SETTING. This is for the safe.</p> <p>KEY LOCK : It is limit key of all operations. It can be set and removed in case of only pressing ACCEPT and CLEAR at the same time. Green LED is lighted when it sets.</p> <p>CLEAR : It is used for clearance of menu when menu is controlled by Knob</p>
<p>④</p>		<p>Encoder : Rotary switch for selecting function on the pop-up menu.</p>

Table 3-1. Keypad

3.1.4 LED Indicator

No	Picture	Description
㉑		<p>AC POWER ON LED : LED is on when machine is connected to AC, otherwise, LED is OFF.</p>
㉒		<p>BATTERY ON LED : Battery is charged as the machine is working when it is connected AC. If it is not, battery ON LED is lighted because the machine is operating by internal battery. If both AC POWER ON LED and BATTERY ON LED are OFF, battery should be charged.</p>
㉓		<p>EXT. BATT : LED is on when the machine is connected to external battery, otherwise, LED is OFF.</p>

3

Operating Instructions

3.2 Operating



! NOTE :

Buzzer will be sounded when main power is on normally. Otherwise, check something wrong on the main power. (power, battery, fuse)

ALARM LAMP (Warning or INOP) is on and buzzer is sounded if ventilator has defect on working .

Connect the ventilator to the power is earthed.

3.2.1 Operating procedure and how to use

No	Picture	Description
1		<p>Connect high pressure hose and accessories.</p>
2		<p>Connect earthed cable to equipotential port and earthed AC power.</p> <p>Connect power cord to an electric outlet and AC inlet.</p> <p>Switch on POWER on the back side.</p>

3

Operating Instructions



No	Picture	Description
3		<p>Switch on START on the Front Panel. System Fail sound will be disappeared in case power is on normally.</p> <p>Wait for the calibration of each sensor. It is presented normal INITIAL, and then normal ventilator READY is appeared.</p>
4		<p>Set the present time and date. (SYSTEM SET - USER SETUP - DATE/TIME)</p> <p>Enter patient's bed number etc. (SYSTEM SET – USER SETUP – BED NUMBER)</p> <p>Select patient's weight. (VENT MODE – B/W)</p> <p>Set up alarm limit. (ALARM SET)</p> <p>Choose mode suitable for the patient. (VENT MODE)</p> <p>Input other value as prescription. Adjust value by pressing and rotating.</p> <p>Select RUN.</p> <p>Check the wave form and figure on the monitor normally.</p>

Figure 3-3. Operating procedure and how to use

3

Operating Instructions

3.3 Display




Figure 3-4. LCD Display

Explanation		Reference
Status information	ICONS of working status, BED NO and Date/Time.	3.3.1
Measurement information	Measurement of ventilator or patient factor.	3.3.2
Waveform and alarm status	Graph and alarm status of ventilator or patient factor.	3.3.3
Setting information	Setting value according to respiratory mode.	3.3.4

Table 3-2. Description of display

3.3.1 Status information

ICONS of working status, BED NO, Date/Time.









Figure 5-5. Explanation of status information

No	Icon	Explanation
(a)	BED-07	BED Number. You can change at SETUP MENU.
(b)		It displays respiratory mode of VENT Mode. (VACV, VSIMV, SPONT, PACV, PSIMV, TBILEVEL, NIV, PRVC, AUTO)
(c)		It display O ₂ 100% Mode. You can ON/OFF to O ₂ 100% KEY.
(d)		It displays FREEZE MODE (Display Paused). You can ON/OFF to FREEZE KEY.
(e)		It displays NEBULIZER MODE. You can ON/OFF to NEBULIZER KEY. Setting time is 30 minute.
(f)		It displays KEY LOCK MODE. You cannot use another key. If you use another Key, push the ACCEPT KEY. If not, push the EXIT KEY.
(g)		It displays CALIBRATION MODE. When it work Calibration to SETUP MENU it will be displayed.
(h)		When you push the MANUAL INSPIRATION KEY, Ventilator supplies breathing to patient at the same time.

No	Icon	Explanation
①		It makes longer to "INSPIRATION TIME" a regular second.
①		It makes longer to "EXPIRATION TIME" a regular second.
Ⓚ		It displays ALARM SOUND OFF or MUTE. You can ON/OFF to ↵ MIN KEY.
①		It displays AC Power. (It maybe changed Battery Status)
Ⓜ		It displays month-day-hour-minute. You can choose setting value to DATE/TIME of SYSTEM SET menu.

3.3.2 Measurement information

It displays measurement information of ventilator or patient factor. It shows different tables of measurement information to each Mode.

VTIDAL ml	Vi TIDAL ml	V MIN LPM		O ₂ %	I-E RATIO	
607	600	8.6		21	1/2.8	
RRESP BPM	PPEAK cmH ₂ O	PMEAN cmH ₂ O	PPAUSE cmH ₂ O	PEEP cmH ₂ O	sVmin	sRR
14	21	7	17	3.0	0.0	0

Figure ㉚-6. Measurement information explanation (VACV)

3.3.3 Waveform and alarm status

It displays WAVE, ALARM STATUS and MESSAGE to ventilator or patient factor.
Refer to Section ୧ for alarm information.




Figure ୧-7. Waveform and alarm status

No	Display	Explanation
(a)	WAVE	It displays waveforms of Pressure, Flow, Volume, V/P and F/V.
(b)	ALARM STATUS	It displays alarm message. When alarm is working, it displays orange color.
(c)	MESSAGE DISPLAY	It displays Software Version and Message.

Table ୧-୩. Waveform and alarm status explanation

3.3.4 Setting information

It displays respiratory mode of VENT Mode. It shows different table of setting table to each Mode.




Figure ㄚ-8. Setting information Display explanation

No	Display	Explanation
㉑		It displays Setting Mode.
㉒		It displays Weight of setting information. (Unit : Kg)
㉓		It displays ratio of I:E


Table ㄚ-4. Setting information Display explanation

3

Operating Instructions

3.4 Menu


If you press the encoder, the Menu screen is shown up at the center left side of main menu as shown below. If you choose menu on the menu tree and press the encoder, menu is closed. Below Figure shows menu tree.



DESCRIPTION		Reference
V _{ENT MODE}	You can choose respiratory mode according to respiratory status of patient respiratory.	3.4.1
A _{LARM SET}	Choosing alarm mode	3.4.2
T _{TREND MODE}	This is function which presents each value of parameter as table and graphic. It is available for 12 hours.	3.4.3
S _{YSTEM SET}	Choosing System set	3.4.4
EXIT	Remove menu on the front screen	3.4.5

3.4.1 VENT Mode

MV2000 has the ventilation mode which are SIMV(VSIMV,PSIMV), ACV(VACV, PACV, other company calls CMV to Mandatory Mode), Spontaneous Mode, tBi-Level, NIV, PRVC and AUTO Mode.



No.	Item	Mode	Explanation	Refer.
Ⓐ	V _{ACV}	V _{ACV}	Volume based Assistance Controlled Ventilation. It controls Patient Ventilation by force from air volume. And it controls expiration time and inspiration time by force.	۳.۴.۱.۲
Ⓑ	V _{SIMV}	V _{SIMV}	Volume based Synchronized Intermittent Mandatory Ventilation. It controls Patient Ventilation to air volume from Synchronized Intermittent Mandatory Ventilation	۳.۴.۱.۳
Ⓒ	SPONT	SPONT	Spontaneous Ventilation. It supports spontaneous breathing.	۳.۴.۱.۴
Ⓓ	P _{ACV}	P _{ACV}	Pressure based Assistance Controlled Ventilation. It controls Patient Ventilation by force from air pressure.	۳.۴.۱.۵
Ⓔ	P _{SIMV}	P _{SIMV}	Pressure based Synchronized Intermittent Mandatory Ventilation. It controls Patient Ventilation to air pressure from Synchronized Intermittent Mandatory Ventilation.	۳.۴.۱.۶
Ⓕ	T _{BILEVEL}	T _{BILEVEL}	It is same to SPONT. And it has Upper PEEP. (Refer to SPONT)	۳.۴.۱.۷
Ⓖ	NIV	NIV	It is mask type, not tubing type. It work PACV allow to leakage. (Refer to PACV)	۳.۴.۱.۸
Ⓗ	PRVC	PRVC	Breath are delivered mandatorily with a constant flow to assure present volumes.	۳.۴.۱.۹
Ⓘ	AUTO	PACV	Quick Start according to the weight of patient. (PACV Controlled)	۳.۴.۱.۱۰
Ⓢ	EXIT	-	Remove menu on the front screen.	-

Table ۳-5. V_{ENT MODE} Menu tree explanation

3.4.1.1 V_{ACV}

V_{ACV} (Volume based Assistance Controlled Ventilation) : It controls Patient Ventilation by force from tidal volume. And it controls expiration time and inspiration time by force. V_{TIDAL} means inspiration volume to supply patient.



3

Operating Instructions


A/C (Volume Control) - VACV

If MEKANT perceives spontaneous breathing it works the Assisted Mandatory Breath. If not it works Controlled Mandatory Breath.

When MEKANT works in Volume Control Ventilation it decides inspiratory time for setting of tidal volume and peak flow. Also inspiratory time is influenced by flow pattern and plateau.

In the case of I:E ratio decides it for factor of inspiratory time and respiratory rate.

When you change the factor of I:E ratio and inspiratory time, you can know cycle time, inspiratory time, expiratory time and I:E ratio through breath timing bar.



No.	Item	Setting Range	Unit	Step	Explanation
㉑	V _{TIDAL}	15 ~ 2500	ml	10	It decides volume of 1 breath. NOTE : It set the volume of patient weight. It is 10ml per 1kg. (You can set the volume per 1kg in [SYSTEM SET - USER SETUP - B/W FACTOR] Menu.)
㉒	R _{RESP}	2 ~ 120	BPM	1	It set the respiratory rate per 1 minute.
㉓	T _{PAUSE}	0.0 ~ 2.0	Sec	0.1	Inspiration Pause Time. It set the inspiration pause time.
㉔	PEEP	0 ~ 45	cmH ₂ O	1	It set the PEEP pressure.(Base pressure of patient)
㉕	O ₂	21 ~ 100	%	1	It set the percent of O ₂ to patient.
㉖	F _{PEAK}	5 ~ 110	LPM	1	It set the flow-peak (maximum velocity) to patient.
㉗	F _{END}	25 ~ 100	%	5	It set percent of Flow (supply velocity) / F _{peak} (supply maximum velocity) to End-tidal.


No.	Item	Setting Range	Unit	Step	Explanation
㊸	E _{SENS}	1. ~ 6.	%	1.	It set the sensitivity of one breath completion of patient during to setting of percent of patient expiration volume.
㊸	TRIGGER Type	Pressure /Flow	cmH ₂ O /LPM		It sets the type of trigger. It is Pressure or Flow type.
㊸	TRIGGER	0.5 ~ 2.	cmH ₂ O /LPM	0.5	Setting to detect spontaneous breathing signal, There are pressure and Flow type. User can select one.
㊸	SIGH MODE	OFF,30,60, 90,120	BPM	30	It sets breathing rate timing in sigh mode. Don't support over 1500ml volume.
㊸	RUN	-	-	-	It works VACV Mode.
㊸	RETURN	-	-	-	Go back to previous menu

3

Operating Instructions

3.4.1.3 V_{SIMV}

V_{SIMV} (Volume based Synchronized Intermittent Mandatory Ventilation) : It controls Patient Ventilation to tidal volume from Synchronized Intermittent Mandatory Ventilation. V_{TIDAL} means inspiration volume to supply patient.



3

Operating Instructions


No.	Item	Setting Range	Unit	STEP	Explanation
Ⓐ	V _{TIDAL}	15~2500	ml	1	It decides volume of 1 breath. NOTE : It sets the volume of patient weight. It is 10ml per 1kg. But you can set the volume per 1kg.
Ⓑ	R _{RESP}	2~120	BPM	1	It sets the respiratory rate per 1 minute.
Ⓒ	T _{PAUSE}	0.1~2.0	Sec	0.1	Inspiration Pause Time It set the inspiration time.
Ⓓ	PEEP	0~45	cmH ₂ O	1	It set the PEEP pressure. (Base pressure of patient)
Ⓔ	O ₂	21~100	%	1	It sets the percent of O ₂ to patient.
Ⓕ	F _{PEAK}	5~110	LPM	1	It sets the flow-peak (maximum velocity) to patient.
Ⓖ	F _{END}	25~100	%	5	It sets percent of Flow (supply velocity) / Fpeak (supply maximum velocity) to End-tidal.
Ⓗ	E _{SENS}	10~60	%	10	It sets the sensitivity of one breath completion of patient during to setting of percent of patient expiration volume.
Ⓘ	TRIGGER Type	Pressure /Flow	cmH ₂ O /LPM		It sets the type of trigger. It is Pressure or Flow type.
Ⓙ	TRIGGER	0.5~2.0	cmH ₂ O /LPM	0.5	Setting to detect spontaneous breathing signal, There are pressure and Flow type. User can select one.
Ⓚ	PPS	0~45	cmH ₂ O	1	It sets the pressure to patient during spontaneous breathing.
Ⓛ	RUN	-	-	-	It works V _{SIMV} Mode.
Ⓜ	RETURN	-	-	-	Go back to previous menu

Table 4-6. V_{SIMV} Menu Tree

SIMV (Volume Control or Pressure Control) only Mandatory Breath

SIMV is a mixed ventilator mode that allows both mandatory and spontaneous breaths. SIMV mode differs to other company. SIMV of MEKANT guarantees one of mandatory breath per SIMV cycle. This mandatory breath works by patient initiated mandatory (PIM or Assisted Mandatory) breath or ventilator initiated mandatory (VIM or Controlled Mandatory).

As the below figure shows, each SIMV breath cycle (T_b) has two parts. One is the mandatory interval (T_m) and the other is Spontaneous Interval (T_s). PIM is delivered, the T_m interval ends and the ventilator switches to the spontaneous interval (T_s). If a PIM is not delivered, the ventilator delivers a VIM at the end of the mandatory interval, and then switches to the spontaneous interval.



In the SIMV mode of MEKANT, SIMV Breath Cycle is decided by respiratory rate of setting. If respiratory rate is \uparrow times, SIMV Breath Cycle is \downarrow seconds. Mandatory Interval is defined as whichever is less $\cdot \cdot \cdot \times$ SIMV Breath Cycle, or $\cdot \cdot \cdot$ seconds. Mandatory interval can be shorter according to PIM breath.

3.4.1.4 SPONT

SPONT (Spontaneous Ventilation) : It supports spontaneous breathing.

Caution

When MEKANT is working SPONT mode to apnea status during setting time, it change mode(VACV) of Apnea Backup setting automatically. It does change back SPONT mode, when 3 times spontaneous breaths within 10seconds in apnea backup mode.



No	Item	Setting Range	Unit	STEP	Explanation
㉑	PEEP	0 ~ 45	cmH ₂ O	1	It sets the PEEP pressure. (Base pressure of patient)
㉒	O ₂	21 ~ 100	%	1	It sets the percent of O ₂ to patient.
㉓	TRIGGER Type	Pressure /Flow	cmH ₂ O /LPM		It sets the type of trigger. It is Pressure or Flow type.
㉔	TRIGGER	0.5 ~ 2.0	cmH ₂ O /LPM	0.5	Setting to detect spontaneous breathing signal, There are pressure and Flow type. User can select one.
㉕	PPS	0 ~ 45	cmH ₂ O	1	It sets the pressure to patient during spontaneous breathing.
㉖	APNEA	OFF ~ 60	sec	10	If patient has not breath of setting time it works Apnea mode. It does not detect breath of patient during "OFF"
㉗	RUN	-	-	-	It works SPONT Mode.
㉘	RETURN	-	-	-	Go back to previous menu

Table ㉗-7. SPONT Menu tree

SPONT mode only

If spontaneous breath reach to sensitivity, it makes pressure between setting of sensitivity and $1\text{ cmH}_2\text{O}$ through flow supply according to patient flow demand. If Pressure is bigger than $1\text{ cmH}_2\text{O}$, it starts expiration to opening Exhalation valve. If patient has apnea breath, MEKANT occurs Apnea alarm and start Apnea ventilation(VACV Mode).

It works setting of Respiratory Rate, Tidal Volume, Inspiratory Pressure, Peak Flow, Inspiratory Time, FiO_2 and so on.

It works same as pressure support ventilation. Also if patient has apnea breath, MEKANT occur Apnea alarm and start Apnea ventilation. But it change previous spontaneous mode if patient has 3 times spontaneous breaths within 10seconds.

3.4.1.5 P_{ACV}

P_{ACV} (Pressure based Assistance Controlled Ventilation) : It controls Patient Ventilation by force from air pressure.




No	Item	Setting Range	Unit	STEP	Explanation
Ⓐ	P _{INSP}	5 ~ 10	cmH ₂ O	1	It sets the pressure to patient during inspiration.
Ⓑ	T _{INSP}	0.2 ~ 9.9	sec	0.1	It sets the inspiration time.
Ⓒ	R _{RESP}	2 ~ 120	BPM	1	It sets the respiratory rate per minute.
Ⓓ	PEEP	0 ~ 45	cmH ₂ O	1	It sets the PEEP pressure. (Base pressure of patient)
Ⓔ	O ₂	21 ~ 100	%	1	It sets the percent of O ₂ to patient.
Ⓕ	TRISE	FAST,MID, SLOW	-	-	It sets the flow speed of gas to patient.
Ⓖ	E _{SENS}	10 ~ 60	%	10	It sets the sensitivity of one breath completion of patient during to setting of percent of patient expiration volume.
Ⓗ	TRIGGER Type	Pressure /Flow	cmH ₂ O /LPM		It sets the type of trigger. It is Pressure or Flow type.
Ⓘ	TRIGGER	0.5 ~ 20	cmH ₂ O /LPM	0.5	Setting to detect spontaneous breathing signal, There are pressure and Flow type. User can select one.
Ⓙ	RUN	-	-	-	It works PACV Mode
Ⓚ	RETURN	-	-	-	Go back to previous menu

Table 8-8. P_{ACV} Menu tree**A/C (Pressure Control) PACV**

When MEKANT provide pressure control ventilation of mandatory breath, it draws inspiratory pressure instead of tidal volume and it sets I:E ratio or inspiratory time instead of peak flow same as other company.


But as the below figure shows, it is differ to slope of increase – baseline (PEEP) to Inspiratory pressure – according to setting of flow acceleration percent (FAP).

According to an increment of FAP setting value, slope of pressure is rapid. If not slope of pressure is gentle. If you show slope of flow, you know that FAP is the change of flow volume to patient. FAP sets PCV and PSV each other.



MEKANT keeps on pressure regularly during inspiration same as other company. So it has Ramp shape of flow pattern and cannot use plateau function.

MEKANT is more activity and elaborateness for control of airway pressure because it has exhalation valve and active exhalation valve. So it controls pressure to exhalation valve in patient's talk, cough. And it allows spontaneous breath. It is special function only MEKANT. It explain figure of next page.



3

Operating Instructions

3.4.1.6 P_{SIMV}

P_{SIMV} (Pressure based Synchronized Intermittent Mandatory Ventilation) : This mode is that MEKANT control artificial respiration by spontaneous breath and pressure.



NOTE :

MEKANT synchronized spontaneous breathing or in case of without spontaneous breathing signal MEKANT ventilate by pressure.



3

Operating Instructions


No	Item	Setting Range	Unit	STEP	Explanation
㉑	P _{INSP}	5 ~ 10	cmH ₂ O	1	Setting the patient pressure to once breath.
㉒	T _{INSP}	0.2 ~ 9.9	sec	0.1	Setting the inspiration Time.
㉓	R _{RESP}	2 ~ 120	BPM	1	Setting number of breath per minute.
㉔	PEEP	0 ~ 45	cmH ₂ O	1	Setting PEEP pressure. (Base Pressure)
㉕	O ₂	21 ~ 100	%	1	Setting oxygen density.
㉖	T _{RISE}	FAST,MID, SLOW	-	-	It sets the flow speed of gas to patient
㉗	E _{SENS}	10 ~ 60	%	10	Setting sensitivity of expiration by ratio of expiration capacity.
㉘	TRIGGER Type	Pressure /Flow	cmH ₂ O /LPM		It set the type of trigger. It is Pressure or Flow type.
㉙	TRIGGER	0.5 ~ 20	cmH ₂ O /LPM	0.5	Setting to detect spontaneous breathing signal, There are pressure and Flow type. User can select one.
㉚	PPS	0 ~ 45	cmH ₂ O	1	Setting pressure of patient when it is working spontaneous breathing.
㉛	RUN	-	-	-	Execute PSIMV Mode.
㉜	RETURN	-	-	-	Go back to previous menu

Table ㉙-9. P_{SIMV} Menu tree

SIMV (Volume Control or Pressure Control) only Mandatory Breath

SIMV mode is able to work Mandatory Breath and Spontaneous Breath. MEKANT's SIMV mode is difference with other equipments. MEKANT's SIMV mode guarantees at least one Mandatory Breath per each cycle. This Mandatory Breath is one of the alternative Patient Initiated Mandatory (PIM or Assisted Mandatory) Breath or Ventilator Initiated Mandatory (VIM or Controlled Mandatory) Breath.


As bellow figure each SIMV Breath Cycle has two parts one is Mandatory Interval (T_m) and other one is Spontaneous Interval (T_s). After receive PIM Breath signal the Mandatory Interval was terminated and than to start Spontaneous Interval (T_s). If PIM Breath does not happen until Mandatory Interval terminate It will finish Mandatory Interval and than change to Spontaneous Interval with VIM Breath.



MEKANT SIMV mode's SIMV Breath Cycle depend on number of breath. For example if set Respiratory Rate is 12 SIMV Breath Cycle is 5 sec. Mandatory Interval is taken smaller value maximum 1.5 sec or 0.25x SIMV Breath Cycle. Mandatory Interval will be shorter by PIM Breath.

3.4.1.7 T_{BILEVEL}

T_{BILEVEL} : It works same as SPONT mode and it has Upper PEEP.



No	Item	Setting Range	Unit	STEP	Explanation
㉑	T/H _{PEEP}	0.2 ~ 24	sec	0.1	Setting High PEEP Time (Upper pressure holding time).
㉒	R _{RESP}	2 ~ 120	BPM	1	Setting patient's number of breath.
㉓	PEEP _{HIGH}	0 ~ 45	cmH ₂ O	1	Setting High PEEP (Upper pressure).
㉔	PEEP _{LOW}	0 ~ 45	cmH ₂ O	1	Setting PEEP (Base Pressure of patient).
㉕	O ₂	21 ~ 100	%	1	Setting oxygen density of air
㉖	TRIGGER Type	Pressure /Flow	cmH ₂ O /LPM		It sets the type of trigger. It is Pressure or Flow type.
㉗	TRIGGER	0.5 ~ 20	cmH ₂ O /LPM	0.5	Setting to detect spontaneous breathing signal, There are pressure and Flow type. User can select one.
㉘	PPS	0 ~ 45	cmH ₂ O	1	Setting pressure of patient when it is working spontaneous breathing.
㉙	RUN	-	-	-	Execute TBILEVEL Mode.
㉚	RETURN	-	-	-	Move to previous menu.


Table 3-10. T_{BILEVEL} Menu tree

T-BiLevel Mandatory Breath


MEKANT's 'BiLevel mode is complex mode it may happened Mandatory Breath and Spontaneous Breath. Because it works as SIMV mode but it difference because Spontaneous range not fixed. So, anytime it is able to Spontaneous Breath.

'BiLevel has two level of PEEP to supply patient. The upper PEEP is $PEEP_H$ and lower PEEP is $PEEP_L$. Mandatory Breath's Inspiration is from $PEEP_L$ to $PEEP_H$ and Expiration is from $PEEP_H$ to $PEEP_L$. So 'BiLevel mode's Mandatory Breath is Pressure Control Ventilation. Each change time of before and after small zone is changing zone, this zone purpose is to detect patient breath if it detect patient breath it would adjust little early or late to start change inspiration and expiration for synchronizing with patient breathing.

Setting parameters are Respiratory Rate, $PEEP_H$, $PEEP_L$, $PEEP_H$ time (or I:E ratio , $PEEP_L$ time), FAP and so on, In case of Pressure Support Ventilation there is more add Support Pressure, FAP, Essens and so on.



MEKANT's 'BiLevel mode supplies $1.5 \text{ cmH}_2\text{O}$ Pressure in Spontaneous Breath. If user want to more high pressure support it will use Pressure Support Ventilation setting. If Support Pressure + $PEEP_L$ value is higher than $PEEP_H$ It will be working pressure support ventilation difference than Spontaneous breath works Pressure Support Ventilation. Refer to the figure below.





'BiLevel mode is able to work Spontaneous Breath any time it reduce fighting between ventilator and patient. It can reduce sedation so patient gets over more sort

All of Spontaneous Breath supplies Pressure Support Ventilation so patient feel more comfortable. If it using opposite I:E ratio it means very sort Expiratory Time it is Airway Pressure Release Ventilation(APRV) status. This APRV is one of best way to apply ARDS (Adult Respiratory Distress Syndrome).



3.4.1.8 NIV

NIV: Mask type ventilation, it works similar PACV mode with air leakage.





No	Item	Setting Range	Unit	STEP	Explanation
㉑	P _{INSP}	5 ~ 10	cmH ₂ O	1	Setting the patient pressure to once breath.
㉒	T _{INSP}	0.2 ~ 9.9	sec	0.1	Setting the inspiration Time
㉓	R _{RESP}	2 ~ 120	BPM	1	Setting number of breath per minute.
㉔	PEEP	0 ~ 45	cmH ₂ O	1	Setting PEEP pressure (Base Pressure)
㉕	O ₂	21 ~ 100	%	1	Setting oxygen density.
㉖	TRIGGER Type	Pressure /Flow	cmH ₂ O /LPM		It set the type of trigger. It is Pressure or Flow type.
㉗	TRIGGER	0.5 ~ 20	cmH ₂ O /LPM	0.5	Setting to detect spontaneous breathing signal, There are pressure and Flow type. User can select one.
㉘	RUN	-	-	-	Execute NIV Mode.
㉙	RETURN	-	-	-	Move to previous menu.

Table 3-11. NIV menu tree

3.4.1.1 PRVC


PRVC (Pressure Regulated Volume Controlled Ventilation) : Breaths are delivered mandatorily to assure present volumes, with a constant inspiratory pressure continuously adapting to the patient's condition. The flow pattern is decelerating.



Initial operation of PRVC (1st Breathing)

3

Operating Instructions



No.	Item	Setting Range	Unit	Step	Explanation
Ⓐ	V _{TIDAL}	15 ~ 2000	ml	10	It decides volume of 1 breath. NOTE : It set the volume of patient weight. It is 10ml per 1kg. (You can set the volume per 1kg in [SYSTEM SET - USER SETUP - B/W FACTOR] Menu.)
Ⓑ	T _{INSP}	0.2 ~ 9.9	sec	0.1	Setting the inspiration Time
Ⓒ	R _{RESP}	2 ~ 120	BPM	1	It sets the respiratory rate per 1 minute.
Ⓓ	PEEP	0 ~ 40	cmH ₂ O	1	It sets the PEEP pressure.(Base pressure of patient)
Ⓔ	O ₂	21 ~ 100	%	1	It sets the percent of O ₂ to patient.
Ⓕ	TRISE	FAST, MID, SLOW	-	-	It sets the flow speed of gas to patient
Ⓖ	E _{SENS}	10 ~ 60	%	10	It sets the sensitivity of one breath completion of patient during to setting of percent of patient expiration volume.
Ⓗ	TRIGGER Type	Pressure /Flow	cmH ₂ O /LPM		It sets the type of trigger. It is Pressure or Flow type.
Ⓘ	TRIGGER	0.5 ~ 20	cmH ₂ O /LPM	0.5	Setting to detect spontaneous breathing signal, There are pressure and Flow type. User can select one.
⓵	PLIMIT	3 ~ 80	cmH ₂ O	30	It sets the patient pressure limit to set over pressure.
Ⓚ	RUN	-	-	-	It works PRVC Mode.
Ⓛ	RETURN	-	-	-	Go back to previous menu

Table 3-12. PRVC menu tree description

3

Operating Instructions

3.4.1.10 AUTO

AUTO : Emergency Quick Start function it requirement patient body weight. It works PACV Mode. This menu purpose is that in case of emergency patient user have to start quickly so it has minimum stability. All parameters are setup default and user can execute with only body weight. But this mode can work one hour every one hour it displayed warning message. The warning message is that “CAUTION : AUTO MODE”.



No	Item	Setting Range	Unit	STEP	Explanation
Ⓐ	RUN	-	-	-	Execute AUTO Mode.
Ⓑ	RETURN	-	-	-	Move to previous menu.

Table 3-13. AUTO menu tree description


3

Operating Instructions

3.4.2 Alarm Set

Each parameter has own alarm limit. User can set alarm limitation through moving the encoder on the key front.

Leak volume setting value affect airway leakage checking. If you set 100(ml) leak volume, ventilation operates. If inspiratory volume is bigger than expiratory volume about 100ml, [AIR LEAK] alarm sounds and message displays.



3

Operating Instructions

No	Item	Setting Range	Unit	STEP	Explanation
㉓	P _{RESS} HIGH	5 ~ 12.	cmH ₂ O	1	Upper limit of patient's pressure.
㉔	P _{RESS} LOW	0 ~ 5.	cmH ₂ O	1	Lower limit of patient's pressure.
㉕	V _{TIDAL} HIGH	20 ~ 2500, OFF	ml	10	Upper limit of patient's capacity.
㉖	V _{TIDAL} LOW	0 ~ 2500	ml	10	Lower limit of patient's capacity.
㉗	O ₂ HIGH	21 ~ 100	%	1	Upper limit of oxygen's density.
㉘	O ₂ LOW	0 ~ 10	%	1	Lower limit of oxygen's density.
㉙	R _{RESP} HIGH	2 ~ 12.	BPM	2	Upper limit of number of breath per minute.
㉚	V _{MIN} LOW	0.0 ~ 99	LPM	0.1	Lower limit of patient's breathe of capacity.
㉛	L _{EAK} VOLUME	10 ~ 500	ml	10	Difference breathe of capacity between Ins. And Exp. Volume
㉜	RETURN	-	-	-	Move to previous menu.

Table 3-14. ALARM SET menu tree


3

Operating Instructions

3.4.3 Trend Mode

TREND MODE display vital sign from the patient last 12 hours. As the following Figure, it displays table and graph, user can move indicator bar (white vertical line) and check each time patient vital sign. The data stored by one minute, it maximum 12 hours can be stored, when ever user press "CLEAR" key user can be returned VENT MODE



In the TREND MODE user press encoder the trend interval will be change as 1 min, 2 min, 3 min, 4 min and 12 min.



3

Operating Instructions

NOTE :

When user goes from VENT MODE to TREND MODE the icon will be changed from  to  on the status bar.




Figure 9-9. TREND MODE screen

3


Operating Instructions

3.4.4 SYSTEM SET

SYSTEM SET menu has BED number, Date/Time, Sound volume, Graph Scale and CALIBRATION menu when you change parameter using encoder.

No.	Item	Mode	Explanation	Refer.
Ⓐ	SOUND	-	Adjust sound volume.	㉓.㉓.4.1
Ⓑ	NEBULIZER	-	Setting operating time of nebulizer(10 ~ 180 minute)	㉓.㉓.4.2
Ⓒ	GRAPH SETUP		Graph setup menu(Auto scale set, Pressure Limit, Flow Limit, Volume Limit, Trigger Display)	3.4.4.3
Ⓓ	USER SETUP	-	User setup menu (B/W Factor, BED Number, O2 Sensor, and so on)	㉓.㉓.4.4
Ⓔ	SYSTEM SETUP	-	System setup menu (System Information, System Self Calibration)	㉓.㉓.4.5
Ⓕ	RETURN	-	Move to previous menu.	

Table ㉓-15. SYSTEM SET menu tree




3

Operating Instructions

3.4.4.1 SOUND

Adjust sound volume.



No	Menu	Range	Unit	STEP	Explanation
Ⓐ	VOLUME	10 ~ 100	%	10	Setting alarm sound volume
Ⓑ	RETURN	-	-	-	Move to previous menu.


Table 3-16. VOLUME menu tree

3

Operating Instructions

3.4.4.2 NEBULIZER

Setting nebulizer operating time in minute.



No	Menu	Range	Unit	STEP	Explanation
Ⓐ	Nebulizer	10 ~ 180	min	10	Setting operating time of nebulizer.
Ⓑ	RETURN	-	-	-	Move to previous menu.

Table 3-17. NEBULIZER menu tree


3

Operating Instructions

3.4.4.3 GRAPH SETUP

It is setting for graph. There are Auto Scale Set, Pressure Limit, Flow Limit, Volume

Limit and Trigger Display.



No	Menu	Range	Unit	STEP	Explanation
Ⓐ	Auto Scale Set	On/Off	-	1	Setting graph auto scale.
Ⓑ	Pressure Limit	10~ 100	cmH2O	5	Setting Pressure Graph Highlimit.
Ⓒ	Flow Limit	10~100	L/min	5	Setting Flow Graph Highlimit.
Ⓓ	Volume Limit	50~3000	ml	50	Setting Volume Graph Highlimit.
Ⓔ	Trigger Display	On/Off	-	1	Setting Trigger Mark Display on graph.
Ⓕ	RETURN	-	-	-	Move to previous menu.

Table ㄿ-18. GRAPH SETUP menu tree

3


Operating Instructions

Auto Scale Set Menu sets highlimit value on graph. If Auto Scale Set is ON, the graph highlimit value will control automatically. If Auto Scale Set is OFF, user sets graph highlimit value.

Pressure, Flow and Volume limit Menu set highlimit value on each parameter.

Trigger Display Menu set trigger mark on graph. If Trigger Display is ON, trigger mark displays on graph at occurring spontaneous breathing.

It is setting for user. There are B/W(Body weight) Factor, BED Number, PATIENT ID, O2 Sensor Enable/Disable, O2 Sensor Calibration, Language, HL7 INTERVAL and Date/Time Setting.



No	Menu	Range	Unit	STEP	Explanation
Ⓐ	B/W Factor	5 ~ 15	MI/Kg	'	Volume / Body weight
Ⓑ	BED Number	1~99		1	Setting BED Number
Ⓒ	PATIENT ID				Setting Patient ID.
Ⓓ	O2 Sensor				Setting O2 Sensor. Enable/Disable
Ⓔ	O2 Calibration				Calibrate O2 Sensor.
Ⓕ	LANGUAGE				Setting Language(English,German,Spanish)
Ⓖ	HL7 INTERVAL				Setting HL7 data sending interval.
Ⓗ	Date/Time				Setting Date and time.
Ⓘ	RETURN	-	-	-	Move to previous menu.


Table 19-19. USER SETUP menu tree

3

Operating Instructions

3.4.4.4.1 B/W Factor

Setting body weight of patient in ml/kg.




No	Menu	Range	STEP	Explanation
Ⓐ	B/W Factor	5 ~ 15	1	Volume / Body weight
Ⓘ	RETURN	-	-	Move to previous menu.

Table 20. B/W Factor menu tree

3.4.4.4.2 BED Number

Setting Bed Number.



No	Menu	Range	Unit	STEP	Explanation
Ⓐ	BED Number	1~99		1	Setting Bed number.
Ⓑ	RETURN	-	-	-	Move to previous menu.

Table ㄿ-21. BED Number menu tree

3

Operating Instructions

3.4.4.4.3 PATIENT ID

Setting patient ID. The patient ID is used for HL7 data format.



No	Menu	Range	Unit	STEP	Explanation
㉑	PATIENT ID	0 ~ 99999999			Setting Patient ID.
㉒	RETURN	-	-	-	Move to previous menu.

Table ㉒-22. PATIENT ID menu tree

3


Operating Instructions

3.4.4.4.4 O2 Sensor

Setting O2 Sensor Enable / Disable.

O₂ Sensor Enable : It is default. It is setting O₂ sensor is normal.

O₂ Sensor Disable : Setting O₂ sensor is abnormal or O₂ Gas is empty.



No	Menu	Range	Unit	STEP	Explanation
Ⓐ	O ₂ Sensor	ENABLE/ DISABLE			Setting O ₂ Sensor.
Ⓑ	RETURN	-	-	-	Move to previous menu.

Table ㄿ-23. O₂ Sensor menu tree


3

Operating Instructions

3.4.4.4.5 O₂ CALIBRATION

Calibrating O₂ Sensor. Push the encoder button and O₂ calibration start.

After few seconds O₂ calibration is stopped.



No	Menu	Range	Unit	STEP	Explanation
Ⓐ	O ₂ Sensor Calibration	START			Calibration O ₂ Sensor
Ⓑ	RETURN	-	-	-	Move to previous menu.

Table Ƴ-24. O₂ Calibration menu tree

3

Operating Instructions

3.4.4.4.6 LANGUAGE

Setting the Language. Supported language is English, German, Spanish,

Chinese, Polish, French and Russian. Default language is English.



No	Menu	Range	STEP	Explanation
①	LANGUAGE	ENGLISH/GERMAN /SPANISH/CHINESE /POLISH/FRENCH/RUSSIAN		Setting the language.
②	RETURN	-	-	Move to previous menu.

Table 2-25. Language menu tree


3

Operating Instructions

3.4.4.4.7 HL7 INTERVAL

Setting HL7 data sending period. The period is 30sec, 1,5,10,30 and 60min.

Default value is 1min.



No	Menu	Range	Unit	STEP	Explanation
㉑	HL7 INTERVAL	30Sec, 1, 5, 10, 30, 60 Min			Setting sending data interval.
㉒	RETURN	-	-	-	Move to previous menu.


Table ㉑-26. HL7 Interval menu tree

3

Operating Instructions

3.4.4.4.8 DATE/TIME

Setting date and time.



No	Menu	Range	Unit	STEP	Explanation
Ⓐ	YEAR	2000~2099			
Ⓑ	MONTH	1 ~ 12			
Ⓒ	DATE	1 ~ 31			
Ⓓ	HOUR	0 ~23			
Ⓔ	MINUTE	0 ~ 59			
Ⓕ	RETURN	-	-	-	Move to previous menu.

Table Ɣ-27. DATE/TIME menu tree


3

Operating Instructions

3.4.4.5 SYSTEM SETUP

Setting system setup menu. There are system information, system

calibration and expiratory gas kind



No	Menu	Range	Unit	STEP	Explanation
①	SYSTEM INFORMATION				Main S/W Version, Module S/W Version, Working Time
②	SYSTEM SELF CALIBRATION	-	-	-	Calibration Flow Sensor.
③	CALIBRATION INSP-EXP				Calibration Inspiration-Expiration Sensor Balance
④	EXP. GAS				Setting expiratory gas kind.(Dry or Humid)
⑤	RETURN	-	-	-	Move to previous menu.

Table 3-28. SYSTEM SETUP menu tree

3

Operating Instructions

3.4.4.5.1 SYSTEM INFORMATION

There are SW version information and working time information.



No	Menu	Range	Unit	STEP	Explanation
Ⓐ	MAIN S/W VERSION				
Ⓑ	MODULE S/W VERSION	-	-	-	
Ⓒ	WORKING TIME				Total working time
Ⓓ	RETURN	-	-	-	Move to previous menu.

Table ㄿ-29. SYSTEM INFORMATION menu tree

3

Operating Instructions

3.4.4.5.2 SYSTEM SELF CALIBRATION

It is calibration mode for flow sensors. It has 2 flow sensor calibrations. One is offset calibration of Inspiration flow sensor and the other is offset

calibration of expiration flow sensor.

CAUTION : Before calibration, MV2000 must be separated from patient.



No	Menu	Range	Unit	STEP	Explanation
①	CALIBRATION				Inspiration and expiration flow sensor offset calibration.
②	RETURN	-	-	-	Move to previous menu.

Table 30. CALIBRATION menu tree

3


Operating Instructions

3.4.4.5.3 CALIBRATION INSP-EXP

It is calibration mode for striking balance between inspiration flow sensor and expiration flow sensor. In F-V graph, when the start point of wave falls in with the end point of wave, the calibration meets the case. BF means Base Flow.

BF is the difference of inspiratory flow and Expiratory flow. It set that [Exp. Gas] is [DRY] status.

CAUTION : Before calibration, the test lung must be connected to ventilator.



3

Operating Instructions

3.4.4.5.4 EXP. GAS

You select the expiratory gas kind of dry gas or humid gas.

If you use humidifier. you select the humid. Default set is dry.



No	Menu	Range	Unit	STEP	Explanation
①	EXP. GAS				Select the respiratory gas kind of dry or humid.
②	RETURN	-	-	-	Move to previous menu.

Table 3-31. EXP. GAS menu tree

3

Operating Instructions

3.4.4.5.5 ALTITUDE

Setting altitude height in meter.



No	Menu	Range	Unit	STEP	Explanation
①	ALTITUDE	-100 ~ 5000	meter	20	Setting height of altitude
②	RETURN	-	-	-	Move to previous menu.

Table 3-32. ALTITUDE menu tree

3.4.4.5.6 AIR TEMPERATURE

Setting air temperature in centigrade.



No	Menu	Range	Unit	STEP	Explanation
①	TEMPERATURE	0 ~ 50	centigrade	1	Setting air temperature of centigrade
②	RETURN	-	-	-	Move to previous menu.

Table ㉙-33. AIR TEMPERATURE menu tree

Alarms

SECTION 4

4.1 General Information

All alarm indicates condition of each situation on LED display window. LCD monitor shows their warning message.




Figure 4-1. LCD Display Window

4

Alarms

4.2 Alarm Message

If alarm message is indicated, check the current mode and solve the problem

according to it.

Item	Alarm range
High pressure	$\Delta \text{ cmH}_2\text{O} \sim 12 \cdot \text{ cmH}_2\text{O}$
Low pressure	$\cdot \text{ cmH}_2\text{O} \sim \Delta \cdot \text{ cmH}_2\text{O}$
High tidal volume	20 ml $\sim 250 \cdot \text{ ml}$, OFF
Low tidal volume	$\cdot \text{ ml} \sim 250 \cdot \text{ ml}$
High minute volume	1 $\sim 99\text{L}$
Low minute volume	0 $\sim 98\text{L}$
High respiration rate	2 $\sim 120 \text{ bpm}$
Leak Volume	10 $\sim 500 \text{ ml}$
Apnea	OFF, 10 $\sim 60 \text{ sec}$
O ₂ % high	21 $\sim 100\%$, OFF
Power loss	AC/DC power loss
Low Air/O ₂ pressure	
Obstructed tube	
Airway leak	
Vent in-operation	
Open Circuit	

Table 4-1. Manual alarm description

1	Low Battery	Discharge the internal battery. Connect the AC power.
2	Low Air/O ₂ Pressure	Supplied pressure of Air and O ₂ is low. Check the Air and O ₂ line.
3	Obstructed tube	It is problem between breathing circuit and patient. Check the breathing circuit.
4	Airway leak	In breathing circuit airway leaks. Check the breathing circuit.
5	Vent in- operation	The ventilator is abnormal status. Check the all ventilator system.
6	Leak volume	Difference volume between inspiration and expiration volume. Affected to airway leakage checking.
7	Open Circuit	It is problem between breathing circuit and patient. Check the breathing circuit.

Operating theory

APPENDIX A

Pneumatic System of MEKANT is composed of 7 pcs of gas circuit. One is for Oxygen,

the other is for Air. Proportional solenoid Valve (PSOLs) is the core of Pneumatic System. The PSOLs is controlled accurately by Flow sensor and Microprocessor so that provide accuracy respiration to patient.

MEKANT's Pneumatic System is divided into inspiratory module, expiratory module and patient circuit system. inspiratory module is subdivided into Gas supply, Flow control, Safety Valve and Inspiration monitoring system.




Figure A-1. PNEUMATIC Block Diagram

Blending Part

Gas supply system provides the equipment high pressure air (120-150 psi) and Oxygen from external source of supply, and sense that supplying pressure is not adequate. The high supplied pressure air and Oxygen is decompressed into 100 psi through each regulator. Adjustment pressure of Regulator should be measured by pressure valve (SV1, SV2).

Appendix A

Operating theory

Moisture which may be occurred in the high supplied pressure is removed by water

trap in Air Inlet. Two peaces of Check Valve (CV₁, CV₂) prevent backflow between air and oxygen.

Inspiration Part

Flow control system is the most important part in Pneumatic System. This is composed of Flow sensor (DP2) and PSOLs (PV₁). This system adjusts Tidal volume, Peak Flow, Flow pattern, Respiratory Rate, Oxygen Percentage. Control of PSOLs is done on measured value by Microprocessor. Flow sensor is Differential Pressure type.

Safety System

This is divided in Safety Valve section and Inspiration monitoring section. Safety valve section release excess pressure in the equipment and make a patient breathe with room air as a passageway.

Safety valve open and release pressure in case interior pressure reach to 114 cmH₂O. And it is a role of passageway to connect to room air in case equipment is defect or no power, patient's expiratory limb is closed.

Inspiration monitoring section monitor gas pressure and oxygen concentration provided the patient. Monitor is measured by Absolute Pressure Transducer (PA) which is attached inspiratory Pressure Transducer (P2). Oxygen is done by Oxygen sensor.

Appendix A

Operating theory

Expiration Part

This section is for Exhalation. PEEP/CPAP is controlled by Exhalation Valve during exhalation. Exhalation Valve of MEKANT is Active type and is similar with PSOL. This active exhalation valve makes BiLevel mode possible and prevent pressure increasing suddenly.

Flow sensor (DP4) measure Flow and Volume of patient's expiration gas, and Expiratory pressure Transducer(P4) measure pressure.

Patient Circuit System


Patient circuit system transits gas which is adjusted in Flow control system to patient, and makes patient's exhalation gas to exhalation system.

There is each Bacteria Filter on the connecting part of patient circuit system and equipment. So they make patient circuit system separate from the ventilator. It makes gas to be provided the patient heat and humidify in case humidifier is used.

Communications

APPENDIX B

VGA



Pin	Name	Direction	Description
①	RED	→	Red Video (V _Δ ohm, .V _{p-p})
②	GREEN	→	Green Video((V _Δ ohm, .V _{p-p})
③	BLUE	→	Blue Video(V _Δ ohm, .V _{p-p})
④	ID _Y	←	Not connected
⑤	GND	—	Ground
⑥	RGND	—	Red Ground
⑦	GGND	—	Green Ground
⑧	BGND	—	Blue Ground
⑨	KEY		Not connected
⑩	SGND	—	Sync Ground
⑪	ID _B	←	Not connected
⑫	ID ₁ or SDA	←	Not connected
⑬	^T HSYNC or CSYNC	→ →	Horizontal Sync (or Composite Sync)
⑭	VS _Y NC	←	Vertical Sync
⑮	ID _R or SCL		Not connected

App
endi


Table B-1. VGA serial port pin out

x B

Communications

RS-232 (serial) port

A 9-pin male connector configured as data terminal equipment (DTE). Allowable current is 0.1A at 5V DC (maximum).




Pin	Name	Direction	Description
①	CD (Data Carrier Detect)		Not connected
②	RxD (Receive Data)	←	Receive data (RxD)
③	TxD (Transmit Data)	→	Transmit data (TxD)
④	DTR (Data Terminal Ready)		Not connected
⑤	SG(Signal Ground)	—	Ground(GND)
⑥	DSR(Data Set Ready)		Not connected
⑦	RTS(Request to Send)		Not connected
⑧	CTS(Clear to Send)		Not connected
⑨	RI(Ring Indicate)		Not connected

Table B-1. RS-232 serial port pin out


Communication capabilities (cont)

NOTE :

The connection of accessories or equipment to the ventilator's RS-232 and remote alarm ports requires electrical separation of conductive earth (ground) connections between different items of equipment of a system. The system cable assembly shielding can interconnect the ventilator and other equipment, increasing the risk of excessive enclosure leakage current from external equipment. For safe connection, the shield conductor must have separation of 4mm between the equipment and ventilator conductors, as show below:



LCD MONITOR POWER SUPPLY Connector













Pin	Name	Direction	Description
①	NC		Not connected
②	GND		GND
③	+VCC(+12V/5A)		12V/ Max 5A, OUTPUT LCD MONITOR POWER SUPPLY
④	NC		Not connected









Table B-۳. LCD Monitor Power Connector port pin out

Parts accessories

APPENDIX C

Standard

Item	Photo	Part Number	Q'ty	Description
Circuit Arm		14.4.35 2.	1	
High Pressure Hose		1414... 2.	1	Air Gas / Yellow
		1414... 3.	1	O ₂ Gas / Green
Display Panel Cable		1424.16 6.	1	
O ₂ Sensor		14.5.01 ..	1	
PEEP Silicone Plate		14.6.16 3.	2	
Test Lung		14.6.16 2.	1	
Operation Manual		1417.03 8.	1	
Power Cord		14.6.05 3.	1	3m, 220VAC
Ground Cable		1413.02 ..	1	

Fuse		۱۴۱۳.۰۱ ۷.	۲	INLET ۵.۰T ۶.۳A ۲۵.۷۲.۰mm
Wrench Driver		۱۴.۸.۰۷ ۲.	2	
Battery		۱۴.۸.۰۷ ۲.	1	
Circuit Arm Fix Knob		۱۴.۴.۱.۰ ..	1	
53010(O-Ring)		۱۴۲۱.۰.۴ ..	12	O-Ring for Peep Module
53030(O-Ring)		۱۴۲۱.۰.۴ ۱.	3	O-Ring for Peep Module
Monitor Fix Bolt		۱۴.۸.۰.۳ ۳.	3	M5x10mm Wrench Bolt(For Monitor)
Spring Washer		۱۴.۸.۰.۷ ..	3	M5 Bolt Spring Washer
Body Cover		۱۵.۰.۰.۱ ۶.	1	

Contact Us

Service calls

Followings are telephone numbers and addresses for contacting various service, product supplies and sales personnel.
Reaching open a service call with MEKICS Co., Ltd., contact the numbers listed below.

MEKICS Co., Ltd.

5F, A-BLDG, WoolimLion, #144-3, Sangdaewon-dong,

Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea 462-725

Tel : + 82-70-7119-2500 Fax : + 82-31-735-2761

URL : www.mek-ics.com

e-mail : service@mek-ics.com

- ※ In the event of a malfunction or failure, contact Service Dept. of MEKICS Co., Ltd. along with the model name, serial number, date of purchase and explanation of failure.

Authorized Representative

Pedro Pestana, Lda

Nif : 507662407

Fax : +351214414138

Telef. : +351914881444

Adress: Rua Bernardo Marques, nº3, 2º Dtº

2770-199 Paco de Arcos