neurotech®
embrace life.

neurotech
Vital™

multipath™
neurotech

+ INSTRUCTION MANUAL
For the treatment of stress Incontinence

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Table of Contents

1.0	Safety Information	3
1.1	Warranty	3
1.2	Contraindications, Warnings & Precautions	4
2.0	NEUROTECH VITAL™ Pack Contents	6
2.1	Literature Pack	7
2.2	Program Information	7
2.3	The NEUROTECH VITAL™ Controller	8
2.4	The NEUROTECH VITAL™ Unit Cradle	8
2.5	Description of the NEUROTECH VITAL™ Controls	9
2.6	Controller Display Icons	9
3.0	Setting Up The Controller	10
3.1	Electrode Fault Detector	11
3.2	Setting the Date and Time	12
4.0	Program 1 - NMES	12
4.1	Before the Treatment	12
4.2	Preparing the Garments	13
4.3	Connecting the Garments to the Controller	14
4.4	Running Program 1	14
4.5	Pausing a Treatment Session	15
4.6	Concluding Program 1	15
5.0	Program 2 - Biofeedback Assisted Exercise	16
5.1	Program 2 Set-Up (Option A and Option B)	16
5.2	Running Program 2	17
5.3	Pausing a Treatment Session	17
5.4	Concluding Program 2	18
6.0	Program 3 - Biofeedback and NMES	19
6.1	Running Program 3	19
6.2	Pausing a Treatment Session	20
6.3	Concluding Program 3	20
7.0	Program 4 - biofeedback Triggered NMES (ETS)	21
7.1	Running Program 4	21
7.2	Pausing a Treatment Session	22
7.3	Concluding Program 4	22
8.0	Problem Solving Guide	23
9.0	System Maintenance	25
9.1	Caring for Your NEUROTECH VITAL™ Garments	26
9.2	Caring for Your NEUROTECH VITAL™ Electrodes	26
10.0	Technical Specifications	27

1.0 Safety Information

Validity

The information and technical data contained in this document relates to the NEUROTECH VITAL™ muscle stimulator and accompanying cradle provided with this manual. Each NEUROTECH VITAL™ controller and cradle is assigned a serial number, which is located on the front underside of the controller and on the underside of the cradle.

The information and technical data disclosed in this document are proprietary to Bio-Medical Research Ltd. and may be used and disseminated only for the purposes of and to the extent specifically authorised in writing by the company.

This manual is intended for the guidance of the patient in the correct application of the product.

Important:

You should not begin using this device until you have received full instructions on its correct usage from your clinician.

Restrictions

The sale and/or operation of this equipment is subject to law in various countries. Compliance with this legislation rests with the importer, dealer or user of the equipment as appropriate.

Intended Use

- NEUROTECH VITAL™ is intended for rehabilitation of the pelvic floor muscles in stress urinary incontinence. Stress urinary incontinence is described as the leakage of urine associated with an increase in intra-abdominal pressure, for example, during a cough, sneeze, laugh or physical exertion.

1.1 Warranty

Bio-Medical Research Limited (BMR Ltd.) warrants (a) that your NEUROTECH VITAL™ controller and cradle will be free from defects in material and workmanship for a period of two years from the date of original purchase and (b) that the NEUROTECH VITAL™ garments will be free from defects in material and workmanship for a period of six months from the date of original purchase. These warranties do not apply to the batteries supplied with your Neurotech Vital Product and exclude consumables such as the electrodes.

Should your controller or cradle, or either garment develop a fault within the warranty period, BMR Ltd. will replace or repair them free of charge, provided the controller, cradle and/or garment:

- have been used for their intended purpose and in the manner described in this instruction manual;
- have not been connected to an unsuitable power source;
- have not been subjected to misuse or neglect;
- have not been modified or repaired by anyone other than an approved **neurotech®** distributor or agent.

This warranty complements existing applicable national guarantee obligations and does not affect a consumers statutory rights.

1.2 Contraindications, Warnings & Precautions

Contraindications:

- Do not use this device if you have a cardiac pacemaker, implanted defibrillator, or other implanted electronic device because this may cause electric shock, burns, electrical interference, or death.
- Do not use this device if you are pregnant. The safety of electrical stimulation during pregnancy has not been established.
- Do not use the device if you have any implanted metal devices in the abdominal or pelvic areas, for example hip prostheses.

Warnings:

- Do not apply stimulation on the front or the sides of the neck, because this could cause severe muscle spasms resulting in closure of the airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- Do not apply stimulation across the chest, because the introduction of electrical current into the chest may cause rhythm disturbances to the heart, which could be lethal.
- Do not apply stimulation over open wounds or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions (e.g. phlebitis, thrombophlebitis, varicose veins).
- Do not apply stimulation over, or in proximity to, cancerous lesions.
- Do not apply stimulation in the presence of electronic monitoring equipment (e.g. cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation device is in use.
- Do not apply stimulation when in the bath or shower.
- Do not apply stimulation while driving, operating machinery, or during any activity in which electrical stimulation can put you at risk of injury.
- Do not use in close proximity to shortwave or microwave therapy equipment which may produce instability in the stimulator output.
- Use this device only with the garments, leads, electrodes, and accessories recommended by the manufacturer.
- Always ensure that the charger/power supply is positioned so that it is easy to unplug the power supply from the socket.
- To isolate this equipment from the mains, the charger/power supply plug should be removed from the socket.
- Apply stimulation only to normal, intact, clean, healthy skin.
- Do not attach any cables or connectors to NEUROTECH VITAL™, other than those supplied by **neurotech®**. Never connect any computer leads to the controller.

Precautions:

- The long-term effects of electrical stimulation are unknown.
- Since the effects of stimulation of the brain are unknown, stimulation should not be applied across the head and electrodes should not be placed on opposite sides of the head.
- Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or the electrical conductive medium.
- Patients with suspected or diagnosed heart disease should follow precautions recommended by their clinician.
- Patients with suspected or diagnosed epilepsy should follow precautions recommended by their clinician.
- Use caution if you have a tendency to bleed internally, such as following an injury or fracture.
- Use caution following recent surgical procedures when stimulation may disrupt the healing process.
- Use caution if stimulation is applied over areas of skin which lack normal sensation.
- Keep this device out of the reach of children.
- Use this device only under the continued supervision of a licensed practitioner.
- The size, shape, and type of electrodes may affect the safety and effectiveness of electrical stimulation and recording.
- Using stimulation electrodes that are too small or incorrectly applied could result in discomfort or skin burns.
- You should contact the manufacturer of this electrical stimulator if you do not know if the electrode you have can be used with this stimulator.
- Electrodes and garments should not be shared with other people. To prevent cross-infection, each person should have their own set of electrodes and garments.
- Replace the electrodes if they no longer stick to your skin properly. See the advice on electrode maintenance (Page 26).

Additional Warnings and Precautions:

- Please wait before using NEUROTECH VITAL™ until:
 - At least six weeks after the birth of your baby (consult a clinician first).
 - At least three months after having a Caesarean section (consult a clinician first).
 - The heavy days of your period have finished.
- Do not use the NEUROTECH VITAL™ stress programs if you have had an IUD fitted, which is made from any electrically conductive material such as copper. Please consult your doctor in this instance.

- Stimulation should not be applied over the carotid sinus nerves, particularly in patients with a known sensitivity to the carotid sinus reflex.
- Although compliant with applicable EMC requirements, this device may still interfere with more sensitive equipment, please move away or switch off this device if necessary.
- Specific patient groups (e.g. the elderly) may require further attention and care.
- Before using the equipment, a visual inspection must be performed (check for a broken enclosure, damaged connector or exposed wires). If there is any doubt, the equipment must not be used.
- Once the treatment has ceased, always remember to switch off the controller.
- If in doubt, always seek medical advice.
- Medical advice must be obtained before use on persons who are insulin dependent diabetics or for persons who are under medical supervision for any cognitive dysfunction.
- When repositioning electrodes during treatment, always pause or stop the controller first.
- Simultaneous connection of a patient to high frequency surgical equipment may result in burns at the site of the stimulator electrodes and possible damage to the stimulator.
- This device can deliver current densities in excess of 2mA/cm² when used at a high intensity with small electrodes.
- Do not use the device with the electrodes positioned on injection sites (of medications/drugs), such as hormone treatment sites.
- Do not touch the electrodes, connectors or lead pins with your fingers while the stimulation session is in progress. Always pause or switch off the controller before handling the electrodes.
- The device must not be used with any other unit which delivers electrical current to the body (e.g. interferential or another muscle stimulator).

Adverse Reactions:

- A small number of isolated skin reactions have been reported by people using muscle stimulation devices, including allergies, a prolonged reddening of the skin and acne.
- On very rare occasions, first-time users of NMES have reported feeling light-headed or faint. We recommend that you have somewhere to sit nearby in case you feel faint.
- Patients should stop using the device and should consult their clinician if they experience adverse reactions from the device.

Note:

- You should always leave at least 1 day between treatment sessions to reduce the potential for muscle fatigue. This restriction may be changed depending on the prescribing clinician.
- An effective treatment should not cause undue discomfort.

Important:

- The safety information provided in this User Manual must be followed.
- Cease using your unit if you are feeling light headed or faint. Consult your doctor if this happens.
- **neurotech®** will not accept responsibility if the guidelines and instructions supplied with this unit are not followed.

IMPORTANT

- **NEVER** connect leads or wires from other equipment such as computers, games, cameras etc. to the controller.
- **NEVER** connect the leads provided with NEUROTECH VITAL™ to any other equipment.
- If in doubt, contact **neurotech®**.

IMPORTANT ELECTRODE INFORMATION

In order to prolong the life of the electrodes, after every few treatment sessions, you should sprinkle a few drops of clean water over the adhesive surface of each electrode with your fingertips and then gently rub it over the complete surface of the electrode. This will clean away debris and help to rehydrate the adhesive surface.

Note: This should only be done with the side of the electrode which adheres to the skin. The side of the electrode which adheres to the garment should be left affixed to the garment.

2.0 NEUROTECH VITAL™ Pack Contents

1. NEUROTECH VITAL™ Controller

A portable neuromuscular electrical stimulator that generates signals which are sent, via electrodes and garments, to activate your muscles. It also records EMG biofeedback data to track your progress. It is rechargeable and will take approximately 5 hours to fully charge.

2. Unit Cradle

The unit cradle is used for charging the controller.

3. Battery Charger

NEUROTECH VITAL™ is supplied with a power supply for charging the battery. The charger will be either a 2-pin or a 3-pin, depending on your country of residence.

4. Left Stress Garment

This conductive garment fastens around the left thigh and buttock area. It has a **blue** connector pin, which is plugged into the stress lead.

5. Right Stress Garment

This conductive garment fastens around the right thigh and buttock area. It has a **red** connector pin, which is plugged into the stress lead.

6. Stress Lead

This connects the controller to the garments.

7. Large NMES Electrodes

These 8 electrodes are for use on the stress garments only.

8. Biofeedback Probe

This is a special electrode for insertion into the vagina to acquire pelvic floor EMG biofeedback data.

9. Surface EMG Biofeedback Electrodes

These are small electrodes which may be used as an alternative to the biofeedback probe during the EMG biofeedback programs.

10. EMG Reference Electrode

This is used during the EMG biofeedback programs to improve the quality of the EMG signal.

11. EMG Biofeedback Lead

This lead is used in conjunction with the EMG reference electrode and either the biofeedback probe or the EMG biofeedback electrodes.

12. Med Gel

Used to lubricate the biofeedback probe prior to insertion.

2.1 NEUROTECH VITAL™ Literature Pack

13. Stress Incontinence Treatment Instructions

These are the instructions specific to the treatment of stress urinary Incontinence.

14. Quick Start Guide

This is an additional guide, which illustrates how to set up and apply the garments for stress incontinence.

2.2 Program Information

NEUROTECH VITAL™ has 4 treatment programs specifically designed for stress urinary incontinence. These can also be combined with a separate urge urinary incontinence kit to treat mixed urinary incontinence.

Stress Urinary Incontinence Programs

Prog.	Treatment	Duration	What you need for this program	Recommended Prescription
1	Stress Incontinence - NMES	30 mins	Controller, stress lead, left & right garments, 8 x NMES electrodes.	1 x 30 min. session per day, 5 days per week, for 12 weeks
2	Biofeedback Assisted Exercise	10 trials*	Controller, EMG biofeedback lead, 1 x reference electrode, Med Gel and biofeedback probe OR 2 small EMG electrodes, if preferred.	1 session per day
3	Biofeedback and NMES	10 trials*	Controller, stress lead, left & right garments, 8 x NMES electrodes, 1 x reference electrode, EMG biofeedback lead, Med Gel and biofeedback probe OR 2 small EMG electrodes, if preferred.	1 session per day
4	EMG Triggered Stimulation (ETS)	30 trials*	Controller, stress lead, left & right garments, 8 x NMES electrodes, 1 x reference electrode, EMG biofeedback lead, Med Gel and biofeedback probe OR 2 small EMG electrodes, if preferred.	1 session per day

* A Trial is defined as one single Work/Rest phase. During the Rest phase, you can relax your muscles and then, during the Work phase, you should actively contract your pelvic floor muscles, with the aim of reaching a pre-defined target EMG level. The duration of each phase is counted down on the controller screen.

2.3 The NEUROTECH VITAL™ Controller

NEUROTECH VITAL™ is a single unit providing Neuromuscular Electrical Stimulation (NMES - also known as Electronic Muscle Stimulation, EMS), Electromyography (EMG) biofeedback and EMG Triggered Stimulation (ETS). Easy-to-follow screens lead you through the set-up process and treatment programs. The controller can be customised to suit your needs and you can bring your controller back to your clinician so they can monitor and record your treatment progress.

The controller has two connector ports.

- Port 1:** is the main output connector, where the stress lead is attached. It also connects the controller to the unit cradle for recharging the battery.
- Port 2:** is where the EMG biofeedback lead is attached to the controller.

NEVER ATTACH COMPUTER LEADS TO THESE PORTS

The NEUROTECH VITAL™ controller has a built-in audio indicator which emits a high-pitched tone when a valid button is pressed and a low-pitched tone when an invalid button is pressed.

Controller (rear view)


2.4 The NEUROTECH VITAL™ Unit Cradle

NEUROTECH VITAL™ is supplied with a unit cradle for charging the controller's battery.

The unit cradle has a battery charger socket and an indicator LED.


Battery Charger Socket: is where the battery charger is attached to the cradle.

LED 1: lights to show that the battery is being charged.

1. To charge the battery, connect **Port 1** of the controller to the cradle connector and plug the charger into the back of the cradle (Fig. 1). Then plug the charger into an electrical socket and switch it on.
2. You should initially charge the battery for approximately 5 hours. It is fully charged when the green battery symbol at the top-right of the screen  stops flashing and remains a solid bar.

Fig. 1

2.5 Description of the NEUROTECH VITAL™ Controls

1. On/Off Switch (I / O)

This button is located on the left-hand side of the controller. Slide the button down to turn the controller on (I). Slide the button up to turn the controller off (O).

2. Screen Button 1 / Screen Button 2 / Screen Button 3 (⊖)

These three buttons are used for various controller commands. The icon for the relevant command will appear on the display above the corresponding button (See below for a complete list of command icons).

3. Up Intensity Button (⬆)

Press this button to increase the stimulation intensity during treatment.

4. Down Intensity Button (⬇)


Press this button to decrease the stimulation intensity during treatment.

5. Navigation Buttons (⊗)

These buttons are used to navigate through the various menus on the controller display. The Up/Down Buttons will move the controller cursor up and down through the menu options on the display. The Left/Right Buttons will decrease (left) or increase (right) any adjustable treatment parameters or controller settings.

6. Pause/Stop Button (⏸)

Press this button to pause a treatment session. Press the button a second time to resume the session. Press and hold this button (for more than 2 seconds) to end a treatment session.

2.6 Controller Display Icons

The following icons appear on the controller display at various times throughout the set-up and treatment.

ICON	DESCRIPTION	ICON	DESCRIPTION
Menu Icons		Screen Button Icons	
	Program		Menu
	History		Select
	Settings		Save
Status Icons			Back
	Program		Cancel
	History		Edit
	Communicating with PC (for use by clinician only)		Key lock
	Therapist Mode (for use by clinician only)		Start EMG Biofeedback Program
	Warning		
	Mute		

3.0 Setting Up The Controller

The following adjustments can be made to the controller settings **BEFORE BEGINNING** a treatment session.

Program: ()

When you switch on the controller, you can select your prescribed program. To do this, press **Screen Button 3** to access the **Menu** (). Then press it again to **Select** () **Programs** (). Use the **Up/Down Navigation Buttons** () to highlight the desired program, and then press **Screen Button 3** again to **Select** that program.

Note: Once selected, the controller automatically changes to that program's start screen.

History: ()

For a quick review of your previous usage, access the **Menu** as above and then use the **Up/Down Navigation Buttons** to highlight the **History** () option. **Select** () that option and use the **Up/Down Navigation Buttons** () to choose one of the following:

- **Last Session**, lets you view details of your last treatment session.
- **Summary**, lets you view a summary of all of your treatments to date.
- **Calendar**, shows you the dates when you performed all of your treatment sessions (in green). The current date is highlighted in red.

Note: Once you have set up the controller, use **Screen Button 1** to go **Back** () to the previous screen. You may have to press this again to get to the program start screen.

Settings: ()

To set the controller to suit your specific needs, you can access the **Menu** as above and then use the **Up/Down Navigation Buttons** to highlight the **Settings** () option. **Select** () that option and use the **Up/Down Navigation Buttons** () to choose one of the following:

- **Volume**, increase or decrease the volume of the beeps emitted by the controller. The range is 0 (silent) to 3 (max). When the Volume option is highlighted, use the **Left/Right Navigation Buttons** () to set the desired level. If '0' is chosen, the Mute icon () will appear at the **top** of the screen.
- **Brightness**, increase or decrease the brightness of the controller screen. The range is 1 (darkest) to 4 (brightest). When the Brightness option is highlighted, use the **Left/Right Navigation Buttons** () to set the desired level.
- **Language**, choose your preferred language by highlighting the language option, pressing **Select** () and then choosing your desired language. The languages available are; German, English, Spanish, French, Italian and Norwegian.

Note: Once you have set up the controller, use **Screen Button 1** to go **Back** () to the previous screen. You may have to press this again to get to the program start screen.

The following adjustments can be made to the controller settings **DURING** a treatment session.

Locking/Unlocking the Intensity:

When you have reached a signal strength you are comfortable with, you can disable the Intensity Buttons, so the strength is not changed accidentally. To do this, press **Screen Button 1** which **Locks** () the Intensity Buttons. You can unlock the buttons by pressing **Screen Button 1** again. The **Lock** symbol will appear at the **top** of the screen when the intensity is locked.

Settings:

You can change the following settings during a treatment:

- **Volume:** As above.
- **Brightness:** As above.
- **Change Program:** there are two reasons why you might change the program during a running treatment; first, if you have started the wrong program and wish to change to the correct one. Second, your clinician may have prescribed that you run part of one program in conjunction with a second program. For example, if they consider 30 minutes of Program 1 is too much for your specific needs, they may recommend that you run 15 minutes of Program 1 and then change to Program 2. In this instance, you must have both the NMES and the biofeedback elements of the treatment set up before starting (i.e. the garments **and** the biofeedback probe or electrodes).

3.1 Electrode Fault Detection

During a treatment session, the controller continuously monitors its connection with the electrodes and can detect if they are in proper connection with your skin. If it detects a poor quality connection, it will alert you that there is a problem. This can be caused by one of the following:

- An electrode is missing from a garment.
- An electrode is not properly covering its silver stud on the garments.
- Poor contact between the electrode(s) and the skin (e.g. liner left on an electrode or the electrode is not flat against the skin).
- Poor electrode condition (i.e. with damaged and/or dry surface following too many sessions).

If any of these issues are detected, the controller will:

- Highlight the specific electrode(s) at fault.
- Stop the stimulation until the problem is fixed.
- Beep repeatedly until the controller is switched off or the problem is fixed.

CHECKING / FIXING ELECTRODE FAULTS

1. If an electrode detection fault occurs during a treatment session, the controller will highlight the faulty electrode(s) in **red** on a diagram, showing where that electrode is positioned on the body (Fig. 2). Please note that there may be more than one electrode concerned (Fig. 3).
2. **Switch off the controller.** Disconnect the stress lead connectors from the garment plugs.
3. Carefully remove the garment(s), as described on Page 15 and move or replace the electrode(s) causing the problem. Make sure none of the liners have been left on the electrodes.
4. Then put the garment(s) back on and press the whole surface of each electrode firmly against the skin. Reconnect the stress lead to the garments and begin the treatment again.
5. If the screen still indicates the same electrode(s) as being faulty, you should try another set of electrodes. If this still doesn't fix the problem, you may have a faulty lead or garment and should contact **neurotech®** or your local distributor for advice.

Important:

The Electrode Fault Detection feature does not work during Program 2 of the Stress Treatments, as there is no NMES stimulation.

NB: NEVER TOUCH THE ELECTRODES WITH YOUR FINGERS WHILE THE TREATMENT IS IN PROGRESS.

Fig. 2



Fig. 3

3.2 Setting the Date and Time

Every controller should be delivered with the battery partially charged. However, if the battery has been fully discharged, you may have to reset the date and time. You will be requested to do this when you switch the unit on.


1. To set the date and time, simply use the **Up/Down Navigation Buttons** (⊗) to select **Date**, **Month**, **Year**, **Hour** or **Minute** and then use the **Left/Right Navigation Buttons** (⊗) to increase or decrease the values as necessary (Fig. 4). You can also decide whether to set the time to a 12h or 24h clock when setting the time.

2. Press **Screen Button 3** to **Save** (💾) your settings.

Note:

- Do not leave the charger plugged into the electrical socket when not in use.

Fig. 4

4.0 Program 1 - NMES

Program 1 is an NMES program. The garments are fitted around the thighs and buttocks and are then connected to the controller so that an electronic signal can be sent via the electrodes to the pelvic floor muscles, strengthening them and helping to treat stress urinary incontinence.

4.1 Before the Treatment

1. Before using NEUROTECH VITAL™, you should ensure the battery is adequately charged. Place the controller into the cradle, connect the charger to the back of the cradle and plug the charger into an electrical socket. You should charge the battery for approximately 5 hours. It is fully charged when the battery symbol at the top-right of the screen (🔋) stops flashing and remains a solid bar.
2. The controller will switch on automatically while charging and a warning will appear on screen saying you should connect an appropriate garment (Fig. 5). As this cannot be done while the unit is charging, simply press any one of the 3 **Screen Buttons** (⊖) beneath the LCD screen to clear this. If it appears again, just press one of the **Screen Buttons** a second time.

Notes:

This warning screen appears every time the controller is switched on except in the following situations:

- a) If the stress lead has already been connected.
- b) When Program 2 has been selected (the garments are not used for Program 2).

Fig. 5

4.2 Preparing The Garments

- Please note that this will only be required for Programs 1, 3 and 4.

N.B. We recommend wearing high-thigh underwear (e.g. a thong), which leaves the buttocks exposed, because the electrodes need to make contact with the lower buttocks and the hips.

1. Clean the skin thoroughly where treatment will be applied. The electrodes will not adhere properly if any dirt, oils, creams or other cosmetics are on the skin. Use warm, soapy water to clean the skin. Dry the skin thoroughly after cleaning.
2. Place the garment for your right-leg on a flat surface with the inner surface facing upwards.
3. Remove the covers from the **patterned side** of the electrodes (Fig. 6). Do not remove the covers from the black side of the electrodes yet.
4. Place the electrodes on the garment in either the A or B positions, as determined in conjunction with your clinician. Press the whole surface of each electrode firmly onto the garment to ensure they adhere properly.
5. Repeat these steps with the garment for your left-leg.
6. Remove the covers from the **black side** of the electrodes. Do not throw these covers away as you have to put them back on the electrodes at the end of the session. With the electrodes facing your body, wrap the right garment firmly around your right leg and fasten. Please refer to the separate **Garment Quick Start Guide** for more detail. Then repeat the same steps, wrapping the left garment around your left leg.
7. When the garments are positioned correctly, the electrodes should lie on the front and rear of your thighs and on the bottom and hips, exactly as depicted in Fig. 7, 8 and 9.

Note:

- The whole surface of each electrode must be in contact with the skin to ensure a successful treatment. Press each electrode firmly onto the skin.

- Only use electrodes specified for use with NEUROTECH VITAL™.

Fig. 6


Fig. 7


Fig. 8



Fig. 9

4.3 Connecting the Garments to the Controller

1. Ensure the electrodes and garments are positioned correctly.
2. Connect the stress lead to port 1 on the controller (Fig. 10) and switch the controller on by sliding the **On/Off Switch** to "I". You will hear three rising tones and the **neurotech®** logo will appear on screen, followed by a screen showing NEUROTECH VITAL™. After that the Program screen will appear (Fig 11).
3. Connect the stress lead to the garments by inserting the **red** plug on the **right** garment to the **red** connector on the stress lead (Fig. 12) and the **blue** plug on the **left** garment into the **blue** connector on the stress lead (Fig. 13).

Fig. 10



Fig. 11

4.4 Running Program 1

1. When you switch on the controller for the first time with the stress lead connected, Program 1 (Stress Incontinence - EMS) will automatically be selected (Fig. 11). If it is not already selected, press **Screen Button 3** to access the **Menu** (⋮), then press **Select** (✓) to choose **Programs**. Use the **Up/Down Navigation Buttons** (⬆/⬇) to highlight Program 1 and then **Select** (✓). The screen will appear as shown in Fig. 11.
2. Start the treatment session by pressing the **Up Intensity Button** (⬆). You will begin to feel stimulation and the intensity bar will rise on-screen (Fig. 14). The Intensity Buttons increase/decrease the intensity in single increments. The intensity range is 0 - 100. The value appears directly above the orange intensity bar.

As the intensity increases, you should begin to feel a strong muscle contraction, first in the buttock area and then, as the intensity increases, a contraction should become apparent between the legs.

3. You should continue to increase the Intensity until you reach the target intensity level advised by your clinician. The intensity should be strong enough to give a good, strong contraction though you should not exceed your comfort limit. It may take you a few minutes to become accustomed to the strong sensation. On screen, the treatment timer counts down from 30 minutes as the program is running (Fig. 14).

During each contraction, the stimulation will ramp-up to the selected intensity level, stay at that intensity through the contraction phase and then ramp-down to the relaxation phase. The treatment runs through pre-set contraction and relaxation phases until the timer counts down to zero, at which point the stimulus is reduced to zero and the program stops.

4. To optimise the results from a NEUROTECH VITAL™ treatment session, it is recommended that you stand slightly bent over, with your feet flat on the ground and your hands resting on a stationary object, such as a table. However, if the standing position is not suitable, you may sit or lie down. A chair or bed should be nearby in case you need to do so.

Fig. 12


Fig. 13



Fig. 14

4.5 Pausing a Treatment Session

- The **Pause Button** () is located above the Navigation Buttons. If you need to pause the treatment for any reason (e.g. if the garments are not positioned correctly), simply press this button. To unpause, just press the Pause Button again. When paused, the stimulation stops and the Intensity Up Button is disabled. The intensity may, however, be decreased while the controller is paused and then, when the controller is unpaused, the stimulation will only reach the new, lower intensity.
- If you press and hold the Pause Button for 2-3 seconds, it will terminate the treatment session and the controller will revert to the program screen.

Important: Terminating a session is not the same as switching off the controller (see boxed text below).

4.6 Concluding Program 1


1. The treatment session is completed when the counter reaches zero and the controller beeps to indicate the end of a session.
2. Disconnect the stress lead from the left and right garments.
3. Switch the controller off by sliding the On/Off Switch to “0”. If you do not switch off the controller, the **Power Off Unit** screen will appear after 2 minutes (Fig. 15).
4. Unfasten and then slowly and carefully remove the garment from your left leg. First, peel the garment from the front of the leg until it is half detached and then peel the other end from the rear. Then unfasten and slowly remove the garment from your right leg in the same way.
5. Lay the garments on a flat surface and replace the covers on the **black side** of the 8 electrodes, then pack everything away until the next session.

Important: You should only fold the garments once in half, which allows the electrodes to remain flat between sessions.

Notes:

- Always switch off the controller after use to ensure that the battery does not run down.
- Always replace the covers on the black side of the electrodes after use. The electrodes can remain affixed to the garment in between treatment sessions.

Fig. 15

Switching the controller off:

At the end of a session, switch off the controller using the On/Off Switch (to “0”). If you don’t do this, after two minutes of inactivity, the controller will display a **Power Off Unit** screen (Fig. 15) and will emit three beeps every 30 seconds to remind you to switch it off. Should you press a screen button to simply clear this from the screen, it will appear again after 2 minutes. Otherwise, it will remain on screen and the controller will continue to beep every 30 seconds until it is switched off.

5.0 Program 2 - Biofeedback Assisted Exercise

Program 2 combines pelvic floor muscle exercises with EMG biofeedback. You are encouraged to reach a target EMG level by contracting your pelvic floor muscles. The EMG level of that contraction is then measured and you are encouraged to reach a percentage of that level with each subsequent contraction. The garments are not required for this program.



There are two ways you can set up program 2:

- **Option A - uses the following:**
 - 1 x Biofeedback probe
 - 1 x Med Gel
 - 1 x EMG biofeedback lead
 - 1 x Square EMG reference electrode.
- **Option B - uses the following:**
 - 2 x Small EMG biofeedback electrodes
 - 1 x EMG biofeedback lead
 - 1 x Square EMG reference electrode.

5.1 Program 2 Set-Up

SET-UP - OPTION A


1. Connect the **red** and **black** lead pins on the EMG biofeedback lead, marked "PLUS" and "MINUS", to the vaginal biofeedback probe. It does not matter which lead pin you insert into which probe socket. The lead pins must be fully inserted into the biofeedback probe.
2. Spread a light coating of the conductive, lubricating Med Gel over the body of the biofeedback probe to ease insertion.
3. Insert the biofeedback probe into the vagina, ensuring an 'East-West' position (i.e. the metal surfaces face towards each hip - Fig. 16).
4. Ensure the biofeedback probe's collar remains outside the body (Fig. 16).
5. Connect one of the square EMG reference electrodes to the **green** lead pin, marked "REFERENCE" on the EMG biofeedback lead and remove the cover from the adhesive side of the electrode.
6. Position the adhesive side of the EMG reference electrode on the lower abdomen (Fig. 17).

SET-UP - OPTION B

You may require assistance for Option B.

1. You should lie down in a comfortable position.
2. Connect the **red** and **black** lead pins on the EMG biofeedback lead, marked "PLUS" and "MINUS", to two small EMG biofeedback electrodes. It does not matter which lead pin you connect to which electrode.
3. Remove the two electrodes from their liner and position them as shown in Fig. 18. Press each electrode firmly onto the skin. The leads should run down from the lowest part of the electrodes to avoid peeling the electrodes from the skin during treatment.
4. Connect one of the square EMG reference electrodes to the **green** lead pin, marked "REFERENCE" on the EMG biofeedback lead and remove the cover from the adhesive side of the electrode.
5. Position the adhesive side of the EMG reference electrode on the lower abdomen (Fig. 17).

5.2 Running Program 2

1. Ensure the controller is adequately charged.
 2. After you have set up either Option A or Option B, switch on the controller by sliding the On/Off Switch to "I". You will hear three rising tones and the **neurotech®** logo will appear on screen, followed by a screen showing NEUROTECH VITAL™. After that the Program screen will appear (Fig 19).
 3. Insert the EMG biofeedback lead into connector port 2 (Fig. 20).
 4. If Program 2 is not already selected, press **Screen Button 3** to access the **Menu** (ⓘ), then **Select** (✓) **Programs**. Use the **Up/Down Navigation Buttons** (⊗) to highlight Program 2 and then **Select** (✓). NEUROTECH VITAL™ will now operate solely as a biofeedback unit and not as a muscle stimulator.
 5. Press either **Screen Button 2 to Start** (👉) or the **Up Intensity Button** (⬆) to begin the treatment session (please note that the intensity does not need to be increased as there is no stimulation in this program).
 6. When started, the display indicates the Rest and Work phases with Stop and Go signs (Fig. 21 and 22). Above this you can see the number of the current Trial as well as the total number of Trials to be completed. Beneath the Stop/Go signs, the number of seconds left in each Work and Rest phase of a Trial is counted down.
- Note:** A Trial is defined as one single Work/Rest phase.
7. Every time the Go sign appears on screen, you should make an effort to contract your pelvic floor muscles. The quality of the effort will be indicated on the EMG biofeedback graph. Your biofeedback level is indicated by the blue bar graph on the left hand side of the display (Fig. 22). When the Stop sign appears, you can relax again until the next Work phase.

5.3 Pausing a Treatment Session

- The **Pause Button** (⏸) is located above the Navigation Buttons. If you need to pause the treatment for any reason, simply press this button and the program will stop. To unpause, just press the Pause Button again.
 - If you press and hold the Pause Button for 2-3 seconds, it will terminate the treatment session and the controller reverts to the program screen (Fig. 19).
- Important:** Terminating a session is not the same as switching off the controller.

Fig. 19


Fig. 20



Fig. 21

Fig. 22

5.4 Concluding Program 2

1. The treatment session is completed when the full number of Trials have been completed. The controller will beep to indicate the end of a session.
2. Disconnect the EMG biofeedback lead from the controller and then remove the biofeedback probe (or EMG biofeedback electrodes) and the reference electrode. Replace the electrodes on their respective covers and if used, clean the biofeedback probe - it can be washed in mild soapy water and rinsed. Do not let the connectors get wet though and always dry the biofeedback probe fully after cleaning.
3. Switch the controller off by sliding the On/Off Switch to “ 0 ”. If you do not switch off the controller, the **Power Off Unit** screen will appear after 2 minutes (see boxed text, Page 15).

Notes:

- Always switch off the controller after use to ensure that the battery does not run down.
- Always replace the covers on the adhesive side of the electrodes after use.

6.0 Program 3 - Biofeedback and NMES

Program 3 combines volitional pelvic floor exercises and EMG biofeedback with NMES. Like Program 2, you must manually contract your pelvic floor muscles. The EMG level of that contraction is then measured and that data is used to set the target level for subsequent work phases. At the end of each Work phase, the controller begins an NMES contraction/relaxation cycle, before once again continuing through another Rest/Work Trial.

6.1 Running Program 3

1. Ensure the controller is adequately charged.
2. Switch on the NEUROTECH VITAL™ controller by sliding the On/Off Switch to "I". You will hear three rising tones and the **neurotech®** logo will appear on screen, followed by a screen showing NEUROTECH VITAL™. After that the Program screen will appear (Fig 23).
3. Ensure the garments are set up and applied correctly (see Page 14 and the Garment Quick Start Guide).
4. Ensure the EMG biofeedback (either Option A or B) is set up correctly (Page 16).
5. If Program 3 is not already selected, press **Screen Button 3** to access the **Menu** (ⓘ), then **Select** (✓) **Programs**. Use the **Up/Down Navigation Buttons** (⊗) to highlight Program 3 and then press **Select**.

6A. If your clinician has recommended that you should run a number of EMG biofeedback Trials before the NMES part of the treatment, you should begin the treatment by pressing **Start** (▶). When started, the display indicates the Rest and Work phases with Stop and Go signs (Fig. 24 and 25). Above this you can see the number of the current Trial as well as the total number of Trials to be completed. Beneath the Stop/Go signs, the number of seconds left in each Work and Rest phase of a Trial is counted down.

Every time the Go sign appears on screen, you should make an effort to contract your pelvic floor muscles. The quality of the effort will be indicated on the EMG biofeedback graph. Your biofeedback level is indicated by the blue bar graph on the left hand side of the display (Fig. 25). You should aim to reach the target EMG level indicated by a red line on the bar graph. When the Stop sign appears, you can relax again until the next Work phase.

To initiate the NMES part of the treatment, you should simply begin increasing the intensity after the prescribed number of Trials. Once you reach a level of stimulation which delivers a good, strong yet still comfortable pelvic floor contraction, you can stop increasing the intensity and the controller will run through an NMES contraction, followed by a relaxation period. The controller then moves into the EMG Rest phase and the EMG Work phase before once again beginning an NMES contraction and relaxation etc.

6B. If your clinician has prescribed that the NMES and EMG phases should run from the outset of the treatment, you should start the treatment by pressing the Intensity Up Button (⬆). You should increase the intensity to a level which delivers a good, strong yet still comfortable pelvic floor contraction and then let the program run as normal.

Fig. 23



Fig. 24

Fig. 25


6.2 Pausing a Treatment Session

- The **Pause Button** () is located above the Navigation Buttons. If you need to pause the treatment for any reason (e.g. if the garments are not positioned correctly), simply press this button. To unpause, just press the Pause Button again. When paused, the stimulation stops and the Intensity Up Button is disabled. The intensity may, however, be decreased while the controller is paused and then, when the controller is unpaused, the stimulation will only reach the new, lower intensity.
- If you press and hold the Pause Button for 2-3 seconds, it will terminate the treatment session and the controller will revert to the program screen (Fig. 26).

Important: Terminating a session is not the same as switching off the controller.

Fig. 26

6.3 Concluding Program 3

1. The treatment session is completed when the full number of Trials and NMES phases have been completed. The controller will beep to indicate the end of a session.
 2. Disconnect the EMG biofeedback lead from the controller and then remove the biofeedback probe (or EMG biofeedback electrodes) and the reference electrode. Replace the electrodes on their respective covers and if used, clean the biofeedback probe - it can be washed in mild soapy water and rinsed. Do not let the connectors get wet though and always dry the biofeedback probe fully after cleaning.
 3. Disconnect the stress lead from the left and right garments.
 4. Unfasten and slowly remove the garment from your left leg. First, peel the garment from the front of the leg until it is half detached and then peel the other end from the rear. Then unfasten and slowly remove the garment from your right leg in the same way.
 5. Lay the garments on a flat surface and replace the covers on the **black side** of the 8 electrodes, then pack everything away until the next session.
- Important:** You should only fold the garments once in half, which allows the electrodes to remain flat between sessions.
6. Switch the controller off by sliding the On/Off Switch to " O ". If you do not switch off the controller, the **Power Off Unit** screen will appear after 2 minutes (see boxed text, Page 15).

Notes:

- Always switch off the controller after use to ensure that the battery does not run down.
- Always replace the covers on the adhesive side of the EMG reference electrode, the EMG biofeedback electrodes and on the **black side** of the garment electrodes after use. The 8 stress electrodes can remain affixed to the garment in between treatment sessions.

7.0 Program 4 - biofeedback Triggered NMES (ETS)

Program 4 combines volitional pelvic floor exercises and EMG biofeedback with NMES. However, unlike Program 3, you must reach a pre-set EMG target level by contracting your pelvic floor muscles in order to trigger the NMES stimulation. Your clinician may set the target to a specific level, which you must reach in every Trial to trigger the NMES stimulation **OR** they may set the controller so that the NMES Trigger level for each Work phase is based on your ability to achieve the previous target level. A single contraction and relaxation is delivered for each successful triggering of the NMES phase of the treatment.

7.1 Running Program 4

1. Ensure the controller is adequately charged.
2. Switch on the NEUROTECH VITAL™ controller by sliding the On/Off Switch to "I". You will hear three rising tones and the **neurotech®** logo will appear on screen, followed by a screen showing NEUROTECH VITAL™. After that the Program screen will appear (Fig 27).
3. Ensure the garments are set up and applied correctly (see Page 14 and the Garment Quick Start Guide).
4. Ensure the EMG biofeedback (either Option A or B) is set up correctly (Page 16).
5. If Program 4 is not already selected, press **Screen Button 3** to access the **Menu** (⋮), then **Select** (✓) **Programs**. Use the **Up/Down Navigation Buttons** (⬆/⬇) to highlight Program 4 and then press **Select** (✓).
6. As the aim of this program is to make you work towards a target in order to trigger an NMES stimulation, you should start the treatment by pressing the Intensity Up Button (⬆). You should increase the intensity to a level which delivers a good, strong yet still comfortable pelvic floor contraction and then let the program run as normal.
7. When started, the display indicates the Rest and Work phases with Stop and Go signs (Fig. 28 and 29). Above this you can see the number of the current Trial as well as the total number of Trials to be completed. Beneath the Stop/Go signs, the number of seconds left in each Work and Rest phase of a Trial is counted down.

Every time the Go sign appears on screen, you should make an effort to contract your pelvic floor muscles. The quality of the effort will be indicated on the EMG biofeedback graph. To initiate the NMES part of the treatment, you must reach the pre-set target level, indicated by the red line on the bar graph on the left-hand side of the screen. As soon as you reach this level the stimulation will begin and the controller will run through an NMES contraction, followed by a relaxation period. The controller then moves into the EMG Rest phase and back to the next EMG Work phase etc.

Note:

- Because the NMES phase begins as soon as you reach the EMG Target level, you may not have to contract your pelvic floor muscles for the full duration of the EMG Work phase.

Fig. 27


Fig. 28

Fig. 29


7.2 Pausing a Treatment Session

- The **Pause Button** () is located above the Navigation Buttons. If you need to pause the treatment for any reason (e.g. if the garments are not positioned correctly), simply press this button. To unpause, just press the Pause Button again. When paused, the stimulation stops and the Intensity Up Button is disabled. The intensity may, however, be decreased while the controller is paused and then, when the controller is unpaused, the stimulation will only reach the new, lower intensity.
- If you press and hold the Pause Button for 2-3 seconds, it will terminate the treatment session and the controller will revert to the program screen (Fig. 30).

Important: Terminating a session is not the same as switching off the controller.

Fig. 30

7.3 Concluding Program 4

1. The treatment session is completed when the full number of Trials and NMES phases have been completed. The controller will beep to indicate the end of a session.
 2. Disconnect the EMG biofeedback lead from the controller and then remove the biofeedback probe (or EMG biofeedback electrodes) and the reference electrode. Replace the electrodes on their respective covers and if used, clean the biofeedback probe - it can be washed in mild soapy water and rinsed. Do not let the connectors get wet though and always dry the biofeedback probe fully after cleaning.
 3. Disconnect the stress lead from the left and right garments.
 4. Unfasten and slowly remove the garment from your left leg. First, peel the garment from the front of the leg until it is half detached and then peel the other end from the rear. Then unfasten and slowly remove the garment from your right leg in the same way.
 5. Lay the garments on a flat surface and replace the covers on the **black side** of the 8 electrodes, then pack everything away until the next session.
- Important:** You should only fold the garments once in half, which allows the electrodes to remain flat between sessions.
6. Switch the controller off by sliding the On/Off Switch to " O ". If you do not switch off the controller, the **Power Off Unit** screen will appear after 2 minutes (see boxed text, Page 15).

Notes:

- Always switch off the controller after use to ensure that the battery does not run down.
- Always replace the covers on the adhesive side of the EMG reference electrode, the EMG biofeedback electrodes and on the **black side** of the garment electrodes after use. The 8 stress electrodes can remain affixed to the garment in between treatment sessions.

8.0 Problem Solving Guide

SECTION 1:

The following is a list of issues which may arise when using NEUROTECH VITAL™ and how to deal with each one.

Problem	Possible Cause	Solution
The display does not come on and there is no signal	Battery discharged	Recharge battery
Error Message appears on screen	See Section 2	See Section 2
Battery symbol flashing	The battery is low	Recharge the battery
Increasing intensity causes unpleasant sensation	Electrodes incorrectly positioned	Reposition the garment to ensure the electrodes are correctly positioned
	Electrode surface is dry	Moisten the surface of the electrodes
	Possible oils, creams or lotions, pigment marks, dry marks or other factors on the skin	Wash any oils, creams or lotions, from the skin
Alarm symbol () on, controller is beeping	Liner left on one or more electrodes	Remove liners from all electrodes
	Poor electrode contact with the skin	Ensure each electrode is in full contact with skin. Press each one firmly against the skin
	Stress lead not properly connected	Ensure lead is properly connected to controller and garments
	Faulty electrode(s)	Replace electrodes
	Faulty lead	Check connections, replace if broken
	Broken wire within garment	Replace the garment

SECTION 2:

The following is a list of warnings which may potentially appear on the controller screen:

Problem	Cause	Solution
"Connect Appropriate Garment"	Garment not connected properly or incorrect lead connected to the controller	Ensure the correct lead is connected for the chosen program. If the warning persists, contact neurotech®
"Power Off Unit"	Reminder when controller left idle for 2 minutes	Switch off the controller
"Battery Low. Charge Unit"	Battery low and could possibly be damaged	Recharge battery. If the message prsists after charging, contact neurotech®
"No Data Available"	Trying to view history when none available to view	Run one or more treatments and try again
"Memory Full"	Number of sessions performed exceeds maximum memory	Visit your clinican to upload sessions to their PC or reset your controller
"Charger Error"	A fault has been detected in the charger	Ensure the charger is connected correctly. If error persists, contact neurotech®
"Internal Software Error"	A fault has been detected in the internal software of the controller	Restart the controller. If error still present after restart contact neurotech®
"Hardware Error"	Error in the electronic circuitry	Restart controller. If the error persists, contact neurotech®
"Incompatible Hardware"	There is a technical fault within the controller itself	The controller should be returned to neurotech® for repair/replacement
"Features Not Supported"	There is a technical fault within the controller itself	The controller should be returned to neurotech® for repair/replacement

9.0 System Maintenance

If you drop your NEUROTECH VITAL™ controller and there is any visible damage or any indication of loose/dislodged internal components as a result, please return the controller to **neurotech®**.

Battery Power

The battery symbol is displayed in the upper-right corner of the screen. When the battery is adequately charged, the battery symbol will be green. When the battery symbol turns orange, it indicates the battery is almost empty and needs recharging.

If the battery depletes completely, the controller will emit a double beep followed by the battery symbol flashing red. The following warning will appear on screen: **Battery Low. Charge Unit.** If this happens, you must connect the controller to the cradle, connect the battery charger and recharge the battery.

If the battery is completely depleted, you will have to reset the date and time on the controller after recharging. The **Date Screen** on the display will automatically appear the next time the controller is switched on. See Page 12 for resetting the **Date and Time**. The expected life of a full battery charge is 4-5 treatment sessions (ca. 120-150 minutes).

Replacing the Battery

After a period of time, the controller may require more frequent charging between sessions. Rechargeable batteries have a limited number of charge cycles and may eventually need to be replaced. If the battery pack needs to be replaced, you should contact **neurotech®** to arrange the return of the controller.

Disposal of Batteries and Electrodes

Used electrodes and battery packs must never be disposed of in a fire, but in accordance with your country's national or state laws governing the disposal/recycling of such items.

Disposal of device

At the end of the product lifecycle, do not throw out this product with the normal household garbage, but bring it to a collection point for the recycling of electronic equipment.

Cleaning instructions

Never let the controller get wet. However, it can be cleaned using a soft cloth, lightly dampened in soapy water.

9.1 Caring for Your NEUROTECH VITAL™ Garments

Washing Instructions

Important: Detach the controller and stress lead from the garments before cleaning. Before washing, remove all eight electrodes. We recommend that it be washed when you are replacing the electrodes.

You should wash the garments in a receptacle which allows you to lay them flat (e.g. a bath or tub). NEVER wring the garments during or after washing. This protects the internal composition of the garment from damage.

Always follow the instructions on the label when washing the garments.

Never machine wash the garments. Always handwash them in lukewarm water.

Do not use any bleach when washing your garments.

Do not dry clean your garments.

Do not tumble dry your garments. The garments should be hung up to dry. Do not dry over anything hot (e.g. a radiator or any other direct heat source). Ensure the garment is completely dry before use. Do not fold while drying.

The garments should not be ironed.

Note: Replace the garment if the material is damaged and exposes the wiring within.

Garment materials

Outer material: 100% Nylon; Inner Material: 100% Polyurethane; Fastening: 100% Nylon

9.2 Caring for Your NEUROTECH VITAL™ Electrodes

IMPORTANT ELECTRODE INFORMATION

In order to prolong the life of the electrodes, after every few treatment sessions, you should sprinkle a few drops of clean water over the adhesive surface of each electrode with your fingertips and then gently rub it over the complete surface of the electrode. This will clean away debris and help to rehydrate the adhesive surface.

Note: This should only be done with the side of the electrode which adheres to the skin. The side of the electrode which adheres to the garment should be left affixed to the garment.

- The electrodes are for single person use only.
- You should always ensure your skin is clean and free of oils, creams and other lotions before use.
- The durability and effectiveness of the electrodes depends entirely on the proper use, storage and care on the part of the user, certain skin types and the type, duration, number of sessions, intensity used and site of stimulation.
- Ensure your hands are clean before handling electrodes. You should avoid touching the skin-side of the electrodes with your fingers as much as possible when applying or removing the garments or electrodes, as this can transfer oils and skin particles to the electrodes' surface.
- After use, replace the liners on the adhesive side of the electrodes and store your unit in a cool, dry place until the next use.
- Electrodes will need to be replaced periodically as the surface picks up skin debris and becomes dry over time.

10.0 Technical Specifications

TREATMENT PROGRAM INFORMATION

Treatment Parameters - Stress Incontinence

Prog	Frequency Hz	Pulse Width μ S	Contraction (Sec.)	Relaxation (Sec.)	Treatment Time (Mins.)	Number of EMG Trials	Rest Time (Sec.)	Work Time (Sec.)
1	50	620	5	5	30			
2					N/A	10	10	5
3	50	620	5	10	N/A	10	10	5
4	50	620	5	10	N/A	30	10	5

GENERAL SPECIFICATIONS

Product Type: 203
Classification: Internally powered equipment, Class II charger, and Type BF applied parts.
Intended use: Electrical muscle stimulator / EMG biofeedback detector.
Waveform: Symmetrical bi-phasic square current waveform.

Environmental Specifications:

Operating Range: Temperature: 10 to 35°C / 55 to 95°F
 Humidity: 20 to 65% RH
Transport & Storage Range: Temperature: 0 to 55°C / 32 to 131°F
 Humidity: 10 to 90% RH

ELECTRICAL SPECIFICATIONS

Rated Outputs - Voltage/Currents for Stress Incontinence

Parameters	500 Ω	1K Ω	1.5K Ω
Output RMSA	35mA	25mA	17mA
Output RMSV	18V	25V	25V
Output Frequency	50HZ	50HZ	50HZ
DC component	0V	0V	0V
Pulse Width	620 μ S	620 μ S	620 μ S
Max Output Current (peak)	140mA	140mA	140mA

Power Supply (2504-0801) - Home Offering: FW75550/12

Input: 100-240V ~ /47-63Hz/400mA
Output: 12V DC / 1.25A
Important: Do not use any other chargers or Power Supplies

DESCRIPTION OF THE CONTROLLER'S SYMBOLS

There are a number of technical markings on your controller. These can be explained as follows:

The controller, cradle and garments are manufactured for Bio-Medical Research Ltd, Parkmore Business Park West, Galway, Ireland.

The controller requires 1 x 7.2 (6x1.2 AA) volt (NiMH) DC battery pack; DC is indicated by the symbol: 

Power (**P**): Maximum power output measured in Watts (W) into a 500Ω load.

Frequency (**F**): Number of pulses output by the controller per second, measured in Hertz (Hz).

Output (RMSA) means the maximum output root mean square current for each channel.

Output (RMSV) means the maximum output root mean square voltage for each channel.

These symbols mean "Attention, read the accompanying documents".

This symbol means type BF applied parts.

This symbol on your unit is to indicate conformity to the requirements of the Medical Device Directive (93/42/EEC). 0366 is the number of the notified body (VDE).

SN stands for 'serial number'. On a label on the underside of the controller and cradle are the serial numbers specific to those items. The letter preceding the number indicates the year of manufacture, where "Q" denotes 2011, "R" denotes 2012 etc.

The garments' batch number is represented on their packaging, by the number corresponding with the LOT symbol.

At the end of the product lifecycle, do not throw this product into the normal household garbage, but bring it to a collection point for the recycling of electronic equipment.

Some product materials can be re-used if you bring them to a recycling point. By re-using some parts or raw materials from used products you make an important contribution to the protection of the environment. Please contact your local authorities if you need more information about collection points in your area.

Waste Electrical and Electronic Equipment can have potentially harmful effects on the environment. Incorrect disposal can cause harmful toxins to build up in the air, water and soil and can be harmful to human health.

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Copyright © 2012 BMR Ltd. All rights reserved.
neurotech and Kneehab are registered trade marks of
BMR Ltd. in the EU. Registration is pending in the U.S.