

USER'S MANUAL

MODEL: 5203

PROJECTED VERTICAL FIELD ELECTROMAGNET

Date Sold:	
Serial number:	

PROPRIETARY

THIS DOCUMENT CONTAINS CONFIDENTIAL INFORMATION PROPRIETARY TO GMW ASSOCIATES. IT MUST NOT BE REPRODUCED OR DISCLOSED TO OTHERS OR USED IN ANY WAY EXCECPT FOR THE INSTALLATION, OPERATION OR MAINTENANCE OF GMW ASSOCIATES PRODUCTS.

File No: 5203_C_Proj_Vert_Field_Magnet_May_21_2015.pdf Revision Date: 21 May, 2015



SPECIFICATIONS 5203 Electromagnet Specifications	Section 1
WARNINGS [Refer to this section before operation of Electromagnet System]	Section 2
INSTALLATION Mounting Position Electrical Connections Power Supply Electrical Interlocks Water Cooling	Section 3
OPERATION	Section 4
MAINTENANCE	Section 5
EXCITATION CURVES Field versus Current Hysteresis	Section 6
TEST DATA Field Uniformity AC Operation (Not Included)	Section 7
DRAWINGS Drawing 11907-0224-0-D-S2 5203 General Assembly	Section 8
Drawing 11907-0285-0_A Cooling System	
Drawing 13907-0025-0-A 5203 & 1x BOP20-50	
Drawing 13907-0026-0-A 5203 & 2x BOP20-50	
Drawing 11907-0281-0-A Electrical Assembly 5203 / BOP20-50 (1 supply)	

Drawing 11907-0282-0-A Electrical Assembly 5203 / BOP20-50 (2 supplies)

SPECIFICATIONS

Model: 5203 Electromagnet Specifications

Projected Field: (at max current of 100A)	$B_Z = \pm 0.50T (5000G)$
	(X,Y,Z=0,0,5mm)

Projected Field Region (for B _z)	X = -6 to +6 mm
	Y = -6 to +6 mm
	Z = 0 to 12mm

Coil:

Coil resistance (20°C) $173 \text{ m}\Omega$ Max. resistance (80°C) $240 \text{ m}\Omega$ Max. continuous power*60A/14.4V (865W)Max. peak power2000W

Max. DC current60AMax. sinusoidal current85AMax. triangle wave current100A

Self Inductance (Low Field):

Self Inductance (High Field):

Approx. 3500μH at 1 Hz

Approx. 895μH at 1 Hz

(The apparent inductance increases with frequency due to eddy currents in the solid poles)

Cooling: (measured at water I/O manifold) up to 8.0 Liter/min, 4.0 bar [2.1 USG/min, 60 psid]

Thermal Interlock: Open circuit above 80° C (176° F)

Dimensions:

Drawing 11907-0224-0_D_S2
74.0 mm W x 74.0 mm D x 123.5 mm H
2.91 inch W x 2.91 inch D x 4.86 inch H

Mass: 2.5 kg (5.5 lb)

*CAUTION - The value of maximum coil power given should not be exceeded.

At this power the coils are at maximum safe temperature for continuous operation.

WARNINGS

REFER TO WARNINGS BELOW BEFORE OPERATING ELECTROMAGNET SYSTEM

1 Personnel Safety

In operation the magnet fringing field in the vicinity of the pole gap is in excess of 0.5mT (5G). This can cause malfunctioning of sensitive electronic and magnetic components. We recommend that warning signs are posted indicating that a magnetic field may be present.

2 Ferromagnetic Objects

During operation the magnet exerts magnetic attraction towards ferromagnetic objects in the near vicinity of its pole faces. Keep ferromagnetic items clear!

3 Arcing

This magnet stores energy in its field during operation. Do not disconnect any current lead while under load or the magnetic field energy will be discharged across the interruption causing arcing and possible damage to electronic circuits.

4 Coil Hot Resistance

Do not exceed the maximum coil hot resistance given in the specifications or coil overheating and possible damage may occur

5 Watches, Credit Cards, and Magnetic Disks

Do not move magnetically sensitive items into the close vicinity of the magnet pole gap. Even some anti-magnetic watches can be damaged when placed in close proximity to the pole gaps during operation. Credit cards, and magnetic disks are affected by magnetic fields as low as 0.5mT (5G). Depending on the previous operating field and the pole gap, the remanent field in the gap can be in excess of 0.5mT (5G) with the magnet power supply off or disconnected.

6 Power Supply

Refer to the power supply manufacturers manual for additional important safety information.

INSTALLATION

Mounting Position

The magnet system can be mounted in any orientation, including being completely inverted. Four M3 clearance holes are provided on the magnet transition plate for mounting the magnet.

Electrical Connections

The magnet system comes with integrated wiring for the magnet. Never connect or remove cables from the magnet system with the DC power energized otherwise damage to the magnet power supply may occur. Follow instruction below for making electrical connections.

Power Supply (Refer to drawings 13907-0025-0-A & 13907-0025-0-A)

- 1. Firstly ensure the power supply is turned off and the AC power cable is disconnected.
- 2. Connect the three wires to the output connector block on the rear of the Kepco BOP power supply as detailed below.
 - Red wire to Output
 - Black wire to Common
 - Green wire to Ground

Note: Reconnect AC power cable to power supply . The magnet system is now ready to use. Do not power up the magnet unless the cooling water is turned on and flowing at 8.0 liters/min.

Electrical Interlocks

The Model 5203 has a single thermostat, Selco part no UP62-080C. It is located between the base yoke plate and the lower cooling plate. The thermostat is normally closed, opening when the coil heatsink temperature exceeds 80° C, +/- 5° C. Connect the temperature interlock to the rear of the Kepco supply.

Water Cooling

The Model 5203 can be operated to a cooling plate temperature of 80° C. The coil thermostat will open when the cooling plate temperature exceeds approximately 80° C. If the temperature switch opens then the magnet power supply circuit breaker will trip to the off position. Clean, cool (16° C - 20° C) water at 8.0 l/min at 4.0 bar (60 psid) should be used to cool the 5203 magnet. This can be provided either by house water supply or from a recirculating chiller. The magnet is shipped with the cable junction box connected. Also 10m of hose is provided to give two 5m meter lengths for connection between the magnet and the recirculating chiller (if ordered), brass fittings on the the chiller are already provided. Note that 4 x hose clamps have also be provided for fitting these water hoses. The power cable that connects the power supply to the magnet (16907-0116-0) is also provided.

OPERATION

Electromagnet System (Kepco Power Supply operating in Current Control)

The power supply is a Kepco 20-50GL is optimized for very low ripple and noise. This model does not have a front panel interface and can be computer controlled using either digital or analogue interface. For assistance communicating with this supply please contact Yuqiang Qin at GMW Associates:

Yuqiang Qin GMW Associates 955 Industiral Road San Carlos 94070 California

Tel: (650) 802 8292 ext.18 Email: <u>yuqiang@gmw.com</u>

The Kepco 20-50MG may also be used and comes with a front user interface. This unit may also be controlled through either digital or analogue interface.


MAINTENANCE


Electrical Connections on the cable junction box and power supply should be checked annually. The electrical connections should be clean and tight. Discoloration is a sign that the connection is overheating and must be rectified before further use of the magnet.

Water Hoses should be checked regularly for water leaks. Any leaks should be rectified before further use of the magnet.


Section 6


EXCITATION CURVES








TEST DATA









































DRAWINGS

