

USER MANUAL

Version 3.0 for Microsoft® Windows

MasterPlex™ QT

Multiplex Data Analysis Software

MiraiBio

A Group of Hitachi Software

For Research Use Only

601 Gateway Blvd.
Suite 100
South San Francisco, CA 94080

TELEPHONE

1.888.615.9600 (toll free)

1.650.615.7600

FACSIMILE

1.650.615.7639

Part no. P-33020-10200

TRADEMARKS

MicroSoft® is a registered trademark of Microsoft Corporation.

COPYRIGHT

© 2001-2007 Hitachi Software Engineering America, Ltd. All Rights Reserved.

LICENSE AGREEMENT

BEFORE OPENING THIS PACKAGE, YOU SHOULD CAREFULLY READ THE FOLLOWING TERMS AND CONDITIONS. BY OPENING THIS PACKAGE YOU AGREE TO BECOME BOUND BY THE TERMS AND CONDITIONS OF THIS AGREEMENT, WHICH INCLUDES THE SOFTWARE LICENSE AND LIMITED WARRANTY. IF YOU DO NOT AGREE WITH THESE TERMS AND CONDITIONS, YOU SHOULD PROMPTLY RETURN THE PACKAGE UNOPENED TO HITACHI SOFTWARE ENGINEERING AMERICA, LTD. ("HISAL") or HISAL Distributor AND YOUR MONEY WILL BE REFUNDED.

The enclosed software is licensed, not sold, to you for use only upon the terms of this Agreement, and HISAL reserves any rights not expressly granted to you. You are responsible for the selection of the Software to achieve your intended results, and for the installation, use and results obtained from the Software. You own the media on which the Software is originally or subsequently recorded or fixed, but HISAL retains ownership of all copies of the Software itself.

LICENSE

You may:

- a. Use the Software on a single machine at any given time.
- b. Obtain limited numbers of Copy Protection Devices. Additional, Copy Protection Devices are provided only as a convenience of running the software.
- c. In no manner engineer or reverse-engineer the copy protection hardware, or whole or part of the software.
- d. Copy the software only for backup provided that you reproduce all copyright and other proprietary notices that are on the original copy of the Software provided to you. Certain Software, however, may include mechanisms to limit or inhibit copying. Such Software is marked copy protected.
- e. Transfer of the Software and all rights under this Agreement to another party together with a copy of this Agreement if the other party agrees to accept the terms and conditions of this Agreement. If you transfer the Software, you must at the same time either transfer all copies

whether in printed or machine-readable form, to the same party or destroy and copies not transferred.

RESTRICTIONS

You may not use, copy, modify, or transfer the Software, or any copy, in whole or in part, except as expressly provided for in this Agreement. Any attempt to transfer any of the rights, duties or obligations hereunder except as expressly provided for in this Agreement is void. **YOU MAY NOT RENT, LEASE, LOAN, RESELL FOR PROFIT, OR DISTRIBUTE.**

TERM

This Agreement is effective until terminated. You may terminate it at any time by destroying the Software together with all copies in any form. This Agreement will immediately and automatically terminate without notice if you fail to comply with any term or condition of this Agreement. You agree upon termination to promptly destroy the Software together with all copies in any form.

LIMITED WARRANTY

HISAL warrants, for the period of ninety (90) days from the date of delivery of the Software to you as evidenced by a copy of your receipt, that:

(1) The Software, unless modified by you, will perform the function described in the documentation provided by HISAL. Your sole remedy under the warranty is that HISAL will undertake to correct within a reasonable period of time any marked Software Error (failure of the Software to perform the functions described in the documentation).

HISAL does not warrant that the Software will meet your requirements, that operation of the Software will be uninterrupted or error-free, or that all Software Errors will be corrected.

(2) The media on which the Software is furnished will be free from defects in materials and workmanship under normal use. HISAL will, at its option, replace or refund the purchase price of the media at no charge

to you, provided you return the faulty media with proof of purchase to HISAL. HISAL will not have any responsibility to replace or refund the purchase price of the media damaged by accident, abuse or misapplication.

THE ABOVE WARRANTIES ARE EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR

PURPOSE. NO ORAL OR WRITTEN INFORMATION OR ADVICE GIVEN BY HISAL, ITS EMPLOYEES, DISTRIBUTORS, OR AGENTS SHALL INCREASE THE SCOPE OF THE ABOVE WARRANTIES OR CREATE ANY NEW WARRANTIES. SOME STATES DO NOT ALLOW THE EXCLUSION OF IMPLIED WARRANTIES, SO THE ABOVE EXCLUSION MAY NOT APPLY TO YOU. IN THAT EVENT, ANY IMPLIED WARRANTIES ARE LIMITED IN DURATION TO NINETY (90) DAYS FROM THE DATE OF DELIVERY OF THE SOFTWARE. THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS. YOU MAY HAVE OTHER RIGHTS, WHICH VARY FROM STATE TO STATE.

LIMITATIONS OF REMEDIES

HISAL's entire liability to you and your exclusive remedy shall be the replacement of the Software media or the refund of your purchase price as set forth above. If HISAL or the HISAL's distributors are unable to deliver replacement media which is free of defects in materials and workmanship, you may terminate this Agreement by returning the Software and your money will be refunded.

REGARDLESS OF WHETHER ANY REMEDY SET FORTH HEREIN FAILS ITS ESSENTIAL PURPOSE, IN NO EVENT WILL HISAL BE LIABLE TO YOU FOR ANY DAMAGES, INCLUDING ANY LOST PROFITS, LOST DATA OR OTHER INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING OUT OF THE USE OR INABILITY OF SUCH DAMAGES, OR FOR ANY CLAIM BY ANY OTHER PARTY. SOME STATES DO NOT ALLOW THE LIMITATION OR EXCLUSION OR LIABILITY FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES TO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

GOVERNMENT LICENSEE

If you are acquiring the Software on behalf of any unit or agency of the United States Government, the following provisions apply:

The Government acknowledges HISAL's representation that the Software and its documentation were developed at private expense and no part of them is in the public domain.

The Government acknowledges HISAL's representation that the Software is Restricted Computer Software as that term is defined in Clause 52.227-19 of the Federal Acquisition

LICENSE AGREEMENT

Regulations (FAR) and is commercial Computer Software as that term is defined in Subpart 227.401 of the Department of Defense Federal Acquisition Regulations supplement (DFARS) The Government agrees that:

If the Software is supplied to the Department of Defense (DOD), the Software is classified as Commercial Computer Software and the Government is acquiring only restricted rights in the Software and its documentation will be as defined in Clause 52.227-19 (c) (2) of the FAR. If the Software is supplied to any unit or agency of the United States Government other than DOD, the Governments rights in Software and its documentation

RESTRICTED RIGHTS LEGEND

Use, duplication, or disclosure by the Government is subject to restrictions as set forth in subparagraph.

(c) (1) (11) of the rights in Technical Data and computer software clause of DFARS 52.227-7013.

Hitachi Software Engineering America, Ltd.
601 Gateway Boulevard, Suite 100
South San Francisco, CA 94080

EXPORT LAW ASSURANCES

You acknowledge and agree that the Software is subject to restrictions and controls imposed by the United States Export Administration Act ("The Act") and the regulations thereunder. You agree and certify that neither the Software nor any direct product thereof is being or will be acquired, shipped, transferred or reexported, directly or indirectly, into any country prohibited by the Act and the regulations thereunder or will be used for any purpose prohibited by the same.

GENERAL

This agreement will be governed by the laws of the State of California, except for that body of law dealing with conflicts of law.

Future updates of the Software will be available for purchase by licensees for a fee provided a registration card has been received by Hitachi Software Engineering America, Ltd.

Should you have any questions concerning this Agreement, you may contact HISAL at <http://www.miraibio.com>.

You acknowledge that you have read this Agreement, understand it and agree to be bound by its terms and conditions. You further agree that it is the complete and exclusive statement of

LICENSE AGREEMENT

the agreement between us which supersedes any proposal or prior agreement, oral or written, and any other communications between us in relation to the subject matter of this Agreement.

MiraiBio

MasterPlex™ QT 3.0

Analysis software for multiplex data from
the Luminex® 100/200 system.

CONTENTS

CHAPTER 1	PAGE
Welcome	
About This Manual	1
Technical Support	2
 CHAPTER 2	
Installing MasterPlex QT	
Requirements	3
Installing MasterPlex QT	4
Installing a License	8
 CHAPTER 3	
Getting Started	
Overview of MasterPlex™ QT Analysis	9
Starting MasterPlex™ QT	10
Importing Luminex® Results	11
Using Windows Explorer to Import .csv, .xls, .lxd or Open .mlx Files	15
Viewing Data in the Plate Window	16
Thresholds	21
The Plate Navigator	23
Saving Plate Data	28
 CHAPTER 4	
Defining a Plate	
Designating Well Type and Group	30

Working With Templates	45
Setting Preferences	50
Saving a Plate	58

CHAPTER 5

Standard Curves & Analyte Concentration

Specifying Standard Data	60
Working With Diluted Unknowns	67
Generating Standard Curves & Computing Analyte Concentrations	70
Printing the Well Grid	74
Working With Standard Curves	75
Importing Standard Curves	79

CHAPTER 6

Virtual Plates

Creating a Virtual Plate	82
Working With the Virtual Analyte Filter	88
Generating a Dose-Response Curve	91
Auto Plate Merging	96

CHAPTER 7

Data Charts

Viewing a Data Chart	100
Data Chart Types	101
Working with Data Charts	112

CHAPTER 8

MasterPlex QT Reports

Generating a Report	137
Working With Reports	140
Custom Report	155

APPENDIX A

Preferences

Application Tab	159
Plate Tab	160
Calculations Tab	164

APPENDIX B

MasterPlex QT Toolbars

Main Toolbar	167
Plate Toolbar	168
Calculation Toolbar	170
Chart Toolbar	171
Report Toolbar	172

APPENDIX C

Model Equations

Four Parameter Logistic Curve	173
Five Parameter Logistic Curve	174
Heteroscedasticity	176
Weighted Nonlinear Least Square	178
Results of Weighting	180

APPENDIX D

Dose-Response Analysis

Dose-Response Analysis	182
------------------------------	-----

APPENDIX E

Relative Quantification Assay

Relative Quantification Assay	184
-------------------------------------	-----

CHAPTER

1

Welcome to the MiraiBio MasterPlex™ QT User Manual. MasterPlex QT software analyzes results files (.csv, *.xls or *.lxd) from the Luminex® 100/200 system.*

1.1

About This Manual

This manual explains how to use the MasterPlex QT 3.0 software to:

- import results files (*.csv, *.xls or *.lxd) from the Luminex system
- designate standard, unknown, control, and background wells
- generate standard curves
- compute analyte concentrations
- generate data charts and reports

What's New in MasterPlex QT 3.0

MasterPlex QT 3.0 offers new features, including the ability to:

- Merge plates using virtual plate feature so that it can analyze beyond 100 panels at one time
- Make a sample marking and groups easily and quickly using Auto-grouping feature or dragging grouping feature
- Calculate a fold change especially for being used relative gene analysis
- Normalize the data so that it can analyze between difference plates
- generate a custom reports using style sheet

Conventions Used in This Manual

This manual describes the steps required to perform the various tasks associated with the MasterPlex QT software. The manual uses a step format to explain the various tasks associated with MasterPlex QT. The symbol may follow a step instruction. It indicates the software response to the action performed by the user.

Screen Captures

Screen captures may accompany the step instructions for further illustration. The screen captures in this manual may not exactly match those displayed on your screen.

1.2

Technical Support

You can contact MiraiBio Technical support at:

Hitachi Software Engineering America, Ltd.

601 Gateway Boulevard, Suite 100

South San Francisco, CA 94080

USA

Tel: +1 (650) 615-7600

Toll Free: +1 (888) 615-9600

Fax: +1 (650) 615-7639

E-mail: support@miraibio.com

www.miraibio.com

CHAPTER 2

This chapter explains the minimum hardware and software requirements needed to install and use MasterPlex™ QT 3.0. It provides installation instructions for a computer connected to the Luminex® 100/200 system.

2.1

Requirements

For optimum performance, MasterPlex QT requires hardware and software that meet or exceed the following specifications. It is also strongly recommended that you use the Luminex XY platform.

Minimum Hardware Requirements

Platform	PC
CPU	Intel Pentium 4 1.5 GHz or equivalent, Intel Pentium4 2 GHz or better recommended
Memory (RAM)	512MB or higher for Windows 98/ME/2000/XP
Storage space (HDD)	30 MB available space for the installation
Input devices	Keyboard and mouse or any other pointing device
Video RAM	16MB or higher
Monitor resolution	XGA (1024x768 pixels or higher; 1280 x1024 recommended)
Monitor color	16-bit color (high color) or higher
CD-ROM drive	Required for CD media version. Not applicable for download version.

Software Requirements

Operating system	Microsoft Windows 98/98SE/Me/NT4.0 SP6/2000/XP (Windows 2000 or XP recommended)
------------------	--

2.2

Installing MasterPlex QT

1. Insert the MasterPlex™ QT CD-ROM in the workstation computer and double-click MasterPlex QT.exe.
⇒ The installation begins and the InstallShield Wizard appears (Figure 2.1).




Figure 2.1 InstallShield Wizard, Welcome screen

2. To continue the installation, click **Next**.
⇒ The Choose Destination Location window appears (Figure 2.2).




Figure 2.2 Install Shield Wizard, Choose Destination Location window

3. To accept the default destination folder, click **Next**.

To specify a different destination folder, click **Browse**, choose the folder, and click **Next**.

⇒ The Start Copying Files window appears (Figure 2.3).




Figure 2.3 InstallShield Wizard, Start Copying Files window

4. To copy the files to the selected directory, click **Next**.
⇒ After the installation is completed, the InstallShield Wizard Complete window appears (Figure 2.4).





Figure 2.4 InstallShield Wizard Complete window

5. Choose the **View release notes** option, and click **Finish**.

2.3

Installing a License

1. Double-click the MasterPlex™ QT icon  on the workstation desktop.
⇒ The License Information dialog box appears (Figure 2.5).




Figure 2.5 License Information dialog box

2. To view instructions on how to obtain a license (*.lic), click **Obtain Product Licenses**.
3. After you have obtained a license, click **Install New License**.
⇒ The Open dialog box appears.
4. Use the Open dialog box to locate the license (*.lic) and double-click the file.
⇒ The license is installed.

CHAPTER 3

This chapter provides a brief overview of data analysis using MasterPlex™ QT 3.0. It also explains how to start the software, import a Luminex® 100/200 results file (.csv, .xls or .lxd), and the user interface components.

3.1

Overview of MasterPlex™ QT Analysis

MasterPlex QT software analyzes results files (.csv, .xls or .lxd) from the Luminex 100/200 system. The analysis steps include:

- Import a Luminex results file (.csv, .xls or .lxd)
- Designate well types (standard, unknown, background, or control) and well groups (identifies members of a standard data set or replicate unknowns)
- Define the standard data set (enter standard concentrations and select a model equation for the standard curve)
- Associate or *link* a standard data set to an unknown group(s)
- Compute the analyte concentrations
- Save the Luminex results file in MasterPlex QT file format (.mlx).



The .mlx file includes information associated with the file (for example, well definitions and interpolated concentrations)

After the concentrations are calculated, you can:

- view the results in graphs or several different report formats
- create a *virtual plate* (a simulated microtiter plate) that contains data from user-selected actual plates (.csv, .xls, .lxd or .mlx)
- generate a Dose-Response curve and determine the Log EC50 value for user-selected data in a virtual plate

3.2

Starting MasterPlex™ QT

- On the desk top, double-click the MasterPlex QT icon . Alternatively, you can click the Windows start menu button  and select **Programs > MasterPlex QT 3.0 > MasterPlex QT 3.0**.

⇒ The MasterPlex QT user interface appears and displays the Plate Wizard and Navigator window (Figure 3.1).

For more information about the Plate Navigator window, see page 3.14.



NOTE: The Plate Wizard appears if the **Enable plate wizard at start up** option is chosen in the Application preferences or the **Display wizard at startup** option is chosen in the Plate Wizard.

The Plate Wizard guides you through the steps to import a Luminex® results file (.csv, .xls or .lxd) or create a virtual plate. For more information on virtual plates, see Chapter 6 on page 6.1.




Figure 3.1 MasterPlex™ QT user interface

3.3


Importing Luminex® Results

To begin a MasterPlex QT analysis, import a .csv, .xls or .lxd file from the Luminex 100/200 system using the Plate Wizard, toolbar, or menu bar commands.



NOTE: The Luminex default directory is named Output.

Importing Luminex Results Using the Plate Wizard

1. If the Plate Wizard is not open, click the **Plate Wizard** button 
 - ⇒ The Plate Wizard appears (Figure 3.1).
2. In the Welcome tab of the Plate Wizard, click **Next**.
 - ⇒ The Select Plate Type tab appears (Figure 3.2).




Figure 3.2 Plate Wizard

Select Plate Type tab

3. Choose the **Import a new plate** option and click **Next**.
 - ⇒ The Import File tab appears (Figure 3.3).




Figure 3.3 Plate Wizard

Import File tab

4. Enter the file path for the .csv, .xls or .lxd that you want to import.

Alternatively, click the  **Browse** button.

⇒ The Open dialog box appears (Figure 3.4).




Figure 3.4 Open dialog box

5. Navigate to the directory of the .csv, .xls or .lxd that you want to import.

6. Select one or more .csv, .xls or .lxd files and click **Open**.

To select adjacent files, press and hold the **Shift** key while you click the first and last file in the selection. To select nonadjacent files, press and hold the **Ctrl** key while you click the files of interest.

7. Click **Finish** in the Plate Wizard.

⇒ The Plate window opens and displays the results data (Figure 3.5).





Figure 3.5 Plate window displaying results data

More than one Plate window can be open at the same time.

8. To import additional Luminex results files using the Plate Wizard, click the Plate Wizard button  and repeat step 1 to step 6. Each set of results data is displayed in a separate Plate window.

Importing Luminex Results Using the Toolbar or Menu Bar

You can import a Luminex results file using the toolbar or menu bar.

1. To import a .csv, .xls or .lxd file, click the **Open** button  or select **File > Open** from the menu bar.

⇒ The Open dialog box appears (Figure 3.6).




Figure 3.6 Open dialog box

2. Navigate to the directory of the .csv, .xls or .lxd that you want to import.



NOTE: The Luminex default directory is named Output.

3. Select one or more .csv, .xls or .lxd files and click **Open**.
To select adjacent files, press and hold the **Shift** key while you click the first and last file in the selection. To select nonadjacent files, press and hold the **Ctrl** key while you click the files of interest.
⇒ The Plate window opens and displays the results data (Figure 3.7).




Figure 3.7 Navigator window, Plate window

More than one Plate window can be open at the same time.

4. To import additional Luminex results files, repeat step 1 to step 3.
Each set of results data is displayed in a separate Plate window.

3.4

Using Windows Explorer to Import .csv, .xls or Open .lxd, .mlx Files

1. Open Windows Explorer and adjust the window size so that you can view both the MasterPlexTM QT and Windows[®] Explorer application windows.
2. Use Windows Explorer to navigate to the .csv, .xls, .lxd or .mlx file(s) that you want to open.
3. Select the file(s) of interest, then click and hold the mouse button while you drag the selected file(s) to the MasterPlex QT application window (Figure 3.8).
To select adjacent files, press and hold the **Shift** key while you click the first and last file in the selection. To select nonadjacent files, press and hold the **Ctrl** key while you click the files of interest.
4. Release the mouse button.
⇒ The file(s) open in MasterPlex QT.







Figure 3.8 MasterPlex QT and Windows® Explorer application windows
Use a drag-and-drop operation to open a .csv, .xls, .lxd or .mlx file(s) in the MasterPlex QT application window


3.5

Viewing Data in the Plate Window

1. If more than one Plate window is open, select the  **Cascade**,  **Tile Horizontally**, or  **Tile Vertically** menu from the window menu bar to arrange the Plate windows for easier viewing.
2. To change the data displayed in the well grid:
 - a. Click an analyte in the Bead Set panel.
 - b. Make a selection from the data type drop-down list.

⇒ The well grid displays the data for the selected analyte.

Figure 3.9 shows the components of the Plate window. Table 3.1 lists the types of data available for display in the well grid.

3. To view background-subtracted data, click the **Subtract background** button .

⇒ The Plate window displays background-subtracted data.

For more information on background calculation options, see *Background Type* on page A.4.)

Data type drop-down list

Edit box




Figure 3.9 Plate window and Bead Set tab



Bead Set tab shows the analytes in the plate.

Plate Window Components

Well grid

microtiter plate that displays the well contents for the analyte selected from the Bead Set panel and data type selected from the data drop-down list.

Data type drop-down list

Shows the types of data available for display in the well grid. Make a selection from this drop-down list to choose the data type displayed in the well grid. Click the dropdown arrow to view the list and select a data type. Alternatively, click the   arrows to scroll forward or backward through the list. (See Table 3.1 for a description of the data types.)

Edit box

A representation of a

active (selected) well. Some data types can be edited (see Table 3.1). Enter a new value in this box to edit well data.

Bead set tab

Displays a list of the analytes (bead sets) in an assay.

Color tab

Shows the color that represents each analyte in the multi-well chart. To change a color for an analyte, right-click the color swatch, and choose or define a color in the color palette that appears.

Standard tab

Displays a list of the local or imported standard curves for the plate.

Displays the selected data type value for the

Table 3.1 Data Types in the well grid

Data Type	Description	Edit Data
Median Fluorescence Intensity (MFI)	The median fluorescence intensity measured by the Luminex [®] 100/200 system for a bead set count.	No
Count	The number of beads (per bead set) detected by the Luminex [®] 100/200 system (specified by the user in the Luminex software).	No
Concentration	The analyte concentration that is computed (interpolated or extrapolated) from the user-selected standard curve.	Yes
Dilution	The dilution factor for the well.	Yes
Standard/ Independent Value	Standard Value: Analyte concentration for a standard.	Yes
	Independent Value: The agent concentration associated with a well that is a member of a regression data set. A Dose-Response curve is generated from a regression data set.	Yes
Sample Names	User-specified name for the well.	Yes
Outlier selection	A check mark indicates the well is outlier and the well data are not included in the calculation of concentrations or a Dose-Response curve.	Yes
Well Colors	The color that represents each analyte in the bar chart.	Yes
Group Numbers	The group number of the well. Wells that belong to the same group have the same group number.	No
Replicate type	User-specified replicate type name for the well.	Yes
Standard/ Control Links	Shows the standard number that is linked to each well or well group.	Yes
Sample Average	Shows the average within the group.	No
Standard Deviation	Shows the std within the group.	No
%CV	Shows the %CV within the group.	No
Normalized Data	Shows the Normalized data within the group.	No




Figure 3.10 Plate window and Color tab

The Color tab shows the color that represents each analyte in the multi-well chart. To change a color for an analyte, right-click the color swatch, and choose or define a color in the color palette that appears.




Figure 3.11 Plate window and Standards tab

Click a tab at the bottom of the Standards tab to display a list of the local or imported standards for the plate. For more information about standards, see Chapter 5 on page 5.1.

3.6

Thresholds

You can enter an MFI, count, or concentration threshold for a plate.

The software can identify wells that contain data less than the user-specified threshold.

To set a threshold(s):

1. Click the **Preferences** button .

⇒ The Preferences dialog box appears (Figure 3.12).





Figure 3.12 Preferences dialog box, Plate tab

2. Click the Plate tab.
3. Enter the concentration, MFI, or bead count threshold.
4. Click **Apply** when you are finished.
5. To return the plate preferences to the factory set defaults, click **Default**.


To identify wells in the Plate window that contain data less than threshold:

1. Make a selection from the Bead Set panel (Figure 3.13).
2. Select the data type (MFI, count, or concentration) from the data-type


drop-down list.

3. To identify wells with MFI data less than threshold, select the  show < MFI Limit from plate menu.

⇒ A red border is placed around the well (Figure 3.13).

To identify wells with bead count data less than threshold, select the  show < Count Limit menu from plate menu.

⇒ A red border is placed around the well.

To identify wells with concentrations less than threshold, select the  Show < CON Limit menu from plate menu

⇒ A red border is placed around the well.







Figure 3.13 Well grid in the Plate window

Red border identifies wells with a MFI value less than threshold for the selected analyte.




NOTE: The threshold menus (,  or ) are mutually exclusive. Only one is active at a time.

3.7

The Plate Navigator

At startup, MasterPlex™ QT displays the Plate Navigator (Figure 3.14). The plate navigator workspace shows a tree of the windows that are open during a session.

The Plate Navigator organizes the open windows that are associated with a plate. It provides a convenient way to change the view in the main display area.

- To view a particular window, click the item in the Plate Navigator window (Figure 3.15).
- To show or hide the Plate Navigator, select the  Navigator menu from Show > show/hide > Navigator.




Figure 3.14 Plate Navigator window at start up

If no windows are open in the main display area, the Plate Navigator window is empty.




Figure 3.15 Plate Navigator window

If multiple windows are open, click an item in the Plate Navigator to view it in the main display area.




Figure 3.16 Plate Navigator

Shows the windows that are open in the main display area.

Table 3.2 shows the types of windows available in the software.

Table 3.2 MasterPlex™ QT windows

Window Name	Displays...	See...
Plate	a .csv, .xls, .lxd or .mlx file	Figure 3.17
Standard Curves	the Standard Curves window for the current Plate window	Figure 3.18
Data Chart	the current Data Chart window for the selected plate	Figure 3.19
Report	a report that was generated for the current plate	Figure 3.20




Figure 3.17 Plate window




Figure 3.18 Standard Curves window



Figure 3.19 Data Chart window

MasterPlex QT Plate Format Report			
Report Date: 4/27/2007	Run Date: 3/20/01		
Report Time: 9:41:02 AM	Run Time: 12:45:55 PM		
Data File: IL5.csv	Hardware Serial No. : LX10000280002		
Plate Name: Untitled Plate 2	Operator: AC		
MasterPlex QT Version: 3.0.1.166	Analyst:		

Data Type : MFI

Analyte : 034 Mouse IL-5								Background : 0.00				
	1	2	3	4	5	6	7	8	9	10	11	12
A	20.00	367.00	806.50	341.00								
B	1308.00	429.00	480.50	392.00								
C	1084.50	82.00	88.00	78.00								
D	525.00	88.50	82.00	96.50								
E	201.00	37.50	38.00	37.00								
F	66.50	32.00	32.00	34.00								
G	39.00	18.00	18.00	23.00								
H	26.50	19.00	17.50	27.00								

Analyte : 038 Mouse IL-10								Background : 0.00				
	1	2	3	4	5	6	7	8	9	10	11	12
A	13.00	710.50	1086.00	833.00								
B	5476.50	758.50	1028.00	1001.00								
C	3242.00	105.50	165.00	124.00								
D	1332.00	110.00	161.00	126.00								
E	402.00	32.50	40.00	33.00								
F	107.00	30.00	40.00	35.00								
G	41.00	13.00	12.00	12.00								
H	21.00	12.00	12.00	14.00								

Analyte : 054 Mouse IL-2						Background : 0.00						
	1	2	3	4	5	6	7	8	9	10	11	12
A	28.00	1550.50	1729.50	1223.00								
B	7855.00	1617.50	1691.00	1298.00								
C	5706.50	106.00	100.00	70.00								
D	2445.00	111.00	95.00	74.00								
E	336.00	42.50	38.00	45.50								
F	57.00	50.00	40.50	47.00								
G	34.00	43.00	32.00	39.00								
H	33.00	46.00	37.00	36.50								

Analyte : 073 Mouse IFN gamma								Background : 0.00				
	1	2	3	4	5	6	7	8	9	10	11	12
A	102.00	1185.00	1674.00	847.00								
B	3141.50	1357.00	1516.00	808.00								
C	2515.00	292.00	293.50	172.50								


Figure 3.20 Report viewer, Plate format report

3.8

Saving Plate Data

After you import a Luminex results file (.csv, .xls or .lxd), the data can be saved to a MasterPlex™ QT file format (.mlx). The .mlx file includes all data associated with a plate such as well definitions and computed (interpolated or extrapolated) concentrations.

To save results data (.csv, .xls or .lxd) to a MasterPlex file (.mlx):

1. Click the **Save** button . Alternatively, select **File > Save** from the main menu.
⇒ The Save As dialog box appears (Figure 3.21).





Figure 3.21 Save As dialog box

2. Confirm the default directory where the file will be saved or choose another directory.
3. Enter a file name and click **Save**.

Opening a MasterPlex File (.mlx)

1. Click the **Open** button . Alternatively, select **File > Open Plate** from the main menu.
⇒ The Open dialog box appears (Figure 3.21).




Figure 3.22 Open dialog box

2. Confirm the default directory or choose another directory.
3. Select a file name (.mlx) and click **Open**.
⇒ A Plate window opens and displays the results data (Figure 3.23).




Figure 3.23 Plate window

CHAPTER 4

After you import a Luminex results file (.csv, .xls, .lxd), your analysis begins by defining a plate. This chapter explains how to define and save a plate. The steps to define a plate include:

- **Designate well type** to identify the standard, unknown, background, and control wells.
- **Designate well groups** to identify replicate unknowns, members of a standard data set, or members of a regression data set.
- **Create a standard data set(s)** by entering the concentration for each well in the standard data set and selecting a model equation for the standard curve. A plate can have more than one standard data set.
- **Link each well group to a standard data set** to specify the standard that is used to compute (interpolate or extrapolate) the analyte concentrations.

The plate definition can be saved as a template that can be applied to other plates. The Template Manager helps you manage your templates. For more information on templates, see *Working With Templates* on page 4.9.

4.1

Designating Well Type and Group

Selecting Wells

To select a well in the Plate window, click the well in the well grid. There are three ways to select multiple wells:

- To select adjacent wells (Figure 4.1), press and hold the mouse key while you drag the mouse cursor over the wells that you want to select. Click and release the mouse button to select the highlighted wells.
- To select adjacent wells, press and hold the **Shift** key while you click the first and last well in the selection.
- To select nonadjacent wells (Figure 4.2), press and hold the **Ctrl** while you click the wells.

	1	2	3	4
A	0.00	347.00	486.50	321.00
B	1288.00	409.00	460.50	372.00
C	1064.50	62.00	68.00	58.00
D	505.00	68.50	62.00	76.50
E	181.00	17.50	18.00	17.00
F	46.50	12.00	12.00	14.00
G	19.00	-2.00	-2.00	3.00
H	6.50	-1.00	-2.50	7.00

Figure 4.1 Well grid

To select adjacent wells, press and hold the **Shift** key while you click the first and last well in the selection. Alternatively, press and hold the mouse key while you drag the mouse over the wells of interest.

	1	2	3	4
A	0.00	347.00	486.50	321.00
B	1288.00	409.00	460.50	372.00
C	1064.50	62.00	68.00	58.00
D	505.00	68.50	62.00	76.50
E	181.00	17.50	18.00	17.00
F	46.50	12.00	12.00	14.00
G	19.00	-2.00	-2.00	3.00
H	6.50	-1.00	-2.50	7.00

Figure 4.2 Well grid

To select nonadjacent wells, press and hold the **Ctrl** key while you click the wells of interest.

Designating Well Type

Table 4.1 shows the types of wells that are available.

1. Select the well(s) that you want to define.
2. To define (or *mark*) the well(s), right-click the selection and choose **Mark Wells** and a well type from the pop-up menu that appears Figure 4.3. You can also use toolbar buttons or menu bar commands to define the wells (Table 4.1).
⇒ The well type is applied to the selected well(s).



NOTE: If the **Automatic Well Grouping** option is selected in the Preferences dialog box, the wells are automatically grouped when you set the well type. For more information on this option, see *Application Tab* on page A.1.












Figure 4.3 Well grid pop-up menu

Right click a well to display the pop-up menu

Table 4.1 Button and context menu to define wells


Well Type	Button	Context menu on the plate window
Background Wells that contain no analytes.		Mark Wells > Absolute Quantification > Background
Standard Wells that contains analyte of known concentration.		Mark Wells > Absolute Quantification > Standard
Unknown Wells that contains analytes of unknown concentration.		Mark Wells > Absolute Quantification > Unknown
Control Wells that contain analytes that function as controls for a particular assay design.		Mark Wells > Absolute Quantification > Control
Regression Wells of a virtual plate that are members of a regression data set. (See <i>Generating a Dose-Response Curve</i> on page 6.11 for more information on regression analysis data sets.)		Mark Wells > Absolute Quantification > Regression
Assay Background Wells that contains analytes of assay background concentration.		Mark Wells > Relative Quantification > Assay background
Treatment Sample Wells that contain analytes that function as treatment (patient sample).		Mark Wells > Relative Quantification > Treatment Sample
Control Sample Wells that contain analytes that function as controls (known as reference sample).		Mark Wells > Relative Quantification > Control Sample



NOTE: To show or hide the color-coded well types, click the **Show well type**



button .

3. To unmark a well, right-click the well and select **Un-Mark Wells** from the pop-up menu. Alternatively, select the well(s) and click the **Unmark wells** button  or select **Template > Unmark well** from the menu bar.



If a well belongs to a group, unmarking the well also removes the well from the group.

4. Repeat step 1 and step 2 to mark and group other well(s).


Designating Well Groups

After you have defined the wells, the wells must be organized into *groups* so that the software can identify:

- Replicate unknowns
- A standard data set
- A regression analysis data set (for a virtual plate only) (For more information on virtual plates, see Chapter 6 on page 6.1.)

MasterPlex™ QT automatically places all background wells into one group. You can define one or more groups of control wells per plate.



NOTE: A group can include nonadjacent wells. A plate can have more than one group of standards or unknowns. A virtual plate can have more than one set of regression data. To show or hide the group borders, click the **Show well grouping** button .



NOTE: Grouping does not set replicate unknowns. Standard, unknowns, control and background groups can be treated as replicates to obtain statistics such as mean, standard deviation and CV% in the final report.

Automatic Well Grouping

If the **Automatic Well Grouping** option is chosen in the Preferences dialog

box (Figure 4.4), the software automatically groups wells when you set the well type. The automatic well grouping is the factory set default.





Figure 4.4 Preferences dialog box

To display the Preference dialog box, click the button  or select **File > Preferences** from the menu bar.

Manually Grouping Wells


You can only place the same type of wells in a group.

1. Select the wells that you want to group.



NOTE: The grouping function is only available if you have selected wells that have been defined.

2. Right-click the selection, and select **Group Wells** from the pop-up menu that appears. Alternatively, click the selection, and click the **Group Wells**

- button  or select **Template > Group Wells** from the menu bar.
- ⇒ A border appears around the grouped wells (Figure 4.5).
3. To display the well grouping (Figure 4.5), click a well.




Figure 4.5 Well groups


A red border identifies a group in the well grid.

Ungrouping Wells

1. Right-click the group.



NOTE: The ungrouping function is only available if you have selected grouped wells.

2. Select **Ungroup wells** from the pop-up menu that appears. Alternatively, click the **Ungroup wells** button  or select **Template > Ungroup wells** from the menu bar.
- ⇒ The grouping is removed.

Grouping Wells by Pattern

The purpose of pattern grouping is to provide users another way to easily and quickly group wells. Pattern here means two things: the group type (e.g.,

standard, unknown...) and the dimensions of the group (i.e., rows and columns). This function acts similarly to the Resizing feature of Microsoft Excel. It is especially useful when the plate has many groups/replicates which follow similar group patterns.

1. Define the group pattern by selecting a group of wells, and marking and grouping them together. We will group other wells into this pattern.
2. Select all wells of the pattern group(Figure 4.6).

	1	2	3	4	5	6
A	20.00	367.00	506.50	341.00		
B	1308.00	429.00	480.50	392.00		
C	1084.50	82.00	88.00	78.00		
D	525.00	88.50	82.00	96.50		
E	201.00	37.50	38.00	37.00		
F	66.50	32.00	32.00	34.00		
G	39.00	18.00	18.00	23.00		
H	26.50	19.00	17.50	27.00		

Figure 4.6 Well groups

2. Point mouse cursor to the bottom-right corner of the selection. When you see the cursor turns into a black cross, push down the left button of mouse and start to drag. During dragging, you will see in real-time that new wells are selected and grouped into the pattern, as indicated by a red-line border (Figure 4.7).


	1	2	3	4	5
A	20.00	367.00	506.50	341.00	
B	1308.00	429.00	480.50	392.00	
C	1084.50	82.00	88.00	78.00	
D	525.00	88.50	82.00	96.50	
E	201.00	37.50	38.00	37.00	
F	66.50	32.00	32.00	34.00	
G	39.00	18.00	18.00	23.00	
H	26.50	19.00	17.50	27.00	

Figure 4.7 Well groups

- Once you are satisfied with the selection, just release the mouse button. The software will automatically finish the grouping(Figure 4.8).

	1	2	3	4	5
A	20.00	367.00	506.50	341.00	
B	1308.00	429.00	480.50	392.00	
C	1084.50	82.00	88.00	78.00	
D	525.00	88.50	82.00	96.50	
E	201.00	37.50	38.00	37.00	
F	66.50	32.00	32.00	34.00	
G	39.00	18.00	18.00	23.00	
H	26.50	19.00	17.50	27.00	

Figure 4.8 Well groups

 **NOTE:** When starting drag, you can move the cursor either downwards or rightwards, which results in different ways to select wells. To switch between the two modes, just drag the cursor back into the pattern group, and then drag out in either direction. So, it is determined by your first move direction when you are dragging the cursor out of the pattern group.

	1	2	3	4	5
A	177.00	2966.00	162.50	609.50	11.00
B	93.00	3083.00	190.00	788.00	12.50
C	247.00	9867.00	188.00	2203.50	27.00
D	212.00	10478.00	206.50	2440.00	26.00
E	356.50	16523.00	199.50	4631.00	80.50
F	352.00	13902.00	207.00	4344.00	79.00
G	848.00	202.00	274.00	190.00	341.00
H	847.50	221.00	320.00	201.00	344.00

	1	2	3	4	5
A	177.00	2966.00	162.50	609.50	11.00
B	93.00	3083.00	190.00	788.00	12.50
C	247.00	9867.00	188.00	2203.50	27.00
D	212.00	10478.00	206.50	2440.00	26.00
E	356.50	16523.00	199.50	4631.00	80.50
F	352.00	13902.00	207.00	4344.00	79.00
G	848.00	202.00	274.00	190.00	341.00
H	847.50	221.00	320.00	201.00	344.00


Figure 4.9 Well groups

Dragging downwards as the first move (above) vs. dragging rightwards as the first move (below)

Grouping Wells by Sample Name

The purpose of grouping wells by sample name is to provide users an alternative way to easily and quickly group wells based on wells' sample names. When using this tool, users will be asked, in three steps, to provide information on what string segment of sample name are related to sample type, and what segment is related to replicate group. The segment of sample type provides information about well's group type (e.g., standard, unknown, background ...); the segment of replicate group provides information about replicates, assuming wells in one replicate group have an identical string segment in their names. Based on the information, this tool will automatically group wells together into different types. The auto grouping tool is especially useful when the data file is

generated strictly following the naming rules for its analytes.

1. Start wizard by click toolbar button  or right-click dropdown button “Auto Group”.
2. Step 1: specify delimiters. Users choose or set delimiters to segment sample names. The segmented results are displayed at the same time (Figure 4.10).




Figure 4.10 Specifying delimiters

3. Step 2: specify sample type / replicate group. Users choose which segmented column represents sample type (e.g., standard, unknown, background ...), and which column represents replicate group (Figure 4.11).




Figure 4.11 Specifying sample types

4. Step 3: specify group type. Based on users' choice on sample type in step 2, they need to choose the segmented string for each group type (Figure 4.12).





Figure 4.12 Specifying group types

5. Finish.



	1	2	3	4	5	6	7
A	bkg	std	ctrl 1	unk 3	unk 7	unk 11	unk 15
B	bkg	std	ctrl 1	unk 3	unk 7	unk 11	unk 15
C	std	std	ctrl 2	unk 4	unk 8	unk 12	unk 16
D	std	std	ctrl 2	unk 4	unk 8	unk 12	unk 16
E	std	std	unk 1	unk 5	unk 9	unk 13	unk 17
F	std	std	unk 1	unk 5	unk 9	unk 13	unk 17
G	std	std	unk 2	unk 6	unk 10	unk 14	unk 18
H	std	std	unk 2	unk 6	unk 10	unk 14	unk 18


Figure 4.13 Autogrouping result



NOTE: User's choices in the wizard will be saved into a file each time. Next time, user can simply choose "Use the last settings" to further speed up the grouping procedure.

Setting Standard Concentrations

After you define and group the standard wells, use the autofill feature to help you automatically enter the standard concentrations. For more information on specifying standard data, see Chapter 5 on page 5.1.

1. Right-click the standard data set and select **Auto Fill** from the popup menu. Alternatively, click the **Auto Fill** button  or select **Template > Auto Fill** from the menu bar.

⇒ The Auto Fill dialog box appears (Figure 4.14).






Figure 4.14 Auto Fill dialog box

2. Make a selection from the Analyte drop-down list.
3. Enter the starting concentration for the standard data set.
4. Enter the dilution factor.
5. Make a selection from the concentration unit drop-down list
6. To select a dilution direction for the standard data set, click a dilution direction arrow.
 - ⇒ The gradient map shows the location and direction of the dilution gradient(s) (Figure 4.15).



This gradient map specifies a separate dilution gradient in each column of the standard data set. The starting concentration is at the top of a column.



This gradient map specifies one dilution gradient per standard data set. The starting concentration is at the upper left well and the end concentration is at the lower right well. Click an arrow to choose a dilution direction.

Figure 4.15 Example dilution gradient maps

Click a dilution direction arrow to choose the dilution gradient configuration for a standard data set

7. To specify the same starting concentration, dilution factor, and concentration units for all analytes in the standard data set, choose the **Fill in for all bead sets** option. To specify a different starting concentration, dilution factor, or concentration unit for a different analyte, repeat step 2 through step 4.
8. Click **Fill** when finished entering the concentration, dilution, and dilution direction for all analytes in the standard data set.
 - ⇒ A confirmation message appears (Figure 4.16).




Figure 4.16 Message box confirming autofill

Linking a Standard Data Set


Background, control, and unknown wells must be associated with or *linked* to the standard data set that will be used to calculate concentrations. By default, the first standard that you define will be linked to the background, control, and unknown well groups.

If there is more than one standard on a plate, you can link a user-selected standard to a user-selected well group(s).

1. To link a well group to a standard data set, press and hold the **Ctrl** key while you click the group and the standard data set that you want to link.



NOTE: A standard data set can be linked to multiple groups of the same well type, but each group can have only one standard.

2. Right-click the selection and choose **Link Standard** from the popup menu that appears. Alternatively, click the **Link** button  or select **Template > Link Standard** from the menu bar.

4.2

Working With Templates

A plate definition includes:

- well types and well groups
- standards (including standard concentrations, associated model equation, and concentration units)

-
- links between the standard(s) and well groups
 - data calculated for the plate (for example, analyte concentrations, standard data curves, or regression analysis curves (for virtual plates only)
 - data manually entered in the plate (for example, dilution factors or independent data values)
 - plate preferences


You can save the plate definition as a template. You can apply a template to an active plate. Templates may also be exported, imported, or deleted.



NOTE: If a plate window is not open, you can still delete, import, or export templates; however, you cannot load, save, or overwrite a template.

Opening the Template Manager

The Template Manager is a tool that helps you manage your templates.

1. Click the **Template Manager** button . Alternatively select **Plate > Template Manager** from the menu bar.
⇒ The Template Manager appears (Figure 4.17).
2. Click a template in the Available Templates list to view information about the template.





Figure 4.17 Template Manager shows available templates

Click a template to view information about the template.

Saving a Template

You can save the current plate definition to a template.

1. After you have finished defining a plate, open the Template Manager and click the **Save** button .

⇒ The Template Name box appears (Figure 4.18).





Figure 4.18 Template Name box

2. Enter a name for the template and click **OK**.

⇒ The new template is added to the Available Template list.


Loading a Template

You can apply or *load* a saved template to the current plate.

1. In the Template Manager, select the template that you want to apply to the plate.
2. Click the **Load** button .
 - ⇒ The template is applied and the well grid shows the new well attributes (well type, well group, and links to standard data sets).

Overwriting a Template

You can overwrite an existing template with the current plate definition.

1. In the Template Manager, select the template that you want to overwrite
2. Click the **Overwrite** button .
 - ⇒ A confirmation box appears (Figure 4.19).





Figure 4.19 Confirmation box

3. Click **OK** to overwrite the selected template with the current plate definition.

Exporting a Template

You can export a template to a user-specified location.

1. In the Template Manager, click the template that you want to export.
2. Click the **Export** button .
 - ⇒ The Save As dialog box appears (Figure 4.20).




Figure 4.20 Save As dialog box

3. Choose the directory for the template that you want to export.
4. Enter a name for the template (*.mxt).



NOTE: A template must have a .mxt file extension. Changing the extension will render the exported template unusable.

Importing a Template

You can import a template (.mxt) from a user-specified location. Templates from MasterPlex QT 1.0 can also be imported.

1. In the Template Manager, click **Import** button.
 - ⇒ The Open dialog box appears (Figure 4.21).





Figure 4.21 Open dialog box

2. Choose the directory with the template that you want to import.
3. Select the template and click **Open**.
 - ⇒ The template name is added to the Template Manager.

Deleting a Template

You can delete a template (.mxt) from the system.

1. In the Template Manager, click the template that you want to delete.
2. Click **Delete** button .
 - ⇒ A confirmation box appears (Figure 4.22).




Figure 4.22 Confirmation box

3. Click **OK** to delete the template.
 - ⇒ The template is removed from the Template Manager.



WARNING: This permanently removes the template from the system.

4.3

Setting Preferences

Preferences are user-modifiable software settings. They are displayed in the Preferences dialog box.

- To open the Preferences dialog box (Figure 4.23), click the **Preferences** button . Alternatively, select **File > Preferences** from the menu bar.




Figure 4.23 Preferences dialog box

Application tab default settings

Application Preferences

The application preferences (Figure 4.23) include:

Default Virtual Plate Dimensions

The default number of rows and columns for a virtual plate (displayed by the plate wizard during virtual plate setup).

Automatic Well Grouping

Choose this option to automatically group wells after they are defined.

Enable plate wizard at startup

Choose this option to display the plate wizard at startup.

1. To set the row and column dimensions of the well grid for a virtual plate displayed by the plate wizard, enter the number of rows and columns.
2. If you do not want to automatically group wells after they are defined, remove the check mark from the **Automatic Well Grouping** option.
3. If you do not want to display the plate wizard at program start up, remove the check mark from the **Enable plate wizard at start up** option.
4. Click **Apply** when you are finished.
5. To return the application preferences to the factory set defaults, click **Default**.

Plate Tab

The Plate preferences (Figure 4.24) specify:

- how to compute background and when to subtract background
- threshold values for concentration, MFI, and bead count
- the plate and analyst name




Figure 4.24 Preferences dialog box

Plate tab default settings



NOTE: If a plate (.csv, .xls, .lxd or .mlx) is not open, the Preferences dialog box does not display the Plate tab.

Background Options

You can specify whether you want to consider background data in the calculation of analyte concentrations. There are two methods of computing background-subtracted analyte concentrations: All Calculations or Regressions Only.

Subtract Background

Choose this option if you want to compute background-subtracted analyte concentrations. If this option is not chosen, the background MFI value is not considered during calculation of the analyte concentrations.

All Calculations

This method of computing background-subtracted analyte concentrations subtracts the background MFI from each member of the standard data set, then fits the standard curve. The method subtracts background MFI from unknown MFI, then interpolates or extrapolates the unknown analyte concentration from the standard curve.

Regressions Only

This method of computing background-subtracted analyte concentrations subtracts the background MFI from each member of the standard data set, then fits the standard curve. The method does not subtract the background MFI from the unknown MFI before interpolating or extrapolating the unknown analyte concentration.



NOTE: The **All Calculations** method is recommended. The **Regressions Only** method provides backward compatibility with data generated in MasterPlex™ QT 1.0/2.0.

Background Type

If there are three or more background wells in the assay, choose one of the following methods for computing background MFI.

Average Background (Bkg) MFI = (Bkg MFI₁ + Bkg MFI₂ +... Bkg MFI_n)/n

Peak Value where n = the number of background wells in the plate
Highest background MFI value.

Lowest Value Lowest background MFI value.

Plate Info

Plate Name Displays the name assigned to the result file in the
Luminex® 100/200 software. To edit the plate name, enter
a new name.

Analyst Name Displays the analyst name entered in the Luminex®
100/200 software. To edit the analyst name, enter a new
name.

Thresholds

You can enter an MFI, count, or concentration threshold for a plate. The software automatically marks wells that contain data less than the user specified threshold with a red border (Figure 4.18).

To set a threshold(s):

1. Enter the MFI, count, or concentration threshold in the Preferences dialog box (Figure 4.25).







Figure 4.25 Preferences dialog box, Plate tab

2. Click **Apply** when you are finished.
3. To return the plate preferences to the factory set defaults, click **Default**.


To identify the wells that contain data less than threshold:

1. Choose an analyte in the Bead Set panel.
2. Select the data type (MFI, count, or concentration) from the datatype drop-down list.
3. To identify wells with:
 - a. MFI less than threshold, click the  button.
 - b. Bead count data less than threshold, click the  button.
 - c. Concentration less than threshold, click the  button.

A red border marks wells that contain data less than threshold for the selected

analyte (Figure 4.26).

 **NOTE:** The threshold buttons () are mutually exclusive. Only one is active at a time.



	1	2	3	4	5
A	9.30	193.57	243.01	183.99	
B	998.69	215.85	233.95	202.64	
C	500.28	66.15	69.94	63.55	
D	249.45	70.26	66.15	75.12	
E	127.80	31.90	32.38	31.41	
F	55.70	26.17	26.17	28.33	
G	33.38	4.41	4.41	14.58	
H	19.57	7.11	2.59	20.23	

Figure 4.26 Well grid

Red border identifies wells with concentration less than user-specified threshold for the selected analyte

Calculations Preferences

The Calculations preferences (Figure 4.27) specify how to:

- fit the standard curve when there are replicate standard data sets
- display concentrations for unknowns that were diluted prior to the assay (the diluted concentration or the original undiluted concentration can be displayed)




Figure 4.27 Preferences dialog box, Calculations tab

Standard Points

If a plate contains replicate standard data sets, there are two ways to fit a standard curve.

Average standards

The replicate standard data points are averaged and the standard curve is fitted to the single set of averaged data. Choose this option if the experimental errors are dependent on each another.

Individual points	The replicate standard data points are not averaged and the standard curve is fitted to all of the data points. For example, if there are three replicates of eight standard wells, the standard curve is fitted using all 24 data points. Choose this option if the experimental errors are independent of each other. This option is set by default. In most cases, if you're not sure about the nature of the distribution of the experimental errors, try this option first.
--------------------------	--

Dilution for Unknowns

Samples can be diluted prior to the assay and analysis. After MasterPlex QT interpolates the diluted unknown analyte concentrations from the standard curve, it can compute and display the original, undiluted concentration in the Plate window.

$$\text{Original concentration} = \text{Diluted concentration} * \text{Dilution Factor}$$

The dilutions factors are manually entered in the Plate window. (For more information see *Working With Diluted Unknowns* on page 5.8.

4. Click **Apply** when you are finished.


To return the plate preferences to the factory set defaults, click **Default**.

4.4

Saving a Plate

The software saves the imported results (.csv, .xls or .lxd) in MasterPlex™ QT native file format (.mlx). The .mlx file also includes the:

- well and group definitions or the associated template
- well and bead set (analyte) color
- data calculated for the plate (for example, analyte concentrations, standard data curves, or regression analysis curves (for virtual plates only)
- data manually entered in the plate (for example, dilution factors or independent data values)
- plate preferences

1. To save a plate, click the **Save** button .
⇒ The Save As dialog box appears (Figure 4.28).




Figure 4.28 Save As dialog box

Save the current plate data in native .mlx format.

2. Enter a name for the file.
3. Click **Save**.

CHAPTER 5

This chapter explains how to generate standard curves and compute (interpolate or extrapolate) analyte concentrations from the standard curves.

5.1

Specifying Standard Data

Each well in a standard data set represents an x,y data point. The MFI value is plotted on the y-axis and the concentration is plotted on the x-axis.

MasterPlex™ QT uses regression analysis to fit a user-specified model equation to the standard data set and generate a standard curve.



NOTE: The standard curve may not pass through each point in the standard data set.

The software computes the R^2 value ($0 \leq R^2 \leq 1$) for the model equation. R^2 measures the goodness of fit of the model equation to the standard data set (where $R^2 = 1$ is the probability that the model predicts the data perfectly).

The steps to create a standard curve include:

1. Mark the standard wells.
2. Group the wells in a standard data set.
3. Link the standard data set to the unknown well group(s) of interest.
(The analyte concentrations are interpolated from the standard curve that is linked to the unknown well group.)
4. Enter the standard concentrations.
5. Select a model equation for the standard data set.
6. Calculate the standard curves.



NOTE: A plate can have more than one standard data set. The standard data sets may have different concentrations or model equations.

Entering Standard Concentrations

Enter the standard concentrations after you mark the standard wells, group them into a standard data set, and link the standard data set to a group(s) of unknown

wells. (For more information on setting well types and groups, see *Designating Well Type and Group* on page 4.1.)

There are two ways to enter standard concentrations:

- Use the autofill feature to automatically enter the analyte concentrations
- Manually enter the analyte concentrations

Using the Autofill Feature

The autofill feature enters the analyte concentrations for selected standard wells based on the user-specified starting concentration and dilution factor.

1. Open the results of interest (.csv, .xls, .lxd or .mlx) and select the standard data set (Figure 5.1).





Figure 5.1 Standard wells selected in the well grid

2. Right-click a standard well and select **Auto Fill** from the pop-up menu that appears. Alternatively, click the **Auto Fill** button 
 - ⇒ The Auto Fill dialog box appears (Figure 5.2).

The Auto Fill dialog box is used for setting up dilution parameters. It includes a 'Dilution parameters' section with fields for Analyte (034 Mouse IL-5), Starting concentration (10000.000), Dilution Factor (2.000), and Concentration Unit (pg/mL). There is also a checkbox for 'Fill in for all bead'. The 'Details' section shows a table with 4 rows (A, B, C, D) and 3 columns (1, 2, 3), all containing downward arrows. A 'Dilution direction' section contains a 3x3 grid of arrows (up, down, left, right, and diagonals) in green, blue, and red. At the bottom are 'Fill' and 'Cancel/Exit' buttons.

	1	2	3
A	↓	↓	↓
B	↓	↓	↓
C	↓	↓	↓
D	↓	↓	↓

Figure 5.2 Auto Fill dialog box

3. Make a selection from the Analyte drop-down list.
4. Enter the starting concentration for the standard data set.
5. Enter the dilution factor.
6. Make a selection from the concentration unit drop-down list
7. To select a dilution direction for the standard data set, click a dilution direction arrow.
 - ⇒ The gradient map shows the location and direction of the dilution gradient(s) (Figure 5.3).





Figure 5.3 Example dilution gradient maps

Click a dilution direction arrow to choose the dilution gradient configuration for a standard data set

8. If you want to specify the same starting concentration, dilution factor, and concentration units for all analytes in the standard data set, choose the **Fill in for all bead sets** option. If you want to specify a different starting concentration, dilution factor, or concentration unit for a different analyte, repeat step 2 through step 4.
9. Click **Fill** when finished entering the concentration, dilution, and dilution direction for all analytes in the standard data set.

⇒ A confirmation message appears (Figure 5.4).




Figure 5.4 Message box confirming autofill

10. Select Standard/Independent Values from the data type drop-down list to view the analyte concentrations for the selected standards (Figure 5.5).

Select Standard/Independent Values from the data type drop-down list.





Figure 5.5 Standard concentrations for the selected analyte

The autofill feature entered the analyte concentrations for the selected standard set based on a user-specified starting concentration of 10,000 and a dilution factor of 3.

Entering Standard Concentrations Manually

1. Open the results of interest (.csv, .xls, .lxd or .mlx).
2. Click the **Well Editor** button . Alternatively, select **Plate > Edit Wells**

from the menu bar.

⇒ The Plate window is in edit mode.



NOTE: You can only edit the standard concentrations when the plate window is in edit mode.

3. Select an analyte from the Bead Set panel (Figure 5.6).
4. Select **Standard/Independent Values** from the data type dropdown list.
5. In the well grid, click the well that you want to edit.
⇒ The concentration value for the selected well is displayed.
6. Enter a concentration value and click **Enter**.
⇒ The well grid displays the new concentration value.
7. To enter other standard concentration values for the same analyte, repeat step 4 and step 5.
8. To enter standard concentration values for a different analyte, repeat step 2 through step 5.

1. Click the **Well Editor** button to put the Plate window in edit mode.

2. Select an analyte.

3. Select **Standard/Independent Values**.

5. Displays the concentration for the selected well. Edit the value and click **Enter**.

The screenshot shows the MasterPlex QT software interface. On the left is the 'Bead Set panel' with a list of analytes: 034 Mouse IL-5, 038 Mouse IL-10, 054 Mouse IL-2, 073 Mouse IFN gamma, and 077 Mouse IL-4. The '034 Mouse IL-5' analyte is selected. In the center is the 'Well grid' table. The table has columns 1 through 7 and rows A through H. The value '0.00' is displayed in cell A1, and it is highlighted with a red border. At the bottom of the window, it says 'Number of Wells: 96' and '034 Mouse IL-5'.


	1	2	3	4	5	6	7
A	0.00	347.00	486.50	321.00			
B	1288.00	409.00	460.50	372.00			
C	1064.50	62.00	68.00	58.00			
D	505.00	68.50	62.00	76.50			
E	181.00	17.50	18				
F	46.50	12.00	12				
G	19.00	-2.00	-2.00	3.00			
H	6.50	-1.00	-2.50	7.00			


Figure 5.6 Plate window

Steps to manually enter a standard concentration

4. Click the well that you want to edit.

Selecting a Model Equation for the Standard Data Set

1. Select an analyte from the Bead Set panel (Figure 5.6).
2. In the well grid, select a standard data set.
3. Right-click the standard data set and select **Assign Model Equation** from the pop-up menu that appears. Alternatively, click the **Select Model** button  or select **Calculations > Model Equations** from the menu bar.
⇒ The Model Equations dialog box appears (Figure 5.7).

 **NOTE:** The Model Equations dialog box only appears if you selected a standard data set.




Figure 5.7 Model Equations dialog box

Model equations available for regression analysis of a standard data set

4. Select a model equation.
5. To apply the selected model to all analytes, choose the **Select this model for all analytes** option.
6. To fix the bottom asymptote of the four parameter logistics curve or the five parameter logistics curve to zero (sets $A = 0$), choose the **Fix bottom at zero** option.
7. To apply weighting during curve fitting, choose the **Use Weighting** option and select a weighting method from the drop-down list.
8. Click **Assign Model**.
9. To choose a model equation for another analyte, repeat step 1 to step 3, and click **Assign Model**.



NOTE: For more information about model equations and weighting methods, see Appendix C.

5.2

Working With Diluted Unknowns


If you need to dilute a sample prior to an assay, you can specify a dilution factor in the well grid. MasterPlex™ QT can compute the diluted or original analyte concentration.

For a diluted unknown:

$$\text{Original concentration} = \text{Dilution factor} * \text{Calculated concentration.}$$

Setting the Concentration Calculation Preference


When you work with diluted unknowns, you must specify whether the software computes the diluted or original analyte concentration.

1. Click the **Preferences** button . Alternatively, select **File > Preferences** from the menu bar.
⇒ The Preferences dialog box appears (Figure 5.8).
2. Click the Calculations tab.
3. To display the diluted analyte concentration for a diluted unknown in Plate window, confirm that the **Calculate dilution factor** option is not chosen.

- To display the original, undiluted analyte concentration for a diluted unknown in the Plate window, choose the **Calculate dilution factor** option.



NOTE: If you change the Undo dilution option, you must recalculate the analyte concentrations.



Choose this option to display the original, undiluted analyte concentration in the well grid.

Figure 5.8 Preferences dialog box, Calculations tab



Editing a Dilution Factor

- Click the **Well Editor** button . Alternatively, select **Plate > Edit Wells** from the menu bar.

⇒ The Plate window is in edit mode.



NOTE: You can only edit a dilution factor when the Plate window is in edit mode.

2. Select **Dilution** from the data type drop-down list.
3. In the well grid, click the dilution factor(s) that you want to edit. To select all of the dilution factors for a group, right-click a dilution factor and choose **Select Entire Group** from the pop-up menu that appears.
4. To edit a single selected dilution factor, use the  arrows. Alternatively, double-click the dilution factor in the well, edit the value in the edit box, and click **Enter**.
5. To edit the dilution factor for several wells simultaneously:
 - a. Click the **Selection Tool** button  to put the plate in selection mode.
 - b. Click the wells that you want to edit,
 - c. Click the edit box and enter the dilution factor.




Figure 5.9 Steps to edit a dilution factor

5.3

Generating Standard Curves & Computing Analyte Concentrations


MasterPlex™ QT carries out a two step calculation sequence when it fits the standard curves. The software:

- fits a standard curve for all defined standard data sets
- interpolates or extrapolates analyte concentrations for the unknown groups that are linked to the standard data set

Standard Curve Preferences

If a standard data set includes replicates, MasterPlex™ QT can fit a standard curve two different ways:

Individual points method (default)	Treats each point in the standard data set individually for curve fitting. For example, if there are three replicates of eight standard wells, the software fits a standard curve using all 24 data points. If the experimental errors are independent, choose this option.
Average standards method	Computes the average of the replicate standard data points and fits the standard curve to the averaged points. If the experiment errors are dependent on each another, choose this option.

1. To view or change the standard curve preference, click the **Preferences** button . Alternatively, select **File > Preferences** from the menu bar.
⇒ The Preferences dialog box appears (Figure 5.10).





Figure 5.10 Preferences dialog box, Calculations tab

2. Click the Calculations tab.
3. Choose a **Standard points** option, click **Apply**, and click **OK**.
4. To return all preference settings to the factory set defaults, click **Default**, and click **OK**.

Generating Standard Curves & Viewing Analyte Concentrations

1. To generate the standard curves and compute (interpolate or extrapolate) the analyte concentrations, click the **Calculate** button . Alternatively, select **Calculations > Calculate STD Curves** from the menu bar.
⇒ A message box confirms the calculations are completed (Figure 5.11).




Figure 5.11 Message box

2. To view analyte concentrations in the Plate window, select **Concentration** from the data type drop-down list, and select an analyte in the Bead Set tab (Figure 5.12).




Figure 5.12 Well grid displays concentration data for the selected analyte

If the standard curve uses a sigmoidal model (for example, the four parameter logistics equation, Figure 5.13), the software interpolates the analyte concentration when:

$$\text{Highest standard MFI} \leq \text{MFI} \leq \text{Lowest standard data point}$$

The software extrapolates the analyte concentration when:

$$A < \text{MFI} < \text{Lowest standard MFI}$$

or

$$\text{Highest standard MFI} < \text{MFI} < D$$

where A is the bottom asymptote and D is the top asymptote of the sigmoidal

curve (Figure 5.13).

A MFI value less than A or greater than D is beyond the range of the standard curve model and the concentration value cannot be extrapolated.

If $MFI < A$, the well grid displays the lowest standard MFI preceded by < (Figure 5.12). If $MFI > D$, the well grid displays the highest standard MFI preceded by >.




Figure 5.13 Four parameter logistics model equation, x-axis log scale


$A = 16.335$ (bottom asymptote), $D = 1373.175$ (top asymptote). MFI values less than A or greater than D are beyond the range of the model equation.

Right-click the graph to select log scale for the x-axis.

5.4

Printing the Well Grid

You can print the contents of the well grid.

1. To set print preferences, click the **Plate Print Settings** button 
⇒ The Print Settings dialog box appears (Figure 5.14).





Figure 5.14 Print Settings dialog box

2. To display a print preview of the well grid, click the **Print** button 
⇒ The Print Preview dialog box displays the currently selected well contents (data type and analyte).




Figure 5.15 Print Preview dialog box

Preview the well grid contents.


3. To print the well grid, click **Print**.

5.5

Working With Standard Curves

After the standard curves are generated, you can view the curves, residual plots, the standard data, and information about the standard data in the Standard Curves window (Figure 5.16).

Opening the Standard Curves Window

1. Click the **Open Standard Curve** button . Alternatively, select **Plate > View Curves** from the menu bar.

⇒ The Standard Curves window opens and displays the Standard Curve Chart tab (Figure 5.16).

The window displays a tree of local standard curves, imported standard curves, and regression curves for the active plate.





Figure 5.16 Standard curves window, Standard Curve Chart tab

2. To display a standard curve, click the curve(s) of interest in the data tree. To select multiple curves, press and hold the **Ctrl** key while you click the curves.



NOTE: If the Standard Curves window shows data points only, you need to generate the standard curves (click the **Calculate** button ).

Viewing the Residual Plot

The residual plot helps you identify outliers and detect non-constant variability in a standard curve. In a residual plot, each point represents the difference between the expected standard concentration and the calculated standard concentration for a selected analyte.

$$\text{Residual} = \text{Observed (or calculated) concentration} - \text{Expected concentration}$$

The observed (calculated) concentration is interpolated or extrapolated from the

standard curve. The expected concentration is the standard data value that you enter using the autofill function (or manually). If the residual is positive, we have detected a concentration that is greater than the expected value. If the residual is negative, the detected concentration is less than the expected value.

1. To view the residual plot, click the Residual Plot tab in the Standard Curves window.
 2. Click the standard curve of interest.
- ⇒ The residual plot is displayed (Figure 5.17).




Figure 5.17 Residual plot

Each point represents the difference between the expected and calculated concentrations of the standard data for the selected analyte.

The residual plot helps you identify non-constant variability (*heteroscedasticity*) in the data. For more information on detecting heteroscedasticity using a residual plot, see Appendix C.

Viewing Standard Data & Specifying Outliers

In the Standard Data tab, you can view information about each standard and

STANDARD CURVES & ANALYTE CONCENTRATION

specify outliers in a standard data set. The software ignores outliers in the calculation of the standard curve.

1. In the Standard Curves window, click the Standard Data tab.

2. Click the standard curve of interest.

⇒ The data for the selected standards are displayed (Figure 5.18).

Table 5.1 explains the information available in the Standard Data tab.






Figure 5.18 Standard Data tab

Table 5.1 Standard Data

Item	Description
Outliers	To make a standard data point an outlier, place a check mark in the check box next to the standard data point. The software ignores outliers during the calculation of the standard curve.
MFI	Mean fluorescence intensity of the data point in the selected analyte standard.
Expected	The standard concentration for the selected analyte that was entered using the autofill function (or manually).
Calculated	The standard concentration for the selected analyte that was interpolated or extrapolated from the standard curve.
Residuals	Calculated - Expected
% Recovery	(Calculated/Expected) x 100
Background	MFI background value for the plate.

3. To specify an outlier, place a check mark next to the standard data point of interest (Figure 5.18).


Alternatively, you can select outliers in the Standard Curve Chart tab:

- a. Click the Standard Curve Chart tab.
 - b. Click the **Chart Selection Tool** button  then click a point in the standard curve that you want to make an outlier.
 - ⇒ The Standard Data tab highlights the selected point.
 - c. Place a check mark in the highlighted outlier check box.
4. Recalculate the standard curve and analyte concentrations. (Click the **Calculate** button  or select **Calculations > Calculate STD Curves** from the menu bar.)

5.6

Importing Standard Curves

You can import standard curves from another plate so that you can easily compare standards from different plates.

1. Open the plate that contains the standard curves for export (source plate).
2. Open the plate that the standard curves will be imported to (destination plate).
3. Select the  **Import standards** menu from plate menu.
 - ⇒ The plate windows are in import/export mode.



NOTE: The plate windows must be in import/export mode to export standard curves from the destination plate to the source plate.

4. Use the click-and-drag method to move the standard data set of interest from the source plate to the well grid of the destination plate.
 - ⇒ The Standard Curves window displays a list of the standard curves imported from the source plate (Figure 5.19).

STANDARD CURVES & ANALYTE CONCENTRATION





Figure 5.19 Standard curves can be exported

Local and imported standard curves are listed in the Standards tab of the destination Plate window.



NOTE: To view and conveniently compare standard curves, open the Standard Curves window (click the  button). For more information, see *Opening the Standard Curves Window* on page 5.16.

CHAPTER 6

A virtual plate is a software simulation of an empty microtiter plate. User-selected unknown and standard data from one or more actual plates (.csv, .xls, .lxd or .mlx) can be added to a virtual plate. In a virtual plate you can:

- Compare and analyze samples from user-selected actual plates (.csv, .xls, .lxd or .mlx)
- Specify custom row and column dimensions to increase analysis throughput
- Generate a Dose-Response curve for user-selected data





Figure 6.1 Add user-selected data from actual plates (.csv, .xls, .lxd or .mlx) to a virtual plate

In a virtual plate, you can perform a regression analysis on user-selected data to generate a Dose-Response curve and determine the LogEC50 value.

6.1

Creating a Virtual Plate

1. Open the Luminex results files (.csv, .xls or .lxd) or MasterPlex™ QT files (.mlx) that are the data sources for the virtual plate.
2. If the Plate Wizard is not open, click the **Plate Wizard** button 
⇒ The Plate Wizard appears (Figure 6.2).



NOTE: If the **Display wizard at start up** option is chosen, the Plate Wizard is displayed when MasterPlex™ QT starts.



Figure 6.2 Plate Wizard, Welcome tab

3. Click **Next**.
⇒ The Select Plate Type tab appears (Figure 6.3).





Figure 6.3 Plate Wizard, Select Plate Type tab

4. Choose the **Virtual Plate** option and click **Next**.
⇒ The Plate Dimensions tab appears (Figure 6.4).




The row and column defaults can be edited in the Preferences dialog box (click the  button).

Figure 6.4 Plate Wizard, Plate Dimensions tab

5. Enter the number of rows and columns for the virtual plate. Click **Finish**.
⇒ A Plate window opens and displays the empty well grid of the virtual plate (Figure 6.5).





Figure 6.5 Virtual plate

Selecting Data from a Source Plate

The virtual pipette copies (*aspirates*) data from user-selected wells in a source plate and pastes (*dispenses*) the data into a virtual plate. The virtual pipette aspirates all of the analyte data in a well, including the computed analyte concentrations. It remains loaded until you dispense or clear the pipette.



NOTE: The data source plates must contain the same type and number of analytes, otherwise concentrations cannot be calculated. If the source plates contain the same number of analytes, but they are named differently, use the virtual analyte filter to rename analytes so that the nomenclature is consistent. (See *Selecting and Renaming Analytes* on page 6.9.)

1. Click the **Virtual Pipette** button .
⇒ The virtual pipette appears when the mouse cursor is positioned over a plate window.

2. In the source plate, select the wells of interest.

To select adjacent wells, press and hold the mouse button while you drag the mouse pointer to select the wells of interest.




NOTE: Selecting non-adjacent wells is not recommended.

3. Right-click the selected wells and select **Aspirate** from the pop-up menu that appears (Figure 6.6).

⇒ The data for the analytes in the selected wells are added to the virtual pipette and is ready to dispense into a virtual plate.



NOTE: If the background is subtracted in the source plate, the virtual pipette aspirates and transfers background-subtracted values. If you do not want to aspirate background-subtracted values, make sure the background subtraction is turned off before you aspirate data into the virtual pipette. (Click the  button to turn background subtraction on or off.)

	1	2	3	4
A	0.00	347.00	486.50	321.00
B	1288.00	409.00	460.50	372.00
C	1064.50	62.00	68.00	58.00
D	50		62.00	76.50
E	18		18.00	17.00
F	4		12.00	14.00
G	19.00	-2.00	-2.00	3.00
H	6.50	-1.00	-2.50	7.00




Figure 6.6 Well grid

Right-click selected wells to display the pop-up menu.

4. To clear the data from the virtual pipette, right-click and select **Clear** from the pop-up menu.

Adding Data to a Virtual Plate

After the virtual pipette aspirates data from the source plate, it is ready to dispense the data into the virtual plate.


1. Position the virtual pipette over the virtual plate.
2. Click the first well to which the data will be added.
3. Right-click the well and select **Dispense** from the pop-up menu that appears.
⇒ The data are added to the virtual plate (in the same configuration as in the source plate) (Figure 6.7).




NOTE: If the number or names of the analytes in the virtual pipette is different from that in the virtual plate, the virtual analyte filter automatically appears. For more information on using the filter, see *Selecting and Renaming Analytes* on page 6.9.




NOTE: Data in a virtual plate cannot be removed, but can be overwritten.

4. After you finish adding data to the virtual plate, click the **Virtual Pipette** button  to turn off the virtual pipette.

1. Open a .mlx or .csv. Click the  button to turn on the virtual pipette.

2. Select the wells of interest in the source plate (.csv or .mlx). Right-click the selected wells and choose aspirate from the pop-up menu.



IL5 Project.mlx <IL5.csv>


H1 20.00

Bead Set Properties Star MFI 1288.00

	1	2	3	4
A	0.00	347.00	486.50	321.00
B	1288.00	409.00	460.50	372.00
C	1064.50	62.00	68.00	58.00
D	505.00		62.00	76.50
E	181.00		18.00	17.00
F	46.00		12.00	14.00
G	19.00	-2.00	-2.00	3.00
H	6.50	-1.00	-2.50	7.00

Number of Wells:96 034 Mouse IL-5

3. In the virtual plate, select the first well where you want to dispense the data.



Untitled Plate 2 <Virtual Plate>


A3 0.00

Bead Set Properties MFI

	1	2	3	4
A				
B				
C				
D				
E				
F				
G				
H				

Number of Wells:96

4. Right-click the well and select **Dispense** from the pop-up menu. The data are added to the virtual plate (starting at the selected well) in the same configuration as in the source plate.



Untitled Plate 3 <Virtual Plate>


A3 0.00

Bead Set Properties MFI

	1	2	3	4
A				
B				
C				
D				
E				
F				
G				
H				

Number of Wells:96

Figure 6.7 Adding data to a virtual plate

Open a source plate (.mlx or .csv, .xls or .lxd) and create a virtual plate (click the  button to start the virtual plate wizard).

6.2

Working With the Virtual Analyte Filter

In a multiplex assay, all of the plate wells must contain:

- the same types of analytes (bead sets) with the same nomenclature
- the same number of analytes

This is true for virtual plates as well. When you add data to a virtual plate, MasterPlex™ QT compares the name and number of the analytes in the virtual pipette to those in the virtual plate. The virtual pipette will not dispense if there are discrepancies between the number or names of analytes in the pipette and the virtual plate. If the number of analytes in the pipette is greater than that in the destination plate, the virtual analyte filter automatically appears (Figure 6.9).

The virtual analyte filter displays a list of the analytes that are present in the virtual pipette. It enables you to choose the analytes that you want to add to the virtual plate and, if necessary, rename them to be consistent with the number and name of analytes in the virtual plate.

If you add data to a virtual plate from source wells that contain different analyte names or a different number of analytes, data holes are created. As a result, a well in the virtual plate appears blank if the analyte selected in the analyte panel is not present in the well. If a plate (.csv, .xls, lxd, .mlx, or virtual) contains data holes, the concentrations cannot be calculated.



NOTE: In order to prevent data holes, if the number of analytes in the virtual pipette is less than the number of analytes in the destination plate, the data cannot be added to the virtual plate.




Figure 6.8 Virtual analyte filter shows the analytes in the virtual pipette

Selecting and Renaming Analytes

If the virtual analyte filter appears, you must select and, if necessary, rename the analytes to match the number and names of the analytes in the virtual plate.

1. In the virtual analyte filter (Figure 6.9), place a check mark next to each analyte that you want to add to the virtual plate. To select all analytes for the virtual plate, click **Check All**.
2. To rename an analyte so that it is consistent with the nomenclature in the virtual plate:
 - a. Click **here to assign** next to the analyte that you want to rename.
 ⇒ A drop-down list shows the names of the analytes in the virtual plate (Figure 6.9).
 - b. Select a name from the drop-down list.
 ⇒ The virtual analyte filter displays the new name for the analyte.




Figure 6.9 Virtual analyte filter

3. To save the renaming assignments for use again with the same source plate (.csv, .xls or .lxd or .mlx) during the current session, choose the **Save this assignment** option.

If you want to aspirate other data from the same source plate, choose the **Use last saved assignments** option in the virtual analyte filter to automatically rename all of the analytes in the filter.

4. Click **OK**.

⇒ The data are added to the virtual plate and the virtual analyte filter closes.


6.3

Generating a Dose-Response Curve


A Dose-Response curve plots the dose of an agent (independent variable) on the x-axis and the response (computed concentration) on the y-axis.


After you have created a virtual plate that contains the concentration data of interest, you can run a regression analysis on the data in the virtual plate to generate a Dose-Response curve (Figure 6.10). (For more information on Dose-Response curves, see *Appendix D, Dose-Response Analysis* on page D.1.) The software determines the LogEC50 value from the Dose-Response curve. The logEC50 value is the log of the agent concentration (for example, a drug or other chemical) that produces a response half way between the baseline and maximum response.


1. Select concentration from the data-type drop-down list.



2. Create a virtual plate.
3. Transfer response data to the virtual plate using the virtual pipette.
4. Mark and group wells into a regression data set.
5. Select Standard/Independent values from the data-type drop down list and enter the agent concentrations (use the Autofill function or enter manually).



6. To select a Dose-Response model, right-click the regression data set and select **Assign Model Equation**.
7. Click the **Calculate** button  , then select regression radio button.

8. To view the Dose-Response curve, click the Standard Curves button .








Figure 6.10 Steps to generate a Dose-Response curve

1. Create a virtual plate and use the virtual pipette to add the data for the Dose-Response curve to the virtual plate. (See “Creating a Virtual Plate” on page 6.2.)



NOTE: After you finish creating the virtual plate, click the  button to turn off the virtual pipette.

2. Select the wells that you want to include in the regression analysis data set. The regression analysis data set includes the wells that contain the response data to different dosages of the agent.
3. Right-click the selection and choose **Mark Wells > Regression** from the pop-up menu that appears. Alternatively, click the **Mark regression wells** button  or select **Template > Mark Wells > Mark Regression Wells** from the menu bar.
4. Right-click the regression data set and select **Group Wells** from the pop-up menu. Alternatively, click the **Group** button  or select **Template > Group Wells** from the menu bar.
5. If the agent was serially diluted, follow these steps to enter the agent concentration using the auto fill function:
 - a. Choose **Standard/Independent Values** from the data type dropdown list.
 - b. Right-click the regression data set and select **Auto Fill** from the pop-up menu that appears. Alternatively, click the **Auto Fill** button 
⇒ The Auto Fill dialog box appears (Figure 6.11).






Figure 6.11 Auto Fill dialog box

Use the auto fill function to set concentrations for a serially diluted agent.

- c. Enter the starting concentration for the analyte displayed in the Analyte drop-down list.
 - d. Enter the dilution factor
 - e. If the starting concentration and dilution factor are the same for all analytes, choose the **Fill in for all bead sets** option. If not, make another selection from the Analyte drop-down list and repeat step b and step c.
 - f. Click **Fill** after the starting concentration and dilution factor are entered for all analytes.
6. If the agent was not serially diluted, follow these steps to manually enter the agent concentration for each analyte:
- a. Click the **Edit mode** button .
 - b. Choose **Standard/Independent** from the data type drop-down list.
 - c. Select an analyte in the Analyte panel.
 - d. Click a well in the regression data set and enter the agent concentration in the edit box.
 - e. Click **Enter**.

⇒ The value is displayed in the well grid.
 - f. Repeat step b to step e for the remaining analytes.

7. Right-click the regression data set and select **Assign Model Equation** from the pop-up menu that appears. Alternatively, click the regression data set and click the **Select Model** button .
- ⇒ The Model Equations dialog box appears (Figure 6.12).






Figure 6.12 Model Equations dialog box

Select a model for the Dose-Response curve

8. Select the Dose-Response (variable slope) model equation and choose the **Select this model for all analytes** option.
9. Click **Assign Model**.
10. Click the **Calculate** button  then select the regression.
 - ⇒ The software computes the Dose-Response curve and displays a confirmation message when the calculation is completed.

11. To view the Dose-Response curve, click the **Standard Curves** button . The Dose-Response curve shows the agent dose (independent variable) on the x-axis and the response (interpolated concentration) on the y-axis (Figure 6.13). It displays the parameters of the model equation and the LogEC50 value.
12. To compute the EC50 value, take the antilog of the LogEC50.




Figure 6.13 Dose-Response curve


Dose (independent variable) on the x-axis and response (computed concentration, dependent variable) on the y-axis.

6.4

Auto Plate Merging

The purpose of auto plate merging is to quickly merge multiple plates, at analytes dimension, into a virtual plate. That is, assuming all plates have different analytes, we want to merge plates together by adding all their analytes as new. So, if a plate with 5 analytes is merged with another plate with 7

analytes, we will get a virtual plate with 12 analytes. This tool is especially useful for miRNA analysis where assay results of the same sample pool are usually distributed in several plates. It will be extremely useful if we can put these results into one place and do analysis together. Compared to use virtual pipette to manually add wells, this tool makes the merging task much faster and seamless.

1. Start the plate merge wizard by clicking the toolbar button  “Merge Plates”.
2. You will see a list of plate names which are currently opened (Figure 6.14). If you want to add more plates to merge, you can click Open button to load more plates. Then, you need to check the plates to be merged from the list, and go to the next page.




Figure 6.14 Selecting plates

3. Now you will see a list of bead (analyte) names which appear in all selected plates(Figure 6.15). If you confirm they are the common house keeping genes and want to exclude them from the virtual plate, you may choose “Automatically” merge, then finish wizard. Otherwise, you may choose to “Manually select house keepers”, then go to the next step.




Figure 6.15 Merge method selection

3. At this step, you can manually select analytes(Figure 6.16), from each plate, to be excluded from the virtual plate. Then, go the next page to finish wizard.




Figure 6.16 Manual analyte selection




NOTE: If there are duplicated analyte names in the final virtual plate, the tool will automatically rename them by adding [n] at the end of the name. “n” indicates the number of duplicates. Finally, a report of the renaming process will be presented in a popup window.

CHAPTER 7

MasterPlex™ GT displays MFI, count, or concentration data in several graph formats in the Data Chart window.

7.1

Viewing a Data Chart

1. In the Plate window, select one or more wells of interest.
To select adjacent wells, press and hold the mouse button while you drag the mouse pointer to select the wells of interest. To select nonadjacent wells, press and hold the **Ctrl** key while you click the wells of interest.
2. Click the **Open Data Charts** button . Alternatively, select **Plate > Data Charts** from the menu bar.
⇒ The Data Chart window opens and displays the well chart(s) (default) for the selected well(s) (Figure 7.1).






Figure 7.1 Multi-well chart

MFI data for the selected analytes in well B1

3. To display another graph format for the selected wells, click the button  and make a selection from the drop-down list.
 - ⇒ The well data are displayed in the selected format.
4. To change the data type in the chart (MFI, count, or concentration), make a selection from the data drop-down list (Figure 7.1).
5. To remove an analyte from a chart, remove the check mark next to the analyte name in the analyte panel.
6. To view other well data in the current chart format:
 - a. Return focus to the Plate window.
 - b. Select the well(s) of interest.
 - c. Click the Data Chart window.
 - ⇒ The Data Chart window is updated and displays the new well data in the selected format.
7. To change the y-axis scale, move the slider up or down the scale bar.

7.2

Data Chart Types

MasterPlex QT provides seven different chart types (Table 7.1). To choose a chart format for the selected well data, click the  button, and make a selection from the drop-down list of data charts.







NOTE: The data chart toolbar is available if a Data Chart window is open.

Table 7.1 MasterPlex™ QT data charts

Data Chart	Displays...	See Page
Multi-Well Chart	Separate bar graphs of analyte data for each user selected well.	7.4
3D Depth Chart	Bar graph of analyte data that includes all selected wells. Analytes from a particular well are organized by row.	7.6
Inverted 3D Depth Chart	Bar graph of analyte data that includes all selected wells. Analytes from a particular well are organized by column.	7.7
Intensity Map	A qualitative map of the intensity data for the analytes in one well. The map colors indicate the relative intensity of an analyte compared to the other analytes in the same well.	7.8
Single Analyte Chart	Bar graph of a single analyte across all user selected wells.	7.10
Stacked Analyte Chart	Bar graph of user selected wells. Each color-coded segment of a bar represents the contribution of an analyte to the total data value for the well.	7.11
Stacked Well Chart	Bar graph of analytes. Each color-coded segment of a bar represents the contribution from a well to the total data value for the analyte.	7.12

Multi-Well Chart

The multi-well chart displays a bar graph of analyte data (MFI, count, or concentration) for a user-selected well(s). Each bar represents an analyte. The multi-well chart is the default format that is displayed when you open the Data Chart window.

1. To minimize a multi-well, click the  arrow (Figure 7.2).
 2. To restore a chart, click the  arrow.
Alternatively, click the  button to minimize all charts (Figure 7.3).
Click the button again to restore all charts.
 3. To normalize the data for all open multi-well charts, click the **Normalize** button .
- ⇒ The well data are plotted using the y-axis scale from the chart with the

largest y-axis value.



NOTE: You can choose different colors for the chart bars in the Color tab of the Plate window (see *Plate Window Components* on page 3.9).

Select the data type for the chart




Figure 7.2 Multi-well chart for C2

Multi-well charts for A2 and well B2 are minimized. Click a bar to display a tool tip (analyte name or y-axis data). Right-click a bar for tool tip options.




Figure 7.3 Multi-well charts are minimized



TIP: You can conveniently reorder the multi-well charts in the Data Chart window by dragging a chart to the desired position in the list.

3D Depth Chart

The 3D depth chart is a convenient way to compare the MFI, count, or concentration data in one well and across different wells (Figure 7.4). Each bar represents an analyte. The analytes of a particular well are identified by color and row.





Figure 7.4 3D depth chart compare analyte data in the same well or between different wells.

Each bar represents an analyte. Each row of color-coded bars represents the analyte data from one well. Click a bar to display a tool tip (analyte name or y-axis data). Right-click a bar for tool tip options.

Inverted 3D Depth Chart

Each bar of the chart represents an analyte in a selected well. The bars are organized into columns and each column represents the analytes in respective selected wells (Figure 7.5).



Each column of bars represents
the analytes in one well

Figure 7.5 Inverted 3D depth chart

Click a bar to display a tool tip (analyte name or y-axis data). Right-click a bar for tool tip options.

Intensity Map

The intensity map is a qualitative picture of the intensity data in the wells selected in the Plate window. The map provides a convenient way to look at a large number of samples and identify analytes that consistently have a high (or low) intensity compared to the other analytes in the same well.



NOTE: The intensity map compares and ranks data in the same well. It does not compare data across wells.




Figure 7.6 Intensity map

Click the map to display a tool tip (analyte name or y-axis data). Right-click the map for tool tip options.

The intensity map is color-coded. Figure 7.6 shows the default colors:

- Red = peak intensity
- Green = mid-range intensity
- Blue = low intensity

The default colors may be changed. (See *Intensity Map Settings* on page 7.35.) Each row in the map represents the relative intensities of the analytes in one well.




Figure 7.7 Intensity map

In this example map, IL-1Beta has the highest intensity and GM-CSF has the lowest intensity in the third well that was selected in the Plate window.

Single Analyte Chart

This chart shows the data (MFI, count, or concentration) for a single analyte across all user-selected wells (Figure 7.8). Each bar represents the selected analyte data from one well.




Figure 7.8 Single analyte chart

Each bar represents data for a single analyte from one well. Click a bar to display a tool tip (analyte name or y-axis data). Right-click a bar for tool tip options.

Stacked Analyte Chart

Each bar represents a single well. The color-coded segments of a bar represent the contribution of each analyte to the total MFI, count, or concentration value in the well (Figure 7.9).




Figure 7.9 Stacked analytes chart

Click a bar to display a tool tip (analyte name or y-axis data). Right-click a bar for tool tip options.

Stacked Well Chart

Each bar in the stacked well chart represents the sum total of the MFI, count, or concentration values of an analyte across the user-selected wells. The color-coded segments of a bar show the contribution of each well to the total value.




Figure 7.10 Stacked well chart












Click a bar to display a tool tip (analyte name or y-axis data). Right-click a bar for tool tip options.

7.3

Working With Data Charts

Table 7.2 explains the toolbar button functions in the Data Chart window. The available buttons depend on the type of data chart that is open.

Table 7.2 Toolbar buttons in the Data Chart window

Button	Click to...
	Choose another data chart format. The software plots the data for the selected wells in the new format.
	Toggle the chart display between a two and three dimensional view.
	Put the chart in rotate mode so that you can rotate the chart around the x or y-axis.
	Put a chart in zoom mode so that you can magnify a user selected area of the chart.
	Opens the Chart Print Preview window so a chart may be printed.
	Normalize the data for the open multi-well charts. Plots the selected well data using the y-axis scale from the chart with the largest y-axis value. Note: This function is only available for the multi-well chart.
	Select all analytes in the analyte panel for display in the chart.
	Removes all analytes from the chart.
	Minimize or restore all well charts in the Data Chart window.
Menu	Click to...
	Hide or show the analyte panel.
	Hide or show the scale bar. To change the y-axis scale, move the scale slider up or down the scale bar.

Copying, Saving, or Sending a Chart

You can copy a chart to the system clipboard or save a chart.

1. Open the chart of interest.
2. Select **Edit > Export Chart** from the menu bar.
⇒ The Export Chart dialog box appears (Figure 7.11).




Figure 7.11 Export Chart dialog box

3. In the Picture tab, choose a file format option:
Metafile saves the chart in a scalable .wmf format. Choose the **Enhanced** option to save the charts in enhanced metafile format (.emf).
Bitmap saves the chart in .bmp format. Depending on the chart dimensions, the .bmp file size can become large.
VML (HTM) saves the chart as an HTML document (.htm) that can be viewed using a web browser.
4. To specify a size option for the chart:
 - a. Click the Size tab.
 - b. Enter the width and height (in pixels) for the chart.
 - c. Choose the **Keep aspect ratio** option to maintain the same dimension ratio at the new chart size.
5. To copy the chart, click **Copy**.
6. To save the chart, click **Save**.
⇒ The Save As dialog box appears (Figure 7.12).




Figure 7.12 Save As dialog box

7. To e-mail the chart:

a. Click **Send**.

⇒ The Choose Profile dialog box appears (Figure 7.13).




Figure 7.13 Choose Profile dialog box

b. Select a recipient for the e-mail and click **OK**.



NOTE: This dialog is only functional with Microsoft Outlook.

Copying or Saving Data

The software can export data from a chart to other applications. The data may be copied to the system clipboard or saved in different file formats.

1. Open the chart of interest.
2. Select **Edit > Export Chart** from the menu bar.
⇒ The Export Chart dialog box appears (Figure 7.11).
4. Click the Data tab.




Figure 7.14 Export Chart dialog box, Data tab

4. Make a selection from the **Series** drop-down list. To choose the data from all the wells in the chart, select **All**.
5. Choose a file format:
 - Text** option saves the data in plain ASCII text format (Figure 7.16).
 - XML** option saves the data in extensible markup language.
 - HTML Table** option saves the data to a table in an HTML file.
 - Excel** option saves the data in .xls file format. Copying to the system clipboard is not available for this option.
6. If you are exporting data to a text file, make a selection from the Delimiter drop-down list. The delimiter specifies how data values are separated.
7. If desired, choose the **Point Index**, **Point Labels**, or **Header** options.

Point Index	Numbers the members of the data series.
Point Labels	Displays the analyte name for each series.
Header	Displays a header for each column of the data series.

8. To copy the data series to the system clipboard, click **Copy**.
9. To save the data to a file:
 - a. Click **Save**.

⇒ The Save As dialog box appears (Figure 7.15).
 - b. Enter a name for the file.
 - c. Select a file format from the Save as type drop-down list.
 - d. Click **Save**.




Figure 7.15 Save As dialog box




Figure 7.16 Example of data (series) for well A3 saved to a text file

Printing a Chart

1. Click the **Print current chart** button .

⇒ The Chart Print Preview window opens (Figure 7.17).






Figure 7.17 Chart Print Preview window

3. To print the chart, click **Print**.

Rotating or Magnifying a Chart

1. To rotate the chart around the x or y-axis:
 - a. Click the **Chart rotate tool** button .
 - b. Press and hold the mouse button while you move the mouse pointer to rotate the graph.
2. To zoom in on a chart area:
 - a. Click the **Chart zoom** button .

- b. Press and hold the mouse button while you draw a rectangle (from the upper left to lower right corner) over the graph area that you want to magnify.
3. To zoom out and display a chart without magnification, press and hold the mouse button while you draw a rectangle (from the lower right to upper left corner) over the graph.

Displaying Chart Tool Tips

1. To show a tool tip that displays information in a bar chart:
 - a. Click the bar of interest.
 - ⇒ A tool tip displays the analyte name.
 - b. To change the tool tip display, right-click a well and select **Mouse Tip Styles > Y Values** from the pop-up menu that appears (Figure 7.19).
 - ⇒ The tool tip displays the MFI, concentration, bead count value for the well (depending on the data type selected in the analyte panel)
2. To change the tool tip display, right-click a well and select **Mouse Tip Styles > Y Values** from the pop-up menu that appears (Figure 7.19).
 - ⇒ The tool tip displays the MFI, concentration, bead count value for the well (depending on the data type selected).




Figure 7.18 Multi-well chart

Click a bar to display a tool tip.

Select the data type for the chart y-axis (MFI, bead count, or concentration).




Figure 7.19 Multi-well chart

Right-click a bar to display a pop-up menu of options.

Modifying Chart Properties

You can modify the appearance of a chart.

- Select **Chart > Chart Properties** from the menu bar.
⇒ The Chart Properties dialog box appears and displays a list of the charts in the Data Chart window (Figure 7.20).




Figure 7.20 Chart Properties dialog box



NOTE: The types of chart property options available in each tab depends on the type of chart that is open.

Charts Tab


Alternative formats are available for the multi-well chart.

1. To view the available multi-well chart formats, click **Modify Series** (Figure 7.20).
⇒ The Chart Gallery appears (Figure 7.21).



NOTE: The **Modify Series** button is only available is a multi-well chart is open in the Data Chart window.

Click the down arrow to display a sub-menu of bar chart formats.



Sub-menu of bar chart formats.




Figure 7.21 Chart gallery of available formats for the multi-well chart

2. To change the chart format, click the format of interest.
3. To view a sub-menu of styles for the selected format, click the down arrow (Figure 7.21).
4. If you want to display the chart format in 3 dimensions, select the **Show 3D** option.
5. Click **OK**.
⇒ The selected format is applied to the multi-well chart.

3D Tab

To view 3D chart properties, click the 3D tab in the Chart Properties dialog box (Figure 7.22).




Figure 7.22 Chart Properties dialog box, 3D tab

User-Modifiable 3D Chart Properties

Enable 3D Choose this option to display a three dimensional view of the chart.

3D% Increases or decreases the three dimensional thickness of the graph bars.

Orthogonal Choose this option to display a front view of the chart (no rotation around the vertical axis).

Angle Enter a value to change the display angle for the graph bars when the orthogonal option is chosen.

Rotation Move the slider to rotate the chart around its vertical axis.

Elevation Move the slider to rotate the chart around its horizontal axis.

Perspective Move the slider to change the perspective angle of the chart. The perspective effect is the dimensional appearance of the chart with respect to distance from the viewer. If this option is set to zero, the perspective option is disabled.

Note: If the orthogonal option is chosen, the chart perspective cannot be changed.

Horiz. Offset Move the slider to adjust the horizontal position of the chart in the Data Chart window.

Vert. Offset Move the slider to adjust the vertical position of the chart in the Data Chart window.

Legends Tab

To view the chart legend properties, click the Chart tab in the Chart Properties dialog box (Figure 7.23). These settings determine the location and appearance of the legend in a data chart (Figure 7.24).




Figure 7.23 Chart Properties dialog box, Legends tab




Figure 7.24 Well chart

Settings in the Legends tab determine the location and appearance of a chart legend.

User-Modifiable Legend Settings

Gradient

Visible Choose this option to display the legend with a gradient background.

No Middle

Direction Make a selection from the drop-down list to set the direction of the gradient displayed in the legend.

Start, Mid, End Click to open a color palette. Make a selection from the color palette to set the background color for the start, middle, or end portion of the gradient.

Positions Choose an option to display the legend to the left, right, top, or bottom of the chart.

Shadows

Color Click to open a color palette. Make a selection from the color palette to set the color for the shadow at the legend perimeter.

Size Select a number to set the thickness of the shadow at the perimeter of the legend.

Text style Make a selection from the drop-down list to select among

display options that include order of items, placement of data, hide or unhide data, display percentage values, and display the x-axis value in the legend.

Show Legends Choose this option to display a legend in a data chart.

Walls Tab

1. To view the chart wall properties, click the Walls tab in the Chart Properties dialog box (Figure 7.25).

These settings determine the appearance of the left, back, and bottom walls in a three dimensional chart (Figure 7.25).

2. To change a setting, click a wall (Left, Back, or Bottom), and modify the settings.




Figure 7.25 Chart Properties dialog box, Walls tab




Figure 7.26 3D depth chart

User-Modifiable Wall Settings

Appearance

- Color** Click to open a color palette. Choose a color for the wall from the color palette.
- Visible** Choose this option to display the wall color.
- Thickness** Change this value to increase or decrease the wall thickness.
- Dark 3D** Choose this option to apply a darker color to the inside surface of the wall.
- Transparent** Choose an option to display the a transparent wall without color.

Gradient

- Visible** Choose this option to display a gradient wall color.
- No Middle** Choose this option if you do not want to display a gradient with only two colors (start and end color). If this option is not chosen, the gradient has three colors (start, middle, and end color).
- Direction** Make a selection from the drop-down list to set the direction of the gradient.
- Start/Mid/End** Click to open a color palette. Choose a color for the start, middle, or end of the gradient. Note: If the **No Middle** option is chosen, the **Mid** button is not displayed.

Series Tab

A series is a group of related bars in a chart. For example, in the 3D depth chart, the bars that represent the analytes in a particular well comprise a series (Figure 7.27). In the series tab, you can modify the display of a series.

1. To view the series properties for an open chart(s), click the Series tab in the Chart Properties dialog box (Figure 7.28).
2. To modify the settings for a series, make a selection from the Series drop-down list. If you want to apply the settings to all series in the graph, choose the **Apply to all series** option.




Figure 7.27 3D depth chart

In the 3D depth chart, the bars that represent the analytes in a well are an example of a series.




Figure 7.28 Chart Properties dialog box, Series tab, Marks settings

User-Modifiable Series Settings

Marks Settings

These settings determine the display of the name labels (*marks*) in a data chart (Figure 7.27).

Gradient

- Visible** Click to display a gradient background in the name label.
- Direction** Make a selection from the drop-down list to set the direction of the gradient.
- Start/Mid/End** Click to open a color palette. Choose a color for the start, middle, or end of the gradient.
- No Middle** Choose this option to display a two color gradient. If this option is not chosen, the gradient has three colors (start, middle, and end color)
- Visible** Choose an option to display the name labels for the selected series.
- Marks Angle** Click the up or down arrows to set the display angle for the name labels.

Shadow

- Color** Click to open a color palette. Choose a color for the shadow of the name label.
- Size** Click the up or down arrows to set the thickness of the name label shadow.

General Settings




Figure 7.29 Chart Properties dialog box, Series tab, General settings

- Series Color** Click to open a color palette. Choose a color for the selected series.
- Horizontal Axis** Make a selection from the drop-down list to display the horizontal axis labels at the bottom, top, or bottom and top of the chart.
- Vertical Axis** Make a selection from the drop-down list to display the vertical axis labels to the left, right, or left and right of the chart.
- Show in Legend** Choose this option to display the series name in the chart legend.
- Color Each Point** Choose this option to display each bar in the series in a different color.
- Sort** Make a selection from the drop-down list to choose a descending or ascending sort order for the series. Note:
After a series is sorted, it cannot be unsorted. To display the series prior to sorting, simply close the chart, then regenerate it.

Bar Settings




Figure 7.30 Chart Properties dialog box, Series tab, Bar settings

- Multi Bars None** Choose this option to display the bars in the default orientation. For example, Figure 7.31 shows the default orientation of the 3D depth chart.
- Side** Choose this option to display the bars side-by-side along the x-axis (Figure 7.31).
- Stacked** Choose this option to combine the separate bars of a series into one vertical column (Figure 7.32).
- Stacked 100%** Choose this option to display the percentage data for a series in one vertical bar. Each segment of the bar represents the percent a particular series value contributes to the total bar value (Figure 7.32).
- Side All** Choose this option to display all bars of the same type in a side-by-side orientation.
- Appearance**
- Bar Style** Make a selection from the drop-down list to choose a bar style.
- % Bar Width** Click the up or down arrows to change the bar width.
- Bar Dark 3D** Choose this option to display three dimensional bar demarcations.
- Border** Click **Border** to open a dialog box of style options for borders of bars.




Figure 7.31 3D depth chart

None option is the default display (left). Side option (right).




Figure 7.32 3D depth chart

Stacked option (left). Stacked 100% option (right)




Figure 7.33 3D depth chart

Side All option.

Intensity Map Settings

If an intensity map is open, the Series tab shows Intensity settings (Figure 7.34).





Figure 7.34 Chart Properties dialog box, Series tab

If an intensity map is open in the data chart window, the intensity map color scheme settings are available.

To change a color (peak, mid-range, or low color) in the intensity map:

1. Click the button for the color that you want to change (for example, Peak ).

⇒ The color palette appears (Figure 7.35).




Figure 7.35 Color palette

2. To select a predefined color, click one of the basic colors.

3. To define a custom color:
 - a. Click **Define Custom Colors**, then use the click-and-drag method to move the cross hairs in the custom color field.
 - b. Adjust the color brightness using the luminosity scale to the right.
 - c. When finished, click **Add to Custom Colors** to apply the color.
4. Click **OK** to close the color palette.

Panel Tab

The panel settings determine the appearance of the display area behind the chart.




Figure 7.36 Chart properties dialog box, Panel tab

Gradient

- | | |
|------------------|--|
| Visible | Choose this option to display a gradient in the area behind the chart. |
| No Middle | Choose this option to display a two color gradient (start and end color) in the area behind the chart. If you do not choose this option, you can specify three colors for the gradient (start, middle, and end). |
| Direction | Make a selection from the drop-down list to set the direction |

of the gradient.

Start/Mid/End Click to open a color palette. Choose a color for the start, middle, or end of the gradient.

Background

Chart Panel Click to open the color palette. Choose a color for the display area behind the chart.
Background-Color

Log scale on X-axis Choose this option to display a log scale on the x-axis for a standard curve only.



NOTE: The Gradient option overrides the Background option.

CHAPTER 8

This chapter explains the types of reports that are available and how to use the report generator. You can preview, print, or save a report.


8.1

Generating a Report

1. Open the plate (.csv, .xls, lxd or .mlx) for the report.



NOTE: The report generator is only available if a plate window is open and active.

2. Click the **Report Generator** button .
⇒ The report generator appears (Figure 8.1).




Figure 8.1 Report generator

3. To generate report for absolute quantification analysis, choose a report category and the **Absolute Quantification** option. To generate a report for relative quantification analysis, choose the report categories of interest and the **Relative Quantification** option.

Table 8.1 and Table 8.2 lists the available types of reports.

Table 8.1 MasterPlex™ QT reports for Absolute Quantification analysis

Report Name	Contents
Report per analytes	MFI, bead count, and concentration data organized by analyte name (Figure 8.7 on page 8.8)
Report per sample well	MFI, bead count, and concentration data organized by sample well (Figure 8.8 on page 8.9).
Plate format report	Organized by data types (MFI, bead count, or concentration) that are displayed in the well grid format (Figure 8.9 on page 8.10)
Report per unknown replicates	Considers each group of unknown wells a replicate. Includes the mean, standard deviation, and CV% for each unknown group.
Report per control replicates	Considers each group of control wells a replicate. Includes the mean, standard deviation, and CV% for each control group.
Report per background replicates	Includes the mean, standard deviation, and CV% for the background wells.
Report per standard replicates	Considers each group of standard wells a replicate. Includes the mean, standard deviation, and CV% for each standard group (Figure 8.11 on page 8.12).
Report per standard data	Standard data with expected and calculated values, residuals, and per cent recovery (Figure 8.12 on page 8.13).
Standard curves report	Local standard curves generated on the current plate (Figure 8.14 on page 8.15).
Tabular MFI report	Displays a table of MFI data per analyte (column) and well (row)
Tabular concentration report	Displays a table of concentration data per analyte (column) and well (row)
Well map report	A layout of the plate that shows the well type, the group number for each well, and the links between standards and unknown groups.

Table 8.2 MasterPlex™ QT reports for Relative Quantification analysis

Report Name	Contents
Report per analytes	MFI, bead count, and concentration data organized by analyte name (Figure 8.7 on page 8.8)
Report per sample well	MFI, bead count, and concentration data organized by sample well (Figure 8.8 on page 8.9).
Plate format report	Organized by data types (MFI, bead count, or concentration) that are displayed in the well grid format (Figure 8.9 on page 8.10)
Report per unknown replicates	Considers each group of unknown wells a replicate. Includes the mean, standard deviation, and CV% for each unknown group.
Report per control replicates	Considers each group of control wells a replicate. Includes the mean, standard deviation, and CV% for each control group.
Report per background replicates	Includes the mean, standard deviation, and CV% for the background wells.
Tabular MFI report	Displays a table of MFI data per analyte (column) and well (row).
Tabular Fold Change report	Displays a table of Fold Change data per analyte (column) and well (row).
Well map report	A layout of the plate that shows the well type, the group number for each well, and the links between standards and unknown groups.

4. To apply a log scale to the x-axis of standard curves in the Standard curves report, choose the **Log on X axis** option.
5. Click **Generate** to display the report.
 - ⇒ The report window opens and displays the report. The report name is added to the plate navigator under the report node for the plate.
6. To display a log scale x-axis in the Standard curves report, choose the **Log on X axis** option.
7. To print the Standard curves report on a monochromatic printer, choose the **Gray Scale** option.

8.2










Working with Reports

If the report window is open, the report toolbar is available (Figure 8.2). You can print report or save a report.



Figure 8.2 Report window toolbar

Table 8.2 Report window toolbar button functions

Toolbar Button	Click to...
	Open files
	Exit report window
	Change the magnification in the Report window. (Make a selection from the drop-down list of view options.)
	Display the first page in the report.
	Display the previous page.
	Display the next page.
	Display the last page.
	Open the Save As dialog box.
	Open the Print dialog box (Figure 8.6).

Saving a Report

1. Click the **Export Report** button .
⇒ The Save As dialog box appears (Figure 8.3).




Figure 8.3 Save As dialog box

2. Enter a file name
3. Choose a file type from the **Save as type** drop-down list.
A report may be saved in .frp, .csv, .txt, .emf, .wmf, .pdf, .bmp, .rtf, or .htm format.
4. Click **Save**.



NOTE: A report saved in the *.frp file format can only be opened in the MasterPlexTM QT application.

Viewing Reports

To view a report generated in the current session, click the report of interest in the Plate Navigator. To view a report saved during a previous session, use the report viewer.

1. To open a saved report, select **File > Report Viewer** from the menu bar.
⇒ The Report Viewer appears (Figure 8.4).




Figure 8.4 Report viewer

2. Click the **Open** button .
⇒ The Open dialog box appears (Figure 8.5).






Figure 8.5 Open dialog box

3. Double-click the report that you want to open.
⇒ The plate viewer displays the report.

4. To close the report viewer, click the **Close** button .

Printing a Report

1. Double-click the report that you want to print.
 2. Click the **Print Report** button .
- ⇒ The Print dialog box appears (Figure 8.6).




Figure 8.6 Print dialog box

3. Choose the print options of interest and click **OK**.

CHAPTER 8

MASTERPLEX QT REPORTS

MasterPlex QT Report Per Analytes	
Report Date: 5/8/2007	Run Date: 3/20/01
Report Time: 12:09:14 AM	Run Time: 12:45:55 PM
Data File: IL5.csv	Hardware Serial No.: LX10000280002
Plate Name: IL5 Project.mlx	Operator: AC
MasterPlex QT Version: 3.0.1.166	Analyst:

Analyte Name: 034 Mouse IL-5				Background: 0.00	
Well	Sample Name	MFI	Concentration	Unit	Count
A1	curve	20	9.30	pg/mL	63
B1	2	1308	998.69	pg/mL	95
C1	3	1084.5	500.28	pg/mL	72
D1	4	525	249.45	pg/mL	79
E1	5	201	127.80	pg/mL	81
F1	6	66.5	55.70	pg/mL	80
G1	7	39	33.36	pg/mL	78
H1	8	26.5	19.57	pg/mL	68
A2	Ind mouse serum	367	193.57	pg/mL	67
B2	10	429	215.85	pg/mL	66
C2		82	66.15	pg/mL	79
D2		88.5	70.26	pg/mL	66
E2		37.5	31.90	pg/mL	82
F2		32	26.17	pg/mL	79
G2		18	4.41	pg/mL	81
H2		19	7.11	pg/mL	69
A3	90-10 buffer	506.5	243.01	pg/mL	76
B3		480.5	233.95	pg/mL	78
C3		88	69.94	pg/mL	75
D3		82	66.15	pg/mL	76
E3		38	32.39	pg/mL	85
F3		32	26.17	pg/mL	66
G3		18	4.41	pg/mL	68
H3		17.5	2.59	pg/mL	80
A4	Sigma pooled	341	183.99	pg/mL	80
B4		392	202.64	pg/mL	79
C4		78	63.55	pg/mL	70
D4		96.5	75.12	pg/mL	74
E4		37	31.41	pg/mL	96
F4		34	28.33	pg/mL	74
G4		23	14.58	pg/mL	69
H4		27	20.22	pg/mL	42

Figure 8.7 Report by analyte
Data are organized by analyte name.

CHAPTER 8 MASTERPLEX QT REPORTS

MasterPlex QT Report Per Sample Well							
Report Date: 5/8/2007				Run Date: 3/20/01			
Report Time: 12:26:05 AM				Run Time: 12:45:55 PM			
Data File: IL5.csv				Hardware Serial No.: LX10000280002			
Plate Name: IL5 Project.mlx				Operator: AC			
MasterPlex QT Version: 3.0.1.166				Analyst:			
Sample Well : A1							
Analyte Name	Sample Name	MFI	Concentration	BKG	Unit	Count	
034 Mouse IL-5	curve	20	9.30	0.00	pg/mL	63	
038 Mouse IL-10	curve	13	<15.62	0.00	pg/mL	138	
054 Mouse IL-2	curve	28	<15.62	0.00	pg/mL	58	
073 Mouse IFN gamma	curve	102	<15.62	0.00	pg/mL	103	
077 Mouse IL-4	curve	68	13.22	0.00	pg/mL	112	
Sample Well : B1							
Analyte Name	Sample Name	MFI	Concentration	BKG	Unit	Count	
034 Mouse IL-5	2	1308	998.69	0.00	pg/mL	95	
038 Mouse IL-10	2	5476.5	1000.00	0.00	pg/mL	120	
054 Mouse IL-2	2	7855	1000.28	0.00	pg/mL	65	
073 Mouse IFN gamma	2	3141.5	1059.74	0.00	pg/mL	128	
077 Mouse IL-4	2	7792	999.79	0.00	pg/mL	92	
Sample Well : C1							
Analyte Name	Sample Name	MFI	Concentration	BKG	Unit	Count	
034 Mouse IL-5	3	1084.5	500.28	0.00	pg/mL	72	
038 Mouse IL-10	3	3242	500.01	0.00	pg/mL	125	
054 Mouse IL-2	3	5706.5	499.82	0.00	pg/mL	60	
073 Mouse IFN gamma	3	2515	470.38	0.00	pg/mL	115	
077 Mouse IL-4	3	5772	500.20	0.00	pg/mL	128	
Sample Well : D1							
Analyte Name	Sample Name	MFI	Concentration	BKG	Unit	Count	
034 Mouse IL-5	4	525	249.45	0.00	pg/mL	79	
038 Mouse IL-10	4	1332	249.97	0.00	pg/mL	125	
054 Mouse IL-2	4	2445	250.11	0.00	pg/mL	77	
073 Mouse IFN gamma	4	1707.5	259.86	0.00	pg/mL	136	
077 Mouse IL-4	4	2980	248.66	0.00	pg/mL	83	

Figure 8.8 Report per sample well

Data are organized by well location.

CHAPTER 8 MASTERPLEX QT REPORTS

MasterPlex QT Plate Format Report

Report Date: 5/8/2007

Report Time: 12:28:57 AM

Data File: IL5.csv

Plate Name: IL5 Project.mlx

MasterPlex QT Version: 3.0.1.166

Run Date: 3/20/01

Run Time: 12:45:55 PM

Hardware Serial No. : LX10000280002

Operator: AC

Analyst:

Data Type : MFI

Analyte : 034 Mouse IL-5

Background : 0.00

	1	2	3	4	5	6	7	8	9	10	11	12
A	20.00	367.00	806.50	341.00								
B	1308.00	429.00	480.50	382.00								
C	1084.50	82.00	88.00	78.00								
D	525.00	88.50	82.00	86.50								
E	201.00	37.50	38.00	37.00								
F	66.50	32.00	32.00	34.00								
G	39.00	18.00	18.00	23.00								
H	26.50	19.00	17.50	27.00								

Analyte : 038 Mouse IL-10

Background : 0.00

	1	2	3	4	5	6	7	8	9	10	11	12
A	13.00	710.50	1086.00	933.00								
B	5476.50	758.50	1028.00	1001.00								
C	3242.00	105.50	165.00	124.00								
D	1332.00	110.00	161.00	126.00								
E	402.00	32.50	40.00	33.00								
F	107.00	30.00	40.00	35.00								
G	41.00	13.00	12.00	12.00								
H	21.00	12.00	12.00	14.00								

Analyte : 054 Mouse IL-2

Background : 0.00

	1	2	3	4	5	6	7	8	9	10	11	12
A	28.00	1550.50	1729.50	1223.00								
B	7855.00	1617.50	1681.00	1286.00								
C	5706.50	106.00	100.00	70.00								
D	2445.00	111.00	95.00	74.00								
E	336.00	42.50	38.00	45.50								
F	57.00	50.00	40.50	47.00								
G	34.00	43.00	32.00	38.00								
H	33.00	46.00	37.00	35.50								

Analyte : 073 Mouse IFN gamma

Background : 0.00

	1	2	3	4	5	6	7	8	9	10	11	12
A	102.00	1185.00	1674.00	847.00								
B	3141.50	1357.00	1516.00	808.00								
C	2515.00	292.00	293.50	172.50								

IL5 Project

Page : 1 of 4

IL5 Project

IL5.csv

Page - 1 of 4

IL5 Project.mlx

Figure 8.9 Plate format report

Each data type is displayed in a well grid format.

CHAPTER 8 MASTERPLEX QT REPORTS

MasterPlex QT Report Per Unknown Replicates					
Report Date: 5/8/2007			Run Date: 3/20/01		
Report Time: 12:30:53 AM			Run Time: 12:45:55 PM		
Data File: IL5.csv			Hardware Serial No.: LX10000280002		
Plate Name: IL5 Project.mlx			Operator: AC		
MasterPlex QT Version: 3.0.1.166			Analyst:		
Unknown Samples					
034 Mouse IL-5			Background: 20.00		
Replicate Set : 0					
Well	Sample Name	MFI	Concentration	Unit	
A2	Ind mouse serum 10	347.0000	1935.7353	pg/mL	
B2		409.0000	2158.5012	pg/mL	
C2		62.0000	661.5248	pg/mL	
D2		68.5000	702.5527	pg/mL	
E2		17.5000	319.0262	pg/mL	
F2		12.0000	261.7078	pg/mL	
G2		-2.0000	44.0515	pg/mL	
H2	90-10 buffer	-1.0000	71.1064	pg/mL	
A3		486.5000	2430.0934	pg/mL	
B3		460.5000	2339.4646	pg/mL	
C3		68.0000	699.4469	pg/mL	
D3		62.0000	661.5248	pg/mL	
E3		18.0000	323.9322	pg/mL	
F3		12.0000	261.7078	pg/mL	
G3	Sigma pooled	-2.0000	44.0515	pg/mL	
H3		-2.5000	25.9035	pg/mL	
A4		321.0000	1839.8562	pg/mL	
B4		372.0000	2026.4169	pg/mL	
C4		58.0000	635.5244	pg/mL	
D4		76.5000	751.2116	pg/mL	
E4		17.0000	314.0754	pg/mL	
F4		14.0000	283.3361	pg/mL	
G4		3.0000	145.8495	pg/mL	
H4		7.0000	202.2244	pg/mL	
Mean		120.1667	797.4510		
Std. Dev		169.2687	819.1842		
CV%		140.8616	102.7253		

Figure 8.10 Report per unknown replicates

CHAPTER 8 MASTERPLEX QT REPORTS

MasterPlex QT Report Per Control Replicates					
Report Date: 5/8/2007			Run Date: 3/20/01		
Report Time: 12:32:30 AM			Run Time: 12:45:55 PM		
Data File: IL5.csv			Hardware Serial No.: LX10000280002		
Plate Name: IL5 Project.mlx			Operator: AC		
MasterPlex QT Version: 3.0.1.166			Analyst:		

Control Samples					
034 Mouse IL-5			Background: 20.00		
Replicate Set : 0					
Well	Sample Name	MFI	Concentration	Unit	
A2	Ind mouse serum 10	347.0000	1935.7353	pg/mL	
B2		409.0000	2158.5012	pg/mL	
C2		62.0000	661.5248	pg/mL	
D2		68.5000	702.5527	pg/mL	
E2		17.5000	319.0262	pg/mL	
F2		12.0000	261.7078	pg/mL	
G2		-2.0000	44.0515	pg/mL	
H2		-1.0000	71.1064	pg/mL	
Mean		114.1250	769.2757		
Std. Dev		165.8249	826.3844		
CV%		145.3012	107.4237		
Replicate Set : 1					
Well	Sample Name	MFI	Concentration	Unit	
A3	90-10 buffer	486.5000	2430.0934	pg/mL	
B3		460.5000	2339.4646	pg/mL	
C3		68.0000	699.4469	pg/mL	
D3		62.0000	661.5248	pg/mL	
E3		18.0000	323.9322	pg/mL	
F3		12.0000	261.7078	pg/mL	
G3		-2.0000	44.0515	pg/mL	
H3		-2.5000	25.9035	pg/mL	
Mean		137.8125	848.2656		
Std. Dev		208.9944	980.1522		
CV%		151.6512	115.5478		

Figure 8.11 Report per control replicates

CHAPTER 8
MASTERPLEX QT REPORTS

MasterPlex QT Report Per Standard Replicates				
Report Date: 5/8/2007		Run Date: 3/20/01		
Report Time: 12:35:53 AM		Run Time: 12:45:55 PM		
Data File: IL5.csv		Hardware Serial No.: LX10000280002		
Plate Name: IL5 Project.mlx		Operator: AC		
MasterPlex QT Version: 3.0.1.166		Analyst:		
Analyte: 034 Mouse IL-5			Background: 20.00	
Standard Data Set : 0				
Model: Five Parameter Logistics				
Replicates Set: 0				
Well	Sample Name	MFI	Concentration	Unit
B1	2	1288	9986.8756	pg/mL
Mean		NA	NA	
Std. Dev		NA	NA	
CV%		NA	NA	
Replicates Set: 1				
Well	Sample Name	MFI	Concentration	Unit
C1	3	1064	5002.7725	pg/mL
Mean		NA	NA	
Std. Dev		NA	NA	
CV%		NA	NA	
Replicates Set: 2				
Well	Sample Name	MFI	Concentration	Unit
D1	4	505	2494.4558	pg/mL
Mean		505.0000	2494.4558	
Std. Dev		505.0000	2494.4558	
CV%		100.0000	100.0000	
Replicates Set: 3				
Well	Sample Name	MFI	Concentration	Unit

Figure 8.12 Report per standard replicates

This partial view of the report shows data for the first three IL-10 standard replicates.

CHAPTER 8


MASTERPLEX QT REPORTS

MasterPlex QT Report Per Standard Data								
Report Date: 5/8/2007				Run Date: 3/20/01				
Report Time: 12:37:16 AM				Run Time: 12:45:55 PM				
Data File: IL5.csv				Hardware Serial No.: LX10000280002				
Plate Name: IL5 Project.mlx				Operator: AC				
MasterPlex QT Version: 3.0.1.166				Analyst:				
Analyte : 034 Mouse IL-5						Background: 20.00		
Standard Data Set: 0 Equation Model: Five Parameter Logistics								
Well	Sample Name	MFI	Calculated	Expected	Unit	Residuals	%Recovery	
H1	8	6.5000	195.7183	156.2500	pg/mL	39.4683	125.260	
G1	7	19.0000	333.6156	312.5000	pg/mL	21.1156	106.757	
F1	6	46.5000	556.9976	625.0000	pg/mL	-68.0024	89.120	
E1	5	181.0000	1277.9756	1250.0000	pg/mL	27.9756	102.238	
D1	4	505.0000	2494.4558	2500.0000	pg/mL	-5.5442	99.778	
C1	3	1064.5000	5002.7725	5000.0000	pg/mL	2.7725	100.055	
B1	2	1288.0000	9986.8758	10000.0000	pg/mL	-13.1244	99.869	
Analyte : 038 Mouse IL-10						Background: 13.00		
Standard Data Set: 0 Equation Model: Five Parameter Logistics								
Well	Sample Name	MFI	Calculated	Expected	Unit	Residuals	%Recovery	
H1	8	8.0000	124.9841	156.2500	pg/mL	-31.2659	79.990	
G1	7	28.0000	333.3337	312.5000	pg/mL	20.8337	106.667	
F1	6	94.0000	618.4314	625.0000	pg/mL	-6.5686	98.949	
E1	5	389.0000	1251.3058	1250.0000	pg/mL	1.3058	100.104	
D1	4	1319.0000	2499.7311	2500.0000	pg/mL	-0.2689	99.969	
C1	3	3229.0000	5000.0816	5000.0000	pg/mL	0.0816	100.002	
B1	2	5483.5000	9999.9688	10000.0000	pg/mL	-0.0312	100.000	
Analyte : 054 Mouse IL-2						Background: 28.00		
Standard Data Set: 0 Equation Model: Five Parameter Logistics								
Well	Sample Name	MFI	Calculated	Expected	Unit	Residuals	%Recovery	
H1	8	5.0000	<156.2500	156.2500	pg/mL	NA	NA	
G1	7	6.0000	<156.2500	312.5000	pg/mL	NA	NA	
F1	6	29.0000	743.5749	625.0000	pg/mL	118.5749	118.972	
E1	5	308.0000	1238.4097	1250.0000	pg/mL	-11.5903	99.073	
D1	4	2417.0000	2503.5060	2500.0000	pg/mL	3.5060	100.140	
IL5.csv			Page : 1 of 3			IL5 Project.mlx		

Figure 8.13 Report per standard data

MasterPlex QT 3.0 Standard Curve Report	
Report Date: 5/8/2007	Run Date: 3/20/01
Report Time: 12:38:02 AM	Run Time: 12:45:55 PM
Data File: IL5.csv	Hardware Serial No.: LX10000280002
Plate Name: IL5 Project.mlx	Operator: AC
MasterPlex QT Version: 3.0.1.166	Analyst:

Analyte : 034 Mouse IL-5



IL5.csv Page - 1 of 5 IL5 Project.mlx

Figure 8.14 Standard curves report
Shows the standard curve(s) for each analyte.

CHAPTER 8
MASTERPLEX QT REPORTS

MasterPlex QT Tabular MFI Report						
Report Date: 5/8/2007			Run Date: 3/20/01			
Report Time: 12:38:42 AM			Run Time: 12:45:55 PM			
Data File: IL5.csv			Hardware Serial No. : LX10000280002			
Plate Name: IL5 Project.mlx			Operator: AC			
MasterPlex QT Version: 3.0.1.166			Analyst:			

Well	Sample Name	034 Mouse IL-5	038 Mouse	054 Mouse IL-2	073 Mouse IFN	077 Mouse IL-4
A1	curve	0.00	0.00	0.00	0.00	0.00
B1	2	1288.00	5483.50	7827.00	3039.50	7724.00
C1	3	1084.50	3229.00	5678.50	2413.00	5704.00
D1	4	505.00	1319.00	2417.00	1605.50	2912.00
E1	5	181.00	389.00	308.00	488.00	1408.50
F1	6	48.50	94.00	29.00	75.00	516.00
G1	7	19.00	28.00	6.00	23.50	197.00
H1	8	6.50	8.00	5.00	8.50	55.00
A2	Ind mouse serum	347.00	897.50	1522.50	1083.00	1331.50
B2	10	409.00	745.50	1589.50	1255.00	1554.00
C2		62.00	92.50	78.00	190.00	332.00
D2		68.50	97.00	83.00	180.00	321.50
E2		17.50	19.50	14.50	38.00	90.00
F2		12.00	17.00	22.00	38.00	95.00
G2		-2.00	0.00	15.00	-1.00	6.50
H2		-1.00	-1.00	18.00	2.00	2.00
A3	90-10 buffer	488.50	1073.00	1701.50	1572.00	2653.50
B3		480.50	1015.00	1653.00	1414.00	2509.00
C3		68.00	152.00	72.00	191.50	717.00
D3		62.00	148.00	67.00	196.00	781.00
E3		18.00	27.00	10.00	22.00	191.50
F3		12.00	27.00	12.50	20.50	204.00
G3		-2.00	-1.00	4.00	7.00	1.00
H3		-2.50	-1.00	9.00	3.50	-1.00
A4	Sigma pooled	321.00	920.00	1195.00	745.00	1478.00
B4		372.00	988.00	1270.00	706.00	1604.00
C4		58.00	111.00	42.00	70.50	382.00
D4		76.50	113.00	46.00	85.00	382.00
E4		17.00	20.00	17.50	38.00	90.00
F4		14.00	22.00	19.00	33.00	93.00
G4		3.00	-1.00	11.00	23.00	-1.00
H4		7.00	1.00	7.50	18.00	8.00

Figure 8.15 Tabular MFI report

CHAPTER 8
MASTERPLEX QT REPORTS

MasterPlex QT Tabular Concentration Report						
Report Date: 5/8/2007			Run Date: 3/20/01			
Report Time: 12:39:22 AM			Run Time: 12:45:55 PM			
Data File: IL5.csv			Hardware Serial No. : LX10000280002			
Plate Name: IL5 Project.mlx			Operator: AC			
MasterPlex QT Version: 3.0.1.166			Analyst:			

Well	Sample Name	034 Mouse IL-5	038 Mouse	054 Mouse IL-2	073 Mouse IFN	077 Mouse IL-4
A1	curve	93.0342	<156.2500	<156.2500	<156.2500	132.1908
B1	2	9986.8756	9999.9688	10009.5370	10635.4325	9997.8745
C1	3	5002.7725	5000.0816	4991.8516	4895.0817	5002.0447
D1	4	2494.4558	2499.7311	2503.5060	2600.6457	2498.6258
E1	5	1277.9758	1251.3058	1238.4097	1208.1002	1295.9702
F1	6	558.9978	618.4314	743.5749	670.5801	573.9481
G1	7	333.6156	333.3337	<156.2500	453.2189	304.8798
H1	8	195.7183	124.9841	<156.2500	<156.2500	181.2247
A2	Ind mouse serum	1935.7353	1710.6917	2023.2823	1876.7282	1236.1352
B2	10	2158.5012	1775.7768	2059.3493	2089.7139	1413.3113
C2		661.5248	613.5848	915.2461	877.1919	419.8751
D2		702.5527	628.0033	927.0284	862.6805	410.9990
E2		319.0282	289.9925	617.0083	549.3872	212.0405
F2		261.7078	247.4404	698.5710	549.3872	216.4221
G2		44.0515	<156.2500	624.2880	<156.2500	138.0313
H2		71.1084	<156.2500	680.4154	<156.2500	133.9892
A3	90-10 buffer	2430.0934	2197.2950	2119.3512	2548.0257	2281.7300
B3		2339.4848	2124.7271	2093.4040	2306.2835	2167.5778
C3		699.4469	782.0496	900.3220	879.3375	739.7345
D3		661.5248	771.9113	887.1300	885.7284	792.0854
E3		323.9322	326.6574	504.1899	435.9786	300.1497
F3		261.7078	326.6574	581.7038	414.5085	310.8942
G3		44.0515	<156.2500	<156.2500	<156.2500	133.0901
H3		25.9035	<156.2500	417.7450	<156.2500	131.2908
A4	Sigma pooled	1839.8562	2004.2695	1843.4487	1500.5032	1352.9028
B4		2026.4189	2090.7112	1885.3335	1459.2362	1452.9992
C4		635.5244	670.7417	805.5859	659.0877	445.1800
D4		751.2116	676.6131	821.0448	694.4052	445.1800
E4		314.0754	274.2278	655.0945	549.3872	212.0405
F4		283.3361	290.4100	670.4287	523.8016	214.6700
G4		145.8495	<156.2500	543.4965	447.8207	131.2908
H4		202.2244	<156.2500	<156.2500	354.1332	139.3772

Figure 8.16 Tabular concentration report

CHAPTER 8 MASTERPLEX QT REPORTS

MasterPlex QT Well Map Report												
Report Date: 5/8/2007						Run Date: 3/20/01						
Report Time: 12:40:31 AM						Run Time: 12:45:55 PM						
Data File: IL5.csv						Hardware Serial No. : LX10000280002						
Plate Name: IL5 Project.mlx						Operator: AC						
MasterPlex QT Version: 3.0.1.166						Analyst:						
	1	2	3	4	5	6	7	8	9	10	11	12
A	Bkg-1	Unk-1	Unk-1	Unk-1								
B	Std-1	Unk-1	Unk-1	Unk-1								
C	Std-1	Unk-1	Unk-1	Unk-1								
D	Std-1	Unk-1	Unk-1	Unk-1								
E	Std-1	Unk-1	Unk-1	Unk-1								
F	Std-1	Unk-1	Unk-1	Unk-1								
G	Std-1	Unk-1	Unk-1	Unk-1								
H	Std-1	Unk-1	Unk-1	Unk-1								

Figure 8.17 Well map report

8.3

Custom Report

Custom report is a powerful and flexible tool for presenting and exporting your data. Compared to regular report, custom report has great flexibility on what and how to present data. As QT stores its analysis results in the format of XML document, it is possible for users to present their data in whatever format they want. The only thing users need to do is to define their presenting formats in XSL files (Extensible Stylesheet Language). QT will apply the stylesheets to transform and present data with certain format.

1. Start Report Generator. Choose “Custom report” as report type.
2. Click “View” button to bring out Style Sheet Manager(Figure 8.18).




Figure 8.18 Custom report viewer

3. Click “Import” to load a stylesheet into manager(Figure 8.19).




Figure 8.19 Import custom report stylesheet

4. Select a stylesheet from the manager window. The “XSL file details” window will display the information of the selected stylesheet, including name, created time, and output type (e.g., text, html, xml) (Figure 8.20).




Figure 8.20 Stylesheet information window

5. Click “Apply” to generate custom report, which is then displayed in a report window(Figure 8.21).

View Save Print												
Report Per Unknown Replicates												
GroupName	WellAddress	Sample Name	034 Mouse IL-5				038 Mouse IL-10				054 Mouse IL-12	
			Concentration	Mean Concentration	Stdev	CV	Concentration	Mean Concentration	Stdev	CV	Concentration	Mean Concentration
Unknown 1	A2	Ind mouse serum	1935.74	797.45	169.27	140.86	1710.69	987.88	389.71	148.89	2023.28	1056.8
	B2		2158.50				1775.78				2059.35	
	C2		661.52				613.58				915.25	
	D2		702.55				628.00				927.03	
	E2		319.03				269.99				617.01	
	F2		261.71				247.44				696.57	
	G2		44.05				NaN				624.29	
	H2		71.11				NaN				660.42	
	A3	90-10 buffer	2430.09				2197.29				2119.35	
	B3		2339.46				2124.73				2093.40	
	C3		699.45				782.05				900.32	
	D3		661.52				771.91				887.13	
	E3		323.93				326.66				504.19	
	F3		261.71				326.66				581.70	
	G3		44.05				NaN				NaN	
	H3		25.90				NaN				417.74	
	A4	Sigma pooled	1839.86				2004.27				1843.45	
	B4		2026.42				2090.71				1885.33	
	C4		635.52				670.74				805.60	
	D4		751.21				676.61				821.04	
	E4		314.08				274.23				655.09	
	F4		283.34				290.41				670.43	
	G4		145.85				NaN				543.50	
	H4		202.22				NaN				NaN	

Figure 8.21 Opened custom report

6. From the report window, user has options to save and print report, as well as to load report previously saved(Figure 8.22).





Figure 8.22 Save and Print icon for custom report



NOTE: When creating custom report in html/xml format, user has an option to export Standard Curve images in base64 string into the report. To convert the base64 string from the custom report back to images, user has to use some tool or utility to do the conversion.

APPENDIX A

Some of the software settings can be modified. This appendix explains the preferences (user-modifiable software settings) that are available to you. These settings are displayed in the Preferences dialog box.

To open the Preferences dialog box (Figure A.1), click the **Preferences** button . Alternatively, select **File > Preferences** from the menu bar.




Figure A.1 Preferences dialog box

Application tab default settings

A.1

Application Tab

The application preferences (Figure A.1) specify:

- the number of rows and columns for a virtual plate displayed in the plate wizard
 - whether the plate wizard appears when the MasterPlex™ QT software starts
1. To change the row and column dimensions of the well grid, enter the number of rows and columns for the well grid.
 2. If you do not want to display the plate wizard at program start up, remove the check mark from the **Enable plate wizard at start up** option.
 3. Click **Apply** when you are finished.
 4. To return the application preferences to the factory set defaults, click **Default**.

A.2

Plate Tab

The Plate preferences specify:

- how to compute background and when to subtract background
- threshold values for concentration, MFI, and bead count
- the plate and analyst name




Figure A.2 Preferences dialog box

Plate tab default settings



NOTE: If a plate (.csv, .xls, .lxd or .mlx) is not open, the Preferences dialog box does not display the Plate tab.

Background Options

You can specify whether you want to consider background data in the calculation of analyte concentrations. There are two methods of computing background-subtracted analyte concentrations: All Calculations or Regressions Only.

Subtract Background

Choose this option if you want to compute background-subtracted analyte concentrations. If this option is not chosen, the background MFI value is not considered during calculation of the analyte concentrations.

All Calculations	This method of computing background-subtracted analyte concentrations subtracts the background MFI from each member of the standard data set, then fits the standard curve. The method subtracts background MFI from unknown MFI, then interpolates the unknown analyte concentration from the standard curve.
Regressions Only	This method of computing background-subtracted analyte concentrations subtracts the background MFI from each member of the standard data set, then fits the standard curve. The method does not subtract the background MFI from the unknown MFI before interpolating the unknown analyte concentration.



NOTE: The **All Calculations** method is recommended. The **Regressions Only** method provides backward compatibility with data generated in MasterPlex™ QT 1.0.

Background Type

If there are three or more background wells in the assay, choose one of the following methods for computing background MFI.

Average Background (Bkg) MFI = (Bkg MFI₁ + Bkg MFI₂ + ... Bkg MFI_n)/n
where n = the number of background wells in the plate

Peak Value Highest background MFI value.

Lowest Value Lowest background MFI value.

Plate Info

Plate Name Displays the name assigned to the result file in the Luminex® 100/200 software. If you want to edit the plate name, enter a new name.

Analyst Name Displays the analyst name entered in the Luminex® 100/200 software. If you want to edit the analyst name, enter a new name.


Thresholds


You can enter an MFI, count, or concentration threshold for a plate.
The software can identify wells that contain data less than the user specified threshold.

To set a threshold(s):


1. Enter the MFI, count, or concentration threshold in the Preferences dialog box (Figure A.2).
2. Click **Apply** when you are finished.
3. To return the plate preferences to the factory set defaults, click **Default**.

To identify the wells in the Plate window that contain data less than threshold:

1. Make a selection from the analyte panel.
2. In the Plate window, select the data type (MFI, count, or concentration) from the data-type drop-down list.
3. To identify wells with MFI data less than threshold, click the  button.
⇒ A red border is placed around the well (Figure A.3).

To identify wells with bead count data less than threshold, click the  button.

⇒ A red border is placed around the well.

To identify wells with concentrations less than threshold, click the  button.

⇒ A red border is placed around the well.




Figure A.3 Plate window
Red border identifies wells with a MFI value less than the user-specified threshold for the selected analyte.

A.3
Calculations Tab

The Calculations preferences (Figure A.4) specify how to:

- fit the standard curve when there are replicate standard data sets
- display concentrations for diluted unknowns that were diluted (the diluted concentration or the original undiluted concentration)




Figure A.4 Preferences dialog box, Calculations tab

Standard Points

If the plate contains replicate standard data sets, there are two ways to fit a standard curve.

Average standards

The replicate standard data points are averaged and the standard curve is fitted to the single set of averaged data.

Choose this option if the experimental errors are not independent of each other

Individual points

Replicate standard data points are not averaged and the standard curve is fitted to all of the data points. For example, if there are three replicates of eight standard wells, the standard curve is fitted using all 24 data points. Choose this option if experimental errors are independent of each other.

Diluted Unknowns

Samples can be diluted prior to the assay and analysis. After MasterPlex QT interpolates the diluted unknown analyte concentrations from the standard curve, it can compute and display the original, undiluted concentration in the Plate window.

Original concentration = Diluted concentration * Dilution Factor

The dilutions factors are manually entered in the Plate window. (For more information see

4. Click **Apply** when you are finished.

To return the plate preferences to the factory set defaults, click **Default**.

APPENDIX **B**

The toolbars that are available depend on the types of windows that are open in the main display area.




B.1





Main Toolbar



Figure B.1 Main toolbar

Table B.1 Main toolbar buttons and functions

Menu Bar Command	Main Toolbar Button	Function
File → Open		Displays the Open dialog box so that a Luminex results file (.csv, .lxd), BioPlex file (.xls) or MasterPlex QT file (.mlx) may be opened.
File → Report Viewer	---	Displays the report viewer.
File → Close	---	Closes the current opened window.
File → Save		Displays the Save As dialog box so that a results file (.csv, .lxd or .xls) may be saved in MasterPlex QT file format (.mlx).
File → Save as		Displays the Save As dialog box so that a results file (.csv, .lxd or .xls) may be saved in MasterPlex QT file format (.mlx).
File → Preferences		Displays the Preferences dialog box.
Edit → Copy	---	Copies the data in the active window to the system clipboard.
Edit → Copy Horizontal List	---	Copies the data horizontally in the active window to the system clipboard.

Edit → Copy Vertical List	---	Copies the data vertically in the active window to the system clipboard.
Edit → Export Chart	---	Displays the Export Chart dialog box.
Charts → Chart Properties	---	Displays the Chart Properties dialog box. Displays the Plate Wizard.
Plate → Template Manager		Opens the Template Manager.
Show → Show/Hide → Navigator	---	Shows or hides the Plate Navigator.
---		Displays the Plate Wizard.
---		Displays the Plate Merge Wizard.
---		Displays the Auto Grouping Wizard.

B.2



Plate Toolbar

The Plate toolbar is available when a Plate window is open. To show or hide the plate toolbar, select **Show > Show/Hide > Plate Tools** from the menu bar.









Figure B.2 Plate toolbar

Table B.2 Plate toolbar buttons and functions

Menu Bar Command	Main Toolbar Button	Function
Plate → Selection Tool		Click this button to enable users to select wells in the Plate window using the mouse cursor.
Plate → Edit Wells		Click this button to enable users to edit user-modifiable data in the Plate window.

APPENDIX B

MASTERPLEX QT TOOLBARS

Plate → Virtual Pipette		Enables the virtual pipette that is used to transfer analyte data from a source plate to a virtual plate.
Plate → Analytes Filter		Opens the virtual analyte filter. Note: This button is only available when the virtual pipette is turned on.
Plate → Import Standards	---	When toggled on, this button enables the user to copy a standard data set from one plate to another using a drag-and-drop operation.
Plate → Data Charts		Opens a Data Chart window and displays the well chart of the selected data type for the selected well(s).
Plate → View Curves		Opens the Standard Curves window for the active Plate window. Displays the Report Generator dialog box.
Show → Show Well Types	---	Displays a well grid using colors to denote the different well types. Blue = standard well, green = unknown well gray = background well, orange = control well.
Show → Show Well Groups	---	In the well grid, displays a border around the wells in a group.
Plate → Show MFI < Limit	---	Displays a red border around wells that contain MFI data less than the user-specified threshold.
Plate → Show CON < Limit	---	Displays a red border around wells that contain MFI data less than the user-specified concentration threshold.
Plate → Show Count < Limit	---	Displays a red border around wells that contain count data less than the user-specified count threshold.
Plate → Subtract Background		Displays background-subtracted data in the Plate window. Displays the Print Settings dialog box.
Plate → Print Plate		Displays a print preview of the well grid for the selected analyte in the active plate window.

B.4



Calculation Toolbar

The calculation toolbar is available when a Plate window is open. To show or hide the calculation toolbar, select **Show > Show/Hide > Calculation Tools** from the menu bar.



Figure B.4 Calculation toolbar

Table B.4 Calculation toolbar buttons and functions

Menu Bar Command	Main Toolbar Button	Function
Calculation → Model Equations		Displays the Model Equations dialog box so that you can select a model for a standard or regression analysis data set.
Calculate → Calculate STD Curves, Fit Regression Curves		Generates the standard curves for the selected standard data set using the selected model and interpolates analyte concentrations from the standard curves. Generates a Dose-Response curve for the selected regression analysis data set using the selected model and determines LogEC50.

B.5






Chart Toolbar

The chart toolbar is available when a Data Chart or Standard Curves window is open. To show or hide the chart toolbar, select **Show > Show/ Hide > Chart Tools** from the menu bar.



Figure B.5 Chart toolbar

Table B.5 Chart toolbar buttons and functions

Menu Bar Command	Main Toolbar Button	Function
Charts → Select Charts		Displays a drop-down list of chart types.
Charts → Chart 3D		Toggles the display of a chart between a two or three dimensional view.
Charts → Rotate Chart		Enables you to rotate a chart using a click-and-drag operation.
Charts → Zoom Chart		Enables you to magnify a user selected area in a chart using a click-and-drag operation. Displays or hides the analyte panel in the Data Chart window. Displays or hides the scale bar in the Data Chart window.
Charts → Print Chart		Opens the Chart Print Preview dialog box.

B.6










Report Toolbar

The report toolbar is available when the Report window is open.



Figure B.6 Report toolbar

Table B.6 Report toolbar buttons and functions

Menu Bar Command	Main Toolbar Button	Function
---		Open saved reports.
---		Close current report window.
---		Displays a drop-down list of magnification options for viewing the report in the Report window.
---		Displays the first page in the report.
---		Displays the previous page in the report.
---		Displays the next page in the report.
---		Displays the last page in the report.
Report → Save Report		Displays the Save As dialog box so that you may save a report (.frp).
Report → Print Report		Displays the Print dialog box so that you may print a report.

APPENDIX **C**

This appendix provides background on the four and five parameter logistic curves. It also explains how weighting methods can improve the fit of nonlinear models to data with non-constant variability (heteroscedasticity).

C.1

Four Parameter Logistic Curve




Figure C.1 Four parameter logistics curve

In Figure C.1, the asymptotes D and A are the upper and lower limits of the model. An asymptote is a value that the function never reaches. Therefore, the value of a function at or beyond an asymptote cannot be

predicted. If $MFI \leq A$ or $MFI \leq D$, the MFI value is out of the calculable range of the four parameter logistics model. It is not possible to mathematically extrapolate values that are equal to or beyond the asymptote values.

For example, imagine the function $F(x) = \text{Log}(x)$. The vertical asymptote for this function is the line $x = 0$. This means the function can never reach the line where $x = 0$. Try calculating $\text{Log}(0)$ to see what happens.

In Figure C.1, B is the slope at the inflection point. It is the speed of the function as it moves away from the inflection point. C is the most interesting parameter since it corresponds to the x value that is associated with the y value at the midpoint between the minimum and maximum limits of the function. In a biochemical assay, C corresponds to the concentration (since it is on the x-axis) that produces a 50% response. In drug discovery terms, it is also known as ED50, the 50% effective dosage. It can also be interpreted as LD50 if this dosage kills 50% of the sample population.

C.2

Five Parameter Logistic Curve




Figure C.2 Five parameter logistic curve

A five parameter logistic curve introduces an additional parameter to the four parameter logistic model to compensate for asymmetric data. In Figure C.2, we can see that the curve is not symmetric because the bottom part of the curve behaves quite differently from the top portion. The additional parameter E compensates for this behavior and corrects for asymmetry.

Fixing Parameter A to Zero

If an assay contains background or blank wells, readings from these wells can be considered noise in the detection. By deducting this noise from the MFI values we can obtain more accurate responses. If we subtract a reasonable background value, then we can further assume that the lowest detectable value should now be zero. Since parameter A represents the lowest detectable concentration value, it can be set to zero.




Figure C.1 Five parameter logistics curve

Log-Log Model

The Log-Log model transforms the data to log scale for the x and y values. It applies linear regression to fit a straight line through these points. This model is appropriate for data that are intrinsically linear.

C.3

Heteroscedasticity

Fitting nonlinear models to observed data is often complicated by non-constant or heterogeneous variability. Heterogeneous variability or *heteroscedasticity* occurs in most types of observed data. This is especially true for biochemical assays where concentration or dose is the predictor and response is often based on count. Therefore, we can expect that measurement error varies with respect to the mean. In the Luminex 100/200 system, MFI (median fluorescence intensity) values are based on bead counts and vary with the concentration. In this case, we expect the error in detecting MFI values to increase as concentration increases. This is best seen in Figure C.3, a residual plot from a Luminex 100/200 cytokines assay.




Figure C.3 Residual plot

A residual plot is a graphical representation of how far away an observed concentration is from its expected value. It plots residuals against observed concentrations. In Figure C.3, we can see that the deviation of the observed concentration from the expected value increases as concentration increases. This means the variability is not constant.

Residual plots help you detect non-constant variability as well as outliers. If a residual plot exhibits data points in a wedge or funnel-shaped pattern, then we can expect the underlying data to have non-constant variability.

Non-constant variability complicates curve fitting because the regression process assumes the errors are constant across all data points. When the data violate this assumption, the resulting curve fit is less than optimal. This is illustrated in Figure C.4.




Figure C.4 Data with non-constant variability

The residual increases as concentration increases.

When we fit a model to data, the curve is applied to all of the data points as closely as possible so that the distances between the predicted and expected concentrations are minimized. In Figure C.4 we can think of the lines that

represent the residuals as ropes and each data point as a wrench. Curve fitting can be imagined as pulling the curve line so that it is as close to each point as possible without snapping the actual curve.

The best curve fit is reached when the curve is pulled as close as possible to each data point without breaking the actual curve model.

The nonlinear least square algorithm accomplishes this task. The nonlinear (or linear) least square algorithm assumes that all points have the same variability, so all points influence the curve fit equally.

However, data that exhibit non-constant variability violate this assumption. As a result, data points with greater variability assert more ‘pull’ on the curve. Data points at higher concentrations have more variability than those at lower concentrations, and have greater influence on the curve fit than the points at a lower concentration. As a result, accuracy or sensitivity at the lower concentrations decreases.

C.4

Weighted Nonlinear Least Square

The weighted nonlinear least square method of curve fitting is one way to correct for non-constant variability. In this method, weights are assigned to each point so that all points have equal influence on the curve. Instead of minimizing the residuals, the method minimizes residuals based on the weight at each point. In mathematical terms, the non-constant variability is made constant again by these additional weighting factors. If the weight for a point is higher, it will influence the curve fit more. The weight is the inverse of the variance, so points with low variability have more influence during curve fitting, which seems logical. As a result, the curve fit represents the data better and the sensitivity often increases at lower concentrations.

Weighting By $1/Y^2$

We almost never know what the weight should be at each point, so it must be estimated. One of the most widely used estimation methods is based on a working hypothesis that variance (variability) depends on the square of the mean (the y value). This implies that the standard deviation of a detected MFI value for a particular concentration is directly proportional to the expected

concentration. In other words, standard deviation is proportional to the mean.

$$\sigma = \mu$$

$$\sigma = [\text{Some Constant}] * \mu$$

$$\sigma / \mu = [\text{Some Constant}]$$

where σ = standard deviation and μ = mean

Since standard deviation/mean is the coefficient of variation (CV), we can now expect CV to be constant for all points. In this scheme, weight (W) can be estimated using the following expressions:

$$W = 1 / \text{Variance}$$

$$W = 1 / \sigma^2$$

where weight is the inverse of the variance and variance = standard deviation²

Since we assume, $\sigma \propto \mu$:

$$W \cong 1 / Y^2$$

(Note that μ = Mean = $F(x)$ = Y)

where W is the estimated weight and Y = F(x) and F(x) is the mean function which is represented by the equation of the model.

In this scheme, we can think of minimizing the residuals based on their relative MFI values. This method is appropriate when higher MFI values have larger residuals than lower MFI values. The wedge shape pattern of data on a residual plot usually confirms this.

Weighting by 1/Y

Another weighting scheme is 1/Y. This method is useful when the MFI values follow a Poisson distribution. In a Poisson distribution, the standard error of the MFI values is equal to the square root of that value.

Since weights are inverses of variances we have:

$$W = 1 / \text{Variance}$$

$$W = 1 / \sigma^2$$

$$W = 1 / (\sqrt{Y})^2$$

$$W \cong 1 / Y$$

where:

weight (W) is the inverse of variance

variance = (standard deviation)²

$\sigma = \sqrt{Y}$

$Y = \text{MFI}$

This scheme works well with most Luminex 100/200 assays. It is a compromise between minimizing absolute residuals and minimizing residuals based on their relative intensity values. If the $1/Y^2$ scheme is too strong for a particular data set, this scheme should be tried even if the underlying data distribution is not understood.

Weighting by $1/X$ and $1/X^2$

These schemes are rarely used. They give more weight to a data point that is closer to the left side of the graph. A data point on the right side of the graph receives less weight.

Weighting by $1/\text{Stdev}^2$

If the distribution of MFI values for a particular assay type is understood, we can derive standard errors for the values. In this case, the weights can be calculated exactly as $1/\text{Stdev}^2$. You can enter the standard deviations into the well grid under the "Standard Deviation" data type and MasterPlex™ QT will automatically calculate the weights.

C.5

Results of Weighting

In addition to the R-Square value, MasterPlex™ QT computes % Recovery values that measure how well the calibration curve (standard curve) fits the observed data.

% Recovery = Observed Concentration/Expected Concentration) * 100

The closer % Recovery is to 100, the better the curve fit at that point. If % Recovery is less than 100, the point is below the curve. If it is greater than 100, the point is above the curve.

Figure C.5 and Figure C.6 show the % Recovery for two different standard curve

fits. The standard data in Figure C.6 is from a four parameter logistic equation with no weighting. Note that the % Recovery for H1 is poor. This means the sensitivity deteriorates significantly at lower concentrations.

The data in Figure C.6 is from the same equation with 1/Y weighting. After we apply weighting, the % Recovery for G1 and H1 improves significantly.

Standard Curve Chart				Residual Plot						Standard Data			
	Well	Sample Name	Outliers	MF1	Calculated	Expected	Residuals	% Recovery	Background				
7	H1	8	<input type="checkbox"/>	55.00	2.12	29.08	-26.96	7.30	68.00				
6	G1	7	<input type="checkbox"/>	197.00	84.81	87.24	-2.44	97.21	68.00				
5	F1	6	<input type="checkbox"/>	516.00	273.63	261.73	11.90	104.55	68.00				
4	E1	5	<input type="checkbox"/>	1406.50	864.88	785.19	79.70	110.15	68.00				
3	D1	4	<input type="checkbox"/>	2912.00	2199.04	2355.56	-156.52	93.36	68.00				
2	C1	3	<input type="checkbox"/>	5704.00	7345.53	7066.67	278.86	103.95	68.00				
1	B1	2	<input type="checkbox"/>	7724.00	20793.86	21200.00	-406.14	98.08	68.00				

Figure C.5 Example standard curve data

Four parameter logistics equation fitted to standard data with no weighting.

Standard Curve Chart				Residual Plot						Standard Data			
	Well	Sample Name	Outliers	MF1	Calculated	Expected	Residuals	% Recovery	Background				
7	H1	8	<input type="checkbox"/>	55.00	27.92	29.08	-1.16	96.00	68.00				
6	G1	7	<input type="checkbox"/>	197.00	92.56	87.24	5.32	106.10	68.00				
5	F1	6	<input type="checkbox"/>	516.00	256.86	261.73	-4.87	98.14	68.00				
4	E1	5	<input type="checkbox"/>	1406.50	822.09	785.19	36.90	104.70	68.00				
3	D1	4	<input type="checkbox"/>	2912.00	2192.21	2355.56	-163.35	93.07	68.00				
2	C1	3	<input type="checkbox"/>	5704.00	7625.68	7066.67	559.01	107.91	68.00				
1	B1	2	<input type="checkbox"/>	7724.00	20322.65	21200.00	-877.35	95.86	68.00				

Figure C.6 Example standard curve data

Fitting the four parameter logistics equation to standard data using 1/Y weighting improves the % Recovery at lower concentrations.

APPENDIX

D

This appendix provides an introduction to Dose-Response analysis.

Dose-Response Analysis

MFI values are raw response values measured by the Luminex® 100/200.

Based on these measured responses, standard curves are used to approximate the actual analyte concentration in a sample. Analyte concentrations interpolated from a standard curve are quantified responses. Quantified responses can be used to study dependent variables inherent to a biological or chemical process.

For example, suppose a chemical agent that is given to mice causes an increase in the level of certain cytokines. The chemical agent dosage is the independent variable, and the cytokine concentration is the dependent variable. To understand the potency of the agent, you can perform a Dose-Response analysis. The Dose-Response analysis provides information on the relationship between the agent and changes in the level of cytokines.

Dosages

To conduct a Dose-Response analysis, different dosages of the chemical, biological, or other agent are prepared. For example, eight serial dilutions of an agent (mg/ml) could be prepared and given to eight different mice.

A blood sample taken from each mouse is then used to perform a cytokine multiplex assay on the Luminex 100/200.

Responses

Using MasterPlex QT, you can interpolate the cytokine concentrations of the blood samples from the standard curves of the reagent kit. Now you have both dose and quantified response data (mg/ml).

Regression Analysis

The software uses regression analysis to fit a dose response model to the data points: x = dose and y = response. After the model is fitted, the software provides the LogEC50. Take the antilog of the logEC50 to determine the EC50, the agent dosage that elicits 50% of the response. If 50% of the test subjects die

at the EC50 value, the EC50 value is the LD50 value. In this way, you can determine the effectiveness or the toxicity of an agent.

APPENDIX E

This appendix provides an introduction to Relative Quantification.

Relative Quantification Assay

Researchers sometimes are interested in studying the quantified relationships between groups. Instead of calculating the absolute concentration values, they may want to know the response ratio (Fold Change) between Control and Treatment groups. The following demo will guide your through setting up and analyzing a single QuantigenePlex plate.

1. First, we will begin by designating the bead regions (House Keeper) that you want to normalize the data with. Right click on a bead region, and select housekeeping from the drop-down menu.
2. Now we will start to group our wells according to their sample type. There are basically three types of the groups in relative quantification: Background, Control, and Treatment. There is at most one background group, while the other two types may have multiple groups. Each treatment group is linked to one control group.
3. Now we are ready to calculate the Fold Change. Click “Calculation” button from the toolbar to open the calculation dialog. On the Calculation dialog, there are three options for relative quantification analysis.
 - 1) Fold Change without normalization – this will only calculate the MFI ratio between the Control and Treatment groups without any data normalization.
 - 2) Fold Change with normalization – this will perform a normalization step first, and then calculate the Fold Change between groups based on the normalized data.
 - 3) Normalization only – this will only perform normalization on plate in one of the three ways.
 - Normalize by the geometric mean of several house keeper genes.
 - Normalize by a single house keeper gene.
 - Normalize by a constant value.

5. For this demo, we will choose Fold Change with Normalization, based on our previously selected house keeper gene.
6. Click “Ok”, and a popup window will display the calculation status.
7. To view the normalized data, select “Normalized data” from the Data Type pull-down menu. To view the final Fold Change results, select “Concentration/Fold Change” from the drop-down menu. To view results for a different bead region, just click its name from the bead set window.