

mobilegs™



Upgrade your **Mobilegs** at: www.mobilegs.com

[facebook.com/mobilegs](https://www.facebook.com/mobilegs)

twitter.com/mobilegs


CAUTION: Please inspect your new Mobilegs for visible defects or damage prior to each use. If present, do not use. Please contact us at 612.605.5898 or info@mobilegs.com to have your product replaced where manufacturing or factory shipping damage is present.

PN:1001 (12.01.11)



mobilegs™

user manual & upgrades



mobi™

perfecting
human mobility

“Our passion is to redefine mobility and bring assisted mobility products to the market that instill dignity and build self-esteem through disruptive design.”

Upgrades & Accessories

- **Style your Mobilegs with Mobiskins**
- **Replacement feet**
- **Replacement handgrips**

Visit Mobilegs.com for more information.



Warning and Instructions:

IMPORTANT RULES FOR SAFETY AND COMFORT

- ★ Follow your physician's or licensed medical professionals instructions for **Mobil**egs use and weight bearing on your affected limb. Call your physician if you have any questions or difficulties. If need help locating a professional in your area, please contact **Mobi**, LLC.
- 👉 Don't look down. Look straight ahead as you normally do when you walk.
- 👉 Don't use **Mobil**egs if you feel dizzy or drowsy.
- 👉 Don't walk on slippery surfaces. Avoid snowy, icy, or rainy conditions.
- 👉 Maximum weight 300 lbs or 136 kg.

WWW.MOBILEGS.COM

👁️ Check for and replace **worn** parts such as Feet prior to each use.

👁️ Verify **Mobilegs** are in good condition prior to use. If you discover a problem with your **Mobilegs**, please discontinue use and notify Mobi, LLC.

👁️ Wear well-fitting, low-heeled shoes (or shoe).

👁️ Position the handgrips correctly and take frequent rests to release pressure.

✖️ Get your **Mobilegs** properly fitted.

CAUTION: Do not bear weight on the saddle through your armpit(s).

CAUTION: Patients with heart disease should use **Mobilegs** only under the supervision of a physician, and with Partial-Weight Bearing gait, and only for short distances with frequent rests.

Ultra Fit Guide



1 Locate left and right-handed **Mobile**legs

2 Squeeze push buttons to adjust

Set arm length



- 3
1. Hold onto grip
 2. Place saddle snug underarm, keeping shoulders relaxed
 3. Arm comfortably bent at 15°-30° angle
 4. Leave a slight space between armpit and Saddle
 5. Adjust to fit

Set leg length

- 4
1. Adjust lower section to 6 inches in front of toe when leg is straight in standing position

5 Confirm final fit

When necessary, please seek assistance to ensure your safety and comfort during the fitting process.



Scan with your smartphone





www.mobilelegs.com
www.facebook.com/mobilelegs

Universal Fitting



- 1 Locate left and right-handed **Mobilegs**
- 2 Squeeze push buttons to adjust
- 3 Set arm length
 1. Hold onto grip
 2. Arm slightly bent at 15°-30° angle
 3. Leave a two finger space between armpit and Saddle
 4. Do not place Saddle tight into armpit
 5. Adjust to fit
- 4 Set leg length
 1. Adjust lower section to 6 inches in front of toe when leg is straight in standing position
- 5 Confirm final fit

When necessary, please seek assistance to ensure your safety and comfort during the fitting process.



Scan with your smartphone



www.mobilegs.com
www.facebook.com/mobilegs

Instructions:

**PLEASE READ ALL
INSTRUCTIONS
CAREFULLY**



Orientation

Mobilegs are specific to the left and right side of the body. They are not interchangeable. Use Mobilegs only as specified. In the correct orientation the handgrips face to the rear and point slightly outward. The bend in the tube will face forward and will naturally curve slightly away from the hips.




The tripod position

The tripod position is the position in which you stand when using Mobilegs. It is also the position in which you begin walking.

TO GET INTO THE TRIPOD POSITION

Place the crutch tips about 4" to 6" to the side and front of each foot. Then stand on your "good" foot (the one that is weight-bearing).

IMPORTANT RULES FOR SAFETY AND COMFORT



Walking with your Mobilelegs

(non-weight bearing)

If your physician has told you to avoid all weight bearing on the affected limb, it is important to follow these instructions carefully. You will need sufficient upper body strength to support your weight with just your arms and shoulders.

TO WALK WITH MOBILEGS:

- 1.** Begin in the tripod position—and remember, keep all of your weight on your “good” (weight-bearing) foot.
- 2.** Lean forward to create momentum.
- 3.** Move the “good” weight-bearing foot/leg forward (beyond the Mobilegs).
- 4.** Advance both Mobilegs and then the affected foot/leg.
- 5.** Repeat.

Product Warranty and Liability

Mobi, LLC warrants its Products for the original purchaser to be free from defects in materials and workmanship for a period of 12 months after the original purchase date. If you believe the product you purchased is defective within the warranty period, return the product to Mobi, LLC or the place of purchase. The product must be received by Company within 12 months after the original purchase date. The product will be inspected by the Company and, if found defective will, at Company's option, be repaired or replaced. If warranty coverage is not applicable, the product will be returned to you at your expense. The warranty is voided if the defect is caused by improper use or repair or if the defect is due to general abuse, neglect, accident, modification, tampering or any other reason not related to defects in materials or workmanship. Please visit mobilegs.com for details on warranty returns.

This warranty does not extend to non-durable components of the product, including but not limited to feet/tips, foam hand grips, foam saddle padding, or other parts subject which are prone to wear under normal use.

PRECAUTIONARY NOTE: You should always consult with a medical professional before using any crutch or other ambulatory aide.

THE SOLE AND EXCLUSIVE REMEDY FOR BREACH OF WARRANTY IS, AT COMPANY'S OPTION, PRODUCT REPLACEMENT OR REPAIR. NO WARRANTY OTHER THAN THIS LIMITED WARRANTY IS MADE AND THE PRODUCT IS OTHERWISE SOLD ON AN "AS-IS" BASIS WITH NO OTHER EXPRESS OR IMPLIED WARRANTIES. ANY IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE AND ANY IMPLIED WARRANTIES ARISING FROM COURSE OF PERFORMANCE, COURSE OF DEALING OR USAGE OF TRADE ARE EXPRESSLY DISCLAIMED.

COMPANY SHALL NOT BE LIABLE UNDER ANY CIRCUMSTANCES FOR ANY INCIDENTAL, CONSEQUENTIAL OR OTHER INDIRECT DAMAGES, INCLUDING WITHOUT LIMITATION, DAMAGES OR INJURY CAUSED BY DEFECT, FAILURE OR OTHERWISE, WHETHER ARISING UNDER WARRANTY, CONTRACT, NEGLIGENCE OR OTHERWISE. IN NO EVENT SHALL COMPANY'S LIABILITY UNDER ANY CAUSE OF ACTION RELATING TO ANY PRODUCT EXCEED THE PURCHASE PRICE OF SUCH PRODUCT. YOUR EXCLUSIVE REMEDY REGARDING ANY GOODS PROVIDED BY US SHALL BE THE REPLACEMENT OR REPAIR OF ANY SUCH GOOD FOUND TO BE DEFECTIVE. BECAUSE SOME STATES OR JURISDICTIONS DO NOT ALLOW THE EXCLUSION OR LIMITATION OF LIABILITY FOR CONSEQUENTIAL OR INCIDENTAL DAMAGES, IN SUCH STATE OR JURISDICTIONS, OUR LIABILITY, AND THE LIABILITY OF OUR SUPPLIERS, SHALL BE LIMITED TO THE FULLEST EXTENT PERMITTED BY LAW.

Features + Benefits

- **Collapsible** for easy storage
- **1 Size** for adults (up to 300lbs) - 4'9" - 6'4"
- **Swivel-Pivot saddle** that moves with you
- **Lightweight** aluminum frame
- **Push button height adjustment**
- Padded **contoured handgrips**
- **Sold as a pair**
- **Non-skid rocker feet** for traction on slippery surfaces

Mobilegs are designed to cause less stress on the body than typical crutches because of patent pending contoured frame design, superior saddles, smaller profile, and light-weight construction.

High impact, lightweight anodized aluminum frame angles in towards the body to provide vertical stability support while keeping a **slim profile** to aid navigation through tight spaces.

US Patent 7,717,123,, 7,926,498, 7,926,499

US and Foreign Patents Pending

Mobilegs is a trademark of Mobi, LLC

Register Your Mobilegs

Please register your product at <http://www.mobilegs.com/register.html>

Comments & Questions

Do you have a product question or want to share your experience using your **Mobilegs**? Your feedback matters to us and will help make **Mobilegs** even better.

www.mobilegs.com
info@mobilegs.com

3500 American Blvd West Suite 640
Bloomington, MN 55431
Phone: 952-562-5580 Fax: 952-562-5581

Join us at facebook and twitter:

www.facebook.com/mobilegs
www.twitter.com/mobilegs