

FoodTech
Calibration Software
for Windows™ 98 and up

Version 1.5

COLE-PARMER INSTRUMENT CO.

625 East Bunker Court
Vernon Hills, IL 60061

Phone: (800)323-4340
Fax: (847)247-2929
E-mail: info@coleparmer.com
Website: <http://www.coleparmer.com>

This manual provides you with all the information to operate the Cole-Parmer FoodTech CP-50

Table of Contents

SECTION 1: GENERAL INFORMATION3
 SYSTEM REQUIREMENTS3
 INSTALLING THE PROGRAM.....3
 RUNNING THE PROGRAM3
SECTION 2: CALIBRATING A FOODTECH INSTRUMENT.....5
 BEFORE BEGINNING6
 STEP 1: CREATE AN OPTICAL DATA FILE AND COLLECT DATA FROM THE ANALYZER.....7
 STEP 2: PLOT CURVES.....9
 STEP 3: OBTAIN LAB ANALYSIS FOR ALL CALIBRATION SAMPLES.....10
 STEP 4: CREATE OR EDIT A CONSTITUENT DATA FILE10
 STEP 5: MERGE CONSTITUENT DATA FILE TO OPTICAL DATA FILE11
 STEP 6: FOODTECH REGRESSION ANALYSIS11
 STEP 7: ENTER FINAL K'S.....13
 STEP 8: CHECK ACCURACY16

Please read this entire manual before beginning a calibration.

Examples are given in this manual for some of the FoodTech analyzers' most common applications. If you have questions about the recommended settings for a product that is not mentioned, contact Cole-Parmer.

Phone: (800)323-4340
Fax: (847)247-2929
E-mail: info@coleparmer.com

A note on the text: Throughout this manual, the Enter/Return key on the computer keyboard is referred to as **Return** to distinguish it from the **Enter** button on the FoodTech analyzer's keypad.

Section 1: General Information

System Requirements

- Windows™ 98 or above
- VGA Color monitor
- 1 com port
- 3-1/2 inch disk drive
- 15 megabytes of hard drive space
- printer

Installing the Program

Before installing the program, make sure you have closed all open Windows™ applications. From the **Start** menu, select **Run**. Insert Disk 1 of the Calibration Software Program into your 3-1/2 inch disk drive. On the **Command Line**, type **a:** or **b:** (as appropriate) **setup.exe**. Press **Return** or select **OK**. The computer will walk you through the setup process, which requires little more than waiting patiently to insert Disk 2.

After installing the program, store the disks in a safe place in case of hard-drive failure or accidental erasure.

Running the Program

Select FoodTech Calibration Software from the Programs menu. The following screen will appear:



As with other Windows™ applications, you can access the menu bar with the mouse or by pressing **Alt** and the underlined letter simultaneously. Shortcut keystrokes for some of the options in each menu are listed to the right of the command when the menu appears. The toolbar of shortcut buttons directly below the menu bar provides quick access to some of the most commonly used commands. The right mouse button accesses the **Data** menu when an optical data file is open.


File names should be no more than 8 characters long and consist only of letters and numbers. The computer automatically assigns the following extensions to the files and stores them in the DATA directory:

- “.dt1” to optical data files
- “.st1” to statistics files
- “.con” to constituent data files
- “.ks1” to K values files
- “.zbk” to instrument settings files
-

Any optical data file can be viewed in two ways: **Plot View** or **Grid View**. The grid view option is designed mainly for use by Cole-Parmer researchers. The plot view is the most useful as it shows the curves of the sample readings. Whenever multiple samples are graphed in plot view, you can click on any curve or point on the graph to determine the sample's ID number.

Section 2: Calibrating a FoodTech Instrument

After you become familiar with the process of calibrating a FoodTech instrument you will be able to alter the order in which you perform the following tasks, but beginners should follow this basic flowchart. The additional information you will need to perform each step is given later in this section.



Before Beginning

Before beginning a calibration, you should have 100 to 300 samples of the product for which you will be calibrating. Label these samples with FoodTech sample ID numbers from the Calibration Numbering Setup Manual.


CUPS AVAILABLE for the FoodTech

The sliding window cup for the FoodTech includes the following:

Item/Stock #	Quantity
Cup Frame/2672	1
Window/A0420-1	2
Window/A0420-2	2
Window/A0420-3	2
Window/A0420-4	2
Window/A0420-5	1

The combination of these windows in cup frame #2672 produces a cup with a sample thickness as indicated in the following table (rounded to nearest mm).

Shaded boxes indicate preferred window combinations.

Window #	A0420-1	A0420-2	A0420-3	A0420-4	A0420-5
A0420-1	15	14	13	12	8
A0420-2	14	13	12	11	7
A0420-3	13	12	10	9	6
A0420-4	12	11	9	8	5
A0420-5	8	7	6	5	

The combination of these windows in cup frame #A0612 produces a cup for liquid and slurries.

Window #	A0612-1	A0612-2	A0612-3	A0612-4
A0612-1	21			
A0612-2	18	15		
A0612-3	17	14	13	
A0612-4	13	10	9	5

See next page for products that require the indicated cup size.


	Stock #
Wheat, Barley	#2670 Standard Black molded cup (16mm)
Corn, Soybean	#A0562 Black molded cup w/windows on outside (24mm)
Milo/Sorghum	Use the chart on previous page to make a 13mm cup with the sliding windows.
Flax, Mustard, Soymeal	Use the chart on previous page to make a 10mm cup with the sliding windows.
Canola, Rapeseed	Use the chart on previous page to make an 8mm cup with the sliding windows.
Clamshell/ Multipurpose cup	#A0581 (Uses: slurry/liquid samples)

Step 1: Create an Optical Data File and Collect Data from the Analyzer

To collect data on calibration samples from the FoodTech analyzer, you will need to

- create a new optical data file (**File, New, Optical Data**) or
- open an existing optical data file to which you wish to add additional samples.

Select **Collect** from the **Data** menu. Answer no to the questions, then the **Collect** screen will appear:



Select the **Options** button. Connect the RS232 cable to one of the computer's com ports and to the instrument's RS232 port (beside the AC adapter port). Select the **OK** button on the **Communication Options** screen.

Ensure that the instrument is ready to communicate with the PC.
(Press SELECT, (down arrow) [Data Link], Enter)

Enter the FoodTech sample ID number from the *Calibration Numbering Setup Manual* and select **Collect** or press **Return**. The computer screen will display “Attempting to Connect” until it is successful in doing so. The screen will then prompt you to read a sample. (For more information on each of these steps, see the *User’s Manual* that came with your instrument.)

After the instrument has read the first sample in the new file, the **Save As** screen will appear. Enter a file name and click **OK**.

The **Collect** screen will reappear. Enter the next sample number and repeat the above steps. (The **Save As** screen will not reappear each time, as the computer automatically saves additional data to the same file name.)

If you make a mistake while running a sample or if the curve of the sample appears questionable (see examples of good and bad curves under Step 2: Plot Curves), highlight the sample ID number under “Last ___ samples collected” and select **Cut**.

Keep in mind: The samples do not all have to be run the same day; the file can be saved and additional samples added later.

Temperature Variation during a Calibration

Because of temperature variations in the field, it is important to vary the instrument temperature during the calibration.

To vary instrument temperature, place heating pads, set on high, around the instrument for about 15 minutes or until the instrument temperature reaches approximately 38°C. (See “*To check instrument temperature*” below.) Remove the heating pads and continue running samples.

At another point in the calibration, wrap the instrument with plastic and surround it with ice packs for 15 to 30 minutes or until the instrument temperature is lowered to 17°C. (See “*To check instrument temperature*” below.) Remove the ice packs and continue running samples.

To check instrument temperature, collect data from a sample and allow the graph of the sample reading to appear. The bottom of the window lists the sample and instrument temperatures respectively in °C.


The sample temperatures must also be varied during calibration. Cole-Parmer recommends using 3 to 5 different temperature ranges of samples during calibration covering the warmest to the coldest sample temperatures that will happen during normal use.

Step 2: Plot Curves


You can view the sample curves a number of ways. With the optical data file (*.dt1) open,

- use **P**lot, **N**ext (**P**age **D**own) and **P**lot, **P**revious (**P**age **U**p) to view one sample graph at a time;
- use **P**lot, **O**verlay, **P**age **D**own or **P**age **U**p to view more than one sample at a time (the current sample is outlined in red, while all others are blue); or
- use **P**lot, **S**how **A**ll to view all the sample graphs at once.

Use the method that works best for you. The reason for viewing the curves is to search for any whose basic shape is radically different from the other curves in the calibration. Sample curves should curve gradually rather than peaking sharply (see the following examples). Such curves are usually a result of operator error and should be deleted from the calibration. To delete the sample highlighted in red, select **D**el~~e~~te from the **E**d~~i~~t menu, press **D**elete on the keyboard, or use the **C**ut button on the toolbar.



Example A: In this example, all the curves shown are similar in shape and are acceptable.



Example B: In this second example, the highlighted curve on the bottom is clearly due to operator error—probably failure to close the light shield properly or to insert the sample.

You may wish to print the optical data gathered from the instrument but be forewarned: if you are printing the data in **Grid View**, this could take awhile. It's best to do this only if you have a good reason for doing so. From the **File** menu, select **Print**, or click on the **Print** button on the shortcut toolbar.

Step 3: Obtain Lab Analysis for all Calibration Samples

Either immediately before or after collecting optical data on the calibration samples, obtain laboratory analysis for the samples. This analysis should be from a lab you trust, as the calibration can never be more accurate than the lab values. Enter the lab analysis information into the *FoodTech Calibration Numbering Setup Manual* that came with your instrument.

Step 4: Create or Edit a Constituent Data File

The computer stores the laboratory values of constituent content for the calibration samples in a constituent data file. The same constituent data file can be used over and over for calibrations of the same product and sample set.

To create a constituent data file, select **New** from the **File** menu or use the shortcut key on the toolbar. Under **Document Type**, select **Constituent Data**. A spreadsheet screen will appear. To enter values, select **Insert** from the **Data** menu, or press the **Insert** key on the computer keyboard. An empty row of cells will appear. Type the FoodTech sample ID number and press **Return**. Use the cursor keys or the mouse to move to the next cell and enter the constituent data. To add the next sample, press **Insert**. When you are finished entering the information, select **Save As** from the **File** menu. Assign a file name of up to 8 characters that will remind you of the product values stored in the file. The software will automatically assign a ".con" ending to the file and place it in the C:\DATA directory.

To edit an existing constituent data file, select **Open** from the **File** menu or use the shortcut key on the toolbar. Under **List Files of Type**, select **Constituent Data (*.con)**. Select or type the correct file name.

To navigate within the constituent data file, use the cursor keys or the mouse. Also, the **Search** option under the **Edit** menu and on the toolbar allows you to move quickly to a particular sample ID number.


After adding or changing the appropriate information, be sure to save your changes before closing the file.

Step 5: Merge Constituent Data File to Optical Data File

With the optical data file open on the screen, from the **Data** menu select **Merge Constituents**. Select the appropriate constituent data file (*.con) for the product you are calibrating the instrument to analyze. Wait while the computer merges the two files. When an optical data file is merged with a constituent data file, the computer adds the constituent data to the optical file. For later regressions, the files should not have to be re-merged (unless changes have been made to the constituent data file).

Step 6: FoodTech Regression Analysis

After merging the files, select **Reduce** from the **Data** menu. The **Regression Parameters** screen will appear:



Tumbles: Refers to the number of samples to be “tumbled” from the regression. Samples which are statistical “outliers” should be eliminated from the regression so that they do not adversely affect the calibration. The recommended number of tumbles will be set automatically.

Type: Refers to the tumbling method.

Fast tumbling predicts the analyzer’s constituent value for each sample and compares it to the lab value. If the analyzer’s predicted value for a particular sample is off by more than double the standard error for the entire calibration set, that sample is “tumbled.” This continues until selected number of samples has been tumbled. The calibration constants (K values) are determined from the remaining samples.

Slow tumbling eliminates one sample at a time and calculates the calibration constants, comparing the standard error each time. It selects the set of K values with the lowest standard error. If more than one tumble was selected, it repeats this process with the remaining samples until the number of samples specified has been tumbled. Usually the results are the same, but as it is more accurate, Cole-Parmer recommends using slow tumbling.

Filters: This option was designed for use by Cole-Parmer researchers only. All the filters listed should have a checkmark.

Constituents: If you would like to perform the regression on more than one constituent at a time, make sure there is a 3 mark next to all the constituents for which you are calibrating. However, you may prefer to perform the regression on one constituent at a time if, for example, you want to exclude certain sample numbers from the calibration for a particular constituent.


Exclude: Click on this button to specify specific sample numbers that should be eliminated from the regression. You might choose to do this for samples that

- have no lab value or a questionable lab value for a particular constituent or
- have a bad curve.

After typing each sample number, click **ADD** to confirm its exclusion. When all the sample numbers to be excluded are listed, click **OK**.

After making all of the appropriate selections on the **Regression Parameters** screen, click **OK**. Be patient while the computer performs the regression. If the file contains a large number of samples, this would be a good time for a coffee break.

When the regression is complete, the following will appear.




Cole-Parmer recommends that you **PRINT** all these results so that they will be readily accessible when you select K values in Step 7.

Step 7: Enter Final K's


The Instrument Settings file stores all of the instrument settings and K values in one place. (Later calibrations should be added to this same file.) The instrument will need a set of K values, or constants, for each constituent within each product. A set of K's consists of a total of 17 K's (K0 through K16). For example, if you are calibrating the instrument to measure Wheat for Protein and Moisture, you will need to enter two sets of K's into the instrument settings file.

Create a new instrument settings file (**File, New, Instrument Settings**). In a cell under **Name**, type the name of the product you are calibrating for.



Click on these gray buttons under each constituent to reach the following screen.

Click on the button under Con 1 to enter the data for the first constituent in that product. The following screen will appear:



Name: Type the constituent name.

Slope & Bias: The slope should be 1.0 and the bias 0.0 (for exceptions, see *Bias and Slope Adjustments* in Step 8.)

Variance: The Variance should be set to 10.0

Too High Value:

Find the Lab Maximum for the constituent on the *FoodTech Regression Analysis* printout. Round it up to the nearest whole number. Enter this value.

Too Low Value:

Find the Lab Minimum for the constituent on the *FoodTech Regression Analysis* printout. Round it down to the nearest whole number. Enter this value.

Click on **Load K's**. From the *FoodTech Regression Analysis* printout, select the recommended set of K values for this constituent. Enter the file name for the selected set of K's (The file name is listed on the printout.). The computer will access the file and automatically enter the K's. Click **OK**. Repeat the above for each constituent in the product.

*When selecting K's for the instrument, you may want to predict the instrument's accuracy with several sets of K's to help you decide which set to use. To predict the accuracy of a given set of K values, open the optical data file. From the **Data** menu, select **Predict**. Enter the appropriate **constituent** number. Enter the file name for the **K's** you are considering. **Exclude** any specific samples which may be outliers affecting the calibration. Click on **OK**. Select **Print**. You can select to print the plot, residuals, or both. The plot provides most of the necessary information; the residuals list all the samples with their lab values, predicted readings on the instrument, and the difference between the two. The printed plot gives information on the standard error with a bias adjustment and with bias and slope adjustments. Select the standard error with the preferred degree of accuracy and use the associated K values.*

Staple the printouts of the selected K values into your *Calibration Numbering Setup Manual*.

Back in the Instrument Settings file . . .

Points 1, 2, & 3: After all the constituent-specific settings have been entered, enter the product-specific information. From the drop-down list, select the following settings.

	Point 1	Point 2	Point 3
GRAIN products	998	1037	950

For other products, examine the curves in the optical data file. **Point 1** should be the point that is always one of the two *highest* points on the curve. **Point 3** should be the point that is always one of the two *lowest* points on the curve. **Point 2** should be a point that always lies between **Point 1** and **Point 3**.

Hi and Lo Temperatures:

The Hi and Lo Ambient Temperatures should be 38.0 and 7.0 respectively.

*To print the Instrument Settings information, select **P**rint from the **F**ile menu. After the file prints, staple it into the *Calibration Numbering Setup Manual*.*

To upload the Instrument Settings data to the instrument: Connect the instrument to the computer with the RS232 cable. Select **C**omm Options from the **I**nstrument Edit menu, and make sure the correct port is selected. Put the instrument in **DATA LINK** mode. From the **I**nstrument Edit menu, select **U**pload to Instrument.

Select Save As, Enter file name. The backup file will be stored with a .ZBK extension


Step 8: Check Accuracy

To check the accuracy of the instrument after entering the final K values, take readings on ten lab-analyzed samples. Read each sample twice, and enter the readings on a photocopy of the following chart. (For instructions on taking readings, consult your *User's Manual*.)

Sample ID #	Con 1 Lab Value	Con 1 Unit Reading	Con 2 Lab Value	Con 2 Unit Reading	Con 3 Lab Value	Con 3 Unit Reading
		1 st		1 st		1 st
		2 nd		2 nd		2 nd
		1 st		1 st		1 st
		2 nd		2 nd		2 nd
		1 st		1 st		1 st
		2 nd		2 nd		2 nd
		1 st		1 st		1 st
		2 nd		2 nd		2 nd
		1 st		1 st		1 st
		2 nd		2 nd		2 nd
		1 st		1 st		1 st
		2 nd		2 nd		2 nd
		1 st		1 st		1 st
		2 nd		2 nd		2 nd
Total:						
Divided by:	10	20	10	20	10	20
Average:						

Average each column at the bottom of the table. Compare the Lab Value Average and the Unit Reading Average for each constituent. If there is a difference, a bias and/or slope adjustment may help.

To determine whether a bias or slope adjustment would correct the problem, select **New** from the **File** menu. Under "Document Type," select **Statistics**. The following screen will appear:



Using the completed accuracy chart from the previous page, enter the ten Lab Values (enter each one twice) in the **Lab** column and twenty Unit Readings in the **FoodTech** for one constituent. The computer will calculate the regression information as each sample is entered. If a particular sample is an outlier, select **Exclude** from the **Data** menu to exclude it. When all the sample readings have been entered, save the file and select **Print** from the **File** menu. The printout will include the graph along with the following information:

SEP After Bias : 0.139
 Only Adjustment
 SEP After Slope & : 0.130
 Bias Adjustment

Correlation : 0.989

* FoodTech recommends that only
 * a bias adjustment be made.

Slope : 0.937
 Slope Bias : 0.681
 Bias (No slope) : -0.165

Recommendation: It is a good idea to follow the FoodTech recommendation.

Slope Bias: This value is entered as the bias *only when* the slope is also adjusted.

Bias (No slope): When the bias is adjusted, the instrument's readings for the constituent are raised or lowered by this amount.

The printed plot gives information on the standard error with a bias adjustment and with bias and slope adjustments. Select the standard error with the preferred degree of accuracy and make adjustments accordingly. Keep in mind:

1. A slope adjustment should be made *only if*
 - the range of the samples used to test the accuracy approximates the range used in the calibration *and*
 - the suggested slope adjustment is less than 0.95 or greater than 1.05.
2. If the printout does not show a bias adjustment of ± 0.15 or more, the adjustment should not be made.

Bias Adjustments

On the completed accuracy chart on the previous page, subtract the Unit Reading Average for the constituent from the Lab Value Average to determine the bias adjustment. For instance, if the Lab Value Average is 23.5 and the Unit Reading Average is 23.8, the bias adjustment would be -0.3. Or, you can use the **Bias (No slope)** from the *Statistical Analysis* printout.

See user manual for instructions to enter the Bias adjustment.

Slope Adjustments

Slope adjustments are generally not recommended. However, if you feel that a slope adjustment may be necessary, gather the following information and contact Cole-Parmer's Calibration Department for advice.

- a completed accuracy chart for the product in question (see previous page)
- the *FoodTech Regression Analysis* printout for the calibration
- a printout of the Instrument Settings file for the product in question
- the *FoodTech Prediction Analysis* or *Statistical Analysis* printout

If it is determined that a slope adjustment is necessary, enter the **Instrument Settings file**. Make sure that you are in the appropriate product and constituent.

When making a slope adjustment, multiply the **Slope** from the prediction printout to the **Slope** in the instrument settings file, and multiply the **Bias** in the instrument settings file with the **Slope** from the prediction printout, and then add the **Slope Bias** from the printout.

After making bias or slope adjustments, be sure to **Upload** the new instrument settings to the instrument.

Appendix 1: Glossary

accuracy - how closely the readings of the calibration samples match the samples' predetermined lab values in a given instrument

bias - the average shift of the regression curve from the ideal laboratory regression curve on all samples

constants (also called K values or K's) - figures determined during calibration of the instrument which are used in the equation:

$$\text{Percentage of Constituent} = \text{Bias} + (\text{Slope})(X)$$

where X is the summation of the optical data multiplied by the constants.

constituent - the component in the sample for whose content the instrument is measuring

diffuse transmittance - light energy is transmitted and distributed evenly through the sample being tested

diodes - near-infrared light sources

K's or K values - see constants

lab analysis - the chemical or physical analysis of constituent(s) in a product

linear regression - a graphical representation of a statistical analysis of the relationship between two or more measured variables

light absorption - the amount of light energy absorbed by various constituents in the product being measured

NIR - abbreviation for near-infrared

near-infrared - the range of light energy near the visible spectrum

OD - optical data expressed in terms of voltage values for each wavelength used in a FoodTech analyzer

regression (also called **reduction**) - mathematical process to fit the optical data to the lab data

sensor - a device used in FoodTech analyzers that converts light energy to electrical signals for processing optical data

sigma - one standard deviation

spectrum - the range of OD energy measured in wavelength increments

transmittance - the ease with which light or other forms of energy can pass through a sample being tested

wavelengths - specific band widths of light energy, usually expressed in nanometers in near-infrared instrumentation